
Medtronic

BIO-MEDICUS®

Femoral Percutaneous Cannula Kit
Kit de canule fémorale percutanée
Perkutanen Femoralkanülen-Set
Equipo de cânula femoral percutânea
Femorale percutane canuleset
Set per cannula femorale percutanea
Perkuntant femoraliskanylesett
Perkutaaninen reiteen asetettava kanyylikpakkaus
Perkutan femoraliskanyl
Κιτ διαδερμικής μηριαίας κάνουλας
Perkutan femoraliskanyle
Zestaw kaniuli do przezskórnego dostępu udowego
Conjunto de cânulas percutâneas femurais
Femoral Perkutan Kanül Kiti
Sada femorální perkutánní kanyly
Perkután femorális kanüláló szett
Комплект бедренной чрескожной канюли

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic Inc.

Hav Aug 3, 2017

*Reg. Affairs Spec.
Gayatri Ghadge*

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, USA
Tel: +1 763 514 4200

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγίες χρήσης • Bruksanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati utasítás
Инструкция по эксплуатации

Caution: Federal law (USA) restricts this
device to sale by or on the order of a physician.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf der Verpackung / Explicación de los símbolos que aparecen en el etiquetado del envase / Verklaring van de symbolen op de verpakkinglabels / Splegazione dei simboli sull'etichetta della confezione / Forklaring av symboler på emballasjen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objaśnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketindeki sembollerin açıklamaları / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Объяснение символов, приведенных на этикетке упаковки

Refer to the outer package label to see which symbols apply to this product / Se referer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Verpackung / Consulte la etiqueta del envase exterior para comprobar qué símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Vedere l'etichetta sulla confezione esterna per sapere quali simboli si applicano al prodotto / Se på etiketten utenpå pakken for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksesta, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Ανατρέξτε στην ετικέτα της εξωτερικής συσκευασίας για να δείτε ποια σύμβολα αφορούν αυτό το προϊόν / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consulte as etiquetas exteriores da embalagem para ver que símbolos se aplicam a este produto / Hangi sembollerin bu ürün için geçerli olduğunu görmek için ambalaj etiketine başvurun / Podle štítků na vnější straně balení zjistíte, které symboly se týkají tohoto produktu / A termekre vonatkozó szimbólumok a külső csomagoláson láthatóak / См. на упаковке, какие символы применимы к данному продукту /

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Femoral Percutaneous Cannula Kit

Product Description

Extracorporeal Circulation Cannula with Percutaneous Insertion Set (all contents are sterile/single-use only).

Contents Include: (one set)

- 1 - Femoral Cannula
- 1 - Hemostasis Cap
- 1 - Introducer Needle - 18 Gauge
- 1 - Vessel Dilator - stepped, 8/10 French
- 1 - Vessel Dilator - stepped, 12/14 French
- 1 - Large Vessel Dilator (multi-stage cannula only)
- 1 - Cannula Introducer
- 1 - Scalpel - #11 Blade
- 1 - Guidewire 0.097 cm (0.038 in) diameter x 100 cm length for cannula tip length up to 18 cm or guidewire 0.097 cm (0.038 in) diameter x 150 cm length for cannula tip length 18 cm to 50 cm or guidewire 0.097 cm (0.038 in) diameter x 180 cm length for cannula tip length over 50 cm.

Indications for Use

This percutaneous cannula is for use by trained physicians only, to cannulate vessels, perfuse vessels or organs in a patient for cardiopulmonary bypass circulation. Standard surgical or percutaneous insertion techniques can be employed. This product is intended for use up to six hours or less.

Contraindications

Alone, this cannula is not a medical treatment device. Selection of the patient as a candidate for such procedures is the physicians' responsibility. The outcome is dependent on many variables, including patient pathology, surgical procedure, and perfusion procedures. Do not use if the patient has severe peripheral atherosclerosis or severe atrial dissection.

Warnings and Precautions

- Instructions for Use Manuals for all related extracorporeal devices should be read prior to use.
- Each cannula and introducer is for single-use only. **DO NOT** reuse or resterilize. Discard after use.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- Do not use alcohol or alcohol-based fluids for lubrication.
- Attach and manipulate cannula in such a manner as to prevent kinks or any restrictions that may alter flow such as false lumen cannulation.
- Caution should be used when attaching the cannula to other devices in the extracorporeal circuit so as not to damage the device being connected or the cannula.
- Clamping too close to connector may damage connector end.
- Do not clamp wire-reinforced sections. If wire-reinforced cannulae are used, clamping on reinforced sections may result in permanent wall distortion and/or lumen collapse.
- Since the cannula wall is opaque due to the presence of the spring, be sure to completely walk the air out of the cannula during the priming procedure.
- A cannula inserted into the left ventricle through the left atrium can perforate the left ventricular apex if undue pressure is exerted by the tip.
- Care must be taken to avoid overinsertion of the cannula and possible impingement of the cannula tip against the vessel/ventricular/atrial walls.
- The use of cold solutions and cardioplegic agents, administered externally or via the coronary circulation, may increase the stiffness of urethane plastic and contribute to ventricular perforation under conditions of cannula tip pressure during manipulation of the heart.
- Forcing the introducer hub or hemostasis cap into the straight connector could result in stress cracking of the connector.
- Do not attempt to cross-clamp the cannula with the introducer in place.
- Carefully monitor for signs of occlusion on both the inlet and outlet cannulae.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.

- Thoroughly inspect the cannula prior to use to verify that the lumen and side holes are patent and that the cannula has not been damaged or kinked.
- As with all medical devices, this device is to be used by a trained physician only.
- Initiate adequate anticoagulation therapy prior to patient cannulation.
- Only physicians trained and experienced in using percutaneous cannulation (Seldinger) technique and cardiopulmonary bypass should use this device.
- Ischemia may occur distal to the cannula placement site.
- Cannula tip placement should be confirmed by using fluoroscopy.
- Extreme care must be exercised to ensure that the arterial cannula and venous cannula are not confused while being connected to the perfusion circuit.
- Do not wipe the introducer, it has been treated to improve lubricity.

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

Complications

Possible complications include those normally associated with cardiopulmonary bypass, anticoagulation, and large bore venous and arterial catheterization.

Note: The benefits of cannulation for extracorporeal circulation must be weighed against the risk of systemic anticoagulation and subsequent propensity for hemorrhage.

Sterilization

All Bio-Medicus® Cannulae are sterilized (EtO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single-use only.

Products with Carmeda® BioActive Surface¹

The blood contact surface of the product is coated with a non-leaching Carmeda® BioActive Surface (heparin) that provides thromboresistant blood contact surfaces. Model numbers for Carmeda® coated products are prefixed with "CB."

Warnings and Precautions

Product coated with a Carmeda® BioActive Surface is intended for single-use only. Resterilization may adversely affect the Carmeda® BioActive Surface. Do not store Carmeda® coated product above 40°C.

A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The Carmeda® BioActive Surface provides a non-leaching, non-thrombogenic coating on all blood contact surfaces. The benefits of extracorporeal support must be weighed against the risk of systemic anticoagulation and must be assessed by the prescribing physician.

Instructions for Use

1. Sterile percutaneous cannula should be placed into the vascular system according to percutaneous techniques. Various components are included in this kit to assist in the access of the vessel. Various sizes of dilators and cannulae are available in order to provide the physician the opportunity of selecting the most appropriate device for each specific procedure.
When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels, and their placement sites are repaired according to standard surgical practice. (Use aseptic techniques.)
2. Carefully open the sterile pouch containing cannula and components tray.
3. Determine the length of cannula to be introduced and note the insertion markings on the cannula if present. Some cannulae are labeled with length markings in ten (10) centimeter increments for depth reference during cannula placement.
4. With the vessel identified, introduce the 18 gauge introducer needle into the desired vessel and aspirate blood.
Note: If no free flow of blood is present, reposition the needle and aspirate until adequate blood flow is established.
Caution: Ensure that the vessel is of adequate size to allow distal perfusion of the limb when the cannula is in place.
5. Insert the "J" tip of the guidewire through the needle into the vessel.
Note: Guidewire should advance freely without resistance.

¹ Manufactured under license from Carmeda AB.

Caution: Do not withdraw spring guidewire back into needle as shearing of the guidewire may result. Always remove the needle first.

6. Advance the guidewire into the vessel until it is positioned at a point beyond the desired cannula tip position.
7. Hold the guidewire in place and remove the introducer needle by sliding it away from the patient over the guidewire.
8. Enlarge the skin puncture site with the #11 scalpel blade.
9. Pass the smaller stepped vessel dilator over the guidewire. Be sure to hold the guidewire firmly in place.
10. Advance the dilator approximately three to five centimeters past the step into the vessel to enlarge the opening.
11. Remove the vessel dilator by sliding it away from the patient over the guidewire. Be sure to apply pressure to the percutaneous site and to hold the guidewire firmly in place.
12. Repeat steps 9 through 11 with the larger vessel dilator(s).
13. Assemble the femoral cannula/introducer.
 - a. For 0.95 cm (3/8 in) barbed connectors, place the hemostasis cap/introducer assembly over the barbed connector and advance the cannula introducer through the cannula.
 - b. For all other connectors, advance the cannula introducer through the cannula until the introducer hub is firmly seated in the cannula connector.

Note: Tighten the luer fitting cap on the cannula connector (if present) to ensure that the luer fitting is completely sealed.
14. Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
15. Insert tip of the cannula introducer into the vessel approximately three to five centimeters. Grasp the guidewire and hold stationary while advancing the cannula to the desired depth.

Note: Advancement of cannula and introducer with a rotating motion will make cannula placement easier.

Note: When the bold band on the introducer is visible outside the connector cap, the cannula tip and the introducer tip are at the same depth.
16. When advancing the multistage cannula into the SVC, verify positioning with tactile feel or indirect visualization method.

Note: Graduated marks on the hub end of introducer allows user to approximate introducer length beyond tip of cannula.
17. Remove the guidewire.
18. Remove the introducer until the second bold band is observed at the hemostasis cap. The introducer is now withdrawn far enough so the cannula can be clamped. Clamp the cannula and remove the introducer completely.

Note: Perfusion tubing clamps must be used. Many clamps are too weak to ensure complete clamping.
19. Remove the hemostasis cap, if used, exposing the 0.95 cm (3/8 in) barbed connector.
20. Slowly fill the remainder of the cannula tube with sterile priming fluid displacing any air out through the tube connector end.
21. Attach the extracorporeal circuit, taking care to connect the venous cannula to the venous return line (inlet to pump) and the arterial cannula to the arterial return line (outlet from membrane oxygenator), ensuring air-free connection.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE FEMORAL PERCUTANEOUS CANNULA KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE FEMORAL PERCUTANEOUS CANNULA KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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USA **Caution:** Federal law (USA) restricts this device to sale by or on the
order of a physician.

Explanation of symbols on package labeling / Explication des symboles des étiquettes sur l'emballage / Erläuterung der Symbole auf dem Verpackungsetikett / Explicación de los símbolos en el etiquetado del envase / Verklaring van de symbolen op de verpakingslabels / Spiegazione dei simboli sulle etichette della confezione / Forklaring av symboler på emballasjen / Pakkauksen symbolien selitykset / Förklaring av symboler på förpackningsetiketten / Επεξήγηση των συμβόλων στην ετικέτα της συσκευασίας / Forklaring af symboler på emballagens mærkater / Objaśnienie symboli zamieszczonych na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Ambalaj etiketi üzerindeki sembollerin açıklaması / Vysvětlení symbolů uvedených na štítcích obalu / A csomagoláson látható szimbólumok jelentése / Объяснение символов, приведенных на этикетке упаковки

Refer to the outer package label to see which symbols apply to this product / Se référer à l'étiquette sur l'emballage extérieur pour savoir quels symboles s'appliquent à ce produit / Welche Symbole für dieses Produkt zutreffen, entnehmen Sie bitte dem Etikett auf der Außenverpackung / Consulte la etiqueta del envase exterior para comprobar que símbolos son aplicables a este producto / Controleer het label op de buitenverpakking om te zien welke symbolen op dit product van toepassing zijn / Vedere l'etichetta sulla confezione esterna per sapere quali simboli si applicano a questo prodotto / Se på etiketten utenpå pakningen for å fastslå hvilke symboler som gjelder for dette produktet / Katso ulkopakkauksesta, mitkä symbolit koskevat tätä tuotetta / Se etiketten på ytterförpackningen för de symboler som gäller denna produkt / Ανατρέξτε στις εξωτερικές ετικέτες της συσκευασίας για να δείτε ποια σύμβολα ισχύουν για το προϊόν αυτό / Se mærkaten på den ydre emballage for de symboler, der gælder for dette produkt / Należy zapoznać się z etykietą na zewnętrznym opakowaniu, aby określić, które symbole dotyczą niniejszego produktu / Consultar as etiquetas exteriores da embalagem para ver que símbolos se aplicam a este produto / Bu üründe hangi sembollerin geçerli olduğunu görmek için dış ambalaj etiketine başvurun / Podle štítků na vnějším obalu zjistíte, které symboly se týkají tohoto produktu / A termékre vonatkozó szimbólumok a külső csomagoláson találhatóak / См. на этикетке внешней упаковки, какие символы применимы к данному продукту

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Femoral Insertion Kit

Product Description

This product contains the necessary components to achieve insertion of a femoral cannula and introducer. The included items are: an 18 ga Seldinger Needle, a 180 cm (71 in) × 0.038 in (0.097 cm) guidewire, a #11 scalpel blade, 2 stepped dilators (8 Fr / 10 Fr and 12 Fr / 14 Fr), and a 60 cc catheter tip syringe.

Sterile, nonpyrogenic, single use.

Indications for Use

This kit is for use by trained physicians only to assist in vessel cannulation for cardiopulmonary bypass circulation. Standard surgical or percutaneous insertion techniques can be employed. This kit is intended for use for up to 6 hours.

Contraindications

Alone, this kit is not a medical treatment device. Selection of patient as a candidate for such procedures is the physician's responsibility. The outcome is dependent on many variables including patient pathology, surgical procedure, and perfusion procedures. Do not use kit if the patient has severe peripheral atherosclerosis or severe atrial dissection.

Warnings and Precautions

- Instructions for Use manuals for all related extracorporeal devices should be read prior to use.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Do not autoclave.
- Do not use if the sterile package has been damaged or if the seal has been broken in any way.
- A guidewire inserted into the left ventricle through the left atrium can perforate the left ventricular apex if undue pressure is exerted by the tip.
- User assumes responsibility for any variations in extracorporeal procedures that could compromise the intended use of this device.
- Thoroughly inspect the components prior to use.
- As with all medical devices, this device is to be used by a trained physician only.
- Only physicians trained and experienced in using percutaneous cannulation (Seldinger) technique and cardiopulmonary bypass should use this device.

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, embolic events and dislocation, infections, blood loss, and thrombus formation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

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Sterilization

All Bio-Medicus® Femoral Insertion Kits are sterilized (EtO) prior to shipment. The lot number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is not recommended under any circumstances as this device is intended for single use only.

Instructions for Use

1. Only sterile cannula should be placed into the vascular system. Various components are included in this kit to assist in the access of the vessel. When it is determined that extracorporeal bypass is no longer required, the cannulae are carefully removed from the vessels, and their placement sites are repaired according to standard surgical practice. (Use aseptic techniques.)
2. Carefully open the sterile insertion kit tray.

3. Determine the length of cannula to be introduced and note the insertion markings on the cannula, if present. Some cannulae are labeled with length markings in 10 centimeter increments for depth reference during cannula placement.
4. With the vessel identified, introduce the 18 ga introducer needle into the desired vessel and aspirate blood.
Note: If no free flow of blood is present, reposition the needle and aspirate until adequate blood flow is established.
Caution: Ensure that the vessel is of adequate size to allow distal perfusion of the limb when the cannula is in place.
5. Insert the "J" tip of the guidewire through the needle into the vessel.
Note: Guidewire should advance freely without resistance.
Caution: Do not withdraw spring guidewire back into needle as shearing of the guidewire may result. Always remove the needle first.
6. Advance the guidewire into the vessel until it is positioned at a point beyond the desired cannula tip position.
7. Hold the guidewire in place and remove the introducer needle by sliding it away from the patient over the guidewire.
8. Enlarge the skin puncture site with the #11 scalpel blade.
9. Pass the smaller, stepped vessel dilator over the guidewire. Be sure to hold the guidewire firmly in place.
10. Advance the dilator an additional 3 to 5 centimeters past the dilator step into the vessel to further enlarge the opening.
11. Remove the vessel dilator by sliding it away from the patient over the guidewire. Be sure to apply pressure to the percutaneous site and hold the guidewire firmly in place.
12. Repeat steps 9 through 11 with the larger vessel dilator.
13. Pass the cannula/introducer assembly over the guidewire, advancing until the guidewire can be grasped beyond the introducer hub.
14. Insert tip of the cannula introducer into the vessel approximately 3 to 5 centimeters. Grasp the guidewire and hold stationary while advancing the cannula to the desired depth.
Note: Advancement of cannula and introducer with a rotating motion will make cannula placement easier.
15. Remove the guidewire.
16. Remove the introducer so the cannula can be clamped.
17. See cannula Instructions for Use for additional instructions.

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Disclaimer of Warranty

ALTHOUGH THE FEMORAL INSERTION KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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Disclaimer of Warranty

ALTHOUGH THE FEMORAL INSERTION KIT, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT.

MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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Declaration Letter


Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation:

- 1- Appendix to IFU - 66 pages;
Are true and unaltered copies of the original document.

Huda Yusuf
Sr. Regulatory Affairs Specialist Coronary & Structural Heart Disease
Management

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Overall length, cm	47,0
Distance from the distal end of the cannula to the first port, cm	1,0
Port diameter, mm	2,25
Number of ports at the distal end	12
Distance from the distal end of the cannula to the first mark	0,5
Depth mark thickness, mm	1,0
Breaking force (min), N	15
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa
Tightness	no air leakage at the junction points of the units at a load of 700 kPa
Introducer characteristics	

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Overall length cm	31,8	31,8	31,8
Distance from the distal end of the cannula to the first port, cm	1,0	1,0	1,0
Port diameter, mm	2,25	2,25	2,25
Number of ports at the distal end	6	6	6
Distance from the distal end of the cannula to the first mark, cm	5	5	5
Depth mark thickness, mm	0,9	0,9	0,9
Breaking force (min), N	15		
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa		
Tightness	no air leakage at the junction points of the units and at the junction of the tap with the Luer cap at a load of 700 kPa		
Connector	0,95 cm (3/8") Connector vented connection		
Branch technical characteristics			
Luer cap diameter, cm	1,20 – 1,22		
Luer cap mass, g	0,85 ± 0,06		
Height of the luer cap, cm	1,35 ± 0,2		
Height of the protruding part of the luer cap, cm	1,14 ± 0,2		
Thread pitch, mm	0,75 ± 0,05		
Thread length, cm	0,8 ± 0,05		
Branch length, cm	0,70 ± 0,05		
Branch outer diameter, cm	0,75 ± 0,05		
Branch inner diameter, cm	0,46 ± 0,05		
luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7)			
Introducer characteristics			

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
– Femoral arterial cannula 96530-121/-123 respectively

Outer diameter, mm	7,0	7,7
Inner diameter, mm	6,0	6,7
Mass, g	11	11
Tube length, cm	8,0	18,0
Overall length cm	31,8	31,8
Distance from the distal end of the cannula to the first port, cm	1	1
Port diameter, mm	2,25	2,25
Number of ports at the distal end	12	12
Distance from the distal end of the cannula to the first mark, cm	10	10
Depth mark thickness, mm	0,9	0,9
Breaking force (min), N	15	
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa	
Tightness	no air leakage at the junction points of the units and at the junction of the tap with the Luer cap at a load of 700 kPa	

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Distance from the distal end of the cannula to the first port, cm	1	1	1
Port diameter, mm	2,25	2,25	2,25
Number of ports at the distal end	14	14	14
Distance from the distal end of the cannula to the first mark, cm	10	10	10
Depth mark thickness, mm	1,0		
Breaking force (min), N	5		
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa		
tightness	no air leakage at the junction points of the units at a load of 700 kPa		
connector	0.95 cm (3/8") Connector non-vented connection		
luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7)			
Introducer			
			
Tube length, mm	37	37	37
Outer diameter, mm	4,0	4,7	5,3
Inner diameter, mm	1	1	1
Distance from distal tip to first radiopaque marker, cm	5,0	5,0	5,0
Mass,g	7	8	21
Breaking force/ min, N	5	5	5

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Example, specification of medical device

Outer diameter, mm	7,0	7,7	8,3	9,0	9,7
Inner diameter, mm	6,0	7,7	7,3	8,0	8,7
Mass, g	17	20	22	23	29
Tube length, cm	55,0	60,0	60,0	60,0	60,0
Overall length, cm	69,9	76,2	76,2	76,2	76,2
Distance from the distal end of the cannula to the first port, cm	1	1	1	1	1
Port diameter, mm	2,25	2,25	2,25	2,25	2,25
Number of ports at the distal end	6	6	6	18	16
Distance from the distal end of the cannula to the first mark, cm	10	10	10	10	10
Breaking force (min), N	15				
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa				
Tightness	no air leakage at the junction points of the units at a load of 700 kPa				
Connector	0,95 cm (3/8") Connector non-vented connection				
luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7)					
Introducer					

Tube length, mm	60	65	65	65	65
Outer diameter, mm	6,0	7,7	7,3	8,0	8,7

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Inner diameter, mm	1,0	1,0	1,0	1,0	1,0
Distance from distal tip to first radiopaque marker, cm	6	6	6	6	6
Mass, g	38	36	43	54	65
Breaking force, N	5				

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			
Outer diameter, mm	6,3	7,0	8,3
Inner diameter, mm	5,3	6,0	7,3
Mass, g	17	18	22
Tube length, cm	60	60	60
Overall length	76,2	76,2	76,2
Distance from the distal end of the cannula to the first port, cm	1,0	1,0	1,0
Port diameter, mm	10	10	10
Number of ports at the distal end	12	14	16

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Distance from the distal end of the cannula to the first mark, cm	6	7	8
Depth mark thickness, mm	1,0	1,0	1,0
Breaking force (min), N	15		
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa		
tightness	no air leakage at the junction points of the units at a load of 700 kPa		
Connector	0,95 cm (3/8") Connector non-vented connection		

luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7)

Guidewire 0.097 cm diameter and 180 cm length

Guidewire 0.097 cm diameter and 180 cm length

Length, cm	180 ± 2
------------	-------------

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Mass, g	$7 \pm 0,3$
Outer diameter, mm	$0,97 \pm 0,05$
Breaking force (min), N	100
Type of the distal tip	J – type
Curvature radius mm	$3 \pm 0,6$
Torque response	5-10 rotations of the proximal end lead to one complete rotation of the J-shaped distal end of the guidewire



Guidewire 0.097 cm diameter and 100 cm length,



Guidewire 0.097 cm diameter and 100 cm length,

Length, cm	100 ± 2
Mass, g	$7 \pm 0,3$
Outer diameter, mm	$0,97 \pm 0,05$
Breaking force (min), N	100
Type of the distal tip	J – type
Curvture radius, mm	$3 \pm 0,6$

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Torque response	5-10 rotations of the proximal lead to one complete rotation of the J-shaped distal end of the guidewire
Guidewire 0.054 cm diameter and 60 cm length,	
	
Length, cm	60 ± 2
Mass, g	$7 \pm 0,3$
Outer diameter, mm	$0,54 \pm 0,05$
Breaking forc (min), N	100
Type o he distal tip	J – type
Curvature radius, mm	$3 \pm 0,6$
Torque response	5-10 rotations of theproximal end lead toone complee rotation of the J-shaped distal end of the guidewire
Hemostatic cap	
	
Length, cm	$3,099 \pm 0,051$
Outer diameter, cm	1,453
Inner diameter, cm	1,140

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


Width, mm	27,91 - 28,47
Diameter , mm	16,104 - 16,307
Inner diameter, mm	14,00
Tip wall thickness, mm	1,0
Tip length, mm	10,0
Tip diameter, mm	1,1
Marking ml	bigger 1 smaller 0,2
Distance between marks, mm	6,0 ±0,1
Mass, g	7

Seldinger needle 18 G

Overall length, mm	94,0
Inner diameter, mm	1,016 +0,064/-0,013
Outer diameter, mm	1,28 +0,01/-0,13
Needle length, mm	69,85 ±1,27
Sharpening angle, °	15 ± 2
Hub width, mm	26,67±0,25
Hub length, mm	19,75
Hub height, mm	9,65 ±0,25
Mass, g	0,30±0,03

Adjustable radiopaque suture ring

Mass, g	0,5 ± 0,05
Outer diameter	7,0 ± 0,05
Inner diameter	3,0 ± 0,2
Thickness	3,0 ± 0,2

		radiopacity	radiopaque
		Others	
		Thickness of marks on all sets, mm	1,0
		The cannula of the femoral arterial / venous cannula Bio-Medicus is radiopaque	
		The bending moment of all cannula species, N·m	<12
		The location of the tip of the cannula is confirmed by fluoroscopy length tube = working length	


The permissible deviations of the technical characteristics are $\pm 5\%$, unless otherwise specified.

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The permissible deviations of the technical characteristics are ± 5%, unless otherwise specified.

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
1. Pressure loss (mmHg)
2. Flow rate (L/min of water)

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

Component	Material
Femoral percutaneous arterial kit Bio-Medicus model 96540-023	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
Cannula introducer, Large Vessel Dilator	- Introducer tube, Large Vessel Dilator: Silicone emulsion, 360 Medical Fluid, 100 CST; Thermoplastic polyurethane based on aliphatic polyester, Tecoflex EG-93A-B40, - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone
Introducer needle, 18 Gauge (1,02 mm) Seldinger needle, 18 Gauge (1,02 mm)	- Needle: Stainless steel, Type 304 (Iron - 70%, chrome - 18.5%, nickel - 9%, manganese - 1%, cobalt - 0.2%, impurities - 1.3%); - Needle pavilion: Pink Polypropylene Polymer 1895C PINK PP.

	Vessel Dilator stepped, 2,3-3,3 mm, Vessel Dilator stepped, 4,0-4,7 mm, Vessel Dilator stepped, 5,3-6,0 mm,	- Dilator tube: EM812AA low density polyethylene; Colorant blue. - Connector: Polyethylene hexene copolymer Marlex HHM 5202BN - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone
	Guidewire 0.097 cm diameter and 100 cm length, Guidewire 0.097 cm diameter and 180 cm length Guidewire 0.054 cm diameter and 60 cm length	Stainless steel, Type 304 (Iron - 70%, chrome - 18.5%, nickel - 9%, manganese - 1%, cobalt - 0.2%, impurities - 1.3%) (Material Safety Data Sheet Stainless steel, Type 304)
	Scalpel № 11 Blade,	- Blade: Chromium stainless steel Sandvik 13C26 - Holder: Low density polyethylene Indothene 16MA400; Colorant blue
	Femoral arterial cannula Bio-Medicus	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%) (Material Safety Data Sheet Stainless steel, Type 304); - Connector: Polyvinyl chloride 2212C2B / II-118 Clear 0003 - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A- 50%; Silicone - MED-4860 Part B- 50% Radiocontrast mark within the adjustable radiopaque suture ring: Barium sulphate - 100; - Luer cap: Polypropylene Pro-fax 6524 - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
	Cannula introducer	- Introducer tube: Polyvinyl chloride 9066M-RAD14, - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10) - Connector: Polyvinyl chloride 2212C2B / II-118 Clear 0003, - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone

	Hemostasis Cap for Femoral arterial cannula Bio-Medicus	Silicone LIM 6050-D2
	Vessel Dilator 7,0 mm Vessel Dilator 8,3 mm	- Dilator tube: Polyvinyl chloride 9066M-RAD14; - Connector: Polyethylene hexene copolymer Marlex HHM 5202BN, - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone
	Femoral percutaneous venous kit Bio-Medicus	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%) (Material Safety Data Sheet Stainless steel, Type 304); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A, - 50%; Silicone - MED-4860 Part B, - 50% Radiocontrast mark within the adjustable radiopaque suture ring: Barium sulphate - 100% ; - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
	Femoral venous cannula Bio-Medicus	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%) (Material Safety Data Sheet Stainless steel, Type 304); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
	Syringe, 60 cc.	- Cylinder: Polypropylene Pro-fax 6524; - Piston: Polypropylene Pro-fax 6524 - Sealant: synthetic resin free rubber
	Cannula retention clip	Polyethylene Marlex

	Protective cover for Guidewire	Polyethylene hexene copolymer Marlex HHM 5202BN
	Protective cover sleeve for Guidewire	Ethylene propylene copolymer Titanpro SM546; Colorant permanent blue VF159
	Protective cover connector for Guidewire	Polypropylene copolymer Titanpro SM546
	Protective cover clamps for Guidewire Straightener for Guidewire	Polyvinyl chloride 3006-77 CLEAR 0003; Colorant permanent blue VF159
	Straightener for Guidewire	Polyvinyl chloride 3006-77 CLEAR 0003; Colorant 8538-E1 red
	Protective cover cap for Guidewire	Silicone LIM 6050-D2
	Protective tube Protective cover for scalpel Protective tube for Introducer needle Protective tube for Seldinger needle	Polymer polypropylene P9G1Z-047
8. Данные о маркировке медицинского изделия и его упаковке/ Labelling and Packaging of medical device	Пример маркировки упаковки медицинского изделия для всех видов / Example labeling of medical device for all types 	

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Cannula primary package

Tyvek pouch made of tyvek/PE, with stucked label.

Pouch dimensions:

Length $109,2 \pm 12,7$ cm

Width $14,6 \pm 6,4$ cm

Seal width $9,5 \text{ mm} \pm 1,5 \text{ mm}$

Weight $33 \text{ g} \pm 2 \text{ g}$

Accessories primary package

Scalpel, syringe and dilators are fixed on a carton lid, packaged in a labeled tyvek/PTFE pouch

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



- for N2: 3 cm2 / s-kgs / cm2
 - Impact resistance (bursting strength):
 At 23 ° C: 725 g $\pm$ 30 g
 At -29 ° C: 550 g $\pm$ 25 g
 12. pH: Tyvek has a neutral pH of 7
 13. Adhesive characteristics:
Peel strength (strength of the adhesive layer): 8.53 N-22.14 N
Weight of the adhesive layer: 4.5 kg $\pm$ 0.7 kg / 279 m2
Adhesion transfer:> 75%
Density of glue application:> 1.03 (g / cm3)
 14. Components of packages do not react, do not contaminate, do not pass, do not affect negatively on products before, during or after sterilization.
 15. The opening of the packages by delamination is continuous and homogeneous, without delamination / rupture of the material.
 16. Used dyes:
 See packaging materials item 8.

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Symbol	Description
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Council Directive 93/42/EEC.
	Nonpyrogenic
Non-toxic	Medical device is not toxic
	Sterilized Using Ethylene Oxide
	Do Not Reuse

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<i>Stainless</i>	Marking on scalpel blade
Rx only	Only doctors agreement
	Use By
	Do Not Resterilize
	Lot Number
	Quantity
	Open Here
	Pull Tab to Open
	Caution, Consult Accompanying Documents
	Consult Instructions for Use
	For US Audiences Only
	Upper Limit of Temperature
	Catalog Number
	Manufacturer
	Date of Manufacture
	Do Not Use if Package is Damaged
	Manufactured In
	Guidewire

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**Sterilization method,
Sterilization Validation,
Packaging Validation for
Sterile devices**

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Table 1 Summary of EO Residuals

EO Exposure	Sample ID	Ethylene Oxide		Ethylene Chlorohydrin	
		EO (mg/device)	EO Limit/device	ECH (mg/device)	ECH Limit/device
1X	1	1.50	20	<0.60	9
	2	1.53		<0.60	
	3	1.61		<0.60	
	Average	1.55		<0.60	

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Disclaimer of Warranty

ALTHOUGH THE PRODUCT HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS.

THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

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- Femoral arterial cannula Bio-Medicus 7,7 mm,
- Cannula introducer,
- Hemostasis Cap,
- Introducer needle, 18 Gauge (1,02 mm),
- Vessel Dilator stepped, 2,3-3,3 mm,
- Vessel Dilator stepped, 4,0-4,7 mm,
- Vessel Dilator stepped, 5,3-6,0 mm,
- Vessel Dilator 7,0 mm
- Guidewire 0.097 cm diameter and 100 cm length,
- Scalpel № 11 Blade,
- Cannula retention clip,
- Instructions for Use.

7. Femoral percutaneous venous kit Bio-Medicus model 96600-115 contains:

- Femoral venous cannula Bio-Medicus, sizes 5.0 mm
- Cannula introducer,
- Hemostasis Cap,
- Introducer needle, 18 Gauge (1,02 mm),
- Vessel Dilator stepped, 2,3-3,3 mm,
- Vessel Dilator stepped, 4,0-4,7 mm,
- Guidewire 0.097 cm diameter and 180 cm length,
- Scalpel № 11 Blade,
- Cannula retention clip,
- Instructions for Use.

8. Femoral percutaneous venous kit Bio-Medicus model 96600-117 contains:

- Femoral venous cannula Bio-Medicus , sizes 5.7 mm
- Cannula introducer,
- Hemostasis Cap,
- Introducer needle, 18 Gauge (1,02 mm),
- Vessel Dilator stepped, 2,3-3,3 mm,
- Vessel Dilator stepped, 4,0-4,7 mm,
- Guidewire 0.097 cm diameter and 180 cm length,
- Scalpel № 11 Blade,
- Cannula retention clip,
- Instructions for Use

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

Overall length, cm	47,0
Distance from the distal end of the cannula to the first port, cm	1,0
Port diameter, mm	2,25
Number of ports at the distal end	12
Distance from the distal end of the cannula to the first depth mark	0,5
Depth mark thickness, mm	1,0
Breaking force (min), N	15
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa
Tightness	no air leakage at the junction points of the units at a load of 700 kPa

Introducer characteristics

Tube length, cm	50
Outer diameter, mm	6,7
Inner diameter, mm	1
Mass, g	7
Distance from distal tip to first radiopaque marker, cm	5

– Femoral arterial cannula, 96530-115/-117/-119 respectively

Main technical characteristics of medical devices (tolerances are $\pm 5\%$ unless otherwise specified)

Outer diameter, mm	5,0	5,7	6,3
Inner diameter, mm	4,0	4,7	6,3
Mass, g	9	9	10
Tube length, cm	18,0	18,0	18,0

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Overall length, cm			
Distance from the distal end of the cannula to the first port, cm	1,0	1,0	1,0
Port diameter, mm	2,25	2,25	2,25
Number of ports at the distal end	6	6	6
Distance from the distal end of the cannula to the first mark, cm	5	5	5
Depth mark thickness, mm	0,9	0,9	0,9
Breaking force (min), N	15		
Fluid leak under the pressure	no fluid leakage from the junctions and the cannula tube at 350 kPa		
Tightness	no air leakage at the junction points of the units and at the junction of the tap with the Luer cap at a load of 700 kPa		
Connector	0,95 cm (3/8") Connector vented connection		
Branch technical characteristics			
Luer cap diameter, cm	1,20 – 1,22		
Luer cap mass, g	0,85 ± 0,06		
Height of the luer cap, cm	1,35 ± 0,2		
Height of the protruding part of the luer cap, cm	1,14 ± 0,2		
Thread pitch, mm	0,75 ± 0,05		
Thread length, cm	0,8 ± 0,05		
Branch length, cm	0,70 ± 0,05		
Branch outer diameter, cm	0,75 ± 0,05		
Branch inner diameter, cm	0,46 ± 0,05		
luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7)			
Introducer characteristics			

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Femoral arterial cannula 96530-121/-123 respectively

Outer diameter, mm	7,0	7,7
Inner diameter, mm	6,0	6,7
Mass, g	11	11
Tube length, cm	8,0	18,0
Overall length cm	31,8	31,8
Distance from the distal end of the cannula to the first port, cm	1	1
Port diameter, mm	2,25	2,25
Number of ports at the distal end	12	12
Distance from the distal end of the cannula to the first mark, cm	10	10
Depth mark thickness, mm	0,9	0,9
Breaking force (min), N	15	

Fluid leak under the pressure
no fluid leakage from the junctions and the cannula tube at 350 kPa

Tightness
no air leakage at the junction points of the units and at the junction of the tap with the Luer cap at a load of 700 kPa

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1.20 - 1.22 Branch technical characteristics Connector vented connection Connector		Отсутствует утечка воздуха в местах соединения узлов и в месте соединения канюли Отсутствует утечка		Герметичность	
Внешний диаметр, мм 5,0 Внутренний диаметр, мм 4,0 Масса, г 15 Длина трубки канюли, см 50,0 Общая длина, см 64,8 Расстояние от дистального конца канюли до первого отверстия, см 1 Диаметр отверстия, мм 2,25 Количество отверстий на одной стороне дистального конца канюли 14 Расстояние от дистального конца канюли до первой метки, см 10 Толщина отметки глубины, мм 1,0 Минимальное усилие на разрыв, Н 5 Утечка жидкости под давлением Отсутствует утечка жидкости из соединительных узлов и трубки канюли под давлением 350 кПа Герметичность Отсутствует утечка воздуха в местах соединения узлов нагрузке 700кПа Соединитель 3/8" (0,95 см) Соединитель без отвода Канюля бедренная артериальная соответствует требованиям ISO 594-1 и ISO 594-2 (ISO 80369-7) Спецификация интубатора Длина трубки, см 37 Внешний диаметр, мм 4,0 Внутренний диаметр, мм 1 Расстояние от дистального конца интубатора до первой метки, см 5,0		Distance from the distal end of the cannula to the first port, cm 1 Port diameter, mm 2,25 Number of ports at the distal end 14 Distance from the distal end of the cannula to the first mark, cm 10 Depth mark thickness, mm 1,0 Breaking force (min), N 5 Fluid leak under the pressure no fluid leakage from the junctions and the cannula tube at 350 kPa tightness no air leakage at the junction points of the units at a load of 700 kPa Connector 0,95 cm (3/8") Connector non-vented connection luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7) Introducer Tube length, cm 37 Outer diameter, mm 4,0 Inner diameter, mm 1 Distance from distal tip to first radiopaque marker, cm 5,0 Mass,g 7 Breaking force/ min, N 5		Distance from the distal end of the cannula to the first port, cm 1 Port diameter, mm 2,25 Number of ports at the distal end 14 Distance from the distal end of the cannula to the first mark, cm 10 Depth mark thickness, mm 1,0 Breaking force (min), N 5 Fluid leak under the pressure no fluid leakage from the junctions and the cannula tube at 350 kPa tightness no air leakage at the junction points of the units at a load of 700 kPa Connector 0,95 cm (3/8") Connector non-vented connection luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7) Introducer Tube length, cm 37 Outer diameter, mm 4,0 Inner diameter, mm 1 Distance from distal tip to first radiopaque marker, cm 5,0 Mass,g 7 Breaking force/ min, N 5	
Длина трубки, см 37 Внешний диаметр, мм 4,0 Внутренний диаметр, мм 1 Расстояние от дистального конца интубатора до первой метки, см 5,0		Distance from the distal end of the cannula to the first port, cm 1 Port diameter, mm 2,25 Number of ports at the distal end 14 Distance from the distal end of the cannula to the first mark, cm 10 Depth mark thickness, mm 1,0 Breaking force (min), N 5 Fluid leak under the pressure no fluid leakage from the junctions and the cannula tube at 350 kPa tightness no air leakage at the junction points of the units at a load of 700 kPa Connector 0,95 cm (3/8") Connector non-vented connection luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7) Introducer Tube length, cm 37 Outer diameter, mm 4,0 Inner diameter, mm 1 Distance from distal tip to first radiopaque marker, cm 5,0 Mass,g 7 Breaking force/ min, N 5		Distance from the distal end of the cannula to the first port, cm 1 Port diameter, mm 2,25 Number of ports at the distal end 14 Distance from the distal end of the cannula to the first mark, cm 10 Depth mark thickness, mm 1,0 Breaking force (min), N 5 Fluid leak under the pressure no fluid leakage from the junctions and the cannula tube at 350 kPa tightness no air leakage at the junction points of the units at a load of 700 kPa Connector 0,95 cm (3/8") Connector non-vented connection luer port corresponds to the requirements of ISO 594-1 и ISO 594-2 (ISO 80369-7) Introducer Tube length, cm 37 Outer diameter, mm 4,0 Inner diameter, mm 1 Distance from distal tip to first radiopaque marker, cm 5,0 Mass,g 7 Breaking force/ min, N 5	

Femoral venous cannula, 96600-121/-123/-125/-127/-129

inner diameter, mm	1.0	1.0	1.0	1.0	1.0
Distance from distal tip to first radiopaque marker, cm	6	6	6	6	6

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Example, specefication of medical device

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Component	Material
Femoral percutaneous arterial kit Bio-Medicus model 96540-023	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
Cannula introducer, Large Vessel Dilator	- Introducer tube, Large Vessel Dilator: Silicone emulsion, 360 Medical Fluid, 100 CST; Thermoplastic polyurethane based on aliphatic polyester, Tecoflex EG-93A-B40, - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone
Introducer needle, 18 Gauge (1,02 mm) Seldinger needle, 18 Gauge (1,02 mm)	- Needle: Stainless steel, Type 304 (Iron - 70%, chrome - 18.5%, nickel - 9%, manganese - 1%, cobalt - 0.2%, impurities - 1.3%); - Needle pavilion: Pink Polypropylene Polymer 1895C PINK PP,
Vessel Dilator stepped, 2,3-3,3 mm, Vessel Dilator stepped, 4,0-4,7 mm, Vessel Dilator stepped, 5,3-6,0 mm,	- Dilator tube: EM812AA low density polyethylene; Colorant blue, - Connector: Polyethylene hexene copolymer Marlex IIIIM 5202BN - Adhesive: PVC Oxy Vinyls suspension resin Tetrahydrofuran, Cyclohexanone

Stainless steel 0.054 cm diameter and 180 cm length
Guidewire 0.054 cm diameter and 60 cm length

Scalpel № 11 Blade,

Femoral arterial cannula Bio-Medicus

Cannula introducer

Hemostasis Cap for Femoral arterial cannula Bio-Medicus

Stainless steel, Type 304 (Iron - 70%, chrome - 18.5%, nickel - 9%,
manganese - 1%, cobalt - 0.2%, impurities - 1.3%) (Material Safety
Data Sheet Stainless steel, Type 304)

- **Blade:** Chromium stainless steel Sandvik 13C26

- **Holder:**

Low density polyethylene Indothene 16MA400; Colorant blue

- **Cannula tube:** Polyurethane - Pelletan 5863-91A;

Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium
18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%)
(Material Safety Data Sheet Stainless steel, Type 304);

- **Connector:** Polyvinyl chloride 2212C2B / H-118 Clear 0003

- **Adjustable radiopaque suture ring:** Silicone - MED-4860 Part
A- 50%;

Silicone - MED-4860 Part B- 50%

Radiocontrast mark within the adjustable radiopaque suture ring:

Barium sulphate - 100;

- **Luer cap:**

Polypropylene Pro-fax 6524

- **Adhesive:** PVC Oxy Vinyls suspension resin,

Tetrahydrofuran, Cyclohexanone

- **Depth marks:** Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl
alcohol - 10%, Cellulose nitrate - 10%)

- **Introducer tube:** Polyvinyl chloride 9066M-RAD14,

- **Depth marks:** Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl
alcohol - 10%, Cellulose nitrate - 10)

- **Connector:**

Polyvinyl chloride 2212C2B / H-118 Clear 0003,

- **Adhesive:** PVC Oxy Vinyls suspension resin,

Tetrahydrofuran, Cyclohexanone

Silicone LIM 6050-D2

Permanent blue VFI59 Polyethylene propylene copolymer Titanium SM546: Colorant	
Vessel Dilator 7,0 mm Vessel Dilator 8,3 mm	- Dilator tube: Polyvinyl chloride 9066M-RAD14; - Connector: Polyethylene hexene copolymer Marlex HHM 5202BN, - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone
Femoral percutaneous venous kit Bio-Medicus	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%) (Material Safety Data Sheet Stainless steel, Type 304); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adjustable radiopaque suture ring: Silicone - MED-4860 Part A, - 50%; Silicone - MED-4860 Part B, - 50% Radiocontrast mark within the adjustable radiopaque suture ring: Barium sulphate - 100% ; - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
Femoral venous cannula Bio-Medicus	- Cannula tube: Polyurethane - Pelletan 5863-91A; Cannula frame: Stainless steel, Type 304 (Iron 70%, chromium 18.5%, nickel 9%, manganese 1%, cobalt 0.2%, impurities 1.3%) (Material Safety Data Sheet Stainless steel, Type 304); - Connector: Polyvinyl chloride 2212C2B / H-118 Clear 0003 - Adhesive: PVC Oxy Vinyls suspension resin, Tetrahydrofuran, Cyclohexanone - Depth marks: Ink L-219 (Methyl ethyl ketone - 80%, Isopropyl alcohol - 10%, Cellulose nitrate - 10%)
Syringe, 60 cc,	- Cylinder: Polypropylene Pro-fax 6524; - Piston: Polypropylene Pro-fax 6524 - Sealant: synthetic resin free rubber
Cannula retention clip	Polyethylene Marlex
Protective cover for Guidewire	Polyethylene hexene copolymer Marlex HHM 5202BN

- Dilator tube: Polyvinyl chloride 9066M-RAID14;
- Connector: Polyethylene hexene copolymer Marlex 5202BN.

Vessel Dilator 7,0 mm
Vessel Dilator 8,3 mm

polyethylene propylene copolymer Titanpro SM546, Colorant permanent blue VF159

Protective cover connector for Guidewire

Polypropylene copolymer Titanpro SM546

Protective cover clamps for Guidewire
Straightener for Guidewire

Polyvinyl chloride 3006-77 CLEAR 0003; Colorant permanent blue VF159

Straightener for Guidewire

Polyvinyl chloride 3006-77 CLEAR 0003; Colorant 8538-EI red

Protective cover cap for Guidewire

Silicone LIM 6050-D2

Protective tube
Protective cover for scalpel
Protective tube for Introducer needle
Protective tube for Seldinger needle

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primary package

Tyvek pouch made of Tyvek/PE, with sticked label.

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
Inner sheet in primary package

High-density polyethylene 0.35 HDPE White

Tray in primary package

- Tray: Thermoformed plastic Tyvek 1059B, coated TPT-
- Cover:
Thermoformed plastic Tyvek 1059B, coated TPT-0265

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Technical characteristics according to ISO 11607-1/ISO 11607-2

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«...ability of substrate: 2% P_{max}

2. Toxicological properties:
 See packaging materials item 8
5 Physical and chemical properties:
 Overall dimensions, see item 8
 Film thickness: 0.25 mm ± 0.05 mm
 Gloss level: 108 gloss units, at an angle of 45°
 Stress at the yield point: 52 MPa
 Young's modulus: 1720 MPa
 Coefficient of friction: > 1
 Coefficient of transmission of water vapor: 6 g / m²
3.1 Characteristics of Tyvek fibrous material:
 Bendtsen air permeability: 572 ml / min
 Permeability rate of water vapor: 1615 g / m²
 Height of hydrostatic pressure: 1442 Pa
 Tensile strength in the longitudinal direction: 203 N
 Tensile strength in transverse direction: 200 N
 Relative elongation in the longitudinal direction: 20%
 Elongation in the transverse direction: 24%
4. Compatibility with the methods of sterilization with which the use of the material is intended:
Methods of sterilization:
 - Ethylene oxide
5. Compatibility with the formation and gluing processes:
 Termoskleivaetsya with a large number of thermoplastics, as PVC, styrene, PETG (PETG), HIPS; with polyethylene some kinds of polypropylene
6. Shelf life after sterilization: 3 years
7. The materials are non-alkaline, non-toxic and odorless.
8. There are no holes, cracks, tears, wrinkles and refinements in the materials.
9. Density: 1.75 g / cm³
10. The level of purity of the material: acceptable.
11. Specific and minimum physical characteristics:
 Tensile strength (tear resistance): 2. 4 Н ± 0.1 Н
 Tensile load (tensile strength): 50 MPa
 Thickness deviation: ± 0.05 mm
Breathability:
 - for CO₂: 20 cm² / s-kgs / cm²
 - for O₂: 8 cm² / s-kgs / cm²

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At -29 ° C: 550 g ± 25 g
At 23 ° C: 725 g ± 30 g
- Impact resistance (bursting strength):
- for N2: 3 cm² / s-kgs / cm²

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Stainless	Marking on scalpel blade
Rx only	Only doctors agreement
	Use By
	Do Not Resterilize
	Lot Number
	Quantity
	Open Here
	Pull Tab to Open
	Caution, Consult Accompanying Documents
	Consult Instructions for Use
	For US Audiences Only
	Catalog Number
	Manufacturer
	Date of Manufacture
	Do Not Use if Package is Damaged
	Manufactured In
	The sign of the appeal, which indicates that the products marked by them have passed all the procedures of assessment established in the technical

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GUIDE WIRE		Guidewire
EC	REP	Authorized Representative in the European Community

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Information in Russian sticker

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable).

Sample of the sticker

 	Name of medical device	
	Manufacturer: Medtronic Inc.	Code UPN
	Country of origin: USA	Code CFN
	Manufacturer's authorized representative: Medtronic LLC	
	Regrestretion Certificate No _____ or _____	
Non toxic!	Date of manuf./Lot/UBD: See on pack	See IFU See packaging symbols

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All Bio-Medicus Cannulae are sterilized (EtO) prior to shipment. The Lot Number is designated on the outer package label. This device should not be used if the sterile package has been damaged or if the seal has been broken. Resterilization is

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not recommended under any circumstances as this device is intended for single-use only.

Cycle characteristics

Sterilizer	
Preconditioning	43±4,3°C, 65% RH
Initial Vacuum	5,1 кПа
Humidity Inject	3,4 кПа
Product Temperature (exposure)	51°C (avg) 41,6 - 61°C
Product Humidity (exposure)	70% (avg) 48 - 100%
Dynamic Steam Conditioning	70 Minutes
EO Inject	32 kPa
EO Concentration	500 mg/L
Exposure Time	60 Minutes
Post Exposure Purges	14
Aeration	48,9±4,3°C

An assessment of the EO residuals on the Femoral Venous Cannula was conducted. Sterilization was conducted using the maximum fully-loaded worst-case load. Three (3) devices were tested after one (1) and two (2) EO exposures. Results of testing demonstrated the Ethylene oxide and Ethylene chlorohydrin (ECH) residuals met requirements of ISO 10993-7:2008. Data are summarized in Table 1.

Table 1 Summary of EO Residuals

EO Exposure	Sample ID	Ethylene Oxide		Ethylene Chlorohydrin	
		EO (mg/device)	EO Limit/device	ECH (mg/device)	ECH Limit/device

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Contact for all questions, including complaints, repair and service
of medical devices the Authorized Representative in Russia:

Medtronic, LLC,
123317, Moscow, Presnenskaya Naberezhnaya., 10
Tel.: (495) 580-73-77, Fax: (495) 580-73-78,
www.medtronic.ru

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© 2003, 2011 Medtronic