State of Minnesota

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SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Documents for 3M Health Care Business

2. has been signed by Irma Vargas

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Irma Vargas, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: July 21, 2017

7. by: Secretary of State, State of Minnesota

8. File No: 9583231-6

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF in English is truthful and accurate.


Regina Feferman
International Regulatory Affairs Associate
3M Health Care


Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 20th day of July 2017.

Irma Vargas, Notary Public
My Commission Expires: January 31, 2020

Instruction for use

3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF

Description: The 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF is equivalent in performance to the Bowie-Dick towel pack described in AAMI ST-79. The test pack consists of lead-free steam-sensitive chemical indicator test sheets positioned in a package of layered porous materials. The primary test sheet consists of a lead-free steam-sensitive chemical indicator ink printed on paper as a yellow-colored diagonal pattern and is positioned near the center of the porous pack. The primary test sheet will turn a uniform dark brown/black color except when air removal failures such as air leaks occur. An air removal failure is indicated by a lighter-colored area in the indicator ink pattern of an otherwise dark-colored test sheet. The test pack also contains an Early Warning sheet which is placed at a lower point in the pack than the primary test sheet. The Early Warning sheet has a lead-free steam-sensitive chemical indicator circle that turns from yellow to dark brown/black for early detection of air removal failures before they appear on the primary test sheet. 3M™ Comply™ Bowie-Dick Test Packs meet ISO 11140-5:2007 and are Type 2 (Category s2) Special Indicators as categorized by ISO 11140-1:2014.

Indications for Use

The 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF is designed for testing air removal efficiency of 132- 134°C (270- 273°F) dynamic-air-removal steam sterilizers.

Contraindications

Do not use the 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF for air removal monitoring of gravity displacement steam, dry heat, ethylene oxide, or other low temperature sterilization processes.

Instructions for Use

1. Test Frequency

The Bowie-Dick test must be carried out each day in an empty sterilizer before the first processed load. A warm-up cycle (dry time may be omitted) should be run first to properly heat the sterilizer prior to the Bowie-Dick test. If the sterilizer is in use for 24 hours each day, the test may be done at any time but preferably at the same specified time each day. The Bowie-Dick test should also be performed as part of qualification testing following sterilizer installation, relocation, malfunction, major repairs, and sterilization process failures.

2. Test Cycle

Before putting the 3M™ Comply™ Bowie-Dick Plus Test Pack 00135LF in the sterilizer, inspect the lead-free indicator on the label of the test pack. If any part of the indicator on the label has changed from yellow to light brown or darker, the test pack should not be used. Place the test pack horizontally, label side up, on the bottom shelf of the rack, over the drain in an otherwise empty, dynamic-air-removal steam sterilizer. Run the sterilizer for 3.5 minutes at 132-134°C (270- 273°F). If timer must be set for full minute intervals, the test may be conducted for 4 minutes. (Do not process for more than 4 minutes; results will be invalid.)

3. Interpretation of Results

Promptly remove the test pack after completion of the cycle. Once processed, the indicator ink on the label will change from yellow to light brown or darker. Allow test pack to cool; retrieve and interpret test sheets.

4. Primary Test Sheet

Pass - A satisfactory test result is indicated by a primary test sheet that shows uniform dark brown/black color development. A uniform color change result indicates rapid steam penetration, adequate air removal, and lack of significant air leaks. (See picture for example)

Fail - An unsatisfactory test result is indicated by a primary test sheet that shows non-uniform color development with a lighter-colored area on the test sheet, usually near the center of the indicator ink pattern. A non-uniform color change result is caused by residual air in the sterilizer chamber and indicates incomplete air removal from the chamber or failure of the sterilizer to hold a vacuum during the test cycle. (See picture for example) If a non-uniform color change is noted, it should be reported immediately to the supervisor and the sterilizer should not be used until it is examined for possible malfunction.

5. Early Warning - If the yellow indicator ink circle on the Early Warning test sheet does not turn to a color that is as dark as or darker than the color standard, impending sterilizer failure is indicated; routine sterilizer maintenance should be scheduled soon.

Note: A color interpretation chart is available upon request from the 3M Health Care Helpline at 1-800-228-3957.

6. Maintain all test sheets in records.

Storage

Best stored under normal room conditions: 59-86°F (15-30°C), <50% relative humidity.

- Protect from direct light. Do not store near strong alkaline or acidic products such as cleaning or disinfecting agents.
- After use the indicator sheets will not change visually within 24 months when stored at above conditions.

Information about repair and maintenance

This medical device isn't subject to maintenance and repair.

Utilization

Dispose of the device in accordance with applicable local/regional/national/international regulations.

Transport conditions

Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Warranty

3M Health Care. warrants this product will be free from defects in material and manufacture.

User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Health Care sole obligation shall be repair or replacement of the 3M Health Care product.

Explanation of Symbols

 REF Catalogue Number

Caution, see instructions for use

Do not reuse

Use by date

Batch code

Manufacture

Date of manufacture

Product is designed for testing air removal efficiency in dynamic air-removal steam cycles.

Product is designed for use with steam sterilization cycles.

Made in U.S.A. by
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957
www.3M.com/C3CD

Authorised representative of the
manufacturer in Russia:
ZAO "3M Russia", Russia 121614,
Moscow, Krylatskaya str.
bldg.17/3
T.+74957847474

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Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Comply™ Bowie-Dick Plus Test Pack 1233LF in English is truthful and accurate.

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

Instruction for use

3M™ Comply™ Bowie-Dick Type Test Pack 1233LF

Description: The 3M™ Comply™ Bowie-Dick Type Test Pack 1233LF is equivalent in performance to the Bowie-Dick towel pack described in AAMI ST-79. The test pack consists of a lead-free steam-sensitive chemical indicator test sheet positioned in a package of layered porous materials. The test sheet consists of a lead-free steam-sensitive chemical indicator ink printed on paper as a yellow-colored diagonal pattern and is positioned near the center of the porous pack. The test sheet will turn a uniform dark brown/black color except when air removal failures such as air leaks occur. An air removal failure is indicated by a lighter-colored area in the indicator ink pattern of an otherwise dark-colored test sheet. 3M™ Comply™ Bowie-Dick Test Packs meet ISO11140-5:2007 and are Type 2 (Category s2) Special Indicators as categorized by ISO11140-1:2014.

Indications for Use

The 3M™ Comply™ Bowie-Dick Type Test Pack 1233LF is designed for testing air renewal efficiency of 132- 134°C (270- 273°F) dynamic-air-removal steam sterilizers.

Contraindications

Do not use the 3M™ Comply™ Bowie-Dick Plus Test Pack 1233LF for air removal monitoring of gravity displacement steam, dry heat, ethylene oxide, or other low temperature sterilization processes.

Instructions for Use

1. Test Frequency

The Bowie-Dick test must be earned out each day in an empty sterilizer before the first processed load. A warm-up cycle (dry time may be omitted) should be run first to properly heat the sterilizer prior to the Bowie-Dick test. If the sterilizer is in use for 24 hours each day, the test may be done at any time but preferably at the same specified time each day. The Bowie-Dick test should also be performed as part of qualification testing following sterilizer installation, relocation, malfunction, major repairs, and sterilization process failures.

2. Test Cycle

Before putting the 3M™ Comply™ Bowie-Dick Type Test Pack 1233LF in the sterilizer inspect the lead-free Indicator on the label of the test pack. If any part of the indicator on the label has changed from yellow to light brown or darker, the test pack should not be used. Place the test pack horizontally, label side up, on the bottom shelf of the rack, over the drain in an otherwise empty dynamic-air-removal steam sterilizer. Run the sterilizer for 3.5 minutes at 132- 134°C (270- 273°F). If timer must be set for full minute intervals the test may be run for 4 minutes (Do not process for more than 4 minutes: results will be invalid.)

3. Interpretation of Results

Promptly remove the test pack from the sterilizer after completion of the cycle. Once processed, the indicator ink on the label will change from yellow to light brown or darker. Allow the test pack to cool: retrieve and examine the test sheet.

Pass - A satisfactory test result is indicated by a test sheet that shows uniform dark brown or black color development. A uniform color change result indicates rapid steam penetration, adequate air removal, and lack of significant air leaks. (See picture for sample)

Fail - An unsatisfactory test result is indicated by a test sheet that shows non-uniform color development with a lighter-colored area in the indicator ink pattern, usually near the center of the

test sheet A non-uniform color change result is caused by residual air in the sterilizer chamber and indicates incomplete air removal from the chamber or failure of the sterilizer to hold a vacuum during the test cycle. (See picture for example) If a non-uniform color change is noted, it should be reported immediately to the supervisor and the sterilizer should not be used until it is examined for possible malfunction.

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Warranty

3M Health Care warrants this product will be free from defects in material and manufacture.

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 Manufacture:

 Date of manufacture:

 AIR REMOVAL Product is designed for testing air removal efficiency in dynamic air-removal steam cycles.

 STEAM Product is designed for use with steam sterilization cycles.

Made in U.S.A. by
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA
1-800-228-3957

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9. Seal/Stamp:

10. Signature:


A handwritten signature in black ink that reads "Steve Simon".

Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Comply™ 00104 Steam Multi-Parameter Chemical Indicator in English is truthful and accurate.

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

Stated Value: 21.1 minutes at 121°C (250°F)

Description: The 3M Comply 00104 Steam Multi-Parameter Chemical Indicator consists of a steam and temperature sensitive indicator ink printed on a 2.2 cm (7/8 in.) wide by 5.7 cm (2 5/16 in.) long plastic laminated paperboard strip. The indicator ink responds to steam, time, and temperature by changing from a purple to a green color. The indicator ink is printed in a flower pattern; two rectangular spots labeled clockwise with number 1 and 2, and a circular spot to the center of the pattern labeled with the number 3. The indicator is designed to respond to exposure of sterilization conditions in a progressive manner – number 1-3 indicator spots turning from purple to green in sequential fashion.

Indications for Use

Use the Comply 00104 steam multi-parameter chemical indicator for inside the pack monitoring of all 134°C (273°F) steam sterilization cycles. The flower indicator pattern provides more information than a “pass or fail” color interpretation; the sequential clockwise pattern of color development also indicates one level of insufficient exposure and one level of overexposure.

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Directions for Use

1. Place a Comply 00104 steam multi-parameter chemical indicator in the center of each pack, peel pouch or tray to be steam sterilized. For rigid containers, place a strip in each corner, in two opposing diagonal corners, or in the area determined by product testing to be the greatest challenge. (When using the plastic holder (product no. 81102), place the indicator in holder at the center of the pack and extend the string to a point slightly beyond the inner pack contents).
2. Process the load according to established procedures.
3. After processing, remove the Comply 00104 steam multi-parameter chemical indicator strip and observe for color change. The top two indicator spots (numbers 1 and 2) should have turned from purple to green indicating that the items in the pack, peel pouch, unwrapped tray, or rigid container were exposed to sufficient steam sterilization conditions.
4. The failure of indicator spot number 2 to change from purple to green indicates insufficient exposure to steam sterilization conditions inside the pack – these items should be returned for reprocessing.
5. Color changes with indicator spots numbers 1 and 3 provide information on the degree of insufficient or sufficient conditions of exposure. Refer to graphics below for a visual

way of indicator color changes and interpretation.

Storage/Shelf Life

Store unopened packages containing Comply 00104 steam multi-parameter chemical indicator in a dry (<50% relative humidity) condition at room temperature 15-30°C (59-86°F). Comply 00104 steam multi-parameter chemical indicator contained in an unopened package have a 18 month shelf life from the date of manufacture when stored at recommended conditions. The expiration date is printed on the package label.

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Warranty

3M Health Care. warrants this product will be free from defects in material and manufacture. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Health Care sole obligation shall be repair or replacement of the 3M Health Care product.

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
Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Comply™ 00103 Steam Multi-Parameter Chemical Indicator in English is truthful and accurate.

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
3M™ Comply™ 00103 Steam Multi-Parameter Chemical Indicators

STEAM

Class D Multi-variable Indicator as classified by EN-867-1.

Class 4 Multi-parameter Indicator as classified by ANSI/AAMI ST60 and ISO 11140-1.

Stated Value: 21.1 minutes at 121°C (250°F)

Description: The 3M Comply 00103 Steam Multi-Parameter Chemical Indicator consists of a steam and temperature sensitive indicator ink printed on a 2.2 cm (7/8 in) wide by 5.7 cm (2 5/16 in) long, plastic-laminated paperboard strip. The indicator ink responds to steam, time, and temperature by changing from a purple to a green color. The indicator ink is printed in a flower pattern; three rectangular spots labeled clockwise with numbers 1, 2, & 3 with a circular spot to the center of the pattern labeled with the number 4. The indicator is designed to respond to exposure of sterilization conditions in a progressive manner- number 1-4 indicator spots turning from purple to green in sequential fashion.

Indications for Use

Use the Comply 00103 steam multi-parameter chemical indicator for inside the pack monitoring of all 121°C (250°F) steam sterilization cycles. The flower indicator pattern provides more information than a “pass or fail” color interpretation; the sequential clockwise pattern of color development also indicates two levels of insufficient exposure and one level of overexposure.

Contraindications

No contraindications have been identified for these products.

Directions for Use

1. Place a Comply 00103 steam multi-parameter chemical indicator in the center of each pack, peel pouch, or tray to be steam sterilized. For rigid containers, place a strip in each corner, in two opposing/diagonal corners, or in the area determined by product testing to be the greatest challenge. (When using with the plastic holder (product no. 81102), place the indicator in holder at the center of the pack and extend the string to a point slightly beyond the inner pack contents.)
2. Process the load according to established procedures.
3. After processing, remove the Comply 00103 steam multi-parameter chemical indicator and observe for color change. All the three top indicator spots (numbers 1, 2, and 3) should have turned from purple to green indicating that the items in the pack, peel pouch, unwrapped tray, or rigid container were exposed to sufficient steam sterilization conditions.
4. The failure of indicator spot number 3 to change from purple to green indicates insufficient exposure to steam sterilization conditions inside the pack-these items should be returned for reprocessing.
5. Color changes with indicator spots numbers 1, 2, and 4 provide information on the degree of insufficient or sufficient conditions of exposure. Refer to graphics below for a visual display of indicator color changes and interpretation.

Storage/Shelf Life

Store unopened and opened packages containing Comply 00103 steam multi-parameter chemical indicator in a dry (<50% relative humidity) condition at room temperature 15-30°C (59-86°F).

Comply 00103 steam multi-parameter chemical indicator contained in an unopened package has an 18 month shelf life from the date of manufacture when stored at recommended conditions. The expiration date is printed on the package label.

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This medical device isn't subject to maintenance and repair

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Authorised representative of the
manufacturer in Russia:
ZAO "3M Russia", Russia 121614,
Moscow, Krylatskaya str. bldg.17/3
T.+74957847474

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Steve Simon
Secretary of State

To whom it may concern

In my capacity as International Regulatory Affairs Associate, 3M Health Care Business, St. Paul, MN, USA, I certify that the attached Instruction for Use (IFU) for 3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A in English is truthful and accurate.

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
Instruction for use

3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A

Description:

3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A are chemical indicators consisting of a paper wick and steam and temperature sensitive chemical pellet contained in a paper/film and paper/film/foil laminate. The chemical pellet melts and migrates as a dark color along the paper wick. The migration is visible through a window marked ACCEPT or REJECT; the extent of migration depends on steam, time, and temperature. 3M™ Comply™ SteriGage™ Steam Chemical Integrators are Type 5 (Category i5) Integrating Indicators as categorized by ISO 11140-1:2014.

Stated Values

See Stated Values listed in the table on the front of the package. The table lists the time in minutes at the corresponding temperature for chemical indicators in this lot to reach their endpoint when tested in a steam resistometer under highly controlled conditions with saturated steam. This compares with a 1×10^6 spore population of *Geobacillus stearothermophilus* with a D-value of 1.6 minutes at 121°C and a z-value of 10.0°C.

This Product Is Not Made With Lead.

Indications for Use

Use 3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A for pack control monitoring of all 118-138°C (245-280°F) steam sterilization cycles.

Contraindications

None.

Precautions

Do not use 3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A to monitor dry heat, ethylene oxide, or other low temperature sterilization processes.

Directions for Use

1. To preserve as much package area as possible for folding and resealing, use scissors on the marked area at the top of the foil package to make the initial opening. Remove only the number of 3M™ Comply™ SteriGage™ Steam Chemical Integrators 1243A needed. Reseal the package by folding the opened end over at least two times.
2. Place a 3M™ Comply™ SteriGage™ Steam Chemical Integrator 1243A in each pack, peel pouch, container system or tray to be steam sterilized in the area determined to be the least accessible to steam penetration.
3. Process the load according to established procedures.
4. After processing, the dark color should have entered the ACCEPT window of the 3M™ Comply™ SteriGage™ Steam Chemical Integrator 1243A. If the dark color has not entered the ACCEPT window, a REJECT result is indicated and the items in the pack, peel pouch, container system, or tray were not exposed to sufficient steam sterilization conditions. These items should be returned for reprocessing.

Storage

- Best stored at normal room temperature conditions, 15-30C (59-86.F) and 40-60% relative humidity. Store away from direct sunlight. Do not store near strong alkaline or acidic products such as cleaning/disinfecting agents.
- After use, the indicator will not change visually within 6 months when stored at above conditions

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 Manufacturer

 Date of manufacture

 Product is designed for use with steam sterilization cycles.

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