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**Gambro Dialysatoren GmbH**

## Theranova dialyzer.

The variants of execution to which this instruction applies:

*Theranova 400*

*Theranova 500*

### Definition of expressions used in this manual

The Theranova dialyzer, in the embodiments in the text of this document, is referred to as a "dialyzer", but to specify the version as "Theranova 400" or "Theranova 500".

**WARNING!** is used to alert the operator not to take a specific action that can cause potential hazard which, if ignored, can result in an adverse reaction, injury or death.

**CAUTION!** is used to alert the operator to take a specific action to protect against a possible hazard which, if ignored, could have an adverse effect on the patient or equipment.

**NOTE!** is used to remind the operator of normal treatment functions and what is the suitable action in a particular situation.

### INDICATIONS FOR USE

Theranova dialyzers are indicated for treatment of chronic and acute renal failure by Hemodialysis.

### PURPOSE OF MEDICAL DEVICE

Theranova dialyzers are indicated for treatment of chronic and acute renal failure by Hemodialysis.

### CONTRAINDICATIONS

There are no known contraindications for the use of the dialyzers.

**CAUTION!** Do not use THERANOVA dialyzers for HDF (hemodiafiltration) or HF (hemofiltration) due to higher permeability of larger molecular weight proteins such as albumin.

### ADVERSE EVENTS

Certain types of adverse reactions may occur due to operational factors associated with the treatment.

**WARNING!** Selection of a suitable dialyzer configuration and of treatment parameters should be considered based on patient characteristics (body size, weight, cardiovascular status, therapy tolerance etc.) and clinical requirements.

**WARNING!** Adherence to set-up and priming procedures, proper management of fluid removal, electrolyte balance, adjustment of the pH-value of the dialysate, anticoagulation, blood and dialysate flow rates as well as monitoring of the overall treatment parameters are essential to avoid side effects which may be associated with Hemodialysis.

**WARNING!** To minimize bacterial and endotoxin contamination, compliance to water and dialysate quality standards such as ANSI/AAMI RD62 or ISO 11663 should be monitored.

In rare cases (e.g. for patients having a history of being highly sensitive to a variety of substances), hypersensitivity reactions may occur during dialysis, particularly at the onset of treat-

ment. In severe cases, dialysis must be discontinued and appropriate medical intervention administered. After rinsing the dialysate compartment the fluid in the blood compartment may contain bicarbonate. During patient connection, the volume and rate of infusion of this fluid needs to be considered, especially for patients with metabolic acidosis. Discard as much as possible of the priming solution.

## CLINICAL SIDE EFFECTS

**WARNING!** If an **internal blood leak** is observed, the operator must stop the treatment session and replace the dialyzer. Do not return the blood to the patient since it may have become contaminated by the dialysis fluid. If necessary, administer adequate replacement solution to the patient to compensate for the blood loss.

**WARNING!** If an **external blood leak** is observed, secure connections or replace the dialyzer. If necessary, administer adequate replacement solution to the patient to compensate for the blood loss.

**WARNING!** If a **dialysate compartment leak** is observed the operator has to check the correct placement of the dialysate connectors or to stop the treatment and replace the dialyzer. If necessary, administer adequate replacement solution to the patient to compensate fluid imbalance.

**WARNING!** If **air enters the extracorporeal blood circuit**, an air embolism may occur. This can be hazardous to the patient. To minimize the risk of air embolism constant monitoring of the extracorporeal blood circuit, both visually and with an air detector, is recommended. Strict adherence to the manufacturer's recommended activities will facilitate the removal and prevent the accumulation of air in the dialyzer before the treatment session. If air enters or is identified in the dialyzer during priming and cannot be removed with the use of additional priming, the dialyzer must be replaced.

**CAUTION!** If **clotting occurs in the dialyzer** both the dialyzer and the blood lines must be replaced. Flush the vascular access devices according to clinic procedure. Discard the dialyzer and the blood lines. Do not return the blood to the patient.

## DESCRIPTION

See Specifications at the end of the instructions.

Device with Sterile and non-pyrogenic fluid pathways. Capillary dialyzer, intended only for single use.

## GENERAL WARNINGS AND CAUTIONS

The dialyzers are for single use only.

The quality of the product is guaranteed for first use only, and only when prepared and used according to the procedures described.

Reprocessing this dialyzer may cause serious damage to the product resulting in patient hazard.

**WARNING!** Treatment parameters should be chosen within the limitations given in the specification table

**WARNING!** Use the dialyzers only with dialysis equipment which can accurately control and monitor the ultrafiltration rate, and blood in the dialysate circuit (blood leak detector).

**WARNING!** Elimination of drugs by the dialyzers requires specific considerations.

**CAUTION!** Federal law (USA) restricts this device to sale by or on the order of a physician.

**CAUTION!** Do not apply isolated/sequential ultrafiltration when using THERANOVA dialyzers, due to higher permeability for larger plasma proteins such as free haemoglobin. In dialysis patients plasma free haemoglobin is usually present in low concentrations (up to 239 mg/L). During isolated/sequential ultrafiltration free haemoglobin is filtered and concentrated in the dialysate compartment. This leads to a reddish coloration of the ultrafiltrate which may trigger the internal blood leak detector.

### **WARRANTY**

The manufacturer warrants that the capillary dialyzer has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the dialyzer becoming apparent before the expiry date.

**Expiry date:** Refer to information on the unit container label

**The guaranteed shelf life** - 3 years when stored under recommended storage conditions. The warranty period of storage is guaranteed provided that the packaging or packaging is not damaged and that the protective caps are in place.

The warranty above is in lieu of, and to the exclusion of, any other warranties, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or fitness or other warranties, which extend beyond those described.

The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the dialyzer and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the dialyzer, whether as a result of any defect therein or otherwise.

The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings, directions and instructions in the labeling, damage arising from events after the manufacturer's release of the dialyzer before use in order to insure that the dialyzer is in proper condition, or any warranty given by independent distributors or dealers.

### **PERFORMANCE DATA**

Performance data obtained in accordance with the ISO 8637 standard.

Information about test methods used to obtain performance data is available from the manufacturer upon request.

## **DIRECTIONS FOR USE**

### **Set-up instructions**

**WARNING!** Do not use if package or product is damaged or if protection caps are not in place.

**WARNING!** Aseptic technique is required throughout dialysis preparations and treatment to avoid contamination.

**WARNING!** For automated priming modes please refer to the clinic's protocol and the operator's manual of the dialysis equipment in use.

Recommendations for manual priming:

- Prepare and connect a bag with isotonic saline (e.g. 0.5 l) to the blood lines.
- Remove the protective caps from the blood ports and use them to close the dialysate ports.
- Position the filter vertically and connect the blood lines to the filter to fill up the blood compartment from bottom to top.
- Start the blood pump at 100 ml/min.

**WARNING!** Eliminate as much air as possible from the blood compartment.

- Attach the dialysate tubing for counter current flow.
- Activate the dialysate flow and remove air from the dialysate compartment of the dialyzer.
- The dialyzer is now ready for connecting to the patient!

### **Treatment instructions**

- Follow instructions in the operator's manual of the dialysis equipment in use.

**CAUTION!** Anticoagulation is recommended to prevent clotting in the extracorporeal circuit. Anticoagulation therapy should be selected adjusted to the needs of the patient and monitored under the direction of a physician.

If Heparin is used, it is recommended to administer a loading dose to the patient 2 to 5 minutes before start of the treatment.

- Connect the arterial blood line to the vascular access of the patient.
- Connect the venous line to the vascular access of the patient.
- Start the blood pump and increase the blood flow rate to the prescribed value.

To avoid haemoconcentration in the dialyzer, do not start ultrafiltration before desired blood flow rate is reached.

### **Termination of Treatment**

- Follow instructions in the operator's manual of the dialysis equipment in use.
- Connect the arterial blood line to a bag of isotonic saline solution and adjust the blood flow rate to return the blood to the patient.

This does not apply if online-fluid is used for rinse back (e.g. with Artis and ULTRA prime line).

**WARNING!** Do not deactivate the air detection system on the machine until the patient is disconnected.

- When the required amount of blood has been returned to the patient, stop the blood pump and disconnect the patient.
- Follow the instructions in the operator's manual of the dialysis equipment in use how to empty and disconnect the dialyzer and blood lines.
- After use, this product may be a potential hazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

• Organization authorized to accept claims on the territory of the Russian Federation –  
An authorized representative of the manufacturer in the Russian Federation:

AO Company "Baxter"  
Leningradskoe shosse, 16A, Building 1  
125171, Moscow, Russia.  
tel .: +7 (495) 647-68-07,  
Fax: +7 (495) 647-68-08  
[russia@baxter.com](mailto:russia@baxter.com)

**Delivery set**

Consumer box contains:

1. Device as (individually in primary packaging), 24 pcs.
2. Instructions for use, 1 pc.

## SPECIFICATIONS

| Model                          | Theranova 400          | Theranova 500          |
|--------------------------------|------------------------|------------------------|
| Surface Area (m <sup>2</sup> ) | 1,7±0,1 m <sup>2</sup> | 2,0±0,1 m <sup>2</sup> |
| Blood compartment volume       | 91±6 ml                | 105±7 ml               |
| Residual blood volume          | <1 ml                  | <1 ml                  |
| Maximum TMP                    | 600 mmHg               | 600 mmHg               |
| Minimum TMP                    | 200 mmHg               | 200 mmHg               |
| Flow rate                      |                        |                        |
| Blood                          | 200–600 ml/min         | 200–600 ml/min         |
| Dialysate                      | 300–800 ml/min         | 300–800 ml/min         |
| Weight                         | 170 g (±10%)           | 185 g (±10%)           |

### Hemodialysis clearance

Permissible deviations from the nominal value:

± 10% for urea, phosphate, creatinine, vitamin B12 and inulin

± 20% for Cytochrome C

± 30% for myoglobin

### Theranova 400

QD = 300 ml/min, UF = 0 ml/min

| Q <sub>B</sub> [ml/min]           | 200 | 300 | 400 | 500 | 600 |
|-----------------------------------|-----|-----|-----|-----|-----|
| Urea (60 Da)                      | 191 | 246 | 272 | 285 | 291 |
| Phosphate (95 Da)                 | 179 | 225 | 250 | 266 | 276 |
| Creatinine (113 Da)               | 184 | 232 | 258 | 273 | 282 |
| Vitamin B <sub>12</sub> (1,4 kDa) | 148 | 178 | 199 | 214 | 226 |
| Inulin (5,2 kDa)                  | 119 | 140 | 156 | 169 | 180 |
| Cytochrome C (12 kDa)             | 109 | 128 | 142 | 153 | 164 |
| Myoglobin (17 kDa)                | 93  | 108 | 119 | 129 | 138 |

QD = 500 ml/min, UF = 0 ml/min

| Q <sub>B</sub> [ml/min]           | 200 | 300 | 400 | 500 | 600 |
|-----------------------------------|-----|-----|-----|-----|-----|
| Urea (60 Da)                      | 198 | 282 | 344 | 388 | 418 |
| Phosphate (95 Da)                 | 192 | 261 | 311 | 348 | 376 |
| Creatinine (113 Da)               | 194 | 269 | 323 | 362 | 391 |
| Vitamin B <sub>12</sub> (1,4 kDa) | 164 | 207 | 239 | 264 | 285 |
| Inulin (5,2 kDa)                  | 133 | 161 | 183 | 200 | 216 |
| Cytochrome C (12 kDa)             | 122 | 146 | 165 | 180 | 194 |
| Myoglobin (17 kDa)                | 104 | 123 | 137 | 150 | 161 |

$Q_D = 800 \text{ ml/min}$ ,  $UF = 0 \text{ ml/min}$

| $Q_B \text{ [ml/min]}$ | 200 | 300 | 400 | 500 | 600 |
|------------------------|-----|-----|-----|-----|-----|
| Urea (60 Da)           | 199 | 293 | 376 | 445 | 502 |
| Phosphate (95 Da)      | 196 | 279 | 345 | 400 | 446 |
| Creatinine (113 Da)    | 198 | 285 | 357 | 416 | 465 |
| Vitamin B12 (1,4 kDa)  | 174 | 227 | 267 | 301 | 329 |
| Inulin (5,2 kDa)       | 144 | 178 | 204 | 225 | 245 |
| Cytochrome C (12 kDa)  | 133 | 161 | 183 | 202 | 219 |
| Myoglobin (17 kDa)     | 114 | 135 | 152 | 166 | 180 |

### Theranova 500

$Q_D = 300 \text{ ml/min}$ ,  $UF = 0 \text{ ml/min}$

| $Q_B \text{ [ml/min]}$ | 200 | 300 | 400 | 500 | 600 |
|------------------------|-----|-----|-----|-----|-----|
| Urea                   | 192 | 250 | 276 | 288 | 294 |
| Phosphate              | 182 | 230 | 256 | 271 | 281 |
| Creatinine             | 186 | 237 | 263 | 278 | 286 |
| Vitamin B12            | 152 | 185 | 206 | 222 | 235 |
| Inulin                 | 124 | 147 | 164 | 178 | 189 |
| Cytochrome C           | 114 | 134 | 150 | 162 | 173 |
| Myoglobin              | 98  | 114 | 127 | 138 | 148 |

$Q_D = 500 \text{ ml/min}$ ,  $UF = 0 \text{ ml/min}$

| $Q_B \text{ [ml/min]}$ | 200 | 300 | 400 | 500 | 600 |
|------------------------|-----|-----|-----|-----|-----|
| Urea                   | 199 | 285 | 351 | 397 | 428 |
| Phosphate              | 194 | 267 | 320 | 358 | 388 |
| Creatinine             | 196 | 274 | 331 | 372 | 402 |
| Vitamin B12            | 169 | 215 | 249 | 277 | 299 |
| Inulin                 | 139 | 170 | 193 | 213 | 230 |
| Cytochrome C           | 128 | 155 | 175 | 192 | 208 |
| Myoglobin              | 110 | 130 | 147 | 161 | 173 |

$Q_D = 800 \text{ ml/min}$ ,  $UF = 0 \text{ ml/min}$

| $Q_B \text{ [ml/min]}$ | 200 | 300 | 400 | 500 | 600 |
|------------------------|-----|-----|-----|-----|-----|
| Urea                   | 200 | 295 | 381 | 454 | 515 |
| Phosphate              | 197 | 283 | 354 | 413 | 462 |
| Creatinine             | 199 | 288 | 365 | 428 | 481 |
| Vitamin B12            | 178 | 236 | 280 | 317 | 348 |
| Inulin                 | 150 | 188 | 216 | 241 | 262 |
| Cytochrome C           | 139 | 171 | 196 | 217 | 236 |
| Myoglobin              | 120 | 144 | 163 | 180 | 195 |

**Compartment pressure drop**

| Model                               | Theranova 400                              | Theranova 500                              |
|-------------------------------------|--------------------------------------------|--------------------------------------------|
| Blood compartment pressure drop     | Reducing the pressure in the blood circuit | Reducing the pressure in the blood circuit |
| QB = 200 ml/min                     | ≤ 90 mmHg                                  | ≤ 80 mmHg                                  |
| QB = 300 ml/min                     | ≤ 130 mmHg                                 | ≤ 120 mmHg                                 |
| QB = 400 ml/min                     | ≤ 170 mmHg                                 | ≤ 160 mmHg                                 |
| QB = 500 ml/min                     | ≤ 210 mmHg                                 | ≤ 200 mmHg                                 |
| QB = 600 ml/min                     | ≤ 250 mmHg                                 | ≤ 240 mmHg                                 |
| Dialysate compartment pressure drop |                                            |                                            |
| QD = 300 ml/min                     | ≤ 20 mmHg                                  | ≤ 15 mmHg                                  |
| QD = 500 ml/min                     | ≤ 30 mmHg                                  | ≤ 25 mmHg                                  |
| QD = 800 ml/min                     | ≤ 50 mmHg                                  | ≤ 40 mmHg                                  |

**Direct or indirect blood contact material**

| Component        | Material                                                 |
|------------------|----------------------------------------------------------|
| Membrane         | Polyarylethersulfone (PAES) / Polyvinylpyrrolidone (PVP) |
| Potting          | Polyurethane (PUR)                                       |
| Housing / header | Polycarbonate (PC)                                       |
| Sealing ring     | Silicone rubber (SIR)                                    |
| Supporting Ring  | Polypropylene (PP)                                       |

**CONDITIONS OF TRANSPORT, STORAGE AND USE****Conditions of operation**

|                     | Conditions of operation                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Temperature         | In accordance with the operating conditions of the dialysis machine (not included in the scope of delivery): from +18 °C till +35 °C |
| Relative Humidity   | In accordance <b>with</b> the operating conditions of the dialysis machine (not included in the scope of delivery): 30-85%           |
| Atmosphere pressure | In accordance with the operating conditions of the dialysis machine (not included in the scope of delivery): 795-1060 hPa            |

**Transportation and storage conditions**

The TheraNova series dialer is transported by all modes of transport in covered vehicles in accordance with the requirements and rules for the transportation of loads acting on each **mode of transport**.

|                     | Storage           | Transportation    |
|---------------------|-------------------|-------------------|
| Temperature         | from 0 to + 30 °C | from 0 to + 30 °C |
| Relative Humidity   | from 30 to 80%    | from 30 to 80%    |
| Atmosphere pressure | 79.5-106.0 kPa    | 79.5-106.0 kPa    |

## STERILIZATION

Dialyzer is supplied sterile in a sealed package. Sterilization is performed by steam. The sterilization process was validated according to the requirements of the standard EN ISO 17665-1: 2006.

Dialyzers are intended for single use only. Re-sterilization and re-use is prohibited.

## SHELF LIFE (MAINTAINING STERILITY)

Term of keeping sterility (shelf life) - 3 years. After the expiry of preservation of sterility, the dialyzer must be disposed of in accordance with the requirements of the local authorities.

A dialyzer with expired and/or damaged packaging is disposed of as medical waste of class A according to SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste".

## DISPOSAL

After use, this product may be a potential hazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations as class a waste, according to SanPiN 2.1.7.2790-10.

Do not dispose of the dialyzer after use with normal household waste, as they pose a potential biological hazard due to blood contamination. Recycling should be carried out in accordance with the norms adopted in the medical practice, as well as the relevant local, state and Federal laws and regulations, as waste class B on SanPiN 2.1.7.2790-10 "Sanitary and epidemiological re-quirements for the treatment of medical waste."


**Manufacturer:** «Gambro Dialysatoren GmbH», Holger-Crafoord-Strasse 26, 72379 Hechingen, Germany

T +49 7471 17-0, F +49 7471 17-1157

[www.gambro.com](http://www.gambro.com)

## EXPLANATION OF SYMBOLS

|                                                                                     |                                                   |
|-------------------------------------------------------------------------------------|---------------------------------------------------|
|  | Manufacturer                                      |
|  | Sterile, sterilized by steam (sterilization date) |
|  | Catalogue number                                  |
|  | Lot number                                        |
|  | Do not reuse                                      |
|  | Do not use with damaged packaging                 |
|  | Use by                                            |
|  | Caution! Refer to accompanying documentation!     |
|  | Upper limit of Temperature                        |
|  | Pyrogen-free                                      |

Do not use dialyzers to HDF (hemodiafiltration) or HF (hemofiltration)

Conformity symbol: The product complies with the essential requirements of the European Medical Product Directive 93/42 / EEC with the identification number of the notified body that has supplied the EC marking

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