ACKNOWLEDGEMENT


State of Utah)

§

County of Salt Lake)

On this 20th day of February, 2018, before me, Jennifer Wojciechowski, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.

A handwritten signature in blue ink, appearing to be "JWojciechowski", written over a horizontal line.

Notary Public

My Commission Expires: March 6, 2019

Dan Clark
Director, Regulatory Affairs
Merit Medical Systems, Inc.
20 February 2018

**EXTRACT FROM TECHNICAL DOCUMENTATION FOR
MEDICAL DEVICE**

« Kits of The ReSolve Biliary Locking Drainage Catheter »

Manufactured by

Merit Medical Systems, Inc., the USA

Redaction 2

Manufactured by:

Merit Medical Systems, Inc.

1600 West Merit Parkway South Jordan, Utah, 84095, USA

The place of production of the medical device

1) Merit Medical Systems, Inc.

1600 West Merit Parkway

South Jordan, Utah, 84095, USA

2) Merit Maquiladora Mexico, S. DE R.L. DE C.V., Avenida Sor Juana Ines de la Cruz
19970 Interior B, Edificio 2, Parque Industrial Frontera Tijuana, Baja California, MEXICO 22630

3) Merit Medical Systems, Inc.

14646 Kirby Drive Houston, TX 77047

The authorized representative of the manufacturer on the territory of the Russian Federation:

"Merit technologies"

Of the Russian Federation, 115054, Moscow Bakhrushina street 32, building 1

phone: +7 (495) 221-89-02

1. DEVICE DESCRIPTION

Kits of The ReSolve Biliary Locking Drainage Catheter designed for biliary drainage, namely procedures in order to ensure drainage of the clogged biliary system through the biliary drainage catheter placed through the skin into the bile duct.

The products are designed for use in clinics, hospitals and other medical institutions.

Kits of The ReSolve Biliary Locking Drainage Catheter is packaged with the following components:

- 1) Kit of the ReSolve RBC Biliary Locking Drainage Catheter, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
 - RBC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
- 2) Kit of the ReSolve RBDC Biliary Locking Drainage Catheter, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
 - RBC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
- 3) Kit of the ReSolve RBC Biliary Locking Drainage Catheter Plus, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
 - RBC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
 - Revolution catheter securement device – 1 pc
- 4) Kit of the ReSolve RBDC Biliary Locking Drainage Catheter Plus, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
 - RBDC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc

- Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
 - StayFix catheter securement device – 1 pc
- 5) Kit of the ReSolve RBC Biliary Locking Drainage Catheter Plus, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
- RBC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
 - StayFix catheter securement device – 1 pc
- 6) Kit of the ReSolve RBDC Biliary Locking Drainage Catheter Plus, the catheter diameter 8,5 F (2,8 mm); 10 F (3,3 mm); 12 F (4,0 mm); 14 F (4,6 mm); length 40 cm, consisting of:
- RBDC catheter ReSolve Biliary Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula – 1 pc
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
 - Revolution catheter securement device – 1 pc

The ReSolve Biliary Locking Drainage Catheter is made of single lumen tubing with 17 or 18 drainage holes depending on French size with one of the drainage holes also serving as a suture hole. The catheter is available in 8.5, 10, 12 and 14 French sizes with depth markers placed every centimeter beginning one (1) cm from the distal catheter tip and ending at the hub. The 8.5 and 10 French catheter configurations have 17 drainage holes and the 12 and 14 French configurations have 18 drainage holes as the latter have larger shaft diameters that allows for an additional drainage hole. The RBC (long) configuration has a lesion gap (distance between the proximal pigtail and the distal shaft hole) of 5cm. The RBDC (standard) has a lesion gap of 2 cm. This provides two (2) possible options of drainage hole configurations for the clinician to choose from based on patient anatomy. Both RBC and RBDC configurations are offered in the same French sizes (8.5, 10, 12, and 14) with a standard 40 cm working length.

A single radiopaque marker band is located adjacent to the most proximal drainage hole and helps in accurate placement of the drainage holes within the biliary duct under direct fluoroscopic guidance.


Figure 1 – RESOLVE BILIARY CATHETER

In all sets, biliary drainage catheters ReSolve includes the following auxiliary components: a metal cannula; a flexible cannula; the device repositioning (Repositioning Tool) to unlock the locking mechanism biliary drainage catheters; rectifier loop of the "pigtail" (Pigtail Straightener), to facilitate the introduction of the rigidity of the cannula, and end cap (Deadend cap) (see figure 1)

Sets, biliary drainage catheters ReSolve can include the locking device of the catheter or Revolution StayFix designed to secure the catheter and prevent movement and accidental removal of the catheter.

The fixing device Revolution catheter includes a filament for stabilizing and securing the catheter in place while leaving the injection site open, and has a special cover to provide additional protection against ingress of water (see figure 2,3).


Figure 2 - StayFix


Figure 3 - Revolution

2. INTENDED USE / INDICATIONS FOR USE, CONTRAINDICATIONS.

Indications:

The ReSolve Biliary Drainage Catheter with locking pigtail and hydrophilic coating is used for drainage of bile within the biliary system.

Contraindications:

The ReSolve Biliary Drainage Catheter is contraindicated for use where percutaneous biliary drainage catheterization is unacceptable. The ReSolve Biliary Drainage Catheter is contraindicated for intravascular use.

Cautions:

Where long-term use is indicated, it is recommended that indwelling time not exceed 90 days. This catheter should be evaluated by the physician on or before 90 days post-placement. To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

Warnings:

- RX Only Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Contents of unopened, undamaged package are sterile and non-pyrogenic
- This device is intended for single use only
- DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating
- DO NOT over tighten the connection between the drainage catheter and drainage tubing or dead end cap • If using alcohol to clean the catheter hub, allow sufficient time for alcohol to dry before connecting the drainage tubing or dead end cap
- Do not modify the catheter as modification may damage the catheter REUSE PRECAUTION STATEMENT For single patient use only. Do not reuse, reprocess or resterilize. Reuse,

reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Potential Complications:

- Hemobilia
- Cholangitis
- Pancreatitis
- Pneumothorax
- Hemothorax
- Bilithorax
- Sepsis
- Biloma
- Fever
- Catheter leak
- Catheter occlusion
- Catheter dislodgement
- Hemorrhage
- Bile peritonitis
- Skin Infection
- Perforation of the bile ducts, liver and/or duodenum

Stay Fix

INTENDED USE

The StayFIX Fixation Device is intended to secure non-vascular catheters and tubes.

INDICATIONS FOR USE The StayFIX Fixation Device is indicated to adhere to the skin and secure a percutaneous non-vascular catheter or tube to help minimize movement and/or accidental removal.

Contraindications

- Use on patients with known tape or adhesive allergies
- Use with vascular catheters

RESTRICTED DEVICE

Rx Only: Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device.

Precautions:

- Read manufacturer's instructions prior to use.
- Contents are sterile.
- Do not use if packaging is opened, damaged, or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

- Do not use after expiration date.
- Store in a cool, dry place.
- After seven (7) days (sooner if the fixation device becomes loose for any reason) remove the old StayFIX Fixation Device and replace with a new StayFIX Fixation Device
- When removing the old StayFIX Fixation Device, use care not to move or accidentally pull out catheter or tube.

REVOLUTION

Catheter securement device

INDICATIONS FOR USE: The Revolution catheter securement device is designed to secure a non-vascular catheter placed in a patient to help minimize catheter movement and accidental removal.

CONTRAINDICATIONS:

The Revolution catheter securement device is not intended to be used on vascular catheters.

CAUTIONS:

- Read instructions prior to use.
- Rx Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Contents of unopened, undamaged package are sterile.

REUSE PRECAUTION STATEMENT For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

3. Classification

Class IIb

Rule 8

According to Annex II of the Medical Device Directive 93/42/EEC

4. Sterilization Method

Merit Medical Systems, Inc. utilizes Ethylene Oxide (EO) sterilization for all ReSolve product configurations. The EO sterilization cycle is validated in accordance to the requirements of The European International Standard ISO 11135. Annual re-validation occurs assessing a worst-case reference load based on density and boxing configuration, which allows routine sterilization of a mixed-pallet load. Merit's products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals.

The protocols and reports validating the sterilization process are maintained in Merit's Document Control system and are referenced on Merit's Sterilization Summary (Table 942).

This device is not intended to be re-sterilized by the user.

ReSolve Catheter product is less complicated than a standard process challenge device (PCD) that represents a worst case configuration and is readily adopted into Merit's validated sterilization cycle.

Revolution and StayFix or are sterilized using gamma radiation and have passed the validation process. Validation was performed according to EN ISO 11137-2 - Sterilization of health care products. Exposure. Part 2. Establishing a sterilizing dose. On the basis of this check the minimum dosage for sterilization is 25 kGy.

The standard was verified dosages and verification of bioburden in accordance with EN ISO 11137-2 to verify that the dose for sterilization continues to achieve the desired SAL (10^{-6}).

5. DEVICE SPECIFICATIONS

5.1 ReSolve Biliary Locking Drainage Catheter pigtail

Characteristic	8.5 F	10 F	12 F	14 F
Catheter OD	0,112" ± 0,002" (2,84 mm ± 0,051 mm)	0,131" ± 0,002" (3,33 mm ± 0,051 mm)	0,158" ± 0,002" (4,01 mm ± 0,051 mm)	0,184" ± 0,002" (4,67mm ± 0,051 mm)
Catheter ID	0,070" - 0,073" (1,78mm - 1,85 mm)	0,089" ± 0,002" (2,26 mm ± 0,051 mm)	0,106" ± 0,002" (2,69 mm ± 0,051 mm)	0,118" ± 0,002" (3 mm ± 0,051 mm)
Catheter working Length (cm)	40±1	40±1	40±1	40±1
The distance between the loop-lock type pigtail (pigtail) and the hole of shaft (cm)	5±0,5 (RBC)	2±0,5 (RBD C)	5±0,5 (RBC)	2±0,5 (RBD C)
Pigtail Loop Size (cm)	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,794±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	3,175±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.
Number of Drainage Holes	17	17	18	18
Distance between handles depth (cm)	1±0,05	1±0,05	1±0,05	1±0,05
The distance between the loop-lock type pigtail (pigtail) and a radiopaque marker band (cm)	14 ± 0,5	14 ± 0,5	14 ± 0,5	14 ± 0,5
The strength of the thread	Not less than 6 N (1,371 lb)			
Suture Latch (Rotary cap to hub)	Not less than 10 N (2.25 lbs)			
The hardness of the catheter along the Shore	100A±0,5A			

Characteristic	8.5 F	10 F	12 F	14 F
Strength of the thread in a simple knot	Not less than 10 N (2.25 lbs)			
Strength of thread fastening	Not less than 10 N (2.25 lbs)			
Radiographic visibility	The catheter should be viewed under fluoroscopy			
The tensile strength of the catheter	Not less than 40 N (9.0 lbs)			
The tensile strength of the connection, "catheter-tip"	Not less than 6.5 N (1.48 lbs)			
The flow rate of the fluid	At least 5 CC /min.			
Compatibility guidewire	guidewire of 0.038 inches (0.97 mm) with an installed cannula needs to pass through the entire structure in the Assembly.			
Weight of the product	6,5 ± 2g	7,0± 2 g	7,2 ± 2 g	8,4± 2g
Compatibility connection, Luer-type	The hub of the catheter must be compatible with swivel and fixed adapters items covered by type (hemostatic valves, Luer lock connection, syringes, catheter valves, etc.).			
A cone-shaped tip	At the distal end of the tip should be tapered.			
Safety of magnetic resonance imaging (MRI)	The catheter should not interfere with the MRI.			
Surface	On the outer surface should not be technology and surface defects. This surface should be minimally traumatic to the vessel during use.			
Outer diameter connection, Luer-type	6,60 mm ± 0,5 mm			
The inner diameter of the connection, Luer-type	4,0 mm ± 0,5 mm			
The taper connection of Luer-type	6%			
Dimensions Repositioning Tool (length x width x height)	(50 mm x 12 mmx 25 mm) ± 3 mm			
Length Pigtail Straightener	100 m m ± 5 mm			
Outer diameter Pigtail Straightener	6 mm ± 2 mm			
The inner diameter Pigtail Straightener	4 mm ± 2 mm			
Length Dead end cap	13 mm± 2 mm			
Outer diameter Dead end cap	11 mm ± 2 mm			
Inner diameter Dead end cap	8 mm ± 2 mm			

5.2 Metal Stiffening Cannula

Characteristic	8.5 F	10 F	12 F	14 F
----------------	-------	------	------	------

Cannula Size	18	15	14	13
Cannula OD	0,050" ± 0,0005" (1,27 mm ±0,013 mm)	0,072" ± 0,0005" (1,83 mm±0,013 mm)	0,083" ± 0,001" (2,11 mm ±0,013 mm)	0,095" ± 0,001" (2,41 mm±0,013 mm)
Cannula ID	0,033" ± 0,0005" (0,84 mm ± 0,013 mm)	0,054" ± 0,0005" (1,372 mm ± 0,013 mm)	0,063" ± 0,0005" (1,60 mm ±0,013 mm)	0,071" ± 0,0005" (1,80 mm ±0,013 mm)
Exposed Length	19,452" ± 0,0395" (494 mm±1 mm)	19,442" ± 0,0395" (493,8 mm±1 mm)	19,362" ± 0,0395" (491,8 mm±1 mm)	19,342" ± 0,0395" (491,2 mm±1 mm)
Joint Strength	Not less than 10 N (2.25 lbs)			
Fracture strength	Not less than 10 N (2.25 lbs)			
Rockwell hardness	53-55 HRC			
Roughness Ra	≤ 0,4 μm			
Allowable angle of deflection, without compromising the integrity	15° ± 1°			
Surface	On the outer surface should not be technology and surface defects.			
Weight	2,05± 1 g	2,7± 1 g	3,5± 1 g	5,05± 1 g

5.3 Flexible (Plastic) Stiffening Cannula

Characteristic	8.5 F	10 F	12 F	14 F
Cannula OD	0,052" ± 0,002" (1,32 mm ±0,05 mm)	0,068" ± 0,002" (1,73 mm ±0,05 mm)	0,080" ± 0,002" (2,03 mm ±0,05 mm)	0,091" ± 0,002" (2,31 mm ±0,05 mm)
Cannula ID	0,040" ± 0,002" (1,02 mm ±0,05 mm)	0,051" ± 0,002" (1,3 mm ±0,05 mm)	0,063" ± 0,002" (1,6 mm ±0,05 mm)	0,071" ± 0,002" (1,8 mm ±0,05 mm)
Exposed Length	19,295" ± 0,0395" (490 mm±1 mm)	19,285" ± 0,0395" (489,8 mm±1 mm)	19,205" ± 0,0395" (487,8 mm±1 mm)	19,185" ± 0,0395" (487,3 mm±1 mm)
The tensile strength of the connection with the hub (strength of fixing the connector of a flexible cannula with a tube of a flexible cannula)	Not less than 10 N (2.25 lbs)			
Tensile strength	Not less than 10 N (2.25 lbs)			
Allowable angle of deflection, without compromising the integrity	15° ± 1°			

Characteristic	8.5 F	10 F	12 F	14 F
Surface	On the outer surface should not be technology and surface defects.			
Weight	1,3± 1 g	1,5± 1 g	1,5± 1 g	1,7± 1 g

5.4 StayFix

StayFix	Small	Large
Total area, cm ²	120 ± 2	120 ± 2
The area of the adhesive layer, cm ²	88± 2	90± 2
The peel force of the adhesive layer	Not less than 2N	Not less than 2N
Hydrocolloid bandage (length x width) in mm	(60 x 50) ± 5	(60 x 50) ± 5
Weight	6.2 ± 0.5 g	6.2 ± 0.5 g
catheter compatibility	For fixing the catheter outer diameter of 5F-14F (1,67 mm - 4.67 mm)	For fixing the catheter outside diameter 12F-22F (4,0 mm - 7.3 mm)

5.5 Revolution

Total area, cm ²	111 ± 2
The area of the adhesive layer, cm ²	95 ± 2
The outer diameter of the protective rings, mm	67 ± 5
The inner diameter of the protective rings, mm	40 ± 5
The maximum diameter of the lid, mm	70 ± 2
Weight	19.9 ± 0.5 g
Locking thread (length x width x height) in mm	(12 x 12 x 4) ± 2
The ratchet (length x width x height) in mm	(60 x 30 x 5) ± 2
Force of separation thread	Not less than 5H
The peel force of the adhesive layer	Not less than 2H
catheter compatibility	For fixing the catheter outer diameter of 6F-24F (2,0 mm - 8.0 mm)

5.6 Requirements for packaging products (main characteristics) that is subjected to finish sterilization.

5.6.1 ReSolve Biliary Locking Drainage Catheter pigtail with a hydrophilic coating

The width of the package, mm	133 mm $\pm$ 5
The Length of the package, mm	606 mm $\pm$ 5
The width of the adhesive layer, mm	8.0 $\pm$ 0.5 mm
The strength of the adhesive layer on the gap	$\geq$ 1 lbf (4.45 N)
The density of packing material	0.92 $\pm$ 0.05 g / cm ³
The tensile strength in the longitudinal direction	> 0.7 kN / cm ²
The tensile strength in the transverse direction	> 0.9 kN / cm ²
Bursting strength	> 650 kPa

5.6.2 StayFix

The width of the package, mm	143 $\pm$ 3 mm
The Length of the package, mm	198 $\pm$ 3 mm
The width of the adhesive layer, mm	7.3 $\pm$ 0.5 mm
The strength of the adhesive layer on the gap	$\geq$ 1 lbf (4.45 N)
The density of packing material	0.92 $\pm$ 0.05 g / cm ³
The tensile strength in the longitudinal direction	> 0.7 kN / cm ²
The tensile strength in the transverse direction	> 0.9 kN / cm ²
Bursting strength	> 650 kPa

5.6.3 Revolution

The width of the package, mm	140 $\pm$ 3 mm
The Length of the package, mm	215 $\pm$ 3 mm
The width of the adhesive layer, mm	9.0 $\pm$ 0.5 mm
The strength of the adhesive layer on the gap	$\geq$ 1 lbf (4.45 N)
The density of packing material	0.92 $\pm$ 0.05 g / cm ³
The tensile strength in the longitudinal direction	> 0.7 kN / cm ²
The tensile strength in the transverse direction	> 0.9 kN / cm ²

Bursting strength	> 650 kPa
-------------------	-----------

6. MATERIALS OF CONSTRUCTION

Component	Material Description
Drainage Catheter	
Hub	Polyester polycarbonate blend Xenoy HX5600HP Dye - solid fraction Xenoy EXN X5600HP titanium oxide, color: bright white, code classification Pantone 11-0601 TCX
Rotary Cap	Polyester polycarbonate blend Xenoy HX5600HP Dye - solid fraction Xenoy EXN X5600HP titanium oxide, color: bright white, code classification Pantone 11-0601 TCX
Rotary Cam Lock	<i>Medium viscosity acetal homopolymer Delrin® 500P NC010 POM DuPont™ (For 8,5F diameter)</i> Pigment - titanium solid fraction Xenoy EXN X5600HP oxide Colorant - the particulate aluminum oxide, color: green, classification 354S code Pantone <i>Medium viscosity acetal homopolymer Delrin® 500P NC010 POM DuPont™ (For 10F diameter)</i> Colorant - the solid fraction Xenoy EXN X5600HP oxide Colorant - alumina powder color: Purple, classification 272S code Pantone <i>Medium viscosity acetal homopolymer Delrin® 500P NC010 POM DuPont™ (For 12F diameter)</i> Colorant - the solid fraction Xenoy EXN X5600HP oxide Pigment - Aluminum oxide powder Color: Yellow, Pantone 109C classification code <i>Medium viscosity acetal homopolymer Delrin® 500P NC010 POM DuPont™ (14F to diameter)</i> Pigment - titanium oxide solid fraction Xenoy EXN X5600HP, Pigment - Aluminum oxide powder Color: Orange, on the classification of Pantone code 021S Red ink INK-PG Pad TYPE UV 2002 ICT Trans Tech Color: Red, Classification Pantone code 485S (for 12F diameter) White ink INK-PG Pad TYPE UV 2002 ICT Trans Tech Color: white, code classification Pantone WHITE C (for 6,5F diameter - 10F, 14F)
	The hydrophilic urethane ester-based HydroMed™ D640
Catheter shaft, all sizes	thermoplastic polyurethane elastomer (TA) Pellethane 2363 80AE + 25% barium sulfate (BaSO4)
Catheter Shaft Tip Plug (8.5F RBC and RBC/MB sizes only)	thermoplastic polyurethane elastomer (TA) Pellethane Shore 65D, RO20%
Marker Band	Platinum 92% / 8% iridium

Depth Markings	Black ink INK-PG Pad PLT15 ICT Trans Tech. Color: black, code classification Pantone Black C
Lubricants	Highly fluoro-silicone grease lubrication MED-6731, Polydimethylsiloxane liquid lubricant DOW CORNING® 360 Medical Fluid, 12500 CST.
Pigtail Straightener	polypropylene 5317 Pigment - titanium oxide solid fraction PP MEVOPUR® Clariant Color: light blue, according to Pantone classification code 13-4804 TPX
Suture Assembly	
Suture	polyamide 6 monofilament, 3/0
Latch Seal	Silicone Brand EXEMPT Parker Seal
Ferrule Crimp	stainless steel alloy AISI 304 (composition: C -0,08%, Mn - 2,0%, Si - 1,0%, Cr - 18,0 / 20,0%, Ni - 8,0 / 10,5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0,08%, Mn - 2,0%, Si - 1,0%, Cr - 16,0 / 18,0%, Ni - 10,0 / 14,0%, Mo - 2.0-3.0 %, Fe - 68.92-61.92%)
Suture Latch	Polyester polycarbonate blend Xenoy HX5600HP Dye - solids fraction of titanium oxide Xenoy EXN X5600HP The ratio of dye is not more than 6.25% Color: bright white, code classification Pantone 11-0601 TCX Red ink INK-PG Pad TYPE UV 2002 ICT Trans Tech Color: Red, Classification Pantone code 485S
Repositioning Tool	polycarbonate CALIBER™ mrad™ 2081 Pigment - titanium oxide solid fraction A32329C PolyOne™ Color: bright white, code classification Pantone 11-0601 TCX Red ink INK-PG Pad TYPE UV 2002 ICT Trans Tech Color: Red, Classification Pantone code 485S
Cap	Thermoplastic copolymer Lustran® 248 NR Color: white, code classification Pantone White
Metal Cannula	
Hub	polycarbonate CALIBER™ mrad™ 2081 Dye - 10200442 TRANS granules BLUE PolyOne™ Color: blue, Pantone 3005U Classification code (for 8,5F diameter) Dye - 10050274 TRANS beads PURPLE PolyOne™ Color: lilac, classification code Pantone 15-3817 (10F to diameter) Dye - polycarbonate granules PC 2081-15 Caliber™ Mrad™ Color: Yellow, Pantone 109C classification code (for 12F diameter)

	Dye - polycarbonate granules PC 020-0102 ORANGE Clariant Color: orange, Pantone 021S Classification code (for 14F diameter)
Cannula	stainless steel alloy AISI 304 (composition: C -0,08%, Mn - 2,0%, Si - 1,0%, Cr - 18,0 / 20,0%, Ni - 8,0 / 10,5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0,08%, Mn - 2,0%, Si - 1,0%, Cr - 16,0 / 18,0%, Ni - 10,0 / 14,0%, Mo - 2.0-3.0 %, Fe - 68.92-61.92%)
Coating	Medical Silicone Dow Corning (R) MDX4-4159 + 50% lubricant medical purpose
Flexible Stiffening Cannula	
Non-Swivel Nut	Copolyester Eastar (TM) DN003
extruded tube	polypropylene 5317
hub	polypropylene 5317 Dye - solids fraction of titanium oxide, aluminum oxide powder PE PE / PP MEVOPUR® Clariant Color: blue, Pantone 293C classification code (for 8,5F diameter) dye - solids fraction of titanium oxide PP PURPLE PolyOne™ Color: purple, Pantone 272S Classification code (for 10F diameter) Dye - solids fraction of titanium oxide, aluminum oxide powder PE PE / PP MEVOPUR® Clariant Color: Yellow, Pantone 108C classification code (for 12F diameter) Dye - solids fraction of titanium oxide, aluminum oxide powder PE PE / PP MEVOPUR® Clariant Color: Orange, classification code Pantone 021C (14F to diameter)
Revolution	
A thread	polyamide 6 monofilament, 3/0
The ring-mount	The CALIBER™ Polycarbonate Mrad™ 2081
Cap	copolymer of styrene with methyl methacrylate ZYLAR® 530
Ratchet	Polyethylene resin Zytel® 101F NC010 Dye - NY Zytel beads AB63631417 Clariant (composition: Nickel antimony titanium - 22.91%, a cobalt compound with titanium - 12.60% kaolin (adjuvant) - 2%). Color: green, the code classification of 101F
the pull ring	The CALIBER™ Polycarbonate Mrad™ 2081 Pigment - titanium oxide solid fraction PC 2081-15 Clariant. Color: light blue, code classification Pantone 292U
adhesive patch	Adhesive non-woven material c silicone coated Release Liners
StayFix	Adhesive non-woven material c silicone coated Release Liners hydrocolloid dressing

	Polyethylene resin Zytel® 101F NC010
--	--------------------------------------

** in the manufacture of this medical product, a barium sulfate chemical element is used. Barium sulphate acts as a constituent element in the manufacture of these products, in this case it is not a dosage form or substance. Barium sulphate is one of the main radiopaque materials that is added to polymer medical products. Barium sulphate does not have a pharmacological and metabolic effect on the body. The chemical formula is BaSO₄. Barium sulphate because of its low solubility is not toxic to the body and therefore its use as an X-ray contrast substance is possible. When using the listed products for the intended purpose, the ingestion of BaSO₄ into the human body is excluded. Barium sulphate is used as part of a medical device for visualizing a medical device during surgical interventions.*

All ReSolve products are packaged and labeled as sterile, single use devices. Sterile barrier systems are specified per ISO11607 to prevent ingress of microorganisms and allow aseptic presentation of product at point of use. Package integrity testing is performed as part of the production process.

Packaging Component	Material
Pouch	Tyvek
Tray	White Plastic Eastman Kodak 6763
Carton	Corrugate
Shipper	Corrugate

7. Compliance with Relevant Regulations and Technical Standards

The product is in compliance with the Council Directive 93/42/EEC, Medical Device Directive of the European Union. The applicable CE certificate is CE 541900.

The product is legally marketed according to the requirements of the United States Food and Drug Administration in compliance with the premarket notification 510(k) process. The 510(k) clearance letter is K141408.

The qualification activities that were performed for the product were based on standards that were in effect at the time of the testing. Standards applied include:

Standard	Title
General standards	
ISO 14971	Medical devices – Application of risk management to medical devices
ISO 10993-1	Biological evaluation of medical devices Part 1: Evaluation and testing within a risk management process
ISO 11135	Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of a sterilization process for medical devices
ISO 10993-3	Biological Evaluation of Medical Devices – Part 3: Tests for Genotoxicity, Carcinogenicity and Reproductive Toxicity.
ISO 10993-4	Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood.
ISO 10993-5	Biological Evaluation of Medical Devices - Part 5: Tests for Cytotoxicity – In Vitro Methods
ISO 10993-6	Biological Evaluation of Medical Devices - Part 6: Tests for Local Effects after Implantation.
ISO 10993-7	Biological Evaluation of Medical Devices - Part 7: Ethylene Oxide Sterilization Residuals

ISO 10993-10	Biological Evaluation of Medical Devices - Part 10: Tests for irritation and sensitization
ISO 10993-11	Biological Evaluation of Medical Devices – Part 11: Tests for systemic toxicity
ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems, and packaging systems.
ISO 11607-2	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ASTM D4169	Standard practice for performance testing of shipping containers and systems
ASTM F1980	Standard guide for accelerated aging of sterile barrier systems for medical devices
BS EN 980	Symbols for use in the labeling of medical devices
BS EN 1041	Information supplied by the manufacturer of medical devices
ISO 15223-1	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements
BS EN 15986	Symbols for use in labeling of medical devices – Requirements for labeling of medical devices containing phthalates
EN556	Sterilization of medical devices – Requirements for terminally sterilized medical devices to be labeled “Sterile”
Product Specific Standards	
ISO 592-2	Conical Fittings with a 6% (Luer) Taper for Syringes, Needles and Certain Other Medical Equipment - Part 2: Lock Fittings
EN 1617	Sterile Drainage Catheters and accessory devices for single use
EN 1618	Catheters other than Intravascular – Test methods for common properties

8. Risk Analysis System

Potential risks are evaluated during the design, product development and manufacturing processes according to ISO 14971. The appropriate risk assessment tool (FMEA, Clinical Risk Assessment - CRA.) is used to identify potential failure modes and potential effect on patients or users. Appropriate control and prevention mechanisms are defined and implemented as part of the risk management process throughout the product lifecycle. Risk Analysis reports are revisited and updated as necessary, and controlled under Merit's Document Control system. A cross-functional team, including persons with clinical experience, reviewed the Risk Analysis reports to ensure that all potential risks were addressed.

Table Risk Assessment

REFERENCE	TITLE
RMR 00004	Risk Management Report for ReSolve Biliary and Locking Drainage Catheter
FMEA 0010	Design FMEA for ReSolve NL Drainage Catheter
FMEA 0011	Process FMEA for ReSolve NL Drainage
FMEA 0074	Design FMEA for ReSolve Locking Catheter
FMEA 0075	Process FMEA for ReSolve Locking Catheter

CRA0045	Clinical Risk Analysis for the ReSolve Biliary Catheter
FMEA0187	Process FMEA for ReSolve Locking Catheter
FMEA 0210	ReSolve Non-Locking and Locking Catheter FMEA
FMEA 0206	ReSolve Biliary Drainage Catheters
QRMS0009	Quality Risk Management Summary for ReSolve Biliary Drainage Catheters

Important Hazards

The conclusions drawn by the cross-functional team are located in the Risk Assessment reports. The most important hazards include exposure to blood borne pathogens, infection, tissue trauma, toxic reaction, foreign body embolization, air ingress, tissue irritation and delayed procedure.

Manufacturer statement regarding safety

Based on the documentation provided and/or referenced in this technical file describing the design, manufacture, packaging and sterilization of the ReSolve products and the assessment of potential risks compared to the benefits of the use of the device through formal risk analysis, clinical evaluation, and product/process qualifications, Merit Medical Systems has determined that the Essential Requirements of the European Community Council Medical Device Directive 93/42/EEC have been satisfied and that the ReSolve Catheters are safe and effective for their intended use as labeled.

9. BIOCOMPATIBILITY

Drainage Catheters as a function of their original CE and FDA clearance.

Categorization of Device(s)

ReSolve Biliary Drainage Catheters are categorized below. Components of these devices with patient contact will be represented in the Biocompatibility Testing Summary below.

Table - Device Biocompatibility Categorization

DEVICE	CATEGORIZATION
ReSolve Biliary Catheter	Externally Communicating Device with Tissue Contact for a Permanent Duration (> 30 days).

A brief report on the results of tests on biocompatibility

Tests Performed to Evaluate the Biological Safety of the ReSolve Biliary Catheter per Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Stainless Steel - All samples were prepared according to ISO 10993-12 and were sterilized by ethylene oxide.		
Cytotoxicity/ ISO 10993-5/ MEM Elution	Non-Cytotoxic	Nelson Laboratories 466472
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizing	Toxikon GLP Report 10-2857-G1
Irritation/ ISO 10993-10/ Intracutaneous Injection	Negligible Irritant	Toxikon GLP Report 10-2857-G2
Purple Delrin – Sample preparation and sterilization information is unavailable.		
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution	Non-Cytotoxic	Nelson Laboratories Report 309477
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	AppTec Laboratories Report 37899
Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Northview Pacific Laboratories Report X5K063G
Systemic Toxicity (Acute)/ ISO 10993-11/ Systemic Injection	Negative	Northview Pacific Laboratories Report X5K062G
Genotoxicity/ No Standard Referenced / Salmonella Typhimurium and Escherichia Coli Reverse Mutation Assay	Non-Genotoxic	Nelson Laboratories Report 309478

Tests Performed to Evaluate the Biological Safety of the ReSolve Biliary Catheter per Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Implantation/ ISO 10993-6/ 7 Day Implantation	Negative	Northview Pacific Laboratories Report X5K060G
Implantation/ ISO 10993-6/ 90 Day Implantation	Negative	Northview Pacific Laboratories Report X5K059G
Pellethane , Suture , Hydrophilic Coating, Lubricant Polyethylene, Polyurethane Solvent, Platinum, Xenoy, White Colorant, White Printing Ink, , Hardener, Orange Delrin, Black Ink, Retarder, Thinner, Hardener, Red Ink– All samples were prepared in accordance with ISO 10993-12 and were sterilized by ethylene oxide prior to testing.		
Cytotoxicity/ ISO 10993-5/ L-929 MEM Elution	Non-Cytotoxic	Nelson Laboratories 841022
Sensitization/ ISO 10993-10/ Kligman Maximization Study	Non-Sensitizer	Wuxi AppTec GLP Report 192776
Irritation/ ISO 10993-10/ Intracutaneous Injection	Negligible Irritant	Wuxi AppTec GLP Report 192777
Acute Systemic Toxicity/ ISO 10993-11/ Acute Systemic Injection	Negative	Wuxi AppTec GLP Report 192778
Pyrogenicity/ USP <151>/ Material Mediated Pyrogenicity	Non-Pyrogenic	Wuxi AppTec GLP Report 192779
Genotoxicity/ ISO 10993-3/ Bacterial Mutagenicity Test (AMES Assay)	Non-Genotoxic	Wuxi AppTec GLP Report 192782
Implantation/ ISO 10993-6/ 2 Week Implantation Study	Negative	Wuxi AppTec GLP Report 192783
Implantation/ ISO 10993-6/ 4 Week Implantation Study	Negative	Wuxi AppTec GLP Report 192784
Implantation/ ISO 10993-6/ 13 Week Implantation Study Including Systemic Effects	Negative	Wuxi AppTec GLP Report 193893
Chemical Characterization/ USP <661>/ Physicochemical Tests for Plastics	Pass	Nelson Laboratories Report 721109
Tubing – All samples were prepared in accordance with ISO 10993-12 and were sterilized by ethylene oxide prior to testing.		

Tests Performed to Evaluate the Biological Safety of the ReSolve Biliary Catheter per Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution Test	Non-Cytotoxic	Toxikon 12-5727-G1
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	Toxikon 12-5727-G2
Irritation or Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Toxikon 12-5727-G3
Acute Systemic Toxicity/ ISO 10993-11/ Systemic Toxicity Test	Negative	Toxikon 12-5727-G4
Acute Systemic Toxicity/ ISO 10993-11, USP <151>/ Rabbit Pyrogen Test (Material Mediated)	Non-Pyrogenic	Toxikon 12-5727-G5
Genotoxicity/ ISO 10993-3/ Salmonella Typhimurium and Escherichia Coli Reverse Mutation Assay	Non-Genotoxic	Toxikon 12-5727-G6
Hemocompatibility/ ISO 10993-4/ Complement Activation	Meets Acceptance Criteria	Nelson Laboratories 666341
Hemocompatibility/ ASTM F756, ISO 10993-4/ Hemolysis – Direct and Indirect Contact	Non-Hemolytic	Toxikon 12-5727-G7
Hemocompatibility/ ISO 10993-4/ Thrombogenicity	Meets Acceptance Criteria	Toxikon 12-5727-G8
Chemical Characterization/ USP <661>/ USP Physicochemical Tests for Plastics	Pass	Nelson Laboratories 666340

Conclusion:

A biological safety assessment has been completed for the ReSolve Non-Locking, Locking, and Biliary Drainage Catheters. After review of the pertinent, objective evidence available, Merit Medical Systems, Inc. concludes that the biological and toxicological risks associated with these devices are low and acceptable.

Material mediated Pyrogenicity testing is performed in conjunction with biocompatibility tests during product design and development.

Bacterial Endotoxin Testing is conducted using the in-vitro Limulus Amebocyte Lysate (LAL) test method and meets the current U.S. Pharmacopeia (USP) limit of 20EU/device. This test is performed to support Merit Medical's label claim of 'non-pyrogenic.'

In addition, samples of assembled, packaged products are pulled weekly for LAL Endotoxin testing. LAL testing is not performed on every lot of product; however, the samples taken represent the product families manufactured during the week the samples were pulled.

10. STABILITY AND DURABILITY

Merit determines shelf life based on accelerated aging studies performed according to the Arrhenius Model for both the packaging and the product. The Arrhenius

model $[Q_{10} = 2^{\frac{T_2 - T_1}{10}}]$ quantifies the doubling of reaction rate for every 10°C increase in temperature; providing the number of weeks of aging needed for a material under specified conditions. It is widely accepted as a standard in American industry.

The Arrhenius model relates the variables of time and temperature to the deterioration of materials and therefore is useful in calculating accelerated thermal aging parameters. It is preferred for this purpose over the n-degree rule in that it has a firmer theoretical foundation and has been verified for many materials. (Reference: Electric Power Research Institute EPRI NP-1558; Project 890-01; Final Report September 1980. "A Review of Equipment Aging Theory and Technology")

Aging studies are followed by real-time aging of product. Merit's internal procedure, STP0002, is derived from the above standard and is utilized in the aging studies of the majority of Merit's products.

The RLC and RNL Drainage Catheters and the accessories are labeled for a 3 year shelf life. ReSolve Biliary catheters have a shelf life of 3 years.

11. MANUFACTURING PROCESS

A flowchart representing the manufacturing process is attached in Appendix A. All ReSolve products are manufactured in a suitably-controlled environment as described in the following standard operation procedures:

- GPS 999.031 Clean-room Environmental Action Requirements
- GPS 999.044 Controlled Environments
- QAP 999-950 Environmental Monitoring Procedure

Processes and equipment used for medical device manufacturing are validated according to established protocols and results are documented and approved. Records are maintained for validated/qualified processes, equipment and personnel, and include the work instructions used, inspection and testing data, the date on which the process was performed, the identity of persons who perform or monitor operations and, where appropriate, the equipment used. The flowchart included in Appendix A provides a visual layout of the manufacturing process. * StayFIX/ Revolution included in the compositions of the kits are produced and Packed at the production site: Merit Maquiladora Mexico S. DE R. L. DE C. V., Avenida Sor Juana Ines de la Cruz B 19970 Interior, Edificio 2, Parque Industrial Frontera Tijuana, Baja California, MEXICO 22630.

* Catheters and other related products are manufactured and packaged at production sites: Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah, 84095, USA; Merit Medical Systems, Inc. 14646 Kirby Drive Houston, TX 77047, the final Assembly of finished kits is carried out there. The country of manufacture is indicated on the marking.

Quality Control

Final Product Release Criteria

All Resolve product assemblies are inspected at various points during the manufacturing process according to parameters defined in Merit procedures and specifications referenced in the Device Master Record. (See Appendix A- Manufacturing Flow Chart(s)).

The ReSolve Catheter inspection requirements include the following:

Visual inspection includes:

- Contaminants and foreign material.
- Damage or defects caused by the manufacturing process.
- Correct product configuration and labelling. Verify all components are in the tray per the work order drawing.
- Kinks, bends or misalignment in the catheter, plastic stiffener, trocar or cannula relative to the hub.
- Cracked, broken, deformed or otherwise damaged components or parts.
- Compare the suture knot to the work order drawing. Reject if smashed, kinked, incomplete, and incorrect or damaged. Correct Merit logo orientation. Suture is wrapped around tip.
- Printed hub information is correct and is not missing incomplete, smeared, illegible or incorrectly placed or oriented.
- Echogenic surface must be present at the end of the metal cannula.
- Trocar and cannula bond for excessive bubbles in the bond area or burnt, melted, or severely discolored adhesive.
- The sharp tip of the trocar must not be damaged or dull.
- Verification that the catheter tip is not deformed, damaged, slanted, incompletely formed or incompletely trimmed.
- The tip of the catheter should not exhibit discoloration, cracks, bubbles, blisters, sinks or pits and the tip must not be plugged.
- Inspect for visible separation of marker band from tubing.
- Inspect printed cm marks for missing marks, illegible print, crooked or out of alignment.
- Verify the hydrophilic coating is uniform along the coated section or if coating is missing flaky, dirty or spotted.
- Reject for evidence of hydrophilic solution (cured/uncured) inside the tubing
- Punched holes for plugged or un-punched holes. Debris from punching that is not removed. Holes skewed or cut off centerline and into tube sidewall. Suture hole punched on hub side of marker band
- Correct drainage hole pattern and number of holes along the tubing shaft.
- Molded hub assembly for short shots, sinks, or voids or sharp edges.

- Molded component critical dimension inspection.
- Component fit.
- ID of the metal cannula at the hub end.
- Metal cannula length.
- Measurement of the trocar protrusion distance from the cannula. Reject if distance is not within specification per work order drawings.
- Verify proper placement and length of 5 cm dots on the catheter.
- Distance from tip edge of marker band to tip of catheter is within specification (Catheters with Marker Band only).
- Curved catheter fits within drawing specification tolerances.
- Location of drainage holes.
- Distance from nearest punched hole to marker band (Catheters with Marker Band Only).
- Plastic Stiffener length and tip inner diameter after placing radius.
- Catheter length.

Functional inspection includes:

- Tensile strength of the trocar to hub bond and tubing to hub connection.
- Pouch burst test.
- Catheter leak test.

Components/assemblies that do not meet specifications identified in procedures are dispositioned for rework or destroyed.

After the finished devices are released from the manufacturing process, they are packaged and sterilized. Sterilized product is placed in a post-sterilization quarantine area until the sterilization documentation is verified. Product is released for distribution by Quality Assurance personnel signature.

12. PACKAGING QUALIFICATION

For products "Merit Medical Systems, Inc." Acceptable standard warehousing, unless special conditions of storage, in accordance with the risk management reports are not listed on the packaging. The product is packaged and labeled as a sterile disposable article. The sterile barrier systems specified in the standard ISO 11607, are designed to prevent the ingress of microorganisms and ensures aseptic preservation of the product to the final application. Checking package integrity is carried out within the production process on the package. Sets biliary drainage catheters Resolve the shelf life of 3 years.

The sets of the Resolve biliary drainage catheter consist of the following components: one (1) ReSolve biliary drainage catheter (RBC or RBDC) together with a pigtail loop straightener and repositioning tool, one (1) metal cannula, one (1) flexible cannula (1) one end deadend cap is placed in an impact-resistant polystyrene tray with a PET lid to support the products, which is then placed in a bag of Tyvek® material sealed with a validated thermal welding process, designed to cover and Protecting the product during sterilization, loading, transportation, manipulation and storage.

The sets of the Resolve Plus biliary drainage catheter consist of the following components: one (1) ReSolve biliary drainage catheter (RBC or RBDC) together with a pigtail loop straightener and repositioning tool, one (1) metal cannula, one (1) flexible cannula, one (1) end deadend cap is

placed in an impact-resistant polystyrene tray with a PET lid, which is then placed in a bag of Tyvek® material sealed with a validated thermal welding process. The fixing device (StayFix or Revolution) is packaged in individual packages in the form of a packet of Tyvek® material, sealed with a validated thermal welding process.

Packages made of Tyvek® material sealed with a validated thermal welding process are designed to cover and protect the product during sterilization, loading, transportation, manipulation and storage.

Individual packages of the biliary drainage catheter and fixation devices are connected by an external fixing carton with a common label on it, on which the composition of the kit is presented.

Instructions for use are included in the group packaging of the kits (in the group pack, five sets of the biliary drainage catheter Resolve are laid). Sets of catheters in a group box are placed in a transport box.

An example of a group package is shown in the figure below.

ReSolve's packaging qualification was conducted to meet the requirements of the following standards:

ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems.

ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation Requirements for forming, sealing and assembly processes

ASTM D4169: Standard Practice for Performance Testing of Shipping Containers and Systems.

13. LABEL

Product label formats have been developed and approved for unit product, sterile packages and transport packaging. The approved labeling is located in the product Device Master Record.

Merit utilizes the following standards for labeling content and symbol generation in order to meet current requirements of the MDD:

EN 980: Symbols for use in the Labeling of Medical Devices.



EN 1041: Information Supplied by the Manufacturer.

ISO 15223 - 1: Medical Devices - Symbols to be used with medical device labels, labeling and information to be supplied Part 1: General Requirements.

ISO 15223 - 2: Medical Devices - Symbols to be used with medical devices, labeling, and information to be supplied Part 2: Symbol Development, Selection and Validation.

- FDA Guidance on Medical Device Patient Labeling; Final Guidance for Industry and FDA Reviews (4/19/2001).

Description of symbols labeling shown in the table.

	Symbol.	Explanation of symbols
1.		EC official
2.		manufacturer
3.		It sold only to doctors and hospitals
4.		For single use only. Do not reuse
5.		catalog number
6.		shelf life
7.		lot number
8.		sterilized ethyleneoxide
9.		Sterilized by gamma radiation
10.		attention
11.		CONDITIONAL compatibility with MR
12.		CE Mark Mark
13.		Pyrogen-free
14.		Protect from moisture
15.		Do not use if package is damaged
16.		This side up

17.		Fragile, handle with care
-----	---	---------------------------

Examples of marking are shown below:

ReSolve RBC Biliary Locking Drainage Cathete

ReSolve RBC Biliary Locking Drainage Cathete Plus

Example of transport packages

14. TRANSPORTATION, STORAGE AND HANDLING

The shelf life of sterile - 36 months from the date of sterilization.

Tools that are exempt from transport packaging, must be kept at an ambient temperature of 10° C to 40° C and a relative humidity of 30% to 75% non-condensing.

Transporting tools made according to the conditions of storage.

15.UTILIZATION

Use a set of catheters disposed of as hazardous waste epidemiologically class B (presence of contact with blood and mucosal tissues of the patient).

After the expiry date the unused sets of catheters not in contact with blood and mucous membranes should be disposed of as hazardous waste epidemiologically, close in composition to the solid waste. Dispose of used conductors shall be in accordance with the laws and regulations on the disposal of products made in the country of use of the medical device.

16. WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered,

defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

17. COMPLAINT


The organization, which takes in the Russian Federation claims from consumers on the quality of the product:

Limited Liability Company "Merit Technologies" (LLC "Merit Technologies"), BIN 1117746679533, TIN 7705960130

Legal address: 115054, Moscow, Bakhrushina Str. 32, p. 1

Tel / fax: +7 (495) 221-89-02

A block diagram of the manufacturing process

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