Apostille

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**Important information
for users of KARL STORZ
devices**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
other products, this product is the result of years
of experience and great care in manufacture. You
and your organization have decided in favor of
a modern, high-quality item of equipment from
KARL STORZ.

This instruction manual is intended to serve as an
aid in the proper setup, installation, and operation
of the KARL STORZ CALCULASE II. All essential
details of the equipment and all actions required
on your part are clearly presented and explained.
We thus ask that you read this manual carefully
before proceeding to work with the equipment.
Insert this manual in its protective wallet and keep
it available for ready reference in a convenient and
conspicuous location near the equipment.

All product illustrations, product descriptions and
texts are the intellectual property of KARL STORZ
GmbH & Co. KG. Their use and reproduction
by third parties require the express approval of
KARL STORZ GmbH & Co. KG.

All rights reserved.

CAUTION: Federal (USA) law restricts this device
to sale by or on the order of a physician.

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



Controls, displays, connectors, and their uses

- ① Keyswitch (0, 1)
- ② Emergency Stop Switch
- ③ ± keys to set aiming beam brightness
- ④ ± keys to set frequency
- ⑤ LED frequency display
- ⑥ ± keys to set pulse energy
- ⑦ LED pulse energy display
- ⑧ Ready/Standby key
- ⑨ LED display for READY mode
- ⑩ LED display for STANDBY mode
- ⑪ 'Red' LED display to indicate no fiber/
'Yellow' LED display to indicate LASER ACTIVE
- ⑫ Footswitch connection socket
- ⑬ Fiber connector/Laser outlet aperture
- ⑭ Potential equalization connector
- ⑮ Power supply cord receptacle
- ⑯ Power switch for turning the unit on and off (fused with 16 A)
- ⑰ Fuses for internal 24 V module
- ⑱ Remote interlock connection
- ⑲ SCB sockets
- ㉑ Service port (for service purposes only)
- ㉒ Label acc. to IEC 60825
- ㉒ Fiber tip
- ㉓ Protective cover around laser fiber
- ㉔ Fiber connector
- ㉕ Sterilization counter
(on reusable fiber)
- ㉖ Footswitch with guard

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The SCB interface, which is based on the CAN
field bus, makes it possible to control unit functions
and display unit parameters remotely.

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	Follow instructions for use
	Laser radiation
	Potential equalization connector
	Remote interlock connection
	Service port
	Footswitch connection
	Fiber (laser fiber) (symbol over LASER ACTIVE LED)
	Ready/Standby button
	Symbol above STANDBY LED
	Symbol above READY LED
	Energy
	Frequency
	Aiming beam dark
	Aiming beam bright
	Laser Emergency stop
	Degree of protection against electric shock: Applied part type B
	This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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Important information for users of KARL STORZ devices	III
Images of the unit	V
Controls, displays, connectors, and their uses	VII
Symbols employed	VIII

System descripton	3
General information	3
Ho:YAG laser theory	4

Safety instructions	5
Warnings and cautions	5
Safety precautions at the site of installation	9
Safety precautions when operating the system	9
Operator's responsibility	11
Intended use	12
Improper use	12
Organizational measures	13
Informal safety measures	13
Laser area	14
Personal safety	14
Patient safety	15
Safety features of the system	16
Keyswitch	16
Emergency stop switch	17
Remote interlock connection	17
Fault detection	18

Installation and operating instructions	19
Unpacking and inspection	19
Basic equipment	19
Installing and connecting up the unit	20
SCB	21
Emergency stop switch	22
Switching the unit on	22
Connecting the fiber	23
Removing the fiber	28
Trimming the fiber	29

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Ho:YAG laser theory

The CALCULASE II is a Ho:YAG laser, which emits a mid-infrared beam.

A laser rod, coated with Holmium, is optically excited by a flash lamp, causing energy to be transmitted into the crystal. This energy excites electrons in the Holmium to a higher level of energy. As soon as these electrons revert to their original level of energy, photons with a precisely defined energy difference (Holmium wavelength 2.1 µm) are released. As a result of a process of 'stimulated emission', these photons in the Ho:YAG lasing medium cause other photons to be released, inducing, in a snowball effect, the formation of a beam of photons of identical direction of travel, wavelength, phase and polarization. To amplify this effect the beam is passed back and forth between mirrors in the laser resonator and a small percentage is released via a fiber in the form of a laser beam.

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Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patients or physician may be involved.
Disregarding a Warning or improper use of the laser could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A Note indicates special information about operating the unit, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using the CALCULASE II system. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

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WARNING: The KARL STORZ CALCULASE II is a Class 4 laser according to IEC 60825-1. The laser may only be operated in designated laser areas on the authority of a surgeon.

WARNING: The laser may only be operated after receiving instruction by KARL STORZ or an authorized representative on the basis of this instruction manual. Further instructions to be followed before operating the unit for the first time are given on p. 13.

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WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1).

WARNING: Risk of explosion if the CALCULASE II system is used in the presence of flammable anesthetics.

WARNING: To avoid the risk of electrical shock, do not open the unit! Refer servicing to qualified personnel.

WARNING: Grounding reliability for the unit can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e. approved for use in an operating room environment). Routinely inspect the electrical plug and cord. Do not use if inspection reveals damage.

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Safety instructions

Warnings and cautions

WARNING: Inadvertent activation of the footswitch in READY Mode can trigger laser pulses. For safety reasons, the laser device should be switched from READY Mode into STANDBY Mode between treatments.

WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards or equivalent national regulations.

WARNING: Keep the unit out of reach of patients.

WARNING: Always unplug the unit before doing any maintenance work on it (e.g. cleaning).

WARNING: Caution – Use of the operator controls or setting options in ways other than those described here can result in dangerous irradiation.

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

CAUTION: Only operate the CALCULASE II with the voltage indicated on the nameplate.

CAUTION: A laser unit, just as any other high performance medical device, requires special knowledge and care with regards to handling and use. The laser may only be operated by physicians who have been thoroughly trained in handling the device in accordance with this operating manual and are familiar with its therapeutic effect and possible risks.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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The system may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

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Safety precautions when operating the system

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CE marked equipment has been tested and found to comply with the EMC limits for the Medical Device Directive 93/42/EEC (EN 55011 Class B and EN 60601-1-2). Please read the Electromagnetic Compatibility Information in the manual.

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The laser radiation of the CALCULASE II has a power of 20 Watts max. and a wavelength of 2080 nm which is invisible to the human eye (the green aiming beam operates on a wavelength of 532 nm and is visible to the human eye). Incorrect use may lead to third degree burns even if the irradiation is unfocussed.

Flammable drapes, surgical gowns, gauze and other ignitable materials must be kept clear of the laser beam path.

The use of non-flammable materials and instruments is advised. Flame retardant surgical drapes, gowns, etc., are recommended.

! WARNING: There is a risk of fire or explosion if the LASER BEAM is used in the presence of flammable materials, solutions or gases or in an oxygenated environment. Some materials which are saturated with oxygen, for example, cotton, may ignite when exposed to the high temperatures produced during the INTENDED USE of the laser device. The solvents in adhesives and flammable solutions, which are used for cleaning and disinfection, should be allowed sufficient time to evaporate before the laser equipment is used. Please note that endogenous gases are flammable.

! WARNING: Laser smoke can contain viable tissue particles. Smoke evacuation is necessary for applications which produce smoke gas.

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Operator's responsibility

The operator and safety officer have responsibility for carrying out all the safety measures to prevent any risk to the patient, attending physician and other people present during use of the laser. In this respect, the German Social Accident Insurance (DGUV) regulation 11 on laser radiation must be observed (Germany only). In addition, the accident prevention rules and regulations applicable to the place of use must be adhered to.

The operator must keep the following documents accessible at all times and must update them regularly (Germany only).

- Inventor list
- Equipment book and/or equipment chart
- Instruction manual

Accidents and personal injury must be reported to the responsible authorities without delay.

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Organizational measures

The required personal protection equipment must be provided by the operator.

All the safety devices fitted must be checked at regular intervals.

When the laser unit is not in use, the key must be removed from the keyswitch ① to prevent unauthorized use.

Before operating for the first time

The initial installation of CALCULASE II should be carried out by a trained, authorized person in compliance with the information contained in the manual entitled 'CALCULASE II Initial Installation'. The system must first be filled with cooling water, otherwise start-up will not be possible (also see 'caution note' on page 19).

The guidelines of the country of use for installing and setting up laser systems must always be followed.

Informal safety measures

WARNING: The instruction manual must be kept with the CALCULASE II at all times. In addition to the instruction manual, the generally applicable and local rules on accident prevention and environmental protection must be available and must be observed.

All the safety and hazard information on the unit must be kept in a legible condition and replaced if necessary.

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Safety features of the system

The safety standards IEC 60601 and IEC 60825-1 determine the general requirements relating to the design and manufacture of the KARL STORZ CALCULASE II. In compliance with international standards, the CALCULASE II is equipped with the following:

- Energy indicator ⑦
- Keyswitch ①
- Emergency stop switch ②
- Remote interlock connection ⑩
- Proper signs and labels (acc. to IEC 60825), see images on left

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Remote interlock connection

To provide added safety, the CALCULASE II laser unit is equipped with a remote interlock connection point ⑩. The connection point is situated on the service panel at the rear of the unit. An external switch can be connected to this point to create a remote interlock system.

This switch should be mounted at the entrance of the treatment room, such that if the door opens, the switch contacts also open and disable the laser radiation.

Once the switch contacts are opened:
the unit stops laser radiation;
laser mode changes to Standby;
Message display indicates:
REMOTE INTERLOCK DISCONNECTED

To resume the operation:
make sure the switch contacts are closed;
return the system to the READY state;
press the footswitch for laser radiation.

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① NOTE: The connection of an external interlock must be performed by an authorized electrician. For the remote interlock, switches with basic insulation suitable for 230 V AC must be used.

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Fault detection

The unit is equipped with a self-check routine that continuously monitors system operation. If a fault condition is detected:

- The LED displays issue a fault message (e.g. E 04 – footswitch pressed when turned on).
- The audible indicator emits an alarm signal.
- The system ceases or reduces laser emission to ensure user and patient safety and to prevent damage to the unit.

The list of fault messages, their meaning and possible remedies is given on p. 40.

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Unpacking and Inspection

Carefully remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can be used if the unit has to be transported.

CAUTION: The initial installation of CALCULASE II may only be carried out by qualified service personnel. To ensure safe operation, CALCULASE II should be filled with cooling water and an ion filter should be installed. Detailed instructions are provided in the 'CALCULASE II Initial Installation' instruction manual available to qualified personnel and on request.

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Only operate the CALCULASE II at the line voltage stated on the manufacturer's identification plate. Connect up the power cord and insert power cord connector fully into the power socket ⑮.

The power cord connector must be pulled out to ensure an all-pole disconnection from the mains power supply. For this reason, the device may not be set up in a manner which makes it difficult to disconnect it from the mains power supply.

! WARNING: The unit must not be operated in potentially explosive atmospheres.

! CAUTION: Do not connect the laser either to a multiple power socket outlet or a socket which can be switched off.

After inserting the power cord connector, turn on the main switch ⑯.

The main switch has a two-pole automatic circuit breaker. This can be reset after it has been tripped due to a current overload.

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

For added protection the CALCULASE II is equipped with a connection for remotely interrupting laser emissions. An external switch is connected to this outlet (18), which enables remote interlocking. The external switch is connected e.g. to the door of the operating room, so that when the door opens the laser beam is interrupted, the unit returns to standby mode and the error message I is issued. To resume operation, the external switch contact must be closed and then the unit turned to ready mode once again (green LED (9) lights up). The remote interlock jumper (27 7502 26) can be inserted instead of the switch.

NOTE: If neither the remote interlock connection nor the remote interlock jumper is connected, it is not possible to operate the laser.

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Emergency stop switch

To begin laser operation the emergency stop switch (2) must be unlocked by turning it.

Push the emergency stop switch (2) in an emergency; the laser is switched off instantly. To resume operation the switch must be unlocked.

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WARNING: For normal system shutdown the keyswitch (1) and the main power switch (16) must be used. Pushing in the EMERGENCY STOP switch (2) stops laser emission instantly. If the EMERGENCY STOP switch (2) has been pushed in, it must be unlocked before the unit is started up again. To do so, turn the switch.

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
Switching the unit on

The unit is switched on (position '1') and off (position '0') using the keyswitch (1) on the front of the unit.

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CAUTION: Never leave the unit unattended in activated condition. Secure the unit against unauthorized use by pulling the key from the keyswitch (1).

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Connecting the fiber

You cannot operate the unit until a fiber is connected. The red LED ⑪ lights up when no fiber is connected.

To insert the fiber, open it, move the lever on the socket ⑬ to the left or right and simultaneously attach the fiber by turning it. Once the fiber is connected, the red LED ⑪ goes out and you can operate the unit.

WARNING: Only clean fibers in perfect condition that have been sterilized according to the instructions may be used for treatment. The reprocessing of the fibers must always be monitored.

NOTE: The fibers are color coded on the anti-kink protection.

Fiber	Color
230 µm	blue
365 µm	green
600 µm	red

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WARNING: Never activate the laser when the fiber is defective.

WARNING: Avoid extreme bending or tight winding of the fiber. This could cause the fiber to break. The minimum bending radius is between 1.5 and 3.5 cm, depending on the fiber

(230 µm fiber: 1.5 cm, 365 µm fiber: 2.0 cm, 600 µm fiber: 3.5 cm).

Ensure that this is also taken into consideration for storage.

WARNING: Insufficient fixation or severe bending of the fiber, or the use of fiber optics, can pose the following risks:

- Injury to the patient, user and/or assistant
- Damage to the fiber.

WARNING: Pushing the laser fibers forward in the instrument under poor visual conditions can result in the perforation of the bladder or ureter as a result of the mechanical force exerted.

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Installation and operating instructions

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The fiber must be checked at the beginning of every laser treatment. To do so, the aiming beam must be rechecked and reset every time you use it. The brightness of the aiming beam can be increased by pressing the '+' key ③ and decreased using the '-' key ③.

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
The frequency is set using the '+' and '-' keys ④ in the frequency field. There are up to five different frequency levels to choose from depending on the fiber used: 4, 6, 8, 10, and 15 Hz.

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
The energy is set using the '+' and '-' keys ⑥ in the energy field. Up to five different energy levels are available here depending on the fiber and specified frequency: 0.5/0.8/1.2/1.7, and 2.0 J.

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The laser is always in standby mode while and after the parameters are set. The STANDBY LED ⑩ (yellow) lights up.

Changing to ready mode

To 'activate' the laser, press the ready/standby key ⑧ once. The READY LED ⑨ (green) lights up in the READY mode.

① **NOTE:** This procedure is accompanied by a three-time change between the yellow and green LED (STANDBY and READY LED).

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① **NOTE:** The laser beam can only be activated once the target energy has been attained (takes approx. 1 second).

The activation/output of laser energy is displayed by the yellow LASER 'ACTIVE' LED ⑪ and by an activation tone.

⚠ **WARNING:** For safety reasons, the laser should be switched to standby mode during extended periods of (foreseen) non-activation.

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The CALCULASE II is equipped with a fiber recognition system. A permissible frequency/energy combination is provided depending on the fiber used, from which the user can select the corresponding capacity (see also p. 23).

CAUTION: The use of fibers which have not been authorized by KARL STORZ for use with the device may damage the fiber recognition system and thus lead to failure of the device.

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Parameter matrix for 230 µm fibers

Energy [J]	Frequency [Hz]				
	4	6	8	10	15
0.5	2.0 W	3.0 W	4.0 W	5.0 W	–
0.8	3.2 W	4.8 W	6.4 W	8.0 W	–
1.2	4.8 W	7.2 W	9.6 W	12.0 W	–
1.7	–	–	–	–	–
2.0	–	–	–	–	–

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
NOTE: Thanks to its desktop design the KARL STORZ CALCULASE II is an extremely mobile system. The laser can be transported using special retractable handles.

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Removing the fiber

For safety reasons, switch the device to standby mode; to do so, press the ready/standby key ⑧. Then remove the fiber.

! WARNING: Disposable fibers must be disposed of and reusable fibers reprocessed.

Seal the reusable fiber with the cover cap for reprocessing.

Inspect the reusable fiber for defects prior to reprocessing. Dispose of defective fibers or trim them afresh (see Section 'Trimming the fiber').

⑧ NOTE: You will find information on the sterilization of reusable fibers in the instruction leaflet accompanying the fibers.

⑧ NOTE: Prior to every sterilization of a reusable fiber, a number must be cut off the sterilization counter (from 10 to 0). The fiber should be changed when you reach ZERO at the latest.

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To break the fiber cleanly, it is scored with a ceramic scoring pen (27750284). It should then be held taut over the silicon pad (27750285) and broken.

NOTE: Complete fiber stripping set 27750280 (see also recommended accessories).

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Laser protective eyewear

Laser protective eyewear prevents the eyes from being harmed by intensive light during a light pulse.

The glasses are equipped with side guards for additional safety against lateral reflections.

WARNING: During treatment, the physician and everyone who is present in the laser area must wear special laser protective eyewear for 2080 nm wavelength (e.g. laser protective eyewear 27 750095).

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NOTE: The frame or bow of the goggles should have a label indicating the specific wavelength spectrum they are designed for. The wavelength on the label must correspond to the laser's wavelength or include it.

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Fuse replacement

! NOTE: The power switch 16 (rocker switch) acts as a line fuse. The fuses 17 fuse the 24 V module.

Procedure (for non-functioning devices)

- Check whether power switch is turned on.
- Check fuses 17.

To do so, turn off the device and disconnect it from the power supply.

Remove the fuse holder 17 with a screwdriver.

! CAUTION: Use only fuses with the correct nominal ratings.

Insert new fuses (see table).
Refit the fuse holder 17.

	100...240 V~
Fuse	2 x T 10AL, 250 V (5x20mm)

Reconnect the power.

Perform function testing.

! NOTE: If the device still does not function correctly after the fuse has been changed, please contact the Service department.

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Reprocessing

General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **WARNING:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **WARNING:** National laws and regulations must be observed.

 **NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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Manual wipe-down disinfection

WARNING: Before cleaning or disinfecting, unplug the power cord. Personnel are only allowed to clean the exterior surfaces. Never open the unit.

CAUTION: When cleaning the unit make sure that no liquid penetrates into the device since this may cause malfunction to the unit. Take special care of the fiber connection part/laser outlet aperture!

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility must be observed. Do not use sprays because no liquid must be allowed to enter the interior of the unit.

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Maintenance and safety inspection

WARNING: At least once per year (within 12 months), maintenance and a safety inspection/repeat inspection of the complete laser system as defined by IEC 62353 must be carried out by a service technician authorized by KARL STORZ.

Maintenance

Annual maintenance is essential. Regular maintenance can contribute to identifying potential problems before they become serious, as well as enhancing the device's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer.

Maintenance of the CALCULASE II must include the following work at least:

- Changing the coolant
- Replacing the particle filter
- Replacing the ion exchanger

Safety inspection

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety inspections/repeat inspections as defined by IEC 62353 must be carried out by a service technician authorized by KARL STORZ at least once a year.

The inspection includes accessories, wearing parts, disposable articles and the interfaces connected to the unit that are used to operate the unit if this accessory equipment has any impact on the unit's operational safety.

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Visual inspection

1. Check the device and accessories for any
mechanical damage which may impair
functionality.
2. Check that inscriptions relevant to safety are
legible.

Electric measurements

- Inspect the device safety fuses
- Protective ground resistance measured accord-
ing to IEC 62353: For limit values, please refer to
the current standard.
- Earth leakage current measured according to
IEC 62353: For limit values, please refer to the
current standard.
- Touch current measured according to
IEC 62353: For limit values, please refer to the
current standard.
- Patient leakage current measured according to
IEC 62353: For limit values, please refer to the
current standard.

Test for proper functioning

The test for proper functioning must include the
following tests at least:

- Emergency Stop
- Door interlock switch (remote interlock)
- Displays and control keys
- Checking all energy/frequency combinations of
the laser parameters

Documentation

The safety inspection/repeat inspection and the
results must be recorded in the medical device
logbook.

WARNING: If any defects or shortcomings
which could endanger patients, personnel
or others are ascertained during the safety
inspection/repeat inspection, the unit must
not be operated until such time as these
defects or shortcomings have been elimi-
nated by qualified technical servicing.

NOTE: For detailed information on the range
and implementation of the safety inspections/
repeat inspections, please refer to the current
version of the relevant Service Manual.

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Servicing and repair

WARNING: Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life*, dispose of the unit as electronic waste.

* expected useful operating life: 7 years

Please ask KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your distributor for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff).

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General

In the event of the unit malfunctioning, the troubleshooting list provided in this chapter will assist the user in troubleshooting. The guide explains how to rectify the various malfunctions.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ Customer Service technician.

WARNING: Maintenance work that is not mentioned in this instruction manual must only be carried out by KARL STORZ authorized technical personnel.

WARNING: Maintenance work or modifications to the laser unit which are carried out contrary to this instruction may result in personal injury, damage to the unit and loss of any guarantee rights.

WARNING: Always unplug the unit from the power supply before carrying out any maintenance work.

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Troubleshooting guides

The troubleshooting tables do not attempt to list all possible system failures. Any fault not listed, should be referred to KARL STORZ customer service.

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Error message	Possible cause	Remedy
No function	Power plug not connected to outlet	Check power plug and, if necessary, connect to power supply.
	Emergency stop switch activated	Unlock emergency stop switch by turning it.
Red LED lights up.	Fiber not connected	Connect fiber or inspect it.
E 01 – 02	Technical error	Note error code and notify Service department.
E 03	Footswitch is not connected or is defective	Check and connect footswitch.
E 04	Footswitch pressed when turned on	Release footswitch.
E 05	Remote interlock active when turned on	Close remote interlock and turn on again.
E 06	Flash lamp aged – standard service life: < 10 %	Confirm with the ready key to continue use, arrange a servicing date, and replace the flash lamp.
E 07	Flash lamp defective	Send unit in for servicing and have flash lamp replaced.
E 08	Key pressed when turned on	Release all keys and turn on again.
E 09 – 13	Technical error	Note error code and notify Service department.
E 14	Water temperature too low	Wait until heating-up phase is complete.
E 15	Water temperature too high	Allow device to cool in standby mode, check room temperature.
E 16 – 47	Technical error	Note error code and notify Service department.
E 28	Laser energy too low	Select another energy/frequency setting to continue use, note the error code, and notify the Service department.
E 29	Laser energy too high	Select another energy/frequency setting to continue use, note the error code, and notify the Service department.
I	Interlock in operation	Interlock in operation

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Technical data

Specifications

Laser type:	Pulsed Ho:YAG laser
Wavelength:	2080 nm
Power:	Po = max. 20 Watts, ($\pm 10\%$) Pp = max. 10 kW
Pulse frequency:	4 / 6 / 8 / 10 / 15 Hz
Pulse duration:	100...500 μ s
Energy:	0.5 / 0.8 / 1.2 / 1.7 / 2.0 J (max. 2 J $\pm 20\%$)
Beam transmission:	Fiber system
Fiber types:	Bare fiber, 230 μ m, 365 μ m, 600 μ m fiber
Fiber connection:	Fiber recognition via modified SMA fiber connector
Aiming beam:	4 mW, 532 nm
Operating voltage:	100...240 VAC ($\pm 10\%$)
Line frequency:	50/60 Hz
Power consumption:	Max. 1800 VA
Cooling:	Integrated water cooling system with water-/air- heat exchanger
Dimensions:	610 x 250 x 435 mm
Weight:	35 kg (incl. 2 kg cooling water)
Storage/transport conditions:	
Temperature:	-10° C ... 60 °C
Air humidity:	5 % ... 95 %
(RH, non-condensing)	
Operating conditions:	
Temperature:	18 °C ... 28 °C
Air humidity:	30 % ... 70 %

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(for 27 7502 20-1)

According to IEC 60601-1, IEC 60601-2-22:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type B ⬆
- Type of protection against moisture: drip-water protection as per IPX 1

According to IEC 60825-1:

- Laser class: 4

IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 49-62).

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

① NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
Spare parts, accessories

Article	Cat. no.
Laser safety goggles Ho:YAG Laser 2080	27.750095
Safety goggles surgical LASER, for users with and without glasses	UL020
Safety goggles surgical LASER, sporty version for users who don't wear glasses	UL021
Reusable fiber sets CALCULASE II Bare Fiber 230 µm, reusable, sterile, package of 6 pieces CALCULASE II Bare Fiber 365 µm, reusable, sterile, package of 6 pieces CALCULASE II Bare Fiber 600 µm, reusable, sterile, package of 6 pieces	27.750271-P6 27.750272-P6 27.750273-P6
Single-use fiber sets CALCULASE II Bare Fiber 230 µm, single use, sterile, package of 6 pieces CALCULASE II Bare Fiber 365 µm, single use, sterile, package of 6 pieces CALCULASE II Bare Fiber 600 µm, single use, sterile, package of 6 pieces	27.750277-P6 27.750278-P6 27.750279-P6
Warning sign 'Laser radiation'	27.750092
Set of keys CALCULASE II	27.750235
Remote interlock jumper CALCULASE II	27.750226
Fiber stripping set	27.750280
Footswitch CALCULASE II	27.750124

CAUTION: When handling accessories and substances that come in contact with patients, the safety regulations of the country of use must be followed to prevent the transmission of infections.

CAUTION: Never operate with defective accessories or a damaged power cord.

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Appendix

Electromagnetic compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The CALCULASE II model 27750220-1 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The CALCULASE II model 27750220-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix**Electromagnetic compatibility (EMC)
information**

WARNING: The CALCULASE II model 27 7502 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the CALCULASE II model 27 7502 20-1. The cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
Footswitch with cable	Yes	3	No	Activating the unit functions
Power Cord 10A with non-heating appliance connector	No	5	No	Power supply

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The CALCULASE II model 27750220-1 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCULASE II model 27750220-1 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using an AC line filter.
Surge IEC 61000-4-5	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	Complies ± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T^*$ (>95 % dip in U_T) for 0.5 cycle $40\% U_T$ (60 % dip in U_T) for 5 cycles $70\% U_T$ (30 % dip in U_T) for 25 cycles $<5\% U_T$ (>95 % dip in U_T) for 5 sec	Complies $<5\% U_T^*$ (>95 % dip in U_T) for 0.5 cycle Complies $40\% U_T$ (60 % dip in U_T) for 5 cycles Complies $70\% U_T$ (30% dip in U_T) for 25 cycles Complies $<5\% U_T$ (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 4

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The CALCULASE II model 27.750220-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the CALCULASE II model 27.750220-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V and 10 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the CALCULASE II model 27.750220-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m and 10 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

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^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds

Table 6

**Recommended separation distances between portable and mobile RF communications equipment and the
CALCULASE II model 27750220-1**

The CALCULASE II model 27750220-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopy México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, Mexico, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F10-Q, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai (Zhangjiang)
Room 701A, 7F Building 5, No. 3000 Longdong
Avenue, Zhangjiang Hi-Tech Park, Pudong District,
Shanghai, 201203, People's Republic of China
Phone: +86 21 61009688, Fax: +86 21 61009686
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: Infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

The laser urological holmium device CALCULASE II in current instruction is mentioned as device, device CALCULASE II, CALCULASE II, device CALCULASE II model 27 7502 20-1 or laser.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the KARL STORZ CALCULASE II. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

All product illustrations, product descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG. Their use and reproduction by third parties require the express approval of KARL STORZ GmbH & Co. KG.

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
- 1 Keyswitch (0, 1)
- 2 Emergency Stop Switch
- 3 ± keys to set aiming beam brightness
- 4 ± keys to set frequency
- 5 LED frequency display
- 6 ± keys to set pulse energy
- 7 LED pulse energy display
- 8 Ready/Standby key
- 9 LED display for READY mode
- 10 LED display for STANDBY mode
- 11 'Red' LED display to indicate no fiber/
'Yellow' LED display to indicate LASER ACTIVE
- 12 Footswitch connection socket
- 13 Fiber connector/Laser outlet aperture
- 14 Potential equalization connector
- 15 Power supply cord receptacle
- 16 Power switch for turning the unit on and off (fused with 16 A)
- 17 Fuses for internal 24 V module
- 18 Remote interlock connection
- 19 SCB sockets
- 20 Service port (for service purposes only)
- 21 Label acc. to IEC 60825
- 22 Fiber tip
- 23 Protective cover around laser fiber
- 24 Fiber connector
- 25 Sterilization counter
(on reusable fiber)
- 26 Footswitch with guard

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The SCB interface, which is based on the CAN field bus, makes it possible to control unit functions and display unit parameters remotely

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General information

The KARL STORZ CALCULASE II (Art. no.: 27 7502 20-1) is a Ho:YAG desktop laser. It emits a mid-infrared beam at a wavelength of 2080 nm.

The laser energy is transmitted into the tissue via a fiberoptic fiber, the maximum laser power being 20 W. The laser energy generated by the CALCULASE II system enables the optimal lithotripsy of all calculus compositions in the human body and soft tissue treatment such as the ablation of strictures, stenoses, and tumors, for example.

Thanks to its desktop design the unit is extremely mobile and is easy to transport using special, retractable handles.

All controls and displays required for operating the unit are located on the front of the CALCULASE II. The laser has a user-friendly control panel and easy-to-read LED displays, which make the unit very easy to operate, and safe.

Property rights

This product is protected in the USA by (at least one of the following) US Patent no(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539

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A laser rod, coated with Holmium, is optically excited by a flash lamp, causing energy to be transmitted into the crystal. This energy excites electrons in the Holmium to a higher level of energy. As soon as these electrons revert to their original level of energy, photons with a precisely defined energy difference (Holmium wavelength 2.1 µm) are released. As a result of a process of «stimulated emission», these photons in the Ho:YAG laser medium cause other photons to be released, inducing, in a snowball effect, the formation of a beam of photons of identical direction of travel, wavelength, phase and polarization. To amplify this effect, the beam is passed back and forth between mirrors in the laser resonator and a small percentage is released via a fiber in the form of a laser beam.

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The KARL STORZ CALCULASE II has been specially designed to minimize accidental exposure to hazardous laser radiation. The Holmium laser beam is mainly absorbed in water.

During treatment, contact lithotripsy is the method used, meaning that there is direct contact between the laser fiber and calculus, and the laser energy is transmitted directly onto the calculus. The method employs short pulse sequences of a duration of approx. 1-2 seconds.

The laser unit is also equipped with a visible, green aiming beam with a wavelength of 532 nm, which helps the surgeon to aim the laser beam.

Reflected and direct eye exposure hazard

During the procedure the laser emits high radiation power in the invisible spectral range of 2080 nm that is hazardous to the human eye and skin.

As a precaution against accidental exposure to the output laser beam or its reflections, all personnel and patients must use appropriate safety eyewear.

Do not stare into the laser beam or allow it to be reflected from any surface. Even rough metal can reflect the laser beam.

WARNING: Do not stare into the tip of the fiber. To prevent skin burns, do not aim the laser beam directly at unprotected skin.

WARNING: Do not stare directly into the aiming laser beam either. Also ensure that the aiming beam is not directed at anyone's eyes. Although the beam is low powered, direct exposure can be hazardous to the eyes.

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Safety precautions at the site of installation

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the system must not be operated inside the demarcated hazard zone. This also applies to highly combustible and potentially explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

Safety precautions when operating the device

It is the user's responsibility to make sure the system is safe and operates properly before using it.

During treatment with the CALCULASE II system the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the unit.

CE marked equipment has been tested and found to comply with the EMC limits for the Medical Device Directive 93/42/EEC (EN 55011 Class B and EN 60601-1-2).

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The laser radiation of the CALCULASE II has a power of 20 Watts max. and a wavelength of 2080 nm which is invisible to the human eye (the green aiming beam operates on a wavelength of 532 nm and is visible to the human eye). Incorrect use may lead to third degree burns even if the irradiation is unfocused.

Flammable drapes, surgical gowns, gauze and other ignitable materials must be kept clear of the laser beam path. The use of non-flammable materials and instruments is advised. Flame retardant surgical drapes, gowns, etc., are recommended

WARNING: There is a risk of fire or explosion if the LASER BEAM is used in the presence of flammable materials, solutions or gases or in an oxygenated environment. Some materials which are saturated with oxygen, for example, cotton, may ignite when exposed to the high temperatures produced during the INTENDED USE of the laser device. The solvents in adhesives and flammable solutions, which are used for cleaning and disinfection, should be allowed sufficient time to evaporate before the laser equipment is used. Please note that endogenous gases are flammable.

WARNING: Laser smoke can contain viable tissue particles. Smoke evacuation is necessary for applications which produce smoke gas.

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Intended use

The device is used to create laser energy for the destruction of calculi and/or for soft tissue treatments such as the cutting of tissue and the opening of stenoses and strictures during endoscopic urological applications.

CAUTION: The device CALCULASE II may only be used with accessories authorized by KARL STORZ.

Incompatible accessories can cause damage to the unit, not all unit parameters can be set, or the treatment is not possible.

Indications

Disintegration of stones in the urinary tract via endoscopic access ports and/or for soft tissue treatments such as the cutting of tissue and the opening of stenoses and strictures during endoscopic urological applications.

Contraindications

The use of the CALCULASE II is contraindicated if calculus treatment or soft tissue treatment with a laser is contraindicated, and if in the opinion of the attending physician such an application would endanger the patient, e.g., due to the general condition of the patient.

Improper use

WARNING: Improper use which does not comply with the contents of the instruction manual may result in dangerous radiation.

Use in fields other than those indicated above is not allowed.

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Laser area

The laser area is the area in which the levels for maximum permissible irradiation may be exceeded. This means the entire operating room unless building measures have been carried out to limit the area. The possibility of unintentional deflection of the laser beam must be taken into account. The access routes to the laser areas must be identified with warning lamps. A flashing light is recommended.

Signs must be placed in the laser area giving an explicit warning about the intensive light and the laser beam used. The illustration opposite shows the sign (27 7500 92) for the laser area, which is optionally available with CALCULASE II.

Warning signs for the laser area

Position the warning signs conspicuously.

Only allow access to the CALCULASE II laser area to those members of staff who are absolutely necessary for the treatment and who have been thoroughly trained in the necessary safety measures.

CAUTION: In the laser area no reflective objects, e.g. surgical instruments, mirrors, watches or jewelry are permitted which might reflect the laser.

Personal safety

- Ensure that all members of staff in the laser area are familiar with controlling the CALCULASE II and know how to shut down the system rapidly in an emergency.
- Never look directly into the light beam coming from the laser fiber, even if you are wearing protective eyewear.
- Never direct the laser fiber around the room.
- When the CALCULASE II is not being used, switch it off immediately and remove the key.

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WARNING: During treatment, the physician and everyone who is present in the laser area must wear special laser protective eyewear for 2080 nm wavelength (e.g. laser protective eyewear 27 7500 95.)

CAUTION: Non-compliance can cause irreversible eye damage! Please verify the perfect condition of the laser protection glasses before use, i.e. the glasses must be free of any mechanical damage. Only the type of protective eyewear indicated below may be used.

NOTE: The frame or bow of the goggles should have a label indicating the specific wavelength spectrum they are designed for. The wave length on the label must correspond to the laser's wavelength or include it.

Patient safety

The safety of the patients is ensured primarily by well-trained staff and a laser area equipped according to the regulations. The patients must be informed about the safety measures which are important for them and about the type of treatment.

- Ensure that the patient's eyes are protected against the laser beam.
- If the patient is not able to wear safety goggles, provide them with opaque eye protection which entirely prevents the light from reaching their eyes.

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
Safety features of the device

The safety standards EC 60601 and IEC 60825-1 determine the general requirements relating to the design and manufacture of the KARL STORZ CALCULASE II.

In compliance with international standards, the CALCULASE II is equipped with the following:

- Energy indicator ⑦
- Keyswitch ①
- Emergency stop switch ②
- Remote interlock connection ⑩
- Proper signs and labels (acc. to EC 60825), see images on left

WARNING: The unit is grounded by a grounded conductor in the power cord. This grounding is necessary for safe operation.

Before switching on the unit, all the safety devices must be properly fitted and fully functional.

Keyswitch

The keyswitch ① is used to turn the unit on and off. The system can only be operated with the key provided by KARL STORZ. This prevents unauthorized use of the unit. Always remove the key after switching the unit off and store it in a safe place.

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Emergency stop switch

The red round button is used for shutting down the unit rapidly in an emergency. If you press the emergency stop switch ② the power supply to the system is immediately disconnected. The complete switch-off procedure using the keyswitch ① is bypassed.

Pull the emergency stop switch ② by turning it before switching the unit back on. The system cannot be switched on if the emergency stop switch is still pressed in.

Remote interlock connection

To provide added safety, the CALCULASE II laser unit is equipped with a remote interlock connection point ④. The connection point is situated on the service panel at the rear of the unit. An external switch can be connected to this point to create a remote interlock system.

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To resume the operation: make sure the switch contacts are closed; return the system to the READY state; press the footswitch for laser radiation.

NOTE: The connection of an external interlock must be performed by an authorized electrician. For the remote interlock, switches with basic insulation suitable for 230 V AC must be used.

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Installing and connecting up the unit

Set up the CALCULASE II on a flat surface as a free-standing unit (ideal for a KARL STORZ equipment trolley).

NOTE: The CALCULASE II may only be used in medical environments set up in accordance with applicable national regulations for laser use.

¹⁴ CALCULASE II is equipped with a connector for attaching a ground line.

The ground line should be installed by a qualified electrician.

Protection Class I device

! WARNING: To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.

Only operate the CALCULASE II at the line voltage stated on the manufacturer's identification plate. Connect up the power cord and insert power cord connector fully into the power socket ¹⁵.

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! CAUTION: Do not connect the laser either to a multiple power socket outlet or a socket which can be switched off. After inserting the power cord connector, turn on the main switch ¹⁶. The main switch has a two-pole automatic circuit breaker. This can be reset after it has been tripped due to a current overload.

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SCB

NOTE: To prevent the SCB connecting cable from being inadvertently pulled out, the SCB connector has a safety device.

Pull back the safety device on the SCB connector and insert the connector into one of the sockets (19). Connect the other end of the cable to a KARL STORZ SCB control unit (KARL STORZ Communication Bus) or other SCB devices (see the instruction manual for the KARL STORZ SCB System in relation to this).


Insert the footswitch connecting cord into socket (12) (in alignment with red dot).

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
Emergency stop switch

To begin laser operation the emergency stop switch W must be unlocked by turning it.

Push the emergency stop switch W in an emergency; the laser is switched off instantly. To resume operation the switch must be unlocked.

NOTE: The emergency stop switch does not necessarily have to be pushed in when operating the unit for the first time. To unlock it, turn the switch.

! WARNING: For normal system shutdown the keyswitch Q and the main power switch r must be used. Pushing in the EMERGENCY STOP switch W stops laser emission instantly.

If the EMERGENCY STOP switch W has been pushed in, it must be unlocked before the unit is started up again. To do so, turn the switch.

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
Switching the unit on

The unit is switched on (position '1') and off (position '0') using the keyswitch Q on the front of the unit.

The device is ready for use when the keyswitch is in the '1' position.

! CAUTION: Never leave the unit unattended in activated condition. Secure the unit against unauthorized use by pulling the key from the keyswitch Q.


NOTE: The key can only be removed when in the '0' position.

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Once the laser is activated it begins a self-test. During the self-test, all display elements (7 segments and LED) are activated several times. Ensure that all the numbers and LEDs illuminate completely. This process takes approx. 10 seconds, during which time it is not possible to operate the unit.

The laser is equipped with a fiber recognition system, which provides a permissible parameter profile according to the employed fiber diameter.

The set of parameters is selected via the polarizing key (see Table 'Parameter matrix' on p. 27).

Connecting the fiber

You cannot operate the unit until a fiber is connected. The red LED { lights up when no fiber is connected.

To insert the fiber, open it, move the lever on the socket q to the left or right and simultaneously attach the fiber by turning it.

Once the fiber is connected, the red LED { goes out and you can operate the unit.

! WARNING: Only clean fibers in perfect condition that have been sterilized according to the instructions may be used for treatment. The reprocessing of the fibers must always be monitored.

i NOTE: The fibers are color coded on the anti-kink protection.

Fiber	Color
230 µm	blue
365 µm	green
600 µm	red

Only KARL STORZ CALCULASE II fibers (p. 47) may be used with this unit.

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The fiber must be checked at the beginning of every laser treatment. To do so, the aiming beam must be rechecked and reset every time you use it. The brightness of the aiming beam can be increased by pressing the '+' key ③ and decreased using the '-' key ④.

WARNING: As the AIMING BEAM follows the same path through the laser transmission system as the WORKING BEAM, it offers a good means of checking the integrity of the laser transmission system. If the AIMING DOT does not appear at the distal end of the laser transmission system, its intensity is too weak, or if it appears diffuse, this is a possible sign that the laser transmission system is damaged or not functioning cleanly.

To check the fiber, the aiming beam is projected onto a flat surface or the unit itself and the beam emitted is assessed.

G: Do not stare into the tip of the fiber. To prevent skin burns, do not aim the laser beam directly at unprotected skin.

CAUTION: Always check the aiming beam in the standby mode and not in the READY mode. The beam emitted by the laser fiber should be round and homogeneous.

WARNING: Do not activate the Ho:YAG radiation while setting the brightness of the aiming beam!

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
The energy is set using the '+' and '-' keys ^⑥ in the energy field. Up to five different energy levels are available here depending on the fiber and specified frequency: 0.5/0.8/1.2/1.7, and 2.0 J.

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The laser is always in standby mode while and after the parameters are set. The STANDBY LED P (yellow) lights up.

Changing to ready mode

To 'activate' the laser, press the ready/standby key ^⑧ once. The READY LED ^⑨ (green) lights up in the READY mode.


! NOTE: This procedure is accompanied by a three-time change between the yellow and green LED (STANDBY and READY LED).

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Activating the laser beam

The unit is in the READY mode. Press the footswitch to activate the laser beam.

! NOTE: The laser beam can only be activated once the target energy has been attained (takes approx. 1 second).

The activation/output of laser energy is displayed by the yellow LASER 'ACTIVE' LED ^⑪ and by an activation tone.

! WARNING: For safety reasons, the laser should be switched to standby mode during extended periods of (foreseen) non-activation.

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The CALCULASE II is equipped with a fiber recognition system. A permissible frequency/ energy combination is provided depending on the fiber used, from which the user can select the corresponding capacity (see also p. 23).

! **CAUTION:** The use of fibers which have not been authorized by KARL STORZ for use with the device may damage the fiber recognition system and thus lead to failure of the device.

Parameter matrix for 230 µm fibers					
Energy [J]	Frequency [Hz]				
	4	6	8	10	15
0.5	2.0 W	3.0 W	4.0 W	5.0 W	–
0.8	3.2 W	4.8 W	6.4 W	8.0 W	–
1.2	4.8 W	7.2 W	9.6 W	12.0 W	–
1.7	–	–	–	–	–
2.0	–	–	–	–	–

Parameter matrix for 365/600 µm fibers					
Energy [J]	Frequency [Hz]				
	4	6	8	10	15
0.5	2.0 W	3.0 W	4.0 W	5.0 W	7.5 W
0.8	3.2 W	4.8 W	6.4 W	8.0 W	12.0 W
1.2	4.8 W	7.2 W	9.6 W	12.0 W	18.0 W
1.7	6.8 W	10.2 W	13.6 W	17.0 W	–
2.0	8.0 W	12.0 W	16.0 W	20.0 W	–

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Removing the fiber

For safety reasons, switch the device to standby mode; to do so, press the ready/standby key .

Then remove the fiber.

WARNING: Disposable fibers must be disposed of and reusable fibers reprocessed.

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NOTE: You will find information on the sterilization of reusable fibers in the instruction leaflet accompanying the fibers.

NOTE: Prior to every sterilization of a reusable fiber, a number must be cut off the sterilization counter (from 10 to 0). The fiber should be changed when you reach ZERO at the latest.

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Fuse replacement

NOTE: The power switch (16) (rocker switch) acts as a line fuse. The fuses (17) fuse the 24 V module. Procedure (for non-functioning devices)
a) Check whether power switch is turned on.
b) Check fuses *.
To do so, turn off the device and disconnect it from the power supply.

Remove the fuse holder with a screwdriver.

CAUTION: Use only fuses with the correct nominal ratings.
Insert new fuses (see table).
Refit the fuse holder (17).

	100 ...240 V~
Fuse	2 x T10 AL, 250 V (5x20mm)

Reconnect the power. Perform function testing.

NOTE: If the device still does not function correctly after the fuse has been changed, please contact the Service department.

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Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Protective ground resistance measured according to EC 62353: For limit values, please refer to the current standard.
- Earth leakage current measured according to EC 62353: For limit values, please refer to the current standard.
- Touch current measured according to EC 62353: For limit values, please refer to the current standard.
- Patient leakage current measured according to EC 62353: For limit values, please refer to the current standard.

Test for proper functioning

The test for proper functioning must include the following tests at least:

- Emergency Stop
- Door interlock switch (remote interlock)
- Displays and control keys
- Checking all energy/frequency combinations of the laser parameters

Documentation

The safety inspection/repeat inspection and the results must be recorded in the medical device logbook.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during the safety inspection/repeat inspection, the unit must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range ^ and implementation of the safety inspections/ repeat inspections, please refer to the current version of the relevant Service Manual.

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We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installations in the room in which the unit is connected and operated meet applicable national and local electrical codes; and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Error message	Possible cause	Remedy
No function	Power plug not connected to outlet	Check power plug and, if necessary, connect to power supply.
	Emergency stop switch activated	Unlock emergency stop switch by turning it.
Red LED lights up.	Fiber not connected	Connect fiber or inspect it.
E 01 - 02	Technical error	Note error code and notify Service department.
E 03	Footswitch is not connected or is defective	Check and connect footswitch
E 04	Footswitch pressed when turned on	Release footswitch
E 05	Remote interlock active when turned on	Close remote interlock and turn on again.
E 06	Flash lamp aged - standard service life < 10 %	Confirm with the ready key to continue use, arrange a servicing date, and replace the flash lamp.
E 07	Flash lamp defective	Send unit in for servicing and have flash lamp replaced.
E 08	Key pressed when turned on	Release all keys and turn on again.
E 09 - 13	Technical error	Note error code and notify Service department.
E 14	Water temperature too low	Wait until heating-up phase is complete.
E 15	Water temperature too high	Allow device to cool in standby mode, check room temperature.
E 16 - 47	Technical error	Note error code and notify Service department.
E 26	Laser energy too low	Select another energy/frequency setting to continue use, note the error code, and notify the Service department.
E 29	Laser energy too high	Select another energy/frequency setting to continue use, note the error code, and notify the Service department.
I	Interlock in operation	Interlock in operation

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Specifications

Laser type:	Pulsed Ho:YAG laser
Wavelength:	2080 nm
Power:	Po = max. 20 Watts, (± 10 %) Pp = max. 10 kW
Pulse frequency:	4 / 6 / 8 / 10 / 15 Hz
Pulse duration:	100...500 µs
Energy:	0.5 / 0.8 / 1.2 / 1.7 / 2.0 J (max. 2 J ± 20 %)
Beam transmission:	Fiber system
Fiber types:	Bare fiber, 230 µm, 365 µm, 600 µm fiber
Fiber connection:	Fiber recognition via modified SMA fiber connector
Aiming beam:	4 mW, 532 nm
Operating voltage:	100...240 VAC (± 10 %)
Line frequency:	50/60 Hz
Power consumption:	Max. 1800 VA
Cooling:	Integrated water cooling system with water-/air- heat exchanger
Dimensions:	610 x 250 x 435 mm
Weight:	35 kg (incl. 2 kg cooling water)
Storage/transport conditions:	
Temperature:	-10 °C ... 60 °C
Air humidity:	5 % ... 95 % (RH, non-condensing)
Operating conditions:	
Temperature:	18 °C ... 28 °C
Air humidity:	30 % ... 70 % (non-condensing)
Atmosph. pressure:	700...1060 hPa

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Standard compliance
(for 27750220-1)

According to IEC 60601-1, IEC 60601-2-22:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type B ⤴
- Type of protection against moisture: drop-water protection as per IPX 1

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 **NOTE:** We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof

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Spare parts, accessories

Article	Cat. no.
Laser safety goggles Ho:YAG Laser 2080	27 750095
Safety goggles surgical LASER, for users with and without glasses	UL020
Safety goggles surgical LASER, sporty version for users who don't wear glasses	UL021
Reusable fiber sets CALCULASE II Bare Fiber 230 µm, reusable, sterile, package of 6 pieces CALCULASE II Bare Fiber 365 µm, reusable, sterile, package of 6 pieces CALCULASE II Bare Fiber 600 µm, reusable, sterile, package of 6 pieces	27 750271-P6 27 750272-P6 27 750273-P6
Single-use fiber sets CALCULASE II Bare Fiber 230 µm, single use, sterile, package of 6 pieces CALCULASE II Bare Fiber 365 µm, single use, sterile, package of 6 pieces CALCULASE II Bare Fiber 600 µm, single use, sterile, package of 6 pieces	27 750277-P6 27 750278-P6 27 750279-P6
Warning sign 'Laser radiation'	27 750092
Set of keys	27 750235
Remote interlock jumper	27 750226
Fiber stripping set	27 750280
Footswitch	27 750124

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The CALCULASE II model 27 7502 20-1 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The CALCULASE II model 27 7502 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
Footswitch with cable	Yes	18	No	Activating the unit functions
Power Cord 10A with non-heating appliance connector	No	1.5	No	Power supply

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The CALCULASE II model 27750220-1 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCULASE II model 27750220-1 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The CALCULASE II model 27750220-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The CALCULASE II model 27750220-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The CALCULASE II model 27 7502 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCULASE II model 27 7502 20-1 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 8 kV air	Complies ± 8 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using an AC line filter.
Surge IEC 61000-4-5	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	Complies ± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T^*$ (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles $<5\% U_T$ (>95 % dip in U_T) for 5 sec	Complies $<5\% U_T^*$ (>95 % dip in U_T) for 0.5 cycle Complies 40 % U_T (60 % dip in U_T) for 5 cycles Complies 70 % U_T (30 % dip in U_T) for 25 cycles Complies $<5\% U_T$ (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The CALCULASE II model 27 7502 20-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the CALCULASE II model 27 7502 20-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V and 10 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the CALCULASE II model 27 7502 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance: $d = [3.5 \cdot 3] \sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5 \cdot 3] \sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7 \cdot 3] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m and 10 V/m	

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 6

**Recommended separation distances between portable and mobile RF communications equipment and the
CALCULASE II model 277502 20-1**

The CALCULASE II model 277502 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Midcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 616-4500, Fax: +1 905 658-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
St. Louis, MO 63103-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 421-1204 (USA only)
E-Mail: info@kasea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 686-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latin-America, Inc.
815 N. W. 57th Avenue, Suite 430
Miami, FL 33126-2042, USA
Phone: +1 305 262-6360, Fax: +1 305 262-5966
E-Mail: info@kasea.com

KARL STORZ Endoscopy
Miramar Trade Center
Edmond Jerusalem, Oficina 108,
La Hadrana, Cuba
Phone: +53 72041097, Fax: +53 72041099

KARL STORZ Endoscopy México S.A. de C.V.
Lago Comanzapa No 326, Col. Chapultepec Morales
D.F.C.P. 1520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

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KARL STORZ GmbH & Co. KG
Resident Representative Office
80/23 (44/19) Dang Van Ngu,
F10-0, Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +84 991 5442, Fax: +84 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 16F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5626186, Fax: +86 10 5628199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1, Grand Gateway
No. 1 Hongqiao Road, Xuhui District
Shanghai, 200030, People's Republic of China
Phone: +86 21 6112-1188, Fax: +86 21 6112-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai (Zhangjiang)
Room 701A, 7F Building 5, No. 3000 Longdong
Avenue, Zhangjiang Hi-Tech Park, Pudong District
Shanghai, 201203, People's Republic of China
Phone: +36 21 61009688, Fax: +86 21 61009686
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610014, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shente District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 028, 25F Team Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 3 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 3 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
1-Perk Bundang 1, 102dong, 8F
239 Jeongjil-ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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Karl Storz GmbH & Co. KG