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Apostille

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**Important information
for users of KARL STORZ
devices**

**1 Important information
for users of KARL STORZ
devices**

Thank you very much for your expression of trust in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the UNIDRIVE® S III ARTHRO motor system. The motor system is used with the following application instruments:

- Standard shaver handpiece (8K) 28 7210 35 or 28 7210 45 (without head buttons)
- Hybrid shaver handpiece (7K) 28 7210 50
- High-speed shaver handpiece (15K) 28 7210 30
- Small joint shaver handpiece (6K) 28 7210 40
- Multifunction handpiece 28 7210 36 for use with optional attachments 28 7210 31/37/38/39/33/34

Please read these instructions carefully. Keep this manual in a convenient and conspicuous location and in its protective casing close to the device.

Users of the various application instruments should familiarize themselves with the operating instructions included with the respective devices.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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


**2 Controls, displays,
connectors, and their uses**

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Controls, displays, connectors, and their uses

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3 Controls, displays, connectors, and their uses

- ① Power switch for turning the unit on and off
- ② TFT touch screen
- ③ Connection socket for four-pedal footswitch
- ④ Connection socket for dual-pedal/triple-pedal footswitch
- ⑤ Connection socket for multifunction handpiece
- ⑥ Connection socket for shaver
- ⑦ Drawer for operator's quick reference guide
- ⑧ USB port (for maintenance only)
- ⑨ Potential equalization connection
- ⑩ Power socket
- ⑪ Fuse holder
- ⑫  connection sockets*
- ⑬ Dual-pedal footswitch
- ⑭ Button for activating the function which enables the speed to be reduced/increased using the footswitch.

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*** NOTE:** The KARL STORZ SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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1	Important information for users of KARL STORZ devices	III
2	Images of the equipment	IV
3	Controls, displays, connectors, and their uses	V
4	Symbols employed	VI
5	General information	4
5.1	Description of the system	4
5.1.1	Combination of UNIDRIVE® S III ARTHRO and ARTHROPUMP® Power	5
5.2	Special features	6
6	Safety instructions	7
6.1	Explanation of warnings and cautions	7
6.2	Intended use	13
6.3	User qualifications	13
6.4	Safety precautions at the site of installation	14
6.5	Safety precautions when operating the system	14
7	Operating instructions	16
7.1	Unpacking	16
7.2	Basic equipment	16
7.3	Installing and connecting the unit	17
7.4	KARL STORZ SCB	18
7.5	Connecting the footswitch	18
7.6	Connecting the multifunction handpiece/shaver handpiece	19
7.7	Attaching the blade to the shaver handpiece (6K, 8K, 15K)	19
7.8	Connecting the suction tube (6K, 7K HYBRID, 8K, 15K)	20
7.9	Attaching the blade to the 7K shaver handpiece	21
7.10	Connecting the suction tube (7K)	21
7.11	Operation of the handpieces (6K, 8K, and 15K)	22
7.12	Handpiece button operation (7K shaver handpiece)	23
7.13	Small joint shaver handpiece (6K)	24
7.14	Shaver handpiece 28721045	24
7.15	Footswitch operation	25
7.15.1	Dual-pedal footswitch	25
7.15.2	Triple-pedal footswitch	26
7.15.3	Four-pedal footswitch	26

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7.16	Multifunction handpiece (drill/saw)	27
7.17	Functional check of the shaver blades before renewed application	28
7.18	Initial operation	29
7.19	UNIDRIVE® S III function	29
7.20	Selecting the language	30
7.21	Selecting the motor connection	31
8	Quick overview of screen functions	32
9	Mode: Settings	34
9.1	System information	34
9.2	Language	35
9.3	Date/Time	35
9.4	Service	35
9.5	Error messages	36
10	Maintenance	37
10.1	Changing the fuse	37
10.2	Cleaning, disinfection and sterilization	38
10.2.1	Cleaning and disinfection instructions for the shaver handpieces and reusable blades	38
10.2.2	Manual cleaning and disinfection	40
10.2.3	Machine cleaning and disinfection	42
10.2.4	Sterilization	43
10.2.5	References for cleaning and sterilization	49
10.3	Maintenance and safety check	50
10.3.1	Maintenance	50
10.3.2	Safety check	50
10.4	Service and repair	52
10.5	Disposal	52
10.6	Repair program	53
10.7	Important information	53
10.8	Limitation of liability	54
10.9	Manufacturer's warranty	54
11	Technical description	55
11.1	General	55
11.2	UNIDRIVE® error table	57
11.3	Troubleshooting	59
11.4	Technical Data	60
11.5	Intended handpieces	61
11.6	Standard conformity	63
11.7	Directive compliance	63
11.8	Technical documentation	64
11.9	Block diagram	65

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12 Spare parts, recommended accessories	66
12.1 Spare parts	66
12.2 Recommended accessories	67
12.3 Recommended shaver blades and burrs	68
12.4 Recommended shaver blades and burrs for small joint handpiece	80
12.5 Accessories for shaver handpieces and SL blades	83
12.6 Accessories for pin driver	83
12.7 Accessories for sagittal saw	83
13 Single-use blades and burrs	84
13.1 Intended use	84
13.2 Use of the single-use products and explanation of the symbols employed	84
14 Electromagnetic Compatibility (EMC) Information	86
15 Subsidiaries	100

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5.1 Description of the system

The UNIDRIVE® S III ARTHRO can be used with:

- Standard shaver handpiece (8K) 28721035
- Standard shaver handpiece without head buttons (8K) 28721045
- Hybrid shaver handpiece (7K) 28721050
- High-speed shaver handpiece (15K) 28721030
- Small joint shaver handpiece (6K) 28721040
- Multifunction handpiece 28721036

Dual-pedal footswitch 20016831 or Three-pedal footswitch 20012832

It is used in conjunction with a shaver blade or burr (see Accessories page 66-81) to perform surgical interventions under endoscopic control in the field of arthroscopy.

With the corresponding blade, the system is designed to provide controlled cutting, shaving, and abrading of bone and tissue during arthroscopic surgical procedures of the knee, shoulder, ankle, elbow, wrist, hip, and TMJ.

The ergonomically designed and lightweight shaver handpiece has a suction channel for connection to a suction system. Using a rocker lever, the suction channel can be opened or closed.

The easily visible LCD/touch screen on the control unit of the UNIDRIVE® S III ARTHRO provides information on the most important operating parameters. The touch screen provides control of the operating parameters.

With the optionally available multifunction handpiece and its accessories, such as the drill attachment and the saw attachment, the UNIDRIVE® S III ARTHRO becomes a multifunctional instrument that can be used in a wide range of areas in orthopedics.

The control unit automatically recognizes the shaver handpieces or the multifunction handpiece as soon as they are connected.

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5.1.1 Combination of UNIDRIVE® S III ARTHRO and ARTHROPUMP® Power

Either a dual-pedal footswitch or the combined suction and irrigation pump ARTHROPUMP® Power 28340720 can be connected to the right, lower connection socket ④ (see page IV).

If an ARTHROPUMP® Power is connected and a shaver is connected to position ⑥ on the UNIDRIVE® S III ARTHRO and turned on, the shaver sends the signal 'Shaving in progress' to the ARTHROPUMP® Power (increased suction). The tube squeeze valve of the ARTHROPUMP® Power then switches from the function 'Sheath' to the function 'Shaver suction'.

Please observe the instruction manual for the ARTHROPUMP® Power 96156010DER.

The four-pedal footswitch 20016834 is used if the suction and irrigation pump 'ARTHROPUMP® Power' is connected to the socket ③ on the UNIDRIVE® S III ARTHRO via a connection cable. It is then inserted on the rear of the ARTHROPUMP® Power.

The two left buttons (gray) perform the same function as the dual-pedal footswitch (counterclockwise/clockwise rotation).

The blue pedal triggers the lavage function. When the inflow is closed, pressing the pedal activates the suction function.

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5.2 Special features

- The shaver handpieces (6K, 7K, 8K and 15K) have two modes of operation: 'Normal' and 'Oscillation'.
- The rotational speed of the small joint shaver handpiece (6K) ranges from 1,000 rpm to 6,000 rpm.
- The rotational speed of the hybrid shaver handpiece (7K) ranges from 1,000 rpm to 7,000 rpm.
- The rotational speed of the standard shaver handpiece (8K) ranges from 1,000 rpm to 8,000 rpm and the high-speed shaver handpiece (15K) can be used from 1,000 rpm to 15,000 rpm*.
- The rotational speed of the standard, high-speed, hybrid and small joint shaver handpiece in oscillation mode ranges from 500 rpm to 3,000 rpm.
- The speed and mode of the shaver handpieces may be controlled by either the touch screen, the handpiece (except 28721040 and 28721045) or the footswitch.
- Fast and simple blade change using a precise, quick-release coupling.
- UNIDRIVE® S III ARTHRO is equipped with the KARL STORZ Communication Bus (SCB).
- The control unit automatically recognizes the shaver handpiece and the multifunction handpiece connected in parallel.

* **Note:** The recommended maximum speed is 4,000 rpm, or 1,500 rpm in oscillation mode.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

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6.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.

CAUTION: A Caution indicates that particular maintenance procedures or safety precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information

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WARNING: Read this instruction manual thoroughly before operating the system. Read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

WARNING: Do not operate the handpiece unless the blade/burr is visible.

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WARNING: Continuous checks must be performed both before and during use to establish whether the required channel and/ or mode is set.

WARNING: The UNIDRIVE® S III ARTHRO may only be used by persons with an appropriate medical qualification and who are acquainted with the technique.

WARNING: When replacing fuses, use only fuses of the same rating.

WARNING: The device may only be operated at the voltage stated on the rating plate.

WARNING: If the motor does not stop as required, switch the unit off.

WARNING: Store/keep the instrument at a safe distance from the patient when not in use.

WARNING: The UNIDRIVE® S III ARTHRO must only be used in medical facilities that conform to applicable national, state and local electrical codes.

WARNING: Keep out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/ or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only safe if:

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

WARNING: Connect the device to a properly grounded 'Hospital grade' or 'Hospital use' electrical outlet. Check grounds regularly. Routinely inspect the power cord and plug and do not use if the plug or the cord is damaged.

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-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.
-  **WARNING:** Do not use in the presence of flammable anesthetics. There is a risk of explosion.
-  **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.
-  **WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional raw surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.
-  **WARNING:** Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication of no preparation or incorrect preparation. Prepare the instruments before initial use and before and after every subsequent use using validated procedures.

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WARNING: Risk of injury: Incorrect application of medical instruments or bending of these instruments by applying excessive side load force poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: When using applied parts outside of the field of vision there is a risk of unintentionally injuring tissue and damaging accessories. Always keep the applied parts of the active electrode and the laser and other instruments which transmit energy pointed at the target in the field of vision during use.

WARNING: The UNIDRIVE® S III ARTHRO, in conjunction with the multifunction handpieces/shavers used, is only designed for intermittent operation, not for continuous operation. (See Technical Data on p. 60). For this reason, always turn the motor off again, as it will otherwise become too hot.

WARNING: Risk of infection: Incorrectly prepared instruments expose patients and users to a risk of infection as well as the danger of an instrument malfunction. Observe the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

NOTE: The 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions can be downloaded or requested at www.karlstorz.com.

WARNING: When using HF surgical devices, the shaver handle and shaver blade must be put out of reach of the patient.

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WARNING: Always unplug the unit before performing any maintenance, including cleaning.

WARNING: To avoid the risk of electrical shock, do not open the device. Refer servicing to the manufacturer or qualified personnel that are authorized by the manufacturer. Removal of the cover by non-authorized personnel will void the warranty.

CAUTION: The UNIDRIVE® S III ARTHRO may only be used with accessories, and disposable items which have been designed for use by KARL STORZ.

CAUTION: Do not store liquids on or above the device. Do not allow fluids to get into the unit.

CAUTION: If the motor becomes blocked for any length of time, it can lead to overheating and cause damage to the motor. If the blockage cannot be disengaged using oscillating operation (if possible), the motor must be shut off and the blockage must be removed manually.

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- ❗ **NOTE:** The USB port fitted with a cover is intended for use by Service.
- ❗ **NOTE:** Before using the unit, a check must be performed to establish whether the touch screen and the unit display are functioning correctly.
- ❗ **NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.
- ❗ **NOTE:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices. Current lists of chemicals can be found on the Internet at www.karlstorz.com under 'Hygiene'.
- ❗ **NOTE:** If the unit is not functioning as required, please consult the display or contact the KARL STORZ Service.
- ❗ **NOTE:** Limits can be adjusted when working with RUI (= the SCB software). When working with RUI, please check the parameters before and after use.
- ❗ **NOTE:** This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your distributor for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Before using the UNIDRIVE® S III ARTHRO motor system on a patient, it is imperative that the instructions in this manual are understood and are followed carefully.

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6.2 Intended use

When used in arthroscopy/orthopedics, the motor control unit serves to regulate and control motors or handpieces and accessories during orthopedic, arthroscopic and minimally invasive interventions for the excision, removal, ablation and transection of soft tissue and/or bones.

Use of the UNIDRIVE® S III ARTHRO in fields other than those indicated above is not permitted.

The UNIDRIVE® S III ARTHRO may only be used with accessories, replacement parts, and disposable items which have been designated as suitable for the unit by KARL STORZ.

Unauthorized conversions or modifications to the system are not permitted.

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6.3 User qualifications

The UNIDRIVE® S III ARTHRO motor system in conjunction with an application instrument may only be used by physicians and/or medical assistants who have been instructed in the use of this device/instrument.

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6.4 Safety precautions at the site of installation

The UNIDRIVE® S III ARTHRO motor system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The UNIDRIVE® S III ARTHRO is equipped with a connector for attaching a ground wire. Connect this in accordance with the state, local and national electrical codes.

The UNIDRIVE® S III ARTHRO motor system is not intended for use in hazardous areas. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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
6.5 Safety precautions when operating the system

WARNUNG: The system must be tested for operational safety and proper condition prior to each surgical procedure.

During surgery with the UNIDRIVE® S III ARTHRO, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

CE marked equipment has been tested and is found to comply with the EMC limits of IEC 60601-1-2:2001 CISPR 11 Class B.

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WARNING: When using HF surgery, the shaver handpiece and shaver blade must not come in contact with the patient.

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7 Operating instructions

7.1 Unpacking

Carefully remove the unit and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can be used if the unit has to be transported.

7.2 Basic equipment

- 1 UNIDRIVE® S III ARTHRO 28723020-1
- 1 Power cord 400 A
- 1 SCB connecting cable, length 100 cm, 200901 70
- 1 Instruction manual 96156013DER

i NOTE: Optional accessories are described in section 'Recommended accessories'.

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7.3 Installing and connecting the unit


NOTE: The UNIDRIVE® S III ARTHRO including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

! WARNING: The UNIDRIVE® S III ARTHRO is not intended for use in hazardous areas. Do not operate the system within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
Place the UNIDRIVE® S III ARTHRO on a level surface.

The UNIDRIVE® S III ARTHRO is equipped with a connector ⑨ for attaching a ground wire. The unit's ground should be installed by a qualified electrician.

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Before connecting the UNIDRIVE® S III ARTHRO to the power supply line, verify that the line voltage stated on the unit's identification plate conforms with that of the power supply lines to be used.


Connect the power, insert the power cord end into the receptacle ⑪ of the unit.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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
7.4 KARL STORZ SCB

ⓘ NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector is equipped with a protection mechanism. Pull the protection mechanism of the SCB connector and insert the plug into one of the SCB sockets ⑫. Connect the other end of the cable to a KARL STORZ SCB control unit (KARL STORZ Communication Bus) or other SCB devices (see instruction manual KARL STORZ SCB System).

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
7.5 Connecting the footswitch

To use the footswitch ⑬, insert the connecting cord* of the footswitch into the socket for the footswitch (③ or ④) on the front of the UNIDRIVE® S III ARTHRO unit. Align the red dot on the connecting cord with the red dot on the footswitch socket.

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7.6 Connecting the multifunction handpiece/shaver handpiece

Prepare the multifunction handpiece/shaver handpiece as described in the instruction manual.

NOTE: Please follow the instruction manual enclosed with the multifunction handpiece/shaver handpiece!

Connect the connecting cord from the multifunction handpiece/shaver handpiece to the connection socket ⑤ (multifunction handpiece) and/or ⑥ (shaver).

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Align the red point on the connecting cord with the red point on the connection socket to guarantee a proper connection.


NOTE: If no 'motor' is connected when the unit is switched on, the message: 'E12: NO MOTOR' flashes in the status line of the display.

NOTE: To unplug the connecting cable, pull out the plug — do not pull directly on the cable!

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7.7 Attaching the blade to the shaver handpiece (6K, 8K, 15K)

NOTE: Insert the blade either with direction setter, and thus with the blade opening facing upwards, or rotated through 180°, with the opening facing downwards.

Insert the blade with the direction setter into one of the two notches on the lock.

Turn the lock on the lever counterclockwise and hold it.

Slide the blade **as far as it will go** (see arrow) into the lock.

Finally, let go of the lever and make sure that the blade is securely attached by pulling it.

NOTE: If disposable products are used, follow the instructions on pages 84 ff.

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

NOTE: Removing the blade from shaver handpiece: Turn the lock on the lever counterclockwise and hold in this position.

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7.8 Connecting the suction tube (6K, 7K HYBRID, 8K, 15K)

Attach the suction tube to the connection piece 15 on the handpiece and connect the tubing to a suction pump (ARTHROPUMP® Power, for example).

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
! CAUTION: During shaver activation the handpiece needs a certain amount of flow for the cooling of the handpiece. Make sure that the irrigation/suction system is properly set up and that tubes are not occluded or kinked.

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Suction is controlled using the rocker lever 16. If the rocker lever is at the rear limit stop, the suction channel is closed; at the front limit stop, the suction channel is open (maximum suction).

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7.9 Attaching the blade to the 7K shaver handpiece

Press the button on the side of the handpiece collar and insert the selected shaver blade into the receptacle on the handpiece until the shaver blade clicks into place. Be sure to align the notch on the blade with the corresponding notch on the handpiece.

i NOTE: Removal of Blade from the 7K Shaver Handpiece: Release the shaver blade by pressing the button on the side of the handpiece collar and simultaneously removing the blade.

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7.10 Connecting the suction tube (7K)

Attach the suction tube to the connection 12 on the handpiece and connect the tubing to a suction pump or wall suction. Suction is controlled using the rocker lever 10. If the rocker lever is at the rear limit stop, minimum suction is provided, at the front limit stop, maximum suction is provided.

! CAUTION: The minimum suction position restricts or eliminates fluid movement. Operating the handpiece with minimal fluid movement can result in excess heat retention and damage to the unit.

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
7.11 Operation of the handpieces (6K, 8K, and 15K)

The Unidrive® S III ARTHRO system features handpieces which can be operated using the LCD/touch screen of the Unidrive® S III ARTHRO, the handpiece buttons (except 28721040 and 28721045) or the footswitch. The handpieces 28721030 and 28721035 are equipped with three buttons, whose function is identical on both types of handpiece. Below is a definition of the button functionality:

Button 1: right (Oscillation) (+)

When the handpiece has stopped, pressing this button will change the handpiece to Oscillation mode. Each time the button is pressed again, the blade will oscillate for one second.

When the handpiece is running, pressing button 1 will increase the speed of the blade.

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Button 2: left (Forward/Reverse) (-)

When the handpiece has stopped, pressing button 2 will change the direction of the blade rotation from forward to reverse.

When the handpiece is running, pressing button 2 will decrease the speed of the blade.

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
7.12 Handpiece button operation (7K shaver handpiece)

The Unidrive® S III ARTHRO system features a fourth shaver handpiece (7K), which also can be controlled via the touch screen/LCD, handpiece buttons, or footswitch. The 7K handpiece is equipped with two sets of buttons. The first set consists of three buttons as depicted on the left. These buttons are located at the top of the handpiece. The second set consists of a single button located at the bottom of the handpiece. This button is used to start and stop the handpiece.

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With the handpiece stopped, pressing this button will change the direction of the blade rotation from forward to reverse. When the handpiece is running, pressing this button will decrease the speed of the blade.

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
When the handpiece is stopped, pressing this button will switch to oscillation mode. Each time the button is pressed again, the blade will oscillate for one second. When the handpiece is running, pressing this button will increase the speed of the blade.

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7.13 Small joint shaver handpiece (6K)

The small joint shaver handpiece is operated using either the LCD/touch screen of the Unidrive® S III ARTHRO system or the footswitch. It only has one button, with the following function:

Button 1:


This button starts and stops the handpiece.

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7.14 Shaver handpiece 28721045

The shaver handpiece 28721045 is operated using either the LCD/touch screen of the Unidrive® S III ARTHRO system or the footswitch.

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
7.15.2 Triple-pedal footswitch

The triple-pedal footswitch 20016832 has three distinct pedals for mode selection and a button for speed control. Below is a definition of the triple-pedal footswitch functionality:

Left pedal:

Reverse (Left) < Counterclockwise

Middle pedal:

Oscillate

Right pedal:

Forward (Right) > Clockwise

Button on top of footswitch:

Pressing the button on the top of the footswitch once switches the footswitch into speed adjust mode; the foot pedals can now adjust the speed up and down.

Pressing the right foot pedal increases the speed (by the programmed rate) and pressing the left foot pedal decreases the speed (by the programmed rate).

By pressing the top button again, the unit returns to normal operating mode with the new speed setting.

7.15.3 Four-pedal footswitch

The four-pedal footswitch 20016834 is used if the suction and irrigation pump 'ARTHROPUMP® Power' 28340720-1 is connected to the socket ④ on the UNIDRIVE® S III ARTHRO via a connection cable. It is then inserted into the rear of the ARTHROPUMP® Power.

The two left buttons (gray) perform the same function as the dual-pedal footswitch (counterclockwise/clockwise rotation). The blue pedal triggers the lavage function. When the inflow is closed, pressing the pedal activates the suction function.

The red pedal triggers the hemostasis mode (temporary increase in pressure). After a predetermined period of time x has elapsed, but after no more than 2 minutes, the pump returns to normal mode. Pressing the red pedal again during hemostasis mode deactivates this mode.

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7.16 Multifunction handpiece (drill/saw)

The Unidrive® S III ARTHRO control unit is capable of driving a multifunction handpiece. This handpiece may be attached to drill chucks, sagittal saws, etc. A list of available accessories is shown in a separate section of this manual.

Each accessory is connected to the handpiece by attaching the accessory and handpiece together until solid contact is made. Press the rear of the button on the top ②④ of the handpiece to remove the accessory.

The motor is cannulated and accepts pins up to 0.125 inches (3.2 mm) in diameter.

The Unidrive® S III ARTHRO allows the speed of the handpiece to be controlled by moving the trigger ②③. The speed (rpm) of the handpiece is proportional to the position of the trigger. Pushing the trigger all the way in toward the handpiece provides the maximum speed.

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








7.17 Functional check of the shaver blades before renewed application

WARNING: Risk of injury to the patient and/or user: Before every renewed application of the shaver blades, it is compulsory to carry out a functional check. This includes the following tests:

- The inner blade should be free to rotate easily within the outer blade and must not show any run-in marks (negative example, see photos 1 and 2).
- The inner blade must lock securely into place in its front end position, in order to ensure that it is not unintentionally lost.
- The distal cutting edges must be free from defects (negative example, see photo 3).
- The outer blade of the burr must not show any points of abrasion (negative example, see photos 3 and 4).

Always operate the shaver/shaver blades with running suction

WARNING: Risk of injury to patient and damage to property:

It is strongly recommended never to operate the shaver/shaver blades without running suction. Suction of the irrigation fluid serves to cool the running surfaces and the shaver motor and prevents injuries/damage caused by overheating and mechanical "seizure". To this end, the rocker lever used for setting the suction power must at least be in a vertical position (middle position) (see photo 5).

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
7.18 Initial operation


Switch the UNIDRIVE® S III ARTHRO on at the power supply ① (power switch lights up).

NOTE: Prior to turning the device on, the user must connect a multifunction handpiece and/or shaver handpiece to the unit. This is checked by the software and if necessary (if nothing is connected) indicated by a corresponding error message in the top image margin.

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The screen shown on the left appears for the time it takes the application to start.

NOTE: In addition, the following pages describe the simplest application case as an example, as well as the adjustment of the speed. A single signal tone can be heard during the switch-on procedure.

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7.19 UNIDRIVE® S III function


A start screen is displayed once the unit has started. The function only becomes active once the multifunction handpiece and/or shaver handpiece has/have been connected to the unit.

CAUTION: In order to avoid damaging the delicate touch screen interface, do not use any objects to touch the touch screen. The screen displayed shows the start screen with a shaver handpiece connected.

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NOTE: Error messages are shown in text in the status line of the operation screen if necessary (flashing text).

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
7.20 Selecting the language

It may be required to define a different language for the operation screen since the device is set to English as the default in the factory. This is done by pressing the 'menu' symbol (wrench) in the bottom left.

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Activate the 'LANGUAGE' button (2nd button from the top).

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X = Aktion abbrechen

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Select the desired language*.

Select one of the ten available languages by pressing the respective language button. The languages currently available are displayed over two screens.

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7.21 Selecting the motor connection

The handpieces described above (which must be operated within various speed ranges) can be connected up to the control unit. To do so, the required instrument must be connected to the corresponding motor connection socket (connect multifunction handpiece up to upper left socket ⑤, connect shaver handpiece up to upper right socket ⑥).

Having detected the handpiece automatically, the UNIDRIVE® S III ARTHRO loads the standard parameters and displays them on the screen. The symbol for the respective instrument is displayed in the bottom right of the screen.

It is only possible to switch between the individual motor connection sockets when both channels are occupied.

The switching is performed via the forward/backward lever on the multifunction handpiece. In the 'OFF' setting, the shaver handpiece is active, otherwise the multifunction handpiece is active.

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
**8 Quick overview of screen
functions**

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1 Name of the currently active instrument (multi-function handpiece/shaver handpiece). Error message E12 appears ('No motor') if neither a multifunction handpiece nor a shaver handpiece is connected.

If both a shaver handpiece and a multifunction handpiece are connected to the unit, the multifunction handpiece is automatically active unless the switch on the lower end of the multifunction handpiece is set to OFF. Then the shaver handpiece is active.

2 Currently set speed limit, i.e., the maximum number of revolutions currently possible.

The speed limit of the currently active instrument is automatically recognized and displayed by the unit.

3 Increase/reduction of speed limit in increments dependent on instrument, e.g., 100 rpm for multifunction handpiece.

4 Display of speed range possible for the currently active instrument.

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Quick overview of screen functions

5 Footswitch symbols

If the multifunction handpiece is active, the + or - sign appears on the footswitch symbols.

This means that the speed limit 2 on the connected footswitch can be increased/reduced. (For shaver handpieces, the yellow button on the dual-pedal footswitch must be pressed first).

6 Symbols for the forward rotation / oscillation / reverse rotation modes (see also pages 22 / 23)

Currently active mode = orange.

The 'Oscillation' symbol only appears if a shaver handpiece is currently activated (e.g., when the name of a shaver handpiece is in the header 1).

It is possible to switch between the forward and reverse rotation modes

- via the forward/reverse head buttons of the shaver handpiece
- via the switch at the lower end of the multifunction handpiece (FWD – forward, REV – reverse)

The shaver handpiece can be switched to the oscillation mode via the 'Oscillation' head button.

It is only possible to switch modes via the instruments, not the touch screen.

7 This button can be used to enter the 'Settings' menu. This area allows you to change the screen language (e.g., from English to German) among other functions.

8 The currently active instrument (multifunction handpiece/shaver handpiece) is displayed here in orange. If both instruments are connected to the unit, the multifunction handpiece is automatically active unless the switch on the lower end of the multifunction handpiece is set to OFF.

9 Activate/Deactivate the head buttons of the shaver handpiece (orange = active). This button can thus also be used to turn off a running shaver handpiece.

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9 Mode: Settings

9.1 System information

To call up the menu functions, press the 'Menu' button (wrench) in the bottom left.

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
This opens the 'MENU' screen (image on the left). Press the 'SYSTEM INFORMATION' button to obtain detailed information on your UNIDRIVE® S III ARTHRO.

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The device number (serial number) is displayed at the top.

The operating hour counter and the software version statuses are displayed in the lower section.

Then click the tick in the bottom right (= 'Confirm action' button; in this case, the 'Confirm action' triggers the operation).

Activating the 'OPERATING MODE' button returns you to the operation screen, from which the MENU was called. This applies for every menu level.

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Siehe Seite 30.

9.2 Language

Here you can set the language to be displayed on the screen.

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9.3 Date/Time

Here you can set the date and time, which appear at the top of the screen.

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
9.4 Service

This menu is reserved for the authorized Service and, for this reason, access to this menu is password-protected.

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
The upper section displays the current date and time; the lower section displays a list of all the error messages that have occurred in order of importance.

E01-E10: Hardware error messages which are potentially caused by a unit fault.

E11-E30: User error messages, which are due to incorrect configuration of the unit (e.g. if neither a multifunction handpiece nor shaver handpiece is connected). See Error table on p. 57.

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
The example on the left displays an error message:

E 12: NO MOTOR
(No motor/handpiece connected)

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

10 Maintenance

10.1 Changing the fuse

Switch the control unit off and disconnect the power.

Remove the power fuse holder ⑪ with a screwdriver or another suitable tool.

! CAUTION: Use only fuses with the correct nominal ratings.

Insert new fuses.
Insert the power fuse holder again.

	100...240V
Line fuse	2xT4.0 AH 250 V

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Reconnect the power.
Perform function testing.

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10.2 Cleaning, disinfection and sterilization

NOTE: The instructions 'Cleaning, disinfection, care and sterilization of KARL STORZ instruments' must be followed. The procedures for cleaning, disinfection and sterilization are explained there in detail.

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UNIDRIVE® SIII ARTHRO

WARNING: Risk of fatal injury from electric shock: Unplug from power supply before every cleaning!

CAUTION: Damage to the control unit by liquid: Avoid allowing liquid to enter the control unit.

The design of the device allows only for mechanical wipe-down disinfection. The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

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10.2.1 Cleaning and disinfection instructions for the shaver handpieces and reusable blades

WARNING: Exposure prophylaxis: Improperly prepared handpieces and reusable blades/burns pose a health risk to both the patient and the staff. Always clean, disinfect and sterilize using validated preparation procedures!

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* The use of a highly alkaline cleaning agent pH > 10.5 (in accordance with the Robert Koch Institute or applicable national requirements) may lead to cosmetic color changes to the coating, which, however, do not influence the functionality of the instruments. It is thus preferable to use mildly alkaline or enzymatic cleaning agents.

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* See far left column

** Shaver blades marked as disposable must not be re-prepared.

*** Optional accessories for validated cleaning: Irrigation adaptor 29205SA for shaver handpieces Irrigation adaptor 28205SAA and 28205SAI for cleaning of reusable blades (Inner and outer blades)

NOTE: Use a chemical agent approved by KARL STORZ for cleaning and disinfection. The current list of KARL STORZ approved materials can be found on the Internet at 'www.karlstorz.com'.

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10.2.2 Manual cleaning and disinfection

10.2.2.1 Manual Precleaning

WARNING: Observe the guidelines of the Employer's Liability Insurance Association and equivalent organizations for all cleaning and disinfection of contaminated instruments.

CAUTION: When preparing and using solutions for manual cleaning and disinfection, follow the chemical manufacturer's instructions closely.

Place the instruments in a container with cleaning solution immediately after use (as per the manufacturer's instructions). This prevents impurities on the instruments from drying out.

Disassembly for preparation/Manual cleaning and disinfection

Shaver handpiece/Multifunction handpiece and accessories

Remove blade/burr from the handpiece. All moving parts must be in the open position. Remove the suction tube from the handpiece. Remove the accessory (e.g. Jacobs chuck) from the multifunction handpiece. Clean the handpiece, including the suction channel and the area where the blades/burrs connect into the handpiece, mechanically with water or a cleaning and disinfectant solution (see note, page 39) and a soft brush (Art. no. 27650 E). Clean the accessory for the multifunction handpiece (e.g. Pin Driver) with a sponge and brush. Do not use hard or sharp cleaning tools. In order to avoid material damage, carry out the final flushing with microbiologically correct/sterile water. Finally dry with filtered compressed air (residual moisture can lead to operational malfunctions).

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Rinsing

If the shaver is given to nursing staff or laid down, the attachment should be briefly rinsed out with sterile water with the cutting opening open in order to prevent the attachment from becoming blocked.

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Shaver blades/burrs (attachments)

Thoroughly rinse the instruments to remove all gross debris. Remove contaminants with brushes, sponges, soft cloths or a cotton applicator. Cleaning accessories are available from KARL STORZ (see catalog).

Cleaning brushes must be cleaned and disinfected immediately.

Remove the insert from the outer sheath of the shaver blade. Thoroughly clean the lumens of the reusable blades/burrs with a cleaning brush. The outer sheath as well as the insert of the shaver blade must then be rinsed out with a cleaning and disinfectant solution until no more residue can be seen. In order to avoid material damage, carry out the final flushing with microbiologically correct/sterile water. Finally dry with filtered compressed air.

Ultrasound cleaning (cleaning)

If required, the blade can be cleaned in an ultrasound bath in the disassembled state with an exposure time of approx. 3-5 min.

WARNING: Risk of injury to the patient: The sheaths and tips of the blades/burrs must not have any grooves, scratches or dents. After cleaning, inspect the handpieces, blades and burrs for cleanliness and damage.

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10.2.3 Machine cleaning and disinfection

Machine preparation

The handpiece (and, if applicable, accessories) as well as the shaver blades/burrs* can then be thermally disinfected at 93°C (device specific) in the washer. For this we recommend using the irrigation adaptor provided.

First photo: To achieve the best possible cleaning result when machine cleaning the shaver handpiece, use the irrigation adaptor 28205 SA. This irrigation adaptor is connected to the shaver handpiece in the same way as a shaver blade (see photo). Afterwards, an irrigation tube is connected to the irrigation adaptor. The lumen of the shaver can then be rinsed through.

Second photo: For the optimal machine cleaning of reusable outer blades, use the cleaning adaptor 28205SAA. This cleaning adaptor is equipped with a LUER lock, which can be connected to the irrigation tube in a washer and disinfector.

Third photo: For the machine cleaning of reusable inner blades, use the cleaning adaptor 28205SAI. Loosen the LUER lock slightly before inserting the proximal (= furthest from the patient) end of the inner blades into the adaptor. Subsequently tightening the LUER lock keeps the inner blade in the intended position. The cleaning adaptor is then connected to the irrigation tube in the washer and disinfector.

*Shaver blades/burrs: Care must be taken to ensure that the sheaths are thoroughly rinsed during machine cleaning and disinfection (special connection in the washer/disinfector).

* Shaver handpiece: To ensure the proper functioning of the handpiece, steam sterilization with fractionated prevacuum should be performed after the washer cycle (+134°C + 3°C), in particular to remove any moisture left in the area around the electric connections.

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10.2.4 Sterilization Care prior to sterilization

NOTE: Lubrication of the inner blade prior to sterilization increases the useful service life of the blade.

1. Apply a small amount of lubricant to the hub of the blade.
2. Insert the inner blade into the outer blade.
3. Turn the inner blade several times.
4. Dismantle the blades again.

NOTE: The new generation can be recognized by the fact that the driver is manufactured from black plastic, thereby rendering the use of O-rings and washers unnecessary.

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10.2.4.1 Safety instructions for sterilization

! WARNING: Successful sterilization is only possible on clean, disinfected surfaces.

! WARNING: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

! WARNING: Any deviations from the recommended parameters for sterilization should be validated by the user.

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WARNING: KARL STORZ recommends that, in order to assure sterility, all blades/burrs should be sterilized while disassembled and then reassembled in a sterile field.

WARNING: Before sterilization, the handpieces and reusable blades/burrs must be thoroughly cleaned and disinfected and all residues completely removed.

WARNING: The handpieces should not be treated with care agents such as instrument spray.

WARNING: Risk of burns! Parts are hot after steam sterilization. Allow them to cool down.

KARL STORZ has validated the following steam sterilization methods:

- Prevacuum
- Fractionated prevacuum

NOTE: Due to the design of the handpieces suitable for the system, we recommend the 'fractionated prevacuum procedure'.

CAUTION: Before putting into operation after completing sterilization, allow the shaver handpiece to cool down to room temperature under ambient conditions.

CAUTION: Intense build-up of heat at the tip of the shaver attachments: Interchanged attachments can lead to an intense build-up of heat. Closely observe, therefore, the printed article numbers during assembly!

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10.2.4.2 Possible sterilization methods

NOTE: Before sterilization, the items must be thoroughly cleaned and disinfected (see pages 40, 41).

Then pack the items into appropriate containers (Art. no. 39315 A).

WARNING: Danger to patients due to contaminated instruments: Subject the instruments to a complete cleaning cycle before the first and every further use. The instruments can also be sterilized with EtO gas. To achieve the desired sterility assurance level (SAL) of 10^{-6} , KARL STORZ recommends only the following gas and steam sterilization methods.

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Steam Sterilization

Prevacuum procedure
Fractionated prevacuum procedure

Gas sterilization

Gas sterilization with ethylene oxide (EtO)

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WARNING: Risk of burns!

Parts are hot after steam sterilization. Allow them to cool down.

10.2.4.3 Steam sterilization

Sterilization trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, and so that the air can escape and the condensate flow away. Load the sterilizer according to its instruction manual. After completion of the steam sterilization cycle, the sterilized instruments must be allowed to cool down slowly.

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Prevacuum method (steam sterilization)

Prevacuum or high vacuum sterilization consists of four phases: a conditioning phase, a sterilization phase, an evacuation phase and a drying phase. In the conditioning phase, a vacuum is generated in the sterilization chamber; then the instruments are heated by injecting steam. Sterilization occurs during the sterilization phase at a temperature of 134°C and a pressure (absolute) of 3 bar. The sterilization phase lasts 4 minutes. The evacuation phase removes the steam from the chamber. In the drying phase, the instruments are dried in a vacuum for approximately 20 minutes. The following sterilization conditions have been validated in a prevacuum sterilizer.

Temperature:	134°C (+ 3°C)
Pressure (p _a):	2 bar
Sterilization time:	4 minutes

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**Fractionated prevacuum method
(steam sterilization)**

The fractionated prevacuum procedure also consists of the four previously mentioned phases. In the conditioning phase, a vacuum is generated up to four times. Steam is then blown in and sterilization takes place at 134°C (+ 3°C) 5 minutes (minimum time). Before the final drying phase, the steam is removed from the chamber during the evacuation phase. Drying takes place under vacuum.

**10.2.4.4 Steam sterilization instructions
for the steam sterilization of the
multifunction handpiece**

WARNING: The multifunction handpiece requires additional exposure to steam for sterilization. The validated sterilization cycles are listed below.

WARNING: In order to guarantee sterility, KARL STORZ recommends removing all attachments from the multifunction handpiece, and sterilizing these separately.

**Prevacuum procedure
for multifunction handpiece:**

The following sterilization conditions have been validated for the multifunction handpiece and its accessories in a prevacuum sterilizer:

Temperature: 134°C (+3°C) (270 to 272°F)
Pressure (absolute): 3 bar
Exposure time: min. 4..., max. 18 minutes

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Ethylenoxid: 4 mg (< 24h)
Ethylenchlorhydrin: 9 mg (< 24h)

10.2.4.5 Ethylene oxide (EtO) gas sterilization

WARNING: KARL STORZ recommends that all blades/burrs be sterilized with ethylene oxide while disassembled and then reassembled in a sterile field.

WARNING: Sterility of KARL STORZ blades/burrs sterilized with ethylene oxide cannot be assured if the instruments are sterilized while assembled.

1. Clean the instrument as described in the preparation section of this instruction manual. Place the instrument securely on a sterilization tray.
2. All sterilization cycles must include a preconditioning cycle:
Temperature: $54 \pm 2^\circ\text{C}$ ($130 \pm 5^\circ\text{F}$) /
Time: 30 minutes, nominal
Humidity: $\geq 70\%$ RH
Vacuum: 1.3 psia
3. KARL STORZ has validated EtO sterilization using the following parameters:
Gas mixture: 100% EtO
Temperature: Set point of 55°C
Relative humidity: $\geq 70\%$ RH
Humidity dwell time: 30-45 minutes
Exposure time: 180 minutes, (full cycle)
EtO concentration: $735 \pm 30 \text{ mg/L}$
4. Aeration of the equipment can be accomplished in any aeration chamber at a temperature of $50 - 54^\circ\text{C}$ ($120 - 130^\circ\text{F}$) for 12 hours.

Maximum acceptable levels of residues following EtO sterilization are as follows (per ANSI/AAMI/ISO 10993-7):

Ethylene oxide: 4 mg (< 24h)
Ethylene chlorohydrin: 9 mg (< 24h)

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- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

Maintenance

Cleaning, disinfection and sterilization

i NOTE: Sterilization methods which deviate from the validated processes described must be validated individually. The country-specific and sterilization device-specific requirements and outline conditions concerning the effectiveness of the sterilization must be taken into account, as must, where applicable, the necessary airing times.

10.2.4.6 Gas sterilization with formaldehyde (FO)

When using formaldehyde the instruments can be used immediately after the sterilization cycle has been completed. An additional airing time is not necessary.

i NOTE: Observe national laws and regulations for the gas sterilization method. Sterilization using formaldehyde is a procedure not validated by KARL STORZ.

10.2.5 References for cleaning and sterilization

- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, No. 1, pp. 23 - 29, 1995.

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10.3 Maintenance and safety check

10.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

10.3.2 Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.

Test for proper functioning

Carry out the functional check as per the instruction manual (see page 28).

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Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/ repeat inspection, please refer to the latest issue of the relevant Service Manual.

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10.4 Service and repair

Repairing defective devices is only to be performed by personnel that have been authorized by KARL STORZ and using original KARL STORZ parts.

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10.5 Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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10.6 Repair program

Usually to bridge the repair period, you will receive a unit on loan, which you will then return to KARL STORZ as soon as you receive the repaired unit.

In Germany, you can refer the repairs directly to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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10.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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10.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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10.9 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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11 Technical description

11.1 General

In the event of the unit malfunctioning, the troubleshooting list provided in this chapter will assist the user in finding the possible cause. The troubleshooting list indicates how to find and rectify the various malfunctions.

WARNING: High voltage is present in the unit's interior. The unit must only be opened by an authorized KARL STORZ customer service technician.

WARNING: Maintenance work that is not mentioned in this chapter must only be performed by technical personnel that are authorized by KARL STORZ.

General information

Error messages are output directly to the screen on the UNIDRIVE® S III ARTHRO. Send the device in for servicing when a hardware error message is displayed (E01-E10). Take appropriate action for user error messages (E12-E30). E.g. E12 'NO MOTOR' - check whether the motor/handpiece is properly connected. Detailed information can be obtained from the error table on page 57 (see section 'Service mode' as well).

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11.2 Error table UNIDRIVE®

Error message	Cause	Remedy
1 = E01: MEMORY ERROR	Parameters could not be saved successfully.	Turn device on/off. Load default values. Or send unit in for repair.
2 = E02: MOTORBOARD DEFECTIVE	Motor control overloaded.	Turn device on/off. Or send unit in for repair.
3 = E03: HARDWARE CONFIGURATION ERROR	Hardware configuration is incorrect.	Send unit in for repair.
4 = E04: SOFTWARE CONFIGURATION ERROR	Software configuration is incorrect.	Send unit in for repair.
5 = E05: -		
6 = E06: -		
7 = E07: -		
8 = E08: -		
9 = E09: MOTOR DEFECTIVE M1	Motor on channel 1 is defective.	Release footswitch and restart. If the error persists, a different handpiece must be used.
10 = E10: MOTOR DEFECTIVE M2	Motor on channel 2 is defective.	See E09
12 = E12: NO MOTOR	No motor/handpiece connected.	Connect motor/handpiece.
13 = E13: TOO HIGH TEMPERATURE MOTOR M1	Motor/handpiece on channel 1 overheated.	Allow the connected handpiece to cool down or operate an additional handpiece on the 2nd channel.
14 = E14: -		
15 = E15: TOO HIGH TEMPERATURE MOTOR M2	Motor/handpiece on channel 2 overheated.	See E13
17 = E17: MOTOR NOT SUPPORTED M1	The motor/handpiece connected on channel 1 is not supported.	Only use handpieces approved by KARL STORZ. A service update may be necessary.
18 = E18: MOTOR NOT SUPPORTED M2	The motor/handpiece connected on channel 2 is not supported.	See E17
19 = E19: MOTOR STOPPED M1	Motor/handpiece on channel 1 blocked in operation.	Release footswitch and restart. If the error persists, clean the handpiece or exchange the attachment.
20 = E20: MOTOR STOPPED M2	Motor/handpiece on channel 2 blocked in operation.	See E19
21 = E21: NO SCB MODULE	Unit is not connected to an active SCB bus.	Insert SCB cable and activate SCB bus.

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Power switch lamp will not light up.

Possible cause

- 1 Power failure.
- 2 Defective power supply fuse.
- 3 Power switch lamp defective.

Recommended remedy

- 1 Check that there is electricity to the wall outlet.
- 2 Change fuses as described in the instruction manual. Make sure you fit the correct type of fuse.
- 3 Notify customer service.

Problem

- 1 Power switch lamp is illuminated.
- 2 No indication of rotation.
- 3 No energy output when operating the foot-switch/buttons.
- 4 Speed cannot be regulated.

Possible cause

- 1 Motor control electronics faulty.
- 2 Connection to motor interrupted.
- 3a Faulty handpiece.
- 3b Handpiece blocked or jammed.

Recommended remedy

- 1, 2, 4 Notify customer service.
- 3a Connect another handpiece.
- 3b Clean the handpiece.

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11.4 Technical Data

UNIDRIVE® S III ARTHRO	28723020-1
Line voltage	100–240 V~
Power frequency	50/60 Hz
Line fuse	2 x T 4.0 AH 250 V, IEC 127 Format
Power consumption	100–200 W
Operating temperature	10–40°C
Relative humidity	30%...75%
Dimensions	305 x 165 x 280 mm (W x H x D)
Weight	5.2 kg
Operating mode (in connection with the shaver multifunction handpiece)	5-10 min with cooling
Relative On-time	max. 25 min
Storage/transport conditions:	
Relative humidity (non-condensing)	20%...95%
Temperature	-10°C...60°C
Atmosph. pressure	500 hPa...1080 hPa

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(für 28723020-1)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock:
Protection class I
- Degree of protection against electric shocks:
Applied part type BF
- Type of moisture protection:
drip-water protection as per IPX 1

Conforms with IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 86-98).

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11.7 Directive compliance

Medical product of class IIa

This medical product bears the CE mark in
accordance with the Medical Device Directive
(MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

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11.8 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
11.9 Blockschaltbild

11.9 Block diagram

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**Spare parts, recommended
accessories**

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**12 Spare parts,
recommended accessories**

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Curved Aggressive Full Radius Resector,

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Curved full radius resector,

working length 180 mm, with keyhole, distal end angled at 15°, for partial, subtotal, and total synovectomy, for removal of degenerative meniscus tissue and adhesion filaments. The curvature enables inaccessible meniscus and cartilage areas, and also synovial areas to be reached.

Item	Cat. no.
<i>Sterile, single use, pack of 6</i>	
Diameter 4.5 mm (violet) color code light blue	28208 EGS

Curved aggressive full radius resector,

working length 180 mm, with keyhole, distal end angled at 15°, for partial, subtotal, and total synovectomy, for removal of degenerative meniscus tissue and adhesion filaments. The curvature enables inaccessible meniscus and cartilage areas, and also synovial areas to be reached.

Item	Cat. no.
<i>Sterile, single use, pack of 6</i>	
Diameter 4.5 mm (violet), color code light green	28208 EHS

Semi-hooded barrel burr,

working length 180 mm, for removing large areas, for removing osteophytes or a notchplasty; for smoothing bone irregularities, for milling the subchondral bone during subchondral abrasion for induction of fibrous cartilage, extensive cartilage ablation.

Item	Cat. no.
<i>Sterile, single use, pack of 6</i>	
Diameter 5.5 mm (dark gray), color code pink	28208 IDS

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Pro Line Shaver Blades,
working length 120 mm. Partial, subtotal and total synovectomy, removal of degenerative meniscus tissue, removal of cicatricial strands and adhesion filaments. Removal of osteophytes or a notchplasty; smoothing of bone irregularities, milling of the subchondral bone during subchondral abrasion for induction of fibrous cartilage, removal of large areas of cartilage, e.g., with subacromial decompression, for safe removal of meniscus tissue and adhesion filaments

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Agressive Pro Line Shaver Blades,
working length 120 mm. Partial, subtotal and total
synovectomy, removal of degenerative meniscus
tissue, removal of cicatricial strands and adhesion
filaments. Removal of osteophytes or a notchplas-
ty; smoothing of bone irregularities, milling of the
subchondral bone during subchondral abrasion
for induction of fibrous cartilage, removal of large
areas of cartilage, e.g., with subacromial decom-
pression, for the safe removal of meniscus tissue
and adhesion filaments

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Aggressive cutter,


working length 120 mm, for removal of degenerative meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.5 mm, color code green	28205 AK
<i>Sterile, single use, pack of 6</i>	
Diameter 4.5 mm (violet), color code green	28205 AKS
Diameter 4.5 mm (violet), color code green	*28208 AKS

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Full Radius Resector,

working length 120 mm. Partial, subtotal and total synovectomy, removal of degenerative meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.5 mm, color code yellow	28205 CK
<i>Sterile, single use, pack of 6</i>	
Diameter 4.5 mm (violet), color code yellow	28205 CKS
Diameter 4.5 mm (violet), color code yellow	*28208 CKS

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* 28208 Blades feature a keyhole on the distal tip
* For US orders, please add a /6 to the catalog number.

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Aggressive Full Radius Resector,
working length 120 mm. Partial, subtotal and total
synovectomy, removal of degenerative meniscus
tissue and adhesion filaments.

Item	Cat. no.
------	----------

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Full radius resector,
working length 120 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 3.5 mm, color code yellow	28205 CB
Diameter 4.2 mm, color code yellow	28205 CC
Diameter 5.5 mm, color code yellow	28205 CD
<i>Sterile, single use, pack of 6</i>	
Diameter 3.5 mm, color code yellow/light gray	28205 CBS
Diameter 4.2 mm, color code yellow/middle gray	28205 CCS
Diameter 5.5 mm, color code yellow/dark gray	28205 CDS

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Aggressive full radius resector,
working length 120 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 3.5 mm, color code blue	28205 DB
Diameter 4.2 mm, color code blue	28205 DC
Diameter 5.5 mm, color code blue	28205 DD

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Round burr,
working length 120 mm, for removal of osteo-
phytes, for smoothing bone irregularities, for milling
the subchondral bone during subchondral abra-
sion, for induction of fibrous cartilage

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.2 mm, color code orange	28205 FC
Diameter 5.5 mm, color code orange	28205 FD
Diameter 6.5 mm, color code orange	28205 FE
<i>Sterile, single use, pack of 6</i>	
Diameter 4.2 mm, color code orange/middle gray	28205 FCS
Diameter 5.5 mm, color code orange/dark gray	28205 FDS
Diameter 6.5 mm, color code orange/black	28205 FES

**Finish barrel burr
(10 flutes, 4 notches),**
working length 120 mm, for removing large areas,
e.g. during subacromial decompression at the shoul-
der; for removing osteophytes or a notchplasty; for
smoothing bone irregularities, for milling the sub-
chondral bone in subchondral abrasion, for induction
of fibrous cartilage; for removing large areas of carti-
lage, e.g. in subacromial decompression

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 4.2 mm, color code purple	28205 GD
Diameter 5.5 mm, color code purple	28205 GE

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28205 HE

**Aggressive barrel burr
(8 flutes, 4 notches),**

working length 120 mm, for removing large areas, e.g. during subacromial decompression at the shoulder; for removing osteophytes or a notch-plasty; for smoothing bone irregularities, for milling the subchondral bone in subchondral abrasion, for induction of fibrous cartilage; for removing large areas of cartilage, e.g. in subacromial decompression

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Slotted Whisker,

working length 120 mm, for safe removal of meniscus tissue and filaments of adhesion.

Item	Cat. no.
------	----------

Reusable, sterilizable

Diameter 4.2 mm, color code light orange	28205 LC
Diameter 5.5 mm, color code light orange	28205 LD

*Sterile, single use,
pack of 6*

Diameter 4.2 mm, color code light orange/middle gray	28205 LCS
Diameter 5.5 mm, color code light orange/dark gray	28205 LDS

Auger,

working length 120 mm, for removing large areas, e.g., during subacromial decompression at the shoulder; for removing osteophytes or a notch-plasty; for smoothing bone irregularities, for milling the subchondral bone in subchondral abrasion to induce fibrous cartilage; for removing large areas of cartilage, e.g. in subacromial decompression; for safe removal of meniscus tissue and filaments of adhesion.

Item	Cat. no.
------	----------

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**12.4 Recommended shaver blades
and burrs for small joint hand-
piece**

Aggressive cutter,
working length 70 mm, for removal of degenerative
meniscus tissue and filaments of adhesion.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code green	28206 AA
Diameter 3.5 mm, color code green	28206 AB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code green/white	28206 AAS
Diameter 3.5 mm, color code green/light gray	28206 ABS

Full radius resector,
working length 70 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code yellow	28206 CA
Diameter 3.5 mm, color code yellow	28206 CB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code yellow/white	28206 CAS
Diameter 3.5 mm, color code yellow/light gray	28206 CBS

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
Aggressive full radius resector,
working length 70 mm, for partial, subtotal, and
total synovectomy, for removal of degenerative
meniscus tissue and adhesion filaments.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code blue	28206 DA
Diameter 3.5 mm, color code blue	28206 DB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code blue/white	28206 DAS
Diameter 3.5 mm, color code blue/light gray	28206 DBS

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Slotted Whisker,
working length 70 mm, for safe removal of men-
iscus tissue and filaments of adhesion.

Item	Cat. no.
<i>Reusable, sterilizable</i>	
Diameter 2.5 mm, color code light orange	28206 LA
Diameter 3.5 mm, color code light orange	28206 LB
<i>Sterile, single use, pack of 6</i>	
Diameter 2.5 mm, color code light orange/white	28206 LAS
Diameter 3.5 mm, color code light orange/light gray	28206 LBS

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Small Joint Burr,
working length 70 mm, for removal of osteo-
phytes, smoothing of bone irregularities.

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Diameter 2.5 mm (white),
color code orange 28206 FAS

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Diameter 3.5 mm (light gray),
color code orange 28206 FBS

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**12.5 Accessories for shaver
handpieces and SL blades**

Irrigation adaptor, for validated cleaning of shaver
handpieces 28205 SA
Irrigation adaptor, for validated cleaning of
SL outer blades 28205 SAA
Irrigation adaptor, for validated cleaning of
SL inner blades 28205 SAI

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12.6 Accessories for pin driver

Drill wire, diameter 2.4 mm,
spiral-ground 28729 D
Kirschner wire, diameter 2.4 mm,
triangular-ground 28729 E

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12.7 Accessories for sagittal saw

Saw blade, blade width 10 mm, blade
thickness 0.4 mm, cutting depth 0.6 mm,
working length 50 mm 28207 SLA
Same, blade width 14 mm 28207 SLB

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Saw blade, blade width 10 mm, blade
thickness 0.4 mm, cutting depth 0.6 mm,
with stopper at 15 mm, 28207 SLC
working length 50 mm
Same, blade width 6 mm 28207 SLD
Same, blade width 8 mm 28207 SLE

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Saw blade, blade width 12 mm, blade
thickness 0.4 mm, cutting depth 0.6 mm,
working length 22 mm 28207 SLF

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Saw blade, blade width 6 mm, blade
thickness 0.4 mm, cutting depth 0.6 mm,
working length 31 mm 28207 SLG
Same, blade width 10 mm 28207 SLH

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13 Single-use blades and burrs

13.1 Intended use

The single-use sterile blades/burrs are intended for use with the UNIDRIVE® S III ARTHRO.

Use of the single-use sterile blades/burrs in any way other than that described is not permitted.

13.2 Use of the single-use products and explanation of the symbols employed

WARNING: Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the content of the pack must not be used under any circumstances.

WARNING: This product is intended for **single use** and must not be resterilized. Dispose of the product after use. Do not allow the blade to run 'dry'. The recommended speed of rotation for the blades is 4,000 rpm.

NOTE: This product is individually wrapped in protective packaging and sterilized with gamma radiation.

Manufacturer:

Manufacturer in accordance with Directive 93/42/EEC:

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Traceability

To ensure traceability, the self-adhesive identification labels enclosed in the original packaging can be kept and stuck in the operation report and in the patient dossier.

The person responsible for stores at the hospital can maintain a record of the products currently in store, including article no. (REF) and lot no. (LOT) so that rapid action can be taken in the event of changes or recalls.

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14 Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The UNIDRIVE® S III ARTHRO model 28723020-1 corresponds to IEC 60601-1-2:2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The UNIDRIVE® S III ARTHRO model 28723020-1 is a Group 1 unit (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

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NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**Electromagnetic
Compatibility (EMC)
Information**

-  **WARNING:** The UNIDRIVE® S III ARTHRO model 28723020-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.
-  **WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.
-  **WARNING:** The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the UNIDRIVE® S III ARTHRO model 28723020-1. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Table 201 Guidance and manufacturer's declaration – electromagnetic emissions		
UNIDRIVE® S III ARTHRO model 28723020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28723020-1 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group1	UNIDRIVE® S III ARTHRO model 28723020-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The UNIDRIVE® S III ARTHRO model 28723020-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class D (see NOTE below)	
Voltage fluctuations/ flicker emissions acc. to IEC 61000-3-3	not applied (see NOTE below)	
 NOTE: The measurement of harmonic emissions is not required (power < 75 W) The measurement of voltage fluctuations/flicker is not required (activated manually, unlikely to produce significant voltage fluctuations or flicker)		

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Table 202 Guidance and manufacturer's declaration – electromagnetic immunity			
UNIDRIVE® S III ARTHRO model 28723020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28723020-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	2 kV/1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	0 V, 10 msec 92 V, 100 msec 168 V, 500 msec 0 V, 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Table 204			
Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
UNIDRIVE® S III ARTHRO model 28723020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the UNIDRIVE® S III ARTHRO model 28723020-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF disturbance acc. to IEC 61000-4-6	3 Vrms 150 kHz to < 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the UNIDRIVE® S III ARTHRO model 28723020-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3,5/3] \cdot \sqrt{P}$ 150 kHz to < 80 MHz $d = [3,5/3] \cdot \sqrt{P}$ 80 MHz to < 800 MHz $d = [3,5/3] \cdot \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbance acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment
and the UNIDRIVE® S III ARTHRO model 28723020-1**

UNIDRIVE® S III ARTHRO model 28723020-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kse.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytten 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231
KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gounari 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

OOO KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 199939*
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Sary-Arka str., 6, BC "Arman", off.430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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/Karl Storz GmbH & Co. KG/