# STORZ

## KARL STORZ — ENDOSKOPE

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## 1 Important information for users of KARL STORZ equipment

Thank you for your confidence in the brand name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the UROMAT E.A.S.I.. All required details and all actions on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

## General warning

Users of the UROMAT E.A.S.I. must be thoroughly trained in its use for the applicable medical procedures.

The device and accessories should not be used until you have familiarized yourself fully with the instruction manual, assembly, and the care instructions. Instruction manuals should be carefully studied and be available to the surgical team during the procedure. It is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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


E.A.S.I. =  
Endoscopic Automatic System for Irrigation

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\*SCB = KARL STORZ Communication Bus

### 3 Controls, displays, connectors and their uses

- ① Power switch, illuminated  
0 = off
- ② Display and touch screen
- ③ Irrigation pump with locking lever
- ④ Suction pump with locking lever
- ⑤ Pump roller
- ⑥ Pressure sensor for the pressure measuring chamber of the irrigation tube
- ⑦ Tube holder for the suction tube
- ⑧ Footswitch connection socket
- ⑨ Connection for external systems
- ⑩ Connections for SCB\*
- ⑪ Potential equalization connector
- ⑫ Power cord socket
- ⑬ Fuse holder

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Potential equalization connector

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Applied part type BF in accordance with IEC 60601-1

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Electronic information product pollution control (China RoHS)

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This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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This material is intended for single use and must not be resterilized.

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|          |                                                                |            |
|----------|----------------------------------------------------------------|------------|
| <b>1</b> | <b>Important information for users of KARL STORZ equipment</b> | <b>III</b> |
| <b>2</b> | <b>Images of the equipment</b>                                 | <b>IV</b>  |
| <b>3</b> | <b>Controls, displays, connectors and their uses</b>           | <b>V</b>   |
| <b>4</b> | <b>Symbols employed</b>                                        | <b>VI</b>  |
| <b>5</b> | <b>General information</b>                                     | <b>4</b>   |
| 5.1      | The UROMAT E.A.S.I. device                                     | 4          |
| 5.2      | Intended use                                                   | 4          |
| 5.3      | User qualification                                             | 5          |
| 5.4      | Operating modes                                                | 5          |
| 5.5      | Operation                                                      | 5          |
| <b>6</b> | <b>Warnings and cautions</b>                                   | <b>6</b>   |
| 6.1      | Definitions                                                    | 6          |
| 6.2      | Training in the operation and function of the device           | 12         |
| 6.3      | Patient profile                                                | 13         |
| 6.4      | Intended region of treatment on the body                       | 13         |
| 6.5      | User profile (specialist physician)                            | 14         |
| 6.6      | User profile (assistant)                                       | 14         |
| 6.7      | Use                                                            | 15         |
| 6.8      | User position                                                  | 16         |
| 6.9      | Safety precautions at the site of installation                 | 16         |
| 6.10     | Safety precautions when operating the device                   | 17         |
| <b>7</b> | <b>Explanation of terms (in alphabetical order)</b>            | <b>18</b>  |
| <b>8</b> | <b>Installation and operating instructions</b>                 | <b>21</b>  |
| 8.1      | Unpacking the equipment                                        | 21         |
| 8.2      | Basic equipment                                                | 21         |
| 8.3      | Installing and connecting up the device (preparations)         | 21         |
| 8.4      | Operating the unit                                             | 24         |
| 8.5      | Preparing connection of the irrigation tube                    | 25         |
| 8.6      | Preparing connection of the suction tube                       | 28         |
| <b>9</b> | <b>Menu description</b>                                        | <b>31</b>  |
| 9.1      | Components of the header/footer                                | 31         |
| 9.1.1    | Components of the display                                      | 31         |
| 9.2      | Start screen                                                   | 32         |
| 9.3      | 'Operating modes' selection menu                               | 32         |
| 9.3.1    | Function and use                                               | 32         |
| 9.4      | Procedure List                                                 | 33         |
| 9.4.1    | Function/use                                                   | 33         |

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|           |                                                        |           |
|-----------|--------------------------------------------------------|-----------|
| 9.4.2     | Buttons                                                | 33        |
| 9.5       | Tube insertion/animation                               | 34        |
| 9.5.1     | Function/use                                           | 34        |
| 9.6       | Normal Operation Continuous Flow                       | 36        |
| 9.6.1     | Function/ use                                          | 36        |
| 9.6.2     | Buttons                                                | 37        |
| 9.6.3     | Displays                                               | 37        |
| 9.7       | Normal Operation Single Flow                           | 38        |
| 9.7.1     | Function/use                                           | 38        |
| 9.7.2     | Buttons                                                | 39        |
| 9.7.3     | Displays                                               | 39        |
| 9.8       | Normal Operation Irrigation/Suction                    | 40        |
| 9.8.1     | Function/use                                           | 40        |
| 9.8.2     | Buttons                                                | 41        |
| 9.8.3     | Displays                                               | 41        |
| 9.9       | Normal Operation Morcellation                          | 42        |
| 9.9.1     | Function/use                                           | 42        |
| 9.9.2     | Buttons                                                | 43        |
| 9.9.3     | Displays                                               | 43        |
| 9.10      | Normal Operation Lithotripsy                           | 44        |
| 9.10.1    | Function/use                                           | 44        |
| 9.10.2    | Buttons                                                | 44        |
| 9.10.3    | Displays                                               | 45        |
| 9.11      | Fluid supply monitoring of the<br>irrigation bag       | 46        |
| 9.11.1    | Function/use                                           | 46        |
| 9.11.2    | Buttons                                                | 47        |
| <b>10</b> | <b>Menu structure for 'Settings'</b>                   | <b>48</b> |
| 10.1      | Basic 'Settings' menu                                  | 48        |
| 10.1.1    | Function/use                                           | 48        |
| 10.1.2    | Buttons                                                | 48        |
| 10.2      | Pressure/flow settings in 'Continuous<br>Flow' mode    | 50        |
| 10.2.1    | Function/use                                           | 50        |
| 10.2.2    | Buttons                                                | 50        |
| 10.3      | Pressure/flow settings in 'Single Flow'<br>mode        | 51        |
| 10.3.1    | Function/use                                           | 51        |
| 10.3.2    | Buttons                                                | 51        |
| 10.4      | Pressure/flow settings in 'Irrigation/Suction'<br>mode | 52        |
| 10.4.1    | Function/use                                           | 52        |
| 10.4.2    | Buttons                                                | 52        |
| 10.5      | Pressure/flow settings in 'Morcellation'<br>mode       | 53        |

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## **5 General information**

### **5.1 The UROMAT E.A.S.I. device**

The UROMAT E.A.S.I. is a combined pressurized suction and irrigation pump with a continuous fluid flow for urological interventions. Alternatively, the two pumps can be operated independently of each other with flow regulation, e.g. for laparoscopic interventions ('Irrigation/Suction').

The device accumulates the irrigation liquid supplied and can inform the user when preset limit values have been exceeded (monitoring of irrigation fluid supply).

The UROMAT E.A.S.I. is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface. The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

### **5.2 Intended use**

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs and operating fields as well as the suctioning off of irrigation and bodily fluids, secretions and tissue during endoscopic interventions.

For safety reasons, the device must only be used as intended. Use of the UROMAT E.A.S.I. is contraindicated if, in the opinion of the responsible physician, the patient's health is endangered through its use.

#### **Contraindication**

Suction/irrigation pumps without a pressure monitoring system must not be used for the dilation of hollow organs and joints. Serious injuries to the patient may result if this use exclusion is not observed.

The UROMAT E.A.S.I. may only be operated with tubing sets, sheaths and accessories which have been designated as suitable for the device by

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### 5.3 User qualification

The UROMAT E.A.S.I. may only be used by physicians and trained medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the device.

### 5.4 Operating modes

The UROMAT E.A.S.I. has 5 operating modes for surgical and diagnostic interventions/operations. The operating modes are 'Continuous Flow', 'Single Flow', 'Irrigation/Suction', 'Morcellation' (operation with the UNIDRIVE® SIII), and 'Lithotripsy' (operation with the CALCUSON).

### 5.5 Operation

The device features a color monitor with a touch-sensitive screen for making entries (touch screen). The information displayed and the possible entries are adapted to the specific work flow status.

Both pump beds feature a rotating tensioning lever with which the user can open or close the pump bed for insertion and removal of the tubes.

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## 6 Warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this instruction manual, they should be carefully reviewed to ensure the safe and effective operation of this device. To make the signal words stand out more clearly, they are accompanied by a pictogram.

### 6.1 Definitions

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** Notes contain additional information on the safe operation of the device.

**WARNING:** Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel or yourself at risk.

**WARNING:** The electrical installations in the operating room in which the device is connected and operated must comply with the applicable IEC standards.

**WARNING:** Keep the unit out of reach of patients.

**WARNING:** For safety reasons, when the system is in use the output sockets of the device and the patient must not be touched at the same time.

**WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

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**WARNING:** Combinations of medical devices are only assured to be safe if:

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

**WARNING:** To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor

**WARNING:** When used in combination with HF surgical devices: Only HF operating modes with a maximum output voltage of 5.5 kVpp may be selected.

**WARNING:** Grounding reliability can only be achieved when the equipment is connected to 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**WARNING:** Test this device prior to each surgical procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

**WARNING:** Always unplug the device before performing any maintenance work on it.

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**WARNING:** Do not open the device! Danger of an electrical shock. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims.

**WARNING:** Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's wellbeing is maintained in the event of such interference.

**WARNING:** Risk of explosion if used in the presence of flammable anesthetics.

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**WARNING:** Reel pump (pressure-controlled pump motor). Do not touch the pump reels 5 when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

**WARNING:** The irrigation tubing set must be filled with irrigation fluid before use on a patient.

**WARNING:** The fluid level in the irrigation bag must be monitored carefully to ensure that no air is pumped into the cavum. Another bag of irrigant should be started prior to emptying the bag being used.

**WARNING:** The intracavitary pressure is always equal to or lower than the pressure at the outlet of the UROMAT E.A.S.I. which is used to regulate the device.

**WARNING:** The user should always use the lowest possible pressure for sufficient distention and/or irrigation.

**WARNING:** Do not reuse disposable tubing sets. Discard after use.

**WARNING:** If a footswitch no longer reacts when pressed, the device must not be operated.

**WARNING:** Before you activate the pump, make sure that the tubing sets have been inserted properly and that the pressure measuring chamber is correctly positioned.

**WARNING:** The time display showing the seconds in the header indicates that the device control is working properly. If the changing seconds are not visible, the device must not be operated.

**WARNING:** Before you start the intervention, make sure that the tubing sets used are not leaking or kinked and that the ends of the tubes are properly secured.

**WARNING:** The device regulates based on the pressure measured at the outlet of the irrigation pump. The pressure in the cavity depends, among other things, on the height difference between the device and cavum. Try to keep the difference as small as possible.

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**WARNING:** Any treatment may only be performed if visual observation of the action of the device is ensured.

**WARNING:** Attaching a pressurized measuring chamber (cambered diaphragm) can result in an erroneous pressure measurement. The irrigation tube should always be unpressurized when installed.

**WARNING:** Observe the applicable local regulations when disposing of disposable tubing.

**WARNING:** Observe the availability signal during power-up. If it can't be heard the device must not be used.

**WARNING:** Screw the patient tube's connector firmly into the pump tube's coupling to avoid leakages.

**WARNING:** Servicing may only be performed by specialists authorized by KARL STORZ.

**WARNING:** The UROMAT E.A.S.I. may be subject to sudden failures. Therefore you should ensure that a replacement unit is available or be prepared to apply alternative procedures if the intervention requires permanent fluid transportation/distension.

**WARNING:** The one-day tubing set for irrigation can remain in the reel pump for a maximum of one surgery day (24 hours) during which it can be used for several operations. Afterwards it must be disposed of.

**WARNING:** The user's position must always be selected so that any information message can be immediately identified. The distance from the device should not be more than 2 m.

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## Safety instructions

### Warnings and cautions

**WARNING:** Valves and tube clamps should be either completely opened or completely closed.

**WARNING:** Make sure that the set volume of the loudspeakers allows for perception of acoustic signals over normal background noise.

**WARNING:** If the touch screen no longer reacts when pressed, the device must not be operated.

**WARNING:** The device must not be operated if the touch screen or display are defective.

**WARNING:** It is not permitted to press more than one touch point simultaneously.

**WARNING:** The device may only be operated under the environmental conditions specified.

**WARNING:** Please ensure that there is sufficient air circulation for your device.

**WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

**WARNING:** The device may only be operated at the voltage stated on the rating plate.

**WARNING:** When cleaning the device, always observe the instructions in the instruction manual (see section »Cleaning and care«).

**WARNING:** When replacing fuses, only use fuses of the correct rating.

**CAUTION:** Avoid allowing liquids to enter the housing. Do not store liquids on or above the unit. If, despite precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies to the formation of condensation).

**CAUTION:** When cleaning the device surface, the area around the pressure sensor ⑥ should be wiped carefully due to its sensitive surface.

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**NOTE:** Observe the applicable local regulations when disposing of disposable tubing.

**NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

**NOTE:** This unit may only be operated with the Irrigation Tubing Set 031717-01 or 031767-01 intended for this purpose and which is indicated by the bar on the packaging.

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Before using the unit on the patient, it is imperative that you be acquainted with how the unit operates and is controlled.

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## **6.2 Training in the operation and function of the device**

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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#### 6.3.3 Urologie (Lithotripsie, Continuous Flow, Single Flow)

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Weight: irrelevant

Nationality: irrelevant

State of health: suitable for treatment in the physician's judgment

#### 6.3.2 Urology (Morcellation)

Gender: male

Age: irrelevant

Weight: irrelevant

Nationality: irrelevant

State of health: suitable for treatment in the physician's judgment

#### 6.3.3 Urology (Lithotripsy, Continuous Flow, Single Flow)

Gender: irrelevant

Age: irrelevant

Weight: irrelevant

Nationality: irrelevant

State of health: suitable for treatment in the physician's judgment

### 6.4 Intended region of treatment on the body

Laparoscopy: abdomen

Urology: lower and upper urinary tract

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**6.5 User profile (specialist physician)**

Recognized medical skills in the area of urology/ laparoscopy (specialist physician).

Adequate powers of comprehension to assess the situation rationally at all times during the application.

Adequate language skills in at least one of the languages used on the device and in the instruction manual.

Be thoroughly trained in the operation and use of the device.

Knowledge of the contents of the instruction manual.

No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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## 6.7 Use

The device may only be used in hospitals and medical practices in the ambient conditions given in the technical data.

- Not in intensive care units.
- Not in mobile facilities (cars or trucks, ambulances, aircraft, ships or trains).
- Position it on a vibration-free, fixed base.
- Visual distance from the physician to the device when operated: max. 2 m.
- Visual distance for operating the device: 30–70 cm.
- Viewing angle to the front of the device: max.  $\pm 45^\circ$ .
- The use of the device requires the division of work.
- Frequency of use: up to several times a day.
- Length of use: from a few minutes to several hours a day.
- Can only be moved by installing it on an equipment trolley.
- Can only be swiveled by installing it on ceiling supply units.
- Patient must not necessarily be under general anesthesia during the intervention (local anesthesia).
- May be used with other devices (light source, HF surgery etc.).
- To be operated in ambient light between 100 Lux and 1500 Lux.
- Ambient temperature in immediate proximity of the device during operation: max. 40°C.
- Sound pressure from background noise < 60 dBA.
- The device is placed outside the sterile area.

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### 6.8 User position

- For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$  at a distance of approx. 30–70 cm to the front panel.
- For observation of the actual values and for operation of the footswitch during the application, a visual distance to the device of 2 m is assumed (the tube length is 2.5 m).

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
### 6.9 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The device is equipped with a connector for potential equalization. It should be connected up in accordance with the applicable national regulations.

**WARNING:** The device must not be used in oxygenated environments. It is not intended for use in hazardous zones.

This means, among other things:

**WARNING:** When using easily combustible and explosive inhalation anesthetics or mixtures thereof, the device may not be operated in the demarcated hazard zone.

This also applies to easily combustible and explosive chemicals, e.g. skin disinfectants and surface disinfectants.

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#### **6.10 Safety precautions when operating the device**

It is the user's responsibility to make sure the device is safe and operates properly before use.

During treatment with the UROMAT E.A.S.I., the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring vital signs and the anesthetic.

Any treatment may only be performed if visual observation of the action of the device is ensured.

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## Explanation of terms in alphabetical order

### 7 Explanation of terms (in alphabetical order)

**i NOTE:** The terms are in alphabetical order in all languages. As such, it is possible that the German text on one page does not correspond to the English or Spanish text in terms of content.

**Activated:** When the line voltage is switched on, the user can move the rollers of the irrigation and suction pump via a button on the screen and, if necessary, via the footswitch of an external system. This is referred to as 'activation of the pumps'.

**i NOTE:** Further preconditions must be satisfied prior to activation. These are described on page 21 (Installation and operating instructions).

**Bolus function:** The 'Single Flow' operating mode → features the special 'Bolus' function, which allows the user to increase the pressure in percentage terms. For example, with a (Bolus) increase of 50%, the pressure increases from 60 to 90 mmHg.

The faster fluid exchange which can be achieved in this way can facilitate improved vision, for example.

**Day set / day tube set:** The day set concerns only the irrigation tube and not the suction tube.

Purpose of the day set: There is an interface in the irrigation tube set; this is kept sterile. Before each intervention a new patient tube (length approx. 2 m) is connected to this interface. After each operation, rather than replacing the entire irrigation tube only the patient tube of the irrigation tube set is replaced.

This means that the costly sterile parts of the irrigation tube set (e.g. pressure measuring chamber) only have to be replaced once a day; only the non-sterile (e.g. contaminated with blood) patient tube is replaced after every intervention.

**Deactivated:** Opposite of → activated.

**Menus:** A distinction is made between

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- **'Settings' menu:** this is used to set all basic settings (e.g. prior to the intervention), e.g. pressure/flow settings, screen language settings etc.

**Operating modes:** There are five different operating modes:

- **Continuous Flow**  
The irrigation pump firstly irrigates in order to build up pressure. The suction pump is inactive until the pressure has been achieved. Once the set pressure has been achieved, the Continuous Flow Mode is activated. Both pumps then irrigate and suck with the same flow rate (e.g., 50 ml/min). As such, there is no fall in pressure. The Continuous Flow Mode offers the advantage that clear fluid is continuously conveyed.
- **Single Flow**  
In this operating mode, the pump irrigates at first in order to build up pressure. Once this has been achieved, it stops. The irrigation pump only reacts to a drop in pressure at the sensor and then pumps accordingly in order to maintain the pressure.
- **Irrigation/Suction**  
The irrigation and suction flow rates can be altered independently of one another and the corresponding pumps activated or deactivated.
- **Morcellation**  
The irrigation pump firstly irrigates in order to build up pressure. Once this has been achieved, it stops. The suction pump is (de)activated via the footswitch of the UNIDRIVE® SIII and conveys the set fluid volume. The irrigation pump only reacts to a drop in pressure at the sensor and then pumps accordingly in order to maintain the pressure.
- **Lithotripsy**  
The suction pump is (de)activated via the footswitch of the CALCUSON and conveys the set fluid volume.

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**Explanation of terms**  
in alphabetical order

**Procedures:** With the UROMAT E.A.S.I., the user works exclusively with procedures. On delivery, the device features predefined procedures for each of the five operating modes.

The user cannot change, rename or delete these predefined procedures.

They can, however, copy a predefined procedure and save this under another name and then change the parameters individually (e.g. the desired pressure or flow).

The following applies: If a given procedure 'X1' in a given operating mode 'Z' is copied and used as the basis for this, only one additional procedure ('X2') can be saved in operating mode 'Z'.

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### 8.3 Installing and connecting up the device (preparations)


**WARNING:** Observe the section 'Safety precautions at the site of installation' on page 16.

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
Place the device on a flat surface.  
The difference in height between the device and the operating field should be as small as possible (see warning on page 8).

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### Potential equalization

Connect the potential equalization cable to the potential equalization connector **9**.

**WARNING:** The device's ground line should be installed by a qualified electrician.

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### Connecting the power cord

Connect the power cord when the device is switched off. Insert the power cord into the socket **10** as far as it will go.

The device may only be operated at the voltage stated on the rating plate.

**WARNING:** The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

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
**WARNING:** Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

**NOTE:** Ensure that access to the power plug is ensured at the site of installation.

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### Connecting the footswitch

Insert the connecting cable of the footswitch into the corresponding socket **7**.

**NOTE:** The footswitch is an optional accessory.

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
### Connecting the external system

Insert the connecting cable of the UNIDRIVE® SIII or the CALCUSON into the corresponding socket **6**.

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### Connecting up the SCB cable (optional)

In order to use this function, insert one end of the SCB cable (both ends are identical) into the connection socket ⑧ of the UROMAT E.A.S.I.. Connect the other end to an available SCB connection socket on a device in the existing SCB network or a device with which you wish to establish a network.

**① NOTE:** To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑧. Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control device or other SCB devices (see 'KARL STORZ SCB control NEO System' Instruction Manual).

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
### 8.4 Operating the unit


**CAUTION:** Before switching the UROMAT E.A.S.I. on, make sure that no tubing set has been inserted into the device and that the cables are properly connected.

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Press the power switch ①.

**CAUTION:** Do not touch the touch screen while you switch the device on.

Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo. After this, a white screen appears.

**WARNING:** Check the white screen for color errors. In the event of any errors the device must not be used.

During this period, the device runs a self test. Once the test has been successfully completed, an acoustic availability signal sounds. An animation appears on the screen which assists you in inserting the tubing set.

If the self test fails, the message 'Safe State' is displayed and an alarm sounds. Verify that the device is properly connected and that the preparations have been made according to the instructions. For further information about troubleshooting please refer to the section 'Technical description: Information signals' from page 69.

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**Animation**


The animation shows the order in which the tubing set should be inserted. The animation runs in an endless loop. This menu is exited when the device detects that:

- the irrigation and/or suction pump levers have been opened and closed and the pressure measuring chamber on the sensor ⑥ has been adapted.

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**8.5 Preparing connection of the irrigation tube**

**WARNING:** The packaging for the tubing set must not be damaged (sterility!). Check the tubing set for damage etc. before use.

Remove the irrigation tube (blue, with pressure measuring chamber) from the sterile packaging and connect it.

Pass the irrigation tube on for installation in the pump while observing aseptic procedures (safety precautions).

**WARNING:** The connection of the irrigation tubing set (sheath) to the patient must remain in the sterile area.

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### Unlocking the irrigation pump

**WARNING:** Do not touch the pump reels ⑤ when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.


Unlock the irrigation pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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
### Positioning the pressure measuring chamber

Position the pressure measuring chamber and turn it clockwise until fixed.

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### Inserting the irrigation tube

Insert the irrigation tube into the pump bed.


**WARNING:** Make sure that the tube is completely inserted into the pump bed.

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### Locking the irrigation pump

Lock the irrigation pump by turning the lever counterclockwise until it clicks into place.

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
Connect the irrigation tube to the irrigation bag(s). Make sure that no liquids are stored on or above the device (see the section 'Safety instructions', p. 10).

**Note:** The tubes for the irrigation bag are fitted with a blue clamp which can be used to pinch off the tubes if necessary.

**WARNING:** The irrigation pump is not self-priming. As such, the irrigation bag must always be located to the side and above the device.

### Day tube sets

If the device is approved for use with day tube sets, the steps for inserting the irrigation tube set 031767-01 must be performed as described on page 25-27. The separately packed patient tube 031162-01 is first connected with the sheath and then with the irrigation tube via the coupler.

**WARNING:** The patient tube must be disposed of after each intervention, while the tube inserted in the device may be reused.

**WARNING:** A new patient tube must be used for each patient.

**WARNING:** To guarantee the sterility of the connected tube set during patient changeover, the patient tubes are equipped with a separately packed, sterile plug. This must be screwed into the tube coupler from removal of the patient tube from the tube set until the next patient tube is connected.

Handling of the day tube set:

- Following the end of the intervention, the patient tube must be removed from the irrigation tube and disposed of correctly.

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## Installation and operating instructions

- The plug is then removed from the small sterile packaging and used to seal off the open end of the pump tube.
- The new patient tube is removed from the tube set for the next patient.
- The plug in the pump's irrigation tube is now removed and disposed of, and the new patient tube is connected to the open pump tube under aseptic conditions (safety measures).
- The new plug in its sterile packaging is then connected from the sterile area.

### 8.6 Preparing connection of the suction tube

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Remove the suction tube (red) from the sterile packaging and connect this to the instrument's sheath.

Pass the suction tube on for installation in the pump while observing aseptic procedures (safety precautions).

**WARNING:** The connection of the suction tubing set to the patient must remain in the sterile area.

**NOTE:** The tube to the sheath is fitted with a red clamp which can be used to pinch it off if necessary.

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#### Unlocking the irrigation pump


**! WARNING:** Do not touch the pump reels 5 when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

Unlock the suction pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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#### Inserting the suction tube


Insert the suction tube into the pump bed. Make sure that the tube is completely inserted into the pump bed.

The tube connector at the patient end should be inserted into the suction tube bed at the top in such a way that the locking lever can be turned counterclockwise to the end position without meeting obstacles. The outflow end of the suction tube is inserted in the tube bed at the bottom.

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#### Locking the suction pump

Lock the suction pump by turning the lever counterclockwise until it clicks into place.

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**Installation and  
operating instructions**

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**Connecting the waste container**

Connect the end of the suction tube which is not  
close to the patient (with the text 'Waste') to the  
waste container.

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- Header
- Footer
- Content area

#### Buttons

① Shortcut to the 'Settings' submenu. Here, amongst other things, the following can be set

- the desired screen language
- the desired volume of the acoustic signals
- the 'Fluid Supply Monitor'

② Switch acoustics signals ON/OFF

Desired volume: The desired signal volume can be set in the 'Settings' menu (see ①)

The state 'switched off' is retained until the user deactivates this function or switches the device off and on again.

③ Fluid Supply Monitor. When this button is pressed, the menu 'Fluid Supply Monitor' (monitoring of the amount of fluid in the irrigation bag) is called up, but only if this function has been activated in the menu 'Extended Settings' (see ①).

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## 9.2 Start screen


After switching on the supply voltage, a KARL STORZ logo with the device name is displayed for approx. 5 s (operating system loading time) followed by a white screen for at least 2 s.

This should allow the user to detect any pixel errors on the display.

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## 9.3 'Operating modes' selection menu

### 9.3.1 Function and use

After starting, the operation type to be performed can be selected.


The operating modes 'Continuous Flow', 'Single Flow', 'Irrigation/Suction', 'Morcellation' (operation with the UNIDRIVE® SIII), 'Lithotripsy' (operation with the CALCUSON) are available, and can be activated using the buttons ①, ②, ③, ④ or ⑤.

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## 9.4 Procedure List

### 9.4.1 Function/use

After selecting the operating mode, the Procedure List with the procedures saved for this kind of operation is displayed. On the first page, the last four procedures to have been used are shown, with the most recently used one first.

As of page 2, the other procedures are listed according to the procedure number.

The user can create new procedures; this is described in section 10.7.

The UROMAT E.A.S.I. contains numerous pre-defined procedures. For Continuous Flow, these begin at number 145, for Single Flow at number 245, for Irrigation/Suction at number 345, for Morcellation at number 445 and for Lithotripsy at number 545.

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
### 9.4.2 Buttons

- ① Jump to the previous page of the Procedure List
- ② Jump to the next page of the Procedure List
- ③ Shows the current page number/total number of pages
- ④ Jump back to 'Operation Selection'
- ⑤ Name of procedures

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## 9.5 Tube insertion/animation

### 9.5.1 Function/use

Animation to aid the assistant with insertion of the tube.

This menu level is exited as soon as it is detected that all of the components necessary for the relevant operation (irrigation function / suction function) have been connected.

Example 1: The Continuous Flow Mode is decisive for the selected procedure 'Resection Standard'. For this mode, both an irrigation and a suction tube must be inserted, each with a pressure dome attached. The menu level is only exited once this has been done.

Example 2: For a Single Flow procedure, only the irrigation tube has to be inserted in order to leave this menu level.

Within a fixed period of 24 hours, operation with a day set is permitted if this has been selected accordingly in the 'Extended Settings'.

This means that the irrigation tube (excluding the patient tube element) can remain in the pump.

After switching on the device, a state change in terms of the position detection of the pressure measuring chamber and irrigation pump is not necessary. The relevant animations cease.

The animation is comprised of the following images (see next page).

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## 9.6.2 Buttons

Explanation of terms: see page 18-20

- ③ By pressing the - button the desired pressure in mmHg is reduced.
- ④ By pressing the + button, the desired pressure in mmHg is increased.
- ⑤ By pressing the - button, the desired flow in ml/min is reduced.
- ⑥ By pressing the + button, the desired flow in ml/min is increased.

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Predefinitions: from number 145 for Continuous Flow, from number 245 for Single Flow, from number 345 for Irrigation/Suction, from number 445 for Morcellation, from number 545 for Lithotripsy.

- ⑦ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

## 9.6.3 Displays

- ⑧ Procedure names are comprised of the following 4 elements:

- Procedure icon
- Procedure number
- Designation 1
- Designation 2

- ⑨ Designation of the pressure/increase in pressure. During Normal Operation: Pressure.
- ⑩ Symbol for pressure
- ⑪ Pressure (three digits during Normal Operation)
- ⑫ Designation for flow during Normal Operation: Flow.
- ⑬ Symbol for flow
- ⑭ Flow (three or four digits during Normal Operation)

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#### 9.7.1 Function/use

**Explanation of terms:** see pages 18–20

**Use:** This working level is selected for URS procedures, for example.

**Setting the desired pressure:** Using the '+/-' buttons, the setpoint of the displayed cavity pressure can be changed. This can be done independently of the pump state (→activated/→ deactivated).

**Flow display ②:** The current actual flow is displayed here.

**Bolus ⑤:** percentage pressure increase for 2 seconds

**Footswitch ⑥:** the footswitch symbol is only displayed if a footswitch is connected. If this is the case, pressing the pedal has the same function as the 'Bolus' button.

**ⓘ NOTE:** The footswitch symbol is only displayed if a footswitch is actually connected and the selected procedure allows for the functionality of the footswitch.

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### 9.7.2 Buttons

Explanation of terms: see page 18-20

- ③ By pressing the – button the desired pressure in mmHg is reduced.
- ④ By pressing the + button, the desired pressure in mmHg is increased.

This predefined parameter can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer.

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- ⑤ The percentage pressure increase (10% to 50%) can be set when the pump is deactivated. The function of the percentage pressure increase is then carried out when the pump is activated.
- ⑦ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

### 9.7.3 Displays

- ⑥ Footswitch symbol, if a footswitch is connected
- ⑧ Procedure names are comprised of the following 4 elements:
  - Procedure icon
  - Procedure number
  - Designation 1
  - Designation 2
- ⑨ Designation of the pressure/increase in pressure. During Normal Operation: Pressure.
- ⑩ Symbol for pressure
- ⑪ Pressure (three digits during Normal Operation)
- ⑫ Designation for flow during Normal Operation: Flow.
- ⑬ Symbol for flow
- ⑭ Flow (four digits during Normal Operation)

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### 9.8 Normal Operation Irrigation/Suction

#### 9.8.1 Function/use

**Explanation of terms:** see pages 18–20

**Use:** This working level is selected for LAP procedures, for example.

**Setting the desired flow ① / ②:** Using the «+/-» buttons, the setpoints of the flow parameters for irrigation\* ① and suction ② can be changed. This can be done independently of the pump state (activated/deactivated).

Each individual pump can be activated and deactivated independently with its own Start/Stop button.

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

**\*NOTE:** In the ENDO Mode, only the setpoint of the flow parameter for suction ② can be changed.

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### 9.8.2 Buttons

Explanation of terms: see page 18-20

- ③ By pressing the – button, the desired irrigation flow rate in ml/min is reduced.
- ④ By pressing the + button, the desired pressure in mmHg is increased.
- ⑤ By pressing this button, the pump is activated. When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.
- ⑥ By pressing the – button, the desired suction flow rate in ml/min is reduced.
- ⑦ By pressing the + button, the desired suction flow rate in ml/min is increased.

These predefined parameters can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer.

Predefinitions: from number 145 for Continuous Flow, from number 245 for Single Flow, from number 345 for Irrigation/Suction, from number 445 for Morcellation, from number 545 for Lithotripsy

- ⑧ By pressing this button, the pump is activated. When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

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

### 9.8.3 Displays

- ⑨ Procedure names are comprised of the following 4 elements:
  - Procedure icon
  - Procedure number
  - Designation 1
  - Designation 2
- ⑩ Text 'Irrigation Pump Flow'
- ⑪ Symbol for (irrigation) flow (is always displayed)
- ⑫ Flow (four digits, in [ml/min])
- ⑬ Text 'Suction Pump Flow'
- ⑭ Symbol for (suction) flow (is always displayed)
- ⑮ Flow (four digits, in [ml/min])

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## 9.9 Normal Operation Morcellation

① **NOTE:** This mode of operating a KARL STORZ pump together with a KARL STORZ morcellator is not authorized for users inside the U.S.

### 9.9.1 Function/use

**Explanation of terms:** see pages 18 – 20

**Use:** This working level is selected for operation with the UNIDRIVE® SIII for morcellation.

**Setting the desired pressure ①:** Using the '+/-' buttons, the setpoint of the displayed cavity pressure can be changed. This can be done independently of the pump state (→activated/→ deactivated).

**Setting the desired (suction) flow ②:** Using the '+/-' buttons, the setpoint of the suction flow rate can be changed. This can be done independently of the pump state (activated/deactivated).

Footswitch ⑦: is displayed if the connecting cable between the UNIDRIVE® S III and the UROMAT E.A.S.I. has been inserted and the UNIDRIVE® S III is switched on. --> If there is no UNIDRIVE® S III connected, a message appears on the display and the pump cannot be started.

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

① **NOTE:** The footswitch symbol ⑦ is only displayed if a footswitch is actually connected.

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
### 9.9.2 Buttons

Explanation of terms: see page 18 – 20

- ③ By pressing the – button the desired pressure in mmHg is reduced.
- ④ By pressing the + button, the desired pressure in mmHg is increased.
- ⑤ By pressing the - button, the desired (suction) flow in ml/min is reduced.
- ⑥ By pressing the + button, the desired (suction) flow in ml/min is increased.

These predefined parameters can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer.  
Predefinitions: from number 145 for Continuous Flow, from number 245 for Single Flow, from number 345 for Irrigation/Suction, from number 445 for Morcellation, from number 545 for Lithotripsy

- ⑥ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

### 9.9.3 Displays

- ⑦ Footswitch symbol, if a footswitch is connected
- ⑨ Procedure names are comprised of the following 4 elements:
  - Procedure icon
  - Procedure number
  - Designation 1
  - Designation 2
- ⑩ Designation of the pressure/increase in pressure. During Normal Operation: Pressure.
- ⑪ Symbol for pressure
- ⑫ Pressure (three digits during Normal Operation)
- ⑬ Designation for (suction) flow during Normal Operation: Suction Pump Flow.
- ⑭ Symbol for (suction) flow
- ⑮ (Suction) flow (four digits during Normal Operation)

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
### 9.10.1 Function/use

**Explanation of terms:** see pages 18 – 20

**Use:** This working level is selected for operation with the CALCUSON for lithotripsy.

**Setting the desired (suction) flow ①:** Using the «+/-» buttons, the setpoint of the suction flow rate can be changed. This can be done independently of the pump state (activated/deactivated).

**NOTE:** This operating mode is activated via the signal from the CALCUSON footswitch if the CALCUSON and the UROMAT E.A.S.I. are connected to one another.

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### 9.11 Fluid supply monitoring of the irrigation bag

#### 9.11.1 Function/use

**Accessing the menu:** This menu is accessed by pressing the button 'Fluid Supply Monitor' (see the orange highlighted irrigation bag at the bottom of the screen).

This menu is only displayed if the 'Fluid Supply Monitor' is activated (enabled) in the menu 'Extended Settings'.

**Use:** This menu monitors the level of fluid in the irrigation bag. This menu is available regardless of whether the pumps are activated or deactivated. For safety reasons, when the pump is activated the display returns to the previous menu if no button is pressed within 3 seconds. Values which were modified will then be discarded.

The selected irrigation bag size and remaining available fluid volume are displayed. The value for the irrigation bag size can be changed with the +/- buttons (buttons ① and ②). Moreover, the remaining volume value can be modified, if required (buttons ③ and ④), for example if the procedure is started with a partly-used irrigation bag. Using the ⑤ button, the available fluid supply can also be reset to the size of the bag, for example when a new full irrigation bag is connected. All functions can be used regardless of the status of the irrigation and suction pumps (activated/deactivated).

- ⑥ Discard the values and exit the menu
- ⑦ Confirm and save the values and exit the menu

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
The modified parameters are retained even if the pumps are deactivated. If the power is switched off and on, however:

- the irrigation bag size is reset to the original value which was set in the menu 'Extended Settings'
- the remaining fluid volume is set to the same value as the irrigation bag size (it is assumed that a new irrigation bag is used for each new intervention)

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### 9.11.2 Buttons

**NOTE:** The values for the irrigation bag size are overwritten by any modifications made in the section 'Irrigation bag size' in the 'Extended Settings' menu.

**NOTE:** Please do not forget to save your settings ⑦

- ①/② Decrease/increase irrigation bag size setpoint. Possible setpoints: 1 – 3 – 5 – 10.  
③/④ Decrease/increase the actual value for available fluid supply.

Range: 'Irrigation bag size' up to 0.0

- ⑤ Reset the available fluid supply to the value of the selected bag size (see input field above) – for example when changing the irrigation bag.  
⑥ Return to the previous menu level, changes are discarded.  
⑦ Return to the previous menu level, changes are saved.

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### 10.1 Basic 'Settings' menu

#### Accessing the menu

This menu is accessed by pressing the 'Settings' button ⑧.

**Explanation of terms:** see pages 18-20

#### 10.1.1 Function/use

This menu is the basic menu (main menu) in the menu structure 'Settings'. It is the starting point for the submenus in which settings for the UROMAT E.A.S.I. can be made and displayed.

#### 10.1.2 Buttons

① Shortcut to the menu 'Pressure/Flow Settings'. This is where you enter the initial values for the 'operating menu'.

② Shortcut to the menu 'Extended Settings'. In this menu, you can:

- activate/deactivate the Fluid Supply Monitor
- set the standard value for the irrigation bag size
- activate the use of → Day sets

③ Shortcut to the 'Device Settings' menu. In this menu, you can:

- set the language
- set the display brightness
- set the volume for the information and alarm signals
- set the volume for button actuation
- set the date format

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④ Shortcut to the menu 'Device Configuration'. In this menu, the following information is displayed:

- serial number
- operating hours
- user interface software version
- software process control

⑤ Shortcut to the menu 'System Log'. This is where the information and alarm signals are recorded.

⑥ Shortcut to the menu 'Service Programs'. Access only with password for authorized service personnel.

⑦ Shortcut to the menu 'Procedure Administration'. In this menu, you can:

- save predefined procedures under another name and set different values for the newly created procedure (e.g. different pressure)
- rename procedures (excluding the procedures predefined on delivery)
- delete procedures (excluding the procedures predefined on delivery)
- overwrite procedures with other values (excluding the procedures predefined on delivery)

⑧ Return to the previous menu level.

⑨ Return to the previous menu level. Saving the settings.

**NOTE:** When you switch to the menu level 'Operation Selection' the buttons ① and ⑦ are not available.

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### 10.2 Pressure/flow settings in 'Continuous Flow' mode

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- Specification of setpoint for irrigation pressure ①. Example: In this menu, a pressure value of 75 mmHg is specified. This setpoint is then saved and displayed in the operating menu 'Continuous Flow'.
- Specification of setpoint for irrigation flow ②. In this menu, a flow rate of 800 ml/min is specified. This setpoint is then saved and displayed in the operating menu 'Continuous Flow'.

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#### 10.2.2 Buttons

- ⑤ Decrease the specified setpoint for the pressure.
- Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) for the pressure/flow which is to be displayed in the 'Pressure'/'Flow' field of the operating menu 'Continuous Flow' is specified.
- ⑥ Increase the specified setpoint for the pressure
- ⑦ Decrease the specified setpoint for the flow
- ⑧ Increase the specified setpoint for the flow

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
## 10.3 Pressure/flow settings in 'Single Flow' mode

### 10.3.1 Function/use

Input menu for the pressure parameters:

Depending on the operating mode selected, different parameters can be selected individually.

- Specification of setpoint for irrigation pressure ①.  
Example: In this menu, a pressure value of 140 mmHg is specified. This setpoint is then saved and displayed in the operating menu 'Single Flow'.
- Specification of setpoint for (Bolus) pressure increase ②.  
Example: In this menu, a Bolus value (percentage pressure increase for 2 seconds) of 20% is specified. This setpoint is then saved and displayed in the operating menu 'Single Flow'.

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#### Overwriting the values of the operating menu:

Any previously changed settings in this operating menu are overwritten upon confirmation of the settings with the button ④.

### 10.3.2 Buttons

Decrease the specified setpoint for the pressure.

Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) for the pressure/Bolus which is to be displayed in the 'Pressure'/'Bolus' field of the operating menu 'Single Flow' is specified.

- ⑥ Increase the specified setpoint for the pressure
- ⑦ Decrease the specified setpoint for the 'Bolus'
- ⑧ Increase the specified setpoint for the 'Bolus'
- ③ Return to the previous menu; the changed values are discarded
- ④ Return to the previous menu; the changed values are saved

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## 10.4 Pressure/flow settings in 'Irrigation/Suction' mode

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- Specification of setpoint for irrigation flow ① (not in the 'ENDO Mode').  
Example: In this menu, a flow rate of 800 ml is specified. This setpoint is then saved and displayed in the operating menu 'Irrigation/Suction'.
- Specification of setpoint for suction flow ②.  
In this menu, a flow rate of 800 ml/min is specified. This setpoint is then saved and displayed in the operating menu 'Irrigation/Suction'.

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### 10.4.2 Buttons

- ⑤ Decrease the specified setpoint for the irrigation flow rate.

Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) for the flow which is to be displayed in the 'Flow' field of the operating menu 'Irrigation/Suction' is specified.

- ⑥ Increase the specified setpoint for the irrigation flow rate

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## 10.5 Pressure/flow settings in 'Morcellation' mode

### 10.5.1 Function/use

Input menu for the pressure and flow parameters:  
Depending on the operating mode selected, different parameters can be selected individually.

- Specification of setpoint for suction pressure ①.  
Example: In this menu, a pressure value of 75 mmHg is specified. This setpoint is then saved and displayed in the operating menu 'Morcellation'.
- Specification of setpoint for suction flow ②.  
In this menu, a flow rate of 1300 ml/min is specified. This setpoint is then saved and displayed in the operating menu 'Morcellation'.

**Saving changes permanently:** Changes which are made here are only retained for as long as no new procedure is selected. If the settings should be retained for longer, the values must be saved in a new procedure (see section 10.7).

**Overwriting the values of the operating menu:** Any previously changed settings in this operating menu are overwritten upon confirmation of the settings with the button ④.

### 10.5.2 Buttons

- ⑤ Decrease the specified setpoint for the pressure.  
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) for the pressure/flow which is to be displayed in the 'Pressure'/'Flow' field of the operating menu 'Morcellation' is specified.
- ⑥ Increase the specified setpoint for the pressure
- ⑦ Decrease the specified setpoint for the flow
- ⑧ Increase the specified setpoint for the flow
- ③ Return to the previous menu; the changed values are discarded
- ④ Return to the previous menu; the changed values are saved

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## 10.6 Pressure/flow settings in 'Lithotripsy' mode

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- Specification of setpoint for suction flow ①.  
In this menu, a flow rate of 800 ml/min is specified. This setpoint is then saved and displayed in the operating menu 'Lithotripsy'.

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### 10.6.2 Buttons

- ⑤ Decrease the specified setpoint for the flow.  
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) for the flow which is to be displayed in the 'Flow' field of the operating menu 'Lithotripsy' is specified.
- ⑥ Increase the specified setpoint for the flow.
- ③ Return to the previous menu level; the changed values are discarded.
- ④ Return to the previous menu level; the changed values are saved.

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

## 10.7 'Extended Settings' menu

### 10.7.1 Function/use

This is an entry menu for further settings.

**NOTE:** Do not forget to confirm your entries by pressing button ⑦.

### 10.7.2 Buttons (orange)

- ① Here you can call up the menu which serves to monitor the irrigation bag fluid supply.

- ②/③ Decrease/increase the initial value of the irrigation bag size.

The irrigation bag size is shown in the 'Fluid Supply Monitor' menu if this function has been selected with the button ①.

**NOTE:** If the initial values are confirmed by pressing the button ⑥, the values set in the 'Fluid Supply Monitor' menu are overwritten!

- ④ Enable the 'Day Set' function.  
⑤ Return to the previous menu level and discard the selected values.  
⑥ Return to the previous menu level and save the selected values.

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### 10.8.1 Function/use

This is the menu for setting

- ①/② the language displayed on-screen,
- ⑦/⑧ monitor brightness,
- ③/④ the volume of information/alarm signals,

**WARNING:** Make sure that the set volume of the loudspeakers allows for perception of the acoustic signals despite normal background noise.

- ⑨/⑩ volume of button actuation,
- ⑤/⑥ date format.

In addition, pressing the ⑪ button allows you to switch to the settings menu for date/time. Any modified values that have not been saved are then discarded.

⑫ Return to the basic 'Settings' menu. The modified values are discarded.

⑬ Return to the basic 'Settings' menu. The modified values are saved.

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## 10.10 'Procedure Administration' menu

### 10.10.1 Function/use

① By pressing this button, you can overwrite an existing procedure that you defined yourself.

Procedures which were defined on delivery cannot be overwritten. This means that the button does not react when you click it for a predefined procedure.

② By pressing this button, you can delete an existing procedure that you defined yourself.

Procedures which were defined on delivery cannot be deleted. This means that the button does not react when you click it for a predefined procedure.

③ Save as New Procedure. If you have changed the values of a predefined procedure and want to have these changed values available permanently, we recommend creating a new procedure. To do so, press this button.

④ Rename Procedure. Self-defined procedures can be renamed.

### 10.10.2 Save as New Procedure

#### Function/use

The current settings are saved as a new procedure together with a designation, name and number. The designation and name are restricted to 20 alpha-numeric characters or a total width of 200 pixels.

#### Buttons

a Tab – switches between the input windows (display elements a, b and c)

b Return to the menu level 'Settings'; the entries are discarded

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### **Example: Defining a new procedure and setting other values yourself**

Our example has the following objectives:

- A procedure with the sample name '156 URS Standard' is to be copied
- The following value should be changed:  
- Pressure: from 75 mmHg to 65 mmHg
- The new procedure should be saved under the name '111 Cont. Flow 65mmHg' under the number '111'.

#### Method for saving a new procedure

- Switch the device on and wait a few seconds until the menu level 'Type of operation' appears.
- Select an operation type
- On the touch screen, click on the procedure '156 URS Standard'.
- Click on the tool symbol in the bottom left of the screen
- (The 'Settings' menu appears)
- Click on the 'Procedure Administration' submenu
- Click on 'Save as New Procedure'
- In the orange highlighted field, enter 'Cont. Flow 65mmHg'
- Press the TAB key (the large button on the far bottom left). The second input field is then highlighted orange
- Leave this field empty
- Press the TAB key again (the large button on the far bottom left). The third input field is then highlighted orange
- Enter the identification number for the new procedure here, i.e., '111'
- Click on the check mark in the far bottom right to save your settings
- (The updated procedure list appears. The new procedure 'Cont. Flow 65mmHg' appears at the top)

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### Method for changing values

- Click on the tool symbol in the bottom left of the screen
- (The 'Settings' menu appears)
- Click on 'Pressure/Flow Settings'
- Under the header 'Pressure', click twice on the minus sign to reduce the value from 75 to 65
- Click on the check mark in the bottom right to save your settings
- (You are once again in the 'Settings' menu)
- (The procedure 'Cont. Flow 65mmHg' has now been permanently changed).
- (The updated procedure list appears. The new procedure '111 Cont. Flow 65mmHg' appears at the top).

### Method for renaming the procedure

- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Rename Procedure'
- The flashing character (cursor) is now on the left, next to the '1' of '160'
- Click on the 'arrow-to-the-right button' in the bottom, right-hand area until the cursor is to the right of '65'
- Enter a space. The procedure is now called '111 Cont. Flow 65 mmHg'
- Save your entries by clicking on the check mark

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
#### Method for deleting the procedure

- (You are now in the Procedure List)
- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Delete Procedure'
- Click on the check mark on the bottom right to delete the procedure permanently

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#### **10.11 'Device Configuration' menu**

##### **Function/use**

The following elements are displayed:

- Device serial number
- Operating hours
- User interface software version
- Process control software version

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
### 10.12 'System Log' menu

#### Function/use

Presentation of the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

The most recent entry is allocated the number 0.

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
### 10.13 'Service Programs' menu

#### Function/use

The service programs are only available to service personnel. Access is password-protected.

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
### 10.14 'Safe State' menu

- The Safe State is set for certain causes of risks. In this state, the device can no longer be operated normally.
- However, in the 'Settings' area you can access another examination.
- This menu level is a special case. Switching off the line voltage ends the Safe State.

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


## 11 Service and repair

### 11.1 Fuse replacement

Turn the device off and remove the power cord from the electrical outlet.

**WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Loosen the line fuse holder with a screwdriver or other suitable tool. Remove the line fuse holder.

**WARNING:** Only use fuses of the correct rating. Replace defective fuses and screw the line fuse holder back on.

Line fuse: 2 x T1.6 AH 250 V

Reconnect the power cord.

Check that the device can be switched on and that the KARL STORZ start screen is displayed.

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### 11.2 Cleaning and disinfection

**! WARNING:** Always pull out power plug before cleaning!

**! CAUTION:** Avoid allowing liquids to enter the housing at all costs.

The design allows only for mechanical wipe-down disinfection. The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

**! CAUTION:** When cleaning the device surface, the area around the pressure sensor (6) should be wiped carefully due to its sensitive surface.

**! WARNING:** Do not reuse disposable tubing sets. Discard after use.

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### 11.3 Maintenance

Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the reliability of the device and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

**WARNING:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, we recommend a functional or safety test of the device at least once a year.

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**FOR DETAILED INSTRUCTIONS PLEASE REFER TO THE LATEST VERSION OF THE SERVICE MANUAL.**

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#### 11.4 Service and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts only.

#### 11.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

**WARNING:** At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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### 11.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

When shipping the device, use either the original packaging or equivalent packaging marked with the symbols illustrated here.

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### 11.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

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### 11.9 Warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot bear the cost of transportation or the risks associated with shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen, Germany

**NOTE:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str.  
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall

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## 12.1 Information signals

## 12.1.1 Configuration

All the information signals created, with the exception of those which display the cause for an inoperable device, (entry 'Safe State' in the 'Reaction' column of the table from page 77) are not self-retaining. This means that the signal is only issued whilst the signal conditions are present (unless something different is defined under 'Reaction').

If the signal conditions are only very brief, there are regulations as to the minimum visual/acoustic signal duration to which the described information signals also have to comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by an alarm signal with a higher priority or a low priority information signal is suppressed so long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized first will be signaled. Should signals deviate from this rule, this is noted in the corresponding entry from page 77.

## 12.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, a minimum display duration of 5s is specified.

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### 12.1.3 Acoustic signals


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. An entire sequence is played through as a minimum.

The acoustic signal can be switched on or off using a button in the footer of the device display.

The acoustic signal, however, can only be switched off temporarily by the user. After switching on the line voltage, the setting is always 'Audio Alarm On'.

#### Information signals (general)

The pitch is modulated by +/- 2 Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics. In the case of messages which are related to the Safe State, the sequence is repeated every 15 s, but otherwise just played through once.

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


### Availability signal


After completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by +/- 2 Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics.


### Bolus function

The following periodic signal is issued on starting the Bolus function.

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## Technical description

### Specification for acoustic information signals

#### Low irrigation fluid supply

If the Fluid Supply Monitor is activated and the irrigation fluid supply is low, a double gong is produced. It lasts for 1.5 seconds.

#### Key click

Pressing a key on the touch screen is confirmed by a short beep. The volume of this beep can be set separately from the volume of all other information signals in the 'Device Settings' menu and can also be permanently deactivated.

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| Brief description (1 line with approx. 40 characters possible)<br>Possible cause                 | Priority | Reaction                                                                                                                     |
|--------------------------------------------------------------------------------------------------|----------|------------------------------------------------------------------------------------------------------------------------------|
| 20: Pump footswitch recognized                                                                   | 3        | None                                                                                                                         |
| Pump footswitch was plugged in                                                                   |          |                                                                                                                              |
| 21: Pump footswitch disconnected                                                                 | 3        | End footswitch special function                                                                                              |
| Pump footswitch was disconnected                                                                 |          |                                                                                                                              |
| 22: Shaver system recognized                                                                     | 3        | None                                                                                                                         |
| The shaver system is connected and switched on                                                   |          |                                                                                                                              |
| 23: Shaver system disconnected                                                                   | 3        | End shaver special function                                                                                                  |
| The shaver system was switched off or disconnected                                               |          |                                                                                                                              |
| 24: Footswitch/Shaver not supported                                                              | 3        | None                                                                                                                         |
| The current operating context has not assigned a function to the actuated foot- or shaver switch |          |                                                                                                                              |
| 26: Low irrigation fluid supply                                                                  | 1        | None. Message remains until the irrigation fluid supply is above the limit value or the fluid supply monitor is deactivated. |
| The irrigation fluid supply is coming to an end (< 200 ml)                                       |          |                                                                                                                              |
| 27: Pressure dome detached                                                                       | 2        | Cancel activation (if activated)                                                                                             |
| Pressure dome was detached                                                                       |          |                                                                                                                              |
| 28: Irrigation pump unlocked                                                                     | 2        | Cancel activation (if activated)                                                                                             |
| Irrigation pump locking lever was opened                                                         |          |                                                                                                                              |
| 29: Suction pump unlocked                                                                        | 2        | Cancel activation<br>(if activated) Single Flow: no reaction                                                                 |
| Suction pump locking lever was opened                                                            |          |                                                                                                                              |
| 2A: Illegal password                                                                             | 3        | None                                                                                                                         |
| Incorrect password entered                                                                       |          |                                                                                                                              |
| 2B: Password entry mismatch                                                                      | 3        | None                                                                                                                         |
| On entering/confirming a new password, different entries were made                               |          |                                                                                                                              |
| 2C: Re-attach pressure dome                                                                      | 2        | Menu change not executed, message remains until the pressure dome is removed                                                 |
| On turning on the line voltage, the pressure dome was attached but the day set was not enabled   |          |                                                                                                                              |
| 2D: Invalid procedure number                                                                     | 3        | None                                                                                                                         |
| The procedure number which has been selected for saving has already been used                    |          |                                                                                                                              |

**Information messages**

| Brief description (1 line with approx. 40 characters possible)<br>Possible cause                                                                                                          | Priority | Reaction                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------------------------------------------------------------------------------------|
| 40: Contact KARL STORZ Service<br>Voltage of the SRAM back-up battery on powering up < 2.7 V                                                                                              | 2        | None                                                                                  |
| 41: Contact KARL STORZ Service<br>Voltage of the SRAM back-up battery on powering up < 2.3 V                                                                                              | 1        | None                                                                                  |
| 42: Reset of settings/data<br>The settings/data were reset manually or due to an error (one or more saving areas)                                                                         | 2        | None                                                                                  |
| 43: Duration of Bolus mode exceeded<br>The current duration of the Bolus Mode exceeds the setpoint by more than 5 s                                                                       | 1        | Cancel activation                                                                     |
| 44: Perform maintenance<br>On switching on the line voltage incompatible software versions were found by the User Interface and Process Control                                           | 1        | Safe State                                                                            |
| 45: Activated ext. switch during Power-up<br>On switching on the line voltage a pedal on the footswitch or a button on the shaver handpiece was actuated                                  | 1        | Switch function is not executed, message remains until the switch element is released |
| 46: Inserted external switch is activated<br>On insertion a pedal on the footswitch or a button on the shaver handpiece was actuated                                                      | 1        | Pedal function is not executed, message remains until the switch element is released  |
| 47: Touch screen short circuit<br>On switching on the line voltage the touch screen was pressed                                                                                           | 1        | Safe State                                                                            |
| 48: Touch screen short circuit<br>Maximum touch screen actuation time of 40 s exceeded                                                                                                    | 1        | Safe State                                                                            |
| 4A: Contact KARL STORZ Service<br>On switching on the line voltage a faulty component in the drive system was detected                                                                    | 1        | Safe State                                                                            |
| 4B: Contact KARL STORZ Service<br>On switching on the line voltage a faulty component in the drive system was detected                                                                    | 1        | Safe State                                                                            |
| 4C: Bolus footswitch depressed<br>The Bolus footswitch remains depressed even though the duration of the Bolus Mode has expired. Notification 10 s after the expiration of the Bolus Mode | 1        | None                                                                                  |

| Brief description (1 line with approx. 40 characters possible)<br>Possible cause                                                                    | Priority | Reaction                                                                                                                               |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|----------|----------------------------------------------------------------------------------------------------------------------------------------|
| 50: Reattach pressure dome<br>The value of the exposed sensor is outside the tolerance of 24.4 mV – 427 mV at the sensor output (approx. 4–70 mmHg) | 1        | Cancel activation. Pressure dome animation<br>Renewed activation is only possible once the pressure dome has been removed and replaced |
| 53: Reattach pressure dome<br>A pressure is measured by the sensor which could result in the latter's destruction                                   | 1        | Cancel activation. Pressure dome animation<br>Renewed activation is only possible once the pressure dome has been removed and replaced |
| 60: Internal temperature increased<br>Temperature inside the device is greater than 55°C                                                            | 3        | None, message remains until the temperature falls below the limit value or is superseded by message 61                                 |
| 61: Critical internal temperature<br>Temperature inside the device is outside the 10°C – 70°C range                                                 | 2        | Cancel activation, message remains until the temperature falls below the limit value                                                   |
| 62: Irrigation motor temperature too high<br>Temperature > 70°C                                                                                     | 3        | None, message remains until the temperature falls below the limit value or is superseded by message 63                                 |
| 63: Irrigation motor temperature critical<br>Temperature > 100°C or disconnection of a sensor line                                                  | 2        | Cancel activation, message remains until the temperature falls below the limit value                                                   |
| 64: Suction motor temperature too high<br>Temperature > 70°C                                                                                        | 3        | None, message remains until the temperature falls below the limit value or is superseded by message 65                                 |
| 65: Suction motor temperature critical<br>Temperature > 100°C or disconnection of a sensor line                                                     | 2        | Cancel activation, message remains until the temperature falls below the limit value                                                   |

**Information messages**

| Brief description (1 line with approx. 40 characters possible)<br>Possible cause                                                                   | Priority | Reaction                                                                                                                                                                                                                                                            |
|----------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 70: Reattach pressure dome in <4h<br>More than 20 hours have elapsed since the last sensor test                                                    | 2        | None, message remains until the pressure dome is detached                                                                                                                                                                                                           |
| 71: Reattach pressure dome<br>The offset value of the pressure sensor is no longer up-to-date (pressure dome must be detached and reattached)      | 1        | <b>Activated:</b><br>Activation remains<br>Message remains<br><b>Not activated</b> (or immediately after deactivation):<br>Pressure dome animation<br>Activation and deletion of the message is only possible once the pressure dome has been removed and replaced. |
| 72: Restart device in <4h<br>Since the last self-test over 20 hours ago, the device has been in continuous operation                               | 2        | None, message remains until the device is switched off and back on again                                                                                                                                                                                            |
| 73: Restart device<br>Since the last self-test over 24 hours ago, the device has been in continuous operation and a new self-test is now necessary | 1        | <b>Activated:</b><br>Activation remains<br>Message remains<br><b>Not activated</b> (or immediately after deactivation):<br>Safe State                                                                                                                               |
| 80: Contact KARL STORZ Service<br>Supply voltage(s) outside the specification                                                                      | 1        | Safe State                                                                                                                                                                                                                                                          |
| 81: Contact KARL STORZ Service<br>Supply voltage below 15 V                                                                                        | 1        | Safe State                                                                                                                                                                                                                                                          |
| 90: Unable to download a procedure<br>An attempt was made to download a procedure from the RUI whilst the device was not ready-to-receive.         | 2        | None                                                                                                                                                                                                                                                                |
| 91: Invalid procedure download data<br>The data sent from the RUI are faulty                                                                       | 2        | None                                                                                                                                                                                                                                                                |
| FD: Internal system error I<br>An unexpected system state is detected                                                                              | 3        | None                                                                                                                                                                                                                                                                |
| FE: Internal system error II<br>An unexpected system state is detected                                                                             | 2        | Cancel activation                                                                                                                                                                                                                                                   |
| FF: Internal system error III                                                                                                                      | 1        | Safe State                                                                                                                                                                                                                                                          |

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- Type of protection against electric shock:  
Protection Class I
- Degree of protection against electric shock:  
Applied part type BF
- Type of moisture protection:  
Drip-proof as per IPX1

#### 12.5 Directive compliance (for UP410)

This medical product belongs to Class IIb.

This medical device bears the CE mark in accor-  
dance with the Medical Device Directive (MDD)  
93/42/EEC.

**NOTE:** The code number after the CE mark  
indicates the responsible notified body.

#### 12.6 Technical documentation

On request, the manufacturer will provide available  
circuit diagrams, detailed parts lists, descriptions,  
instructions for adjustments, and other documenta-  
tion deemed useful for suitably qualified user per-  
sonnel authorized by the manufacturer when repair-  
ing device components that have been designated  
as reparable by the manufacturer.

Supply of such technical documentation relating to  
the device shall not be construed as constituting  
manufacturer's authorization of user's personnel,  
regardless of their level of technical training, to  
open or repair the device.

Explicitly exempted herefrom are those main-  
tenance and repair operations described in this  
manual.

**NOTE:** We reserve the right to make engi-  
neering modifications in the interest of pro-  
moting technological progress and generating  
performance improvements without obligation  
on our part to submit prior notice thereof.

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## Electromagnetic compatibility (EMC) information

### 14 Electromagnetic compatibility (EMC) information

**WARNING:** Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The UROMAT E.A.S.I. model UP410 corresponds to EN/IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The UROMAT E.A.S.I. model UP410 is a Group 1 device (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the devices to different electrical circuits.

If you have any further questions, please contact

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 UP410 should not be used adjacent to or stacked with other equipment. If it must be operated next to or on top of other devices, you should observe the device or the system to verify its operation according to its designated use in the configuration in which it is used.

 **WARNING:** The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

 **WARNING:** The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the electronic UROMAT E.A.S.I. model UP410. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:

| Cable type   | Shielded | Length [m] | Ferrite | Used for               |
|--------------|----------|------------|---------|------------------------|
| MKF 4 E-4396 | Yes      | 4          | No      | Footswitch             |
| Power Cord   | No       | 2          | No      | Power Supply           |
| SCB Cord     | Yes      | 1.0 m      | No      | Connecting cable       |
| PE           | No       | > 3.0 m    | No      | Potential equalization |

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| Table 1<br>Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                                  |            |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>UROMAT E.A.S.I. model UP410 is intended for use in the electromagnetic environment specified below.<br/>The customer or user of the UROMAT E.A.S.I. model UP410 should make sure that it is used in such an environment.</p> |            |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Emissions test                                                                                                                                                                                                                  | Compliance | Electromagnetic environment – guidelines                                                                                                                                                                                                                                                                                                                                                                                                                   |
| RF emissions CISPR 11                                                                                                                                                                                                           | Group 1    | <p>The UROMAT E.A.S.I. model UP410 uses only RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.</p> <p>The UROMAT E.A.S.I. model UP410 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p> |
| RF emissions CISPR 11                                                                                                                                                                                                           | Class B    |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Harmonic emissions IEC 61000-3-2                                                                                                                                                                                                | Class A    |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Voltage fluctuations/flicker emissions IEC 61000-3-3                                                                                                                                                                            | Complies   |                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

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**Table 2**

**Guidance and manufacturer's declaration – electromagnetic immunity**

UROMAT E.A.S.I. model UP410 is intended for use in the electromagnetic environment specified below.  
The customer or user of the UROMAT E.A.S.I. model UP410 should ensure that it is used in such an environment.

| Immunity test                                                                                                  | EN/IEC 60601 test level                                                                                                                                                                                         | Compliance level                                                                                                                                                                                                | Electromagnetic environment – guidelines                                                                                                                                                                                                                                                         |
|----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic Discharge (ESD)<br>acc. to IEC 61000-4-2                                                         | ± 6 kV contact<br>± 8 kV air                                                                                                                                                                                    | ± 6 kV contact<br>± 8 kV air                                                                                                                                                                                    | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                                    |
| Electrical fast transient/burst<br>acc. to IEC 61000-4-4                                                       | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                                                  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                                                  | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                              |
| Surge<br>acc. to IEC 61000-4-5                                                                                 | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth                                                                                                                                                            | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth                                                                                                                                                            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                              |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>acc. to IEC 61000-4-11 | <5% $U_T$ *<br>(>95% dip in $U_T$ )<br>for 0.5 cycle<br>40% $U_T$<br>(60% dip in $U_T$ )<br>for 5 cycles<br>70% $U_T$<br>(30% dip in $U_T$ )<br>for 25 cycles<br><5% $U_T$<br>(>95% dip in $U_T$ )<br>for 5 sec | <5% $U_T$ *<br>(>95% dip in $U_T$ )<br>for 0.5 cycle<br>40% $U_T$<br>(60% dip in $U_T$ )<br>for 5 cycles<br>70% $U_T$<br>(30% dip in $U_T$ )<br>for 25 cycles<br><5% $U_T$<br>(>95% dip in $U_T$ )<br>for 5 sec | Mains power quality should be that of a typical commercial or hospital environment.<br><br>If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply. |
| Magnetic field at the supply frequency (50/60 Hz)<br>acc. to IEC 61000-4-8                                     | 3 A/m                                                                                                                                                                                                           | 3 A/m                                                                                                                                                                                                           | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                        |

\* Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

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**Table 4**  
**Guidance and manufacturer's declaration – electromagnetic immunity –**  
**for equipment and systems that are not life-supporting**

UROMAT E.A.S.I. model UP410 is intended for use in the electromagnetic environment specified below.  
The customer or user of the UROMAT E.A.S.I. model UP410 should ensure that it is used in such an environment.

| Immunity test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | EN/IEC 60601 test level                 | Compliance level   | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                         |                    | Portable and mobile RF communications equipment should be used no closer to any part of the UROMAT E.A.S.I. model UP410, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br>Recommended separation distance:                                                                                                                                                                                                                                                          |
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 3 V <sub>rms</sub><br>150 kHz to 80 MHz | 3 V <sub>rms</sub> | $d = [3.5/3]\sqrt{P}$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Radiated RF<br>IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 3 V/m<br>80 MHz to 2.5 GHz              | 3 V/m              | $d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz<br>$d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                         |                    | <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p>Note 1: At 80MHz and 800MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.</p> <p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> |                                         |                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

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**Table 6**

**Recommended separation distances between portable and mobile HF communications equipment and the UROMAT E.A.S.I. model UP410**

The UROMAT E.A.S.I. model UP410 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum<br>output power of<br>transmitter<br>[W] | Separation distance <b>d</b> [m]<br>according to frequency of transmitter |                                                                  |                                                                 |
|--------------------------------------------------------|---------------------------------------------------------------------------|------------------------------------------------------------------|-----------------------------------------------------------------|
|                                                        | 150 kHz to 80 MHz<br>$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$          | 80 MHz to 800 MHz<br>$d = \left[ \frac{3.5}{3} \right] \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = \left[ \frac{7}{3} \right] \sqrt{P}$ |
| 0.01                                                   | 0.12                                                                      | 0.12                                                             | 0.23                                                            |
| 0.1                                                    | 0.38                                                                      | 0.38                                                             | 0.73                                                            |
| 1                                                      | 1.2                                                                       | 1.2                                                              | 2.3                                                             |
| 10                                                     | 3.8                                                                       | 3.8                                                              | 7.3                                                             |
| 100                                                    | 12                                                                        | 12                                                               | 23                                                              |

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80MHz and 800MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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# WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:  
Company stamp/signature:

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To be filled out by device owner:  
Return address/company stamp:

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga  
Ontario, L5N 3R3, Canada  
Phone: +1 905 816-4500, Fax: +1 905 858-4599  
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100,  
Toll free: 800 421-0837 (USA only)  
Fax: +1 424 218-8526  
Toll free: 800 321-1304 (USA only)  
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.  
175 Cremona Drive  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@kseas.com

KARL STORZ Endoscopia  
Miramar Trade Center  
Edificio Jerusalem, Oficina 108,  
La Habana, Cuba  
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.  
Lago Constanza No 326, Col. Chapultepec Morales  
D.F.C.P. 11520, Mexico, México  
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso  
B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
P.O. Box 153, Rolf Olsensvei 28  
N-2007 Kjeller, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14,  
12724 Skärholmens, Sweden  
Postal address: Po Box 8013,  
14108 Kungens Kurva, Sweden  
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Valimotie 1 B, 4. kerros, 00380 Helsinki, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG  
Representation Office  
Kestucio g. 59  
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
392 Edinburgh Avenue, Slough  
Berkshire, SL1 4UF, Great Britain  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscope Nederland B. V.  
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.  
Phone: +32 473 810 451, E-Mail: info@stople.be

KARL STORZ Endoscopy France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.  
Zadarska 80  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Verovškova c. 60A  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si  
KARL STORZ Endoscopia Romania srl  
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4  
041393 Bukarest, Romania  
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801  
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.  
Gounari 40,  
15343 Agia Paraskevi/Athens, Greece  
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392  
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.\*  
Ipsilantou Str. 32  
54248 Thessaloniki, Greece  
Phone: +30 2310 304868 Fax: +30 2310 304862

\*Repair & Service Subsidiary

KARL STORZ Industrial\*\*  
Gedik Is Merkezi B Blok  
Kat 5, D 38-39, Bagdat Cad. No: 162  
Maltepe Istanbul, Turkey  
Phone: +90 216 442 9500, Fax: +90 216 442 9030

\*\*Sales for Industrial Endoscopy

OOO KARL STORZ Endoscopy – WOSTOK  
Derbenyevskaya nab. 7, building 4  
115114 Moscow, Russia  
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41  
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine  
Obolonska naberezhna, 15  
building 3, office 3  
04210 Kiev, Ukraine  
Phone: +380 44 42668-14, -15, -19  
Fax: +380 44 42668-20

KARL STORZ ENDOSKOPE  
East Mediterranean and Gulf S.A.L.  
Block M, 3rd Floor  
Beirut Souks, Weygand Street  
2012 3301 Beirut, Lebanon  
Phone: +961 1 1999390, Fax: +961 1 1999391  
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.  
P.O. 6061, Roggebaai 8012, South Africa  
Phone: +27 21 417 2600, Fax: +27 21 421 5103  
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan  
Sary-Arka str., 6, BC "Arman", off. 430  
010000 Astana, Republic of Kazakhstan  
Phone: +7 7172 552549, 552788  
Fax: +7 7172 552444  
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf  
Building 25 (Al-Hayyan) – Units 205 – 206  
Dubai Health Care City, P.O.Box 118069  
Dubai - United Arab Emirates  
Phone: +971 (0)4 4471230  
Fax: +971 (0)4 4471231  
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.  
D-181, Okhla Industrial Area, Phase-1  
New Delhi 110020, India  
Phone: +91 11 26815445-51,  
Fax: +91 11 268129 86  
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG  
Resident Representative Office  
80/33 (44/19) Dang Van Ngu,  
F.10-Q. Phu Nhuan  
Ho Chi Minh City, Vietnam  
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.  
Hong Kong Representative Office  
Unit 1601, Chinachem Exchange Square  
1 Hoi Wan Street, Quarry Bay,  
Hong Kong, People's Republic of China  
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch Company,  
Room 610, China Life Tower No. 6, Chaowai Street,  
Beijing, 100020, People's Republic of China  
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Unit 3901-3904, Tower 1 Grand Gateway  
No. 1 Hong Qiao Road  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch Company  
F-5, 24/F., Chuanxing Mansion,  
No. 18 Renming Road South  
Chengdu, Sichuan, 610016, People's Rep. of China  
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch Company,  
Rm 2225, Tower B, City Plaza,  
No. 83, Zhongshan Road, Heping District,  
Shenyang, Liaoning, 110001,  
People's Republic of China  
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch Company  
Room 1119-20, Dongshan Plaza  
69 Xianlie Road Middle, Dongshan District,  
Guangzhou, Guangdong, 510095,  
People's Rep. of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.  
Jeilong Tower 1, Suite 905  
16 Neuti-ro, Bundang-Gu  
Seongnam, 463-847, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
6F-1, No. 10, Sec. 1, Beisin Rd.  
Sindian, Taipei County 231, Taiwan (R.O.C.)  
Phone: +886 933 014 160  
Fax: +886 2 8672 6399  
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**

**KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen  
Postfach 230  
78503 Tuttlingen  
Germany

Tel.: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)