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Apostille

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Table of Contents

About This Manual	4
Safety Considerations	5
Intended Use and Contraindications	6
Declaration of conformity	7
Dimensions	10
Monitor housing	10
Mounting	11
Cable entry point	11
Connection field	12
Monitor Mounting to Stand	12
Video and Graphics Inputs/Outputs	13
Quick Start	14
Standard Timings RGB, DVI	15
Output Gamma	15
Monitor classification	15
Aspect Ratio	15
User Interface	15
Control Interface	16
OSM (On-Screen-Manager)	17
Unlocking the on-screen display (OSD)	20
Overview-Menu-Structure	21
Rotation and Mirroring	26
Single Source	26
Dual Display	27
Dual Source	28
Remote Control	29
Boot Time & Splash screen	29
Power Save Mode	29
BIOS (Update)	29
Languages	29
Cleaning and care of the HD monitor	30
Maintenance	31
Servicing and repair	31
Repair program	32
Limitation of liability	33
Manufacturer's warranty	33
Technical Specifications	34
General circuit diagram	35
Electromagnetic Compatibility (EMC) Information	36
Warranty card	52
Subsidiaries	54

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About This Manual

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This manual is designed to assist the user with proper installation, setup and operation of the 9619NB -19" and 9626NB (-2) -26" LCD flat screen monitor. Depending on the model and options that were purchased, some of the features and options in this manual may not apply to the monitor you are using.

The functional descriptions in this manual apply to:

Part Number: 9619NB KARL STORZ HD Monitor 19", 9626NB (-2) KARL STORZ HD Monitor 26" and 9627NB(-2) KARL STORZ HD Monitor 27"

Firmware BIOS: 1.0 and later.

All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Keep this manual available for reference in a convenient and conspicuous location near the device.

Proprietary Notice and Liability Disclaimer

The information disclosed in this document, including all designs and related materials, is the valuable property of KARL STORZ and / or its licensors and, as appropriate, they reserve all patent, copyright, and other proprietary rights to this document, including all design, manufacturing reproduction, use, and sales rights thereto, except to the extent said rights are expressly granted to others.

All product names and trademarks or registered trademarks are the property of their respective holders.

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Safety Considerations

CAUTION

This symbol alerts the user that important literature concerning the operation of this unit has been included. Therefore, it should be read carefully in order to avoid potential problems.

This symbol warns user that un-insulated voltage within the unit may have sufficient magnitude to cause electrical shock. Therefore, it is dangerous to make contact with any part inside the unit. To reduce the risk of electric shock, DO NOT remove cover (or back). There are no user serviceable parts inside. Refer servicing to qualified service personnel.

To prevent fire or shock hazards, do not expose this device to rain or moisture. Also, do not use this device's polarized plug with an extension cord receptacle or other receptacles unless the prongs can be fully inserted. The display is designed to meet the medical safety requirements for devices in the vicinity of patients according to EN60601-1 Ed.3 2007.

Underwriters Laboratories (UL) Classification:

UL Safety Compliance:

This LCD monitor is U.L. Classified as medical - GENERAL MEDICAL EQUIPMENT as to ELECTRIC SHOCK, FIRE and MECHANICAL HAZARDS only in accordance with ANSI/AAMI ES60601-1 (2005), CAN/CSA-C22.2 No. 60601-1 (2008).

EEC Safety Compliance:

This display meets the requirements of EN-60601-1 so as to conform to the Medical Device Directive 93/42/EEC (general safety information).

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**9619 NB / 9626 NB(-2) /
9627 NB(-2) – HD Monitor**

This monitor complies to the above standards only when used with the supplied medical grade power supply.

24 Volt DC: Order no. 9626 PS
(Typ: DT150-24 MOPP)

The monitor should be powered from a center tapped circuit when used in the US at voltages over 120 volts. The monitor is intended for continuous operation.

This display is energized from an external electrical power source for class 1 equipment. It is the responsibility of the installer to test the display's grounding and to verify that it complies with the hospital, local and national impedance requirements.

A ground post, located on the back of the display, may be used for the purpose of potential equalization and grounding of the display's housing. Any such ground must be installed in accordance with applicable electrical regulations. The ground post is shown on the drawing found on page 11 (cable entry point).

Intended Use and Contraindications

Intended Use

This monitor is intended for use in a medical environment to display high-quality video and graphic images, and complies with the safety requirements for devices used in the vicinity of patients.

Contraindications

Sterility: The device can only be sterile when used with a plastic sheet containing a cut-out for the active area of the display and the ring slider. No air vent in which dirt can collect. The smooth vitreous surface is significantly bigger than the active display area, meaning it is easily possible to mask the glass without restricting the field of vision. The cable cover conceals the connectors, which are impossible to disinfect, and protects the sheet from damage

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**9619 NB / 9626 NB(-2) /
 9627 NB(-2) – HD Monitor**

The monitor may not be used in the presence of flammable anesthetics mixture with air, oxygen or nitrous oxide.

For mission critical applications, we strongly recommend that a replacement unit be immediately available.

Recycling

The lamps in this equipment could contain a very small amount of Mercury. Follow local governing ordinances and recycling plans regarding the recycling or disposal of this equipment, or contact the Electronic Industries Alliance (www.eiae.org).

Declaration of conformity
FCC and Council Directives of European Standards

This device complies with Part 15 of FCC rules and 93/42/EEC of the Council Directives of European Standards. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesirable results.

1. Use the attached specified cables with the color monitor so as not to interfere with radio and television reception. Use of other cable and adapters may cause interference with other electronic equipment.

2. This equipment has been tested and found to comply with the limits pursuant to FCC Part 15 and CISPR 11 (Comité International Spécial des Perturbations Radioélectriques). This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

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IEC

This device has been tested and found to comply with the limits for medical devices pursuant to EN 60601-1-2 :2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity.

FCC, Council Directives of European Standards and IEC

There is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult your dealer or an experienced radio/TV technician for help.

Accessory equipment connected to this monitor must be certified according to the respective IEC standards (i.e., IEC 60950-1 for data processing equipment and EN 60601-1 Ed. 3 2007 for medical equipment). Any person who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements according to EN 60601-1 Ed. 3 MOPP. Any person who is responsible for mounting the monitor to a system must ensure that the mounting equipment used with this monitor complies to the standard EN 60601-1. If in doubt, consult the medical technical services department or your local representative.

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


Dimensions



Monitor housing

- Special housing made of magnesium alloy/glass
(9626NB(-2) and 9627NB(-2)) or
aluminum/glass (9619NB)
- No fans

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


Mounting

The monitors are equipped with a VESA 100 mounting compatible to most boom arms or cart mounts. Please check whether the designated extension or ceiling arm or the equipment trolley is able to take the weight of the monitor.

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Cable entry point

Hidden cable entry, through VESA 100 mounting or the lower side of the cable cover.

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



Monitor Mounting to Stand

Attach the stand to the monitor using four supplied
torx screws and special bit.

The stand supports rotation ($\pm 60^\circ$) and tilt ($\pm 25^\circ$)
adjustment.

Note: The monitor stand is an optional accessory.
Order no. 9626 SF

Caution: Risk of tipping! Ensure that the
monitor is stable. Ensure that the stand
with the monitor is screwed correctly to
the equipment trolley, ceiling arm or the
suspended supply unit.

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Quick Start

Switching on the Device

Remove the cable cover at the back of the monitor and connect the power supply to the display using the mains cable. Insert the power plug. Connect the video source to the KARL STORZ monitor. Switch on the peripheral device and then the display. The KARL STORZ logo is displayed, followed by the active video signal.

The electronics, designed by KARL STORZ, incorporate 'state of the art' technology which, at initialization, examines the incoming signal and automatically displays the video image in its proper format. This eliminates the need to make adjustments to most video sources. To fine tune the image, please refer to how to use the OSM on page 17.

ON/OFF-Function

The monitor can be turned off via the OSD (see page 17) or by holding down the middle of the control wheel for more than 5 seconds.

To power up the monitor again from this state, the middle of the control wheel must be held down for more than 2 seconds. All of the usual functions can now be used again.

The maximum power consumption when switched off is:
monitor 0.2 W and power supply 0.5 W.

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**9619 NB / 9626 NB(-2) /
9627 NB(-2) – HD Monitor**

Standard Timings RGB, DVI

VGA, SVGA, XGA, WXGA, SXGA, SXGA+, UXGA, WUXGA will be detected with nearly any timing. Depending on the image present at detection it may be necessary to adjust size, position, (frequency and phase only RGB) once.

Output Gamma

Various gamma settings are supported and detected by the monitors. Gamma cannot be set or adjusted manually.

Monitor classification

Tolerances according to EBU Tech 3320 Grade 3 Monitor

Aspect Ratio

The aspect ratio of the input signal is retained as standard. All dual input modes retain the aspect ratio of the respective source.

User Interface

The individual settings can be made using a control with visual feedback via an on-screen display (OSD).

The following functions are comprised in the relevant menus and can be set via the user interface:

- Input Selection // Automatic Input Search
- Brightness, Contrast, Saturation, Hue, Sharpness, Noise reduction level
- H-Mirror, V-Mirror, H&V-Mirror, 90° Rotation, 270° Rotation
- Auto Setup (Full, Phase, Geometry)
- Scaling for main channel (Native, Best Fit, Fit Height, Fit Width)
- Underscan/Overscan to expand main image continuously from -50% to 150%

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
**9619NB / 9626NB(-2) /
9627NB(-2) – HD Monitor**

- The aspect ratio is set automatically
- PiP (fixed PiP arrangements all four corners, three PiP sizes, horizontal and vertical split screen)
- Power Safe configuration
- OSD Language Options
- OSD time out
- Control interface lock down mode – by remote command
- Info (Serial Number, Run Time, Backlight Run Time, SW Revision, HW Revision, Signal detection status and format at various inputs)
- Locking the OSD (see page 20)

Most of these settings are specific to the particular input selected, i.e. input specific color, scaling, brightness, etc. settings.

Control Interface

The control interface (located at the lower right corner on the monitor) for all user interface functions, works similar to an control wheel. The round out circle is engraved into the glass and etched over. A sliding circular action or even a direct touch in this area serves as the menu decide function.

Touching the "ok"-button of the element corresponds to a select action. It serves to select menu items and to drill down or move up in the menu hierarchy.

The LED does not light up if the control is inactive. When the control is activated the LED is permanently illuminated (see Unlocking the OSD, page 20).

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
OSM (On-Screen-Manager)

The following images roughly describe "how to use the OSM" implementation.

Touching the center of the "control wheel" will open up the main menu

Note: Pressing the control wheel for 12 seconds will reset the monitor.

Note: It is not allowed to touch the control wheel of the monitor during the boot-up process.

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You could slide to or even direct touch the function (i.e. Picture) you want.

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After this you have to touch and release the center ("ok"-button) again.

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Now you are one step down in menu hierarchy.

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
Again you could slide to, or direct touch the function (i.e. Position) you want to select.

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
After this you have to touch and release the center ("ok"-button) again

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
Now you are down in menu hierarchy to the menu "Position".

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To change the Picture Position you only have to touch several times to move the Picture into the desired direction (i.e. right) or leave your Finger at the desired area to start a repeat funktion until you remove your Finger.

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
Unlocking the on-screen display (OSD)

The OSD is automatically locked after 120 seconds in order to prevent the settings being accidentally changed during cleaning, for example. The default setting of 120 seconds can be set to a time between 30 and 600 seconds depending on your requirements ('MainMenu'>>Setup – Menu Lock').

Proceed as follows to unlock the OSD:

1. Touch the button in the middle of the on-screen display.
2. Move your finger clockwise on the control wheel through half a circle to the bottom of the wheel.
3. Move your finger counterclockwise on the control wheel through half a circle to the bottom of the wheel.
4. Repeat step 2.

The OSD is then unlocked automatically.

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
Horizontal and vertical mirroring

The source can be mirrored horizontally as well as vertically and/or rotated 180° (= mirrored horizontally and vertically at the same time).

The following rotation and mirror modes can be selected via the user interface:

- Mirror: mirror horizontally
- Flip: mirror vertically
- Rotate 180°: mirror horizontally and vertically


Note: The 'Mirror and Flip function' is only available when PiP, PbP and PoP are deactivated (see Menu, p. 24).

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Dual Display Modes

Dual Display refers to the simultaneous display of two input sources in various configurations. The size and position of the two input sources can be configured via the user interface. The functions for switching between primary and secondary sources can also be selected via RS232. The following table illustrates the various input modes and their respective status as a primary or secondary source.

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Dual Sources

Dual source refers to the simultaneous display of two input sources according to the above table.

The PIP positions below serve as one example. All four corners can be selected via the user interface.

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Remote Control

An RS232 interface and Ethernet 100Mbit/s connection (RJ-45) allows remote control of all functions of the monitor through a suitable driver software from a PC or other device.

Boot Time & Welcome Screen

The short startup time (from power on to display of a source without input search) is optimized for fast response.

During the startup time the KARL STORZ logo is displayed.

Signal detect and lock times are kept to the absolute minimum.

Power Save Mode

Two states – back to KARL STORZ logo, sleep with blank screen.

BIOS (Update)

The BIOS can be updated on site using a conventional notebook with an Ethernet 100Mbit/s connection (RJ-45) or RS232.

Languages

The following languages are supported in the user interface:

English
German
French
Italian
Swedish
Spanish
Dutch
Polish
Russian
Danish
Japanese
Chinese (Mandarin)
Korean

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Cleaning and care of the HD monitor

Warning: Always disconnect the power cord before cleaning the HD monitor.

- Wet a disposable cloth with a mild cleaning solution (wring out excess well) before wiping the screen.
- Avoid applying heavy pressure to the LCD monitor to prevent permanent damage to the display.
- Do not expose the LCD monitor directly to rain, water, moisture or sunlight.

Wipe-down disinfection

The design allows only for mechanical wipe-down disinfection. The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

Surfaces and objects should be wiped until fully coated. Leave for the specified exposure time and then wipe with a disposable cloth if necessary.

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Maintenance

Performance of preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the camera's reliability and extending its useful operating life. Maintenance services and information on maintenance contracts can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

- Do not attempt to carry out maintenance work yourself. Incorrect maintenance may invalidate the manufacturer's warranty. Only allow qualified maintenance personnel to carry out maintenance work.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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**9619 NB / 9626 NB(-2) /
9627 NB(-2) – HD Monitor**
Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
 Abt. Reparaturservice
 Dr. Karl-Storz-Str. 34
 78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**Electromagnetic Compatibility
 (EMC) Information**

Warning: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The 9619NB/9626NB(-2)/9627NB(-2) corresponds to EN 60601-1-2: 2007 and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The 9619NB/9626NB(-2)/9627NB(-2) generates or uses RF energy only for its internal functioning.

Note: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment or systems or non-medical electric equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Guidance and manufacturer's declaration – electromagnetic emissions

9619 NB / 9626 NB(-2) / 9627 NB(-2) is intended for use in the electromagnetic environment specified below. The customer or user of the 9619 NB / 9626 NB(-2) / 9627 NB(-2) should ensure that it is used in such an environment.

KARL STORZ touch screen monitor model 20090437 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.

Emissions test	Compliance	Electromagnetic environment – guidance
Radiated emission CISPR 22 EN 55022: 2006 +A1:2007 30 MHz ... 1 GHz - Class B 1 GHz ... 6 GHz - Class B	Class B 30 MHz ... 6 GHz	9619 NB / 9626 NB(-2) / 9627 NB(-2) is suitable for use in all establishments (e.g. hospitals and medical practices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Conducted at mains terminal EN 55022: 2006 +A1:2007 - 0.15 MHz ... 30 MHz Class B	Class B 0,15 MHz ... 30 MHz	
Harmonic emissions IEC 61000-3-2: 2005 +A1:2008 +A2:2009	100 Hz ... 2 kHz	
Voltage fluctuations/flicker emission IEC 61000-3-3: 2008	$P_{st} \leq 1,0$; $P_{et} \leq 1,0$; $d_c \leq 3,3\%$; $d_{max} \leq 4\%$; $d_{(t)} \leq 3,3\%$ (t=500 ms)	

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Guidance and manufacturer's declaration – electromagnetic immunity

9619NB/9626NB(-2)/9627NB(-2) is intended for use in the electromagnetic environment specified below. The customer or user of the **9619NB/9626NB(-2)/9627NB(-2)** should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2: 2008	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4: 2004	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5 : 2005	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle <40% U_T (60% dip in U_T) for 5 cycles <70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies <40% U_T (>60% dip in U_T) for 5 cycles Complies <70% U_T (>30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8: 2004	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

9619NB/9626NB(-2)/9627NB(-2) is intended for use in the electromagnetic environment specified below.
The customer or user of the 9619NB/9626NB(-2)/9627NB(-2) should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the 9619NB/9626NB(-2)/9627NB(-2) including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Recommended separation distances between portable and mobile RF communications equipment and the 9619NB/9626NB(-2)

The **9619NB/9626NB(-2)/9627NB(-2)** is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the **9619NB/9626NB(-2)/9627NB(-2)** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **9619NB/9626NB(-2)/9627NB(-2)** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.2	0.2	0.4
0.1	0.4	0.4	0.8
1	1.2	1.2	2.4
10	3.7	3.7	7.4
100	11.7	11.7	23.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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Firmenstempel/Unterschrift:

WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:

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To be filled out by instrument owner:

Return address/company stamp:

Field of application:

Type of instrument:

Serial no.:

Purchase date:

Signature/Date:

Please attach
sufficient
postage

Reply card

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen
Germany

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Absender/Firmenstempel:

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Geräte-Typ:

Serien-Nr.:

Kauf-Dat.:

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Antwortkarte

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 7461 708-0, Fax: +49 7461 708-105
E-Mail: info@karlstorz.de

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 30 30 69090, Fax: +49 30 30 19452

KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseia.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America Do Sul Ltda.
R. Joaquim Floriano, 413 – 20º andar –
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsatragrand 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Vallmotie 1 B, 4. kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231
KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopleer.be

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskope Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopiija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 31 4250800, Fax: +40 31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Gounari 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

OOD KARL STORZ Endoscopy – WOSTOK
Darbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 1999390, Fax: +961 1 1999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061
Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
Sary-Arka str., 6, BC "Arman", off.430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 4 4471230
Fax: +971 4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeijon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Belsin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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/Karl Storz GmbH & Co. KG/