2

STORZ

KARL STORZ — ENDOSKOPE

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20 1337 20-1 D-LIGHT P

INSTRUCTION MANUAL

20 1337 20-1 D-LIGHT P

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Important information for users of KARL STORZ Instruments

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the D-LIGHT P system. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

© All product images and descriptive texts are the intellectual property of KARL STORZ GmbH & Co. KG

Utilization and reproduction by third parties requires express authorization by KARL STORZ GmbH & Co. KG
All rights reserved.

Trademarks

TRICAM® and **SCB**® are registered trademarks of KARL STORZ.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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Controls, displays, connectors, and their uses

- ① Power switch ("O" = off)
- ② Lamp replacement indicator
- ③ Warning light: lamp service life
- ④ Standby button, with pilot lamp
- ⑤ Bar graph of brightness
- ⑥ ± buttons for brightness control
- ⑦ Digital readout of brightness
- ⑧ M button for changing the application/
operating mode (Service Mode) or for Manual
Mode/Autovideo changeover
(*Std' application)
- ⑨ White light/ICG fluorescence light mode
(color light) changeover button
(Standard application only has white light mode)
- ⑩ Footswitch connector
- ⑪ Light exit, unscrewable
- ⑫ Label strip (replaceable)
- ⑬ Cooling vents
- ⑭ Potential equalization connector
- ⑮ Power fuse holder
- ⑯ Power supply cord receptacle
- ⑰ * SCB connectors
- ⑱ Insulation foils
- ⑲ Lamp module (has 'coding' to avoid confusion
with lamp modules from other light sources;
see Section 'Lamp module replacement')

* The KARL STORZ SCB interface (KARL STORZ
Communication Bus), based on the CAN field bus,
permits remote control of equipment functions, as
well as remote display of equipment parameters.

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ICG = Indocyaningrün
IR = Infrarot
NIR = Nahes Infrarot

Abbreviations used

ICG = IndocyanineGreen
IR = Infrared
NIR = Near Infrared

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	Read the instructions carefully before operating the equipment!
	Standby mode
	Automatic brightness setting (Δ)
	Manual brightness setting
	Operating mode with filter
	Light exit
	Lamp brightness (%)
	Footswitch connector
	Lamp service life warning light (lamp operating time > 350 hours)
	Lamp replacement indicator (lamp operating time > 400 hours)
	Potential equalization connector
	Applied part of type CF
	Use of electronic devices to avoid environmental pollution (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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Images of the unit	IV
Controls, displays, connectors, and their uses	V
Symbols employed	VI

General information

Description of the system	2
---------------------------------	---

Safety instructions

Warnings and cautions	4
Intended use	7
User qualification	9
Safety precautions at the site of installation	9
Safety precautions when using the unit	9

Installation and operating instructions

Unpacking the equipment	10
Basic equipment	10
Installing and connecting up the unit	10
KARL STORZ SCB fluorescence system	13
Operating the unit	15
Manual brightness setting	16
Automatic brightness setting	16
Changing the light mode	17
Standby operation	18
Function of the warning and indicator lights	18
Service mode	19

Maintenance

Lamp module replacement	21
Fuse replacement	23
Reprocessing	24
Maintenance and safety check	25
Servicing and repair	26
Disposal	26
Repair program	27
Limitation of liability	28
Manufacturer's warranty	29

Technical description

Troubleshooting	30
Technical data	32
Technical documentation	33
Block diagram	34

Spare parts,

recommended accessories

Spare parts list	35
Recommended accessories	35

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Description of the system

The D-LIGHT P system is a high-intensity light system which provides up to 2 different types of application (see table).

NOTE: The abbreviations used are explained on page V.

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For the D-LIGHT P system for fluorescence diagnostics, there are rigid and flexible endoscopes as well as the KARL STORZ ENDOVISION® TRICAM® SL II camera system and IMAGE 1 SPIES™, especially modified for fluorescence diagnostics. The KARL STORZ ENDOVISION® TELECAM SL II FI camera control unit is used in combination with fluorescence videoscopes.

The D-LIGHT P system may be used only with recommended accessories.

Refer to the system description for KARL STORZ Fluorescence Imaging (FI) Systems (see instruction manual 96206144DER) for the recommended configuration.

As with virtually all diagnostic methods, false positive and false negative results can also occur in ICG-supported visualization. Validation using other methods is necessary.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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D-LIGHT P system

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel. If unauthorized persons take off covers which are not explicitly indicated as being removable, the guarantee will be invalidated.

WARNING: Danger of burns. The lamp and lamp module may be very hot.

WARNING: Grounding reliability can only be achieved when the equipment is connected to a properly installed "Hospital Only" or "Hospital Grade" receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: The electrical installation of the operating room must comply with the applicable IEC requirements.

WARNING: Keep out of reach of patients.

WARNING: Always unplug the unit before all maintenance work (e.g. cleaning).

WARNING: As with virtually all diagnostic methods, false positive and false negative results can also occur in ICG-supported visualization. Validation using other methods is necessary.

WARNING: The use of endoscopes and light

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Safety instructions

Warnings and cautions

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

CAUTION: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

CAUTION: Do not block any air vents. Danger of overheating. Do not place the unit directly against the wall since the air vents are on the left-hand side and rear of the unit (light source viewed from front) (see also page 11).

CAUTION: When replacing fuses, use only fuses of the correct rating.

CAUTION: Do not touch lamp ceramic or window with bare fingers; use gloves or cloth or tissue to handle lamp. Touching the lamp with bare fingers may result in starting difficulty.

CAUTION: Federal law restricts this device to sale by or on the order of a physician (USA only).

Light cables

WARNING: Danger of glare. Never look into the open end of a light cable.

WARNING: Danger of glare! Do not operate the light source without a light cable adaptor and never look directly into the light outlet point.

WARNING: Danger of burns. Contact with the open end of a light cable connected to the D-LIGHT P system can possibly lead to burns.

WARNING: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and at the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera (standard mode only).

WARNING: Never leave an illuminated light cable with or without an endoscope connected, on or near a patient; the risk of patient burns can result.

CAUTION: Do not place connected light cables on combustible materials such as textiles (drapes)

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WARNING: Only the fluid light cables developed specially for the fluorescence endoscopes may be used in fluorescence applications. It is of the utmost importance to observe the correct assignment of the fluid light cable to the corresponding fluorescence endoscope in order to prevent any possible tissue burning. The respective combination of fluid light cable to fluorescence endoscope can be found in the KARL STORZ catalog for the corresponding medical specialty.

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NOTE: Please read the Electromagnetic Compatibility Information in the appendix (p. 37-49).

NOTE: The ready-reference guide is not intended to replace the unit's instruction manual. Refer to the unit's instruction manual for more detailed operating instructions.

NOTE: KARL STORZ accepts no responsibility whatsoever for the administration and effect of suitable marker substances (ICG; refer to application mode 2).

NOTE: When disposing of parts, equipment components or an entire instrument, the respective national regulations must be observed.

NOTE: The UV radiation in the 350-380 nm range complies with the threshold limit values specified by the ACGIH 'Ultraviolet Radiation' standard*.

*The standard pertains to the skin and eyes.

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Intended use

The D-LIGHT P system is used as a light source for endoscopic procedures, in particular:

- as a standard high-intensity light system for use with rigid and flexible endoscopes (Standard – Std application)
- as a fluorescent excitation light system for exciting the fluorescence of ICG for visualizing the superficial circulation of soft tissue and/or organs. (ICG application)

Use of the unit in fields other than those indicated

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The D-LIGHT P system may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

As a marker substance for the visualization of circulation in tissue, for example in the recognition of tissue areas with no or poor circulation, only use pharmaceuticals specially approved for this purpose according to the manufacturer's specifications (e.g., ICG).

KARL STORZ accepts no responsibility whatsoever for the administration and effect of suitable marker substances.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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Contraindication

When using the D-LIGHT P system for the visualization of peripheral circulation, the medical contraindications apply which are listed for the marker substance ICG, and use is contraindicated if endoscopic examination techniques contradict the current state of academic knowledge.

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User qualification

The D-LIGHT P system may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit. Defective color vision hinders normal application of the system.

NOTE: KARL STORZ strongly recommends participation in advanced clinical training courses. Detailed information regarding the date and location of these courses can be requested from the KARL STORZ company.

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Safety precautions at the site of installation

The unit may only be used in medical rooms where the installations comply with the applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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
Safety precautions when operating the unit

The user must be convinced of the operational safety and proper condition of the unit, as well as of the relevant accessories, before using. During a procedure using the D-LIGHT P system the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as

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Unpacking the equipment

Carefully unpack the D-LIGHT P system and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Basic equipment

- 1 D-LIGHT P system 20 1337 20-1
- 1 Footswitch including connecting cable, 200141 30
- 1 SCB connecting cable 200901 70
- 1 Power cord 400 A (240V~) Intl. or Power cord 400 B (110V~) U.S.
- 1 Instruction manual
- 1 Cover for instruction manuals
- 1 Operator's quick reference guide

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Installing and connecting up the unit

NOTE: The D-LIGHT P system including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

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WARNING: The D-LIGHT P system is not intended for use in hazardous zones. Do not operate the unit inside the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
Fitting a label strip to the front panel


i NOTE: The appropriate controls/indicators on the unit are designated on the label strip.

Undo the left-hand front side cover with a screwdriver.

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Remove any unmarked label strips. Select a label strip (currently available in 6 languages) and slide it between the plastic and the metal front.

Fold the strip along the broken line to an angle of 90° and approx. 2 mm further on (toward the end of the strip) fold it back to an angle of 90° so that the side cover can be refitted without any obstructions.

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Set the unit on a flat surface.

i CAUTION: Do not block any air vents. Danger of overheating. Do not place the unit directly against the wall since the air vents are on the left-hand side and rear of the unit (light source viewed from front). Ensure that the left hand side of the unit is located at least 10 cm away from the wall. The rear of the light source must remain uncovered in order to ensure optimal air circulation (do not use a closed equipment trolley).

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The unit is equipped with a connector ⑭ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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WARNING: The device may only be operated at the voltage stated on the rating plate.

Connect the power cord.

WARNING: Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

WARNING: The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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
Insert light cable until it clicks into socket ⑪.

WARNING: Danger of glare! Never look directly into the light outlet point. Do not operate the light source without a light cable adaptor (see arrow)

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Connect the fluorescence endoscope with the light cable.

NOTE: Please draw the suitable light cables and endoscopes from the KARL STORZ catalog.

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
Connect the cable from the footswitch to its mating connector ⑩.

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Connect SCB cable to one of the sockets ⑪ (the two SCB sockets are interchangeable). Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ SCB System Instruction Manual).

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KARL STORZ SCB fluorescence system

**D-LIGHT P – TRICAM® SL II/ IMAGE 1
HUB™ HD/ TELECAM SL II FI**

Communication/synchronization between the SCB versions of TRICAM® SL II 20223020-1 and/or TELECAM SL II FI 20213020-1 Camera Control Units or IMAGE 1 HUB™ HD 22201020-1xx / 22202020-1xx with (light source) D-LIGHT P 20133720-1 is carried out via the KARL STORZ SCB System (KARL STORZ Communication Bus). The KARL STORZ fluorescence system thus always operates in sync.



This means that if the filter mode is changed on the D-LIGHT P, the camera control unit adapts automatically, and vice versa.

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Connecting the SCB D-LIGHT P to the SCB camera control unit 202230 20-1, 202130 20-1, TC200/TC300 or 22201020-1xx/22202020-1xx.

At the rear of each SCB unit there are two SCB sockets. With the enclosed SCB cable the D-LIGHT P system is connected to the camera control unit (CCU) (the two SCB sockets are identical).

! WARNING: Ensure that there is no SCB cable which is only connected at one end, or else there may be transmission interference within the system. In case of malfunctions, please contact the Service Department.

! WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely. Please read the instruction manual for the KARL STORZ SCB system carefully and follow the instructions exactly.

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Putting the SCB system into operation

It does not matter in what sequence the units of the SCB system are switched on. When the camera control unit is switched on, the following SCB message appears: New Device ...

As standard, the system starts in the camera mode last selected.

! NOTE: A camera connection is required to display the ICG signal; consequently, it must be ensured that the camera and light source are in ICG mode

! NOTE: Changing the application is described in the section 'Service Mode – changing application' on page 19.

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When a fluorescence mode is selected, e.g., 'ICG (W)' is displayed on the monitor. 'W' is the abbreviation for 'white'. That is, the D-LIGHT P first switches to the reduced white light mode if an ICG mode is selected.

Pressing the 'blue' function key or Ctrl-1 on the keyboard or footswitch or changeover button ⑨ switches over to the ICG fluorescence light mode. 'ICG (C)' is displayed on the monitor. 'C' is the abbreviation for 'Color'.

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Operating the unit

Activate power switch ①.

The self-test starts (all LEDs light up). The lamp ignites. The unit version B is displayed briefly. The application that was active when the unit was switched off is reactivated (e.g.: ICG). The lamp operating time is displayed briefly.

① NOTE: With a lamp operating time > 350 h, a signal sounds for 2 s and the lamp service life warning light ③ glows. With a lamp operating time > 400 h, an intermittent signal sounds for 10 s and the lamp replacement indicator ② glows (also refer to the section, 'Function of the warning and indicator lights on the front panel' on page 18).

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① HINWEIS: Hat das NOVRAM (Speicher)

The D-LIGHT P system always initializes to 100% brightness in the reduced white light mode* of the application or to 25% in 'Std' mode.

① *NOTE: To guarantee efficient fluorescence excitation, it is generally not recommended to reduce the brightness from 100%, as the fluorescence excitation modes only emit light in a small wavelength spectrum.

① NOTE: If the NOVRAM (memory) has a fault, the error message 'E03' (see section 'Troubleshooting' on p. 31) appears instead of the application and operating hours display.

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Manual brightness setting*

With the $\pm$ buttons (6), the brightness can be adjusted from 5% to 100%. The current brightness level is shown on the bar display (5) and the digital display (7).

- * NOTE:** To guarantee efficient fluorescence excitation, it is generally not recommended to reduce the brightness from 100%, as the fluorescence excitation modes only emit light in a small wavelength spectrum.
- NOTE:** In the white light mode of an application, the digital display (7) indicates the brightness in percent. In the color mode of an application, the abbreviation of the application set (e.g. 'ICG') is displayed. The brightness emitted in percent appears only after pressing one of the $\pm$ buttons (6). Approx. 1.5 seconds after last pressing the button, the abbreviation of the application reappears on the display.
- NOTE:** The brightness can be adjusted separately for white light and color light of the active application. The Standard application has no color mode.

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Automatic brightness setting (only with Standard application)


By pressing the mode (M) button (8) either the MANUAL or AUTO(MATIC) operating mode may be selected, as shown by the symbols.

- NOTE:** The 'AUTO' operating mode (Autovideo) only functions together with a camera by communication via SCB.

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With the $\pm$ buttons (6), the video level (Autovideo Mode) can be adjusted from -10 to +10. The setting is shown on the bar display (5) and the digital display (7).

- NOTE:** In the 'Std' (= Standard) application, the D-LIGHT P system can be used as a standard light source for endoscopy.

It is only possible to change the set operation type/application in the service mode (see Section 'Service Mode – Changing application' on page

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Changing light mode (switching between white and color light mode)

Press the changeover button ⑨ or footswitch to change the light mode (white light <-> color light) within the active application. This change is indicated by a short acoustic signal. If the color mode is active, the digital display ⑦ shows the abbreviation of the application (see page V). In addition, the LED in the changeover button lights up.

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① NOTE: In the color mode of an application set 'ICG' is displayed. The brightness emitted in percent appears only after pressing one of the ± buttons ⑥. Approx. 1.5 seconds after last pressing the key, the abbreviation of the application reappears on the display.

① NOTE: The brightness can be adjusted separately for white light and color light of the active application. The Standard application has no color mode.

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The white light mode (reduced white light) is used for normal viewing during an ICG session.

① NOTE: In the 'Std' (= Standard) application, the D-LIGHT P system can be used as a standard light source for endoscopy.

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Changing application


Switching into a different application is described in the Service Mode section (see page 19).

Selectable applications:
Std and ICG

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Standby operation

By pressing the Standby button ④, light output can be reduced to a minimum value. After pressing the Standby button ④ again, light output returns to the previously selected value.

ⓘ NOTE: Constant switching on and off of the D-LIGHT P system leads to increased wear of the lamp. For brief interruptions in operation, the unit should therefore be switched to standby operation.

Pressing the + button ⑥ also ends standby mode. The light intensity can then be newly adjusted as from 5%.

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Function of the warning and indicator lights on the front panel

(also refer to note 'Operating the unit' on p.15)

With a lamp operating time of 350–399 hours, the lamp service life warning light ③ is lit on the front panel. This warning light only goes out when the operating hour counter of the XENON lamp is set to 0.

With a XENON lamp operating time of 400 hours or more, the lamp replacement indicator ② is lit on the front panel. This indicator light only goes out when the operating hour counter of the XENON lamp is set to 0 (see section 'Lamp module replacement' p. 21/22).

The indicator light ② is also lit in the event of a fault, e.g. if the XENON lamp has not ignited. The digital display ⑦ additionally displays an error code (E01). If the lamp still does not ignite after switching on again, change the lamp. The error codes are explained in 'Troubleshooting' on p. 30, 31).

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
Service mode

The Service mode allows setting of the various applications and various information concerning the operational status of the D-LIGHT P system to be requested.

To activate the 'Service Mode', the Standby button (4) must be held down for 8 seconds. The application set then appears on the digital display (7).

Changing application

To change the set application, press the M button (8) repeatedly until the application required is displayed.

Selectable applications:
Std or ICG

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
Display device information

If the + button (Ⓢ) is pressed repeatedly, the following information can be called up:

- Active application
- Operating hours of the lamp
- Number of lamp starts (ignitions)/10
- Operating hours of the unit /10
- Current operating temperature in °C
- Highest operating temperature reached to date in °C
- Error message codes 01 ... 10, 15 (only for service personnel)
- ROM version number
- Serial number (position 1... 5)
- Device type
- SCB device type
- PLD version number

By pressing the – button (Ⓢ) the displays can be called up again in reverse sequence.

(PLD version number ... SCB device type ... Lamp operating hours ... Active application)

Leaving Service Mode

To leave Service Mode, hold down the Standby button (Ⓢ) again for 8 seconds or switch off the main power switch (1).

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Lamp module replacement

WARNING: Switch off the unit. Disconnect the power plug.

Unlatch the quarter-turn fastener counterclockwise and pull the lamp door out.

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WARNING: Danger of explosion. The lamp is filled with xenon gas at very high pressure. It must therefore be handled with appropriate care and must only be replaced by qualified service personnel. High mechanical forces could cause the lamp to break.

Wear safety goggles! We urgently recommend wearing protective clothing (face protection, gloves, etc).

Fold insulation foils to the side (see arrow).

Pull both card extractors on the lamp module simultaneously outwards to unlock the lamp module. Pull the lamp module assembly straight out of the unit.

NOTE: The lamp module 20 1337 25 for the D-LIGHT P has a 'coding' (see arrow) to prevent confusion with lamp modules of other light sources due to its special spectral emission characteristic.

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Lamp module replacement

Replace insulation foils correctly over the lamp module.

Reattach the lamp door and fasten by a quarter turn of the fastener clockwise.

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
Write down the number of the XENON lamp operating hours indicated on the counter and keep a record as to the lamp's operating life (see guarantee information). Then reset the operating hour counter to zero.

To do so, change to Service Mode (→ press Standby button ④ for 8 s, see p.19) and press the + button ⑥ until the lamp operating hours are displayed.

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Then press the Standby button ④ and the – button ⑥ simultaneously for 8 seconds. This will reset the operating hour counter to zero.

NOTE: Counting of the hours run by the unit as a whole will not be affected by this.

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
NOTE: The defective lamp module can be sent back to KARL STORZ for disposal or for repair. It must be shipped in the original packaging (see also p. 26, "Disposal").

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
Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Remove the power fuse holder ⑮ with a screwdriver or other tool.

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! **CAUTION:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

	100...240 VAC
Power fuse	2 x T 8.0 AL250V

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Refit the power fuse holder ⑮.
Reconnect the power cord.
Test the unit for proper operation.

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Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

WARNING: Always pull out power plug before cleaning! Make sure the lamp has cooled down. Fluid must never be allowed to enter the unit.

Wipe-down disinfection of device and footswitch

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not

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Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety checks/repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section "Operation").

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/ repeat inspection, please refer to the latest issue of the relevant Service Manual.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

D-LIGHT P requires custom repair. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices.

All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Warranty for XENON lamp or lamp module

Independent of the general warranty conditions for our equipment, we guarantee 500 hours of life for our XENON lamps or lamp modules. Lamps or lamp modules that burn out within the first 50 hours will be replaced free of charge; for lamps or lamp modules that burn out after the first 50 hours, the price of the replacement lamp or module will be calculated according to the number of hours the original lamp was used.

This is applicable only if the defective lamp/lamp module, the serial number of the D-LIGHT P system, and a record of the lamp hours are presented.

NOTE: To ensure a consistently high light quality, replace the XENON lamp or the lamp module after 400 lamp hours ('Replace lamp' indicator ② is lit).

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall release us from any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.*

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Symptom:

No light is emitted. Lamp replacement indicator ② lights up.

Possible causes:

- 1 Lamp is defective (digital display ⑦ shows 'E01').
- 2 Defect in electronic circuitry.

Remedy:

- 1 Have a qualified electrician change the lamp.
- 2 Notify Service.

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Possible causes:

- 1 The NOVRAM (memory) is defective (E03), fault in the positioning of the filter disk (E04) or the attenuation filter (E07).
- 2 'E10': Overheating of the D-LIGHT P due to air vents being covered (lamp has shut down due to the unit overheating).

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Troubleshooting

WARNING: Always unplug the unit before all maintenance work.

Symptom:

Too little light.

Possible causes:

- 1 Lamp is defective.
- 2 Ends of light cable and/or endoscope are soiled.
- 3 Light cable is defective.

Remedy:

- 1 Have the lamp exchanged.
- 2 Clean ends of light cables, and light entrances and exits on the endoscope.
- 3 Replace light cable.

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Error codes

The digital display ⑦ shows the error code: This code is for the user's information or provides Service with preliminary information.

- 'E' flashes: current error
- 'E' glows continuously: stored error (occurred 'at some time or other')
- No 'E' before the number (e.g. -01): there is no error

- E01 Lamp does not operate
- E02 Not assigned
- E03 NOVRAM (memory) fault (affects other parameters and error messages)
- E04 Filter disk fault
- E05 Lamp service life warning
- E06 Unit overheating
- E07 Attenuation filter fault
- E08 Not assigned

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Standard compliance (for 20 1337 20-1)

**According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:**

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type CF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 37-49).

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Nach Medical Device Directive (MDD):

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Class II Medical Device in accordance with MDD 93/42/EEC

Directive compliance (for 20 1337 20-1)

According to Medical Device Directive (MDD):

This medical product belongs to Class II a

This medical product bears the CE mark in

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

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
Block diagram

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Spare parts, recommended accessories

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Recommended accessories

Item	Cat. no.
Fluorescence endoscopes	(see catalog; corresponding specialties)
Fluid light cable, for fluorescence endoscopes	(see catalog; corresponding specialties)
TRICAM® SL II Camera Control Unit	202230 20-1
TRICAM® PDD camera head with Parfocal-Zoom lens, PAL	202210 37
Same, NTSC	202211 37

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**Spare parts,
recommended accessories**
Recommended accessories

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Item	Cat. no.
TELECAM SL II FI Camera Control Unit	20213020-1
Compatible fluorescence videoscopes	(see catalog; corresponding specialties)
IMAGE 1 HUB™ HD Camera Control Unit	22201020-1xx/ 22202020-1xx
IMAGE 1 H3-Z FI Three-Chip FULL HD camera head	22220085-3
KARL STORZ video endoscope adaptor	487 VE

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**Recommended combination of light cable with
endoscope**

For safety reasons, the diameters of fiberoptic/fluid
light cables and endoscopes should be used in the
combination shown in the table opposite.

*** NOTE:** Special endoscopes may deviate
slightly.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The D-LIGHT P model 20 1337 20-1 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The D-LIGHT P model 20 1337 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment to an outlet on a circuit different from that to which the other

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Appendix

Electromagnetic Compatibility
(EMC) Information

 **WARNING:** The D-LIGHT P model 20 1337 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

 **WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

 **WARNING:** The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the D-LIGHT P model 20 1337 20-1. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
				Footswitch model 20014130

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Table 1

Guidance and manufacturer's declaration – electromagnetic emissions

D-LIGHT P model 20 1337 20-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the D-LIGHT P model 20 1337 20-1 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidelines
RF emissions CISPR 11	Group 1	D-LIGHT P model 20 1337 20-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The D-LIGHT P model 20 1337 20-1 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2

Guidance and manufacturer's declaration – electromagnetic immunity

D-LIGHT P model 20133720-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the D-LIGHT P model 20133720-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – earth	± 1 kV voltage outer conductor – outer conductor ± 2 kV voltage outer conductor – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T^* (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T^* (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

D-LIGHT P model **20 1337 20-1** is intended for use in the electromagnetic environment specified below.
The customer or user of the D-LIGHT P model **20 1337 20-1** should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the D-LIGHT P model 20 1337 20-1 , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional

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Table 6

**Recommended separation distances between portable and mobile RF communications equipment and the
D-LIGHT P model 20 1337 20-1**

The D-LIGHT P model 20 1337 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseaa.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufrieren 627 6º Piso, B1638 CAA - Vicente Lopez
Buenos Aires, Argentina
Phone: +54 11 4222000, Fax: +54 11 4222001
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai (Zhangjiang)
Room 701A, 7F Building 5, No. 3000 Longdong
Avenue, Zhangjiang Hi-Tech Park, Pudong District,
Shanghai, 201203, People's Republic of China
Phone: +86 21 61009688, Fax: +86 21 61009686
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

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