Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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**Important information for
users of KARL STORZ
devices**

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**1 Important information for
users of KARL STORZ
devices**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
products, this product too is the result of years
of experience and careful manufacturing. You
and your organization have decided in favor of a
modern high quality product from KARL STORZ.

This instruction manual is intended to serve
as an aid in the proper installation, connection
and operation of the ENDOFLATOR® 50. All of
the necessary details and actions are clearly
explained.

Please read these instructions carefully. Keep this
manual in a convenient and conspicuous location
and in its protective casing close to the device.

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


CAUTION: Federal (USA) law restricts this device
to sale by or on the order of a physician.

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Controls, displays, connectors, and their uses

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3 Controls, displays, connectors, and their uses

- ① Power switch ("O" = off)
- ② Touch screen
- ③ Insufflation connection to patient
- ④ Connection for patient tube heating
- ⑤ Room temperature sensor
- ⑥ Gas connection (American connection)
- ⑦ SCB connectors
- ⑧ Service interface (concealed)
- ⑨ Ethernet interface (concealed)
- ⑩ Potential equalization connector
- ⑪ Power cord socket
- ⑫ Line fuse holder
- ⑬ Holder for CO₂ bottle
(optional – cat. no. UI005)

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4 Symbols employed

4.1 Symbols on the device

	Follow instructions for use
	ON
○	OFF
↓	Potential equalization connector
	Applied part of type CF
~	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

Symbols on label and packaging of the device

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. This can be downloaded from www.karlstorz.com.

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4.2 User interface symbols

	Settings
	Procedure list
	Working area
	Start/Stop
	Sensitive Mode
	High Flow Mode
	Gas bottle
	House connection
	Flow
	Pressure
	Audio signals on
	Audio signals off
	House connection, input pressure acceptable
	House connection, input pressure too high or too low
	Gas bottle; input pressure above 40 bar
	Gas bottle; input pressure between 40 and 30 bar
	Gas bottle; input pressure between 30 and 10 bar
	Gas bottle; input pressure below 10 bar
	Gas bottle; input pressure too high
	Display of gas consumption/ Reset the display by clicking on the -0- button
	Patient heating plugged in and in operation
	Defective patient heating/heating control

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	Manufacturer
	Article catalog no.
	Product must not be reused
	Product must not be resterilized
	CE mark
	Do not use the product if the packaging has been damaged
	Consult instructions for use
	In the USA, federal law restricts this device to sale or use by or on the order of a physician
	Store in a dry place
	Store out of direct sunlight
	Sterilized with ethylene oxide
	Use by date of the product
	Production batch label

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1 Important information for users of KARL STORZ devices	III
2 Images of the equipment	IV
3 Controls, displays, connectors, and their uses	V
4 Symbols employed	VI
4.1 Symbols on the device	VI
4.2 User interface symbols	VII
4.3 Gas filter label symbols	VIII
5 General information	4
5.1 Description of the device	4
5.1.1 Property rights	5
6 Safety instructions	6
6.1 Explanation of warnings and cautions	6
6.2 Intended use	11
6.2.1 Indication	11
6.2.2 Contraindications	12
6.3 User qualification	14
6.4 User profile of physician and assistants	14
6.5 Intended conditions of use	14
6.5.1 Use	14
6.5.2 Other intended conditions	14
6.5.3 User position	15
6.6 Safety precautions at the installation site	16
6.7 Safety precautions when operating the device	16
6.8 Safety features	17
6.8.1 Self-test	17
6.8.2 Monitoring during operation	17
7 Installation and operating instructions ..	19
7.1 Unpacking the equipment	19
7.2 Basic equipment	19
7.3 Installing and connecting the device	19
7.3.1 Connecting the ground line	20
7.3.2 Connecting the power cord	20
7.3.3 KARL STORZ SCB	21
7.3.4 Affixing the CO ₂ bottle to the device (optional)	22
7.3.5 Connecting the CO ₂ bottle	22

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8.3	Maintenance and safety check.....	47
8.3.1	Maintenance.....	47
8.3.2	Safety check/repeat inspection according to IEC 62353.....	47
8.4	Servicing and repair.....	48
8.5	Disposal.....	48
8.6	Repair program.....	49
8.7	Important information.....	49
8.8	Limitation of liability.....	50
8.9	Warranty.....	50
9	Technical description.....	51
9.1	Signals.....	51
9.1.1	Visual signals.....	51
9.1.2	Acoustic signals.....	52
9.2	Troubleshooting.....	54
9.3	Technical data.....	57
9.3.1	Standard compliance (for UI500).....	58
9.3.2	Directive compliance (for UI500).....	58
9.4	Technical documentation.....	59
10	Spare parts, recommended accessories.....	60
10.1	Spare parts/accessories.....	60
10.2	Accessories.....	61
11	Electromagnetic Compatibility (EMC) information.....	62
12	Subsidiaries.....	78

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5 General information

5.1 Description of the device

The ENDOFLATOR® 50 is an insufflation device for universal application in laparoscopic and thoracoscopic examinations and operations as well as endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy). In addition, the ENDOFLATOR® 50 is designed to expulse ambient air in open and endoscopically assisted cardiac surgery as well as endoscopic vessel harvesting. The very latest technologies for pressure and flow measurement and control allow for various operating modes, including modes tailored to specific situations.

The 'SENSITIVE' and 'HIGH FLOW' modes enable various applications to be accommodated. Easy handling and clearly arranged controls together with several safety circuits ensure maximum patient safety.

For particularly sensitive areas such as pediatric applications, the device features a 'SENSITIVE' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min). This ensures considerably enhanced patient safety during operations on neonates and infants, among others.

To quickly compensate for the considerable loss of gas occurring in complex laparoscopic operations, the ENDOFLATOR® 50 has been designed for a high flow rate of up to 50 l/min.

Use of an insufflation tube with integrated gas heater can reduce the risk of hypothermia.

The ENDOFLATOR® 50 is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface.

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The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

Thus the current patient pressure, current flow and the amount of gas used during the current intervention are displayed simultaneously during operation.

The user is able to save his standard settings as a procedure (max. 30).

5.1.1 Property rights

This product is protected in the USA by (at least one of the following) US Patent(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this instruction manual through carefully and follow its instructions exactly. The words Warning, Caution, and Note convey special meanings. Wherever they are used in this instruction manual, the following text must be carefully reviewed to ensure the safe and effective operation of this device. In addition, to make the signal words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note provides special information regarding the operation of the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel, or yourself at risk.

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Safety instructions

Warnings and cautions

-  **WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.
-  **WARNING:** Set up the device out of reach of patients.
-  **WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.
-  **WARNING:** Combinations of medical devices are only assured to be safe if
 - they are identified as such in the respective instruction manuals or
 - the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1).
-  **WARNING:** Grounding reliability can only be achieved when the device is connected to a 'Hospital only' or 'Hospital Grade' outlet. Check the plugs and cables regularly and do not use if damaged.
-  **WARNING:** To avoid the risk of an electric shock, this device may only be connected to a power supply network with a protective conductor.
-  **WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.
-  **WARNING:** For safety reasons, do not work with a leaky device. There is a risk that the intraabdominal pressure will increase in an uncontrollable manner.
-  **WARNING:** Always rinse out the system with 1l CO₂ before use.

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Safety instructions

Warnings and cautions

-  **WARNING:** Make sure there are no external forces acting on the cavity. This can increase the intracavitary pressure and/or cause pressure fluctuations.
-  **WARNING:** No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.
-  **WARNING:** Always use the lowest possible pressure and flow for insufflation.
-  **WARNING:** Always check which operating mode is set. For safety reasons, use the 'Sensitive' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.
-  **WARNING:** Any treatment may only be performed if visual observation of the action of the device is ensured.
-  **WARNING:** After 24 h of continuous operation, a self-test must be performed. Take note of the warnings to ensure that the device is fully functional.
-  **WARNING:** The device must not be used if either the touch screen or display is defective.
-  **WARNING:** If the touch screen does not react when pressed, you must turn the device off.
-  **WARNING:** Do not press several points on the touch screen simultaneously.
-  **WARNING:** Only sterile CO₂ gas may be used.
-  **WARNING:** You must secure CO₂ bottles connected up to the device to prevent them from falling over.
-  **WARNING:** Follow the safety instructions when changing the CO₂ bottle.
-  **WARNING:** The CO₂ bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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Safety instructions

Warnings and cautions

-  **WARNING:** Always unplug the device from the power supply before carrying out any cleaning and maintenance work.
-  **WARNING:** Danger of electric shock! Do not open the device. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unauthorized persons will void warranty claims.
-  **WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.
-  **WARNING:** Only operate the device under the specified environmental conditions.
-  **WARNING:** Always keep a spare device to hand in case the first device should fail.
-  **WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.
-  **WARNING:** The device may only be operated with patient tube heating (heatable insufflation tube sets) supplied by KARL STORZ. The connection socket ④ on the front panel is intended solely for connection of the patient tube heating.
-  **WARNING:** If the patient heating symbol lights up red, this indicates that the control unit of the patient tube heating is defective. In this case, do **not** connect the cable for the patient tube heating to the device.
-  **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

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Safety instructions

Warnings and cautions

-  **WARNING:** Always use a hydrophobic bacterial filter to prevent cross-contamination of patients.
-  **WARNING:** You must insert a sterile CO₂ filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. Replace the CO₂ filter after every use.
-  **WARNING:** Ensure that liquid cannot flow back into the device. Position the device higher than the patient and use a hydrophobic filter between the insufflation connection and insufflation tube.
-  **WARNING:** Check the packaging of the tube set for damage! If the packaging is damaged, the sterility of the tube set may be impaired and the tube set must not be used.
-  **WARNING:** Follow the country-specific regulations/laws for the disposal of accessories.
-  **WARNING:** Avoid allowing fluids to enter the housing. Do not store liquids on or directly above the device. If, despite the precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies for the formation of condensation).
-  **WARNING:** The device may only be operated at the voltage stated on the rating plate.
-  **WARNING:** When replacing fuses, use only fuses of the correct rating.

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Before using the device on the patient for the first time it is imperative that you be acquainted with how the device operates and is controlled.

6.2 Intended use

CO₂ insufflators and their accessories are used to create and maintain a cavity during diagnostic or therapeutic endoscopic interventions as well as to expulse ambient air during open, endoscopically assisted or endoscopic interventions.

6.2.1 Indication

The insufflator is used to create a cavity in the following diagnostic and therapeutic interventions:

- Laparoscopy
- Thoracoscopy
- Endoscopy of the upper and lower gastrointestinal tract: (e.g., TEO®, colonoscopy)

In addition, the ENDOFLATOR® 50 is also designed to expulse ambient air in open and endoscopically assisted cardiac surgery and endoscopic vessel harvesting.

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Use of the device in fields other than those indicated above is not allowed for safety reasons. The ENDOFLATOR® 50 may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

WARNING: No responsibility for the safe operation of the device can be accepted if any tubing systems other than those specified are used.

WARNING: Unauthorized conversions or modifications to the device are not allowed for safety reasons.

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WARNING: Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see IEC 60601-1-1 or para 16 of the 3rd edition of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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6.2.2 Contraindications

No contraindications relating directly to the products used for CO₂ insufflation are currently known.

The use of an insufflator is contraindicated if, in the opinion of a qualified physician, such an application would pose risks to the patient, e.g. due to the patient's general condition.

The device is contraindicated for hysteroscopic applications.

CO₂ absorption

CO₂ absorption occurs during insufflation. The body absorbs a portion of the CO₂ used for the insufflation. A CO₂ concentration in the blood or in the respiratory system which is too high can lead to the death of the patient in isolated cases. Observe vitals throughout the entire procedure and ensure sufficient patient respiration. High pressure or high gas flow encourage CO₂ absorption.

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Changes to circulatory, metabolic and respiratory functions

The following symptoms can occur during CO₂ insufflation:

- Hypercapnia
- Acidosis
- Decrease in cardiac output
- Decrease in venous return

Embolism

Incorrect positioning of the insufflation instrument or a high intraabdominal pressure can lead to a CO₂ embolism. Avoid high pressures and verify the correct positioning of the insufflation instrument.

Hypothermia

The CO₂ insufflation can lead to a decrease in body temperature. Monitor the body temperature throughout the intervention.

Idiosyncratic reactions

The risk of a metabolic imbalance due to increased CO₂ absorption is increased for patients with sickle cell anemia or pulmonary insufficiency.

Dehydration

For longer operations with high gas levels, insufflation can lead to a drying out of the tissues or cause tissue damage. Unnecessary leaks must be avoided for this reason.

The decisive factor for the incidence of these complications is the length of the pneumoperitoneum, the level of the intraabdominal pressure and the CO₂ absorption associated with it. Increased pressure and flow values must thus be avoided. Pressures in excess of 15 mmHg are only required in isolated cases.

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6.3 User qualification

The ENDOFLATOR® 50 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

6.4 User profile of physician and assistants

Recognized medical skills in the field of application (specialist physician, qualified medical staff).

Adequate powers of comprehension to rationally assess the situation in hand.

Familiarity with the use of technical apparatus
Suitable command of one of the languages used by the device and in the instruction manual.

Must have been trained thoroughly at least once in the use of the device

No physical impairments which limit or prevent perception of acoustic or visual alarm signals.

6.5 Intended conditions of use

6.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

6.5.2 Other intended conditions

- Frequency of use: one or more times a day
- Length of use: from a few minutes to several hours a day.
- Place of installation: positioning on a level, vibration-free surface
- Mobility: can be moved if positioned on a cart
- Combination: can be used on the patient at the same time as other equipment required for the operations

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- Control: can be controlled via the KARL STORZ SCB.
- The device is placed outside the sterile area.

6.5.3 User position

- To operate the device, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$, at a distance of approx. 30 – 70 cm to the front panel.
- For observation of the actual values during the application, a visual distance from the device of 2 m is assumed.

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6.6 Safety precautions at the installation site

The device may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

WARNING: DO NOT use the device in the presence of flammable anesthetics.

WARNING: The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

WARNING: The device must not be used in oxygenated environments.

This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with the applicable national regulations.

6.7 Safety precautions when operating the device

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

This device has been tested and found to comply with the EMC limits of the Medical Device Directive 93/42/EEC (EN 55011 Class B and IEC 60601-1-2).

Please observe the Electromagnetic Compatibility information in section 11.

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6.8 Safety features

The ENDOFLATOR® 50 has the following safety features:

6.8.1 Self-test

A self-check procedure is performed each time the device is switched on (for this the device must be supplied with gas). If the test is negative, a corresponding error message is shown and the device can no longer be used as normal.

The following components are tested:

- Leakproofness of the system (valves)
- Checking pressure measurement
- Checking flow measurement
- Functionality of the patient tube heating control

6.8.2 Monitoring during operation

- **Overpressure:** If, during insufflation, the intracorporeal pressure rises above the set setpoint value by a limit dependent on the setpoint value for more than 3 seconds, an acoustic information signal sounds (see chapter 9.1.2) and the gas supply is stopped. After a time period which can be set (5...99 seconds), the overpressure is actively reduced via a blow-off valve.
- **Mode-specific pressure monitoring**
- **Checking gas supply**
If the system recognizes a lack of CO₂ gas, an optical signal is displayed.
- **Checking gas preheating**
The internal gas preheating is monitored by a second circuit, independently of the control. If a fault arises the heating is switched off.

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- **Checking the gas temperature to the patient**
The gas temperature and the energy output are recorded in the device. If a system error occurs, the heating will be switched off. The color of the symbol for the patient tube heating changes from white to red. Insufflation with (cold) CO₂ is, however, still possible.

- **Passive safety measures**
 - The pressure and flow measurement systems have a redundant stand-by.
 - Both pressure control stages and the flow measurement are equipped with safety pressure release valves.
 - The function of the mechanical high pressure regulator is monitored.

i NOTE: If the monitoring system recognizes an error which renders safe functioning seemingly impossible, the device is transferred into a "safe state". From this state, the device can only be reset by turning off the power supply.

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7 Installation and operating instructions

7.1 Unpacking the equipment

Carefully remove the ENDOFLATOR® 50 and the accessories from their packaging. Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately.

If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

7.2 Basic equipment KARL STORZ Set no.: UI500S1

- 1 ENDOFLATOR® 50 UI500
- 1 Power cord 400A
- max 5 Insufflation tube with gas filter for single use, sterile
- 1 Universal wrench 20 4000 30
- 1 SCB connecting cable 20 0901 70, length 100 cm
- 1 Instruction manual
- 1 CO₂ supply tube

7.3 Installing and connecting the device

 **NOTE:** The ENDOFLATOR® 50 including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

 **WARNING:** Position the device higher than the patient.

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Installation and operating instructions

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1. Place the device on a flat surface.
! WARNING: Ensure that there is sufficient air circulation for the device.
! WARNING: The device must be set up in such a way that the power cord can be simply pulled out at any time.

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- 7.3.1 Connecting the ground line**
! WARNING: The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.
2. Connect the potential equalization cable to the potential equalization connector.
3. Where necessary, have grounding performed by qualified personnel.

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- 7.3.2 Connecting the power cord**
! WARNING: The device may only be operated at the voltage stated on the rating plate.
4. Push the power cord into the power socket as far as it will go.
! WARNING: Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.
! WARNING: The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.
! WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

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7.3.3 KARL STORZ SCB


NOTE: To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets. Connect the other end of the cord to the KARL STORZ SCB (KARL STORZ Communication Bus) control unit or other SCB devices (see KARL STORZ SCB control NEO System Instruction Manual).

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7.3.4 Affixing the CO₂ bottle to the device (optional)

Attach the optionally available holder to the rear of the device. Use the enclosed screws of length 6 mm.

Insert the CO₂ bottle vertically into the holder and close the locking clip.

WARNING: Only sterile CO₂ gas may be used.

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
7.3.5 Connecting the CO₂ bottle

Connect the gas outlet of the CO₂ bottle to the device's gas connection using the high-pressure tube supplied.

WARNING: CO₂ bottles connected up to the device must be secured to prevent them from falling over.

WARNING: Follow the safety instructions when changing the CO₂ bottle.

WARNING: The CO₂ bottle must be placed vertically while operating the device, otherwise the functioning of the device cannot be guaranteed.

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7.3.6 CO₂ bottles equipped with German-standard or ISO connection

Place the high-pressure tube (e.g. 20 4002 22; see the section 'Accessories') onto the outlet of the CO₂ bottle and screw tight using the universal wrench provided.

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7.3.7 CO₂ bottles equipped with PIN-Index connection


Place the connecting piece of the high-pressure tube (e.g., 20 4000 22; see the section 'Accessories') onto the CO₂ bottle and screw this tight.

Use the wrench provided (see the illustration) to open the CO₂ bottle.

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
7.3.8 Opening the valve of the CO₂ bottle, if necessary.

Slowly open the valve stopcock of the CO₂ bottle by turning it approx. one-half turn counterclockwise.

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7.3.9 Connecting to the central gas supply

Connect the low-pressure tube (e.g., UI001; see the section 'Accessories') to the device's gas connection and connect the other end directly to the wall outlet of the central gas supply

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7.4 Commissioning

NOTE: The gas supply must be connected when the device is switched on, otherwise an error message is emitted.

7.4.1 Commissioning for the first time

The ENDOFLATOR® 50 is delivered as follows:

- (Screen) Language: English
- Start screen: Working screen
- High Flow mode (pressure in mmHg)

NOTE: Changing/adjustment of the settings for language, start screen, audio and the gas source are described on the following pages.

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1. Press the power switch ①.
Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo. After this, a white screen appears.
The device is now running a self-test. Once the test has been successfully completed, an acoustic availability signal sounds and a touch screen appears.

2. Confirm the message.

If no gas is connected or if the input pressure is too low, the warning 'Input pressure too low' appears.
Now check whether the bottle is open.

- ① **NOTE:** Normal use is only possible after successful conclusion of the tests.

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If the self-test fails, the screen **Safe State** appears with an information message and an acoustic information signal sounds.

NOTE: If restarting (turning the device off and on again) fails, verify that the device is properly connected and the preparations have been made according to the instructions.

Further information on troubleshooting can be found in the chapter 'Technical description - Troubleshooting'.

The switch-on message is displayed after the **Working screen**.

NOTE: In routine operation, insufflation can be started, the device settings changed or the procedure information called up from here.

3. Press the button ① (wrench) to change to the settings mode.

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

Changing the language and other presets

NOTE: The screen language is preset as English.

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The **Device Settings** screen appears with the list of available user settings.

- Language
- Start screen
- Audio
- Gas Source

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
5. Select your preferred language under **Language**. Confirm the selection with 15 or reject the changes with 14.

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6. Under **Start screen**, select the desired option.

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
7. Select the desired volume for information signals and key tones (from 0 ... 4) under **Audio**.

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8. Select the desired gas supply source under **Source of gas supply**.

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7.4.2 Normal commissioning


After successfully starting (see section 'Commissioning for the first time'), the **Working screen** or **Procedure list** screen appears depending on the setting.

NOTE: In routine operation, insufflation can be started, device settings changed or the procedure information called up from the **Working screen**.

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If the **Procedure list** working screen is not set as standard (see section 'Commissioning for the first time'), the **Procedure list** screen will appear after pressing the button ②.

For every unit 30 user procedures can be saved.

NOTE: The user procedure list is arranged according to when the procedures were last used. The most recently used procedure is shown at the top of the list.

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1. Attach the sterile CO₂ gas filter to the insufflation connection ③.

WARNING: You must insert a sterile CO₂ filter between the insufflation connection and the insufflation tube to prevent the risk of contamination of the device due to the reverse flow of gas or body fluid. The CO₂ filter must be replaced after each use.

2. Attach the insufflation tube to the CO₂ gas filter and connect the other end of the insufflation tube to the VERESS needle.

NOTE: In the case of the disposable insufflation tube set, the filter is firmly connected to the tube.

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7.4.3 Heating up the CO₂ gas

The ENDOFLATOR® 50 enables the CO₂ gas to be preheated to body temperature. The risk of hypothermia during the insufflation of sizeable quantities of CO₂ can thus be reduced.

Use the heatable insufflation tube, Art. no. 031210-10; pack of 10.

Connecting the heatable patient tube

Connect the connecting cable of the heatable patient tube to the socket ④ on the ENDOFLATOR® 50.

The heating function is activated by plugging in the connecting cable.

NOTE: The patient tube heating symbol ⑦ is displayed in the **Working screen**.

The CO₂ gas flowing through is heated to a temperature of approx. 37 °C. The unit adapts to various flow rates automatically.

NOTE: If the connecting cable is disconnected, the ENDOFLATOR® 50 can be operated without heating.

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7.4.4 Test for proper functioning

WARNING: Test the device prior to each surgical procedure to ensure that it functions correctly. The device should not be used if any damage is evident.

Pressure stability and leakproofness

1. Set the mode to 'Sensitive'.
2. Use 15 mmHg and 2.0 l/min as the setpoint values.
3. Check the connection of the patient tube to the device.
4. Close the stopcock on the instrument.
5. Start insufflation.

■ Pressure then builds up in the tube which increases above the setpoint value and which, after a few seconds, is released via the blow off valve. This is repeated cyclically.

6. Observe the instrument's conduct over a 30 second period.

NOTE: The flow display may only show values greater than zero for short periods. If the set pressure is only reached slowly even when the instrument's stopcock is closed or if you note stationary flow rates of 0.2 l/min or more, there is a leak in the system filter – patient tube – instrument.

7. In this case, check the connections and, if necessary, change the tube set.

WARNING: Do not work with a leaky device. There is a risk that the intraabdominal pressure will increase in an uncontrollable manner.

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7.4.5 Preparing insufflation

You have two possibilities for entering the desired operating parameters for insufflation mode*, pressure and flow.

a) Manual selection

1. In the **Working screen**, click on the High-Flow or Sensitive button to select the mode.
2. Set the setpoint values for pressure and flow using the $\pm$ buttons.

b) Selecting a saved procedure

1. Click on the button ② in the **Working screen**. ➤ The **Procedure list** screen appears.
2. Select the desired procedure by clicking (e.g. ③). (By clicking on the button ④, the device jumps back to the **Working screen** without any changes being made). ➤ The **Working screen** appears.

*The ENDOFLATOR® 50 is equipped with two different insufflation modes:

- High Flow ('man' symbol)
- Sensitive ('baby' symbol)

NOTE: To compensate for the considerable loss of gas occurring in complex operations, the High-Flow mode has been fitted with a flow rate of up to 50 l/min. For particularly sensitive areas such as pediatric applications, the device features a 'Sensitive' mode which is characterized by low insufflation pressure as well as special safety limits in the pressure (up to max. 15 mmHg) and flow ranges (up to max. 15 l/min).

WARNING: Always check the mode set. For safety reasons, use the 'Sensitive' mode if the required cavity pressure is smaller or equal to 15 mmHg and a maximum flow of 15 l/min is sufficient.

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7.4.6 Performing CO₂ insufflation

 **WARNING:** Before starting insufflation rinse the entire system with 1 l CO₂.

 **WARNING:** Always use the lowest possible pressure and flow for insufflation.

3. Reset the volume display using the '0' button.
4. Insert the VERESS needle and open the locking lever. Press the START/STOP button to initiate the insufflation process.

 **NOTE:** In order to use the maximum device flow in the High Flow mode, HiCap® instruments with a suitable connection must be used. The LUER-lock connection fitted to the end of the insufflation tube must be removed in this case.

 **NOTE:** In the case of disposable tubing sets, cut the tube off behind the LUER-lock.

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
Displays

The insufflation pressure, flow and volume of the gas used are continuously displayed in the *Working screen*.

If the patient pressure rises above the preselected value, an audible warning signal will be sounded. After the set delay (see *Settings* – Administration screen), the pressure is released via the overpressure valve.

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7.4.7 Decommissioning

Close the valve on the CO₂ bottle on ending the intervention (when using a high-pressure gas supply) and remove the insufflation tube (incl. filter) used from the device.

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7.5 Menu explanations

7.5.1 Procedure list

Presaved values for mode, pressure and flow can be called up directly via the procedure list.

Press the button ② to call up the **Procedure list** screen (text is language-dependent) (if not set as the start screen) (see Fig. B).

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Up to 30 procedures can be saved. Four procedures are displayed on each page. If there are more than four procedures, the buttons ⑥ are displayed for turning the pages.

Saved procedures can be deleted or modified, or saved as a new procedure.

NOTE: By clicking on the button ④, the device jumps back to the Working screen without any changes being made.

Creating a new procedure

1. Press the button ① (see Fig. A) to call up the **Procedure list** screen if necessary (also see Section 7.4.2).

The **Procedure list** screen shows any saved procedures (see Fig. B).

2. Press the button ③ (Create a new procedure).

NOTE: If 30 procedures have already been created, the following information dialog appears when you press the button ③: "A maximum of 30 procedures can be saved in the procedure list. All memory locations are currently occupied. Before a new procedure can be saved, an existing procedure must be deleted from the list."

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The *Set procedure values* screen with the current values for mode, pressure and flow appears (see Fig. C).

3. Press the button 15 to confirm the set values. Alternatively, the desired values for mode, pressure and flow can be changed with the corresponding buttons 8, 9, 10, 11, 12 or 13.

NOTE: Press the button 4 or 14 to reject the changes and return to the Working screen or Procedure list screen.

The *Assign procedure name* screen appears (see Fig. D).

4. Enter a name for the procedure.

NOTE: The on-screen keyboard can type lower-case or capital letters, numbers and special characters depending on the status of the buttons 16 and 17. By pressing the arrow buttons 18 or 19, the cursor moves one character to the left or right. Use the button 20 to delete the character to the left of the cursor.

5. Press the button 15 to confirm the input. When a new procedure is created, the device returns to the Procedure list screen (see Fig. E).

Editing an existing procedure

1. Press the button 8 in the *Procedure list* screen (see Fig. E).

The *Edit procedure list* screen appears (see Fig. F).

2. If necessary, select the procedure using the buttons 6 (procedure is marked orange).

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
- The *Set procedure values* screen with the current values appears (see Fig. C).
- 3. Set the desired values for mode, pressure and flow using the corresponding buttons (8, 9, 10, 11, 12 or 13).
- 4. Confirm the set values by pressing the button (15) or reject the changed settings using (4) or (14).

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
- The *Assign procedure name* screen appears (see Fig. D).
- 5. Enter a new name for the procedure or overwrite the old one.
- After pressing the button (15) to confirm, or canceling using the button (14), the *Edit procedure list* screen appears again (see Fig. F).

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
- ### Deleting an existing procedure
1. Call up the *Edit procedure list* screen in the same way as described for editing an existing procedure (see Fig. F, page 34).
 2. Select the procedure using the buttons (6) if necessary.
 3. Press the trash can button next to the procedure to be deleted.

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
- The *Delete procedure* screen appears (see Fig. G).
- 4. Confirm the deletion of the procedure by pressing the button (22) or cancel the process using the button (23).
- The *Edit procedure list* screen appears again (see Fig. F).

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
7.6 Settings

- Device Settings
- Service
- System Log
- Administration
- Device information

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7.6.1 Device settings

The Device Settings are described in the section on commissioning for the first time.

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7.6.2 Service


The **Service** screen is protected by a password.

NOTE: The service programs are only available to service personnel. The settings are described in the Service Manual.

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7.6.3 System Log

The System Log contains the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

The most recent entry is allocated the number 000.

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
7.6.4 Administration

- Export system log
- Unit pressure
- Date and Time
- Overpressure valve time delay
- Change administration password


NOTE: The Administration screen is protected by a password. On delivery the password is '2132'.

It is the responsibility of the user to change this if necessary.

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- Export system log
- For service purposes, data can be exported via the service interface ⑧ using the USB stick.

1. Open the cover on the rear of the device with a Torx driver (size T10).
2. Insert a USB stick into the USB interface.
3. Select Export system log on the **Administration** screen.


Exporting of the system log starts.

NOTE: The current status is shown via a progress bar. This progress bar shows that the export is underway, however it does not show how far it has advanced.

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
4. Screw the cover back on once export is complete.

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- Unit pressure


Select the desired unit (mmHg or hPa) under **Unit pressure**.

Confirm the selection with 15 or reject the changes with 14.

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
- Date and Time
The setting of date and time is self-explanatory.
Confirm the selection with (15) or reject the changes
with (14).

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
- Overpressure valve time delay
Select the desired time (5...99s) under
Overpressure valve time delay.
Confirm the selection with (15) or reject the changes
with (14).

NOTE: If there is overpressure for longer
than the set time, this is reduced. If the
overpressure cannot be reduced, the device
emits a warning.

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- 7.6.5 Device information**
The following elements are displayed:
- Device serial number
 - Operating hours
 - Software versions
 - Software license information

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


8 Maintenance

8.1 Fuse replacement

1. Turn the device off and disconnect from the power supply.

2. Remove the screw bits on the line fuse holder ⑫ using a screwdriver or another suitable tool.

! WARNING: Only use fuses of the correct rating.

3. Insert new fuses.

	100...240 V
Line fuse	2 x T2.5 AH250V

4. Now put the line fuse holder ⑫ on again.
5. Restore the power supply.
6. Switch the device on and ensure that the start screen appears.

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8.2 Reprocessing

- UI 500 – ENDOFLATOR® 50
- Heatable insufflation tubing set with gas filter, single use, sterile

General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: The reusable insufflation tube is not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect the reusable insufflation tube for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the reusable insufflation tube before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: Risk of infection/product damage: KARL STORZ has validated the procedures given in the instruction manual. Only use these procedures. The reprocessing method must be validated by the user on site.

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Maintenance

Reprocessing

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Risk of damaging the medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

WARNING: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

8.2.1 Reprocessing of the ENDOFLATOR® 50

WARNING: Always disconnect the device from the mains before cleaning!

WARNING: Liquid must not be allowed to enter the housing under any circumstances.

Manual wipe-down disinfection

The design of the device and supply tube allows only for wipe-down disinfection. The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility must be observed.

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8.2.2 Reprocessing of heatable/non-heatable insufflation tubing set with gas filter, single use, sterile

WARNING: Disposable (heatable) insufflation tubing sets with gas filter, locks and hydrophobic bacteria filters must not be reprocessed.

Sterile disposable products are not approved for reprocessing.

8.2.3 Reprocessing of reusable insufflation tube

WARNING: The reusable insufflation tube must be cleaned, disinfected and sterilized prior to use and subsequent reuse using validated reprocessing procedures.

Disassembly of the connectors

As far as possible, separate the reusable insufflation tube into its individual components prior to cleaning and disinfection.

8.2.3.1 Preparation for cleaning and disinfection

The reusable insufflation tube must be kept moist for transportation to the reprocessing station in order to prevent contaminations from drying on or reduce this as far as possible.

**8.2.3.2 Manual precleaning
Immersion in cold water for 5 minutes**

To dissolve dried organic residues, immerse the reusable insufflation tube in cold water (10-20 °C) for 5 minutes.

Brushing the surfaces

Any visible contamination or heavy soiling must be completely removed from all external surfaces by cleaning under cold running water with the aid of a brush (Art. no. 27652) or a sponge.

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Brushing the lumina/hollow spaces

In order to remove all visible contamination, the lumina of the connectors and the tube ends of the reusable insufflation tube must be precleaned under cold running water with the aid of a suitable brush (Art. no. 27650A/27650B).

Rinsing out with a pressure water gun

The lumina of the reusable insufflation tube must be rinsed out using a pressure water gun and a suitable rinsing attachment. For rinsing, use pulses of at least 3 x 10 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) with the special attachment (Art. no. 27660F) is ideal for this purpose.

8.2.3.3 Manual cleaning

The reusable insufflation tube must be completely immersed in a cleaning solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically using a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

8.2.3.4 Manual disinfection

The reusable insufflation tube must be completely immersed in a disinfectant solution. Through disassembling the instrument and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. Following the exposure time, the reusable insufflation tube and connectors must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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8.2.3.5 Machine cleaning and disinfection

Prior to machine cleaning and disinfection, carry out the steps of manual precleaning. The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method. The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

Connecting

In order to guarantee effective machine cleaning and disinfection, connect the LUER-Lock to the washer and disinfectant using the LUER connector provided. Place the reusable insufflation tube on the rinsing lance in order to ensure that it is rinsed out completely. In addition, put the plastic connector with the tube adapter side to the top on the rinsing lance.

NOTE: If necessary, the reusable insufflation tube must be dried off afterwards by hand.

8.2.3.6 Assembly, inspection and care

WARNING: A visual and functional check of the reusable insufflation tube must be carried out!

WARNING: Check the reusable insufflation tube for leaks after cleaning and disinfection. By doing this it is possible to detect a leak early, through which gas could escape.

NOTE: Discolorations impede visual inspections for possible contamination.

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
The cleaned and disinfected reusable insufflation tube must be visually inspected to ensure that it is clean, complete, dry and has not been damaged:

- If residues or contamination are still present, the reusable insufflation tube must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged medical devices must be withdrawn from use.
- There should be no holes, cuts, dents or scratches on the tube surface. If the reusable insufflation tube is damaged as described above, it must be replaced. If the reusable insufflation tube is discolored as a result of the reprocessing process, it must be removed immediately from the sterile supply cycle.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

The integrity of the reusable insufflation tube can be tested in the following way:

Attach a syringe to one end of the tube and kink the other end. Place the insufflation tube in completely demineralized water or microbiologically pure, sterile water. Force air into the insufflation tube using the syringe. If any bubbles appear in the water while forcing air through the reusable insufflation tube, the tube must not be reused.

ⓘ NOTE: The leakage tester 13242 XL may also be used to check for leaks.

To check for leakage use leak tester 13242 XL, for example. To do so, connect one tube end to the leakage tester and block off the other end, for example, with your finger. If there is a leak in the tube set, a pressure drop will be observed on the pressure gauge of the leakage tester.

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8.2.3.8 Sterilization

WARNING: Prior to the sterilization of the reusable insufflation tube, particularly in cases of manual reprocessing, it must be ensured that all chemical residue has been rinsed off or removed.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method is to be selected taking into account the specific national requirements and in consultation with the device manufacturers.

The following sterilization method has been validated and approved by KARL STORZ for this medical device:

Steam sterilization in the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 – 137 °C, 0,20 MPa with a minimum exposure time of 4 minutes and maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

8.2.3.9 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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8.3 Maintenance and safety check

8.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.3.2 Safety check/repeat inspection according to IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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Test for proper functioning

A functional inspection must be carried out in accordance with the instruction manual (see section 7.4.4).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/ repeat inspection, please refer to the latest issue of the relevant Service Manual.

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8.4 Servicing and repair

WARNING: Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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8.5 Disposal

This device has been marked in accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE).

WARNING: At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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8.6 Repair program

During the repair period, you will generally be loaned a device, which must then be returned to KARL STORZ as soon as the repaired device is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

8.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

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8.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the device has been used in accordance with its operating instructions at all times.

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8.9 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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9 Technical description

9.1 Signals

All of the signals created, with the exception of those which display the cause of an inoperable device, are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by an alarm signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized last will be signaled, i.e., the previous signal is overwritten.

9.1.1 Visual signals

Signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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9.1.2 Acoustic signals


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

i NOTE: An acoustic signal is always issued as a complete sequence (regardless of how long the signal conditions apply).

The acoustic signal can be switched on or off in the 'Device settings' menu.

The acoustic signal can, however, only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

• Information signal (general)

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15 s, but otherwise it is just played once.

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
• **Availability signal**

On completion of the successful self-test after switching on the device, an availability signal with the following characteristics sounds:

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

• **Low gas supply**

When using a high-pressure supply, the drop in supply pressure shows that there is no more liquid CO₂ in the reservoir.

In addition to the visual messages described in the 'Visual signals' section, when the pressure levels 40 bar – 30 bar – 20 bar are fallen short of, a double gong is produced as an acoustic signal.

This is a sampled sequence. It lasts for 1.5 seconds.

If a minimum supply pressure is fallen short of, the information signal is output as specified both with a high-pressure supply and low pressure supply.

• **Key click**

Clicking on a key on the touch screen is indicated by a short beep. The volume of this beep can be set in the 'Device Settings' menu independently of all other information signals. The key click can also be permanently deactivated.

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Device defective.

Measure:

Switch unit off and back on. If the warning message still appears after restarting, contact Service. Prior to this, open the System Log under Settings and note the ID, Word1, Word2 and information text of the last 5 entries.

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No gas flow.

Measure:

Observe the warnings (device temperature exceeded, temperature of high pressure regulator too low, gas supply not sufficient).

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No pressure build-up.

Possible causes:

1. Excessively large leak.
2. Excessively high instrument resistance.
3. Leaky insufflation system.

Measure:

1. Increase the flow setpoint value
(for information: if the gas outflow from the cavity with the desired intraabdominal pressure is larger than the set flow setpoint value, the insufflation regulation goes to flow limitation and the cavity pressure setpoint value is not achieved).
2. Check whether it is possible to change the instrument (for information: there are situations in which the insufflation pressure provided is not sufficient to build up the desired cavity pressure due to a high level of instrument resistance and simultaneous gas outflow. Depending on the instrument and leak conditions, the current values for pressure and flow are below the setpoint values).
3. a) Check the connections of the accessories and the leakproofness of the trocar. If the test for proper functioning is passed, the likelihood of a loose connection is very high.
b) Change the tube/tube set and perform a test for proper functioning.
c) If necessary, change the filter and repeat the test for proper functioning.
d) Contact service if the leak remains.

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
Error description:

After switching on the device without the heatable patient tube plugged in, the patient heating symbol lights up red and a message is displayed.

Possible causes:

The device is defective and the patient heating is permanently switched off.

Action:

Work can be continued without the patient tube heating.

Contact Service.

Error description:

After plugging in the heatable patient tube, the patient heating symbol lights up red and the message 'Defective patient tube heating' is displayed.

Possible causes:

The patient heating is defective.

Action:

Replace the heatable patient tube and put the patient heating into operation by switching the device off and on again.

Error description:

After plugging in a heatable patient tube, the patient heating symbol lights up white and then changes its color from white to red, and the message 'Heated tube defective, please replace' is displayed.

Possible causes:

The patient heating is defective.

Action:

Replace the heatable patient tube. If the patient heating symbol disappears after the heatable patient tube is unplugged, a new heatable patient tube can be used immediately.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

- Type of IP protection: IP21

According to IEC 60601-1-2, FCC-B:

Observe the information on electromagnetic compatibility in chapter 11 (pages 62-75).

**9.3.2 Directive compliance
(for UI500)**

According to the Medical Device Directive (MDD):

This medical product belongs to Class IIb.

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

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9.4 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as reparable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

① NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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**Spare parts,
recommended accessories**

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**10 Spare parts,
recommended accessories**

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Article	Cat. no.
High-pressure tube, American connector/German connector, length 55 cm	20 4000 21
Same, length 102 cm	20 4000 27
High-pressure tube, American connector/ISO connector, length 102 cm	20 4002 22
CO ₂ bottle, 1 l, empty, with German connector	26 4000 90
CO ₂ bottle, 1 l, empty, with PIN-Index connector	26 4000 91
Low-pressure tube, for the central CO ₂ gas supply, length 150 cm	U1001
Same, length 300 cm	U1002
Same, length 600 cm	U1003
HD-Inline gas filter	20 4000 32
Line fuses T2.0AH 250V, pack of 10	1858390
Power cord (grounded)	400 A

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10.2 Accessories

Article	Cat. no.
Insufflation tube, sterilizable, length 2.5 m, ID 9 mm	UI004
Gas filter, disposable, sterile, pack of 25	031122-25
Insufflation tube set with gas filter for single use, sterile, pack of 10	031200-10
Heatable insufflation tubing set with gas filter, single use, sterile, pack of 10	031210-10
VERESS pneumoperitoneum needle with spring-loaded blunt stylet, LUER lock, length 10 cm	26120 J
Same, length 12 cm	26120 JL
Same, length 7 cm	26120 JK
Same, length 15 cm for creating pneumoperitonea by means of DOUGLAS puncture with adipose patients	26120 JF
Bottle holder for the KARL STORZ CO ₂ bottle	UI005
Leakage tester	13242 XL

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Electromagnetic Compatibility (EMC) information

11 Electromagnetic Compatibility (EMC) information

WARNING: Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The KARL STORZ ENDOFLATOR® 50 Model UI500 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

The limits used are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The KARL STORZ ENDOFLATOR® 50 Model UI500 is a Group 1 device (as per CISPR 11).

Group 1 contains all of the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the individual devices
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic Compatibility (EMC) information

WARNING: The device should not be used adjacent to or stacked with other equipment. If use adjacent to or stacked with other devices is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

WARNING: The use of accessories / transducers and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the device. The accessories / transducers and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories / transducers and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Cables which have been shown to comply with IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power Cord	No	3	No	Power supply

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions

The ENDOFLATOR® 50 model UI500 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 50 model UI500 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ENDOFLATOR® 50 model UI500 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The ENDOFLATOR® 50 model UI500 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2
Guidance and manufacturer's declaration – electromagnetic immunity

The ENDOFLATOR® 50 model UI500 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 50 model UI500 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The ENDOFLATOR® 50 model UI500 is intended for use in the electromagnetic environment specified below. The customer or user of the ENDOFLATOR® 50 model UI500 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	$3 V_{rms}$ 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	$3 V_{rms}$ 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the ENDOFLATOR® 50 model UI500, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 6
**Recommended separation distances between portable and mobile RF communications equipment and the
ENDOFLATOR® 50 model UI500**

The ENDOFLATOR® 50 model UI500 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseia.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latina-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseia.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 158
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14,
12724 Skärholmén, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros, 00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlyst 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Suite 2503-5, 25th Floor, AXA Tower, Landmark East,
100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

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KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016
People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseca.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América do Sul Ltda.
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olssensvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopie Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopie France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell' Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adna Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) - Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 - Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F. Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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Karl Storz GmbH & Co. KG