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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to aid you with proper connection and operation of the 3D TIPCAM*1 in conjunction with the instruction manual for the 3D CCU. All required details and all actions on your part are clearly presented and explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment. Insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the 3D CCU instruction manual for the individual functions of the CCU that can be used with the 3D TIPCAM*1 camera.

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
- ① Inserter (sheath) with CCD sensor
- ② Handle
- ③ Control button 1: Scroll up in menu / programmed functions
- ④ Control button 2: Selection / programmed functions
- ⑤ Control button 3: Scroll down in menu / programmed functions
- ⑥ Connection for camera head cable
- ⑦ Connection for fiberoptic light cable
- ⑧ Camera head cable connector to CCU
- ⑨ Fiberoptic light cable (495TIP)
- ⑩ Light cable connector to KARL STORZ light source
- ⑪ Camera Cable

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1 General information	2
1.1 Device description	2
2 Safety instructions	3
2.1 Warnings and cautions	3
2.2 Definitions	3
2.3 Intended use	10
2.4 User qualification	11
2.5 Safety precautions when operating the instrument	11
3 Installation and operating instructions	12
3.1 Unpacking the equipment	12
3.2 Basic equipment	12
3.3 Compatible camera control unit (CCU)	12
3.4 Before every use	13
3.5 Start-up of the 3D TIPCAM*1	13
3.6 Connecting the camera head cable to the 3D TIPCAM*1	13
3.7 Connecting the light cable to the 3D TIPCAM*1	13
3.8 Connecting the 3D TIPCAM*1 to the light source	14
3.9 Connecting the 3D TIPCAM*1 to the CCU	14
3.10 Switching on the CCU, light source and monitor	14
3.11 Optimum light efficiency	15
4 Maintenance	16
4.1 General care	17
4.2 Care and handling	17
4.3 (Pre-) cleaning	21
4.4 Manual cleaning and disinfection of the 3D TIPCAM*1	21
4.5 Machine preparation of the 3D TIPCAM*1	24
4.6 Sterilization instructions for the 3D TIPCAM*1 (autoclavable)	27
4.7 Steam sterilization	27
4.8 Gas sterilization with ethylene oxide (EO)	28
4.9 Sterilization with the STERRAD* 100NX™ system	29
4.10 Sterilization with the STERRAD* 100S system	31
4.11 Sterilization with the STERIS System1™	33
4.12 Sterilization with the AMSCO V-PRO™ 1 (STERIS®)	35
4.13 Use of sterile covers	37
4.14 References for cleaning and sterilization	37
5 Maintenance	38
5.1 Servicing and repair	38
5.2 Disposal	39
5.3 Repair program	39
5.4 Limitation of liability	40
5.5 Manufacturer's warranty	40
6 Technical description	41
6.1 Technical data 3D TIPCAM*1	41
6.2 Directive compliance (for 3D TIPCAM*1)	41
7 Recommended accessories	42

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1 General information

1.1 Device description

The 3D TIPCAM®1 is a high-resolution video laparoscope with two distally integrated video chips.

Since the 3D TIPCAM®1 does not require an optical imaging system, it is very sturdy.

The 3D TIPCAM®1 is compatible with the 3D CCU.

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2 Safety instructions

2.1 Warnings and cautions

Please read this manual and follow the instructions carefully. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device.

2.2 Definitions

 **WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

 **CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

 **NOTE:** A Note indicates special information about operating the unit.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

WARNING: Also pay particular attention to the instruction manual for the KARL STORZ 3D CCU no. 96206204D.

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WARNING: Risk of injury/material damage if the instruction manual/warnings are not observed: Read and observe the instruction manual and warnings for the products and the products used in combination.

WARNING: Risk of injury to the patient/user and/or damage to the products: Combinations of medical products are only assured to be safe if they are identified as such in the respective instruction manuals or the intended purpose and interface specifications of the products used in combination permit this.

WARNING: Risk of injury to the patient/user: Damaged or incomplete instruments/accessories may cause malfunctions. Before and after every use, always check the entire system (e.g., the 3D TIPCAM® 1, all cables, the 3D CCU, monitor and accessories used) to ensure that everything is configured properly, in perfect working order and nothing is missing/damaged. Replace incomplete/damaged instruments/accessories and do not continue to use them.

• Before every use, check that all the instruments required for the operation are functioning correctly.

• Be certain to check the cover glass of the 3D TIPCAM® 1. In addition, the camera cable and the light cable need to be inspected for damage/cracks.

• Do not use wet/moist cables.

• When attaching components to an instrument, the secure position of the components must be checked by gently pulling and/or turning them in both directions before use.

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WARNING: Check the instrument and accessories used in combination before every application for unintentional rough surfaces, sharp corners, burred edges, and projecting parts, which can endanger personal safety. Damaged, defective, bent, or incomplete instruments or accessories must not be used further under any circumstances.

WARNING: Risk of injury to the patient/user: If you apply or exert too much force on the instrument this may cause the instrument to break, bend, or malfunction, which, in turn, can lead to injury to patients and users. Do not place too much pressure on the instrument. Do not attempt to bend bent instruments back to their original position (tensile cracks lead to corrosion).

WARNING: When using with other equipment, e.g. shavers, make certain not to damage the 3D TIP-CAM® 1 mechanically. Only activate equipment under vision.

WARNING: Risk of injury: Only employ the instruments if you possess an appropriate medical qualification and are acquainted with the application.

WARNING: Risk of burns for patient/user! The high level of light energy means that the distal end, the light port, and adjacent components may become so hot that they can cause burns to the patient/user and/or to the surgical accessories (e.g., to sterile cover cloths)! If the distance between the light outlet and the tissue is too small this can result in tissue damage. Avoid tissue coming into direct contact with the distal end and with the light port of the light carrier and the light cable. Never place the 3D TIP-CAM® 1 and the light cable on the patient or in direct contact with the surgical accessories. Select the lowest possible lighting power of the cold light source that still allows optimal illumination of the operating field. Please avoid using the light for long periods at a high intensity.

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WARNING: Danger of glare! Never look into the open end of a connected light carrier or the 3D TIPCAM® 1. The radiant energy being emitted can cause irreversible damage to eyes.

WARNING: The 3D TIPCAM® 1 must be replaced if the image is cloudy or if no image or only parts thereof are visible.

WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the 3D TIPCAM® 1 is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the handle, which can cause cardiac pacemakers or other electronic devices to malfunction.

WARNING: The video system can malfunction during use. An adequate replacement system must be kept available.

WARNING: In the event of image failure immediately remove 3D TIPCAM® 1 from the patient.

WARNING: The use of HF probes, HF electrodes, HF instruments or other energy emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.

WARNING: Patient leakage currents
If various power-operated systems are used on the patient at the same time, patient leakage currents could add up. To reduce patient leakage current in such cases, combining Type CF applied parts is recommended.

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WARNING: Interaction with HF
If the 3D TiPCAM[®] 1 is used in combination with endoscopic accessories for HF equipment, the following must be observed:
Avoid direct contact of the insertion part with an activated HF electrode. This can lead to irreversible insulation defects in the instrument or the CCU, and thus to risk for the patient or user by electrical shock. Activate the electrode only under vision and with sufficient distance from the insertion parts on the endoscope.

WARNING: Live parts or instruments may not come into contact with 3D TiPCAM[®] 1.

WARNING: Unintentional grounding of the patient: The metal part of the camera connector is electrically connected to the sheath of the endoscope. During use, make certain that these metal parts do not come into contact with ground potential.

WARNING: HF probes, HF electrodes and other energy emitting probes may only be used under continuous visual endoscopic control.

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WARNING: Risk of electric shock: Do not attempt to open the 3D TIPCAM®1. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

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CAUTION: Periodically check the video system.

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Moisture which has penetrated the system due to damage can lead to faulty insulation paths within the video laparoscope. Rough surfaces or sharp edges can lead to patient injury.

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NOTE: Follow the applicable national ordinances and recycling plans regarding disposal or recycling of device components.

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Acquaint yourself with the operation and control of the instrument before initial use of the 3D TIPCAM®1.

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2.3 Intended use

The 3D TIPCAM®1 combines camera and endoscope and is used to generate and transmit video and control signals to the camera control unit in medical video endoscopy.

Use of the 3D TIPCAM®1 in fields other than those indicated above is not allowed for safety reasons.

The 3D TIPCAM®1 must only be used with equipment designated by KARL STORZ as suitable for the video laparoscope, or the safe and compatible usage of which has been documented.

Unauthorized conversions or modifications to the 3D TIPCAM®1 are prohibited for reasons of safety.

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2.4 User qualification

The 3D TIPCAM®1 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the 3D TIPCAM®1.

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2.5 Safety precautions when operating the instrument

It is the user's responsibility to make sure the instrument is safe and operates properly before using the instrument.

During treatment in conjunction with the 3D TIPCAM®1, the patient must be treated and observed with the standard medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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3 Installation and operating instructions

3.1 Unpacking the equipment

Carefully unpack the 3D TIPCAM[®]1 and remove the unit from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

WARNING: Risk of infection for patient/personnel: The 3D TIPCAM[®]1 is **not sterile** when delivered. The use of non-sterile instruments poses a risk of infection for patients and users. Prepare the instrument before initial use and before every subsequent use using validated procedures. See point 4 Maintenance.

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3.3 Kompatible Kamerakontrolleinheiten (CCU)

CCU
3D CCU 22204020-114

3.3 Compatible camera control unit (CCU)

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3.4 Before every use

Before every use, the instrument must be checked for damage, especially around the distal end (cover glass, light outlet, ...)

Before each use, the camera cable must be checked for damage, especially breaks and cuts.

Instruments with damaged cables must not be used further under any circumstances. Cables must be replaced (light cable to be replaced by the user or the Service Dept., and the camera head cable by the Service Dept.).

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3.5 Start-up of the 3D TIPCAM®1

NOTE: The 3D TIPCAM®1 including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

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3.6 Connecting the camera head cable to the 3D TIPCAM®1

The data cable is connected by plugging it into the cable gland provided. The red dot on the data cable (see arrow in photo) marks the top side of the data cable. The cable clicks and locks into place when connected correctly. Check that it is firmly connected!

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3.7 Connecting the light cable to the 3D TIPCAM®1

The light cable is mounted by fitting to the light port provided and turning clockwise. Check that it is firmly connected!

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! CAUTION: Please note that the circular symbol on the 3D TIPCAM®1 and the circular symbol on the light cable must correspond (see table at the top of the next page).

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Overview: Recommended combinations*

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* Special endoscopes may deviate slightly

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3.8 Connecting the 3D TIPCAM®1 to the light source

Insert the light cable into the socket on the light source until it clicks into place. Only hold the light cable by the handle piece. Never pull on the cable.

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3.9 Connecting the 3D TIPCAM®1 to the CCU

Insert the 3D TIPCAM®1 camera connector into the connection socket of the 3D CCU.

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3.10 Switch on the CCU, light source and monitor

NOTE: Read the instruction manuals for the Camera Control Unit (CCU), the light source, the light cable and the monitor.

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3.11 Optimum light efficiency

NOTE: The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light transmission is to be maintained. Therefore remove light cable (photo 1).
Following this, clean the light inlet and outlet surfaces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70% isopropyl alcohol. Please note that the glass surface inside the light port of the 3D TIPCAM1™ (see arrow in photo 2) and the light cable itself (see photo 3) must be cleaned in this way. Please pay attention to page 26:

CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

CAUTION: Please note that when attaching the light cable, the circular symbol on the 3D TIPCAM1™ and the circular symbol on the light cable must correspond (see table at the top of page 14).

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4 Maintenance

CAUTION: The manual „Cleaning, disinfection, care and sterilization of KARL STORZ instruments“ (96216003D) must be followed. The cleaning, disinfection, care and sterilization procedures are explained therein in detail. Only approved preparation procedures and cleaning agents may be used. The up-to-date list of KARL STORZ approved materials can be found on the internet at „www.karlstorz.com“ or in extracts in the manual referred to above.

WARNING: Risk of infection for patient/staff: The 3D TIPCAM®1 is **not sterile** when delivered. Clean, disinfect, and sterilize the instruments before initial use and before every subsequent use using validated procedures.

WARNING: Risk of infection for patient: Prior to sterilization inspect the instrument for visible soiling (blood, tissue remnants), which indicates that cleaning and disinfection have not been performed or have been performed incorrectly. Repeat cleaning and disinfection if necessary.

WARNING: Ensure that the instruments and surroundings are sterile!

WARNING: The 3D TIPCAM®1 must be dismantled for preparation. That means: The data cable and the light cable must be removed. To do so, pull the data cable out by the ribbing near the red dot and turn the light cable to the left.

WARNING: The product described here must be used in accordance with the medical rules and procedures of endoscopy recognized for this method.

WARNING: Before sterilization, the 3D TIPCAM®1 must be freed of all organic materials and cleaning residues. Sterilization is only effective on clean, disinfected surfaces.

WARNING: Any deviation from the recommended preparation parameters must be validated by the user.

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4.1 General care

WARNING: Always pull out power plug before cleaning!

Store the 3D TiPCAM[®]1 in the storage/sterilization box available as an optional extra. Do not coil the 3D TiPCAM[®]1 into tight loops or cause it to kink while storing. Do not store the 3D TiPCAM[®]1 in direct sunlight or excessive heat. When returning the 3D TiPCAM[®]1 for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

CAUTION: Chemical residue on the card-edge connector can cause camera failure.

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4.2 Care and handling

To maximize the life of your 3D TiPCAM[®]1, please follow the guidelines given below for proper handling and care of the camera head.

DO NOT PLUG THE CAMERA CONNECTING CORD INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. Where necessary, remove remaining moisture with a low lint, disposable cloth.

The connector of the data cable and of the fiberoptic light cable must also be completely dry before connecting!

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NEVER PICK UP THE 3D TIPCAM*1 BY ITS CABLE! Always pick up the set by the camera head first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE 3D TIPCAM*1 FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to device. When transporting, use the original shipping box to prevent damage to the device.

DO NOT STORE THE 3D TIPCAM*1 WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Always arrange the cable in loose coils of 6 inches (15 cm) or more in diameter.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful "kinks" can result. Always allow the cable to uncoil gently and loosely.

NOTE: The 3D TIPCAM*1 is autoclavable.

NOTE: We strongly recommend drying using sterile compressed air following manual preparation and, if required, after machine preparation.

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CAUTION: Inspect the camera connecting cord and fiberoptic light cable for breaks or cracks before and after disinfection or sterilization, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cable and cause damage to it. Where a 3D TIPCAM® 1 has a damaged fiberoptic light cable or damaged data cable, the relevant cable must be replaced.

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WARNING: The following procedures alone do not guarantee decontamination or sterility. Decontamination or sterility can only be determined by the staff when practicing recognized procedures for disinfection and sterilization.

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NOTE: The use of other disinfection/sterilization methods than those specified in the sections above may cause adverse effects on the materials of the camera and instrument coupling. Resulting damage will not be covered by the warranty.

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4.3 (Preliminary) cleaning

Immediately after use, clean the 3D TIPCAM®1 with camera connecting cord and fiberoptic light cable with a moistened, disposable wipe-down cloth (alcohol-free), in order to prevent contaminants from coagulating and drying out on the surface of the instruments.

Water quality requirements

Tap water can be used for preparing the cleaning and disinfectant solution; however, completely demineralized water is preferred.

Microbiologically pure/sterile water must be used to rinse (final rinsing) after exposure.

4.4 Manual cleaning and disinfection of the 3D TIPCAM®1

WARNING: When cleaning and disinfecting contaminated instruments, the local guidelines of the Employer's Liability Insurance Association and equivalent organizations are always to be observed.

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1. Completely immerse the 3D TIPCAM®1 (including cables) in a cleaning solution (according to the manufacturer's instructions). KARL STORZ advises you not to use any chemicals other than those approved for preparation.
2. Remove soiling from the 3D TIPCAM®1 using a disinfected soft brush (e.g., Art. no. 27562) or a disinfected sponge while taking care to avoid spattering.
3. Cleaning accessories which have been designed for the materials and sensitive surfaces are available from KARL STORZ (see catalog).
4. Rinse the 3D TIPCAM®1 (including cables) and inspect for residual dirt. Carefully check that the connectors of the camera connecting cord and fiberoptic light cable are clean.
5. Clean and disinfect any cleaning accessories (brushes, etc.) after manual preparation.
6. Immerse the 3D TIPCAM®1 (including cables) in a suitable basin with disinfectant (concentration and exposure time according to the chemical manufacturer's instructions), making sure that there are no air bubbles. Use plastic basins so as to avoid scratching the instruments and to prevent electrolytic corrosion.
7. Following exposure, rinse the instrument with microbiologically pure/completely demineralized water. Use unsoiled, sterile water for each rinse. Rinse the 3D TIPCAM®1 and cables thoroughly several times to remove any residues which might interfere with sterilization etc.
8. Dry carefully with sterile compressed air (max. pressure 0.5 bar). Clean the glass windows with a cotton swab moistened with 70% isopropyl alcohol to remove smears and residual moisture.
9. After preparation, inspect the 3D TIPCAM®1, the camera connecting cord and the fiberoptic light cable for damage and residual soiling.

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CAUTION: Inspect the camera connecting cord and fiberoptic light cable for breaks or cracks. Breaks or cracks will allow fluids to enter the cables and cause damage. 3D TIPCAM®1 Instruments with damaged cables may only be wipe-down disinfected. They should be sent to KARL STORZ for repair.

10. Before sterilizing, disconnect the camera head cable and the fiberoptic light cable from the 3D TIPCAM®1. Coil the camera connecting cable and fiberoptic light cable into loops measuring at least 15 cm in diameter. Do not kink or twist cables. Place the cables in the optionally available storage/sterilization container.

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4.5 Machine preparation of the 3D TIPCAM®1

If possible, cleaning and thermal disinfection should be carried out by machine (thermal disinfection at 93°C, device-specific), followed by steam sterilization at 134°C. The use of a machine for preparation offers the advantage of reducing the risks to personnel and standardizing the cleaning process (processes which can be validated or which are validated are required).

NOTE: As there is sometimes strong movement of the water, it is necessary to secure all parts in holding devices to avoid damage.

NOTE: It is not advisable to use lubricants etc. in the machine because they may cause deposits to form on the instruments and may damage plastics.

Use a cleaning/disinfecting agent approved by KARL STORZ for machine preparation (see extract from the chemicals list on page 44).

! WARNING: The 3D TIPCAM®1, which is designated as autoclavable, can be prepared in machines for semi or fully automatic cleaning and disinfection. Preliminary cleaning as described above must also be carried out for machine preparation. The choice of rinsing programs must be agreed with the machine manufacturer. In most cases, special programs for endoscopic instruments are available which are suitable for cleaning without damaging the surface.

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-  **WARNING:** Before the 3D TIPCAM®1 is placed in the cleaner/disinfector, all parts should be manually pretreated according to the cleaning instructions regardless of whether they were used or not.
-  **WARNING:** The cleaner/disinfector must fulfill all the demands made of it in terms of process validation, cleaning performance, disinfection and material compatibility.
-  **WARNING:** Ensure that the required environmental conditions with regard to maximum temperature (60 °C) are observed.
-  **WARNING:** Switching between manual and automatic preparation is not recommended.
-  **CAUTION:** Always check the 3D TIPCAM®1, the camera connecting cord and fiberoptic light cable for external damage prior to cleaning in semi- and fully automatic machines.
Procedure for cleaning, disinfection and sterilization: after the operation, rinse the 3D TIPCAM®1, immerse it in the optionally available sterilization tray, then perform machine cleaning, disinfection and sterilization and store it without further repacking.
-  **CAUTION:** After the program run, the 3D TIPCAM®1 is greatly at risk of corrosion due to the ambient conditions in the preparation machine. Preferably, remove the 3D TIPCAM®1 from the cleaner/disinfector immediately after the program has finished.
-  **CAUTION:** After machine preparation, clean the distal lens and connector of the camera connecting cable and fiberoptic light cable with a lens cleaner or a cotton tip applicator moistened in 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement. Please pay attention to page 16 (Optimum light efficiency).

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CAUTION: After machine preparation, the connection for the fiberoptic light cable must be removed in order to clean the optical fiber area with a cotton tip applicator moistened with 70% isopropyl alcohol. Then the connector needs to be screwed on again. The connector 6947095 can be removed by turning it to the left (counterclockwise). The connector can be fastened by turning it to the right (clockwise).

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4.6 Sterilization instructions for autoclavable 3D TIPCAM®1

The 3D TIPCAM®1 is autoclavable. Steam sterilization using the fractionated prevacuum procedure is recommended for initial and subsequent sterilization of the 3D TIPCAM®1. The 3D TIPCAM®1 can also be sterilized using the ethylene oxide (EtO) or the STERIS® or STERRAD method. To achieve the desired sterility assurance level (SAL) of 10, KARL STORZ recommends the following steam (or EtO, STERIS® or STERRAD) sterilization methods:

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4.7 Steam sterilization

1. Clean the 3D TIPCAM®1 designated as autoclavable as described in the preparation section of this instruction manual.

NOTE: The 3D TIPCAM®1 is in its storage/sterilization container.

2. KARL STORZ has approved the following steam sterilization methods:

- Fractionated prevacuum procedure:
Temperature: 134°C (+ 3°C)
Exposure Time: min. 4 minutes

3. Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

The recommended drying time for the 3D TIPCAM®1 is at least 4 minutes to ensure that the connector of the camera connecting cord and the fiberoptic light cable are completely dry before being plugged in.* Before connecting the light cable with the camera connector, the 3D TIPCAM®1 must first cool down to room temperature.

4. After steam sterilization do not touch the 3D TIPCAM®1 until it has cooled down.

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Maintenance

Cleaning, disinfection and sterilization

Caution: Do not soak the 3D TIPCAM®1 in cold liquid after removing it from the steam sterilizer. This may lead to a 'temperature shock' which can cause damage to the camera lens, and to rendering the warranty void.

- For steam sterilizers which do not have an integrated drying cycle, the connectors of the camera connecting cord and the light cable may be dried with a lint-free cloth once they have cooled down. All moisture must be removed from the connectors before insertion.

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4.8 Ethylene Oxide (EO) gas sterilization

- Clean the camera head unit as described in the preparation section of this instruction manual. Place the camera head unit in a sterilization tray. Do not place any instruments on top of the camera.
- All sterilization cycles must include a pre-conditioning cycle:
- Temperature: $54 \pm 2^\circ\text{C}$ /
Time: 30 minutes
Humidity: $\geq 70\%$ (relative humidity)
Vacuum: 1.3 psia

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Maintenance

Cleaning, disinfection and sterilization

4. KARL STORZ has validated EtO sterilization using the following parameters:
- | | |
|----------------------|------------------------------|
| Gas mixture: | 100% EtO |
| Temperature: | Set point of 55°C |
| Relative humidity: | ≥ 70% RH |
| Humidity dwell time: | 30-45 minutes |
| Exposure time: | 180 minutes,
(full cycle) |
| EtO concentration: | 735 ± 30 mg/l |
5. Airing can be accomplished in an airing chamber at a temperature of 50 - 54°C for 12 hours.
- Maximum permissible levels of residue following EtO sterilization (according to ANSI/AAMI/ISO 10993-7):
- | | |
|------------------------|--------------|
| Ethylene oxide: | 4 mg (< 24h) |
| Ethylene chlorohydrin: | 9 mg (< 24h) |

NOTE: Sterilization methods which deviate from the validated processes described must be validated individually. The country-specific and sterilization device-specific requirements and outline conditions concerning the effectiveness of the sterilization must be taken into account, as must, where applicable, the necessary airing times.

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4.9 Sterilization with the STERRAD® 100NX™ system

Two sterilization cycles are available for selection with the STERRAD® 100NX™ system. The STANDARD cycle time is approximately 47 minutes. The FLEX cycle time is approximately 42 minutes.

CAUTION: STERRAD® sterilization may cause cosmetic changes to the 3D TIPCAM® 1 that do not usually impact the camera's functionality.

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Maintenance

Cleaning, disinfection and sterilization

CAUTION: The 3D TIPCAM®1 must be thoroughly dried before loading into the STERRAD® sterilization chamber. Loads containing moisture may cause a cycle cancellation.

1. Clean and prepare the camera head unit as described in the preparation section of this instruction manual. It must be completely dry.
2. Place the camera head unit in the STERRAD® sterilizer so that the hydrogen peroxide can flow around it. Do not allow the camera head unit to touch the wall of the sterilizer.
3. Press the START button. The door closes automatically.
4. Detailed instructions can be found in the instruction manual for the STERRAD® 100NX™ sterilization system. The latest information on sterilization with the STERRAD® 100NX™ system is available from Advanced Sterilization Products (ASP) (www.aspij.com).

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**4.10 Sterilization with the
STERRAD 100S system**

CAUTION: STERRAD® sterilization may cause cosmetic changes to the 3D TIPCAM® 1 that do not necessarily impact the camera's functionality.

CAUTION: The 3D TIPCAM® 1 must be thoroughly dried before loading into the STERRAD® sterilization chamber. Loads containing moisture may cause a cycle cancellation.

1. Clean and prepare the camera head unit as described in the preparation section of this instruction manual. Be sure that the camera head assembly is completely dry.
2. Load the camera head assembly into the STERRAD® sterilizer, arranging it so that the hydrogen plasma can surround it. Do not allow the camera head assembly to touch the wall of the sterilizer.

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Maintenance

Cleaning, disinfection and sterilization

1. Please consult the STERRAD[®] 100S Sterilization System operator's manual for detailed instructions for use. Please contact Advanced Sterilization Products (ASP) for the most up to date information regarding sterilization with the STERRAD[®] 100S Sterilization System.

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4.11 Sterilization with the STERIS System1®

WARNING: Failure to thoroughly clean and properly position devices so that all surfaces will be exposed to the sterilant or overloading of the processing container may result in an ineffective sterile process and/or damage to the devices.

WARNING: Use only STERIS® processing containers to ensure proper sterilization of the devices. Do not use the STERIS® system for devices that cannot be immersed in liquid.

WARNING: Always verify that the sterilant container is empty after the cycle is complete. If the sterilant container is not empty, then the load cannot be considered sterile.

WARNING: Devices are not sterilized and/or adequately rinsed when a sterile processing cycle is cancelled.

1. Make sure that all instruments/devices have been completely cleaned in accordance with the methods recommended in the preparation section of this instruction manual.
2. Visually inspect and test the devices to ensure that they are working properly before processing in the STERIS® processor.

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3. Place water resistant caps on cameras and other devices as required according to manufacturers' instructions. Failure to properly place the caps may result in damage to the camera or device.
4. Position the device in the STERIS® tray to ensure that all surfaces will be exposed to the sterilant. All devices should be disassembled to the smallest components possible.
5. The equipment is ready to use immediately on completion of the STERIS® cycle.
6. Please consult the STERIS® operator's manual for more detailed instructions for use of the STERIS® processor. Please contact STERIS Corporation for the most up to date information regarding processing of cameras in the STERIS SYSTEM1® Processor.

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Maintenance

**Cleaning, disinfection and
sterilization**

**4.12 Sterilization with the
AMSCO® V-PRO™ 1
(from STERIS®)**

WARNING: If the units are not carefully cleaned and positioned so that all surfaces come into contact with the sterilization medium, or if the container is overloaded, this can lead to an ineffective sterilization process or cause damage to the units.

WARNING: Only use STERIS® containers to ensure correct sterilization of the units.

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WARNING: If a sterilization cycle is canceled, the devices are not properly sterilized and/or rinsed and dried. The procedure must be repeated under all circumstances (check the cause of the error!).

1. Ensure that all instruments/units are fully cleaned according to the recommended methods given in the preparation section of this instruction manual.
2. Perform a visual check of the instruments/units before preparation in the STERIS® unit and check them for correct functioning.

Observe the detailed descriptions pertaining to the Amsco® V-PRO™ 1 which can be found at www.steris.com.

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- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- The Difficulty of Reprocessing Reusable Rigid Laparoscopic Forceps and other Endoscopic Accessories: Are Disposables the Answer? Health Devices, Vol. 23, Nos. 1-2, pp. 57-58, 1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

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4.13 Use of sterile covers

1. Clean and prepare the 3D TIPCAM*1 as described in the preparation section of this instruction manual.
2. Follow the instructions provided by the manufacturer when using sterile covers.

4.14 References for cleaning and sterilization

- Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.
- Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.
- The Difficulty of Reprocessing Reusable Rigid Laparoscopic Forceps and other Endoscopic Accessories: Are Disposables the Answer? Health Devices, Vol. 23, Nos. 1-2, pp. 57-58, 1994.
- Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, No. 1, pp. 23 - 29, 1995.

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5 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid.

NOTE: To ensure safe operation of the camera control unit (CCU), you should perform the following procedures every 12 months:

- Check leakage current from chassis < 100 µA.
- Check ground impedance < 0.1 ohms.
- Check power consumption in rated power.

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5.1 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

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5.2 Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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5.3 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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5.4 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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5.5 Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

! NOTE: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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6 Technical description

6.1 Technical data 3D TIPCAM®1

- Image sensor: 1/6" CCD
- Depth of focus: 15 mm ... 100 mm
- Outside diameter: 10 mm
- Effective length: 310 mm
- Overall length: approx. 550 mm
- Direction of vision: 0°
- Operating temperature: 10°C ... 40°C

Storage:

- Ambient temperature: -10°C ... +60°C
- Humidity: 5 ... 95% (non-condensing)
- Atmospheric pressure: 500 ... 1080 hPa

In conjunction with a KARL STORZ

- 3D CCU the 3D TIPCAM®1 satisfies the IEC 60601-1:
- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF 

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets

Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOUR[®] telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflexors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Exelabo GmbH & Co. OHG	Seucon Plus	1,2,3
	Seucon P. bei Camph	1,2,3
	Ultrasein A	1,2,3
Exelabo Pharma GmbH	Ultrasein Alpha	1,2,3
	Ultrasein - Wasser	1,2,3
	Acid	2
Fresenius AG	A-1-Job	2
	Synco 4+	2
Hofa Pharma AG, a.s.	Synta-Hofa	1,2
Forbis Industrie - Chemische LTDA	Lebichoco	2
	Paroche 4000 - gelbe	2
Jose Colado S.A.	Paroche 500G	2
	Paroche - Silber (2000)	2
	Aldalen 2000	1,2
Lysocom	Allyval	1,2
Dr. Rosemann GmbH	Deodorant	1,2,3
	Lysoform 500G	1,2
Medichem International	Medi 10	2
	Medi 1000	2
Marti Hygiene GmbH	Exaform 5000	1,2
	Macaroni 5	2
	Macaroni - Granulat	1,2
	Macaroni 1 neu	1,2,3
MEI/AS Medientechnik GmbH (Internat.) L. 102/21	D 50	2
Pharmatec	4 20	1,2,3
Promocel Ltd	Washite	2
	Leucoform 1000	1,2
	Leucoform 1000	2
Schumacher DL GmbH	Leucoform 1000	2,3
	Leucoform 1000	2
	Leucoform 1000	1,2,3
	Leucoform 1000	2
Schulke & Nime GmbH	Glaxoform AF 1000	1,2
	Glaxoform AF 1000	1,2
	Glaxoform AF 1000	1,2,3
	Glaxoform AF 1000	1,2
	Glaxoform AF 1000	1,2
Siera	Proclara 2	1,2
	Proclara 100 2	1,2
	Proclara 1000	1,2
	Proclara 1000	1,2
Wegert, Dr. GmbH & Co	Proclara 1000	1,2
	Proclara 1000	1,2
	Proclara 1000	1,2
	Proclara 1000	1,2
	Proclara 1000	1,2
Wiley Industries PTY Ltd	Proclara 1000	1,2
	Proclara 1000	1,2

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination shaft, **not** the housing/handle part.

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Typ

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes
(according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schulze & Wey, GmbH	Thermosol alkaline	1
	Thermosol CA	2
	Thermosol ED	2
	Thermosol ER	1
	Thermosol PR2-Sym	1
Storz	InstumClean	1, 3
	InstumClean	1, 3
Weigert Dr. GmbH & Co.	Neosol AF	1
	Neosol AF forte	1, 3
	Neosol AF forte	1
	Neosol AF forte	1, 3
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
Weigert Dr. GmbH & Co.	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1
	Neosol AF forte	1

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Firmenstempel/Unterschrift:

WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be completed by the supplier/importer:

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 720 41097, Fax: +53 720 41098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 5250-5607, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
R. Joaquim Floriano, 413 - 20º andar -
Itaim Bibi, São Paulo - SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskopi Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: P.O. Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4th floor, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Høje, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopier.be

KARL STORZ Endoscopia France S.A.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colanin, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 31 4250800, Fax: +40 31 4250801
E-Mail: office@karlstorz.ro

KARL STORZ Endoskope Greece M.E.P.E.
Geunari 40,
15343 Agia Paraskevi/Athens, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
"Repair & Service Subsidiary"

KARL STORZ Industrial**
Gedik is Merkezi B Blok
Kat 5 D. 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

OOO KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
Solidere - Beirut Souks
Block M. 3rd Floor
2012 3301 Beirut, Lebanon
Phone: +961 1 1993390, Fax: +961 1 1993391

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061
Roggebaai 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasakhstan
Bokeykhan Ulica 8/2, VP - 1
010000 Astana, Respublika Kasakhstan
Phone: +7 7172 57-52 16 / 57 28 49 / 57 09 34
Fax: +7 7172 43 96 96
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Phone: +971 4 4471230
Fax: +971 4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911
Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325548, Fax: +65 65323832
E-Mail: info@karlstorz.com.sg
E-Mail: servicesg@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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

**Important information
for users of KARL STORZ
devices**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern high quality piece of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid to the proper setup, installation, and operation of the 3D CCU camera system. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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
General warning

Users of KARL STORZ devices should be thoroughly trained in their use and the applicable medical procedures. Use of these devices and accessories should not be undertaken until the user has fully familiarized himself/herself with the instructions for use, assembly and care. Instruction manuals should be carefully studied and be available to the surgical (endoscopic) team during the procedure; it is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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**Controls, displays,
connectors
and their uses**

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
① **NOTE:** 3D is only on the DVI connection ⑨ possible.

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Read the instructions carefully before operating the
equipment!

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Aus

Power off

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Potential equalization conductor

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Type CF equipment, defibrillation proof

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DANGER: Risk of explosion if used in the
presence of flammable anesthetics.

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CAUTION: To reduce the risk of electrical
shock, do not remove cover
Refer servicing to qualified service personnel.
Do not store liquids on or above the unit.
Keep out of reach of patients.

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Device is subject to the requirements of the WEEE Directive, 2002/96/EC

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4.8.3	Entering patient information (USB storage device and keyboard required)	31
4.8.4	Image Capture Menu.....	32
4.8.5	Capturing Still images.....	34
4.8.6	Storing still images to a USB storage device.....	35
4.8.7	Printing Still Images	36
4.8.8	Print During Capture Feature	37
4.8.9	Capturing Video.....	38
4.8.10	Ending a Capture Session	40
4.8.11	Low USB Memory Warnings.....	40
4.8.12	Reviewing Video and Still Images on a Computer	41
4.9	Image	42
5	Camera Setup.....	43
5.1	SCB.....	43
6	Maintenance.....	47
6.1	Cleaning, disinfection and sterilization.....	47
6.1.1	Care and handling	47
6.1.2	Cleaning instructions.....	48
6.2	Fuse replacement.....	49
6.3	Maintenance and safety check	50
6.4	Servicing and repair	50
6.5	Disposal.....	51
6.6	Repair program.....	51
6.7	Limitation of liability.....	52
6.8	Manufacturer's warranty	52
7	Technical description	55
7.1	Troubleshooting	55
7.2	Technical specifications	58
7.3	Equipment classification	60
7.4	Equipment test certificates.....	60
7.5	Technical documentation.....	60
7.6	Software User Rights.....	61
7.7	KARL STORZ Communication Bus (SCB).....	62
8	Appendix	63
8.1	Electromagnetic Compatibility (EMC) Information.....	63

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1 General information**1.1 Device description**

The 3D CCU is a revolutionary endoscopic 3D video camera system. It is used together with the KARL STORZ 3D TIPCAM1[®].

A keyboard is included with the 3D CCU to access camera functions and set displays such as time/date and patient information.

The 3D CCU incorporates the KARL STORZ Communication Bus (SCB) for integration into a complete operating room system. The SCB system allows communication between the KARL STORZ devices connected to the system. This communication enables the parameters of the SCB devices to be displayed centrally. Further information on the use and function of the SCB system can be found in the instruction manual for the KARL STORZ SCB system which can be requested free of charge from KARL STORZ or downloaded by registered persons as a pdf from www.karlstorz.com.

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Safety instructions Warnings and cautions

WARNINGS:

WARNING: The instructions and interface specifications for medical devices used in combination must be observed precisely.

WARNING: The product described here must be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

WARNING: Test this equipment prior to each surgical use. In the event that the image becomes unusable during surgery, it is left to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

WARNING: Do not use in the presence of flammable anesthetics. There is a risk of explosion.

WARNING: Danger of electric shock. Do not open the unit. Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

WARNING: Grounding reliability can only be achieved when the equipment is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC, CEC and NEC standards. Use only in hospital grade receptacles.

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WARNING: Always unplug the unit before performing any maintenance including cleaning.

WARNING: Keep out of reach of patients.

WARNING: For reasons of safety, do not simultaneously touch the device output connectors and the patient.

WARNING: It is essential to follow the instruction manual closely when cleaning the unit and preparing/sterilizing unit accessories.

WARNING: Any deviations from the recommended parameters for sterilization should be validated by the user.

WARNING: Accessory equipment connected to the analog and digital interfaces (signal inputs and outputs) must be set up and tested according to IEC 60601-1 for medical electrical equipment and additionally to IEC 60950 for data processing equipment. Furthermore, all configurations shall comply with the system standard IEC 60601-1-1. Any person who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical service department or your local representative.

WARNING: To avoid burns during endoscopic procedures, operators must use caution when employing non-BF/CF type equipment while using electrosurgical devices.

WARNING: The high-intensity light that radiates through the 3D TIPCAM® can lead to high temperatures at the light outlet, at the end faces of the light cable and at the tip of the 3D TIPCAM®. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera. This does not apply for fluorescence applications; see instruction manual for D-Light C/AF.

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WARNING: When endoscopes are used with energized endoscopically-used accessories, the patient leakage currents may be additive. This is particularly important if a Type CF endoscope is used, in which case a Type CF endoscopically-used accessory should be used in order to minimize the total patient leakage current.

WARNING: Before each use, the outer surfaces of the portions of the 3D TIPCAM® and the endoscope accessories which are to be inserted into the patient should be checked to ensure there are no unintended rough surfaces, sharp edges or protrusions which may cause a safety hazard.

WARNING: When using the camera system in a surgical discipline employing a laser, caution should be taken to prevent eye damage to the operator. The operator should wear protective glasses whenever viewing the surgical site directly through the endoscope.

WARNING: Perceived shocks may be experienced when using the camera system in conjunction with electrosurgical units.

WARNING: Electrostatic discharge (ESD) events may be perceived as shocks when initially touching the camera system.

WARNING: During a defibrillation discharge, the display image may be interrupted by ≤ 1 second.

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CAUTION: This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.

CAUTION: Do not store liquids on or above the device. Avoid allowing fluids to enter the unit.

CAUTION: Only operate the unit with the voltage stated on the unit's identification plate.

CAUTION: When replacing fuses, use only fuses of the same rating.

CAUTION: Federal law restricts this device to sale to or on the order of a physician.

NOTE: Disposal of these system products, at the end of their useful life, must be in accordance with local regulations. **CCU component contains lithium battery.**

NOTE: Do not discard as unsorted municipal waste.

NOTE: Discard as electrical/electronic waste; recycle or reuse accordingly.

NOTE: Consult local authorities for reuse/recycle instructions.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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* CCU = Kamerakontrolleinheit (camera control unit)

2.3 Intended use

The 3D CCU is connected to a KARL STORZ 3D TIPCAM1[®] and used for general endoscopic procedures. Use of the camera in other applications is not allowed for safety reasons. The 3D System is intended to be used by physician/endoscopist in the visual examination and treatment. The 3D CCU may only be used with accessories, wearing parts, and disposable items which are designated by KARL STORZ as suitable for the camera or the safe use of which is proven. For safety reasons, do not perform unauthorized conversions or modifications to the camera.

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
2.4 User qualification

The 3D CCU may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the equipment.

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2.5 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The unit is equipped with a connector for attaching a ground line. It should be connected before power is applied to the unit in compliance with current national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

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2.6 Safety precautions when operating the equipment (electromagnetic compatibility)

It is the user's responsibility to make sure the equipment is safe and operates properly before use.

CE marked equipment has been tested and found to comply with the EMC limits for the Medical Device Directive 93/42/EEC (EN 55011 Class A and EN 60601-1-2). These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. The equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference with other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected
- Consult the manufacturer or field service technician for help

NOTE: A KARL STORZ camera system comprises a camera processor and a 3D TIPCAM1[®], both of which must bear the CE mark for the system to meet the requirements of the Medical Device Directive. If either the 3D TIPCAM1[®] or processor are not CE marked, this combination of components does not meet the requirements of the Directive.

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3.3 Installing the 3D CCU


WARNING: The 3D CCU may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

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
WARNING: This CCU is not intended for use in hazardous zones. Do not operate the CCU within the hazard zone shown in the diagram while explosive anesthetic gases are in use. The CCU is equipped with a connector for attaching an equipotential wire (redundant grounding wire) line. The CCU's equipotential line should be installed by a qualified electrician.

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3.3.1 Connecting power

Set the CCU on a flat surface. Make sure there is sufficient distance on all sides to other equipment (especially radio frequency surgical equipment) and objects. Before plugging in the CCU, ensure that the voltage on the nameplate corresponds to the voltage of the local power line.

Connect power cord. Insert power cord into power cord receptacle as far as it will go. (This should only be done outside potentially explosive locations.)

WARNING: Always use a hospital grade power cord with this camera.

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
NOTE: The image Capture System is compatible with a USB printer only.

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3.3.3 Powering up CCU

! WARNING: Test the 3D TIPCAM1* and video camera before each procedure. Ensure that the proper video image appears on all video monitors before beginning each procedure. Switch on the camera control unit and the video monitor. A test pattern will appear on screen.

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3.3.4 Monitor test-pattern

1. A test pattern appears when the CCU is turned on without a 3D TIPCAM1* connected.

The test pattern (picture) displays the letters 'L' (red) and 'R' (blue).


- Put the 3D glasses on
- When you close your left eye, you should only be able to see blue Rs
- When you close your right eye, you should only be able to see red Ls
- If this is the case, the system has been set correctly.

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3.3.5 Connecting the light cable to the 3D TIPCAM1*

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3.3.6 Connecting the 3D TIPCAM1* to the KARL STORZ light source

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4.1 Basics

4.1.1 Connecting the 3D TIPCAM1® to the CCU

Plug the cable connector into the cable receptacle on the CCU's front panel. If the connector is not plugged into the CCU correctly and completely, the color test pattern will remain on the monitor.

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
4.2 3D TIPCAM1* button basics

The 3D TIPCAM1* has the following 3 buttons for function access and menu control:

- ① - **scroll up in menu / programmed functions:** scrolls up in the menus / performs the function programmed for this button.
- ② - **select / function:** selects the highlighted option / programmed functions.
- ③ - **scroll down in menu / programmed functions:** scrolls down in the menus / performs the function programmed for this button.

① NOTE: Buttons can be programmed. The steps required to this end are: select the main menu (press the button ■ on the 3D TIPCAM1* or the ENTER button on the keyboard) → Settings → Button assignment. For details see page 25.

For simplicity, this manual provides keyboard instructions for programming (i.e., “Press ENTER to select”). However, camera head buttons can also be used (i.e., “Press n to select”).

The buttons of the 3D TIPCAM1* can, however, also be programmed to operate two functions each (e.g., freeze and capture).

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
4.2.1 Keyboard basics

A keyboard may be used to access, program, and control camera features and functions. The buttons for controlling the menus correspond to the 3D TIPCAM1[®] buttons:

- ➡ scroll up
- ENTER menu/select
- ➡ scroll down

The keyboard is also used to input patient information. When inputting patient information, all text keys are functional and specific control keys are available. For more information, see "Operating instructions – Patient Information section of this manual.

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- Einstellungen
- Patientendaten
- Dokumentation
- Bild
- SCB
- Beenden

User Menu

- Setup
- Patient Information
- Capture
- Image
- SCB
- Exit

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4.2.2 Menu basics

Software menus are provided to program and control the camera features/functions.

The User Menu is accessed by pressing the ③ button on the 3D TIPCAM1[®] ▼ or ENTER on the keyboard.

The symbols used in the menus include:

- ▶ indicates that options are available with that particular menu item. To access these options, highlight the menu item and then press ENTER on the keyboard or the ■ button on the 3D TIPCAM1[®] (long press).
- • • indicates that there are sub-menus within that menu item. To access the sub-menus, highlight the menu item and then press ENTER on the keyboard or the ■ button on the 3D TIPCAM1[®].
- ▼ indicates that more menu items are available below those displayed.
- ▲ indicates that more menu items are available above those displayed.
- () indicates a selected value within the parenthesis.

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4.3 Camera functions

The 3D video system has the following functions:

Accessory 1: For use with a KARL STORZ AIDA image and video capture system.

Accessory 2: For use with a KARL STORZ AIDA image and video capture system.

Accessory Display Messages: Displays accessory in use (i.e. "Printer") or no accessory ("None").

Automatic Exposure Control: Patented system automatically sets exposure depending on the available light.

Brightness: Selects the brightness of the image, including *Low, Medium and High* settings.

Color: Allows adjustment of red and blue color values.

Color bars: Turns test pattern On or Off.

Date and Time Setup: Allows entry of date and time.

Enhancement: Electronically controls the level of contrast and definition of the image, including *Off, Low and High* settings.

Image Capture: Allows still image and video capture to a connected USB storage device and/or USB printer.

3D Mode: Allows the DVI signal to be switched between 2D and 3D.

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Patient Information: Using keyboard, patient information is entered and may be displayed on the monitor and/or captured videos and prints.

Button Programming: Allows functions/options to be activated with the buttons of the 3D TIPCAM1[®] (▲, ▼ or ■).

SCB: Displays up to 4 SCB devices for on-screen monitoring.

Setup: Configures and stores settings in camera system.

Shutter (exposure control): This camera features a patented auto exposure mode which selects exposure automatically depending on the available light. The automatic exposure system may be overridden to select one of seventeen manual shutter speeds ranging from 1/60 to 1/17,000 second. Exposure settings can be changed through the "Shutter" camera function. This function can be programmed to allow it to be controlled via the 3D TIPCAM1[®] or the User Menu.

White Balance: Determines the camera system's color controls based on the light source in use.

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4.4 White Balance

The white balance determines the camera system's color controls based on the light source in use.

To perform this adjustment:

1. Turn on the light source and point the 3D TIPCAM1[®] at a white surface, with no colored objects visible within the camera's field of view. The view should fill at least 70% of the area of the screen.
2. Start the white balance using the white balance button on the CCU (see photo on the left). Once activated, 'White balancing' appears on the monitor.
3. The white balance has been successfully performed when 'White balance O.K.' appears on the monitor. If white balancing has not been performed correctly, one of the following messages will appear:

'White balance fail – Dark' Indicates that the field of view is too dark. In this case, shorten the distance between the 3D TIPCAM1[®] and the white surface and perform the white balance again.

'White balance fail – Bright' Indicates that the field of view is too bright. In this case, lengthen the distance between the 3D TIPCAM1[®] and the white surface and perform the white balance again.

'White balance fail – Color' Indicates that the color temperature exceeds the white balance range.

NOTE: The last white balance setting remains stored in the 3D TIPCAM1[®] even after the camera has been switched off. When the camera is next used, a new white balance should be performed if a different light source or 3D TIPCAM1[®] is used.

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

Notes:

- ^{A2} SCB Menu and commands are only visible in KARL STORZ cameras, if at least one SCB device is identified including camera.
- ^{A4} SCB Devices is a dynamic list, based on the devices found on the bus
- ^{A5} SCB Parameters list is a function of the device that has been selected from the SCB Devices menu.
- ^{A6} Only if 3D TIPCAM1* is connected.

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PROGRAM HEAD BUTTONS

Button Programming

1 Primary (Accessory 1 •••)
1 Secondary (None •••)
2 Primary (User Menu •••)
2 Secondary (None •••)
3 Primary (Accessory 2 •••)
3 Secondary (None •••)
Back
Exit

DATE AND TIME SETUP – To program the date/time:

1. Select a category (Year, Month, Day, Hour, or Minute) using the ◀ or ▶ keys. Press ENTER.
2. The ◀ or ▶ keys will now allow you to scroll through options for each category (i.e. "Month" will have values from 1-12). Once the desired value is on-screen, press ENTER. The highlight will now return to the list of categories for further changes.

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CAPTURE SETUP

Capture Setup

Print During Capture	[OFF]
Prints Per Page	[6]
Print Copies	[1]
Paper Size	[Not Available]
Counter Display	[ON]
Still Counter	[ON]
Print Counter	[ON]
Video Counter	[ON]
Back	
Exit	

IMAGE SETUP

Image Setup

Brightness Default	[Medium]
Enhancement Default	[OFF]
Back	
Exit	

DISPLAY PATIENT SETUP

Display

Name	[OFF]
Identification	[OFF]
Date of Birth	[OFF]
Institution	[OFF]
Display Patient Record	[OFF]
Still Counter	[OFF]
Back	
Exit	

CCU SETUP

CCU Setup

Language Display	[English]
Keyboard	[English]
Institution	[• • •]
Accessory 1	[None]
Accessory 2	[None]
Back	
Exit	

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Language Selection:

Change the language of the on-screen display and the keyboard language display (in the menu 'CCU Settings'. This can be accessed from the main menu by selecting the menu items 'Settings' - 'CCU'.)

1. Using the  or  buttons highlight the words **Language Display**.
2. Press ENTER to select.
3. Cursor will move to the right. Scroll through the list of languages using the  or  key (English, Espanol, Francais, Italiano, Portugues, Deutsch). When the desired language is on-screen, press ENTER. The highlighted area will move back to the menu on the left.
4. To indicate the language of the keyboard being used, highlight the word **Keyboard** and press ENTER.
5. Scroll through the list of languages (English, Espanol, Francais, Italiano, Portugues, Deutsch) using the  or  key. When the appropriate language appears on-screen, press ENTER.

To proceed to the next screen, highlight **Next** and press ENTER.

Institution

1. To enter Institution Name, use the arrow keys to highlight the word Institution.
2. Press ENTER to select.
3. Enter Institution name, which will automatically appear in Patient Information Records.

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
Controlling accessory devices (Acc 1 and Acc 2)

If controlling accessory devices such as an AIDA capture device, video printer, or VCR with the 3D TIPCAM1[®]'s user control buttons, the accessory cables must be installed. Plug one end of the cable into the chosen accessory output on the CCU and the other end to the footswitch input on the accessory device. It is recommended that the ACC1 output be used to control the still capture of the AIDA or video printer, and the ACC2 output be used to control the AIDA video capture or the VCR. If a video recorder, printer, or similar accessory device is connected to outputs ACC1 or ACC2, (Note: A video recorder may be connected to ACC2 only), the following functions can be activated via the video endoscope buttons, User Menu.


Printer: Capture and print an image.

AIDA (digital image capture device): Capture an image (one camera head button); video recording **Record and Stop** (the other button).

The options for the accessory display message, which will appear on the main monitor, include: None, Accessory 1/2, VCR, Printer, and Acc1-Still/Acc2-Video.

NOTE: Use "Accessory1" and "Accessory2" to control the AIDA capture functions. The "Capture Still" and "Capture Video" functions in the Image Capture menu only control the internal Image Capture in the 3D CCU.

Example (diagram below): If you select the ACC1 function with your 3D TIPCAM1[®], a control signal is, for example, sent from the ACC1 of the 3D CCU to an AIDA and consequently the required function (e.g., capture still image) is performed.

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4.7 Patient Information

NOTE: Patient Information option is not available unless a keyboard is plugged into the CCU.

When selecting **Patient Information**, the following menu will appear:

1. Highlight a function from the following list:

Display patient information

Display	
Name	[OFF]
Identification	[OFF]
Date of Birth	[OFF]
Institution	[OFF]
Display Patient Record	[OFF]
Back	
Exit	

Select Patient Record: To select one of ten pages available for patient information, scroll through the list of pages and select the desired page.

Select Patient Record	^v
Name	<patient name>
Identification	<patient ID>
Date of Birth	YYYY-MM-DD

Clear Patient Information: Clears all stored text from patient information pages.

Clear All Patient Records?	
Yes	
No	

Edit Patient Record: When inputting or editing patient information, the keyboard becomes completely functional, including the following key assignments:

Edit Patient Record	^v
Name	<patient name>
Identification	<patient ID>
Date of Birth	YYYY-MM-DD
Physician	<dr name>
Procedure	<procedure>
Institution	<institution>
Next	
Select	
Clear Record	

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Capture

Capture Still
Capture Video
Print Now
Patient Information * * *
Capture Setup * * *
Exit

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Capture Setup

Print During Capture Off/On
Prints Per Page 2 4 6 8
Print Copies 1 2 3
Paper Size
Counter Display Off/On
Still Counter Off/On
Printer Counter Off/On
Video Counter Off/On
Back
Exit

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4.8 Image and Video Capture

① NOTE: Images and videos are always captured in 2D.

The 3D CCU is equipped with a built-in image Capture Module that allows Still Capture and Video Capture of images taken with the 3D CCU. Stills and videos may be captured onto a USB storage device, and still images may be printed via a USB printer.

4.8.1 USB storage device and printer connection

① NOTE: A USB storage device must be formatted to FAT 32 in order to store still images and/or video.

A USB storage device and/or a USB printer may be inserted/connected directly into the USB output on the front panel of the 3D CCU or the USB output on the back of the CCU at any time during a procedure. On-screen messages will temporarily appear on the video monitor indicating connection or disconnection. Only one storage device and one printer may be connected at the same time. If two USB storage devices or printers are connected at the same time, an error message will appear. If two printers are connected at the same time, both printers must be disconnected and one printer then reconnected in order to print.

① NOTE: It is recommended to always have a USB storage device connected to ensure captured stills are not "lost" in the event of accidental powering off of the CCU prior to the completion of printing.

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4.8.2 Starting a Capture Session

To start a capture session, switch on the CCU and connect the 3D TIPCAM1* and an USB storage device and/or a USB printer. If a keyboard is attached, the patient information menu will appear.

i NOTE: To start a capture session using a USB printer only, with no USB storage device connected, "Print During Capture" must be set to "On" in the Image Capture Menu. See "Print During Capture Feature" below section on the next pages.

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4.8.3 Entering patient information (USB storage device and keyboard required)

Patient data may be entered before a 3D TIPCAM1* is connected. Patient data records may also be added after connection of a 3D TIPCAM1*. 10 patient records may be entered and available as long as the CCU is not turned off. If the CCU is turned off, patient information will be lost. A capture session may be started with or without first entering patient information. Entered patient information may be selected, at which time it will be added to the folder name (stored to the USB storage device) that already includes the date and time of the capture session.

After turning on the 3D CCU and connection of a USB storage device, if a keyboard is connected the patient information menu will appear. If a keyboard is not connected, this menu will not appear. Using a keyboard input the desired patient information. Use the Tab, Enter or Down/Up Arrow keys to move to the next field and hit Enter when "Exit" is highlighted to exit the menu.

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4.8.4 Image Capture Menu

The Image Capture Menu provides access to the capture and print functions and to the Image Capture Setup menu.

- **Capture Still:** Activates a still image capture. See “Capturing Still Images” section below for more information.
- **Capture Video:** Activates and deactivates recording of standard definition video. See “Recording Video” section below for more information.
- **Print Now:** Instantly prints all images in the print queue. The print queue size is based upon the selected number of prints per page (see below “Prints Per Page” section on the next page) in the Image Capture Setup menu.

NOTE: When the print queue is full and the next still is captured, the previous still images in the queue are cleared and are no longer available for printing. For example, if the prints per page is set to 4, and the 5th still is captured, only the 5th still will print as it is the only image now in the queue.

NOTE: To print all captured still images, “Print During Capture” must be set to “On”. See “Image Capture Setup” and “Print During Capture Feature” sections on the next pages.

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Image Capture Setup

Print During Capture (On)/Off
Prints Per Page ►
Print Copies ►
Paper Size ►
Counter Display ►
Print Counter (On)/Off
Video Counter On/Off
Back
Exit

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- Patient Information
- Capture Setup
- Image Capture Setup:

Before beginning a capture session, check the settings to insure that the system is set to the desired preferences. The Image Capture Setup menu allows selection of the following settings:

- **Print During Capture:** Selects the ability to print still images immediately after they are captured, and the selected prints per page layout is full. When set to "On", this feature also enables still image capture with only a USB printer (no USB storage device) connected. See "Print During Capture Feature" below section on the next page.
- **Prints Per Page:** Selects 2, 4, 6 or 8 images to print on one page.
- **Print Copies:** Selects 1, 2 or 3 copies of each printed page.
- **Paper Size:** Selects paper size.
- **Counter Display:** Selects the display of the Print Counter and Video Counter below. The display options are Off (never display), Temporary (temporary display upon each still image capture and video capture), or Permanent (always display).
- **Print Counter:** Displays the number of captured stills currently in the print queue based on the selected prints per page (above) in the lower left corner of the surgical monitor. Options are On or Off.
- **Video Counter:** Displays the remaining minutes of video that will fit onto the connected USB storage device in the lower right corner of the surgical monitor. Options are On or Off.

NOTE: Messages generated by the Image Capture system could interfere with SCB messages that are displayed in the same location on the surgical monitor.

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4.8.5 Capturing Still Images

The Image Capture system captures still images. To capture still images, a USB storage device or a USB printer must be connected to one of the USB outputs on the 3D CCU. A USB storage device and a USB printer may also be connected at the same time.

ⓘ NOTE: if using a printer only, with no USB storage device connected, "Print During Capture" must be set to "On" in the Image Capture Menu. See "Print During Capture Feature" under "Printing Still Images" below section on the next page.

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Still images can be 'captured' using an specifically programmed button on the 3D TIPCAM1® (further information on the programming of buttons on the 3D TIPCAM1® can be found in the section 'Programming of camera head buttons' in this manual). Still images may also be captured via a keyboard using the 'Capture Still' command in the Image Capture Menu or via the buttons on the 3D TIPCAM1® as described on the previous page.

4.8.6 Storing still images to a USB storage device

To store captured still images, a USB storage device must be connected to one of the USB outputs on the 3D CCU. A printer is not required to store images to a USB device.

All still images are stored to a USB storage device in .jpeg format. Captured still images may then be viewed on, and saved to, a computer.

ⓘ NOTE: Only one USB storage device may be connected at one time.

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4.8.7 Printing Still Images

To capture and print still images, a USB printer must be connected to one of the USB outputs on the 3D CCU. A USB storage device is not required to only capture and print (not store) still images.

Still images may be automatically printed as they are captured (see 'Print During Capture Feature' below) or manually printed via a correspondingly programmed 3D TIPCAM1[®] button or via a button on the 3D TIPCAM1[®] using the 'Print Now' command in the Image Capture Menu as described above. The number of prints per page and the 1, 2 or 3 copies of each page can be set in the Image Capture Setup menu as described above in the section on the previous page.

ⓘ NOTE: The print queue is erased when the CCU is powered off. Be sure that all prints have been printed prior to powering off the CCU. It is recommended to always have a USB storage device connected to ensure captured stills are not "lost" in the event of accidental powering off of the CCU prior to the completion of printing.

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4.8.8 Print During Capture Feature

If "Print During Capture" is set to "On", all captured images will print immediately once a page is full ("full" is based on the selection of images per page set in the Prints Per Page option in the Image Capture Setup menu). For example, if Prints Per Page is set to 4, printing will automatically begin after the 4th still capture. Setting "Print During Capture" to "On" also enables still image capture with only a USB printer (no USB storage device) connected.

If "Print During Capture" is set to "On", any printed pages which are not full when capturing is completed (and therefore have not printed) will print automatically upon removing the USB storage device. Setting "Print During Capture" to "On" also enables still image capture with only a USB printer (no USB storage device) connected.

If "Print During Capture" is set to "Off", captured images will not print unless the manual "Print Now" command is used (see earlier Image Capture Menu section). If "Print During Capture" is set to "Off", all previous still images are cleared from the print queue once the page is full and the next still is captured. For example, If Prints Per Page is set to 4, still captures 1 through 4 are cleared from the print queue (and may not be printed) once the 5th still is captured.

ⓘ NOTE: The print queue is erased when the CCU is powered off. Be sure that all prints have been printed prior to powering off the CCU. It is recommended to always have a USB storage device connected to ensure captured stills are not "lost" in the event of accidental powering off of the CCU prior to the completion of printing.

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4.8.9 Capturing Video

The image capture system records standard definition video. To capture a video, a USB storage device must be connected to one of the USB outputs on the 3D CCU. Long video captures are stored as separate 12 minute video files in mpeg4 format.

Video capturing may be started and stopped via a specifically programmed button on the 3D TIPCAM1* (further information on the programming of buttons on the TIPCAM1* can be found in the section 'Programming of camera head buttons' in this manual).

Video capturing may also be activated via a keyboard or the buttons on the 3D TIPCAM1* using the 'Capture Video' command in the Image Capture Menu as described above.

When video capturing is activated, a green circle icon will display on the bottom right of the screen. The green circle icon will remain on-screen until video capturing is stopped.

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4.8.10 Ending a Capture Session

To end a capture session, turn off the CCU, remove the 3D TIPCAM1[®] or the USB storage device. Disconnecting a USB printer will not end a capture session.

ⓘ NOTE: The print queue is erased when the CCU is powered off. Be sure that all prints have been printed prior to powering off the CCU. It is recommended to always have a USB storage device connected to ensure captured stills are not “lost” in the event of accidental powering off of the CCU prior to the completion of printing.

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4.8.11 Low USB Memory Warnings

Warning messages will be displayed at the bottom of the video surgical monitor when a connected USB storage device has remaining memory space for less than 5 minutes of new video or less than 100 new still image captures. Both the video and still warning messages will display a countdown until there is no remaining memory space on the connected USB storage device, displaying the message, "USB Memory Full". At that time, no additional video or still images can be stored to the USB Storage device.

ⓘ NOTE: If the USB storage device memory is "Full", and Print During Capture is set to "On", still images may still be captured, but they will be sent to the connected USB printer only; they will not be stored on the USB storage device.

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4.8.12 Reviewing Video and Still Images on a Computer

To view video or captured images on a computer, insert the USB storage device into the computer's USB drive. Go to My Computer to access the USB drive. Open the "KARL STORZ" folder to access the capture session(s) stored on the USB storage device. Please contact KARL STORZ Technical Support at 800-421-0837 for assistance.

The "Temp" file is stored on the USB storage device at the completion of a capture session. This is used by the system to end one capture session and set the system to begin another capture session. The "Temp" file contains no capture session data and may be safely deleted.

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4.9 Image

To change the image settings of the camera, select the desired feature, and scroll through the options using the up and down arrows. When the selection is highlighted, the camera will make the change and the menu will disappear.

Image	
Brightness	[Low]
Enhancement	[Off]
Zoom	[A]
Shutter	[1/60]
Color • • •	
Exit	

Zoom - Allows image size adjustment

Shutter - Auto or selectable shutter speeds from 1/60 to 1/17,000 second.

Color (After white balance, red and blue may be adjusted according to user preference.)

Color	
Red	-1, -2, -3, +1, +2, +3
Blue	-1, -2, -3, +1, +2, +3
Reset	
Exit	

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SCB Display – Image 1

Shutter Off
Brightness Off
Enhancement Off
Devices
Exit

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5 Camera Setup

5.1 SCB

NOTE: SCB data is saved in the CCU and not in the 3D TIPCAM1[®].

When selecting **SCB**, the following sub-menu will appear:

SCB
Exit
SCB Text (On)/Off
Configure SCB Display ►
SCB Warning/Status (On)/Off ►

SCB Text: Turns the on-screen text ON or OFF. This function overrides the 'Configure SCB Display' option but not the 'SCB Warning/Status Display.'

i NOTE: To set camera default for SCB Text, please see "Setup Wizard (CCU)".

Configure SCB Display: Allows selection of up to eight SCB devices/settings for display (i.e. light source, gain, etc). To select the device parameters for display:

1. Highlight the 'Configure SCB Display' option and select by pressing ENTER on the keyboard or the center button on the 3D TIPCAM1[®]. The devices connected to the SCB system will be displayed on-screen, for example:

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3. As an example, to display the Brightness setting on the monitor, select Brightness by highlighting and then pressing on the keyboard or the center button on the 3D TIPCAM1[®]. The cursor will move to the options on the right.

SCB Display (3D CCU)

Brightness	(Off, 1-4)	→	Off
Enhancement	(Off, 1-4)		1
Devices...			2
Exit			3
			4

Use the  or  to scroll through the following list of options: Off, 1, 2, 3, 4, 5, 6, 7, 8. The Off option turns the Brightness display off, and the options 1 – 4 represent the following regions on-screen:

Region 1	Region 2
Region 3	Region 4

4. Use  or  to scroll through the list and select the desired option/region by pressing the ENTER button on the keyboard or the center button on the 3D TIPCAM1[®].

 **NOTE:** If a device parameter is already displayed in the selected region, that parameter will be automatically forwarded to a blank region.

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5. Once selected, the highlighted area will move back to the sub-menu on the left side to continue configuring the display for that device. After completing the display for that device, select Devices... to return to previous screen (SCB Display – Devices) to set display for remaining devices or select Exit to return to live image.

ⓘ NOTE: 'SCB Text On/(Off)' command overrides the text set in this 'Configure SCB Display' submenu.

SCB Warning/Status Display: Turns this on-screen Warning/Status display On or Off.

SCB

Exit
SCB Text (On)/Off
Configure SCB Display ►
SCB Warning / Status (On)/Off

To change this option from its current status (indicated in parenthesis on-screen):

1. Highlight SCB Warning/Status Display.
2. Press ENTER to change status. A message will appear on-screen to confirm your selection.

ⓘ NOTE: "SCB Text On/(Off)" command does not override this function.

3. Once confirmed, the image will return to live.

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6 Maintenance

6.1 Cleaning, disinfection and sterilization

6.1.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for proper handling and care of the camera head.

DO NOT PLUG THE 3D TIPCAM® CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Wipe off any moisture with a lint-free disposable cloth before plugging it in.

PROTECT THE CAMERA SYSTEM FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

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6.1.2 Cleaning instructions

CAUTION: Wear protective gloves, clothing and a face mask for cleaning of contaminated equipment.

Camera control unit (CCU) or other box

WARNING: Always disconnect the power cord before cleaning the camera control unit.

Wipe-clean the exterior surfaces of the unit with a low lint disposable cloth moistened with disinfectant. For wipe-down disinfection, we recommend using only surface disinfectants based on quarternary ammonium compounds (no solutions predominantly containing alcohol such as rapid-action disinfectants) and "disposable disinfectant methods" based on ammonium chloride...

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WARNING: Please read the 3D TIPCAM1® Manual GA 96206259D.

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6.2 Fuse replacement

Switch off the CCU and remove the power plug from the electrical outlet.


CAUTION: Only use fuses approved to IEC 127 Sheet III, Type T for replacement.

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Loosen the line fuse holder with a screwdriver or other suitable tool. Remove the line fuse holder.

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CAUTION: Only use fuses of the correct rating.

	22 204020 100 – 240 V-
Power fuse	2 x T 1.6AL250V

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Reconnect the power cord.
Test the camera system for proper operation.

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6.3 Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspecting the device safety fuses
- Protective ground resistance measured according to IEC 62353 ed 1.0. Limit value: 0.2 ohms
- Earth leakage current measured according to IEC 62353 ed 1.0. Limit value: NC 500 µA, SFC: 1000 µA
- Patient leakage current measured according to IEC 62353 ed 1.0. Limit value: 100 µA (BF), 10 µA (CF)
- Leakage currents must not exceed the limit values specified above.

Test for proper functioning

A function check must be carried out in accordance with the instruction manual. The safety check must be entered in the unit logbook and the results of the check must be documented. If the device does not function reliably and/or safely, it must be repaired.

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6.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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6.5 Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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6.6 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit. In Germany you can refer repairs direct to
KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.
In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important Information

WARNING: In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

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6.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

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KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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7 Technical description

7.1 Troubleshooting

NOTE: The picture quality of the camera is dependent upon optimal light conduction. To ensure the best possible light conduction, make sure that the light entrances and exits of the light cable and 3D TIPCAM1® are always kept clean.

Symptom	Possible causes	Remedy
Camera system will not switch on	No power from power line. Defective power fuse. Power cord is not properly connected to power cord connector	Check that there is electricity to the wall outlet. Exchange fuses. Push power cord firmly into power cord connector.
No picture Exchange fuses.	No connection between video-out of CCU and monitor. Cable connector is not inserted completely into receptacle. Defect in camera electronics. Defective monitor.	Check video cable connection, and if necessary, change cables. Reinsert cable connector into receptacle. Send camera/CCU in for repair. Check monitor.
Dark picture	Too little light due to: – soiled lens – defective light cable broken fibers in video endoscope light cable.	Clean the endoscope. Exchange endoscope.
Blurred picture, streaks, smears, etc.	Telescope or lens of the 3D TIPCAM1® or camera soiled.	Clean with cotton swab and alcohol solution
Inaccurate color reproduction	White balance adjustment faulty. Color adjuster on monitor is not adjusted properly. Defective monitor or camera.	Repeat white balance adjustment. Inspect monitor using instruction manual. Send equipment in for repair.
Color reproduction changes	Endoscope cable broken.	Send equipment in for repair to have camera head cable replaced.

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7.3 Equipment classification

- Type of protection against moisture: Ordinary
- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: CF-type instrument , defibrillation proof
- Mode of operation: Continuous

7.4 Equipment test certificates

Type 22204020 equipment is certified and manufactured in accordance with EN 60601-1, EN 60601-2-18, UL 2601-1, CSA 22.2 No. 601.

CE marked equipment has been tested and found to comply with the EMC limits for Medical Devices for the Medical Device Directive 93/42/EEC (EN 60601-1-2: 2001 Class A).

According to Medical Device Directive (MDD):
Class I medical product.
CE mark in accordance with MDD 93/42/EEC

7.5 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their levels of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
7.6 Software User Rights

1. With respect to KARL STORZ products containing software components, the Customer has a non-exclusive License, which is limited by the meaning of the following terms and conditions, to use the programmed logic, computer programs and/or software supplied by KARL STORZ and incorporated in KARL STORZ products („Software“).
2. This License applies for such products sold to the Customer only to the form delivered to Customer and for the sole purpose of operating such products in accordance with written instructions by KARL STORZ. Use of the Software in connection with other products bought by the Customer from KARL STORZ or from third parties is prohibited.
3. A transfer of this License by the Customer is permissible only if the Customer provides the buyer's name to KARL STORZ, and said buyer agrees to comply with the terms and conditions regarding the usage of the software as set forth in this License. In the case of a transfer, the Customer's License shall terminate.
4. All copyrights to the Software, as well as all modifications, enhancements, and upgrades thereto, will remain the property of KARL STORZ. In particular, the Customer may not duplicate, copy, reverse-engineer, decompile, or disassemble the Software or in any way modify the same. The Customer has no right to create derivatives of the software, and the Customer may not copy, decipher or re-create the source code of the Software. § 69g paragraph 2 of the Copyright remains unaffected by this. (This refers to German Copyright Law)
5. Any modification or enhancement of the Software by the Customer in violation of this Agreement is considered a violation of this License; KARL STORZ immediately shall acquire all property rights and copyrights to such modifications or enhancements. The Customer hereby in advance transfers an exclusive, free, and unlimited license to use such modifications and enhancements to KARL STORZ.
6. KARL STORZ does not recommend use of the Software in connection with other software not supplied by KARL STORZ for the respective products. KARL STORZ is not liable for software malfunctions resulting from the use of software not supplied or recommended by KARL STORZ in connection with products supplied by KARL STORZ.
7. In the event of a failure of the Customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the License will, without termination or any further action by KARL STORZ or any other party, immediately terminate.

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7.7 KARL STORZ Communication Bus (SCB)

General information

The 3D CCU has been designed to incorporate a KARL STORZ Communication Bus (SCB) system. This SCB system is based on a serial communication interface between the KARL STORZ devices. This SCB interface, in conjunction with SCB input/output devices, permits the parameters of the SCB devices connected to be displayed centrally and remotely controlled, even from the sterile area.

3D CCUs are shipped with a 2', 6 pin mini DIN cable, and a "KARL STORZ Communication Bus (SCB) - System Manual" (IM-SCB-XX).

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Warnings and cautions

In addition to the Warnings and Cautions in this manual, refer also to all Warnings and Cautions in the "KARL STORZ Communication Bus (SCB) - System Manual" (IM-SCB-XX). In order to assure the safest, most effective use of this KARL STORZ product, it is essential that you read the entire "KARL STORZ Communication Bus (SCB) - System Manual" before attempting to utilize the SCB.

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Installation and operating instructions

Refer to the "KARL STORZ Communication Bus (SCB) - System Manual" for installation and operation of your 3D CCU camera with the KARL STORZ Communication Bus.

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8 Appendix

8.1 Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The 3D CCU model 222040 20 corresponds to EN/IEC 60601-1-2:2001 [CISPR 11 Class A] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The 3D CCU model 222040 20 is a Group 1 unit (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The 3D CCU model 222040 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The 3D CCU model 222040 20 is intended for use by healthcare professionals only. This is a CISPR 11 Class A medical system (RF emissions) and is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes. In a domestic environment the 3D CCU model 222040 20 may cause EMC interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the 3D CCU model 22203020 or filtering the connection to the public mains network.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the 3D CCU model 222040 20. The cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Cables which have been shown to comply with EN/IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
S-Video [Y/C]	Yes	1.8	No	Connection of S-Video signals
BNC to BNC	Yes	1.8	No	Connection of FBAS or SDI video signals
BNC to RGB Video	Yes	1.8	No	Connection of RGB signals
4 Pin Mini DIN	Yes	1.8	No	Connection of SCB signals
DVI	Yes	3	Yes	Connection of DVI signals
AC Power Cord	No	3	No	Power supply
USB	No	1.8	No	Connection to USB printer or storage device.

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Table 201

Guidance and manufacturer's declaration – electromagnetic emissions

3D CCU model 222040 20 is intended for use in the electromagnetic environment specified below. The customer or user of the 3D CCU model 222040 20 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	3D CCU model 222040 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	3D CCU model 222040 20 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202 Guidance and manufacturer's declaration – electromagnetic immunity			
3D CCU model 22204020 is intended for use in the electromagnetic environment specified below. The customer or user of the 3D CCU model 22204020 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery. If image distortion occurs it may be necessary to further filter the AC input to the 3D CCU model 22204020 by using an external AC line filter.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

3D CCU model 22204020 is intended for use in the electromagnetic environment specified below. The customer or user of the 3D CCU model 22204020 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the 3D CCU model 22204020, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
3D CCU model 22204020**

The 3D CCU model 22204020 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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To be filled out by instrument owner:

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanning Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
68 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg


KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325548, Fax: +65 65323832
E-Mail: info@karlstorz.com.sg
E-Mail: service@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisun Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw


KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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