STORZ

KARL STORZ — ENDOSKOPE

40 1500 20

Stroboskop

STORZ- μ PULSAR
KARL STORZ — ENDOSKOPE

40 1500 20

Stroboscope

STORZ- μ PULSAR
KARL STORZ — ENDOSKOPE

40 1500 20

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**Important information
for users of KARL STORZ
instruments**

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended as an aid to correct usage, connection and operation of the KARL STORZ μ -PULSAR. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.


© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG. Utilization and reproduction by third parties require express authorization by KARL STORZ GmbH & Co. KG. All rights reserved.







Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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**Controls, displays,
connectors
and their uses**

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Displays

- ⑥ Operating mode
(stroboscopy mode or continuous light mode)
- ⑦ Battery status
- ⑨ Current frequency
- ⑩ Saved frequency

Controls

- ② Tensioning lever
- ③ ON button/Mode change
- ④ OFF button

Controls

- ① Endoscope connection
- ⑤ Microphone
- ⑧ Power supply jack

Charger

- ⑪ LED display
- ⑫ Primary adapter connection
- ⑬ Primary adapter

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

External system components

Charger

Using the medically approved charger, the Li-Ion battery of the μ -PULSAR stroboscope can be recharged as quickly as possible.

Moreover, this charger enables the μ -PULSAR stroboscope to be operated with a discharged Li-Ion battery.

The charger has a green LED display (11), which indicates that the output voltage is correct.

Country specific configuration

The charger features interchangeable primary adapters (13), which allow the system to be adapted to the country-specific power supply outlets.

Primary adapter for Europe: 1530290
Primary adapter for UK: 1530390
Primary adapter for US/Japan: 1530490

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-  Akkuzustandsanzeige
-  Stroboskopiemodus
-  Dauerlichtmodus

Symbols employed

-  Read the instructions carefully before operating the equipment!
-  Applied part type BF
-  OFF button (No system decoupling when charger is connected to the power supply!)
- Note:** To decouple from the mains supply, pull the charger out of the power supply outlet.*
-  ON button/Mode button
-  Connection for the power supply secondary connector
-  Indication of the power supply to be used
-  Tensioning lever 'open'
-  Tensioning lever 'locked/closed'
-  Electronic information product pollution control (China RoHS)

Display symbols

-  Battery status indicator
-  Stroboscopy mode
-  Continuous light mode

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Garantiekarte

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Photographs of the unit	IV
Controls, displays, connectors and their uses	V
Symbols employed	VII

General information

Device description	2
--------------------	---

Safety instructions

Warnings and cautions	3
Intended use	6
User qualification	6
Safety precautions at the site of installation	7
Safety precautions when operating the unit	7
Battery safety instructions	8

Usage and operating instructions

Unpacking the equipment	11
Basic equipment	11
Using and connecting up the unit	11
Charging the battery	12
Operating the unit	13
Description of operation	14
Initializing the unit	14
Automatic shut-down	14
Stroboscopy mode	14
Slow motion function	15
Continuous light mode	15
Video function	16

Maintenance

Changing the battery	17
Cleaning	17
Maintenance	17
Servicing and repair	17
Disposal	18
Repair program	18
Important information	18
Limitation of liability	19
Manufacturer's warranty	19

Technical description

Technical data	20
Technical documentation	21
Block diagram	22

Appendix

Cleaning agents and disinfectants	23
Electromagnetic Compatibility (EMC)	24
Information	24

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Device description

The µ-PULSAR stroboscope is a mobile device that generates light with an LED lamp for stroboscopic examination (flashing light) and also for normal viewing (continuous light). To be more precise, the continuous light is actually a pulsating (flashing) light that, due to its high frequency (3000 Hz), is perceived as continuous light. In spite of this minor technical difference, for the purpose of simplicity, conventional terminology will be used (flashing light/continuous light) in this manual.

Thanks to its compact and light-weight design, without any cumbersome cables, the µ-PULSAR stroboscope can be directly and safely connected to a laryngoscope via the quick-release coupling, thus allowing it be used straight away.

The light source is provided by a highly efficient LED, which offers numerous advantages:

- As the service life is almost unlimited, the lamp does not have to be changed
- Completely silent operation
- Bright white light without any infrared or UV light.

The unit features a modern Li-Ion battery, and can thus be operated independently of the mains supply.

The µ-PULSAR stroboscope can be quickly recharged using the enclosed charger (401500 30). Once recharged the unit can be subsequently continuously operated for approx. 1 hour.

If examinations are planned when the battery is discharged, work can be continued straight away by connecting the charger to the mains supply*.

***Note:** The battery is not simultaneously charged during this time.

Note: All three country-specific* adapters are included in the charger's scope of supply.

* Primary adapter for Europe: 1530290
Primary adapter for UK: 1530390
Primary adapter for US/Japan: 1530490

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Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device.

Definitions

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

Note: A Note indicates special information about operating the unit.

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Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Warning: The electrical installations at the site of installation must comply with the requirements of the applicable national standards.

Warning: The unit may NOT be used in operating rooms in which combustible mixtures of anesthetics and air, oxygen or laughing gas can form. Risk of explosion if used in the presence of flammable anesthetics.

Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: Danger of glare! Do not look directly into the light outlet when the unit is switched on.

Warning: To avoid the risk of electrical shock, do not remove cover. Refer servicing to the manufacturer or personnel approved by the manufacturer. Removal of covers by nonauthorized personnel will void the equipment's warranty.

Warning: Patient leakage currents can accumulate if powered accessories are used.

Warning: Do not touch the cable end of the charger (secondary connector) and the patient simultaneously or rest the cable end on the patient. Always connect or disconnect the cable from the μ -Pulsar outside of the patient area!

Warning: Risk of burns at high ambient temperature! The surface of the μ -PULSAR can reach temperatures in excess of 41 °C.

Warning: Risk of burns! Higher temperatures can occur around the light outlet window of the endoscope.

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Caution: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.
Caution: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

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Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Note: This unit has been tested and complies with EMC limit values as per the Medical Device Directive 93/42/EEC. Please read the Electromagnetic Compatibility Information (see pages 24 ff.).

Note: This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Note: Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The μ -PULSAR stroboscope is designed for larynx diagnostics in both hospitals and doctors' practices.

The use of this instrument in any way other than indicated above is not permitted for safety reasons.

Video function

The μ -PULSAR stroboscope is only conditionally suitable for operation together with a video camera. Since the stroboscope and camera system are not synchronized, the video image may be subject to flickering.

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The μ -PULSAR stroboscope may only be operated with accessories which have been designated by KARL STORZ as suitable for the unit. Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
User qualification

The μ -PULSAR stroboscope may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

Warning: The unit may NOT be used in operating rooms.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

Warning: The unit is exclusively designed for indoor use and must not come into contact with water or dust.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before use.

When using the μ -PULSAR stroboscope, the patient must be treated and kept under observation with the usual medical care.

Any treatment may only be performed if there is visual observation of the action of the unit/instrument.

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Safety Instructions
Battery


The μ -PULSAR stroboscope contains a modern Li-Ion battery. Before using the unit, please read the following safety instructions on use of the Li-Ion battery.

- Do not disassemble or modify the battery pack. The battery pack is equipped with built-in safety/protect-ion features. Should these features be disabled, the battery pack can leak acid, over-heat, emit smoke, burst and/or ignite.
- Do not use or leave battery near a heat source such as fire or a heater (80°C or higher). If the resin separator should be damaged owing to overheating, internal short-circuiting may occur to the battery pack, possibly leading to acid leakage, overheating, smoke emission, bursting and/or ignition or the battery pack.
- Do not immerse the battery pack in water or sea-water and do not allow it to get wet. Otherwise, the protective features in it can be damaged, it could be charged with extremely high current and voltage, abnormal chemical reactions may occur in it, possibly leading to acid leakage, overheating, smoke emission, bursting and/or ignition.
- Do not recharge the battery pack near a naked flame or in extremely hot weather. High temperature can trigger its built-in protective features, inhibiting recharging, or can damage the built-in protective features, causing it to be charged with an extremely high current and voltage, and, as a result, abnormal chemical reactions can occur in it, possibly leading to acid leakage, overheating, smoke emission, bursting and/or ignition.

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- The guaranteed recharging temperature range is 0 to 40°C. A recharging operation outside this temperature range can lead to acid leakage and/or overheating of the battery pack and may cause damage to it.
- If acid leaking from the battery pack comes into contact with your skin or clothing, wash it off with running water immediately. Otherwise, skin inflammation can occur.
- Before use, carefully study the Operation Manual and Precautions. Also, keep the manual safe for future reference.
- For recharging procedures, refer to the Operation Manual of your battery pack charger.
- If you find rust, a bad odor, overheating and/or other irregularities when using the battery pack for the first time, return it to your supplier or vendor.

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Unpacking the equipment

Carefully remove the μ -PULSAR stroboscope and accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Using and connecting up the unit

Note: The stroboscope μ -PULSAR including the accessories connected may be used only in medical rooms with electrical installations conforming to the applicable national, state, and local electrical regulations.

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Country-specific configuration

The charger features interchangeable primary adapters, which allow the system to be adapted to the country-specific power supply outlets. Primary adapter for Europe: 1530290 Primary adapter for UK: 1530390 Primary adapter for US/Japan: 1530490 To connect the primary adapter, push it as far as it will go into the charger's connection point.

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

Warning: Do not operate the unit within the hazard zone shown in the diagram while explosive, anesthetic gases are in use.

Note: To decouple from the mains supply, pull the charger out of the power supply outlet.

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Charging the battery

Before using for the first time, fully charge the internal Li-Ion battery once!

Warning: Do not touch the cable end of the charger (secondary connector) and the patient simultaneously or rest the cable end on the patient. Always connect or disconnect the cable from the μ-Pulsar outside of the patient area!

To do this, firstly insert the charger's secondary adapter into the power supply jack ⑧ on the bottom of the μ-PULSAR stroboscope.

Then connect the charger to the power supply outlet. The LED display ⑪ lights up green.

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Ensure that the stroboscope is not operating by pressing the OFF button ④.

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Use the battery status indicator ⑦ to check that the battery is being charged (the inside of the battery symbol is "filled" repeatedly).

Once the battery is completely charged, a full battery symbol is shown. As a rule, this takes approx. 2 hours (depending on the battery status at the start of charging).

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Operating the unit

Only operate the charger using AC voltage 100...240 VAC, 50-60 Hz.

Caution: Only use the unit in line operation with the charger specified by KARL STORZ!

Ensure that the μ-PULSAR stroboscope is switched off by pressing the OFF button ④. This will prevent you (and others) from being dazzled by the LED.

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
Connect the μ-PULSAR stroboscope to the laryngoscope by firstly turning the tensioning lever ② to the left as far as it will go and then inserting the unit into the endoscope connection ① as far as it will go via the light connection/Storz connection.

Note: When the tensioning lever is 'open', an 'open' lock is displayed (see arrow).

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Now align the microphone ⑤ to the distal end of the laryngoscope or the patient's mouth.

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
Secure the connection by turning the tensioning lever ② to the right as far as it will go. Check that the connection is correctly locked in place. The μ-PULSAR stroboscope is now securely connected to the laryngoscope and can be operated by pressing the ON button/mode button ③.

Note: When the tensioning lever is 'locked' (closed), a closed lock is displayed (see arrow).

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Description of operation

Initializing the unit

On pressing the ON button ③, the unit is switched on and is initialized. The unit name and firmware version are briefly shown on the display. The unit is now in stroboscopy mode. Pressing the ON button ③ again serves to switch between continuous light mode and stroboscopy mode.

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Automatic shut-down

The μ-PULSAR stroboscope has an automatic shut-down feature, which shuts down the unit 10 minutes after the last phonation.

Note: This feature only works in stroboscopy mode. In continuous light mode, the unit remains switched on until the battery is fully discharged. Even in stroboscopy mode, you should ensure that the unit is switched off following the examination using the OFF button ④.

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Stroboscopy mode


In stroboscopy mode, the μ-PULSAR stroboscope is triggered by audio sources (speech, patient tones), which are recorded by the inbuilt microphone. The trigger feature, which works via modern digital modules, enables the light pulses to be controlled in such a way that movement of the patient's vocal cords can be clearly and easily monitored. The use of an LED together with a specially developed telescope allows a brightness to be achieved the like of which has not been seen before with LEDs.

In display position ⑥ the symbol for stroboscopy mode (below the battery symbol) is shown. The currently recorded frequency (right) is shown in the first line of the display – ⑨. This display always shows the currently recorded pitch.

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If an output frequency remains constant over a period of time, this is shown in the second line of the display (right) – ⑩. This value is saved until there is a new stable pitch. This makes it possible to still read off and document the frequency after the examination.

Slow motion function

The μ-PULSAR stroboscope is set to a slow motion frequency of approx. 1 Hz. This allows easy monitoring of vocal cord movement over their entire range.


Note: The slow motion frequency cannot be altered.

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Continuous light mode

In the continuous light mode, the μ-PULSAR stroboscope is not triggered by an audio signal but rather outputs continuous light with a set frequency.

Due to the extremely high frequency of the light pulses (3000 Hz), the output light is perceived as continuous light.

Note: In continuous light mode, the automatic shut-down function is not activated and it must therefore be ensured that the unit is switched off using the OFF button ④ following examinations, so as not to waste the battery.

The continuous light mode is shown by the continuous light symbol in display position ⑥. No frequency is shown.

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Video function

The μ-PULSAR stroboscope is only conditionally suitable for operation together with a video camera. Since the stroboscope and camera system are not synchronized, the video image may be subject to flickering.

In order to achieve the best possible results, we recommend that examination documentation only be performed with our TELECAM® camera system, C-Mount camera head and C-Mount lens (f=12 mm).

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Warning: Risk of burns at high ambient temperature!
The surface of the μ-PULSAR can reach temperatures in excess of 41 °C.

Warning: Risk of burns!
Higher temperatures can occur around the light outlet windows of the endoscope.

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Changing the battery

For safety reasons, the battery inside the unit (typical service life of approx. 2 years) may only be changed by service staff authorized by KARL STORZ. Further information is available from our service hotline (Tel. +49 (0)7461 708-490 / Fax: +49 (0)7461 708-404).

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
Cleaning

Warning: Always pull out power plug before cleaning!
Caution: Avoid allowing fluids to enter the unit.

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The exterior of the unit can be wipe-down disinfected using a disposable cloth moistened with alcohol-free disinfectant. Reclean the light outlet opening with a cotton tip applicator soaked in 70% isopropyl alcohol.

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Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit. In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 TUTTLINGEN, GERMANY.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender. Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations in the room where the unit is connected and operated complies with the applicable laws and standards, and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty*

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or work-manship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 TUTTLINGEN, GERMANY

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 TUTTLINGEN, GERMANY

***Note:** The battery is not included in the warranty as this is a wearing part (cp. page 17).

Opening the equipment or performance of any repairs or modifications to the equipment by unau-thorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modi-fication performed during the warranty period shall void all warranty.

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- Type of protection against electric shocks: Protection Class II
- Degree of protection against electric shocks: Applied part of type BF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the Appendix (p. 24-36).

Directive compliance (for 40 1500 20)
According to Medical Device Directive (MDD):

This medical product belongs to Class I

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated as reparable by their respective manufacturers.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
Type

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
A = Not suitable for flexible endoscopes
(according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets
Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Pulver Classic	1, 2, 3
	Sekusept Pulver Classic + Aktivator	1, 2, 3
	Ultrasemit AF	1, 2, 3
Esteer Pharma GmbH	Ultrassept Aktiv	1, 2, 3
	Ultrassept Classic	1, 2, 3
	Afid	1, 2
Fresenius AG	Afid plus (neu)	1, 2
	Sporcid FF	1, 2
Indeba Industria E Comercio LTDA	Letahdeldio	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
Lonza AG	ID 50	2
	Aldasan 2000	1, 2
Lysoform Dr. Rosemann GmbH	Almyrol	1, 2
	Desoform	1, 2
	Lysoformin 3000	1, 2
Medichem International	MedDis	2
	MediZyme	1
Merz Consumer Care GmbH	Edisonite Super	3
	Mucadont-IS	2, 3
	Mucadont-Zymaktiv	1, 3
	Mucocit-T	1, 2, 3
Orochemie	A 10	2
	ID 210	2
Promagent AB	Wavacide	2
	Gigasept	2, 3
	Gigasept FF	2, 3
Schülke & Mayr GmbH	Gigasept Instru AF	1, 2, 3, A
	Lysetol V	2
Schumacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1
Weigert, Dr. GmbH & Co.	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
Whiteley Industries PTY Ltd	Aidal Plus	2

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination shaft, **not** the housing/handle part.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix during installation and commissioning.

The μ -PULSAR model 401500 20 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The μ -PULSAR model 401500 20 is a Group 1 unit (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the μ -PULSAR, model 401500 20. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200

Accessories which have been shown to comply with EN/IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
Charger 401500 30	No	2	No	For charging the internal battery

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Table 201

Guidance and manufacturer's declaration – electromagnetic emissions

The μ -PULSAR model 40 1500 20 is intended for use in the electromagnetic environment specified below. The customer or user of the μ -PULSAR model 40 1500 20 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	THE μ -PULSAR model 40 1500 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The μ -PULSAR model 40 1500 20 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202

Guidance and manufacturer's declaration – electromagnetic immunity

μ-PULSAR model 40 1500 20 is intended for use in the electromagnetic environment specified below. The customer or user of the μ-PULSAR model 40 1500 20 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact air (not applicable)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 204

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

μ-PULSAR model **40 1500 20** is intended for use in the electromagnetic environment specified below. The customer or user of the μ-PULSAR model **40 1500 20** should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the μ-PULSAR model 40 1500 20, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
μ-PULSAR model 40 1500 20**

The μ-PULSAR model 40 1500 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

To be completed by the supplier/importer:

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To be filled out by instrument owner:

Return address/company stamp:

Please attach
sufficient
postage

Field of application:

Type of instrument:

Serial no.:

Purchase date:

Signature/Date:

REPLY CARD

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

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Absender/Firmenstempel:

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Geräte-Typ:

Serien-Nr.:

Kauf-Dat.:

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ANTWORTKARTE

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0,
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 (0)30 30 69090
Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
2345 Argenta Road, Suite 100
Mississauga, ON, L5N 8K4, Canada
Phone: +1 905 816-8100, Fax: +1 905 858-0933
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, +1 800 421-0837
Fax: +1 424 218-8526
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, Mexico
Phone: +52 5552 5056 07, Fax: +52 55 5545 0174

KARL STORZ Endoscopia Argentina S.A.
Cervino 4449 Piso 10°
1425 Buenos Aires C. F., Argentina
Phone: +54 11 4772 4545, Fax: +54 11 4772 4433
E-Mail: info@karlstorz.com.ar

KARL STORZ ENDOSKOPI NORGE AS
P.O.Box 153
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14, 12739 Skärholmen, Sweden
Postal address: P.O. Box 8013, 14108 Kungens Kurva,
Sweden
Phone: +46 8 50 56 480, Fax: +46 8 50 56 4848
E-Mail: kundservice@karl-storz.se

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopie Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopie France S. A.
12, rue Georges Guynemer
Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Francia - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 071, Fax: +385 1 6406 077
mail: b.vlajo@karlstorz.hr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868 Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030
**Sales for Industrial Endoscopy

KARL STORZ Endoscopia ROMANIA srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

TOV KARL STORZ Ukraine
18b Geroev Stalingrada avenu
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19, -20
Fax: +380 44 42668-14

000 KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenicia Street, Mina Elhosn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax: +961 1 365151

KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kasachstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@it.te.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51, Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 - 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

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