Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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**1. Important information for
users of KARL STORZ
instruments**

**It is recommended to check the suitability of
the product for the intended procedure prior
to use.**

Thank you for your expression of confidence in the
KARL STORZ brand name. Like all of our other
products, this product is the result of years of expe-
rience and great care in manufacture. You and your
organization have decided in favor of a modern,
high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an
aid in the proper setup, installation, and operation
of the cold light fountain CO₂mbi LED SCB. All
essential details of the equipment and all actions
required on your part are clearly presented and
explained. We thus ask that you read this manual
carefully before proceeding to work with the equip-
ment. Insert this manual in its protective wallet and
keep it available for ready reference in a convenient
and conspicuous location near the equipment.

KARL STORZ is constantly working on the further
development of all products. Please appreciate that
changes to the scope of supply in form, equipment
and technology may occur for this reason. There-
fore, no claims may be derived from the informa-
tion, illustrations, and descriptions in this manual.

© All product pictures, photos and descriptions and
texts are the intellectual property of KARL STORZ
GmbH & Co. KG

Utilization and reproduction by third parties require
express authorization by KARL STORZ GmbH &
Co. KG

All rights reserved.

CAUTION: Federal (USA) law restricts this device
to sale by or on the order of a physician.

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3. Controls, displays, connectors, and their uses

- ① Power switch
- ② TFT touch screen
- ③ Connection for light cable adaptor and insufflation (video endoscopes for gastroenterology)
- ④ Gas connection (standard American type)
- ⑤ Air filter
- ⑥ SCB connectors*
- ⑦ Potential equalization connector
- ⑧ Power cord socket
- ⑨ Line fuses (fuse holder)
- ⑩ Mount for the CO₂ bottle holder (optional)

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* The KARL STORZ-SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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4. Symbols employed

4.1 Symbols on the device

	Follow instructions for use
	ON
○	OFF
⇩	Potential equalization connector
	Applied part of type CF
~	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

Symbols on label and packaging of the device

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. This can be downloaded from www.karlstorz.com.

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

	Menu ①
	Start/Stop (Off) ② or ⑤
	Start/Stop (On) ② or ⑤
	Minus ③
	Plus ④
	A Low ⑥
	A High ⑦
	Medium Air ⑧ in CO ₂ /Air Mode
	Medium 'CO ₂ /Air' ⑨ in CO ₂ /Air Mode or CO ₂ /Air Mode
	CO ₂ 100 % Mode
	CO ₂ available ⑩ (high pressure supply)
	CO ₂ supply low (high pressure supply)
	CO ₂ supply empty (high pressure supply)
	House connection ⑬ (low pressure supply)
	Automatic light regulation (on) ⑪
	Reject
	Accept
	Remove
	Back
	Continue
	Audio signals off

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8.5	Disposal	34
8.6	Repair program	35
8.7	Limitation of liability	36
8.8	Manufacturer's warranty	36
9.	Technical description	37
9.1	Information signals	37
9.1.1	Configuration	37
9.1.2	Visual signals	37
9.1.3	Acoustic signals	38
9.1.4	Information signals (general)	38
9.1.5	Availability signal	39
9.2	Troubleshooting	40
9.3	Technical data	41
9.4	Technical documentation	42
10.	Spare parts, recommended accessories	43
10.1	List of Spare parts	43
10.2	Recommended accessories	44
11.	Electromagnetic Compatibility (EMC) Information	45
12.	Subsidiaries	59

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5. General information

5.1 Device description

The CO₂mbi LED cold light fountain combines a high-performance LED cold light fountain and an insufflation pump for diagnostic and therapeutic endoscopic applications.

The CO₂mbi LED cold light fountain is compatible with all gastroenterological KARL STORZ videoscopes. In the configuration with a light source adaptor (Art. no. 495KS), it can be used with several KARL STORZ light cables.

The device is specially designed for use during gastroenterological interventions, and offers a wide range of functions. These include:

- Efficiency: High light intensity at with low power consumption
- Comfortable working environment: Intelligent cooling concept for minimal noise emission
- Durable: No costs or effort required to change the lamp if used properly
- Unique: Insufflation can be set using either air or CO₂
- Intelligent: CO₂ is only supplied when the user actuates the Insufflation button via the videoscope
- Flexible: There is the option of switching the CO₂ insufflation to pure CO₂ supply
- Clear: Direct operation via a color screen with touch-sensitive surface (touchscreen)
- Integration: The SCB function enables the light source to be controlled remotely and the system status information to be retrieved, while providing simple configuration during initialization
- Up-to-date: Software updates can be installed easily and conveniently via an interface

5.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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6. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly before putting the device into operation. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

NOTE: Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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WARNING: The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep device out of reach of patients.

WARNING: For safety reasons, do not simultaneously touch the device output sockets and the patient.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 or Chap. 16, IEC 60601-1:2005).

WARNING: Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

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WARNING: At the end of its useful operating life, dispose of the unit as electronic scrap.

WARNING: Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

WARNING: Ensure sufficient air circulation for the device.

WARNING: Pressing more than one touch point simultaneously can lead to incorrect settings.

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WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

WARNING: Only use fuses of the correct rating. Replace defective fuses and screw the line fuse holder back on.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: Test this device prior to each procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

WARNING: Do not use the device if the input controls or the displays are defective.

WARNING: The device must not be used in oxygenated environments.

WARNING: The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind.

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WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

WARNING: Always use the lowest possible settings for your intervention.

WARNING: Danger of glare! Never look into the light outlet when the light system is switched on.

WARNING: Insufficient ventilation of the device can cause an internal build-up of heat which can lead to a safety switch-off of the lamp.

WARNING: The light source transmits high energy light which leads to high temperatures at the light window of the endoscope. Tissue damage is possible.

WARNING: To prevent water from flowing back into the device, the water bottle must be placed in the holder intended for this purpose. The filled water bottle must not be turned upside down or positioned higher than the device.

WARNING: Only use sterile CO₂.

WARNING: The correct configuration of the device must be checked after each software update.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

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6.2 Intended use

Cold light fountains with insufflation pumps are used for illumination in diagnostic and therapeutic applications in endoscopy. The integrated insufflation pump is specially designed for use with flexible endoscopes for insufflation and irrigation as well as dilation in open systems for gastroendoscopic applications.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

Indication

Illumination in diagnostic and therapeutic applications in endoscopy for use in insufflation pumps for insufflation and irrigation as well as dilation.

Contraindication

Use is contraindicated, if, in the opinion of the attending physician, the general condition of the patient does not permit the use of the device.

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6.3 User qualification

The cold light fountain CO₂mbi LED may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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6.4 Training in the operation and function of the device

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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6.5 Patient profile

The use of the CO₂mbi LED is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

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6.6 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been thoroughly trained at least once in the use of the device
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

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6.7 Intended conditions of use

6.7.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

6.7.2 Other intended conditions

Frequency of use: one or more times a day
Length of use: from a few minutes to several hours a day.
Place of installation: Positioning on a level, vibration-free surface
Mobility: can be moved if positioned on a cart
Combination: can be used on the patient at the same time as other equipment required for the operations
Control: can be controlled via the KARL STORZ-SCB.

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
6.8 User position

For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30 – 70 cm to the front panel.

6.9 Safety precautions at the site of installation

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations. It is not intended for use in hazardous zones.

! WARNING: The device may not be operated in the demarcated hazard zone wherever highly inflammable inhalation anesthetics and their mixtures are used.

This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

! WARNING: Ensure that sufficient air circulation is guaranteed for the device.

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6.10 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operates properly before using it. During treatment using the CO₂mbi LED cold light fountain, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

! WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

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6.11 Safety features

The cold light fountain CO₂mbi LED is equipped with the following safety features:

6.11.1 Self-test

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

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- Adapter 495KS (optional)
- KARL STORZ Lichtkabel 495NE (optional)

7. Installation and operating instructions

7.1 Unpacking the equipment

Carefully unpack the cold light fountain CO₂mbi LED and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

7.2 Basic equipment

KARL STORZ Set No.: TL100 S1

- 1 CO₂mbi LED cold light fountain TL100
- 1 Bottle holder TL001
- 1 Bottle ring for KARL STORZ irrigation bottles TL002
- 1 Water bottle 13992BS (incl. lid)
- 1 400 A power cord, length 300 cm
- 1 SCB connecting cable 20 0901 70, length 100 cm
- 1 Universal wrench 20 4000 30
- 1 Instruction manual 96206538 DER

Optional:

- 1 Bottle ring for Ampuwa Plastipur irrigation bottles TL003
- 1 CO₂ bottle holder TL004

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7.3 Installation

NOTE: The cold light fountain CO₂mbi LED including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

WARNING: The device is not intended for use in hazardous zones. Do not operate the device within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Installation and operating instructions

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Set the device on a horizontal, flat surface.

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7.3.1 Connecting the ground line


The device is equipped with a connector ⑦ for attaching a ground line.

The device's ground line should be installed by a qualified electrician.

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7.3.2 Connecting the power cord

Connect power cord and insert the power plug fully into the power cord receptacle ⑧.

! WARNING: Only operate the device with the voltage stated on the manufacturer's identification plate.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

i NOTE: Ensure that access to the power plug is ensured at the site of installation.

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
7.3.3 KARL STORZ-SCB


① NOTE: To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑥. Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ-SCB control NEO System Instruction Manual).

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7.3.4 Connecting the CO₂ bottle


① NOTE: The CO₂ bottle must be connected before switching on the device. If the CO₂ bottle is only connected during operation, it will no longer be possible to switch to CO₂ mode, as the CO₂ channel can no longer be tested.

Connect the gas outlet of the CO₂ bottle to the device's gas connection ④ using the high pressure tube supplied.

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7.3.5 CO₂ bottles equipped with German-standard or ISO connection

Place the high pressure tube (see the section 'Accessories') onto the outlet of the CO₂ bottle and screw tight using the universal wrench provided.

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7.3.6 CO₂ bottles equipped with PIN-Index connection

Place the connecting piece of the high pressure tube (see the section 'Accessories') onto the CO₂ bottle and screw this tight.

Use the wrench provided to open the CO₂ bottle.

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
7.3.7 Opening the valve of the CO₂ bottle

Slowly open the valve stopcock on the CO₂ bottle by turning it approx. one-half turn counterclockwise.

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
7.3.8 Connecting to the central gas supply

Connect the low pressure tube (see the section 'Accessories') to the device's gas connection ④ and connect the other end directly to the wall outlet of the central gas supply (3.3...7 bar).

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7.3.9 Attaching the CO₂ bottle to the device (optional)

To attach the optionally available CO₂ bottle holder TL004 (see the section 'Accessories'), position the device vertically on its right-hand side (viewed from the front).

Remove the rear screw from the bottom of the device using the Allen wrench with star-shaped head supplied.

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
First attach the bearer plate for the actual CO₂ bottle holder (see images on the left). Use the enclosed screws of length 6 mm.

Attach the bearer plate using two screws on the bottom of the device and one screw on the back of the device.

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Screw the bottle holder onto the bearer plate. Use the enclosed screws of length 6 mm.
Insert the CO₂ bottle in the holder vertically and secure using the locking clip.

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7.3.10 Attaching the bottle holder TL001 to the device

Attach the water bottle holder (see the section 'Accessories') to the side of the device (left-hand side of the device when viewed from the front). Use the enclosed screws of length 6 mm.

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Insert the water bottle for irrigation/insufflation in the bottle holder.

NOTE: The bottle holder is available in two different versions (see the section 'Accessories'). The standard size (TL002) can be used for KARL STORZ water bottles, and the second size (TL003) can be used for disposable water bottles.

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To fill the water bottle, remove it from the bottle holder.

The water bottle should not be filled in the vicinity of electrical devices.

NOTE: The water bottle should not be filled to more than 2/3 of its capacity. The remaining third is required to cushion air pressure for efficient irrigation/insufflation. For further information on irrigation and insufflation, please refer to the instruction manual of the videoscope to be used.

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
The supply tube for the water bottle can be hung up in the designated holder.

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7.3.11 Connecting the videoscope

Hold the videoscope firmly by the handle and insert the connector (light and insufflation) into the light outlet point ③ as far as it will go.

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
Insert the supply tube, which is attached to the water bottle lid, into the endoscope supply connector.

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
7.3.12 Connecting KARL STORZ standard endoscopes

In combination with a KARL STORZ light cable, KARL STORZ standard endoscopes can be connected to the CO₂mbi LED light source using the 495KS adaptor. This adaptor is connected to the light outlet point in place of the video endoscope.

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- NOTE:** The light cables 495NCS, 495NCSC, 495NAC, 495TIP and 495NB cannot be connected to the CO₂mbi LED light source using the 495KS adaptor.
- NOTE:** If the CO₂mbi LED is only used for illumination purposes in combination with the 495KS adaptor, only the light intensity should be set. For additional information, see Section 7.4.1.

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General applications

The CO₂mbi LED is suitable for use in a wide range of applications. The correct light intensity depends on the particular application, the instrument being used, and the connected light cable.

If the camera control unit used is not SCB-compatible, the intensity should always be kept to a minimum.

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7.4 Commissioning the CO₂mbi LED light source

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
Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo.

The device is now running a self test. Once the test has been successfully completed, an acoustic availability signal sounds and a touch screen appears.

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
Confirm the message to complete the startup process.

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If the self test fails for critical components, the device goes into 'Malfunction' mode (Safe State) and an acoustic information signal sounds. Verify that the device is properly connected and that the preparations have been made according to the instructions.

For further information on troubleshooting, please refer to the section 'Technical description – Troubleshooting'.

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In routine operation, all functions for the use of the CO₂mbi LED light source can be operated and parameters can be adjusted via the menu ①.


NOTE: The Settings menu can be called up when the light and insufflation are switched off (the icon on the Start/Stop button is 'blue').

Press button ① to switch to the Settings menu. Further details can be found in the section 'Settings'.

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7.4.1 Light

Using the automatic light source control

If the CO₂mbi LED light source is connected to the camera control unit via an SCB connection cable, automatic light source control can be used. Confirm the information message via the camera control unit. The 'Automatic Light Regulation (On)' icon ⑪ will appear next to the intensity display.

Pressing the Start/Stop button ② switches the light on (the icon on the Start/Stop button ② is orange). The light intensity is then automatically controlled via the camera control unit, ensuring optimal illumination in every situation.

① NOTE: The automatic light regulation is deactivated as soon as the button ③ or ④ is pressed during operation. The automatic light source control can be switched on again via the camera control unit. Further information can be found in the instruction manual IMAGE1 SPIES™.

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
Setting the brightness manually without an SCB connection to the camera control unit

The CO₂mbi LED starts automatically when the light and insufflation are switched off.

The brightness can be modified using the buttons ③ and ④, both when the light is switched on and off.

Press the button ④ to increase the brightness.

Press the button ③ to reduce the brightness.

The current set brightness is shown in % on the display ⑫.

① NOTE: The orange triangles indicate the set value in relation to the possible setting range.

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Switching the light on and off


The light is switched on and off using the Start/Stop button ②. The current status of the light (on/off) is indicated by the color of the icon on the Start/Stop button (orange = on, blue = off).

① NOTE: The intensity can be set between 5 % and 100 %. This can either be done in stages by briefly pressing (and releasing) the button ④ or ③, or continuously by pressing the corresponding button for an extended period. Pressing and holding the button ④ or ③ changes the value in the corresponding direction until the maximum or minimum setting is reached.

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7.4.2 Insufflation


Switching Insufflation on and off

The internal pump for insufflation/irrigation is switched on and off using the Start/Stop button ⑤. The current status of insufflation (on/off) is indicated by the color of the icon on the Start/Stop button (orange = on, blue = off).

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Setting the insufflation level (flow)

The flow can be set to one of two levels. Pressing the button ⑥ sets the lowest flow level. Pressing the button ⑦ sets the highest flow level.

The flow levels (in CO₂/Air mode and CO₂ 100 % mode) and the insufflation medium (only in CO₂/Air mode) can be changed, both when insufflation is switched on or off.

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NOTE: The selected insufflation level may be affected by system resistance (instrument, tubes, connections), the included volumes (water bottle, gastrointestinal tract) and losses resulting from any leaks.

NOTE: The water bottle should not be filled to more than $\frac{2}{3}$ of its capacity.
The remaining third is required to cushion air pressure for efficient irrigation/insufflation.
For further information on irrigation and insufflation, please refer to the instruction manual of the videoscope to be used.

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7.4.3 CO₂ Insufflation modes

The CO₂mbi LED light source can be operated in one of two CO₂ insufflation modes. The device can be operated in CO₂/Air mode or CO₂ 100 % mode.

CO₂/Air mode

In this mode, CO₂ is only supplied when the user actuates the Insufflation button via the videoscope. If insufflation is not required, the light source automatically switches to air supply. This reduces the CO₂ consumption.

NOTE: In CO₂/Air mode, a very low air content may be included in the insufflation, depending on the system.

In CO₂/Air mode, it is possible to switch between the insufflation media CO₂ and air (see Section 7.4.5).

CO₂ 100 % mode

In this mode, CO₂ is supplied constantly. The supplied amount can be controlled by adjusting the flow level.

NOTE: If there is sufficient CO₂ supply pressure when the device has been started, the CO₂mbi LED light source sets the CO₂ mode and the insufflation medium CO₂ automatically as the standard insufflation settings.

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7.4.4 Selecting the CO₂ mode

Click on the 'Menu' icon ① to open the Settings menu (see page 22 and Section 7.5 'Settings').

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ⓘ NOTE: If CO₂ 100 % mode is selected, air insufflation is not available.

7.4.5 Selecting the insufflation medium in CO₂/Air mode

In CO₂/Air mode, the CO₂mbi LED light source can supply either CO₂ or air as the insufflation medium. The insufflation medium can be changed by pressing the ⑧ or ⑨ button.

The medium can be changed regardless of the flow level (⑥ or ⑦) or the position of the Start/Stop button ⑤ (on/off).

If there is sufficient pressure for the CO₂ supply, the CO₂mbi LED light source automatically sets CO₂ insufflation as the standard insufflation medium after the device has been started.

ⓘ NOTE: If the CO₂ supply cannot guarantee the minimum pressure when the device has been started, only the air insufflation button ⑧ is displayed, and the CO₂ insufflation button ⑨ is hidden.

A more detailed description of the device behavior is provided in Sections 7.4.8 and 7.4.9.

The icon ⑩ indicates the currently selected supply system for CO₂ insufflation. The gas supply can be operated via the central supply (low pressure) or via a CO₂ compressed air bottle (high pressure).

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7.4.6 CO₂ gas supply in CO₂/Air mode

If CO₂/Air mode ⑨ is selected, intelligent CO₂ supply is activated automatically. This CO₂ regulating mechanism optimizes CO₂ consumption during use of the CO₂mbi LED light source.

The CO₂ supply is only activated if the user actuates the insufflation valve on the handle of the video-scope. The light source then switches automatically from air supply to the selected medium CO₂. Once the insufflation valve has been actuated, the CO₂mbi LED starts to continuously supply CO₂. If the insufflation valve on the handle of the video endoscope is not actuated after a certain time delay, the CO₂mbi LED automatically switches the insufflation medium from CO₂ back to air.

This efficiently reduces CO₂ consumption to a minimum during the examination.

NOTE: When CO₂/Air mode ⑨ is selected, the automatic switch from air to CO₂ supply is not shown on the display.

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
7.4.7 Selecting the gas supply

Click on the Settings icon ① to open the Settings menu (see Section 7.5 'Settings').

NOTE: The Settings menu can be called up when the light and insufflation are switched off. (Icons on the Start/Stop buttons ② and ⑤ are blue)

Select the menu item 'CO₂ settings' in the Settings window.

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
Select the connected gas supply in the Settings menu. When using a CO₂ compressed air bottle (high pressure), select button 10. When using a CO₂ central supply (low pressure), select button 13 accordingly.

NOTE: Please make sure that the gas supply is properly connected. Further information can be found in the section 'Setup and Operating Instructions' (see Sections 7.3.3 to 7.3.9).

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7.4.8 CO₂ information message when input pressure is too low

The CO₂mbi LED light source automatically recognizes when the connected CO₂ supply falls below a critical value. The user is informed of this via an information message in the top screen bar and an icon (bottom right).

NOTE: When using a CO₂ bottle, different bottle icons are used to indicate the fill level. Only one icon is used when the gas is connected via the central supply.

NOTE: If the fill level of the CO₂ bottle is low, please replace the CO₂ bottle as quickly as possible in order to ensure that the device runs smoothly.

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


7.4.9 Changing CO₂ to CO₂/Air mode when input pressure is too low

If the CO₂ supply does not guarantee sufficient insufflation via the connected supply system (CO₂ bottle/central supply), this is indicated as an information message in the top screen bar and as an icon (bottom right). The user can choose whether to continue insufflation with air or to stop insufflation.

Select 'Continue insufflation with air' if there is no longer sufficient CO₂ supply during the examination and if you wish to end the procedure using air as an insufflation medium.

NOTE: Once 'Continue insufflation with air' has been selected, the insufflation medium cannot be changed.

Select 'Stop insufflation' to replace the CO₂ bottle during operation.

Select 'Change to CO₂' once the CO₂ bottle has been replaced.

7.4.10 Decommissioning

Close the valve on the CO₂ bottle or the valve on the central supply at the end of the procedure.

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7.5 Settings

Click on the menu icon ❶ to call up the Settings menu.

❶ **NOTE:** The Settings menu can be called up when the light ❷ and insufflation ❸ are switched off (icons on the Start/Stop buttons ❷ and ❸ are blue).

In the Settings menu, you can select various setting and management fields.

- Service (only authorized service technicians)
- Language
- Event protocol
- Device info
- Audio (on/off)
- CO₂ settings (for information on selecting the CO₂ mode, see Section 7.4.4; for information on selecting the gas supply, see Section 7.4.7)

7.5.1 Service

The 'Service' window is opened by pressing the 'Service' button in the Settings menu.

❶ **NOTE:** This area is password protected.

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
7.5.2 Language

The 'Language' window is opened by pressing the 'Language' button in the Settings menu.

The following screen languages are available: Danish, English, Finnish, French, German, Italian, Polish, Portuguese, Russian and Spanish (displayed in the appropriate language).

Select a new language if required, and then confirm.

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
7.5.3 Event protocol

The corresponding window is opened by pressing the 'Event protocol' button in the Settings menu.

The Event protocol contains the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

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7.5.4 Device Info


The 'Settings' window is opened with device-specific data by pressing the 'Device Info' button in the Settings menu.

In this window, you will find information about your device (e.g., the serial number, software version, and operating hours).

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7.5.5 Audio settings


The acoustic signals can be switched on or off (Audio on/Audio off) directly in the Settings menu.

NOTE: When the acoustic signals are switched off, this is indicated by means of an icon in the top row of the main menu.

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

8. Maintenance

8.1 Fuse replacement

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Switch off the device and remove the power plug from the electrical outlet.

Remove the fuse holder ⑧ with a screwdriver or other suitable tool.

WARNING: Only use fuses of the correct rating.

Insert new fuses.

	100...240 V~
Line fuse	2 x T 1.6 A H 250 V

Refit the fuse holder ⑧.
Reconnect the power cord.
Test the device for proper operation.

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8.2 Reprocessing

8.2.1 General warnings for medical devices

WARNING: Risk of infection: Incorrectly prepared medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: Country-specific laws and regulations must be observed.

NOTE: The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting www.karlstorz.com.

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8.2.2 Wipe-down disinfection of device

WARNING: Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

CAUTION: Make absolutely sure that no liquid can get into the device housing. When reprocessing the videoscope or its channels, it must always be disconnected from the CO₂mbi LED.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. Remove any disinfectant with a low lint disposable cloth after the exposure time.

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8.2.3 Reprocessing of the water bottle 13992BS

NOTE: For information on reprocessing, please refer to the respective instruction manual (water bottle model 13992BS: 96096310 DF).

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8.3 Maintenance and safety check

8.3.1 Maintenance

It is recommended to service the insufflation pump every 2.5 years or after 5000 operating hours. Regular maintenance can contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

8.3.2 Safety checks/ repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

Carry out the functional check as per the instruction manual (see section 7.4).

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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8.4 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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8.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

WARNING: At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this device.

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8.6 Repair program

For videoscopes and devices, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany you can refer repairs directly to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 NEUHAUSEN

Servicehotline: +49 7461/708 980
E-Mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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8.7 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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8.8 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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9.1 Information signals

9.1.1 Configuration

All of the information signals created, with the exception of those which display the cause of an inoperable device ('Malfunction', see e.g. page 19 Self test), are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration, with which the described information signals also comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by a signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, only the condition most recently detected is signaled.

9.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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9.1.3 Acoustic signals

The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. But at least one complete signal sequence is played.

The acoustic signal can be switched on or off using the button in the 'Menu'.

But the acoustic signal can only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

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


9.1.4 Information signals (general)

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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

9.1.5 Availability signal


On completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

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9.2 Troubleshooting

WARNING: Always unplug the device before all maintenance work.

Symptom:

- Complete failure of the device.

Possible causes:

- Power supply failure.
- Defective power fuse.
- Power cord connector is not properly connected to device socket.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord connector firmly into device socket.

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Symptom:

- Pump on – no pressure/too low pressure.

Possible causes:

3. Videoscope not properly inserted.
4. Insufflation pump worn out.
5. Internal/external device cable kinked.
6. Air filter ⑤ blocked/dirty.

Remedy:

1. Fully insert the videoscope.
2. Notify Service.
3. Rectify the kink (have this performed by qualified service personnel if the kink is in an internal device cable).
4. Unscrew both air filters ⑤ and replace with two new air filters (Art. no. 2786890).

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
No light is emitted.

Possible causes:

1. Electronics are defective.
2. Overheating of the CO₂mbi LED (LED has switched off due to overheating of the unit caused by covered air vent).

Remedy:

1. Notify Service.
2. Remove objects from in front of the air vents. Switch off the CO₂mbi LED, let it cool down (approx. 10–15 minutes) and switch it back on. Ensure that it is sufficiently ventilated.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

According to IEC 60601-1-2:

Observe the information on electromagnetic compatibility in chapter 11 (p.45-58).

Directive compliance (for TL 100)

According to Medical Device Directive (MDD):

This medical device belongs to Class IIa

This medical device bears the CE mark in accordance with MDD 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

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9.4 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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**10. Spare parts, recommended
accessories**

10.1 List of Spare parts

Item	Cat. no.
Water bottle	13992 BS
Irrigation bottle holder	TL001
Bottle ring for KARL STORZ water bottle	TL002
Bottle ring for use with Ampuwa Plastipur irrigation bottles	TL003
Bottle holder for KARL STORZ CO ₂ container	TL004
Low pressure tube, for central CO ₂ gas supply, length 150 cm	UI001
Same, length 300 cm	UI002
Same, length 600 cm	UI003
High pressure tube, American/ German connection, length 55 cm	20 4000 21
Same, length 102 cm	20 4000 27
High pressure tube, American connection/PIN-Index connector, length 55 cm	20 4000 22
Same, length 102 cm	20 4000 28
High pressure tube, American connection, length 102 cm	20 4000 25
Same, length 55 cm	20 4000 29
High pressure tube, American connection/ISO connection, length 102 cm	20 4002 22
Air filter	2786890
SCB connecting cable, length 100 cm	20 0901 70
Line fuses T1,6AH1250V, pack of 10	1973290
Power cord (Schuko)	400 A
Power cord "Hospital Grade" (USA)	400 B
Instruction manual	96206538 DER

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10.2 Recommended accessories

Item	Cat. no.
Leakage tester	13242LX
Adapter for KARL STORZ Standard fiberoptic and fluid light cables	495KS
Videoscope *All KARL STORZ endoscopes are suitable for use.	*
Light cable *All KARL STORZ light cables are suitable for use (with the exception of 495TIP, 495NCS, 495NCSC, 495NAC, 495ND)	*

For detailed information please contact:
KARL STORZ GmbH & Co. KG
Postfach 230
78503 TUTTLINGEN/GERMANY

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11. Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The CO₂mbi LED model TL100 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The CO₂mbi LED model TL100 is a Group 1 device (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The device should not be used adjacent to or stacked with other equipment. If it must be operated next to or on top of other devices, you should observe the device or the system to verify its operation according to its designated use in the configuration in which it is used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ Instruction manual may result in increased emissions or decreased immunity of the device. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
CO ₂ mbi LED model TL100 is intended for use in the electromagnetic environment specified below. The customer or user of the CO ₂ mbi LED model TL100 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	CO ₂ mbi LED model TL100 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The CO ₂ mbi LED model TL100 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions acc. to CISPR 11	Class B	
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	

Table 2
Guidance and manufacturer's declaration – electromagnetic immunity

CO₂mbi LED model TL100 is intended for use in the electromagnetic environment specified below.
The customer or user of the CO₂mbi LED model TL100 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrical interferences/bursts acc. to IEC 61000-4-4	± 2 kV for mains lines ± 1 kV for input/output lines	± 2 kV for mains lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<5 % U_T^* (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T^* (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	The magnetic fields at the mains frequency should be at the values typical of business or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
CO₂mbi LED model TL100 is intended for use in the electromagnetic environment specified below. The customer or user of the CO₂mbi LED model TL100 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6 Radiated RF disturbances acc. to IEC 61000-4-3	3 V_{rms} 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V_{rms} 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the CO₂mbi LED model TL100, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 6
Recommended separation distances between portable and mobile RF communications equipment and the CO₂mbi model TL100

The CO₂mbi LED model TL100 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvets.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar — Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA — Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) — Units 205 — 206
Dubai Health Care City, P.O.Box 118069
Dubai — United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguai Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhua District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjall-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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