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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.




This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the Stroboscope PULSAR II. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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- ① Power on/off
- ② Display for:
 - Tone frequency (in Hz)
 - Pitch (in musical notation)
- ③ Display for:
 - Volume (in dB)
 - Brightness (in % of max. brightness)
- ④ ± buttons for setting and/or selection of:
 - brightness
 - frequency
 - display values
 - speed of the slow motion
- ⑤ Drawer for ready-reference guide
- ⑥ Slow motion
- ⑦ Freeze frame
- ⑧ LED Light source on/off
- ⑨ Automatic brightness control on/off
- ⑩ Impact sound microphone connection
- ⑪ LED light source connection (LED stroboscope and microphone)
- ⑫ Footpedal connection
- ⑬ ACC1 socket (communication with 'documentation')
- ⑭ ACC1 socket (ACC connection to external units)
- ⑮ S-Video video input from CCU (Camera Control Unit)
- ⑯ S-Video video output to monitor
- ⑰ Audio Line Out
- ⑱ Potential equalization connection
- ⑲ Power cord socket with integrated fuse
- ㉑ Impact sound microphone 40 1600 30
- ㉒ Footpedal 40 1600 32
- ㉒ LED light source 40 1600 21

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Photographs of the unit	V
Controls, displays, connectors, and their uses	VII
Symbols employed	VIII
General information	
Description of the unit	2
Safety instructions	
Warnings and cautions	6
Intended use	9
User qualification	9
Safety precautions at the site of installation	10
Safety precautions when operating the unit	10
Installation and operating instructions	
Unpacking the unit	11
Basic equipment	11
Installing and connecting up the unit	11
Operation	13
Description of operation	16
Standby mode	17
Continuous light mode	19
Stroboscopy mode	20
Freeze frame function	21
Slow motion function	21
Setting the brightness	22
Automatic brightness control	22
Video functions	23
Maintenance	
Fuse replacement	26
Cleaning	27
Maintenance and safety check	28
Servicing and repair	29
Repair program	30
Limitation of liability	31
Manufacturer's warranty	31
Technical description	
Troubleshooting	32
Technical data	35
Technical documentation	36
Block diagram	37
Spare parts, accessories	
List of spare parts	38
Recommended accessories	39
Appendix	
Table of musical notes and frequencies	40
Electromagnetic Compatibility (EMC) Information	41

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Description of the unit

The PULSAR II stroboscope is a modern device that generates light with a high performance LED (light emitting diode) lamp for stroboscopic examination (stroboscopy mode) and also for normal viewing (continuous light mode). To be more precise, the continuous light is actually a pulsating (flashing) light that is perceived as continuous light due to its high frequency. Even though there is no obvious difference in technical terms, for the sake of simplicity the conventional terminology (flashing light/continuous light) is still used below. The advantages of this system concept are:

Compact design of unit and footpedal

- Low weight
- Latest LED technology translated to low energy consumption and no need to replace lamp
- No difference in brightness or color temperature for flashing and continuous light
- Flash performance also suitable for thin and flexible endoscopes
- No noise generation
- Possibility to superimpose text in the video image (overlay)
- Display in musical notation

Compared with uniform illumination, pulsating light does not have any disadvantages if, as practically obligatory in endoscopy, a video system is used in conjunction for visualization and documentation, as the video system also functions intermittently (50/60 Hz). In contrast, the PULSAR II stroboscope is adapted to a video system in a special way. By means of special trigger switching, disruptive effects are avoided (bright-dark modulation, drop outs, beats, etc.), whose cause can be found in the differing frequencies of the voice and video signal. In order to make good quality video recordings possible with small instruments (fiber instruments), a particularly powerful high-performance LED lamp is integrated.

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Other advantages of this concept are:

- no flicker or interruptions in the video image; uniform brightness
- no reduction ratio dependent purely on the voice signal; smooth transitions
- high flash capacity for video recordings even with thin-caliber, flexible endoscopes

Other video functions available as options are:

- automatic brightness control
- superimposing texts on the video image

The unit's functions can be further expanded with the following options:

- Connection to KARL STORZ documentation units
- Connection to KARL STORZ video systems

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
External system components

LED light source with integrated microphone

The LED light source is composed of a compact light source in multichip SMT technology with extremely high brightness and light density and a highly sensitive airborne sound microphone. The audio signal is captured by the microphone integrated in the handle, amplified and analyzed to trigger the stroboscope with voice synchronization. The displays ② and ③ on the unit show the frequency in Hz and the volume in dB. You have the choice to switch between the frequency and the musical notation.


For sound recording, the audio signal is available on the unit's rear panel, according to standard.

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
Impact sound microphone

The functions described in the 'LED light source with integrated microphone' section (see above) can alternatively be realized using the impact sound microphone which can be found in the scope of supply.

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Footpedal

The following functions are controlled with the footpedal:

- tuning the phase in freeze mode
- switching between freeze frame mode and slow motion mode
- starting and stopping the documentation unit (optional)
- Setting the slow motion frequency

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Other equipment

Other equipment can be connected to the stroboscope:

- Camera control units,
e.g., KARL STORZ TELECAM/TRICAM®
- Documentation, e.g., KARL STORZ AIDA®
- Audio equipment for amplification and/or
recording
- KARL STORZ monitors

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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WARNING: The electrical installation of the operating room must comply with the applicable IEC requirements.

WARNING: The unit must NOT be used in operating rooms.

WARNING: Keep the device out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (MDD: Appendix 1: 13.6.c, IEC 60601-1: 6.8.2.c, 19.2.b, 19.2.c).

WARNING: Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: The light outlet in the LED light source must not be coated with fluid or immersed at all.

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

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Safety instructions

Warnings and cautions

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Risk of burns! The surface temperature of the LED light source can exceed 41°C. Leave the LED light source to cool before touching it.

WARNING: Danger of glare! Never look directly into the light outlet when the light source is switched on. Under real conditions (for exposure time, eye pupil, viewing distance), the LED light source does not pose any risk whatsoever to eyes. However, intensive light sources present a high secondary hazard potential from their glare. As when having looked into other bright light sources (e.g. car headlights) vision may be temporarily compromised and afterimages may be seen which, depending on the situation, can result in irritation, annoyance, impairment or even accidents.

WARNING: Always unplug the unit before doing any maintenance work on it (e.g. cleaning).

WARNING: Danger of electric shock! Do not open the LED light source! Refer servicing to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims.

CAUTION: Only operate the unit with the voltage stated on the identification plate.

CAUTION: When replacing fuses, use only fuses of the correct rating.

CAUTION: The stroboscope PULSAR II may only be used with accessories, wearing parts, and disposable items which have been designated as suitable for the stroboscope by KARL STORZ.

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

CAUTION: No liquids or dirt particles may come into direct contact with the LED light outlet. This could cause malfunctions or defects at the light source.

① **NOTE:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The Stroboscope PULSAR II is a unit for diagnosis of larynx disorders in hospitals and medical practices.

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Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The Stroboscope PULSAR II may only be operated with accessories which have been designated by KARLSTORZ as suitable for the unit.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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User qualification

The Stroboscope PULSAR II may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

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It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The Stroboscope PULSAR II is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

WARNING: The unit is exclusively designed for indoor use and must not come into contact with water or dust.

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
Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it (see the section 'Test for proper operation').

During treatment using the Stroboscope PULSAR II the patient must be treated and kept under observation with the usual medical care.

Any treatment may only be performed if there is visual observation of the action of the unit.

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Unpacking the unit

Carefully unpack the Stroboscope PULSAR II and remove the unit and its accessories from their packaging.

Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Basic equipment

- 1 Stroboscope PULSAR II
- 1 LED light source with integrated microphone
- 1 Footpedal
- 1 Impact sound microphone with holder
- 1 Power cord
- 1 Instruction manual

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Installing and connecting up the unit


ⓘ NOTE: The Stroboscope PULSAR II, including any accessories connected, may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

ⓘ NOTE: If a video camera with shutter function is used, set it to 1/50 s (PAL) or 1/60 s (NTSC).

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⚠ WARNING: The unit is not intended for use in operating rooms. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
The unit is equipped with a connector ⑮ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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Operation

Only operate the instrument with an AC voltage of 100...240 V, 50/60 Hz.


CAUTION: Only operate the unit at the voltage printed on it!

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
Fit the unit with fuses corresponding to the line voltage (fuse holder is located in the unit's power cord socket ⑯).

Line voltage	Fuse rating
100...240 V~	T 2.0 AL250V

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
! WARNING: The unit should only be connected to or disconnected from the power supply outside of areas subject to explosion hazards.

Connect the power cord by inserting the connector into the unit's power socket ⑯ as far as it will go.

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Connect the LED light source up to the PULSAR II (connection ⑪).

! NOTE: The functions of the PULSAR II stroboscope are blocked when there is no LED light source inserted.

! WARNING: Do not remove the LED light source when the light is activated. This could otherwise result in damage to the LED light source or the PULSAR II stroboscope.

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Installation and operating instructions

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
Connect the LED light source with the desired telescope (using the telescope's fine threaded connection).

Ensure that the microphone integrated in the LED light source is pointing in the direction of the patient.

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
The turning mechanism integrated in the LED light source makes it possible to turn the LED light source to the desired position even after adaptation to the telescope.

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Connect the corresponding foot pedal to the PULSAR II stroboscope (connection ⑫).

① **NOTE:** No foot pedal is required for the activation of the continuous light.

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optional:


If you are using the impact sound microphone contained in the scope of supply, connect it to the PULSAR II (connection ⑩).

① **NOTE:** If you use the impact sound microphone, the microphone in the LED light source is turned off.

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
optional:
Connect the camera control unit (CCU) up to the
PULSAR II for recording videos.
Connect the S-Video output on the CCU up to the
S-Video input on the PULSAR II ⑮.

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
optional:
Connect the PULSAR II up to the monitor for
recording videos.
Connect the S-Video output of the PULSAR II ⑯
with the video input on the monitor used.

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optional:
Connect the ACC2 interface ⑭ with the ACC
interface of the respective documentation unit
for controlling the video documentation using the
PULSAR II footpedal.

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① **NOTE:** The function of the footpedal needs
to be set in the respective KARL STORZ
documentation unit.

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The unit is now ready for use and can be switched
on.
Switch on the PULSAR II using the power switch
①.

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Description of operation

Initializing the unit

The unit is initialized after switching on. The following functions are carried out:

- All display elements light up for approx. 2 s.
- Internal and external system components are checked and if necessary, an error message appears.
After activation, the PULSAR II checks whether the footpedal and microphone have been connected and are functioning correctly. In the event of a malfunction, a corresponding error message is displayed.
- The last set brightness of the LED light source or the brightness of the capture screen in videostroboscopy is displayed respectively.
- The upper display also displays the 'PULSAR II' ticker, which signifies that the unit is ready for use.
- The system settings valid before the unit was last switched off are activated. These are:
LED light source:
- display value (frequency, musical notation)
- brightness setting (percentage brightness)
- Switchover to standby mode

NOTE: The properties of the LED light source give it a very high service life. To avoid unnecessary heating up of the LED light source, the PULSAR II has an automatic turn-off function for the LED light source in continuous light mode.

NOTE: If the temperature of the LED light source exceeds 50°C, the LED light source turns off automatically in order to protect the safety of the user. The LED light source can then only be reactivated if the temperature of the LED light source falls back below 41°C.

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Standby mode

After the unit has initialized, the unit switches to standby mode. The brightness display appears. The standard setting is 50%. Simultaneously, the 'PULSAR II' ticker is shown in the upper display, which signifies that the unit is ready for use. The LED light source is turned off in standby mode.

Standby mode is automatically activated if the following conditions are fulfilled:

- the LED light source was in continuous light mode last
- 10 minutes have elapsed since the last user action

The **P** button is used to exit standby mode and switch to continuous light mode. The ticker in the upper display disappears.

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Settings

Standby mode allows you to make settings concerning the brightness of the LED light source, the later display and the slow motion function. In detail, these are:

- switching of the display between frequency and musical notation
- brightness of the LED light source
- speed of the slow motion (in slow motion function)


Further details about these settings are given in the sections dealing with the respective modes and functions.

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
Switching frequency to musical 'notation'

To switch between the two specified indicators, press the **±** buttons at the same time and hold for around 3 seconds (in standby mode). The symbol off the respective displays illuminates shortly after activation as confirmation.

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Setting the slow motion frequency


The speed of the visible vocal cord movement can be set infinitely between 0.2 and 2.5 Hz. This switching is only possible in standby mode. It is done by pressing the foot pedal in standby mode and at the same time pressing the + button ④ or - button ④ to select the speed (frequency) of the slow motion.

The respective setting is shown in the display ② (in Hz). When the button ④ is released, the previously shown display reappears. The set frequency is also displayed in the video image (overlay).

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Setting the brightness

The brightness is displayed in the range from 1% to 100% (max. brightness) on display ③. It can be changed using the +/- buttons ④.

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Continuous light mode

In continuous light mode the LED light source is activated and the stroboscope is ready for examinations.

The brightness of the LED light sources can either be set manually or, if a valid video signal is connected to the video input, controlled automatically. The automatic button (8) is used to switch between manual and automatic.

NOTE: Audio signals are not evaluated in continuous light mode. The audio signal is, however, present at the Audio Line Out (17).

NOTE: Manual setting of the brightness is described in detail on p. 22.

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Stroboscopy mode

In stroboscopy mode, the PULSAR II stroboscope is triggered by audio sources (speech, notes from the patient), which are captured by the microphone integrated in the LED light source or the impact sound microphone. The LED light source is triggered at the same time as the audio signal of the microphone.

If the footpedal is (fully) pressed in continuous light mode and a valid audio signal* is present (sufficient volume and periodicity), the unit switches to stroboscopy mode. Triggering is then (quasi-) synchronous with the audio signal. If either of the above conditions is no longer fulfilled, the unit switches back to the continuous light mode.

* If the unit cannot detect the correct trigger frequency, the footswitch must be released and pressed again. The unit can then carry out triggering again.

When successfully triggered by the audio signal, one of the following indicators appears in the upper display ② based on the current settings (see Standby mode):

- tone frequency in Hz
- pitch in musical notation

① **NOTE:** To switch between these indicators see 'Standby mode' section.

The lower display ③ shows the following indicator:

- volume in dB

Screen display

The indicators specified above are additionally overlayed in the video image when a video system is connected.

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Freeze frame function

In the upper section of the pedal movement range, freeze frame mode is active, i.e. a stationary image of the vibrating vocal cords appears. The freeze frame symbol ⑦ lights up. The respective phase, i.e., the position of the (apparently) stationary vocal cords, can be adjusted using the foot pedal within a range from 0 to 360 degrees.

❗ **NOTE:** The display of the volume (in dB) is only correct if the microphone integrated in the LED handle is not covered or the optional impact sound microphone is correctly affixed to the patient and connected properly to the PULSAR II.

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Slow motion function


If the pedal is pressed further down beyond the pressure point (as far as the bottom limit stop), slow motion function is activated. The vocal cords then appear to move more slowly. The slow motion symbol ⑥ lights up (⑥ and ⑦ are only active alternately).

The speed of the (apparent) vocal cord movement can be set infinitely between 0.2 and 2.5 Hz. This setting is only possible in standby mode. It is done by pressing the foot pedal in standby mode and at the same time pressing the + button ④ or – button ⑤ to select the speed (frequency) of the slow motion.

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Setting the brightness

The brightness is displayed in the range from 1% to 100% (max. brightness) on display ③. It can be changed using the ± buttons ④.

① **NOTE:** Setting the brightness affects both the manually set brightness and the brightness in automatic brightness control.

① **NOTE:** Increased brightness raises the power dissipation of the LED and thus also the surface temperature of the LED light source. For continuous operation, the brightness should be reduced as far as possible so that the LED light source does not become too warm.

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Automatic brightness control

Pressing the automatic button ② allows you to switch between the manual brightness setting and automatic brightness control, provided the following conditions are fulfilled:

- video camera connected to endoscope
- video signal (from camera) is present at the video input ⑤


⚠ **WARNING:** Observe the safety instructions on page 7.

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The level of the video signal is kept constant for controlling the brightness. This avoids brightness fluctuations or blooming of the video image as the result of varying object distances. The basic brightness of the video image can be changed by pressing the ± buttons ④. The current value is shown on the display ③ as long as automatic brightness control is active.

① **NOTE:** This value does not reflect the actual brightness of the light source, since that depends on many parameters, such as object distance, endoscope, camera etc. The task of automatic brightness control is to compensate for the influence of these parameters on the overall result.

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NOTE: If video cameras have a shutter function (automatic brightness control), it must be turned off using the camera control unit or the shutter speed defined at 1/50 s or 1/60 s. Otherwise the image will flicker.

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Video functions

There are also additional display functions available in addition to the described automatic brightness control. Parallel to the display in the displays ②/③, the status information is displayed via overlaying in the video image. This information is particularly valuable for archiving.

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• Status information

The pitch and volume values are displayed continuously or temporarily along the bottom edge of the video image. The pitch is given either in Hz or in musical notation, depending on how the unit is set. The volume is displayed in dB. In addition, the set slow motion frequency in slow motion mode is also displayed in the bottom right corner of the screen as 'Δ Hertz'.

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A KARL STORZ documentation unit is normally used for documenting video images and/or sequences.


! WARNING: Observe the safety instructions on page 7.

The PULSAR II stroboscope (in conjunction with the KARL STORZ documentation units) is suitable for controlling the recorder functions that are important for the examination, namely start of recording, pause and stop. Control is by means of the foot pedal.

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Starting, interrupting and ending a recording

Recording can now be started by briefly operating the foot pedal twice. It does not matter how far the foot pedal is pressed down, provided that it is back in its upper home position within half a second. If the foot pedal is operated a second time in the same way, the KARL STORZ documentation unit will interrupt the recording with the Pause function. The recording can either be resumed as described above or finished completely by manually stopping the KARL STORZ documentation unit.

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! NOTE: Briefly pressing the foot pedal is not typical of normal use during laryngostroboscopy. That is why this action was chosen for controlling the recording. It means that the operator does not have to move his foot, as would be the case with a separate button for the remote control function.

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! NOTE: To be able to use this function, the corresponding documentation unit is connected with the PULSAR II via the ACC2 interface. The function of the foot pedal needs to be set in the respective KARL STORZ documentation unit.

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
**Other functions****Sound recording**

The signal picked up by the microphone can be recorded with a sound recording device.

WARNING: Observe the safety instructions on page 7.

This signal, amplified to normal level, is available at output ⑰. Even if a video recorder is used, this outlet can be used for the additional sound recording (onto the video tape). As the signal is not two-channel, the recording device must be switched manually to 'mono' if it is a stereo recorder (be sure to use the correct input). Alternatively, a special Y cable (available from retailers) must be used, which connects the two stereo inputs of the recorder to the single audio output of the stroboscope.

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

Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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
Remove the power fuse holder with a screwdriver or other suitable tool.

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! **CAUTION:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

	100...240 V~
Power fuse	2 x T 2.0 AL250V

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Refit the power fuse holder.
Reconnect the power cord.
Test the unit for proper operation.

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Cleaning

WARNING: Always pull out power plug before cleaning!

CAUTION: It is essential to prevent fluids from entering the unit.

CAUTION: It is essential to prevent fluids from entering the light outlet on the LED light source.

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Cleaning the PULSAR II and accessories (wipe-down disinfection only)

In detail these are:

- PULSAR II basic unit
- LED light source
- Impact sound microphone with holder
- Footpedal

NOTE: For wipe-down disinfection, we only recommend surface disinfectants (no solutions predominantly containing alcohol such as rapid-action disinfectants) and 'disposable disinfectant methods' based on ammonium chloride.

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Maintenance and safety check

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspecting the device safety fuses
- Protective ground resistance measured according to IEC 62353 ed1.0: Limit value: 0.2 ohms
- Earth leakage current measured according to IEC 62353 ed1.0: Limit value: NC: 500 µA, SFC: 1000 µA
- Patient leakage current measured according to IEC 62353 ed1.0: Limit value: 100 µA (BF), 10 µA (CF)
- Leakage currents must not exceed the limit values specified above.

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Test for proper functioning

A functional check must be carried out in accordance with the instruction manual.

The safety check must be entered in the unit logbook and the results of the check must be documented.

If the device does not function reliably and/or safely, it must be repaired.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

NOTE: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Troubleshooting

Symptom:
Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Be sure to fit the correct type of fuse.

WARNING: Always unplug the unit before carrying out any maintenance work!

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Display: **E01**

Possible causes:

- Foot pedal not correctly plugged in or defective.

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Display: **E02**

Possible causes:

- LED light source is not connected correctly or is defective.

Remedy:

- Switch the unit off.
- Check the connector.
- Switch the unit on again.

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Display: **E03**

Possible causes:

- LED light source is too hot (warmer than 50°C).
- Temperature sensor is defective.

Remedy:

- Allow the LED light source to cool down.
- Notify Service.

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Display: **E04**

Possible causes:

- Defective memory module in stroboscope for user data in PULSAR II.

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Display: **E05**

Possible causes:

- Defective memory module in LED light source.

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Display: **E06**

Possible causes:

- Defective components in stroboscope.

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Display: **E07**

Possible causes:

- Defective components in the LED light source or incorrect transmission.

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- ACC interface error.

Remedy:

- Check connection to external units.
- Notify Service.

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- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type BF 

According to IEC 60601-1-2:2007:

Please read the Electromagnetic Compatibility
Information in the appendix (p. 41-53).

Directive compliance (for 40160120)

According to Medical Device Directive (MDD):
This medical product belongs to Class IIa

This medical product bears the CE mark in
accordance with MDD 93/42/EEC.

NOTE: The code number after the CE mark
indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those
circuit diagrams, itemized parts listings, de-
scriptions, sets of adjustment instructions, and
other items of available documentation to suit-
ably qualified user personnel duly authorized by
the manufacturer for their use in repairing those
components of the unit that have been designated
as repairable by their respective manufacturers.

Supply of such technical documentation relating
to the unit shall not be construed as constituting
manufacturer's authorization of user's personnel,
regardless of their level of technical training, to
open or repair the instrument.

Explicitly exempted herefrom are those
maintenance and repair operations described in
this manual.

We reserve the right to make engineering modi-
fications in the interest of promoting technological
progress and generating performance improve-
ments without obligation on our part to submit
prior notice thereof.

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
List of spare parts

Item	Cat. no.
Footpedal	40 1600 32
Impact sound microphone	40 1600 30
LED light source	40 1600 21
Power fuses T 2.0 AL250V [100...240 V~]	2027590
Power cord (grounded)	400 A
'Hospital Grade' Power cord (USA)	400 B
Instruction manual	96056023 D
Operating instruction quick reference guide	96056023 Z

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Recommended accessories

We recommend the following telescopes for use with the KARL STORZ PULSAR II:

Item	Cat. no.
HOPKINS® lateral telescope 90° for laryngoscopy, autoclavable, dia. 5.8 mm, length 20 cm, color code: blue	8700 DKA
BENJAMIN tele-laryngoscope with integrated HOPKINS® lateral telescope 70°, viewing angle 50°, dia. 5.8 mm, length 19 cm, autoclavable, with fiberoptic light transmission, color code: yellow	8700 CKA
BENJAMIN slim line tele-laryngoscope with integrated HOPKINS® lateral telescope 70°, viewing angle 50°, dia. 4 mm, length 18 cm, autoclavable, with fiberoptic light transmission, color code: yellow	8705 CKA
Strobo-laryngoscope with integrated HOPKINS® lateral telescope 70°, working length 17 cm, autoclavable, oval shaft, 7.2 mm x 9.3 mm, with fiberoptic light transmission, color code: yellow	8706 CA
Tele-laryngo-pharyngoscope with HOPKINS® lateral telescope 90°, 4-fold magnification, can be focused, autoclavable, with fiberoptic light transmission	8707 DA

Further devices can also be connected to the stroboscope, e.g.:

- KARL STORZ AIDA® compact NEO SD documentation system for the digital storage of still images, video sequences and audio files
- Camera control unit TELECAM SL II with, for example, C-mount camera head and matching C-mount lens

For detailed information please contact:
 KARL STORZ GmbH & Co. KG
 Postfach 230
 78503 Tuttlingen/Germany

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Table of musical notes and frequencies

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The PULSAR II model 40 1601 20 corresponds to EN/IEC 60601-1-2 :2007 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive 2007/47/EEC. These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The PULSAR II model 40 1601 20 is a Group 1 unit (as per CISPR 11).

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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Appendix

Electromagnetic Compatibility (EMC) Information

! WARNING: The PULSAR II model 40 1601 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

! WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

! WARNING: The use of accessories, transducers and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the PULSAR II model 40 1601 20. The accessories, transducers and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories, transducers and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200

Accessories, transducers and cables which have been shown to comply with EN/IEC 60601-1-2:

Type	Shielded	Length [m]	Ferrite	Used for
40 1600 32	No	–	No	Footpedal
40 1600 30	Yes	–	No	Impact sound microphone
40 1600 21	Yes	–	No	LED light source (handpiece)
Power cord	No	3	No	Power supply

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Table 201

Guidance and manufacturer's declaration – electromagnetic emissions

PULSAR II model **40 1601 20** is intended for use in the electromagnetic environment specified below.
The customer or user of the PULSAR II model **40 1601 20** should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	PULSAR II model 40 1601 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PULSAR II model 40 1601 20 is suitable for use in all establishments (e.g. hospitals and medical practices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class D	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

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Table 202
Guidance and manufacturer's declaration – electromagnetic immunity

PULSAR II model 40160120 is intended for use in the electromagnetic environment specified below.
 The customer or user of the PULSAR II model 40160120 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 204
**Guidance and manufacturer's declaration – electromagnetic immunity –
 for equipment and systems that are not life-supporting**

PULSAR II model 40 1601 20 is intended for use in the electromagnetic environment specified below.
 The user of the PULSAR II model 40 1601 20 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to < 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the PULSAR II model 40 1601 20, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz $d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^a . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
PULSAR II model 40 1601 20**

The PULSAR II model 40 1601 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment.

The warranty referred to in our Standard Conditions of Business shall apply.

To be completed by the supplier/importer:

Company stamp/signature:

Complete warranty card upon purchase/delivery and return immediately to:

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Please attach
sufficient
postage

To be filled out by instrument owner:
Return address/company stamp:

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KARL STORZ Endoscopy — Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O.Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegne@karlstorz-rce.com

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0176, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 (0)24 6258 9911
Fax: +86 (0)24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
#06-07 e-Centre @ Redhill
Singapore 159471
Phone: +65 65325548, Fax: +65 65323832
E-Mail: info@karlstorz.com.sg
E-Mail: servicesg@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 (0)3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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