STORZ
KARL STORZ — ENDOSKOPE

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This instruction manual is intended to serve as an aid in the proper connection and operation of the C-CAM™ camera head. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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


**3 Controls, displays,
connectors, and
their uses**

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①	Coupling for fiberscope/telescope
②	Focus
③	C-MAC® (7" monitor)
④	C-HUB
⑤	Example of fiberscope
⑥	Example of telescope

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1	Important information for users of KARL STORZ instruments	4
2	Photographs of the unit	5
3	Controls, displays, connectors, and their uses	6
4	General information	9
4.1	Description of the unit	9
4.2	Compatibility	9
4.3	Explanations of connections	9
5	Safety instructions	10
5.1	Explanation of warnings and cautions	10
5.2	Warnings	10
5.3	Intended use	11
5.4	User qualification	11
6	Indications and Contraindications	12
6.1	Indications for use	12
7	Care and Cleaning	13
7.1	Care and handling	13
7.2	Approved Sterilization and Disinfection Methods Reference Chart	15
7.3	Cleaning, disinfection and sterilization guidelines	16
7.4	Preparation for cleaning and sterilization	17
7.5	Cleaning instructions	17
7.6	Camera Head Assembly	17
7.7	Sterilization instructions	19
7.7.1	Ethylene Oxide (EtO) gas sterilization ..	20
7.7.2	STERRAD® Sterilization Systems (STERRAD® 100S, NX, and 100NX)	21
7.7.3	STERIS Amsco V-PRO Sterilization System	22
7.8	High Level Disinfection Instructions	24
7.9	Use of sterile drapes	27
7.10	References for cleaning and sterilization ..	27

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8 Repairs and Warranty.....	28
8.1 Maintenance	28
8.2 Servicing and repair	29
8.3 Repair program.....	29
8.4 Limitation of liability	30
8.5 Manufacturer's warranty	30
8.6 Disposal.....	31
9 Technical Specifications.....	32
9.1 General technical data	32
9.2 Operating and storage conditions.....	32
9.3 Standard conformity	32
9.4 Directive compliance.....	32

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4 General information

4.1 Description of the unit

The C-CAM™ 20290131/20290132 is a camera head for the transmission of image and data signals in medical video endoscopy.

The C-CAM™ camera head is available in two versions:

20290131 (6-pin head connector); for use with the C-MAC 8401ZX

20290132 (8-pin head connector); for use with the C-MAC 8402ZX and C-HUB 20290120

4.2 Compatibility

The C-CAM™ camera head has a standard KARL STORZ grasping mechanism. All standard KARL STORZ endoscopes and fiberscopes can be used in combination with the camera head.

The C-CAM™ head connectors are color coded to denote 6-pin or 8-pin. The C-CAM™ camera heads must be used with the C-HUB or one of the following C-MAC™ monitors.

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4.3 Explanations of connections

NOTE: In chapter 3 on page 6 you will find an overview with images of the following explanations.

- The cable ① of the C-CAM™ camera head is connected to either a C-MAC or a C-HUB.
- The C-MAC® is a 7 inch monitor
- The C-HUB is a device for processing and transmitting image and control signals to external devices in medical video endoscopy.

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5 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

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5.1 Explanation of warnings and cautions

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NOTE: A Note indicates special information about operating the device or clarifies important information

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5.2 Warnings

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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 **WARNING:** Test this equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

 **WARNING:** Instruments and accessories used in combination with this camera head must be checked for perfect condition, functionality and completeness immediately before and after every application. Always have a spare available.

 **WARNING:** This camera head may only be used with accessories, wearing parts, and disposable items which are designated by KARL STORZ as suitable for the camera or the safe use of which is proven. For safety reasons, do not perform unauthorized conversions or modifications to the camera.

5.3 Intended use

The C-CAM™ camera head (20290131 / 20290132) is intended to serve as a mount for optical systems and for the creation and transmission of image and data signals to the camera control unit in medical video endoscopy.

The use of the camera head in any way other than described above is not permitted for safety reasons.

5.4 User qualification

The C-HUB may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the instrument.

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7 Care and Cleaning

7.1 Care and handling

① NOTE: The cap with which the camera head is delivered is designed only to protect the device during transport. Afterwards, it is not reused (no preparation).

To maximize the life of your camera, please follow the guidelines given below for proper handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE CONNECTION SOCKET ON THE MONITOR IF IT IS WET OR NOT CLEAN. Always ensure the cord connector is free of debris and completely clean. If moisture is present, dry thoroughly with low lint disposable cloth prior to insertion.

① NOTE: The electrical contacts are fragile and must be dried carefully.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, then provide secondary support for the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

STORE CAMERA SYSTEM AWAY FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

DO NOT STORE THE CAMERA WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Always arrange the cable in loose coils of 6 inches (15 cm) or more in diameter.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful "kinks" can result. Always allow the cable to uncoil gently and loosely.

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DO NOT AUTOCLAVE THE CAMERA HEAD. Safely sterilize by ETO gas with temperature not exceeding 140°F (60°C) or by the STERIS® or STERRAD® methods, according to manufacturer's instructions.

WARNING: The following alternative is only permitted in the USA: Alternatively, the camera head may be safely soaked in a no greater than 2.4% Glutaraldehyde disinfection solution for no longer than 45 minutes.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to device. When transporting, use the original shipping box to prevent damage to the device.

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7.2 Approved Sterilization and Disinfection Methods Reference Chart

CAUTION: The C-CAM™ camera head is Non-Autoclavable. DO NOT steam sterilize (autoclave) any part of the camera head or optical adaptor; irreparable damage will occur.

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WARNING: The "High Level Disinfection Methods" detailed below are only permitted in the USA.

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7.3 Cleaning, disinfection and sterilization guidelines

The latest country-specific cleaning, disinfection and sterilization regulations must be strictly adhered to.

To ensure optimal reproduction of the endoscopic image, the light entrances and exits of the camera head assembly, light cable and endoscope must always be kept clean. Clean exposed faces with a cotton swab dipped in a neutral soap solution and wipe with a cotton swab soaked in alcohol.

! CAUTION: Before and after disinfection or sterilization, inspect the camera head cable for breaks or cuts. Camera heads with damaged cables should NOT be disinfected or sterilized. A break or cut in the cable will allow fluid to enter the camera cable and cause damage. Return the damaged camera head to KARL STORZ for repair.

! CAUTION: Allowing build-up of disinfectant or sterilant to occur on the camera head card-edge connector can cause camera failure.

! CAUTION: The camera head assembly should be thoroughly cleaned and disinfected/sterilized according to validated infection control procedures prior to use and subsequent reuse. Any deviations from the recommended parameters for cleaning, disinfection and sterilization should be validated by the user.

! CAUTION: Use of sterilization/disinfection methods other than those recommended by KARL STORZ may cause adverse effects to the materials of the camera and endoscopic adaptor. Resulting damage will not be covered by the warranty.

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7.4 Preparation for cleaning and sterilization

Place the camera head assembly in containers and soak with a neutral pH (pH 6.0 - 8.0) enzymatic cleaning solution* immediately after use to prevent blood, protein and other contaminants from drying onto the equipment.

*** Examples for the US market:** Enzol, Metrizyme or equivalent diluted to proper concentrations per manufacturer's instructions)

7.5 Cleaning instructions

CAUTION: Wear protective gloves, clothing and a face mask for cleaning of contaminated equipment.

7.6 Camera Head Assembly

1. Thoroughly rinse the camera head assembly to remove all gross debris. Distilled water is recommended for cleaning and rinsing of the camera head assembly.
2. Completely immerse the camera head assembly and camera head cable in a neutral pH enzymatic cleaning solution* and distilled water. KARL STORZ does not recommend the use of detergents alone, as they contain high concentrations of surfactants which can leave a film on the equipment.
*** Examples for the US market:** Enzol, Metrizyme or equivalent per manufacturer's instructions)
3. Remove any residual debris and contaminants from the camera head assembly with a soft brush. Cleaning accessories are available from KARL STORZ.
4. Be sure that the card-edge connector is thoroughly clean.
5. **The following only applies to the US market:** Cleaning brushes should be cleaned and high level disinfected or sterilized daily.
The following applies to the non-US market: Cleaning brushes should be cleaned, disinfected and sterilized daily.

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6. Triple-rinse the camera head assembly with distilled water, for a minimum of one minute for each rinse. The rinse water should be discarded at the end of each rinse, as it will be contaminated with the cleaning solution. Thorough rinsing of the camera head assembly is necessary for removing any debris or detergent which could interfere with sterilization.
7. Dry the equipment with a lint-free soft cloth. Wipe the card-edge connector with 70% isopropyl alcohol to remove any residual detergent.
8. Do not allow exposed glass windows to air dry. 70% isopropyl alcohol may be applied to glass surfaces with a soft cotton applicator to prevent streaks and spots. Dry the surfaces thoroughly with a cotton applicator after applying the alcohol.
9. After cleaning, inspect the camera head assembly and camera head cable for cleanliness and damage.

! **CAUTION:** Inspect the camera head cable for breaks and cuts. Camera heads with damaged cables should not be sterilized or disinfected. Return camera heads with damaged cables to KARL STORZ for repair.

10. Before sterilization and/or disinfection, coil the camera head cable into loops at least six inches (15 cm) in diameter. Do not kink or twist the cable.

i **NOTE:** To prolong the life of the camera, KARL STORZ recommends the use of sterile drapes to cover the camera head assembly during surgical procedures.

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7.7 Sterilization instructions

Routine ethylene oxide (EtO) sterilization is recommended for initial and subsequent sterilization of all non-autoclavable camera head assemblies. They may also be Steris V-PRO 1™ or STERRAD® Sterilization System sterilized. To achieve the desired sterility assurance level (SAL) of 10⁻⁶, KARL STORZ recommends the following EtO, Steris V-PRO 1™ and STERRAD® Sterilization System sterilization methods.

CAUTION: DO NOT steam sterilize (autoclave) any part of the camera head or optical adaptor; irreparable damage will occur.

CAUTION: Before sterilization, the camera head assembly must be thoroughly cleaned and all organic material, blood and cleaning solution completely removed.

CAUTION: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and calibrated.

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7.7.1 Ethylene Oxide (EtO) gas sterilization

1. Clean the camera head assembly as described in the cleaning section of this manual. Place the camera head assembly in a sterilization tray. Do not place any instruments on top of the camera.

2. All sterilization cycles must include a preconditioning cycle:

Temperature: $130 \pm 5^\circ\text{F}$ ($54 \pm 2^\circ\text{C}$)
Time: 30 minutes nominal
Humidity: $\geq 70\%$ RH
Vacuum: 1.3 psia

3. KARL STORZ has validated Ethylene Oxide (EtO) sterilization using the following parameters:

Gas mixture: 100% EtO
Temperature: Set point of 55°C
Relative humidity: $\geq 70\%$ RH
Humidity Dwell Time: 30-45 minutes
Exposure time: 180 minutes, (full cycle)
EtO concentration: $735 \pm 30\text{ mg/L}$

4. Aeration of the equipment can be accomplished in any aeration chamber at a temperature of $120 - 130^\circ\text{F}$ ($50 - 54^\circ\text{C}$) for 12 hours.

Maximum permissible levels of residue following EtO sterilization (according to ANSI/AAMI/ISO 10993-7):

Ethylene oxide: 4 mg (< 24h)
Ethylene chlorhydrin: 9 mg (< 24h)

ⓘ NOTE: Sterilization methods which deviate from the validated processes described must be validated individually. The country-specific and sterilization device-specific requirements and outline conditions concerning the effectiveness of the sterilization must be taken into account, as must, where applicable, the necessary airtight times.

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7.7.2 STERRAD® Sterilization Systems (STERRAD® 100S, NX, and 100NX)

The STERRAD® Sterilization Systems, manufactured by Advanced Sterilization Products (ASP), utilize a synergism between hydrogen peroxide and low temperature gas plasma to produce a rapid, low temperature, low moisture inactivation of microorganisms. It is intended for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices. STERRAD® NX™ is a table top sterilizer with a rectangular chamber size of 51.3L / 26L (usable). It utilizes the same sterilant (hydrogen peroxide) as the other STERRAD® unit, the STERRAD® 100S. The major difference is the vaporization system, which renders the hydrogen peroxide sterilant more concentrated and improves diffusion into lumens. This combined vaporization system combined with a shorter cycle time results in a similar exposure to the sterilant. The STERRAD® 100NX™ System is very similar to the current STERRAD® NX System. The STERRAD® 100NX™ System is a large sterilizer with a rectangular chamber size of 152L / 93.4L (usable).

! **CAUTION:** STERRAD® sterilization may cause cosmetic changes to the camera head assembly which do not necessarily impact the camera's functionality.

! **CAUTION:** The camera head assembly must be thoroughly dried before loading into the STERRAD® sterilization chamber. Loads containing moisture may cause a cycle cancellation.

! **CAUTION:** Use only STERRAD®-compatible instrument trays in the sterilization chamber. These trays are specially designed to allow the plasma to surround the items.

! **CAUTION:** Use only polypropylene sterilization wrap and polyolefin pouches*. Do not use paper pouches or sterilization wraps containing wood pulp or cotton.

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1. Clean and prepare the camera head assembly as recommended in the cleaning section of this manual. Be sure that the camera head assembly is completely dry.
2. Load the camera head assembly into the STERRAD® sterilizer, arranging it so that the hydrogen peroxide vapor can surround it. Do not allow the camera head assembly to touch the wall of the sterilizer.
3. Please consult the specific STERRAD® Sterilization System Operator's Manual for detailed instructions for use.
4. Please contact Advanced Sterilization Products [ASP] for the most up to date information regarding sterilization with the STERRAD® Sterilization Systems.

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7.7.3 STERIS Amsco V-PRO Sterilization System

The Amsco V-PRO Low Temperature Sterilization System (i.e., V-PRO 1 and V-PRO 1 Plus) from STERIS® utilizes vaporized hydrogen peroxide (no gas plasma) to inactivate microorganisms. The V-PRO™ 1 Sterilization System is a large sterilizer with a rectangular chamber size of 136L. It is intended for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.

-  **CAUTION:** V-PRO™ sterilization may cause cosmetic changes to the camera head assembly which do not necessarily impact the camera's functionality.
-  **CAUTION:** The camera head assembly must be thoroughly dried before loading into the V-PRO™ sterilization chamber. Loads containing moisture may cause a cycle cancellation.

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CAUTION: Use only V-PRO™-compatible instrument trays in the sterilization chamber. These trays are specially designed to allow the plasma to surround the items.

CAUTION: Use only polypropylene sterilization wrap and polyolefin pouches*. Do not use paper pouches or sterilization wraps containing wood pulp or cotton.

* for US-market: FDA-cleared

1. Clean and prepare the camera head assembly as recommended in the cleaning section of this manual. Be sure that the camera head assembly is completely dry.
2. Load the camera head assembly into the V-PRO™ sterilizer, arranging it so that the hydrogen peroxide vapor can surround it. Do not allow the camera head assembly to touch the wall of the sterilizer.
3. Please consult the specific Amsco V-PRO™ Sterilization System Operator's Manual for detailed instructions for use.
4. Please contact STERIS® for the most up to date information regarding sterilization with the V-PRO™ Sterilization Systems.

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7.8 High Level Disinfection Instructions

! WARNING: The "High Level Disinfection Methods" detailed below are only permitted in the USA.

! WARNING: High level disinfection is recommended ONLY for instruments which come into contact only with intact mucous membranes.

! CAUTION: Any deviations from the recommended disinfection parameters must be validated by the user.

! CAUTION: Before disinfection, the instruments must be thoroughly cleaned, rinsed, and dried.

! CAUTION: To avoid damage to the instruments, do not immerse the devices in disinfectant solution for longer than one hour.

KARL STORZ camera head assemblies can be chemically disinfected using high-level disinfectant solutions containing a 2.4% concentration of glutaraldehyde (e.g. CIDEX, a 14-day glutaraldehyde solution) or 2.0% concentration of hydrogen peroxide (e.g. Resert XL HLD). KARL STORZ does not recommend the use of CIDEX PLUS or other 28-day glutaraldehyde solutions, as they contain high concentrations of surfactants, which may dry and crystallize on the endoscopes if they are not thoroughly rinsed. The crystalline form of the surfactant can become conductive to electricity providing a pathway for arcing. Glutaraldehyde solutions with concentrations greater than 2.4% should be avoided, as a higher percentage of glutaraldehyde may damage the instruments.

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1. Prepare the camera head assembly for disinfection as described in the cleaning section of this manual.
2. Once the camera head assembly has been properly cleaned, prepare the disinfecting solution for use:

CIDEX (2,4% Glutaraldehyde Concentration Solution)

Activate the glutaraldehyde solution by adding the entire contents of activator vial to the solution in the container. Shake well. Activated solution immediately changes color to green, thereby indicating the solution is ready to use. Check the pH of the glutaraldehyde solution to be sure that it is between 8.2 and 8.9 for optimal antimicrobial activity. Do not use activated solution beyond stated 14 day reuse life. Record the date of activation and the expiration date on the container.

Resert XL HLD (2.0% Hydrogen Peroxide Concentration Solution)

No activation is necessary. Use Verify® Chemical Monitoring Strip for Resert XL HLD Solution to confirm hydrogen peroxide concentration before each use. Record the date the solution was poured out of the original container.

3. Pour the desired amount of HLD solution from its original container into a secondary container. Approved plastic containers should be used to avoid scratching of the instruments and to eliminate electrolytic corrosion, which may occur when dissimilar metals are soaked in the same solution.
4. Place the pre-cleaned camera head assembly into approved plastic containers. Do not soak camera head assembly with other instruments to prevent potential damage. Completely immerse the camera head assembly; care must be taken to remove any air bubbles adhering to the surface of the immersed device.

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5. Utilize the following disinfection conditions to achieve manual high-level disinfection:

CIDEX (2.4% Glutaraldehyde Concentration Solution)

Completely immerse the camera head assembly in the undiluted 2.4% glutaraldehyde solution for a minimum of 45 minutes at 25°C (77°F).

CAUTION: DO NOT soak the camera head assembly longer than 45 minutes.

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Completely immerse the camera head assembly in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20°C (68°F).

6. After immersion in the high level disinfection solution, thoroughly rinse the camera head assembly by immersing it completely in a large volume of sterile water (e.g., two gallons). Rinse (immersion) the camera head assembly three times. Each immersion rinse should last for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfecting solution. Use fresh sterile water for each rinse. Thorough rinsing of the camera head assembly with sterile water is essential for removing any residual disinfection solution.
7. Please refer to the disinfectant manufacturer's instructions for use for more detailed information regarding the use of the disinfectant solution, including proper rinsing techniques.
8. Dry the camera head assembly with a lint-free sterile cloth.
9. Special care must be taken to thoroughly rinse the cord connector and carefully dry it with a lint-free sterile cloth. Be sure that the cord connector is completely clean and free of moisture before inserting into the connection socket on the monitor.

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1. The Difficulty of Reprocessing Reusable Rigid Laparoscopic Forceps and other Endoscopic Accessories: Are Disposables the Answer? Health Devices, Vol. 23, Nos. 1-2, pp. 57-58, 1994.
2. Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, No. 1 pp. 23-29, 1995.
3. Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 12th edition. The C.V. Mosby Company, St. Louis, MO 2003.
4. Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-2004.

7.9 Use of sterile drapes

1. Clean and prepare the camera head assembly as recommended in the cleaning section of this manual.
2. Follow the instructions provided by the drape manufacturer when using sterile drapes.

7.10 References for cleaning and sterilization

1. The Difficulty of Reprocessing Reusable Rigid Laparoscopic Forceps and other Endoscopic Accessories: Are Disposables the Answer? Health Devices, Vol. 23, Nos. 1-2, pp. 57-58, 1994.
2. Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, No. 1 pp. 23-29, 1995.
3. Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 12th edition. The C.V. Mosby Company, St. Louis, MO 2003.
4. Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-2004.

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8 Repairs and Warranty**8.1 Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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8.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

8.3 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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8.4 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

8.5 Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

ⓘ NOTE: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen/Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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8.6 Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

- ① **NOTE:** Do not discard as unsorted municipal waste.
- ① **NOTE:** Discard as electrical/electronic waste; recycle or reuse accordingly.
- ① **NOTE:** Consult local authorities for reuse/recycle instructions

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9 Technical Specifications
(for 20290131 / 20290132)

9.1 General technical data

- Camera Head: C-CAM™
- Image Sensor: 1/3"
- Aspect Ratio: 4:3
- Sensor Resolution: 640 x 480
- Min. Sensitivity: 3 lux
- Focal Length: 20 mm fixed
- Dimensions (Length): 3.20" (8.15 cm)
- Dimensions (Height): 1.10" (2.80 cm)
- Dimensions (Width): 1.10" (2.80 cm)
- Weight: 4.16 oz (117.9 g)

9.2 Operating and storage conditions

- Operating temperature: 5°C to 35°C (41°F to 95°F)
- Storage/transport conditions:**
- Ambient Temperature: 10°C to 60°C (14°F to 140°F)
- Humidity 5% to 95% (non-condensing)
- Atmospheric pressure: 50 kPa...+106 kPa

9.3 Standard conformity
(for 20290131 / 20290132)

According to IEC 60601-1, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock: Protection class I
- Degree of protection against electric shocks: Applied part type BF 

9.4 Directive compliance
(for 20290131 / 20290132)

Medical product of class I
This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment.

The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

To be completed by the supplier/importer:
Company stamp/signature:

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 5250-5607, Fax: +52 55 55450174

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KARL STORZ Endoscopy - Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O.Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.kse@karlstorz-rca.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911
Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianle Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 63761066, Fax: +65 63761068
E-Mail: infoasia@karlstorz.com.sg
E-Mail: serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
06-07 e-Centre @ Redhill
Singapore 159471, Singapore
Phone: +65 65325548, Fax: +65 65323832
E-Mail: info@karlstorz.com.sg
E-Mail: servicesg@karlstorz.com.sg

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Baisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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
**Important information for
users of KARL STORZ
instruments**

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**1 Important information for
users of KARL STORZ
instruments**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to help setting up, connecting and operating the C-HUB properly. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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
①	Camera Connector
②	White Balance Button
③	Power Supply Connector
④	Service Button
⑤	USB 2.0
⑥	S-Video

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4 Symbols employed

	USB 2.0 port
	Video output
	Power supply connector
	White Balance
	Caution
	Observe the instruction manual
	Type BF-applied part according to IEC 60601-1
	At the end of its useful operating life, dispose of the unit as electronic scrap.

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1 Important information for users of KARL STORZ instruments	4
2 Photographs of the unit	5
3 Controls, displays, connectors, and their uses	6
4 Symbols employed	7
5 General information	9
5.1 Description of the unit	9
5.2 System Chart	9
6 Safety instructions	10
6.1 Explanation of warnings and cautions	10
6.2 Warnings	11
6.3 Intended use	11
6.4 User qualification	11
7 Operating instructions	12
7.1 Basic equipment	12
7.2 Installation/Interfaces	12
7.2.1 Connect Power Supply	12
7.2.2 Connect S-Video	12
7.2.3 Connect USB 2.0	12
7.2.4 Connect Camera	13
7.2.5 Perform White Balance	13
7.2.6 Service Button	13
8 Maintenance	14
8.1 Cleaning, disinfection and sterilization	14
8.2 Maintenance	14
8.3 Servicing and repair	15
8.4 Repair program	15
8.5 Limitation of liability	16
8.6 Manufacturer's warranty	16
8.7 Disposal	17
9 Technical description	18
9.1 Technical data	18
9.2 Storage/transport conditions	18
9.3 Standard conformity	18
9.4 Directive compliance	18
10 Electromagnetic compatibility (EMC) information	19

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5 General information

5.1 Description of the unit

The C-HUB (20290120) is a device to process and transmit image- and control signals to external devices in medical video endoscopy.

The C-HUB is a plug & play device.

5.2 System Chart

The following system charts shows how to setup the system correctly.

If the interfaces are used other than described in this manual the user/system manufacturer is responsible for the proper and safe operation.

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6.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.

CAUTION: A Caution indicates that particular maintenance procedures or safety precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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WARNING: Read this instruction manual thoroughly before operating the system. Read the following safety instructions carefully to avoid putting your patients, personal or yourself at risk.

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6.2 Warnings

 **WARNING:** The instructions and interface specifications for medical devices used in combination must be observed precisely.

 **WARNING:** Instruments and accessories used in combination with them must be checked for perfect condition, functionality and completeness immediately before and after every application. Always have a spare available.

 **WARNING:** The C-HUB may only be used with accessories, wearing parts, and disposable items which are designated by KARL STORZ as suitable for the camera or the safe use of which is proven. For safety reasons, do not perform unauthorized conversions or modifications to the camera.

 **WARNING:** A reliable, low latency image via USB can not be guaranteed. Therefore the USB port is for documentation purposes only.

 **WARNING:** The auxiliary devices which can be connected up to the analog and digital interfaces (signal inputs and outputs) must be assembled and tested in accordance with the IEC 60601-1 for medical electrical equipment or the IEC 60950 for data processing equipment. In addition, all combinations must comply with the IEC 60601-1-1 system standard. Any person who connects additional equipment to the signal inputs or signal outputs configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements according to IEC 60601-1-1.

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7 Operating instructions

7.1 Basic equipment

- 1 C-HUB
- 1 Power Supply including clips (US, GB)
- 1 USB Cable
- 1 S-Video Cable
- 1 Instruction manual

7.2 Installation/Interfaces

7.2.1 Connect Power Supply

Connect to power supply (20290120-PS) to the unit's receptacle, pushing it firmly into place. The C-HUB will not work without the power supply connected.

7.2.2 Connect S-Video

Connect the monitor (and if applicable additional equipment) to the C-HUB using the S-Video cable.

NOTE: The S-Video delivers a NTSC signal. The monitor needs to support the NTSC format to show a proper image.

7.2.3 Connect USB 2.0

The USB port is suitable for connecting up USB hosts (optional). The C-HUB is identified by the host as a standard USB video device. No drivers need to be installed using Windows XP SP2 or higher.

NOTE: There is different freeware like AMCAP from Microsoft available to visualize and capture the video stream.

Due to countless hardware/software combinations which can not be controlled by KARL STORZ the USB port is designated for documentation purposes only.

Minimal System requirements

Operating System: Windows XP SP2 (or higher)
Processor: Pentium II, 2.4GHz
RAM: 1GB
USB Ports: USB 2.0 HighSpeed (480 Mbps)

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8 Maintenance

8.1 Cleaning, disinfection and sterilization

WARNING: Always pull out power plug before cleaning!

CAUTION: Avoid allowing fluids to enter the unit.

Wipe clean the exterior surfaces of the unit with a low-lint disposable cloth moistened with disinfectant. For wipe-down disinfection, we recommend using only surface disinfectants (no solutions predominantly containing alcohol such as rapid-action disinfectants) and 'disposable disinfectant methods' based on ammonium chloride.

8.2 Maintenance and safety check

8.2.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

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8.2.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements:

- Inspecting the device safety fuses
- Protective ground resistance measured according to IEC 62353 ed1.0: Limit value: 0.2 ohms
- Patient leakage current measured according to IEC 62353 ed1.0: Limit value: 100 µA (BF), 10 µA (CF)
- Leakage currents must not exceed the limit values specified above.

Test for proper functioning:

A functional check must be carried out in accordance with the instruction manual.

The safety check must be entered in the unit log-book and the results of the check must be documented.

If the device does not function reliably and/or safely, it must be repaired.

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9.2 Storage/transport conditions (for 20290120)

- Operating temperature: 10°C ... 40°C
- Storage/transport conditions:**
- Air humidity 5% ... 95% (RH, non-condensing)
 - Temperature -10°C ... 60°C
 - Abs. air pressure +860 hPa ... +1060 hPa

9.3 Standard conformity (for 20290120)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electrical shock: Protection class I
- Degree of protection against electric shocks: Applied part type BF [A]
- Type of moisture protection: IP 20

9.4 Directive compliance (for 20290120)

In accordance with the Medical Device Directive (MDD):

Medical product of class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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10 Electromagnetic compatibility (EMC) information

 **WARNING:** Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

 **WARNING:** The C-HUB model 20290120 is intended for use by healthcare professionals only. In residential areas, the C-HUB model 20290120 may cause EMC interference, in which case it may be necessary to take adequate mitigation measures, such as re-orienting, relocating, or shielding the C-HUB model 20290120 or filtering the connection to the public mains network.

 **WARNING:** The C-HUB model 20290120 should not be used adjacent to or stacked on other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

 **WARNING:** The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

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WARNING: The use of accessories, transducers, and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the C-HUB model 20290120. The cables listed in the following table (Table 200) have been shown to comply with the requirements of EN IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN IEC 60601-1-2.

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NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems, or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device;
- increase the distance between the equipment;
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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Table 200 Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
Power supply AD 5V/2F (PSM11R 050)	No	2	No	Supply
4P mini-DIN cable	Yes	2,5	No	Connecting S video monitor
BB-8003-06	Yes	1,8	No	Connecting PC USB port

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Table 201 Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The C-HUB is intended for use in the electromagnetic environment specified below. The customer or the user of the C-HUB should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidelines
RF emissions CISPR 11	Group 1	The C-HUB uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The C-HUB is suitable for use in domestic establishments.
Harmonic emissions IEC 61000-3-2	(N/A)	not directly connected to the mains
Voltage fluctuations/flicker emissions IEC 61000-3-3	(N/A)	not directly connected to the mains

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Tabelle 202 Guidance and manufacturer's declaration – electromagnetic immunity			
The C-HUB is intended for use in the electromagnetic environment specified below. The customer or the user of the C-HUB should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level, considering IEC 60601-2-18	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (E-SD) IEC 61000-4-2	8 kV contact 15 kV air		Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be > 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines		typical commercial or hospital environment; no mains (approved power supply module) signal/telecomm. lines considered
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth		N.A. no mains (approved power supply module) typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_i (0 V) (>95 % dip in U_i) for 0,5 cycle (10 ms) 40 % U_i (92 V) (60 % dip in U_i) for 5 cycles (100 ms) 70 % U_i (168 V) (30 % dip in U_i) for 25 cycles (500 ms) <5 % U_i (0 V) (>95 % dip in U_i) for 5 s		N.A. no mains
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m		no components sensitive to magnetic fields

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Table 204
Guidelines and manufacturer's declaration – electromagnetic interference immunity for medical electrical devices that are not life-supporting

 The C-HUB is intended for use in the electromagnetic environment specified below.
 The customer or the user of the C-HUB should assure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbance acc. to IEC 61000-4-6 Radiated RF disturbance acc. to IEC 61000-4-3	3 V_{rms} 150 kHz to < 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V_{rms} 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the C-HUB, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5 / 3] \sqrt{P}$ 150 kHz to < 80 MHz $d = [3.5 / 3] \sqrt{P}$ 80 MHz to < 800 MHz $d = [7 / 3] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the nominal power of the transmitter in Watt [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. The field strength of stationary transmitters at all frequencies on site should be, according to a study, less than the conformity level. Interferences may occur in the vicinity of devices that carry the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

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To be filled out by instrument owner:

Return address/company stamp:

Please attach
sufficient
postage

Field of application:

Type of instrument:

Serial no.:

Purchase date:

Signature/Date:

Reply card

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen
Germany

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Absender/Firmenstempel:

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Geräte-Typ:

Serien-Nr.:

Kauf-Dat.:

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Antwortkarte

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen

KARL STORZ GmbH & Co. KG
Mittelstraße 8, 78532 Tuttlingen, Germany
Postfach 230, 78503 Tuttlingen, Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Ohlauer Straße 43
10999 Berlin, Germany
Phone: +49 (0)30 30 69090
Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
2345 Argenta Road, Suite 100
Mississauga, ON, L5N 8K4, Canada
Phone: +1 905 858-0933
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, +1 800 421-0837
Fax: +1 424 218-8526
E-Mail: info@kseca.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2642, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No. 326
Col. Chapultepec Morales
D.F.C.P. 11520, México, Mexico
Phone: +52 5552 5056 07, Fax: +52 55 5545 0174

KARL STORZ Marketing América Do Sul Ltda.
R. Joaquim Floriano, 413 - 20º andar -
Itaim Bibi, São Paulo SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Cervino 4449 Piso 10º
1425 Buenos Aires C. F., Argentina
Phone: +54 11 4772 4545, Fax: +54 11 4772 4433
E-Mail: info@karlstorz.com.ar

KARL STORZ ENDOSKOP NORGE AS
P.O. Box 153
N-2007 Kjeller, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14,
12739 Skärholmens, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Itälähenkatu 23a, 00210 Helsinki, Finland
Phone: 0035-896824774, Fax: 0035-896824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytten 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopia (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customer.service@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B. V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopia Belgium N. V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S. A.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Francia - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 071, Fax: +385 1 6406 077
E-Mail: b.vloj@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskope Greece M.E.P.E.
Sokratous & Kyprou 2, 151 27 Melissia, Greece
Phone: +30 210 61 31 386, Fax: +30 210 61 31 392
E-Mail: sales@karlstorz.gr

KARL STORZ Endoskope Greece Ltd.*
Ipsilantou Str. 32
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

KARL STORZ Endoscopia ROMANIA srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: office@karlstorz.ro

TOV KARL STORZ Ukraine
Prospekt Geroev Stalingrada, 18b
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19, -20
Fax: +380 44 42668-14

ООО KARL STORZ Endoscopy - WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: kste-wostok@karlstorz.ru

KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenicia Street, Mina Elhosn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax: +961 1 365151

KARL STORZ Endoscopy South Africa (Pty) Ltd.
P.O. 6061,
Roggebaai
Cape Town 8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSKOPY Kazakhstan
Korgalynskoe shosse, 2A
Hotel „Duman“
Astana 010000, Kazakhstan
Phone/Fax: +7 7172 79 15 62, 79 15 63
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy - Gulf & Near East
Villa # 7, Mushrif Business Park
Mushrif District, P.O. Box: 30635
Abu Dhabi, United Arab Emirates
Phone: +971-2-4477593, Fax: +971-2-4477594
E-Mail: salesadmin.ksegne@karlstorz-rce.com

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 - 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595
Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 (0)24 6258 9911
Fax: +86 (0)24 6258 9922
E-Mail: info@karlstorz.com.cn

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