Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 • 78532 Tuttlingen
Tel.: 07461/98372 • Fax: 07461/98358

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**1 Important information for
users of KARL STORZ
devices**

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the CALCUSPLIT®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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


CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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**3 Controls, displays,
connectors, and
their uses**

- ① Power switch ("O" = off)
- ② Pressure gauge for displaying the preselected compressed air force
- ③ Function display for single-shot operation
- ④ Function display for continuous-shot operation
- ⑤ Socket for footswitch
- ⑥ Socket for handle (compressed air outlet)
- ⑦ M button for preselecting the compressed air pulse frequency
- ⑧ Knob for preselecting the compressed air force
- ⑨ Nozzle for compressed air feed
- ⑩ Potential equalization connector
- ⑪ Power fuse holder
- ⑫ Power supply cord receptacle
- ⑬ Handle
- ⑭ Connection for compressed air
- ⑮ Lithotripsy probe
- ⑯ Adaptor for handle

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4 Symbols employed

	Follow instructions for use
	Read the instructions carefully before operating the equipment!
	Output pressure in bar
	Power off
	Change in outlet pressure
	Socket for handpiece
	Footswitch connection point
	Single shot operation
	Continuous shot operation
	Potential equalization connector
	Applied part of type BF
	Alternating current
	Environment protection use period of 50 years (China RoHS)
	Disposal – This device has been labelled in Accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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1 Important information for users of KARL STORZ devices	III
2 Images of the unit	IV
3 Controls, displays, connectors, and their uses	V
4 Symbols employed	VI
5 General information	3
6 Safety instructions	4
6.1 Explanation of warnings and cautions	4
6.2 Intended use	8
6.3 Normal use	8
6.4 User qualification	9
6.5 Patient profile	9
6.6 Use	9
6.7 Safety precautions at the site of installation ..	10
6.8 Safety precautions when operating the unit ..	10
7 Installation and operating instructions	11
7.1 Unpacking the equipment	11
7.2 Basic equipment	11
7.3 Installing and connecting up the unit	11
7.4 Setting of the system voltage	13
7.5 Compressed air supply	13
7.6 Lithotripsy probes	15
7.7 Operating the unit	17
7.8 Test for proper operation	18
7.9 How to use	19
7.10 Decommissioning	20
8 Maintenance	21
8.1 Fuse replacement	21
8.2 Replacing the sealing rings or the lever set on the handle	22
8.3 Testing the system	25
8.4 Cleaning, disinfection and sterilization	26
8.5 Table of preparation procedures	27
8.6 Cleaning, disinfection and sterilization of reusable instruments/tubing	27
8.7 Preparation	28
8.8 Sterilization	31
8.9 Maintenance	33
8.10 Safety checks and inspections	33
8.11 Servicing and repair	34

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5 General information**Device description**

The CALCUSPLIT® is a pneumatic lithotripter for intracorporeal lithotripsy. The kinetic energy of a projectile driven by compressed air is transmitted by an elastic shock to the proximal end of a probe via a special transformation element. Due to the shock motion of the distal probe part, the stone with which it is in contact is fragmented.

The CALCUSPLIT® is suitable for the lithotripsy of cystic, urinary, and renal calculi.

For the application, probes with various lengths and in various diameters are available, which permits adaptation to the stone size as well as to the available lumen of the working channel.

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6 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

6.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the signal words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

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WARNING: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

WARNING: Keep out of reach of patients.

WARNING: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1 and Chap. 16, IEC 60601-1).

WARNING: Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor.

WARNING: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

WARNING: The device may be subject to sudden failure during use. Remove the probe immediately! In case this happens, you should plan to have a replacement device at the ready or employ alternative procedures.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

WARNING: Danger of injury! Do not connect the handle to the control unit if the handle has not been completely assembled. In case of inadvertent triggering of a shock wave, the projectile can cause injuries when it leaves the guiding tube. Before every disassembly of the handle, remove the silicone tube from its connection on the control unit.

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WARNING: Risk of injury! Tube system is under pressure. Hold the tube still when disconnecting it from the connecting piece ⑨.

WARNING: Always unplug the unit before doing any maintenance work on it (e.g. cleaning).

WARNING: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to qualified personnel. Removal of covers by nonauthorized personnel will void the equipment's warranty.

WARNING: To ensure protection against electrical shock, only original KARL STORZ parts may be used for connecting the handle control unit and compressed air system control unit.

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
WARNING: Danger of injury! Although the probe only exerts blunt force, there is a certain risk that particularly the ureter may be punctured in the event of improper handling. The tip of the probe should therefore not extend further than 10–20 mm past the distal end of the tube. During the lithotripsy and when moving the unit, always ensure that there is sufficient visual control of the probe tip.

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
WARNING: Before and after each application the lithotripsy wire probes used must be inspected for possible damage, particularly deformation of the collar edge and any formation of burr or chips. If damage is found, the probe must not be used again. Generally speaking, the lithotripsy wire probes with silicone damper have to be changed after a maximum of 5 applications or after damage.

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WARNING: Never connect the pressure tube of the CALCUSPLIT® control unit to the suction nozzle handpiece (see image) or of any other instrument used. For handpieces with stopcock, the latter **must** be closed, i.e. the stopcock must point sideways.

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CAUTION: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

CAUTION: When replacing fuses, use only fuses of the correct rating.

CAUTION: The CALCUSPLIT® may only be used with accessories, wearing parts, and disposable items which have been designated as suitable for the unit by KARL STORZ.

CAUTION: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

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NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

NOTE: This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

NOTE: Make certain to observe national regulations/laws regarding disposal.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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6.2 Intended use

The CALCUSPLIT® lithotripsy control unit is used to create energy for the destruction of stones during endoscopic urological applications.

The lithotripsy handpiece is used to transfer energy to the lithotripsy probe for destroying stones during endoscopic urological applications.

The probe is used to transfer energy from the lithotripsy handpiece onto the stone during endoscopic urological applications.

6.2.1 Indication

Disintegration of stones in the urinary tract via endoscopic access ports.

6.2.2 Contraindication

Use of the CALCUSPLIT® is contraindicated, if, in the opinion of the attending physician, the general condition of the patient does not permit the use of the CALCUSPLIT® for lithotripsy.

6.2.3 Possible side effects

- The occurrence of abnormal openings and narrowing of the ureter after surgery;
- Separation of the ureter;
- The movement of kidney stones;
- Damage to the mucous membrane;
- Piercing the ureter;
- The need for open surgical intervention;
- Pain in the side;
- Pelvic discomfort;
- Macroscopic hematuria true;
- Urinary tract infection.

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6.3 Normal use

The CALCUSPLIT® intracorporeal lithotripter, which is comprised of a control unit, handle and appropriate lithotripsy probes, is suited for the endoscopic lithotripsy of calculi (cystic, urinary, and renal calculi). The CALCUSPLIT® may only be used in conjunction with the appropriate endoscopes*.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

6.4 User qualification

The CALCUSPLIT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

6.5 Patient profile

Gender: irrelevant
Age: irrelevant
Weight: irrelevant
Nationality: irrelevant
State of health: suitable for treatment in the physician's judgment

6.6 Use

The device may only be used in hospitals and doctors' offices in the ambient conditions given in the technical data (page 38).

Other conditions of use:

- Frequency of use up to several times a day
- Length of use from a few minutes to several hours a day
- Position it on a vibration-free, fixed base
- Can be moved only by installing it on an equipment trolley
- Can be swiveled only by installing it on ceiling supply units

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6.7 Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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6.8. Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment with the CALCUSPLIT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the unit.

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7 Installation and operating instructions

7.1 Unpacking the equipment

Carefully unpack the CALCUSPLIT® and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

7.2 Basic equipment

- 1 Calcusplit® 27 6300 20
- 1 Handle 27 6300 38
- 1 Fabric reinforced tube for connecting CALCUSPLIT® to the central compressed air supply, 27 6300 21
- 1 Silicone tube for connecting the handle to CALCUSPLIT®, 27 6300 40
- 1 Footswitch 20 0142 30
- 1 Package with 5 sets of sealing rings, 27 6301 80
- 20 Replacement damping elements, 27 6300 35
- 1 Cleaning brush 27 650 C
- 1 Instrument oil 27 656 B, 50 ml
- 1 Power cord 400 A
- 1 Instruction manual

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
7.3 Installing and connecting up the unit


NOTE: The CALCUSPLIT® including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

WARNING: This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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
The unit is equipped with a connector ⑩ for attaching a ground line.

The unit's ground line should be installed by a qualified electrician.

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7.4 Setting of the system voltage

Remove the power fuse holder (11) and take out the voltage selection module.

Insert the voltage selection module into the power fuse holder in such a way that the system voltage displayed in the viewing window on the back of the power fuse holder corresponds to the local power supply.


Make certain that the safety fuses which are used are appropriate for the selected system voltage.

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Connect power cord. Insert power cord into power cord receptacle (12) as far as it will go.

! WARNING: Operate the unit only with the system voltage stated on the rating plate.

! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

! WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

i NOTE: Please ensure that access to the power plug is ensured at the site of installation.

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7.5 Compressed air supply

! WARNING: The pressure of the compressed air used at the unit input must be between 3.5 bar (minimum pressure) and 5.0 bar (maximum pressure).

This pressure must be checked at the compressed air feed system used! If the pressure is higher, an excess pressure safety device is triggered in the unit. If necessary, a pressure reducer will have to be used.


! WARNING: Only dry, oil-free compressed air filtered to $\leq 5 \mu\text{m}$ may be used. Under no circumstances should CO_2 gas be used for the compressed air supply.

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
Attach the connecting nozzle for compressed air supply ⑦ to the compressed air wall outlet (3.5 bar – 5 bar) using the fabric reinforced tube.

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
Alternatively, a compressor system can be used to supply compressed air. In such a case we recommend the compressor system 20031001 (see accessories). Output pressure must be regulated on the compressor system so that it cannot exceed 5 bar.

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
Plug the connecting cable for the footswitch into the socket ⑤.

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7.6 Lithotripsy probes

WARNING: Never connect the pressure tube of the CALCUSPLIT® control unit to the suction nozzle handpiece (see image) or of any other instrument used. For handpieces with stopcock, the latter must be closed, i.e. the stopcock must point sideways.

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Preassembly of the wire probe

Always slide the silicone damper from the distal direction up to the limit point of the wire probes.

WARNING: In the event of damage or after the probe has been used for a maximum of five times, it must be replaced with a new one.

WARNING: In the event of damage or after the preassembled wire probe with distally assembled silicone damper (Art. no. 27 6300 35) has been used for a maximum of five times, it must be replaced with a new one.

WARNING: Never use the wire probe without the silicone damper!

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Inserting wire probe into handpiece

Screw the conical cap off the handpiece. Insert the preassembled wire probe with distally assembled silicone damper (see image) through the cap. Screw the cap back onto the handpiece and, at the same time, ensure that the proximal end of the wire probe passes into the hole of the handpiece.

WARNING: Never use the handpiece without the wire probe!

WARNING: Normal functioning of the wire probe can only be guaranteed in conjunction with the distally assembled silicone damper.

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Attach the silicone tube to the connector of the handle all the way and connect the other end to the compressed air outlet ⑥ on the CALCUSPLIT®.

WARNING: Danger of injury! Do not connect the handle to the control unit if the handle has not been completely assembled. In case of inadvertent triggering of a shock wave, the projectile can cause injury when it leaves the guiding tube.

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(Only applicable to rigid KARL STORZ uretero-renaloscopes)

Push the handpiece adaptor over the 'instrument inlet' of the uretero-renaloscope until it clicks into place.

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Then lock the adaptor using the slide.

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The lithotripsy probe is now guided through the working channel of the endoscope. When the adapter is completely fitted into the handle, the probe should extend a maximum of 10–20 mm past the distal end of the working channel.

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
Preselect the force of the compressed air pulses using the knob ⑧. The preselected value is displayed on the pressure gauge ②.

① NOTE: The maximum pressure setting for compressed air is limited to 2.0 bar. To guarantee normal operation a pressure of at least 1 bar is necessary.

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With the M button ⑦, select the desired frequency of the compressed air pulses. According to the preselected frequency, either Display ③ (single-shot) or Display ④ (continuous-shot) lights up.

① NOTE: In the continuous-shot operating mode, approx. 12 shots per second are emitted.

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
The compressed air pulses can be triggered by pressing the footswitch.

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7.8 Test for proper operation

With the handle disconnected from the unit and the unit switched on, press the footswitch. From socket ⑥, a compressed air pulse (single-shot) or a series of 12 pulses (continuous-shot) should now be emitted.

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Testing the system

Before each treatment the CALCUSPLIT® system must be tested for proper functioning. To do this, press the footswitch while CALCUSPLIT® is switched on. The wire probe must make an axial beat movement. **This must be verified!** If the wire probe does not move with a jerk, the system must be checked, also for correct assembly. **A system that is not fully functional must not be used.**

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7.9 How to use

WARNING: Each treatment with the CALCUSPLIT® system and wire probe should be controlled endoscopically. In the event of a malfunction, the treatment should be stopped and the CALCUSPLIT® system inspected.

The lithotripsy probe is guided through the working channel of a rigid endoscope to the stone which is to be crushed. In order to fragment the stone, the probe must be set directly upon it. It is advantageous to hold the stone fast with a grasping forceps during the procedure.

The lithotripsy probe should, when it is completely inserted, extend a maximum of 10 to 20 mm past the distal end of the working channel. This makes it easier for the physician to determine the position of the probe, guarantees sufficient visual control of the probe tip and reduces the danger of possible perforation due to uncontrolled movement.

Lithotripsy probes with a small diameter (0.8 mm) can be inserted into a slightly curved working channel. However it must be taken into account that, in bending the probe, energy is lost. Nevertheless this energy loss can virtually always be compensated by an increase in the strength of the compressed air pulses.

It is recommended that the lithotripsy begin in single-shot mode and with a compressed air pulse strength of approx. 1 bar. According to the consistency and resistance of the stone, the pressure can then be successively increased or switched over to continuous-shot operation.

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Although the probe only exerts a blunt force upon the mucous membrane and there is no fear of heat development, constant visual control of the probe tip must be ensured. Particularly with small stone fragments, proceed with increased care. It is recommended that the lithotripsy be ended when the stone fragments have reached a size which facilitates manual extraction or makes natural passage probable. With increasing fragmentation, reduce the irrigation pressure as much as necessary, so that no unwanted movements of the calculus fragments occur.

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

WARNING: When working in continuous-shot operation, the progress of the lithotripsy and position of the stone and of the probe tip should be checked regularly every 1–2 seconds (i.e. after 12–24 pulses) for safety reasons.

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7.10 Decommissioning

Switch off the CALCUSPLIT® at power switch ①.

Remove the instrument tube from the front of the device and carry out preparation on the handpiece as well as the tube.

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
WARNING: Risk of injury! Tube system is under pressure (even after switching off the compressed air supply). Hold compressed air hose still when removing it from connecting piece ②.

Remove compressed air hose from the back of the device.

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8 Maintenance


8.1 Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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
Remove the power fuse holder ⑪ with a screwdriver or other tool.

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! CAUTION: Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

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
Replace the power fuse holder ⑪.
Reconnect the power cord.
Test the unit for proper operation.

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8.2 Replacing the sealing rings or the lever set on the handle

The four sealing rings in the handle must be replaced regularly (after approx. 30 sterilizations). To do this, the handle has to be disassembled.

In performing this procedure the transmission lever must be examined. If there are any signs of wear the complete lever set has to be replaced (see page 24).

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
! WARNING: Danger of injury! Never disassemble the handle while it is still connected via the silicone tube to the control unit. By inadvertent triggering of a shock wave, the projectile can cause injuries when it leaves the guide tube.

First disconnect the silicone tube from the handle connection, then take off the front cap and remove the lithotripsy probe.

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Unscrew the rear cap (with tube nozzle) and pull it, with the guide tube, out of the handle sleeve.

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
Pull the guide tube out of the cap and take out the projectile.

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i NOTE: If the projectile is sluggish (this would be noticeable during operation by a decline in performance), the guide tube and projectile must then be cleaned or changed.

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Unscrew the handle sleeve from the handle head (see arrow).


WARNING: The metal plate (lever bearing) fastened with screws in the handle head may only be removed by trained personnel.

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(see also exploded drawing on page 24)

Rear cap: 2 sealing rings

- 1 sealing ring on the connector to the guiding tube and

- 1 sealing ring on the exterior threading.

Handle head: 2 sealing rings

- 1 sealing ring on the threading for the front cap
- 1 sealing ring on the threading for the handle sleeve.

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
CAUTION: When replacing the sealing ring on the attachment for the guiding tube in the rear cap, make certain that the surface does not become scratched.

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a-d Set O-Ringe, Packung zu 5 Sets, 27630180

8.2.2 Checking the transmission lever/ changing the lever set 276300 83

Unscrew the metal plate (lever bearing) fastened with flat head screws. Remove the plate (lever bearing) (see exploded drawing).

Check the lever. If there are traces of wear in evidence, change the complete lever set 27 6300 83.

WARNING: This work must only be done by trained personnel.

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
After replacing the sealing rings, reassembly is performed in reverse order.

CAUTION: Lightly oil the projectile, guide tube and metal plate (lever bearing) with instrument oil 27656 B.

WARNING: Make absolutely certain that the plate has been remounted correctly.

When inserting the guide tube through the handle sleeve, make certain that the groove on the guide tube is locked into the transmission lever. To do this, move the guide tube lightly to and fro while pushing forward until it engages and can no longer be turned. Then insert the projectile into the guide tube (the projectile must not stick).

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
WARNING: In the event of damage or after the probe has been used for a maximum of five times, it **must** be replaced with a new one.

WARNING: In the event of damage or after the damper (Art. no. 27630035) has been used for a maximum of five times, it **must** be replaced with a new one.

NOTE: Regardless of this, always use a new damping piece whenever a new probe is inserted.

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
8.3 Testing the system

After any maintenance work and before each treatment, the proper functioning of the CALCUSPLIT® System must be tested. To do this, press the footswitch while CALCUSPLIT® is switched on. The wire probe must make an axial beat movement. **This must be verified!** If the wire probe does not move with a jerk, the system must be checked, also for correct assembly. **A system that is not fully functional must not be used.**

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! WARNING: The CALCUSPLIT® handpiece is pressure-driven. Before each use the system (consisting of a CALCUSPLIT® handpiece and instrument, such as a ureterorenoscope) must be checked for leakage when assembled. The system is to be immersed in a container with sterile liquid (left-hand image). When the system is activated by pressing the footswitch in multiple-pulse mode, no air bubbles must exit the instrument channel from distal (see drawing). **If air can escape from the instrument channel, the system must not be used** as air insufflation via the instrument (e.g. a ureterorenoscope) can lead to safety problems.

! WARNING: Always ensure that there is a seal on the proximal instrument port (e.g. on the ureterorenoscope). Each treatment with the CALCUSPLIT® system and wire probe should be **controlled endoscopically**. In the event of a malfunction, the treatment should be stopped and the CALCUSPLIT® system inspected.

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8.4 Cleaning, disinfection and sterilization

i NOTE: The instructions for "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

! WARNING: Always unplug power plug before cleaning!

! CAUTION: Avoid allowing liquids to enter the unit.

Wipe-down disinfection of device and footswitch

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility must be observed.

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WARNING: Risk of infection/Product damage: KARL STORZ has validated the procedures given in the instruction manual. Only use these procedures. The preparation process must be validated by the on-site operator.

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8.6 Cleaning, disinfection and sterilization of reusable instruments/tubing

WARNING: The local guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed when performing all cleaning and disinfection procedures on contaminated instruments/tubing.

WARNING: Reusable instruments/tubing must be cleaned, disinfected and sterilized using validated preparation procedures before use and subsequent reuse.

WARNING: Prior to the sterilization of reusable instruments/tubing, particularly in cases of manual preparation, it must be ensured that all chemical residue has been rinsed off or removed (see also the preparation steps).

WARNING: Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

NOTE: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices.

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CAUTION: When preparing and using the solutions, follow the chemical manufacturers' instructions, paying close attention to proper concentration and exposure times. Prolonged immersion can lead to material changes.

8.6.1 Water quality requirements

Tap water can be used to prepare the cleaning and disinfectant solutions. When rinsing for the last time for example to remove any chemical residue, microbiologically pure/sterile water must be used.

8.6.2 Preliminary cleaning of reusable tubing

Disconnect the silicone tube (with quick-release coupling) from the unit after use.

Place the reusable tubing in a container filled with a cleaning solution (as per the manufacturer's instructions) immediately after use in order to prevent contaminants from drying on the surface of the tubing.

8.7 Preparation

8.7.1 Manual preparation of reusable tubing

WARNING: The local guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed when performing all cleaning and disinfection procedures on contaminated instruments.

1. Immerse reusable tubing in cleaning solution (as per the chemical manufacturer's instructions).
2. Using a syringe, rinse out the tubing underneath the (water) surface. If necessary, repeat the procedure several times. Remove soiling with disinfected brushes/sponges, soft, low-lint disposable cloths or a cotton applicator.
Cleaning accessories are available from KARL STORZ (see catalog).

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3. After cleaning, rinse the tubing and place it in a suitable container with disinfectant solution, making sure that all air bubbles have been removed. Ensure that the tubing is completely covered in the solution for the duration of the exposure.
4. Cleaning brushes must be cleaned and disinfected immediately.
5. After disinfecting the tubing, remove it and rinse with microbiologically pure/sterile water. Ensure that all chemical residue in and on the tubing is removed.
6. Dry the tubing with sterile compressed air.
7. After preparation, inspect the tubes for cleanliness and damage. If soiling can still be detected, the preparation procedure must be repeated.

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8.7.2 Machine preparation of reusable tubing

For machine preparation (up to 93°C), connect the reusable tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed out.

Prior to machine preparation, a (manual) preliminary cleaning is generally required.

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! WARNING: Check the silicone tube for leaks after preparation. This makes it possible to detect a leak early, through which air could escape.

To check for leaks, use leakage tester 13242 XL for example. For this purpose, connect one end to the leakage tester and block off the other end, with your finger for example. If there is leak in the tubing, a pressure drop will be observed on the leakage tester.

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8.7.3 Getting handpiece/probe ready for preparation

Take off the front cap (probe holder) and remove the lithotripsy probe.

For preparation, disassemble the handpiece as far as it is possible to do so without tools (see section 'Replacing the sealing rings ... on the handpiece', page 22).

NOTE: The metal plate (lever bearing) is not removed for preparation.

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8.7.4 Manual cleaning of handpiece and probe

To clean the suction channel (if present), use the supplied cleaning brush (27650 C). Also clean the guide tube with a soft disinfected cleaning brush (Art. no. 27650 B).

CAUTION: Make certain that the inner surface of the guide tube is not scratched.

NOTE: If the projectile is hard to move, the guide tube and projectile must be replaced.

After cleaning (pay attention to note on page 22), rinse out the probe and the handpiece's individual parts thoroughly using microbiologically pure/sterile water, in order to remove all chemical residue. Then dry the parts carefully using sterile compressed air.

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8.7.5 Machine preparation of handpiece/probe

Due to the partially strong motion of the water, it is necessary that all parts be fixated in holders, which prevent damage. For sufficient rinsing out of the instruments, instrument channels must be correctly connected.

If necessary, dry the handpiece/probe carefully using sterile compressed air.

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8.7.6 After-treatment following manual/ machine preparation – handpiece

Oil the projectile and guide tube with the supplied instrument oil 27656 B. Also apply a drop of the instrument oil to the metal plate (lever bearing).

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8.8 Sterilization

NOTE: Before sterilization, the components must be cleaned and disinfected (see table, page 27):

- mechanically (manual) and/or
- washer/disinfector (thermal disinfection at 93 °C, device-specific)

Then pack the components into appropriate containers (packaging according to EN 868, part 1).

NOTE: The handpiece is steam sterilized as an entire part but without the front cap (probe holder) and with the stopcock open (if present).

The probe and handpiece are steam sterilized separately at 134 °C (+ 3 °C).

NOTE: Due to the design of the handpiece, we recommend the 'fractionated prevacuum method'.

NOTE: Use the fractionated prevacuum method (according to ISO 17665-1), taking into account the applicable national requirements.

Steam sterilize the silicone tube (with quick-release coupling) at 134 °C (+3 °C) using the fractionated prevacuum method.

WARNING: Successful sterilization is only possible on clean, disinfected surfaces.

WARNING: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

WARNING: Any deviations from the recommended parameters for sterilization should be validated by the user.

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8.8.1 Steam sterilization

KARL STORZ has validated the following sterilization method:

Fractionated prevacuum method (steam sterilization)

! WARNING: Tubing sets may only be sterilized using the fractionated prevacuum method.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

On completion of the steam sterilization cycle, all sterilized parts should be allowed to cool down slowly.

! WARNING: Risk of burns! Parts are hot after steam sterilization. Allow them to cool down!

8.8.2 Fractionated prevacuum method

The fractionated prevacuum method comprises four phases.

- **Conditioning phase**
The conditioning phase removes air from the chamber by creating a vacuum up to four times. Then steam is blown into the chamber.
- **Sterilization phase**
Sterilization takes place at a temperature of 134 °C (+3 °C) and a pressure (absolute) of 300 kPa (3 bar) over a period of 5 minutes (minimum time).
- **Evacuation phase**
The steam is released.
- **Drying phase**
Drying takes place under renewed application of a vacuum over a period of approx. 5–8 minutes.

Validated parameters for the fractionated prevacuum method:

Temperature:	134 °C (+3 °C)
Pressure (p _{absolute}):	300 kPa (3 bar)
Exposure time:	5 minutes for all instruments.

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8.9 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, a safety check must be carried out on the unit at least once a year.

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8.10 Safety checks and inspections

KARL STORZ stipulates that safety checks must be carried out once a year, in order to prevent accidents due to aging, wear or failure of medical devices.

If any defects or shortcomings which could endanger patients, operating personnel or others are ascertained during these safety checks, the unit must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

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In accordance with the Medical Device Directive (MDD) this instruction manual contains information on the type, extent and frequency of safety checks.

When carrying out safety checks, the following items must be tested:

8.10.1 Visual inspection

- Housing and accessories
- Labelling, Identification plate data
- CE symbol, KARL STORZ test label
- Instruction manual
- Line fuses
- Fuse label
- Voltage preselection

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8.10.2 Performance check

- Mains power switch
- Buttons
- Pressure regulation
- Output pressure display
- Output pressure
- Footswitch
- Footswitch socket
- Drawer for the short instruction manual

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8.10.3 Electrical safety measurements (IEC 62353, ed. 1.0)

- Protective ground resistance: $\leq 0,3 \text{ Ohm}$
(with power cord)
- Earth leakage current: $\leq 1,0 \text{ mA}$
- Housing leakage current: $\leq 0,5 \text{ mA}$
- Patient leakage current: $\leq 0,5 \text{ mA}$
- Patient leakage current: $\leq 5,0 \text{ mA}$
(line voltage on the applied part)

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NOTE: For detailed information on the range and implementation of the safety checks please refer to the Service Manual SV 3345 (English version).

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8.11 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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8.12 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as ~~electronic scrap~~.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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8.13 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

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In other countries please contact your local KARL STORZ branch or authorized dealer.

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8.14 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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8.15 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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8.16 Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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9 Technical description

9.1 Troubleshooting

WARNING: Always unplug the unit before all maintenance work.

Symptom:

Complete failure of the unit.
LED for operating mode does not light up.

Possible causes:

No power from the power line.
Defective power fuse.

Remedy:

Check that there is electricity to the wall outlet.
Change fuses as described in the instruction man-
ual. Make sure you fit the correct type of fuse.

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9.2.1 Standard compliance
(for 27 6300 20)

According to IEC 60601-1, UL 60601-1,
CAN/CSA C22.2 601.1-M90:

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type BF 

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility
Information in the appendix (p. 44-57).

9.2.2 Directive compliance
(for 27 6300 20)

According to Medical Device Directive (MDD):

This medical product belongs to Class IIb

This medical product bears the CE mark in
accordance with MDD 93/42/EEC.

 **NOTE:** The code number after the CE mark
indicates the responsible notified body.

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9.3 Technical documentation

On request, the manufacturer will provide
those circuit diagrams, itemized parts listings,
descriptions, sets of adjustment instructions, and
other items of available documentation to suitably
qualified user personnel duly authorized by the
manufacturer for their use in repairing those com-
ponents of the unit that have been designated as
reparable by their respective manufacturers.

Supply of such technical documentation relating
to the unit shall not be construed as constituting
manufacturer's authorization of user's personnel,
regardless of their levels of technical training, to
open or repair the unit.

Explicitly exempted herefrom are those mainte-
nance and repair operations described in this
manual.

*We reserve the right to make engineering modifi-
cations in the interest of promoting technological
progress and generating performance improve-
ments without obligation on our part to submit prior
notice thereof.*

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
**10 Spare parts,
recommended accessories**

10.1 Spare parts list

Item	Cat. no.
Fabric reinforced tube for connecting CALCUSPLIT® to the central compressed air supply, length 400 cm	276300 21
Handle, autoclavable	276300 38
Sealing rings, package with 5 sets	276301 80
Lever set	276300 83
Cleaning brush, diameter 2.5 mm, length 35 cm	27650 C
Instrument oil, 50 ml	27656 B
Replacement damping elements (20 pcs), autoclavable	276300 35
Silicone tube, for connecting the handle 276300 38 to CALCUSPLIT®, length 200 cm, autoclavable	276300 40
Footswitch	200142 30
Power fuse T 0.125 AL250V (230/240 V~), pack of 10	1078290
Power fuse T 0.2 AL250V (100/120 V~), pack of 10	1067800
Power cord (grounded)	400 A
'Hospital Grade' Power cord (USA)	400 B
Instruction manual	96136004 DER

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Item	Cat. no.
Adaptor for KARL STORZ Ureterorenoscopes and Ureterosopes, for use with CALCUSPLIT®, autoclavable	27 6300 36
Compressor, 3,5-5 bar, power supply: 230 V~, 50/60 Hz, consisting of: 20 0310 20 Compressor 27 6300 22 Fabric reinforced Tube, for connecting CALCUSPLIT® to compressor	20 0310 01
Same, power supply: 115 V~, 50/60 Hz	20 0310 01 C

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11 Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The CALCUSPLIT® model 27 6300 20 corresponds to IEC 60601-1-2:2007 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The CALCUSPLIT® model 27 6300 20 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the transmitting device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The CALCUSPLIT® model 27 6300 20 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the CALCUSPLIT® model 27 6300 20. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:

Type	Shielded	Length [m]	Ferrite	Used for
20010230/20014230	Yes	3	No	Footswitch
Power cord	No	3	No	Power supply

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions

CALCUSPLIT® model 27 6300 20 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCUSPLIT® model 27 6300 20 should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	CALCUSPLIT® model 27 6300 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The CALCUSPLIT® model 27 6300 20 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity

CALCUSPLIT® model 27 6300 20 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCUSPLIT® model 27 6300 20 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

Table 4

 Guidance and manufacturer's declaration – electromagnetic immunity –
 for equipment and systems that are not life-supporting

CALCUSPLIT® model 27630020 is intended for use in the electromagnetic environment specified below. The customer or user of the CALCUSPLIT® model 27630020 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the CALCUSPLIT® model 27630020, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ $d = [3.5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 6

**Recommended separation distances between portable and mobile RF communications equipment and the
CALCUSPLIT® model 27 630020**

The CALCUSPLIT® model 27 630020 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m], according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326, Col. Chapultepec Morales
D.F.C.P. 11520, México, México
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012, Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

KARL STORZ Endoscopy Kasachstan LLP
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) - Units 205 - 206
Dubai Health Care City, P.O.Box 118069
Dubai - United Arab Emirates
Phone: +971 (0)4 4471230
Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty Ltd
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang 1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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