Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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This instruction manual is intended to aid you with proper connection and operation of the TIPCAM®1 S 3D in conjunction with the instruction manual for the IMAGE1 S. All required details and all actions on your part are clearly presented and explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment. Insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the IMAGE1 S instruction manual for the individual functions of the CCU that can be used with the TIPCAM®1 S 3D camera.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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
3 Controls, displays, connectors, and their uses

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Read accompanying documents before use

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8	Reprocessing not validated for US customers	19
8.1	General care.....	20
8.1.1	Care and handling	20
8.2	General warnings	23
8.3	Accessories.....	24
8.4	Disassembly.....	24
8.5	Preparation for cleaning and disinfection	24
8.6	Manual precleaning	24
8.7	Manual cleaning	25
8.8	Manual disinfection.....	25
8.9	Machine cleaning and disinfection.....	25
8.10	Assembly, inspection and care.....	27
8.11	Packaging systems	27
8.12	Sterilization.....	27
8.12.1	Accessories	29
8.13	Reprocessing limits	29
9	Reprocessing for US customers only	30
9.1	General care.....	30
9.1.1	Care and handling	30
9.2	General warnings	32
9.3	Disassembly	33
9.4	Pre-cleaning	33
9.5	Manual cleaning	34
9.6	Inspection after cleaning.....	35
9.7	Sterilization.....	36
9.8	Limits of reprocessing	36
10	Maintenance	37
10.1	Service and maintenance	37
10.2	Servicing and repair.....	38
10.3	Disposal	38
10.4	Repair program	38
10.5	Limitation of liability	39
10.6	Warranty.....	40

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5 General information**5.1 Device description**

The TIPCAM®1 S 3D is a high-resolution videoendoscope with two distally integrated video chips.

Since the TIPCAM®1 S 3D does not require an optical imaging system, it is very sturdy.

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6 Safety instructions

6.1 Warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of the device.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A Note indicates special information about operating the unit.

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WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

WARNING: Also pay particular attention to the instruction manual for the IMAGE1 S CCU no. 96206286DER.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

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WARNING: Risk of injury to the patient/user: Damaged or incomplete instruments/accessories may cause malfunctions. Before and after every use, always check the entire system (e.g., the TIPCAM®1 S 3D, all cables, the CCU, monitor and accessories used) to ensure that everything is configured properly, in perfect working order and nothing is missing/damaged. Replace incomplete/ damaged instruments/accessories and do not continue to use them.

- Before every use, check that all the instruments required for the operation are functioning correctly.
- Be certain to check the cover glass of the TIPCAM®1 S 3D. In addition, the camera cable and the light cable need to be inspected for damage cracks.
- Do not use wet/moist cables.
- When attaching components to an instrument, the secure position of the components must be checked by gently pulling and/or turning them in both directions before use.

WARNING: Check the instrument and accessories used in combination before every application for unintentional rough surfaces, sharp corners, burred edges, and projecting parts, which can endanger personal safety. Damaged, defective, bent, or incomplete instruments or accessories must not be used further under any circumstances.

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WARNING: Risk of injury: Overloading the instrument by exerting too much force on it may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original positions.

WARNING: Risk of injury and risk of damaging the products: When using applied parts outside of the field of vision there is a risk of unintentionally injuring tissue and damaging accessories. Always keep the applied parts of the active electrode and the laser and other instruments which transmit energy pointed at the target in the field of vision during use.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

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WARNING: Risk of injury and risk of damaging products: The high level of light energy means that the distal end, the light port, adjacent components and tissue in front of the light outlet window to heat up. This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Avoid excessively long (> 10 min) continuous direct tissue contact with the shaft location.
- Never allow the endoscope to rest on the patient or come into direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected fiberoptic light cable or an endoscope. Always wear standardized protective equipment when using lasers.

WARNING: The TIPCAM®1 S 3D must be replaced if the image is cloudy or if no image or only parts thereof are visible.

WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the TIPCAM®1 S 3D is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the handle, which can cause cardiac pacemakers or other electronic devices to malfunction.

WARNING: The video system can malfunction during use. An adequate replacement system must be kept available.

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-  **WARNING:** In the event of image failure immediately remove TIPCAM®1 S 3D from the patient.
-  **WARNING:** The use of HF probes, HF electrodes, HF instruments or other energy emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.
-  **WARNING:** Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope applied part type CF is used. In this case an endotherapy device with a type CF applied part must be used in order to reduce the total patient leakage current.
-  **WARNING:** Live parts or instruments may not come into contact with TIPCAM®1 S 3D.
-  **WARNING:** Unintentional grounding of the patient: The metal part of the camera connector is electrically connected to the sheath of the endoscope. During use, make certain that these metal parts do not come into contact with ground potential.
-  **WARNING:** Risk of electric shock: Do not attempt to open the TIPCAM®1 S 3D. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

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WARNING: For use only by clinicians trained in laparoscopy!

CAUTION: When in contact with or in close proximity of an active electrosurgical instrument temporary image interference or loss can occur. An upgrade of the TC.302 hardware is recommended if the system reports this.

CAUTION: Ultrasonography before laparoscopy may identify clinical conditions that will alter patient management.

CAUTION: Abdominal puncture sites 7 mm or greater (for the introduction of auxiliary instrumentation) may be a source of herniation.

CAUTION: Lasers and electrosurgical probes should not be activated simultaneously. Further, when one is activated, the tip of the other should be completely retracted. This avoids deflection of the energy to the other tip.

CAUTION: During monopolar electrosurgery, inadvertent burns can occur when the appropriate patient return path is obstructed or as a result of capacitive coupling.

CAUTION: Avoid combinations of devices made of both conduction and insulating materials, especially metal trocars with nonconductive trocar anchors, screws, etc.

CAUTION: Periodically check the video system.

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Moisture which has penetrated the system due to damage can lead to faulty insulation paths within the videoendoscope. Rough surfaces or sharp edges can lead to patient injury.

ⓘ NOTE: Follow the applicable national ordinances and recycling plans regarding disposal or recycling of device components.

Acquaint yourself with the operation and control of the instrument before initial use of the TIPCAM®1 S 3D.

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6.2 Indications for use

6.2.1 Intended use

The Rigid Videoendoscope used together with the camera control unit is for use during diagnostic and/or surgical procedures when endoscopic video assistance is required.

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6.3 Contraindications

This instrument must not be used for observation or treatment of the heart and must not contact the heart or any area near the heart. In addition, this instrument must not come into contact with any device or therapeutic accessory that contacts the heart or any area near the heart.

Laparoscopy may be contraindicated by pregnancy greater than 16 weeks gestation or abdominal mass of comparable size.

Laparoscopy may also be contraindicated by the following conditions, depending on their severity or extent:

- Class IV cardiac decompensation
- Bowel obstruction
- Ileus
- Cardiac disease
- Intraperitoneal hemorrhage
- Diaphragmatic hernia
- Infection with acute peritonitis
- Previous abdominal surgery
- Obesity
- Thin nulliparous patient
- Chronic obstructive lung disease
- Liver failure with established collateral vessels

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6.4 User qualification

The TIPCAM®1 S 3D may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the TIPCAM®1 S 3D.

6.5 Safety precautions when operating the instrument

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During treatment in conjunction with the TIPCAM®1 S 3D, the patient must be treated and observed with the standard medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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7.5 Start-up of the TIPCAM®1 S 3D

! **NOTE:** The TIPCAM®1 S 3D including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

7.6 Connecting the light cable to the TIPCAM®1 S 3D

The light cable is mounted by fitting to the light port provided and turning clockwise. Check that it is firmly connected (see fig. A).

! **CAUTION:** Please note that the circular symbol on the TIPCAM®1 S 3D and the circular symbol on the light cable must correspond (see following table).

Overview: Recommended combinations*

Diameter light cable	
4.8 – 5.0 mm	○
3.0 – 3.5 mm	○
2.0 – 2.5 mm	●
Diameter endoscope	
6.6 – 12.0 mm	○
3.0 – 6.5 mm	○
0.8 – 2.9 mm	●

*Special endoscopes may deviate slightly

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

7.9 Switching on the CCU, light source and monitor

NOTE: Read the instruction manuals for the IMAGE1 S, the light source, the light cable and the monitor.

7.10 Optimum light efficiency

NOTE: The entrance and exit faces of the light cable and TIPCAM®1 S 3D must be kept clean if optimal light transmission is to be maintained. The light cable must first be unscrewed (see Fig. C).
Following this, clean the light entrance and exit faces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70 % isopropyl alcohol. Please note that both the glass surface inside the light port of the TIPCAM®1 S 3D (see arrow, Fig. D.) and the light cable itself (see Fig. E) must be cleaned in this way. Please note the instructions on cleaning the light cable port on page 24.

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**8 Reprocessing
not validated for
US customers**

WARNING: The light cable must be removed for reprocessing the TIPCAM®1 S 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM®1 S 3D and must not be removed for reprocessing.

NOTE: For US instructions, please refer to chapter 9.

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Store the TIPCAM®1 S 3D in the storage/sterilization box available as an optional extra. Do not coil the TIPCAM®1 S 3D video connecting cable into tight loops or cause it to kink while storing. Do not store the TIPCAM®1 S 3D in direct sunlight or excessive heat.

When returning the TIPCAM®1 S 3D for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

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The connectors of the video connecting cable and of the fiberoptic light cable must also be completely dry before connecting!

NEVER PICK UP THE TIPCAM®1 S 3D BY ITS CABLE! Always pick up the set by the camera head first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE TIPCAM®1 S 3D FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to device. When transporting, use the original shipping box to prevent damage to the device.

DO NOT STORE THE TIPCAM®1 S 3D WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Always arrange the cable in loose coils of 6 inches (15 cm) or more in diameter.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful "kinks" can result. Always allow the cable to uncoil gently and loosely.

NOTE: We strongly recommend drying using sterile compressed air following manual reprocessing and, if required, after machine reprocessing.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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8.2 General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

 **WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

 **CAUTION:** Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

8.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

8.4 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible. The light cable must be removed for reprocessing the TIPCAM*1 S 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM*1 S 3D and must not be removed for reprocessing.

8.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

8.6 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning them under cold running water with the aid of a brush or sponge.

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8.7 Manual cleaning

NOTE: In order to protect the contacts of the video connecting cable connector, manual cleaning and disinfection is only possible if the protective cap is fitted (see Fig. F).

The medical device must be rinsed under cold running water for 2 minutes. Then completely immerse the medical device in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of the brush for 3 minutes. Finally, it must be rinsed with cold water for 2 minutes to ensure neutralization.

8.8 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations.

Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

8.9 Machine cleaning and disinfection

NOTE: In order to protect the contacts of the video connecting cable connector, machine cleaning and disinfection is only possible if the protective cap is fitted (see Fig. F).

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

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Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

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CAUTION: After the program has run, the TIPCAM®1 S 3D is at a high risk of corrosion due to the ambient conditions in the reprocessing machine. Once the program has ended, remove the TIPCAM®1 S 3D from the washer and disinfectant as soon as possible.

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CAUTION: After machine reprocessing, the connection for the fiberoptic light cable must be removed in order to clean the optical fiber area with a cotton tip applicator moistened with 70 % isopropyl alcohol. Then the connector needs to be screwed on again. The connector can be removed by turning it to the left counterclockwise. The connector can be fastened by turning it to the right clockwise (see fig. G – I).

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8.10 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

8.11 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

8.12 Sterilization

NOTE: The protective cap of the video connecting cable connector must be disassembled for sterilization (see Fig. J).

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and following consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

NOTE: The 3D TIPCAM® 1 is located in its storage/sterilization container.

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Reprocessing (worldwide except US)

For sterilization of the medical device in its fully assembled state, a fractionated prevacuum procedure (DIN EN ISO 17665-1) at 134 °C – 137 °C with a minimum exposure time of 4 minutes and a maximum of 18 minutes must be used. This procedure is only suitable for thermostable instruments.

! **CAUTION:** Do not immerse the TIPCAM® 1 S 3D in cold liquid after removing it from the steam sterilizer. Due to the 'temperature shock', this may result in damage to the camera lens and invalidation of the warranty.

i **NOTE:** Greased components must be dismantled for sterilization so that steam penetration can be ensured.

Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

! **WARNING:** Medical devices must be dismantled for sterilization.

! **WARNING:** Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

! **CAUTION:** For more detailed information, please consult the user handbook of the respective device.

i **NOTE:** Please refer to the 'STERRAD® Sterility Guide' to verify whether the respective medical device can be sterilized using the various STERRAD® devices.

i **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S without booster
- STERRAD® NX® Standard Cycle
- STERRAD® 100NX® Standard Cycle

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Reprocessing (worldwide except US)

8.12.1 Accessories

Accessories recommended for reprocessing, transport and storage:

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CAUTION: Inspect the video connecting cable and fiberoptic light cable for breaks or cracks before and after disinfection or sterilization, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cable and cause damage to it. Where a TIPCAM®1 S 3D has a damaged fiberoptic light cable or damaged video connecting cable, the relevant cable must be replaced.

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Reprocessing for US customers only

9 Reprocessing for US customers only

9.1 General care

WARNING: Always pull out power plug before cleaning.

Store the TIPCAM®1 S 3D in the sterilization tray. The tray is an optional accessory that may be purchased from KARL STORZ. Do not coil the TIPCAM®1 S 3D video connecting cable into tight loops or cause it to kink while storing.

Do not store the TIPCAM®1 S 3D in direct sunlight or excessive heat.

Prior to returning for service, the TIPCAM®1 S must be thoroughly cleaned and sterilized per validated procedures. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

CAUTION: Chemical residue on the cardedge connector can cause camera failure.

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9.1.1 Care and handling

To maximize the life of your TIPCAM®1 S 3D please follow the guidelines given below for proper handling and care of the camera head.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining mois-

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9.2 General warnings

 **WARNING:** Risk of infection: Improperly reprocessed medical devices may pose a risk of infection for patients, users, and third parties and may cause malfunction of the medical device.

 **WARNING:** Risk of infection: These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users, and third parties. Check medical devices for visible soiling. Visible soiling indicates missed or improper cleaning. Reprocess medical devices using validated procedures before initial use and after each subsequent use.

 **WARNING:** During all work on contaminated medical devices, observe the occupational safety guidelines.

 **CAUTION:** When preparing the cleaning solution, follow the manufacturer's instructions for dilution to proper concentration and temperature.

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CAUTION: Country-specific laws and regulations must be observed.

9.3 Disassembly

Before reprocessing, disconnect the light cable from the endoscope and disconnect the endoscope connecting cable from the camera control unit.

9.4 Pre-cleaning

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the scope with a soft, lint-free, disposable cloth moistened with distilled/demineralized or enzymatic cleaning solution to remove all heavy soiling.
2. Contain and transport the scope to the decontamination area for cleaning.
3. If immediate cleaning is not feasible, place the scope in a container and soak for no longer than 1 hour in a neutral/mild pH enzymatic cleaning solution (e.g. Enzol), diluted with water to proper concentration per manufacturer's instructions.

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9.5 Manual cleaning

1. Thoroughly rinse the scope with distilled/ demineralized water at 20 °C, using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
2. Completely immerse the scope in a mild/ neutral pH enzymatic cleaning solution (e.g. Enzol) at 21 °C. Keep immersed for a minimum of 5 minutes.
3. While immersed, using a soft bristle brush (P/N: 27652/3), brush all surfaces of the scope, including the silicone cable, card edge connector and operating buttons.

i NOTE: KARL STORZ cleaning brushes are provided non-sterile, single-use and must be discarded after each use. Additional brushes may be purchased from KARL STORZ.

4. After immersion, rinse off the detergent with distilled/demineralized water at 20 °C, using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
5. Dry the scope with a lint-free soft cloth or clean, dry, compressed air (<5 psi).

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9.6 Inspection after cleaning

The cleaned medical device must be visually inspected to verify that it is clean, complete, damage-free, and dry.

- If soiling or residues remain, the medical device must be cleaned manually and undergo another complete cleaning cycle.
- The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.

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9.7 Sterilization

WARNING: Fully clean the endoscope prior to sterilization.

KARL STORZ recommends to sterilize the 3D TIPCAM® via steam sterilization.

1. Place the endoscope into a 3D TIPCAM® SteriTite container. Follow the container's manufacturer's instructions for details on how to properly set up the container.
2. Load the container into the sterilizer.
3. Use the following cycle parameters for steam sterilization:
 - Temperature: 132 °C (270 °F)
 - Exposure Time: 4 minutes
 - Drying Time: 20 minutes
4. Following cycle completion, crack open the sterilizer door with the container still within the sterilization chamber. Allow the load to cool for 15 minutes with the door cracked open to eliminate any condensation.
5. Remove the container from the sterilization chamber and allow to cool on a stainless steel wire rack for an additional 25 minutes.

9.8 Limits of reprocessing

The end of the product's lifespan is decisively determined by wear, reprocessing methods, employed chemicals, and any use-related damage.

Always perform a pre-inspection functional test with the full 3D system prior to using the device. The system should not be used if any of the accessories are damaged or non-functional.

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10 Maintenance

10.1 Service and maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

- ⓘ NOTE:** To ensure safe operation of the camera control unit (CCU), you should perform the following procedures every 12 months:
- Check leakage current from chassis < 100 µA.
 - Check ground impedance < 0.1 ohms.
 - Check power consumption rated power.

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10.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

10.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

10.4 Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to:

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

10.5 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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10.6 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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11.1.1 Directive compliance (for TIPCAM®1 S 3D)

According to Medical Device Directive (MDD):
This medical product belongs to Class IIa
This medical product bears the CE mark in accordance with MDD 93/42/EEC.

 **NOTE:** The code number after the CE mark indicates the responsible notified body.

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12 Recommended accessories

39501 XTC*

Tray for cleaning, sterilization and storage of the TIPCAM®1 S 3D and a fiberoptic light cable.

External dimensions (w x d x h):
640 x 150 x 87 mm

495 TIP

Fiberoptic light cable, diameter 4,8 mm, length 300 cm, resistant to very high temperatures, for use with TIPCAM®1 S 3D

495 NCSC*

Fiber optic light cable, extremely heat-resistant, safety lock, diameter 4,8 mm, length 250 cm, for use with TIPCAM®1 S 3D

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
KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33129-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseala.com

KARL STORZ Endoscopia
Miramar Trade Center,
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P., 11520 Mexico D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itam Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 158
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14,
12724 Skärholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Valimotie 1 B, 4. kerros
00380 Helsinki, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ GmbH & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TU, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 430
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552549, 552788
Fax: +7 7172 552444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf
Building 25 (Al-Hayyan) – Units 205 – 206
Dubai Health Care City, P.O.Box 118069
Dubai – United Arab Emirates
Phone: +971 (0)4 4471230, Fax: +971 (0)4 4471231
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 43743000-30,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu, F10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +84 981 8442, Fax: +84 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610044, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongail-Ro, Bundang-Gu
Seongnam, 463-859, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.jp

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: mfo@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kscla.com

KARL STORZ Endoscopy
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col Chapulterec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda.
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, São Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsatragrand 14
12724 Skarholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4 Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopie Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopie France S.A.S
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) - Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai - United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd
D-181, Okhla Industrial Area, Phase-I
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 - Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Important Information
for users of KARL STORZ
devices**

**1 Important information for
users of KARL STORZ
devices**

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to aid you with proper connection and operation of the TIPCAM® 1 SPIES™ 3D in conjunction with the instruction manual for the IMAGE 1 SPIES™. All required details and all actions on your part are clearly presented and explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment; insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the IMAGE 1 SPIES™ instruction manual for the individual functions of the CCU that can be used with the TIPCAM® 1 SPIES™ 3D camera.

! **CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

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
**3 Controls, displays,
connectors, and their uses**

- ① Inserter (sheath) with sensors
- ② Handle
- ③ Control button 1: Scroll up in menu/ programmed functions
- ④ Control button 2: Selection/programmed functions
- ⑤ Control button 3: Scroll down in menu/ programmed functions
- ⑥ Connection for fiberoptic light cable
- ⑦ Video connecting cable connector to CCU
- ⑧ Fiberoptic light cable (495 TIP)
- ⑨ Light cable connector to KARL STORZ light source
- ⑩ Video connecting cable

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4 General information

4.1 Device description

The TIPCAM® 1 SPIES™ 3D is a high-resolution video laparoscope with two distally integrated video chips.

Since the TIPCAM® 1 SPIES™ 3D does not require an optical imaging system, it is very sturdy.

The TIPCAM® 1 SPIES™ 3D is compatible over D3-LINK with the IMAGE 1 SPIES™.

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WARNING: Also pay particular attention to the instruction manual for the IMAGE 1 SPIES™ CCU no. 96206286DER.

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product! Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

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- ! WARNING:** Risk of injury to the patient/user. Damaged or incomplete instruments/accessories may cause malfunctions. Before and after every use, always check the entire system (e.g., the TIPCAM® 1 SPIEST™ 3D, all cables, the CCU, monitor and accessories used) to ensure that everything is configured properly, in perfect working order and nothing is missing/damaged. Replace incomplete/damaged instruments/accessories and do not continue to use them.
- Before every use, check that all the instruments required for the operation are functioning correctly.
 - Be certain to check the cover glass of the TIPCAM® 1 SPIEST™ 3D. In addition, the camera cable and the light cable need to be inspected for damage cracks.
 - Do not use wet/moist cables.
 - When attaching components to an instrument, the secure position of the components must be checked by gently pulling and/or turning them in both directions before use.

- ! WARNING:** Check the instrument and accessories used in combination before every application for unintentional rough surfaces, sharp corners, burred edges, and projecting parts, which can endanger personal safety. Damaged, defective, bent, or incomplete instruments or accessories must not be used further under any circumstances.

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- ⚠️ WARNING:** Risk of injury and risk of damaging products: The high level of light energy means that the distal end, the light port, adjacent components and tissue in front of the light outlet window to heat up. This can cause burns to the patient, user and surgical accessories.
- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
 - Avoid excessively long (> 10 min) continuous direct tissue contact with the shaft location.
 - Never allow the endoscope to rest on the patient or come into direct contact with surgical accessories.
 - Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.
- ⚠️ WARNING:** Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected fiberoptic light cable or an endoscope. Always wear standardized protective equipment when using lasers.
- ⚠️ WARNING:** The TIPCAM® 1 SPIES™ 3D must be replaced if the image is cloudy or if no image or only parts thereof are visible.
- ⚠️ WARNING:** Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the TIPCAM® 1 SPIES™ 3D is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the handle, which can cause cardiac pacemakers or other electronic devices to malfunction.
- ⚠️ WARNING:** The video system can malfunction during use. An adequate replacement system must be kept available.

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- ⚠️ WARNING:** In the event of image failure immediately remove TIPCAMP® 1 SPIESTM 3D from the patient.
- ⚠️ WARNING:** The use of HF probes, HF electrodes, HF instruments or other energy emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.
- ⚠️ WARNING:** Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope applied part type CF is used. In this case an endotherapy device with a type CF applied part must be used in order to reduce the total patient leakage current.
- ⚠️ WARNING:** Live parts or instruments may not come into contact with TIPCAMP® 1 SPIESTM 3D.
- ⚠️ WARNING:** Unintentional grounding of the patient: The metal part of the camera connector is electrically connected to the sheath of the endoscope. During use, make certain that these metal parts do not come into contact with ground potential.
- ⚠️ WARNING:** Risk of electric shock: Do not attempt to open the TIPCAMP® 1 SPIESTM 3D. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.
- ⚠️ CAUTION:** Periodically check the video system

Moisture which has penetrated the system due to damage can lead to faulty insulation paths within the video laparoscope. Rough surfaces or sharp edges can lead to patient injury.

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NOTE: Follow the applicable national ordinances and recycling plans regarding disposal or recycling of device components. Acquaint yourself with the operation and control of the instrument before initial use of the TIPCAM® 1 SPIES™ 3D.

5.2 Intended use

The Rigid Videendoscope used together with the camera control unit is for use during diagnostic and/or surgical procedures when endoscopic video assistance is required.

5.3 Indication

For use in all endoscopy and endoscopic surgery within the peritoneal and thoracic cavity.

5.4 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

5.5 User qualification

The TIPCAM® 1 SPIES™ 3D may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the TIPCAM® 1 SPIES™ 3D.

5.6 Safety precautions when operating the instrument

It is the user's responsibility to make sure the instrument is safe and operates properly before using the instrument.

During treatment in conjunction with the TIPCAM® 1 SPIES™ 3D, the patient must be treated and observed with the standard medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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Installation and operating instructions

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6 Installation and operating instructions

6.1 Unpacking the equipment

Carefully unpack the TIPCAM® 1 SPIES™ 3D and remove the unit from its packaging. Check for missing items and evidence of shipping damage.

File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

6.2 Basic equipment

- 1 KARL STORZ TIPCAM® 1 SPIES™ 3D
- 1 Instruction manual 96206364DER

6.3 Compatible camera control unit (CCU)

CCU

IMAGE 1 SPIES™
IMAGE 1 CONNECT (TC 200) in combination with IMAGE 1 D3-LINK (TC 302)

6.4 Before every use

Before every use, the instrument must be checked for damage, especially around the distal end (cover glass, light outlet, ...).

Before each use, the video connecting cable and light cable must be checked for damage, especially breaks and cuts.

Instruments with damaged cables must not be used further under any circumstances. Cables must be replaced (light cable to be replaced by the user or the Service Dept., and the video connecting cable by the Service Dept.).

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6.5 Start-up of the TIPCAM® 1 SPIES™ 3D

① **NOTE:** The TIPCAM® 1 SPIES™ 3D including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

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6.6 Connecting the light cable to the TIPCAM® 1 SPIES™ 3D

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⚠ **CAUTION:** Please note that the circular symbol on the TIPCAM® 1 SPIES™ 3D and the circular symbol on the light cable must correspond (see following table).

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Diameter light cable	
4.8 – 5.0 mm	○
3.0 – 3.5 mm	●
2.0 – 2.5 mm	●
Diameter endoscope	
6.6 – 12.0 mm	○
3.0 – 6.5 mm	●
0.8 – 2.9 mm	●

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* Special endoscopes may deviate slightly

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6.9 Switching on the CCU, light source and monitor

- ① **NOTE:** Read the instruction manuals for the IMAGE 1 SPIES™, the light source, the light cable and the monitor.

6.10 Optimum light efficiency

- ① **NOTE:** The entrance and exit faces of the light cable and TIPCAMP® 1 SPIES™ 3D must be kept clean if optimal light transmission is to be maintained. The light cable must first be unscrewed (see Fig. C). Following this, clean the light entrance and exit faces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70 % isopropyl alcohol. Please note that both the glass surface inside the light port of the TIPCAMP® 1 SPIES™ 3D (see arrow, Fig. D.) and the light cable itself (see Fig. E) must be cleaned in this way. Please note the instructions on cleaning the light cable port on page 33.

- ⚠ **CAUTION:** Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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**Reprocessing
not validated for
US customers**

7.1 General care
! WARNING: Always pull out power plug before cleaning.

Store the TIPCAM® 1 SPIES™ 3D in the storage/sterilization box available as an optional extra. Do not coil the TIPCAM® 1 SPIES™ 3D video connecting cable into tight loops or cause it to kink while storing. Do not store the TIPCAM® 1 SPIES™ 3D in direct sunlight or excessive heat.

When returning the TIPCAM® 1 SPIES™ 3D for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

! CAUTION: Chemical residue on the cardedge connector can cause camera failure.

! CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

7.1.1 Care and handling

To maximize the life of your TIPCAM® 1 SPIES™ 3D please follow the guidelines given below for proper handling and care of the camera head.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. Where necessary, remove remaining moisture with a low lint, disposable cloth.

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① **NOTE:** We strongly recommend drying using sterile compressed air following manual reprocessing and, if required, after machine reprocessing.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- ✗ On the grounds of potential damage, this reprocessing step must not be performed.

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7.2 General warnings

! WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the „Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments“ instructions and the accompanying documentation.

! WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

! WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

! CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

7.3 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible. The light cable must be removed for reprocessing the TIPCAM® 1 SPIES™ 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAM® 1 SPIES™ 3D and must not be removed for reprocessing.

7.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device by wiping down and rinsing, for example.

As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

**7.5 Manual pre-cleaning
Brushing the surfaces**

Any heavy soiling of the medical device surfaces can be completely removed by cleaning them under cold running water with the aid of a brush (Art. no. 27652) or sponge.

Auxiliary pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

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7.6 Manual cleaning

① **NOTE:** In order to protect the contacts of the video connecting cable connector, manual cleaning and disinfection is only possible if the protective cap is fitted.

The medical device must be completely immersed in a cleaning solution. Through disassembling the instrument, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with the aid of a brush (Art. no. 27652) or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

7.7 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. Through disassembling the instrument, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. Following the exposure time, the medical device must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces are dried completely with sterile compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

7.8 Machine cleaning and disinfection

① **NOTE:** In order to protect the contacts of the video connecting cable connector, machine cleaning and disinfection is only possible if the protective cap is fitted.

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments".

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Reprocessing not validated for US customers

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_{v} value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

CAUTION: After the program has run, the TIPCAM® 1 SPIES™ 3D is at a high risk of corrosion due to the ambient conditions in the reprocessing machine. Once the program has ended, remove the TIPCAM® 1 SPIES™ 3D from the washer and disinfectant as soon as possible.

CAUTION: After machine reprocessing, clean the distal lens and connector of the video connecting cable and fiberoptic light cable with a lens cleaner or a cotton tip applicator moistened in 70 % isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement. Please pay attention to chapter "Optimum light efficiency".

CAUTION: After machine reprocessing, the connection for the fiberoptic light cable must be removed in order to clean the optical fiber area with a cotton tip applicator moistened with 70 % isopropyl alcohol.

Then the connector needs to be screwed on again. The connector can be removed by turning it to the left counterclockwise. The connector can be fastened by turning it to the right clockwise (see fig. F–H).

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7.9 Assembly, inspection and care
The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness.

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Threads and glide surfaces must be lubricated locally with instrument oil.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

① **NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

① **NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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7.11 Sterilization

- ① **NOTE:** The protective cap of the video connecting cable connector must be disassembled for sterilization.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual „Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments“. The method must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device.

Steam sterilization in the fractionated prevacuum procedure

- ① **NOTE:** The 3D TIPCAM® 1 is located in its storage/sterilization container.

The medical device must be sterilized in its fully assembled state using a fractionated prevacuum procedure (ISO 17665-1) at 134-137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

- ① **NOTE:** Greased components must be dismantled for sterilization so that steam penetration can be ensured.

- ⚠ **CAUTION:** Do not immerse the TIPCAM® 1 SPIES™ 3D in cold liquid after removing it from the steam sterilizer. Due to the „temperature shock“, this may result in damage to the camera lens and invalidation of the warranty.

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The following procedures are approved in respect of material compatibility, but must be validated by the on-site operator in respect of efficacy.

Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

⚠️ WARNING: Medical devices must be dismantled for sterilization.

⚠️ WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

⚠️ CAUTION: For more detailed information, please consult the user handbook of the respective device.

ⓘ NOTE: Refer to the „STERRAD® Sterility Guide“ to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

ⓘ NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization procedures are approved in respect of material compatibility, but must be validated by the on-site operator in respect of efficacy.

- STERRAD® 50, 100S, 200
- STERRAD® NX®
- STERRAD® 100NX®
- STERRAD® 100NX® DUO Cycle

Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™ 1

Detailed information on the selection of the appropriate cycle for the various device generations is available from the manufacturer STERIS®.

ⓘ NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

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7.12 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

CAUTION: Inspect the video connecting cable and fiberoptic light cable for breaks or cracks before and after disinfection or sterilization, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cable and cause damage to it. Where a TIPCAM® 1 SPIES™ 3D has a damaged fiberoptic light cable or damaged video connecting cable, the relevant cable must be replaced.

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**8 Reprocessing
for US customers only**

**8.1 Cleaning, disinfection
and sterilization
for US customers only**

WARNING: The light cable must be removed for reprocessing the TIPCAMP® 1 SPIESTM 3D. To do this, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the TIPCAMP® 1 SPIESTM 3D and must not be removed for reprocessing.

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8.2 General care

WARNING: Always pull out power plug before cleaning

Store the TIPCAM® 1 SPIES™ 3D in the storage/sterilization box available as an optional extra. Do not coil the TIPCAM® 1 SPIES™ 3D video connecting cable into light loops or cause it to kink while storing.

Do not store the TIPCAM® 1 SPIES™ 3D in direct sunlight or excessive heat.

Prior to returning for service, the TIPCAM® 1 SPIES™ must be thoroughly cleaned and sterilized/high-level disinfected per validated procedures. The light entrances and exit faces of the light cable and the endoscope must always be kept clean (please refer to the corresponding instruction manual) to ensure optimal reproduction of the endoscopic image.

CAUTION: Chemical residue on the cardedge connector can cause camera failure

CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the endoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET: Always ensure the connector is completely dry and free of debris. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. Where necessary, remove remaining moisture with a low lint, disposable cloth.

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NEVER PICK UP THE TIPCAM® 1 SPIES™ 3D BY ITS CABLE! Always pick up the set by the camera head first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE TIPCAM® 1 SPIES™ 3D FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

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DO NOT STORE THE TIPCAM® 1 SPIES™ 3D WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Always arrange the cable in loose coils of 6 inches (15 cm) or more in diameter.

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- This reprocessing step is part of the validated procedure.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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8.3 General warnings

! WARNING: Risk of infection: Improperly reprocessed medical devices may pose a risk of infection for patients, users, and third parties and may cause malfunction of the medical device.

! WARNING: Risk of infection: These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users, and third parties. Check medical devices for visible soiling. Visible soiling indicates missed or improper cleaning. Reprocess medical devices using validated procedures before initial use and after each subsequent use.

! WARNING: During all work on contaminated medical devices, observe the occupational safety guidelines.

! CAUTION: When preparing and using solutions, strictly observe the chemical manufacturer's instructions on concentration, exposure time, and shelf life. Excessively long immersion or incorrect concentration may cause damage. Note the employed chemicals' spectrum of antimicrobial activity.

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CAUTION: Country-specific laws and regulations must be observed

8.4 Disassembly

Before cleaning and disinfection, the medical device must be opened or disassembled into its individual component parts to the fullest extent possible.

8.5 Pre-cleaning

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the scope with a soft, lint-free, disposable cloth moistened with water or enzymatic cleaning solution to remove all heavy soiling.
2. Contain and transport the scope to the decontamination area for cleaning.
3. If immediate cleaning is not feasible, place the scope in a container and soak for no longer than 1 hour in a neutral/mild pH enzymatic cleaning solution (e.g. EnzoI), diluted with water to proper concentration per manufacturer's instructions

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8.6 Manual cleaning

1. Thoroughly rinse the scope with distilled/demineralized water at 17-19 °C, using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
2. Completely immerse the scope in a mild/neutral pH enzymatic cleaning solution (e.g. Enzol) at 17.5 – 20.5 °C. Keep immersed for a minimum of 3 minutes.
3. While immersed, using a soft bristle brush (P/N: 27652/3), brush all surfaces of the scope, including the silicone cable, card edge connector and operating buttons.
4. After immersion, rinse off the detergent with distilled/demineralized water at 17-19 °C, using a cleaning pistol. Rinse the scope head for 1.5 minutes and the silicone cable/card edge connector for 30 seconds, for a total rinse of 2 minutes.
5. Dry the scope with a lint-free soft cloth or clean, dry, compressed air (<5 psi).

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**Reprocessing
for US customers only**

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- If soiling or residues remain, the medical device must be cleaned manually and undergo another complete cleaning cycle.
- Damaged or corroded medical devices must be taken out of circulation.
- This must be followed by functional testing.

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8.8 Sterilization

The instructions "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" describe in detail the workflows and process-relevant parameters of the individual validated procedures. The procedure must be chosen taking into account applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved for this medical device by KARL STORZ:

**Steam sterilization using the fractionated
pre-vacuum procedure**

The medical device must be sterilized in assembled condition using the fractionated pre-vacuum procedure (DIN EN ISO 17665-1) at 134 °C to 137 °C, for a minimum exposure time of 4 minutes and a maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

① **NOTE:** Lubricated components must be sterilized in disassembled condition to ensure steam penetration.

The following procedures have been approved in terms of material compatibility, but their effectiveness must be validated by the on-site operator.

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Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

! WARNING: Medical devices must be sterilized in disassembled condition.

! WARNING: Please note that materials and lumen dimensions determine which products can be sterilized using the various STERRAD® sterilization systems.

! CAUTION: For detailed information, please refer to the user manual of the respective device.

① NOTE: Please refer to the STERRAD® Sterility Guide to verify whether the respective medical device can be sterilized using the various STERRAD® devices.

① NOTE: Lubricated and oiled surfaces cannot be sterilized.

The following STERRAD® sterilization procedures have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 50, 100S, 200
- STERRAD® NX®
- STERRAD® 100NX®
- STERRAD® 100NX® DUO Cycle

Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™ I

Detailed information on choosing the appropriate cycles in the various device generations is available from the manufacturer STERIS®.

① NOTE: Lubricated and oiled surfaces cannot be sterilized.

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9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ, all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's sealability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

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① **NOTE:** The code number after the CE mark indicates the responsible notified body

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11 Recommended accessories

39501 XTC*

Tray for cleaning, sterilization and storage of the TIPCAM® 1 SPIES™ 3D and a fiberoptic light cable.

External dimensions (w x d x h):
640 x 150 x 87 mm

495 TIP

Fiberoptic light cable, diameter 4,8 mm, length 300 cm, resistant to very high temperatures, for use with TIPCAM® 1 SPIES™ 3D

495 NCSC*

Fiber optic light cable, extremely heat-resistant, safety lock, diameter 4,8 mm, length 250 cm, for use with TIPCAM® 1 SPIES™ 3D

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
KARL STORZ ENDOSKOPE

East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Barid Souks, Weyland Street
2012 3301 Barid, Lebanon
Phone: +961 1 9063390; Fax: +961 1 9063391
E-Mail: info@kharizor-erng.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 8001, Roggebaai 8012, Cape Town, South Africa
Phone: +27 21 417 2900; Fax: +27 21 421 5103
E-Mail: info@kharizor.co.za

KARL STORZ ENDOSCOPE Kazakhstan LLP
c/o Sanyerka Ltd. BC "Arman", office-420
010000 Astana, Republic of Kazakhstan
Phone +7 7172 552549, 552788
Fax +7 7172 552444
E-mail info@karstorz.kz

KARL STORZ ENDOSCOPE East Mediterranean and Gulf
Building 25 (A) "Hermes" - Unit 205 - 208
110000 Beirut, Lebanon
Phone +961 3 333 118606
Fax +961 3 333 118606
Doha - United Arab Emirates
Phone +971 04 4471230, Fax +971 04 4471231
E-Mail info.gulf@karstorz-ameg.com

KARL STORZ ENDOSCOPE India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone +91 11 427 43000-30;
+91 11 265122 86
E-Mail info@karstorz-india.com

KARL STORZ GmbH & Co. KG
89073 447169 Dengen Veng. F10-0, Ph. Nhusan
Chi Mai City, Vietnam

KARL STORZ Endoscopy Australia Pty Ltd
Phone +61 81 901 8442; Fax +61 81 844 0320
E-Mail: info@karl-storz.com.au

KARL STORZ Endoscopy Australia Pty Ltd
15 Oron Road, Lane Cove NSW 2060
P O Box 50 Lane Cove NSW 1595, Australia
Phone +61 (0)2 9480 8700
Fax +61 (0)2 9480 8702
Toll free 1 800 965 5402 (Australia only)
Phone +61 (0)2 9420 0905
E-Mail: info@karl-storz.com

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503/68, 23F AIA Tower, Landmark East,
No. 100 Hing Shing Street, Tuen Tin, Kowloon,
Hong Kong
Phone +852 2381 0000
Fax +852 2381 0001
E-Mail: karl@karl-storz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jangongnawar Street, Chaoyang District,
Beijing, 100022, Capital of the People's Republic of China
Phone: +86 10 56338198, Fax: +86 10 56381900
E-Mail: info@karlstorz.com on

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongcun Road, Xuhui District

Shanghai, 200030, People's Republic of China
Phone: +86 21 611131188. Fax: +86 21 611131190
E-Mail: info@kltor.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Jiantang Road, Wuxiao District,
Chengdu, 610014, People's Republic of China
Phone: +86 28 68587977. Fax: +86 28 68587975
E-Mail: info@kltor.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Sanyang Branch
Room 2001-2005, 20F N-MEDIA International Centre,
No. 167 North Sichuan Road, Shanghai District,
Shanghai 200031, People's Republic of China
E-Mail: info@karlstorz.com.cn
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 2001-2005, 20F N-MEDIA International Centre,
No. 167 North Sichuan Road, Shanghai District,
Shanghai 200031, People's Republic of China
E-Mail: info@karlstorz.com.cn
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 3 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69220130 Fax: +65 69622015
E-Mail: info@karlstorz.com.sg
E-Mail: info@karlstorz.com.sg

KARL STORZ Endoscopy Sales Pte Ltd
No. 3 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69220130 Fax: +65 69622015
E-Mail: info@karlstorz.com.sg
E-Mail: info@karlstorz.com.sg

Singapore 148855, Singapore
Phone +85 90229150, Fax +85 90229155
E-Mail: info@karl-storz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
1-Park Bundang 1, 102Bong, Bf
239 Jeongnail-ro, Bundang-gu
Seongnam, 461-700, Korea
Phone +82 31-827-3200
Fax +82 31-827-3290
E-Mail: asptel-70-kar@karl-storz.com

KARL STORZ Endoscopy Taiwan Ltd.
125 No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone +886 933 014 160

KARL STORZ Endoscopy Japan K. K.
 Stage Bldg. 8F, 2-7-2 Fujimi
 Chiyoda-ku, Tokyo 102-0077, Japan
 Phone: +81 3 380-8022 Fax: +81 3 6340-8633
 E-Mail: info@karlstorz.co.jp
 www.karlstorz.com

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) – Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai – United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-I
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 – Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beisin Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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Karl Storz GmbH & Co. KG