INSTRUCTIONS FOR USE OF FRAMES FOR CORRECTIVE EYEGLASSES with transport demo lenses: metal, plastic and combined

1. PURPOSE:

Frames for corrective eyeglasses with transport demo lenses: metal, plastic and combined, produced by **OPTA FILAO SARL** have been designed for fixing the lens in a proper position relative to the eye.

2. SPECIFICATIONS AND DESIGN

Frames for corrective eyeglasses produced by **OPTA FILAO SARL** consist of a rim, temples, nose pads, and hinges. Hinges can be simple and tension.

Frames for corrective eyeglasses are characterized by three basic geometric parameters: the width (a), the width of the bridge (b) and the length of the temple (L/L1).

The dimensions of the **OPTA FILAO SARL** frames used are given in the table below.

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Limit deviations dimensions α and b must be: $\pm 0,2$ mm for frames such as OM; $\pm 0,3$ mm for frames type OP; ± 1 mm in size L.

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Fig.1 Scheme of frames

The angle γ and groove depth s should be from 80° to 120° and from 0.60 to 1.2 mm, accordingly (Fig. a) for type FP and FC frames.

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The angle γ_1 and the groove depth s_1 should be from 90° to 110° and from 0.30 to 1.0 mm, accordingly (Fig.2 b) for metal frames.

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Hard temples should be bent in accordance with Figure 5b. Numeric values should be as set out in Table. 4.

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Bend of temple is fitted on the anthropometric indicators of the user.

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The temple of the frames with a sprung-loaded hinge must have two fixed positions: horizontal and vertical. Tension springs rejecting temples at the same angle must be the same for both the temples.

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The temple of the frames without a sprung- loaded hinge montage at an angle less than 90 ° to a horizontal front should not drop under its own weight.

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Temples shall turn smoothly without sticking along the axis.

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Basic dimensions of the surface of the nose support are listed in Fig, with a range of numerical values of these dimensions should be as given in the table 5 .

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
Fig.6 "The dimensions of the support surface for the nose".

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On the surface of the frame, except faceted grooves should not be sharp edges and burrs. On the surface of the frames should not be visible to the naked eye shells, scratches, sagging, blistering and other defects which adversely affects their appearance. On the surface of the metal frame is allowed points presence diameter of not more than 0.3 mm in an amount not more than 5, on the inner surface - Hairline is not more than 10 mm in an amount not more than 2 pieces.

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Compounds frames parts must be durable and withstand the following forces:

- permanent connections - $(5 \pm 0,5) \text{ H } (0,5 \pm 0,05) \text{ kgf}$;
- tips temples with the metal part - $(2 \pm 0,2) \text{ H } (0,2 \pm 0,02) \text{ kgf}$.

The frame must not have damage and / or to be deformed after removal effort .

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The weight of frame without lenses should not exceed 38 g.

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Decorative coatings must not be applied.

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Frames for corrective glasses must meet the standards: ISO EN 12870.

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3. INDICATION FOR USE:

Correction of defects in the optics of the eyes:

- myopia (myopia) - short-sightedness,
- Presbyopia (presbyopia) - hyperopia
- astigmatism (astigmatismus)
- Hypermetropia (hypermetropia), a weak type of clinical refraction
- anisometropia (difference in refractions on the eyes)
- keratoconus.

CONTRAINDICATIONS TO THE USE can be the following:

- Idiosyncrasy of vision correction when wearing glasses;
- Impossibility of vision correction when wearing glasses.

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4. MATERIALS USED:

Plastic frames are made from material: celluloid acetate, MAZZUCHELLI code GB-023574.
Metal, combined and rimless frames are made from stainless steel H12/HV330-350 from Sonia B Design Arbent 01100; Titanium (titanium 90%, aluminium 6%, vanadium 4%), from Panavision Optical ; tips of the temples - from the celluloid acetate MAZZUCHELLI code GB-023574; nose pads - from the silicone **LSR50 shore** (toughness scale);
For production of hinges, screws, nuts and rivets, use: stainless steel, code: H12/HV330-350
Titanium Titanium (titanium 90%, aluminium 6%, vanadium 4%) **JIS TP270H grade 1**
Transport demo lenses are made from acryl Polycryl Clear (antiblick). The materials meet the requirements for biological safety

5. MARKING AND PACKING:

Frames for corrective eyeglasses made by OPTA FILAO SARL have their markings on the temples and includes:

- Code or model name
- Technical dimensions
- Color Code

Inside the left temple there is the name of the manufacturer and the mark of conformity.

Inside the right temple there are the technical parameters: horizontal width of the lens installed (a), the width of the bridge (b) and the total length of the temples (L); the code or model name; and color code.

The name of the model, colour and technical parameters (a, b and L) are labelled on the card in the packaging plastic bag

In the polyethylene packaging of each frame or on the label, deposited in polyethylene packaging, should be specified:

- The nominal value of the main dimensions;
- model number;
- Color of the model;
- barcode
- lot
- manufacturer
- release date
- expiry date

Due to the diversity of models and different sizes of temples are

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
not may be point out the manufacturer on the temples at the specified model code (in accordance with EU rules).

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Frames, packed in polyethylene packaging should be tightly packed in group packaging boxes.

In the Collective packing box should be stacked frames of the same type with the same basic dimensions, in the amount of 10 pieces.

For each Collective packing box must be glued label enterprise

- the manufacturer, which should be the inscription:
- model number;
- Color of the model;
- The nominal value of the main dimensions;
- The number of frames.

Dimensions of the shipping box in cm (length x width x height): 23 x 15 x 5.8.

6. PREPARING THE PRODUCT FOR USE:

After the insertion of lenses the frame must not lose the original shape. Rims must be symmetrical. Chips, scratches, gaps between the lenses and frame rims, and other defects in the assembly points are not allowed.

Lenses inserted in the frame must not shift and fall out of rims during normal operating.

7. STORAGE AND USE:

- It is not allowed to apply considerable mechanical efforts to the frames;
- Frames must be stored in a case or in any other packaging, which excludes the possibility of mechanical damage;
- The full average service life of the rim is at least 3 years.
- The storage conditions: in heated storage at ambient temperature from + 5 °C to + 50 °C, and a relative humidity of 80%.
- During storage do not expose to corrosive substances (acids, alkalis) and organic solvents or their vapors.

8. CONDITIONS OF USE:

Spectacle frames for corrective glasses are used in all climatic zones, at any time of the day.

9. PRECAUTIONS OF USE:

When using eyeglass frames one must be careful when doing sports, and carrying out works involving sudden movements of the head, slopes, jumps and bumps

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