CERTIFICATION DE SIGNATURE

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 10/04/2018

Luc DEMANGEON

Development, Regulatory & Quality Director

Instructions for use

Y-connectors for endovascular procedures in kits

Super Ketch™ - Speed Ketch™

1 Description

1.1 Y connectors



The **Super Ketch™** and **Speed Ketch™** are Y connectors comprising a tricuspid star shaped valve, designed to provide haemostatis during endovascular procedures. The **Super Ketch™** and **Speed Ketch™** include a straight and a side access lumen.

The haemostatic valve, located at the proximal end of the straight lumen, can be opened and closed by operation of the “push-pull” button valve. The distal end has a rotating male luer lock.

The side access lumen allows pressure monitoring and perfusion.

The side access lumen of the **Speed Ketch™** ends with a female luer lock,

The side access lumen of the **Super Ketch™** with a side arm tubing ending by a standard three-way stopcock.

1.2 Metallic Insertion Tool

The metallic insertion tool is an accessory designed to facilitate the introduction of guidewires (up to 0.018”) used in conjunction with interventional devices during vascular interventional procedures.


1.3 Plastic Insertion Tool

This accessory is designed to facilitate the introduction of (up to 3 French) catheters through the haemostatic valve.


1.4 Torquer

This device is designed to facilitate steering of 0.010” to 0.035” (0.89mm) guidewires.

**Supplied with KESU0002, KESU0005 and KESU0006.

1.5 Product specifications

Components Specifications		
Speed Ketch™ / Super Ketch™	Inner Lumen Compatibility	7F / 0.092" / 2.33 mm maximum
	Pressure Resistance ⁽¹⁾ after introduction of a device	8 atm. maximum
Metallic Insertion Tool	External Diameter	2.7 F / 0.90 mm
	Inner Lumen Diameter	1.9 F / 0.025" / 0.64 mm
	Usable Length	10 cm
Plastic Peelable Insertion Tool	Length	9 cm
	Inner Diameter	3.4 F / 0.044" / 1.14 mm
Torquer	Guidewire Compatibility	From 0.010" minimum to 0.035" (0.89mm) maximum

⁽¹⁾ Haemostatic valve resistance measured in vitro, by applying injection pressure through the side port.

2 Indications

The **Super Ketch™** and **Speed Ketch™** Y connectors are intended to maintain haemostatis before and during the introduction / the use and withdrawal of diagnostic / interventional devices with outside diameter $\leq$ 7F (2.33 mm, 0.092") during interventional procedures.

3 Warnings

- The **Super Ketch™** or **Speed Ketch™** devices should exclusively be used by physicians trained in interventional procedures. It is recommended to the physician to consult current peer-reviewed publications on the interventional techniques.
- After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.

4 Precautions

- This device is designed and intended for single use only. In case of re-use, device sterility shall be lost and performances altered. Do not resterilise or reuse it. Use prior to the "use before" date noted on the packaging.
- Inspect the **Super Ketch™** or the **Speed Ketch™** prior to use. Do not use opened or damaged packages.
- Ensure that the **Super Ketch™** or **Speed Ketch™** is flushed correctly and all air removed before use (refer to Chapter 7 Instructions for use). Any aspiration through the side access of the **Super Ketch™** or **Speed Ketch™** should be effectuated slowly and carefully to avoid the creation of a void in the proximal part of the connector, leading to an effect of suction on the haemostatic valve.

- Always open the haemostatic valve using the push-pull button valve opener before inserting any device with a fragile distal portion, in order to avoid damage of the product when passing through the valve.
- Always retrieve the metallic insertion tool before retrieval of the guidewire, to avoid any damage of the guidewire.
- The metallic insertion tool should be inserted carefully through the valve in order not to damage or dislodge it.
- Do not attempt to force a device with an external diameter larger than 7 French (2.33 mm, 0.092") through the haemostatic valve.

5 Contra-indications

- **Super Ketch™** or **Speed Ketch™** are not intended to be used with an automatic injector.
- **Super Ketch™** or **Speed Ketch™** valves are not intended to be used for injections with a pressure > 8 atm (118 psi).
- Do not inject liquid if air bubbles are visible inside the **Super Ketch™** or **Speed Ketch™**.

6 Possible complications

Potential complications include, but are not limited to:

- Contrast crystal formation
- Clotting
- Emboli

7 Instructions for use

- Connect the **Super Ketch™** or **Speed Ketch™** side access to the perfusion (manifold) device using a stopcock if needed (**Speed Ketch™**).
- Set up continuous flush with a normal saline solution to remove trapped air and connect the male rotating luer lock of the **Super Ketch™** or **Speed Ketch™** to the guiding catheter or any other device chosen by the operator.

- To flush the connector segment comprising the valve, open the haemostatic valve by pushing the push-pull button valve opener and continue to fill the assembled devices. To close the valve, pull-out the push-pull button valve opener into its original position.

Make sure by visually checking that no air bubble is trapped into the **Super Ketch™** or **Speed Ketch™** device.

- Aspirate any air that may have entered the system before the introduction of the devices, to verify that no air bubble has been captured.
- Establish pressure above arterial pressure in order to prevent any retrograde flow of blood into guiding catheter or any other device chosen by the operator. Check that all connections are secure so that air is not introduced during continuous flush.
- Introduce the guiding catheter or any other device chosen by the physician following the recommendations of respective device manufacturers.
- Insert the guidewire, or the guidewire and dilatation catheter together, or any chosen device into the **Super Ketch™** or **Speed Ketch™** straight lumen.

The valve shall be opened using the push-pull button valve opener, when any device is inserted through the valve into the straight lumen of the **Super Ketch™** or **Speed Ketch™**.

- The metallic insertion tool supplied with the **Super Ketch™** and the **Speed Ketch™** must be used when the guidewire alone is inserted through the valve, in order to protect the guidewire distal tip.
- The plastic insertion tool provided with the **Speed Ketch™** must be used when a microcatheter alone is introduced through the valve, in order not to damage the device when passing through the closed valve.
- Perform the rest of the procedure following the recommendations of respective device manufacturers.

Note: When introducing or retrieving devices, it is important to manipulate the devices slowly and with care. Rapid accelerations can provoke a suction phenomenon at the level of the valve, which could introduce air into the system.

Rapid accelerations can also provoke coating damage.

Torquer instructions for use

- The maneuverability of the guidewires when introduced through the inner lumen of **Super Ketch™** can be improved by the use of the Torquer supplied as an accessory (**Super Ketch™ Plus**).
- Torquer instruction of use:
 - Unscrew the cap of the Torquer.
 - Insert the proximal end of the guidewire into the funnel-shaped hole on the distal end of the Torquer. Once the Torquer is positioned at the desired location, tighten the Torquer to secure the Torquer to the guidewire.
 - To move the Torquer to a new position, loosen the proximal end of the Torquer, slide the device along the guidewire to the desired location, and tighten the Torquer.

8 Packaging

- **Super Ketch™** and **Speed Ketch™** Y connectors are delivered in an individual peelable pouch. 10 units are delivered per cardboard box
- STERILE. Sterilised with Ethylene Oxide Gas. Non pyrogenic.

9 Storage requirements

- Use before the expiry date clearly indicated on the label.
- Store at room temperature below 40°C, in a dry place, protected from light.

10 Liability

Minvasys has endeavoured to ensure that the products comply with all relevant standards and regulations currently in force and to ensure that the quality of the products meets the requirements of the above mentioned standards and regulations for a period ending upon the indicated expiry date. The above statement does not apply when the products are used for a purpose other than its intended purpose. Where any loss or damage is caused (other than death or personal injury) due to a defective product, Minvasys shall not be liable for such loss or damage.

11 Conversion Chart

1 French	0.0131"	0.33 mm	
1 bar	0.99 atm.	14.5 PSI	10 ⁵ Pa

12 Shelf life

4 years.

13 Symbols

Quantity per
package

Maximum inner
diameter

Do not use
damaged
packaging

Do not
resterilize

Year CE marking obtained: 2004

7, rue du Fossé Blanc, Bâtiment C1
92230 Gennevilliers – France
Tel: +33 (0)1 47 90 70 30
Fax: +33 (0)1 47 91 05 85

Legal representative in Russia:

Delrus
No. 23, Posadskaya Str.,
Office 810
Yekaterinburg 620086, Russia

Registration number: XXX

Y connectors information

A. Description of the device


Figure 1: Y connectors with side arm


Figure 2: Y connectors without side arm

Measures are in millimetres.

1. Valve opener: Transparent
2. Bride-valve: Transparent
3. Valve: Translucent
4. Y hub: Transparent
5. Intermediate ring: Transparent
6. Joint: Translucent
7. Luer Lock rotating: Transparent
8. Side arm (20cm for KESU0001 and KESU0002 and 50 cm for KESU0003 and KESU0005): Transparent
9. Stopcock:
 - Transparent blue
 - White
 - Blue

B. Technical information

The side arm tubing (8) has a tear force higher than 20N.

The Y connectors pressure resistance is as follow:

- Without material inserted:
 - Air leak test:
 - no leakage at 0.3bar
 - no leakage at 5bar
 - Liquid leak test: 0 drop per 30s at 2bar
- With 0.014" guidewire + PTCA hypotube inserted
 - Liquid leak test:
 - ≤ 30 drops per 30s at 0.3bar
 - ≤ 30 drops per 30s at 2bar
 - ≤ 30 drops per 30s at 8bar

To prevent any deterioration of the valve (3) that could lead to embolia, the valve thickness dimension is monitored. Furthermore, a bride-valve (2) is used to maintain the valve in place on the Y connectors. The bride valve is screwed with a controlled force on the Y connector hub in addition of being paste up to it.

Figure 3 describes Y connectors when connected with other devices.


Figure 3: Y connectors connected with other devices

Super Ketch Plus Inflation Kit

Y- connectors and accessories

Sterile (sterilized with ethylene oxide). For single use. Do not re-use. Not to use if packaging is opened or damaged. Use before the expiry date mentioned on the packaging.
Store away from light in a cool and dry place.



Caution

Reuse of the Y-Connector is strictly prohibited. If a product is reused its functionality would be entirely altered. In addition, the iatrogenic risk from the inter-patient cross-transmission would be high. Please refer to the indications provided by the manufacturers of devices for use with the Connector. After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.

Description

Y connectors

The Y-Connectors comprise a tricuspid star shaped valve, designed to provide haemostatis during endovascular procedures. The Y-Connectors include a straight and a side access lumen. The haemostatic valve, located at the proximal end of the straight lumen, can be opened and closed by operation of the "push-pull" button valve. The distal end has a rotating male Luer lock. The side access lumen allows pressure monitoring and perfusion. The side access lumen of the Y-connectors with a side arm tubing ending by a standard three-way stopcock.

Metallic Insertion Tool

The metallic insertion tool is an accessory designed to facilitate the introduction of guidewires (up to 0.018") used in conjunction with interventional devices during vascular interventional procedures.

Torquer

This device is designed to facilitate steering of 0.010" to 0.035" (0.89mm) guidewires.

Indication

The Y-Connectors are intended to maintain haemostatis before and during the introduction / the use and withdrawal of diagnostic / interventional devices with outside diameter $\leq 7F$ (2.33 mm, 0.092") during interventional procedures.

Contra-indications

- The Y connectors are not intended to be used with an automatic injector.
- The Y connectors' valves are not intended to be used for injections with a pressure > 8 atm (118 psi).
- Do not inject liquid if air bubbles are visible inside the Y connectors

Precaution for use

Any aspiration through the side access of the Y-Connector should be effectuated slowly and carefully to avoid the creation of a void in the proximal part of the connector, leading to an effect of suction on the haemostatic valve.

Always open the haemostatic valve using the push-pull button valve opener before inserting any device with a fragile distal portion, in order to avoid damage of the product when passing through the valve.

Always retrieve the metallic insertion tool before retrieval of the guidewire, to avoid any damage of the guidewire.

The metallic insertion tool should be inserted carefully through the valve in order not to damage or dislodge it.

The Y-Connector devices should exclusively be used by physicians trained in interventional procedures. Do not attempt to force an external diameter device larger than 7 French (2.33 mm, 0.092") through the haemostatic valve.

Instructions for use

1. Connect the Y-Connector side access to the perfusion (manifold).
2. Set up continuous flush with a normal saline solution to remove trapped air and connect the male rotating Luer lock of the Y-Connector to the guiding catheter or any other device chosen by the operator.
3. To flush the connector segment comprising the valve, open the haemostatic valve by pushing the push-pull button valve opener and continue to fill the assembled devices.
To close the valve, pull-out the push-pull button valve opener into its original position. Make sure by visually checking that no air bubble is trapped into the Y-Connector device.
4. Aspirate any air that may have entered the system before the introduction of the devices, to verify that no air bubble has been captured.
5. Establish pressure above arterial pressure in order to prevent any retrograde flow of blood into guiding catheter or any other device chosen by the operator. Check that all connections are secure so that air is not introduced during continuous flush.
6. Introduce the guiding catheter or any other device chosen by the physician following the recommendations of respective device manufacturers.
7. Insert the guidewire, or guidewire and dilatation catheter together, or any chosen device into the Y-Connector straight lumen.

Warning: The valve shall be opened using the push-pull button valve opener, when any device is inserted through the valve into the straight lumen of the Y-Connector.

- a. The metallic insertion tool supplied with the Y-connector must be used when the guidewire alone is inserted through the valve, in order to protect the guidewire distal tip.
- b. Perform the rest of the procedure following the recommendations of respective device manufacturers

Note: When introducing or retrieving devices, it is important to manipulate the devices slowly and with care. Rapid accelerations can provoke a suction phenomenon at the level of the valve, which could introduce air into the system.

Torquer instructions for use

1. The maneuverability of the guidewires when introduced through the inner lumen of Y-Connector can be improved by the use of the Torquer supplied as an accessory.
2. Torquer instruction of use:
 - Unscrew the cap of the Torquer.
 - Insert the proximal end of the guidewire into the funnel-shaped hole on the distal end of the Torquer. Once the Torquer is positioned at the desired location, tighten the Torquer to secure the Torquer to the guidewire.
 - To move the Torquer to a new position, loosen the proximal end of the Torquer, slide the device along the guidewire to the desired location, and tighten the Torquer.

Specifications

Components Specifications		
Y connector	Inner Lumen Compatibility:	7F / 0.092" / 2.33 mm maximum
	Pressure Resistance ⁽¹⁾ after introduction of a device:	8 atm. maximum
Metallic Insertion Tool	External Diameter	2.7 F / 0.90 mm
	Inner Lumen Diameter	1.9 F / 0.025" / 0.64 mm
	Usable Length	10 cm
Plastic Peelable Insertion Tool	Length	9 cm
	Inner Diameter	3.4 F / 0.044" / 1.14 mm
Torquer	Guidewire Compatibility	From 0.010" minimum to 0.035" (0.89mm) maximum

(1) Haemostatic valve resistance measured in vitro, by applying injection pressure through the side port.

Inflation Device

Sterile (sterilised with ethylene oxide). Do not re-sterilise. For single use. Do not use if the packaging is opened or damaged. Use before the expiry date shown on the packaging. Use with balloon dilatation catheters. Read the instructions before use. Store away from light in a cool and dry place.

Caution

Reuse of the Inflation Device is strictly prohibited. If a product is reused its functionality would be entirely altered. In addition, the iatrogenic risk from inter-patient cross-transmission would be high.

After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.

Description:

The Inflation Device is for single use. It's consisting of a 20 cc syringe body with a plunger which may be engaged by a control button, a rotatory handle, a manometer and a high pressure connector with a rotating male luer lock. The manometer is graduated between vacuum and 30 atm with an equivalent scale in PSI. There are two gripping choices by changing manometer's position.

Indications:

The Inflation Device is intended for use during cardiovascular procedures to create, maintain and monitor pressure in the balloon catheter.

Operation:

The plunger is free to move when the control button is rocked towards the back of the device (opposite point of manometer).

To lock the plunger, push completely on the control button towards the manometer. Turn the handle clockwise to inject the mixture of contrast media. The estimated pressure is indicated on the manometer.

Preparation:

1. Prepare the mixture of contrast media for balloon catheter following the manufacturer's recommendations for the contrast media.
2. Check the plunger is unlocked.
3. Hold the Inflation Device with its manometer downwards and insert the end of connector into the contrast solution to fill the syringe.
4. Hold the device with the manometer upwards to purge the air contained in the syringe and in the connector.
5. Connect the Inflation Device to the balloon catheter.
6. To inflate the balloon catheter, point the manometer downwards and turn **slowly** the handle clockwise up to posting of desired pressure. This pressure is kept thanks to the locking system.
7. To deflate the balloon catheter, unlock the plunger and pull completely the handle.

Precautions for use:

The Inflation Device must be used by doctors. Before each inflation, confirm the needle of the manometer is set to the zero position +/- 1 atm. When releasing the deflation, it is recommended to go with the handle. If the control knob can't return in unlocked position, turn the handle anticlockwise to deflate the balloon catheter. Refer to the balloon catheter instructions for use (precautions for use, maximum pressure, etc).

If the key is difficult to lock, do not force it as this could damage the plunger. Turn the handle gently and re-engage the key.

Shelf life:

4 years.

Year CE marking obtained: 2009

7, rue du Fossé Blanc, Bâtiment C1
92230 Gennevilliers – France
Tel: +33 (0)1 47 90 70 30
Fax: +33 (0)1 47 91 05 85

Legal representative in Russia:

Delrus

No. 23, Posadskaya Str.,

Office 810

Yekaterinburg 620086, Russia

Registration number: XXX

SuperKetch Plus inflation kit information

A. SuperKetch


Figure 1: SuperKetch (dimension in millimeters)

1. Valve opener: Transparent
2. Bride-valve: Transparent
3. Valve: Translucent
4. Y hub: Transparent
5. Intermediate ring: Transparent
6. Joint: Translucent
7. Luer Lock rotating: Transparent
8. Side arm (20cm for KESU0001 and KESU0002 and 50 cm for KESU0003 and KESU0005): Transparent
9. Stopcock:
 - Transparent blue
 - White
 - Blue

The side arm tubing (8) has a tear force higher than 20N.

The Y connectors pressure resistance is as follow:

- Without material inserted:
 - Air leak test:
 - no leakage at 0.3bar
 - no leakage at 5bar
 - Liquid leak test: 0 drop per 30s at 2bar
- With 0.014" guidewire + PTCA hypotube inserted
 - Liquid leak test:
 - ≤ 30 drops per 30s at 0.3bar
 - ≤ 30 drops per 30s at 2bar
 - ≤ 30 drops per 30s at 8bar

To prevent any deterioration of the valve (3) that could lead to embolia, the valve thickness dimension is monitored. Furthermore, a bride-valve (2) is used to maintain the valve in place on the Y connectors. The bride valve is screwed with a controlled force on the Y connector hub in addition of being paste up to it.

Figure 2 describes Y connectors when connected with other devices.


Figure 2: Y connectors connected with other devices

B. Inflation device



Figure 3: Inflation device (dimension in millimeters)

1. Handle
2. Button
3. Body
4. Manometer
5. Connection line

The Inflation device consists in a 20 mL syringe body with a plunger which may be engaged by a control button, a rotatory handle, a manometer and a high-pressure connector with a rotating male luer lock. The manometer is graduated between vacuum and 30 atm with an equivalent scale in PSI. There are two gripping choices by changing manometer's position.

The document contains 44 bounds, numbered, signed and stamped pages.

Luc DEMANGEON
Development, Regulatory & Quality Director
MINVASYS

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