1001

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IA 12 Sabine Schmidt

VB Sabine Schmidt

Instructions for use of the medical device
Synovial fluid substitute «SYNOCROM®», sterile 1% solution for intra-articular injection, in 2 ml syringes

1. Brief description

Synovial fluid substitute «SYNOCROM®», sterile 1% solution for intra-articular injection, in 2 ml syringes, (hereinafter SYNOCROM®) is a sterile ready to use 1 % aqueous solution of sodium hyaluronate with a molecular weight of ~ 1.6 million daltons in syringes of 2 ml. Used as a substitute for synovial fluid to restore the natural viscoelastic properties of the synovial fluid when injected into synovial joints (viscosupplementation therapy) and prescribed for pain and limitation of mobility due to degenerative or post traumatic changes in the joints that contains synovial fluid.

Method of sterilization: steam (high-pressure syringe in blister).

SYNOCROM® is a sterile non-pyrogenic isotonic viscoelastic solution for intra-articular injection.

2. Composition

1 syringe (2 ml) SYNOCROM® contains the active ingredient sodium hyaluronate; excipients: sodium chloride, disodium hydrogen phosphate dodecahydrate, citric acid, water for injections.

2.1. Information on medicines contained in the medical product

The medical device contains Sodium hyaluronate.

2.2. Information on materials of animal or human origin contained in the medical product

Materials of animal or human origin are not used in the production process of a medical product or as a raw material.

3. Indications for use

The pains and movement restriction as a result of degenerative or post-traumatic changes in the joints, containing synovial fluid.

4. Principles of Operation of the Medical Device

All synovial joints contain viscoelastic sodium hyaluronate.

This substance has lubricating and absorbing properties, which provide normal, painless movement of these joints.

Viscoelasticity of synovial fluid significantly decreases at degenerative diseases (osteoarthrosis) or post-traumatic changes in the joints, as a result of which mechanical loading on the joint increases and the destruction of articular cartilage enhances (that manifests by restriction and pains on movement in the joints). Intra-articular administration of highly purified sodium hyaluronate with optimized viscoelastic properties allows to improve the quality of joint lubrication.

Lubricating and absorbing properties of SYNOCROM® allow to decrease the pain and to increase joint mobility.

This effect may remain within several months after performing the course of application, consisting of five intra-articular injections.

5. Method of Use

Take the ready-to-use syringe out of the blister, remove the rubber cap from the Luer-Lock adapter by turning it gently, attach an appropriate sterile needles and lock it into place by turning carefully.

The injection site is treated with a disinfectant substance containing alcohol. SYNOCROM® is injected accurately into the joint cavity by standard methods, given the anatomical features. Depending on the size of the joint, up to 2 ml may be administered intra-articularly.

In the knee, three to five injections of the product at weekly intervals are recommended. In hip, ankle and shoulder, the number of injections may range from one to five and should be chosen based on the clinical experience and the relief obtained by the patients. Several joints may be treated simultaneously.

Course of 5 injections maybe effective for more than 6 months (depending on disease severity degree). Repeated courses of application are acceptable. SYNOCROM® administration in joint cavity should be stopped at pain appearance during injection.

The unconsumed solution is not applicable for storage.

6. Contraindications

SYNOCROM® should not be applied in patients with increased sensitivity to one of medical product's component.

The patients with bacterial arthritis should refrain from application this product in order to avoid possible complications because of intra-articular administration of SYNOCROM®.

Doctor in charge should also be familiar with all immunological and other potential risks associated with the use of biological material.

7. Possible Adverse Effects of Medical Device Use

In patients, applying SYNOCROM®, a manifestation of local symptoms (pain, feeling of heat, redness and swelling) may occur. After administration of the product, it is recommended that an ice-pack be placed on the treated join for 5-10 minutes in order to prevent pain and swelling.

It was reported the following adverse events that occur in the treatment of similar products: mild or moderate arthralgia, in rare cases skin rashes, aseptic joint effusion, pruritus and muscle cramps. In very rare cases was observed following adverse events: allergic reaction, anaphylactic shock, hemarthrosis, phlebitis, pseudosepsis, severe acute inflammatory reaction, nasopharyngitis, joint stiffness, tendonitis, bursitis, fever, and myalgia.

In the case of pronounced manifestations, accompanied by severe pain, the solution should be removed from the affected joint.

8. Safety Requirements

SYNOCROM® assigned only for intra-articular administration. During the first 24 hours after the procedure the patient is advised to adhere to a sparing mode, and to avoid an excessive load on the joint.

Application of SYNOCROM® was not studied in pregnant women and children younger 18 years.

SYNOCROM® should be kept out of reach of young children.



SYNOCROM® assigned only for single use. Repeated sterilization it is not allowed. Repeated use of product is associated with a potential risk of infection to the patient or physician.

Do not use the syringe from open and/or damaged sterile blister.

Do not use the syringe with open and/or damaged cap.

Do not apply after the expiration date.

It is necessary to hold the syringe with SYNOCROM® during injection as is shown in Pic.1.

Picture 1

9. **Possible Combinations with Other Equipment**

No information available at present time about incompatibility of SYNOCROM® with other solutions for intra-articular administration.

One should not apply disinfectants, containing quaternary ammonium salts, since hyaluronic acid sediments in the presence of these substances.

SYNOCROM® syringe is compatible with any injection needles that have a Luer Lock connector. The product must be used with an appropriate needles (considering the body type of the patient and the type of joint, in which the injection is administered) and used only by physicians trained in the technique of intra-articular injections

10. **How supplied and packing**

2 ml of a 1% solution for intra-articular injection in a syringe made of glass.

One syringe without a needle, in a contour acheikova packing (blister). Needle for injection not supplied.

One contour acheikova package together with the instructions for use in a cardboard pack. The pack has a sticker from self-adhesive material on both sides for control opening.

11. **Sterility**

Sterilization method: by steam of the high pressure (syringe in a blister) in accordance with EN ISO 17665-1 for 15 minutes at 121°C.

Sterility of the gel is guaranteed inside the undamaged syringe-stopper system.

12. **Methods and Detergents for Disinfection and Pre-cleaning**

Not Applicable.

13. **Storage condition**

Store in a dry place at temperatures from 2 to 25 °C.

Do not heat, do not freeze, protect from direct sunlight.

14. **Transportation conditions**

Products should be transported by all kinds of roofed transport in accordance with cargo carriage regulations applicable to each mode of transport at the temperature from 2 to 25 °C.

15. **Shelf life**

3 years.

16. **Maintenance**

Not applicable.

17. **Repairs**

Not applicable.

18. **Medical Device Recycling or Disposal**

According to the degree of danger SYNOCROM® belongs to a class B (epidemiologically hazardous waste) and must be disposed of in accordance with SanPiN 2.1.7.2790-10. There are no specific disposal instructions for the gel or the glass syringe.

19. **This product is designed, manufactured, tested and packaged in compliance with all relevant requirements.**

The manufacturer guarantees the quality of the medical device in the blister before the expiry of its validity, provided that the packaging is not damaged and the storage conditions are met.

Warranty period corresponds to the period of validity is 3 years

20. **Registration certificate**

№ ΦC3 2010/08100 from October 15, 2010

21. **Manufacturer**

Croma-Pharma GmbH,

Industriezeile 6, 2100, Leobendorf, Austria

www.croma.at

22. The authorized representative of the manufacturer/Packer

ZAO «PharmFirm «Sotex»

bid. 11, village Belikova, rural settlement Bereznikovskoe, Sergiev-Posad municipal district,
Moscow Region, 141345, Russia

Tel/fax: (495) 956-29-30

www.sotex.ru

Customer complaints to be send to the address of ZAO «PharmFirm «Sotex».

23. International List of symbols

Symbol	Meaning
	Prohibition on repeated administration
	Caution! Refer to administration instructions
	Temperature range
	Sterilization using steam
	Prevent the action of sunlight
	Keep dry
	Do not use if package is damaged
	Manufacturer.
	Use by
	Batch number
	Indication of the location of the package opening

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