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Description:

ESTELINE Vita 18 mg/ml hyaluronic acid-based viscoelastic intradermal implant [1.8% Sodium Hyaluronate with 2% Glycerol i.d.] is a sterile, biodegradable, viscoelastic, clean, transparent, isotonic and homogenized injectable gel implant. Medical device contains hyaluronic acid produced by *Streptococcus equi*, normalized to the total concentration of 18 mg/ml in physiological buffer, and glycerol for the regulation of the osmolality of the product and purging radicals.

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intra dermal implant is intended to be used for adult human beings, to replenish the loss of natural dermal hyaluronic acid that may occur due to aging. Thereby it may improve moisture, tone and elasticity of the skin and/or to act as filler for small lines such as crow's feet, smile lines or smoke lines surrounding the mouth. The device is designed to be injected into the superficial dermal tissue where it is totally resorbed. The device is totally resorbed.

The device is intended for single use and is sterile.

The medical device is non-active and has no measuring function.

No materials of animal origin are used during manufacturing or as raw materials.

No medicinal products are integrated into the device.


Scope of Delivery of the Medical Device:

ESTELINE Vita 18 mg/ml hyaluronic acid-based viscoelastic intra dermal implant consists of:

1. Hyaluronic acid 18 mg/ml (sodium hyaluronate (raw material) manufactured by HTL S.A.S.)
2. 1.0 ml syringe with a backstop
3. 0.3 mm x 13 mm sterile single use needle (30G) - 2 pc.
4. Instructions for use
5. Label – 3 pc.

A syringe filled with gel implant and two 30G sterile needles are packaged into a blister sealed with a blister cover. The blister is packaged into a labelled carton box containing an instruction for use. The carton box has three labels with a batch number. One label (identification clamp) remains on the package, the second label shall be attached to the patient chart and the third label shall be given to the patient for control implementation.

Picture of the Medical Device General View:

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0.3 mm x 13 mm sterile single use needle (30G)

No.	Item	Material
7	Cannula	Stainless steel AISI 304
8	Luer hub	Polypropylene Yellow pigment PMS 108
9	Tip cap	Polypropylene unpainted

ATTENTION!

Read the instruction before use.

Indications:

- replenishment of the loss of natural dermal hyaluronic acid
- improvement of facial hydration
- improvement of skin tone and elasticity
- correction of small lines
- correction of crow's feet, smile lines and smoke lines surrounding the mouth.

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Target patient population:

Adults with facial dermal and soft tissue defects.

Contraindications:

The medical device shall not be used in:

- patients who tend to develop hypertrophic scarring, pigmented lesions or sensitivity to keloid formations;
- patients with a history of autoimmune disease or who are receiving immune therapy;
- patients who are known to be hypersensitive to hyaluronic acid;
- patients with previously implanted permanent fillers in the treated area;
- pregnant or breastfeeding women;
- patients under 18 years of age.

Patients receiving anticoagulants or patients receiving blood platelet aggregation inhibitors (for example, acetylsalicylic acid) shall not be treated with this medical device without consulting a physician.

The medical device shall not be used in areas with skin inflammatory and/or infectious processes (such as acne, herpes ...). The device shall not be used simultaneously with laser therapy, chemical peeling, dermabrasion or mesotherapy.

Precautions for Use:

The device shall be used only by physicians who have completed special training in wrinkle filling injections in sterile premises.

Sensitive skin can be pretreated with a topical anesthetic patch or cream. The syringe with remaining gel implant and the needle shall be disposed of in a special container immediately after use.

Undesired Side Effects:

Physicians must inform the patient that there are potential side effects and/or

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incompatibilities associated with implantation of this device, which may occur immediately or may be delayed.

The use of this device or similar devices has been associated with the following events and reactions:

abscess, angioedema, bacterial infections, nodules/induration, bleeding, ecchymoma, burning pain, discoloration (bleaching), discomfort, oedema, erythema, callosity, granuloma, hematoma, herpes reactivation, hypersensitivity, accumulation of implant at the injection site, implant visibility, hardening, inflammation, irritation, itching, marble pattern on skin after inadvertent arterial embolisation, jugal oedema, migration, necrosis due to vascular disorders, nodes (inflated and non-inflated), numbness, pain, paresthesia, puncture marks, redness, retinal artery occlusion, scar sarcoidosis, scleromyxedema (generalized), swelling, telangiectasia, painfulness, vasculitis, angiospasm, vasovagal reaction during injection, loss of vision due to retinal artery occlusion.


Injectations of other fillers for soft tissues into glabellar region and nose alae in very rare instances may lead to retinal artery occlusion. In rare instances, the herpes virus reaction is possible due to direct damage of neural axons with a needle or manipulations with tissues and inflammatory reaction after injection. Dissection of subepidermal plane due to fan-shaped needle introduction may increase the risk of local adverse events. Too superficial placement of the filler or uneven distribution of the injected implant may lead to visible anaemic nodes in skin. It is therefore important to take these possible complications into account. Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other adverse events. The physician should treat these as appropriate. Any other undesirable effects associated with injection of gel implant must be reported to the authorized representative in Russia and/or manufacturer.

ATTENTION!

The device shall not be used by inexperienced medical staff.

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

Method of Correct Tip Cap Removal:

Hold the syringe at the adapter as shown (1). Rotate the tip cap counterclockwise with the other hand (2).

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

ATTENTION!


Do not use the syringe, if the tip cap is open or displaced in the package.

Method of Needle Attachment:

Hold the syringe as shown (4). Attach the supplied 0.3 mm x 13 mm needle (30G) as shown (5) (do not use other needles).

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


To attach, rotate the needle clockwise firmly holding it until a slight click (6).


Method of Syringe Fixation in the Hand:

When injecting gel implant, hold the syringe as shown in picture 1.

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

CORRECT! (a slot in the backstop is directed backwards, towards the palm)

INCORRECT! (a slot in the backstop is directed forward, away from the palm)

Method of Use:

The device must be injected into non-inflamed, disinfected, healthy skin. Using the described injection technique has a great importance for obtaining

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a successful result. Use the 30G needle which is provided with the syringe and inject slowly by applying the appropriate injection technique. Consider that the device will attract water and will increase its volume in situ. Avoid injecting large volumes in the periorbital area. The amount injected will depend on the skin's condition and need. Initial treatment can be supplemented with one or two additional sessions with a three-week interval between sessions, depending on the skin's condition. After the injection, doctors may apply a light massage in order to distribute the product uniformly.

Warnings:

- Check the syringe integrity and shelf-life before use.
- Do not use the syringe if its tip is open or displaced inside the protective package.
- Do not use other needles or syringes instead of those supplied by the manufacturer.
- Do not modify or bend the needle.
- Do not reuse; the quality and sterility can be guaranteed only for syringes in the original package. Reuse of the device may result in infection of patients or users.
- Should 30G needle be blocked, do not increase pressure on the barrel, stop injection and replace the needle.
- There are no clinical data (effectiveness and tolerability) on gel implant injection into the area, which has previously been treated with other gel implant.
- Do not inject into blood vessels, bones, tendons, ligaments, nerves and muscles.
- Do not inject into the eye contour area (including tear trough, upper periorbit and eyelids), glabellar region and nose alae.
- Do not inject gel implant into birthmarks.
- Avoid excessive correction.
- Patients shall be recommended to refrain from makeup within 12 hours after the injection and to avoid long exposure to sunlight and ultraviolet rays and visits to saunas and baths within one week after the injection.
- In order to avoid possible migration of gel implant, the patient shall be

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
implant medical device – steam sterilization
0.3 mm x 13 mm sterile single use needle (30G) – gas sterilization (Ethylene oxide).

Shelf life:
3 years with observation of transportation and storage conditions.

Disposal:
Medical device waste class of hazard is B (epidemiologically hazardous waste). Used devices, devices with the expired shelf life or damaged package shall be disposed of according to the established procedure in compliance with the sanitary epidemiological requirements for medical waste treatment SanPiN 2.1.7.2790-10.
Collection and processing of medical waste at healthcare facilities shall be controlled by the representatives of the disposal company. For this purpose, special rigid containers with vacuum covers shall be used. All containers shall be painted yellow indicating that the waste is dangerous.
After that, waste of the same type is disposed of at incinerators or disinfected and stored for further dumping at storage premises.

Label Symbols:

No.	Symbol	Description	No.	Symbol	Description
1		Manufacturer	2		Use by
3		Manufacturing date	4		Do not resterilize
5		Do not reuse	6		Caution!
7		Lot number	8		Do not use if package is damaged
9		Consult instructions for use	10		Contains no latex

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

Эл. почта: medmelis@yandex.ru

11		Keep away from direct sunlight	12		Keep dry
13		Sterilized using ethylene oxide	14		Sterilized using steam
15		Recommended temperature	16		Recommended humidity
17		Sterile single use needle	18		Sterile single use syringe
19		Authorized representative in the Russian Federation	20		Device registration certificate in Russia
21		CE mark	22		Needle manufacturer

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ESTELINE Voluma 23 mg/ml Hyaluronic Acid-Based Viscoelastic
Intradermal Implant

1ml of gel implant contains:

Substance	Function	Quantity
Sodium hyaluronate	Active ingredient	23 mg
Sodium chloride	Buffer salt	8.3 mg
Sodium dihydrogen phosphate × 12H ₂ O	Buffer salt	0.95 mg
Sodium dihydrogen phosphate × 2H ₂ O	Buffer salt	0.4 mg
Water for injections	Solvent	980.6 mg

Description:

ESTELINE Voluma 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant [2.3% Highly Crosslinked Sodium Hyaluronate i.d.] is a sterile, biodegradable, viscoelastic, clean, transparent, isotonic and homogenized injectable gel implant. The hyaluronic acid used in this device is of a bioenzymatic origin (*Streptococcus equi*) and contains a certain part of linked hyaluronic acid and a certain part of unlinked hyaluronic acid, which are in equal amount produced by streptococci and normalized to a total concentration of 23 mg/ml in a physiological buffer.

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Intended Use and Application of the Medical Device:

ESTELINE Voluma 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant is intended to be used for adult human beings to correct wrinkles and folds, increase or restore volume of facial tissues and to remodel facial contours. The product is designed to be injected into the dermis or subcutis. The device is completely absorbed.

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ESTELINE Voluma 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant consists of:

1. Hyaluronic acid 23 mg/ml (sodium hyaluronate (raw material) manufactured by HTL S.A.S.)
2. 1.0 ml syringe with a backstop
3. 0.4 mm x 12.7 mm sterile single use needle (27G) - 2 pc.
4. Instructions for use
5. Label – 3 pc.

A syringe filled with gel implant and two 27G sterile needles are packaged into a blister sealed with a blister cover. The blister is packaged into a labelled carton box containing an instruction for use. The carton box has three labels with a batch number. One label (identification clamp) remains on the package, the second label shall be attached to the patient chart and the third label shall be given to the patient for control implementation.

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Read the instruction before use.

Indications:

- correction of deep facial wrinkles and folds
- increase of facial tissues' volume
- restoration of facial contour
- treatment of facial loss
- remodeling of face line
- treatment of destructive scars
- treatment of morphological asymmetry

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Undesired Side Effects:

Physicians must inform the patient that there are potential side effects

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and/or incompatibilities associated with implantation of this device, which may occur immediately or may be delayed.

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
Injecting other fillers for soft tissues into glabellar region and nose alae in very rare instances may lead to retinal artery occlusion. In rare instances, the herpes virus reaction is possible due to direct damage of neural axons with a needle or manipulations with tissues and inflammatory reaction after injection. Dissection of subepidermal plane due to fan-shaped needle introduction may increase the risk of local adverse events. Too superficial placement of the filler or uneven distribution of the injected implant may lead to visible anaemic nodes in skin. It is therefore important to take these possible complications into account. Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other adverse events. The physician should treat these as appropriate. Any other undesirable effects associated with injection of gel implant must be reported to the authorized representative in Russia and/or manufacturer.

ATTENTION!

The device shall not be used by inexperienced medical personnel.

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



Method of Needle Attachment:

Hold the syringe as shown (4). Attach the supplied 0.4 mm x 12.7 mm needle (27G) as shown (5) (do not use other needles).

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

Method of Use:

Gel implant shall be injected into uninflamed, disinfected healthy skin. Using the described injection technique has a great importance for

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obtaining a successful result. Use 27G needle, which is supplied with a syringe, and inject slowly using the appropriate technique. Inject small volumes in two or several steps but not the whole volume at once. The injection volume depends upon the type of wrinkles to be corrected. Following the injection, the physician may gently massage the area for uniform distribution of gel implant. Leave the patient in the room for several minutes after the injection so that blanching due to artery embolism could be detected.

Warnings:

- Check the syringe integrity and shelf-life before use.
- Do not use the syringe if its tip is open or displaced inside the protective package.
- Do not use other needles or syringes instead of those supplied by the manufacturer.
- Do not modify or bend the needle.
- Do not reuse; the quality and sterility can be guaranteed only for syringes in the original package. Reuse of the device may result in infection of patients or users.
- Should 27G needle be blocked, do not increase pressure on the barrel, stop injection and replace the needle.
- There are no clinical data (effectiveness and tolerability) on gel implant injection into the area, which has previously been treated with other gel implant.
- Do not inject into blood vessels, bones, tendons, ligaments, nerves and muscles.
- Do not inject into the eye contour area (including tear trough, upper periorbit and eyelids), glabellar region and nose alae.
- Do not inject gel implant into birthmarks.
- Avoid excessive correction.
- Patients shall be recommended to refrain from makeup within 12 hours after the injection and to avoid long exposure to sunlight and ultraviolet rays and visits to saunas and baths within one week after the injection.
- In order to avoid possible migration of gel implant, the patient shall be

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recommended not to massage the treated area within several days after the injection.

- Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solution. For this reason, the gel implant shall never contact these substances or other medical surgical tool, which have previously contacted these substances.

Result Preservation Duration:

The good cosmetic results and high degree of patient satisfaction are observed after applying this medical device. The successful result is preserved for 24 months.

Removal of Gel Implant in Case of Unsuccessful Implantation.

Should it be necessary to remove gel implant from tissues, a medicinal product based on hyaluronidase active substance is injected. Inject locally into the implant area.

Storage:

The medical device shall be stored at $+2^{\circ}\text{C} - +25^{\circ}\text{C}$ in dry premises in an original package and protected against light, high and low temperatures. Handle with care.

Transportation:

The medical device can be transported in covered vehicles of all types depending upon the rules of carriage for such types of vehicles at temperatures $+2^{\circ}\text{C} - +25^{\circ}\text{C}$ and humidity of up to 40 %.

Individual Package of the Device:

- 1 unit of the medical device

Transportation Box:

- 168 individual packages of products in a single transport packaging products.

Sterilization method:

ESTELINE Voluma 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant medical device – steam sterilization

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TEMPLATE OF THE MEDICAL DEVICE INSTRUCTIONS FOR
USE ON THE TERRITORY OF THE RUSSIAN FEDERATION
(initial version)

Medical Device Instructions For Use

ESTELINE Basic 23 mg/ml Hyaluronic Acid-Based Viscoelastic
Intradermal Implant

1ml of gel implant contains:

Substance	Function	Quantity
Sodium hyaluronate	Active ingredient	23 mg
Sodium chloride	Buffer salt	8.3 mg
Sodium dihydrogen phosphate × 12H ₂ O	Buffer salt	0.95 mg
Sodium dihydrogen phosphate × 2H ₂ O	Buffer salt	0.4 mg
Water for injections	Solvent	979.8 mg

Description:

ESTELINE Basic 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant [Crosslinked Hyaluronic Acid 2.3% i.d.] is a sterile, biodegradable, viscoelastic, clean, transparent, isotonic and homogenized injectable gel implant. The hyaluronic acid used in this device is of a bioenzymatic origin (*Streptococcus equi*) and contains a certain part of linked hyaluronic acid and a certain part of unlinked hyaluronic acid, which are in equal amount produced by streptococci and normalized to a total concentration of 23 mg/ml in a physiological buffer.

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ESTELINE Basic 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant is intended to be used for adult human beings, as a soft tissue filler in order to correct moderate to severe wrinkles and folds and to increase lip volume. The device is designed to be injected into the mid to deep dermis.

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ESTELINE Basic 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant consists of:

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Read the instruction before use.

Indications:

- correction of moderate to severe facial wrinkles and folds,
- increase of lip volume,
- restoration of facial volume,
- treatment of facial loss,
- remodelling of face line,
- treatment of destructive scars,
- treatment of morphological asymmetry.

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Method of Use:

Gel implant shall be injected into uninflamed, disinfected healthy skin. Using the described injection technique has a great importance for obtaining

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a successful result. Use 27G needle, which is supplied with a syringe, and inject slowly using the appropriate technique. Inject small volumes in two or several steps but not the whole volume at once. The injection volume depends upon the type of wrinkles to be corrected. Following the injection, the physician may gently massage the area for uniform distribution of gel implant. Leave the patient in the room for several minutes after the injection so that blanching due to artery embolism could be detected.

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ESTELINE Basic 23 mg/ml hyaluronic acid-based viscoelastic intradermal implant medical device – steam sterilization
0.4 mm x 12.7 mm sterile single use needle (27G) – gas sterilization (Ethylene oxide).

Shelf life:

3 years with observation of transportation and storage conditions.

Disposal:

Medical device waste class of hazard is B (epidemiologically hazardous waste). Used devices, devices with the expired shelf life or damaged package shall be disposed of according to the established procedure in compliance with the sanitary epidemiological requirements for medical waste treatment SanPiN 2.1.7.2790-10.

Collection and processing of medical waste at healthcare facilities shall be controlled by the representatives of the disposal company. For this purpose, special rigid containers with vacuum covers shall be used. All containers shall be painted yellow indicating that the waste is dangerous.

After that, waste of the same type is disposed of at incinerators or disinfected and stored for further dumping at storage premises.

Label Symbols:

No.	Symbol	Description	No.	Symbol	Description
1		Manufacturer	2		Use by
3		Manufacturing date	4		Do not resterilize
5		Do not reuse	6		Caution!
7		Lot number	8		Do not use if package is damaged

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Эл. почта: medmelis@yandex.ru

9		Consult instructions for use	10		Contains no latex
11		Keep away from direct sunlight	12		Keep dry
13		Sterilized using ethylene oxide	14		Sterilized using steam
15		Recommended temperature	16		Recommended humidity
17		Sterile single use needle	18		Sterile single use syringe
19		Authorized representative in the Russian Federation	20		Device registration certificate in Russia
21		CE mark	22		Needle manufacturer

Manufacturer:

CROMA GmbH

Address: Industriezeile 6, 2100, Leobendorf, Austria

Tel.: +43 (0) 2262 684 68-0 / + 43 (0) 2262 684 68-165

Email: office@croma.at

Manufacturing site: Industriezeile 6, 2100, Leobendorf, Austria

Authorized Representative in the Russian Federation:

MELIS, LLC

Address: 1, Pokryshkina str., building 1, 119602 Moscow

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