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Notary Public duly Admitted, Authorised and Sworn
Practising at the Hong Kong Special Administrative Region of the People's Republic of China
DO HEREBY CERTIFY that the document hereunto annexed is the original of the document.

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
- Proximal aortic diameter ≥ 15 mm.
 - Distal iliac artery anchorage zone ≥ 15 mm.
 - Distal iliac artery diameter of 10–16 mm.
 - Morphology suitable for endovascular repair.
- Ankura AUI Stent Graft with delivery system is indicated for endovascular repair of patients with abdominal aortic aneurysms, requiring:
- Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - Non-aneurysmal proximal aortic neck length ≥ 15 mm.
 - Non-aneurysmal proximal aortic neck diameter of 18–32 mm.
 - Proximal aortic neck angulation $\leq 60^\circ$.
 - Distal iliac artery anchorage zone ≥ 15 mm.
 - Distal iliac artery diameter of 10–16 mm.
 - Morphology suitable for endovascular repair.
- Ankura TAA Stent Graft System is indicated for endovascular repair of patients with thoracic aortic aneurysms or thoracic aortic dissection, requiring:
- Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - Non-aneurysmal aortic proximal neck length ≥ 1.5 cm.
 - Non-aneurysmal aortic distal anchorage zone ≥ 1.5 cm.
 - Non-aneurysmal aortic diameter in the range of 18–42 mm.
 - Morphology suitable for endovascular repair.
- Additional considerations for patient selection include:
- Patient's age and life expectancy;
 - Co-morbidities (e.g., cardiac, pulmonary or renal

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N anev...
 or > 12 mm.
 proximal aortic neck an...
 iliac artery anchorage...
 iliac artery diameter < 8 mm or...
 Stent Graft is contraindicated in:
 with acute systemic infection;
 with mesenteric blood flow main...
 accessory renal artery o...
 argic re...

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- ♦ Non-aneurysmal proximal aortic neck length <18mm or >32mm.
 - ♦ Proximal aortic neck angulation > 60°.
 - ♦ Distal iliac artery anchorage zone <15mm.
 - ♦ Distal iliac artery diameter <10mm or >16mm.
- The Ankura AUI Stent Graft is contraindicated in:
- ♦ Patients with acute systemic infection;
 - ♦ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ♦ Patients with accessory renal artery original from abdominal artery;
 - ♦ Patients who have allergic reaction to the device;
 - ♦ Patients who are not suitable for endovascular repair in vascular morphology;
 - ♦ Patients who can not tolerate contrast agents due to renal insufficiency;
 - ♦ Patients who are allergic to contrast agents;
 - ♦ Aneurysms neck with thrombus.
 - ♦ Non-aneurysmal proximal aortic neck length <15mm.
 - ♦ Non-aneurysmal proximal aortic neck diameter <18mm or >32mm.
 - ♦ Proximal aortic neck angulation > 60°.
 - ♦ Distal iliac artery anchorage zone <15mm.
 - ♦ Distal iliac artery diameter <10mm or >16mm.
- The Ankura TAA Stent Graft is contraindicated in:
- ♦ Patients with acute systemic infection.
 - ♦ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ♦ Patients who have allergic reaction to the device;
 - ♦ Patients who are not suitable for endovascular repair in vascular morphology;
 - ♦ Patients who can not tolerate contrast agents due to

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inner diameter.

- ♦ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, check the specification to confirm that it meet the patient' anatomy.

- **Patient preparation**

- ♦ Measure proximal aortic neck length and vessel inner diameter.

- ♦ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.

- ♦ Evaluate the quality of aortic neck.

- ♦ Access the morphology of the vessel compatible with endovascular aortic repair.

- ♦ Undergo any necessary common examination, and monitoring of vital signs.


- ♦ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.

- ♦ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.

- ♦ Adequate iliac/femoral arteries access.


- ♦ Expose both common femoral arteries access using standard surgical technique.

- ♦ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, alternative anticoagulant should be considered.

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


(a) Ankura AAA Stent Graft

(b) Ankura cuff Stent Graft


The stent grafts used in the Ankura Stent Graft System is modular

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
graft components, one Ankura CUFF Stent Graft, and two additional iliac components. Each component is introduced separately into the patient's vascular system. The Ankura Stent Graft is a tube composed of ePTFE film supported by a metal mesh called a stent. The mark points are fixed on the stent. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.

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

The Ankura AUI stent graft is a tube composed of ePTFE film supported by a metal mesh called a **stent**. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.

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
(д) Ankura TAA Stent Graft


The Ankura TAA stent graft is a tube composed of ePTFE film supported by a metal mesh called a **stent**. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.

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

to pass through the
can provide stable and
open the graft.

Ankura Delivery System (AAA, AUI, TAA, CUFF)


(a) Ankura Delivery System (AAA, AUI, TAA Stent Graft)

(b) Ankura Delivery System (CUFF Stent Graft)

1. Tapered Tip
2. Markband
3. Proximal Anchor
4. Stent Graft
5. Inner Tube
6. immobile Pole
7. Sheath
8. Front Grip

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


The Ankura Stent Graft is a tube composed of ePTFE film supported by a metal mesh called a stent. The mark points are fixed on the stent. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery.

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
HLK China. Delivery System - Polytetrafluoroethylene Polyflon F-205, supplier ZETUS, Switzerland; Stainless steel spring - 304 (C ≤0.08%; Si ≤1.00%; Mn ≤2.00%; P ≤0.045%; S ≤0.03%; Cr 18-20%; Ni 8-11%; Fe – balance), supplier Motion, USA; marker band - Pt 90 / Ir 10 MDL-PtIr#1, supplier Medelec, Switzerland; outer layer of sheath - Polyether block amides Arkema Pebax 5533, supplier ZEUS, USA; Hydrophilic coating, supplier Coatings2GO, USA; coating 100% PVP (polyvinylpyrrolidone), which is produced through polymerization reaction of three raw materials: 1) Urethane Coating – Formula B (MSDS # 1006); 2) Primer 5-017 (MSDS # 1030); and 3) Crosslinker A (MSDS # 1007).. Shell handle – Acrylonitrile Butadiene Styrene Sabic HMG94MD, Sabic HMG47MD nature, supplier Xin Hua Jia, China; stainless steel pipe - 304 (C ≤0.08%; Si ≤1.00%; Mn ≤2.00%; P ≤0.045%; S ≤0.03%; Cr 18-20%; Ni 8-11%; Fe – balance), supplier Medelec, Швейцария; Silicon gasket –Dow Corning qp1, supplier Xin Yuan, China; hose - Thermoplastic polyurethane TPU58315, supplier Kist, China. Haemostatic Vale: Polycarbonate (PC) RX1805 (BAYER MAKROLON), supplier Qiao Pai, China; Stopcock (three-way) -

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
Systemic Toxicity (Intraperitoneal injection)	ISO10993-11:2006	There was not acute toxicity for the article	Pass
Systemic Toxicity (Intravenous injection)	ISO10993-11:2006	There were not acute toxicity for the article	Pass
Intracutaneous Reactivity	ISO10993-10:2002/Amd.1:2006	No obvious erythema and oedema was observed on the skin or the rabbits.	Pass
Pyrogen test	ISO10993-11:2006	Febrile response of the test rabbits was not observed.	Pass
Maximization sensitization test (polar extraction)	ISO10993-10:2002/Amd.1:2006	There was not delayed-type hypersensitivity response for the test article	Pass
Maximization sensitization test (non-polar extraction)	ISO10993-10:2002/Amd.1:2006	There was not delayed-type hypersensitivity response for the test article.	Pass
Hemolysis test	ISO10993-4:2002/Amd.1:2006	The hemolytic grade of the test article	Pass

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Implantation Study in Rabbit with Histopathology	ISO10993-6:2007	The biological responses of the test article were similar with the control articles	Pass
Ankura Delivery System			
Cytotoxicity	ISO10993.5:2009	There was no cytotoxicity for the test article	Pass
Maximization sensitization test (polar extraction)	ISO10993-10:2002/Amd.1 :2006	There was not delayed-type hypersensitivity response for the test article.	Pass
Maximization sensitization test (non-polar extraction)	ISO10993-10:2002/Amd.1 :2006	There was not delayed-type hypersensitivity response for the test article.	Pass
Systemic Toxicity (Intraperitoneal injection)	ISO10993-11:2006	There were not acute toxicity for the article.	Pass
Systemic Toxicity (Intravenous injection)	ISO10993-11:2006	There were not acute toxicity for the article.	Pass

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Intracutaneous Reactivity	ISO10993-10:2002/Amd.1 :2006	No obvious erythema and oedema was observed on the skin or the rabbits.	Pass
Hemolysis test	ISO10993-4:2002/Amd.1 :2006	The hemolytic grade of the test article was nonhemolytic	Pass
Pyrogen test	ISO10993-11:2006	Febrile response of the test rabbits was not observed.	Pass
Thrombosis Test	ISO10993-4:2002/Amd.1 :2006	There was no significant difference of thrombosis characteristics between test articles and the control articles.	Pass

Information Search for the Ankura Stent*Graft Carcinogenic Toxicity

The main materials of Ankura Stent Graft are nitinol and e-PTFE (e-Polytetrafluoroethylene) membrane. We used the PubMed searching service as the main way to search the information of the Ankura Stent Graft carcinogenic toxicity. PubMed is a service of the U.S. National Library of Medicine and the National Institutes of Health. The PubMed can be considered complete.

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selected following key words
make the search.

Key words	Limits	
Nitinol/carcinogenic	title/abstract	0
Nitinol/carcinogenicity	title/abstract	50
Nitinol/cancer	title/abstract	196
Nitinol/cancer	-	20
Nitinol/toxicity	-	4
Polytetrafluoroethylene /carcinogenic	title/abstract	3
Polytetrafluoroethylene /carcinogenicity	title/abstract	54
Polytetrafluoroethylene /cancer	title/abstract	497
Polytetrafluoroethylene /cancer	-	119
Polytetrafluoroethylene /toxicity	-	

The present report gives a summary of biocompatibility of Nitinol and e-PTFE.

Nitinol There is a substantial experience accumulated on the use of Nitinol in clinical practice. Nitinol have been shown to be biocompatible through extensive animal and human. Metallic implants represent a pool of trace elements in the body which by corrosion are poured out and may disturb the delicate balance of trace elements in the surrounding tissues as well as possibly of the entire organism, or even cause allergic, toxic, and carcinogenic actions. As for the Ankura Stent Graft, Ni release is responsible for the adverse effects.

The majority of observations based on short and long term implantation agree that nitinol implants do not show any corrosion signs in the body [19-20]. Though occasionally slight corrosion is noticed on the surface of retrieved implants, the corrosion rate falls into the range of "excellent" according to Fonana's classification [19]. Impaired cell behavior occurs at $4 \times 10^3 \sim 9 \times 10^3 \mu\text{g/l}$ Ni concentration. The higher the Ni

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concentration achieves the more severe situation the cell or even the entire organism is. Estimation of the Ni concentrations that could be induced upon exposure to a Nitinol implant in the body ranges from 1 to 100 $\mu\text{g/l}$ according to various authors [21]. That is to say, there is no adverse effect result from the release Ni of Nitinol in the body. Some did the relevant research thus proved that exposure to nickel from Nitinol filter is well within accepted limits. An effective mechanism of their elimination is a factor in metal biocompatibility that is difficult to underestimate. Studies on the distribution of metals that have followed dermal and systemic (through injection) administration in vivo show that Ni is efficiently eliminated from the body through urine and feces.


Another important component of Nitinol is Titanium. Traditionally Ti was considered as a perfectly biocompatible metal because of its native “bioinert” surface oxide film. However, in the last decade, the attitude to Ti and its alloys has been changed. Several reports have warned different toxic effects of Ti on the tissues. So the renewal of investigations on Ti and Ti-base alloys is observed presently. It can not be concluded that Nitinol is absolute non-carcinogenicity, but compared with other alloys such as stainless steel, 316L used in medical practice, Nitinol is superior to them in use and whose biocompatibility is similar to or in certain respect even better, than the biocompatibility of pure titanium.

Take into account the fact that Nitinol is widely used in clinical practice without being reported carcinogenic and the researches suggesting the safety of Nitinol implant in vivo, we think that the risk on carcinogenicity is negligible.

e-PTFE (e-Polytetrafluoroethylene)

The ePTFE membrane, which used in the Ankura Stent Graft, can immediately block the irregular blood flow and prevent the hemolysis. The ePTFE has good biocompatibility and has been used in the clinical practice without being reported carcinogenic.

With the analysis above, we think that the material is safe and have a good biocompatibility. Of course, it is seriously considered for further investigation under review by the relevant authority relating to this topic.

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package which includes several pieces of products. Details are shown in **Table 1**.

Table 1 Package Material

Layer	Material	Unit Amount
Primary Package (Inner)	Two Dialyzing pouches: Tyvek 1073B sheet+PE58 /PET12 film, manufacturer Bao Bo, China; Preformed tray: PETG6763, manufacturer Sen Hua, China	One set
Secondary Package (Outer)	Color box:400g art paper; Scale board: single-wall corrugated board	One set
Transportation Package	Carton: double-wall corrugated board	Several

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
Keep away from sunlight

Non-sterile

LifeTech Scientific (Shenzhen) Co., Ltd.
Floor 1-5, Cybio-Eurochem Building, Longshan Road, North Area of High-Tech Park,
Nanshan District, Shenzhen 518057, P.R.China
Tel: +86-755-86028250 Fax: +86-755-86028251

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

Non-sterile

Non-sterile

Lifetech Scientific (Shenzhen) Co., Ltd.
Floor 15, Ciptech Science Building, Langshan 2nd Street, North Area of High-tech Park,
Nanshan District, Shenzhen 518057, P.R. China
Tel: +86-755-88282550 Fax: +86-755-88282551

Lifetech Scientific (Europe) Corporation U.A.
Kruisstraat 64, 6222 PP Maastricht, the Netherlands

21-04-013 000000

	Label of Ankura TAA Stent Graft
--	---------------------------------

Ankura™

RFF TAA2420B120

LOT 100901

TAA Stent Graft System 2013-08-31

Ankura™

TAA2420B120

TAA Stent Graft System

EN TAA Stent Graft System IT Sistema dell'Innesto dello
ES Sistema de endoprótesis cubierta TAA Stent per TAA
PT Sistema de Enxerto de Stent TAA DE TAA-Stentgraft
FR Système d'endoprothèse TAA TH TAA Stent Graft System
SV TAA stentgraftsystem NL TAA-stentgraftsystem
RU Система стент-графта Ankura TAA

REF TAA2420B120 LOT 100901 SN 001


2010-09-01 2013-08-31

Proximal Dia.(D₁) 24mm Distal Dia.(D₂) 20mm

Bare Stent Length(L₁) 14mm







Covered Length(L₁) 120mm

Sheath Length(L₂) 780mm Sheath Dia. 21F

Consult instructions for use Do not reuse
Sterilized using ethylene oxide No hot water
Non-proprietary Do not use if damaged
Keep away from sunlight
Non-toxic
CE 0344
Lifetech Biomedical (Shenzhen) Co., Ltd.
Floor 1-5, Cyber-Electronics Building, Longtan 2nd Street, North Area of High-Tech Park,
Nanshan District, Shenzhen 518067, P.R.China
Tel: +86-755-85232550 Fax: +86-755-85022251
Lifetech Biomedical (Shenzhen) Co., Ltd.
Lifetech Biomedical (Shenzhen) Co., Ltd.

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design and material characteristics. The sterilization process is validated and subject to routine control.

The sterilization process is carried out according to ISO 11135-2014 and ISO 10993-7

EXECUTIVE SUMMARY

1. At this time, it confirms that the sterility assurance level is

$$SAL = \frac{1}{1.4 \cdot 10^6} = 7.14 \cdot 10^{-7} \leq 10^{-6}.$$

2. The transportation products can continuously, stably produce qualified sterile product through the sterilization in figure 6.

3. After confirming the completion of report approval, it should revise LT/TS/00T-02 Transportation Product EO Sterilization Operation Instruction.

Figure 6.0: sterilization parameter chart of transportation product family

Sterilization temperature/°C	50±3
Heat preservation period/min	60±2
Fore-vacuum/ Kpa	-65±2
Dwell time/min	10±2
Humidity of sterilization	Pressure difference before and after humidification: 2±1Kpa (humidity: 30 %RH ~ 80 %RH)
Injection rate of sterilizing agent	Pressure difference before and after medication: 91±5 Kpa (humidity: 4.8% ~ 8.8%)
Sterilizing time/min	5 30±5 120 hours 33-40

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
Ankura AAA/AUI/TAA Stent Graft System:1490±200 g

Ankura CUFF Stent Graft System:980±200g

⇒ **The size and location of the radiopaque markers parameters and characteristics of the standard:**

The markers, made from platinum-iridium alloy, are embedded between the inner and outer layer of ePTFE film, and will not contact with blood directly. The location is shown in Figure.

A view
 o-shaped mark point
 s-shaped mark point
 Connecting bar
 B view
 v-shaped mark point

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а) **The possibility of oxidation : Challenge account contact corrosion on the relevant parts of the system ;**

- The marker is embedded between the inner and outer layer of the porous PTFE film, and will not contact with blood, thus it has no possibility of corrosion. Moreover, the markers are made from corrosion-resistant material.

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3) The force transmitted by the surrounding tissue;

The vessel wall will exert a force on the device after deployment. Considering this force, an over-sized device (10-20%) is recommended during selecting a device for implantation. This recommendation has been addressed in the instructions for use. Also, the radial outward force and migration testing have been performed.

4) The effort required to deform the implant;

Evaluations are performed for radial outward force, migration and, simulation deployment to prevent the deformation of the implants.

Table 1 The requirement of radial outward force

Diameter of Stent Graft (mm)	Radial Outward Force (N)
≤20	>2
20-30	>3
≥30	>4

Test the migration of the device was performed under the pressure of 180mmHg and 60mmHg respectively. The test results demonstrated that the device did not migrate during testing.

Pushability of Delivery System

The pushing force should be less than 80N.

d) the expected ingrowth into the fabric, tissue perforation, and displacement;

There will be endothelialization of the device after implantation. The test results demonstrated that the device did not migrate during testing. Strength of Covered Film is larger than 25N. The results demonstrated that the device did not migrate during testing.

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Lifetech evaluated the burst strength of the device after implantation and stent-graft migration, which indicated that the device did not migrate during testing. The test result indicates that the Burst Strength of Covered Stent is larger than 25N. Test the migration of the device was performed under the pressure of 180mmHg and 60mmHg respectively. The test results demonstrated that the device did not migrate during testing.

е) administration and delivery systems.

а) the ability of the implant to unfold properly with the necessary precision and safety;

- The correct and safety release can be guaranteed by fine tuning of screwing the delivery system and the structure design of proximal release after: By assist of B- ultrasound or CT, the doctor send the leading end of delivery system (TIP and the leading end of sheath) to the lesion tissue of the patient through the vessel. The leading end of delivery system can be smoothly reach the desired location by the flexible & buckling resistance feature of sheath and the hydrophilic coating on it. The platinum markerband at the leading end of sheath can help the accuracy of stent position for the doctor. Secondly, release slowly by fine tuning of screwing the handle of delivery system: Grip the front handle, contrarotate sliding by right hand, to slow release the bare part and forepart of the stent graft, the action is slow and steady until at least one circle of stent contact closely with the vessel wall, so the safety of release can be guaranteed.

Thirdly, into the process of quick release: Grip the front handle by left hand, steady the delivery system, keep pressing the trigger, relieve the lock of slider, pull back the sliding handle, make the stent graft quick release, until the stent graft completely spread, the process can be guaranteed by the self-expanding feature of stent.

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Fourthly, through the near-end releaser: Grip the front handle by left hand, remove the safety catcher, retreat the near-end releaser, the near-end bare part of stent has released and close contact with the vessel wall at this moment, the releasing process is completed.

The pushability of the delivery system has been evaluated.

b) The ability of the implant is securely fastened in the bloodstream;

- During release, the bleeding valve should be toward the patient, so that the accuracy of stent position can be guaranteed. When the sheath arrive to the lesion tissue, and during the slow release by fine tuning of screwing, we can confirm if the keel of stent is toward the greater curvature of vessel by the assist of imaging equipment. The accurate release of stent can be guaranteed by fine tuning of screwing the handle of delivery system. During the quickly release of stent by tensioning the handle, the close contact between the released part of stent and vessel wall should be check after the retreat of sheath, by fine tuning of screwing the delivery system. The stent will stay at the accurate position after completely correct release.

When the stent is completely released, the TIP should be back to the entrance of sheath by screwing the handle. The smoothly retreat of TIP can be guaranteed by the rigidity and slim feature of the inner sheath core, and will not affect the position of stent.

Evaluations are performed for radial outward force and migration.

The requirement of radial outward force

Diameter of Stent Graft (mm)	Radial Outward Force (N)
≤20	>2
20-30	>3
≥30	>4

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Test the migration of the device at 180mmHg and 60mmHg respectively. The test device did not migrate during testing.

c) the capacity to maintain adequate implant integrity;

- First, the stent assemble into the delivery system by its contractility, and the stent can be successfully carrying with the sheath by the lock of inner sheath and proximal anchor.

Secondly, the innermost layer of the inner sheath, which is contact with the stent directly, the material of it is pebax, pebax has excellent flexibility, the relative motion with PTFE (the material of the film covered on the stent) is very smooth, never feeling block, and never damage the stent. It is free contact between near-end anchor and bare part of the stent, not close contact. The release would be very smooth, and never damage the stent. The rest part of the delivery system does not touch the stent, so it is impossible for them to damage the stent.

Moreover, after release, the retreat process of TIP is visible, it is impossible for the TIP to damage the stent in the correct operation, not even touch the stent again, so the integrity of the stent can be guaranteed. The integrity of the stent can be guaranteed

d) capability of the implant when used properly prevent blood flow through the walls of the implant; Should take into account possible changes in the permeability of the vessel wall after implantation.

- Evaluate the leakage of the ePTFE film.

The permeability of the Ankura Stent Graft are tallyed with the requirements of designed, and all of them are less than 100mL/cm²/min.

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e) the interaction between modules (components) endovascular system;

- Performed connection strength testing on the delivery system. Connect force between sheath and Grip shall be larger than 60N. Connect force between Tip and Inner tubing shall be larger than 15N. Connect force between Tip and Grip shall be larger than 15N. The force of sheath shall be larger than 60N.

f) matching the size and shape of the implant from the viewpoint of compatibility when used in certain vascular diameters;

The stent graft has good flexibility and suitable radial outward force which can make the device complete apposition to the vessel wall.

As addressed in the Instructions for Use, the size of the device should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversized device by 10-20%.

g) the ability of the implant to maintain adequate blood flow through the lumen (cross-country);

The stent graft has good flexibility and suitable radial outward force which can make the device complete apposition to the vessel wall.

As addressed in the Instructions for Use, the size of the device should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversized device by 10-20%.

h) the compatibility of the implant with the patient's anatomy;

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in magnetic resonance imagin

Evaluate the compatibility of the implants which indicate design specification and suggest no to check for MRI after accordingly.

The angle is larger than 5° and the temperature is higher than 1.5°C. This isn't safe to human for MRI. So we suggest the patient not to check for MRI after implantation

i) relevance of the requirements of ISO 10993-1 and other parts of ISO 10993;

- Yes, comply with the relevant ISO 10993 series standard.

j) distinguishability implant fluoroscopy or other methods.

- During the operation, doctor judges the direction and position of the implant by markers under X-ray observation, to facilitate accurate positioning of the implants into the vessel, and then inject contrast agent to confirm the implant shape and operation effect.

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⇒ The delivery system - joints strength - bending strength - strength of the connection torsional - The force required to disclose stent - Breaking Local compressive resistance to the movement of the radial recoil force tensile strength of internal / external Corrosion strength anastomotic durability Evaluation performance, and release shall be larger than 60N. ing shall be larger than larger than 1.5N. 's no crack

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- For delivery system, tests are performed between the components, torsion performance, Connect force between sheath and Grip shall be larger than 15N. Connect force between Tip and Inner tubing shall be larger than 15N. Connect force between Tip and Grip shall be larger than 15N. The force of sheath shall be larger than 60N.

- Fatigue strength.

After 380 million fatigue cycles of each sample, there is no crack or rupture observed.

⇒ **Strength in multiple puncture (for arteriovenous shunts during intravascular access) The strength of the stent / fastening systems when mounting graft (eg, worn on the adhesive mounting using sutures) Permeability The bending / twisting resistance migration test gap for modular components free surface area of the stent**

- The stent is intended to be implanted into aorta, not involving the arteriovenous shunt operation.

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		<table><tr><td>24 B 06 0</td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td>T A A 46 38 B 12 0</td><td>5. 3 9</td><td>5. 0 8</td><td>5. 0 2</td><td>5. 1 0</td><td>5. 2 3</td><td>5. 1± 0. 2</td></tr><tr><td>T A A 30 26 B 06 0</td><td>5. 3 6</td><td>5. 4 8</td><td>5. 9 6</td><td>5. 6 8</td><td>5. 7 7</td><td>5. 6± 0. 2</td></tr></table>	24 B 06 0							T A A 46 38 B 12 0	5. 3 9	5. 0 8	5. 0 2	5. 1 0	5. 2 3	5. 1± 0. 2	T A A 30 26 B 06 0	5. 3 6	5. 4 8	5. 9 6	5. 6 8	5. 7 7	5. 6± 0. 2			
24 B 06 0																										
T A A 46 38 B 12 0	5. 3 9	5. 0 8	5. 0 2	5. 1 0	5. 2 3	5. 1± 0. 2																				
T A A 30 26 B 06 0	5. 3 6	5. 4 8	5. 9 6	5. 6 8	5. 7 7	5. 6± 0. 2																				
Corrosion Test of Fatigue Tested Stent Grafts	Measure the polarized curve with the PS-268A electrochemical measurement system, and calculate the corrosion current density with polarized curve, then measure the Ni2+ concentration by ultraviolet spectrograph. The critical pitting potential must be larger than 600mV.	The occurrence of corrosion current at high voltage indicated that only very high voltage (higher than 600mV) could produce corrosion spots on the surface of Ankura Stent Graft. Ankura Stent Graft met the requirements, and all of critical pitting potential is higher than 600mV.																								
Permeability of Covered Film Testing Report	Select one stent, Seal one end of the stent, and through it's other end inject distilled water into the stent by using a force pump, until the liquid pressure is 16.0±0.3KPa, and then record the volume of the permeability with 10min. The leakage should be less than 100mL/cm2/min.	The permeability of the Ankura Stent Graft are tallyed with the requirements of designed, and all of them are less than 100mL/cm2/min. Results: <table><tr><th>Sample</th><th>Area/cm²</th><th>Volume of the permeability N.O.</th></tr><tr><td>TAA2020 B040</td><td>25.12</td><td>N.O.</td></tr><tr><td>TAA4646 B200</td><td>288.88</td><td>N.O.</td></tr><tr><td>AUI20120 120</td><td>120.58</td><td>N.O.</td></tr><tr><td>AAA2010 B100</td><td>120.58</td><td>N.O.</td></tr><tr><td>AAA3414 B160</td><td>12.56</td><td>N.O.</td></tr><tr><td>cuff1010C 040</td><td>71.59</td><td>N.O.</td></tr><tr><td>cuff1622C 120</td><td></td><td>N.O.</td></tr></table> Note: N.O.= No observed	Sample	Area/cm²	Volume of the permeability N.O.	TAA2020 B040	25.12	N.O.	TAA4646 B200	288.88	N.O.	AUI20120 120	120.58	N.O.	AAA2010 B100	120.58	N.O.	AAA3414 B160	12.56	N.O.	cuff1010C 040	71.59	N.O.	cuff1622C 120		N.O.
Sample	Area/cm²	Volume of the permeability N.O.																								
TAA2020 B040	25.12	N.O.																								
TAA4646 B200	288.88	N.O.																								
AUI20120 120	120.58	N.O.																								
AAA2010 B100	120.58	N.O.																								
AAA3414 B160	12.56	N.O.																								
cuff1010C 040	71.59	N.O.																								
cuff1622C 120		N.O.																								

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Stent Grafts Connection Strength Testing Report	Stent fracture strength of each connection point shall be larger than 5N. Connect force between the film and Stent shall be larger than 7N.	4	Y	52,7	Good	060915
		5	Y	48,6	Good	060915
		6	Y	49,2	Good	060915
		7	Y	53,4	Good	060915
		8	Y	56,2	Good	060915
		9	Y	57,8	Good	060915
		10	Y	56,2	Good	060915
		X± S		53,4 ±3,2		
		Note: Y says Ankura Stent Graft System can be easily delivery through the blood vessel model.				
		The test result indicates that the Stent Grafts Connection Strength meets the product technical requirements. Those forces meet the requirement of criterion. Results: The mean value value the maximum burst force of Connect force between the film and Stent is 20.9±1.9 N The mean value the maximum burst force of Stent fracture strength is 40.9±2.0 N				
		The results show that the Ankura Delivery System doesn't have blood leakage and bubble aspiration during operation. Durind the tests were no leakage and aspiration observed. Delivery System observed requirements were met.				
		The results show that the Ankura Delivery System has good resistance. Delivery System's				

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...the date of ... use the expired ...
 ... instructions carefully. Failure to ...
 ... consequences or injury to the patient.
 ... Graft can not be used after ...
 ... is observed prior to use, do ...
 ... return it to Lifetech.
 ... only be used by ...
 ... r interventic

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a minimum, annual X-ray, contrast or non-contrast CT imaging is required.

- ◆ Patients with specific clinical findings, e.g, enlarging aneurysms or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms size and/or persistent endoleak may lead to aneurysms rupture.

- ◆ Patients with infection after the device placement should be treated with oral antibiotics.

- ◆ A vascular surgery team should always be available during the implantation, or re-intervention procedure in which situation that open surgical repair is needed.

- ◆ Do not re-use or re-sterilize the Ankura Stent Graft System. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or cause device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also cause a risk of contamination of the device and/or patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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- ◆ The safety and effectiveness of the Ankura AAA Stent Graft has not yet been evaluated in following populations:

- ◆ Traumatic aortic injury
- ◆ Uncorrectable coagulopathy
- ◆ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos syndrome)
- ◆ Concomitant thoracic aortic aneurysms
- ◆ Patients with active aneurysms
- ◆ Pregnant or nursing females
- ◆ Child obesity
- ◆ Less than 18 years of age
- ◆ Life expectancy less than one year
- ◆ Diameter (measured inner wall to inner wall) (minimal tortuosity, occlusive anatomy) should be compatible with the profile of the device
- ◆ Previous vasculature malformations
- ◆ Previous deployment

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- ♦ Patients with
- ♦ Pregnant or nursing
- ♦ Morbid obesity
- ♦ Less than 18 years of age
- ♦ Patient's life expectancy less than one year
- ♦ Access vessel diameter (measured inner wall to outer wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and the profile of the delivery system. Too tortuous vasculature may increase the risk of failure in stent graft deployment or proximal end release.
- ♦ The Ankura Stent Graft System is not recommended in patients who weigh greater than 350 Pounds (150kg) or can not undergo accurate fluoroscopy examination due to obesity.
- ♦ Patients with a systemic infection may be at increased risk of endovascular stent-graft infection.
- ♦ Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications.
- Implant Procedure
 - ♦ Systemic anticoagulation should be used during the procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
 - ♦ Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
 - ♦ Maintain the guide wire position during delivery system insertion.
 - ♦ Do not bend or kink the delivery system. Doing so may

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cause damage to the delivery system and the stent graft.

- ◆ Always use fluoroscopy for guidance, delivery and observation of the device within the vasculature.
- ◆ The use of the Ankura Stent Graft System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast agent used during the procedure.
- ◆ To avoid any twist in the endovascular stent-graft, during any rotation of the delivery system, rotate all of the components of the system together carefully.
- ◆ Inaccurate placement of the stent graft within the vessel may result in increased risk of endoleak, migration or embolization.
- ◆ Inadequate fixation of the stent graft may result in increased risk of migration of the stent graft, and incorrect deployment and migration of the stent graft requires surgical intervention.
- ◆ Do not continue advancing if resistance is felt during advancement of the guide wire or delivery system. Stop and assess the cause of resistance.
- ◆ Unless medically indicated, do not deploy the stent graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities.

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- ♦ Attach the syringe with heparinized saline to the luer hub and hemostasis valve and flush separately to evacuate air from the inner tube and the outer sheath. (Fig.3)

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- ♦ Monitor important vital signs.
- ♦ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution. Sheath and guidewire should be flushed before every insertion.
- ♦ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)

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
- ◆ Expose bilateral femoral arteries by standard surgery.
- ◆ Insert one 0.035"/260cm super stiff guidewire in the femoral artery, and advance it until above the renal arteries.
- ◆ Insert marker sheath, and double check the stent specification by radiography.
- Deploy the stent graft
 - ◆ Advance the delivery system to the abdominal aneurysm along the guidewire, and always maintain the hemostasis valve towards the other side of the patient.(Fig. 5)

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- ◆ Inject contrast agents into the abdominal aorta; mark the position of the target location, either on the imaging screen or with radiopaque marker or angiographic ruler on the patient's body. Adjust the position of the Stent Graft through the proximal radiopaque markers so that the top edge of the graft is just below the lowest renal artery.
- ◆ Ensure that the proximal radiopaque

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- ♦ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 7)



Fig. 7 Remove the Secure Buckle

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**Fig. 11 Deploy the cuff Stent Graft in the short branch of
AAA Stent Graft**

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- Insert another cuff stent graft system over the long branch of the AAA Stent Graft.

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radiopaque marker on the
Graft, with at least two waves over

(Fig. 13)


Fig. 13 Introduction of another cuff Stent Graft

- Deploy the cuff stent graft under monitoring, and withdraw the delivery system and guidewire. (Fig. 14)

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The Ankura

- **Products preparation**

- ◆ The Stent Graft diameter should be selected to vessel inner diameter not outer diameter. An uncorrectable size of the device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times larger than the landing zone vessel inner diameter.
- ◆ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, re-review the specification conforming to the patient required.

- **Patient preparation**

- ◆ Measure proximal aortic neck length and vessel inner diameter.
- ◆ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
- ◆ Evaluate the quality of aortic neck.
- ◆ Assess the morphology of the vessel compatible with endovascular aortic repair.
- ◆ Undergo any necessary common examination, and monitoring of vital signs.

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
- ◆ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs
- ◆ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
- ◆ Adequate iliac/femoral arteries access and expose both common femoral arteries assess using standard surgical technique.
- ◆ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

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- ◆ Monitor important vital signs.
- ◆ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution, sheath and guidewire should be flushed after each exchange.
- ◆ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)

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- ◆ Expose bilateral femoral arteries by standard surgery.
- ◆ Insert one 0.035"/260cm super stiff guidewire in the femoral artery, and advance it above the ranal arteries.

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- ◆ Insert marker sheath, and determine again the stent specification using radiography.

- Deploy the stent graft
 - ◆ Advance the delivery system to the abdominal aneurysm along the guidewire, and always maintain haemostasis valve in opposite side of the patient. (Fig. 5)

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Note, the connecting bar of the AUI Stent Graft must be in the opposite side of the patient. And the connecting bar has the same side with the hemostatic valve of the delivery system when the AUI Stent Graft is preloaded into the delivery system.

- ◆ Inject contrast agents into the abdominal aorta, mark the position of the target location, either with a radiographic ruler on the stent graft such that the stent graft is just below the lowest renal artery, stationary, then rotate the slider slowly deploy the stent graft to the vascular wall. Hold the front of the stent graft.

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
Fig.6 Withdraw outer sheath

- ♦ Slowly deploy the stent graft. Hold the delivery system stationary, then rotate the grip counterclockwise to slowly deploy the stent graft until the covered stent adhere to the vascular wall.
- ♦ Rapidly deploy the stent graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 6)

- Proximal end release
 - ♦ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system.

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(Fig. 7)


Fig. 7 Remove the Secure Buckle

- ♦ Pull the Proximal Releaser to release the proximal spring of the stent graft. (Fig. 8)


Fig. 8 Pull the Proximal Releaser

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- ◆ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
- ◆ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
- ◆ Adequate iliac/femoral artery access and expose one common femoral artery using standard surgical technique.
- ◆ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

Operation Procedure

- Preparation before implant
 - ◆ Attach the syringe with heparinized saline to the luer hub and hemostatic valve and flush separately to evacuate air from the inner tube and the outer sheath.(Fig.3)

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- ♦ Expose one side femoral arteries by standard surgery.
- ♦ Place a 0.035"/260cm super stiff guidewire in the femoral artery and advance it until to the ascending artery. (Fig. 5)

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- ◆ Ensure the Stent Graft positioning at the aneurysms or dissection sites. Ensure the connecting bar of the TAA Stent Graft in the outer curved side of aortic arch.
- ◆ Repeat the angiography to verify the delivery system locating at the aneurysms or dissection sites.
- ◆ Decrease the systolic blood pressure (SBP) to 10.5 to 12kpa.
- ◆ Slowly deploy the stent graft. Hold the front grip of the

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delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the proximal part of the stent graft until the covered stent adhere to the vascular wall. (Fig. 7)


Fig. 7 Slowly deploy the stent graft

- ♦ Rapidly deploy the stent graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 8)

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

Fig. 8 Rapidly deploy the stent graft

- Proximal end release
 - ♦ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 9)



Fig. 9 Remove the Secure Buckle

- ♦ Pull the Proximal Releaser to release the proximal end

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

process
the main step
with our quali
fertilization, and the

Manufacturing processes of the Ankur
of Lifetech Scientific (Shenzhen)
raw materials were purchased exactly
Cleaning should be performed carefully i
Steps from Cleaning to QC-07 before
performed in cleaning room.

A close-up photograph of a red wax seal on a white envelope. The seal is a thick, irregular, red wax-like substance that has been pressed onto the envelope, creating a textured, slightly raised surface. The envelope is white and has some faint, dark markings on it, including a small cross-like symbol. The background is a plain, light-colored surface.

-
- A close-up photograph of a red wax seal on a white envelope. The seal is a thick, irregular, reddish-brown shape, likely made of wax, used to secure the envelope. It is positioned over the top edge of the envelope, which has a small black mark or hole visible above it.

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Manufacturing process


The manufacturing process shows the main steps of manufacturing of the Ankura Stent Graft with Delivery System in accordance with our quality management system. The special procedures are last and sterilization, and the key procedures include stent shapsetting and stent

The manufacturing processes of the Ankura Stent Graft with Delivery System were performed in the factory of Lifetech Scientific (Shenzhen) Co., Ltd.

All raw materials were purchased exactly according to requirements and drawings.

Cleaning should be performed carefully in water for injection or alcohol solution.

- Steps from Cleaning to QC-07 before EO sterilization (the area inside the dashed lines) were performed in cleaning room.

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TO ALL TO WHOM THESE PRESENTS SHALL COME : I, LAU TAI CHIM

Notary Public duly Admitted, Authorised and Sworn

practising at the Hong Kong Special Administrative Region of the People's Republic of China

DO HEREBY CERTIFY that the documents hereunto annexed are the originals of the

documents.

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued.

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Reference Code 參考編號 9AFF9DC3

IN TESTIMONY whereof I have
hereunto subscribed my name and affixed
my Seal of Office this 13th day of
December in the year of Our Lord Two
thousand and seventeen.

LAU TAI CHIM
NOTARY PUBLIC,
HONG KONG SAR.

Ankura AAA and Cuff Stent Graft Systems Instructions for Use

Regulatory Affairs Director
Lifetech Scientific (Shenzhen) Co., Ltd

Lily Shi Name

(sign and seal)

Dec.12,2017

Lifetech Scientific (Shenzhen) Co., Ltd

1 / 27

Contents

1	Device Description	3
2	Indications for Use.....	8
3	Contraindications.....	8
4	Pre-implantation Preparation	9
5	Implantation Procedure.....	11
6	Potential Adverse Events	17
7	Storage and Transport Requirements	18
8	Warnings	18
9	Follow Up	22
10	Shelf Life	22
11	Package and Label.....	23
12	Warranty	23
13	Operating conditions	23
14	Utilization.....	23
15	Authorized Representative of Manufacturer in the territory of Russia.....	24
	Company name, address, telephone, website, e-mailThe symbols used are as the following:	24
	The list of national and international regulatory documents/standards, which matches medical device	25

Read all instructions carefully. Failure to follow the instructions, warnings and precautions may lead to serious complications or injury to the patients.

1 Device Description

● Ankura Stent Graft System

✧ The Ankura AAA Stent Graft System includes:


- The Ankura AAA Stent Grafts
- The Ankura Delivery System

✧ The Ankura cuff Stent Graft System includes:


- The Ankura cuff Stent Grafts
- The Ankura Delivery System

✧ The Ankura Stent Graft is pre-loaded into the Delivery System. The loaded delivery system is inserted endoluminally via the femoral or iliac artery and tracked through the patient's vasculature to deliver the stent graft to the target site.

● Ankura Stent Grafts

(a) Ankura AAA Stent Graft


(b) Ankura cuff Stent Graft


Fig. 1 Ankura AAA Stent Graft and cuff Stent Graft

- ✧ The stent grafts used in the Ankura AAA Stent Graft System are modular and comprises of three stent graft components, one Ankura AAA Stent Graft and two Ankura cuff Stent Grafts. The modular system consists of one AAA Stent Graft and two cuff Stent Grafts (Fig.1). Each component is introduced separately into the patient's vascular system. The Ankura Stent Graft is a tube composed of ePTFE membrane supported by a metal wireframe called a stent. The markers are fixed on the stent. The ePTFE membrane is sintered on to the wireframe without suture, it isolates the aneurysm from blood pressure. The stent is expandable and provides radial force to keep the graft open.

● Ankura Delivery System

(a) Ankura Delivery System (AAA Stent Graft)

(b) Ankura Delivery System (cuff Stent Graft)

1	Tapered Tip	2	Marker band	3	Proximal Anchor
4	Stent Graft	5	Inner Tube	6	immobile Pole
7	Sheath	8	Front Grip	9	Slider Grip
10	Trigger	11	Screw Gear	12	hemostasis valve
13	Rear Grip	14	Proximal Releaser	15	Secure Buckle
16	Luer Connector				

Fig.2 Ankura Delivery Systems

✧ The Stent Graft is pre-loaded inside the Delivery System. The Delivery System facilitates the placement of the stent graft via the arterial vasculature (e.g., femoral arteries). Under the fluoroscopic guidance, the Ankura Delivery System is properly positioned within the patient's vasculature and the stent graft is deployed from the delivery system. The AAA Stent Graft with Delivery System is shown in **Figure 2(a)**, and the cuff Stent Graft with Delivery System is shown in **Figure 2(b)**.

● Device specifications


✧ The specification of Ankura AAA and cuff Stent Graft System are indicated as follow and the details are shown respectively in table 1 and Table 2.

AAA 20 10 B 100

_____ Covered Length (mm)
 _____ Distal Diameter (mm)
 _____ Proximal Diameter (mm)
 _____ Stent Graft Type

Proximal Diameter (D0/mm)	10	12
---------------------------	----	----

(a) AAA Stent Graft System

(b) cuff Stent Graft System

Table 1 Summary of Ankura AAA Stent Graft System

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Covered Length (L2/mm)	Covered Length (L3/mm)	Delivery System
20	10	100,110,120,130,140,160	50	30	21F
22		100,110,120,130,140,160			21F
24		100,110,120,130,140,160			21F
26	12	100,110,120,130,140,160			21F
28		100,110,120,130,140,160			23F
30		100,110,120,130,140,160			23F
32		100,110,120,130,140,160			23F
34		100,110,120,130,140,160			23F
36		100,110,120,130,140,160			23F

The materials from Ankura AAA Stent Graft System

Stent-Graft – Nitinol; Covered film – ePTFE LQ-ePTFE; Steel tube –SUS316L; Platinum marker – Platinum 0213; Suture – Polypropylene. Delivery System: The sheath: Sheath - Polytetrafluoroethylene (PTFE) liner Polyflon F-205; Stainless steel spring –304; marker band –Pt 90/ Ir 10; outer layer of sheath – Polyether block amides Arkema Pebax; Hydrophilic coating Coatings2GO. Shell handle – Acrylonitrile Butadiene Styrene Sabc; stainless steel pipe - 304; Silicon gasket –Dow Corning; hose - Thermoplastic polyurethane TPU58315. Haemostatic Vale: Polycarbonate (PC) RX1805; Stopcock (three-way) - Acrylonitrile Butadiene Styrene 757K; handle of Stopcock (three-way) – polyethylene PE2911. Inner tubing: Sleeve – stainless steel 304; Tapered Tip – Polyamide Pebax + BaSO₄ + TiO₂; Proximal anchor – stainless steel 304; PI tubing – Polyimides. Immobile pole: Immobile pole – polyethylene + BaSO₄; metal tube - stainless steel 304; Locked connector - Acrylonitrile Butadiene Styrene Sabc; Gripping handle - stainless steel 304.

Table 2 Summary of Ankura cuff Stent Graft System

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
10	10	40,60,80	18F
12		40,60,80,90,100, 120	18F
14		40,60,80, 90,100, 120	18F
12	12	40,60,80, 90,100, 120	18F
14		40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F
12	14	40,60,80, 90,100, 120	18F
14		40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F
12	16	40,60,80, 90,100, 120	18F
14		40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F
12	18	40,60,80, 90,100, 120	18F
14		40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F
18		40,60,80, 90,100, 120	18F
14	20	40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F
14	22	40,60,80, 90,100, 120	18F
16		40,60,80, 90,100, 120	18F

The materials from Ankura cuff Stent Graft System

Stent-Graft – Nitinol; Covered film – ePTFE LQ-ePTFE; Steel tube –SUS316L; Platinum marker – Platinum 0213; Suture – Polypropylene. **Delivery System:** The sheath: Sheath - Polytetrafluoroethylene (PTFE) liner Polyflon F-205; Stainless steel spring –304; marker band –Pt 90 / Ir 10; outer layer of sheath – Polyether block amides Arkema Pebax; Hydrophilic coating Coatings2GO. Front grip – Polyoxymethylene; Silicon gasket – Dow Corning; Hose - Thermoplastic polyurethane (TPU) TPU58315. Haemostatic valve : Polycarbonate (PC) RX1805; Stopcock (three-way) - Acrylonitrile Butadiene Styrene AEC 757K; handle of Stopcock (three-way) – polyethylene PE2911. Inner tubing : tubing – stainless steel 304; Tapered tip – polyamide Pebax + BaSO₄ + TiO₂; Luer connector – Polycarbonate (PC). Immobile pole: Immobile pole – Polyethylene + BaSO₄; Metal tube - stainless steel 304; Rear grip –

Polyoxymethylene (POM); Gripping handle - stainless steel 304; Fix
Polyoxymethylene (POM); Immobile pole connector - stainless steel 304; Silo
gasket - Dow Corning.

Physicians should strictly adhere to indications and contraindications.

2 Indications for Use

- Ankura AAA Stent Graft System is indicated for endovascular repair of patients with abdominal aortic aneurysms, requiring:
 - ✧ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - ✧ Non-aneurysmal proximal aortic neck length $\geq 15\text{mm}$.
 - ✧ Non-aneurysmal proximal aortic neck diameter of 18–32mm.
 - ✧ Proximal aortic neck angulation $\leq 60^\circ$.
 - ✧ Distal iliac artery anchorage zone $\geq 15\text{mm}$.
 - ✧ Distal iliac artery diameter of 8–22mm.
 - ✧ Morphology suitable for endovascular repair.
- Additional considerations for patient selection include:
 - ✧ Patient's age and life expectancy;
 - ✧ Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity);
 - ✧ Toleration to the general, regional or local anesthesia;

The final decision is at the discretion of the physician and patient.

● Device Selection

- ✧ The size of the Ankura Stent Graft should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversized Ankura Stent Graft by 10%-20%.

3 Contraindications

- The Ankura AAA Stent Graft is contraindicated in:
 - ✧ Patients with acute systemic infection;
 - ✧ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ✧ Patients with accessory renal artery original from abdominal artery;

- ✧ Patients who have allergic reaction to the device;
- ✧ Patients who are not suitable for endovascular repair in vascular morphology;
- ✧ Patients who can not tolerate contrast agents due to renal insufficiency;
- ✧ Patients who are allergic to contrast agents;
- ✧ Aneurysms neck with thrombus.
- ✧ Non-aneurysmal proximal aortic neck length <15mm.
- ✧ Non-aneurysmal proximal aortic neck diameter <18mm or >32mm.
- ✧ Proximal aortic neck angulation > 60°.
- ✧ Distal iliac artery anchorage zone <15mm.
- ✧ Distal iliac artery diameter <8mm or >22mm.

4 Pre-implantation Preparation

Read the instructions carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician rather than replacing the physician judgment.

- Materials required
 - ✧ 0.035"/260cm Super Stiff Guidewire. For example, Surpass Super Stiff Guidewire, Amplatz Super Stiff Guidewire or comparable.
 - ✧ Multi-plane fluoroscopic imaging system;
 - ✧ Angiographic sheath and marker sheath;
 - ✧ Contrast agent;
 - ✧ High pressure syringe;
 - ✧ Heparinized saline solution;
 - ✧ Any other devices required by the physician for the procedure (including devices used for femoral artery dissection).
- A multidisciplinary team that has combined procedural experience with:
 - ✧ Local anesthesia or general anesthesia;
 - ✧ Femoral artery exposure, arteriotomy and repair;
 - ✧ Percutaneous access and closure techniques;

- ✧ Fluoroscopic and angiographic image interpretation;
 - ✧ Endovascular stent placement;
 - ✧ Appropriate use of radiographic contrast agent;
 - ✧ Techniques to minimize radiation exposure;
 - ✧ Expertise in necessary patient follow-up modalities;
 - ✧ Experience in traditional vascular surgery for any possibility of failure in EVAR.
- Products preparation
 - ✧ The Stent Graft diameter should be selected by the vessel inner diameter ,and not outer diameter. Selecting an improper size of device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times larger than the landing zone vessel inner diameter.
 - ✧ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, check the specification to confirm that it meet the patient' anatomy.
 - Patient preparation
 - ✧ Measure proximal aortic neck length and vessel inner diameter.
 - ✧ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
 - ✧ Evaluate the quality of aortic neck.
 - ✧ Access the morphology of the vessel compatible with endovascular aortic repair.
 - ✧ Undergo any necessary common examination, and monitoring of vital signs.
 - ✧ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
 - ✧ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
 - ✧ Adequate iliac/femoral arteries access.
 - ✧ Expose both common femoral arteries access using standard surgical technique.
 - ✧ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative

anticoagulant should be considered.

5 Implantation Procedure

- Preparation before implantation

- ✧ Attach the syringe with heparinized saline to the luer hub and hemostasis valve and flush separately to evacuate air from the inner tube and the outer sheath. (Fig.3)


Fig 3 Flush the inner tube and the outer sheath

- ✧ Monitor important vital signs.
- ✧ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution. Sheath and guidewire should be flushed before every insertion.
- ✧ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)


Fig. 4 Activate the hydrophilic coating

- ✧ Expose bilateral femoral arteries by standard surgery.
 - ✧ Insert one 0.035"/260cm super stiff guidewire in the femoral artery, and advance it until above the renal arteries.
 - ✧ Insert marker sheath, and double check the stent specification by radiography.
- Deploy the stent graft

- ✧ Advance the delivery system to the abdominal aneurysm along the guidewire, always maintain the hemostasis valve towards the other side of the patient. (Fig. 5)


Fig. 5 Introduction of the AAA Stent Graft with Delivery System

Note: The Short branch of the AAA Stent Graft must be in the same side with the hemostasis valve of the delivery system when the AAA Stent Graft is preloaded into the delivery system.

- ✧ Inject contrast agents into the abdominal aorta; mark the position of the target location, either on the imaging screen or with radiopaque marker or angiographic ruler on the patient's body. Adjust the position of the Stent Graft through the proximal radiopaque markers so that the top edge of the graft is just below the lowest renal artery.
- ✧ Ensure that the proximal radiopaque "O-shaped" marker is on the same side with the "V-shaped" marker on the long branch of the AAA stent graft: Ensure that the proximal radiopaque "8-shaped" marker is on the same side with the "e-shaped" marker, the "O-shaped" marker and the "✓-shaped" marker on the short branch of the AAA stent graft.
- ✧ Ensure that the distal portion of the short branch is above the aortic bifurcation and within the aneurismal sac.
- ✧ Slowly deploy the stent graft. Hold the front grip of the delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the stent graft until the short branch is deployed. (Fig. 6)


Fig. 6 Deploy the short branch of AAA Stent Graft

- Proximal end release
 - ✧ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. **(Fig. 7)**



Fig. 7 Remove the Secure Buckle

- ✧ Pull the Proximal Releaser to release the proximal spring of the stent graft. **(Fig. 8)**


Fig. 8 Pull the Proximal Releaser

- Insert another guidewire through the other side of the femoral artery into the short branch of placed AAA Stent Graft. (Fig. 9)


Fig. 9 Introduction of another guidewire

- Insert the cuff stent graft with delivery system over guidewire into the short branch of the AAA stent graft. Align the radiopaque marker on the cuff stent with the radiopaque marker on the short branch of the AAA Stent Graft. Ensure that the radiopaque “O-shaped” marker on the cuff stent coincide with the “e-shaped” marker on the short branch of the AAA Stent Graft (Fig. 10)


Fig. 10 Introduction of the cuff Stent Graft with Delivery System

- Withdraw outer sheath to deploy the extension branch cuff under monitoring. After the stent is deployed, withdraw the delivery system and guidewire. (Fig. 11)


Fig. 11 Deploy the cuff Stent Graft in the short branch of AAA Stent Graft

- Rapidly deploy the long branch of the AAA Stent Graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 12)
- Hold the slider grip of the delivery system stationary, unlock the trigger and pull back the front grip to make the tip back to the distal end of the out sheath. Then remove the delivery system.


Fig. 12 Deploy the long branch of the AAA Stent Graft

- Insert another cuff stent graft system over guidewire into the long branch of the AAA Stent Graft. Align the radiopaque marker on the cuff Stent Graft with the radiopaque marker on the long branch of the AAA Stent Graft, with at least two waves overlapped between them.

(Fig. 13)


Fig. 13 Introduction of another cuff Stent Graft

- Deploy the cuff stent graft under monitoring, and withdraw the delivery system and guidewire. **(Fig. 14)**


Fig. 14 Deploy the cuff Stent Graft in the long branch of AAA Stent Graft

6 Potential Adverse Events

Adverse events associated with use of the Ankura AAA Stent Graft System include, but are not limited to the following:

- Aneurysms enlargement
- Aneurysms/vessel rupture
- Arterial or venous thrombosis (including thrombosis and embolization)
- Arterial damage
- Arterial perforation
- Arteriovenous fistula or pseudoaneurysm
- Bleeding/Hematoma
- Bowel complications
- Cardiac failure/Infarction
- Conversion to open surgery
- Death
- Edema
- Erectile dysfunction
- Infection
- Neurological damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- Pulmonary complications
- Renal function insufficiency/kidney failure
- Wound complications and subsequent attendant problems

- Device complications due to:
 - ✧ Insertion and removal difficulties
 - ✧ Release difficulties
 - ✧ Failure to deploy
 - ✧ Inaccurate placement
 - ✧ Endoleak
 - ✧ Stent migration
 - ✧ Stent fracture
 - ✧ Rupture of graft material
 - ✧ Erosion
 - ✧ Stent thrombosis or embolization

7 Storage and Transport Requirements

- The Ankura AAA Stent Graft and cuff Stent Graft are pre-loaded in delivery system, and stored in sterile package.
- Store in a well-ventilated and clean space with non-corrosive gas..
- During transportation, no extrusion, no wet exposure; handle with care.
- No exposure to organic solvent.
- Recommendations on storage: temperature from 15° to 30° C, humidity 30-85%, atmospheric pressure is not controlled.
- Recommendations on transportation: temperature from -30° to 55° C, humidity 30-85%, atmospheric pressure is not controlled.

8 Warnings

- General
 - ✧ Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
 - ✧ The Ankura AAA Stent Graft and cuff Stent Graft can not be used after resterilization.
 - ✧ If any sign of damage is observed prior to use, do not use the product and return it to Lifetech.
 - ✧ This device should only be used by physicians and teams trained in vascular

interventional techniques and the use of this device.

- ✧ The long term performance of endovascular grafts has not yet been confirmed. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft (e.g., endoleak, migration, distortion, aneurysms enlarging etc).
- ✧ Clinical studies show that incorrect patient selection may result in poor treatment effect. In addition, unsuccessful endovascular aortic repair which requires open surgery may result in an increased risk of complications related to both procedures.
- ✧ Patient's condition needs to be evaluated before implantation on the following aspects, cardiac failure/infarction, pulmonary complications, edema, and renal function insufficiency.
- ✧ After the device placement, patients should be regularly monitored for the state of the device, including endoleak, migration, change of the Stent Graft structure or position, aneurysms enlarging etc. At a minimum, annual X-ray, contrast or non-contrast CT imaging is required.
- ✧ Patients with specific clinical findings, e.g, enlarging aneurysms or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms size and/or persistent endoleak may lead to aneurysms rupture.
- ✧ Patients with infection after the device placement should be treated with oral antibiotics.
- ✧ A vascular surgery team should always be available during the implantation, or re-intervention procedure in which situation that open surgical repair is needed.
- ✧ Do not re-use or re-sterilize the Ankura AAA Stent Graft System and the Ankura cuff Stent Graft System. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or cause device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also cause a risk of contamination of the device and/or patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

- Patients selection, treatment and follow-up

- ✧ The safety and effectiveness of the Ankura AAA Stent Graft has not yet been evaluated in following populations:

- ◆ Traumatic aortic injury
- ◆ Uncorrectable coagulopathy
- ◆ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
- ◆ Concomitant thoracic aortic or thoraco abdominal aneurysms
- ◆ Patients with active systemic infections
- ◆ Pregnant or nursing females
- ◆ Morbid obesity
- ◆ Less than 18 years of age
- ◆ Patient's life expectancy less than one year

- ✧ Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and the profile of the delivery system. Too tortuous vasculature may increase the risk of failure in stent graft deployment or proximal end release.

- ✧ The Ankura AAA Stent Graft System is not recommended in patients who weigh greater than 350 Pounds (150kg) or can not undergo accurate fluoroscopy examination due to obesity.

- ✧ Patients with a systemic infection may be at increased risk of endovascular stent-graft infection.

- ✧ Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications.

- Implant Procedure

- ✧ Systemic anticoagulation should be used during the procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

- ✧ Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.

- ✧ Maintain the guide wire position during delivery system insertion.

ot yet been evaluated

- ✧ Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the stent graft.

- ✧ Always use fluoroscopy for guidance, delivery and observation of the device within the vasculature.

- ✧ The use of the Ankura Stent Graft System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast agent used during the procedure.

- ✧ To avoid any twist in the endovascular stent-graft, during any rotation of the delivery system, rotate all of the components of the system together carefully.

- ✧ Inaccurate placement of the stent graft within the vessel may result in increased risk of endoleak, migration or embolization.

- ✧ Inadequate fixation of the stent graft may result in increased risk of migration of the stent graft, and incorrect deployment and migration of the stent graft requires surgical intervention.

- ✧ Do not continue advancing if resistance is felt during advancement of the guide wire or delivery system. Stop and assess the cause of resistance.

- ✧ Unless medically indicated, do not deploy the stent graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities.

● MRI Information

- ✧ Nonclinical testing has demonstrated that the Ankura AAA Stent Graft is MR Conditional. A patient with this device can be scanned safely in 3.0 Tesla under the following conditions:

- Static magnetic field of 3.0 Tesla only

- Spatial gradient field of ≤ 720 Gauss/cm

- Maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning

Based on non-clinical testing, the device was determined to produce a temperature rise of less than 5°C at a maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15

minutes of MR scanning. The maximum whole-body-averaged specific absorption rate was derived by calculation and verified by calorimetry.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. The image artifact extends approximately 20 mm outside the device lumen when scanned in non-clinical testing using the spin echo sequence and gradient echo sequence, respectively in a GE Healthcare SIGNA HD 3.0 T Scanner with a whole body coil..

9 Follow Up

- Periods of follow up
 - ✧ Clinical follow up should be performed at 24 hours, 1, 3, 6, 12 and 24 months after the procedure. The follow up period can be adjusted at any time according to the circumstances of each individual patient.
- Methods used during follow up
 - ✧ Haematological examination, X-ray, B-US, CTA, MRI, DSA etc.
- Follow up items:
 - ◆ Size of aneurysms;
 - ◆ Embolization;
 - ◆ Changing of aneurysms pulsatility;
 - ◆ Migration;
 - ◆ Endoleak;
 - ◆ Stent Graft distortion.
- ✧ In addition, patients with following specific findings during following up, endovascular intervention or surgical treatment should be considered:
 - ◆ Aneurysms enlarged 5mm and above;
 - ◆ Changing of aneurysms pulsatility with or without aneurysms enlarging or endoleak;
 - ◆ Continual endoleak with or without aneurysms enlarging;
 - ◆ Endoleak caused by Stent Graft migration.

10 Shelf Life

The Ankura AAA Stent Graft System is sterilized with ethylene oxide. Shelf life: 36 month.
Do not use the expired product.

11 Package and Label

The Ankura AAA Stent Graft and cuff Stent Graft System are supplied sterile. Each Ankura Stent Graft is individually contained within a Delivery System. The Ankura AAA Stent Graft with Delivery System are immobilized and protected in a PETG box respectively, then sealed with two Tyvek1073B dialyzing pouches, on which a primary label is stuck. The product is sterilized and put in a box with IFU, patient card, customer feedback form and certificate of compliance. A label is stuck on the box.

12 Warranty

The Company Lifetech Scientific (Shenzhen) Co., Ltd. China, guarantees the systems stent-grafts the requirements of normative documents at observance by the consumer of conditions of storage, transportation and operation, set in the normative document.

13 Operating conditions

The Ankura stent-graft system is operated under operating conditions by specially trained personnel.

Ambient temperature $25 \pm 10^{\circ} \text{C}$; Relative humidity 30-75%; Atmospheric pressure 40-106 kPa.

The stent-graft is implanted into the internal environment of the body and is operated at a temperature of $32-42^{\circ} \text{C}$, humidity 100%.

14 Utilization

After use, the product must be disposed of in accordance with the disposal instructions in force in the treatment facility.

The class of danger of medical waste, in accordance with SanPiN 2.1.7.2790-10: B (dangerous, potentially infected medical waste).

15 Authorized Representative of Manufacturer in the territory of Russia

Company BW, LLC



129085, Moscow, Prospect Mira 101, bld. 1 space 17 (office 306a),

Tel: +7(499)281-67-68,

Website: beawire.com,

E-mail: info@beawire.com

The symbols used are as the following:

	Date of manufacture
	Use by
	Serial number
	Sterilized using ethylene oxide
	batch code
	Do not reuse
	Consult Instructions for Use
	Keep away from the sunlight
	Keep dry
	Catalogue number
	Do not use if package is damaged
	Do not resterilize
	Non-pyrogenic
Non-toxic	Non-toxic
	Manufacturer
	Authorised representative in the European Community
	Store at room temperature
	MR conditional

The list of national and international regulatory documents/standards,

matches medical device

ISO 15223-1 Medical devices - Symbols to be used with medical product labels, labeling, and provide information - Part 1: General requirements

EN 1041 Information supplied by the manufacturer with medical devices

ISO 10555-1 and amendments: Sterile, disposable vascular catheters - Part 1: General requirements

ISO 10993-1 + Assessment of biological effects of medical devices: Evaluation and testing

ISO 10993-3 Biological Evaluation of medical treatments products: Tests for genotoxicity,

Carcinogenicity and reproductive toxicity

ISO 10993-4 Assessment of biological effects of medical devices: Selection of tests for interactions with blood

ISO 10993-5 Biological Evaluation of medical treatments products: Test for cytotoxicity in vitro

ISO 10993-6 Biological Evaluation of medical treatments products: test for the evaluation of local effects after implantation

ISO 10993-7 Biological Evaluation of medical treatments products: ethylene oxide (EtO)

Remains after sterilization

ISO 10993-10 Assessment of biological effects of medical products: Test for irritation and sensitization skin

ISO 10993-11 Biological Evaluation of medical treatments products: Tests for General toxicity

ISO 10993-12 Assessment of biological effects of medical products: sample Preparation and reference materials

ISO 11135-1 Test cycles sterilization with ethylene oxide

ISO 13485 Medical devices - quality management System - Requirements for regulatory

ISO 14155 Clinical investigation of medical devices for patient people - good clinical practice

ISO 11607-1 Packaging for final sterilization of medical products - Part 1: Requirements for materials, sterile barrier and packaging systems

ISO 14630 General requirements to the non-active surgical implants

ISO 14971 Medical devices - Application of risk management to medical devices.


ISO 25539-2 Cardiovascular implants - Endovascular devices - Part 2: Vascular stents

MDD 93/42/EEC Basic requirements of the Directives on medical devices, as amended in Annex 1 2007/47/EEC

Manufacturer:

Lifetech Scientific (Shenzhen) Co., Ltd
Floor 1-5, Cybio Electronic Building,
Langshan 2nd Street, North Area of High-tech Park,
Nanshan District, Shenzhen 518057
P. R. China
Tel: +86 755 86026250
Fax: +86 755 86026251
Email: cs@lifetechmed.com


EU authorized representative

LifeTech Scientific (Europe) Coöperatief U.A
Kruisdonk 64, 6222 PH Maastricht, the Netherlands

Document No.LT/TS/121-02RU V1.0 2017-08-08

CE 0344

Ankura TAA Stent Graft System

Instructions for Use

Regulatory Affairs Director
Lifetech Scientific (Shenzhen) Co., Ltd

Lily Shi Name

(sign and seal)

Dec.12,2017

Lifetech Scientific (Shenzhen) Co., Ltd

Contents

1	Device Description	3
2	Indications for Use.....	7
3	Contraindications.....	8
4	Pre-operative Preparation	8
5	Operation Procedure.....	10
6	Potential Adverse Events	14
7	Storage and Transport Requirements	15
8	Warnings	15
9	Follow Up	18
10	Shelf Life	19
11	Package and Label.....	19
12	Warranty	19
13	Operating conditions	19
14	Utilization.....	20
15	Authorized Representative of Manufacturer in the territory of Russia.....	20
	The symbols used are as the following:	21
	The list of national and international regulatory documents/standards, which matches medical device	22


Read all instructions carefully. Failure to follow the instructions, warnings and precautions may lead to serious complications or injury to the patients.

1 Device Description

● Ankura TAA Stent Graft System

- ♦ The Ankura TAA Stent Graft System includes:
 - ♦ The Ankura TAA Stent Graft
 - ♦ The Ankura Delivery System
- ♦ The Ankura TAA Stent Graft is pre-loaded into the Delivery System. The loaded delivery system is inserted endoluminally via the femoral or iliac artery and tracked through the patient's vasculature to deliver the stent graft to the target site

● Ankura TAA Stent Graft



Fig.1 Ankura TAA Stent Graft

- The Ankura stent graft is a tube composed of ePTFE film supported by a metal mesh called a stent. The mark points are fixed on the stent. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.

● Ankura Delivery System

1	Tapered Tip	2	Markband	3	Proximal Anchor
4	Stent Graft	5	Inner Tube	6	Immobile Pole
7	Sheath	8	Front Grip	9	Slider Grip
10	Trigger	11	Screw Gear	12	Haemostatic Vale
13	Rear Grip	14	Proximal Releaser	15	Secure Buckle
16	Luer Connector				

Fig.2 Ankura Delivery System

- The Ankura stent graft is loaded inside the Ankura Delivery System. The Ankura Delivery System facilitates the placement of the stent graft via the arterial vasculature (e.g., femoral arteries). Using fluoroscopic guidance, the Ankura Delivery System is properly positioned within the patient's vasculature and the stent graft is deployed from the delivery system.

- Device specifications

- The specification of Ankura TAA Stent Graft System is indicated as follow and the details are shown in table 1.

TAA 26 22 B 080


Table 1 Summary of Ankura TAA Stent Graft System

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
20	20	40,60,80,100,120	21F
22	22	40,60,80,100,120,	21F
24	20	60,80,100,120,140,160	21F
	24	40,60,80,100,120,140,160	
26	22	60,80,100,120,140,160	21F
	26	40,60,80,100,120,140,160	
28	22	120,140,160,180,200	21F
	24	60,80,100,120,140,160,180,200	
	28	40,60,80,100,120,140,160,180,200	
30	22	120,140,160,180,200	21F
	24	120,140,160,180,200	
	26	60,80,100,120,140,160,180,200	
	30	40,60,80,100,120,140,160,180,200	
32	24	120,140,160,180,200	21F
	26	120,140,160,180,200	

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
	28	120,140,160,180,200	
	32	40,60,80,100,120,140,160,180,200	
34	26	120,140,160,180,200	23F
	28	120,140,160,180,200	
	30	120,140,160,180,200	
	34	40,60,80,100,120,140,160,180,200	
36	28	120,140,160,180,200	23F
	30	120,140,160,180,200	
	32	120,140,160,180,200	
	36	40,60,80,100,120,140,160,180,200	
38	30	120,140,160,180,200	23F
	36	120,140,160,180,200	
	34	120,140,160,180,200	
	38	40,60,80,100,120,140,160,180,200	
40	32	120,140,160,180,200	24F
	34	120,140,160,180,200	
	36	120,140,160,180,200	
	40	40,60,80,100,120,140,160,180,200	
42	34	120,140,160,180,200	24F
	36	120,140,160,180,200	
	38	120,140,160,180,200	
	42	40,60,80,100,120,140,160,180,200	
44	36	120,140,160,180,200	24F
	38	120,140,160,180,200	
	40	120,140,160,180,200	
	44	40,60,80,100,120,140,160,180,200	
46	38	120,140,160,180,200	24F
	40	120,140,160,180,200	
	42	120,140,160,180,200	
	46	40,60,80,100,120,140,160,180,200	

The materials from Ankura TAA Stent Graft System

Stent-Graft – Nitinol; Covered film – ePTFE LQ-ePTFE; Steel tube –SUS316L;

Platinum marker – Platinum 0213; Suture – Polypropylene. Delivery System: The sheath: Sheath - Polytetrafluoroethylene (PTFE) liner Polyflon F-205; Stainless steel spring –304; marker band –Pt 90 / Ir 10; outer layer of sheath – Polyether block amides Arkema Pebax; Hydrophilic coating Coatings2GO. Shell handle – Acrylonitrile Butadiene Styrene Sabic; stainless steel pipe - 304; Silicon gasket –Dow Corning; Hose - Thermoplastic polyurethane TPU58315. Haemostatic Valve: Polycarbonate (PC) RX1805; Stopcock (three-way) - Acrylonitrile Butadiene Styrene 757K; Handle of Stopcock (three-way) – polyethylene PE2911. Inner tubing: Sleeve – stainless steel 304; Tapered Tip –Polyamide Pebax + BaSO4 + TiO2; Proximal anchor – stainless steel 304; PI tubing – Polyimides. Immobile pole: Immobile pole – polyethylene + BaSO4; metal tube - stainless steel 304; Locked connector - Acrylonitrile Butadiene Styrene Sabic; Gripping handle - stainless steel 304.

Physicians should strictly adhere to indications and contraindications.

2 Indications for Use

- Ankura TAA Stent Graft System is indicated for endovascular repair of patients with thoracic aortic aneurysms or thoracic aortic dissection, requiring:
 - ◆ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - ◆ Non-aneurysmal aortic proximal neck length ≥ 1.5 cm.
 - ◆ Non-aneurysmal aortic distal anchorage zone ≥ 1.5 cm.
 - ◆ Non-aneurysmal aortic diameter in the range of 18-42mm.
 - ◆ Morphology suitable for endovascular repair.
- Additional considerations for patient selection include:
 - ◆ Patient's age and life expectancy;
 - ◆ Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity);
 - ◆ Toleration to the general, regional or local anesthesia;

The final decision is at the discretion of the physician and patient.
- Device Selection
 - ◆ The size of the Ankura Stent Graft should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversizing of the Ankura Stent Graft is 10%-20% to the landing zone vessel inner diameter..

3 Contraindications

- The Ankura TAA Stent Graft is contraindicated in:
 - ◆ Patients with acute systemic infection.
 - ◆ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ◆ Patients who have allergic reaction to the device;
 - ◆ Patients who are not suitable for endovascular repair in vascular morphology;
 - ◆ Patients who can not tolerate contrast agents due to renal insufficiency;
 - ◆ Patients who are allergic to contrast agents;
 - ◆ Aneurysms neck with thrombus.
 - ◆ Non-aneurysmal aortic proximal neck length <1.5cm.
 - ◆ Non-aneurysmal aortic distal anchorage zone <1.5cm
 - ◆ Non-aneurysmal aortic diameter <18mm or >42mm.

4 Pre-operative Preparation

Read the instructions carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help to guide the physician and do not replace the physician judgment.

- Materials required
 - ◆ 0.035"/260cm Super Stiff Guidewire. For example, Surpass Super Stiff Guidewire, Amplatz Super Stiff Guidewire or comparable.
 - ◆ Fluoroscope with digital angiography capabilities;
 - ◆ Angiographic catheter, and marker catheter;
 - ◆ Contrast agent;
 - ◆ High pressure Syringe;
 - ◆ Heparinized saline solution;
 - ◆ Any other devices required by the physician for the operation (including devices used for femoral artery dissection).
- A multidisciplinary team that has combined procedural experience with:
 - ◆ Local anesthesia or general anesthesia;

- Femoral artery exposure, arteriotomy and repair;
- Percutaneous access and closure techniques;
- Fluoroscopic and angiographic image interpretation;
- Endovascular stent placement;
- Appropriate use of radiographic contrast agent;
- Techniques to minimize radiation exposure;
- Expertise in necessary patient follow-up modalities;
- Experience in traditional vascular surgery for any possibility of failure in EVAR.

- Products preparation

- The Stent Graft diameter should be selected preferred to vessel inner diameter not outer diameter. The uncorrectable size of the device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times larger than the landing zone vessel inner diameter.
- Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, re-review the specification conforming to the patient required .

- Patient preparation

- Measure proximal aortic neck length and vessel inner diameter.
- Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
- Evaluate the quality of aortic neck.
- Assess the morphology of the vessel compatible with endovascular aortic repair.
- Undergo any necessary common examination, and monitoring of vital signs.
- Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
- Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
- Adequate iliac/femoral artery access and expose one common femoral artery using standard surgical technique.
- Systemic anticoagulation should be used during the deployment procedure based on

hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

5 Operation Procedure

- Preparation before implant
 - ◆ Attach the syringe with heparinized saline to the luer hub and hemostatic valve and flush separately to evacuate air from the inner tube and the outer sheath. (Fig.3)


Fig 3 Flush the inner tube and the outer sheath

- ◆ Monitor important vital signs.
- ◆ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution; sheath and guidewire should be flushed after each exchange.
- ◆ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)


Fig. 4 Activate the hydrophilic coating

- ◆ Expose one side femoral arteries by standard surgery.
- ◆ Place a 0.035"/260cm super stiff guidewire in the femoral artery and advance it until to

the ascending artery. (Fig. 5)

Fig. 5 Introduction of super stiff guidewire

- Deploy the stent graft
 - ♦ Advance the delivery system to the aneurysms or dissection sites along the guidewire, and always maintain the hemostatic valve toward the left side of the patient (Fig. 6)


Fig. 6 Introduction of the delivery system

Note: The connecting bar of the TAA Stent Graft must be in the outer curved side of the aortic arch. And the connecting bar has the same side with the hemostatic valve of the delivery system when the TAA Stent Graft is preloaded into the delivery system.

- ♦ Ensure the Stent Graft positioning at the aneurysms or dissection sites. Ensure the connecting bar of the TAA Stent Graft in the outer curved side of aortic arch.
- ♦ Repeat the angiography to verify the delivery system locating at the aneurysms or

dissection sites.

- Decrease the systolic blood pressure (SBP) to 10.5 to 12 kPa.
- Slowly **deploy** the stent graft. Hold the front grip of the delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the proximal part of the stent graft until the covered stent adheres to the vascular wall. (Fig. 7)


Fig. 7 Slowly deploy the stent graft

- Rapidly deploy the stent graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 8)


Fig. 8 Rapidly deploy the stent graft

- Proximal end release

- Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 9)


Fig. 9 Remove the Secure Buckle

- Pull the Proximal Releaser to release the proximal end of the stent graft. (Fig. 10)


Fig. 10 Pull the Proximal Releaser

- Verify the position and the state of the completely deployed stent graft under fluoroscopy.
- Remove the deliver system
 - Hold the slider grip of the delivery system stationary, unlock the trigger and pull back the front grip to make the tip back to the distal end of the out sheath. (Fig.11)


Fig. 11 Pull back the tip

Note: Pull back the front grip should be gently and under fluoroscopy guidance to avoid stent graft migration.

- ♦ Remove the delivery system and guidewire completely and close the puncture site with appropriate techniques.

6 Potential Adverse Events

Adverse events associated with use of the Ankura TAA Stent Graft System include, but are not limited to the following:

- Aneurysms enlargement
- Aneurysms/vessel rupture
- Arterial or venous thrombosis (including thrombosis and embolization)
- Arterial damage
- Arterial perforation
- Arteriovenous fistula or pseudoaneurysm
- Bleeding/Hematoma
- Bowel complications
- Cardiac failure/Infarction
- Conversion to open surgery
- Death
- Edema
- Infection
- Neurological damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- Pulmonary complications
- Renal function insufficiency/kidney failure
- Wound complications and subsequent attendant problems
- Device complications due to:
 - ♦ Insertion and removal difficulties
 - ♦ Release difficulties
 - ♦ Failure to deploy
 - ♦ Inaccurate placement
 - ♦ Endoleak

- Stent migration
- Stent fracture
- Rupture of graft material
- Erosion
- Stent thrombosis or embolization

7 Storage and Transport Requirements

- The Ankura TAA Stent Graft is pre-loaded in delivery system, and stored in sterile package.
- Store in a non-corrosive gas, well-ventilated, clean room.
- During transportation, no extrusion, no wet exposure, handle with care.
- No exposure to organic solvent.
- Recommendations on storage: temperature from 0° to 30° C, humidity 30-85%, atmospheric pressure is not controlled.
- Recommendations on transportation: temperature from -30° to 55° C, humidity 30-85%, atmospheric pressure is not controlled.

8 Warnings

- General
 - ♦ Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
 - ♦ The Ankura TAA Stent Graft can not be used after resterilization.
 - ♦ If any sign of damage is observed prior to use, do not use the product and return it to Lifetech.
 - ♦ This device should only be used by physicians and teams trained in vascular interventional techniques and in the use of device.
 - ♦ The long term performance of endovascular grafts has not yet been confirmed, all patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft (e.g., endoleak, migration, distortion, aneurysms enlarging, dissection tearing etc.).
 - ♦ Clinical studies show that uncorrectable patient selection may result in poor curative effect. In addition, if unsuccessful endovascular aortic repair require open surgery may

result in an increased risk of complications relative to the two procedures.

- ◆ Patient's condition needs to be evaluated before operation on the following aspects:
 - cardiac failure/infarction, pulmonary complications, edema, renal function insufficiency.
 - ◆ After the device placement, patients should be regularly monitored for the circumstances of the device (including endoleak, migration, aneurysms enlarging, dissection tearing etc.). At a minimum, annual X-ray, contrast or non-contrast CT imaging is required.
 - ◆ Patients with specific clinical findings, e.g, enlarging aneurysms or tearing dissection or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms or dissection size and/or persistent endoleak may lead to aneurysms or dissection rupture.
 - ◆ Patients with infection after the device placement should administer oral antibiotics.
 - ◆ Always have a vascular surgery team available during stent-graft implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.
 - ◆ Do not re-use or re-sterilize the Ankura TAA Stent Graft System.
- Patients selection, treatment and follow-up
 - ◆ The safety and effectiveness of the Ankura TAA Stent Graft has not yet been evaluated in following populations:
 - ◆ Traumatic aortic injury
 - ◆ Uncorrectable coagulopathy
 - ◆ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
 - ◆ Patients with active systemic infections
 - ◆ Pregnant or nursing females
 - ◆ Morbid obesity
 - ◆ Less than 18 years of age
 - ◆ Patient's life expectancy less than one year
 - ◆ Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular

access techniques and the profile of the delivery system. Too tortuous vasculature may increase the failure in stent graft deployment or proximal end release.

- The Ankura TAA Stent Graft System is not recommended in patients who weight greater than 350 Pounds (150kg) or can not undergo accurate fluoroscopy examination due to obesity
- Patients with a systemic infection may be at increased risk of endovascular stent-graft infection
- Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications.

● Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain the guide wire position during delivery system insertion.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the stent graft.
- Always use fluoroscopy for guidance, delivery and observation of the device within the vasculature.
- The use of the Ankura TAA Stent Graft System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast agent used during the procedure.
- To avoid any twist in the endovascular stent graft, during any rotation of the delivery system, be careful to rotate all the components of the system together.
- Inaccuracy placement of the stent graft within the vessel may result in increased risk of endoleak, migration, or embolization.
- Inadequate fixation of the stent graft may result in increased risk of migration of the stent graft, and incorrect deployment of the stent graft requires surgical intervention.

- ♦ Do not continue advancing if resistance is felt during advancement of the guide wire delivery system. Stop and assess the cause of resistance.
- ♦ Unless medically indicated, do not deploy the stent graft in a location that will occlude arteries necessary to supply blood flow to organs or extremities.

● MRI Information

Nonclinical testing has demonstrated that the Ankura TAA Stent Graft is MR Conditional. A patient with this device can be scanned safely in 3.0 Tesla under the following conditions:

- ♦ Static magnetic field of 3.0 Tesla only
- ♦ Spatial gradient field of ≤ 720 Gauss/cm
- ♦ Maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning

Based on non-clinical testing, the device was determined to produce a temperature rise of less than 5°C at a maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of MR scanning. The maximum whole-body-averaged specific absorption rate (SAR) was derived by calculation and verified by calorimetry.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. The image artifact extends approximately 20 mm outside the device lumen when scanned in non-clinical testing using the spin echo sequence and gradient echo sequence, respectively in a GE Healthcare SIGNA HD 3.0 T Scanner with a whole body coil.

9 Follow Up

● Periods of follow up

- ♦ Clinical follow up should be performed at 24 hours, 1, 3, 6, 12 and 24 months after the procedure. The follow up period can be adjusted at anytime according to the circumstances of each individual patient.

● Methods used during follow up

- ♦ Haematological examination, X-ray, B-US, CTA, MRI, DSA etc.

● Follow up items:

- ♦ Size of aneurysms or dissection;

- ♦ Embolization;
- ♦ Changing of aneurysms pulsatility or dissection;
- ♦ Migration;
- ♦ Endoleak;
- ♦ Stent Graft distortion.
- ♦ In addition, patients with following specific findings during following up, endovascular intervention or surgical treatment should be considered:
 - ♦ Aneurysms enlarge 5mm and above;
 - ♦ Changing of aneurysms or dissection pulsatility with or without aneurysms enlarging or endoleak or dissection tearing;
 - ♦ Continual endoleak with or without aneurysms enlarging or dissection tearing;
 - ♦ Stent Graft migration results in endoleak.

10 Shelf Life

- The Ankura TAA Stent Graft System is sterilized with ethylene oxide. The shelf life is marked on the label. Do not use the expired product.
- Shelf life 36 month.

11 Package and Label

- The Ankura TAA Stent Graft with Delivery System is immobilized and protected in a PETG box, then sealed with two Tyvek1073B dialyzing pouches, on which a primary label is stuck. The product is sterilized and put in a box with IFU, patient card, customer feedback form and certificate of compliance. A label is stuck on the box.

12 Warranty

The Company Lifetech Scientific (Shenzhen) Co., Ltd. China, guarantees the systems stent-grafts the requirements of normative documents at observance by the consumer of conditions of storage, transportation and operation, set in the normative document.

13 Operating conditions

The Ankura stent-graft system is operated under operating conditions by specially trained personnel.

Ambient temperature $25 \pm 10^{\circ} \text{C}$; Relative humidity 30-75%; Atmospheric pressure 40-106 kPa.

The stent-graft is implanted into the internal environment of the body and is operated at a temperature of $32-42^{\circ} \text{C}$, humidity 100%.

14 Utilization

After use, the product must be disposed of in accordance with the disposal instructions in force in the treatment facility.

The class of danger of medical waste, in accordance with SanPiN 2.1.7.2790-10: B (dangerous, potentially infected medical waste).

15 Authorized Representative of Manufacturer in the territory of Russia

Company BW, LLC

129085, Moscow, Prospect Mira 101, bld. 1 space 17 (office 306a),

Tel: +7(499)281-67-68,

Website: beawire.com,

E-mail: info@beawire.com

... specially trained
... 40-106 kPa.
... at a
... symbols used are as the following:

	Date of manufacture
	Use by
SN	Serial number
STERILE E0	Sterilized using ethylene oxide
LOT	batch code
	Do not reuse
	Consult Instructions for Use
	Keep away from the sunlight
	Keep dry
REF	Catalogue number
	Do not use if package is damaged
	Do not resterilize
	Non-pyrogenic
Non-toxic	Non-toxic
	Manufacturer
EC REP	Authorised representative in the European Community
	Store at room temperature
	MR conditional


The list of national and international regulatory documents/standards, which matches medical device

ISO 15223-1 Medical devices - Symbols to be used with medical product labels, labeling, and provide information - Part 1: General requirements
EN 1041 Information supplied by the manufacturer with medical devices
ISO 10555-1 and amendments: Sterile, disposable vascular catheters - Part 1: General requirements
ISO 10993-1 + Assessment of biological effects of medical devices: Evaluation and testing
ISO 10993-3 Biological Evaluation of medical treatments products: Tests for genotoxicity, Carcinogenicity and reproductive toxicity
ISO 10993-4 Assessment of biological effects of medical devices: Selection of tests for interactions with blood
ISO 10993-5 Biological Evaluation of medical treatments products: Test for cytotoxicity in vitro
ISO 10993-6 Biological Evaluation of medical treatments products: test for the evaluation of local effects after implantation
ISO 10993-7 Biological Evaluation of medical treatments products: ethylene oxide (EtO) Remains after sterilization
ISO 10993-10 Assessment of biological effects of medical products: Test for irritation and sensitization skin
ISO 10993-11 Biological Evaluation of medical treatments products: Tests for General toxicity
ISO 10993-12 Assessment of biological effects of medical products: sample Preparation and reference materials
ISO 11135-1 Test cycles sterilization with ethylene oxide
ISO 13485 Medical devices - quality management System - Requirements for regulatory
ISO 14155 Clinical investigation of medical devices for patient people - good clinical practice
ISO 11607-1 Packaging for final sterilization of medical products - Part 1: Requirements for materials, sterile barrier and packaging systems
ISO 14630 General requirements to the non-active surgical implants
ISO 14971 Medical devices - Application of risk management to medical devices.
ISO 25539-2 Cardiovascular implants - Endovascular devices - Part 2: Vascular stents
MDD 93/42/EEC Basic requirements of the Directives on medical devices, as amended in Annex 1 2007/47/EEC

Manufacturer:

Lifetech Scientific (Shenzhen) Co., Ltd
Floor 1-5, Cybio Electronic Building,
Langshan 2nd Street, North Area of High-tech Park,
Nanshan District, Shenzhen 518057
P. R. China
Tel: +86 755 86026250
Fax: +86 755 86026251
Email: cs@lifetechmed.com

EU authorized representative:

LifeTech Scientific (Europe) Coöperatief U.A
Kruisdonk 64, 6222 PH Maastricht, the Netherlands

Document No. LT/TS/12I-01E V1.7A 2017-08-08

Ankura AUI Stent Graft System

Instructions for Use

Regulatory Affairs Director
Lifetech Scientific (Shenzhen) Co., Ltd

City Shi Name

(sign and seal)

Dec.12,2017

Lifetech Scientific (Shenzhen) Co., Ltd

Contents

1	Device Description	3
2	Indications for Use.....	6
3	Contraindications	6
4	Pre-operative Preparation	7
5	Operation Procedure.....	8
6	Potential Adverse Events	12
7	Storage and Transport Requirements	13
8	Warnings	14
9	Follow Up	17
10	Shelf Life	18
11	Package and Label.....	18
12	Warranty	18
13	Operating conditions	18
14	Utilization.....	18
15	Authorized Representative of Manufacturer in the territory of Russia.....	19
	The symbols used are as the following:	20
	The list of national and international regulatory documents/standards, which matches medical device	21

Read all instructions carefully. Failure to follow the instructions, warnings and precautions may lead to serious complications or injury to the patients.

1 Device Description

- Ankura AUI Stent Graft System

- ✧ The Ankura AUI Stent Graft System includes:

- ◆ The Ankura AUI Stent Graft
- ◆ The Ankura Delivery System
- ◆ The Surpass Super Stiff Guidewire (optional)

- ✧ The Ankura AUI Stent Graft is pre-loaded into the Delivery System. The loaded delivery system is inserted endoluminally via the femoral or iliac artery and tracked through the patient's vasculature to deliver the stent graft to the target site.

- Ankura AUI Stent Graft



Fig.1 Ankura AUI Stent Graft

The Ankura stent graft is a tube composed of ePTFE film supported by a metal mesh called a stent (**Fig.1**). The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge.

The stent can provide stable and expansionary force necessary to open the graft.

- **Ankura Delivery System**

1	Tapered Tip	2	Markband	3	Proximal Anchor
4	Stent Graft	5	Inner Tube	6	Immobile Pole
7	Sheath	8	Front Grip	9	Slider Grip
10	Trigger	11	Screw Gear	12	Haemostatic Vale
13	Rear Grip	14	Proximal Releaser	15	Secure Buckle
16	Luer Connector				

Fig.2 Ankura Delivery System

✧ The Ankura AUI stent graft is loaded inside the Ankura Delivery System. The Ankura Delivery System facilitates the placement of the stent graft via the arterial vasculature (e.g., femoral arteries). Using fluoroscopic guidance, the Ankura Delivery System is properly positioned within the patient's vasculature and the stent graft is deployed from the delivery system.

- **Device specifications**

✧ The specification of Ankura AUI Stent Graft Delivery is indicated as follow and the details are shown in table 1.


Table 1 Summary of Ankura AUI Stent Graft System

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
20	12	120,140,160,180	21F

22	12	120,140,160,180	21F
	14	120,140,160,180	21F
24	12	120,140,160,180	21F
	14	120,140,160,180	21F
	16	120,140,160,180	21F
26	12	120,140,160,180	21F
	14	120,140,160,180	21F
	16	120,140,160,180	21F
28	12	120,140,160,180	21F
	14	120,140,160,180	21F
	16	120,140,160,180	21F
30	12	120,140,160,180	21F
	14	120,140,160,180	21F
	16	120,140,160,180	21F
32	12	120,140,160,180	21F
	14	120,140,160,180	21F
	16	120,140,160,180	21F
34	12	120,140,160,180	23F
	14	120,140,160,180	23F
	16	120,140,160,180	23F
36	12	120,140,160,180	23F
	14	120,140,160,180	23F
	16	120,140,160,180	23F

The materials from Ankura AUI Stent Graft System

Stent-Graft – Nitinol; Covered film – ePTFE LQ-ePTFE; Steel tube –SUS316L; Platinum marker – Platinum 0213; Suture – Polypropylene. Delivery System: The sheath: Sheath - Polytetrafluoroethylene (PTFE) liner Polyflon F-205; Stainless steel spring –304; marker band –Pt 90 / Ir 10; outer layer of sheath – Polyether block amides Arkema Pebax; Hydrophilic coating Coatings2GO. Shell handle – Acrylonitrile Butadiene Styrene Sabc; stainless steel pipe - 304; Silicon gasket –Dow Corning; hose - Thermoplastic polyurethane TPU58315. Haemostatic Vale: Polycarbonate (PC) RX1805; Stopcock (three-way) - Acrylonitrile Butadiene Styrene 757K; handle of Stopcock (three-way) – polyethylene PE2911. Inner tubing: Sleeve – stainless steel 304; Tapered Tip – Polyamide Pebax + BaSO₄ + TiO₂; Proximal anchor – stainless steel 304; PI tubing – Polyimides. Immobile pole : Immobile pole –polyethylene + BaSO₄; metal tube - stainless steel 304; Locked connector - Acrylonitrile Butadiene Styrene Sabc; Gripping handle - stainless steel 304.

Physicians should strictly adhere to indications and contraindications.

2 Indications for Use

- Ankura AUI Stent Graft with delivery system is indicated for endovascular repair of patients with abdominal aortic aneurysms, requiring:
 - ✧ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
 - ✧ Non-aneurysmal proximal aortic neck length $\geq 15\text{mm}$.
 - ✧ Non-aneurysmal proximal aortic neck diameter of 18–32mm.
 - ✧ Proximal aortic neck angulation $\leq 60^\circ$.
 - ✧ Distal iliac artery anchorage zone $\geq 15\text{mm}$.
 - ✧ Distal iliac artery diameter of 10–16mm.
 - ✧ Morphology suitable for endovascular repair.
- Additional considerations for patient selection:
 - ✧ Patient's age and life expectancy;
 - ✧ Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity);
 - ✧ Toleration to the general, regional or local anesthesia;

The final decision is at the discretion of the physician and patient.
- Device Selection
 - ✧ The size of the Ankura Stent Graft should be appropriate to fit the patient's aortic anatomy. Lifetech recommends an oversizing of the Ankura Stent Graft by 10%-20%.

3 Contraindications

- The Ankura AUI Stent Graft is contraindicated in:
 - ✧ Patients with acute systemic infection;
 - ✧ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
 - ✧ Patients with accessory renal artery original from abdominal artery;
 - ✧ Patients who have allergic reaction to the device;
 - ✧ Patients who are not suitable for endovascular repair in vascular morphology;
 - ✧ Patients who can not tolerate contrast agents due to renal insufficiency;

- ✧ Patients who are allergic to contrast agents;
- ✧ Aneurysms neck with thrombus.
- ✧ Non-aneurysmal proximal aortic neck length <15mm.
- ✧ Non-aneurysmal proximal aortic neck diameter <18mm or >32mm.
- ✧ Proximal aortic neck angulation > 60°.
- ✧ Distal iliac artery anchorage zone <15mm.
- ✧ Distal iliac artery diameter <10mm or >16mm.

4 Pre-operative Preparation

Read the instructions carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help to guide the physician and do not replace the physician judgment.

- Materials required
 - ✧ 0.035"/260cm Super Stiff Guidewire. For example, Surpass Super Stiff Guidewire, Amplatz Super Stiff Guidewire or comparable.
 - ✧ Fluoroscope with digital angiography capabilities;
 - ✧ Angiographic catheter, and marker catheter;
 - ✧ Contrast agent;
 - ✧ High pressure syringe;
 - ✧ Heparinized saline solution;
 - ✧ Any other devices required by the physician for the operation (including devices used for femoral artery dissection).
- A multidisciplinary team that has combined procedural experience with:
 - ✧ Local anesthesia or general anesthesia;
 - ✧ Femoral artery exposure, arteriotomy and repair;
 - ✧ Percutaneous access and closure techniques;
 - ✧ Fluoroscopic and angiographic image interpretation;
 - ✧ Endovascular stent placement;
 - ✧ Appropriate use of radiographic contrast agent;

- ✧ Techniques to minimize radiation exposure;
 - ✧ Expertise in necessary patient follow-up modalities;
 - ✧ Experience in traditional vascular surgery for any possibility of failure in EVAR.
- **Products preparation**
 - ✧ The Stent Graft diameter should be selected preferred to vessel inner diameter not outer diameter. The uncorrectable size of the device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times larger than the landing zone vessel inner diameter.
 - ✧ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, re-review the specification conforming to the patient required .
- **Patient preparation**
 - ✧ Measure proximal aortic neck length and vessel inner diameter.
 - ✧ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
 - ✧ Evaluate the quality of aortic neck.
 - ✧ Assess the morphology of the vessel compatible with endovascular aortic repair.
 - ✧ Undergo any necessary common examination, and monitoring of vital signs.
 - ✧ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs
 - ✧ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
 - ✧ Adequate iliac/femoral arteries access and expose both common femoral arteries assess using standard surgical technique.
 - ✧ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.

5 Operation Procedure

- **Preparation before implant**

- ✧ Attach the syringe with heparinized saline to the luer hub and hemostatic valve and flush separately to evacuate air from the inner tube and the outer sheath. (Fig.3)


Fig 3 Flush the inner tube and the outer sheath

- ✧ Monitor important vital signs.
- ✧ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution, sheath and guidewire should be flushed after each exchange.
- ✧ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)


Fig. 4 Activate the hydrophilic coating

- ✧ Expose bilateral femoral arteries by standard surgery.
- ✧ Insert one 0.035"/260cm super stiff guidewire in the femoral artery, and advance it above the renal arteries.
- ✧ Insert marker sheath, and determine again the stent specification using radiography.
- Deploy the stent graft
 - ✧ Advance the delivery system to the abdominal aneurysm along the guidewire, and always maintain haemostasis valve in opposite side of the patient. (Fig. 5)


Fig. 5 Introduction of the delivery system

Note, The connecting bar of the AUI Stent Graft must be in the opposite side of the patient. And the connecting bar has the same side with the hemostatic valve of the delivery system when the AUI Stent Graft is preloaded into the delivery system.

- ✧ Inject contrast agents into the abdominal aorta; mark the position of the target location, either on the imaging screen or with radiopaque marker or angiographic ruler on the patient's body. Adjust the position of the stent graft such that the top edge of the graft fabric is just below the lowest renal artery.
- ✧ Slowly deploy the stent graft. Hold the front grip of the delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the stent graft until the covered stent adhere to the vascular wall.
- ✧ Rapidly deploy the stent graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. **(Fig. 6)**


Fig.6 Withdraw outer sheath

- Proximal end release
 - ✧ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 7)


Fig. 7 Remove the Secure Buckle

- ✧ Pull the Proximal Releaser to release the proximal spring of the stent graft. (Fig. 8)


Fig. 8 Pull the Proximal Releaser

- Hold the slider grip of the delivery system stationary, unlock the trigger and pull back the front grip to make the tip back to the distal end of the out sheath. And then remove the delivery system.(Fig. 9)


Fig. 9 Remove the delivery system

- Place a plug in the other side of iliac artery. (Fig. 10)


Fig. 10 Plug placement

- Build a bypass between the two femoral arteries.

6 Potential Adverse Events

Adverse events associated with use of the Ankura AUI Stent Graft System include, but are not limited to the following:

- Aneurysms enlargement
- Aneurysms/vessel rupture
- Arterial or venous thrombosis (including thrombosis and embolization)

- Arterial damage
- Arterial perforation
- Arteriovenous fistula or pseudoaneurysm
- Bleeding/Hematoma
- Bowel complications
- Cardiac failure/Infarction
- Conversion to open surgery
- Death
- Edema
- Erectile dysfunction
- Infection
- Neurological damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- Pulmonary complications
- Renal function insufficiency/kidney failure
- Wound complications and subsequent attendant problems
- Device complications due to:
 - ✧ Insertion and removal difficulties
 - ✧ Release difficulties
 - ✧ Failure to deploy
 - ✧ Inaccurate placement
 - ✧ Endoleak
 - ✧ Stent migration
 - ✧ Stent fracture
 - ✧ Rupture of graft material
 - ✧ Erosion
 - ✧ Stent thrombosis or embolization

7 Storage and Transport Requirements

- The Ankura AUI Stent Graft is pre-loaded in delivery system, and stored in sterile package.
- Store in a non-corrosive gas, well-ventilated, clean room..
- During transportation, no extrusion, no wet exposure, handle with care.

- No exposure to organic solvent.
- Recommendations on storage: temperature from 0° to 30° C, humidity 30-85%, atmospheric pressure is not controlled.
- Recommendations on transportation: temperature from -30° to 55° C, humidity 30-85%, atmospheric pressure is not controlled.

8 Warnings

- General
 - ✧ Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
 - ✧ This Ankura AUI Stent Graft should only be used by physicians and teams trained in vascular interventional techniques and in the use of device.
 - ✧ The long term performance of endovascular grafts has not yet been confirmed, all patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft (e.g., endoleak, migration, distortion, aneurysms enlarging etc.).
 - ✧ Clinical studies show that uncorrectable patient selection may result in poor curative effect. In addition, if unsuccessful endovascular aortic repair require open surgery may result in an increased risk of complications relative to the two procedures.
 - ✧ Patient's condition needs to be evaluated before operation on the following aspects, cardiac failure/infarction, pulmonary complications, edema, renal function insufficiency.
 - ✧ After the device placement, patients should be regularly monitored for the circumstances of the device (including endoleak, migration, change of the Stent Graft structure or position, aneurysms enlarging etc.). At a minimum, annual X-ray, contrast or non-contrast CT imaging is required.
 - ✧ Patients with specific clinical findings, e.g, enlarging aneurysms or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms size and/or persistent endoleak may lead to aneurysms rupture.
 - ✧ Patients with infection after the device placement should administer oral antibiotics.
 - ✧ Always have a vascular surgery team available during stent-graft implantation or

re-intervention procedures in the event that conversion to open surgical repair is necessary.

❖ Do not re-use or re-sterilize the Ankura AUI Stent Graft System.

- Patients selection, treatment and follow-up

❖ The safety and effectiveness of the Ankura AUI Stent Graft has not yet been evaluated in following populations:

- ♦ Traumatic aortic injury
- ♦ Uncorrectable coagulopathy
- ♦ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
- ♦ Concomitant thoracic aortic or thoraco abdominal aneurysms
- ♦ Patients with active systemic infections
- ♦ Pregnant or nursing females
- ♦ Morbid obesity
- ♦ Less than 18 years of age
- ♦ Patient's life expectancy less than one year

❖ Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and the profile of the delivery system. Too tortuous vasculature may increase the failure in stent graft deployment or proximal end release.

❖ The Ankura AUI Stent Graft System is not recommended in patients who weight greater than 350 Pounds (150kg) or can not undergo accurate fluoroscopy examination due to obesity

❖ Patients with a systemic infection may be at increased risk of endovascular stent-graft infection

❖ Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications

- Implant Procedure

❖ Systemic anticoagulation should be used during the implantation procedure based on

- ♦ Changing of aneurysms pulsatility with or without aneurysms enlarging or endoleak;
- ♦ Continual endoleak with or without aneurysms enlarging;
- ♦ Stent Graft migration results in endoleak.

10 Shelf Life

- The Ankura AUI Stent Graft System is sterilized with ethylene oxide. The shelf life is marked on the label. Do not use the expired product.
- Shelf life 36 month.

11 Package and Label

- The Ankura AUI Stent Graft with Delivery System is immobilized and protected in a PETG box, then sealed with two Tyvek1073B dialyzing pouches, on which a primary label is stuck. The product is sterilized and put in a box with IFU, patient card, customer feedback form and certificate of compliance. A label is stuck on the box.

12 Warranty

The Company Lifetech Scientific (Shenzhen) Co., Ltd. China, guarantees the systems stent-grafts the requirements of normative documents at observance by the consumer of conditions of storage, transportation and operation, set in the normative document.

13 Operating conditions

The Ankura stent-graft system is operated under operating conditions by specially trained personnel.

Ambient temperature $25 \pm 10^{\circ} \text{C}$; Relative humidity 30-75%; Atmospheric pressure 40-106 kPa.

The stent-graft is implanted into the internal environment of the body and is operated at a temperature of $32-42^{\circ} \text{C}$, humidity 100%.

14 Utilization

After use, the product must be disposed of in accordance with the disposal instructions in force in the treatment facility.

The class of danger of medical waste, in accordance with SanPiN 2.1.7.2790-10: B (dangerous,

enlarging

potentially infected medical waste).

15 Authorized Representative of Manufacturer in the territory of Russia

Company BW, LLC

129085, Moscow, Prospect Mira 101, bld. 1 space 17 (office 306a),

Tel: +7(499)281-67-68,

Website: beawire.com,

E-mail: info@beawire.com

The symbols used are as the following:

	Date of manufacture
	Use by
	Serial number
	Sterilized using ethylene oxide
	batch code
	Do not reuse
	Consult Instructions for Use
	Keep away from the sunlight
	Keep dry
	Catalogue number
	Do not use if package is damaged
	Do not resterilize
	Non-pyrogenic
Non-toxic	Non-toxic
	Manufacturer
	Authorised representative in the European Community
	Store at room temperature
	MR conditional


the list of national and international regulatory documents/standards, which matches medical device

- ISO 15223-1 Medical devices - Symbols to be used with medical product labels, labeling, and provide information - Part 1: General requirements
- EN 1041 Information supplied by the manufacturer with medical devices
- ISO 10555-1 and amendments: Sterile, disposable vascular catheters - Part 1: General requirements
- ISO 10993-1 + Assessment of biological effects of medical devices: Evaluation and testing
- ISO 10993-3 Biological Evaluation of medical treatments products: Tests for genotoxicity, Carcinogenicity and reproductive toxicity
- ISO 10993-4 Assessment of biological effects of medical devices: Selection of tests for interactions with blood
- ISO 10993-5 Biological Evaluation of medical treatments products: Test for cytotoxicity in vitro
- ISO 10993-6 Biological Evaluation of medical treatments products: test for the evaluation of local effects after implantation
- ISO 10993-7 Biological Evaluation of medical treatments products: ethylene oxide (EtO) Remains after sterilization
- ISO 10993-10 Assessment of biological effects of medical products: Test for irritation and sensitization skin
- ISO 10993-11 Biological Evaluation of medical treatments products: Tests for General toxicity
- ISO 10993-12 Assessment of biological effects of medical products: sample Preparation and reference materials
- ISO 11135-1 Test cycles sterilization with ethylene oxide
- ISO 13485 Medical devices - quality management System - Requirements for regulatory
- ISO 14155 Clinical investigation of medical devices for patient people - good clinical practice
- ISO 11607-1 Packaging for final sterilization of medical products - Part 1: Requirements for materials, sterile barrier and packaging systems
- ISO 14630 General requirements to the non-active surgical implants
- ISO 14971 Medical devices - Application of risk management to medical devices.
- ISO 25539-2 Cardiovascular implants - Endovascular devices - Part 2: Vascular stents
- MDD 93/42/EEC Basic requirements of the Directives on medical devices, as amended in Annex 1 2007/47/EEC

Manufacturer:

Lifetech Scientific (Shenzhen) Co., Ltd
Floor 1-5, Cybio Electronic Building,
Langshan 2nd Street, North Area of High-tech Park,
Nanshan District, Shenzhen 518057
P. R. China
Tel: +86 755 86026250
Fax: +86 755 86026251
Email: cs@lifetechmed.com

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