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2017/MAR/16 at 11:56 a.m. Doc number: M967612A001 Rev. 1A [global]
Printspec w - 3.5 x 3.0 inches


Medtronic

Melody™
Ensemble™
Ensemble™ II

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 **Caution:** Federal law (USA) restricts this
device to sale by or on the order of a physician.

Melody™ Transcatheter Pulmonary Valve • Katetrizačné implantovaná pulmonální chlopeč Melody™ • Melody™ transkateter lungeklap • Melody™ Transkatheter-Pulmonalklappe • Διακαθετηριακή πνευμονική βαλβίδα Melody™ • Válvula pulmonar transcatéter Melody™ • Katetrilla implantoitava Melody™-keuhkovahtimoläppä • Valve pulmonaire transcathéter Melody™ • Melody™ transzkatéteres pulmonalis billentyű • Valvola polmonare transcatetere Melody™ • Melody™-transkatheterpulmonalklep • Melody™ transkateter pulmonalklaff • Przezcewnikowa zastawka pnia płucnego Melody™ • Válvula pulmonar transcateter Melody™ • Клапан легочной артерии транскатетерный Melody™ • Pulmonálna chlopňa pre transkatétrové zavedenie Melody™ • Kateterburen pulmonalklaff Melody™ • Melody™ Transkateter Pulmoner Kapak

Ensemble™ Transcatheter Valve Delivery System • Systém pro katetrizační zavedení chlopeč Ensemble™ • Ensemble™ transkateter klapindføringsystem • Ensemble™ Transkatheter-Klappeneinführsystem • Διακαθετηριακό σύστημα προώθησης και τοποθέτησης βαλβίδας Ensemble™ • Sistema dispensador de la válvula transcatéter Ensemble™ • Katetrilla implantoitavan läpän Ensemble™-asetin • Système d'implantation pour valve transcathéter Ensemble™ • Ensemble™ transzkatéteres billentyűbevezető rendszer • Sistema di posizionamento per valvola transcatetere Ensemble™ • Ensemble™-plaatsingssysteem voor transkatheterkleppen • Ensemble™ transkateter klaffinnføringsystem • Przezcewnikowy system wprowadzający zastawki Ensemble™ • Sistema de colocação da válvula transcateter Ensemble™ • Устройство доставки клапана транскатетерного Ensemble™ • System pre

transkatétrové zavedenie chlopne Ensemble™ • Kateterburet
klaffinföringssystem Ensemble™ • Ensemble™ Transkateter Kapak Taşıma
Sistemi

Ensemble™ II Transcatheter Valve Delivery System • Systém pro katetrizační
zavedení chlopně Ensemble™ II • Ensemble™ II transkateter
klapinföringssystem • Ensemble™ II Transkatheter-Klappeneinführsystem •
Διακαθετηριακό σύστημα προώθησης και τοποθέτησης βαλβίδας Ensemble™
II • Sistema dispensador de la válvula transcatéter Ensemble™ II • Katetrilla
implantoitavan läpän Ensemble™ II -asetin • Système d'implantation pour valve
transcathéter Ensemble™ II • Ensemble™ II transzkatéteres billentyűbevezető
rendszer • Sistema di posizionamento per valvola transcatetere Ensemble™ II •
Ensemble™ II-plaatsingssysteem voor transkatheterkleppen • Ensemble™ II
transkateter klaffinnföringssystem • Przezcewnikowy system wprowadzający
zastawki Ensemble™ II • Sistema de colocação da válvula transcateter
Ensemble™ II • Устройство доставки клапана транскатетерного Ensemble™
II • Systém pre transkatétrové zavedenie chlopne Ensemble™ II • Kateterburet
klaffinföringssystem Ensemble™ II • Ensemble™ II Transkateter Kapak Taşıma
Sistemi

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1. Device description

The implant system consists of the Melody™ transcatheter pulmonary valve, model PB10 (stented bovine jugular vein valve) and either the Ensemble™ transcatheter valve delivery system, model NU10 or the Ensemble™ II transcatheter valve delivery system, model ENS10.

1.1. Melody™ transcatheter pulmonary valve (TPV)

The TPV consists of a heterologous (bovine) jugular vein valve sutured within a laser-welded, platinum-iridium stent with gold brazing of the welds. The TPV is available in a 16 mm bovine jugular vein (nominal length of 30 mm) and an 18 mm bovine jugular vein (nominal length of 28 mm). The 16 mm TPV (catalog number PB1016) can be expanded to ≤ 20 mm. The 18 mm TPV (catalog number PB1018) can be expanded to ≤ 22 mm. A final sterilization step is performed using a proprietary sterilant that contains 1% glutaraldehyde and 20% isopropyl alcohol, and in which the TPV is preserved and packaged until used. Adequate rinsing with isotonic saline solution must be performed before implantation to reduce the glutaraldehyde concentration.

1.2. Ensemble™ and Ensemble™ II transcatheter valve delivery system (delivery system)

Each delivery system consists of a balloon-in-balloon catheter with a retractable polytetrafluoroethylene (PTFE) sheath large enough to cover the TPV after crimping. The delivery system has a 22 Fr crossing profile. At inflation, the inner

20 Instructions for Use English

balloon is half the diameter of the outer balloon. Both balloons are made of nylon. The Ensemble™ delivery system comes in outer balloon sizes of 18 mm (catalog number NU1018), 20 mm (catalog number NU1020), and 22 mm (catalog number NU1022). The Ensemble™ II delivery system comes in outer balloon sizes of 18 mm (catalog number ENS1018), 20 mm (catalog number ENS1020), and 22 mm (catalog number ENS1022). The delivery system's sheath has a side-port used to flush the system and a hemostatic sleeve over the sheath to minimize bleeding at the insertion site. The delivery system has a polyether block amide (Pebax®, Arkema Corporation) distal obturator that is conical in shape. The delivery system is compatible with a 0.889 mm (0.035 in.) guidewire. The Ensemble™ II delivery system also includes radiopaque marker bands on the delivery system beneath the balloons. The marker bands are intended to aid in visibility of balloon location under fluoroscopy.

2. Indications

The Melody™ transcatheter pulmonary valve is indicated for use in patients with the following clinical conditions:

- patients with regurgitant prosthetic Right Ventricular Outflow Tract (RVOT) conduits with a clinical indication for invasive or surgical intervention, OR
- patients with stenotic prosthetic RVOT conduits where the risk of worsening regurgitation is a relative contraindication to balloon dilatation or stenting
- existence of a full (circumferential) RVOT conduit that was equal to or greater than 16 mm in diameter when originally implanted

Instructions for Use English 21

The intended lifetime for the Melody™ device is 2 years.

3. Contraindications

- venous anatomy unable to accommodate a 22 Fr size introducer sheath
- implantation in left heart
- unfavorable right ventricular outflow tract for good stent anchorage
- severe right ventricular outflow obstruction, which cannot be dilated by balloon
- obstruction of the central veins
- clinical or biological signs of infection
- active endocarditis
- known allergy to aspirin or heparin
- pregnancy

4. Warnings and precautions

Physician training: Physicians are required to have successfully completed the Physician/Center training program provided by Medtronic before using the Melody™ TPV and either the Ensemble™ or Ensemble™ II delivery system.

4.1. Warnings

4.1.1. TPV

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise

22 Instructions for Use English

the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Do not resterilize the TPV by any method.

Exposure of the device and container to irradiation, steam, ethylene oxide, or other chemical sterilants will render the device unfit for use.

Do not use the device if:

- the device has been dropped, damaged, or mishandled in any way
- the Use By date has elapsed
- each tamper-evident seal is broken
- the serial number tag on the TPV does not match the serial number on the container label
- the shipping temperature indicator window inside the shelf carton is black. If the shipping temperature indicator is black, the TPV is not suitable for clinical use.
- the storage solution does not completely cover the device or there is evidence of leakage

Do not expose the device to solutions other than the storage and rinsing solutions.

Do not add antibiotics to either the storage or the rinse solution. Do not apply antibiotics to the device.

Do not allow the device to dry. Maintain tissue moisture with irrigation or immersion.

Instructions for Use English 23

Do not attempt to repair a damaged device.

Do not handle or use forceps to manipulate the valve leaflet tissue.

Do not use forceps to manipulate the stent as this can damage the surface of the stent leading to stent fracture.

Do not overexpand the device beyond the maximum recommended size, which is 20 mm for Melody™ TPV catalog number PB1016 and 22 mm for Melody™ TPV catalog number PB1018, as this may result in a regurgitant TPV.

4.1.2. Delivery system

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Do not dilate the balloon appreciably beyond the diameter of the obstructive vessel and the intended implant site.

Do not use air or any gaseous substance as a balloon-inflation medium.

Do not advance the guidewire, balloon-dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.

Do not remove the guidewire from the delivery system at any time during the procedure.

4.2. Precautions

- Rinsing procedures of the TPV must be strictly followed to avoid necrosis of the tissue at the site of implantation.
- Exposure to glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to or breathing of the chemical vapor. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water for a minimum of 15 minutes. In the event of eye contact, flush with water for a minimum of 15 minutes and seek medical attention immediately.
- The sealed delivery system packaging should be inspected prior to opening. If the seal is broken or the packaging has been damaged, sterility cannot be assured.
- Proper functioning of the delivery system depends on its integrity. Use caution when handling the delivery system. Damage may result from kinking, stretching, or forceful wiping of the delivery system.
- This delivery system is not recommended for pressure measurement or delivery of fluids.
- Maintain tight catheter connections and use aspiration before proceeding to avoid air introduction into the system.
- The delivery system must be carefully flushed to avoid the introduction of air bubbles.
- Before crimping (reducing) the size of the TPV on the balloon, the orientation and function of the valve must be tested.

Instructions for Use English 25

- Do not take off the labels. These determine the orientation of the TPV before the crimping procedure on the balloon is completed.
- The crimping procedure must be carried out carefully. While crimping, the orientation of the valve must be known at all times. No change of orientation should occur as the TPV is mounted on the balloon.
- The crimping procedure must be carried out carefully. Do not place excessive pressure on the device during crimping as this can overstress the stent and lead to stent fracture.
- The inflation diameter of the balloon used during TPV delivery should approximate the diameter of the obstructive vessel and the intended implant site.
- Use of 2 appropriate size inflation devices (1 for each balloon) with pressure gauges is highly recommended during this procedure when inflating the balloon to deliver the TPV.
- The TPV is rigid and may make navigation through vessels difficult.
- Do not, under any circumstances, advance any portion of the delivery system against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- Balloon deployment should be conducted under fluoroscopic guidance with appropriate X-ray equipment.
- It is very important to ensure the balloons are completely deflated before pulling the delivery system back into the sheath.

- If resistance is felt upon removal, the delivery system and guidewire should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the guidewire and delivery system as a unit and withdrawing both together, using a gentle twisting motion combined with traction.

5. Patient information

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on each patient's condition.

5.1. Traceability labels

Product traceability labels are provided for hospital use only and are not to be returned to Medtronic.

5.2. Registration information

A patient registration form is included in TPV packages for countries that require patient registration.

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU.

After implantation, please complete all requested information. The serial number is located on both the package and the identification tag attached to the TPV.

Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

Medtronic will provide an Implanted Device Identification Card to the patient. The card contains the name and telephone number of the patient's physician as well as information that medical personnel would require in the event of an emergency. Patients should be encouraged to carry this card with them at all times.

6. How supplied

6.1. Packaging

The TPV is chemically sterilized and provided **sterile** and **nonpyrogenic** in a sealed glass container with a screw cap. Sterility is compromised if each tamper-evident seal is broken, the container is damaged, or leakage is evident. The outside of the container is **nonsterile** and should not be placed in the sterile field.

The delivery system is **nonpyrogenic** and is sterilized with ethylene oxide gas and packaged in a double pouch. The delivery system is sterile if the pouches are undamaged and unopened. Do not use the delivery system if the outer pouch is damaged. The delivery system should never be stored in only the inner pouch. The inner pouch does not provide a sterile barrier. The outer surfaces of the outer pouch are **nonsterile** and must not be placed in the sterile field.

6.2. Storage

Store the TPV at 15°C to 25°C (59°F to 77°F). Store the delivery system at room temperature and away from direct sunlight. The delivery system should never be stored in only the inner pouch. The inner pouch does not provide a sterile barrier. Appropriate inventory control should be maintained so that devices with the earliest Use By dates are preferentially used to avoid expiration. Devices must be used by the Use By date shown on the product labels.

7. Instructions for use

The following is a sequential outline of the catheterization/implant procedure. The type of diagnostic catheters, guidewires, dilation balloons, sizing balloons, or other tools needed will be at the discretion of the operator.

7.1. Access site preparation and preimplant diagnostics

1. Perform sterile preparation and draping of the femoral site.
2. Gain arterial and venous femoral access.
3. Administer heparin to achieve a target ACT of >250 seconds.
4. Introduce a catheter into the arterial sheath and advance it into the ascending aorta. Perform an aortogram to demonstrate the coronary arteries are not adjacent to the RVOT conduit, and that there is no risk of coronary compression when a bare metal stent or a TPV is implanted in the conduit.

If the coronaries appear to be in close proximity to the implant site, and coronary compression appears to be possible, further diagnostic investigations will be required at a later stage of the procedure.

5. Insert a catheter through the venous introducer and advance it to the right heart. Record pressure measurements through this catheter in the Right Ventricle (RV) and the Pulmonary Artery (PA). Record RV pressure simultaneously with the aortic pressure so that the systolic RV to aortic pressure ratio can be established.
6. Advance an angiographic catheter to the RV or proximal part of the RVOT conduit for angiography. The angiographic projections obtained will be based on the relative position of the conduit.
7. Obtain angiographic measurements of the intended implantation site to assess suitability of the conduit for TPV implantation. If the narrowest conduit dimension (usually in the lateral projection) is ≤ 18 mm and less than the nominal (implanted) conduit diameter, the site of conduit obstruction (the implantation site) should be predilated in order to facilitate optimal relief of conduit obstruction. The size of the predilation balloon should be at least 2 mm greater than the narrowest diameter of the conduit in any projection, no more than 110% of the nominal conduit diameter, and ≤ 20 mm for Melody™ TPV catalog number PB1016 and ≤ 22 mm for Melody™ TPV catalog number PB1018. Perform another conduit angiogram to ensure the conduit is intact, and measure the narrowest dimension.
8. Perform final conduit sizing using a balloon diameter that clearly demonstrates the narrowest portion of the conduit (at a pressure of

30 *Instructions for Use English*

<811 kPa [<8 atm]). If the prepared conduit waist is ≥ 14 mm and ≤ 20 mm, the RVOT conduit anatomy is considered suitable for Melody™ TPV catalog number PB1016 implantation. If the prepared conduit waist is ≥ 14 mm and ≤ 22 mm, the anatomy is considered suitable for Melody™ TPV catalog number PB1018 implantation.

If questions about the coronary anatomy remain, a selective coronary angiogram of the coronary artery close to the conduit should be carried out simultaneously with inflation of the sizing balloon. Should a coronary artery appear to be compressed by the balloon, the subject is considered to have unsuitable anatomy for the TPV, and the procedure should be abandoned.

9. Remove the diagnostic catheter, and place a guidewire across the RVOT with its tip located as far distal in the pulmonary arterial bed as possible. Gradually dilate the femoral vein up to 22 Fr. Leave the 22 Fr venous dilator in place while the TPV and the delivery system are being prepared.

7.2. Preparation of the TPV

Open the jar and carefully rinse the TPV as follows:

1. Before opening, carefully examine the jar and lid for damage, leakage, or broken seals. **The jar should contain enough sterilant to cover the TPV. The TPV should be rinsed for a minimum of 2 minutes to reduce the glutaraldehyde concentration from the TPV as directed in the following steps.**

- a. Using aseptic technique, prepare 3 sterile bowls; 1 remaining empty and 2 containing isotonic saline solution (500 mL) for rinsing.
 - b. Remove the TPV by grasping the serial number tag with atraumatic forceps and lifting it from the jar. The outside of the jar is **nonsterile**. Do not allow the TPV to come into contact with the outside of the jar.
 - c. A serial number tag is sutured to the outflow aspect of the TPV. Verify that the serial number on the tag matches the jar label serial number and the serial number on the patient registration form. If any differences are noted, the TPV should not be used. This tag should not be detached from the TPV until implantation is imminent.
 - d. Drain the residual storage solution from the valve into the empty discard bowl (bowl 1) by holding the TPV with the serial tag (outflow end) downward.
 - e. Transfer the empty TPV to the rinse bowl (bowl 2). Fill the TPV with rinse solution and alternately empty and fill by inverting and swirling, emptying and filling the valve for 1 minute, then empty the solution from the valve into the bowl.
 - f. Transfer the empty TPV to the rinse bowl (bowl 3) and repeat step e for a minimum of 1 minute. Leave the TPV in the rinse bowl until implantation is imminent to prevent the tissue from drying.
 - g. Empty the rinse solution from the TPV before loading the TPV onto the delivery system.
2. Remove the serial number tag by cutting the suture, which attaches the tag to the TPV.

7.3. Delivery system loading and placement of the TPV

1. Carefully examine the delivery system to confirm it was not damaged in shipment, and the balloon size is suitable for the intended procedure. The size of the delivery system to be used is based on the prepared conduit waist, as described in step 8 of Section 7.1, and should measure ≥ 14 mm and ≤ 22 mm. The delivery system used should contain a TPV delivery balloon that is 0 to 2 mm larger than the prepared conduit waist, but no more than 110% of the nominal conduit diameter.
2. Flush the delivery system guidewire lumen and the side port. Prepare the balloons by completely deflating with a syringe. Connect inflation syringes to the inner and outer balloon lumens. Inflation devices with pressure monitoring capabilities are recommended. The solution in the inflation syringes should be a one-third to two-thirds mixture of contrast to saline.
3. Reduce the size of the TPV while crimping it using mandrels of decreasing size. It is recommended to use a 2.5 mL syringe for crimping to an intermediate size prior to final crimping onto the balloon catheter. When reduced to the intermediate size, slide the TPV over the tip of the delivery system and center over the balloons.
4. Verify that the outflow end of the TPV is oriented toward the distal end of the delivery system. The outflow end of the TPV can be detected by the blue suture used to attach the valve at this location.

5. Gently crimp the TPV onto the balloon using finger pressure and a rolling action to exert equal pressure on all sides of the TPV to elongate the stent. Crimp only until no movement is felt on the delivery system.
Avoid bending or twisting the TPV.
Ensure the blue suture is adjacent to the blue tip of the delivery system to ensure proper orientation of the TPV for delivery.
6. Carefully slide the sheath over the TPV and balloons, ensuring that the crowns at the inflow side of the TPV do not get caught on the sheath as it is advanced over the TPV. Flush the sheath using the sidearm to remove air from the delivery system. Continue advancing the sheath and flushing until the sheath fits snugly over the proximal end of the blue tip.
Note: Check to ensure that the blue suture is adjacent to the blue tip of the delivery system.
7. Remove the venous dilator from the femoral vein, and advance the delivery system over the guidewire. Carefully introduce the delivery system through the skin and advance it toward the implantation site. Should bleeding be noted at the femoral venous site, the sleeve on the proximal end of the shaft can be advanced into the femoral vein to stop the bleeding. Advance the delivery system to the site where the TPV is intended to be deployed. This requires manipulations of the delivery system and the guidewire. **Pay attention to maintain adequate guidewire position at all times.**
8. Once the TPV has reached the level of the implantation site, uncover it. Hold the shaft of the delivery system in place while pulling back the outer sheath.

The TPV will be fully uncovered when the flushing port of the outer sheath is aligned with the proximal marker on the delivery system shaft. While uncovering, care must be taken that there is no inadvertent movement of the TPV, which could possibly lead to the loss of the appropriate TPV position.

9. Once the TPV is fully uncovered, connect a syringe with pure contrast to the flush port of the delivery system and inject a small amount of contrast to confirm good position of the TPV and full retraction of the sheath. Further minor adjustments of the position of the TPV are still possible at this point.
10. Inflate the inner balloon first, and then inflate the outer balloon. Once the TPV is properly deployed against the conduit wall, deflate the balloons and carefully remove the delivery system.
11. Repeat pressure measurements. If a significant RV-PA gradient persists, a balloon dilatation with a high-pressure balloon can be performed in order to minimize the gradient.
12. Repeat pressure measurements and once again compare the RV pressure with the systemic pressure measured through the femoral arterial approach.
13. Perform a contrast injection into the proximal main PA (or distal conduit) to demonstrate valve function and position. Take care to ensure that the valve is not held open by the guidewire and the delivery system as this would give the false impression of pulmonary regurgitation.
14. Remove the catheters, guidewires, and sheaths while maintaining hemostasis.

Note: The delivery system is for single use and should only be inflated once. The delivery system should not be used to expand the TPV after initial implantation. The TPV may be expanded after placement by using a different balloon catheter. Do not exceed the maximum recommended size of 20 mm for Melody™ TPV catalog number PB1016 and 22 mm for Melody™ TPV catalog number PB1018 to minimize the risk of creating regurgitation.

Note: Circumferential tearing of the delivery balloon catheter prior to complete expansion of the TPV may cause the balloon to become tethered to the TPV, requiring surgical removal. In case of rupture of an adequately sized balloon after TPV expansion, the balloon can be withdrawn and a new balloon catheter can be exchanged over a guidewire to complete expansion of the TPV.

Table 1. Melody™ TPV sizing chart

Delivery System Size – Inner Bal- loon / Outer Balloon	Inner Balloon Maximum Applied Pressure RBP		Outer Balloon Applied Pres- sures		Correspond- ing TPV Out- side Diame- ter (mm) (bal- loon inflated)
	atm	kPa	atm	kPa	
Size 18 mm Inner: 9 mm × 3.5 cm / Outer: 18 mm × 4 cm	5	506	1	101	17.9
			2	203	18.6
			3	304	19.4
			4 RBP	405	20.1
Size 20 mm Inner: 10 mm × 3.5 cm / Outer: 20 mm × 4 cm	5	506	1	101	19.7
			2	203	20.7
			3	304	21.7
			4 RBP	405	22.4
Size 22 mm Inner: 11 mm × 3.5 cm / Outer: 22 mm × 4 cm	4.5	456	1	101	21.8
			2	203	22.8
			3 RBP	304	24.1

Instructions for Use English 37

Delivery System Size – Inner Bal- loon / Outer Balloon	Inner Balloon Maximum Applied Pressure		Outer Balloon Applied Pres- sures		Correspond- ing TPV Out- side Diame- ter (mm) (bal- loon inflated)
	RBP				
	atm	kPa	atm	kPa	
Notes:					
1. Do not exceed bolded pressure values for either the inner or outer balloon of the delivery system size.					
2. RBP = rated burst pressure = maximum applied pressure					
3. atm = atmosphere					
4. kPa = kilopascal					

8. Potential complications/adverse events

The Melody™ TPV is a new device. Most of the reference literature addressing tissue or mechanical valved conduits describes hybrid devices composed of woven polyester tubing material and either a mechanical or biological valve sewn onto the lumen of the tube. The Melody™ TPV is a totally integral native valved conduit sewn into a stent. Some of the complications reported with currently commercially available conduits may not be entirely applicable.

Prosthetic heart valves have been associated with serious complications, sometimes leading to reoperation and/or death. This may also be the case with the Melody™ TPV. In addition, complications caused by immunogenic response

to the device, or to physical, chemical, or biological changes, may occur at undetermined intervals, and may require reoperation and replacement of the device.

General complications reported with valved conduits and biological tissue valves implanted in the heart include the following:

- hemorrhage
- bleeding diathesis due to anticoagulant use
- residual or increasing transvalvular gradients
- progressive neointimal thickening and peeling
- progressive stenosis and obstruction of the implant
- progressive pulmonary hypertension
- graft infection, endocarditis
- regurgitation through an incompetent valve
- hemolysis
- valve malfunction
- implant deterioration, physical or chemical
- pulmonary embolism
- thrombus

Types of tissue deterioration include:

- calcification, intrinsic and extrinsic
- infection
- thickening

- perforation
- myxomatous degeneration
- valve leaflet deterioration

Cardiac catheterization carries certain risks. Potential complications and related adverse events include, but are not limited to, the following:

- device embolization or migration
- laceration or rupture of blood vessels
- device misorientation or misplacement
- bleeding
- stent fracture
- loss of structural integrity due to stent fracture
- hematoma
- endocarditis
- arrhythmias
- thrombosis
- pain at the catheterization site
- distal thromboembolism
- cerebrovascular incident
- valve dysfunction
- allergic reaction to contrast agents
- paravalvular leak
- perforation of cardiac chambers

40 *Instructions for Use* English

- pseudoaneurysm
- rupture of the prosthetic RVOT conduit
- sepsis/infection
- heart failure
- fever
- compression of the coronary arteries
- death

These complications may present clinically with abnormal heart murmur, shortness of breath, exercise intolerance, dyspnea, orthopnea, anemia, fever, arrhythmia, hemorrhage, low cardiac output, pulmonary edema, congestive heart failure, and death.

9. Postimplantation information

9.1. Identification and management of stent fractures

The potential for stent fracture should be considered in all patients who undergo TPV placement, regardless of the previously discussed or subsequently characterized risk factors.

9.1.1. Identification

For patients found to have a substantial increase in the degree of RVOT obstruction, the possibility of an associated stent fracture should be considered and evaluated. Radiographic assessment of the stent with chest radiography or

Instructions for Use English 41

fluoroscopy should be included in the routine postprocedural evaluation of patients who receive a TPV.

9.1.2. Management

If a stent fracture is detected, continued monitoring of the stent should be performed in conjunction with clinically appropriate hemodynamic assessment. In patients with stent fracture and significant associated RVOT obstruction or regurgitation, reintervention should be considered in accordance with usual clinical practice.

Reintervention may include implantation of an additional Melody™ TPV or surgical conduit replacement. Note that limited data are available in the Melody™ US Clinical Study on reimplantation of another Melody™ TPV within the original Melody™ TPV.

9.2. Endocarditis

Endocarditis is a potential adverse event associated with all bioprosthetic valves (Section 8). Patients should make their health care providers aware that they have a bioprosthetic valve before any procedure.

A low incidence of suspected endocarditis has been reported in patients implanted with the Melody™ TPV. Unexplained, prolonged fever may be an indication of infection, and patients with these conditions should be advised to seek medical attention.

9.3. Magnetic resonance imaging (MRI) compatibility

Nonclinical testing has demonstrated that the Melody™ transcatheter pulmonary valve is MR (magnetic resonance) Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of ≤ 3.0 T
- Spatial gradient field of ≤ 2500 gauss/cm
- Normal operating mode only with a maximum whole body specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning

Based on nonclinical testing and modeling, the device was calculated to produce a maximum in vivo temperature rise of less than 2.0°C at 64 MHz (1.5 T) and less than 2.1°C at 128 MHz (3.0 T) for a whole body averaged SAR of 2.0 W/kg for 15 minutes of MR scanning.

MR image quality may be compromised if the area of interest is in the same area, or relatively close to the position of the device. It may be necessary to optimize MR imaging parameters for the presence of this implant. When tested at 3.0 T, the device exhibited slight distortion, with approximate extent of distortion from the device of 3 mm for spin echo images and 6 mm for the gradient echo images. The lumen of the device was obscured.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

10. Return of explanted bioprostheses

Medtronic is interested in obtaining recovered TPVs. Specific pathological studies of the explanted TPV will be conducted under the direction of a consulting pathologist. A written summary of the findings will be returned to the physician. To obtain a product return kit, contact a Medtronic distribution center or a Medtronic Sales Representative. If a kit is not available, place the explanted TPV in a container of glutaraldehyde or 10% buffered formalin immediately after excision. For further instructions on the return of an explanted device, contact a Medtronic Sales Representative.

11. Disclaimer of warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

ALTHOUGH THE MELODY™ TRANSCATHETER PULMONARY VALVE MODEL PB10, ENSEMBLE™ TRANSCATHETER VALVE DELIVERY SYSTEM MODEL NU10, AND ENSEMBLE™ II TRANSCATHETER VALVE DELIVERY SYSTEM MODEL ENS10, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC,

44 *Instructions for Use* English

THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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Declaration Letter

Medtronic Heart Valves Division, a corporation organized under the laws of the United States, having its principal place of business at 1851 East Deere Ave., Santa Ana, CA 92705, USA hereby confirm that the attached documentation:

1- Appendix to the IFU 27 pages;


Are true and unaltered copies of the original document.

Candace Johnson

Candace Johnson
Principal Regulatory Affairs
Specialist

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Fig. 1. Appearance of Ensemble delivery system with primary

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Fig. 2. Appearance of Ensemble II delivery system with primary
package

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Fig. 4. Delivery system schema (dimensions of different parts of device are provided in Fig 4A Fig 4B and Fig 4C)

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Fig. 5. Schema and dimensions of the delivery system luer hubs (dimensions are provided in cm)

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Fig. 6. Schema and dimensions of the delivery system tip (dimensions are provided in cm)

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Fig. 7. Schema and dimensions of the Ensemble delivery system balloons (dimensions are provided in cm)

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Fig. 8. Schema and dimensions of the Ensemble delivery system balloons (dimensions are provided in cm)

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
Fig. 9. main functional part of medical device

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Heat-sealing with a large number of thermoplastics have been PVC, styrene, PETG (PETG), HIPS; with polyethylene some kinds of polypropylene

6. Shelf life after sterilization: no more than 2 years

7. The materials are non-alkaline, non-toxic and odorless.

8. There are no holes, cracks, tears, wrinkles and refinements in the materials.

9. The purity level of the material: Acceptable.

10. Specific and minimum physical characteristics:

Tensile strength (tear resistance): $2.4 \text{ N} \pm 0.1 \text{ N}$

Stretching load (tensile strength): 50 MPa

The deviation in thickness: $\pm 0,05 \text{ mm}$

11. The value of pH: Tyvek material component has a neutral $\text{pH} = 7$

12. Adhesive characteristics:

Peel strength (strength of the adhesive layer): $\geq 3,3 \text{ Н}$

Adhesion transfer: $> 75\%$

Glue Density: $> 1.03 \text{ (g/cm}^3\text{)}$

13. Components of packages do not react, do not contaminate, do not pass, do not affect negatively on products before, during or after sterilization.

14. Width of the adhesive layer: $8,5 \text{ mm} \pm 0,5 \text{ mm}$

15. The opening of the packages by delamination is continuous and homogeneous, without delamination / rupture of the material.

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Fig. 13. marking on the Ensemble delivery system outer balloon hub (balloon size and hub functional designation is provided)

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Fig. 14. marking on the Ensemble delivery system sleeve (balloon size and hub is provided)

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Fig. 15. marking on the Ensemble II delivery system outer balloon hub (hub functional designation is provided)

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

Effect	severity of effect	severity rank
No effect	No effect on performance or safety of the device or process. Most unlikely to be noticed by user.	1
Very slight (pyformance)	Effect or defect only noticed by very discriminating user. Customer not annoyed. Will not influence device or process safety. Very slight effect on device performance. A portion (less than 100%) of the product may require minor rework.	2
Slight (performance)	Effect or defect noticed by average user. Fit & finish minor annoyance. Slight effect on device performance. Customer slightly annoyed. Will not influence device or process safety. The product may have to be sorted and	3

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© 2016 Medtronic, B.V.
Trademarks may be registered and are the property of their respective owners.