USER INSTRUCTIONS

SURGISPON Surgical hemostatic resorbable material (powder)

CAUTION: This product is restricted to sale by or on the order of a physician.

INDICATION

SURGISPON is used to stop bleeding in various surgical procedures.

DESCRIPTION

SURGISPON is a sterile, absorbable, water-insoluble, odorless product manufactured from gelatin intended for hemostatic use by applying to a bleeding surface.

SPECIFICATION

No.	Product ID	Form	Weight of an individual package	Diameter, mm	Error, mm
1	SSP-000001	powder	1 g	30 -70 μ m	$\pm$ 3%
2	SSP-000002	powder	2 g	30 -70 μ m	$\pm$ 3%

TRUE COPY

[Signature]
J. P. PATEL
NOTARY
GOVT. OF GUJARAT

[Signature]
TUSHAR PANDYA
ASST. SECRETARY SERVICE CENTER
GUJARAT CHAMBER OF COMMERCE & INDUSTRY
AHMEDABAD

 भारत सरकार GOVERNMENT OF INDIA
अपोस्टिल / APOSTILLE
(Convention de La Haye du 5 octobre 1961)
INDIA

This public document of the type
COMMERCIAL DOCUMENT

is issued to AGIS LIFESCIENCES

has been signed by TUSHAR PANDYA


with the seal, stamp or ASST SECRETARY, GUJARAT
CHAMBER OF COMMERCE & INDUSTRY, AHMEDABAD


Certified by
Section Officer(OI) MINISTRY OF EXTERNAL AFFAIRS
on 04-Sep-2018 at NEW DELHI INDIA

with reference to GJAH0028609818

Seal & Stamp

Signature
(DIPANKAR PAUL)
COMMANDER PAUL
Ministry of External Affairs
New Delhi

No.	Product ID	Image	Product shape
20	SSP-000001		Powder
21	SSP-000002		Powder

INDICATION/INTENDED USE

SURGISPON Absorbable Haemostatic Surgical Sponge Powder is effectively used in various surgical procedures for local haemostasis, when control of capillary, venous and arterial bleeding by pressure, ligature and other conventional methods is either ineffective or dangerous.

ACTION MECHANISM

Surgical hemostatic absorbable material SURGISPON has a porous structure and absorbs the blood that came into the wound from the damaged vessels as a result of surgical intervention. Activated platelets contained in the blood from the damaged blood vessels contain negatively charged phospholipids (PF₃) on the surface and emit a factor V, thus, the internal pathway of blood coagulation is activated. The VA factor binds to the formation of PF₃ together with the XA factor in the presence of Ca²⁺, a prothrombinase complex that converts prothrombin into thrombin, which stimulates the phospholipase of the platelet membrane, releasing arachidonic acid from its phospholipid layer. Arachidonic acid by platelet cyclooxygenase (COX-1) is converted into prostaglandin H₂ and G₂ (PGG₂ and PGH₂), of which the formation of thromboxane A₂ (TxA₂) and B₂ (TxB₂). These prostaglandins, TxA₂, together with the released platelet granule adenosine diphosphate (ADP) and platelet activation factor (FAT) initiate irreversible platelet aggregation inside the porous structure of the material SURGISPON, forming a combined clot.

BIODEGRADATION AND BIOCOMPATIBILITY

Gelatin, which is the material of surgical hemostatic absorbable SURGISPON, is a denatured collagen, homologous to human collagen. During the healing period postoperative polymorphonuclear leukocytes, fibroblasts, and macrophages secrete the enzyme collagenase, which cleaves the collagen molecule at the 1 position from the C-terminal end of the molecule. The resulting fragments are available to the action of other proteolytic enzymes and subsequently metabolized by the body naturally.

SURGISPON material is completely resorbed over a period of three to four weeks, without causing scarring. The product is biocompatible.

BIOCOMPATIBILITY TEST

Product	Test
SURGISPON Absorbable Haemostatic Gelatin Sponge	Cytotoxicity
	Skin irritation / Intracutaneous reaction

CONTRAINDICATIONS

SURGISPON Absorbable Haemostatic Surgical Sponge Powder is not recommended:

- For skin incisions when it may become a mechanical obstacle and prevent healing of skin edges.
- For patients with allergies to collagen
- In intravascular compartments because of the risk of embolization
- In conjunction with methyl methacrylate adhesives
- In conjunction with autologous blood salvage circles
- As a primary treatment of coagulation disorders
- For intraoperative haemostasis when the presence of infection has been confirmed
- For control of postpartum hemorrhage or menorrhagia
- For pulsating arterial bleeding control
- In large pools of blood or other liquids and when the source of bleeding is not clearly identifiable due to excessive liquid.

SIDE EFFECTS

- Hearing loss has been reported in association with the use of the product in tympanoplasty.
- The cases of local fibrosis requiring prolonged fixation and rehabilitation after severed tendon reconstruction have been reported.

WARNINGS


- In urologic surgical procedures the product should not be left in renal pelvis, calyces, bladder, urethra, ureter in order to minimize the risk of lithogenesis. The safety and efficacy of the gelatin sponge in urology has not been evaluated in a randomized clinical trial.
- The safety and efficacy of SURGISPON sponge powder in ophthalmology has not been evaluated.
- The safety and efficacy of SURGISPON sponge powder for children and pregnant women has not been evaluated.
- Although the safety and efficacy of the haemostatic absorbable gelatin product used in conjunction with other agents such as topical thrombine, antibiotic solutions and powders has not been evaluated in controlled clinical trials, this combination may be possible after careful study of indications of each agent.
- The haemostatic absorbable gelatin sponge should not be excessively pressed before wadding cavities or intertissue spaces. Application of excess amount of the product should be avoided, since swelling to original size may interfere with normal function and/or possibly result in compression syndrom of surrounding organs, tissues and structures.
- Although some surgical procedures may require tight packing of a wound, the haemostatic absorbable gelatin product should not be used in this manner unless the excess product not needed for haemostasis is removed after bleeding has been stopped.
- When the haemostatic absorbable gelatin product is used as a haemostasis agent in brain and maxillo-facial surgery in the immediate proximity to brain cavities, optic nerve and chiasm, other cranial and spinal nerves, the sponge and fragments thereof should always be removed after the bleeding has been stopped. Otherwise, this may result in unintentional compression of tissues, which can potentially lead to severed innervation and other neurological disorders.
- The haemostatic absorbable gelatin sponge powder is not intended as a substitute for meticulous surgical technique and proper application of ligatures or other conventional methods of haemostasis.
- The haemostatic absorbable gelatin-based product should be used with extreme caution in contaminated areas of the body. If signs of infection or abscess develop, reoperation may be necessary to remove the infected material and allow drainage.

DIRECTIONS FOR USE

Open the package immediately before use by pulling to sides as indicated on the package.

SURGISPON sponge powder is applied in a surgical procedure directly on a wound immediately after opening (when used dry) or it may be first placed in a plastic container and mixed with saline solution (not included).

To get the proper product consistency, mix 1 g of the sponge powder with 3-4 ml of sterile saline solution (for wet doughy consistency) or 4-5 ml of saline solution (for liquid paste consistency) in a plastic container; mix 2 g of the product with 6-8 ml of saline solution (for wet doughy consistency) or 8-10 ml of saline solution (for liquid paste consistency), respectively in a container. After that close the container with a lid and shake 3-4 times. Open

the container and use the product directly on a bleeding surface, applying it manually or with a sterile spatula (not included) or a syringe (not included).

To get a more liquid consistency ("mousse") add to a container with 1 g of the product 6-7 ml of sterile saline solution (not included); respectively, to a container with 2 g of the product 12-14 ml of sterile saline solution (not included) should be added. After that close the container with a lid and shake 3-4 times. Open the container and use the product directly on a bleeding surface, applying it manually or with a sterile spatula (not included) or a syringe (not included).

When mixed with saline solution SURGISPON sponge powder does not lose its haemostatic properties, but acquires plasticity which enable to apply it to hard-to-access areas.

SURGISPON applied to a bleeding spot in the form of paste or mousse should be held for 1-5 minutes on the bleeding site with manual pressure with a sterile pledget of cotton (not included) or a sterile gauze (not included) until the haemostatic effect has been achieved. A sterile gauze can be removed easier by wetting it with sterile saline solution.


The first correct application of SURGISPON usually arrests bleeding, but if not, additional applications may be needed. Use a new single piece or a fragment of SURGISPON prepared according to given instructions for each subsequent application.

PRECAUTIONS

- When SURGISPON is manually pressed against the bleeding site, special caution should be taken in order to avoid displacement of SURGISPON, which may have already transformed into a clot, which may result in recurrent intraoperative bleeding, when displaced.
- Use of suction applied over the pledget of cotton or gauze to draw blood into the SURGISPON is unnecessary, as the spongy texture will absorb sufficient blood by capillary action. Unnecessary mechanical suction may result in displacement of a clot and result in recurrent intraoperation bleeding.
- Use the minimum required amount of SURGISPON, which should be removed when the haemostatic effect has been achieved. After haemostasis the excess amount of gelatin sponge, if any, should be carefully removed to avoid the risk of subsequent displacement or compression of surrounding organs and tissues.
- The sponge is intended for single use and must not be sterilized!

SPECIAL INSTRUCTIONS

- The persons using SURGISPON should be well familiar with surgical techniques and methods of application of absorbable haemostatic gelatin sponge.
- When ordinary treatment is performed in healthcare institutions, the peculiarities of each patient should be taken into consideration.
- SURGISPON is supplied in protected sterile package. Sterility is assured until the package has been opened. It is recommended to use SURGISPON immediately

after opening the package and dispose of the unused amount of product according to applicable national regulations and standards of healthcare institutions.

- Opening of the package by pulling to the sides as indicated on the package should be made with proper care and in accordance with applicable aseptic practices.
- The product is intended for single use! The product must not be reused or resterilized.
- The product in the damaged protective package must not be used! In this case SURGISPON should be disposed of according to applicable national regulations and standards of healthcare institutions.

When used in adequate amount, SURGISPON is absorbed completely within 3-4 weeks.

HANDLING:

- SURGISPON is resistant to vibration load within 10-55 Hz with displacement amplitude of 0.15mm.
- Transportation package of SURGISPON protects the device against vibration load of 10-55 Hz with displacement amplitude of 0.35 mm and shock pulse with peak acceleration of 100 m/s² (10 g) and duration of shock acceleration 16 ms.
- SURGISPON in transportation package can be stored at +5 to +30°C.
- SURGISPON should not be subject to drastic temperature swings.

Operating conditions:

SURGISPON can be stored at +10 to +35°C and relative humidity of 30 to 90%.

The requirements for the containers of the product powder (SP000001 SSP-000002):

- The powder is delivered in a container with a volume of 5 m³, weighing 150 g.
- Containers are white, surface smooth.
- Integrity: hermetically sealed
- Mechanical strength: withstands 98 N compression force
- Resistance to hot water: when immersed in hot water temperature (70±5)°C-stable
- Durability of marking, stickers pasted on the container: must be at least two points in accordance with GOST 33756-2016
- Heat resistance: after holding the product for 2 hours at a temperature of (40±3) °C and subsequent exposure to normal conditions at a temperature of (22±4) °C for 30 min. the package is not deformed.
- Frost resistance: after soaking the product for 2 hours at a temperature of minus (25±2) °C and subsequent soaking at a temperature of (22±4) °C for 30 min. the package is not deformed.

Requirements for the container cap:

- appearance: screw cap without o-ring
- geometric dimensions: diameter 28 mm
- warping: 0.32%
- organoleptic indicator:
 - the smell of water extract - no more than 1 point
 - taste of water extract - not allowed;
 - changing the color and transparency of the water extract is not allowed.

COMPLIANCE WITH APPLICABLE REGULATIONS AND TECHNICAL STANDARDS

- SURGISPON is manufactured by Aegis Lifesciences in compliance with the Medical Device Directive 93/42/EEC.
- SURGISPON is manufactured according to the regulations of the Food and Drugs Administration (India). The production facilities undergo regular control and are certified in compliance with the regulations for manufacture and quality control of drugs (Manufacturing licence No. G/25/1809).
- Aegis Lifesciences has been certified according to ISO 13485 and GMP.
- Aegis Lifesciences has adopted the quality control system according to ISO 13485 (EN ISO 13485:2012/ISO 13485:2003) and the system of quality assurance for manufacturers of medical products in order to comply with the requirements of supervising authorities and ISO 9001:2008. The company regularly updates the quality control system following the applicable standards and requirements of regulation authorities.

WARRANTY

Manufacturer's authorized representative in the Russian Federation is BVS Ltd.
Legal address: 4th Tverskaya-Yamskaya St., 24, floor 1, placement 1, room 7, RMIL,
Moscow, RF, 125047
Phone: +7 499 553 03 53
Fax: +7 499 553 03 53
E-mail: info@b-v-s.net
Web: www.b-v-s.net

The quality of sponges is at all times guaranteed by the manufacturer's warranty.

SHELF LIFE


Shelf life is indicated by special symbols on the label. Shelf life is a period between the manufacturing date (manufacturer's mark) and expiration date (hourglass symbol). The length of time between the two dates is when the device is fit for use. The shelf life of SURGISPON is 3 (three) years.

STERILIZATION

SURGISPON is supplied sterile in a hermetically sealed envelope and gamma sterilized. It has been established that gamma rays have no impact on the product or package. The gamma radiation sterilization equipment conforms to ISO 11137-1. To ensure maximum sterility, the sterility assurance level applied is 10^{-6} . SURGISPON is intended for single use. Sterility of the product is assured within the entire shelf life period unless the package has been damaged. Do not use SURGISPON sponge powder if there is any doubt as to its sterility!

STORAGE

SURGISPON in protective package should be stored in dark, dry and cool places; avoid direct sunlight. SURGISPON should not be stored after the original package has been opened!

SAFE USE

SURGISPON should be used only in healthcare institutions and by medical professionals trained to perform surgical procedures and apply absorbable haemostatic gelatin products.


PACKAGE

SURGISPON product is supplied in primary protective hermetical envelope which assures sterility and protects the product from mechanical and environmental impacts during the entire shelf life. The envelope contains marking detailing the information on the product. The packaging process is fully automated and validated which ensures steady high quality.


Package marking is made according to ISO 15223-1. The symbols are printed on primary and secondary packaging. The symbols provide information on the product, possibility to reuse the product, sterility, recommendation not to use the product if package is damaged, storage and transportation, shelf life and manufacturer. The meaning of symbols is explained in user instructions.

GRAPHIC SYMBOLS ON THE PROTECTIVE ENVELOPE

	Caution! See user instructions
	See user instructions
	Temperature range
	Keep dry
	Keep away from sunlight
	Single use
	Do not re-sterilize
	Do not use if package is damaged
	Date of manufacture: month and year
	Use by: month and year
	Batch number
	Reference number
	Gamma sterilized
	CE-Mark and identification number of the notified body. The product complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC
	Manufactured by
	Authorized representative in the European Community

SURGISPON products with IDs SSP-000001, SSP-000002 include the following on the label:

- CE-Mark and identification number of the notified body. The product complies with the essential health and safety requirements set out in the Medical Device Directive
- Product ID
- Product name, weight
- Keep dry
- Do not re-sterilize
- Do not use if package is damaged
- Temperature range 5-30°C
- Keep away from sunlight
- See user instructions
- Single use
- Gamma Sterilized
- Caution! See user instructions
- Keep dry
- Batch number
- Manufacturing date: month and year
- Expiration date: month and year
- Manufacturing licence number
- Manufactured by
- Authorized representative in the EU

PACKAGE

SURGISPON products in primary envelopes are packaged into carton boxes (secondary package). Each carton box contain one or several sponges. Each carton box contains an enclosed package insert with information on the product and a user instruction.

The number of secondarily Packed products placed in corrugated cardboard boxes (tertiary packaging), as well as the gross weight of corrugated cardboard boxes depend on the individual requests of the quantities of products by the customer.

PACKAGE (carton box)

Product ID	Weight (g)	Height (cm)	Width (cm)	Length (cm)	Error (mm)
SSP-000001	55	17	8.5	24	±2
SSP-000002	14	7.5	6.5	6.5	±2

PACKAGE

SURGISPON products in primary envelopes are packaged into carton boxes. Each carton box contains several pieces. Each carton box contains an enclosed package insert with information on the product and a user instruction.


The secondary package is a standard carton box protecting the product from mechanical and environmental impacts during the entire shelf life. The secondary package protects the product during transportation and storage and contains marking with the product information.

Product ID	Packing	Primary packaging type
SSP-000001	6 pcs	Aluminum bag
SSP-000002	1 pcs	Aluminum bag

The secondary package (carton box) includes the following data:

- CE-Mark and identification number of the notified body. The product complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC 93/42/EEC
- Product ID
- Product name, dimensions
- Consult user instruction
- Do not re-sterilize
- Do not use if package is damaged
- Temperature range 5-30°C
- Keep away from sunlight
- Do not re-use
- Keep dry
- Gamma Sterilized
- Batch number
- Manufacturing date: month and year
- Expiration date: month and year
- Manufacturing licence number
- Manufactured by
- Open here
- Number of pieces in a box
- Barcode
- Authorised representative in the EC

SURGISPON products in secondary packing are further placed into big corrugated boxes (tertiary packing) and sealed across the sides and along the perimeter with control tape. The packing ensures protection against mechanical and environmental impacts during transportation and storage as well as enables to efficiently use the capacity of transportation vehicles and perform loading operations.

TRANSPORTATION

When arranging a transportation of the product the type of transportation vehicle, point of destination, transportation provider as well as periodicity and eventual overloading of vehicles should be taken into consideration.

DISPOSAL

Discarded SURGISPON absorbable haemostatic products should be disposed of according to national regulations and standards of healthcare institutions.



For the EU countries the Directive 93/42 EEC should be applied.

For the Russian Federation the following national regulations should be used: Sanitary regulations and Standards SanPiN 2.1.7.72-99 "Collection, storage and disposal of waste in healthcare institutions"; SanPiN 2.1.7.1322-03 "Sanitary rules for disposal and treatment of industrial and consumer waste"; Federal Law "On industrial and consumer waste" No. 89 dd June 24, 1998.

Excessive SURGISPON product removed from a wound should be discarded as Class B medical waste.


Unused and expired SURGISPON product or SURGISPON in damaged package should be discarded as Class D medical waste.

PREPARED AND APPROVED BY:

	Full name and title	Date	Signature
Reviewed / Prepared by	Bhavin Trivedi (State Registration Services Expert)	18/08/2018	
Controlled by	Deepak Patel (Representative of the Board of Directors)	18/08/2018	
Approved by	Neelam Gupta (Managing Partner)	18/08/2018	

TUSHAR PANDYA
ASST. SECRETARY SERVICE CENTER
GUJARAT CHAMBER OF COMMERCE & INDUSTRY
AHMEDABAD

TRUE COPY

J. P. PATEL
NOTARY
GOVT. OF GUJARAT

भारत सरकार GOVERNMENT OF INDIA
अपोस्टिल / APOSTILLE
(Convention de La Haye du 5 octobre 1961)
INDIA

 Country: India

This public document of the type
COMMERCIAL DOCUMENT

is issued to AGIS LIFESCIENCES

has been signed by TUSHAR PANDYA

with their seal / stamp of ASSTT. SECRETARY, GUJARAT
CHAMBER OF COMMERCE & INDUSTRY, AHMEDABAD

Certified by
Section Officer(OI) MINISTRY OF EXTERNAL AFFAIRS
on 04-Sep-2018 at NEW DELHI, INDIA

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USER INSTRUCTIONS

SURGISPON Surgical hemostatic resorbable material

CAUTION: This product is restricted to sale to or on the order of a physician.

INDICATION

SURGISPON sponge is used to stop bleeding in various surgical procedures.

DESCRIPTION

SURGISPON is a **sterile, absorbable, water-insoluble, odorless** sponge manufactured from **gelatin** intended for haemostatic use by applying to a **bleeding surface**.

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is issued to **AGIS LIFESCIENCES**


has been signed by **TUSHAR PANDYA**

with the seal / stamp of **ASSTT. SECRETARY, GUJARAT**
CHAMBER OF COMMERCE & INDUSTRY, AHMEDABAD


Certified by
Section Officer(OI) MINISTRY OF EXTERNAL AFFAIRS
 on **04-Sep-2018** at **NEW DELHI, INDIA**

With reference no. **GJAH0028610718**

Seal / Stamp

 Signature









TUSHAR PANDYA
 ASSTT. SECRETARY, GUJARAT
 CHAMBER OF COMMERCE & INDUSTRY
 AHMEDABAD







No.	Product ID	Length, mm	Width, mm	Height, mm	Diameter, mm	Error, mm
17	SSP-805001	80	50	1	-	±5%
18	SSP-200705	20	7	5	-	±5%
19	SSP-1208	12	-	-	8	±5%
20	SSP-253005	25	30	5	-	±5%
21	SSP-954805	95	48	5	-	±5%
22	SSP-484805	48	48	5	-	±5%

No.	Product ID	Image	Product shape
1	SSP-805010		Rectangular slab
2	SSP-705010		Rectangular slab
3	SSP-8030		Cylindrical tampon
4	SSP-101010		Cube







No.	Product ID	Image	Product shape
5	SSP-202007		Quadrate slab
6	SSP-801507		Rectangular strip
7	SSP-8007		Cylindrical tampon
8	SSP-303010		Quadrate slab
9	SSP - 954805		Rectangular slab



No.	Product ID	Image	Product shape
10	SSP-804010		Rectangular slab
11	SSP-802010		Rectangular strip
12	SSP-802507		Rectangular strip
13	SSP-602007		Rectangular strip
14	SSP - 484805		Quadrate slab

No.	Product ID	Image	Product shape
15	SSP-100		Rectangular slab
16	SSP-120		Rectangular slab
17	SSP-200		Rectangular slab
18	SSP-705001		Rectangular slab
19	SSP-805001		Rectangular slab


No.	Product ID	Image	Product shape
20	SSP-200705		Rectangular slab
21	SSP -1208		Cylindrical tampon
22	SSP - 253005		Rectangular slab

INDICATION/INTENDED USE

Surgical hemostatic absorbable material SURGISPON is effectively used in various surgical procedures for local haemostasis, when control of capillary, venous and arterial bleeding by pressure, ligature and other conventional methods is either ineffective or dangerous.

ACTION MECHANISM

Surgical hemostatic absorbable material SURGISPON has a porous structure and absorbs the blood that came into the wound from the damaged vessels as a result of surgical intervention. Activated platelets contained in the blood from the damaged blood vessels contain negatively charged phospholipids (PF3) on the surface and emit a factor V. thus, the internal pathway of blood coagulation is activated. The VA factor binds to the formation of PF3 together with the XA factor in the presence of Ca^{2+} , a prothrombinase complex that converts prothrombin into thrombin, which stimulates the phospholipase of the platelet membrane, releasing arachidonic acid from its phospholipid layer. Arachidonic acid by platelet cyclooxygenase (COX-1) is converted into prostaglandin H 2 and G 2 (PGG 2 and PGH 2), of which the formation of

thromboxane A₂ (TxA₂) and B₂ (TxB₂). These prostaglandins, TxA₂, together with the released platelet granule adenosine diphosphate (ADP) and platelet activation factor (FAT) initiate irreversible platelet aggregation inside the porous structure of the material SURGISPON, forming a combined clot.

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
CONTRAINDICATIONS

SURGISPON Absorbable Haemostatic Surgical Sponge is not recommended:

- for skin incisions when it may become a mechanical obstacle and prevent healing of skin edges.
- for patients with allergies to collagen
- in intravascular compartments because of the risk of embolization
- in conjunction with methyl methacrylate adhesives
- in conjunction with autologous blood salvage circles
- as a primary treatment of coagulation disorders
- for intraoperative haemostasis when the presence of infection has been confirmed
- for control of postpartum hemorrhage or menorrhagia
- for control of pulsating arterial bleeding
- in large pools of blood or other liquids and when the source of bleeding is not clearly identifiable due to excessive liquid.

SIDE EFFECTS:

- Hearing loss has been reported in association with the use of the sponge in tympanoplasty.
- The cases of local fibrosis requiring prolonged fixation and rehabilitation after severed tendon reconstruction have been reported.

WARNINGS:

- In urologic surgical procedures the sponge should not be left in renal pelvis, calyces, bladder, urethra, ureter in order to minimize the risk of lithogenesis. The safety and efficacy of the gelatin sponge in urology has not been evaluated in a randomized clinical trial.
- The safety and efficacy of the gelatin sponge in ophthalmology has not been evaluated.
- The safety and efficacy of the gelatin sponge for children and pregnant women has not been evaluated.
- Although the safety and efficacy of the haemostatic absorbable gelatin sponge used in conjunction with other agents such as topical thrombine, antibiotic solutions and powders has not been evaluated in controlled clinical trials; this combination may be possible after careful study of indications of each agent.
- The haemostatic absorbable gelatin sponge should not be excessively pressed before wadding cavities or intertissue spaces. Application of excess amount of the sponge should be avoided, since swelling to original size may interfere with normal function and/or possibly result in compression syndrom of surrounding organs, tissues and structures.
- Although some surgical procedures may require tight packing of a wound, the haemostatic absorbable gelatin device should not be used in this manner unless excess product not needed for haemostasis is removed after bleeding has been stopped.
- When the haemostatic absorbable gelatin device is used as a haemostasis agent in brain and maxillo-facial surgery in the immediate proximity to brain cavities, optic nerve and chiasm, other cranial and spinal nerves, the sponge and fragments thereof should always be removed after the bleeding has been stopped. Otherwise, this may result in unintentional compression of tissues, which can potentially lead to severed innervation and other neurological disorders.

- The haemostatic absorbable gelatin sponge is not intended as a substitute for meticulous surgical technique and proper application of ligatures or other conventional methods of haemostasis.
- The haemostatic absorbable gelatin sponge should be used with extreme caution in contaminated areas of the body. If signs of infection or abscess develop, reoperation may be necessary to remove the infected material and allow drainage.

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SURGISPON Sponge can be used either dry or saturated with sterile isotonic solution (not included) as a single piece or cut to fragments of desired size.

When used dry a single piece or a fragment of SURGISPON should be manually adhered to the bleeding site and held in place until the haemostatic effect is achieved.

When used with sterile isotonic solution, a single piece or a fragment of SURGISPON should be first immersed in the solution for 5-10 sec and then withdrawn, squeezed between gloved fingers (gloves not included) to expel air bubbles, and then, if needed, replaced in the solution for complete saturation. SURGISPON should promptly return to its original size (small increase in size due to saturation is allowed). If it has not returned to its original shape, repeat the procedure. When the saturated SURGISPON is used, it can be wetted with a gauze (not included) before applying to the bleeding spot.

SURGISPON should be held on the bleeding site with manual pressure with a sterile pledget of cotton (not included) or a sterile gauze (not included) until the haemostatic effect has been achieved. A sterile gauze can be removed easier by wetting it with sterile saline solution.

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PRECAUTIONS:

- When SURGISPON is manually pressed against the bleeding site, special caution should be taken in order to avoid displacement of SURGISPON, which may have already transformed into a clot, which may result in recurrent intraoperative bleeding, when displaced.
- Use of suction applied over the pledget of cotton or gauze to draw blood into the SURGISPON is unnecessary, as the spongy texture will absorb sufficient blood by capillary action. Unnecessary mechanical suction may result in displacement of a clot and result in recurrent intraoperative bleeding.
- Use the minimum required amount of SURGISPON which should be removed when the haemostatic effect has been achieved. After haemostasis the excess amount of

gelatin sponge, if any, should be carefully removed to avoid the risk of subsequent displacement or compression of surrounding organs and tissues.

- The sponge is intended for single use and must not be sterilized!

SPECIAL INSTRUCTIONS

- The persons using SURGISPON should be well familiar with surgical techniques and methods of application of absorbable haemostatic gelatin sponge.
- When ordinary treatment is performed in healthcare institutions, the peculiarities of each patient should be taken into consideration.
- SURGISPON is supplied in protected sterile package. Sterility is assured until the package has been opened. It is recommended to use SURGISPON immediately after opening the package. The unused amount of product should be disposed of according to applicable national regulations and standards of healthcare institutions.
- Opening of the package by pulling to the sides as indicated on the package should be made with proper care and in accordance with applicable aseptic practices.
- The product is intended for single use! The product must not be reused or resterilized.
- The product in the damaged protective package must not be used! In this case SURGISPON should be disposed of according to applicable national regulations and standards of healthcare institutions.

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
- SURGISPON is manufactured by Aegis Lifesciences in compliance with the Medical Device Directive 93/42/EEC.
- SURGISPON is manufactured according to the regulations of the Food and Drugs Administration (India). The production facilities undergo regular control and are

certified in compliance with the regulations for manufacture and quality control of drugs (Manufacturing licence No. G/25/1809).

- Aegis Lifesciences has been certified according to ISO 13485 and GMP.
- Aegis Lifesciences has adopted the quality control system according to ISO 13485 (EN ISO 13485:2012/ISO 13485:2003) and the system of quality assurance for manufacturers of medical products in order to comply with the requirements of supervising authorities and ISO 9001:2008. The company regularly updates the quality control system following the applicable standards and requirements of regulation authorities.

WARRANTY

Manufacturer's authorized representative in the Russian Federation is BVS Ltd.
Legal address: 4th Tverskaya-Yamskaya St., 24, floor 1, placement 1, room 7, RM1L,
Moscow, RF, 125047
Phone: +7 499 553 03 53
Fax: +7 499 553 03 53
E-mail: info@b-v-s.net
Web: www.b-v-s.net

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STERILIZATION

SURGISPON Sponge is supplied sterile in a hermetically sealed envelope. The sponge is gamma sterilized. It has been established that gamma rays have no impact on the product or package. The gamma radiation sterilization equipment conforms to ISO 11137-1. To ensure maximum sterility, the sterility assurance level applied is 10^{-6} . SURGISPON is intended for single use. The product sterility is assured within the entire shelf life period unless the package has been damaged. Do not use the product if there is any doubt as to its sterility!

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SAFE USE

SURGISPON Absorbable Haemostatic Gelatin Sponge should be used only in healthcare institutions and by medical staff trained to perform surgical procedures and apply absorbable haemostatic gelatin products.


PACKAGE


SURGISPON Sponge is supplied in primary protective hermetical envelope which assures sterility and protects the product from mechanical and environmental impacts during the entire shelf life. The envelope contains marking detailing the information on the product. The packaging process is fully automated and validated which ensures steady high quality.

GRAPHIC SYMBOLS ON THE PROTECTIVE EVELOPE

	Attention! See instructions for use
	Consult instructions for use
	Temperature limitations
	Keep dry
	Keep away from sunlight
	Single use
	Do not re-sterilize
	Do not use if package is damaged
	Date of manufacture: month and year
	Use by: month and year
	Batch number
	Reference number
	Gamma Sterilized
	CE-Mark and identification number of the notified body. The device complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC
	Manufactured by
	Authorized representative in the European Community


SURGISPON products with IDs SSP-805010, SSP-705010, SSP-8007, SSP-802507, SSP-602007, SSP-8030, SSP-100, SSP-120, SSP-705001, SSP-805001, SSP-200705, SSP-200, SSP-954805 include the following data on the label:

- CE-Mark and identification number of the notified body. The product complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC
- Product ID
- Product name, dimensions
- Keep dry
- Attention! See user instructions enclosed
- Do not re-sterilize
- Do not use if package is damaged
- Temperature range 5-30°C
- Keep away from sunlight
- Single use
- See user instructions enclosed
- Gamma Sterilized
- Batch number
- Manufacturing date: month and year
- Expiration date: month and year
- Manufacturing licence number
- Manufactured by
- Barcode
- Peel here


SURGISPON products with IDs SSP-801507, SSP-303010, SSP-802010, SSP-804010, SSP-253005, SSP-484805, SSP-202007 include the following data on the label:

- CE-Mark and identification number of the notified body. The product complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC 93/42/EEC
- Product ID
- Product name, dimensions
- Keep dry
- Single use
- Do not use if package is damaged
- Attention! See user instructions enclosed
- Gamma Sterilized
- Batch number
- Manufacturing date: month and year
- Expiration date: month and year
- Manufacturing licence number
- Manufactured by
- Peel here

SURGISPON products with IDs SSP-101010, SSP-1208 include the following data on the label:

- CE-Mark and identification number of the notified body. The device complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC 93/42/EEC
- Product ID
- Attention! See user instructions enclosed
- Single use
- Keep dry
- Gamma Sterilized
- Temperature range 5-30°C
- Batch number
- Expiration date: month and year
- Manufacturing licence number

PACKAGE

SURGISPON sponges in primary envelopes are packaged into carton boxes. Each carton box contain one or several sponges. Each carton box contains an enclosed package insert with information on the product and a user instruction.

PACKAGE (carton box)

Product ID	Weight (g)	Height (cm)	Width (cm)	Length (cm)	Error (mm)
SSP-805010	25	8	12	12	±2
SSP-705010	25	8	12	12	±2
SSP-8030	24	23.5	10.5	3,5	±2
SSP-101010	14	15.5	9.5	2,5	±2
SSP-202007	25	8	12	12	±2
SSP-801507	25	8	12	12	±2
SSP-8007	25	8	12	12	±2
SSP-303010	25	8	12	12	±2
SSP-804010	25	8	12	12	±2
SSP-802010	25	8	12	12	±2
SSP-802507	25	8	12	12	±2
SSP-602007	16.5	8	12	12	±2
SSP-100	28	24.5	15.5	2,5	±2
SSP-170	28	24.5	15.5	2,5	±2
SSP-200	28	28	19	2,5	±2
SSP-705001	16.5	16	12	2,5	±2
SSP-805001	16.5	16	12	2,5	±2
SSP-200705	16.5	16	12	2,5	±2
SSP-1208	14	15.5	9.5	2,5	±2
SSP-253005	25	8	12	12	±2
SSP-954805	28	24.5	15.5	2,5	±2
SSP-484805	25	8	12	12	±2

PACKAGE


SURGISPON sponges in primary envelopes are packaged into carton boxes. Each carton box contains one or several sponges. Each carton box contains an enclosed package insert with product information and user instruction.

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Product ID	Number of SURGISPON products in one box	Primary packaging tipe
SSP 805010	10 pcs	Blister tray, covered with medical paper
SSP 705010	10 pcs	Blister tray, covered with medical paper
SSP 801507	10 pcs	Blister tray, covered with medical paper
SSP 8007	10 pcs	Blister tray, covered with medical paper
SSP 303010	10 pcs	Blister tray, covered with medical paper
SSP 804010	10 pcs	Blister tray, covered with medical paper
SSP 802010	10 pcs	Blister tray, covered with medical paper
SSP 802507	10 pcs	Blister tray, covered with medical paper
SSP 602007	10 pcs	Blister tray, covered with medical paper
SSP 253005	10 pcs	Blister tray, covered with medical paper
SSP 484805	10 pcs	Blister tray, covered with medical paper
SSP 8030	6 pcs	Blister tray, covered with medical paper
SSP 101010	32 pcs	Blister tray, covered with medical paper
SSP 1208	32 pcs	Blister tray, covered with medical paper
SSP 100	2 pcs	Blister tray, covered with medical paper
SSP 120	2 pcs	Blister tray, covered with medical paper
SSP 954805	4 pcs	Blister tray, covered with medical paper
SSP 705001	20 pcs	Blister tray, covered with medical paper
SSP 805001	20 pcs	Blister tray, covered with medical paper
SSP 200705	10 pcs	Blister tray, covered with medical paper
SSP 202007	10 pcs	Blister tray, covered with medical paper
SSP 200	2 pcs	Blister tray, covered with medical paper

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- CE-Mark and identification number of the notified body. The device complies with the essential health and safety requirements set out in the Medical Device Directive 93/42/EEC 93/42/EEC
- Product ID
- Product name, dimensions
- See user instructions enclosed
- Do not re-sterilize
- Do not use if package is damaged
- Temperature range 5-30°C
- Keep away from sunlight
- Do not re-use
- Keep dry
- Gamma Sterilized
- Batch number
- Manufacturing date: month and year
- Expiration date: month and year
- Manufacturing licence number
- Manufactured by
- Open here
- Number of pieces in a box
- Barcode
- Authorized representative in the EC

SURGISPON sponges in secondary packing are further placed into big corrugated boxes (tertiary packing) and sealed across the sides and along the perimeter with control tape. The packing ensures protection against mechanical and environmental impacts during transportation and storage as well as enables to efficiently use the capacity of transportation vehicles and perform loading operations.

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DISPOSAL

Discarded SURGISPON absorbable haemostatic sponges should be disposed of according to national regulations and standards of healthcare institutions.

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

Excessive SURGISPON sponge removed from a wound should be discarded as Class B medical waste.
Unused and expired SURGISPON product or SURGISPON in damaged package should be discarded as Class D medical waste.

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№	Номер по каталогу	Изображение	Форма изделия
15	SSP-100	 The image shows the packaging for Surgispon SSP-100, including a box and a manual, with a white rectangular sponge pad placed in front.	Прямоугольная пластина
16	SSP-120	 The image shows the packaging for Surgispon SSP-120, including a box and a manual, with a white rectangular sponge pad placed in front.	Прямоугольная пластина
17	SSP-200	 The image shows the packaging for Surgispon SSP-200, including a box and a manual, with a white rectangular sponge pad placed in front.	Прямоугольная пластина
18	SSP-705001	 The image shows the packaging for Surgispon SSP-705001, including a box and a manual, with two white rectangular sponge pads placed in front.	Прямоугольная пластина
19	SSP-805001	 The image shows the packaging for Surgispon SSP-805001, including a box and a manual, with two white rectangular sponge pads placed in front.	Прямоугольная пластина

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