
APPROVED BY


Silvana BROWN
Sr Director, International Regulatory Affairs

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M. DAVID MALEDON

(seen exclusively to certify the above signature)

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USER MANUAL

MAGIC Microcatheter for Injection and for Occlusion

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This user manual has been prepared to provide information for purposes of registration of a medical device in the territory of the Russian Federation. The information is provided in the amount required in accordance with the Rules for Registration of Medical Devices (Russian Government Decree No. 1416 dated 27.12.2012).

1. Name of the device:

MAGIC microcatheter for injection or occlusion

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Development, production and sale of the device complies with the requirements of EN ISO 13485:2012 + AC:2012 (Certificate No. 513975 MP2012, effective date: 06.02.2018, expiry date: 31.03.2019, notified body: DQS Medizinprodukte GmbH); the device is marked with CE (EC certificate No. 513975 MR2 annex II excluding section 4, effective date: 12.05.2018, expiry date: 31.07.2022, EC certificate No. 534993 MRA annex II section 4, effective date: 01.08.2017, expiry date: 31.07.2022, notified body: DQS Medizinprodukte GmbH).

Device code according to the All-Russian classifier of products OK-005-93: 94 3630;

OKPD2 - 32.50.13.110 - Syringes, needles, catheters, cannulas and similar instruments;

Class depending on the degree of potential risk of medical use in accordance with the medical devices nomenclature classifier approved by Order of the Ministry of Health of the Russian Federation No. 4n-3 dated June 6, 2012;

Device classification in accordance with the medical devices nomenclature classifier approved by Order of the Ministry of Health of the Russian Federation No. 4n-131460 dated June 6, 2012;

Type of contact with the human body: short-term (less than 24 hours) contact with circulating blood.

This device is sterile. For single use.

Conditions of use: medical and medical-preventive health care institutions, only as prescribed by a medical specialist.

2. Purpose

The device is intended for use in peripheral and intracranial vessels, in selective and hyper-selective vascular catheterization, for diagnostic or therapeutic purposes. A sterile flexible tube of small diameter is designed for controlled injection of substances/fluids.

3. Device special properties

MAGIC catheters are designed for general intravascular use. They can be used for controlled selective local infusion of drug products or embolizing materials into vessels. MAGIC catheters are designed to facilitate access through the distal winding vascular tree. Progressive elasticity in the range from superelastic distal shaft to rigid proximal shaft allows the physician to

advance the catheter. A rigid proximal shaft enables to control the torque in order to facilitate the advancement of the catheter. The tip of MAGIC catheter (ring) and the shaft are radiopaque to facilitate radiosopic imaging. The external surface of MAGIC microcatheter is hydrophilic. Catheter 1.2F has a fourth distal part of the distiller. Diameter of MAGIC microcatheter is reduced from the proximal part with a diameter of 2.7F, the diameter of which is 1.2F, 1.5F or 1.8F, depending on the catheter size. Technical characteristics of the devices are given in clause 4. Each device is packed in an individual package, each individual package is packed in a cardboard box, which also contains user manual. For more information on packaging and labeling, see clause 5. Appearance of the device is shown in Figures 1-4.


Figure 1: Appearance of the device without package


Figure 2: Appearance of the device in the final package for sterilization


Figure 3: Device structure

Figure 4: Magic1.2F device structure

4. Technical specification

Table 1 – Materials of devices

Component	Raw material	MAGIC1.2F	MAGIC1.2FM	MAGIC1.2FM180	MA1.2FOLIVE	MAGIC1.5F	MA1.5FMP	MA1.5FOLIVEMP	MAGICSTD	MAGICMP
Connector (Cannula)	Polyamide, Grilamid TR90	+	+	+	+	+	+	+	+	+
	Green pigment (green C), CAS-No. 1328-53-6	+	+	+	+	+	+	+	+	+
Proximal tube	Polyamide, Grilamid TR90	+	+	+	+	+	+	+	+	+
	Green pigment (green C), CAS-No. 1328-53-6	+	+	+	+	+	+	+	+	+
Medial (central) tube	Polyamide, Grilamid L20G	+	+	+	+	+	+	+	+	+
	Thermoplastic polyurethane, Estane 58563	+	+	+	+	+	+	+	+	+
	Green pigment (green C), CAS-No. 1328-53-6	+	+	+	+	+	+	+	+	+
Pink pursil tube 1.8 F	PVC VC 633/1	+	+	+	+	+	+	+	+	+
	Bismuth carbonate (III) basic, CAS # 5892-10-4	+	+	+	+	+	+	+	+	+
	Red pigment (Pigment Red 254) - CAS # 84632-65-5	+	+	+	+	+	+	+	+	+
White pursil tube 2.4 F	PVC VC 611/3	+	+	+	+	+	+	+	+	+
	Bismuth carbonate (III) basic, CAS # 5892-10-4	+	+	+	+	+	+	+	+	+
	White pigment (WHITE TiO2) CAS # 13463-67-7	+	+	+	+	+	+	+	+	+
Turquoise pursil tube 1.2 F	PVC VC 611/3	+	+	+	+					
	Bismuth carbonate (III) basic, CAS # 5892-10-4	+	+	+	+					
	Blue pigment CAS # 57455-37-5	+	+	+	+					
Adhesive (glue)	Loctite 4014 Composition: epoxy resin CAS # 25068-38-6 – 35-50% Butanedioldiglycidyl ether CAS # 2425-79-8 – 25-35% sophorone diamine CAS # 2855-13-2 – 7-10% trimethylhexamethylenediamine CAS # 25620-58-0 -5-7%	+	+	+	+	+	+	+	+	+

	Araldite 2020 Composition: Ethyl-2-cyanoacrylate CAS # 7085-85-0 – 99,9% - 100%	+	+	+	+	+	+	+	+	+
	Hydroquinone CAS # 123-31-9 - 0,01 - < 0,1 %									
Ring	Platinum 90% CAS #7440-06-4 Iridium 10% CAS# 7439-88-5	+	+	+	+	+	+	+	+	+
Pad (PVC/PVP)	Thermoplastic polyurethane, Estane 58277	+	+	+	+	+	+	+	+	+
Hydrophilic coating (P3)	Polyvinyl pyrrolidone (PVP) CAS #9003-39-8	+	+	+	+	+	+	+	+	+
Mandren	Stainless steel (AISI SUS 304L)	+	+	+	+	+	+	+	+	+
Protective rim /Clamps	Polyethylene (Rigidex 5502 GA) / Blue pigment CAS # 57455-37-5	+	+	+	+	+	+	+	+	+

(†): The components are not intended for direct or indirect body contact.

Note: The pigment ratio is generally from 1 g per 5 kg of polymer.

Table 2: Manufacturers of device materials

Raw material	Supplier	Certification
Grilamid TR90 / Grilamid L20G	EMS-Chimie Ag, Switzerland	ISO 16949
Stainless steel (AISI SUS 304L)	Fort Wayne Metals, USA	/
Estane 58563	Lubrizol, Germany	ISO 9001 / ISO 14001
Estane 58277	Lubrizol, Germany	ISO 9001 / ISO 14001
PVC VC 633/1	Ineos Newton Aycliffe Limited, UK	ISO 9001
PVC VC 611/3	Ineos Newton Aycliffe Limited, UK	ISO 9001
Green pigment	Plasticouleurs, France	/
Red pigment	Plasticouleurs, France	/
White pigment	Plasticouleurs, France	/
Blue pigment	Plasticouleurs, France	/
Bismuth carbonate	5N Plus France, France	ISO 9001
Cyanolite 4014	Henkel, Germany	/
Araldite 2020	Huntsman Corporation, USA	/
Platinum-iridium	Johnson Matthey, UK	/
Polyvinyl pyrrolidone (PVP)	Acros Organics, Belgium	ISO 9001
Polyethylene (Rigidex 5502 GA)	INEOS Olefins & Polymers, UK	

Table 3: Device technical characteristics

Characteristic	MAGIC1.2F	MAGIC 1.2FM	MAGIC 1.2FM180	MA1.2F OLIVE	MAGIC 1.5F	MA1.5 FMP	MA1.5F OLIVEMP	MAGICSTD	MAGICMP
Total length, mm	1,650 (+40/-25)	1,650 (+35/-30)	1,800 (+35/-25)	1,650 (+40/-25)	1,600 (+33/-20)	1,650 (+33/-20)	1,600 (+33/-20)	1,550 (+33/-20)	1,650 (+33/-20)
Effective length, mm	1,580 (+40/-25)	1,585 (+35/-30)	1,735 (+35/-25)	1,585 (+40/-25)	1,535 (+33/-20)	1,585 (+33/-20)	1,585 (+33/-20)	1,485 (+28/-25)	1,585 (+33/-20)
Length of turquoise distal part, mm (for MAGIC1.2F catheters)	120 (+10/-5)	30 (+10/-5)	30 (+10/-5)	120 (+10/-5)	-	-	-	-	-
Length of white distal part, mm	30 (+7/-5)	120 (+7/-5)	170 (+7/-5)	30 (+7/-5)	150 (+10/-5)	200 (+10/-5)	200 (+10/-5)	100 (+10/-5)	200 (+10/-5)
Length of pink medial part, mm	300 (+15/-10)	300 (+15/-10)	400 (+15/-10)	300 (+15/-10)	250 (+15/-10)	250 (+15/-10)	250 (+15/-10)	250 (+15/-10)	250 (+15/-10)
Length of green proximal part, mm	1,140 (+8/-5)	1,140 (+8/-5)	1,140 (+8/-5)	1,140 (+8/-5)	1,140 (+8/-5)	1,140 (+8/-5)	1,140 (+8/-5)	1,200 (+8/-5)	1,200 (+8/-5)
Outer diameter of turquoise distal part, mm (for MAGIC1.2F catheters)	0.44 (+0.02/-0.01)	0.44 (+0.02/-0.01)	0.44 (+0.02/-0.01)	0.44 (+0.02/-0.01)	-	-	-	-	-
Outer diameter of white distal part, mm	0.42 (+0.02/-0.01)	0.42 (+0.02/-0.01)	0.42 (+0.02/-0.01)	0.42 (+0.02/-0.01)	0.52 (+0.02/-0.01)	0.52 (+0.02/-0.01)	0.52 (+0.02/-0.01)	0.62 (+0.02/-0.01)	0.62 (+0.02/-0.01)
Outer diameter of pink medial part, mm	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)	0.83 (+0.02/-)

	0.01)	0.01)	0.01)	0.01)	0.01)	0.01)	0.01)	0.01)	0.01)
Outer diameter of green proximal part, mm	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)	0.95 (+0.02/-0.01)
Maximum allowable pressure in use, bar	7	7	7	7	7	7	7	7	7
Dead space, ml	0.38	0.38	0.40	0.38	0.37	0.38	0.38	0.37	0.38
Minimum internal clearance, mm	0.23	0.23	0.23	0.23	0.28	0.28	0.28	0.32	0.32
Weight, g	1.55	1.55	1.55	1.55	1.59	1.60	1.61	1.63	1.65
Connector diameter, mm	0.80	0.80	0.80	0.80	1.02	1.02	1.02	1.12	1.12
Maximum volume used with a syringe, ml	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5	2.5
Compatible with the guiding catheter in diameter, mm	1.2	1.2	1.2	1.2	1.2	1.2	1.2	1.2	1.2
Compatible with the guiding guidewire in diameter, mm	0.21	0.21	0.21	0.21	0.21	0.21	0.21	0.21	0.21
Minimum breaking force of turquoise distal part, N	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Minimum breaking force of white distal part, N	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Minimum breaking force of pink medial part and green proximal part, N	3	3	3	3	3	3	3	3	3

Unspecified permissible values: $\pm 1\%$

Appearance of the device: The device surface has no foreign inclusions, bubbles or streaks. The surface is smooth, has no technological and surface defects and provides minimal trauma to vessels when being used. The catheter is biologically safe and compatible with biological body fluids. Cannula and other parts of the catheter do not allow fluid leakage. The device is not implantable and is removed from the human body after being used.

4.1.Construction of the product body and peculiarities of use:

Figure 5: Construction of medical Devices *MAGIC1,2F*, *MAGIC1,2FM*, *MAGIC1,2FM180*, *MA1,2FOLIVE*

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

Image 5 / Picture 5

Device features: Presence of a super-flexible turquoise distal part. For *MA1,2FOLIVE* model, an olive-shaped distal tip. Straight tip, constant angle.

Device application features: Microcatheter should be individually chosen by a physician performing a selective and/or hyperselective vascular catheterization procedure on the basis of the patient's vessel anatomy. The vessel diameter should exceed the used catheter diameter.

Figure 6: Construction of medical Devices *MAGIC1,5F*, *MA1,5FMP*, *MA1,5FOLIVEMP*

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Image 4 / Picture 4



Image 5 / Picture 5

Device features: For *MA1.5FOLIVEMP* model, an olive-shaped distal tip, and possibility of changing a distal tip angle as required. For *MA1.5FMP* model, the possibility of changing a distal tip angle as required.

Device application features: Microcatheter should be individually chosen by a physician performing a selective and/or hyperselective vascular catheterization procedure on the basis of the patient's vessel anatomy. The vessel diameter should exceed the used catheter diameter.

Figure 7: Construction of medical Devices *MAGICSTD*, *MAGICMP*

Image 1 / Picture 1



Image 2 / Picture 2

Product reference	A (mm)	B (mm)	C (mm)	D +33/-20 (mm)	E +33/-20 (mm)	Marquage embase / Hub marking
MAGICSTD	100	250	1200	1550	1490	MAGIC1.8FSTD
MAGICMP	200	250	1200	1660	1590	MAGIC1.8FMP

Image 3 / Picture 3

Device features: For *MAGICMP* model, the possibility of changing a distal tip angle as required. For *MAGICSTD* model, a straight distal tip, constant angle.

Device application features: Microcatheter should be individually chosen by a physician performing a selective and/or hyperselective vascular catheterization procedure on the basis of the patient's vessel anatomy. The vessel diameter should exceed the used catheter diameter.

4.2. Device range:

MAGIC microcatheter for injection or occlusion, versions:

- MAGIC1.2F
- MAGIC1.2FM
- MAGIC1.2FM180
- MA1.2FOLIVE
- MAGIC1.5F
- MA1.5FMP
- MA1.5FOLIVEMP
- MAGICSTD
- MAGICMP
- user manual (to be supplied with each device) - 1 pc.

5. Marking and packing

1 MAGIC microcatheter for injection or occlusion is supplied packaged in an individual Tyvek pack which is the final package for sterilization and in an individual cardboard package. User manual is attached to each device.

5.1. Technical characteristics of the final package for sterilization:

On one side packs are made of Tyvek® 2 FS UNCo 7502918M coating and on one side they are made of PET/PE 65 in accordance with standard NF EN ISO11607-1 (2009). The final package for sterilization is glued together.

Table 4: Pack technical characteristics

Basic properties	Results
Pack overall dimensions (LxWxT), mm	290±5x255±5x0.3
Width of adhesive layer on the long side, mm	8
Width of adhesive layer on the wide side, mm	8
Thickness of adhesive layer, mm	0.3
Tyvek paper - Unit weight	59.5 g/m ² (56.5-62.5)
Tyvek paper - Thickness	155 µm (75- 235)
Tyvek paper - Porosity	22 s (9-35)
Tyvek paper - Porosity	529 ml/min (155-885)
Tyvek paper - Mullenburst test	925 kPa (760-1,110)
Tyvek paper - MD breaking strength test	156 N/2.54 cm (135 to 177)
Tyvek paper - XD breaking strength test	157 N/2.54 cm (127 - 187)
Tyvek paper - Elmendorf MD breaking test	2.8 N (1.7-4.0)
Tyvek paper - Elmendorf XD breaking test	3.7 N (2.8-4.6)
Tyvek paper - Stratification	2.70 N/2.54 cm (2.05-3.25)
Tyvek paper - Hydrostatic flow pressure	145 cm H ₂ O (125- 165)
PET film - Thickness	65 µm ± 5%
PET film - Total weight	63g/m ²
PET film - Breaking test	30- 55 N/mm ²
PET film - Elongation at break	50- 100 %
PET film - Elongation at break	80- 160 %
PET film - O ₂ permeability	90m ³ /m ² /bar
PET film - CO ₂ permeability	345 m ³ /m ² /bar
PET film - N ₂ permeability	50 m ³ /m ² /bar

PET film - Water vapor permeability	< 3.0 g/m ²
PET film - Sealing parameters	115-150°C
Sealing	> 1.5 N / 15 mm
Manufacturer of final package for sterilization	SUDPACK MEDICA, France

Unspecified permissible values: ±1%

Table 5: Technical characteristics of an individual cardboard package

Basic properties	Results
Box overall dimensions (LxWxT), mm	(260x300x14) ±5
Cardboard density	365 g/m ²
Material	Cardboard
Manufacturer	IMPRIMERIE LLP/BUIGNET Address: 29, avenue des Chataigniers – 95150 Taverny

Unspecified permissible values: ±1%

5.2. Appearance of the final package marking:


Figure 8: MAGIC1.2F package and box marking

Figure 9: MAGIC1.2FM package and box marking

Figure 10: MAGIC1.2FM180 package and box marking

Figure 11: MA1.2FOLIVE package and box marking

Figure 14: MA1.5FOLIVEMP package and box marking

Figure 15: MAGICSTD package and box marking

Figure 16: MAGICMP package and box marking


Figure 14: Box layout

Marking on the final sterilization package is identical to marking on an individual cardboard package.

5.3. Marking contains:

- Device name and model;
- Information on the manufacturer;
- Schematic drawing of the device with specification of the following information:
 - a) Outer diameter
 - b) Effective length
 - c) Total length
 - d) Length of each device part
 - e) Tip diameter
 - f) Tip shape
- Quantity;
- Compatible catheter diameter;
- Compatible guide wire diameter;
- Indication: The MAGIC catheter is not compatible with an immobilizing agent containing dimethyl sulfoxide or methacrylate;
- Lot number;
- Batch number;
- Date of manufacture;
- Expiration date;
- Storage conditions and instructions, including:
 - a) The "Keep away from sunlight" symbol
 - b) The "Keep dry" symbol
 - c) Storage temperature from +15°C to +27°C, humidity of 35% to 80%, atmospheric pressure of 500 to 1,060 hPa.
- The "Do not use if package is damaged" symbol;
- The "Do not reuse" symbol;
- The "Consult instructions for use" symbol;
- The "Caution! Consult accompanying documents" symbol;
- Information on the sterilization method.

5.4. Model code explanation:

Model	Explanation
MAGIC1.2F	MAGIC is the model 1.2F is the outer diameter of the distal tip (F is the French scale of catheter diameters)
MAGIC1.2FM	MAGIC is the model 1.2F is the outer diameter of the distal tip (F is the French scale of catheter diameters) M stands for MORET (name of the neuroradiologist involved in the development of this catheter with different white and turquoise lengths)
MAGIC1.2FM180	MAGIC is the model 1.2F is the outer diameter of the distal tip (F is the French scale of catheter diameters)

	M stands for MORET (name of the neuroradiologist involved in the development of this catheter with different white and turquoise lengths) 180 is the total length of the catheter
MA1.2FOLIVE	MAGIC is the model 1.2F is the outer diameter of the distal tip (F is the French scale of catheter diameters) OLIVE is the olive-shaped distal tip (See 4.1. Device Design and Application Features)
MAGIC1.5F	MAGIC is the model 1.5F is the outer diameter of the distal tip (F is the French scale of catheter diameters)
MA1.5FMP	MAGIC is the model 1.5F is the outer diameter of the distal tip (F is the French scale of catheter diameters) MP stands for multipurpose tip - possibility of changing a distal tip angle
MA1.5FOLIVEMP	MAGIC is the model 1.5F is the outer diameter of the distal tip (F is the French scale of catheter diameters) OLIVE is the olive-shaped distal tip (See 4.1. Device Design and Application Features) MP stands for multipurpose tip - possibility of changing a distal tip angle
MAGICSTD	MAGIC is the model STD stands for the straight distal tip
MAGICMP	MAGIC is the model MP stands for multipurpose tip - possibility of changing a distal tip angle

5.5. List of reference designations

	Catalog number
	Serial number
	Use before
	Date of manufacture
	EU marking
STERILE EO	Sterilized with EO
	Attention! Refer to user manual
	Prohibition on reuse
	Do not use if package is damaged
	Refer to user manual
	Do not expose to sunlight

	Protect from moisture!
	Does not contain diethylhexyl phthalate, may contain particles of polyhexahydrotriazine
EAH13  3 700481 315538	European article number
	Manufacturer
	Temperature limitation

Cardboard packages with the device are packed in a corrugated cardboard transport box.
One transport box contains – no more than 50 cardboard packages.

Transport box marking contains:

- Device name and model;
- Information on the manufacturer;
- Quantity;
- Lot number;
- Serial number;
- Date of manufacture;
- Expiration date;
- "Do not expose to sunlight" symbol;
- "Protect from moisture!" symbol;

6. Compliance with the Russian Federation standards

GOST ISO 10555-1-2011 Intravascular sterile single-use catheters. Part 1. General technical requirements

GOST R 50444-92 "Medical devices, apparatuses and equipment. General technical specifications";

GOST 31814-2012. "Conformity assessment. General sampling rules for products testing upon confirmation of conformity".

GOST 31214-2016 "Medical devices. Requirements for samples and documentation provided for toxicological, sanitary and chemical tests, sterility and pyrogenicity tests".

GOST ISO 10993-1-2011. "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and research".

GOST ISO 10993-4-2011. "Medical devices. Biological evaluation of medical devices. Part 4. Study of devices interacting with blood".

GOST ISO 10993-5-2011 Medical devices. Biological evaluation of medical devices. Part 5. Cytotoxicity studies: in vitro methods

GOST ISO 10993-7-2016. Medical devices. Biological evaluation of medical devices. Part 7. Ethylene oxide sterilization residuals.

GOST ISO 10993-10-2011. Medical devices. Biological evaluation of medical devices. Part 10. Studies of irritating and sensitizing action.

GOST ISO 10993-11-2011. "Medical devices. Biological evaluation of medical devices. Part 11. Study of general toxic action".

GOST ISO 10993-12-2015. "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and control samples".

GOST ISO 10993-17-2011. Medical devices. Biological evaluation of medical devices. Part 17. Establishment of threshold values for leachables

GOST ISO 10993-18-2011. Medical devices. Biological evaluation of medical devices. Part 18. Study of materials chemical properties.

GOST R ISO 15223-1-2014 "Medical devices. Symbols used for medical devices marking, labels and accompanying documentation. Part 1. General requirements";

GOST ISO 11607-2011 "Package for medical devices subject to final sterilization";

GOST ISO 11135-2012 Medical devices. Validation and current control of sterilization with ethylene oxide.

GOST EN 556-1-2011. Sterilization of medical devices. Requirements for sterile medical devices. Part 1. Requirements for medical devices subject to final sterilization.

GOST ISO 11138-2-2012 Sterilization of medical products. Biological indicators. Part 2. Biological indicators for sterilization with ethylene oxide.

GOST ISO 11737-1-2012 Sterilization of medical devices. Microbiological methods. Part 1. Assessment of the population of microorganisms on products.

GOST R 52770-2016. Medical devices. Safety requirements. Methods of sanitary chemical and toxicological tests.

GOST 25375-82 "Methods, means and regimes of sterilization and disinfection of medical devices. Terms and definitions".

GOST ISO 14971-2011 Medical devices. Application of risk management to medical devices.

GOST R IEC 62366-2013 Medical devices. Designing of medical devices subject to operational suitability.

Control methods in the Russian Federation standards

GOST ISO 10555-1-2011 Intravascular sterile single-use catheters. Part 1. General technical requirements

GOST R 50444-92 "Medical devices, apparatuses and equipment. General technical specifications";

GOST 31814-2012. "Conformity assessment. General sampling rules for products testing upon confirmation of conformity".

GOST 31214-2016 "Medical devices. Requirements for samples and documentation provided for toxicological, sanitary and chemical tests, sterility and pyrogenicity tests".

GOST ISO 10993-1-2011. "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and research".

GOST ISO 10993-4-2011. "Medical devices. Biological evaluation of medical devices. Part 4. Study of devices interacting with blood".

GOST ISO 10993-5-2011 Medical devices. Biological evaluation of medical devices. Part 5. Cytotoxicity studies: in vitro methods

GOST ISO 10993-7-2016. Medical devices. Biological evaluation of medical devices. Part 7. Ethylene oxide sterilization residuals.

GOST ISO 10993-10-2011. Medical devices. Biological evaluation of medical devices. Part 10. Studies of irritating and sensitizing action.

GOST ISO 10993-11-2011. "Medical devices. Biological evaluation of medical devices. Part 11. Study of general toxic action".

GOST ISO 10993-12-2015. "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and control samples".

GOST ISO 10993-17-2011. Medical devices. Biological evaluation of medical devices. Part 17. Establishment of threshold values for leachables

GOST ISO 10993-18-2011. Medical devices. Biological evaluation of medical devices. Part 18. Study of materials chemical properties.

GOST R ISO 15223-1-2014 "Medical devices. Symbols used for medical devices marking, labels and accompanying documentation. Part 1. General requirements";

GOST ISO 11607-2011 "Package for medical devices subject to final sterilization";

GOST ISO 11135-2012 Medical devices. Validation and current control of sterilization with ethylene oxide.

GOST EN 556-1-2011. Sterilization of medical devices. Requirements for sterile medical devices. Part 1. Requirements for medical devices subject to final sterilization.
GOST ISO 11138-2-2012 Sterilization of medical products. Biological indicators. Part 2. Biological indicators for sterilization with ethylene oxide.
GOST ISO 11737-1-2012 Sterilization of medical devices. Microbiological methods. Part 1. Assessment of the population of microorganisms on products.
GOST R 52770-2016. Medical devices. Safety requirements. Methods of sanitary chemical and toxicological tests.
GOST 25375-82 "Methods, means and regimes of sterilization and disinfection of medical devices. Terms and definitions".
GOST ISO 14971-2011 Medical devices. Application of risk management to medical devices.
GOST R IEC 62366-2013 Medical devices. Designing of medical devices subject to operational suitability.

7. Transportation and storage

Transportation conditions: air temperature from -10°C to +35°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 - 1,060 hPa far from direct sunlight.

Storage conditions: store the device in the manufacturer's package in a dry, clean place, away from heat sources and direct sunlight, at controlled temperature. Shelf life is indicated on the label of each device. Air temperature from +15°C to +27°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 – 1,060 hPa.

Operating conditions: air temperature from +20°C to +25°C at relative humidity from 35% to 65% (non-condensing) and atmospheric pressure of 500 – 1,060 hPa far from direct sunlight. The device is temporarily operated in the human body at a temperature of +36°C to +42°C. The device is operated in the human body for less than 24 hours.

8. Information on the device sterility

The device is supplied sterile, method of gas sterilization with Ethylene Oxide.

9. Disposal information

The device shall be disposed in accordance with SanPiN 2.1.7.2790-10 for class B waste (epidemiologically hazardous waste). After the apparatus decontamination methods using physical methods and changing the waste appearance, excluding the possibility of their reuse, class B waste may be accumulated, temporarily stored, transported, destroyed and buried together with class A waste (epidemiologically safe waste approximated to solid household waste). A package of disinfected class B medical wastes shall have a mark indicating the completed waste disinfection.

10. General indications, contraindications and side effects

CAUTION: MAGIC microcatheters are not compatible with embolizers containing DMSO or methyl methacrylate (such as EUDRAGIT®)

Indications:

MAGIC microcatheters are designed for selective and hyper-selective vascular catheterization for diagnostic and therapeutic purposes.

Contraindications:

- MAGIC microcatheter is contraindicated if, in the physician's opinion, the use of such catheter may jeopardize the patient's condition.
- MAGIC microcatheters are contraindicated for newborns, prematurely born, infants, children, pregnant and nursing women.

Potential complications:

Possible complications include, but are not limited to,

- Hematoma at the puncture area
- Vessel perforation or dissection
- Vasospasm
- Bleeding
- Pain and weakness
- Infectious disease
- Thromboembolic episodes
- Neurological disorders, including blood stroke
- Vascular thrombosis
- Death

Administration.

Microcatheter preparation.

Carefully inspect the sterile pack and make sure it is not opened or damaged. Carefully open the sterile pack, remove microcatheter from the protective dispenser and inspect microcatheter before use to make sure it is not damaged. After straightening the microcatheter, remove the mandrel. Control the flow through the microcatheter using a sterile saline solution and rinse the microcatheter. Return the mandrel to its original position in the microcatheter to act as a support for flexible middle and super-flexible distal shafts.

Mandrel shall be used to facilitate the insertion of MAGIC catheters for injection into guiding catheters. This device is not a guidewire and shall never extend beyond the tip of MAGIC catheter for injection.

Mandrel

After a long storage time, a mandrel may stick to the inside of the catheter and it will be difficult to remove it. In this case, you need to straighten the catheter, and then, firmly holding both ends, pull out the mandrel.

Insertion of a microcatheter

It is recommended to use MAGIC microcatheter in conjunction with Balt 6F guiding catheter and an introducer for insertion of BALT 6F for angiographic control, although MAGIC microcatheter is compatible with BALT 5F guiding catheter and an introducer for insertion of BALT 5F envelope, note that the used guiding catheter shall provide an internal clearance of at least 1.20 mm / 0.47 inches. After connecting the rotating haemostatic valve to Luer junction of the guiding catheter:

Open the thumbscrew and carefully insert MAGIC microcatheter and mandrel through the rotating haemostatic valve. Tighten O-ring valve around the microcatheter to avoid backflow, but not as tightly to obstruct the advancement of microcatheter.

Push MAGIC microcatheter into the guiding catheter until the junction between the rigid green proximal shaft and the elastic pink mid shaft passes through the rotating haemostatic valve. At this point, take the mandrel out of the catheter.

The rigidity of green proximal shaft allows advancing the catheter into the guiding catheter.

The mandrel must be removed when the junction between the semi-rigid proximal and elastic mid shaft passes through the rotating haemostatic valve. MAGIC catheters for injection are easily advanced further thanks to the rigidity of green proximal shaft, and the flow between the guiding catheter and MAGIC catheter for injection.

Rotating haemostatic valve

Guiding catheter

Proximal shaft TORVER

Mid shaft junction PURSIL

The guiding catheter helps to advance MAGIC microcatheter to the aortic arch. Super-flexible distal shaft of MAGIC microcatheter may be easily advanced at the tip of the guiding catheter due to the arterial flow.

Embolization

Check the microcatheter integrity before repeated placing the guidewire or injecting the embolization drug to prevent vascular injury or unintentional embolization. The microcatheter integrity is angiographically confirmed by the fact that the contrast medium comes out only from the tip of the microcatheter, enabling to view the entire distal area of the microcatheter.

Microcatheter removal

When the injection of the embolization drug is completed, slightly pull in the syringe, and then carefully pull the microcatheter. If it becomes difficult to remove a microcatheter:

Carefully pull the microcatheter to remove "sagging" in the microcatheter.

Gently apply the traction force to microcatheter (approximately 5-10 cm of the microcatheter extension). Do not stretch the microcatheter more than 20 cm to avoid the risk of rupture.

In some complex clinical situations, it may be safer to leave the microcatheter in the vascular system rather than risk to destroy the malformation and therefore hemorrhage, affecting too much on the clamped microcatheter.

This is achieved by stretching the microcatheter and cutting the shaft near the entry point of the vascular access, allowing the microcatheter to remain in the artery.

Precautions for use.

- Do not use if the package is opened or damaged. If the package is intact, they are sterile.
- MAGIC microcatheters are non-reusable. Do not reuse. Any reuse of the device results in a high risk of infection of the patient and deterioration of the device performance.

Do not resterilize these devices.

Keep in a dry, dark place at room temperature.

Observe the shelf life.

These devices shall only be used by medical specialists in the sphere of interventional neuroradiology and/or interventional radiology.

- DEHP Presence of phthalate.

- Follow the devices manual that shall be washed through MAGIC microcatheters.

- The pressure used shall not exceed the maximum value of 7 bar/100 pounds per square inch, as indicated on the label.

- Do not use a syringe with a volume of less than 2.5 ml (with a piston diameter of less than 10 mm). The procedure shall be stopped and the microcatheter shall be immediately removed if there is any resistance during the injection.

- Never use MAGIC microcatheters with embolizers containing DMSO or methyl methacrylate (such as EUDRAGIT®).

- MAGIC microcatheters differ from other microcatheters by their super-flexible distal shaft. Their breaking strength is limited and inversely proportional to the microcatheters diameter.

- Mandrel shall not be used as a guidewire. Incorrect use of mandrel as a guidewire may cause perforation of microcatheter and/or vessel. No actions may be made with a mandrel inside the microcatheter until it is completely removed.

- **Never apply** force to resistance of an intravascular device without an angiographic control. Applying force to resistance may damage the microcatheter or harm the patient.

Warranties

Shelf life - 5 years. Sterilization validity - 5 years.

Address for complaints acceptance:

Perinatal Medicine Center Limited Liability Company

Legal/Actual address: 109004, Moscow, Tovarishchesky per., 20, bldg. 4

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+33 (0)1 39 89 46 41

+33 (0)1 34 17 03 46

Data for development and production

Development, production and sale of the device complies with the requirements of EN ISO 13485:2012 + AC:2012 (Certificate No. 513975 MP2012, effective date: 06.02.2018, expiry date: 31.03.2019, notified body: DQS Medizinprodukte GmbH); the device is marked with CE (EC certificate No. 513975 MR2 annex II excluding section 4, effective date: 12.05.2018, expiry date: 31.07.2022, EC certificate No. 534993 MRA annex II section 4, effective date: 01.08.2017, expiry date: 31.07.2022, notified body: DQS Medizinprodukte GmbH).

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