INSTRUCTION FOR USE


SU-POR Surgical Implant

PORIFEROUS, LLC


535 Pine Road, Suite 206
Newnan, GA 30263
United States


Phone +1 770-683-3855, +1 877-631-1594
Fax +1 877-683-7459
www.su-porior.com

I, Aaron Noble, do hereby state that the attached document is a true copy of the original document.

Signature

Sworn to (or affirmed) and subscribed before me
this 24th day of May, 2018 in the presence of:

Notary Public

05/22/2018, version 1

MEDICAL DEVICE NAME

Su-Por Surgical Implants

MANUFACTURER

Poriferous, LLC

535 Pine Road Suite 206

Newnan, GA 30263, USA

Phone +1 770-683-3855, +1 877-631-1594

Fax +1 877683-7459

DESCRIPTION

SU-POR Surgical Implants are manufactured using porous high-density polyethylene, a biomaterial that is readily modified and shaped to suit the functional and anatomical requirements of the patient. The interconnecting pores of the Porous HDPE material permits fibrovascular ingrowth into the implant. The implant will not show through overlying tissue due to its light color. In vitro and in vivo biocompatibility studies have shown SU-POR Surgical Implants to be free from any observable systemic or cytotoxic effects. SU-POR Surgical Implants are provided STERILE and NONPYROGENIC. SU-POR Surgical Implants are MR Safe and pose no known hazards in all magnetic resonance imaging environments. SU-POR Surgical Implants are not made with natural rubber latex.

TECHNICAL CHARACTERISTICS

1. **Sheet**
The SU-POR (biomaterial) Sheet provides surgeons with an excellent option for craniofacial reconstruction and augmentation

Weight	A	B	C
.34gr	38mm	50mm	0.25mm
.64gr	50mm	76mm	0.25mm
.41gr	38mm	50mm	0.35mm
.87gr	50mm	76mm	0.35mm
.43gr	30mm	50mm	0.40mm
.56gr	38mm	50mm	0.45mm
1.1gr	50mm	76mm	0.45mm
.90gr	38mm	50mm	0.85mm
2.0gr	50mm	76mm	0.85mm
1.5gr	38mm	50mm	1.5mm
3.2gr	50mm	76mm	1.5mm
2.9gr	38mm	50mm	3.0mm

A 3D perspective diagram of a rectangular sheet. Dimension A is the length, dimension B is the width, and dimension C is the thickness, indicated by arrows pointing to the respective edges.

2. **Membrane Sheet**
The SU-POR Membrane Sheet is designed to selectively prevent tissue attachment to one side of the implant surface. The membrane layer is comprised of solid polyethylene and heat bonded to the porous layer without adhesives or additives. Tissue integration occurs in the porous layer just as with the fully porous Sheet implant.


Weight	A	B	C
.80gr	38mm	50mm	0.60mm
1.1gr	38mm	50mm	1.0mm
1.7gr	38mm	50mm	1.6mm
2.3gr	50mm	76mm	1.0mm
3.7gr	50mm	76mm	1.6mm

A 3D perspective diagram of a rectangular membrane sheet. Dimension A is the length, dimension B is the width, and dimension C is the thickness, indicated by arrows pointing to the respective edges.

3. Flor-Tec Composite Orbital Floor

The Flor-Tec is a composite structure of SU-POR biomaterial comprised of a solid primary section and a leading porous strip. Solid high-density polyethylene acts to prevent tissue ingrowth while the porous strip may help to facilitate implant attachment.


Weight	A	B	C
.60gr	38mm	50mm	0.30mm
.63gr	38mm	50mm	0.40mm
.78gr	38mm	50mm	0.50mm

4. Enophthalmos Wedge

The Enophthalmos Wedge is designed to mimic the shape of the orbital floor and designed to restore the shape of the orbit.


Weight	Description	A	B	C
1.3gr	Regular - Left	24mm	33.5mm	7mm
1.3gr	Regular - Right	24mm	33.5mm	7mm
2.1gr	Large - Left	28mm	40mm	7.5mm
2.1gr	Large - Right	28mm	40mm	7.5mm

5. Microplate Single Channel Sheet

The Microplate Single Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support. The Microplate Single Channel Sheet accepts plates 1.0mm wide and smaller. Canal diameter is 2.2 mm.


Weight	A	B	C
1.8gr	38mm	50mm	0.85mm

6. Membrane Microplate Single Channel Sheet


The Membrane Microplate Single Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support. The Membrane Microplate Single Channel Sheet accepts plates 1.0 mm wide and smaller. Canal diameter is 2.2 mm.

Weight	A	B	C
1.9gr	38mm	50mm	0.85mm


7.	Miniplate Channel Sheet			
	The Miniplate Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support.			
	Canal diameter is 2.2__ mm. Number of canals - __5__. Distance between canals is 1.1 mm.			
Weight	A	B	C	
2.0gr	40mm	52mm	2.3mm	

8.	Microplate Channel Sheet			
	The Microplate Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support.			
	Canal diameter is 2.2 mm. Number of canals - 5 . Distance between canals is 1.1 mm.			
Weight	A	B	C	
2.0gr	40mm	52mm	2.3mm	

9.	Membrane Miniplate Channel Sheet			
	The Membrane Miniplate Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support.			
	The Membrane Layer acts to inhibit tissue ingrowth.			
Canal diameter is 2.2 mm. Number of canals - 5 . Distance between canals is 1.1 mm.				
Weight	A	B	C	
2.4gr	40mm	52mm	2.3mm	


10.	Membrane Microplate Channel Sheet			
	The Membrane Microplate Channel Sheet is designed for repair of significant orbital floor and wall trauma where the addition of a rigid fixation plate provides structural support.			
	The Membrane Layer acts to inhibit tissue ingrowth.			
Canal diameter is 2.2 mm. Number of canals - 5 . Distance between canals is 1.1 mm.				
Weight	A	B	C	
2.4gr	40mm	52mm	2.3mm	


11.	Sphere	
	The SU-POR Sphere provides surgeons with an excellent fully porous option for enucleation and evisceration procedures. The Sphere is available in multiple sizes to ensure the best possible fit.	
	Weight	Diameter
	.65gr	14mm
	1.0gr	16mm
	1.4gr	18mm
	1.7gr	19mm
	2.0gr	20mm
	2.3gr	21mm
	2.7gr	22mm
3.0gr	23mm	
		

12. **Cor-Tec Sphere**


The Cor-Tec Sphere is an ocular implant with a thin porous layer surrounding a solid core. Full tissue integration is achieved within the thin outer porous layer.


Weight	Diameter
1.3gr	16mm
2.0gr	18mm
2.9gr	20mm
3.8gr	22mm



13.	Quadro-Port Tunnel Orbital Sphere		
	The Quadro-Port Tunnel Orbital Sphere is an ocular implant with a smooth anterior surface and pre-fabricated tunnels for sutures to pass through when attaching extraocular muscles. Opening diameter is 2.2 mm. Number of openings - 4. Distance between canals is 4.1 mm.		
	Weight	Diameter	
	1.1gr	16mm	
	1.4gr	18mm	
	2.1gr	20mm	
	2.8gr	22mm	


14. Conical Orbital Implant (COI)	The Conical Orbital Implant (COI) provides an excellent option for enucleation and evisceration procedures where more volume is required. The conical shape provides more volume posteriorly.		
Weight	Volume	Diameter	
1.4gr	3.0 ml	16mm	
2.0gr	4.2 ml	18mm	
2.7gr	5.6 ml	20mm	
3.5gr	6.0 ml	22mm	


15. Quadro-Port Tunnel Conical Orbital Implant (COI)	The Quadro-Port Tunnel Conical Orbital Implant (COI) is a conical shaped ocular implant with a smooth anterior surface and pre-fabricated tunnels for sutures to pass through when attaching extraocular muscles. Opening diameter is 2.2 mm. Number of openings - 4. Distance between canals is 4.1mm.		
Weight	Volume	Diameter	
1.6gr	3.0 ml	16mm	
1.9gr	4.2 ml	18mm	
2.8gr	5.6 ml	20mm	
3.6gr	6.0 ml	22mm	


16.	Inferior Orbital Rim				
	The Inferior Orbital Rim can provide up to 5mm of anterior projection and is designed to be trimmed to meet the needs of the individual patient.				
	Weight	Description	A	B	C
	.59gr	Left	43mm	18mm	3.2mm
	.59gr	Right	43mm	18mm	3.2mm
					


17.	Inferior Medial Orbital Rim				
	The Inferior Medial Orbital Rim is designed to be placed over the inferior orbital rim and extend superiorly and inferiorly medial to the inferior orbital nerve.				
	Weight	Description	A	B	C
	.39gr	Left	25mm	26mm	2.5mm
	.39gr	Right	25mm	26mm	2.5mm
					


18.	Superior Lateral Orbital Rim			
	The Superior Lateral Orbital Rim is designed to augment the lateral and superior orbital rims, and is designed to be trimmed to meet the needs of the individual patient.			
	Weight	Description	A	B
	1.0gr	Left	33mm	45mm
	1.0gr	Right	33mm	45mm
				


19.	Extended Orbital Rim				
	The Extended Orbital Rim provides surgeons with a large amount of implant to work with, allowing for trimming as required, to fit a large array of orbital rim defects.				
	Weight	Description	A	B	C
	1.5gr	Left	47mm	40mm	6.33mm
	1.5gr	Right	47mm	40mm	6.33mm
					


20.	Orbital Rim Onlay				
	The Orbital Rim Onlay is designed to augment the lateral and inferior orbital rims and subtly increase the anterior projection.				
	Weight	Description	A	B	C
	1.6gr	Left	40mm	40mm	8.45mm
	1.6gr	Right	40mm	40mm	8.45mm
					

21.	Midface Contour						
	The Midface Contour is an onlay designed to augment or repair the non-load bearing bony structures of the midface.						
	Weight	Description	A	B	C		
	2.5gr	Left	60mm	40mm	4mm		
	2.5gr	Right	60mm	40mm	4mm		

22.	Midface Rim						
	The Midface Rim is designed to augment and repair non-load bearing bony structures of the midface and the inferior orbital rim.						
	Weight	Description	A	B	C		
	1.3gr	Left	47mm	28mm	3mm		
	1.3gr	Right	47mm	28mm	3mm		

23.	Inferior 2/3 Orbit						
	The Inferior 2/3 Orbit is designed to replace non-load bearing bony structures of the inferior orbit.						
	Weight	Description	A	B	C		
	18.2gr	Left	61mm	97mm	60mm		
	18.2gr	Right	61mm	97mm	60mm		


24.	Complete Orbit						
	The Complete Orbit is designed to replace non-load bearing bony structures of the orbit.						
	Weight	Description	A	B	C		
	22.5gr	Left	77mm	97mm	63mm		
	22.5gr	Right	77mm	97mm	63mm		

25.	Extended Malar						
	The Extended Malar is designed to add more volume to the malar or to rebuild the contour of the bony structure.						
	Weight	Description	A	B	C		
	1.2gr	Small - Left	45mm	24mm	3mm		
	1.2gr	Small - Right	45mm	24mm	3mm		
	2.1gr	Medium - Left	50mm	26mm	4mm		
	2.1gr	Medium - Right	50mm	26mm	4mm		
	2.6gr	Large - Left	55mm	27mm	5mm		
	2.6gr	Large - Right	55mm	27mm	5mm		

26. SP1 Malar

The SP1 Malar is designed to subtly augment the malar bone.


Weight	Description	A	B	C
1.1gr	Small - Left	50mm	19mm	3mm
1.1gr	Small - Right	50mm	19mm	3mm
1.6gr	Regular - Left	50mm	19mm	5mm
1.6gr	Regular - Right	50mm	19mm	5mm

27. SP2 Malar

The SP2 Malar is designed to subtly augment the malar bone.


Weight	Description	A	B	C
1.5gr	Small - Left	64mm	19mm	3mm
1.5gr	Small - Right	64mm	19mm	3mm
2.1gr	Medium - Left	64mm	19mm	4.5mm
2.1gr	Medium - Right	64mm	19mm	4.5mm
2.9gr	Large - Left	64mm	19mm	7mm
2.9gr	Large - Right	64mm	19mm	7mm

28. Orbito-Zygomatic

The Orbito-Zygomatic is designed for reconstruction of the superior and lateral surfaces of the orbital roof.


Weight	Description	A	B	C
1.3gr	Left	33mm	38mm	0.8mm
1.3gr	Right	33mm	38mm	0.8mm

29. Contoured Two-Piece Chin

The Contoured Two-Piece Chin is designed with a gradual taper and concave posterior surface to provide an excellent anatomical fit to the bony anatomy.


Weight	Description	A	B	C
.80gr	Small	72mm	42mm	3mm
1.6gr	Medium	74mm	42mm	5mm
2.4gr	Large	78mm	50mm	7mm
3.5gr	Extra Large	80mm	55mm	9mm



30. Two-Piece Chin


The Two-Piece Chin is designed to allow for easy insertion and placement of the implant. The surgeon can then attach the components together for proper alignment.


Weight	Description	A	B	C
1.9gr	Small	56mm	33mm	5mm
2.9gr	Medium	56mm	36mm	7mm
3.7gr	Large	57mm	38mm	9mm


31.	Button Chin				
	The Button Chin is designed for subtle augmentation to the medial anterior point of the chin. Having multiple sizes allows the surgeon to select the most natural looking configuration.				
	Weight	Description	A	B	C
	1.3gr	Small	40mm	25mm	4mm
	1.8gr	Medium	47.5mm	37.5mm	5.5mm
2.2gr	Large	48.5mm	38mm	7mm	
					

32.	Round Extended Chin				
	The Round Extended Chin is designed to provide tri-dimensional projection (anterior, lateral and inferior).				
	Weight	Description	A	B	C
	3.5gr	Small	45mm	47mm	3mm
	5.1gr	Medium	45mm	47mm	5mm
6.4gr	Large	45mm	47mm	7mm	
					

33.	Square Extended Chin				
	The Square Extended Chin is designed to provide tri-dimensional projection (anterior, lateral and inferior).				
	Weight	Description	A	B	C
	3.8gr	Small	45mm	47mm	3mm
	5.7gr	Medium	45mm	47mm	5mm
7.2gr	Large	45mm	47mm	7mm	
					

34.	Geniomandibular Groove				
	The Geniomandibular Groove is designed to augment the geniomandibular groove. The implant is divided medially for separate insertion of the left and right pieces.				
	Weight	A	B	C	
	2.8gr	45mm	41mm	4mm	


35.	Contoured Mandibular Angle				
	The Contoured Mandibular Angle is designed to subtly augment the mandible.				
	Weight	Description	A	B	C
	3.1gr	Left	59mm	29mm	7mm
	3.1gr	Right	59mm	29mm	7mm
					

36.

SP Mandibular Angle

The SP Mandibular Angle is designed to conform to the posterior and inferior borders of the mandibular angle.

Weight	Description	A	B	C
2.4gr	Small - Left	65mm	35mm	3mm
2.4gr	Small - Right	65mm	35mm	3mm
3.8gr	Medium - Left	65mm	35mm	7mm
3.8gr	Medium - Right	65mm	35mm	7mm
5.1gr	Large - Left	65mm	35mm	11mm
5.1gr	Large - Right	65mm	35mm	11mm


37.

Lateral Augmentation Mandible

The Lateral Augmentation Mandible is designed to subtly augment the lateral projection of the mandible.

Weight	Description	A	B	C	D	E
3.3gr	Regular - Left	47mm	38mm	3mm	3mm	6.5mm
3.3gr	Regular - Right	47mm	38mm	3mm	3mm	6.5mm
5.6gr	Large - Left	57mm	40mm	4mm	3mm	10mm
5.6gr	Large - Right	57mm	40mm	4mm	3mm	10mm





38.


Ramus of the Mandible


The Ramus of the Mandible is designed for augmentation of the ramus of the mandible.


Weight	Description	A	B	C	D	E
5.6gr	Large 5mm - Left	79mm	32mm	5mm	10mm	7mm
5.6gr	Large 5mm - Right	79mm	32mm	5mm	10mm	7mm
6.2gr	Large 10mm - Left	79mm	32mm	10mm	10mm	7mm
6.2gr	Large 10mm - Right	79mm	32mm	10mm	10mm	7mm
4.1gr	Regular 5mm - Left	79mm	32mm	5mm	4mm	5mm
4.1gr	Regular 5mm - Right	79mm	32mm	5mm	4mm	5mm
4.5gr	Regular 10mm -	79mm	32mm	10mm	4mm	5mm

		Left				
						
39.	Two-Piece Auricular Implant The design of the SU-POR Auricular implants allow for surgeons to custom shape the height and projection of the helical rim to match the contralateral ear of the patient. The porous material provides a structural base for a temporal parietal fascia flap and skin grafts. The success of the implant depends on the technique of the surgeon; the porous structure requires a vascular tissue flap such as a temporal parietal fascia flap and skin graft, to prevent late exposure of the implant.					
	Weight	Description	A	B		
	3.0gr	Ear Base - Left	32mm	63mm		
	3.0gr	Ear Base - Right	32mm	63mm		
	1.0gr	Helical Rim - Left	39mm	63mm		
	1.0gr	Helical Rim - Right	39mm	63mm		
40.	Ear Wedge The Ear Wedge is designed to enhance the projection of a surgically reconstructed ear. The Ear Wedge can be trimmed to match the projection of the contralateral ear.					
	Weight	Description	A	B	C	
	1.2gr	Right	15mm	44mm	11mm	
	1.2gr	Left	15mm	44mm	11mm	
41.	Nasal Dorsal Shell The Nasal Dorsal Shell provides surgeons with an excellent option for augmenting or correcting deformities of the nose.					
	Weight	A	B	C		
	.35gr	43mm	16mm	22mm		






42.	Nasal Batten			
	The Nasal Batten is designed for nasal reconstruction procedures involving the external nasal valve. Nasal Battens are packaged 2 per package.			
	Weight	A	B	C
	.09gr	25mm	12.5mm	0.60mm
				

43.	Nasal Sheet			
	The Nasal Sheet provides surgeons with a solution for when nasal tip projection is required. The Nasal Sheet can be used to support the tip by implantation between the medial crura of the alar cartilage.			
	Weight	A	B	C
	.21gr	30mm	9mm	1.1mm
				

44.	Nasal Radix			
	The Nasal Radix offers surgeons an excellent option to augment a low nasal radix.			
	Weight	A	B	C
	.30gr	24mm	3mm	10mm
				

45.	Nasal Arch				
	The Nasal Arch can be used effectively to create a nasal onlay where subtle augmentation of the dorsum is required. Care should be taken to place the Nasal Arch appropriately in the dorsum area and to avoid extending the Nasal Arch proximally into the soft nasal cartilage area of the tip. The edge of the Nasal Arch should be feathered to promote a smooth transition from the implant to the patient's natural contour.				
	Weight	Description	A	B	C
	1.2gr	Small	70mm	13mm	2mm
	1.3gr	Medium	70mm	15mm	2mm
	1.5gr	Large	70mm	17mm	2mm
					


46.	Nasal Dorsum				
	The Nasal Dorsum is designed to provide subtle augmentation to the dorsum.				
	Weight	Description	A	B	C
	.75gr	Tapered	55mm	6mm	9mm
	.77gr	Small	54mm	6mm	11mm
	1.9gr	Large	67mm	9mm	14mm

					
47.	Paranasal Paranasal implants are designed for augmentation and restoration of the midface in patients who have midfacial deficiency.				
	Weight	Description	A	B	C
	.63gr	Small - Left	28mm	26mm	4.5mm
	.63gr	Small - Right	28mm	26mm	4.5mm
	1.2gr	Large - Left	30mm	28mm	7mm
	1.2gr	Large - Right	30mm	28mm	7mm
					
48.	Pterional The Pterional is designed to correct temporal hollowing in patients who have had surgery involving the pterional approach to the brain. The implant is placed deep to the temporalis during closure. Number of protrusions – 7. Distance between protrusions is 3.0 mm.				
	Weight	Description	A	B	C
	1.7gr	Left	6mm	43mm	44mm
	1.7gr	Right	6mm	43mm	44mm
					
49.	Mastoid The Mastoid provides surgeons with an excellent option for repair of mastoid defects. The regular Mastoid is available in left and right configurations while the small Mastoid provides a universal fit.				
	Weight	Description	A	B	C
	2.3gr	Small	36mm	45mm	1mm
	2.3gr	Left	56mm	53.3mm	1.5mm
	1.1gr	Right	56mm	53.3mm	1.5mm
					

50. Cranial Flex Grid

The Cranial Flex Grid is designed to fill full thickness cranial defects as an option to calvarial bone grafts. The Cranial Flex Grid has a design that is strong and flexible and allows for the implant to be cut to the desired shape. Number of protrusions – 144. Distance between protrusions is 3.0 mm.


Weight	A	B	C
24.1gr	97mm	106mm	6mm

**51. Temporal Flex Grid**

The Temporal Flex Grid is designed to augment deficient soft tissue in the temporal region. The Temporal Flex Grid has a thin contoured temporal surface designed for a proper anatomical fit and a more natural result.

Number of protrusions – 54 . Distance between protrusions is 2.0 mm.





Weight	Description	A	B	C
6.8gr	Small - Left	61mm	78mm	18mm
6.8gr	Small - Right	61mm	78mm	18mm
9.1gr	Medium - Left	74mm	93mm	20mm
9.1gr	Medium - Right	74mm	93mm	20mm
10.4gr	Large - Left	82mm	105mm	20mm
10.4gr	Large - Right	82mm	105mm	20mm
4.3gr	SP - Small	70mm	70mm	10mm
9.9gr	SP - Medium	88mm	86mm	15mm
13.8gr	SP - Large	98mm	95mm	18mm

**52. Sellar Floor**

The Sellar Floor is designed to repair the sellar floor. It is available in two sizes and configurations. The larger Sellar Floor is designed with three small tabs to facilitate handling and placement while the regular has a single tab. The Large Sellar Floor also has a membrane layer to aid in preventing tissue ingrowth.

Number of protrusions – 1 (Regular) and 3 (Large + Membrane). Distance between protrusions is 7.0 mm.

Weight	Description	A	B	C	D
.13gr	Regular	20mm	20mm	2.5mm	0.45mm
.82gr	Large + Membrane	40mm	40mm	2.5mm	0.73mm

						
53.	Cranial Hemisphere The Cranial Hemisphere is designed to provide surgeons with a reconstructive option for large cranial defects. The Cranial Hemisphere provides alternatives to customized implants, grafts, and other implant materials.					
	Weight	Description	A	B	C	
	58.0gr	Left	124mm	170mm	4.5mm	
	58.0gr	Right	124mm	170mm	4.5mm	
	78.8gr	Left	124mm	170mm	6mm	
	78.8gr	Right	124mm	170mm	6mm	
54.	Cranial Dome The Cranial Dome is designed to provide surgeons with a reconstructive option for large cranial defects. The Cranial Dome approximates the contour of the superior 1/3 of the cranium and is available in two thicknesses.					
	Weight	A	B	C		
	59.1gr	150mm	180mm	4mm		
	79.0gr	150mm	180mm	6mm		
55.	Flex Sheet The Flex Sheet is designed for small to medium sized cranial defects and deformities. It has a smooth exterior surface and a series of pedicles on the interior surface that are designed to provide volume and flexibility. Number of protrusions – 104 . Distance between protrusions is 1.5 mm.					
	Weight	Description	A	B	C	
	6.6gr	Regular	56mm	91mm	4.5mm	
	7.4gr	Regular + Membrane	56mm	91mm	4.5mm	

56. Block
 The SU-POR Block gives surgeons excellent options for craniofacial reconstruction and augmentation. SU-POR Block implants are available in multiple sizes and thicknesses to ensure that the surgeon has the proper amount of biomaterial for almost any possible contour or shape required.

Weight	A	B	C
.83gr	13mm	38mm	3mm
2.0gr	25mm	50mm	3mm
2.8gr	38mm	63mm	3mm
1.4gr	13mm	38mm	6mm
4.0gr	25mm	50mm	6mm
8.1gr	38mm	63mm	6mm
2.5gr	13mm	38mm	9.5mm
6.6gr	25mm	50mm	9.5mm
12.4gr	38mm	63mm	9.5mm

57. Burr Hole Cover
 The Burr Hole Cover is designed to fit into and over holes made by a cranial perforator. The large Burr Hole Cover has a 14mm diameter with a stem that is easily modified. The small Burr Hole Cover is designed for 5mm diameter holes. Provided 3 pcs per package.

Weight	Description	A	B	C
.16gr	Small	15mm	3mm	5mm
.82gr	Large	29mm	7mm	14mm

All the characteristics of the Implants have a permissible deviation of $\pm 7\%$.

Average stress 5% strain with stress tolerance based on maximum and minimum stresses: 220psi (1517 kPa) + 15psi (103.4 kPa) / - 15psi (103.4 kPa).

Average stress 5% yield strength with strength tolerance based on maximum and minimum strengths: 680psi (4688.435 kPa) + 15psi (103.4 kPa) / - 15psi (103.4 kPa).

Implant do not need to be connected to other medical devices or equipment in order to operate as required for its intended use.

The device does not have a measuring function.

INDICATIONS

SU-POR Surgical Implants in block, sheet and anatomical shapes are intended for non-weight bearing applications of craniofacial reconstruction/cosmetic surgery and repair of craniofacial trauma. SU-POR Surgical Implants are also intended for the augmentation or restoration of contour in the craniomaxillofacial skeleton.

CONTRAINDICATIONS

1. Active infection.
2. In areas where there are prolonged, significant, compressive stresses.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate testing should be made to rule out this possibility prior to implantation.
4. Patients who are unwilling or unable to follow postoperative instructions due to accompanying conditions (mental or physical).
5. Limitations in blood supply and/or systemic disorders that may slow healing and increase the possibility of infection and/or rejection of the implants.
6. Use in tissue that has been compromised by cancer therapies.
7. Inadequate coverage with healthy tissue.

8. Any degenerative disease process that would adversely affect proper placement of the implants.
9. Procedures in which a non-sterile environment exists, such as sinuses.
10. Use in weight bearing area such as the temporal mandible joint.

POSSIBLE ADVERSE EFFECTS

1. Vascular changes.
2. Allergic reaction to the implant.
3. Deformation or fracture of the implant.
4. Surgical trauma related nerve damage.
5. Pain, discomfort, and /or abnormal sensation due to presence of the implant.
6. Superficial and/or deep infection.
7. Loosening or migration of the implant.

WARNINGS

1. SU-POR Surgical Implants must never be re-used or re-sterilized. Re-use may increase the risk of infection.
2. Sufficient patient instruction is a crucial factor in the success of the surgical procedure. Postoperative follow-up and care are essential. The patient must be made aware that excessive physical activity or load bearing may cause implant migration, loosening, deformation, or fracture. The patient must be informed of the risks of using the implants, including possible adverse effects.
3. The success of any implant is dependent upon careful handling and proper surgical technique. When shaping and contouring the implant, sharp edges should be avoided to minimize trauma to the surrounding tissue.
4. The implant must never be shaped by use of intense heat producing devices such as electrosurgical devices, heated knives, or lasers.
5. Because of the porous nature of the implant material, there is an increased possibility of contamination. Therefore, it is most important, when handling and preparing porous polyethylene implants that the highest level of aseptic care be used. Do not store, place, or carve the implant on surfaces which may contaminate it with lint or other particulate matter such as: surgical drapes or cloth.
6. The surgeon must be thoroughly familiar with the implants, the method of application, the instrumentation, and the surgical procedure. These implants are not designed for applications with weight or load bearing requirements.
7. Responsibility for completion of technical training, appropriate patient selection, and choice and placement of the implants is at the discretion of the physician.

INSTRUCTIONS FOR USE

Only a physician, having the necessary qualification, should use SU-POR Surgical Implants. This document is in no way intended to be a comprehensive manual of surgical techniques. Surgeons should utilize their surgical training, clinical experience, and applicable surgical techniques to determine appropriate surgical procedures. Experienced surgical judgment should be utilized in the selection and use of SU-POR implants.

The success of an implant is dependent upon surgical technique. The incision should be made such that the closure does not occur directly over the implant whenever possible. The implant shape may be modified with the techniques outlined in the following sections.

SHAPING THE IMPLANT

SU-POR Surgical Implants are provided as block, sheets, spheres, and basic anatomical shapes, which may be modified by the surgeon to produce the desired result. The material may be carved using scalpel, burr, or cut with scissors, care should be taken to smooth the edges of the implant where it transitions to bone. Trim any poorly attached implant material from the edges. After carving and sculpting is complete, wash the implant with sterile saline to remove any loose

particles from its surface and edges. To help the shaping process, the implant can be submerged in a hot, sterile saline bath (82° - 100°C, 180° - 212°F). Having been heated the implant will be more flexible allowing for modification to the shape. If the implant cools and becomes difficult to bend, it should be returned to the hot saline. Once the desired shape is obtained, the implant should be held in that shape and allowed to cool. A cold, sterile bath can be used to accelerate the cooling process. Repeat the above steps if further molding is required.

PLACING THE IMPLANT

As a prophylactic measure, the implants should be rinsed in a sterile antibiotic solution prior to implantation. Surgeons should remove all carved debris from the surgical site with care. A sufficient tissue pocket should be prepared to ensure that tissue is not pulled taut over the implant. Incisions located exactly over the implant site should be avoided. If the potential for the implant to move before tissue ingrowth exists, fixation can be accomplished with sutures, surgical fixation screws, or K-wire. Fixation screws should be tightened until they are flush with the surface of the implant. Multiple pieces of implant material can be stacked and joined by suturing. A cutting needle can easily penetrate the surface of the implant enabling the surgeon to re-suspend tissue or muscle by attaching it to the implant. For implants placed using the intraoral route, a rubber dam or impervious plastic film should be used to isolate the incision to reduce the risk of contaminating the implant.

SURGICAL REVISIONS

The porous nature of SU-POR Surgical Implants allows for soft tissue ingrowth and vascularization of the implant. In patients that may require later surgical revision, the surgeon should be aware of this vascular and soft tissue ingrowth. In the event revision or removal of the implant is required after ingrowth has occurred, the surrounding soft tissue may be raised with a surgical instrument and the implant dissected out with a scalpel or surgical scissors.

REMOVING OR REPLACING

SU-POR Surgical Implants when being implanted are not intended to be removed or replaced. Implants can be removed in case of complications during 2-3 months without any traumatizing. If removal is needed in the later period scar capsule should be cutted.

AVAILABILITY OF MEDICINAL PRODUCTS, MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN

SU-POR Surgical Implants do not contain any medicinal substances, materials of animal and human origin.

STERILIZATION CONDITIONS

The implants are provided sterile by ethylene oxide (EO) terminal sterilization. The devices are non-pyrogenic. Note the product "Use By" date specified on the package. DO NOT use the device if the "Use by" date has passed.

The implants are intended for single-use only; do not reuse. Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing. Inspect all product prior to use. Do not use if the package is open or damaged, or if product is damaged.

TRANSPORTATION, STORAGE, USE TERMS

Handle with care. Store in original packaging. Keep dry. Keep away from sunlight.

Devices packaged by manufacturer may be transported by all modes of transport in covered vehicles in accordance with the rules of transportation of goods for each type of transport.

The temperature at transportation should be between -50°C to +50°C with a range of $\pm 2^{\circ}\text{C}$ and Relative Humidity from 30% to 75% with a range of $\pm 5\%$.

Products should be stored at $-40^{\circ}\text{C} + 50^{\circ}\text{C}$ in an enclosed space in the original manufacturer's packaging for the labeled shelf life.

Implants are serviceable in the process of operation at nominal temperature values from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$. Implants are resistant in biological fluids and secretions of the body tissues.

SHELF LIFE

The shelf life of SU-POR Surgical Implants is 24 months.

ENVIRONMENT PROTECTION REQUIREMENTS

SU-POR Surgical Implants do not adversely impact the environment.

DISPOSAL

Products, which had contact with blood and/or other biological liquids of patients should be disposed and eliminated in accordance with the requirements for handling medical waste, established by SanPiN 2.1.7.2790-2010 for class B and other normative legal acts of the Russian Federation.

Unused products, including the ones which expired (which had no contact with the biological liquids of patients) should be disposed in accordance with the procedure of disposal of solid domestic waste established by the legislation of the Russian Federation for class A by SanPiN 2.1.7.2790.

LABELING SYMBOLS

	Keep dry
	Do not use if packaging is damaged
	MR safe
	Do not resterilize
	For single use only
	Sterilized with ethylene oxide
	Non-pyrogenic
	Dimensions
	Catalog number


	Caution, consult accompanying documents
	Manufacturer: Poriferous, LLC 535 Pine Road Suite 206, Newnan, GA 30263, USA

WARRANTY

All products are warranted to be free from defects in material and workmanship. No warranty is made for any purpose other than in the product specifications and labeling. Additional information pertaining to specific devices or issues may be obtained from Poriferous, LLC.

AUTHORIZED REPRESENTATIVE IN RUSSIA

LLC "Arion",
Repischeva str., bld.4-2A, office 1,
Saint-Petersburg, 197349 Russia

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