Cardiac Multipack
Instruction for Use
Dated: June 28, 2018


ACKNOWLEDGEMENT

State of Utah)

County of Salt Lake)

On this 28th day of June, 2018, before me, Jennifer Fordham, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.

Notary Public

My Commission Expires: June 27, 2021

Regulatory Affairs Director
Merit Medical Systems, Inc.
USA

Dan Clark

(sign and seal)

June 28 2018

Cardiac Multipack

INSTRUCTIONS FOR USE

1600 West Merit Parkway
South Jordan, Utah 84095
Phone: 801-253-1600
Fax: 801-253-1684 / 6918


INSTRUCTIONS FOR USE

Cardiac Multipack

DEVICE DESCRIPTION


1.1. Cardiac Multipack is composed of:

1. Performa angiographic catheter, configuration: JR3.5, JR3.5 ST, JR4.0, JR5.0, PRC (Progressive Right Coronary), dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm.
2. Performa angiographic catheter, configuration: JL4.0, JL3.5; dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm.
3. Performa angiographic catheter, configuration: Straight Pigtail, 145 ° Pigtail, 155 ° Pigtail, MOD AMP Right (AMPLATZ RIGHT MODIFIED), MIV RADIAL PIG ST, Ultimate 1, Ultimate 2, Ultimate 4, Ultimate 4-T40, Ultimate 4 - T45; dimensions: diameter : 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm, 110cm.

1.2. Cardiac Multipack Plus is composed of:

1. Performa angiographic catheter, dimensions: diameter 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm, 110cm; configuration: Straight Pigtail, 145 ° Pigtail, Ultimate 1.
2. Performa angiographic catheter, dimensions: diameter 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm; configuration: JR4.0, JR5.0, PRC (Progressive Right Coronary).
3. Performa angiographic catheter, dimensions: diameter 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm; configuration JL4.0, JL3.5.
4. Diagnostic Guide wire InQwire, diameter 0.035 "(0.89 mm); length: 150cm, 180cm; tip configuration J3, Fixed Core.
5. Prelude or Prelude PRO Sheath Introducer, diameter 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 7cm, 11cm, consisting of:
 - Prelude or Prelude PRO Sheath Introducer;
 - guide wire, diameter: 0.018 "(0.46mm), 0.035" (0.89mm); 0.038 "(0.97mm); length: 40cm, 50cm;
 - Vascular dilator, dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length 12 cm, 16cm.
 - needle, sizes: 18G, length 7cm; 21G length 4cm.



1.3. Radial Multipack Plus

I. Radial Multipack Plus is composed of:

1. Performa angiographic catheter, dimensions: diameter 5F (1.7mm); length: 100cm; configuration: Ultimate 1, Ultimate4-T40.
2. Diagnostic Guide wire InQwire, diameter 0.035 "(0.89 mm); length: 180cm; tip configuration J3, Fixed Core.
3. Prelude Sheath Introducer, diameter 5F (1.7mm); length 7cm, 11cm; consisting of:
 - Prelude Sheath Introducer;
 - guide wire, dimensions: diameter: 0.018 "(0.46mm), 0.035" (0.89mm); 0.038 "(0.97mm); length: 40cm, 50cm;
 - Vascular dilator, dimensions: diameter 5F (1.7mm); length 12 cm, 16cm.
 - needle, sizes: 21G, length 4cm;

II. Radial Multipack Plus is composed of:

1. Performa angiographic catheter, dimensions. diameter 5F (1.7mm); length: 100cm. configuration: JR4.0;
2. Performa angiographic catheter, dimensions: diameter 5F (1.7mm); length: 100cm. configuration: JL3.5;
3. Diagnostic guide wire InQwire, outer diameter 0.035 "(0.89 mm); length: 180 cm; tip configuration J3, Fixed Core.
4. Prelude Sheath Introducer, diameter 5F (1.7mm); length 7cm, 11cm; consisting of:
 - Prelude Sheath Introducer;
 - guide wire, dimensions: diameter: 0.018 "(0.46mm); length: 40cm, 50cm;
 - Vascular dilator, dimensions: diameter 5F (1.7mm); length 12 cm, 16cm.
 - needle, sizes: 21G, length 4cm;

INTENDED USE / INDICATION FOR USE

Cardiac MultiPack is developed for coronary angiography (coronography) procedures. Coronography is indicated for treatment of, but not limited to, patients with acute and chronic coronary heart disease (CHD), and, in some cases, for patients with suspected CHD. Coronography is also performed before heart surgery (including patients without CHD). The PIGTAIL catheter being a part of Multipack is used in endovascular operations on large vessels (aorta, pulmonary artery, etc.) and heart, for angiography of large vessels (aorta, lower and upper hollow veins, pulmonary artery), and cardiac cavities (most of all, for the left ventricle angiography).

CONTRAINDICATIONS

Contraindicated for use in the coronary and central vasculature.

WARNINGS:

- Angiographic catheters are single use.
- Recommended guide wire size: See the individual catheter label.
- Guide wire use has been associated with greater incidence of thrombus formation.
- Do not advance the introducer and/or guide wire if resistance is met
- Do not leave the introducer in place for extended periods of time without a catheter or an obturator to support the cannula wall.
- Do not use device with a power injector.

PRECAUTIONS:

- Angiography should be undertaken only by an experienced angiographer.
- Read instructions prior to use.
- RX Only • Store in a cool d r y place.
- This device is intended for single use only.
- Do not reuse or resterilize.
- Do not use if package is opened or damaged.
- Do not use after expiration date.
- Do not autoclave.
- Carefully remove the catheter from the package by grasping the hub and slowly withdrawing it from the package.
- If resistance is felt when removing the guide wire from the catheter, remove the guide wire and catheter as a unit to prevent potential damage to the vessel wall.
- Flush all devices entering a blood vessel with sterile heparinized saline or similar isotonic solution before use.
- Guide wires will collect blood and other foreign material in their lumens; neither autoclaving nor ultrasonic cleaning will completely remove foreign material, therefore guide wires are recommended for one time use.

ADVERSE EVENTS:

Merit Medical Multipack Plus should be used by physicians familiar with percutaneous catheter introduction. Complications which may be associated with the use of Multipack Plus include, but are not limited to, the following.


- | | |
|----------------------------------|---------------------------------------|
| • Thrombus formation/emboli | • Air embolism |
| • Hematoma at the puncture site | • Arterial wall damage |
| • Plaque dislodgement | • Myocardial infarction |
| • Cardiac arrhythmias | • Infection |
| • Perforation of the vessel wall | • Allergic reaction to contrast media |
| • Pseudo aneurysm formation | • Vascular spasm |
| • Vascular occlusion | • Bleeding |
| • Stroke and death | • Guide wire embolization |

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Inspect all guide wires prior to use. Do not use any unit if the package is open or damaged.
Employ an aseptic technique during removal from the package and during use.

All guide wires are secured in the hoop dispenser by the locking J-tip straightener. To avoid damaging the guide wire during removal from the flush hoop, grasp the J-tip straightener near the base and slide it forward approximately 5mm or until the J-tip straightener is no longer attached to the flush hoop adaptor. Holding both the guide wire and J-tip straightener, continue to dispense guide wire from the hoop.

Do not use excessive force to advance the moveable core while the guide wire is in a vessel. Advancement with excessive force may cause coil penetration and vessel damage.

Avoid withdrawing PTFE coated guide wires back through a metal needle. The sharp edge of the needle may scrape the coating. It is suggested that a catheter or PTFE vessel dilator replace the needle as soon as the guide wire has reached the appropriate position.

During advancement of the catheter and guide wire within the aorta, it is recommended that the guide wire be removed at the appropriate level of the aorta.

Care should be taken when manipulating a catheter during placement and withdrawal to prevent possible intravascular tissue damage. If resistance is felt when removing a guide wire from a catheter, the guide wire and catheter should be removed as a unit to prevent potential damage to the vessel wall.

A guide wire is a delicate instrument and remains the most fallible instrument used in a percutaneous procedure. Any time that a guide wire is used there is a possibility of thrombus formation/emboli, vessel wall damage, and plaque dislodgement, which could result in myocardial infarction, cardiac arrhythmia, stroke or death. The physician should be familiar with the use of angiography products and the literature concerning the complications of angiography.

INSTRUCTIONS FOR USE

Angiographic catheters

1. A coated guide wire is recommended for use with Merit Angiographic Catheters.
 2. Introduce the catheter into the blood vessel using vascular entry technique of choice.
 3. PIGTAIL catheters. Pigtail catheters are fitted with a straightener. It is split for easy removal after the catheter tip is on the guide wire.
 - Slide straightener forward until pigtail is straight.
 - Place catheter tip on guide wire and remove straightener by pulling on the straightener tab.
 4. MARKER BANDS. Anatomical measurement accuracy is $\pm 4\%$.
- Sterile and non-pyrogenic in unopened and undamaged inner package. Store in a cool dry place.

Prelude Sheath Introducer

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
 2. Remove the Prelude Sheath Introducer components from package using proper aseptic technique.
 3. Flush all components with saline or suitable isotonic solution. After flushing side port, turn stopcock to off position to maintain flush in side port and prevent bleed back upon insertion into the vessel.
 4. Insert vessel dilator into Prelude Sheath Introducer through hemostasis valve and snap into place. Dilator must be securely snapped into place to avoid damage to the vessel.
 5. Insert appropriate access needle into vessel. While holding the access needle, place the flexible end or J end of the guide wire through access needle into vessel. Note – Refer to product labeling for the guide wire size that is compatible with the system components.
- Warning: Do not advance the guide wire if resistance is met. Determine the cause of resistance before proceeding.
6. Hold guide wire in place and remove access needle. Hold pressure at the site until the introducer/dilator assembly is placed. Warning: If a needle with a metal cannula is used, do not withdraw the guide wire after it has been inserted because it may damage the guide wire.
 7. Insert the introducer/dilator assembly over the guide wire into the vessel. Using a rotating motion, advance the introducer/dilator assembly through the tissue into the vessel. Grasp the assembly close to the skin as it is being placed into the vessel to avoid buckling.

8. After introducer/dilator assembly has been placed into vessel, detach the dilator from the introducer by bending the dilator hub down slightly (this will un-snap the dilator hub from the introducer cap). While holding the sheath, carefully remove the dilator and guide wire together, leaving the sheath introducer in the vessel.

9. Aspirate from the side port extension to remove any potential air. After aspiration, flush the side port with a suitable solution. Stopcock should be turned off to maintain flush in side port.

10. Insert selected device(s) (wires, catheters, etc.) into Prelude Sheath Introducer.

Note: Hold the sheath in place when inserting, positioning, or removing the devices. Always exchange or remove devices slowly through the sheath.

11. To temporarily suture the sheath in place, use the rotating suture ring.

12. REMOVAL:

- The sheath may be removed when clinically indicated. Compression on the vessel, above the puncture site, should be started as the sheath is slowly removed. Note. Collected fibrin at the tip of the sheath may be aspirated via the side arm tubing prior to removal of the sheath. Discard the sheath appropriately.

- If the sheath is to be left in place, an obturator of appropriate size should be placed into sheath. After the sheath is flushed, place the obturator through the sheath and snap into place. A suitable solution should be flushed through the side arm after obturator is placed. Note. An obturator that is one French size smaller than the sheath introducer should be used to allow flushing, infusion and pressure monitoring. When clinically indicated, the sheath and obturator may be removed (see above).

InQwire guide wire

Flush Hoop Instruction Guide:

1. Attach flush filed syringe to flush port luer
2. Rotate syringe clockwise (as pictured)
3. Inject saline into hoop
4. Detach syringe from hoop
5. Dispense guide wire

Note: In order to reduce the potential of clot formation, it

is recommended that the guide wire be flushed with saline

or heparinized saline prior to use. Attach a filed syringe to



the luer hub located at the end of the plastic hoop, flush

several times. After flushing, remove guide wire from hoop

and use as described above.

Directions for use: The following schematic shows a typical procedure for percutaneous entry utilizing the Seldinger technique. Variations in individual patient anatomy may preclude the utilization of this technique.

1. Vessel puncture with a two part needle
2. When using a two part needle remove stylet leaving cannula in place, insert flexible (distal) end of guide wire through cannula and into vessel.


3. "J" guide wires are shipped with "J" straightener to aid in the insertion of the wire into the puncture needle. Advance the straightener until 2-3 mm of the tip extends from the tip. Insert wire into hub and through needle. Remove "J" straightener proximally and discard.

4. Remove needle cannula leaving the guide wire within the lumen of the vessel.

5. Pass dilator or catheter over the guide wire directly into the vessel.

6. Carefully remove the guide wire leaving the catheter in place.

Note: Device is sterilized as stated on the package label.

How Supplied: Individually packaged 5-10 per box, refer to catalog for ordering information.

DEVICE SPECIFICATION

PERFORMA ANGIOGRAPHIC CATHETER

Description: Angiographic catheters are single-use, sterile, non-pyrogenic, disposable intravascular catheters

The following product configurations are available in the Performa product family for Cardiac Multipack:

1. Performa angiographic catheter, configuration: JR3.5, JR3.5 ST, JR4.0, JR5.0, PRC (Progressive Right Coronary), dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm.
2. Performa angiographic catheter, configuration: JL4.0, JL3.5; dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm.
3. Performa angiographic catheter, configuration: Straight Pigtail, 145 ° Pigtail, 155 ° Pigtail, MOD AMP Right (AMPLATZ RIGHT MODIFIED), MIV RADIAL PIG ST, Ultimate 1, Ultimate 2, Ultimate 4, Ultimate 4-T40, Ultimate 4 - T45; dimensions: diameter : 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 100cm, 110cm.

Note: Pigtail catheters have a spiral design for uniform contrast delivery: 2mm apart for 10mm segment.
Side holes: 5 spiral holes within 10mm from end of pigtail, allowing contrast to remain within patient ventricle.

Specification compliance table

	Diameter		
4	.042"/1.07 mm	.054"/1.37 mm	0.8mm
5	.046"/1.17 mm	.066"/1.68 mm	0.9mm
5	.052"/1.32 mm	.068"/1.73 mm	0.9mm
6	.054"/1.37 mm	.0775"/1.97 mm	1mm
6	.059"/1.49 mm	.079"/2.01 mm	1mm

The main materials used in the manufacture of the Performa Catheters are identified below.


Wire Braid	304V Stainless Steel
Body (Braided)	Pebax
Body (Nonbraided)	Nylon
Tip (4F, 5F, 6F)	Pebax
Bumper Tip (4F, 5F)	Pebax
Tip (6F) (Nonbraided)	Pebax/Nylon
Sleeve	Pebax
Pigtail Straightener	Nylon
Hub	Polycarbonate

PRELUDE AND PRELUDE SHEATH INTRODUCER

The Merit Prelude Sheath Introducer is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries while maintaining hemostasis for a variety of diagnostic and therapeutic procedures.

The following product configurations are available in the Prelude Sheath Introducer product family for Cardiac Multipack: Prelude or Prelude PRO Sheath Introducer, diameter 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length: 7cm, 11cm, consisting of:

- Prelude or Prelude PRO Sheath Introducer;
- guide wire, diameter: 0.018 "(0.46mm), 0.035" (0.89mm); 0.038 "(0.97mm); length: 40cm, 50cm;
- Vascular dilator, dimensions: diameter: 4F (1.3mm), 5F (1.7mm), 6F (2.0mm); length 12 cm, 16cm.
- needle, sizes: 18G, length 7cm; 21G length 4cm.

Sheath Introducer

Vascular dilator

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The main materials used in the manufacture of the Prelude Sheath Introducer are identified below.

Sheath Hub, all sizes	Medical Grade Resin
Sheath Tubing - 4F through 6F, all lengths	Polyethylene Resin / USP Grade Barium Sulfate

Component	Composition
(PRO) Sheath Tubing - 4F through 6F, all lengths	Polypropylene
Sheath Cap, all sizes	Polycarbonate Resin
Suture Ring, all sizes	Polyethylene Resin
Valve, all sizes	Silicone Rubber
Valve Lubricant	Silicone Fluid
Side-Arm Tubing, all sizes	Polyurethane Resin
Stopcock Handle, all sizes	Acetal Resin
Stopcock Body, all sizes	Polycarbonate Resin
Stopcock Lubricant	Silicone Fluid
Stopcock Cap, all sizes	ABS Resin
Dilator Hub, all sizes	Polyethylene Resin
Dilator Tubing (Polypropylene) - 4F and 5F, all lengths	Polypropylene Resin/ USP Grade Barium Sulfate
Dilator Tubing (Polyethylene)- 6F	High Density Polyethylene Resin/USP Grade Barium Sulfate
Dilator Lubricant, all sizes	Silicone Fluid

INQWIRE GUIDE WIRE

The Merit Diagnostic Guide Wires are intended to facilitate the placement of devices during diagnostic and interventional procedures.

Description: Merit Medical guide wires are fabricated from high quality stainless steel utilizing a sophisticated construction process and are available with or without a PTFE coating. Guide wires are supplied sterile, non-pyrogenic, and are intended for single use only.

Merit Medical guide wires are packaged in a plastic hoop, which is fitted with a luer hub. This packaging is provided to facilitate compliance with the manufacturer recommended guidelines that the wire be flushed with saline or heparinized saline prior to use (See directions for use - Note).

The following product configurations are available in the InQwire product family for Cardiac and Radial Multipack Plus:

- Outer Diameter - 0.035" (.89mm);
- Length – 150 and 180 cm;
- Curved Distal Tip Configurations - J tip 3mm.

The materials used in the manufacture of the InQwire guide wire are identified below.

Component Material Specifications	
Component	Material
Core Wire	304V or 304 Stainless Steel
Wire Coil	304 Stainless Steel
Wire Coil Coating	Polytetrafluoroethylene (PTFE)
Safety Wire	304V or 304 Stainless Steel

Compatibility: 0.035" (.89mm) InQwire guide wire may be used with any device with an inner diameter appropriately sized to accommodate a 0.035 inch (.89mm) guide wire.

Applicable symbols for Cardiac Multipack	
	This symbol means that the device fully complies with European Directive 93/42/EEC.
	Manufacturer
	Authorized Representative
	Single use
	Catalog Number
	Lot Number
	Sterilization: Ethylene-Oxide Gas
	Caution: U. S. federal law restricts this device to sale by or on the order of a physician
	Use by
	Manufacturing date
	Caution: Consult Accompanying Documents
	Do not use if package is damaged.
	Non-pyrogenic
	Нетоксично
	FRENCH SIZE
	RECOM.PRESSURE LIMIT
	RECOM.GUIDEWIRE
	CONTENTS
	SIDE HOLES
	FLOW RATES and PRESSURES. FLOW RATES ACHIEVED AT 10.4cP. INDIVIDUAL RESULTS MAY VARY.

PACKAGING

Packaging of Merit Medical products is conducted with reference to the following standards:

- ISO 11607: 2003: Packaging for terminally sterilized products. All components contained inside MultiPACK packaging. MultiPACK packaging consist of Peelpouch (121x13.5cm) made of Beacon Converters uncoated Tyvek sealed to 2 mil PE/pet) and tray (112x12cm) made of Solid Bleached Sulfated paperboard (SBS). The design of tray help to hold the sheath bag in place.

Sheath and sheath introducer wire are placed in small bag. Wire contained along side the catheters in curved tube. Shelf Transport Carton made of Solid Bleached Sulfated paperboard (SBS).

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS

Transportation temperature: -35 °C to +55 °C

Storage temperature: -35 °C to +55 °C

Usage: 15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

Contact with the internal environment of the body at a temperature of 32 ° C - 42 ° C

Use Conditions - in medical and preventive treatment facility.

SHELF LIFE

Shelf life - 3 year.

STERILIZATION

Ethylene oxide sterilization

DISPOSAL

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

Customer Service

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
+ 1 801 253 1600

Merit Technologies LLC
115054, Moscow, Bakhrushina Str., 32, bldg. 1,
flr/room 4/25A, Russia
+ 7 495 221 89 02

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