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
09,JAN.2019

We (applicant) : Healcerion Co., Ltd.

do hereby solemnly and sincerely declare :

1. That the attached document(s) :

CE 0123(User Manual, Rev.2.5)/Document No. HC-IFU-001 Rev.2.5)
SONON300C ULTRASOUND IMAGING SYSTEM, With accessories

2.

☒

is(are) true and correct.

☐

is(are) true translation for the text(s) originally written in
the Korean language conscientiously believing the same to
be true and correct.

Signature


Healcerion Co., Ltd.

Jeongwon Ryu / CEO

SONON 300C ULTRASOUND IMAGING SYSTEM with accessories

USER MANUAL REV. 2.5

Document No. HC-IFU-001 (Rev. 2.5)

1

The SONON 300C ultrasound imaging system with accessories is a diagnostic ultrasound equipment which transmits ultrasound waves into body tissues and forms images from the received echoes.

The SONON 300C ultrasound imaging system with accessories are Class IIa Active Diagnostic Medical Products, according to MDD 93/42/EEC regulations for use on human patients.

The SONON300C ultrasound imaging system with accessories was developed and manufactured by HEALCERION. For more information, please contact HEALCERION.

HEALCERION

804ho, 38-21, Digital-ro 31-gil, Guro-gu, Seoul, 08376, Korea

Tel.: +(82) 70.7582.6326

E-mail: info@healcerion.com

EC REP

Obeliss.a

Bd. General Wahis 53

1030 Brussels, BELGIUM

Tel.: +(32) 2. 732.59.54

Fax: +(32) 2.732.60.03

E-Mail: mail@obelis.net

For USA only:

Caution: Federal law restricts this device by sale by or on the order of a physician.

Revision History

Revision	Date
User Manual Revision 0 - Initial release	2014-02-01
User Manual Revision 1 - Amendment : 9.2 Acoustic Output & Index Determination - Amendment : 9.3 Acoustic Output	2014-05-28
User Manual Revision 2 - Amendment : 6 Operating the Device - Amendment : 9.2 Acoustic Output Reporting Table - Amendment : Appendix A Usage Agreement	2014-10-27
User Manual Revision 2.1 - Amendment : 3.1 Symbol about RF equipment - Amendment : 9.1.3 Marking Plate - Amendment : 9.1.4 Guidance and Manufacturer's Declaration	2015-01-19
User Manual Revision 2.2 - Amendment : 6 Operating the Device	2015-04-06
User Manual Revision 2.3 - Amendment : 4.3 Battery - Amendment : Company Address - Amendment : 9.1.3 Marking Plate	2016-03-07
User Manual Revision 2.4 - Amendment : Typographical Error	2016-05-07
User Manual Revision 2.5 - Amendment : Typographical Error	2016-07-25
User Manual Revision 2.5, version 2 - Revised : 2.1 Intended Use - Revised : 4.1 Device description - Revised : 5 Installation of the Mobile App - Revised : 5.1 Mobile Device Requirement - Revised : 5.2 Installation of the Mobile App - Revised : 6.5 Patient Information	2018-12-14

Chapter 1 General

1.1. About This User Manual

Product name	SONON 300C Ultrasound Imaging System with accessories
Model	SONON300C

- Please read and understand all instructions in this User Manual before attempting to operate the SOA/OW300C equipment.
- Keep this User Manual with the product for future reference.
- Some options or features may not be available in some countries.
- The screen graphics and illustrations in this User Manual are for illustrative purposes only, and may differ from what is displayed on the screen.
- All references to standards and regulations, and their revisions are valid as of the date of publication of this User Manual.

1.2. Text Conventions

The manual uses the following conventions for Warnings, Cautions and Notes.

Warning

Brings attention to a condition or possible situation that could cause injury to the user and/or patient.

General precautions necessary to protect the patient/user's health and equipment.

Notice

Important information that must be read before proceeding.

Chapter 2 Introduction

100-100

2.1. Intended Use

The user should be careful to manage patient information and key security information because SONON 300C is used with personal mobile device.

Destination: The SONON Ultrasound Imaging System (Model: SONON 300C) is indicated for ultrasound echo imaging, measurement, and analysis of the human body for general clinical applications including obstetrics(OB), gynecology (GY) and general(abdominal) imaging.

Scope of application: Its user interface and portability enable it for use in primary care and special care areas, such as medical office settings, clinics and emergency room in hospitals.

This device is intended for passive observations and early disease screenings. If a more thorough examination of a pregnant woman or the judgment of the doctor is required, we recommended that you send the patient to a hospital for a thorough examination.

If any abnormalities are detected, a more thorough examination and a consultation with a specialist are recommended.

SONON Ultrasound Imaging System is intended for the purposes of ultrasound examination within screening procedure (special safe studies to identify the risk of any medical condition) in various areas of medicine including in external environment and some specific examinations and manipulations in emergency medicine, anesthesiology, trauma surgery, pulmonology, surgery and cardiology.

Special examinations and manipulations should be conducted by qualified professionals (having a legal proof of their specialization) and performed according to occupational standards implying presence of specific medical workers, relevant equipment and required medical supplies, which can only be performed in appropriate general health institutions.

Use of SONON Ultrasound Imaging System in screening is intended for examination of stabilized patients without making any diagnosis in order to determine further thorough specific examinations.

Use of the abovementioned Ultrasound Imaging System for non-stabilized patients in external environment is not recommended, as the examination of non-stabilized patients requires certain conditions regulated by occupational standards, and the information acquired by SONON Ultrasound Imaging System could be insufficient to choose the right therapeutic strategy

2.1.1. Intended Patient Population

Item	Description
Age	All ages
Weight	N/A
Gender	Males and females
Health	Do not use in patient who will be harmed caused by applying ultrasound (e.g., implanted pace-maker)
Nationality	Multiple
Patient state	Patient is not user(not relevant), only used by related experts

2.1.2. Intended User Profile

Item	Description
Education	Professionals who completed related training.
Knowledge	having basic knowledge on ultrasound
	having understanding on physiological effects of ultrasound.
	able to intuitively recognize and understand SWUI and HWUI
	able to distinguish mobile devices to which SW is applied
	able to understand updates.
Language understanding	able to understand terms used for the product and its manual
	able to understand the methodology in the manual (English version)
Experience	familiar to the ways to use mobile devices.
	able to set communication between the device and mobile devices.
	able to keep track of communication status.
	familiar to use apps of mobile devices.
	able to set and install SW updates.
	able to understand and use contents of SW UI.
	well trained about existing ultrasound.
Vision	corrected visual acuity 1.0 or better
Memory	read through the user manual and keep in mind precautions and key functions of the product
Permissible impairments	N/A

2.1.3. Operation principle

The SONON 300C utilizes pulsed-echo technology to determine the depth and location of tissue interfaces and measure the duration of an acoustic pulse from the transmitter to the tissue interface and back to the receiver. Ultrasound waves are emitted from the transducer, propagate through various tissues, and return to the transducer as reflected echoes.

The returned echoes are then converted into electrical impulses by transducer crystals then further processed in order to form the ultrasound image presented on the screen.

The echo signals are amplified and processed by several analog and digital circuits with filters for many frequency and time response options, transforming the high-frequency electrical signals into a series of digital image signals which are stored in memory.

The ultrasound image is based on the mechanical oscillations of a crystal excited by electrical pulses (the piezoelectric effect). Numerous crystals are assembled to form a transducer. A transducer converts one type of energy into another. Based on the pulse-echo principle occurring with ultrasound piezoelectric crystals, ultrasound transducers convert electricity into sound.

The SONON 300C Ultrasound Imaging System is single mode (B Mode), single output (3,5 MHz – low ultrasonic powered) device. The SONON 300L Ultrasound Imaging System a device with operating B mode and CF mode (5 MHz / 7,5 MHz / 10 MHz – low ultrasonic powered).

B Mode is a two dimensional mode (2D by sectors) which allows at a certain depth to get a planar image of tissue interfaces and their movement within a timespan. It is the easiest mode to comprehend, as it helps to display anatomical formations like on a cross-sectional cut. This is the main mode to use and the first one to start with during any ultrasonography. The image creation in B mode is based on reflection of ultrasound waves by the tissue interfaces of the examined area: the Ultrasound Imaging Device calculate the amplitude and the phase of the returned signal, which information is used to build an image (that will be displayed in the screen) and determines the pixel brightness and its position on the diagnosis image.

CF mode (duplex Doppler ultrasound) is designed for:

-measuring the quality of blood circulation – to define the blood behavior: laminar (even flow) or turbulent (multiple swirling)

-measuring the blood quantity – to define the blood speed in the vessel. The blood flow characteristics (laminar or turbulent) are encoded by different colors which intensity is determined by the flow speed. Color Doppler ultrasound examination facilitates the diagnosis and helps to reduce the risk of error, especially for cardiac anomaly: the colorful flow is highly illustrative.

2.2. Contraindications

Warning

The SONON 300C is not intended for:


- ophthalmic use or any use causing the acoustic beam to pass through the eye;
 - intra-operative use
- (e.g., defined as introducing a probe into a surgical incision or burr hole);
- endocavity use; or
 - Imaging an open wound.

Possible side effects


There were no side effects.

2.3. Terms of use

The SONON 300C scanner is suitable for use in all medical institutions (i.e. hospitals, clinics, etc.). The scanner SONON 300C is intended for use by specialists only.

Chapter 3 Safety

3.1. Symbols

Some symbols used with electrical medical equipment have been accepted as standard by the IEC. They serve for marking connections, accessories, and also serve as warnings.

Symbol	Description/Function
	Insulated patient application (Type BF)
	Warning: This symbol indicates hazard about Operation condition For more information, see on Section 3.2 "General Safety"
	Consult accompanying documents: This symbol advises the reader to consult the accompanying documents.
	Freeze mode: This symbol means that image is frozen.
	Stand-by: This symbol's function is to turn the device on or put it in standby mode.
	Wi-Fi: This symbol means wireless communication.
	Manufacturer information: This symbol is followed by the name and address of the device manufacturer.
	Manufacture Date: This symbol is followed by the device manufacture date in the form YYYY-MM
	Serial Number: This symbol is followed by the device serial number.
	Model name: This symbol represents the model name of product
	Representatives information: The manufacturer's EU representative information shows with this symbol
	The CE marking is a marking of the conformity of products to the quality requirements in the territory of the European Union. The CE mark ensures that the products manufactured by the manufacturer fully comply with the safety requirements for humans and the environment
	This way up, Indicates correct upright position of the transport package
	Fragile, handle with care
	Keep dry
	Indicates the need for separate collection for electrical and electronic equipment in compliance with the Waste Electrical and Electronic Equipment (WEEE) Directive.
	This symbol means "Caution: Federal law restricts this device to sale by or on the order of a physician"
	Electromagnetic radiation. There may be interference near equipment marked with this icon
	Temperature limitation

Humidity limitation

Atmospheric pressure limitation

3.2. General Safety

Warning

- Only authorized personnel shall perform any type of repair on the SONON 300C. Never attempt to open a transducer or transducer connector. This will void the warranty!
- Probes are not delivered sterile! Before the first use, it is MANDATORY to clean and disinfect probes to avoid infections or disease transmissions!
- Probes must be cleaned and disinfected before they are replaced or disposed of.
- Do not modify this device without authorization of HEALCERION.
- Do not use the probe with high frequency surgical equipment. Doing so may damage the equipment.
- If you lose your tablet PC or smart phone, the stored data cannot be recovered. Please take care not to lose these items.
- Do not touch the patient with iOS or Android mobile device when you use the 300C
- Must comply with Operation condition; **max. 10 min. with 10 min. resting time.**

Caution

- Do not touch the patient and the charging connectors simultaneously.
- This equipment should be used in compliance with applicable laws. Some jurisdictions restrict certain uses, such as gender determination.
- Changing the display settings can affect image quality and compromise diagnostic quality. The user is responsible for using appropriate display settings for achieving appropriate image quality.
- Ultrasound procedures should be performed prudently using the principle of ALARA (As Low As Reasonably Achievable). It is strongly recommended to consider ALARA when conducting ultrasound scans. See Section 3.6 "Bioeffects and Safety of Ultrasound Scans" for additional information.
- Features that facilitate measurements must be used with extreme care. Such measurements are suggestions of the system. If in doubt, verify the measurement results with manual measurement methods. The user is responsible for the diagnostic interpretations of the measurements.

Notice

- Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use care when handling and protect the equipment from damage when not in use.
- DO NOT use a damaged or defective probe. Failure to follow these precautions can result in serious injury and equipment damage.
- Transducer damage can result from contact with inappropriate couplings or cleaning agents.
- Do not soak or saturate the transducer with solutions containing alcohol, bleach, ammonium chloride compounds, hydrogen peroxide, or any incompatible solutions.
- Inspect the probe prior to use for damage or degeneration to the housing, strain relief, lens, or seal. A thorough inspection should be conducted during the cleaning process.
- If a probe has been dropped on the floor or on any other hard surface, do not use the probe any more. That may increase the risk of electric shock due to damaged electrical insulation.

110157-1

107

3.3. Electric Safety

The probe is driven by electrical energy that can harm patients and users if live internal parts come in contact with conductive solutions:

Warning

- DO NOT immerse the probe into any liquid beyond the immersion level. Never immerse the probe connector into any liquid.
- DO NOT drop the probe or subject them to other types of mechanical shock or impact. Degraded performance or damage such as cracks or chips in the housing may result.
- Electrical leakage checks should be performed on a routine basis by qualified hospital personnel.

3.3.1. Immersion Level

Caution

- Probes labeled "IPX7" are watertight up to a maximum of 4cm below the probe.

3.4. Electromagnetic Compatibility

- The system must remain the essential performance of the ultrasound system in EMC environment according to IEC 60601-1-2. The essential performance refers to section 9.1.4 "Guidance and Manufacturer's Declaration" of this manual

3.5. Coupling Gels

Caution

- Do not use non-recommended gels (lubricants). They may damage the probe and void the warranty.

Applying In order to ensure the optimal transmission of energy between the patient and probe, a conductive gel or couplant must be applied liberally to the patient where the scanning will be performed.

Precautions Coupling gels should NOT contain any of the following ingredients, as they are known to cause probe damage:

- Methanol, ethanol, isopropanol alcohol, or any other alcohol-based products
- Mineral oils
- Iodine
- Lotions
- Lanolin
- Aloe vera
- Olive oil
- Methyl or ethyl parabens (para hydroxybenzoic acid)
- Dimethylsilicone

The coupling gels are recommended that use the followings.

Name	Manufacturer
Aquasonics 100	Parker Laboratory Inc.
Clear Image	Sonotech Inc.
Scan	Parker Laboratory Inc.
Sonogel	Sonogel Vertriebs

3.6. Bioeffects and Safety of Ultrasound Scans

Thermal Safety

Maintaining a safe thermal environment for the patient has always been a design priority at Healcerion. The operating temperature of the ultrasound probe must remain below 43°C.

Whenever ultrasound waves travel through tissue, there always is a certain risk of damage. There has been a great deal of research on the impact that high frequency waves can have on different kinds of tissues under defined conditions and, "There is, to date, no evidence that diagnostic ultrasound has produced any harm to humans - including the developing fetus." (Guidelines for the Safe Use of Diagnostic Ultrasound Equipment, Safety Group of the British Medical Ultrasound Society 2010).

Physiological effects due to ultrasound are generally assumed to be deterministic and only occur above a certain threshold, in contrast to ionizing radiation, which causes effects accidentally. Thus, ultrasound examinations can be conducted very safely if certain procedures are followed. Therefore, it is recommended to read the following sections and study the cited literature.

3.6.1. Prudent Use

Despite the relatively low risks of ultrasound scans, as compared to other imaging techniques, the operator must choose exposure levels with caution to minimize the risks of bio-effects.

"A fundamental approach to the safe use of diagnostic ultrasound is to use the lowest output power and the shortest scan time consistent with acquiring the required diagnostic information. This is the ALARA principle (i.e. As Low As Reasonably Achievable). It is acknowledged that in some situations it is reasonable to use higher output or longer examination times than in others: for example, the risk of missing a fetal anomaly must be weighed against the risks of harm from potential bio-effects. Consequently, it is essential for operators of ultrasound scanners to be properly trained and fully informed when making decisions of this nature." (Guidelines for the Safe Use of Diagnostic Ultrasound Equipment, Safety Group of the British Medical Ultrasound Society 2010)

Special care regarding ALARA should be taken with obstetric examinations, as any potential bio-effects are likely to be of greatest significance in the embryo or fetus.

It is strongly recommended to consider ALARA when conducting ultrasound scans. (See Appendix B.)

The SONON Ultrasound Imaging System (Model: 300C) is B Mode and CF Mode, output (5MHz / 7.5MHz/10MHz- low ultrasonic powered) device.

The factors that operator can control are below and the output value is less than 1. (MI and TI)

- TGC (Time Gain Compensation)
- DR (Dynamic Range)
- FL(Focal Length)
- Depth
- Image Filter
- Line Density

Operated B-mode, The SONON300C Ultrasound Imaging System is safe because exposure time to patient applied part is not concentrated in one place. In Addition, Healcerion provide the generic content of the education program with an AIUM "Medical Ultrasound Safety" (AIUM 2014 - Refer to appendix.) and additional education. For example, the doctor who used the SONON Ultrasound Imaging System to verify the status of the baby prohibit inordinate inspection in order to show the family or other meaningless situation.

3.6.2. Bioeffects

• Thermal effects refer to the heating of soft tissue and bone

The thermal indices "TIs" (for soft tissue), "Tib" (for bone near the focus) and "Tic" (for bone near the surface) were introduced to provide the operator the relative potentials for a tissue temperature rise. According to the Standard for Real-time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment (2004) those thermal indices shall be displayed by ultrasound consoles. It should be noted that a TI of 1 does not necessarily mean that tissues being scanned will increase in temperature by 1°C. Almost every scanning situation departs from the assumed model conditions, such as tissue types, blood perfusions, modes of operation and actual exposure times of the scanned area. However, the thermal indices provide information regarding the possible increase in the risks of potential thermal bioeffects and it provides a relative magnitude that can be used to implement ALARA. In addition to tissue heating due to the generated ultrasound field, the temperature of the probe head itself can also increase during the examination. The operator should be aware, that in the tissue region near the ultrasonic transducer, there will be a superposition with the heating due to the ultrasound field, which is not considered by the TI values.

• Non-thermal effects refer to mechanical phenomena such as cavitation

Non-thermal bio-effects are caused by the interaction of ultrasound fields with very small pockets of gas (stabilized gas bodies), i.e., the generation, growth, vibration and possible collapse of microbubbles within the tissue. This behavior is referred to as cavitation (Medical Ultrasound Safety, 2nd Edition, AIUM 2009/American Institute of Ultrasound in Medicine Consensus Report on Potential Bio-effects of Diagnostic Ultrasound, AIUM 2008/Guidelines for the safe use of diagnostic ultrasound equipment, Safety Group of the British Medical Ultrasound Society 2010). The potential of cavitation increases with the rarefactional peak pressure, but decreases with the pulse frequency. Therefore the Mechanical Index (MI) was introduced to take account of both pressure and frequency. The higher the MI, the greater the potential for non-thermal bio-effects.

3.7. Environmental Conditions

Probes must be operated, stored, or transported within the parameters outlined below.

Item	Operational	Storage&Transport
Temperature	+18°C to +35°C (+64.4°F to +86°F)	-10°C to +50°C (+14°F to +122°F)
Humidity	30% to 75% non-condensing	10% to 85% non-condensing
Pressure	700 hPa (3000 m) to 1060 hPa	700 hPa (3000 m) to 1060 hPa

3.8. Other Recommendations

Like most high frequency computing devices, the electronic components of the 300C will generate some heat while operating normally and as intended. The 300C is equipped with safety mechanisms which will automatically reduce computing speed (frame rate), and ultimately shut down the device, before any risk of overheating occurs. The 300C is verified to comply with harmonized safety standards, under any operating condition described in this user manual (see Section 3.7 "Environmental Conditions"). To help keep the 300C operating temperature at the optimal level, and to ensure longer scanning times with the maximum frame rate, it is recommended to hold the 300C so that there is good contact between the device and the hand. The image display of the mobile device is dependent on the ambient light; avoid direct sunlight on the display when scanning and reviewing images.

3.9. Disposal

This symbol indicates that electrical and electronic equipment waste must not be disposed as unsorted municipal waste and must be collected separately. Please contact the manufacturer or an authorized disposal company to decommission your equipment according to local regulations.

Battery disposal

Caution

Lithium batteries are included with this device. Do not puncture, mutilate or dispose of batteries in fire. Replace only with the same type, as recommended by the manufacturer. Dispose of used batteries in accordance with manufacturer instructions and in accordance with your local regulations.

A separate collection symbol is affixed to batteries, or its packaging, to advise you that it must be recycled or disposed of in accordance with local and national laws. To minimize potential impacts to the environment and human health, it is important that all marked batteries that you remove from the product are properly recycled or disposed of. For information on how batteries may be safely removed from the device, please consult the equipment instructions, or your local authorities.

3.10. Transportation

The SONON 300C ultrasound imaging system in the package is transported by all modes of transport, in covered vehicles in accordance with the rules of transportation of goods, operating on each mode of transport.

Transportation and storage of the SONON 300C ultrasound imaging system without the packaging of the manufacturer does not guarantee its safety. Damage to SONON 300C as a result of transportation or storage without packaging of the manufacturer is eliminated at the expense of the consumer.

3.11. Warranty

The manufacturer warrants that there is no rejection of materials and the quality of the SONON ultrasound scanning system when properly used and maintained for 12 months (excluding the battery) from the date of purchase.


3.12. Packaging

The SONON 300C ultrasound imaging system is packed in a cardboard box.

Box dimensions: 323 × 255 × 71,5 mm (W × L × H).

Cardboard label: GZ (Coated SBB).

The box contains: scanner SONON 300C, battery charger SONON, two accumulator batteries (one inserted in the scanner, the other is spare), network adapter (for SONON) and power cord, placed into a separate cardboard box, and AC adapter for the SONON battery charger with the power card, also placed into a separate cardboard box, and user manual.

Box, top view


Box, bottom view

Box, inside

Then ten such boxes are placed into a big corrugated box for transportation.

3.13. Service Life

Expected service life: 3 years

Chapter 4 Device Description

2019


4.1 Product Description

4.1.1. Package Contents:

Scanner SONON 300C	User Manual SONON 300C
	
Accumulator battery - 2 pcs.	Battery Charger SONON
	
Adapter (for Battery Charger) SONON	Power Cord - 2 pcs.
	
Network adapter (for SONON 300C)	Brief manual SONON 300C
	
Software Usage Agreement	
Accessories:	
Protective cap for SONON	Hand strap for SONON
	
Mounting for SONON	
	

4.2 Screen Layout


4.2.1. Outline of apparatus

No.	Name	Function
1	Probe Head	To be applied to the patient
2	Function button	Blue LED is turned on during the freeze screen. Blue LED is flashing during the live screen.
3	Power LED	Turning on the power supply (White) Turning off the power supply (Green LED is displayed for 2 seconds, and then, the LED is turned off.) Charging in progress (Flickering with green color)
4	Power button	Power On /Off
5	Wi-Fi button	Wi-Fi Indicator (Blue)
6	Charging device connected	This is the section where AC adapter is connected for charging the battery of apparatus.

25

4.2.2. Basic Layout

No.	Name	Function
1	Indicator Area	Mode and Major Function Indicator
2	View Area	Display Diagnostic Image and Information
3	Menu Area	Function Menu Button

4.3. Battery

Caution

- Do not place the battery near a heat source or expose it to direct flames. Such exposure may lead to corrosive liquid leakage, electrical shock or fire.
- If any liquid from the battery should come in contact with an eye, immediately wash the eye with plenty of water and seek medical assistance as soon as possible.
- Do not immerse or expose the battery to water.
- The AC adapter must be kept outside the patient environment, (refer to IEC 60601-1)


Figure 1 Patient environment

The 300C is powered by a Lithium Ion battery. The battery is not fully charged prior to shipment. To maximize battery life, it is recommended that the battery be fully charged before initial use. Use only the AC Adapter provided with the 300C.


4.3.1. Battery Level Indicator

Notice

- The battery should not be charged in the Power On state. To charge the battery, power should be turned off.
- We recommend to use battery charger and adaptor because charging by adaptor for 300C do not fully charge to protect the 300C device.
(About 80% charging by adaptor for 300C)
 - Condition Operation: 10 minutes
 - Time to Fully Charge: 3 Hours

The battery level indicator is displayed by color and flash rate.


Battery Level	Description
Charging is going on	Battery level indicator lights green
Fully charged	Battery level indicator is off
Battery low	Battery level indicator flashes


4.3.2. Removing and Inserting the Battery

4.3.2.1. Battery changing Process

(1) Slide down the lock button


(2) Slide the battery cover toward the lock button


(3) Remove the battery cover


(4) Remove the battery

(5) Insert the new battery

(6) Push the battery cover to the inside while lowering the locking button.

4.3.2.2. Battery charging process


(1) Plug the adapter cable into the adapter

(2) Turned the green light after connected the cable

(3) Plug the cable into the Battery Dock

(4) Place the battery on battery dock and push to insert the battery into battery dock.


(5) Slide the battery into the battery dock.

(6) During the charging, red indicator lamp is lighted, and when completed, green indicator lamp is turned on.

(7) Press the button of battery charger

(8) Separate the battery pack


4.3.2.3. Device charging Process


(1) Open the side cover on the device

(2) Plug the USB cable into the device

(3) Confirm whether the green indicator lamp is flickered, and when the charging is completed, the indicator lamp is turned off.

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103101

4.3.3. Battery Specifications

Item	Specification
Description	Rechargeable Li-ion Battery Pack
Capacity	2600 mAh (7.4 V)
Lifetime Guarantee	6 months
Model	MBP-2S1PSD26
Cell Type	ICR18650
Dimensions	L*W*H (L)±0.2* (W)±0.2 * (H)±0.2mm
Safety	PCM Logic

Notice

For information about purchasing a battery, please contact HEALCERION.

Chapter 5 Installation of the Mobile App

30

To operate SONON 300C ultrasound imaging system you need to install SONON Ultrasound App. It is impossible to log into the mobile application without SONON 300C devices (just as the scanner cannot work without the mobile application installed)! The access code / access key (or wifi connection between the mobile device with pre-installed software and the SONON 300C scanner) is a serial number of the scanner indicated on the identification label located on the back of the product. To activate the program, enter the login - the serial number (SN) of the scanner SONON 300C, company name, email and set your own password. Thus, you create your account which you will use as the only way to access the software.

Software Features

- displaying on your phone or tablet the ultrasound images made by SONON 300C scanner;
- processing (converting) the images to be displayed
- displaying the processed (converted) image on the phone/tablet
- managing diagnosis features
- sensor-based control of the main program settings
- saving the diagnosis information (image, measures, user notes/comments and patient information) and managing previous records (edit, copy, delete, send by email)
- measuring of the received ultrasonic data.
- managing scanner settings: wifi connection settings with possibility to modify «default password», upgrading the scanner firmware
- managing transfer and backup of patient information.

5.1. Mobile Device Requirement

Warning

- If you remove and reinstall the app, all stored information will be deleted and cannot be recovered.
- Intended for use with standard Android (v. 4.1 or later) and iOS operating systems (v.7.0 or later). Android Rooting and iOS Jailbreaking removes the limitations on the manufactured devices running the standard OS, and thus the 300C is not intended for use with these devices. Healcerion is not responsible for any error in usage.
- Do not use anything other than a mobile device that recommends.

Caution

After installing the 300C App from the App Store / Google Play, the new version is recognized and it is automatically updated. At this time, firmware also is automatically upgraded.

Notice

When you operate the 300C, HELCERION needs your personal information. Please read 'Appendix A. Usage Agreement'.

Personal information protection and protection of the device from an unauthorized access is performed as follows:

1. The application by Healcerion downloaded from Apple Store or Google Play could not be used without registering it.
2. The registration of the software by Healcerion implies the following:
 - 2.1. Entering of the access code. The access code is the scanner serial number (SN) indicated on the identification label on the back of the product. Thus, one cannot use the application if he/she does not have any device.

2.2. In addition to the access code the user (doctor) should enter the company name, email, user password. Thus, the user account will be created, and further access to the software will only be available after user authorization.

HEALCERION

Fill in your information

Organization


Serial Number

Password

Password

Confirm password

Email address

3. Once the application is registered you should connect the SONON device to the installed software to SONON scanner using Wifi interface, primary password indicated in the User manual. The connection allows using only one device (application) with each SONON scanner simultaneously. The attempt to connect another device (application) at the same time will display a "connection error" message and the application on the device will be unavailable.

The minimum requirements for the mobile device are listed in the table below.

Item	Minimum requirements
Target Device	iPad mini Series / iPad Air Series iPhone 5,5S, 6,6 plus Galaxy Note 10.1 Wi-Fi/LTE Galaxy S5, Galaxy Note 3 Galaxy Tab 8.4 Pro
Mobile OS Version	iOS 8.0 or later Android 4.4 or later
Minimum device storage available for mobile app installation	8 Gb

5.2. Installation of the Mobile App

Notice

For iPhone, installation of the 300C mobile app requires an iTunes registration and password.

5.2.1. iOS


- (1) Launch the "App Store" App
- (2) Search for *SONON 300C*
- (3) Select the *SONON 300C* App and Install
- (4) Enter the user's Apple iTunes password
- (5) Automatic installation

5.2.2. Android

- (1) Launch the "Google Play" App
- (2) Search for *SONON 300C*
- (3) Select the *SONON 300C* App and Install
- (4) Automatic installation

Software version is M1.XX.XX. After installing the SONON Ultrasound App from the App Store / Google Play, a new version will be found, which will be upgraded automatically. It will also upgrade the firmware automatically.

The SONON 300C system is designed to work on standard Android OS (v. 4.1 or higher) and iOS (v. 7.0 or higher).

Chapter 6 Operating the Device

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- Please check for any physical damage of the exterior, including cracks and foreign materials before using the product.
- The system should be cleaned/disinfected per instructions in section 7.3 prior to use.

6.1. Device Start-Up

6.1.1. Power On


- (1) Press the power button for 1-2 seconds.
- (2) Hold the power button until the green light is switched to white color.
- (3) When the LED of Wi-Fi button is turned on to green color, it represents that the device is ready.

6.1.2. Power Off

Press the power button for 1-2 seconds (until the LED of power button is lighted in white color). The power LED is change from white to green. When the power LED is changed from green to turning off, it means that the ending of apparatus operation is completed.


6.2. Network Interface Technical Specifications (Wi-Fi)

- 2.4 GHz
- SoftAP Mode
- SSID: SONON300C- Serial Number (XXXXXX)
- Password: 1234567890

6.3. Probe orientation (Line 1 Marker)

The probe is provided with an orientation marking. This mark is used to identify the end of the probe corresponding to the side of the image having the orientation mark on the scanning screen.

※ The basis of line 1 Marker is based on the protruding part of 300C product's probe head section.

- (1) Select the line 1 Marker.
- (2) The line 1 Marker is displayed at the opposite side.
- (3) The screen is reversed.


6.4. Operating Procedure

6.4.1. User Interaction

Action	Interaction	Icon
Single Tap	Single	
Double Tap	Double Tap	
Drag	Drag (Any Direction)	
Pinch in	Zoom Out	
Pinch out	Zoom In	

6.4.2. Registration in the App

6.4.2.1. Primary adjustment screen

- (1) When pressing the next button under the state of not filling the input values of all fields, popup 1 is outputted. Confirm the changed details.

Please complete all fields.

OK


Popup 1

- (2) If the PW value is not matched with confirm PW value, print Pop-up below.

Please re-input your password and re-confirm the password.

OK

Popup 2


6.4.2.2. App Log-in

- (1) Select OK after inputting the password set up during the subscription.
- (2) Moved to the APP standby screen.

6.4.2.3. Find Password


The value is not correct. Please try again.

OK

Popup 4

- (1) Select Find password.
- (2) If the serial number or E-mail address is different from the saved values, popup 4 is outputted.
- (3) When the serial number and E-mail address established during the member subscription are inputted, the password can be confirmed.

6.4.3. Start Scan

※ * When 300C is connected to the mobile application and if the user taps the start button, it is available to start the scan.

- (1) Select the start scan button.
- (2) The progress icon is displayed.
- (3) After connecting the 300C and mobile application, the information of 300C is normally outputted.
- (4) When 300C is not connected to the mobile application, the user can't use the scan functions.

Not connected with SONON.
Please try again.

OK

6.4.4. Live Screen


6.4.4.1. Live screen description- Phone (Live Default Screen)

No.	Name	Function
1	TGC button	Adjusting the TGC(Time Gain Compensation) value for each depth by utilizing the TGC and 4 slides
2	DR button	Adjusting the DR (Dynamic Range) value by utilizing the DR and 4 slides
3	Freeze switching button	Stopping the image during the scanning, or activating again the stopped image
4	Image save button	Saving the image to the play list (database) during the scanning
5	Patient information button	Selecting the patient in the list or inputting the new patient information
6	Page move button	Moved to the previous page


No	Name	Function
7	File list button	Moved to the saved file list (Play List) screen
8	Setup button	Setting menu is displayed.
9	Probe direction display	Displaying the reversed image
10	Grey bar	Displaying the brightness of image
11	FL display	Displaying the current FL (Focal length) value by locating at the side scale area, and setting up the FL value to the target section to be observed
12	Side scale area	Adjusting the value of depth by dragging this area, or setting up the FL value by selecting the FL indicator.
13	Network connection state display	Green: Displaying that the network connection state is stable Yellow: Displaying that the network connection state is unstable
14	Battery remains display	Displaying the battery remains of <i>SONON 300C</i> main body
15	Acoustic output power display	MI: Displaying the mechanical index value TIB: Displaying the bone thermal index value TIS: Displaying the soft-tissue thermal index value TIC: Displaying the cranial-bone thermal index value
16	Frame rate display	Displaying the number of frames to be transmitted
17	Patient information area	Displaying the patient information saved in the patient list Displaying with the 'Name'-'ID'

6.4.4.2. Zoom Action

- ※ Pinch Out (Zoom in) Action (max 400%)
- ※ Pinch in (Zoom out) Action (min 100%)

6.4.4.3. FL Control


- ※ In Zoom Mode, FL Control function does not work.

- (1) Single tap FL position in side ruler area
- (2) Display New FL marker
- (3) The MI /TI value in accordance with the change of FL value is displayed.

- ※ FL value cannot exceed the value of depth, and maximum FL value shall be set up as Depth-1.
- ※ The range of FL change is 2 cm~16 cm.

6.4.4.4. Depth Control

※ In Zoom Mode, Depth Control function does not work.


(1) Drag the side ruler area.

(2) Display New Depth

- ※ Drag up - increase Depth / Drag down - decrease Depth
- ※ If the depth is reduced by user to less than the FL value, the FL is automatically changed to the finally changed value of depth -1.
- ※ e.g When the depth before the change is 10 cm, and FL is 7 cm, if the depth is changed to 5 cm, the FL is automatically changed to 4 cm.
- ※ The range of depth change is 5 cm~20 cm.

2021

6.4.4.45. TGC Control

※ During Zoom in action, the TGC control function does not work.


- (1) Tap TGC Button.
- (2) Display TGC slide bar.

- ※ Divide the entire length into four sections.
- ※ Drag slide bar to left: Decreases the brightness of the section.
- ※ Drag slide bar to right: Increases the brightness of the section.

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

6.4.4.46. DR Control

- (1) Tap DR Control Button.
- (2) Display DR slide bar.
- (3) Drag slide bar by user action..
 - ※ Drag Right to increase DR XX dB (Max 96 dB).
 - ※ Drag Left to decrease DR YY dB (Min 30 dB).
- (4) Display DR value (unit dB).
- (5) Tap DR Control Button.
- (6) Disappear DR slide bar..


6.4.4.47. Recording


- (1) Select the recording button.
- (2) During the recording, the 'Recording...' message is printed onscreen.
- (3) Select the recording button.
- (4) The recording end message, 'Recording file is saved' is displayed, and the diagnosis data (Recording file) is saved.

6.4.5. Freeze Screen


- ※ While scanning, press Function button or tap freeze button to freeze the image.

6.4.5.1. Freeze screen description (Freeze Screen Layout)

No	Name	Function
1	Clear button	Delete object of measurement
2	Measurement button	Measurement length, ellipse and add comment
3	Live button	Change to Live mode
4	Still-Shot button	Save displayed Image (sill-shot)
5	Patient List button	Select patient of list and add new patient information
6	Rewind Bar	Check previous diagnosis frame

6.4.5.1.1. Length


- ① Single tap the Measurement Button.
- ② The Sub Menu will be displayed.
- ③ Tap the Length Button.
- ④ Draw line where to measure length and yellow line is displayed.
- ⑤ Display Length Value on the screen.

- ※ User can move the position icons with the drag action.
- ※ The maximum number of measurement objects is 10.
- ※ Below the pop-up will be displayed once the object limit has been reached.

Maximum Measurement object limit is 10. You can not add any more objects.


OK

6.4.5.1.2. Ellipse


- ① Single tap the Measurement Button.
- ② The sub menu icon will be displayed.
- ③ Tap the Ellipse Button.
- ④ Drag line where to measure ellipse and displayed Ellipse frame (Yellow).
- ⑤ Double tap the middle of yellow point of ellipse frame, it changes to Green and displayed position icon.
- ⑥ Displayed Ellipse value on ellipse frame.
- ⑦ Displayed Ellipse value on the screen.

6.4.5.1.3. Comment


- ① Single tap the Measurement Button.
- ② Display the sub menu.
- ③ Tap Comment Button.
- ④ The virtual keyboards will be displayed.
- ⑤ Input comment and tap the send of virtual keyboard.
- ⑥ Display user's input note.
- ⑦ Double tap the comment and change to green and revision is available.

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
6.4.5.1.4. Fetal Biometry

- ① Select the Measurement button.
- ② Fetal Biometry sub-menu item is displayed.
- ③ Select the CRL button among the Fetal Biometry menus.
 - CRL: Length (The length from head to hips)
 - AC: Ellipse (The abdomen circumference of fetus)
 - BPD: Length (The head diameter of fetus)
 - HC: Ellipse (The head circumference of fetus)
 - FL: Length (The length of thigh bone)
- ④ After choosing the corresponding item, draw the line by selecting the measuring section.
- ⑤ The length currently measured is displayed at the top left section.
 - ※ The length of item and the predicted number of column, Example_ CRL: 9.0cm 14w6d.
 - ※ This menu is shown only OB mode.
 - ※ GA calculation algorithm is 'Hadlock et al' method.

6.4.5.6. Delete


6.4.5.6.1. Delete Measure

- ① Single tap the Clear icon button.
 - ② Display submenu.
 - ③ Single tap the Delete Measure Button.
 - ④ The measurement values of Length, Ellipse, and Fetal Biometry which were recorded are disappeared.
- ※ The latest measurement object is deleted.
 - ※ Deleting each object is not available.
 - ※ When selecting the Delete Measure, the Comment is not deleted.

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
6.4.5.6.2. Delete Comment

- ① Single tap the Clear icon button.
- ② Display submenu.
- ③ Single tap the clear comment Button.
- ④ The Comment which was recorded is disappeared.

※ When selecting the Delete Comment, the Measure is not deleted.


6.4.5.6.3. Delete All

- ① Single tap the Clear icon button.
- ② Display submenu.
- ③ Select the Delete All.
- ④ All measurement and comment objects are deleted.

6.4.5.55. Saving Diagnosis File

6.4.5.7.1. Saving Diagnosis File (Previous Frame)


(1) Single tap the saving icon button.

(2) The message 'Still-shot image is saved' is displayed for 2 seconds, and the image is saved.

6.5. Patient Information


6.5.1. Patient List

The location of Patient List can be accessed by selecting the 'Patient List' from the Start Screen, or selecting the 'Patient List Icon (Human shape)' from the Live, Freeze screen.


- (1) Tap the patient info icon- Same as the Freeze Mode Patient Info. Button.
- (2) Move to Patient list Screen.

6.5.1.1. Add new patient


- ① Tap Add new patient button.
- ② Input each item at the Patient Info screen.
- ③ Select the Save.
- ④ Popup window is displayed, and select OK.
- ⑤ When selecting the OK from popup window, the patient information inputted is saved.

※ Essential input items for patient information: ID and Name.

6.5.1.2. Patient List Edit

- ① Tap edit button.
- ② Tap selected patient on patient list.
- ③ Tap the activated Delete.
- ④ Print popup for caution about deleting patient.
- ⑤ Selected patient and related diagnosis files are deleted on patient list screen.

Caution!
Selected patient information and all
diagnosis data cannot be recovered
after deletion. Do you still want to
delete the patient information?

OK Cancel

1610101


1610101

6.5.1.3. Playlist


- (1) From the patient list screen, select the patient whose diagnosis data is to be reviewed.
- (2) The patient information screen of the patient chosen is displayed, and select the Playlist.
- (3) From the Playlist screen, the diagnosis data of the patient selected is displayed.

6.5.1.4. Playlist Edit-Delete


- (1) From the Playlist screen, the diagnosis data of the patient selected is displayed.
- (2) After choosing the image whose diagnosis file is to be deleted, select the delete.
- (3) Warning popup window is appeared. Select OK. ⑧ Guide popup window is appeared. Select OK.
- (4) From the Playlist screen, the diagnosed data of patient is displayed. Deleted diagnosis data is not displayed.

6.5.1.5. Playlist Edit-Move

- ① From the Playlist screen, the diagnosed images of patient are displayed.
- ② Select the diagnosis data which is to be moved to another patient.
- ③ Select the Move.
- ④ From the Patient List screen, select the patient whose diagnosis data is to be moved.
- ⑤ Guide popup window is appeared. Select OK.
- ⑥ The Playlist screen of the patient to be moved is displayed.
The diagnosed image list is displayed.
The diagnosis data for which the moving is selected is added.

6.5.1.6. Playlist Edit-Rename

- ① From the Playlist screen, the diagnosed images of patient are displayed. Select the Edit.
- ② Select the Rename of the diagnosis file whose name is to be reset.
- ③ Input the name of the file to be changed. After the inputting, select OK button.
- ④ From the Playlist screen, the diagnosis data of patient is displayed, and the changed name is reflected to be seen.

6.5.1.7. Playlist Edit –Email

From the Play list screen, select the file to be transmitted via e-mail.

- ① When selecting the sharing icon from the regeneration screen of corresponding file, the sub-menu is appeared. Select E-mail.
- ② When connected with SONON AP, No.1 popup window is displayed, and select OK.
- ③ Diagnosis data is automatically attached. When mail is prepared, select the Send button.

Caution


When connected to the **SONON 300CAP**, the user can't use the E-mail Client App.
If the user wants to send an E-mail, Go to setting app and disconnect with **SONON 300CAP**.
Use a cellular network or another Wi-Fi network.

6.5.2. Live Screen as followed in Playlist

- (1) From Patient info screen, select the Start Scan button.
- (2) After moving to the Freeze screen, the Initializing popup window is displayed.
- (3) The Live screen is displayed, and the state of scanning is ready.

6.5.2.1. Diagnosis File Viewer


- (1) Tap the Diagnosis Image.
- (2) The Image Viewer will be displayed.
- (3) Tap the Diagnosis Cine File.
- (4) The Cine File Viewer will be displayed.

After examination, the image (photo or video) saved on smart phone or tablet can be sent to another telephone (or PC) in order to archive it (storage will allow its further analysis to establish diagnosis or monitor patient's condition). This information will be analyzed by ultrasonography professional or consulting physician by visual estimation of the image and recorded measures. The report will be made by the ultrasonography professional based on the image analysis. He/she will provide a specific image description and conclusion (following the ultrasound examination procedure) that will be added to the medical history sheet.


6.6. Settings

6.6.1. Settings Menu

- ✱ Setting menu can be selected from the Live, Freeze, and start screen.


- (1) From the Start screen, select the setting icon.
- (2) The Settings screen is appeared, and the 5 items, i.e., Preferences, Export, Preset, Backup, and Import are displayed.

- (1) From the Live, Freeze screen, select settings. icon
- (2) As the sub-menu, and when the Settings button is selected, the 5 items, i.e., Device Settings, Image Option, Preferences, Backup and Import are displayed.

6.6.1.1. Device Settings

- ✳ Display Firmware Version and WIFI - WIFI SSID, Password.
- ✳ If firmware update is needed, firmware update button will be displayed automatically.


- ① From the Settings screen, select the Device Settings button.
- ② From the Device settings screen, select the Change for changing the Wi-Fi password of **SONON 300C**.
- ③ Select the SSID / Key field. Keyboard is appeared, and it is activated so that the SSID and Password can be changed and set up.
- ④ After the changing, select the 'Save'.
- ⑤ On the Device Settings screen, the changed details are displayed.

61


6.6.1.2. Image Option

- ※ Frame Average: Random Noise alleviating function.
- ※ Image Enhancement: Speckle Noise alleviating filter.

6.6.1.3. Preference


- ※ Thermal Index - Setting up the type of TI value to be displayed on the screen (For MI value, the basic value is displayed).
- ※ Recording File Size - Recording file save unit (Capacity).

- ① Input the desired size and select the Preferences button from the Settings screen.
- ② Select the capacity of the recording file to be saved to the Recording File Size.
- ③ On the Device Settings screen, the changed details are displayed.


6.6.1.4. Backup

- ✱ When selecting the backup after choosing the patient to be backed up, the data of corresponding patient is backed up.
- ✱ The file backed up is saved to the 'Backup Data' folder of Healcerion sub-folder.


6.6.1.5. Import Backup File

- ✱ When selecting the patient with backed up file, the data of the selected patient can be imported.

- ① From the Settings screen, select the Import backup file.
- ② On the Import screen, total capacity and the backed up file are displayed.
- ③ Select the backed up file to be imported again.
- ④ Popup window is displayed, and select OK.

6.6.1.5.1. Import Backup File-Delete

- ① From the Settings screen, select the Import backup file.
- ② After selecting the patient to be deleted, select the delete.
- ③ Popup window is displayed, and select OK.
- ④ Among the Import patient list, the backup file selected is not displayed.

Chapter 7 Maintenance

11111111

Caution

- The user must ensure that safety inspections are performed at least every 12 months according to the requirements of the patient safety standard IEC 60601-1 / UL60601-1.
- Only trained persons are allowed to perform the safety inspections mentioned above.
- The SONON 300C requires regular care and maintenance to function safely and properly.
- To ensure that the SONON 300C operates continuously at maximum efficiency, we recommend that the following procedures be observed as part of the customer's internal routine maintenance program.

7.1. Inspecting Probes

After each use, inspect the probe's lens and casing. Look for any damage that would allow liquid to enter the probe. If any damage is found, the probe must not be placed into any liquid (e.g., for disinfection) and must not be used until it has been inspected and repaired/replaced by a HEALCERION Service Representative.

7.2. Planned Maintenance

Do the Following	Daily	After/Before Each Use	As Necessary
Inspect the Probe head	✓	✓	✓
Clean the Probe head		✓	✓

Caution

- Probes labeled "IPX7" are watertight up to a maximum of 4cm below the probe.

7.3 Cleaning & Disinfecting

Recommendations for cleaning the ultrasound probe:

- Remove all coupling gel and other visible substances from the probe by wiping with a soft dry cloth. If necessary to remove material dried to the surface, the cloth can be moistened with lukewarm water.
- After each use, inspect the probe's lens and casing. Look for any damage that would allow liquid to enter the probe. If any damage is found, the probe must not be placed into any liquid (e.g., for disinfection) and must not be used until it has been inspected and repaired/replaced by a HEALCERION Service Representative.

Recommendations for disinfecting the ultrasound probe (After cleaning):

- (1) Spray disinfect onto the surface of probe head.
- (2) Repeat step one for two or three times.
- (3) Let stand for about 1 minute.
- (4) Wipe out the disinfectant with a clean cloth.

Recommended disinfectants


- In order to provide users with options in choosing a disinfectant, HEALCERION routinely reviews new medical disinfectants for compatibility with the system. Although a necessary step in protecting patients and employees from disease transmission, liquid chemical disinfectants must also be selected to minimize potential damage to the transducer.

The following disinfectants can be used on the system.

Name	Manufacturer
Bleach Germicidal Cleaner	Clorox
Transeptic Cleansing Sol.	Parker

Chapter 8 Error Messages

8.1. Device Message

8.1.1. FAN Error Message

- If there is a fan error, the below pop-up will be displayed.

Caution!
SONON device FAN control is error.
Please reboot SONON device

OK


8.1.2. Battery Charging Message

- If the user connects the power cable when the 300C power is on, the below pop-up will be displayed.

Caution!
Battery charging is activated only
power-off situation. SONON device
will turn off automatically.

OK

8.1.3. Temperature Message

- If the 300C handgrip's temperature is over 40 °C, the below pop-up will be displayed.

Caution!
SONON device temperature is
overheated. Please turn off the
power and re-use after cooling.

OK

- If the 300C handgrip's temperature is over 48 °C, the below pop-up will be displayed.


Caution!
SONON device temperature is
overheated. SONON device turn off
automatically.

OK

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8.1.4. Low Battery Message


- If the SONON 300C battery level is under 30%, the below pop-up will be displayed.


8.2. Connection Error

8.2.1. Firmware Upgrade Error Message


- If the Firmware upgrade has failed or another error has occurred, the below pop-up will be displayed.

8.2.2. Wi-Fi AP missing Error Message

- If the Wi-Fi AP signal is missing, the below pop-up will be displayed.

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OK

- Document No. HC-IFU-001 (Rev. 2.5)

Chapter 9 Technical Data/Information

9.1. Safety Conformance

9.1.1. Conforms to the following safety standards



- IEC 60601-1 Electrical Medical Equipment
- IEC 60601-1-2 Electromagnetic Compatibility
- IEC 60601-1-6 Usability
- IEC 62304 Software Life Cycle Processes
- IEC 60601-2-37 Particular Requirements for the Safety of Ultrasound Medical Diagnostic and Monitoring Equipment
- IEC 62359 Ultrasonics - Field Characterization - Test Methods for the Determination of Thermal and Mechanical Indices Related to Medical Diagnostic Ultrasonic Fields
- ISO 10993 Biological Evaluation of Medical Devices

9.1.2. Essential Performance

- Acquisition of ultrasound images
- Display of ultrasound images on main display
- Measurement of ultrasound images

9.1.3. Marking plate

9.1.3.1. ID Label


Model: SONON300C


- Position: Bottom surface of the device

9.1.3.2. Side Label Marking


⚠ Don't connect micro USB connector to SONON during the use.
only use that for battery charging


- Model: 300C
- Position: the side surface of the device

9.1.3.3. Adapter

Adapter (for 300C)
(for charging device)

AC adapter

9.1.3.4. Battery Charger


9.1.3.5. Package Label

The label on the package contains the following information

- commercial name of the device.
- catalogue number (REF).
- serial number (SN).
- symbol "Manufacturing Date" (year-month)
- symbol "Manufacturer Data". This symbol is followed by the name and address of the device manufacturer.
- symbol "Manufacturer's representative in the EU". This symbol is followed by the name and address of the manufacturer's representative in the EU.
- battery type.
- direct current.
- battery capacity.
- operation mode.
- symbol "Refer to the accompanying documents".
- symbol "Attention! Dangerous working conditions". Please see the "Safety instructions".
- symbol "Electromagnetic Radiation". Interference may occur near the equipment having this sign.
- symbol "Caution: Federal law restricts this device to sale by or on the order of a physician"
- symbol of the conformity of products to the requirements for quality in the territory of the European Union.

Label design for Russia:

	This way up, Indicates correct upright position of the transport package.
	Fragile, handle with care
	Keep dry
	Temperature limitation for storage and transportation -10°C to +50°C (+14°F to +122°F)
	Humidity limitation for storage and transportation 10% to 85% non-condensing
	Atmospheric pressure limitation for storage and transportation 700 hPa (3000 m) to 1060 hPa

9.1.4. Guidance and Manufacturer's Declaration


Warning

- SONON 300C requires special precautions regarding EMC.
- SONON 300C should not be used adjacent to or stacked with other equipment
- Using the wrong cable and accessories may adversely affect the EMC performance

Electromagnetic emissions

The 300C is intended for use in electromagnetic environments, as specified below. The customer or the user of the 300C should ensure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions - CISPR11	Group 1	The 300C uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference to nearby electronic equipment. The 300C is suitable for use in all medical establishments (i.e., hospitals, clinics, etc.). 300C is intended for professional use only.
RF emissions-CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC6100-3-3	Complies	

Electromagnetic Immunity

The 300C is intended for use in electromagnetic environments, as specified below. The customer or the user of the 300C should ensure that it is used in such an environment

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment- Guidance
Electrostatic discharge (ESD) IEC 610004-2	± 6kV contact ±8kV air	±6kV ±8kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 6100044	± 2 kV for power supply lines ± 1 kV for input/output Lines	±2kV ±1kV	The main power quality should be that of a typical commercial or hospital environment
Surge IEC 610004-5	± 1 kV line(s) to line(s) ±2 kV line(s) to earth	±1kV ±2kV	The main power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 610004-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30 % dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Compliance for all test levels. Controlled shutdown with return to pre-disturbance condition after operator's intervention. (Power-on switch)	The main power quality should be that of a typical commercial or hospital environment. If the user of the ME SYSTEM requires continuous operation, and in which the main power is interrupted, it is recommended that the SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 6100048	3A/m	3A/m 50 and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

NOTE: UT is the AC mains voltage prior to application of the test level.

Electromagnetic immunity

The 300C is intended for use in electromagnetic environments, as specified below. The customer or the user of the 300C should ensure that it is used in such an environment

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
Conducted RF EC 61000-4-6 Radiated RF EC 610004-3	3 V rms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3Vrms[V] 3V/m[EI]	Portable and mobile RF communications equipment should be used no closer to any part of the 300C, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 80 MHz–800 MHz $d = 1.2\sqrt{P}$ 800 MHz–2,5 MHz $d = 2.3\sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Held strengths from fixed RF transmitters, as determined by an electromagnetic site surveys, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

SONON Recommended separation distances between portable and mobile RF communications equipment and the SONON

The 300C is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 300C can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 300C as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz – 80 MHz $d = 1.2\sqrt{P}$	80 MHz -800 MHz $d = 1.2\sqrt{P}$	150 kHz – 80 MHz $d = 1.2\sqrt{P}$
0,01	0,12	0,01	0,12
0,1	0,38	0,1	0,38
1	1,2	1	1,2
10	3,8	10	3,8
100	12	100	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance, d, in meters

(m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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9.2. Acoustic Output Reporting Tables

Index label		MI	TIS			TIB	TIC
			Scan	Aaprt ≤ 1 cm ²	Aaprt > 1 cm ²	Non- scan	
Maximum index value		0.7861	0.2535			0.6403	0.3571
Associated acoustic parameters	Pr, α z = 6.3cm(Mpa)	1.4172					
	P (W)	0.01633				0.01633	
	Min of [Pa(zs), I _{ta} , α (zs)]						
	Zs		1.7144				
	Zbp	1.7144					
	Zb (m)						
	z at max lpi α (m)	0.0375					
	deq (zb)	1.143					
	Fawf (MHz)	3.25					
	Dim of Aaprt						
	X	1.14					
	Y	0.9					
Other information	td	1.46					
	Prr(PRF)	0.472ms					
	pr at max lpi(Pa) (Peak rarefactional)	1.417					
	deq at max lpi						
	lpi at max MI	0.16990					
	Focal Length	FLx	0.08				
		FLy	1.04				
Operating control conditions	Control 1	✓		✓		✓	
	Control 2	✓		✓		✓	
	Control 3	✓		✓		✓	
Track3 (For FDA)							
Transducer		ISPTA.3	TIS	MI	IPA.3@M _{imax}		
Convex		0.0627 [W/cm ²]	0.2535	0.7861	0.1699 [W/cm ²]		

9.3. Acoustic output

9.3.1. Definition of the acoustic output parameters

9.3.1.1. Thermal Index

TI is an estimate of the temperature increase of soft tissue or bone. There are three thermal index categories:

- TIS: Soft tissue thermal index, the main TI category. Used for applications that do not image bone.
- TIB: Bone thermal index (for bone located in a focal region). Used for fetal applications.
- TIC: Cranial bone thermal index (for bone located close to the surface). Used for transcranial application.

Reference to calculation of TI can be found in:

- NEMA Standards Publication UD 3: "Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment", Revision 2
- IEC 60601-2-37. Medical Electrical Equipment. Part 2-37: Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment

9.3.1.2. Mechanical Index

MI is the estimated likelihood of tissue damage due to cavitation. The absolute maximum limit of the MI is 1.9, as set by the U.S. Food and Drug Administration (FDA) 510(k) guidance on ultrasound systems issued September 9, 2008 ("Guidance for Industry and FDA Staff, Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers").

9.3.1.3. Ispta

The Ispta is the Spatial Peak Temporal Average Intensity. The absolute maximum limit of Ispta is 720 mW/cm² as set by the FDA 510(k) guidance of September 9, 2008.

9.3.2. Acoustic output and display on the SONON 300C

MI and TI values are displayed on the scanning screen. For all imaging modes of 300C, TIS equals TIB and is displayed as TI.

The 300C chooses the correct category based on mode of operation and chosen application, and presents only one TI to the operator. It is therefore important that the operator chooses the right application.

The maximum possible MI and Ispta on the 300C is within the limits set in Track 3 in the FDA 510(k) guidance of September 9, 2008, MI < 1.9 and Ispta < 720 mW/cm².

9.3.2.1. Display Accuracy and Acoustic Measurement Uncertainty

The display accuracy and measurement precision of the output display are summarized in the table below. Accuracy of the output display parameters depends on the measurement system precision, the acoustic model used to calculate the parameters in the acoustic output of systems. The measurement precision and overall accuracy of the measurements have been assessed by determining both the random and the systematic uncertainties and given in percent at 95% confidence level.

Parameter	Uncertainty
Power	$\pm 5,396\%$

9.3.2.2. Measurement accuracy

The measurement accuracy of the system is shown in the following tables.

Parameter	Accuracy
Length	$\pm 7\%$
Ellipse	$\pm 10\%$

9.3.3. System controls affecting acoustic output


The operator controls that directly affect the acoustic output are discussed in the Acoustic Output Data Tables (see 9.2). These tables show the highest possible acoustic intensity for a given mode, obtainable only when the maximum combination of control settings is selected.

Most settings result in a much lower output. It is important to note the following:

- The duration of an ultrasound examination is as important as the acoustic output, since patient exposure to output is directly related to the exposure time.
- Better image quality yields faster clinical results, making it possible to complete the relevant ultrasound examination in a shorter period of time. Therefore, any control that improves the quality of the examination can help to reduce patient exposure, even though it may not directly affect acoustic output.

9.3.3.1. Clinical application selection

Selecting the application appropriate to a particular ultrasound examination automatically provides acoustic output limits within FDA guidelines for that application. Other parameters which optimize performance for the selected application are also set automatically, and should assist in reducing the patient exposure time.

9.4. Specifications

Parameter	Value
Dimensions (mm)	77.6(W)x218.5(L)x 39.2(H) (± 1 mm)
Weight (g)	390 (with battery) (±10 g)
Material	1. Casing surface – polyphenylsulfone Radel R-5100 NT 15 2. Side plates-Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015 3. Functional button-Polyetherane PU20 4. Scanning surface- organosilicone filling material RTV 627 A+B 5. Power button - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015 6. Wi-Fi button - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015
Battery	Type: rechargeable lithium-ion battery Capacity: 2600 mAh (7.4 V)
Operating time	In scan mode: 2.5-3 h In standby mode: 10-12 h
Wireless communication	Type: Soft AP. Frequency: 2.4 GHz Communication standard IEEE 802.11
Mobile App	OS: iOS/ Android
IEC 60601-1	Type BF working part Short mode of operation (max. 10 min. work and 10 min. rest) Equipment with internal power supply Probe head: IPX7
Application	General (abdominal region) / OB / GY
Ultrasound	Frequency: 3.5 MHz Probe type: convex Scan depth: 5 cm to 20cm
Software version	SONON 300C: SW version M 1.03.01 Mobile device: SONON 300C App version M1.03.01
Battery Charger SONON	Used to charge the lithium-ion battery. Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Dimensions: 115mm*65mm*27mm(± 0.2 mm)
AC adapter for the battery charger SONON	Used to connect to the charger of SONON system. Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Cable length 120 cm (± 1 cm). Dimensions: 90mm*50 mm *30 mm (± 0.2 mm)
AC adapter (for SONON 300C)	Used to connect to the charger of SONON 300C Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Cable length 150 cm (± 1 cm). Dimensions: 90mm*43 mm *30 mm (± 0.2 mm)
Power cord	Used to connect the AC adapter for the battery charger SONON and network adapter for SONON to the power source. 100-240 V, 50/60 Hz. PVC WSP-300. Cable length: 200cm (± 1 cm)
Protective cap for SONON	Is used to protect the probe head of SONON 300C Material: polyoxymethylene POM-C Dimensions: 85 mm *26 mm *26 mm (± 1 mm).
Hand strip for SONON	Is used to carry the SONON 300C (if necessary). Length: 30cm (± 1 cm) Material: polyamide PA66-GF15, polyetherane YU-50D

Fastener for SONON	Ensures fixation of the scanner SONON 300C (if necessary). Dimensions, mm: 550 mm *145 mm *75 mm (± 1 mm). Materials: Scanner fastener - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABSNH-1015 Bracket - aluminum ENAW 7005 Support - steel 1.0143 Handle - polyamide P A6 brand CB-30, stainless steel AISI 201
---------------------------	---

9.5. Reclamations

Complaints about inadequate quality should be addressed to:
 HEALCERION Co., Ltd
 804ho, 38-21, Digital-ro 31-gil, Guro-gu, Seoul, 08376, Korea
 Telephone: +(82) 70-7582-6326
 E-mail: info@healcerion.com

Authorized representative of the manufacturing company HEALCERION Co., Ltd in Russia is OOO Firm "Solo"
 Address: 620137, Russia, Sverdlovsk region, Ekaterinburg, Blyukhera St., house 75, building1, room № 3-17
 Telephone: (343) 351-76-16
 E-mail: mail@solo-trade.ru

Chapter 10 Glossary-Abbreviations

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IEC	Meaning
a	Acoustic Attenuation Coefficient
Aaprt	-12db Output Beam Area
CMI	Normalizing Coefficient
Deq	Equivalent Aperture Diameter
d-6	Pulse Beam Width
deq	Equivalent Beam Diameter
awf	Acoustic Working Frequency
Ipa	Pulse-Average Intensity
Ipa,a	Attenuated Pulse-Average Intensity
Ipi	Pulse-Intensity Integral
Ipi,a	Attenuated Pulse-Intensity Integral
I _{ta} (z)	Temporal-Average Intensity
I _{ta,a} (z)	Attenuated Temporal-Average Intensity
I _{zpta} (z)	Spatial-Peak Temporal-Average Intensity
I _{zpta,a} (z)	Attenuated Spatial-Peak Temporal-Average Intensity
MI	Mechanical Index
P	Output Power
Pa	Attenuated Output Power
P1	Bounded Output Power
pi	Pulse Pressure Squared Integral
pr	Peak-Rarefactional Acoustic Pressure
Pra	Attenuated Peak-Rarefactional Acoustic Pressure
prr	Pulse Repetition Rate
TI	Thermal Index
TIB	Bone Thermal Index
TIC	Cranial-Bone Thermal Index
TIS	Soft-Tissue Thermal Index
td	Pulse Duration / (same)
X, Y	-12 dB Output Beam Dimensions / (same)
z	Distance from the Source to a Specified Point / (same)
zbp	Depth for TIB / Depth at which the relevant index is maximum
zbp	Break-Point Depth / (same)
zs	Depth for TIS / Depth at which the relevant index is maximum

FDA	Meaning
MI	the Mechanical Index.
TISscan	the Soft Tissue Thermal Index in an auto-scanning mode.
TISnon-scan	the Soft Tissue Thermal Index in a non-autoscanning mode.
TIB	the Bone Thermal Index.
TIC	the Cranial Thermal Index.
Aaprt	the area of the active aperture (square centimeters).
pr.3	the derated peak rarefactional pressure associated with the transmit pattern giving rise to the value reported under MI (megapascals).
Wo	the ultrasonic power, except for TISscan, in which case it is the ultrasonic power passing through a one centimeter window (milliwatts).
W.3(z1)	the derated ultrasonic power at axial distance z1 (milliwatts).
ITA.3(z1)	the derated spatial-peak temporal-average intensity at axial distance z1 (milliwatts per square centimeter).
z1	the axial distance corresponding to the location of $\max[\min(W.3(z), ITA.3(z) \times 1 \text{ cm}^2)]$, where $z \geq z_{bp}$ (centimeters).
zbp	$1.69 \sqrt{\frac{f}{f_0}}$ (centimeters).
zsp	the axial distance at which TIB is a global maximum (i.e., $z_{sp} = z_{B.3}$) (centimeters).
z@PII.3max	the axial distance corresponding to the maximum of the derated spatial-peak pulse intensity integral (megapascals).
deq(z)	the equivalent beam diameter as a function of axial distance z. It is equal to $[(4/\pi)(W_o/ITA(z))]^{0.5}$ where $ITA(z)$ is the temporal-average intensity as a function of z (centimeters).
fc	the center frequency (MHz). For MI, fc is the center frequency associated with the transmit pattern giving rise to the global maximum reported value of MI. For TI, for combined modes involving transmit patterns of unequal center frequency, fc is defined as the overall range of center frequencies of the respective transmit patterns.
Dim. of Aaprt	the active aperture dimensions for the azimuthal (x) and elevational (y) planes (centimeters).
PD	the pulse duration (microseconds) associated with the transmit pattern giving rise to the reported value of MI.
PRF	the pulse repetition frequency associated with the transmit pattern giving rise to the reported value of MI (Hz).

FDA	Meaning
pr@PII _{max}	the peak rarefactional pressure at the point where the free-field, spatial-peak pulse intensity integral is a maximum (megapascals). See Section 6.2.4.1 of the Output Display Standard, entitled "Measurement Methodology for Mechanical and Thermal Indices".
deq@PII _{max}	the equivalent beam diameter at the point where the free-field, spatial-peak pulse intensity integral is a maximum (centimeters). See Section 6.2.5.1 of the Output Display Standard, entitled "Measurement Methodology for Mechanical and Thermal Indices".
FL	the focal length, or azimuthal (x) and elevational (y) lengths, if different (centimeters).
IPA.3@MI _{max}	the derated pulse-average intensity at the point of global maximum reported MI (watts per square centimeter).

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
2, Jul. 2018 Under Law No.9416.

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Healcerion Co., Ltd.

Jeongwon Ryu / CEO

SONON300L ULTRASOUND IMAGING SYSTEM with accessories

USER MANUAL, REV. 2

Copyright© All right reserved by Healcerion Co., Ltd.

The SONON 300L with accessories is a diagnostic ultrasound equipment which transmits ultrasound waves into body tissues and forms images from the received echoes.

The SONON 300L with accessories are Class Ma Active Diagnostic Medical Products, according to MDD 93/42/EEC regulations for use on human patients.

The Ultrasound Imaging System SONON 300L with accessories was developed and manufactured by HEALCERION. For more information, please contact HEALCERION.

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
Revision History

Revision	Date
User Manual Revision 0 - First issue	2016-10-25
User Manual Revision 0 - Amendment: Typographical Error	2016-10-25
User Manual Revision 1 - Amendment: 5.1 Mobile Device Requirement - Amendment: 9.2 Acoustic Output Reporting Tables	2017-01-12
User Manual Revision 2 - Amendment: 2.2 Contraindication - Amendment: 3.1 Symbols - Amendment: 5 Installation of the mobile app - Amendment: 6 Operating the device	2017-02-22
User Manual Revision 2, version 2 - Revised: 2.1 Intended Use - Revised: 4.1 Device description - Revised: 5 Mobile App Installation - Revised: 5.1 Mobile Device Requirement - Revised: 5.2 Mobile App Installation - Revised: 6.5 Patient Information	2018-12-14

Chapter 1 General

1.1. About This User Manual

Product name	SONON 300 L Ultrasound Imaging System with accessories
Model	300L

- Please read and understand all instructions in this User Manual before attempting to operate the SONON 300L equipment.
- Keep this User Manual with the product for future reference.
- Some options or features may not be available in some countries.
- The screen graphics and illustrations in this User Manual are for illustrative purposes only, and may differ from what is displayed on the screen.
- All references to standards and regulations, and their revisions are valid as of the date of publication of this User Manual.

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2.1. Intended Use

Caution

The user should be careful to manage patient information and key security information because SONON 300L is used with personal mobile device.

Intended use: The SONON 300L is indicated for ultrasound echo imaging, measurement, and analysis of the human body for general clinical applications including musculoskeletal (MSK), vascular, small parts (breast, thyroid) and lung imaging.

Scope of application: The SONON 300L is indicated for primary and special medical help, e.g. in medical offices, clinics and hospital departments of emergency surgery.

This device is designed for passive observation and early detection of diseases. If you need a more complete examination or doctor's opinion, we recommend that you send the patient to the hospital for a more thorough examination

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2.1.2. Intended User Profile

Item	Description
Education	Professionals who completed related training.
Knowledge	having basic knowledge on ultrasound
	having understanding on physiological effects of ultrasound.
	able to intuitively recognize and understand SWUI and HWUI
	able to distinguish mobile devices to which SW is applied
	able to understand updates.
Language understanding	able to understand terms used for the product and its manual
	able to understand the methodology in the manual (English version)
Experience	familiar to the ways to use mobile devices.
	able to set communication between the device and mobile devices.
	able to keep track of communication status.
	familiar to use apps of mobile devices.
	able to set and install SW updates.
	able to understand and use contents of SW UI.
	well trained about existing ultrasound.
Vision	fully understands the manual accompanying the product.
	corrected visual acuity 1.0 or better
Memory	read through the user manual and keep in mind precautions and key functions of the product
Permissible impairments	N/A

2.1.3. Operation principle

The 300L utilizes pulsed-echo technology to determine the depth and location of tissue interfaces and measure the duration of an acoustic pulse from the transmitter to the tissue interface and back to the receiver. Ultrasound waves are emitted from the transducer, propagate through various tissues, and return to the transducer as reflected echoes.

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The SONON 300L Ultrasound Imaging System is single mode (B Mode), single output (3,5 MHz – low ultrasonic powered) device. The SONON 300L Ultrasound Imaging System a device with operating B mode and CF mode (5 MHz / 7,5 MHz / 10 MHz – low ultrasonic powered).

B mode is a two dimensional mode (2D by sectors) which allows at a certain depth to get a planar image of tissue interfaces and their movement within a timespan. It is the easiest mode to comprehend, as it helps to display anatomical formations like on a cross-sectional cut. This is the main mode to use and the first one start with during any ultrasonography. The image creation in B mode is based on reflection of ultrasound waves by the tissue interfaces of the examined area: the Ultrasound Imaging Device calculate the amplitude and the phase of the returned signal, which information is used to build an image (that will be displayed in the screen) and determines the pixel brightness and its position on the diagnosis image.

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Warning

The 300L is not intended for:

- ophthalmic use or any use causing the acoustic beam to pass through the eye;
- use during surgery (e.g., when a sensor is inserted into a surgical incision or a trephine hole);
- intracavitary application;
- visualization of an open wound.

Possible side effects

There were no side effects.

2.3. Terms of use

The SONON 300L scanner is suitable for use in all medical institutions (i.e. hospitals, clinics, etc.). The scanner SONON 300L is intended for use by specialists only.

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Symbol	Description/Function
	Electrical protection: Insulated patient application (Type BF)
	Warning: This symbol indicates hazard about Operation condition For more information, see on Section 3.2 "General Safety"
	Consult accompanying documents: This symbol advises the reader to consult the accompanying documents.
	Freeze mode: This symbol means that image is frozen.
	Stand-by: This symbol's function is to turn the device on or put it in standby mode.
	Wi-Fi: This symbol means wireless communication.
	Manufacturer information: This symbol is followed by the name and address of the device manufacturer.
	Manufacture Date: This symbol is followed by the device manufacture date in the form YYYY-MM
	Serial Number: This symbol is followed by the device serial number.
	Model name: This symbol represents the model name of product
	Representatives information: The manufacturer's EU representative information shows with this symbol
	The CE marking is a marking of the conformity of products to the quality requirements in the territory of the European Union. The CE mark ensures that the products manufactured by the manufacturer fully comply with the safety requirements for humans and the environment
	This way up, Indicates correct upright position of the transport package.
	Fragile, handle with care
	Keep dry
	Indicates the need for separate collection for electrical and electronic equipment in compliance with the Waste Electrical and Electronic Equipment (WEEE) Directive.
	This symbol means "Caution: Federal law restricts this device to sale by or on the order of a physician"
	Electromagnetic radiation. There may be interference near equipment marked with this icon
	Temperature limit

Humidity limitation

Atmospheric pressure limitation

3.2. General Safety

Warning

- Only authorized personnel shall perform any type of repair on the 300L. Never attempt to open a transducer or transducer connector. This will void the warranty!
- Probes are not delivered sterile! Before the first use, it is MANDATORY to clean and disinfect probes to avoid infections or disease transmissions!
- Probes must be cleaned and disinfected before they are replaced or disposed of.
- Do not modify this device without authorization of HEALCERION.
- Do not use the probe with high frequency surgical equipment. Doing so may damage the equipment.
- If you lose your tablet PC or smart phone, the stored data cannot be recovered. Please take care not to lose these items.
- Do not touch the patient with iOS or Android mobile device when you use the 300L
- Must comply with Operation condition; **max. 10 min. with 10 min. resting time.**

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
3.3.1. Immersion Level

Caution

- Probes labeled "IPX7" are watertight up to a maximum of 4 cm below the probe.

max. 4 cm

3.4. Electromagnetic Compatibility

- The system must remain the essential performance of the ultrasound system in EMC environment according to IEC 60601-1-2. The essential performance refers to section 9.1.4 "Guidance and Manufacturer's Declaration" of this manual.

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- Methanol, ethanol, isopropanol alcohol, or any other alcohol-based products

- Mineral oils
- Iodine
- Lotions
- Lanolin
- Aloe vera
- Olive oil
- Methyl or ethyl parabens (para hydroxybenzoic acid)
- Dimethylsilicone

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"A fundamental approach to the safe use of diagnostic ultrasound is to use the lowest output power and the shortest scan time consistent with acquiring the required diagnostic information. This is the ALARA principle (i.e. As Low As Reasonably Achievable). It is acknowledged that in some situations it is reasonable to use higher output or longer examination times than in others: for example, the risks of missing a fetal anomaly must be weighed against the risks of harm from potential bio-effects. Consequently, it is essential for operators of ultrasound scanners to be properly trained and fully informed when making decisions of this nature." (Guidelines for the Safe Use of Diagnostic Ultrasound Equipment, Safety Group of the British Medical Ultrasound Society 2010)

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The SONON Ultrasound Imaging System (Model: 300L) is B Mode and CF Mode, output (5MHz / 7.5MHz/10MHz- low ultrasonic powered) device.

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Operated B-mode and CF-mode, The SONON Ultrasound Imaging System is safe because exposure time to patient applied part is not concentrated in one place. In Addition, Healcerion provide the generic content of the education program with an AIUM "Medical Ultrasound Safety" (AIUM 2014 - Refer to appendix.) and additional education. For example, the doctor who used the SONON 300L Ultrasound Imaging System to verify the status of the baby prohibit inordinate inspection in order to show the family or other meaningless situation.

3.6.2. Bioeffects

- Thermal effects refer to the heating of soft tissue and bone

The thermal indices "TIs" (for soft tissue), "Tib" (for bone near the focus) and "Tic" (for bone near the surface) were introduced to provide the operator the relative potentials for a tissue temperature rise. According to the Standard for Real-time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment (2004) those thermal indices shall be displayed by ultrasound consoles. It should be noted that a TI of 1 does not necessarily mean that tissues being scanned will increase in temperature by 1°C. Almost every scanning situation departs from the assumed model conditions, such as tissue types, blood perfusions, modes of operation and actual exposure times of the scanned area. However, the thermal indices provide information regarding the possible increase in the risks of potential thermal bioeffects and it provides a relative magnitude that can be used to implement ALARA. In addition to tissue heating due to the generated ultrasound field, the temperature of the probe head itself can also increase during the examination. The operator should be aware, that in the tissue region near the transducer, there will be a superposition with the heating due to the ultrasound field, which is not considered by the TI values.

- Non-thermal effects refer to mechanical phenomena such as cavitation

Non-thermal bio-effects are caused by the interaction of ultrasound fields with very small pockets of gas (stabilized gas bodies), i.e., the generation, growth, vibration and possible collapse of microbubbles within the tissue. This behavior is referred to as cavitation (Medical Ultrasound Safety, 2nd Edition, AIUM 2009/American Institute of Ultrasound in Medicine Consensus Report on Potential Bio-effects of Diagnostic Ultrasound, AIUM 2008/Guidelines for the safe use of diagnostic ultrasound equipment, Safety Group of the British Medical Ultrasound Society 2010). The potential of cavitation increases with the rarefactional peak pressure, but decreases with the pulse frequency. Therefore the Mechanical Index (MI) was introduced to take account of both pressure and frequency. The higher the MI, the greater the risk of non-thermal bio-effects.

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Item	Operational	Storage&Transport
Temperature	+18°C to +35°C (+64.4°F to +86°F)	-10°C to +50°C (+14°F to +122°F)
Humidity	30% to 75% noncondensing	10% to 85% noncondensing
Pressure	700 hPa (3000m) to 1060 hPa	700 hPa (3000m) to 1060 hPa

3.8. Other Recommendations

Like most high frequency computing devices, the electronic components of the SONON 300L will generate some heat while operating normally and as intended. The SONON 300L is equipped with safety mechanisms which will automatically reduce computing speed (frame rate), and ultimately shut down the device, before any risk of overheating occurs. The SONON 300L is verified to comply with harmonized safety standards, under any operating condition described in this user manual (see Section 3.7 "Environmental Conditions"). To help keep the SONON 300L operating temperature at the optimal level, and to ensure longer scanning times with the maximum frame rate, it is recommended to hold the SONON 300L so that there is good contact between the device and the hand. The image display of the mobile device is dependent on the ambient light; avoid direct sunlight on the display when scanning and reviewing images.

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3.10. Transportation

The SONON 300L ultrasound imaging system in the package is transported by all modes of transport, in covered vehicles in accordance with the rules of transportation of goods, operating on each mode of transport.

Transportation and storage of the SONON 300L ultrasound imaging system without the packaging of the manufacturer does not guarantee its safety. Damage to SONON 300L as a result of transportation or storage without packaging of the manufacturer is eliminated at the expense of the consumer.

3.11. Warranty

The manufacturer warrants that there is no rejection of materials and the quality of the SONON ultrasound scanning system when properly used and maintained for 12 months (excluding the battery) from the date of purchase


3.12. Packaging

The SONON 300L ultrasound imaging system is packed in a cardboard box.

Box dimensions: 323 × 255 × 71,5 mm (W × L × H).

Cardboard label: GZ (Coated SBB).

The box contains: scanner SONON 300L, battery charger SONON, two accumulator batteries (one inserted in the scanner, the other is spare), network adapter (for SONON) and power cord, placed into a separate cardboard box, and AC adapter for the SONON battery charger with the power card, also placed into a separate cardboard box, and user manual.

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

4.1 Product Description

4.1.1. Package Contents:

Scanner SONON 300L	User Manual SONON 300L
	
Accumulator battery - 2 pcs.	Battery Charger SONON
	
Adapter (for Battery Charger) SONON	Power Cord - 2 pcs.
	
Network adapter (for SONON 300L)	Brief manual SONON 300L
	
Accessories:	Software Usage Agreement
Protective cap for SONON	Hand strap for SONON
	
Mounting for SONON	
	


4.2 Screen Layout

4.2.1. Outline of apparatus


No.	Name	Function
1	Probe Head	To be applied to the patient
2	Function button	Blue LED is turned off during the freeze screen. Blue LED is turned off during the live screen.
3	Power LED	Turning on the power supply (White) Turning off the power supply (Green LED is displayed for 2 seconds, and then, the LED is turned off.) Charging in progress (Flickering with green color)
4	Power button	Power On /Off
5	Wi-Fi button	Wi-Fi Indicator (Blue)
6	Charging device connected	This is the section where AC adapter is connected for charging the battery of apparatus.

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
Figure 1 Patient environment

The SONON 300L is powered by a Lithium Ion battery. The battery is not fully charged prior to shipment. To maximize battery life, it is recommended that the battery be fully charged before initial use. Use only the AC Adapter provided with the SONON 300L.

28

4.3.1. Battery Level Indicator

Notice

- The battery should not be charged in the Power On state. To charge the battery, power should be turned off.
- We recommend to use battery charger and adaptor because charging by adaptor for SONON 300L do not fully charge to protect the SONON 300L device.
(About 80% charging by adaptor for SONON 300L)
 - Condition Operation: 10 minutes
 - Time to Fully Charge: 3 Hours

The battery level indicator is displayed by color and flash rate. The following icons are displayed.

※ 300L

Battery Level	Description
The battery doesn't charge, or Charging error	Battery level indicator flashes
Charging is going on	Battery level indicator lights green (doesn't flash)
Fully charged	Battery level indicator is off
Battery low	If there is a connection to the application, a pop-up window will appear on the screen indicating that the battery is low. If the scanner is not connected to the application, the function button is disabled (white indicator)

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


4.3.3. Battery Specifications

Item	Specification
Description	Rechargeable Li-ion Battery Pack
Capacity	2600 mAh (7.4 V)
Lifetime Guarantee	6 months
Model	MBP-2S1PSD26
Cell Type	ICR18650
Dimensions	L*W*H 70 mm (± 0.2 mm) * 37 mm (± 0.2 mm) * 20 mm (± 0.2 mm)
Safety	PCM Logic

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Chapter 6

To operate SONON 300L ultrasound imaging system you need to install SONON Ultrasound App. It is impossible to log into the mobile application without SONON 300L devices (just as the scanner cannot work without the mobile application installed)! The access code / access key (or wifi connection between the mobile device with pre-installed software and the SONON 300L scanner) is a serial number of the scanner indicated on the identification label located on the back of the product. To activate the program, enter the login - the serial number (SN) of the scanner SONON 300L, company name, email and set your own password. Thus, you create your account which you will use as the only way to access the software.

Software Features

- displaying on your phone or tablet the ultrasound images made by SONON 300L scanner;
- processing (converting) the images to be displayed
- displaying the processed (converted) image on the phone/tablet
- managing diagnosis features
- sensor-based control of the main program settings
- saving the diagnosis information (image, measures, user notes/comments and patient information) and managing previous records (edit, copy, delete, send by email)
- measuring of the received ultrasonic data.
- managing scanner settings: wifi connection settings with possibility to modify «default password», upgrading the scanner firmware
- managing transfer and backup of patient information

5.1. Mobile Device Requirement

Warning

- If you remove and reinstall the app, all stored information will be deleted and cannot be recovered.
- Intended for use with standard Android (v. 4.4 or later) and iOS operating systems (v.9.0 or later). Android Rooting and iOS Jailbreaking removes the limitations on the manufactured devices running the standard OS, and thus the SONON 300L is not intended for use with these devices. Healcerion is not responsible for any error in usage.
- Do not use anything other than a mobile device that recommends.

Caution

After installing the SONON 300L App from the App Store / Google Play, the new version is recognized and it is automatically updated. At this time, firmware also is automatically upgraded.

Notice

When you operate the SONON 300L, HELCERION needs your personal information. Please read 'Appendix A. Usage Agreement'.

Personal information protection and protection of the device from the unauthorized access is performed as follows:

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
Item	Minimum requirements
Target Device	iPad Air 2/ iPad Mini 4 iPhone 6, 6 plus, 6S, 6S plus Galaxy S5 / S6 / Note 3, Note 4 or later GalaxyNoteTablet версии 10.12013 or later GalaxyTabPro версии 8.4 2014 or later
Mobile OS Version	iOS 9.0 or later Android 4.4 or later
Minimum device storage available for mobile app installation	8 Gb

5.2. Installation of the Mobile App

Notice

For iPhone, installation of the SONON 300L mobile app requires an iTunes registration and password.

5.2.1. iOS


- (1) Launch the "App Store"  App
- (2) Search for "300L"
- (3) Select the SONON 300L App and Install
- (4) Enter the user's Apple iTunes password
- (5) Automatic installation

5.2.2. Android

- (1) Launch the "Play Store"  App
- (2) Search for "300L"
- (3) Select the SONON 300L App and Install
- (4) Automatic installation

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The SONON 300L system is designed to work on standard Android OS (v. 4.4 or higher) and iOS (v. 9.0 or higher).

Chapter 6 Operating the Device

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1010101

- Please check for any physical damage of the exterior, including cracks and foreign materials before using the product.
- **The system** should be cleaned/disinfected per instructions in section 7.3 "Cleaning & Disinfecting" prior to use.

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


6.2.1. Wi-Fi Connection

- (1) Turn on the power button and Wi-Fi LED light should turn to green.
- (2) Go to settings menu at your tablet or mobile phone and search for SONON SSID. SSID: SONON300L - Serial number.
- (3) Input initial password and complete the connection.
Password: 1234567890

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※ The basis of line 1 Marker is based on the protruding part of SONON 300L product's probe head section.

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
6.4.2. Registration in the App

6.4.2.1. Primary adjustment screen

※ Must fill in the registration information only once for the first time.

(1) Select the next button after opening the 300L app.

- (2) Agree with the terms and conditions.
- (3) Fill out your information.
- (4) Select the OK button.


* Be careful not to forget your password.

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
6.4.3.1. Find Password

- (1) Fill out the serial number and e-mail address you filled in the first time.
- (2) If the user fills out the correct information, the password will pop-up.

6.4.4 Start Scan


※ When SONON 300L is connected to the mobile application and if the user taps the start-button, it is available to start the scan.


- (1) Select the start scan button.
- (2) The progress icon is displayed.
- (3) After connecting the SONON 300L and mobile application, the information of SONON 300L is normally outputted.
- (4) When SONON 300L is not connected to the mobile application, the user can't use the scan functions.

6.4.5. Live Screen

6.4.5.1. B mode Default Screen


6.4.5.1.1. Indicator


No.	Name	Function
1	Device type Indicator	Available Device type 300C/310C/300L
2	Mobile device battery capacity Indicator	Mobile device (Phone or Rad) battery capacity
3	End Examination	Return to initial screen (Patient management/ Start Scan selection)
4	Scan status Indicator	Live/Freeze Status
5	Additional Action	Patient Management, Preset, Quick viewer, Setting

6.4.5.1.2 Menu Area

No.	Name	Function
1	B Mode Selection	Default B mode, All device type support
2	B mode image Parameter	Gain / DR / TGC, All device type support
3	Gain	Overall Gain control
4	DR	Dynamic Range control
5	TGC	Time Gain Compensation Control
6	Filter image parameter	Frame Average/SRI /Gray Map, All device type support
7	Frame Average	4 Step (0,1, 2,3)
8	SRI	6 Step (0,1, 2,3,4,5)
9	Gray Map	8 type gray map support
10	ETC image parameter	Multi Focus / Frequency, Convex model type is not supported
11	Multi Focus	On / Off, Max 2 Focus selection
12	Frequency	5,7.5,10 MHz Center Frequency selection
13	CF Mode Selection	B+CF mode, ROI box control, 300C is not supported
14	CF image parameter	In CF mode selection, automatically changed to CF image parameter
15	CF control panel	Invert/Steer/ Box size/ Flow Speed /C Gain / Rejection Control

6.4.5.1.3. Menu Area

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
No.	Name	Function
1	Patient Information	Patient Name/ID
2	Operation Indicator	Network Status / SONON device battery capacity checking
3	Line 1 Marker button	Display Scan Line 1 (Reversal screen is displayed.)
4	Acoustic Output Power Level Indicator	MI, TI value is displayed with Operator & Hospital information
5	Ultrasound Scan Information	Preset / FPS / Frequency information
6	Focus control	Focus value
7	Depth control	Depth Value

6.4.5.2. CF mode Default Screen


No.	Name	Function
1	Invert	Color Invert (On/Off)
2	Angle steer	ROI box steer (-15/ 0/15)
3	Color box size	ROI box size (Normal / Large)
4	Flow speed	Color scale (Slow/ Med / Fast)
5	C gain	Color gain control
6	Rejection	Low signal rejection

1013101


575757

6.4.5.3. Zoom Action


- ※ Pinch Out (Zoom in) Action (max 400%)
- ※ Pinch in (Zoom out) Action (min 100%)


6.4.5.4. FL Control

- (1) Single tap FL position in side ruler area
 - (2) Display New FL marker
 - (3) The MI /TI value in accordance with the change of FL value is displayed.
- ※ FL value cannot exceed the value of depth, and maximum FL value shall be set up as Depth-1.
 - ※ The range of FL change is 0.5 cm~6 cm.


6.4.5.5. Depth Control


- (1) Drag the side ruler area.
- (2) Display New Depth.
 - ※ Drag up - increase Depth / Drag down - decrease Depth
 - ※ If the depth is reduced by user to less than the FL value, the FL is automatically changed to the finally changed value of depth -1.
 - ※ e.g When the depth before the change is 10 cm, and FL is 6 cm, if the depth is changed to 5 cm, the FL is automatically changed to 4 cm.
 - ※ The range of depth change is 3 cm~10 cm.

6.4.5.6. Recording

6.4.5.6.1. Cine File


- (1) Select the recording button.
- (2) During the recording a message «Recoding» appears on screen.


- (3) Press the recording button
- (4) The recording end message, 'Cine file saved' is displayed, and the diagnosis data (Recording file) is saved


6.4.5.6.2. Still-Shot Image

- (1) Select the recording button.
- (2) 'Still shot image is saved' is displayed, and the diagnosis data is saved.

6.4.5.7. Quick viewer

- (1) Select the **More** button in the indicator area, extension menu will show.
- (2) Select Quick viewer in the indicator area.
- (3) Playlist will be displayed.
- (4) Select a video or picture to display a saved diagnostic file.

6.4.5.8. Select and Save Preset

- (1) Select the **More** button in indicator area, extension menu will show.
 - (2) Select preset menu in preset list, then change the setting value.
 - (3) To change the preset values, change the settings in the live screen.
 - (4) The user can overwrite or save as a new name. If it is overwritten, the icon, **More**, will be displayed at the end of the name.
- ※ Go to settings to reset or delete a preset.


6.4.6. Freeze Screen

✖ While scanning, press Function button or tap freeze button to freeze the image.

No.	Name	Function
1	Scan – status Indicator	Display Live/Freeze.
2	Image Parameter	Confirm settings value of B-mode, CF-mode, Filter and ETC.
3	Function button	Select Live / Freeze.
4	Rewind Bar	Drag the rewind bar and confirm the saved diagnosis image. (Saved up to 150 frames)


6.4.6.1. Length

- (1) Single tap the Measurement button.
- (2) The Sub Menu will be displayed.

- (3) Tap the Length button.
 - (4) Draw line where to measure length and yellow line is displayed.
 - (5) Display Length Value on the screen.
- ※ User can move the position icons with the drag action.


6.4.6.2. Ellipse

- (1) Single tap the Measurement button.
- (2) The sub menu icon will be displayed.
- (3) Tap the Ellipse button.
- (4) Drag line where to measure ellipse and displayed Ellipse frame (Yellow).
- (5) Double tap the middle of yellow point of ellipse frame, it changes to Green and displayed position icon.
- (6) Displayed Ellipse value on ellipse frame.
- (7) Displayed Ellipse value on the screen.


6.4.6.3. Comment


- (1) Single tap the Measurement button.
- (2) Display the sub menu.
- (3) Tap Comment button.
- (4) The virtual keyboards will be displayed.
- (5) Input comment and tap the send of virtual keyboard.
- (6) Display user's input note.

6.4.6.4. Delete

- (1) Delete individual.
 - Select the measured data and drag to the trashcan icon at the right bottom of the screen.
- (2) Delete All.
 - Select trashcan icon and delete all measured data.

6.4.6.5. Saving Diagnosis File

- (1) Single tap the saving icon button.
- (2) The message 'Still-shot image is saved' is displayed for 2 seconds, and the image is saved.

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

- (1) Tap the patient info icon- Same as the Freeze Mode Patient Info, button.
- (2) Move to Patient list Screen.

6.5.1.1. Add new patient

- (1) Tap Add new patient button.
- (2) Input each item at the Patient Info screen.
- (3) Select the Save.
- (4) Popup window is displayed, and select OK.

6.5.1.2. Patient List Edit

- (1) Tap edit button
- (2) Tap selected patient on patient list.
- (3) Tap the activated Delete.
- (4) Print popup for caution about deleting patient.


6.5.1.3. Patient Information Edit

- (1) Tap edit button.
- (2) Input each item at the Patient Info screen.
- (3) Select the Save.
- (4) Popup window is displayed, and select OK.


6.5.2. Play List


- (1) From the patient list screen, select the patient whose diagnosis data is to be reviewed..
- (2) The patient information screen of the patient chosen is displayed, and select the Playlist.
- (3) From the Playlist screen, the diagnosis data of the patient selected is displayed.

6.5.2.1. Edit Play list


6.5.2.1.1. Delete

- (1) Select Edit button.
- (2) Select a file you want to delete, and select delete button.
- (3) To delete the file, select OK button if a pop-up window appears


6.5.2.1.2. Move


- (1) Select a file you want to move, and select move button.

- (2) Select the patient name.
- (3) Confirm the process with a completion pop-up.

6.5.2.2. Diagnosis File Viewer


- (1) Tap the Diagnosis Cine File.
- (2) The Cine File Viewer will be displayed.

- (1) Tap the Diagnosis Image.
- (2) The Image Viewer will be displayed.

After examination, the image (photo or video) saved on smart phone or tablet can be sent to another telephone (or PC) in order to archive it (storage will allow its further analysis to establish diagnosis or monitor patient's condition). This information will be analyzed by ultrasonography professional or consulting physician by visual estimation of the image and recorded measures. The report will be provided by the ultrasonography professional based on the image analysis. He/she will provide a specific image description and conclusion (following the ultrasound examination procedure) that will be added to the medical history sheet.

6.5.2.3. Share

- (1) Select **Share** button in the viewer
- (2) Select **Share** button.
- (3) Select **Create new**, and the transformation process will begin.
- (4) Select the file sending method.


6.6. Settings

6.6.1. Settings Menu


- ✖ Setting menu can be selected from the Live, Freeze, and start screen.

- (1) From the Start screen, select the  setting icon.
- (2) The Settings screen is appeared, and the 5 items, i.e., Preferences, Export, Preset, Backup, and Import are displayed.


- (1) From the Live, Freeze screen, select  icon
- (2) As the sub-menu, and when the Settings button is selected, the 6 items, i.e., Device settings, Preferences, Export, Preset, Backup, and Import are displayed.

6.6.2. Device Settings

6.6.2.1. Firmware Update


※ If firmware update is needed, firmware update button will be displayed automatically.

- (1) From the Settings screen, select the Device Settings button.
- (2) From the Device settings screen, select the Change for changing the Wi-Fi password of SONON 300L

6.6.2.2. Change Wi-Fi password


※ Disconnect the network first then, reconnect to the device through Wi-Fi again.

- (1) From the Settings screen, select the Device Settings button.
- (2) From the Device settings screen, select change to reset the Wi-Fi password of SONON 300L.
- (3) Select the SSID / Key field. Keyboard is appeared, and it is activated so that the SSID and Password can be changed and set up.
- (4) After changing the password, select 'Save'.
- (5) On the Device Settings screen, the changed details are displayed.

※ Forgot Wi-Fi password.

- (1) Press the function button and Wi-Fi button together for about 5 seconds
- (2) SONON device will be turned off and SSID / Password is initialized when the SONON device is turned on.
- (3) Password is initialized
 - SSID: SONON300L-H-xxxxxxxx (Serial number of SONON)
 - Password: 1234567890


6.6.3. Preference

- ※ Thermal Index - Setting up
- ※ type of TI value to be displayed on the screen (For MI value, the basic value is displayed).
- ※ Recording File Size - Recording file save unit (Capacity).


- (1) Input the desired size and select the Preferences button from the Settings screen.
- (2) Select the capacity of the recording file to be saved to the Recording File Size.
- (3) On the Device Settings screen, the changed details are displayed.

6.6.4. Export


- (1) Before you share a file, Set types of file when a file is shared in Export menu.
- (2) Set compression ratio of file in Export menu.

6.6.5. Preset

- (1) Select Reset to default the preset value.
- (2) Select delete to delete the saved preset (Saved as New Name)


6.6.6.Backup

- * When selecting the backup after choosing the patient to be backed up, the data of corresponding patient is backed up.
- * The file backed up is saved to the 'BackupFiles' folder of SONONApp.

- (1) Select the Backup.
- (2) Select the patient to be backed up.
- (3) Select the Backup button.
- (4) Pop up window is displayed, and select OK.
- (5) After completing the backup, the data memory used is increased.


6.6.7.Import Backup File

※ When selecting the patient with backed up file, the data of the selected patient can be imported.

- (1) From the Settings screen, select the Import backup file.
- (2) On the Import screen, total capacity and the backed up file are displayed.
- (3) Select the backed up file to be imported again.
- (4) Popup window is displayed, and select OK.

6.6.7.1. Import Backup File-Delete

- (1) From the Settings screen, select the Import backup file.
- (2) After selecting the patient to be deleted, select the delete.
- (3) Popup window is displayed, and select OK.
- (4) Among the Import patient list, the backup file selected is not displayed.

Chapter 7 Maintenance

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1012101

Caution

- The user must ensure that safety inspections are performed at least every 12 months according to the requirements of the patient safety standard IEC 60601-1 / UL60601-1.
- Only trained persons are allowed to perform the safety inspections mentioned above.
- The 300L requires regular care and maintenance to function safely and properly.
- To ensure that the 300L operates continuously at maximum efficiency, we recommend that the following procedures be observed as part of the customer's internal routine maintenance program.

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OK

8.1.2. Battery Charging Message

- If the user connects the power cable when the SONON 300L power is on, the below pop-up will be displayed.

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OK

8.1.3. Temperature Message

- If the 300L handgrip's temperature is over 40 °C, the below pop-up will be displayed.

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OK


- If the SONON 300L handgrip's temperature is over 48 °C, the below pop-up will be displayed.

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OK

8.1.4. Low Battery Message

- If the SONON 300L battery level is under 30%, the below pop-up will be displayed.

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Caution!
Firmware upgrade is failed. Please
retry to upgrade firmware.

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Connection problem has occurred or
already connected other client
Please check your SONON device
network status.

OK

8.2.3. Network Error Message

- When the HW device connection is not ready and if the user taps the start button, the below pop-up will be displayed.

Connection with SONON is Not working.
Please retry to connect with SONON

OK

- Confirm whether the network (Wi-Fi) of the terminal in use is connected with SONON 300L.

If error message is generated continuously when the network is normally connected, please contact with the Healcerion Co., Ltd.

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

9.1.3. Marking plate

9.1.3.1. ID Label

- Model: 300L
- Position: bottom surface of the device

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
⚠ CE 0123 Made in Korea

- Model: 300L
- Position: the side surface of the device

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
Label design for Russia:

Transport packaging contains the following symbols:

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Warning

- SONON 300L requires special precautions regarding EMC.
- SONON 300L should not be used adjacent to or stacked with other equipment
- Using the wrong cable and accessories **may adversely affect the EMC performance**

Electromagnetic emissions

The SONON 300L is intended for use in electromagnetic environments, as specified below. The customer or the user of the SONON 300L should ensure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions - CISPR11	Group 1	The SONON 300L uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference to nearby electronic equipment.
RF emissions-CISPR11	Class B	The SONON 300L is suitable for use in all medical establishments (i.e., hospitals, clinics, etc.). SONON 300L is intended for professional use only.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC6100-3-3	Complies	

Electromagnetic immunity

The SONON 300L is intended for use in electromagnetic environments, as specified below. The customer or the user of the SONON 300L should ensure that it is used in such an environment

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment- Guidance
Electrostatic discharge (ESD) EC 610004-2	± 6 kV contact ± 8 kV air	± 6 kV ± 8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 6100044	± 2 kV for power supply lines ± 1 kV for input/output Lines	± 2 kV ± 1 kV	The main power quality should be that of a typical commercial or hospital environment
Surge IEC 610004-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV ± 2 kV	The main power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 610004-11	$< 5\% UT$ ($> 95\%$ dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles $< 5\% UT$ ($> 95\%$ dip in UT) for 5 s	Compliance for all test levels. Controlled shutdown with return to pre-disturbance condition after operator's intervention. (Power-on switch)	The main power quality should be that of a typical commercial or hospital environment. If the user of the ME SYSTEM requires continuous operation, and in which the main power is interrupted, it is recommended that the SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 6100048	3 A/m	3 A/m 50 and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

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Electromagnetic immunity

The SONON 300L is intended for use in electromagnetic environments, as specified below. The customer or the user of the SONON 300L should ensure that it is used in such an environment

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment-Guidance
Conducted RF EC 61000-4-6 Radiated RF EC 61000-4-3	3 V rms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3Vrms[V] 3V/m[EI]	Portable and mobile RF communications equipment should be used no closer to any part of the SONON 300L, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 80 MHz–800 MHz $d = 1.2\sqrt{P}$ 800 MHz–2,5 MHz $d = 2.3\sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site surveys, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

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The SONON 300L is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SONON 300L can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SONON 300L as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz – 80 MHz $d = 1.2\sqrt{P}$	80 MHz -800 MHz $d = 1.2\sqrt{P}$	150 kHz – 80 MHz $d = 1.2\sqrt{P}$
0,01	0,12	0,01	0,12
0,1	0,38	0,1	0,38
1	1,2	1	1,2
10	3,8	10	3,8
100	12	100	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance, d, in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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
9.3.1.2. Mechanical Index

MI is the estimated likelihood of tissue damage due to cavitation. The absolute maximum limit of the MI is 1.9, as set by the U.S. Food and Drug Administration (FDA) 510(k) guidance on ultrasound systems issued September 9, 2008 ("Guidance for Industry and FDA Staff, Information for Manufacturers Submitting Marketing Clearance of Diagnostic Ultrasound Systems and Transducers").

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9.3.2. Acoustic output and display on the SONON 300L

MI and TI values are displayed on the scanning screen. For all imaging modes of SONON 300L, TIS equals TIB and is displayed as TI.

The SONON 300L chooses the correct category based on mode of operation and chosen application, and presents only one TI to the operator. It is therefore important that the operator chooses the right application. The maximum possible MI and Ispta on the SONON 300L is within the limits set in Track3 in the FDA510(k) guidance of September 9, 2008, MI < 1.9 and Ispta < 720 mW/cm².

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1105101

9.4. Specifications

Parameter	Value
Dimensions (mm)	78(W)x229(L)x 38(H)(± 1 mm)
Weight (g)	370 (with battery) (±10 g)
Material	1. Casing surface – polyphenylsulfone Radel R-5100 NT 15 2. Side plates-Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015 3. Functional button-Polyetherane PU20 4. Scanning surface- organosilicone filling material RTV 627 A+B 5. Power button - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015 6. Wi-Fi button - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015
Battery	Type: rechargeable lithium-ion battery Capacity: 2600 mAh (7.4 V)
Operating time	In scan mode: 2.5-3 h In standby mode: 10-12 h
Wireless communication	Type: Soft AP. Frequency: 2.4 GHz and 5 GHz, double-band Wi-Fi Communication standard IEEE 802.11
Mobile App	OS: iOS/ Android
Ultrasound	Frequency: 5 MHz / 7.5 MHz / 10 MHz Probe type: linear Scan depth: max. 10cm
Application	musculoskeletal system (MSK), vascular, small parts (breast, thyroid) and lungs
IEC 60601-1	Type BF working part Short mode of operation (max. 10 min. work and 10 min. rest) Equipment with internal power supply Probe head: IPX7
Software version	300L: SW version M 1.02 Mobile device: 300L App version M1.02
Battery Charger SONON	Used to charge the lithium-ion battery. Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Dimensions: 115mm*65mm*27mm(± 0.2 mm)
AC adapter for the battery charger SONON	Used to connect to the charger of SONON system. input: (100 ♦ 240) V, 0,5 A, (50 ♦ 60) Hz output (USB port): 8,4 V dc; 1,5 A; Class of protection II Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Cable length 120 cm (± 1 cm). Dimensions: 90mm*50 mm *30 mm (± 0.2 mm)
AC adapter (for SONON 300L)	Used to connect to the charger of SONON 300L system. Material: Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABS NH-1015. Cable length 150 cm (± 1 cm). Dimensions: 90mm*43 mm *30 mm (± 0.2 mm)
Power cord	Used to connect the AC adapter for the battery charger SONON and network adapter for SONON to the power source. 100-240 V, 50/60 Hz. WSP-300. Cable length: 200cm (± 1 cm)
Protective cap for SONON	Is used to protect the probe head of SONON 300L Material: polyoxymethylene POM-C


	Dimensions: 85 mm *35 mm *26 mm (± 1 mm).
Hand strip for SONON	Is used to carry the SONON 300L (if necessary). Length: 30cm(± 1 cm) Material: polyamide PA66-GF15, polyetherane YU-50D
Fastener for SONON	Ensures fixation of the scanner SONON 300L (if necessary). Dimensions, mm: 550 mm *145 mm *75 mm (± 1 mm). Materials: Scanner fastener - Mixture of polycarbonate and acrylonitrile butadiene styrene (PC+ABS)PC-ABSNH-1015 Bracket - aluminum ENAW 7005 Support - steel 1.0143 Handle - polyamide P A6 brand CB-30, stainless steel AISI 201

9.5. Reclamations

Complaints about inadequate quality should be addressed to:
HEALCERION Co., Ltd
804ho, 38-21, Digital-ro 31-gil, Guro-gu, Seoul, 08376, Korea
Telephone: +(82) 70-7582-6326
E-mail: info@healcerion.com

Authorized representative of the manufacturing company HEALCERION Co., Ltd in Russia is OOO Firm "Solo"
Address: 620137, Russia, Sverdlovsk region, Ekaterinburg, Blyukhera St., house 75, building1, room № 3-17
Telephone: (343) 351-76-16
E-mail: mail@solo-trade.ru

Chapter 10 Glossary-Abbreviations

Chapter 10

IEC	Meaning
a	Acoustic Attenuation Coefficient
Aaprt	-12db Output Beam Area
CMI	Normalizing Coefficient
Deq	Equivalent Aperture Diameter
d-6	Pulse Beam Width
deq	Equivalent Beam Diameter
awf	Acoustic Working Frequency
Ipa	Pulse-Average Intensity
Ipa,a	Attenuated Pulse-Average Intensity
Ipi	Pulse-Intensity Integral
Ipi,a	Attenuated Pulse-Intensity Integral
I _{ta} (z)	Temporal-Average Intensity
I _{ta,a} (z)	Attenuated Temporal-Average Intensity
I _{zpta} (z)	Spatial-Peak Temporal-Average Intensity
I _{zpta,a} (z)	Attenuated Spatial-Peak Temporal-Average Intensity
MI	Mechanical Index
P	Output Power
P _a	Attenuated Output Power
P ₁	Bounded Output Power
p _i	Pulse Pressure Squared Integral
p _r	Peak-Rarefactional Acoustic Pressure
P _{ra}	Attenuated Peak-Rarefactional Acoustic Pressure
prr	Pulse Repetition Rate
TI	Thermal Index
TIB	Bone Thermal Index
TIC	Cranial-Bone Thermal Index
TIS	Soft-Tissue Thermal Index
td	Pulse Duration / (same)
X, Y	-12 dB Output Beam Dimensions / (same)
z	Distance from the Source to a Specified Point / (same)
zbp	Depth for TIB / Depth at which the relevant index is maximum
zbp	Break-Point Depth / (same)
zs	Depth for TIS / Depth at which the relevant index is maximum

FDA	Meaning
MI	the Mechanical Index.
TISscan	the Soft Tissue Thermal Index in an auto-scanning mode.
TISnon-scan	the Soft Tissue Thermal Index in a non-autoscanning mode.
TIB	the Bone Thermal Index.
TIC	the Cranial Thermal Index.
Aaprt	the area of the active aperture (square centimeters).
pr.3	the derated peak rarefactional pressure associated with the transmit pattern giving rise to the value reported under MI (megapascals).
Wo	the ultrasonic power, except for TISscan, in which case it is the ultrasonic power passing through a one centimeter window (milliwatts).
W.3(z1)	the derated ultrasonic power at axial distance z1 (milliwatts).
ITA.3(z1)	the derated spatial-peak temporal-average intensity at axial distance z1 (milliwatts per square centimeter).
z1	the axial distance corresponding to the location of $\max[\min(W.3(z), ITA.3(z) \times 1 \text{ cm}^2)]$, where $z \geq z_{bp}$ (centimeters).
zbp	$1.69 \sqrt{A_{aprt}}$ (centimeters).
zsp	the axial distance at which TIB is a global maximum (i.e., $z_{sp} = z_{B.3}$) (centimeters).
z@PII.3max	the axial distance corresponding to the maximum of the derated spatial-peak pulse intensity integral (megapascals).
deq(z)	the equivalent beam diameter as a function of axial distance z. It is equal to $[(4/\pi)(W_o/ITA(z))]^{0.5}$ where $ITA(z)$ is the temporal-average intensity as a function of z (centimeters).
fc	the center frequency (MHz). For MI, f_c is the center frequency associated with the transmit pattern giving rise to the global maximum reported value of MI. For TI, for combined modes involving transmit patterns of unequal center frequency, f_c is defined as the overall range of center frequencies of the respective transmit patterns.
Dim. of Aaprt	the active aperture dimensions for the azimuthal (x) and elevational (y) planes (centimeters).
PD	the pulse duration (microseconds) associated with the transmit pattern giving rise to the reported value of MI.
PRF	the pulse repetition frequency associated with the transmit pattern giving rise to the reported value of MI (Hz).
pr@PIImax	the peak rarefactional pressure at the point where the free-field, spatial-peak pulse intensity integral is a maximum (megapascals). See Section 6.2.4.1 of the Output Display Standard, entitled "Measurement Methodology for Mechanical and Thermal Indices".
deq@PIImax	the equivalent beam diameter at the point where the free-field, spatial-peak pulse intensity integral is a maximum (centimeters). See Section 6.2.5.1 of the Output Display Standard, entitled "Measurement Methodology for Mechanical and Thermal Indices".
FL	the focal length, or azimuthal (x) and elevational (y) lengths, if different (centimeters).
IPA.3@MImax	the derated pulse-average intensity at the point of global maximum reported MI (watts per square centimeter).

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Copyright© Все права принадлежат Healcerion Co., Ltd.