Manufacturer's Statement

To whom this may concern:

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Technical Manual (IFU) for Iforia ICDs

Is an original document issued by BIOTRONIK SE & Co. KG.
The content of English and Russian document versions is identical.


Name:
Position:

Elena Polikarpova
Associate Regulatory Affairs

Company:
Street:
City / Country:

BIOTRONIK SE & Co. KG
Woermannstraße 1
12359 Berlin - Germany

Date and Signature:

12. Feb. 2018 

p.p. Elena Polikarpova

BIOTRONIK
excellence for life

Russian Federation

Application for registration

Approved

BIOTRONIK SE & Co. KG
Woermannkehre 1, 12359 Berlin, GERMANY

Associate Regulatory Affairs

Elena Polikarpova

 **BIOTRONIK**
excellence for life
BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin
(signature, stamp)

«12» February 2018

Implantable Cardioverter-Defibrillator Single-chamber Iforia
Implantable Cardioverter-Defibrillator Dual-chamber Iforia
Implantable Cardioverter-Defibrillator Triple-chamber Iforia

Technical Manual (IFU)


Table of Contents

Product Description	1
Intended Medical Use	1
System Overview	2
Therapeutic and Diagnostic Functions	5
General Safety Instructions	6
Operating Conditions	6
Possible Complications	6
Possible Risks	7
Implantation	8
Implantation Procedure	8
Precautionary Measures while Programming	11
Magnet Response	12
Follow-up	13
• Patient Information	13
Replacement Indications	14
Explantation and Device Replacement	14
Parameters	18
Bradycardia / CRT	15
Tachycardia	18
Sensing	19
Diagnostics	19
Home Monitoring	20
Technical Data	20
Mechanical Characteristics	20
Electrical Characteristics	20
Battery Data	22
Legend for the Label	24

Product Description

Intended Medical Use

Intended use

Iforia 3/5/7 is part of a family of implantable cardioverter-defibrillators (ICDs). The primary objective of the therapy is to prevent sudden cardiac death. Furthermore, the device is capable of treating bradycardia arrhythmias and cardiac resynchronization therapy with multisite ventricular pacing.

An ICD provides symptomatic therapy with the following objectives:

- Termination of spontaneous ventricular fibrillation (VF) through shock delivery
- Termination of spontaneous ventricular tachycardia (VT) through antitachycardia pacing (ATP); in case of ineffective ATP or hemodynamically not tolerated VT, with shock delivery
- Cardiac resynchronization through multisite ventricular pacing (triple-chamber devices)
- Compensation of bradycardia through ventricular (single-chamber devices) or AV sequential pacing (DX, dual- and triple-chamber devices)

Diagnosis and therapy forms

The device monitors the heart rhythm and automatically detects and terminates cardiac arrest resulting from ventricular tachyarrhythmia. All major therapeutic approaches from the field of cardiology and electrophysiology are included. BIOTRONIK Home Monitoring® enables physicians to perform therapy management at any time.

Required expertise

In addition to having basic medical knowledge, the user must be thoroughly familiar with the operation and the operation conditions of a device system.

- Only qualified medical specialists having this required special knowledge are permitted to use implantable devices.
- If users do not possess this knowledge, they must be trained accordingly.

Indications

Iforia can treat life-threatening ventricular arrhythmias with antitachycardia pacing and defibrillation.

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We recommend observing the Indications published by the German Cardiac Society (Deutsche Gesellschaft für Kardiologie, Herz- und Kreislaufforschung) and the ESC (European Society of Cardiology). This also applies to the guidelines published by the Heart Rhythm Society (HRS), the American College of Cardiology (ACC), the American Heart Association (AHA), and other national cardiology associations.

Single-chamber and dual-chamber

Single-chamber and dual-chamber ICDs are indicated for patients with the following risk:

- Sudden cardiac death caused by ventricular arrhythmias

Triple-chamber

Triple-chamber ICDs are indicated for patients with the following risks:

- Sudden cardiac death caused by ventricular arrhythmias
- Congestive heart failure with ventricular asynchrony

Also indicated for primary prophylaxis in congestive heart failure patients is Iforia.

Contraindications

Known contraindications:

- Tachyarrhythmia caused by temporary or reversible irritation, e.g. poisoning, electrolyte imbalance, hypoxia, sepsis or acute myocardial infarction
- Such frequent VT or VF that the therapies would cause an unacceptably rapid depletion of the device batteries
- VT with few or without clinically relevant symptoms
- VT or VF treatable by surgery
- Concomitant diseases that would substantially limit a positive prognosis
- Accelerated intrinsic rhythm

System Overview

Device family

The complete Iforia 3/5/7 device family consists of several device types with a DF-1/IS-1 or DF4/IS-1 connection.

Single-chamber: VR-T and VR-T DX (only devices with a DF-1/IS-1 connection); dual-chamber: DR-T; triple-chamber: HF-T. Not all device types are available in every country.

Device

The device's housing is made of biocompatible titanium, welded from outside and thus hermetically sealed. The ellipsoid shape facilitates implantation in the pectoral muscle area. The connections for bipolar pacing and sensing (and unipolar connections for the triple-chamber device) as well as for shock delivery are found in the device header. The housing serves as a potential antipole during shock delivery or in the case of unipolar lead configuration.

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Note: The device type DX can only be connected using a DF-1/IS-1 connector.

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
Note: A device's DF4 connector port may only be used for connecting leads with a DF4 connector that conform to ISO 27186.

Leads

BIOTRONIK leads are sheathed with biocompatible silicone. They can be flexibly maneuvered, are stable long-term, and are equipped for active or passive fixation. They are implanted using a lead introducer set. Some leads are coated with polyurethane which is known to increase the sliding properties for the lead. Leads with steroids reduce inflammatory processes. The fractal design of the electrodes provides for low pacing thresholds. BIOTRONIK provides adapters to connect already implanted leads to new devices.

Telemetry

Telemetric communication between the device and the programmer can be carried out following initialization either by applying the programming head (PGH) to the device or by using wandless telemetry in the programmer. BIOTRONIK calls this function SafeSync®.

Programmer

Implantation and follow-up are performed with BIOTRONIK's portable programmer. Software PSW as of version 1204.A.

There are programmers with integrated or external SafeSync Module for wandless telemetry.

Leadless ECG, IEGM, markers and functions are displayed simultaneously on the color display.

The programmer allows you to determine the thresholds and to perform all tests during an in-office follow-up; in addition, you can change the permanent program and send it to the implanted device.

In addition to this, the programmer is used to set mode and parameter combinations, as well as for interrogation and saving of data from the device.

Modes

The mode setting depends on the individual diagnosis:

Device type	Modes
VR	WI; WIR; V00; OFF
VR DX	VDD; VDDR; VDI; VDIR; WI; WIR; V00; OFF
DR, HF	DDD; DDDR; DDI; DDIR; VDD; VDDR; VDI; VDIR WI; WIR; AAI; AAIR; V00; D00; OFF

NBD and NBG codes

VE is the NBD code used for the antitachycardia modes of the single-chamber, dual-chamber and triple-chamber devices:

V	Shock in the ventricle
V	Antitachycardia pacing (ATP) in the ventricle
E	Detection via IEGM analysis

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D	Pacing in the atrium and ventricle
D	Sensing in the atrium and ventricle
D	Pulse inhibition and pulse triggering
R	Rate adaptation

D	Pacing in the atrium and ventricle
D	Sensing in the atrium and ventricle
D	Pulse inhibition and pulse triggering
R	Rate adaptation
V	Multisite pacing in both ventricles

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V	Pacing in the ventricle
D	Sensing in the atrium and ventricle
D	Pulse inhibition and pulse triggering
R	Rate adaptation

WVIR is the NBG code for the antibradycardia pacing modes of the single-chamber device:

V	Pacing in the ventricle
V	Sensing in the ventricle
I	Pulse inhibition in the ventricle
R	Rate adaptation

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- With Home Monitoring, diagnostic and therapeutic information as well as technical data are automatically sent to a stationary or mobile transmitter via an antenna in the device header. The data are encrypted and sent from the transmitter to the BIOTRONIK Service Center via the cellular phone network.
- The received data are deciphered and evaluated. Each physician can set the criteria for evaluation to be used for each patient and can configure the time of notification via e-mail, SMS or fax.
- A clear overview of the results of this analysis is displayed for the attending physicians on the protected Internet platform Home Monitoring Service Center (HMSC).

- Data transmission from the device is performed with a daily device message.
- Device messages which indicate special events in the heart or in the device are forwarded immediately.
- A test message can be initiated at any time using the programmer to immediately check the Home Monitoring function.

Order numbers for Moria with DF-1/IS-1 or DF4/IS-1 connection

Not all device types are available in every country:

	Moria 3	
	DF-1	DF4
VR-T	383586	383588
DR-T	383570	383572
HF-T	383554	383556

Order numbers for Moria ProMRI with DF-1/IS-1 or DF4/IS-1 connection

Not all device types are available in every country:

	Moria 5 ProMRI		Moria 7 ProMRI	
	DF-1	DF4	DF-1	DF4
VR-T	390119	390121	390083	390089
VR-T DX	390123	—	390095	—
DR-T	390115	390117	390069	390075
HF-T	390111	390113	390056	390062

Package contents

The storage package includes the following:

- Sterile packaging with device
- Serial number label
- Patient ID card
- Warranty booklet

Note: The technical manual pertaining to the device is either included in hard copy form in the storage package or in digital form on the internet.

The sterile packaging includes the following:

- Device, blind plug DP-1 (if applicable) and blind plug IS-1 for device type HF
- Screwdriver

Therapeutic and Diagnostic Functions

Diagnostic functions

- Data from implantation and the most recent interrogations and follow-ups are recorded as well as arrhythmia episodes; they are stored together with other data to assess patients and the state of the device at any time.
- To check the lead for proper functioning, an automatic impedance measurement using subthreshold pacing pulses is performed in the device.
- Leadless ECG function: for all device types, far-field derivation can be measured without external leads between the right ventricular distal shock coil and housing, which, depending on the implantation site, corresponds to ECG derivation II or III (Einthoven).
- Once a telemetry connection has been established during a test procedure in ag in-office follow-up, the leadless ECG and the IEGM are displayed with markers.

Antitachycardia pacing

- The ICD can treat ventricular tachycardia with antitachycardia pacing (ATP); ATP can also be delivered in the VF zone (ATP One Shot) when the stability criterion indicating that this will be effective before shock delivery (monomorphic rapid VTs) is met.
- Depending on the device type, the software contains not only the ICD functions but also all pacemaker functions for 1, 2, or 3 chambers. The heart rhythm is continuously monitored; each arrhythmia is classified according to the heart rate and the adjustable detection criteria. Depending on the preset values, antibradycardia as well as antitachycardia therapy is inhibited or delivered.

Cardioversion, defibrillation

- The ICD can treat ventricular tachyarrhythmia with cardioversion and/or defibrillation. Shock polarity and energy can be programmed individually. Shock energies between 2.0 J and 40 J are possible. Before delivery of the shock, the ICD can be set to only deliver a shock when ongoing tachyarrhythmia is confirmed; during this time period the device can identify spontaneous conversion of the tachyarrhythmia and cancel the charging process if necessary.
- The shock paths can be set between the different shock coils (SVC/AA) and/or the housing.

Antibradycardia pacing and CRT

- Innovative rate hystereses, automatic sensor functions, and a night program promote the patient's intrinsic rhythm, avoid overdrive pacing, and facilitate adaptation of the device to the individual needs of the patient.
- Setting an upper tracking rate for the atrium prevents unspecific atrial pacing, thus reducing the risk of pacemaker-mediated tachycardia.
- Positive AV hystereses support intrinsic conduction and thus the natural contraction sequence. Negative AV hystereses support the cardiac resynchronization therapy by maintaining pacing in stress situations.
- For resynchronization of the ventricles, triple-chamber devices have functions for multisite ventricular pacing with possible VV delays in either direction.
- To ensure that no additional surgery is necessary in case of a left-sided increase of the pacing threshold or undesired phrenic nerve stimulation, different pacing polarities can be set for the left ventricular lead with a triple-chamber device.
- Automatic active capture control is available for the right and left ventricle with automated tracking of the pacing threshold or automatic threshold monitoring (ATM) for trend analysis.

With the Iforia 3, the pacing threshold is only measured.

Storing programs

The parameter settings can be saved in 3 individual therapy programs.

Home Monitoring functions

- The device automatically sends information to the transmitter once a day. It also sends messages related to events, which are immediately forwarded to the Service Center. In addition to this, test messages can be initiated using the programmer.
- Appointments for Home Monitoring-supported follow-ups can be scheduled via the HMSC. This applies to Iforia 5/7.
- Important medical information in the device messages include the following:
 - Atrial and ventricular arrhythmias
 - Parameters relevant to leads in the atrium and ventricle: pacing thresholds, sensing amplitudes, impedances
 - Current statistics
 - IEGM online RD® with up to 3 high definition channels

Operating Conditions

Technical manuals

The following technical manuals provide information about usage of the device systems:

- Technical manual for the device
- Technical manual for the HMSC
- Technical manuals for leads
- Technical manuals for the programmer and its accessories
- Technical manuals for the user interface
- Technical manuals for cables, adapters and accessories
- Technical manuals are either included in hard copy form in the storage package or in digital form on the internet: manuals.biotronik.com.
- Follow all relevant technical manuals.
- Keep technical manuals for later use.

Care during shipping and storage

- Devices must not be stored or transported close to magnets or sources of electromagnetic interference.
- Note the effects of the storage period; see Battery Data.

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Temperature

Extremely low and high temperatures affect the service time of the battery in the device.

- Permitted for shipping, storage and operation are -15°C to $+45^{\circ}\text{C}$.

The device and its accessories have been guaranteed for sterility if the following conditions are met: the blister and quality control seal have not been damaged.

Sterile packaging

The device and screwdriver are packaged in two separately sealed blisters. The inner blister is also sterile on the outside so that it can be transferred in a sterile state during implantation.

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- Do not use the device if the package is damaged.
- The device must not be resterilized and reused.

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Skeletal myopotentials

Bipolar sensing and control of sensitivity are adapted by the device to the rate spectrum of intrinsic events so that skeletal myopotentials are usually not recorded. Skeletal myopotentials can nonetheless be classified as intrinsic events especially at very high sensing sensitivity and, depending on the interference, may cause inhibition or anti-arrhythmia therapy.

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- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimal.
- Due to the intensity and variety of EMI, there is no guarantee for safety. It is generally assumed that EMI produces only minor symptoms, if any, in patients.
- Depending on the pacing mode and the type of interference, sources of interference may lead to pulse inhibition or triggering, an increase in the sensor-dependent pacing rate or asynchronous pacing.
- Under unfavorable conditions, for example during therapeutic or diagnostic procedures, interference sources may induce such a high level of energy into the pacing system that the cardiac tissue surrounding the lead tip is damaged.

Device behavior in case of EMI

In case of electromagnetic interference, the device switches to asynchronous pacing for as long as the interference rate is exceeded.

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- Therapeutic ultrasound: Harm to the patient via excess warming of body tissue near the device system
- Transcutaneous electrical nerve stimulation
- Hyperbaric oxygen therapy
- Applied pressures higher than normal pressure

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If electrical current from an external source is conducted through the body for diagnostic or therapeutic purposes, then the device can be subjected to interference, which can place the patient at risk:

Arrhythmia or ventricular fibrillation can be induced during diathermic procedures such as electrocautery, HF ablation or HF surgery. For example, damaging heat can result during lithotripsy. Influences on the device are not always immediately clear. If risky procedures cannot be avoided, the following should be observed at all times:

- Electrically insulate the patient.
- Switch off the ICD's detection function. The pacemaker function can remain active. The device may need to be switched to asynchronous modes for this.
- Do not introduce energy near the device system.
- Additionally check the peripheral pulse of the patient.
- Monitor the patient during and after every intervention.

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. Specifically, the current induced in the implanted leads may result in necrotic tissue formation close to the el-+electrode/tissue interface. As a result, sensing properties and pacing thresholds may change.

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- Adhere to instructions for risky therapy and diagnosis procedures.
- Shield device against radiation.
- After applying radiation, double-check the device system to make sure it is functioning properly.

Note: Please contact BIOTRONIK with questions during the risk/benefit analysis.

Magnetic resonance imaging

For implantable devices without ProMRI labeling, Magnetic resonance imaging (MRI) is contraindicated due to the high frequency fields and the associated magnetic flux density: damage or destruction of the device system by strong magnetic interaction and damage to the patient by excessive warming of the body tissue in the area surrounding the device system.

Under certain conditions and when maintaining mandatory measures to protect the patient and device system, magnetic resonance imaging can be performed. BIOTRONIK devices with the "MR conditional" function bear the identification ProMRI.

- The ProMRI manual – MR conditional device systems – contains detailed information on safely conducting an MRI.
 - Download the digital manual from the web site:
www.biotronik.com/manuals/manualselection
 - Order the printed manual from BIOTRONIK.
- Does approval as "MR-Conditional" apply in your country or region?
Request current information from BIOTRONIK.

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- BIOTRONIK device with blind plug and screwdriver
- BIOTRONIK leads and lead introducer set
 - Single-chamber device: one bipolar ICD lead with 1 or 2 shock coils for the ventricle
 - Dual-chamber device: one bipolar lead for the atrium and one bipolar ICD lead for the ventricle with 1 or 2 shock coils
 - Triple-chamber device: an additional unipolar or bipolar LV lead
- Approved connections are DF-1, DF4 and IS-1. Use only adapters approved by BIOTRONIK for leads with different connections or leads from other manufacturers.
- BIOTRONIK programmer (with integrated SafeSync RF telemetry or with separate SafeSync Module) and approved cables
- External multi-channel ECG device
- Keep spare parts for all sterile components.

Keeping an external defibrillator ready

To be able to respond to unforeseeable emergencies or possible technical failures of the device:

- Keep an external defibrillator and paddles or adhesive electrodes ready.

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- Peel the sealing paper off of the outer blister at the marked position in the direction indicated by the arrow. The inner blister must not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments!
- Take hold of the inner blister by the gripping tab and take it out of the outer blister.
- Peel the sealing paper off of the sterile inner blister at the marked position in the direction indicated by the arrow.

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- Check for damage before and after unpacking all parts.
- Replace damaged parts.
- The ICD is shipped with tachyarrhythmia therapy deactivated and is only to be connected and implanted in this state.
- Leads must not be shortened.

Implantation site

- Depending on lead configuration and the patient's anatomy, the ICD is generally implanted subpectorally on the left side.

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Avoiding damage to the header

There is a blind plug for DF-1 and IS-1 connections in the header. The provided set screws must be carefully loosened or tightened.

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- Either leave unused ports closed with the prebunched blind plugs, or close them using the supplied blind plugs.

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- Clean lead connectors with a sterile cloth.
- Rinse connection only with sterile water.

Overview: Implanting

1	Prepare the vein.
2	Implant the leads, perform the measurements, and fixate the leads.
3	Form the device pocket.
4	Connect the lead connector to the device.
5	Insert the device.
6	Guide the fixation suture through the opening in the header and fixate the device in the prepared device pocket.
7	Close the device pocket.
8	Check the device with standard tests.

Connecting the device

The lead connectors are connected to the ports in the header of the device.

1	Disconnect stylets and stylet guides.
2	DF-1/IS-1 connection: • Connect the DF-1 connector for the right-ventricular shock coil to RV. • Connect the DF-1 connector for the supraventricular shock coil to SVC. Or connect a subcutaneous array to SVC. DF4/IS-1 connection: • Connect the DF4 connector to RV
3	DF-1/IS-1 connection: • Connect the bipolar IS-1 connector for the atrium to RA. • Connect the bipolar IS-1 connector for the right ventricle to RV. • Connect the unipolar or the bipolar IS-1 connector for the left ventricle to LV. DF4/IS-1 connection: • Connect the bipolar IS-1 connector for the atrium to RA. • Connect the unipolar or the bipolar IS-1 connector for the left ventricle to LV.

4	Push the lead connector into the header without twisting or bending the connector or conductor until the connector tip (on the DF-1 connector) or the insertion indicator (on the DF4 connector) becomes visible behind the set screw block. This indicator can vary depending on the manufacturer of the lead used.
5	If you cannot easily plug the lead connector into the connection: • Use only sterile water as lubricant.
6	If the lead connector cannot be inserted completely, the set screw may be protruding into the drill hole of the set screw block. • Use the screwdriver to perpendicularly pierce through the slitted point in the center of the silicone plug until it reaches the set screw. • Carefully loosen the set screw without completely unscrewing it, so that it does not become tilted upon retightening.
7	Turn the set screw clockwise until torque control starts (you will hear a clicking sound).
8	Carefully withdraw the screwdriver without retracting the set screw. • In case of IS-1 connections with 2 set screws, tighten both screws! • When the screwdriver is withdrawn, the silicone plug automatically seals the lead connection safely.

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- The distance between 2 shock coils must be greater than 6 cm.
- Tip and ring electrodes must not have contact with each other.

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- Switch on RF telemetry on the programmer.
- Apply the programming head for about 2 s until successful initialization is displayed on the programmer:

The SafeSync symbol is displayed in the navigator and the signal strength is displayed in the status line.

- Remove the programming head.

Activating ICD therapy

- Load the device program that is suitable for the device type in the programmer.
- Activate ICD therapy.
- Shipment mode is permanently deactivated once the leads have been connected and initial measurement of the pacing impedance has been performed. The device data are saved.
- Take precautionary measures while programming.
- If the device induces tachycardia while programming ATPs or does not deliver adequate therapy in the DET test: use emergency shock or an external defibrillator.

Precautionary Measures while Programming

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Cancelling telemetry

Programmer interference or interrupted telemetry during performance of temporary programs (follow-up tests) can result in inadequate pacing of the patient. This is the case if the programmer can no longer be operated due to a program error or a defective touch screen and therefore the temporary program cannot be terminated.

Under these circumstances, it is helpful to cancel telemetry, in which case the device automatically switches to the permanent program.

- In the case of telemetry with programming head: tilt the PGM by at least 80 cm.
- In the case of RF telemetry: switch off and reposition the programmer.
- Turn off possible sources of interference.

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Checking for electrodes suitable for the shock path

Three different shock paths can be set. Two of these form an electrical path to the housing of the implanted device.

- For the RV → SVC shock path, a second shock coil must be available (dual shock coil).

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The asynchronous modes VOO and DOO can only be set if tachyarrhythmia sensing is deactivated. This would leave the patient without sensing and therefore without ICD therapy.

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BIOCRONIK devices feature several functions to prevent conduction of atrial tachycardia to the ventricle(s):

- Set mode switching for indicated patients.
- Set the upper rate and the refractory periods to prevent abrupt ventricular rate switching.
- Give preference to Wenckebach response and avoid 2:1 behavior.
- Set all parameters so as to prevent constant changing between atrial and ventricular-controlled modes.

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- The following post-shock program parameters can be adjusted: post-shock duration, basic rate, rate hysteresis, ventricular pacing, LV-T-wave protection, triggering, AV delay (fixed, not dynamic).
- The default settings for the post-shock program are as follows:
A and RV: 7.5 V, 1.5 ms
LV: settings from the permanent program

Phrenic nerve stimulation that cannot be terminated

In rare cases, chronic phrenic nerve stimulation cannot be terminated by reprogramming the available left ventricular pacing configuration or using other measures.

- Set a right ventricular mode both in the permanent program as well as the ATP, in the post-shock program and for mode switching if need be.

Avoiding risks in the case of exclusive LV pacing

Lead dislodgement in the case of exclusive left ventricular pacing could pose the following risks: loss of ventricular pacing and ATP therapy, induction of atrial arrhythmias.

- Consider sensing and pacing parameters with reference to loss of therapy.
- Exclusive LV pacing is not recommended for patients who depend on the device.
- Please note that active capture control is not available.
- In the case of follow-ups and threshold tests, take loss of synchronized ventricular pacing into consideration.
- Mode switching and post-shock do not allow for exclusive LV pacing. Please note the effects when programming mode switching and the post shock parameters.

Recognizing lead failure

Automatic impedance measurement is always switched on.

- Impedance values that indicate technical failure of a lead are documented in the event list.

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- Do not establish unnecessary RF telemetry.
- After 5 minutes without input, SafeSync switches to the economy mode.
- Check the battery capacity of the device at regular intervals.

Magnet Response

Application of the programming head when ICD therapy is set

If a connected programming head is applied and is communicating with the programmer and ICD therapy is permanently set, detection and therapy remain intact except during the diagnostic tests. If ICD therapy is not set as permanent, no therapy is delivered when the programming head is applied.

Programming head application

When the programming head is applied, time remains for device interrogation and for reprogramming or deactivation of the therapy before the device switches back to the previously set permanent therapy mode. The same applies to programming head application to establish RF telemetry contact.

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- If detection interruptions of longer than 8 hours are required, the magnet has to be briefly removed from the device. The 8 hour countdown restarts when the magnet is applied again.
- Use BIOTRONIK magnets: type M-50 permanent magnets.

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- The first follow-up should be carried out by the physician using the programmer [in-office follow-up] approximately 3 months after implantation following the lead ingrowth phase.
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

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Monitoring using the Home Monitoring function does not serve to replace regular in-office appointments with the physician required for other medical reasons. Follow-up supported by Home Monitoring can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted despite use of the Home Monitoring function if symptoms worsen or if new symptoms arise.
- Device messages are transmitted regularly.
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device system are sufficient. If not, an in-office follow-up needs to be carried out.

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1	Record and evaluate the ECG.
2	Interrogate the device.
3	Evaluate the status and automatically measured follow-up data.
4	Check the sensing and pacing functions.
5	Possibly evaluate statistics and IEGM recordings.
6	Manually perform standard tests if necessary.
7	Possibly customize program functions and parameters.
8	Transmit the program permanently to the device.
9	Print and document follow-up data (print report).
10	Finish the follow-up for this patient.

Patient Information

Patient ID card

A patient ID card is included in delivery.

- Provide the patient with the patient ID.
- Request that patients contact the physician in case of uncertainties.

Prohibitory signs

Places with prohibitory signs must be avoided.

- Draw the patient's attention to prohibitory signs.

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Possible battery levels

- BOS: Beginning of Service: > 70% charge
- MOS 1: Middle of Service: 70% to 40% residual charge
- MOS 2: Middle of Service: < 40% residual charge
- ERI: Elective Replacement Indication, (i.e. RRT: Recommended Replacement Time)
- EOS: End of Service

Elective Replacement Indication (ERI)

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- Replace device soon.
- The device can monitor the heart rhythm for at least 3 more months.
- At least 6 maximum energy shocks can be delivered until EOS occurs.
- The selected parameters in the device program do not change.

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- Replace device immediately.
- Monitor patient constantly until immediate replacement of the device!

- VT and VF detection and all therapies are deactivated!
- The antibradycardia function remains active in the VVI mode:
 - Ventricular pacing: RV; basic rate 50 bpm; without special pacemaker functions such as hysteresis, etc.
 - Pulse amplitude of 6 V; pulse width of 1.5 ms
 - Time of transmission for Home Monitoring: 90 days

Explantation and Device Replacement

Explantation

- Interrogate the device status.
- Deactivate VT and VF therapies prior to explantation.
- Remove the leads from the header. Do not simply cut them loose.
- Use state-of-the-art techniques to remove the device and, if necessary, the leads.

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- Deactivate VT and VF therapies prior to device replacement.
- Insulate connections that are not used.

Basic principles:

- The device must not be resterilized and reused.

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Disposal

BIOTRONIK takes back used products for the purpose of environmentally safe disposal.

- Clean the explant with an at least 1% sodium hypochlorite solution.
- Rinse off with water.

Bradycardia / CRT

General ICD therapy

Parameter	Range of values	Standard	VR	DX	DR	HF
ICD therapy	OFF; ON	ON	x	x	x	x
Programs	Display standard program; Display safe program; Display first interrogated program; Individual 1,2,3	-	x	x	x	x

Timing: Basic rate day/night and rate hystereses

Parameter	Range of values	Standard	VR	DX	DR	HF
Basic rate	30 ... (5) ... 100 ... (10) ... 160 bpm	40 bpm	x	x		
		60 bpm			x	x
Night rate	OFF; 30 ... (5) ... 100 bpm	OFF	x	x	x	x
Night begins	00:00 ... (00:01) ... 23:59 hh:mm	22:00 hh:mm	x	x	x	x
Night ends		06:00 hh:mm				
Rate hysteresis	OFF; -5 ... (-5) ... -25 ... (-20) ... -65 bpm	OFF	x	x	x	x
Scan/repetitive	OFF; ON	ON	x	x	x	x

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Parameter	Range of values	Standard	VR	DX	DR	HF
AV dynamics	Low; Medium; High; Fixed; (Individual)	Low		x	x	x
AV delay (1 or 2) after:						
- Pacing	15; 40 ... (5) ... 350 ms	-			x	x

- Sensing	Either automatic: AV delay after pacing + sense compensation Or: 40 ... (5) ... 350 ms	-	x	x	x	
- At rate 1	50 ... (10) ... 130 bpm	60 bpm				
- At rate 2	60 ... (10) ... 140 bpm	130 bpm				
Sense compensation	OFF; -5 ... (-5) ... -120 ms	-40 ms			x	x
AV hysteresis mode	OFF; Positive; Negative; IRSplus	OFF	x	x		
	OFF; Positive; Negative	OFF				x
AV hysteresis (positive)	70; 110; 150; 200 ms	70 ms	x	x	x	
AV hysteresis (negative)	10 ... (10) ... 150 ms	50 ms	x	x	x	
AV scan/repetitive (positive)	OFF; ON	ON	x	x	x	

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Parameter	Range of values	Standard	VR	DX	DR	HF
Maximum sensor rate	80 ... (10) ... 160 bpm	160 bpm	x	x	x	x
Sensor gain	AUTO; Very low; Low; Medium; High; Very high	Medium	x	x	x	x
Sensor threshold	Very low; Low; Medium; High; Very high	Medium	x	x	x	x
Rate increase	1; 2; 4; 8 bpm/cycle	2 bpm	x	x	x	x
Rate decrease	0.1; 0.2; 0.5; 1.0 bpm/cycle	0.5 bpm	x	x	x	x

Parameter	Range of values	Standard	VR	DX	DR	HF
Pulse amplitude A	0.5 ... (0.25) ... 4.0 ... (0.5)	2.5 V			x	x
Pulse amplitude V/RV	... 6.0; 7.5 V		x	x	x	x
Pulse amplitude LV						x
Pulse width A	0.4; 0.5 ... (0.25) ... 1.5 ms	0.4 ms			x	x
Pulse width V/RV			x	x	x	x
Pulse width LV		0.4 ms				x

Pacing: Ventricular capture control

The following parameters apply to Iforia 5/7:

Parameter	Range of values	Standard	VR	DX	DR	HF
Ventricular capture control	OFF; ATM; ON	ATM	x	x	x	x
Threshold test start	2.5 ... (0.5) ... 5.0 V	ATM: 2.5 V ON: 3.5 V	x	x	x	x
Minimum amplitude	1.0 ... (0.25) ... 4.0 V	1.0 V	x	x	x	x
Safety margin	1.0; 1.2 V	1.0 V	x	x	x	x

The following parameters apply to the Iforia 3:

Parameter	Range of values	Standard	VR	DX	DR	HF
Capture control	OFF; ATM	ATM	x	x	x	x

Parameter	Range of values	Standard	VR	DX	DR	HF
Atrial capture control	OFF; ATM	ATM			x	x

LV lead configuration

Parameter	Range of values	Standard	VR	DX	DR	HF
Pacing polarity LV	LV tip -> LV ring; LV tip -> RV ring; LV ring -> LV tip; LV ring -> RV ring; UNIP	LV-Tip -> RV-Ring				x
Sensing polarity LV	UNIP; BIPL	UNIP				x

MRI program

Valid for devices with ProMRI:

Parameter	Range of values	Standard	VR	DX	DR	HF
Mode	V00; OFF	OFF	x	x		
	V00; D00; OFF				x	x
Basic rate	70 ... (10) ... 160 bpm	90 bpm			x	x
Ventricular pacing	RV; BV	RV				x

Detection

Parameter	Range of values	Standard	VR	DX	DR	HF
Interval AT/AF	280 ... 600 ms	300 ms		x	x	x
Interval VT1	OFF; 270 ... (10) ... 600 ms	OFF	x	x	x	x
Interval VT2	OFF; 270 ... (10) ... 500 ms					
Interval VF	OFF; 240 ... (10) ... 400 ms	300 ms				
Detection counter VT1	10 ... (2) ... 60	26	x	x	x	x
Detection counter VT2	10 ... (2) ... 40	16				
Detection counter VF	6 out of 8; 8 out of 12; 10 out of 14; 12 out of 16; 6 out of 20; 18 out of 24; 20 out of 26; 22 out of 30; 24 out of 30	8 out of 12				
Redetection counter VT1	10 ... (2) ... 30	20	x	x	x	x
Redetection counter VT2		14				
SMART detection VT1/VT2	OFF; ON	ON		x	x	x
SMART detection ON:						
- Onset VT1/VT2	4 ... (4) ... 32 %	20 %		x	x	x
- Stability VT1/VT2	8 ... (4) ... 48 %	12 %				
SMART detection OFF:						
- Onset VT1/VT2	OFF; 4 ... (4) ... 32 %	20 %	x	x	x	x
- Stability VT1/VT2	OFF; 8 ... (4) ... 48 ms	24 ms				
Sustained VF	OFF; 1; 2; 3; 5; 10; 20; 30 min	OFF	x	x	x	x
Forced termination	OFF; 1 ... (1) ... 10 min	1 min		x	x	x

Parameter	Range of values	Standard	VR	DX	DR	HF
Attempt	OFF; 1 ... (1) ... 10	OFF	x	x	x	x
ATP type for VT1/VT2	Burst; Ramp	Burst	x	x	x	x
ATP type for VF	OFF; Burst; Ramp	Burst	x	x	x	x
ATP optimization	OFF; ON	OFF	x	x	x	x
Number S1 for VT1/VT2	1 ... (1) ... 10	5	x	x	x	x
Number S1 for VF		8				
S1 decrement for VT1/VT2 and for VF	5 ... (5) ... 40 ms	10 ms	x	x	x	x
Scan decrement	OFF; 5 ... (5) ... 40 ms	OFF	x	x	x	x
Add. S1 for VT1/VT2	OFF; ON	ON	x	x	x	x
Ventricular pacing for VT1/VT2	RV; LV; BV	RV				x
Ventricular pacing for VF		RV				
R-S1 interval for VT1/VT2	70 ... (5) ... 95%	80%	x	x	x	x
R-S1-Intervall for VF		85%				

Therapy: Shock

Parameter	Range of values	Standard	VR	DX	DR	HF
Number of shocks VT1/VT2	0; 1; 2; 6; 8	8	x	x	x	x
Number of shocks VF	6; 8	8	x	x	x	x
1. shock for VT1/VT2	OFF; 2 ... (2) ... 20 ... (5) ... 40 J	40 J	x	x	x	x
2. shock for VT1/VT2	OFF; 4 ... (2) ... 20 ... (5) ... 40 J	40 J	x	x	x	x
3. - nth shock for VT1/VT2	OFF; 4*40 J; 6*40 J	6*40 J	x	x	x	x
1. shock for VF	2 ... (2) ... 20 ... (5) ... 40 J	40 J		x	x	x
2. shock for VF	4 ... (2) ... 20 ... (5) ... 40 J	40 J		x	x	x
3. - nth shock for VF	4*40 J; 6*40 J	6*40 J	x	x	x	x

For AT/AF	OFF; ON	ON			x	
- Confirmation	Normal; Reverse; Alternating	Normal				
- Shock waveform	Biphasic; Biphasic 2	Biphasic				
- Shock path	RV -> ICD+SVC RV -> ICD RV -> SVC	RV -> ICD+SVC RV -> ICD	x		x	x
		RV -> ICD		x		

Sensing

Sensitivity and thresholds

Parameter	Range of values	Standard	VR	DX	DR	HF
Sensing A	STD; OFF	Std.		x	x	x
Sensing RV	STD; TWS; VFS; IND	Std.	x	x	x	x
Sensing LV	Std.; OFF; IND	Std.				x
Upper threshold RV	50; 75%	50%	x	x	x	x
Upper threshold LV	50; 75%	50%				x
Upper threshold duration after detection	110; 150 ... [50] ... 500 ms VFS: 110 ms	350 ms	x	x	x	x
Upper threshold duration after pacing		400 ms				
Lower threshold RV	25; 50%	25%	x	x	x	x
T-wave suppression after pacing	OFF; ON	OFF	x	x	x	x
Minimum threshold A	0.2 ... [0.1] ... 2.0 mV	0.4 mV		x	x	x
Minimum threshold RV	0.5 ... [0.1] ... 2.5 mV	0.8 mV	x	x	x	x
Minimum threshold LV	0.5 ... [0.1] ... 2.5 ... [0.5] ... 5.0 mV	1.6 mV				x

Parameter	Range of values	Standard	VR	DX	DR	HF
For AT/AF	OFF; ON For Iltia 7: Extended DN	ON		x	x	x
For SVT	OFF; ON	ON		x	x	x
Periodic recording	When Home Monitoring OFF: OFF: 30 ... (30) ... 180 days	90 days	x	x	x	x
IEGM configuration	RA, RV, LV RA, RV, FF FF, RV, LV	RA, RV, LV				x

For Iltia 7 the following can be set:

Parameter	Range of values	Standard	VR	DX	DR	HF
Start resting period	00:00 ... [1:00 AM] ... 23:00 hh:mm	2:00 hh:mm	x	x	x	x
Duration of resting period	0.5 ... (0.5) ... 12 h	4 h	x	x	x	x
AV delay adjustment sensing test	OFF; 300 ms	300 ms		x	x	x
Thoracic Impedance (TI)	OFF; ON	OFF	x	x	x	x

Parameter	Range of values	Standard	VR	DX	DR	HF
Home Monitoring	OFF; ON	OFF	x	x	x	x
Time of transmission	STD; 00:00 ... (01:00) ... 23:00 hh:mm	STD	x	x	x	x
IEGM for:						
- therapy episodes	OFF; ON	ON	x	x	x	x
- Monitoring episodes						
Sustained atrial episode	OFF; 6; 12; 18 h	12 h		x	x	x

Only applies to Iforia 3

Parameter	Range of values	Standard	VR	DX	DR	HF
Cycle duration	OFF; 30; 60; 90; 120; 180 days	30 days	x	x	x	x
Transmission on:	XX.XX.XXXX	Follow-up + 7 days	x	x	x	x

Mechanical Characteristics

Housing

Devices with a DF-1/IS-1 header:

Type	Connector port	W x H x D in mm	Volume cm ³	Mass g
VR	DF-1	65 x 55 x 11	33	80
DX, DR	DF-1	65 x 55 x 11	33	81
HF	DF-1	65 x 58.5 x 11	34	83

Devices with a DF4/IS-1 header:

Type	Connector port	W x H x D in mm	Volume cm ³	Mass g
VR	DF4	65 x 52 x 11	31	80
DR	DF4	65 x 56 x 11	32	82
HF	DF4	65 x 56 x 11	33	82

Materials in contact with body tissue

- Housing: Titanium
- Header: Epoxy resin
- Blind plug and silicone plug: Silopren or silastic
- DF4 seal: Silastic

X-ray identification

NT

Electrical Characteristics

Standards

The specifications are made according to EN 45502-2-2:2008.

Measuring conditions

If not indicated otherwise, all specifications refer to the following conditions:

- Ambient temperature: 37°C ± 2°C
- Pacing/sensing: 500 Ω ± 1%
- Shock: 50 Ω ± 1%


Factory settings

- Arrhythmia zones VT1, VT2, VF: OFF
- Antibradycardia pacing: OFF
- Home Monitoring: OFF

Telemetry data

- Nominal carrier frequency: 403.6 MHz
- Maximum power of transmission: < 25 µW [-16 dBm]

Pulse waveform

The pacing pulse has the following form:
The pulse amplitude reaches its maximum value at the beginning of the pulse (U_a). With increasing pacing duration (t_p), the pulse amplitude is reduced dependent on the pacing impedance.

Resistance to interference

- Note on device type DX (only devices with a DF-1/IS-1 connection): The EMC requirements are met as long as atrial sensitivity is set to 1.0 mV (factory settings) or values > 1.0 mV. Measures must be taken to assure interference-free therapy if more sensitive values are set.
- Note on device type HF: In the case of unipolar sensing, the requirement for interference voltages of < 0.3 mV (peak to peak) is met.

Common mode rejection ratio

Rate	Common mode rejection ratio			
	Atrium: DX*	Atrium: DR, HF	V right: VR, DR, HF	V left: HF
16.6 Hz	58 dB	53 dB	64 dB	66 dB
50 Hz	55 dB	55 dB	64 dB	66 dB
60 Hz	56 dB	56 dB	64 dB	68 dB

*Only devices with a DF-1/IS-1 connection.

ATP amplitude

A burst was measured at 500 D, an amplitude of 7.5 V (tolerance ± 1.5 V), pulse width of 1.5 ms, R-S1 interval of 300 ms and a number S1 of 5:

ATP amplitude	Measured minimum	Measured maximum	Mean value
RV	7.67 V	7.67 V	5.00 V
LV	7.67 V	7.67 V	4.99 V

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Manufacturer:	GREATBATCH, INC. Clarence, NY 14031	LITRONIK GmbH & Co 01796 Pirna, Germany
Battery type	GB 2992	LiS 3410 RA/RR
Battery ID number shown on the programmer	3	4
Device type	VR, (DX), DR, HF	
Battery voltage at ERI	2.5 V	2.85 V
Charge time at BOS	8 s	8 s
Charge time at ERI	10 s	10 s
Usable capacity until ERI	Iforia 3/5: 1390 mAh Iforia 7: 1600 mAh	1390 mAh
Usable capacity until EOS	1730 mAh	1520 mAh

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- Devices should be implanted within 19 months between the manufacturing date and the use by date (indicated on the package).
- If the ICD is implanted shortly before the use by date, the expected service time may be reduced by up to 16 months.

Calculation of service times

- The services times have been calculated as follows – in all chambers depending on the device type:
 - Pulse amplitude: 2.5 V
 - Pulse width: 0.4 ms
 - Pacing impedance: 500 Ω
 - Basic rate: 60 bpm
 - Home Monitoring: ON, 1 device message each day and 12 IEGM online HD transmissions per year
 - Diagnostic functions and recordings: permanently set
- Capacitor reforming is performed 4 times per year and therefore at least 4 maximum charges for shocks have to be assumed per year even if less than 4 are delivered.

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Calculation of the number of shocks:

Service time (in years) x number of shocks per year

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	Manufacturing date		Use by
	Storage temperature	REF	Order number
SN	Serial number	PID	Product identification number
	Dangerous voltages	CE	CE mark
	Contents		Follow the instructions for use!
manuals.biotronik.com 		Follow the electronically available instructions for use!	

STERILE EO Sterilized with ethylene oxide			
	Do not resterilize		Single use only. Do not re-use!
	Do not use if packaging is damaged		Non-sterile

	Transmitter with non-ionizing radiation at designated frequency
	MR conditional: patients having a device system implanted whose components are labeled with this symbol on the packaging, can be examined using an MR scan under precisely defined conditions.

	Device: NBG code and compatible leads
	Factory settings for therapy: OFF
	Screwdriver

	Example of DF-1/IS-1 header
	Examples of DF-1/IS-1 or DF4/IS-1 header
	Bipolar IS-1 connector
	Unipolar IS-1 connector
	Unipolar DF-1 connector
	DF4 connector

© BIOTRONIK SE & Co. KG
All rights reserved. Specifications subject to modification, revision and improvement.

© All product names in use may be trademarks or registered trademarks held by BIOTRONIK or the respective owner.

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
Manufacturer's Statement

To whom this may concern:

We hereby guarantee that the enclosed document:


Addition to the Instruction for Use (Iforia)

is an original document issued by BIOTRONIK SE & Co. KG.
The content of English and Russian document versions is identical.

Name:	Elena Polikarpova
Position:	Associate Regulatory Affairs
Company:	BIOTRONIK SE & Co. KG
Street:	Woermannkehre 1
City / Country:	12359 Berlin – Germany

Date and Signature:

04. July 2018

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin

p.p. Elena Polikarpova

Approved

BIOTRONIK SE & Co. KG
Woermannkehe 1, 12359 Berlin, GERMANY

BIOTRONIK
excellence for life

Associate Regulatory Affairs

BIOTRONIK SE & Co. KG
Woermannkehe 1
12359 Berlin


(signature, stamp)

Elena Polikarpova

« 04 » July 2018

Implantable Cardioverter-Defibrillator Single-chamber Iforia
Implantable Cardioverter-Defibrillator Dual-chamber Iforia
Implantable Cardioverter-Defibrillator Triple-chamber Iforia

Addition for the Instruction For Use (IFU)

1. Name of the medical device

- Implantable Single-chamber Cardioverter / Defibrillator Iforia, models:

1. Implantable Single-chamber Cardioverter / Defibrillator Iforia 3 VR-T

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 3 VR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

2. Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

3. Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T DX

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T DX;
- Screwdriver TW;
- Instruction for Use
- Addition to the Instruction for Use

4. Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

5. Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T DX

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T DX;
- Screwdriver TW;
- Instruction for Use
- Addition to the Instruction for Use

- Implantable Dual-chamber Cardioverter / Defibrillator Iforia, models:

1. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 3 DR-T

Basic scope of delivery:

- Implantable Dual -chamber Cardioverter / Defibrillator Iforia 3 DR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

2. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 5 DR-T

Basic scope of delivery:

- Implantable Dual -chamber Cardioverter / Defibrillator Iforia 5 DR-T

- Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
3. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 7 DR-T
Basic scope of delivery:
 - Implantable Dual -chamber Cardioverter / Defibrillator Iforia 7 DR-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
- Implantable Triple-chamber Cardioverter / Defibrillator Iforia, models:
1. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 3 HF-T
Basic scope of delivery:
 - Implantable Triple -chamber Cardioverter / Defibrillator Iforia 3 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
 2. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 5 HF-T
Basic scope of delivery:
 - Implantable Triple -chamber Cardioverter / Defibrillator Iforia 5 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
 3. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 7 HF-T
Basic scope of delivery:
 - Implantable Triple -chamber Cardioverter / Defibrillator Iforia 7 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use

(Below – ICD, ICD Iforia 3/5/7, implant, device).

2. Information about the manufacturer

Manufacturer:

BIOTRONIK SE & Co. KG
 Woermannkehre 1, 12359 Berlin, Germany.
 Tel.: + 49 (0) 30 68905 – 1221.
 E-mail: tomas.zeiske@biotronik.com

Production Site Contact

BIOTRONIK SE & Co. KG,
 Woermannkehre 1, 12359 Berlin,
 Germany;

Official representative in Russia

LLC

"Biotronik",

Moscow,

Russia,

109028, Moscow, Zemlyanoy Val, 50A/8,

bld. 2. Phone: +7 (495) 789-68-31.

Fax: +7 (495) 789-68-32.

Web: www.biotronik.ru

3. Intended Use

Iforia 3/5/7' is a family of implantable cardioverters/defibrillators (ICDs). Primary objective of these ICDs is the invasive cardiac arrhythmias treatment, brady- and congestive heart failure therapy performing, prevention of sudden cardiac death by termination of ventricular tachyarrhythmias.

4. Indication for use

ICD Iforia can treat life-threatening ventricular arrhythmias and congestive heart failure with cardiac resynchronization therapy with multisite ventricular pacing. The implantation of an ICD is a symptomatic therapy with the following objectives:

- Termination of spontaneous ventricular fibrillation (VF) through shock delivery;
- Termination of spontaneous ventricular tachycardia (VT) through antitachycardia pacing (ATP); in case of ineffective ATP or hemodynamically not tolerated VT, with shock delivery;
- Cardiac resynchronization through multisite ventricular pacing (triple-chamber devices);
- Compensation of bradycardia through ventricular (single-chamber devices) or AV sequential pacing.

Generally approved differential diagnostics methods, indications, and recommendations for ICD therapy apply to BIOTRONIK devices. See the guidelines of cardiology associations for guidance.

We recommend observing the indications published by the:

- German Cardiac Society (Deutsche Gesellschaft für Kardiologie, Herz- und Kreislaufforschung);
- ESC (European Society of Cardiology);
- Heart Rhythm Society (HRS);
- American College of Cardiology (ACC);
- American Heart Association (AHA),

as well as other national cardiology associations.

Single- and Dual-chamber ICD Iforia 3/5/7' are indicated for patients with the following risk:

- Sudden cardiac death caused by ventricular arrhythmias;

Triple-chamber ICD Iforia 3/5/7' are indicated for patients with the following risk:

- Sudden cardiac death caused by ventricular arrhythmias;
- Congestive heart failure with ventricular asynchrony

Also, ICD Iforia 3/5/7' are indicated for prevention of sudden cardiac death in patients with congestive heart failure

5. Contraindications

Known contraindications:

- Tachyarrhythmia caused by temporary or reversible irritation, e.g., intoxication, electrolyte imbalance, hypoxia, sepsis, or acute myocardial infarction;
- VT or VF of such frequent occurrence that the therapies would cause an unacceptably rapid depletion of the ICD batteries;
- VT with few or without clinically relevant symptoms;
- VT or VF treatable by surgery;
- Concomitant diseases that would substantially limit a positive prognosis after implantation;
- Patients with accelerated idioventricular rhythms.

6. Possible complications and side effects while using the product

Possible complications and side effects while using the product can be:

- Tachyarrhythmia acceleration caused by the device;
- Bleeding;
- Air embolism;
- Phenomenon of rejection by the body, including local tissue reaction;
- Heart perforation;
- Cardiac tamponade;
- Death;
- Pressure sores at the site of implantation;
- The growth of the connective tissue;
- Atrial fibrillation or other;
- Formation of hematomas or cysts;
- Rupture of the heart wall or vein wall;
- Infection;
- Abrasion and breakage of the electrode;
- Myocardial and / or neurostimulation;
- The excitability of the myocardium;
- Pericardial effusion;
- Pneumothorax;
- Rejection phenomena (local tissue reaction, connective tissue formation, device migration);
- Thromboembolism;
- Thrombosis;
- Valve damage (particularly in underdeveloped heart with anatomical structures);
- Perforation of a vein or the heart;
- Chronic nerve injury;
- Endocarditis;
- Bedsores in the projection area;
- Extrusion;
- Accumulation of fluid;
- Heart block;
- Hematoma / seroma;
- Formation of keloid scar;

- Migration / electrode displacement;
- Injury infarction;
- Detection of myopotentials;
- Noise pericardial friction;
- The possibility of death due to the inactivated bradyfunction (stimulation performing);
- Increase the stimulus threshold on the electrodes;
- Thrombolytic and air embolism;
- Transvenous thrombosis due to the electrode;
- Blockage of the vein.

IN CASE OF ARISING OF ONE OR MORE OF THE ABOVE-DESCRIBED COMPLICATIONS, IT REQUIRES AN IMMEDIATE PERFORMING OF CORRESPONDING MEDICAL INTERVENTIONS THAT MUST LEAD TO COMPLICATION /-S CORRECTION.

Additional possible complications associated with the use of ICD systems are listed below:

- Inappropriate ICD shock therapy;
- The possibility of death due to the impossibility of defibrillation;
- Resetting electrical energy or isolating the myocardium during defibrillation.

In patients carrying frequent shocks, regardless of medical indications, psychological intolerance to the ICD system may develop, among which are:

- Addiction;
- Depression;
- Fear of premature battery discharge;
- Fear of shock, despite the realization of its necessity;
- Fear of the likelihood of the device performing shock therapy;
- Imaginary shock.

7. Precautions while using the product

Must be excluded the following diagnoses of the patient when making decisions about implantation of this type of device:

- Chronic renal failure;
- Allergic reactions to contrast agent and intolerance to iodine;
- Severe coagulopathy, severe anemia;
- Uncontrolled arterial hypertension;
- Fever and acute infectious diseases;
- Endocarditis;
- Severe underlying non-cardiological disease;
- Decompensation of heart failure and pulmonary edema;
- Intoxication with glycosides, expressed hypokalemia.

8. Possible complications

General information on medical complications

Complications for patients and device systems generally recognized among practitioners also apply to BIOTRONIK devices.

- Normal complications may include fluid accumulation within the device pocket, infections, or tissue reactions. Primary sources of complication information include current scientific and technological knowledge.
- It is impossible to guarantee the efficacy of antitachycardia therapy, even if the programs have proven successful during tests or subsequent electrophysiological examinations. In rare cases the set parameters may become ineffective. It is possible for therapies to induce or accelerate tachycardia and cause sustained ventricular flutter or fibrillation.

Skeletal myopotentials

Bipolar sensing and control of sensitivity are adapted by the device to the rate spectrum of intrinsic events so that skeletal myopotentials are usually not recorded. Skeletal myopotentials can nonetheless be classified as intrinsic events especially at very high sensing sensitivity and, depending on the interference, may cause inhibition or antiarrhythmia therapy.

In the case of undesired myopotentials, the device switches to asynchronous pacing if the interference rate is exceeded.

Possible technical failures

Technical failure of a device system cannot be entirely ruled out. Possible causes can include the following:

- Lead dislodgement, lead fracture
- Insulation defects
- Device component failures
- Battery depletion
- Interrupted telemetry

Electromagnetic interference (EMI)

Any device can be sensitive to interference if external signals are sensed as intrinsic rhythm or if measurements prevent rate adaptation.

- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimal;
- Due to the intensity and variety of EMI, there is no guarantee for safety. It is generally assumed that EMI produces only minor symptoms, if any, in patients;
- Depending on the pacing mode and the type of interference, sources of interference may lead to pulse inhibition or triggering, an increase in the sensor-dependent pacing rate or asynchronous pacing;
- Under unfavorable conditions, for example during therapeutic or diagnostic procedures, interference sources may induce such a high level of energy into the pacing system that the cardiac tissue surrounding the lead tip is damaged.

Device behavior in case of EMI

In the case of electromagnetic interference, the device switches to asynchronous pacing for as long as the interference rate is exceeded.

Static magnetic fields

The reed switch in the device closes starting at a field strength of 1.8 mT. The reed switch opens if the magnetic field falls below 1 mT.

9. Possible risks

Methods to avoid:

The following procedures are contraindicated as they may cause harm to the patient or damage the device and, as a result, put the system functionality at risk:

- Therapeutic ultrasound: Harm to the patient via excess warming of body tissue near the device system;
- Transcutaneous electrical nerve stimulation;
- Hyperbaric oxygen therapy;
- Applied pressures higher than normal pressure.

Risky therapeutic and diagnostic procedures

If electrical current from an external source is conducted through the body for diagnostic or therapeutic purposes, then the device can be subjected to interference, which can place the patient at risk.

Arrhythmia or ventricular fibrillation can be induced during diathermic procedures such as electrocautery, HF ablation or HF surgery. For example, damaging heat can result during lithotripsy. Influences on the device are not always immediately clear.

If risky procedures cannot be avoided, the following should be observed at all times:

- Electrically insulate the patient;
- Switch off the ICD's detection function. The pacemaker function can remain active. The device may need to be switched to asynchronous modes for this;
- Do not introduce energy near the device system;
- Additionally, check the peripheral pulse of the patient;
- Monitor the patient during and after every intervention.

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. Specifically, the current induced in the implanted leads may result in necrotic tissue formation close to the electrode/tissue interface. As a result, sensing properties and pacing thresholds may change.

- Place adhesive electrodes anterior-posterior or perpendicular to the axis formed by the device to the heart at least 10 cm away from the device and from implanted leads.

Radiation therapy

The use of radiation therapy is contraindicated due to possible damage to the device and the resulting impaired functional safety. If this type of therapy is to be used anyway, prior risk/benefit analysis is absolutely necessary. The complexity of influencing factors such as different sources of radiation, a variety of devices and therapy conditions makes it impossible to issue directives that guarantee radiation therapy without an impact on the device. The EN 45502 standard pertaining to active implantable medical devices requires the following measures during the administration of therapeutic ionizing radiation:

- Adhere to instructions for risky therapy and diagnosis procedures;
- Shield device against radiation;
- After applying radiation, double-check the device system to make sure it is functioning properly.

Note: If you have questions regarding the evaluation of risks and benefits, please contact

BIOTRONIK. MRI - Magnetic resonance imaging¹

For implantable devices without ProMRI labelling Magnetic Resonance Imaging (MRI¹) is contraindicated due to the high frequency fields and the associated magnetic flux density: damage or destruction of the device system by strong magnetic interaction and damage to the patient by excessive warming of the body tissue in the area surrounding the device system.

For BIOTRONIK devices with the "MR conditional" function that have an identification ProMRI, MRI can be performed under certain conditions and when maintaining mandatory measures to protect the patient and device system

The ProMRI manual – MR conditional device systems – contains detailed information on safely conducting an MRI.

- Download the digital manual from the web site: manuals.biotronik.com
- Order the printed manual from BIOTRONIK.

10. Information about the impact on the ability to control the vehicle

For all patients with ICD, self-management of the vehicle is a relative contraindication, since a device shocked device to a patient after the detection of life-threatening tachycardia can cause pain, which will lead to loss of control over the road.

Product application features for different categories of patients

Pregnant - Implanted cardioverter-defibrillator has no effect to the pregnancy and lactation. If there is a need for a x-ray investigation performing of the device, care should be taken to minimize exposure to the fetus and the mother;

Nursing mothers - although the respective tests were performed on the biocompatibility of ICD, there is no any quantitative assessment of the presence of leachables in breast milk;

Geriatric patients - most of patients (about 71%), involved in receiving ICDs in clinical studies, were older than 60 years;

¹ For information concerning MRI compatibility of different ICD variants please see Chapter „MRI-compatibility“ on page 14 of this document

Patients with disabilities and people with disabilities - special caution is required when using this device for patients in wheelchairs with the electric drive or for the patients who are using other electrical (external or implanted) devices.

11. Application methods

Required professional experience

All actions with ICD 'Iforia 3/5/7', associated with implantation procedure, the observation of the patient in the postoperative period, remote monitoring system data analysis, device permanent program parameters setting, revision can be performed only by skilled health personnel that must be thoroughly familiar with the operation and the operation conditions of a device system.

- Only qualified medical specialists having this required special knowledge are permitted to use implantable devices;
- If users do not possess this knowledge, they must be trained accordingly.

Storage conditions

Temperature	5 - 45 °C
Relative humidity	30 - 75 %
Pressure	375 - 795 mm Hg

Operating conditions

Product	Temperature	Humidity or other special conditions
ICD Iforia 3/5/7	32 - 42 °C	Resistant to the effects of biological fluids (blood)
The screwdriver (torque wrench) TW	5 - 45 °C	30 - 75 %

12. Classification of the product

Type of contact with the human body

#	Product	Type of contact with the human body
1	ICD Iforia 3/5/7	Implantable device with permanent contact with the internal environment and tissues of the human body
2	The screwdriver (torque wrench) TW	A product that has a short-term mediated contact with the human body

Class according to the potential risk of using a medical device	Third risk class
By multiplicity of application	Single-use only
Sterility	The medical device is sterile.
Protection class	CF

13. Description of principles that underpin the work of the medical device

Iforia 3/5/7 ICDs are generally implanted subpectorally on the left side for the entire lifetime of the device. Before the operation, there is follow-up of the patient: examination and consultation attending the physician, surface ECG registration at rest and during exercise, daily monitoring of Holter ECG. These diagnostic procedures can accurately determine the future parameters settings of the implantable device permanent program.

Implantation procedure is performed under normal environmental conditions.

The implantation is a standard procedure. The preparatory stage involves treatment of the surgical field, analgesia using local anesthetics. Drug solution infiltrate into the skin and deeply lying tissue. General anesthesia is used in exceptional cases.

Next, the physician made the incision 5-6 cm long in the subclavian area.

Then performed a puncture of the subclavian vein in an amount corresponding to the number of implanted endocardial lead / leads. They are introduced into the lumen of the subclavian vein and under X-ray control are positioned in the respective chambers of the heart.

After fixing the lead / leads, the final stage of ICD Iforia 3/5/7 implantation is to create a pocket for the device. Programming the device is individually tailored to the patient's needs, device model features and arrhythmology status of the patient.

Forms of diagnosis and therapy

The ICD shall control the heart rhythm, so it will automatically recognize a cardiac arrest because of ventricular tachyarrhythmias and restore the heart. All significant reasons for undertaking cardiac and electrophysiological therapy are available. BIOTRONIK Home Monitoring® allows physicians to provide around the clock monitoring therapy.

14. Scope of delivery

- Implantable Single-chamber Cardioverter / Defibrillator Iforia, models:

1. Implantable Single-chamber Cardioverter / Defibrillator Iforia 3 VR-T

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 3 VR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

2. Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T
- Screwdriver TW
- Instruction for Use
- Addition to the Instruction for Use

3. Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T DX

Basic scope of delivery:

- Implantable Single-chamber Cardioverter / Defibrillator Iforia 5 VR-T DX;
- Screwdriver TW;
- Instruction for Use
- Addition to the Instruction for Use

4. Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T
Basic scope of delivery:
 - Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
5. Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T DX
Basic scope of delivery:
 - Implantable Single-chamber Cardioverter / Defibrillator Iforia 7 VR-T DX;
 - Screwdriver TW;
 - Instruction for Use
 - Addition to the Instruction for Use
-
- Implantable Dual -chamber Cardioverter / Defibrillator Iforia, models:
 1. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 3 DR-T
Basic scope of delivery:
 - Implantable Dual-chamber Cardioverter / Defibrillator Iforia 3 DR-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
 2. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 5 DR-T
Basic scope of delivery:
 - Implantable Dual-chamber Cardioverter / Defibrillator Iforia 5 DR-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
 3. Implantable Dual-chamber Cardioverter / Defibrillator Iforia 7 DR-T
Basic scope of delivery:
 - Implantable Dual-chamber Cardioverter / Defibrillator Iforia 7 DR-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
- Implantable Triple-chamber Cardioverter / Defibrillator Iforia, models:
 1. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 3 HF-T
Basic scope of delivery:
 - Implantable Triple-chamber Cardioverter / Defibrillator Iforia 3 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use
 2. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 5 HF-T
Basic scope of delivery:
 - Implantable Triple -chamber Cardioverter / Defibrillator Iforia 5 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use

3. Implantable Triple-chamber Cardioverter / Defibrillator Iforia 7 HF-T
- Basic scope of delivery:
- Implantable Triple-chamber Cardioverter / Defibrillator Iforia 7 HF-T
 - Screwdriver TW
 - Instruction for Use
 - Addition to the Instruction for Use

15. Delivery package

Delivery package consists of:


- The implant and screwdriver TW in a sterile package;
- Sticker with the serial number;
- Patient ID-card;
- The warranty booklet;
- Instruction for Use
- Addition to the Instruction for Use

Instruction for use (Technical manual) and Addition to the Instruction for Use are provided in printed versions on every shipment of ICD


16. ICD Family. Model versions

ICD Iforia 3/5/7 family includes the models with different types of connectors: DF-1/IS-1 and DF4/IS-1.


Single-chamber ICDs: VR-T and VR-T DX

Dual-chamber ICDs: DR-T

Triple-chamber ICDs: HF-T

MRI-compatibility

Iforia ICD family includes nonMRI and MRI-compatible devices:

Below is presented the list of models with- or without ProMRI function:


Single-chamber ICD		Dual-chamber ICD		Triple-chamber ICD	
Iforia 3 VR-T	Iforia 5 VR-T Iforia 5 VR-T DX Iforia 7 VR-T Iforia 7 VR-T DX	Iforia 3 DR-T	Iforia 5 DR-T Iforia 7 DR-T	Iforia 3 HF-T	Iforia 5 HF-T Iforia 7 HF-T
None	ProMRI	None	ProMRI	None	ProMRI

17. General ICD description

The device's housing is made of biocompatible titanium, welded from outside and thus hermetically sealed. The ellipsoid shape facilitates implantation in the pectoral muscle area. The connections for bipolar pacing and sensing (and unipolar connections for the triple-chamber device) as well as for shock delivery are found in the **device header**. **The housing serves as** a potential antipole during shock delivery or in the case of unipolar lead configuration.

18. Electrical Functional Description

The main electrical components of the 'Iforia 3/5/7' ICDs are feedthroughs, high voltage capacitors, electronic module, RF-antenna and battery.

Battery: left (blue)

Electronic Module: middle (green)

High Voltage Capacitors: right (brown)

Feedthroughs: top (red)

Connector	Lead Connector	Configuration	Implantation site	Device type
RV	DF-1	Shock coil	Right Ventricle	VR, VR DX, DR, HF
SVC	DF-1	Shock coil	Superior vena cava	VR, VR DX, DR, HF
RA	IS-1	Bipolar	Atrium	VR DX, DR, HF
(R)V	IS-1	Bipolar	(Right) ventricle	VR, VR DX, DR, HF
LV	IS-1	Unipolar, Bipolar	Left Ventricle	HF




BIOTRONIK provides ICDs with headers for different standardized lead connections: DF-1/IS-1 and DF4/IS-1.

Note: Device and leads must be compatible.

- The device type DX can only be connected using a DF-1/IS-1 connector.

Lead connection DF-1/IS-1

The device labeling provides information pertaining to possible lead connections depending on the device type and pertaining to connection assignment:

VR	DX	DR	HF
DF-1 SVC DF-1 RV  IS-1 RV	DF-1 SVC DF-1 RV  IS-1 RA IS-1 RV	DF-1 SVC DF-1 RV  IS-1 RA IS-1 RV	DF-1 SVC DF-1 RV  IS-1 LV IS-1 RA IS-1 RV

Lead connection DF4/IS-1

The device labeling provides information pertaining to possible lead connections depending on the device type and pertaining to connection assignment:

VR	DR	HF
 DF4-LLHH RV	IS-1 RA  DF4-LLHH RV	IS-1 RA  IS-1 LV DF4-LLHH RV

Connector	Lead connector	Configuration	Implantation site	Device type
RA	IS-1	Bipolar	Atrium	DR, HF
LV	IS-1	Unipolar, bipolar	Left ventricle	HF
RV, SVC	DF4	Bipolar and shock	Right Ventricle	VR, DR, HF

Note: Suitable leads must meet the standards.

- The device's DF4 connector port may only be used for connecting leads with a DF4 connector that conform to ISO 27186.

Technical specifications of 'Iforia 3/5/7' ICDs portfolio models

Function	Iforia 3	Iforia 5	Iforia 7
Bradycardia therapy (pacing, sensing, postshock pacing, rate adaptation, etc)	✓	✓	✓
Tachycardia therapy (ATP and shock therapy)	✓	✓	✓
Home Monitoring	✓	✓	✓
RF telemetry (wandless telemetry)	✓	✓	✓
DF-1 connector	✓	✓	✓
DF4 connector	✓	✓	✓
Atrial automatic threshold measurement	--	✓	✓
Right / left ventricular automatic threshold tracking	--	✓	✓
Atrial NIPS	---	---	✓
Intrathoracic Impedance measurement	---	✓	✓
Programmer based scheduling of Home Monitoring- supported follow ups	✓	---	---
Remote scheduling / HMSC based scheduling of Home Monitoring-supported follow ups	--	✓	✓

19. Description of accessories and devices combinations

Firmware

The functionality of the 'Iforia 3/5/7' ICDs is controlled by embedded software. The same firmware (FW 2.0.0) is used in all 'Iforia 3/5/7' ICDs. The embedded software is stored in the ROM and RAM of the IC 'CIC_1'.

The ROM FW application of CIC_1 used in all 'Iforia 3/5/7' devices implements the back-up functionality (back up pacing, VF detection, delivery of defibrillation shocks with maximum energy, communication of the implant with the programmer and with the CardioMessenger II/II-S Home Monitoring patient device) and manual emergency shocks of the ICDs.

The RAM firmware serves to control the clinical functionality of the 'Iforia 3/5/7' ICDs like bradycardia therapy, VT/VF detection, atrial detection, ventricular therapy sequences, tachycardia therapy (ATP and shock), EP test functionality, holter and statistics, Home Monitoring, remote scheduling, wireless wand, EOS detection.

Screwdriver (torque wrench)

The screwdriver (torque wrench) TW is supplied sterile, complete with 'Iforia 3/5/7' ICDs and is used to tighten the screws of the connectors of the leads for a secure fit in the corresponding connectors on the housing of the ICD.

Technical Specifications:

Length, mm	65
Width, mm	13
Height, mm	13
Mass, g	5,13

Leads

BIOTRONIK leads are sheathed with biocompatible silicone. They can be flexibly maneuvered, are stable long-term, and are equipped for active or passive fixation. They are implanted using a lead introducer set. Some leads are coated with polyurethane which is known to increase the sliding properties for the lead. Leads with steroids reduce inflammatory processes. The fractal design of the electrodes provides for low pacing thresholds. BIOTRONIK provides adapters to connect already implanted leads to new devices. The leads are supplied separately.

20. Programming

'Iforia 3/5/7' ICDs can be interrogated and programmed with BIOTRONIK programmer, models: ICS 3000 and Renamic. The data transmission between the ICDs 'Iforia 3/5/7' and the BIOTRONIK programmers is performing using the programming head/wand or via wireless RF-telemetry.

Function	Purpose
Programming and testing functions	Program BIOTRONIK pacemakers and ICDs during the implantation procedure or follow-ups. The follow-up can be performed using the PGH or via wireless telemetry with the implants that have a corresponding RF-module on-board.
ECG recording and ECG monitoring	Display and printout of up to three leads of surface ECGs, as well as up to three intracardiac derivations (IEGMs) and the corresponding event markers, in real-time.
Data management	Store parameter values and ECG/IEGM recordings for computer-aided archiving and evaluation.
Data export	Print out follow-up reports using the internal and/or external printer

21. Programmer

Programmer is a device that is used for programming the permanent program of ICD.

Programming during the implantation and the follow-ups after that can be performed using ICS3000 and Renamic programmers of BIOTRONIK SE & CO. Wireless telemetry can be supported initially or can be realized using the separate SafeSync module. The programmer is used during implantation to send the program to the device. You can define thresholds values and to perform all tests during the follow-up. In addition, the programmer is used as a data manager of all the exported data from the implantable devices. saving e data in the memory of the implant. The programmer displays wireless ECG and IEGM (intracardiac electrogram), marker annotation, and functions.

The programmer is supplied separately.

22. Modes

The mode setting depends on the individual diagnosis.

Device type	Modes
VR	VVI; VVIR; V00; OFF
VR DX	VDD; VDDR; VDI; VDIR; VVI; VVIR; V00; OFF
DR, HF	DDD; DDDR; DDI; DDIR; VDD; VDDR; VDI; VDIR; VVI; VVIR; AAI; AAIR; V00; D00; OFF

23. NBD and NBG codes

VVE is the NBD code used for the antitachycardia modes of the single-chamber, dual-chamber and triple-chamber devices:

V	Shock in the ventricle
V	Antitachycardia pacing (ATP) in the ventricle
E	Detection via IEGM analysis

DDDR is the NBG code for the antibradycardia mode of the dual-chamber devices:

D	Pacing in the atrium and ventricle
D	Sensing in the atrium and ventricle
D	Pulse Inhibition and pulse triggering
R	Rate Adaptation

DDDRV is the NBG code for the antibradycardia mode of the triple-chamber devices:

D	Pacing in the atrium and ventricle
D	Sensing in the atrium and ventricle
D	Pulse Inhibition and pulse triggering
R	Rate Adaptation
V	Multisite pacing in both ventricles

VDDR is the NBG code for the antibradycardia mode of the single-chamber DX devices:

V	Pacing in the ventricle
D	Sensing in the atrium and ventricle
D	Pulse Inhibition and pulse triggering
R	Rate Adaptation

VVIR is the NBG code for the antibradycardia pacing modes of the single-chamber device:

V	Ventricular in the ventricle
V	Sensing in the ventricle
I	Pulse inhibition in the ventricle
R	Rate adaptation

24. BIOTRONIK Home Monitoring

In addition to effective pacing therapy, BIOTRONIK provides a complete therapy management system:

- With Home Monitoring, diagnostic and therapeutic information as well as technical data are automatically sent to a stationary or mobile transmitter via an antenna in the device header. The data are encrypted and sent from the transmitter to the BIOTRONIK Service Center via the cellular phone network.
- The received data are deciphered and evaluated. Each physician can set the criteria for evaluation to be used for each patient and can configure the time of notification via e-mail, SMS or fax.
- A clear overview of the results of this analysis is displayed for the attending physicians on the protected Internet platform Home Monitoring Service Center (HMSC).
- Data transmission from the device is performed with a daily device message.
- Device messages which indicate special events in the heart or in the device are forwarded immediately.
- A test message can be initiated at any time using the programmer to immediately check the Home Monitoring function.

25. Method of product use

Therapy and diagnostic functions

Diagnostic functions

- Data from implantation and the most recent interrogations and follow-ups are recorded as well as arrhythmia episodes; they are stored together with other data to assess patients and the state of the device at any time;
- To check the lead for proper functioning, an automatic impedance measurement using subthreshold pacing pulses is performed in the device;
- Leadless ECG function: For all device types, far-field derivation can be measured without external leads between the right ventricular shock coil and housing, which, depending on the implantation site, corresponds to ECG derivation II or III (Einthoven);
- Once a telemetry connection has been established during a test procedure in an in-office follow-up, the leadless ECG and the IEGM are displayed with markers.

Antitachycardia pacing (ATP)

- The ICD can treat ventricular tachycardia with antitachycardia pacing (ATP); ATP can also be delivered in the VF zone (ATP One Shot) when the stability criterion indicating that this will be effective before shock delivery (monomorphic rapid VTs) is met;
- Depending on the device type, the device program contains not only the ICD functions but also all pacemaker functions for 1, 2, or 3 chambers. The heart rhythm is continuously monitored; each arrhythmia is classified according to the heart rate and the adjustable detection criteria. Depending on the preset values, antibradycardia as well as antitachycardia therapy is inhibited or delivered.

Cardioversion, defibrillation

- The ICD can treat ventricular tachyarrhythmia with cardioversion and/or defibrillation. Shock polarity and energy can be programmed individually. Shock energies between 2.0 and 40 J are possible. Before delivery of the shock, the ICD can be set to only deliver a shock when ongoing tachyarrhythmia is confirmed; during this time period the device can identify spontaneous conversion of the tachyarrhythmia and cancel the charging process if necessary;
- The shock paths can be set between the different shock coils (SVC/RV) and/or the housing.

Antibradycardia pacing and CRT

- Innovative rate hystereses, automatic sensor functions, and a night program promote the patient's intrinsic rhythm, avoid overdrive pacing, and facilitate adaptation of the device to the individual needs of the patient;
- Setting an upper tracking rate for the atrium prevents unspecific atrial pacing, thus reducing the risk of pacemaker-mediated tachycardia;
- Positive AV hysteresis functions support the intrinsic conduction and thus the natural contraction sequence. Negative AV hysteresis functions support the cardiac resynchronization therapy by maintaining pacing in stressful situations;
- For resynchronization of the ventricles, triple-chamber devices have functions for multisite ventricular pacing with possible VV delays in either direction;
- To ensure that no additional surgery is necessary in case of a left-sided increase of pacing threshold or undesired phrenic nerve stimulation, different pacing polarities can be set for the left ventricular lead with a triple-chamber device;
- Automatic active capture control is available for the right and left ventricle with automated tracking of the pacing threshold or automatic threshold monitoring (ATM) for trend analysis.

With the Lforia 3, the pacing threshold is only measured.

Program presets

The parameter settings can be saved in 3 individual therapy programs.

Home Monitoring functions

- The device automatically sends information to the transmitter once a day. It also sends messages related to events, which are immediately forwarded to the Service Center. In addition to this, test messages can be initiated using the programmer.
- Appointments for Home Monitoring-supported follow-ups can be scheduled via the HMSC (Home Monitoring Service Center). This applies to Lforia 5/7 only.
- Important medical information in the device messages include the following:
 - atrial and ventricular arrhythmias;
 - parameters relevant to leads in the atrium and ventricle: pacing thresholds, sensing amplitudes,

- impedances;
 - current statistics;
- IEGM-Online HD with up to three high definition (high resolution) channels.

26. The selection and treatment of the patients

Individualization of treatment

- Identify outweigh the expected benefits of the device, the possibility of early device replacement for patients whose ventricular tachyarrhythmia require frequent shocks;
- To determine whether the implant and programmable options are appropriate for patients with drug-resistant supraventricular tachyarrhythmia (SVT), as drug resistance SVT can trigger unwanted therapy implant.

Changes in anti-arrhythmic drugs or other medicinal drugs of the patient that affect the normal heart rhythm of the patient or the conductivity, can influence the rate of tachyarrhythmia and/or the effectiveness of therapy.

If there is a need in another cardiac surgery procedure performing before the implantation of ICD, perhaps preferable to implant an ICD during this procedure. This can prevent the need for an additional intervention.

Patient consulting

A patient guide is available to the patient, to his relatives and other interested people. It's obligatorily to discuss the content of book/guide with the patient and other interested people before and after ICD implantation to provide a complete understanding of the material.

Evaluation of potential patients for ICD

Estimated candidate for ICD implantation must pass the assessment of the heart in order to classify any and all tachyarrhythmias. In addition, other specific cardiac data of the patient can help in selecting the optimal device settings.

This assessment may include, but is not limited to:

- Estimated rate of the tachycardia;
- Confirmation and/or evaluation of any supraventricular arrhythmias and bradyarrhythmias;
- Evaluation of various ATP and pacing functions;
- The presence of any post-shock arrhythmias;
- Evaluation of the maximum sinus rate during exercise.

If you have changed the mode of drug treatment or adjusted, while the ICD is implanted, you may need additional EF testing to determine whether the parameter settings of the detection or treatment necessary and appropriate.

Empirical changes in the parameters of detection or treatment should be evaluated on the basis of patient safety.

27. Implantation

Parts preparation

The following parts that correspond to the requirements of the EC Directive 90/385/EEC are required:

- The implant with the plug cap and BIOTRONIK screwdriver;
- BIOTRONIK leads and the lead introducer set:
 - Single-chamber device: One bipolar ICD lead with 1 or 2 shock coils for the ventricle;
 - Dual-chamber device: One bipolar lead for the atrium and one bipolar ICD lead for the ventricle with 1 or 2 shock coils;
 - Triple-chamber device: one bipolar ICD lead with 1 or 2 shock coils for the ventricle, one bipolar pacing lead for the right atrium and one bipolar pacing lead for the left ventricle;
- DF-1, DF4 and IS-1 connections are approved. For leads with a different connection or leads from other manufacturers, use adapters approved by BIOTRONIK only;
- BIOTRONIK programmer and approved cable;
- External multi-channel ECG device.

Keep spare parts for all sterile components.

External Defibrillator

In order to be able to respond to unforeseeable emergencies or possible technical failures of the device:

- Keep an external defibrillator and paddles or patch electrodes ready.

Unpacking the device

WARNING

Inadequate therapy due to defective device

If an unpacked device is dropped on a hard surface during handling, electronic parts could be damaged.

- Use a replacement device.
- Return the damaged device to BIOTRONIK.

- Peel the sealing paper off of the outer blister at the marked position in the direction indicated by the arrow. The inner blister must not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments;
- Take hold of the inner blister by the gripping tab and take it out of the outer blister;
- Peel the sealing paper off of the sterile inner blister at the marked position in the direction indicated by the arrow.

Checking parts

Damage to any of the parts can result in complications or technical failures.

- Check for damage before and after unpacking all parts;
- Replace damaged parts;
- The ICD is shipped with tachyarrhythmia therapy deactivated and is only to be connected and implanted in this state;
- Leads must not be shortened.

Implantation site

Depending on lead configuration and the patient's anatomy, the ICD is generally implanted subpectorally on the left side.

Preventing leakage currents

Leakage currents between the tools and the device must be prevented during implantation.

- Electrically insulate the patient.

Preventing unintentional shock delivery

WARNING

Shock delivery with activated ICD

There is a risk of unintended shock delivery when handling an activated ICD.

- Deactivate ICD therapy before touching the device during implantation, device replacement and explantation.

Avoid damage to the header

The provided set screws must be carefully loosened or tightened.

- Loosen set screws with the supplied screwdriver. Use only BIOTRONIK screwdrivers with torque control!
- Do not forcibly pull out the blind plug!
- If lead revision is necessary, re-order sterile screwdrivers from BIOTRONIK.

Preventing short circuits in the header

WARNING

Short circuit due to open lead connector ports

Connector ports in the header which are open and thus not electrolyte-proof may cause undesired current flows to the body and penetration of body fluid into the device.

- Either leave unused ports closed with the premounted blind plugs, or close them using the supplied blind plugs.

Ensure that connections are clean

In case of contamination during implantation:

- Clean lead connectors with a sterile cloth;
- Rinse connection only with sterile water.

Connecting the lead connector to the device

1 Disconnect stylets and stylet guides.

2 **DF-1 connection:**

- Connect the DF-1 connector for the right-ventricular shock coil to RV.
- Connect the DF-1 connector for the supraventricular shock coil to SVC. Or connect a subcutaneous array to SVC.

DF4 connection:

- Connect the DF4 connector to RV

3 **IS-1 connection:**

- Connect the bipolar IS-1 lead connector for the atrium to RA.
- Connect the bipolar IS-1 lead connector for the DX to RA.
- Connect the IS-1 lead connector for the right ventricle to RV.
- Connect the unipolar or the bipolar IS-1 lead connector for the left ventricle to LV.

4 Push the lead connector into the header without twisting or bending the connector or conductor until the connector tip (on the DF-1 connector) or the insertion indicator (on the DF4 connector) becomes visible behind the set screw

block. This indicator can vary depending on the manufacturer of the lead used.

5 If you cannot easily plug the lead connector into the connection:

- Use only sterile water as lubricant.

6 If the lead connector cannot be inserted completely, the set screw may be protruding into the drill hole of the set screw block.

- Use the screwdriver to perpendicularly pierce through the slitted point in the center of the silicone plug until it reaches the set screw.

- Carefully loosen the set screw without completely unscrewing it, so that it does not become tilted upon retightening.

7 Turn the set screw clockwise until torque control starts (you will hear a clicking sound).

8 Carefully withdraw the screwdriver without retracting the set screw.

- In case of IS-1 connections with two set screws, tighten both screws!

- When the screwdriver is withdrawn, the silicone plug automatically safely seals the lead connector port.

Keeping distance between leads

WARNING

Inadequate therapy

When leads are not spaced sufficiently apart or are positioned inappropriately, this can lead to far-field sensing or insufficient defibrillation.

- The distance between 2 shock coils must be greater than 6 cm;
- Tip and ring electrodes must not have contact with each other.

Implanting

1	Prepare the vein.
2	Implant the leads, perform the measurements, and fix
3	Make the device pocket.
4	Connect the lead connector to the device.
5	Insert the device.
6	Guide the fixation suture through the opening in the header and fixate the device in the prepared device pocket.
7	Close the device pocket.
8	Check the device with standard tests.

28. Applying the programming head

The programming head (PGH) features a diagram of the device. This is used to assist in positioning the head to ensure proper telemetry.

- Make sure the PGH is positioned correctly.

29. Establishing telemetry contact

The programmer (or the SafeSync Module) can be no more than 3 m from the device; ideally there should be no hindrances between the patient and the programmer.

- Switch on RF telemetry on the programmer.
- Apply the programming head for about 2 s until successful initialization is displayed on the programmer:
- The SafeSync symbol is displayed in the navigator and the signal strength is displayed in the status line.
- Remove the programming head.

30. Activating ICD therapy

- Load the device program that is suitable for the device type in the programmer;
- Activate ICD therapy;
- Shipment mode is permanently deactivated once the leads have been connected and initial measurement of the pacing impedance has been performed. The device data are saved;
- Take precautionary measures while programming;
- If the device induces tachycardia while programming ATPs or does not deliver adequate therapy in the DFT test: use emergency shock or an external defibrillator.

31. Precautionary Measures while Programming

Performing standard tests and monitoring the patient

Critical conditions can occur for the patient even during standard tests due to inadequate parameter settings or interrupted telemetry.

- Ensure sufficient patient care even during tests.
- After the threshold test, check to determine whether the threshold is clinically and technically justifiable.
- Continuously monitor the ECG and the patient's condition.
- Cancel testing if necessary.

Cancelling telemetry

Programmer interference or interrupted telemetry during performance of temporary programs (follow-up tests) can result in inadequate pacing of the patient. This is the case if the programmer can no longer be operated due to a program error or a defective touch screen and therefore the temporary program cannot be terminated. Under these circumstances, it is helpful to cancel telemetry, in which case the device automatically switches to the permanent program.

- In the case of telemetry with programming head: lift the PGH by at least 30 cm;
- In the case of RF telemetry: switch off and reposition the programmer;
- Turn off possible sources of interference.

Avoiding critical parameter settings

No modes and parameter combinations that pose a risk to the patient should be set.

- Prior to setting rate adaptation, determine the patient's capacity for strain.
- Check compatibility and effectiveness of parameter combinations after making settings.

Check for leads suitable for shock path

Three shock paths can be set, two of which form an electrical path to the device housing.

- A second shock coil (dual shock coil) must be available for the shock path RV -> SVC.

Monitoring the patient when setting asynchronous modes

The asynchronous modes V00 and D00 can only be set if tachyarrhythmia sensing is deactivated. This would leave the patient without sensing and therefore without ICD therapy.

- Continually monitor the patient.
- Keep an external defibrillator ready.

Setting sensing

Manually set parameters can be unsafe. For example, unsuitable far-field protection may impede sensing of intrinsic pulses.

- Note automatic sensitivity control.

Preventing device-induced complications

BIOTRONIK devices feature several functions to prevent device-induced complications to the greatest extent possible:

- Measure the retrograde conduction time.
- Set PMT protection.
- Set the VA criterion.

Preventing conduction of atrial tachycardia

BIOTRONIK devices feature several functions to prevent conduction of atrial tachycardia to the ventricle(s):

- Set mode switching for indicated patients;
- Set the upper rate and the refractory periods to prevent abrupt ventricular rate switching;
- Give preference to Wenckebach response and avoid 2:1 behavior;
- Set all parameters so as to prevent constant changing between atrial and ventricular-controlled modes.

Observing the shock impedance limit

The implanted device could be damaged if the shock impedance is too low.

- The shock impedance must be $> 25 \Omega$.

Preventing recurrence after therapy shock

After a therapy shock, pacing can be performed with a post-shock program if there is no intrinsic rhythm.

Permanent program	Post-shock program
DDD, DDI, AAI	DDI
VDD, VDI	VDI
VVI and OFF	VVI

- The following post-shock program parameters can be adjusted: post-shock duration, basic rate, rate hysteresis, ventricular pacing, AV delay (fixed, not dynamic), LV-T-wave protection;
 - The default settings for the post-shock program are as follows:
A and RV: 7.5 V, 1.5 ms,
LV – the same value as for permanent program

Considering power consumption and service time

RF telemetry requires somewhat more power: Consumption during implantation corresponds to approximately 10 days of service time and consumption during a 20-minute follow-up corresponds to approximately 3 days.

- Do not establish unnecessary RF telemetry;
- After 5 minutes without input, SafeSync switches to the economy mode;
- Check the battery capacity of the device at regular intervals.

Application of a permanent magnet

Applying a permanent magnet interrupts detection and therapy of tachycardia events. After 8 hours of this type of deactivation, the device automatically reactivates the therapy functions to prevent accidental permanent deactivation.

- If detection interruptions of longer than 8 hours are required, the magnet has to be briefly removed from the device. The 8 hour countdown restarts when the magnet is applied again;
- Use BIOTRONIK magnets: type M-50 permanent magnets.

32. Follow-up

Follow-up intervals

Follow-ups must be performed at regular, agreed intervals.

- The first follow-up should be carried out by the physician using the programmer (in-office follow-up) approximately 3 months after implantation following the lead ingrowth phase;
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

Follow-up with BIOTRONIK Home Monitoring®

Monitoring using the Home Monitoring function does not serve to replace regular inoffice appointments with the physician required for other medical reasons. Follow-up supported by Home Monitoring can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted despite use of the Home Monitoring function if symptoms worsen or if new symptoms arise;
- Device messages are transmitted regularly;
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device system are sufficient. If not, an in-office follow-up has to be carried out.

Possible early detection due to information gained via Home Monitoring may necessitate an additional in-office follow-up. For example, the data may indicate at an early stage lead problems or a foreseeable end of service time (ERI). Furthermore, the data could provide indications of previously unrecognized arrhythmias or modification of the therapy by reprogramming the device.

Follow-up with the programmer

Use the following procedure for in-office follow-up:

1	Record and evaluate the ECG.
2	Interrogate the device.
3	Evaluate the status and automatically measured follow-up data.
4	Check the sensing and pacing functions.
5	Possibly evaluate statistics and IEGM recording.
6	Manually perform standard tests if necessary.
7	Possibly customize program functions and parameters.
8	Transmit the program permanently to the device.
9	Print and document follow-up data (print report).
10	Finish the follow-up for this patient.

33. Patient consulting

Patient ID card

A patient ID card is included in delivery packaging;

- Provide the patient with the patient ID.
- Advise patients to consult a doctor if they have any questions.

Prohibitory sign

The patient should avoid areas marked with the prohibitory sign.

- Please attract patient's attention to prohibitory signs.

Possible sources of interference

Electromagnetic interference should be avoided in daily activities. Sources of interference should not be brought into close proximity with the device.

- Draw the patient's attention to special household appliances, security checkpoints, anti-theft alarm systems, strong electromagnetic fields, cell phones, and transmitters among other things.
- Request patients to do the following:
 - Use cell phones on the side of their body that is opposite of the device.
 - Keep the cell phone at least 15 cm away from the device both during use and when stowing.

34. Replacement Indications

Possible battery levels

BOS	Beginning of Service: > 70% charge
MOS 1	Middle of Service: 70% to 40% residual charge
MOS 2	Middle of Service: < 40% residual charge
ERI	Elective Replacement Indication, (i.e. RRT: Recommended Replacement Time)
EOS	End of Service

Elective Replacement Indicator

Elective Replacement Indication can be detected by Home Monitoring.

CAUTION

Temporary limited therapy

If ERI occurs shortly after follow-up and is only detected during the subsequent follow-up, then the remaining service time can be much less than 3 months.

- Replace device soon.
- The device can monitor the heart rhythm for at least 3 more months;
- At least 6 maximum energy shocks can be delivered until EOS occurs;
- The selected parameters in the device program do not change.

EOS replacement indication

End of Service can be detected by Home Monitoring.

WARNING

Patient at risk of death

If EOS replacement indication occurs before replacement of the device, then the patient is without therapy.

- Replace device immediately;
- Monitor patient constantly until immediate replacement of the device!
- VT and VF detection and all therapies are deactivated!
- The antibradycardia function remains active in the VVI mode:
 - Ventricular pacing: RV; basic rate 50 bpm; without special pacemaker functions such as hysteresis, etc;
 - Pulse amplitude of 6 V; pulse width of 1.5 ms;
 - Time of transmission for Home Monitoring: 90 days.

35. Explantation and Device Replacement

Explantation

- Interrogate the device status;
- Deactivate VT and VF therapies prior to explantation;
- Remove the leads from the header. Do not simply cut them loose;
- Use state-of-the-art techniques to remove the device and, if necessary, the leads.

Note: Normal oxidation processes may cause ICD housing discolorations. This is neither a device defect nor does it influence device functionality.

- Explants are biologically contaminated and must be disposed of safely due to risk of infection.

Device replacement

If, upon replacing the device, already implanted leads are no longer used but left in the patient, then an additional uncontrolled current path to the heart can result.

- Deactivate VT and VF therapies prior to device replacement;
- Insulate connections that are not used.
- Basic principles:
- The device must not be resterilized and reused.

Cremation

Devices should not be cremated.

- Explant the device before the cremation of a deceased patient.


36. Technical Specifications

Box specifications

Model	Net weight, g	Gross weight, g
Iforia 3 VR-T	117,2	266,8
Iforia 5 VR-T	116,8	266,4
Iforia 7 VR-T	114,1	263,7
Iforia 5 VR-T DX	117,4	267
Iforia 7 VR-T DX	115	264,6
Iforia 3 DR-T	117,2	266,8
Iforia 5 DR-T	116,7	266,3
Iforia 7 DR-T	114,5	264,1
Iforia 3 HF-T	119,3	268,9
Iforia 5 HF-T	119,1	268,7
Iforia 7 HF-T	117,3	266,9

The measured values are nominal $\pm 10\%$

Box photo




Box sizes for all Iforia ICD models: 120 x 160 x 54 mm.

Inner packaging photo

37. Labeling description


The label icons symbolize the following:

	Manufacturing date		Use by
	Temperature limit		Order number
	Serial number		Product identification number
	Dangerous voltages!		CE mark
	Contents		Follow the instructions for use

  Sterilized with ethylene oxide			
	Do not resterilize		Single use only. Do not re-use!
	Do not use if packaging is damaged		Non-sterile

	Transmitter with non-ionizing radiation at designated frequency
Label icon on devices with ProMRI: 	MR conditional: patients having a device system implanted whose components are labeled with this symbol on the packaging can be examined using an MR scan under precisely defined conditions

 Example	Device: NBG code and compatible leads
 Example	Factory settings for therapy: OFF
	Screwdriver

 Examples of DF-1/IS-1 or DF4/IS-1 header	
	Bipolar IS-1 connector
	Unipolar IS-1 connector
	Unipolar DF-1 connector
	DF4 connector

Labelling also contains:

- ICD model;
- Name and address of the manufacturer

During the realization on the territory of the Russian Federation, medical device will be marked by the local sticker.

Local sticker will contain the following additional information:

- ☐ Name of device (in Russian)
- ☐ Model
- ☐ Type of working part - CF
- ☐ Manufacturer (name)
- ☐ Country of origin
- ☐ Reference to Instruction for use
- ☐ Registration certificate № _____ date of issue: _____

38. Parameters

Bradycardia / CRT

General ICD therapy

Parameter	Range Of Values	Standard	VR	DX	DR	HF
ICD therapy	OFF; ON	ON	x	x	x	x
Programs	Display standard program; Display safe program; Display first interrogated program; Individual 1,2,3	—	x	x	x	x

Timing: Basic rate day/night and rate hysteresis

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Basic rate	30 ... (5) ... 100 ... (10) ... 160 bpm	40 bpm	x	x		
		60 bpm			x	x
Night rate	OFF; 30 ... (5) ... 100 bpm	OFF	x	x	x	x
Night begins	00:00 ... (00:01) ... 23:59 hh:mm	22:00 hh:mm	x	x	x	x
Night ends		06:00 hh:mm				
Rate hysteresis	OFF; -5 ... (-5) ... -25 ... (-20) ... -65 bpm	OFF	x	x	x	x

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Parameter	Range Of Values	Standard	VR	DX	DR	HF
AV dynamics	Low; Medium; High; Fixed; (Individual)	Low		x	x	x
AV delay (1 or 2) after:						
– pacing	15; 40 ... (5) ... 350 ms	—			x	x
– sensing	Either automatic: AV delay after pacing + sense compensation Or: 40 ... (5) ... 350 ms	—		x	x	x
– At rate 1	50 ... (10) ... 130 bpm	60 bpm		x	x	x
– At rate 2	60 ... (10) ... 140 bpm	130 bpm		x	x	x

Sense compensation	OFF; -5 ... (-5) ... -120 ms	-40 ms			x	x
AV hysteresis mode	OFF; Positive; Negative; IRSplus	OFF		x	x	
	OFF; Positive; Negative	OFF				x
AV hysteresis (positive)	70; 110; 150; 200 ms	70 ms		x	x	x
AV hysteresis (negative)	10 ... (10) ... 150 ms	50 ms		x	x	x
AV scan and repetitive (positive)	OFF; ON	ON		x	x	x

Timing: Post-shock pacing

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Post-shock duration	OFF; 10 s; 30 s; 1 min; 2 min; 5 min; 10 min	10 s	x	x	x	x
Post-shock basic rate	30 ... (5) ... 100 ... (10) ... 160 bpm	60 bpm	x	x	x	x
Post-shock AV-delay	50 ... (10) ... 350 ms	140 ms			x	x
Ventricular post-shock pacing	RV; BiV	RV				x

Timing: Upper rate

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Upper rate	90 ... (10) ... 160 bpm	130 bpm		x	x	x
Atrial upper rate	OFF; 175; 200; 240 bpm	200 bpm			x	x

Timing: Mode switching

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Intervention rate	OFF; 120 ... (10) ... 200 bpm	160 bpm		x	x	x
Onset criterion	3 ... (1) ... 8 (из 8)	5		x	x	x
Resolution criterion						
Modification of basic rate	OFF; 5 ... (5) ... 30 bpm	10 bpm		x	x	x
Mode	VDI(R); VDD(R)	VDI		x	x	x
	DDI(R); DDD(R)	DDI			x	x
After mode switching:						
– Rate	OFF; 5 ... (5) ... 50 bpm	10 bpm		x	x	x
– Duration	1 ... (1) ... 30 min	1 min		x	x	x

Timing: Ventricular Stimulation

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Permanent	RV; BiV; LV	BiV				x
Triggering	OFF; RVs; RVs+PVC	RV				x
LV-T wave protection	OFF; ON	ON				x
Maximum trigger rate:						
– DDD(R) and VDD(R)	UTR + 20; 90 ... (10) ... 160 bpm	UTR + 20				x
– DDI(R), VDI(R) and VVI(R)	90 ... (10) ... 160 bpm	130 bpm				x
Initially paced chamber	RV; LV	LV				x
VV-delay after Vp	0 ... (5) ... 100 ms	5 ms				x

Timing: Refractory periods and blanking periods

Parameter	Range Of Values	Standard	VR	DX	DR	HF
PVARP	AUTO; 175 ... (25) ... 600 ms	225 ms		x	x	x
Blanking after atrial pacing	50 ... (10) ... 100 ms	50 ms			x	x
LV blanking after RV pacing		80 ms				x
RV blanking after LV pacing						
Far-field protection after Vs	AUTO; 25 ... (25) ... 225 ms	75 ms		x	x	x

Far-field protection after Vp	50 ... (25) ... 225 ms	75 ms		x	x	x
-------------------------------	------------------------	-------	--	---	---	---

Timing: PMT protection

Parameter	Range Of Values	Standard	VR	DX	DR	HF
PMT detection/termination	OFF; ON	ON		x	x	x
VA criterion	250 ... (10) ... 500 ms	350 ms		x	x	x

Timing: Rate adaptation via accelerometer

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Maximum sensor rate	80 ... (10) ... 160 bpm	120 bpm	x	x	x	x
Sensor gain	AUTO; Very low; Low; Medium; High; Very high	Medium	x	x	x	x
Sensor threshold	Very low; Low; Medium; High; Very high	Medium	x	x	x	x
Rate increase	1; 2; 4; 8 bpm/cycle	2 bpm/cycle	x	x	x	x
Rate decrease	0.1; 0.2; 0.5; 1.0 bpm/cycle	0,5 bpm/cycle	x	x	x	x

Pacing: Pulse amplitude and pulse width

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Pulse amplitude A	0,5 ... (0,25) ... 4,0 ... (0,5) ... 6,0; 7,5 V	2,5			x	x
Pulse amplitude V/RV			x	x	x	x
Pulse amplitude LV						x
Pulse width A	0,4; 0,5 ... (0,25) ... 1,5 ms	0,4 ms			x	x
Pulse width V/RV			x	x	x	x
Pulse width LV		0,4 ms				x

Pacing: Ventricular capture control

(The following parameters apply to Iforia 5/7)

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Ventricular capture control	OFF; ATM; ON	ATM	x	x	x	x
Threshold test start	2,5 ... (0,5) ... 5,0 V	ATM: 2,5 V ON: 3,5 V	x	x	x	x
Minimum amplitude	1,0 ... (0,25) ... 4,0 V	1,0 V	x	x	x	x
Safety margin	1,0; 1,2 V	1,0 V	x	x	x	x

(The following parameters apply to Iforia 3)

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Capture control	OFF; ATM	ATM	x	x	x	x

Pacing: Atrial capture control

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Atrial capture control	OFF; ATM	ATM			x	x

LV lead configuration

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Pacing polarity LV	LV tip -> LV ring; LV tip -> RV ring; LV ring -> LV tip; LV ring -> RV ring; UNIP	LV tip -> RV ring				x
Sensing polarity LV	UNIP; BIPL	UNIP				x

MRI program

(Valid for MR-compatible devices with ProMRI²)

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Mode	V00; OFF	OFF	x	x		
	V00; D00; OFF				x	x
Basic Rate	70... (10)... 160 bpm	90 bpm			x	x
Ventricular Stimulation	RV, BiV	RV				x

Tachycardia

Detection

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Interval AT/AF	240...600 ms	300 ms		x	x	x
Interval VT1	OFF; 270 ... (10) ... 600 ms	OFF	x	x	x	x
Interval VT2	OFF; 270 ... (10) ... 500 ms	OFF	x	x	x	x
Interval VF	OFF; 240 ... (10) ... 400 ms	300 ms	x	x	x	x
Detection counter VT1	10 ... (2) ... 60	26	x	x	x	x
Detection counter VT2	10 ... (2) ... 40	16	x	x	x	x
Detection counter VF	6 out-of 8; 8 out-of 12; 10 out-of 14; 12 out-of 16; 6 out-of 20; 18 out-of 24; 20 out-of 26; 22 out-of 30; 24 out-of 30	8 out of 12	x	x	x	x
Redetection counter VT1	10 ... (2) ... 30	20	x	x	x	x
Redetection counter VT2	10 ... (2) ... 30	14	x	x	x	x
SMART detection VT1/VT2	OFF; ON	ON		x	x	x
SMART detection ON:						
– Onset VT1/VT2	4 ... (4) ... 32 %	20 %		x	x	x
– Stability VT1/VT2	8 ... (4) ... 48 %	12 %		x	x	x
SMART detection OFF:						
— Onset VT1/VT2	OFF, 4 ... (4) ... 32 %	20 %	x	x	x	x
– Stability VT1/VT2	OFF; 8 ... (4) ... 48 ms	24 ms	x	x	x	x
Sustained VT	OFF; 1; 2; 3; 5; 10; 20; 30 min	OFF	x	x	x	x
Forced termination	OFF; 1 ... (1) ... 10 min	1 min		x	x	x

² For information concerning MRI compatibility of different ICD variants please see Chapter „MRI-compatibility“ on page 14 of this document

Therapy: ATP

Parameter	Range Of Values	Standard	VR	DX	DR	HF
ATP type for VT1/VT2	Burst; Ramp	OFF	x	x	x	x
ATP type for VF	OFF; Burst; Ramp	Burst	x	x	x	x
ATP optimization	OFF; ON	OFF	x	x	x	x
Attempts	OFF; 1 ... (1) ... 10	OFF	x	x	x	x
Number S1 for VT1/VT2	1 ... (1) ... 10	5	x	x	x	x
Number S1 for VF		8	x	x	x	x
S1 decrement for VT1/VT2 and for VF	5 ... (5) ... 40 ms	10 ms	x	x	x	x
Scan decrement	OFF; 5 ... (5) ... 40 ms	OFF	x	x	x	x
Additional S1 for VT1/VT2	OFF; ON	ON	x	x	x	x
Ventricular pacing for VT1/VT2	RV, LV, BiV	RV				x
Ventricular pacing for VF		RV				
R-S1 interval for VT1/VT2	70 ... (5) ... 95 %	80 %	x	x	x	x
R-S1 interval for VF		85 %	x	x	x	x

Therapy: Shock

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Number of shocks VT1/VT2	0; 1; 2; 6; 8	8	x	x	x	x
Number of shocks VF	6; 8	8	x	x	x	x
1st Shock for VT1/VT2	OFF; 2 ... (2) ... 20 ... (5)... 40 J	40 J	x	x	x	x
2nd Shock for VT1/VT2	OFF; 4 ... (2) ... 20 ... (5)... 40 J	40 J	x	x	x	x
3rd - nth shock for VT1/VT2	OFF; 4*40 J; 6*40 J	6*40 J	x	x	x	x
1st Shock for VF	2 ... (2) ... 20 ... (5) ... 40 J	40 J	x	x	x	x
2nd Shock for VF	4 ... (2) ... 20 ... (5) ... 40 J	40 J	x	x	x	x
3rd - nth Shock for VF	4*40 J; 6*40 J	6*40 J	x	x	x	x

For shock in VT1/VT2 and VF:						
Confirmation	OFF; ON	ON	x	x	x	x
Polarity	Normal; Reverse; Alternating	Normal	x	x	x	x
Shock form	Biphasic; Biphasic 2	Biphasic	x	x	x	x
Shock path	RV -> ICD+SVC	RV-> ICD+SVC	x		x	x
	RV -> ICD	RV -> ICD				
	RV -> SVC	RV -> SVC		x		

Sensing

Sensitivity and thresholds

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Sensing A	STD; OFF	STD		x	x	x
Sensing RV	STD; TWS; VFS; IND	STD	x	x	x	x
Sensing LV	STD; OFF; IND	STD				x
Upper threshold RV	50; 75 %	50 %	x	x	x	x
Upper threshold LV	50; 75 %	50 %				x
Upper threshold duration after detection	110; 150 ... (50) ... 500 ms VFS: 110 ms	350 ms	x	x	x	x
Upper threshold duration after pacing		400 ms	x	x	x	x
Lower threshold RV	25; 50 %	25 %	x	x	x	x
T-wave suppression after pacing	OFF; ON	OFF	x	x	x	x
Minimum threshold A	0,2 ... (0,1) ... 2,0 mV	0,4 mV		x	x	x
Minimum threshold RV	0,5 ... (0,1) ... 2,5 mV	0,8 mV	x	x	x	x
Minimum threshold LV	0,5 ... (0,1) ... 2,5 ... (0,5) ... 5,0 mV	1,6 mV				x

Diagnostics

The following can be set:

Parameter	Range Of Values	Standard	VR	DX	DR	HF
For AT/AF	OFF; ON For Iforia 7: Extended ON	ON		x	x	x
For SVT	OFF; ON	ON		x	x	x
Periodic recording	When Home Monitoring OFF: OFF; 30 ... (30) ... 180 days	90 days	x	x	x	x
IEGM configuration	RA, RV, LV RA, RV, FF FF; RV; LV	RA, RV, LV				x

The following can additionally be set for Iforia 7:

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Start resting period	0:00 ... (1:00) ... 23:00 hh:mm	02:00 hh:mm	x	x	x	x
Duration of resting period	0,5 ... (0,5) ... 12 h	4 h	x	x	x	x
AV delay modification in sensing test	OFF; 300 ms	300 ms		x	x	x
Thoracic impedance (TI)	OFF; ON	OFF	x	x	x	x

Home Monitoring

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Home Monitoring	OFF; ON	OFF	x	x	x	x
Time of transmission	STD; 00:00 ... (01:00) ... 23:00 hh:mm	STD	x	x	x	x
IEGM for:						
– Therapy episodes	OFF; ON	ON	x	x	x	x
– Monitoring episodes						
Ongoing atrial episode	OFF; 6; 12; 18 h	12 h		x	x	x

Only applies to Iforia 3:

Parameter	Range Of Values	Standard	VR	DX	DR	HF
Cycle duration	OFF; 30; 60; 90; 120; 180 days	30 days	x	x	x	x
Transmission date:	XX.XX.XXXX	Follow-up + 7 days	x	x	x	x

39. Mechanical Characteristics

Housing

Devices with a DF-1/IS-1 header:

Type	Connector port	W x H x D in mm	Volume cm ³	Mass g
VR	DF-1	65 x 55 x 11	33	80
DX, DR	DF-1	65 x 55 x 11	33	81
HF	DF-1	65 x 58.5 x 11	34	83

Devices with a DF4/IS-1 header:

Type	Connector port	W x H x D in mm	Volume cm ³	Mass g
VR	DF4	65 x 52 x 11	31	80
DR	DF4	65 x 56 x 11	32	82
HF	DF4	65 x 56 x 11	33	82

Iforia 3/5/7 materials that are in permanent contact with internal environment and body tissue:

- Housing: Titanium;
- Header: Epoxy resin;
- Blind plug and silicone plug: Silopren or silastic;
- DF4 seal: Silastic

Screwdriver materials that have short-term mediated contact with the human body:


- Stainless chromium steel
- Polyoxymethylene

Inner Packaging materials:

- PETG
- High density polyethylene

40. Technical Characteristics

Pulse waveform

The pacing pulse has the following form: the pulse amplitude reaches its maximum value at the beginning of the pulse (U_a). With increasing pacing duration (t_b), the pulse amplitude is reduced dependent on the pacing impedance.

X-ray marker

«NT»

Resistance to interference

- Note on device type DX (only devices with a DF-1/IS-1 connection): The EMC requirements are met as long as atrial sensitivity is set to 1.0 mV (factory settings) or values ≥ 1.0 mV. Measures must be taken to assure interference-free therapy if more sensitive values are set.
- Note on device type HF: In the case of unipolar sensing, the requirement for interference voltages of ≤ 0.3 mV (peak to peak) is met.

Common mode rejection ratio

Rate	Common mode rejection ratio			
	Atrium: DX*	Atrium: DR, HF	V right: VR, DR, HF	V left: HF
16.6 Hz	58 dB	53 dB	64 dB	66 dB
50 Hz	55 dB	55 dB	64 dB	66 dB
60 Hz	56 dB	56 dB	64 dB	68 dB
*Only devices with a DF-1/IS-1 connection.				

ATP Amplitude

A burst was measured at 500 Ω , an amplitude of 7.5 V (tolerance ± 1.5 V), pulse width of 1.5 ms, R-S1 interval of 300 ms and a number S1 of 5:

ATP Amplitude	Measured minimum	Measured maximum	Mean value
RV	7,67 V	7,67 V	5,00 V
LV	7,67 V	7,67 V	4,99 V

Automatic sensitivity setting

Measurement of actual values and test signal wave shape: standard triangle. For the device type DX, the programmed atrial sensitivity is intensified by a factor of 4.

Sensitivity	Value	Tolerance	Measured value
A: positive	0.2 mV	0.2 ... 0.5	0.24 mV
A: negative			0.24 mV
DX: A: positive	0.2 mV	0.2 ... 0.52 [0.05 to 0.13]	0.05 mV
DX: A: negative			0.05 mV
RV: positive	0.5 mV	0.3 ... 0.7	0.48 mV
RV: negative			0.40 mV
LV: positive	0.5 mV	0.3 ... 0.7	0.48 mV
LV: negative			0.56 mV

Shock energy / peak voltage

With shock path: RV to housing + SVC

Shock energy (Tolerance)	Tolerance peak voltage	Measured value Shock energy	Measured value Peak voltage
1 J (0.7 ... 1.18)	90 ... 120 V	0.84 J	100 V
20 J (16.9 ... 20.9)	440 ... 480 V	18.1 J	469 V
40 J (33.8 ... 41.4)	620 ... 690 V	36.9 J	667 V

Guidance and manufacturer's declaration – electromagnetic emissions

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes

Guidance and manufacturer's declaration – electromagnetic immunity

Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic environment – guidance
Electrostatic discharge (ESD)	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Radiated RF	3 V/m 80 MHz to 2,5 GHz	10 V/m	$d = 1,2\sqrt{P}$ (80 MHz to 800 MHz) $d = 2,3\sqrt{P}$ (800 MHz to 2,5 GHz) where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a) should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE: U_T is the a.c. mains voltage prior to application of the test level. NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

41. Battery data

Battery characteristics

The following data is provided by the manufacturers:

Manufacturer	GREATBATCH, INC. Clarence, NY 14031	LITRONIK GmbH & Co 01796 Pirna, Germany
Battery type	GB 2992	LiS 3410 RA/RR
Battery ID number shown on the programmer	3	4
Device type	VR, (DX), DR, HF	
Battery voltage at ERI	2.5 V	2.85 V
Charge time at BOS	8 s	8 s
Charge time at ERI	10 s	10 s
Usable capacity until ERI	Iforia 3/5: 1390 mAh Iforia 7: 1600 mAh	1390 mAh
Usable capacity until EOS	1730 mAh	1520 mAh

42. Compliance to Recognized Standards

The 'Iforia 3/5/7' ICDs are in compliance with the following international standards:

1.	EN 45502-1:1998	Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer
2.	EN 45502-2-2:2008	Particular requirements for active implantable medical devices intended to treat tachyarrhythmia
3.	EN 980: 2008 A1: 1999	Graphical symbols for use in the labelling of medical devices
4.	ISO 5841-3:2000 EN 50077:1994	Cardiac pacemakers Low profile connectors for implantable pacemakers
5.	ISO 11318:2002	Cardiac defibrillators – Connector assembly for implantable defibrillators – Dimensional and test requirements
6.	ISO 27186:2010	Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices - Dimensional and test requirements
7.	EN ISO 14971:2012	Medical products – Risk management – Risk analysis
8.	IEC 60601-1:2005: Chapter 14	PEMS: Programmable electrical medical systems

9.	ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
10.	EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices
11.	ISO 11135-1:2007	Sterilization of health care products – Ethylene oxide Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
12.	EN 556:2002	Sterilization of medical devices – Requirements for terminally sterilized medical devices to be labelled „Sterile”
13.	EN 62304: 2006	Medical device software -- Software life cycle processes

43. Biocompatibility

ICD 'Iforia 3/5/7' biocompatibility tests

Test details	Compliance to ISO standard
Biological evaluation of medical devices	EN ISO 10993-1:2009/AC:2010
Tests for cytotoxicity	EN ISO 10993-5:2009
	EN ISO 10993-12:2009
The residual amount of ethylene oxide	EN ISO 10993-7:2008
Solid particles	EN 45502-2-2:2008

ICD 'Iforia 3/5/7' biosecurity ICD 'Iforia 3/5/7' is confirmed by the positive results of tests for biocompatibility of materials used in the implant design.

Materials of animal origin

There are no materials of animal or human origin in 'Iforia 3/5/7'

Drug Substances

ICDs 'Iforia 3/5/7' don't contain any drug substances that incorporate into the human blood. ICDs 'Iforia 3/5/7' don't contain other substances that can be considered as a drug in a separate status.

44. General Safety instructions

Transportation conditions

Care during shipping and storage

- Devices must not be stored or transported close to magnets or sources of electromagnetic interference;
- Note the effects of the storage period; see Battery Data.

Temperature	From 5 °C to 45 °C
Humidity	From 30 to 75%
Pressure	From 375 to 795 mm Hg

Delivery in shipment mode

The device is delivered in shipment mode to protect the battery; capacitor reforming required during storage could result in controlled extended charge times of the shock capacitors.

- The shipment mode is displayed on the programmer after loading the device program (it is deactivated during implantation on initial measurement of the pacing impedance).

Sterile Shipping

The device and the screwdriver are delivered after their sterilization with ethylene oxide in the factory. Sterilization is guaranteed if blister packaging and quality control print are not damaged.

Sterile boxing

The device and screwdriver are packaged in two separately sealed blisters. The inner blister is also sterile on the outside so that it can be transferred in a sterile state during implantation;

Single use only

The device and screwdriver are intended for single use only.

- Do not use the device if the package is damaged;
- All ICD family 'Iforia 3/5/7' are intended for use only in one patient. DO NOT REUSE OR EXPOSE FOR SECONDARY STERILIZATION;
- Before use, carefully check 'Iforia 3/5/7' ICD components and packaging for damage during transportation. DO NOT USE in case of any components damage;
- Do not use after UBD expiring.
- The device must not be resterilized and reused .

45. Maintenance and repair requirements

The schedule of planned on hospital device follow-ups

After performing the ICD implantation procedure for individual device programming it's obligatorily to use the programmer to interrogate the device and set the permanent program after this.

The programmer can not only set the permanent program for the device, but also control the patient rhythm in time of performing the follow-ups.

The schedule of the follow-ups is planned by the physician, according to the point 25.1. The patient must adhere to this schedule.

All the follow-ups include the evaluation of the device functionality, viewing the accumulated statistics and the patient rhythm data, performing the leads tests.

The product cannot be repaired.

46. Disposal and special precautions

- Prior to explantation/replacement of the device: interrogate the device status and deactivate VT and VF therapies. It will exclude the risks for unexpected shock therapy;
- All the devices must be explanted before the cremation of a deceased patient;
- Implantable devices of BIOTRONIK SE & Co. KG are intended to single-use only;
- Implantable devices manufactured by BIOTRONIK CE & Co. KG are intended for single use only;

Explanted devices are the subject to mandatory return to the regional representative of the BIOTRONIK CE & Co. KG company for performing the analysis and disposal.

Disposal is made in accordance with the legislation effective in the territory of the Russian

Federation. Waste class: B

47. Environmental protection requirements in time of Medical Device use

Explanted ICDs are biologically contaminated and must be disposed of safely due to risk of infection .
On the expiration date or after explantation the device should be returned to the manufacturer for environmentally safe disposal.

48. BIOTRONIK SE & Co. KG Warranty conditions

The Cardiac Implantable Electronic Devices of BIOTRONIK SE & Co. KG (hereinafter "BIOTRONIK") have been developed, manufactured, and tested in compliance with the highest standards based on current technology.

However, premature battery failure, unpredictable random component failure, and other similar events may cause a Cardiac Implantable Electronic Device (hereinafter "Device") to function improperly or completely fail to function, even before the warranty period has expired. In such a case, BIOTRONIK provides the following limited warranty without limiting applicable statutory rights:

Warranty period

The warranty is valid for the periods indicated on the warranty chart and applies to the indicated Device type. The periods start with the initial implantation date of the Device.

Scope of warranty

If the implanted Device fails to function during the warranty period as indicated on the warranty chart – either due to component failure or battery depletion, despite the Expected Operating Specifications having been met – the original purchaser is entitled to receive, at BIOTRONIK's discretion, either: a replacement Device of an identical or equivalent model type from BIOTRONIK's current line of products or a replacement credit for the value of an identical or equivalent type of Device at the time the warranty was granted.

Warranty conditions

Claims listed under item 2 are valid only under this limited warranty for the initial implantation of the Device and provided that the following conditions have been met:

- a. The Device was implanted not later than the "use by" date shown on the Device package.
- b. The implantation report was sent back to BIOTRONIK or its distributor within 30 days after implantation.
- c. The Device was returned to BIOTRONIK or its distributor within 30 days after explantation, along with a written report stating the reasons for explantation. With this return, the Device shall become BIOTRONIK's property.
- d. The Device was not tampered with in any way, either by the purchaser or any third party.
- e. The failure of the explanted Device was confirmed by BIOTRONIK.

Range of warranty

The warranty does not cover the following items:

- a. The costs incurred by the purchaser of the Device or by any third party in connection with its replacement or any other claims that might have arisen in connection with a failure of the implanted Device.
- b. The costs of explanting the Device.
- c. The implantation of the replacement Device.
- d. Any expenses incurred by the original purchaser or by any third party by reason of such a failure – in

- particular the costs for surgical, medical, or incidental expenses.
- e. Any damage resulting from the use, malfunction, or failure of the Device.

There is no warranty, expressed or implied, for leads, adapters, or any other accessories used in conjunction with the Device. Furthermore, any malfunctions caused by defective leads are specifically excluded.

Statutory Legal Rights

This warranty is an additional guarantee granted by BIOTRONIK and does not affect any mandatory statutory legal rights based on applicable national legislation. This warranty is exclusively governed by German law.

Disclaimer

The manufacturer makes no guarantee that the Device is suitable for the proposed treatment of the individual patient for whom the Device is intended or that an adverse reaction will not take place between the implanted Device and the patient.

Neither BIOTRONIK SE & Co. KG, nor its distributors shall be liable for any loss, damage, or injury of any nature - whether direct or indirect, arising out of or resulting from the use of the product indirect, or consequential - in connection with, or resulting from, the use of the product. No person has any authority to bind BIOTRONIK SE & Co. KG or its distributor to any agreement or warranty other than those listed above.

49. Storage period and service time

The storage period affects the battery service time.

- Devices should be implanted within 19 months between the date of manufacture and the use by date (indicated on the package);
- If the ICD is implanted shortly before the use by date, the expected service time may be reduced by up to 16 months.

50. Calculation of service times

The services times have been calculated as follows:

Pulse amplitude	2.5 V
Pulse width	0.4 ms
Pacing impedance	500 Ω
Basic rate	60 bpm
Home Monitoring	ON, 1 device message each day and 12 IEGM online HD transmissions per year
Diagnostic functions and recordings	Permanently set

Calculation of the number of shocks

Capacitor reforming is performed 4 times per year and therefore at least 4 maximum charges for shocks have to be assumed per year even if less than 4 are delivered.

Iforia 3/5 VR-T

Service times with GB 2992 or LiS 3410 RA/RR battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	10.42	8.39	7.01	6.03	5.28
15%	10.14	8.20	6.89	5.93	5.21
50%	9.55	7.81	6.60	5.72	5.05
100%	8.81	7.31	6.24	5.45	4.83

Iforia 5 VR-T DX

Service times with GB 2992 or LiS 3410 RA/RR battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	9.48	7.76	6.57	5.70	5.03
15%	9.24	7.61	6.46	5.61	4.96
50%	8.75	7.26	6.21	5.42	4.81
100%	8.12	6.83	5.89	5.17	4.62

Iforia 3/5 DR-T

Service times with GB 2992 or LiS 3410 RA/RR battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	9.48	7.76	6.57	5.70	5.03
15%	9.02	7.45	6.35	5.53	4.89
50%	8.10	6.81	5.88	5.17	4.61
100%	7.08	6.07	5.32	4.73	4.26

Iforia 3/5 HF-T

Service times with GB 2992 or LiS 3410 RA/RR battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	8.78	7.29	6.23	5.44	4.82
15%	8.21	6.89	5.94	5.21	4.65
50%	7.14	6.12	5.35	4.76	4.28
100%	6.01	5.27	4.69	4.23	3.85

Iforia 7 VR-T

Service times with GB 2992 battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	11.78	9.52	7.98	6.87	6.03
15%	11.48	9.32	7.84	6.76	5.95
50%	10.81	8.87	7.52	6.53	5.76
100%	9.99	8.31	7.11	6.21	5.52

Iforia 7 VR-T DX

Service times with GB 2992 battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	10.73	8.82	7.48	6.50	5.74
15%	10.48	8.65	7.36	6.40	5.66
50%	9.92	8.26	7.08	6.19	5.50
100%	9.22	7.77	6.71	5.91	5.27

Iforia 7 DR-T

Service times with GB 2992 battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	10.73	8.82	7.48	6.50	5.74
15%	10.22	8.47	7.23	6.31	5.59
50%	9.20	7.76	6.70	5.90	5.27
100%	8.05	6.92	6.07	5.40	4.87

Iforia 7 HF-T

Service times with GB 2992 battery:

Pacing	Service time [in years] at number of shocks per year				
	4	8	12	16	20
0%	9.96	8.29	7.10	6.20	5.51
15%	9.33	7.85	6.77	5.95	5.31
50%	8.12	6.97	6.11	5.43	4.89
100%	6.85	6.01	5.36	4.83	4.40

It is forbidden to use this medical device after the expiration of its service life, since in this case it becomes impossible to perform the prevention of sudden cardiac death, as well as bradytherapy performing.

The remaining battery level is displayed on the screen of the programmer in time of performed follow-up, based on the schedule of planned visits agreed upon between the doctor and the patient.

Since ERI (Elective Replacement Indicator) status will be reached - the doctor informs the patient about the need of device changing procedure performing.

51. Manufacturer information

BIOTRONIK (currently BIOTRONIK SE & Co. KG) was established in 1963 in West Berlin by Professor Max SCHALDACH. The merits of Professor Schaldach in the field of Biophysics have been recognized by leading universities in America, Asia and Europe. More than 20 years ago Professor Schaldach has established extensive cooperation with leading clinical and academic institutions and was elected an honorary doctor of Moscow state University (1994), foreign member of the Ural branch of the Russian Academy of Sciences (1998), honorary doctor of the Russian Military-Medical Academy in St. Petersburg (1999).

During its half-century history BIOTRONIK has become a world leader in the design and manufacturing of implants for interventional cardiology and arrhythmology, equipment and consumables for electrophysiological studies and pacing of the heart. Currently, the centers for development and research and production units BIOTRONIK located on two continents in Germany, USA and Switzerland.

The company produces and supplies the world market implantable devices for the treatment of cardiac arrhythmias, external pacemakers and temporary electrodes, catheters for electrophysiological diagnostic and radiofrequency ablation of the heart, coronary and peripheral stents packaged with other products for angioplasty of blood vessels.

BIOTRONIK carries out successful research and clinical collaboration with leading medical centers in Russia for the past 20 years.

Contacts of official representative of BIOTRONIK:

LLC "Biotronik",
Moscow, Russia,
109028, Moscow, Zemlyanoy Val, 50A/8, bld. 2.
Phone: +7 (495) 789-68-31.
Fax: +7 (495) 789-68-32.
Web: www.biotronik.ru

52. Table of Contents

1.	Name of the medical device	2
2.	Information about the manufacturer	3
3.	Intended Use	4
4.	Indication for use	4
5.	Contraindications	5
6.	Possible complications and side effects while using the product.....	5
7.	Precautions while using the product.....	6
8.	Possible complications	7
9.	Possible risks	8
10.	Information about the impact on the ability to control the vehicle	9
11.	Application methods.....	10
12.	Classification of the product	10
13.	Description of principles that underpin the work of the medical device.....	11
14.	Scope of delivery.....	11
15.	Delivery package.....	13
16.	ICD Family, Model versions.....	13
17.	General ICD description	14
18.	Electrical Functional Description	15
19.	Description of accessories and devices combinations.....	17
20.	Programming.....	18
21.	Programmer	19
22.	Modes	19
23.	NBD and NBG codes	19
24.	BIOTRONIK Home Monitoring	20
25.	Method of product use.....	20
26.	The selection and treatment of the patients.....	22
27.	Implantation.....	23
28.	Applying the programming head.....	26
29.	Establishing telemetry contact.....	26
30.	Activating ICD therapy.....	26
31.	Precautionary Measures while Programming	27
32.	Follow-up.....	29
33.	Patient consulting	30

34.	Replacement Indications	31
35.	Explantation and Device Replacement.....	32
36.	Technical Specifications	32
37.	Labeling description	34
38.	Parameters.....	36
39.	Mechanical Characteristics.....	45
40.	Technical Characteristics	46
41.	Battery data	49
42.	Compliance to Recognized Standards	49
43.	Biocompatibility	50
44.	General Safety instructions	51
45.	Maintenance and repair requirements	52
46.	Disposal and special precautions	52
47.	Environmental protection requirements in time of Medical Device use.....	53
48.	BIOTRONIK SE & Co. KG Warranty conditions	53
49.	Storage period and service time	54
50.	Calculation of service times.....	55
51.	Manufacturer information.....	58
52.	Table of Contents	59

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