NOTARIAL CERTIFICATE

TO ALL TO WHOM these presents shall come. I, **MOHD YACOB KARIM**, a Notary Public, officiating at Lot 10, 3rd Floor, Pusat Dagangan UMNO Shan Alam, Persiaran Damai, Seksyen 11, 40100 Shah Alam, Selangor Darul Ehsan, Malaysia duly authorised and appointed, practicing in the State of Selangor Darul Ehsan, Malaysia DO HEREBY CERTIFY AND ATTEST THAT I was present on this 18th Day of May, 2018 and did see **MR. LEE KIM MEOW**, the Managing Director of **TOP GLOVE SDN. BHD.** (Company No. 220483-T) sign and execute the **OPERATIONAL DOCUMENTATION FOR RUSSIA.**

IN FAITH AND TESTMONY, where of I have hereunto subscribed my name and affixed my Seal of Office Shah Alam, Selangor Darul Ehsan, Malaysia this 18th Day of May, 2018 in this year of our Lord Two Thousand Eighteen (2018).

MOHD. YACOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR
MALAYSIA
Lot 10, 3rd Floor,
Pusat Dagangan UMNO Shah Alam
Persiaran Damai, Seksyen 11
40100 Shah Alam, Selangor
MALAYSIA
Tel : 603-5513 5610
Fax : 603-5513 5609

MOHD YACOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR DARUL EHSAN
MALAYSIA

18 MAY 2018


25 MAY 2018

This is to certify that the signature
appears on this document/Certificate/ Marriage
Certificate/Birth/Death Certificate is that of
MOHD YACOB KARIM who
is Notary Public

The Ministry of Foreign Affairs, Malaysia is not
responsible of the accuracy of the information
contained therein.

ZULKIPLY ENDUT
Consular Officer
Consular Division
Ministry of Foreign Affairs
Putrajaya Malaysia

25 MAY 2018

TOP GLOVE SDN. BHD.


Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia
Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

OPERATIONAL DOCUMENTATION FOR RUSSIA

Gloves examination latex powder free non-sterile DeMAX
produced by

Top Glove Sdn. Bhd.

"Approved"
"Top Glove Sdn. Bhd."
Managing Director
Lee Kim-Meow

/Signature/

Stamp *Top Glove Sdn. Bhd. *
TOP GLOVE SDN. BHD
(Company No. 220483-T)
A member of Top Glove Corporation Bhd
LOT 4969, JALAN TERATAI, BATU 6, Date
OFF JALAN MERU, 41050 KLANG 17.05.2018
SELANGOR D. E. MALAYSIA.
TEL : 603-3392 1992/ 1905
FAX : 603-3392 1291/ 8410

L.S.

MOHD. YACOB KARIM
NOTARIUM
NOTARY PUBLIC
SELANGOR
MALAYSIA
Lot 10, 3rd Floor.
Pusat Dagangan UMNO Shah Alam
Persiaran Damai, Seksyen 11
40100 Shah Alam, Selangor
MALAYSIA
Tel : 603-5513 5610
Fax : 603-5513 5609

18 MAY 2018

25 MAY 2018

This is to certify that the signature
appears on this document/Certificate/ Marriage
Certificate/Birth/Death Certificate is that of
MOHD YACOB KARIM who
is Notary Public

The Ministry of Foreign Affairs, Malaysia is not
responsible of the accuracy of the information
contained therein.

Zulkipli Endut
ZULKIPLY ENDUT
Consular Officer
Consular Division
Ministry of Foreign Affairs
Putrajaya Malaysia

25 MAY 2018

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia
Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

Section	Contents	Page
1	Product information	3
1.1	MD name and designation	3
1.2	Medical device manufacturer	3
2	Description and operation	3
2.1	Brief information about medical product, including its intended use	3
2.2	Technical parameters and characteristics guaranteed by the manufacturer	4
2.3	Package contents (scope of delivery) of medical product	5
2.4	Warnings, precautions and restrictions when using MD	5
2.5	Description of packaging, labeling and warning signs on the product	5
3	MP operation	7
3.1	Installation and commissioning	7
3.2	Working procedure (use of the product for its intended purpose)	7
4	Maintenance and repair	7
5	Disposal	7
6	Guarantees	8
7	Indications and contraindications for use, possible side effects and risks of use.	8
8	List of standards used by the manufacturer (producer) of the medical product	9

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia
Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

Product information

1.1. MD name and designation

Gloves examination latex powder free non-sterile DeMAX, execution options:

- Gloves examination latex powder free non-sterile DeMAX textured, sizes XS, S, M, L, XL;
- Gloves examination latex powder free non-sterile DeMAX textured polymer coated, sizes XS, S, M, L, XL.

1.2. Medical device manufacturer

Manufacturer data:

Name: Top Glove Sdn. Bhd.

Address: Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

E-mail: sales@topglove.com.my

Web-site: www. topglove.com.my

Data on the authorized representative of the manufacturer in the Russian Federation:

Limited Liability Company "Ardeil-Impeks" (LLC "Ardeil-Impeks")

15201, Moscow, Kashirskiy proezd, 23, b.5, office 3, Russia

Phone: 8 495 777 1170, 8 495 777 1171, e-mail: abx@ardl.ru, TIN 7724781425

2. Description and operation

2.1. Brief information about medical product, including its intended use

The product is non-sterile gloves made of natural latex intended exclusively for single use in medical examination, diagnostic and therapeutic procedures to protect a patient and an examiner from mutual contamination.

If applicable, the gloves have a polymer coating inside to facilitate the process of putting on, including on wet hands.

The working surface of the gloves has a pronounced texture that provides a reliable grip and retention of the tool, prevents its sliding in the hands when working in the conditions of high humidity.

Classification of the manufactured gloves:

- 1) depending on the material - natural rubber latex;
- 2) depending on the finish - textured, the texture is applied to the entire surface to ensure a reliable grip and retention of the tool, prevents its sliding in the hands when working in the conditions of high humidity.
- 3) depending on the finish - powder free;
- 4) depending on the edge of the cuff - with the roller at the edge;
- 5) depending on the color - light yellow color;
- 6) depending on the sizes - XS, S, M, L, XL;
- 7) depending on the form - non-anatomical.

The gloves can be supplied with or without polymer coating.

The gloves are universal (for both hands).

Application scope: all medical professions.

Operation conditions and potential consumers: medical personnel within the healthcare facilities.

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

2.2. Technical parameters and characteristics guaranteed by the manufacturer



Material, used for manufacturing

Ingredient	Chemical name	Function	%
Natural latex CAS 9003-31-0		Main component	59.3
Calcium carbonate CAS 1317-65-3	Calcium carbonate	Filler	29.7
Zinc diethyldithiocarbamate CAS 136-23-2	Zinc diethyldithiocarbamate	Accelerator B	1.06
Zinc dibutyldithiocarbamate CAS 14324-55-1	Zinc dibutyldithiocarbamate	Accelerator A	1.15
Calcium carboxylate CAS 471-34-1	Calcium carboxylate	Release agent	2.09
Calcium nitrate CAS 10124-37	Calcium nitrate	Activator	6.7

Without coloring agent

For execution option with acrylic polymer coating: Acrylic polymer coating Nobel 3009 CAS 68163-66-2 1 micron (1.6% by weight) manufactured by Nobel Synthex (Malaysia).

Main technical characteristics

The points for measuring the width and length of the gloves are shown in the figure below.	The points for measuring the thickness of the glove wall are shown in the figure below.
	
	The thickness of the glove wall at the cuff level is measured at a distance of 25 ± 5 from the edge of the cuff

Physical characteristics of gloves, for execution options:

- Gloves examination latex powder free non-sterile DeMAX textured, sizes XS, S, M, L, XL;
- Gloves examination latex powder free non-sterile DeMAX textured polymer coated, sizes XS, S, M, L, XL.

Size	Abbreviated size name (size code)	Length (mm)	Thickness, one layer (mm)			Width (mm)
			At the level of finger tips	At the palm level	At the cuff level	
X-Small	XS	240±5	0.14±0.01	0.11±0.01	0.09±0.01	75±5
Small	S					85±5
Medium	M					95±5
Large	L					105±5
X-Large	XL					115±5

The weight of one glove is not more than 7 g (regardless of the size).

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

Physical properties of gloves, for execution options:

-Gloves examination latex powder free non-sterile DeMAX textured, sizes XS, S, M, L, XL;

-Gloves examination latex powder free non-sterile DeMAX textured polymer coated, sizes XS, S, M, L, XL.

Before accelerated aging		After accelerated aging	
Breaking force, N, minimum	Elongation at break, %, minimum	Breaking force, N, minimum	Elongation at break, %, minimum
7.0	650	6.0	500

Powder free gloves, residual amount of dusting substance < 2.0 mg per glove in accordance with ASTM standard 3578

Gloves are airtight, airtightness quality level (AQL) is 2.5 maximum

Gloves with a circular cross-section roller on the edge of the cuff, the diameter of the roller is 1.8-2.2 mm

2.3. Package contents (scope of delivery) of medical product

One transport package contains 10 consumer packages and one operational document. There are 100 (50 pairs) pieces of gloves in one consumer box.

2.4. Warnings, precautions and restrictions when using MD

Precautions

- 1) The hand skin shall be in good condition, without wounds, micro-cracks and other damages.
- 2) Nails shall be short and neatly cut. Do not use artificial nails.
- 3) Jewelry and wrist watches during operation of gloves shall be removed.
- 4) The gloves are intended exclusively for single use in medical examination, diagnostic and therapeutic procedures to protect a patient and an examiner from mutual contamination.
- 5) May cause allergic reactions to rubber latex in individuals.
- 6) Do not treat with alcohol-containing solutions, since these solutions destroy the protective layer of gloves.

2.5. Description of packaging, labeling and warning signs on the product

Packaging

The sizes of a consumer cardboard box (Length x Width x Height, in cm, maximum): 30x20x10.

The weight of the package with gloves (100 pcs.) ss not more than 900 g (regardless of the execution option)

Group packing of coated double-ply cardboard on the back side of brown DB 500 grade manufactured by Siam Kraft Industry Co., Ltd. (Thailand).

Transport packaging from corrugated cardboard of CORRUGATING MEDIUM grade manufactured by MUDA Paper Mills Snd.Bhd. (Malaysia).

Labeling

Labeling is reliable, verifiable and readable. Labeling is applied to consumer packaging in any way (printing, embossing, stamping) ensuring its clarity, legibility and readability. When using the printed method of labeling, scuffing is not allowed.

The labeling on the consumer package of the gloves contains:

- name and location of the manufacturer;
- name, location and address of the official site of the authorized representative of the manufacturer;
- sign: "Natural latex that is part of the product can cause allergic reaction in individuals" or an equivalent sign and/or pattern that clearly reflects these indications;
- signs: "DeMAX";
- signs describing the abbreviated and full name of the product: "GLOVES EXAMINATION LATEX", "GLOVES EXAMINATION LATEX POWDER FREE NON-STERILE TEXTURED" or "GLOVES EXAMINATION LATEX POWDER FREE NON-STERILE TEXTURED POLYMER-COATED";
- sign: "NON-STERILE" and/or pattern, clearly reflecting this indication;

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

- sign: "SINGLE-USE";
- sign: "ARCHDALE";
- sign: "EXAMINATION" and/or pattern, clearly reflecting this indication;
- sign: "LATEX" and/or pattern, clearly reflecting this indication;
- sign: "NON-ANATOMICAL SHAPE" and/or pattern, clearly reflecting this indication;
- sign: "TEXTURED" and/or pattern, clearly reflecting this indication;
- sign: "POWDER FREE" and/or pattern, clearly reflecting this indication;
- sign (if appropriate): "POLYMER-COATED" and/or pattern, clearly reflecting this indication;
- information on the number and date of the registration certificate (indicated after the fact of state registration);
- information about the number of the party;
- information on the date of production;
- shelf life;
- size and number of pairs of gloves;
- information about the conditions of storage and transportation.

The marking shall be applied in the country of circulation of the product, excluding the proper names and titles.

Symbol	Symbol explanation
	Examination
	Latex
	Non-sterile
	Non-anatomical shape
	Textured
	Powder free
	Polymer-coated
	Batch
	Manufacturing date
	Shelf life

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

3. MP operation

3.1. Installation and commissioning

Open the pack of gloves, remove one pair of gloves and put the gloves on the right, then on the left hand in turn.

3.2. Working procedure (use of the product for its intended purpose)

Turn the bent edge of the left and right gloves in turn so that the glove completely covers the cuff of the dressing gown.

Taking off the gloves

Grip the glove in the palm base area with the fingers of the other hand.

Pull off the glove and hold it with the other hand.

Place the fingers of the free hand in the glove cuff on the other hand and pull it off like the first one.

Turning the glove inside out, use it as a cover for the other glove.

Place the gloves in a designated container.

After taking off the gloves, wash your hands with soap and put alcohol-containing cutaneous antiseptic on them.

Change gloves:

- every two hours of continuous operation;
- after detecting a visible defect;
- after a possible or actual puncture;
- if the glove begins to swell, stretch, stick or is loose on the hand;
- in accordance with the rules adopted at the place of application of these gloves.

Do not use the same pair of gloves when switching from one patient to another or when switching from one manipulation to another.

4. Maintenance, current repair, climatic conditions of operation/transportation/storage

The gloves are for single use and are not subject to maintenance or repair.

Gloves are operated at a temperature (from +10 to +30) °C and humidity (80-100)%.

Store in the original packaging in a dry and well ventilated place, excluding the presence of dust, exposure to sunlight, moisture, X-rays at a temperature not lower than - 50 °C and not more than 40 °C and humidity (60-75)%.

Protect from the effects of oils, gasoline and other substances destroying the rubber.

The product is transported by all types of covered vehicles, in accordance with the transport rules applicable to the mode of transportation at a temperature not lower than - 50 °C and not higher than 50 °C and humidity (60-75)%.

After transportation at minus temperature, unpack, condition for 24 hours at room temperature.

5. Disposal

The gloves shall be disposed of and/or destroyed safely. At the same time, all local, regional and federal regulations shall be complied with.

This product shall be disposed of in accordance with regulatory requirements adopted in the handling area. Information on proper disposal of the gloves can be obtained from the operational documentation for the area where the medical device is handled.

Disposal and destruction of used products in Russia shall be carried out in accordance with the requirements of SanPiN 2.1.7.2790-10, the used gloves are classified as class B waste (epidemiologically hazardous waste). Collection of used gloves shall be carried out in sealed bags of yellow color with the labeling (inscription).

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia
Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

6. Guarantees

The gloves belong to single-use medical products.

The shelf life of gloves is 4 years from the date of manufacture.

Guaranteed shelf life of the gloves is 4 years.

7. Indications and contraindications for use, possible side effects and risks of use.

Indications: for single use in medical examinations, diagnostic and therapeutic procedures to protect a patient and an examiner from mutual contamination.

Contraindications to use: the possibility of allergic reactions to natural latex in individuals.

Possible side effects: the possibility of allergic reaction to natural latex in individuals in the form of contact dermatitis.

Application risks:

Hazards	Risks	Dangerous situations	Harm
A. Biological and chemical risks	Insufficient control of the adjusted production schedule for the leach tank.	An inadequate leaching process to wash out excess protein and chemical that contacts sensitive skin. The strength of gloves is reduced, and the glove can tear when the user uses it.	Skin allergy. Viral or bacterial infection.
B. Operational risks	Insufficient control of the production process. Use after the shelf life expires. Poor control of the main door (metal door), closed near production and storage facilities.	Parameter from the specification, which can lead to a point defect on the glove. Reducing the physical characteristics of the gloves. The glove can tear when the user uses it. The front door has been opened, which can lead to ingress (insect) during the immersion and storage process.	Viral or bacterial infection. Infection of the user. Glove with an insect.
C. Information risks	Insufficient warning of hazards associated with the re-use of disposable gloves. Insufficient warning of side effects.	Re-use of non-sterile disposable gloves Insufficient warning on the package about use in the product.	Use of an infected glove and a decrease in the strength of gloves. A user with increased sensitivity may accidentally use gloves

TOP GLOVE SDN. BHD.

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor Darul Ehsan, Malaysia

Tel: +603-3392 1992 / 1905 E-mail: sales@topglove.com.my

List of standards used by the manufacturer (producer) of the medical product

No.	Standard	Description	Publishing date
1	ASTM 3578-05(2015)	Rubber examination gloves	November 2015
2	ASTM 5151-06(2015)	Detection of holes in medical gloves	November 2015
3	ASTM 6124-06(2017)	Residual powder on medical gloves	May 2017
4	EN 455-1:2000	Disposable medical gloves. Part 1: Requirements and detection of integrity.	October 2000
5	EN 455-2:2015	Disposable medical gloves. Part 2: Methods of detecting physical and mechanical properties.	April 2015
6	EN 455-3:2015	Disposable medical gloves. Part 3: Requirements and testing for biological evaluation.	April 2015
7	EN 455-4:2009	Disposable medical gloves. Part 4: Requirements and testing for determining shelf life.	October 2009
8	ISO 11193-1:2008	Disposable medical gloves. Part 1. Requirements for the gloves made of rubber latex or rubber-based mixtures.	June 2012
9	ISO 14971:2012	Medical devices. Application of risk management to medical devices.	July 2012
10	ISO 2859-1 :2011	Alternative signs selective control procedures	June 2011
11	ISO 10993-1:2003	Biological evaluation of medical devices. Part 1. Evaluation and testing	August 2003
12	ISO 10993-5:2009	Biological evaluation of medical devices. Part 5. In vitro cytotoxicity testing	June 2009
13	ISO 10993-10:2010	Evaluation of biological medical devices. Part 10. Tests for irritation and skin allergy reactions	August 2010
14	ISO 10993-12:2012	Biological evaluation of medical devices. Part 12. Preparation of testing samples and standard samples.	June 2012
15	ISO 14155-1:2011	Clinical testing of medical devices for human. Part 1. General requirements	January 2011
16	ISO 14155-2:2011	Clinical testing of medical devices for human. Part 2. Clinical testing schemes	January 2011
17	ISO 15223-1:2012	Medical devices. Symbols used on the labels of medical devices, when labeling and in the information provided. Part 1. General requirements	July 2012

The following standards are applied to the product in the territory of the Russian Federation:

GOST R 50444-92 "Devices, apparatuses and equipment for medical use. General technical conditions",

GOST R 52239-2004 "Medical diagnostic disposable gloves. Part 1. Specification for gloves made of rubber latex or solution"

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