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
Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Aoibheann Gill, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Aoibheann Gill

Date: 25th Sept 2018

Aoibheann Gill
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co Longford
Ireland

ABBOTT DIAGNOSTICS DIVISION LONGFORD

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Gerardine Courtenay ^{26/09/18} Authorised Official

Alinity i

TSH Calibrators

en

TSH

REF 07P4801

G72192R02

S7P480

Read Highlighted Changes: Revised December 2017.

NAME

Alinity i TSH Calibrators (also referred to as TSH Cals)

INTENDED USE

The Alinity i TSH Calibrators are for the calibration of the Alinity i analyzer when used for the quantitative determination of human Thyroid Stimulating Hormone (TSH) in human serum and plasma. For additional information, refer to the Alinity i TSH reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CAL 1 contains TRIS buffer with protein (bovine) stabilizers.

CAL 2 contains TSH (recombinant) in TRIS buffer with protein (bovine) stabilizers. Preservative: sodium azide.

The calibrators are at the following concentrations:

Calibrator	Quantity	TSH CONC (μ IU/mL)	TSH CONC (mIU/L)
CAL 1	1 x 3.0 mL	0	0
CAL 2	1 x 3.0 mL	40	40

MATERIALS REQUIRED BUT NOT PROVIDED

- 04R1001 Alinity ci-series Calibrator/Control Replacement Caps

STANDARDIZATION

The calibrators are manufactured by addition of Recombinant Human TSH of known concentration to obtain a target concentration. The concentration is referenced against WHO TSH 80/558.

PRECAUTIONS

- IVD**
- For *In Vitro* Diagnostic Use
- Rx ONLY**

The following warnings and precautions apply to: CAL 1 and CAL 2	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped with new replacement cap. Return to refrigerated storage after use.

The analyzer will track In-use Stability, which is the time the calibrator is outside of refrigerated storage while on the analyzer. The analyzer will not allow the use of the calibrator if the In-use Stability has been exceeded. Maximum In-use Stability can be found in the Assay Parameter Report. For additional information on calibrator In-use Stability, refer to the Alinity ci-series Operations Manual, Section 5.

For additional information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

INSTRUMENT PROCEDURE

- Calibrator lots may be configured using the bar code label on the calibrator carton.
- For information on configuring calibrator data, refer to the Alinity ci-series Operations Manual, Section 2.
- For instructions on ordering and loading calibrators on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

QUALITY CONTROL PROCEDURES

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.

For information on ordering controls, refer to the Alinity ci-series Operations Manual, Section 5.

Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of the associated reagent package insert.
- If statistically-based quality control limits are not available then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

For additional information, refer to the assay reagent package insert and the Alinity ci-series Operations Manual.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, if calibration does not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria, or if controls do not meet the appropriate criteria.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	<i>In Vitro</i> Diagnostic Medical Device
	Lot Number
	List Number
	Serial number

Other Symbols	
	Calibrator (1,2)
	Control Number
	Concentration
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Information needed for United States of America only
	Product of Ireland
	For use by or on the order of a physician only (applicable to USA classification only).

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised December 2017.

©2016, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i TSH Calibrators (for Russia)

SHIPPING

Recommended shipping conditions for device is ambient. Device is stable until the expiration date indicated on the package. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 13 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP) as amended, 1907/2006/EC (REACH) as amended, 80/181/EEC, 2009/3/EC, Regulation (EC) 1069/2009, Regulation (EC) 142/2011, CLASS 2013 (as amended), GHS (as amended), WHSR Part 7.1 (as amended), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980, EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23640, EN ISO 20776-1, EN 13975, EN ISO 17511, ISO 80000-1, EN ISO 15193, EN ISO 15194

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbott diagnostics.com or contact your local representative.

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