Instructions for Use
Kits of The ReSolve Locking Drainage Catheter
Dated June 11, 2018

ACKNOWLEDGEMENT


State of Utah)


§

County of Salt Lake)

On this 11th day of June, 2018, before me, Jennifer Wojciechowski, a notary public, personally appeared Dan Clark, Director, Regulatory Affairs, proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand and official seal.

Notary Public

My Commission Expires: March 6, 2019

Dan Clark
Director, Regulatory Affairs
Merit Medical Systems, Inc.
11 June 2018

**MEDICAL DEVICE
INSTRUCTIONS FOR USE**

« Kits of The ReSolve Locking Drainage Catheter »

**Manufactured by
Merit Medical Systems, Inc., the USA**

DESCRIPTION OF THE PRODUCT

The ReSolve Locking Drainage Catheter with locking pigtail and hydrophilic coating is a radiopaque catheter with multiple side holes used for percutaneous drainage of fluids from body cavities with abscesses, abdominal ascites, nephrostomy. The components of the catheter allow for introduction and placement using a trocar stylette or over-the-wire method.

The ReSolve Locking Drainage Catheter may be packaged with the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
- One (1) metal stiffening cannula
- One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
- One (1) trocar
- One (1) Repositioning Tool
- One (1) Pigtail Straightener
- One (1) Dead end cap


Figure 1 ReSolve Locking Drainage Catheters

The second packaging configuration (herein referred to as the “ReSolve NV”) is specifically designed to include components that accommodate the Seldinger’s technique. The ReSolve NV kit contains the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
 - One (1) metal stiffening cannula
 - One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
 - One (1) Repositioning Tool
 - One (1) Pigtail Straightener
 - One (1) Dead end cap
 - One (1) MAK NV system consisting of coaxial introducer with dilator, stiffening cannula, 21 GA introducer needle and 0.018 and 0.038 guide wires (ref. TEC0018) (not shown in figure below)
 - One (1) 600ml MDD with PVC bag (not shown in figure below)
 - One (1) Dilator (8F, 10F, 12F or 14F) (not shown in figure below)
- *6.5F and 7.5F catheter does-not require a dilator


Figure 2 ReResolve Locking Drainage Catheters NV

All configurations of the ReResolve Locking Drainage Catheter catheter kit may also include StayFix or Revolution catheter fixation devices designed to secure the catheter securely, as well as to prevent movement and accidental catheter removal (see Figure 2, 3, 4). If there is a StayFix or Revolution catheter in the set of devices, the kit names are as follows: ReResolve Locking Drainage Catheter Plus or ReResolve Locking Drainage Catheters NV.

INDICATIONS FOR USE

The ReResolve Locking Drainage Catheter with locking pigtail and hydrophilic coating is intended for percutaneous drainage of fluid:

- when removing the contents of an abscess in different localisations;
- in the abdominal cavity during the formation of serous fluid after a hemodynamic disorders such as right heart failure, portal hypertension (liver disease);
- abdominal irritations due to infection;
- in the abdominal cavity due to some tumor metastasis;
- upon termination of the flow of urine in the kidney or ureter proximal portion;
- to access pyelocaliceal system, ureter or bladder for diagnostic or therapeutic procedure;
- the formation of a stoma when ureteral injury.

CONTRAINDICATIONS

The ReResolve Locking Drainage Catheter is contraindicated for use where percutaneous drainage catheterization is unacceptable.

The ReResolve Locking Drainage Catheter is contraindicated for intravascular use.

RESTRICTED DEVICE

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device.

PRECAUTIONS

- Read manufacturer's instructions prior to use.
- Contents are sterile (via ethylene oxide) and non-pyrogenic.

- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/ or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Follow universal precautions when inserting and maintaining this device. Due to the risk of bloodborne pathogens, healthcare professionals should always use standard blood and body fluid precautions in the care of all patients. Sterile technique should always be followed.
- Do not use after expiration date.
- Store in a cool, dry place.

WARNING:

The ReSolve Locking Drainage Catheter is not to be used to deliver nutritional supplements. The ReSolve Locking Drainage Catheter is not intended for use in the biliary system.

ADVERSE REACTIONS*

- Septic shock
- Bacteremia
- Hemorrhage
- Superinfection
- Bowel transgression
- Pleural transgression
- Vascular injury
- Catheter dislodgement
- Catheter occlusion

INSTRUCTIONS FOR USE OPTION 1:

DIRECT PLACEMENT USING TROCAR STYLET

1. Wet the distal portion of the ReSolve Locking Drainage catheter prior to use with sterile water or saline. Keep the distal portion of the catheter wet during placement.

WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.

2. Slide pigtail straightener along distal portion of catheter to straighten curve prior to placing metal stiffening cannula into the catheter. Place the metal stiffening cannula into the catheter and tighten the Luer lock fittings. See Figure 3.

Figure 3

3. Remove the paper spacer from the trocar stylette. Advance the trocar stylette through the metal stiffening cannula and tighten the Luer lock fittings. See Figure 4.


Figure 4

4. Remove pigtail straightener from catheter prior to insertion.

5. Place the catheter/cannula/trocar assembly into the fluid collection site using standard insertion technique. Note: Placement should be confirmed with diagnostic imaging.

6. After placement is confirmed, remove the trocar stylette and stiffening cannula.

7. To engage the suture locking mechanism: Pull the suture until desired pigtail is formed. Rotate the suture locking mechanism distally to hold the suture in place. See Figure 5.


Figure 5

Note: If the catheter needs to be repositioned, unlock the suture locking mechanism by rotating the arm proximally to the point of resistance.

Precaution: Do not rotate the suture locking mechanism beyond the point of resistance. Rotating the suture locking mechanism beyond the point of resistance will not release the suture to allow the pigtail to straighten upon removal.

8. Once placement is confirmed, and the suture locking mechanism has been rotated to the most distal position, press the suture locking mechanism into the hub to secure it. The suture locking mechanism is now locked into position. See Figure 6.

Figure 6

9. The ReSolve Locking Drainage Catheter is now ready to be connected to appropriate drainage bag, tubing or dead end cap.

WARNING: If using alcohol to clean the catheter hub, allow sufficient time for alcohol to dry before connecting the drainage tubing or dead end cap.

WARNING: DO NOT over tighten the connection between the drainage catheter and drainage tubing or dead end cap. Note: A flush regimen should be designed for the circumstances of each patient and the protocol of the physician. Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

INSTRUCTIONS FOR USE OPTION 2: SELDINGER ENTRY TECHNIQUE OR GUIDE WIRE EXCHANGE

1. Remove the stiffening cannula and trocar stylette assembly from catheter.
2. Ensure that the distal portion of the catheter is wet prior to placement. Wet the distal portion of the ReSolve Locking Drainage catheter prior to use with sterile water or saline. Keep the distal portion of the catheter wet during placement. WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.
3. Slide pigtail straightener along distal portion of catheter to straighten curve prior to placing the stiffening cannula into the catheter. Place the stiffening cannula into the catheter and tighten the luer lock fittings. See Figure 7.

Figure 7

4. Remove pigtail straightener from catheter prior to insertion.
5. Place catheter/cannula assembly over appropriate guide wire and advance into the fluid collection site. The catheter accommodates a 0.038" (0.97 mm) wire. See Figure 8. Note: Placement should be confirmed with diagnostic imaging. Figure 8

Figure 8

6. After placement is confirmed, remove the stiffening cannula and guide wire.
7. To engage the suture locking mechanism: Pull the suture until desired pigtail is formed. Rotate the suture locking mechanism distally to hold the suture in place. See Figure 9.


Figure 9

Note: If the catheter needs to be repositioned, unlock the suture locking mechanism by rotating the arm proximally to the point of resistance.

Precaution: Do not rotate the suture locking mechanism beyond the point of resistance. Rotating the suture locking mechanism beyond the point of resistance will not release the suture to allow the pigtail to straighten upon removal.

8. Once placement is confirmed, and the suture locking mechanism has been rotated to the most distal position, press the suture locking mechanism into the hub to secure it. The suture locking mechanism is now locked into position. See Figure 10.

Figure 10

9. The ReSolve Locking Drainage Catheter is now ready to be connected to appropriate drainage bag, tubing or dead end cap.

WARNING: If using alcohol to clean the catheter hub, allow sufficient time for alcohol to dry before connecting the drainage tubing or dead end cap.

WARNING: DO NOT over tighten the connection between the drainage catheter and drainage tubing or dead end cap. Note: A flush regimen should be designed for the circumstances of each patient and the protocol of the physician. Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

MDD drainage bag.

Drainage bags MDD bags are made of translucent PVC, c feed liquid with a label for measuring fluid flow in "mL" on one side. The bag contains a built antireflux valve which prevents reverse fluid flow and disabling of the rotary valve, which facilitates the removal of fluid. The drainage bag is supplied with a strap for fixing to the patient and designed to meet the patient's mobility.


MDD drainage bag 600 ml volume It designed for percutaneous catch drainage when used in conjunction with a suction catheter and supplied in sets of drainage catheter "ReSolve Locking Drainage CatheterNV".

Dosing drainage bag:

1. Close empty port prior to use.


2. Thread belt as shown. Straps for securing the bag used for fixing the drainage bag on the leg with adhesive tape.

3. Connect drainage catheter to bag tubing as shown.

4. The bag is located below the drainage portion with a strap for fixing as shown below

Contact your healthcare professional if your drainage becomes inadequate, or if there is leakage around the catheter or drainage bag.

Disconnecting the bag from the tube occurs through the Luer connector. Emptying the drainage bag going through at the bottom of the drain hole with the cap, it is necessary to open the drain valve by turning it. The connecting tube is not twisted and does not prevent the free flow of fluid.

Represented MAK-NV was designed to facilitate the visualization and easy access to the drainage. The system increases support and ensures a smooth introduction. Represented system MAK-NV provides access to the extravascular system and facilitates placement of wires. MAK-NV system to minimize bleeding and trauma to the surrounding tissues.

Mode of application:

1. Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
2. Insert the 21 gauge introducer needle using standard technique.
3. Carefully advance the flexible end of the 0.018" (0.46mm) guide wire through needle. Advance the guide wire as far as appropriate. Verify correct positioning. **WARNING:** Do not advance guide wire if resistance is met.
4. Remove needle while maintaining the 0.018" (0.46mm) guide wire in position. **WARNING:** Withdrawal, pull back, or manipulation of the guide wire distal tip through

the needle tip may result in breakage or embolization. To avoid guide wire damage during manipulation, remove the introducer needle and proceed to step 5.

5. Insert the co-axial introducer with dilator, with or without stiffening cannula, over the 0.018" (0.46mm) guide wire and advance to the desired position.

6. Remove the dilator/stiffening cannula, leaving the introducer in position. The 0.018" (0.46mm) guide wire may be removed or left in position.

7. Insert up to 0.038" (0.97mm) guide wire through introducer.

8. Remove introducer leaving guide wire(s) in place.

9. Next, install the drainage catheter as described above (method 2).


CATHETER EXCHANGE, REPOSITIONING OR REMOVAL

1. Disconnect catheter from drainage bag, tubing or dead end cap.

2. To release the pigtail loop choose one of the following options:

Option 1:


- Using the Repositioning Tool align the opening of the round section of the Repositioning Tool in line with the handle of the suture locking mechanism.
- Bring the flat back of the Repositioning Tool around the catheter hub.
- Gently squeeze together.
- Remove the Repositioning Tool and rotate the suture locking mechanism to the most proximal position.

Precaution: Do not rotate the suture locking mechanism beyond the point of resistance. Rotating the suture locking mechanism beyond the point of resistance will not release the suture to allow the pigtail to straighten upon removal.

Option 2: For exchange or removal only, cut the hub off the drainage catheter and sever suture. This will release the suture and the pigtail loop.

WARNING: The suture will no longer be secured to the catheter. Take care to remove both the suture and catheter.

3. For catheter exchange or if access is to be maintained, advance appropriate guide wire through catheter; using diagnostic imaging to confirm guide wire placement. Guide wire will maintain access to drainage site. To ease guide wire placement, a stiffening cannula may be used.

WARNING: When long-term use is indicated, it is recommended that indwelling time not exceed 90 days. The ReSolve Locking Drainage Catheter should be evaluated by a physician on or before 90 days post placement.

4. Carefully remove the ReSolve Locking Drainage Catheter. Proceed with either catheter exchange or skin closure.

ATTENTION ATTENDING PHYSICIAN: IF PATIENT WILL NOT BE FOLLOWED UP BY YOU, IT IS RECOMMENDED THAT THE "INSTRUCTIONS FOR USE" OR THE SECTION ON HOW TO REMOVE THE CATHETER BE ATTACHED TO THE PATIENT'S CHART.

Introducer system MAK-NV (mini access for extravascular intervention) was developed to facilitate the visualization and quick access into the drainage area. To improve the visibility of the device has an echogenic needle and the radiopaque tip of the sheath. System increases support and provides a smooth introduction.

Stay Fix

Flaxation Device

DESCRIPTION OF THE PRODUCT

The StayFIX Fixation Device is a sterile adhesive fastener.

INTENDED USE

The StayFIX Fixation Device is intended to secure non-vascular catheters and tubes.

INDICATIONS FOR USE The StayFIX Fixation Device is indicated to adhere to the skin and secure a percutaneous non-vascular catheter or tube to help minimize movement and/or accidental removal.

Contraindications

- Use on patients with known tape or adhesive allergies
- Use with vascular catheters

RESTRICTED DEVICE

Rx Only: Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device.

Precautions:

- Read manufacturer's instructions prior to use.
- Contents are sterile.
- Do not use if packaging is opened, damaged, or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Do not use after expiration date.
- Store in a cool, dry place.
- After seven (7) days (sooner if the fixation device becomes loose for any reason) remove the old StayFIX Fixation Device and replace with a new StayFIX Fixation Device
- When removing the old StayFIX Fixation Device, use care not to move or accidentally pull out catheter or tube.

INSTRUCTIONS FOR USE

To adhere effectively, the StayFIX Fixation Device must be applied to a clean, dry skin surface. Be sure to follow your healthcare provider's instructions for cleaning and managing the area where the catheter or tube exits your body. Application instructions are illustrated below (Reverse instructions for removal).

1. Remove the StayFIX Fixation Device from sterile packaging.
2. Align the StayFIX Fixation Device slit with catheter or tube. (Figure 11 B, Figure 12 B)
3. Slide the StayFIX Fixation Device slit around catheter or tube until catheter or tube is positioned in hole against the blue foam. (Figure 11 C, Figure 12 C)
4. Peel off paper backing and adhere the StayFIX Fixation Device to skin. (Figure 11 B & C, Figure 12 B & C)
5. Peel off paper backing from the blue foam. (Figure 11 D, Figure 12 D)
6. Gently bend catheter or tube into the blue foam channel. (Figure 11 D, Figure 12 D)
7. Peel off red tab to remove paper backing (685ME Large only). (Figure 12 E)
8. Fold adhesive tab over catheter or tube to secure catheter or tube (685ME Large only). (Figure 12 E)
9. Peel off large paper backing from second tab. (Figure 11 E, Figure 12 E)
10. Fold second tab over the catheter or tube. (Figure 11 F, Figure 12 F)

Daily management of your catheter and StayFIX Fixation Device

Your healthcare provider will tell you how to inspect your catheter or tube and exit site for possible problems. Contact your physician if you notice the following:

- excess leakage
- soreness around the catheter's exit site
- movement of the catheter or tube

Inspect the StayFIX Fixation Device daily to make sure that it:

- is still holding the catheter or tube securely
- continues to stick to your skin
- maintains its overall integrity

You may shower with the StayFIX Fixation Device on your body. Avoid soaking for long periods of time and refrain from touching the StayFIX Fixation Device until it has dried completely. Again, if you notice the StayFIX Fixation Device becoming loose, replace it with a new one.

Figure 11



Figure 12

REVOLUTION

Catheter securement device

Description: Revolution catheter securement device.

INDICATIONS FOR USE: The Revolution catheter securement device is designed to secure a non-vascular catheter placed in a patient to help minimize catheter movement and accidental removal.

CONTRAINDICATIONS:

The Revolution catheter securement device is not intended to be used on vascular catheters.

CAUTIONS:

- Read instructions prior to use.
- Rx Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Contents of unopened, undamaged package are sterile.


REUSE PRECAUTION STATEMENT For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

INSTRUCTIONS FOR USE: Inspect the Revolution catheter securement device prior to use to verify that no damage has occurred during shipping.

A. Remove the shower lid and Revolution catheter securement device from the packaging. Use caution to ensure that you do not pull the finger rings prematurely and activate strings.

B. Before placing the Revolution catheter securement device, ensure that the skin where the device will be placed is clean and dry.

C. Position Revolution catheter securement device & secure catheter.

D. If desired, place the supplied lid over the Revolution catheter securement device to cover the catheter and puncture site while showering.

The Revolution catheter securement device is designed to secure catheters from 6F - 24F and to remain secure up to 7 days.

CAUTIONS: In order to secure the catheter, verify that the blue ring is rotated sufficiently. Under-tightening may allow the catheter to move. Over-tightening may compress the catheter and restrict or stop fluid flow.

CAUTIONS: Verify that the catheter is not kinked after deployment.

CAUTIONS: The Revolution catheter securement device can get wet and can be worn while showering. Merit Medical recommends having the device lid in place when showering. Avoid soaking the device for long periods of time.

CAUTIONS: Avoid using alcohol or acetone based solutions when cleaning around or in the center ring of the Revolution catheter securement device. These substances may weaken the device's adhesive. To remove the Revolution catheter securement device, remove the lid if applicable, cut the strings holding the catheter, and carefully peel the Revolution catheter securement device from the patient's skin.

The FastBreak connector.

The FastBreak Breakaway Connector is a two-piece coupling with a patient side connecting to the hub of the drainage catheter, and a second side connecting to the drainage bag during placement. In use, the two (2) connected components allow fluids to pass freely. Upon application of variable strength against the catheter tubing, the coupling will detach. The FastBreak connector lessens the chance of the drainage catheter being inadvertently detached. The FastBreak connector should not be used with biliary catheters.

FASTBREAK INSTALLATION

1. Connect the male luer end of the FastBreak connector to the open end of the allpurpose drainage catheter and twist in a clockwise direction. (See Step 2)
2. Connect the open female end of the FastBreak connector to the male luer end of the external fluid reservoir and twist in a clockwise direction. (See Step 3)
3. Connect the FastBreak male to female connectors together until you hear a "click" noise. (See Steps 4 and 5)

INDICATIONS FOR USE


The FastBreak connector is intended to serve as a conduit between a general drainage catheter and an external fluid reservoir.

PRECAUTIONS


- Follow instructions for use and institutional guidelines for disposal.
- Do not store at extreme conditions.
- Do not use in biliary application **REUSE PRECAUTION STATEMENT** For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

	Do not use if package is damaged				
	Do not re-use. Do not re-sterilize				
	Consult instructions for use				
	Not for individual sale.				
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.				
	Orange (light)		Green (medium)		Yellow (strong)

STEP 1


STEP 2 Male Connector

STEP 3 Female Connector

STEP 4

STEP 5

**MEDICAL DEVICE
ATTACHMENT INSTRUCTIONS FOR USE**

Kits of The ReSolve Locking Drainage Catheter

**Manufactured by
Merit Medical Systems, Inc., the USA**

1. DEVICE DESCRIPTION

Kits of The ReSolve Locking Catheter designed for percutaneous drainage of fluids from body cavities with abscesses, abdominal ascites, nephrostomy. The products are intended for use in clinics, hospitals and other medical institutions.

Kits of The ReSolve Locking Cathete is packaged with the following components:

- 1) Kit of the ReSolve Locking Drainage Catheter, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:
 - ReSolve Locking Drainage Catheter of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Trocar - 1 pc.
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
- 2) Kit of the ReSolve Locking Drainage Catheter, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:
 - ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Trocar - 1 pc.
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
 - FastBreak connector - not more than 6 pcs. / packaging
- 3) Kit of the ReSolve Locking Drainage Catheter, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:
 - ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
 - Flexible cannula
 - Trocar - 1 pc.
 - Repositioning Tool – 1 pc
 - Pigtail Straightener – 1 pc
 - Deadend cap -1 pc.
- 4) Kit of the ReSolve Locking Drainage Catheter, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:
 - ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc

- Flexible cannula
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc.
- FastBreak connector - not more than 6 pcs. / packaging

5) Kit of the ReSolve Locking Drainage Catheter Plus catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc
- StayFix catheter securement device – 1 pc

6) Kit of the ReSolve Locking Drainage Catheter Plus catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc
- StayFix catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

7) Kit of the ReSolve Locking Drainage Catheter Plus, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Flexible cannula
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc.
- StayFix catheter securement device – 1 pc

8) Kit of the ReSolve Locking Drainage Catheter Plus, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Flexible cannula

- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc.
- StayFix catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

9) Kit of the ReSolve Locking Drainage Catheter Plus catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc
- Revolution catheter securement device – 1 pc

10) Kit of the ReSolve Locking Drainage Catheter Plus catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc
- Revolution catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

11) Kit of the ReSolve Locking Drainage Catheter Plus, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Flexible cannula
- Trocar - 1 pc.
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Deadend cap -1 pc.
- Revolution catheter securement device – 1 pc

12) Kit of the ReSolve Locking Drainage Catheter Plus, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Flexible cannula
- Trocar - 1 pc.
- Repositioning Tool – 1 pc

- Pigtail Straightener – 1 pc
- Deadend cap - 1 pc.
- Revolution catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

13) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,

14) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.
- FastBreak connector - not more than 6 pcs. / packaging

15) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc

- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,

16) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
- FastBreak connector - not more than 6 pcs. / packaging

17) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
- StayFix catheter securement device – 1 pc

18) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with

- hydrophilic coating– 1 pc
- The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018” (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038” (0,97 mm), length 150 cm, ‘J’ tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
 - Merit Drainage Depot, 600 ml – 1 pc.,
 - Belt – 1 pc.
 - Luer Lock - 1 pc.,
 - StayFix catheter securement device – 1 pc
 - FastBreak connector - not more than 6 pcs. / packaging

19) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018” (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038” (0,97 mm), length 150 cm, ‘J’ tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
 - Merit Drainage Depot, 600 ml – 1 pc.,
 - Belt – 1 pc.
 - Luer Lock - 1 pc.,
 - StayFix catheter securement device – 1 pc

20) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018” (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038” (0,97 mm), length 150 cm, ‘J’ tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:

- dilator, diameter 6F – 1 pc.,
- the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
- StayFix catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

21) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
 - Revolution catheter securement device – 1 pc

22) Kit of the ReSolve Locking Drainage Catheter NV, catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); length 15 cm, 25 cm, 40 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
 - The metal cannula - 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038" (0,97 mm), length 150 cm, 'J' tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
 - Revolution catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging

23) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc

- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018” (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038” (0,97 mm), length 150 cm, ‘J’ tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
- Revolution catheter securement device – 1 pc

24) Kit of the ReSolve Locking Drainage Catheters NV, catheter diameter 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); 14F (4,6 mm); length 15 cm, 25 cm, 40 cm, dilator, diameter 8F (2,7 mm), 10F (3,3 mm), 12F(4,0 mm); 14F (4,6 mm), length 20 cm, consisting of:

- ReSolve Locking Drainage Catheter (RLC) of the "pigtail" (pigtail) with hydrophilic coating– 1 pc
- The metal cannula - 1 pc
 - Flexible cannula – 1 pc
- Repositioning Tool – 1 pc
- Pigtail Straightener – 1 pc
- Guide wire, diameter 0,018” (0,46 mm), length 60 cm - 1 pc.,
- Guide wire, diameter 0,038” (0,97 mm), length 150 cm, ‘J’ tip - 1 pc
- Introducer needle 21G (0,81mm) with the tip-trokar - 1 pc.
- Introduce:
 - dilator, diameter 6F – 1 pc.,
 - the metal cannula – 1 pc.,
- Merit Drainage Depot, 600 ml – 1 pc.,
- Belt – 1 pc.
- Luer Lock - 1 pc.,
- Revolution catheter securement device – 1 pc
- FastBreak connector - not more than 6 pcs. / packaging


ReSolve Locking Drainage Catheter (RLC)

The RLC is intended for percutaneous drainage of fluids from body cavities for up to 90 days.

The RLC is offered in 6.5, 7.5, 8.5, 10, 12 and 14 French sizes with a working length of 15, 25 and 40 cm with depth markers placed every centimeter beginning one (1) cm from the distal catheter tip and ending at the hub. The catheter shaft is made of radiopaque polyurethane and contains 8 or 9 drainage holes depending on French size with one of the drainage holes also serving as a suture hole. The 6.5F, 7.5F, 8.5F and 10F catheters have 8 drainage holes and the 12F and 14F catheters have 9 drainage holes. All RLC catheters are compatible with a 0.038” guide wire.

The RLC Catheter is packaged in several packaging configurations. The first configuration provides the clinicians with the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
- One (1) metal stiffening cannula
- One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
- One (1) trocar
- One (1) Repositioning Tool
- One (1) Pigtail Straightener
- One (1) Dead end cap

METAL STIFFENING CANNULA
 TROCAR STYLET
 PAPER SPACER
 HUB
 RESOLVE DRAINAGE CATHETER
 WITH LOCKING PIGTAIL
 SUTURE LOCKING MECHANISM
 REPOSITIONING TOOL
 PIGTAIL STRAIGHTENER
 FLEXIBLE STIFFENING CANNULA

Figure 1 ReSolve Locking Drainage Catheters

The second packaging configuration (herein referred to as the “ReSolve NV”) is specifically designed to include components that accommodate the Seldinger’s technique. The ReSolve NV kit contains the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
 - One (1) metal stiffening cannula
 - One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
 - One (1) Repositioning Tool
 - One (1) Pigtail Straightener
 - One (1) Dead end cap
 - One (1) MAK NV system consisting of coaxial introducer with dilator, stiffening cannula, 21 GA introducer needle and 0.018 and 0.038 guide wires (ref. TEC0018) (not shown in figure below)
 - One (1) 600ml MDD with PVC bag (not shown in figure below)
 - One (1) Dilator (8F, 10F, 12F or 14F) (not shown in figure below)
- *6.5F and 7.5F catheter does-not require a dilator


Figure 2 ReResolve Locking Drainage Catheters NV

METAL STIFFENING CANNULA
 TROCAR STYLET
 PAPER SPACER
 HUB
 RESOLVE DRAINAGE CATHETER
 WITH LOCKING PIGTAIL
 SUTURE LOCKING MECHANISM
 REPOSITIONING TOOL
 PIGTAIL STRAIGHTENER
 FLEXIBLE STIFFENING CANNULA

All configurations of the ReResolve Locking Drainage Catheter catheter kit may also include StayFix or Revolution catheter fixation devices designed to secure the catheter securely, as well as to prevent movement and accidental catheter removal (see Figure 2, 3, 4). If there is a StayFix or Revolution catheter in the set of devices, the kit names are as follows: ReResolve Locking Drainage Catheter Plus or ReResolve Locking Drainage Catheters NV.

The fixing device Revolution catheter includes a filament for stabilizing and securing the catheter in place while leaving the injection site open, and has a special cover to provide additional protection against ingress of water (see figure 3,5).

Figure 3 - StayFix

Figure 4 - Revolution

The ReSolve Locking Drainage Catheter Plus kit (or ReSolve Locking Drainage Catheters NV) includes all configuration components 1 (or configuration 2) with the addition of the StayFix or Revolution catheter fixation device (Figure 2,3,4).

The ReSolve Locking Drainage Catheter Plus configuration contains the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
- One (1) metal stiffening cannula
- One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
- One (1) trocar
- One (1) Repositioning Tool
- One (1) Pigtail Straightener
- One (1) Dead end cap
- One (1) StayFix or Revolution catheter fixation device

Configuration ReSolve Locking Drainage Catheters NV contains the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
- One (1) metal stiffening cannula
- One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
- One (1) Repositioning Tool
- One (1) Pigtail Straightener
- One (1) Dead end cap
- One (1) MAK NV system consisting of coaxial introducer with dilator, stiffening cannula, 21 GA introducer needle and 0.018 and 0.038 guide wires (ref. TEC0018) (not shown in figure below)
- One (1) 600ml MDD with PVC bag (ref. TEC0044) (not shown in figure below)
- One (1) Dilator (8F, 10F, 12F or 14F) (not shown in figure below)
- One (1) StayFix or Revolution catheter fixation device

* For a catheter with a diameter of 6.5 F (2.1 mm); A 7.5 F (2.4 mm) dilator is not required.

Also, all sets can be equipped with a FastBreak connector (see Figure 5). The FastBreak Breakaway Connector is a two-piece coupling with a patient side connecting to the hub of the drainage catheter, and a second side connecting to the drainage bag during placement. In use, the two (2) connected components allow fluids to pass freely. Upon application of variable strength against the catheter tubing, the coupling will detach. The FastBreak connector lessens the chance of the drainage catheter being inadvertently detached. The FastBreak connector should not be used with biliary catheters. Connectors are available in the following colors - yellow, orange, green, depending on the degree of fixation.

Figure 5

2. STERILIZATION METHOD

Merit Medical Systems, Inc. utilizes Ethylene Oxide (EO) sterilization for all ReSolve product configurations. The EO sterilization cycle is validated in accordance to the requirements of The European International Standard ISO 11135. Annual re-validation occurs assessing a worst-case reference load based on density and boxing configuration, which allows routine sterilization of a mixed-pallet load. Merit's products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals. The residual ethylene oxide content does not exceed 10^{-6}

The protocols and reports validating the sterilization process are maintained in Merit's Document Control system and are referenced on Merit's Sterilization Summary (Table 942).

This device is not intended to be re-sterilized by the user.

ReSolve Catheter product is less complicated than a standard process challenge device (PCD) that represents a worst case configuration and is readily adopted into Merit's validated sterilization cycle.

Revolution and StayFix or are sterilized using gamma radiation and have passed the validation process. Validation was performed according to EN ISO 11137-2 - Sterilization of health care products. Exposure. Part 2. Establishing a sterilizing dose. On the basis of this check the minimum dosage for sterilization is 25 kGy.

The standard was verified dosages and verification of bioburden in accordance with EN ISO 11137-2 to verify that the dose for sterilization continues to achieve the desired SAL (10^{-6}).

3. TECHNICAL CHARACTERISTICS

3.1 ReSolve Locking Drainage Catheter

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F	14 F
Catheter OD	0,085" ± 0,002" (2,16 MM ± 0,051 MM)	0,098" ± 0,002" (2,49 MM ± 0,051 MM)	0,112" ± 0,002" (2,84 MM ± 0,051 MM)	0,131" ± 0,002" (3,33 MM ± 0,051 MM)	0,158" ± 0,002" (4,01 MM ± 0,051 MM)	0,184" ± 0,002" (4,67 MM ± 0,051 MM)
Catheter ID	0,070" - 0,073" (1,78 MM - 1,85 MM)	0,070" - 0,073" (1,78 MM - 1,85 MM)	0,070" - 0,073" (1,78 MM - 1,85 MM)	0,089" ± 0,002" (2,26 MM ± 0,051 MM)	0,106" ± 0,002" (2,69 MM ± 0,051 MM)	0,118" ± 0,002" (3 MM ± 0,051 MM)
Catheter working Length (+/-1,5 cm)	15, 25,40	15, 25,40	15, 25,40	15, 25,40	15, 25,40	15, 25,40
Number of Drainage Holes	8	8	8	8	9	9
Pigtail Loop Size (cm)	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,54±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	2,794±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.	3,175±0,1 The return of the tip of the tendril to the stem no more than 30 seconds after the complete straightening of the catheter.
Distance between handles depth (cm)	1±0,05	1±0,05	1±0,05	1±0,05	1±0,05	1±0,05
Radiographic visibility	The catheter should be viewed under fluoroscopy					
Compatibility guidewire	guidewire of 0.038 inches (0.97 mm) with an installed cannula needs to pass through the entire structure in the Assembly.					
Compatibility connection, Luer-type	The hub of the catheter must be compatible with swivel and fixed adapters items covered by type (hemostatic valves, Luer lock connection, syringes, catheter valves, etc.).					
A cone-shaped tip	At the distal end of the tip should be tapered.					
Safety of magnetic resonance imaging (MRI)	The catheter should not interfere with the MRI.					

Surface	On the outer surface should not be technology and surface defects. This surface should be minimally traumatic to the vessel during use.
---------	---

3.2 Metal Stiffening Cannula

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F	14 F
Cannula Size	18	18	18	15	14	13
Cannula OD	0,050" ± 0,0005" (1,27 mm ± 0,013 mm)	0,050" ± 0,0005" (1,27 mm ± 0,013 mm)	0,050" ± 0,0005" (1,27 mm ± 0,013 mm)	0,072" ± 0,0005" (1,83 mm ± 0,013 mm)	0,083" ± 0,001" (2,11 mm ± 0,013 mm)	0,095" ± 0,001" (2,41 mm ± 0,013 mm)
Cannula ID	0,033" ± 0,0005" (0,84 mm ± 0,013 mm)	0,033" ± 0,0005" (0,84 mm ± 0,013 mm)	0,033" ± 0,0005" (0,84 mm ± 0,013 mm)	0,054" ± 0,0005" (1,372 mm ± 0,013 mm)	0,063" ± 0,0005" (1,60 mm ± 0,013 mm)	0,071" ± 0,0005" (1,80 mm ± 0,013 mm)
Exposed Length	19,452" ± 0,0395" (494 mm ± 1 mm)	19,452" ± 0,0395" (494 mm ± 1 mm)	19,452" ± 0,0395" (494 mm ± 1 mm)	19,442" ± 0,0395" (493,8 mm ± 1 mm)	19,362" ± 0,0395" (491,8 mm ± 1 mm)	19,342" ± 0,0395" (491,2 mm ± 1 mm)
The length of the exposed portion of the cannula subassembly	0,388" ± 0,040" (9,85 mm ± 1 mm)	0,360" ± 0,040" (9,14 mm ± 1 mm)	0,360" ± 0,040" (9,14 mm ± 1 mm)	0,390" ± 0,040" (9,9 mm ± 1 mm)	0,459" ± 0,040" (11,6 mm ± 1 mm)	0,490" ± 0,040" (12,4 mm ± 1 mm)
Allowable angle of deflection, without compromising the integrity	15° ± 1°					
Surface	On the outer surface should not be technology and surface defects.					

3.3 Flexible (Plastic) Stiffening Cannula

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F	14 F
Cannula OD	N / A *	N / A *	0.056 " ± 0,002 " (1.42 mm ± 0,05 mm)	0.072 " ± 0,002 " (1.82 mm ± 0,05 mm)	0.082 " ± 0,002 " (2.08 mm ± 0,05 mm)	0.095 " ± 0,002 " (4.41 mm ± 0,05 mm)
Cannula ID	N / A *	N / A *	0.042 " ±	0.051 " ±	0.063 " ±	0.071 " ±

			0,002 " (1.06 mm ± 0,05 mm)	0,002 " (1.29 mm ± 0,05 mm)	0,002 " (1.60 mm ± 0,05 mm)	0,002 " (1.80 mm ± 0,05 mm)
Exposed Length for 15, 25, 40 cm	N / A *	N / A *	9.452 " ± 0,0395 " (240 mm ± 1 mm); 13,395 " ± 0,0395 " (340 mm ± 1 mm) or 19.295 " ± 0,0395 " (490 mm ± 1 mm)	9.438 " ± 0,0395 " (239.8 mm), 13.375 " ± 0,0395 " (339.7 mm), or 19.285 " ± 0,0395 " (489.8 mm)	9.761 " ± 0,0395 " (248 mm), 13.305 " ± 0,0395 " (337.9 mm), or 19.205 " ± 0,0395 " (487.8 mm)	9.741 " ± 0,0395 " (247.4 mm), 13.285 " ± 0,0395 " (337.4 mm), or 19.185 " ± 0,0395 " (487.2 mm)
Allowable angle of deflection, without compromising the integrity	N / A *	N / A *	15° ± 1°			
Surface	N / A *	N / A *	On the outer surface should not be technology and surface defects.			

* Drainage catheters ReSolve size 6,5F and 7,5F It does not include a flexible rigid cannula.

3.4 Trocar

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F	Д 14 F
The diameter of the trocar for the catheter length 15, 25 or 40 cm	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)	0,038" ± 0,0001" (0.97 mm ± 0,0025 mm)
Length of Trocar	10,438 " ± 0,125 "(265 mm ± 3,17 mm) 14,375 " ± 0,125 "(365.1 mm ± 3,17 mm) or 20,281 " ± 0,125 "(515.1	10,410 " ± 0,125 "(264,4 mm ± 3,17 mm), 14,347 " ± 0,125 "(364,4 mm ± 3,17 mm) or 20,252 " ± 0,125	10,410 " ± 0,125 "(264,4 mm ± 3,17 mm), 14,347 " ± 0,125 "(364,4 mm ± 3,17 mm) or 20,252 " ± 0,125	10,420 " ± 0,125 "(264,6 mm ± 3,17 mm), 14,357 " ± 0,125 "(364,6 mm ± 3,17 mm) or 20,262 " ± 0,125	10,813 " ± 0,125 "(274,6 mm ± 3,17 mm), 14,356 " ± 0,125 "(364,6 mm ± 3,17 mm) or 20,261 " ± 0,125	10,824 " ± 0,125 "(274,9 mm ± 3,17 mm), 14,367 " ± 0,125 "(364,9 mm ± 3,17 mm) or 20,272 " ± 0,125

	mm ± 3,17 mm)	"(514,4 mm ± 3,17 mm)	"(514,4 mm ± 3,17 mm)	"(514,6 mm ± 3,17 mm)	"(514,6 mm ± 3,17 mm)	"(514,9 mm ± 3,17 mm)
Surface	On the outer surface should not be technology and surface defects.					

3.5 StayFix

StayFix	Small	Large
Total area, cm²	120 ± 2	120 ± 2
The area of the adhesive layer, cm2	88 ± 2	90 ± 2
catheter compatibility	For fixing the catheter outer diameter of 5F-14F (1,67 mm - 4.7 mm)	For fixing the catheter outside diameter 12F-22F (4,0 mm - 7.3 mm)

3.6 Revolution

Total area, cm2	111 ± 2
The area of the adhesive layer, cm2	95 ± 2
catheter compatibility	For fixing the catheter outer diameter of 6F-24F (2,0 mm - 8.0 mm)

3.7 Guide wire, diameter 0,018" (0,46 mm), length 60 cm

diameter	0,018 " (+0005"/ -0.0010") (0.46 mm, +0.0127 mm / mm -0.0254)
length	600 mm ± 25 mm
The diameter of the guard	4 mm ± 1 mm
The radius of the distal tip	3 mm ± 1 mm

3.8 Guide wire, diameter 0,038" (0,97 mm), length 150 cm, J-tip,

diameter	0.038 " (+0005"/ -0.0010") (0.97 mm, +0.0127 mm / mm -0.0254)
length	1500 mm ± 125 mm
The diameter of the guard	7 mm ± 1 mm
The radius of the distal tip	3 mm ± 1 mm

3.9 Introducer needle 21G

The diameter of the needle	21 G (0,81 mm)
Working length of the stylet	170 mm + 5 mm The stylet must be longer than the cannula, when the two components are assembled into a single structure.
The outer diameter of the stylet	0.5 mm ± 0,05 mm
Tip shape stylet	Triangular pointed tip for atraumatic placement device.
The shape of the cannula tip	At the distal end of the cannula tip should be tapered.

The outer diameter of the cannula	1.25 mm $\pm$ 0,05 mm
Surface	When carrying out a visual check on the external surface of the working length of the device no foreign substance or surface defects. At no point of the needle protrusions growths, hooks etc.
The protruding portion of the needle tip	In the assembled structure cannula does not extend beyond the tip of the stylet, or its length does not exceed the length of the stylet more than 2 mm.
Compatibility guidewire	guidewire of 0.018 inches (0.46 mm) with an installed cannula needs to pass through the entire structure in the Assembly.

3.10 Introduce:

-- dilator diameter 6F,

External diameter	The length of the conical narrowing	Length products
6F $\pm$ 1F (2,0 mm $\pm$ 0,3 mm)	6 mm $\pm$ 1 mm	250 mm $\pm$ 10 mm
The inner diameter of the distal (narrowed) end of dilator	For 6F diameter	0,041 " $\pm$ 0,002 " (1.04 mm $\pm$ 0,05 mm)
Compatibility guidewire	guidewire of 0.038 inches (0.97 mm) with an installed cannula needs to pass through the entire structure in the Assembly.	
The geometric shape of the distal end of the dilator	The presence of tapering, polished the edges. Edge should be sharp. Without cracks, scratches.	
Radiographic	Necessarily	

- Metal Stiffening Cannula,

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F	14 F
Cannula OD	0,81 mm +0,0127 mm / -0,0254 mm					
Cannula ID	0,56 mm +0,0127 mm / -0,0254 mm					
The length of the open part (mm + 1 / -2 mm)	238					
Allowable angle of deflection, without compromising the integrity	15 ⁰ $\pm$ 1 ⁰					
The shape of the cannula tip	At the distal end of the cannula tip should be tapered					
Surface	On the outer surface should not be technology and surface defects.					

3.11 MDD bag 600 ml volume

Length	265 mm $\pm$ 5 mm
Width	118 mm $\pm$ 5 mm
volume	600 ml $\pm$ 5 ml

Weight bag tube	46,9+ 2 g
Connecting dimensions ensuring tightness of the connection of the tube of the connecting bag to the catheter	Outer diameter 10 mm $\pm$ 0,5 mm Internal diameter 7 mm $\pm$ 0,5 mm The diameter of the central hole is 3 $\pm$ 0,5 mm
The connecting dimensions ensuring the tightness of the connection of the bag with the connecting pipe	Outer diameter 10 mm $\pm$ 0,5 mm Internal diameter 6 mm $\pm$ 0,3 mm
Tightness	All connectors must not allow fluid leakage
Graduating	The bag scale should be graduated by divisions. The unit of capacity should be indicated at the top of the bag. Graduation lines should have numerical values corresponding to the capacities (100 ml, 200 ml, 300 ml, 400 ml, 500 ml, 600 ml)
Diameter of the drain hole	9 mm $\pm$ 2 mm
Length of the drain cap, mm	27 mm $\pm$ 2 mm
Diameter of the drain cap, mm	16 mm $\pm$ 2 mm
The rate of filling the bag at a pressure of 1 mv.	Not less than 8 ml / s
The speed of emptying the bag at a pressure of 1 m.v.	Not less than 20 ml / s

3.12 Connecting tube for drainage bag

Diameter	6,0 mm $\pm$ 0,3 mm
tube length	615 mm $\pm$ 5 mm
The allowable deviation angle, without compromising the integrity of	900 $\pm$ 10

3.13 dilator

External diameter	The length of the conical narrowing	Weight	Length products
8F $\pm$ 1F (2,7 mm $\pm$ 0,3 mm)	8 mm $\pm$ 1 mm	1.55 g $\pm$ 0,50 g	233 mm $\pm$ 10 mm
10F $\pm$ 1F (3,4 mm $\pm$ 0,3 mm)	13 mm $\pm$ 1 mm	1.75 g $\pm$ 0,50 g	233 mm $\pm$ 10 mm
12F $\pm$ 1F (4,1 mm $\pm$ 0,3 mm)	13 mm $\pm$ 1 mm	2.55 g $\pm$ 0,50 g	233 mm $\pm$ 10 mm
14F $\pm$ 1F (4,7 mm $\pm$ 0,3 mm)	15 mm $\pm$ 1 mm	4.0 g $\pm$ 0,50 g	233 mm $\pm$ 10 mm
The inner diameter of the distal end of dilator	For diameters of 8F-12F		0,041 " $\pm$ 0,002 " (1.04 mm $\pm$ 0,05 mm)
	For diameters 14F		0,098 " $\pm$ 0,002 " (2,49 $\pm$ 0,05 mm)
Compatibility guidewire	guidewire of 0.038 inches (0.97 mm) with an installed cannula needs to pass through the entire structure in the Assembly.		
The geometric shape of the distal end of the dilator	The presence of tapering, polished the edges. Edge should be sharp. Without cracks, scratches.		
Radiographic	Necessarily		

3.14 belt

Length	935 mm $\pm$ 15 mm
Width	25 mm $\pm$ 5 mm

3.15 Tip type Luer-lok

Length	20 mm $\pm$ 2mm
External diameter	10,0 mm $\pm$ 2 mm
Inner diameter	3,0 mm $\pm$ 0,2 mm

3.16 FastBreak

length of the connector, type "male"	20,82 mm $\pm$ 0,13 mm
length of female connector	20,32 mm $\pm$ 0,13 mm
outer diameter of Luer coupling	6,85 mm $\pm$ 0,05mm, thread pitch - 1,4mm $\pm$ 0,025mm
inner diameter of Luer's standard plug connection,	7,23 $\pm$ 0,05 mm thread pitch - 1,4mm $\pm$ 0,025mm
tightness	the connector must be sealed

4. PACKAGING QUALIFICATION

For products "Merit Medical Systems, Inc." Acceptable standard warehousing, unless special conditions of storage, in accordance with the risk management reports are not listed on the packaging. The product is packaged and labeled as a sterile disposable article. The sterile barrier systems specified in the standard ISO 11607, are designed to prevent the ingress of microorganisms and ensures aseptic preservation of the product to the final application. Checking package integrity is carried out within the production process on the package.

The sets of the drainage catheter Resolve Locking Drainage Catheter consist of the following components: one (1) ReSolve Locking Drainage Catheter drainage catheter together with a pigtail loop straightener and repositioner, one (1) metal cannula, one (1) flexible cannula with catheter, diameter 8.5 F (2.8 mm), 10F (3.3 mm), 12F (4.0 mm), 14F (4.6 mm)), (1) one trocar, (1) one the end cap is placed in a shockproof polystyrene tray with a PET lid to support the products, which is then placed in a bag of Tyvek® material, sealed from oschyu validated thermal welding process intended to cover and protect the product during sterilization, handling, transportation, handling and storage.

Also, Resolve Locking Drainage Catheter catheter kits can be shipped with the FastBreak connector. The FastBreak connector is packaged in an individual package in the form of a packet of Tyvek® material, sealed with a validated thermal welding process.

Individual packings of the Resolve Locking Drainage Catheter drainage catheter and the FastBreak connector are connected by an external fixing package with a printed general label on which the kit is presented.

The kit of the Resolve Locking Drainage Catheter Plus drainage catheter consists of the following components: one (1) ReSolve Locking Drainage Catheter drainage catheter with a pigtail loop straightener and repositioner, one (1) metal cannula, one (1) flexible cannula sets with catheter, diameter 8.5 F (2.8 mm), 10F (3.3 mm), 12F (4.0 mm), 14F (4.6 mm)), (1) one trocar, (1) one end cap is placed in an impact-resistant polystyrene tray with a PET lid to support

the products, which is then placed in a bag of Tyvek® material sealed with a validated thermal welding process and a separately packaged fixation device (and / or a separately packaged FastBreak connector). The fixing device (StayFix or Revolution) and / or FastBreak connector are packaged in individual packages in the form of a packet of Tyvek® material, sealed with a validated thermal welding process.

Packages made of Tyvek® material sealed with a validated thermal welding process are designed to cover and protect the product during sterilization, loading, transportation, manipulation and storage. Individual packaging of the drainage catheter and the fixation device (and / or FastBreak connector) are connected by an external fixing package with a common label.

Configuration ReSolve Locking Drainage Catheters NV contains the following components:

- One (1) drainage catheter with locking pigtail and hydrophilic coating
 - One (1) metal stiffening cannula
 - One (1) flexible stiffening cannula (not included with the 6.5F and 7.5F)
 - One (1) Repositioning Tool
 - One (1) Pigtail Straightener
 - One (1) Dead end cap
 - One (1) MAK NV system consisting of coaxial introducer with dilator, stiffening cannula, 21 GA introducer needle and 0.018 and 0.038 guide wires (ref. TEC0018) (not shown in figure below)
 - One (1) 600ml MDD with PVC bag (ref. TEC0044) (not shown in figure below)
 - One (1) Dilator (8F, 10F, 12F or 14F) (not shown in figure below)
 - One (1) StayFix or Revolution catheter fixation device
- * For a catheter with a diameter of 6.5 F (2.1 mm); A 7.5 F (2.4 mm) dilator is not required.

All listed products are placed in a shockproof polystyrene tray with a PET lid to support products, which are then placed in a bag of Tyvek® material sealed with a validated thermal welding process. Also, these sets can be completed with a fixing device (StayFix or Revolution) and / or a FastBreak connector.

The fixing devices (StayFix or Revolution) and the FastBreak connector are packed in individual packages in the form of a packet of Tyvek® material sealed with a validated thermal welding process.

Packages made of Tyvek® material sealed with a validated thermal welding process are designed to cover and protect the product during sterilization, loading, transportation, manipulation and storage.

Individual packaging of the drainage catheter and the fixation device (and / or FastBreak connector) are connected by an external fixing package with a common label applied.

Instructions for use are included in the group packaging of the sets. Sets of catheters in a group box are placed in a transport box.

ReSolve's packaging qualification was conducted to meet the requirements of the following standards:

ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems.

ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation Requirements for forming, sealing and assembly processes

5. LABEL

Product label formats have been developed and approved for unit product, sterile packages and transport packaging. The approved labeling is located in the product Device Master Record.

Merit utilizes the following standards for labeling content and symbol generation in order to meet current requirements of the MDD:

EN 980: Symbols for use in the Labeling of Medical Devices.

EN 1041: Information Supplied by the Manufacturer.

ISO 15223 - 1: Medical Devices - Symbols to be used with medical device labels, labeling and information to be supplied Part 1: General Requirements.

ISO 15223 - 2: Medical Devices - Symbols to be used with medical devices, labeling, and information to be supplied Part 2: Symbol Development, Selection and Validation.

- FDA Guidance on Medical Device Patient Labeling; Final Guidance for Industry and FDA Reviews (4/19/2001).

Description of symbols labeling shown in the table.

	Symbol.	Explanation of symbols
1.		EC official
2.		manufacturer
3.	R _x ONLY	It sold only to doctors and hospitals
4.		For single use only. Do not reuse
5.		catalog number
6.		shelf life
7.		lot number
8.		sterilized ethyleneoxide
9.		Sterilized by gamma radiation
10.		attention
11.		CONDITIONAL compatibility with MR
12.		CE Mark
13.		Pyrogen-free
14.		Compatibility with the conductor

15.		Not for individual sale
16.		Protect from moisture
17.		This side up
18.		Fragile, handle with care
19.		Do not use if package is damaged

6. TRANSPORTATION, STORAGE AND HANDLING

The shelf life of sterile - 36 months from the date of sterilization.

Tools that are exempt from transport packaging, must be kept at an ambient temperature of 10° C to 40° C and a relative humidity of 30% to 75% non-condensing.

Transporting tools made according to the conditions of storage

7. UTILIZATION

Use a set of catheters disposed of as hazardous waste epidemiologically class B (presence of contact with blood and mucosal tissues of the patient).

After the expiry date the unused sets of catheters not in contact with blood and mucous membranes should be disposed of as hazardous waste epidemiologically, close in composition to the solid waste. Dispose of used conductors shall be in accordance with the laws and regulations on the disposal of products made in the country of use of the medical device.

8. WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life ReSolve lockable catheter sets 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided

without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

9. COMPLAINT

The organization, which takes in the Russian Federation claims from consumers on the quality of the product:

Limited Liability Company "Merit Technologies" (LLC "Merit Technologies"), BIN 1117746679533, TIN 7705960130

Legal address: 115054, Moscow, Bakhrushina Str. 32, p. 1, floor/space 4/25A

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