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Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: 
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
REF 02R5101
G72100R01
S2R510

Created November 2016.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

Alinity c Transferrin Calibrator Kit (also referred to as Transf Cals)

INTENDED USE

For use in the calibration of the Alinity c Transferrin assay on the Alinity c analyzer.

For additional information, refer to the Alinity c Transferrin reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CAL 1 through CAL 5 are prepared from transferrin fractions, human IgA, IgG, IgM, C3, C4 and haptoglobin in human serum.

Preservative: sodium azide.

Lot-specific calibrator values are listed in the Alinity c Transferrin Calibrator Kit value sheet, packaged with the calibrator. Verify that the lot number listed on each calibrator carton agrees with the lot number printed on the value sheet. The last two digits of the lot numbers can vary.

Calibrator	Quantity
CAL 1-5	1 x 1.0 mL

STANDARDIZATION

Each bottle is prepared from transferrin fractions in human serum. The calibrator value assignment for each calibrator level is determined using a reference method calibrated against ERM-DA470.

PRECAUTIONS

- IVD
- For In Vitro Diagnostic Use
- Rx ONLY

Safety Precautions

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that this product and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in Levels 1-5 of the Transferrin Calibrators is nonreactive for HBsAg, HIV-1 RNA, HCV RNA, HBV RNA, anti-HIV-1/HIV-2, anti-HCV, anti-HTLV-I/II and anti-HBc.

The following warnings and precautions apply to: CAL 1-5

Contains sodium azide.

EUH032 Contact with acids liberates very toxic gas.

P501 Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	30 days after opening	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- Calibrator values may be configured using e-files accessed and imported from www.abbottdiagnostics.com, or from Abbott Mail.
- Verify that the correct calibrator values have been entered into the calibration file.
- Place an appropriate amount of the calibrator into each sample cup in the assigned position.
- Calibration is performed by running a water blank and Transferrin Calibrator set.
- Water for the blank is provided by the instrument.
- Verify control results are within acceptable limits before reporting patient results.
- For information on configuring calibrator data, refer to the Alinity ci-series Operations Manual, Section 2.
- For instructions on ordering and loading calibrators on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, if calibration does not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria, or if controls do not meet the appropriate criteria.

LIMITATIONS OF THE PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments and reagents, storage of product as directed, and good laboratory technique.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
ALTERNATE UNITS	Alternate Units
ASSAY NAME	Assay Name
ASSAY NO.	Assay Number
CAL 1	Calibrator (1, 2, 3, 4, 5)
CAL 1-5	Calibrators 1 through 5
CN	Control Number
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
DEFAULT UNITS	Default Units
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF JAPAN	Product of Japan
REF	List Number
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).
UNIT	Unit

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Created November 2016.
©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity c Transferrin Calibrator Kit (for Russia)

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

Medical device must be stored at a temperature of 2 °C to 8 °C in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 12 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required

REPAIRS

Not required

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001, ISO 13485, ISO 13485 (including CMDCAS)

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

INDICATIONS AND CONTRA-INDICATIONS FOR USE

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