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Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Nora Fitzgerald, Regulatory Affairs, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed:

Nora Fitzgerald

Nora Fitzgerald
Regulatory Affairs
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co Longford
Ireland

Date: 23 / Nov / 2018

ABBOTT DIAGNOSTICS DIVISION LONGFORD

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Mary Kelly 27/11/2018 Authorised Official

Alinity i

25-OH Vitamin D Controls

en

Vit D 25OH

REF 08P4510

G71165R01

C8P450

Created September 2016.

NAME

Alinity i 25-OH Vitamin D Controls (also referred to as Vit D 25OH Ctrls)

INTENDED USE

The Alinity i 25-OH Vitamin D Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer when used for the quantitative determination of 25-hydroxyvitamin D (25-OH vitamin D) in human serum and plasma. For additional information, refer to the Alinity i 25-OH Vitamin D reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CONTROL L, **CONTROL M**, and **CONTROL H** contain 25-OH vitamin D in PBS buffer with human serum. Preservatives: ProClin 950 and sodium azide.

The controls are at the following ranges and target concentrations:

Control	Quantity	25-OH vitamin D		25-OH vitamin D	
		CONC (ng/mL)	RANGE (ng/mL)	CONC (nmol/L)	RANGE (nmol/L)
CONTROL L	1 x 8.0 mL	20.0	14.0 - 26.0	50.0	35.0 - 65.0
CONTROL M	1 x 8.0 mL	40.0	28.0 - 52.0	100.0	70.0 - 130.0
CONTROL H	1 x 8.0 mL	75.0	52.5 - 97.5	187.5	131.3 - 243.8

NOTE: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Control lot
- Reagent lot
- Calibrator lot
- Instrument

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in the controls is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to:

CONTROL L, **CONTROL M**, and **CONTROL H**

WARNING Contains methylisothiazolone and sodium azide.

H317 May cause an allergic skin reaction.

EUH032 Contact with acids liberates very toxic gas.

Prevention

P261 Avoid breathing mist / vapors / spray.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves / protective clothing / eye protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of water.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

P362+P364 Take off contaminated clothing and wash it before reuse.

Disposal

P501 Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Store in an upright position. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 4 drops of the low control, 4 drops of the medium control, and 4 drops of the high control into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

Abbott

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.


BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Control Number
	Concentration
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Control (Low, Medium, High)
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Ireland
	Range
	List Number

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Created September 2016.

©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i 25-OH Vitamin D Controls for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device are (2-8°C).

Medical device must be stored in freezers to ensure regulated temperature conditions with a daily temperature record.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. The maximum shelf life of the medical device is 15 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

EN ISO 9001
ISO 13485
EN ISO 14971
EN 13641
EN 13612
EN ISO 15223-1
EN 980
EN ISO 18113-1
EN ISO 18113-2
EN ISO 18113-3
EN ISO 23640
EN ISO 20776-1
EN 13975
EN ISO 17511
ISO 80000-1
EN ISO 15193
EN ISO 15194

POTENTIAL USER

Professional technicians.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs.

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage.

TECHNICAL CHARACTERISTICS

- Controls: pH-value at 20 °C: 7,2

COMPATABILITY WITH OTHER MEDICAL DEVICES

No available data.

MATERIALS OF ANIMAL/HUMAN ORIGIN

Materials of animal/human origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbott diagnostics.com or contact your local representative.

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