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
As a registered Notary Public for the State of Minnesota,  
I hereby certify that this is a true and accurate copy of the  
original document. The original is filed at Medtronic Inc.

2012/OCT/10 at 9:12 a.m. Doc number: MB50761A001 Rev. 1B [Multi26]  
Printspec full-b - 8.5 x 11 inches

**Medtronic**

## Cardioblate<sup>®</sup> CryoFlex<sup>®</sup>

Surgical Ablation Console  
Console d'ablation chirurgicale  
Konsole für die chirurgische Ablation  
Consola de ablación quirúrgica  
Chirurgische ablatieconsole  
Console per ablazione chirurgica  
Konsoll for kirurgisk ablasjon  
Kirurginen ablaatiokonsoli  
Kirurgisk ablationskonsol  
Κονσόλα χειρουργικής κατάρυσης  
Kirurgisk ablationskonsol  
Konsola do ablacji chirurgicznej  
Consola de ablação cirúrgica  
Cerrahi Ablasyon Konsolu  
Konzole pro chirurgickou ablaci  
Sebészeti ablációs konzol  
Konzola na chirurgickú abláciu  
Kirurgilise ablaatsiooni konsool  
Конзола за хирургична аблация  
Konzola za kiruršku ablaciju  
Kirurgiskās ablācijas konsole  
Chirurginės ablacijos pultas  
Consolă pentru ablație chirurgicală  
Консоль для хирургической аблации  
Konzola za hiruršku ablaciju  
Konzola za kirurško ablaciju

Operator's Manual • Manuel de l'utilisateur  
Gebrauchsanweisung • Manual del usuario  
Gebruikershandleiding • Manuale dell'operatore  
Brukerhåndbok • Käyttöopas  
Användarmanual • Εγχειρίδιο χειριστή  
Brugsanvisning • Podręcznik użytkownika  
Manual do operador • Teknik Kullanicı El Kitabı  
Uživatelská příručka • Kezelési kézikönyv  
Priručka operatora • Kasutamisinjuhend  
Наръчник на оператора • Příručník za rukovateľa  
Operatora rokasprēmata • Naudojimo instrukcija  
Manual de utilizare • Руководство оператора  
Priručnik za rukovaoca • Priročnik za operaterje

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Cardioblate® and CryoFlex® are trademarks of Medtronic, Inc.  
Cardioblate® et CryoFlex® sont des marques commerciales de Medtronic, Inc.  
Cardioblate® und CryoFlex® sind Marken von Medtronic, Inc.  
Cardioblate® y CryoFlex® son marcas comerciales de Medtronic, Inc.  
Cardioblate® en CryoFlex® zijn handelsmerken van Medtronic, Inc.  
Cardioblate® e CryoFlex® sono marchi di fabbrica della Medtronic, Inc.  
Cardioblate® og CryoFlex® er varemerker for Medtronic, Inc.  
Cardioblate® ja CryoFlex® ovat Medtronic, Inc:n tavaramerkkejä.  
Cardioblate® och CryoFlex® är varumärken som tillhör Medtronic, Inc.  
Οι ονομασίες Cardioblate® και CryoFlex® είναι εμπορικά σήματα της Medtronic, Inc.  
Cardioblate® og CryoFlex® er registrerede varemærker tilhørende Medtronic, Inc.  
Cardioblate® i CryoFlex® są znakami towarowymi firmy Medtronic, Inc.  
Cardioblate® e CryoFlex® são marcas comerciais da Medtronic, Inc.  
Cardioblate® ve CryoFlex® Medtronic, Inc. şirketinin ticari markalarıdır.  
Cardioblate® a CryoFlex® jsou ochranné známky společnosti Medtronic, Inc.  
A Cardioblate® és a CryoFlex® a Medtronic, Inc. védjegye.  
Cardioblate® a CryoFlex® su ochranné známky spoločnosti Medtronic, Inc.  
Cardioblate® ja CryoFlex® on ettevõtte Medtronic, Inc. kaubamargid.  
Cardioblate® и CryoFlex® са търговски марки на Medtronic, Inc.  
Cardioblate® i CryoFlex® zaštitni su znakovi tvrtke Medtronic, Inc.  
Cardioblate® un CryoFlex® ir Medtronic, Inc. reģistrētas preču zīmes.  
„Cardioblate®” ir „CryoFlex®” yra „Medtronic, Inc.” prekių ženklai.  
Cardioblate® și CryoFlex® sunt mărci comerciale ale Medtronic, Inc.  
Cardioblate® и CryoFlex® являются товарными знаками Medtronic, Inc.  
Cardioblate® i CryoFlex® su žigovi preduzeca Medtronic, Inc.  
Cardioblate® in CryoFlex® sta blagovni znamki družbe Medtronic, Inc.

Explanation of Symbols / Explication des symboles / Erläuterung der Symbole / Explicación de los símbolos / Verklaring van de symbolen / Spiegazione dei simboli / Symbolforklaring / Symbolien merkitykset / Förklaring av symboler / Επεξήγηση των συμβόλων / Symbolforklaring / Objaśnienie symboli / Explicação dos símbolos / Sembollerin Açıklaması / Vysvětlení symbolů / Jelmagyarázás / Vysvetlivky k symbolom / Sümbolite seletus / Объяснение на символите / Objašnjenja simbola / Simbolu skaidrojums / Simbolių paaiškinimas / Explicarea simbolurilor / Объяснение символов / Objašnjenje simbola / Razlaga simbolov

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Table of Contents

Table of Contents

|                                                                             |           |
|-----------------------------------------------------------------------------|-----------|
| <b>1 Introduction</b>                                                       | <b>2</b>  |
| Device Description                                                          | 2         |
| Indications                                                                 | 2         |
| Contraindications                                                           | 2         |
| Warnings                                                                    | 2         |
| Precautions                                                                 | 3         |
| Adverse Effects                                                             | 3         |
| <b>2 System Setup</b>                                                       | <b>4</b>  |
| Console                                                                     | 4         |
| Control Panel                                                               | 5         |
| Installation                                                                | 5         |
| <b>3 Pre-Op Setup</b>                                                       | <b>7</b>  |
| <b>4 Performing Cryoablation</b>                                            | <b>11</b> |
| <b>5 Post-Op Take Down</b>                                                  | <b>13</b> |
| <b>6 Maintenance</b>                                                        | <b>14</b> |
| Replacing Fuses                                                             | 14        |
| Refrigerant                                                                 | 14        |
| Refrigerant Usage                                                           | 14        |
| Shipping the Refrigerant Tank                                               | 15        |
| Changing the Refrigerant Tank                                               | 15        |
| Cleaning the Cardioblade CryoFlex Surgical Ablation Console and Accessories | 17        |
| Internal Maintenance                                                        | 18        |
| <b>7 Troubleshooting</b>                                                    | <b>19</b> |
| <b>8 Technical Specifications</b>                                           | <b>20</b> |
| Surgical Cryoablation System Specifications                                 | 20        |
| How Supplied                                                                | 20        |
| Electromagnetic Emissions and Immunity Declarations                         | 20        |
| <b>9 Warranties</b>                                                         | <b>22</b> |
| Limited Warranty for the Cardioblade® CryoFlex® Surgical Ablation Console   | 22        |
| Disclaimer of Warranty for the Accessories                                  | 22        |

## Chapter 1

# Introduction

1

## Device Description

The Cardioblate CryoFlex Surgical Ablation Console is a cryoablative surgical device that is used in the cryosurgical treatment of cardiac arrhythmias. The system utilizes an argon-based cryogen for fast, controlled freezing.

The system consists of the Cardioblate CryoFlex Surgical Ablation Console, connection components and accessories, and disposable cryoablation probes.

### Cryoablation

During a cryoablation procedure, pressurized argon gas flows from the gas tank through the regulator, which maintains the gas flowing through the control panel at a constant pressure. The pressurized argon gas inside the control panel flows through a valve, which controls gas delivery to the cryoablation probe through the high-pressure connector on the front of the control panel. The console is equipped with a vent-valve that relieves the pressure for storage between uses. As the pressurized argon gas reaches the ablation segment of the probe, it passes through microholes that cause the pressure to drop dramatically. This pressure drop causes the temperature of the argon gas to decrease and results in freezing temperatures at the ablation segment of the probe.

## Indications

The Cardioblate CryoFlex Surgical Ablation System is intended for minimally invasive cardiac surgical procedures, including the treatment of cardiac arrhythmias. The Cardioblate CryoFlex 7 cm, 10 cm, and 10-S probes plus the Cardioblate CryoFlex Clamp and Cardioblate CryoFlex Surgical Ablation Console freeze target tissue and block the electrical conduction pathways by creating an inflammatory response and cryonecrosis.

## Contraindications

The Cardioblate CryoFlex Surgical Ablation Console, Clamps, and Probes are not designed for use inside a beating heart.

## Warnings

- Refer to product literature included for each disposable device, including cryoablation probes and clamps, for further instructions and warnings specific to those devices.
- Cryoablation probes contain pressurized refrigerant during operation. Discontinue treatment immediately if a breach in the probe is suspected, as this may result in release of pressurized gas and injury to the patient or the user.
- Avoid contact between the cryoablation probe and the phrenic nerve to avoid injury.
- If an unanticipated event occurs, the injection of refrigerant can be stopped at any time by pushing the Stop Freeze button on the control panel.
- If you are unable to stop the injection, immediately turn the system off by using the power switch located on the side of the control panel or by disconnecting the power cord.
- Do not attempt to disconnect the probe or gas line during operation in the freeze mode as the sudden release of pressure may cause the probe or gas line to recoil, which may injure the operator. The system must be vented prior to disconnecting the probe and gas line.
- If the injection of refrigerant spontaneously stops before the preset time on the control panel has elapsed, replace the probe.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- The power cord may be used to disconnect the device from electrical power. Do not position the control panel so that it is difficult to remove the power cord.
- Do not modify this equipment.
- To avoid the risk of electric shock, do not touch the patient when installing or handling the Pressure Sensor Cable.
- Use caution when transporting the CryoFlex Console over rough or uneven surfaces with changes in height greater than 10 mm.

## Introduction

### Precautions

- Only appropriately trained personnel in a fully equipped facility should perform cardiac cryoablation procedures.
- Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- This device is designed for use in an electromagnetic environment and may generate or be susceptible to electromagnetic interference. Refer to the electromagnetic compatibility section for the recommended environment for electromagnetic emissions and immunity.
- Ensure that equipment used with the system is electrically safe and does not expose the patient to hazards.
- Ensure that all pressure and regulator lines are securely connected.
- Always turn off the main gas cylinder valve prior to venting the inlet pressure line to ensure that excess pressure is safely vented.
- Perform a test injection before placing a probe on tissue to be ablated.
- Follow standard guidelines for the safe handling and storage of high-pressure gas cylinders.
- To prevent unwanted movement of the Surgical Ablation Console on an incline, position the locking casters on the low side of the incline.
- Do not disconnect the probe from the control panel until the procedure is finished, unless the probe needs to be replaced.
- Discard all used probes and clamps in accordance with hospital procedures.

### Adverse Effects

Potential adverse effects with this device are similar to other cardiac surgery procedures and may include the following:

- bleeding
- reoperation
- extension of extracorporeal bypass
- heart rhythm disturbances (atrial and/or ventricular)
- pericardial effusion
- pericarditis
- cardiac tamponade
- pleural effusion
- mediastinitis
- conduction disturbances (SA/AV node)
- acute ischemic myocardial event
- thrombus formation
- low cardiac output
- stroke
- renal, gastrointestinal, or respiratory complications
- sepsis
- adjacent structural damage
- death

Cryoablation involving coronary vessels has been associated with subsequent clinically significant arterial stenosis. It is unknown whether cryoablation with the Cardioblate CryoFlex Surgical Ablation Probe will have such an effect, but as in all such procedures, care should be taken to minimize unnecessary contact with coronary vessels during cryoablation.

## Chapter 2

# System Setup

2

## Console


Figure 1.

Surgical cryoablation probes and clamps are used to apply freezing temperatures to the heart. Refer to product literature included with cryoablation probes for further information. Products that are compatible with the Cardioblade CryoFlex Surgical Ablation Console include:

- Cardioblade CryoFlex Surgical Ablation Probe, 7 cm
- Cardioblade CryoFlex Surgical Ablation Probe, 10 cm
- Cardioblade CryoFlex 10-S Surgical Ablation Probe, 10 cm
- Cardioblade CryoFlex Clamp and Surgical Ablation Probe, 10 cm

**Note:** The probes and clamps are Type CF applied parts and are designed for single patient use. The control panel can be locked to the tank carrier or used separately. A 3 m (10 ft) power cord is shipped separately.

The tank carrier holds the tank and allows for safe transport and storage of the argon tank. It consists of two pieces: the front white section holds the control panel, and the back blue section holds the argon gas tank. The two pieces should be locked together during use and transport.

The regulator attaches the gas tank to the control panel and regulates the pressure of argon gas supplied to the control panel. The 2.4 m (8 ft) gas hose (shipped separately) connects the gas tank and regulator to the control panel. The 2 m (6 ft) pressure sensor cable is shipped with the regulator.

## System Setup

### Control Panel


Figure 2. Control Panel

1. **Display Screen:** Shows freeze time duration, temperature of ablation segment, and user messages.
2. **Freeze Timer Window:** Use the black arrow keys to change the duration (in seconds) of the freeze. The timer on the Display Screen counts up from zero.
3. **Target Temperature Window:** Use the red and blue arrow keys to increase or decrease the target temperature.
4. **Start Freeze Button:** Press and hold for 2 seconds to start the injection of refrigerant into the probe. The green light above the button remains lit while refrigerant is flowing.
5. **Message Indicator:** Lights up when there is a message on the Display Screen for the user to read.
6. **Stop Freeze Button:** Press at any time to stop the injection of refrigerant.
7. **Vent Button:** Press to depressurize the gas line.

### Installation

**Caution:** Follow standard guidelines for the safe handling and storage of high-pressure gas cylinders.

**Caution:** The Cardioblade CryoFlex Surgical Ablation Console meets the requirements of Underwriters Laboratories (UL). It is the user's responsibility after installation to verify and ensure that the system meets the applicable electrical safety requirements with respect to electric shock, fire, and mechanical hazards only in accordance with UL 60601-1, CAN/CSA C22.2 No. 601.1-1, ANSI / AAMI ES60601-1, and CAN/CSA C22.2 No. 60601-1.

#### Unpacking

The system and its accessories are shipped in boxes and packed for protection. Upon receipt, open all the boxes to ensure that all equipment and accessories are included and that there is no evident shipping damage. Manage disposal in accordance with local statutory regulations.

Gather the following items to assemble the Console:

- Control Panel, shipped with the Operator's Manual
- Tank Carrier (blue and white sections)
- Power Cord appropriate for geography
- Gas Tank Regulator (supplied with wrench)
- Pressure Sensor Cable (provided with regulator or available separately)
- Argon Gas Tank (not supplied by Medtronic; obtain from hospital supplier)
- Gas Hose

#### Verify Voltage Setting

This control panel is capable of running 120V or 230V AC power. Inspect the fuse block next to the power switch to verify that the setting is correct for your location. If the setting is not correct, refer to "Replacing Fuses" in Chapter 6 for instructions on how to change the voltage setting.


## Chapter 2

Figure 3. Fuse Block in 120VAC Configuration


### Selecting the Language for the Control Panel

The control panel displays English by default. To change the displayed language, power off the control panel and set the DIP switches to the desired language as shown in Figure 5.

1. DIP Switches

Figure 4.

1. English
2. French
3. Italian
4. German
5. Danish
6. Spanish
7. Swedish
8. Portuguese
9. Dutch

Figure 5.

## Pre-Op Setup

### Pre-Op Setup

3

**Caution:** Only appropriately trained personnel in a fully equipped facility should perform cardiac cryoablation procedures.

**Caution:** Ensure that equipment used with the console is electrically safe and does not expose the patient to hazards.

1. Position the Surgical Ablation Console in the desired location. Push down the locks on the front wheels.

Figure 6. Front Wheels

**Caution:** To prevent unwanted movement of the Surgical Ablation Console on an incline, position the locking casters on the low side of the incline.

2. Remove the control panel from the console by grasping the handles and lifting straight up until the control panel clears the hooks on the tank carrier, then tilt the control panel away from the tank carrier and toward the user. Set the control panel on a table. (The control panel must be removed from the tank carrier to access the back connections.)

**Note:** The control panel can be operated while mounted on the tank carrier or without a tank carrier.

1. Handle


Figure 7. Control Panel Handle

3. Connect the regulator to the gas tank by attaching the nut on the regulator (Figure 8, number 7) to the port on the tank. Use the wrench to tighten the nut.

**Note:** The nut on the North American regulator turns counterclockwise to tighten.

**Caution:** Ensure that all pressure and regulator lines are securely connected.

4. To connect the gas hose to the regulator, pull back on the outer collar of the gas hose to expose the inner collar. Place the inner collar over the gas hose stem (Figure 8, number 5) on the regulator and release the outer collar. The connector locks into place.

1. Regulator
2. Pressure Sensor Cable
3. Pressure Sensor Stem
4. Gas Hose
5. Gas Hose Stem
6. Auxiliary Vent Valve
7. Nut

Figure 8. Regulator and Connections

5. Connect the other end of the gas hose to the metal gas inlet port on the back of the control panel (Figure 9, number 2). (Both ends of the gas hose are the same.)

### Chapter 3

6. The pressure sensor cable comes attached to the regulator. If it is not attached, line up the 4-pin connector to the pressure sensor stem (Figure 8, number 3) on the top of regulator and rotate the collar to tighten.
7. To connect the pressure sensor cable to the control panel, push the white plastic connector on the end of the pressure sensor cable onto the pressure transducer connection (Figure 9, number 3). It will click into place.


Figure 9. Control Panel Connections

8. Place the control panel back onto the tank carrier or position it on a table in the desired location.
9. Connect the power cable to the control panel, plug the power cord into an outlet, and turn on the power switch.  
**Note:** An optional user-provided equipotentialization conductor may be attached to the equipotentialization connector. The equipotentialization connector is used to connect the ground of the control panel to the ground of other instruments in a clinical setting.


Figure 10. Power Connections

10. Open the tank valve to begin the flow of gas and begin the self-test.


Figure 11.

- Note:** Before opening the gas tank, verify that the auxiliary vent valve on the regulator is closed by turning the handle clockwise (Figure 8). If the auxiliary vent valve is not fully closed, gas will escape when the gas tank is opened.
11. Open the gas tank by turning the knob on the top of the tank counterclockwise.
  12. The system performs a self-test then asks the user to attach the probe.  
**Note:** The system will not perform the self-test unless the gas tank is open.

## Pre-Op Setup

Figure 12.

13. To connect the probe to the control panel, orient the probe connector with the green arrow up (Figure 14) and the red stop sign down.
14. Push the probe connector into the control panel connector until the grips snap into place. Hold the connector as shown in Figure 13. Do not squeeze the grips.


Figure 13. Attach Probe

**Note:** Do not use excessive force to connect the probe. If resistance is felt, manually vent the system as described in Chapter 6, Maintenance, and try connecting the probe again.

**Note:** Do not press the top and bottom grips while inserting the probe connector into the control panel. Only one hand is required to make the connection.

15. Gently tug on the probe connector to ensure it is securely connected.
16. Pull out the single-use black pin on the left side of the probe connector and discard.

Figure 14. Remove Pin

### Chapter 3

**Caution:** Do not pull the pin until the probe is securely connected.

**Caution:** Do not disconnect the probe from the control panel until the procedure is finished, unless the probe needs to be replaced, as the probe will not function if reconnected.

17. Once the probe is attached, the control panel will display the time (defaults to zero), the temperature in Celsius (starts at ambient temperature), and an estimated number of minutes of cryoablation that the gas remaining in the tank will permit.

**Note:** Change the tank if there is not enough gas to complete at least one full procedure. (See Changing the Refrigerant Tank in Chapter 6 for details.)

## Performing Cryoablation

### Performing Cryoablation

4

**Caution:** Always perform a test injection before placing the probe on tissue to be ablated.

1. Perform a test injection with the probe in the air. Press the Start Freeze button for at least 2 seconds to start the injection of refrigerant.
2. The green light above the Start Freeze button lights up to indicate that refrigerant is flowing.


Figure 15. Start Freeze Button

3. If performing a test injection with the probe alone, the temperature reading on the control panel should drop to  $-110^{\circ}\text{C}$  within 20 seconds if the equipment is operating properly. At the end of the default time (60 seconds), the injection stops automatically and the temperature rises.

**Warning:** If an unanticipated event occurs, the injection of refrigerant can be stopped at any time by pushing the Stop Freeze button on the console control panel.

**Warning:** If you are unable to stop the injection, immediately turn the system off by using the power switch located on the side of the control panel or by disconnecting the power cord.

**Warning:** Do not attempt to disconnect the probe or gas line during operation in the freeze mode. Always vent the system at the end of the procedure before removing the gas line.

**Note:** If the probe needs to be replaced during a procedure, it is not necessary to vent the system. Refer to Step 6 in Chapter 5 for instructions on removing the probe.

4. Shape the probe and adjust the length of the freeze zone (if using the probe alone) by moving the white insulating sleeve up or down.

Figure 16. Surgical Ablation Probe

For instructions on use of specific probes or clamp, refer to the Instructions for Use supplied with those devices.

5. Select the target time and target temperature for cryoablation.
  - a. The default time for cryoablation is 1 minute. Press the black arrow keys to increase or decrease the ablation time. The default time can be adjusted before or during an injection.

#### Chapter 4

Figure 17.

- b. The default temperature for cryoablation is set to full power (coldest possible temperature), which is displayed as a dashed line. The lowest targeted temperature that will display in the Target Temperature window is  $-95^{\circ}\text{C}$ . Press the red and blue arrow keys to increase or decrease the targeted temperature.

Figure 18.

6. Position the ablation segment of the probe onto the tissue to be ablated, ensuring all adjacent tissue is protected to avoid inadvertent freezing.  
**Warning:** Refer to product literature included for each disposable device, including cryoablation probes and clamps, for further instructions and warnings specific to those devices.
7. Press and hold the Start Freeze button for 2 seconds to start cryoablation.  
**Note:** The timer on the control panel does not start until the temperature on the control panel drops below  $-40^{\circ}\text{C}$ .  
**Warning:** Probes contain pressurized refrigerant during operation. Discontinue treatment immediately if a breach in the probe is suspected, as this may result in release of pressurized gas and possible injury to the patient.  
**Warning:** If the injection of refrigerant spontaneously stops before the default time on the control panel has elapsed, replace the probe.  
**Warning:** Do not attempt to disconnect the probe or gas line during operation in the freeze mode. Always vent the system at the end of the procedure before removing the gas line.
8. Freezing stops at the end of the preset time unless it is terminated earlier by pressing the Stop Freeze button.

Figure 19. Stop Freeze Button

**Warning:** If an unanticipated event occurs, the injection of refrigerant can be stopped at any time by pushing the Stop Freeze button on the control panel.

**Warning:** If you are unable to stop the injection, immediately turn the system off by using the power switch located on the side of the control panel or by disconnecting the power cord.

9. The system automatically vents after each injection. Venting is signaled by a loud hissing noise.

## Post-Op Take Down

### Post-Op Take Down

5

1. Close the gas tank by turning the knob on the top of the tank clockwise.  
**Important:** The tank must always be closed before venting the system to avoid venting all the gas in the tank.
2. Press the Vent button on the control panel.
3. The display screen will read "Close Gas Tank Valve and Press Vent Button." Ensure that the gas tank is closed and press the Vent button again.


Figure 20.

4. The light above the button turns on to indicate that venting has been selected. The display screen will read "System Venting." A loud hissing noise will be heard as the gas is vented from the gas line between the tank and the regulator.
5. The auxiliary vent valve on the regulator assembly can also be used to vent any remaining gas in the gas line. Open the auxiliary vent valve by turning the handle counterclockwise.  
**Caution:** Always turn off the main gas cylinder valve prior to venting the inlet pressure line to ensure that excess pressure is safely vented.
6. Disconnect the probe from the system.

Figure 21.

- a. With one hand, press the top and bottom grips of the probe connector at the same time.
  - b. With the other hand, push the connector flange ring in slightly.
  - c. Pull back to disconnect the probe from the control panel connector.
  - d. Dispose of single-use accessories and probes according to hospital regulations.
7. Turn off the power at the switch on the side of the control panel.
  8. Optional: Disconnect the gas hose from the regulator and control panel.
  9. Optional: Disconnect the pressure sensor cable from the control panel.
  10. Carefully store the reusable accessories.

## Chapter 6

# Maintenance


6

This section describes how to replace the fuses, change the refrigerant tank, and clean the outer surfaces of the console and its nonsterile accessories.

Caution: Follow standard guidelines for the safe handling and storage of high-pressure gas cylinders.

## Replacing Fuses

The fuse block is part of the power switch. Use a standard screwdriver to pop out the fuse block insert. Replace blown fuses using the same amperage rating and physical geometry as the existing fuses. Reinstall the fuse block insert so that the arrow for the correct voltage lines up with the arrow on the housing. The voltage setting must match the supplied voltage for your location. Refer to Chapter 8 for correct fuse type.

1. Voltage Setting Arrows
2. Arrow on Housing

Figure 22.

## Refrigerant

The grade 4.8 argon gas refrigerant can be obtained from local gas suppliers. It must be pressurized as follows and of a purity of at least 99.998%:

| Europe<br>(50L tank)   | USA<br>(type BX, approximately 50L) |
|------------------------|-------------------------------------|
| 31,000 kPa             | 41,000 kPa                          |
| (4,500 psi or 310 bar) | (6,000 psi or 410 bar)              |

A tank of argon gas at these pressures can provide approximately 35 minutes (Europe) or 70 minutes (USA) of freeze time. The control panel sounds an alarm when 2 minutes of gas time remain. At least one spare tank of argon gas should be readily available for each procedure.

When you receive the gas tanks, check that the labels on the tank identify the gas and the purity. If the labels are damaged or missing, call the gas supplier for a replacement.

Always move gas tanks carefully. Be sure that the outlet valve is fully closed and the valve protection cap is properly secured. Always move gas tanks with carts or hand trucks designed for that purpose. Never drop the tanks or allow them to strike another tank. After moving a tank to its point of use, secure it using stands, clamps, or other securing devices recommended by the gas supplier.

## Refrigerant Usage

- The console operates with a refrigerant tank that contains high-pressure argon gas. Store and handle the argon refrigerant tanks according to the hospital gas management policy.
- Do NOT use refrigerant tanks on anesthetic equipment.
- Do NOT open the valve on the argon tank until it is properly connected to the regulator and control panel.

## Shipping the Refrigerant Tank

- Follow standard guidelines for the safe handling and storage of high-pressure gas cylinders.

## Changing the Refrigerant Tank

The tank carrier holds one tank of refrigerant at a time. To change the tank, you will need:

- A full refrigerant tank (available from your hospital supplier)
- An open-ended wrench (supplied with the system)

**Caution:** Only use Argon gas tanks that meet the requirements of local regulations. In Canada, only use Argon gas tanks that meet CSA B51 requirements.

**Caution:** Follow standard guidelines for the safe handling and storage of high-pressure gas cylinders.

The blue section of the tank carrier can function as a dolly for loading and unloading the gas tank.

1. Close the tank by turning the knob on the top of the tank clockwise.
2. Press the Vent button on the control panel. The display screen will read "Close Gas Tank Valve and Press Vent Button."
3. Ensure the gas tank is closed and press the Vent button on the control panel a second time. The light above the button lights up to indicate that venting has been selected. The following message appears on the display screen: "System Venting"
4. A loud hissing noise will be heard as the gas is vented from the gas line between the tank and the regulator.
5. The auxiliary vent valve on the regulator assembly can also be used to vent any remaining gas in the gas line. Open the auxiliary vent valve by turning the handle counterclockwise.  
**Caution:** Always turn off the main gas cylinder valve prior to venting the inlet pressure line to ensure that excess pressure is safely vented.
6. Use the power switch to turn off the control panel and disconnect the power cable.
7. Remove the control panel from the console by grasping the handles and lifting straight up until the control panel clears the hooks on the tank carrier, then tilt the control panel toward the user. Disconnect the pressure sensor cable and gas hose from the back of the control panel before moving it to a table.  
**Warning:** ALWAYS remove the control panel from the tank carrier before changing the tank to prevent user harm. If not, the tank carrier may fall forward when removing the tank.
8. Disconnect the regulator from the tank.
  - a. Use the supplied wrench to loosen the nut.  
**Note:** The nut on the North American regulator turns clockwise to loosen.
  - b. Unscrew the nut until the regulator assembly can be separated from the tank.
  - c. Set the regulator assembly aside.
  - d. Fit a cap onto the tank valve.
9. Push down the levers on the front wheels of the console to lock the wheels.
10. Unlock the supporting legs of the tank carrier.
  - a. Pull on the knob at the back of the tank carrier and turn it 90 degrees to lock the knob open (Figure 23).

## Chapter 6

Figure 23.

- b. Pull the legs away from the tank carrier to extend them.
- c. Release and turn the knob so it snaps into the lower hole on the legs and locks the legs in the extended position (Figure 24).

Figure 24.

11. On the top of the tank carrier, a toggle clasp attaches to a post to lock the two pieces of the tank carrier together. Unlock the tank carrier by releasing the black toggle clasp then lifting the hook away from the post (Figure 25). Repeat on the other side.


Figure 25.

**Warning:** The tank carrier dolly is heavy. Be ready to take its weight.

**Warning:** Only use the tank carrier dolly for loading and unloading refrigerant gas tanks. Do not use the tank carrier to transport refrigerant gas tanks.


12. Separate the dolly from the tank carrier by lowering the dolly until all four wheels are on the ground.

## Maintenance

Figure 26.

13. Wheel the dolly away from the front section of the tank carrier.
14. Wheel the dolly to the location of the new refrigerant tank.
15. Stand the dolly up in a vertical position. (The back wheels will now be off the ground.)
16. Unbuckle the black straps holding the empty tank in the dolly.
17. Remove the empty tank from the dolly and set it aside.
18. Set the full tank into the dolly and buckle the straps.
19. Wheel the dolly back to the tank carrier and dock it with the front section. Line up the dolly so the black pegs on both sides of the dolly fit into the hooks on the front section (Figure 27). Lift the dolly into a vertical position.

1. Hooks
2. Peg

Figure 27.

20. Lock the tank carrier together by placing the hooks on the top of the dolly over the posts on the front section. Press the black toggle clasps down to lock the two sections together (Figure 28).

Figure 28.

21. To fold and lock the supporting legs of the dolly, pull and turn the knob on the back of the dolly to retract the legs back into the tank carrier. Release the knob so it snaps back into the upper hole and locks the legs into place.
22. Connect the regulator assembly to the refrigerant tank.
  - a. Remove the cap from the full argon tank.
  - b. Connect the regulator assembly using the supplied wrench.
  - c. Replace the control panel on the tank carrier.

## Cleaning the Cardioblate CryoFlex Surgical Ablation Console and Accessories

Ensure the console power is off and unplugged before continuing.

## Chapter 6

Use a mild solution of soap and water or other suitable nonabrasive hospital cleaning agent to wipe the exterior surfaces of the console and tank carrier.

Use care when cleaning keypad and display screen. If the control panel is deteriorated or worn, contact Medtronic Customer Service immediately.

**Caution:** Do not expose the console, accessories, or electrical connectors to harsh cleaning fluids or solvents.

## Internal Maintenance

The Cardioblade CryoFlex Surgical Ablation Console has no internal user serviceable parts or operator adjustments. Contact Medtronic Customer Service for service contract information.

## Troubleshooting

## Troubleshooting

7

**Warning:** If an unanticipated event occurs, the injection of refrigerant can be stopped at any time by pushing the Stop Freeze button on the control panel.

**Warning:** If you are unable to stop the injection, immediately turn the system off by using the power switch located on the side of the control panel or by disconnecting the power cord.

If the symptom cannot be resolved, contact Medtronic Technical Support at the phone number on the back cover.

| Symptom                                                                                                        | Possible Causes                                                                                                                                                                                         | Suggested Actions                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No power                                                                                                       | <ul style="list-style-type: none"> <li>Power cord is not plugged in correctly</li> <li>Cables and hoses are not securely attached</li> </ul>                                                            | <ul style="list-style-type: none"> <li>Check that the power cord is securely plugged into the side of the control panel.</li> <li>Check that the power cord is plugged into a hospital-grade outlet.</li> <li>Check that all cables and hoses are securely attached.</li> </ul>                                                                                                                                |
| No time and temperature readings on the Display Screen                                                         | <ul style="list-style-type: none"> <li>Gas tank not open</li> <li>Cables and hoses are not securely attached</li> <li>Probe or clamp is not securely attached</li> </ul>                                | <ul style="list-style-type: none"> <li>Check that the argon tank is open and properly connected to the regulator.</li> <li>Check that the gas hose and the pressure sensor cable are connected to the back of the control panel.</li> <li>Check that the probe/clamp is connected to the control panel. Tug on the cable to ensure that it is securely attached. Ensure that the pin is pulled out.</li> </ul> |
| Probe is difficult to connect                                                                                  | <ul style="list-style-type: none"> <li>System is pressurized</li> <li>Improper connection</li> <li>Attempting to connect the probe improperly</li> </ul>                                                | <ul style="list-style-type: none"> <li>Manually vent the system (see Chapter 5).</li> <li>Gently wiggle connector to connect.</li> <li>Do not press the top and bottom grips while connecting.</li> <li>Do not use excessive force to connect.</li> </ul>                                                                                                                                                      |
| Probe is difficult to disconnect                                                                               | <ul style="list-style-type: none"> <li>Attempting to disconnect the probe incorrectly</li> <li>System is pressurized</li> </ul>                                                                         | <ul style="list-style-type: none"> <li>See Chapter 5, Figure 21</li> <li>Manually vent the system (see Chapter 5).</li> </ul>                                                                                                                                                                                                                                                                                  |
| Probe is connected but Display Screen says "0 minutes of gas remain" or number of minutes remaining fluctuates | <ul style="list-style-type: none"> <li>Argon tank is not open</li> <li>Argon tank is empty</li> <li>Regulator/pressure sensor is damaged</li> </ul>                                                     | <ul style="list-style-type: none"> <li>Ensure tank is open.</li> <li>Replace argon tank if empty.</li> <li>Verify regulator and connections.</li> </ul>                                                                                                                                                                                                                                                        |
| Probe does not freeze                                                                                          | <ul style="list-style-type: none"> <li>Start Freeze button not pressed long enough</li> <li>Wrong kind of gas</li> <li>Gas is not flowing</li> <li>Cable and hose connections are not secure</li> </ul> | <ul style="list-style-type: none"> <li>Press the Start Freeze button for 2 seconds to start a freeze.</li> <li>Check the tank label to ensure that the argon gas being used is of the proper purity and pressure.</li> <li>Check that the gas hose and the pressure sensor cable are connected to the back of the control panel.</li> </ul>                                                                    |
| Probe is connected but is unable to inject, or the system does not recognize the probe                         | <ul style="list-style-type: none"> <li>Improper probe connection; electrical connection is loose/intermittent</li> <li>Pressure sensor cable is loose</li> </ul>                                        | <ul style="list-style-type: none"> <li>Ensure black pin is pulled out.</li> <li>Gently wiggle the connector.</li> <li>If it still will not connect, remove probe and check connector site for damage.</li> <li>Try connecting a second probe. If a second probe will not connect, contact Medtronic Technical Support.</li> <li>Check pressure sensor cable connection.</li> </ul>                             |
| Cannot achieve therapeutic temperatures and handle is cold                                                     | <ul style="list-style-type: none"> <li>Argon tank is closed</li> <li>Injection line is blocked or broken</li> <li>Too much blood/fluid in the operating site</li> </ul>                                 | <ul style="list-style-type: none"> <li>Ensure tank is open.</li> <li>If frost is forming on handle, connect a second probe.</li> <li>Remove excess fluid.</li> </ul>                                                                                                                                                                                                                                           |
| Cannot achieve therapeutic temperatures and device temperature rises                                           | <ul style="list-style-type: none"> <li>Using wrong gas (often helium)</li> </ul>                                                                                                                        | <ul style="list-style-type: none"> <li>Turn Console off.</li> <li>Change tank to Argon.</li> <li>Verify Argon specs via certificate on tank.</li> </ul>                                                                                                                                                                                                                                                        |
| System beeps                                                                                                   | <ul style="list-style-type: none"> <li>2 minutes of gas or less remain in the tank</li> </ul>                                                                                                           | <ul style="list-style-type: none"> <li>Replace Argon tank.</li> </ul>                                                                                                                                                                                                                                                                                                                                          |

Chapter 8

## Technical Specifications

8

### Surgical Cryoablation System Specifications

|                                         |                                                                                                                                       |
|-----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Storage and Shipping Temperature:       | -40°C to 66°C (-40°F to 150°F)                                                                                                        |
| Operating Temperature:                  | 10°C to 44°C (50°F to 110°F) at altitudes ranging from sea level to 3000 meters (10,000 feet)                                         |
| Storage and Operating Humidity:         | 10% to 95%, noncondensing                                                                                                             |
| Electrical:                             | 120 VAC / 230 VAC, 50/60Hz, Max Current: 2A                                                                                           |
| Fuses:                                  | 2A, 250VAC, Type T, Breaking capacity: 1500A at 250V~                                                                                 |
| Type of Protection from Electric Shock: | Class 1                                                                                                                               |
| Mode of Operation:                      | Continuous                                                                                                                            |
| Power Disconnect:                       | Power cord                                                                                                                            |
| Target Freeze Time:                     | 5 seconds to 5 minutes, selectable in 5-second increments                                                                             |
| Target Temperature:                     | -30°C to -95°C, selectable in 5°C increments. Defaults to continuous gas flow (maximum cold mode, down to -150°C), indicated by "ice" |
| Probe Tip Temperature:                  | controlled to ±10°C of Target Temperature                                                                                             |
| Displayed Temperature:                  | ±3°C                                                                                                                                  |

### How Supplied

The Cardioblate CryoFlex Surgical Ablation Console is supplied in a customized crate.

### Electromagnetic Emissions and Immunity Declarations

| IEC 60601-1-2 Table 1, Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                        |            |                                                                                                                                                                                                                                                                                                                              |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The CryoFlex system is intended for use in the electromagnetic environment specified below. The customer or the user of the CryoFlex system should assure that it is used in such an environment. |            |                                                                                                                                                                                                                                                                                                                              |
| Emissions test                                                                                                                                                                                    | Compliance | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                       |
| RF Emissions<br>CISPR 11                                                                                                                                                                          | Group 1    | The CryoFlex system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                                                                                                                     |
| RF Emissions<br>CISPR 11                                                                                                                                                                          | Class A    | The CryoFlex System is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded.                      |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                               | Class A    | Warning: This equipment/system is intended for use by healthcare professionals only. This equipment <b>may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the CryoFlex System or shielding the location.</b> |
| Voltage fluctuations/ flicker emissions<br>IEC 61000-3-3                                                                                                                                          | Complies   |                                                                                                                                                                                                                                                                                                                              |

| IEC 60601-1-2 Table 2, Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                         |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The CryoFlex system is intended for use in the electromagnetic environment specified below. The customer or the user of the CryoFlex system should assure that it is used in such an environment. |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                  |
| Immunity test                                                                                                                                                                                     | IEC 60601 test level                                                                                                                                                                                                                                              | Compliance level                                                                                                                                                                                                                                                  | Electromagnetic environment – guidance                                                                                                                                                                                                                                                           |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                                                                                                    | ±6 kV contact<br>±8 kV air                                                                                                                                                                                                                                        | ±6 kV contact<br>±8 kV air                                                                                                                                                                                                                                        | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                                    |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                                                                                                                  | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                                                                                                                      | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                                                                                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                              |
| Surge<br>IEC 61000-4-5                                                                                                                                                                            | ±1 kV line(s) to line(s)<br>±2 kV line(s) to earth                                                                                                                                                                                                                | ±1 kV differential mode<br>±2 kV common mode                                                                                                                                                                                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                              |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11                                                                                            | <5% U <sub>n</sub> (>95% dip in U <sub>n</sub> ) for 0.5 cycle<br>40% U <sub>n</sub> (60% dip in U <sub>n</sub> ) for 5 cycles<br>70% U <sub>n</sub> (30% dip in U <sub>n</sub> ) for 25 cycles<br><5% U <sub>n</sub> (>95% dip in U <sub>n</sub> ) for 5 seconds | <5% U <sub>n</sub> (>95% dip in U <sub>n</sub> ) for 0.5 cycle<br>40% U <sub>n</sub> (60% dip in U <sub>n</sub> ) for 5 cycles<br>70% U <sub>n</sub> (30% dip in U <sub>n</sub> ) for 25 cycles<br><5% U <sub>n</sub> (>95% dip in U <sub>n</sub> ) for 5 seconds | Mains power quality should be that of a typical commercial or hospital environment.<br>If the user of the CryoFlex system requires continued operation during power mains interruptions, it is recommended that the CryoFlex system be powered from an uninterruptible power supply or a battery |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                                                                                                                        | 3 A/m                                                                                                                                                                                                                                                             | 3 A/m                                                                                                                                                                                                                                                             | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                        |
| Note: U <sub>n</sub> is the a.c. mains voltage prior to application of the test level.                                                                                                            |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                  |

## Technical Specifications

| IEC 60601-1-2 Table 4, Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                         |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The CryoFlex system is intended for use in the electromagnetic environment specified below. The customer or the user of the CryoFlex system should assure that it is used in such an environment. |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Immunity test                                                                                                                                                                                     | IEC 60601 test level        | Compliance level | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                                     | 3 Vrms<br>150 kHz to 80 MHz | 3 Vrms           | Portable and mobile RF communications equipment should be used no closer to any part of the CryoFlex system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><b>Recommended separation distance</b><br>$d = 1.2\sqrt{P}$<br>$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz<br>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. <sup>1, 2</sup><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br> |
| Radiated RF<br>IEC 61000-4-3                                                                                                                                                                      | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

| IEC 60601-1-2 Table 5, Recommended separation distances between portable and mobile RF communications equipment and the CryoFlex system                                                                                                                                                                                                                                                                                                                 |                                                                |                                        |                                         |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------|-----------------------------------------|
| The CryoFlex system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CryoFlex system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CryoFlex system as recommended below, according to the maximum output power of the communications equipment. |                                                                |                                        |                                         |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                                          | Separation distance according to frequency of transmitter<br>m |                                        |                                         |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 150 kHz to 80 MHz<br>$d = 1.2\sqrt{P}$                         | 80 MHz to 800 MHz<br>$d = 1.2\sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.3\sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.12                                                           | 0.12                                   | 0.33                                    |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 0.38                                                           | 0.38                                   | 0.75                                    |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 1.2                                                            | 1.2                                    | 2.3                                     |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 3.8                                                            | 3.8                                    | 7.3                                     |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 12                                                             | 12                                     | 23                                      |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Cautions:**

- Do not use non-Medtronic components with Medtronic in-line-powered external devices. The use of non-Medtronic components may result in damage to Medtronic components, increased emissions, or decreased electromagnetic immunity of the Medtronic devices or systems.
- Do not use Medtronic in-line-powered external devices adjacent to, or stacked with, other electronic devices. Using Medtronic devices in these configurations may result in decreased electromagnetic immunity of the Medtronic devices or systems.

- ↑ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CryoFlex system is used exceeds the applicable RF compliance level above, the CryoFlex system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CryoFlex system.
- ↓ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Chapter 9

## Warranties

9

### Limited Warranty for the Cardioblade® CryoFlex® Surgical Ablation Console

- A. This LIMITED WARRANTY provides the following assurance to the purchaser of the Cardioblade® CryoFlex® Surgical Ablation Console, hereafter referred to as the "Equipment," that should the Equipment fail to function within normal tolerances due to a defect in materials or workmanship within a period of one year, commencing with the delivery of the Equipment to the purchaser, Medtronic, Inc. will at its option: (a) repair or replace any defective part or parts of the Equipment; (b) issue to the purchaser a credit equal to the original Equipment purchase price (but not to exceed the value of the replacement Equipment), against the purchase of replacement Equipment, or (c) provide functionally comparable replacement Equipment at no charge.
- B. To qualify for this repair, replacement or credit, the following conditions must be met:
- (1) The Equipment must be returned to Medtronic within sixty (60) days after discovery of the defect (Medtronic may, at its option, repair the Equipment on site).
  - (2) The Equipment must not have been repaired or altered by someone other than Medtronic in any way that, in the judgment of Medtronic, affects its stability and reliability.
  - (3) The Equipment must not have been subjected to misuse, abuse or accident.
- C. This LIMITED WARRANTY is limited to its express term. In particular, Medtronic is not responsible for any incidental or consequential damages based on any use, defect or failure of the Equipment, whether the claim is based on warranty, contract, tort or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

### Disclaimer of Warranty for the Accessories

Although the Cardioblade® CryoFlex® Surgical Ablation Regulators, pressure hoses, pressure sensor cables, and cord sets, hereafter referred to as "Product," have been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of the DISCLAIMER OF WARRANTY. Medtronic, Inc. therefore disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of the DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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