INSTRUCTION MANUAL

Digital Behind-the-Ear Hearing Instrument ENYA with Accessories

“APPROVE”

GN ReSound A/S, Denmark

Teamlead Partnerships & Quality

(position)

Carin Engelbrecht Petersen

(name)

Laura E. Petersen

(signature)

“ “ 22 10 2018

«day» month (numerals)

GN ReSound A/S

Lautrupbjerg 7

2750 Ballerup - Denmark


The Danish Chamber of Commerce do hereby confirm
that the company is a member of our organization and
known to us as worthy of confidence.

Pia Lind

23 OCT 2018

2018

Danish Chamber of Commerce
Secretary: Pia Lind

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Land:		Denmark Danmark	
This public document Dette offentlige dokument			
2. has been signed by er underskrevet af		Pia Lind	
3. acting in the capacity of i egenskab af		Secretary Sekretær	
4. bears the seal/stamp of er forsynet med segl/stempel af		Danish Chamber of Commerce Dansk Erhverv	
Certified Attesteret			
5. at i	Copenhagen København	6. the den	23 Oct 2018 23 okt 2018
7. by af	Ministry of Foreign Affairs of Denmark Udenrigsministeriet		
8. No nr.	D9786ACA		
9. Seal/stamp: Segl/stempel:		10. Signature: Underskrift:	 Pia Ejdal Frandsen

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, scan the QR code or visit the following website:

<https://e-register.um.dk>

Name of Medical Device: Digital behind-the-ear hearing instrument ENYA with accessories

Description of Medical Device:

Digital Behind-the-Ear Hearing Instrument ENYA, Models:

- High-power digital behind-the-ear hearing instrument ENYA:

I. HI,EY488-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging

20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

II. HI,EY388-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"

15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

III. HI,EY288-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. **Audio shoe**

11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

- Medium-power digital behind-the-ear hearing instrument ENYA:

I. HI, EY477-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left

5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

II. HI,EY377-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

III. HI,EY277-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments

12. Electroacoustic component comprised of:

- Microphone;
- Amplifier;
- Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

- Low-power digital behind-the-ear hearing instrument ENYA:

I. HLEY467-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)
9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)

- Mid-size powerful ear mold, no more than 20 units (if necessary)
- Large powerful ear mold, no more than 20 units (if necessary)
- Small open ear mold, no more than 20 units (if necessary)
- Mid-size open ear mold, no more than 20 units (if necessary)
- Large open ear mold, no more than 20 units (if necessary)
- 11. Instruction for the behind-the-ear hearing instruments
- 12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet “Autophone”
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

II. HI,EY367-DW comprised of:

1. Hook for adolescents (if necessary)
2. Hook for adults (if necessary)
3. Hook for children (if necessary)
4. Hook for children with a filter (if necessary)
5. Hook for adolescents with a filter (if necessary)
6. Hook for adults with a filter (if necessary)
7. Metal hook for children (if necessary)
8. Metal hook for adults (if necessary)

9. Thin tubes (various sizes) for connection with ear molds with a cleaning line, no more than 50 units (if necessary)
10. Ear molds (various sizes) for thin tubes, no more than 20 units:
 - Tulip ear mold, no more than 20 units (if necessary)
 - Small powerful ear mold, no more than 20 units (if necessary)
 - Mid-size powerful ear mold, no more than 20 units (if necessary)
 - Large powerful ear mold, no more than 20 units (if necessary)
 - Small open ear mold, no more than 20 units (if necessary)
 - Mid-size open ear mold, no more than 20 units (if necessary)
 - Large open ear mold, no more than 20 units (if necessary)
11. Instruction for the behind-the-ear hearing instruments
12. Electroacoustic component comprised of:
 - Microphone;
 - Amplifier;
 - Receiver.

Accessories:

1. Programming cable CS 44, right
2. Programming cable CS 44, left
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units)
4. Programming adapter cable 500 mm, left
5. Programming adapter cable 500 mm, right
6. Programming adapter cable 2500 mm, left
7. Programming adapter cable 2500 mm, right
8. Integrated FM shoe NA
9. Integrated FM shoe NB
10. Audio shoe
11. Meter
12. Remover of the battery compartment
13. Hook filters, 4 units
14. Magnet "Autophone"
15. Remover of the housing
16. Case for thin hooks, empty
17. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging
18. Multi microphone MULTI MIC, RS (EU) c/w packaging
19. Micro microphone MICRO MIS, RS (EU) c/w packaging
20. Thin tube adapter
21. Red marker, no more than 5 units
22. Blue marker, no more than 5 units
23. Volume control plug
24. Hearing instrument packaging:
 - Carton;
 - Plastic case;
 - Pouch;
 - Napkin;
 - Brush.

1. Device class

Class IIa medical devices (see the list above and the attached approved product list) according to R&TTE Directive 1999/5/EC (radio and telecommunication terminal equipment) and Medical Devices Directive 93/42/EEC (amendment 2007/47/EC), including the corresponding accessories.

2. Compliance assessment

Range of design, production and final inspection certification of the given class IIa devices and corresponding accessories. Quality standard compliance certificate ISO 9001:2008, certificate No. 099446 QM08, and ISO 13485:2012 + AC:2012, certificate No. 099446 MP2012, issued by DQS Medizinprodukte GmbH, as well as CE certificate No. 099446 MR2 certifying compliance with requirements of Annex II, excluding Article 4 of Council Directive 93/42/EEC concerning medical devices, issued by DQS Medizinprodukte GmbH.

3. Production

Primary manufacturer

GN ReSound A/S, Lautrupbjerg 7, DK-2750 Ballerup, Denmark

Addresses of the places of medical device production

1. GN ReSound A/S, Lautrupbjerg 7, 2750 Ballerup, Denmark;
2. GN ReSound China Ltd., 15, Chuang Xin Road, CTP Area, Xiamen 361006, China;
3. GN ReSound (Malaysia) Sdn Bhd, №.5, Jalan i-Park 1/1, Kawasan Perindustrian i-Park, 81000 Bandar Indahpura, Kulai, Johor, Malaysia.




Digital behind-the-ear hearing instrument ENYA with accessories manufactured by GN ReSound A/S, Denmark) must comply with the internal technical regulation on design and production of digital behind-the-ear hearing instruments ENYA with accessories (confidential commercial information due to the description of technological know-hows) of GN ReSound A/S, Denmark.

GN ReSound A/S, Denmark, is entitled to design and develop these medical devices.

The quality control department performs final inspection in compliance with the requirements approved by GN ReSound A/S.

4. Operating Principle

The operation of digital hearing instrument is based on multiple gain of sound by means of digital and analog components of the product. Digital signal processor enables additional signal processing, including digital small-signal gain, background noise correction and automatic adjustment to the sound environmental conditions.

1. Program selector button
2. On/Off with the battery compartment
3. Output sound hole
4. Front microphone
5. Back microphone
6. Volume control
7. Direct audio input
8. Device labeling (left – blue, right - red)
9. Battery door lock
10. Model
11. Manufacturer
12. Serial No.

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5. Destination of medical device

For electroacoustic correction of mild to severe (deafness) hearing loss.

Indications for Use:

- Mild, moderate and severe hearing loss
- Monaural or binaural hearing loss

Contraindications:

- Any conditions that require further diagnosis: sudden hypoacusis, suspected retrocochlear pathology, etc.
- Surgically correctable forms of conductive hypoacusis

Potential Side Effects:

Allergic reactions while using the medical device are unlikely. Please consult your hearing care professional or doctor if you have itching, redness, soreness, inflammation or burning sensation in the ears or near them.

Method of Application:

The medical device is worn on the upper part of the ear; the ear mold is inserted into the external auditory canal.

Application Conditions:

The medical device is intended only for individual use and can be used with constant wear.

Precautionary Information:

1. Do not leave hearing instruments in the sun, near the open fire, or in a hot parked car.

2. Do not wear hearing instruments while taking a shower, swimming in the pool, in the sauna, in heavy rain or snowfall.
3. If your hearing instrument becomes wet, remove the battery and place it with an open battery compartment in a sealed container with an absorbent. Your hearing care professional will help you in choosing the care products for your hearing instrument.
4. Remove hearing instruments for a period of use (application) of cosmetic products, such as perfume, aftershave, hair spray, sunblock, etc.
5. Low-power digital coded transmission is used during a wireless signal transfer that can affect the performance of nearby electronic equipment. To avoid this, do not use the instruments near devices that can be influenced by this type of transmission.
6. If an electromagnetic interference with other electronic devices occurs during the wireless communication, increase the distance between them.
7. Use only original ReSound accessories, for example, thin tubes, receivers and ear molds. Never attempt to modify the shape of the hearing instrument and the ear mold. This will lead to the withdrawal of the warranty and can cause damage to the device.
8. Connect ReSound devices only with the accessories that are supplied by ReSound and specially designed for this purpose.

Contact with human body

Components	Materials coming into contact with human body
Hearing instrument housing	ABS Terluran GP-22 (housing bottom, housing top, battery door, program selector button, volume control) Stainless steel AISI 316L C – 0.08 Cr – 23.00-26.00 Ni – 10.00-14.00 Mo – 2.00-3.00 Mn – 2.00 P – 0.05 Si – 1.00 S – 0.03 (connector [pin]) Silver ink: TPC 218/79-AB-NT Hardener HN Teca (labeling element – model, manufacturer, serial No.)
Hook: - for adolescents - for children - for adults	ABS (MABS) Terluxe 2802 TR
Hook for children with a filter Hook for adolescents with a filter Hook for adults with a filter	ABS (MABS) Terluxe 2802 TR (hook) Hostaform POM C52021 (filter)
Hook filters	Hostaform POM C52021
Metal hook for children Metal hook for adults	Stainless steel AISI 316L C – 0.08 Cr – 23.00-26.00 Ni – 10.00-14.00 Mo – 2.00-3.00 Mn – 2.00


	P – 0.05 Si – 1.00 S – 0.03
Integrated FM shoe: - NA - NB	Polyamide RILSAN BESNO
Audio shoe	ABS (MABS) Terluxe 2802 TR
Thin hook adapter	ABS (MABS) Terluxe 2802 TR
Meter	Thermoplastic polyurethane Estane 58245 TPU
Color mark, red	ABS Terluran GP-22 (blue)
Color mark, blue	ABS Terluran GP-22 (red)
Thin hooks	Thermoplastic elastomer PEBA 5533 SA 01 MED
Ear mold: - tulip	Silicone ELASTOSIL R401/40
Powerful ear mold: small, mid-sized, large	Silicone ELASTOSIL R401/40
Open ear mold: small, mid-sized, large	Silicone ELASTOSIL LR 3003/20 TR A/B
Volume control plug	ABS Terluran GP-22

Minimum Scope of Delivery:

- hearing instrument ENYA
- instruction manual
- case

6. Technical specification

High-power programmable digital behind-the-ear hearing instrument: HI,EY488-DW, HI,EY388-DW, HI,EY288-DW

Hearing correction level: hypoacusis levels III-IV

Model Overview	HI,EY488-DW	HI,EY388-DW	HI,EY288-DW
Model description			
Battery size		13	
Functional technologies			
Fully configurable programs	4	4	3

Synchronized Program Selector Button	•		
Synchronized Volume Control	•		
SmartStart	•	•	•
PhoneNow	•		
Comfort Phone	•	•	•
2.4 GHz Binaural Connection	•		
ReSound Unite TV Streamer	•	•	
ReSound Unite TV Streamer 2	•	•	
ReSound Unite Remote Control	•	•	•
ReSound Unite Remote Control 2	•	•	•
ReSound Unite Phone Clip+	•	•	
ReSound Unite Midget Microphone	•	•	
Multi Microphone MULTI MIC	•	•	
Micro Microphone MICRO MIC	•	•	
ReSound Control Application (requires Phone Clip+)	•	•	
Acoustical functions			
WARP Signal Processing Channels	10	8	6
SoftSwitching	•	•	
Adaptive Directionality	•	•	•
Fixed Directionality	•	•	•
NoiseTracker II	•	•	•
Expansion	•	•	•
WindGuard	•	•	
DFS Ultra II	•	•	•
Auto DFS	•	•	•
Tinnitus Generator	•	•	•
Setting functions			
Aventa 3.9 or newer	•	•	•
Adjustment channels	No more than 10	No more than 8	No more than 6
Onboard Analyzer II	•	•	•
Safe Setup	•	•	•
In-situ Audiometry	•	•	•
Wireless Setup with Airlink2	•	•	•

* Programmer (programming interface: AIRLINK2) – registration certificate No. RZN 2015/2727 dated June 04, 2015.

Aventa 3.9 or newer (Aventa 3.x, where x is the version No.)

SOFTWARE SAFETY CLASS – class B

The hearing instrument firmware can be accessed only at the manufacturing site. Firmware version No. 3.7.10.

Smart Start – this function delays activation of the hearing instrument by 5 or 10 seconds after closing the battery compartment.

PhoneNow – automatic switch of the hearing instrument to the telephone program when a telephone receiver with PhoneNow magnet is raised to the ear.

Comfort Phone – used in case of binaural hearing loss correction. When a telephone receiver with PhoneNow magnet is raised to one ear and the hearing instrument is automatically switched to the telephone program, the gain of the second instrument goes down by 6 dB to avoid excessive environmental noise.

Softswitching – automatic switch from omnidirectional mode to directional mode and vice versa depending on the presence/absence of speech and environmental noise.

Noise Tracker II – multilevel noise reduction system based on noise recognition and subtraction from the general signal.

Expansion – low-level noise (such as noises made by refrigerators, computer coolers etc.) reduction technology.

Windguard – wind noise reduction technology.

Sound Shaper – high-to-medium frequency sound compression technology. Used if the patient does not react to high frequencies.

DFS Ultra II with Music Mode – feedback (hissing noise) reduction technology with a special mode for listening to music characterized by the absence of reduction of high tones of musical compositions.

Auto DFS – this function activates if feedback calibration proves infeasible. It will automatically determine critical frequencies and gain capable of producing feedback while wearing the hearing instrument.

Tinnitus Generator – tinnitus (ringing in the ears) therapy technology. The user receives different sound signals, e.g. relaxing sounds of nature or constant/modulated white noise, distracting from the basic noise.


Onboard Analyzer II – option of latching of statistical data on the conditions of wearing of hearing instruments by the patient: average wearing time per day, volume control and program selection activity, period of return to different acoustic situations etc.

In-situ Audiometry – possibility of audiometry by means of a hearing instrument. Sound tones used to register response are produced by a hearing instrument. Used if no audiometer is available.

Wireless Setup with Airlink2 – possibility of configuring hearing instruments without using additional shoes and wires by means of an Airlink 2 programmer connected to the specialist's computer.

* All types of directionality work on the principle of mixed directionality: low frequencies are processed in the omnidirectional mode, medium and high frequencies – in the directional mode. The frequencies processed in the omnidirectional mode are determined by the configuration program for each patient individually.

Adaptive directionality – automatic change of the pattern of directionality from behind depending on intensity and location of the noise source.

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









1. Adaptive directionality











2. Fixed directionality


Low-power programmable digital behind-the-ear hearing instrument: HI,EY467-DW, HI,EY367-DW

Hearing correction level: hypoacusis levels II-III

Model Overview		HI,EY467-DW	HI,EY367-DW
Model description			
Battery size	13		
Available colors	5		
Functional technologies			
Fully configurable programs	4	4	
Synchronized Program Button	•		
Synchronized Volume Control	•		
SmartStart	•	•	
PhoneNow	•	•	
Comfort Phone	•		
2.4 GHz Binaural Connection	•		
ReSound Unite TV Streamer	•	•	
ReSound Unite TV Streamer 2	•	•	
ReSound Unite Remote Control	•	•	
ReSound Unite Remote Control 2	•	•	
ReSound Unite Phone Clip	•	•	

ReSound Unite Phone Clip+	•	•
ReSound Unite Midget Microphone	•	•
Multi Microphone MULTI MIC	•	•
Micro Microphone MICRO MIC	•	•
ReSound Control Application (requires Phone Clip+)	•	•
Audiological functions		
WARP Signal Processing Channels	10	8
SoftSwitching	•	•
Adaptive Directionality	•	•
Fixed Directionality	•	•
NoiseTracker II	•	•
Expansion	•	•
WindGuard	•	•
DFS Ultra II	•	•
Auto DFS	•	•
Tinnitus Generator	•	•
Setting functions		
Aventa 3.9 or newer	•	•
Adjustment channels	No more than 10	No more than 8
Onboard Analyzer II	•	•
Safe Setup	•	•
In-situ Audiometry	•	•
Wireless Setup with Airlink2	•	•

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Maximum OSPL (90 dB SPL at the input)	Max. 1,600 Hz/HFA	132 125	122 117	dB SPL
Total harmonic distortion	500 Hz 800 Hz 1,600 Hz	2.2 2.3 0.7	1.8 1.5 0.25	%
Telecoil sensitivity (1 mA/m at the input) HFA-SPLIV @ 31.6 mA/m (ANSI) Full sensitivity of the telecoil @ 1mA/m	Max. HFA 1,600 Hz/HFA	94 82	98 77	dB SPL
Noise ratio at the input		24	22	dB SPL
Noise ratio of 1/3 octave at the input (without noise reduction)		11		dB SPL
Frequency range (DIN 45605/ANSI)		100-7,150	100-7,110	Hz
Current Drain		1.1/1.2	1.1/1.2	mA

Parameter	
Compression coefficient	1:0 3:0
AGC attack time (ms)	10
AGC recovery time (ms)	50
Electrical input sensitivity to 70 dB at 1,600 Hz (mV)	2.6

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


Main Frequency Characteristics


Frequency characteristics at different positions of tone controls


Variation of OSPL90 with the change of supply voltage


Variation of full-on acoustic gain with the change of supply voltage


Variation of harmonic factor with the change of supply voltage


Changing OSPL90 with the change of internal resistance of the power supply

Variation of full-on acoustic gain with the change of internal resistance of the power supply





Variation of harmonic factor with the change of internal resistance of the power supply




Pattern of directionality


1. Adaptive directionality



6300 Hz



7100Hz






8000Hz




2. Fixed directionality

Charts of the change of frequency characteristics of hearing instruments, when external microphones are used.







External microphones are digital devices designed to receive an audiosignal, subtract a useful signal from the surrounding noise (changing the SNR) (SNR — Signal to Noise Ratio), transform it into a digital form and transmit it to a user's hearing aid by a radio channel. External microphones can be used in two modes: in the first case a user will hear the sound coming only from the external microphone, in the second case – both from the external microphone and the hearing aids. Adjustment of the modes can be performed by a hearing care professional in the software Aventa 3. It is important to note that an external microphone can be paired with a hearing aid without the software, although in this case a user will hear only a signal coming from the external microphone.





The signal coming from the external microphone is given at the input of a digital path, therefore frequency characteristics of hearing aids do not change, when external microphones are used.

Directional characteristics (the pattern of directionality) of external microphones.

1. Micro microphone works only in a directional mode.





Patterns of directionality:






2. Multi microphone has two modes: directional, when being worn on the body of a person, and multidirectional, when placed on a flat surface horizontally.

Patterns of directionality:

IP58 Under the influence of dust and moisture, the apparatus does not change its characteristics. Water can penetrate into the battery compartment, but after replacing the battery with a new one, the characteristics of the apparatus become similar to the original characteristics. Since zinc-air batteries used in hearing aids require air for operation, the design of the battery compartment is not sealed and therefore not fully protected from water penetration.

Technical specifications of a magnet for automatic switching of the hearing instrument into phone conversation mode.

Magnet "Autophone" power – 1,370-1,436 mT

1. Programming cable, right 2. Programming cable, left Intended use: used to connect the hearing instrument to a computer. Weight: 24.70 g Length: 2,500.00 mm	
3. Programming adapter (2 units) with a fixture (red – 2 units, blue – 2 units) Intended use: used to connect the hearing instrument to a programming cable. Weight: 1.00 g Dimensions: 15.0 x 11.0 x 7 mm	

4. Programming adapter cable 500 mm, left

Intended use: used to connect the hearing instrument to a computer.

Weight: 7.5 g

Length: 500.0 mm

5. Programming adapter cable 500 mm, right

Intended use: used to connect the hearing instrument to a computer.

Weight: 7.5 g

Length: 500.0 mm

6. Programming adapter cable 2,500 mm, left

Intended use: used to connect the hearing instrument to a computer.

Weight: 13.5 g

Length: 2,500.0 mm


7. Programming adapter cable 2500 mm, right

Intended use: used to connect the hearing instrument to a computer.




Weight: 13.5 g

Length: 2,500.0 mm

8. Integrated FM shoe NA 9. Integrated FM shoe NB Intended use: an adapter for connecting to FM systems. Weight: 1.10 g Dimensions: 19.0 x 12.0 x 7.0 mm	
10. Audio shoe Intended use: an adapter for connecting to audio systems. Weight: 0.8 g Dimensions: 19.0 x 16.0 x 7.0 mm	
11. Meter Intended use: a rounded measuring ruler intended for determining the required tube size for connecting to ear molds (an accessory required for open prostheses). Weight: 1.7 g Dimensions: 60.0 x 40.0 x 1.0 mm	
12. Remover of the battery compartment Intended use: a tool for opening the battery door. Weight: 1.90 g Dimensions: 5.6 x 0.9 x 0.9 mm	
13. Hook for children 14. Hook for adolescents 15. Hook for adults Intended use: a part of a hearing instrument conducting sound from the instrument to the ear mold and the ear. Weight: 0.20 g, 0.22 g and 0.24 g Dimensions: 17.0 x 16.0 x 5.0 mm, 19.0 x 16.0 x 5.0 mm, 21.0 x 15.0 x 5.0 mm	

16. Hook filters, 4 units Intended use: used to prevent ingress of moisture and dirt into the hearing instrument; an acoustic filter forming correct frequency response. Weight: 1.10 g Dimensions: 5.6 x 0.9 x 0.9 mm	
17. Hook for children with a filter 18. Hook for adolescents with a filter 19. Hook for adults with a filter Intended use: a part of a hearing instrument conducting sound from the instrument to the ear mold and the ear. Weight: 0.21 g, 0.23 g and 0.25 g Dimensions: 17.0 x 16.0 x 5.0 mm, 19.0 x 16.0 x 5.0 mm, 21.0 x 15.0 x 5.0 mm	
20. Metal hook for children 21. Metal hook for adults Intended use: a set of acoustic hooks of various sizes without a filter for hearing instruments. Weight: 0.9 g, 1.2 g Dimensions: 17.0 x 15.0 x 5.0 mm, 20.0 x 15.0 x 5.0 mm	
22. Magnet "Autophone" Intended use: used for automatic switching of the hearing instrument into phone conversation mode. Weight: 2.50 g Dimensions: 12.0 x 5.0	
23. Remover of the housing Intended use: used to remove metal pins holding housing elements. Weight: 1.90 g Dimensions: 5.6 x 0.9 x 0.9 mm	

24. Case for thin hooks Intended use: a case for storing thin hooks. Weight: 495.0 g Dimensions: 360.0 x 225.0 x 45.0 mm	
25. Thin hooks (various sizes) for connection with ear molds with a cleaning line, no more than 50 units Configuration: thin hooks with sports locks: 0 (right and left), 1 (right and left), 2 (right and left), 3 (right and left). Hook cleaning line Intended use: thin hooks for open prostheses. Conduct sound from a hearing instrument to the ear. Weight: 0 – 0.11 g, 1 – 0.13 g, 2 – 0.15 g, 3 – 0.17 g Dimensions: 0 – 59.0 mm, 1 – 62.0 mm, 2 – 68.0 mm, 3 – 74.0 mm Sports lock length: 40.0 mm Hook cleaning line: Weight: 0.035 g Dimensions: 100.0 x 0.6 mm	 <i>Спортивная застежка Sports lock</i> 
26. Ear molds (various sizes) for thin tubes, no more than 20 units: - Tulip ear mold, no more than 20 units Intended use: an ear mold for connecting to a hearing instrument and use by hearing-impaired people. Weight: 0.15 g Dimensions: 10.0 x 12.0 x 4.0 mm	
- Small powerful ear mold, no more than 20 units - Mid-size powerful ear mold, no more than 20 units - Large powerful ear mold, no more than 20 units Intended use: an ear mold for connecting to a hearing instrument and use by hearing-impaired people. Weight: 0.07 g, 0.09 g and 0.1 g Dimensions: 7.0 x 7.0 x 7.0 mm, 10.0 x 9.0 x 9.0 mm and 12.0 x 10.0 x 12.0 mm	

- Small open ear mold, no more than 20 units
- Mid-size open ear mold, no more than 20 units
- Large open ear mold, no more than 20 units

Intended use: an ear mold for connecting to a hearing instrument and use by hearing-impaired people.

Weight: 0.01 g, 0.01 g and 0.01 g
 Dimensions: 5.0 x 5.0 x 2.0 mm,
 6.0 x 8.0 x 2.0 mm and 8.0 x 10.0 x 2.0 mm

27. Remote control device RC-2, PLUS, UNITE, RS (EU) c/w packaging

Composition: USB cable, AC/DC adapter, carrying case, case belt.

Intended use: a remote control for a hearing instrument with accessories operating in a 2.4 GHz wireless network.

Adapter:
 Input voltage: 100-240 V, current: 350 mA
 Output voltage: ± 5 V, current: 1,000 mA
 Discharge time: up to 80 hours (depending on the intensity of use)
 Capacity: 240 mA/h

Commercial set package weight: 210.0 g
 Remote control weight: 32.0 g
 Commercial set package dimensions:
 160.0 x 120.0 x 45.0 mm
 Remote control dimensions: 95 x 41 x 12 mm

28. Multi microphone MULTI MIC, RS (EU) c/w packaging


Composition: USB cable, AC/DC adapter, carrying case, case belt.



Intended use: an external wireless microphone for hearing instruments with a wireless connectivity range up to 25 m.

Adapter:
 Input voltage: 100-240 V, current: 200 mA
 Output voltage: ± 5 V, current: 1,000 mA
 Discharge time: 6 hours in the conversation mode,
 80 hours in the sleeping mode
 Capacity: 240 mA/h

Commercial set package weight: 195.0 g
 Microphone weight: 21.8 g

Commercial set package dimensions: 160 x 120 x 45 mm Microphone dimensions: 60.0 x 28.0 x 21.0 mm 29. Micro microphone MICRO MIC, RS (EU) c/w packaging Composition: USB cable, AC/DC adapter, interface cable (3.5 mm jack), carrying case, case belt. Intended use: an external wireless microphone for hearing instruments with a wireless connectivity range up to 25 m. Adapter: Input voltage: 100-240 V, current: 200 mA Output voltage: ± 5 V, current: 1,000 mA Discharge time: 6 hours in the conversation mode, 80 hours in the sleeping mode Capacity: 240 mA/h Commercial set package weight: 195.0 g Microphone weight: 17.4 g Commercial set package dimensions: 173 x 110 x 58 mm Microphone dimensions: 53.0 x 25.0 x 18.0 mm	
30. Thin hook adapter Intended use: used to connect a thin hook. Weight: 0.20 g Dimensions: 9.0 x 7.0 x 8.0 mm	
31. Red marker, no more than 5 units 32. Blue marker, no more than 5 units Intended use: red and blue plastic labeling plugs for hearing instruments (left – blue, right – red). Weight: 0.09 g Dimensions: 32.0 x 5.0 x 2.0 mm	
33. Volume control plug Intended use: closing the gap left after removing the volume control from the housing. Weight: 0.05 g Dimensions: 8.0 x 2.0 x 2.0 mm	

34. Hearing instrument packaging: - Carton - Plastic case - Pouch - Napkin - Brush Intended use: storage of a hearing instrument. Weight: 85.0 g Dimensions: 160.0 x 120.0 x 45.0 48.0 g 90.0 x 43.0 x 30.0 0.50 g 90.0 x 70.0 x 15.0 0.13 g 145.0 x 145.0 x 1.0 0.8 g 62.0 x 9.0 x 4.0	    
35. Instruction for the behind-the-ear hearing instruments Intended use: contain description of the hearing instrument, recommendations for use and a warranty card.	

Note:

Dimensions are within $\pm 2\%$; maximum weight.

See the following links:

<https://www.resound.com/en-us/help/accessories#anchor3>

containing the instructions for use of:

- Multi microphone MULTI MIC, RS (EU)
- Micro microphone MICRO MIS, RS (EU)

RESISTANCE TO WEAR



Battery door – at least 5200 cycles of operation.

- Operating switches and switches with handles (program selector button and volume control) – at least 5200 cycles of operation.
- Ports – 300 cycles of operation.

Operation of the hearing instrument**On/Off function**

1. When the battery door is closed, the hearing instrument is turned on, and the default program is activated.
2. To turn on the hearing instrument, open the battery door.

Always keep the device in the off position when you are not using it.

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
“Smart Start” function

The hearing instrument can be turned on once you have placed it on your ear. If you prefer to turn on the device in advance, just before putting on, then you can use the special function “Smart Start”. This function delays activation of the device after closing the battery compartment.

When “Smart Start” turn-on delay is on, you will hear beep tones at 1 second intervals. Delay activation and setting is performed by an expert.

Battery replacement

1. Open the battery door.
2. Remove the discharged battery and replace it with a new one. The battery is inserted into the device with the “+” sign up according to the designation on the door.
3. Carefully close the battery door.

Useful tips

1. Always use zinc-air batteries with the remaining shelf life of at least 1 year.
2. When not using the device, keep the battery door open. This will prevent corrosion and prolong battery life.
3. Before inserting the battery remove the sticker from it and expose it to the air for about 2 minutes.
4. If your hearing instrument often loses contact with ReSound Unite wireless accessories, contact your specialist to rule out low resistance of batteries.

Low battery indicator

Your hearing care professional can activate the low battery indicator. The device will reduce gain and emit a warning signal when the battery charge is running low. It will be repeated every 15 minutes until the device automatically turns off. The term of device deactivation after actuation of the indicator may vary depending on the type of battery. We recommend to always keep spare batteries.

Low battery indicator in the devices paired with ReSound Unite accessories

Active use of ReSound Unite accessories (remote control device, TV streamer, Phone Clip Bluetooth hands-free kit, midjet microphone) or FM receiver results in a quicker battery discharge of devices. This means that the battery life directly depends on the amount of wireless accessories connected to the hearing instrument. When the battery power drops to the limits when the device will no longer maintain a connection with accessories you will hear two descending signals. After that your remote control device will operate normally, but you will not be able to use a hands-free kit, midjet microphone and TV streamer.

After some time the battery power will decrease even more. You will hear two descending signals again and you will cease control over device settings using the remote control device. Hearing instruments will continue to function normally. Once the device batteries are replaced by new ones function of the accessories will be restored.

How to put on/off your hearing instrument

How to put on a hearing instrument with custom ear molds

1. Hold the ear mold between your thumb and index finger directing the output sound hole in your ear canal.
2. Insert the ear mold into your ear canal by means of careful rotational movements.
3. Turn the ear mold first up, then down and gently press it down. Opening and closing your mouth will facilitate insertion.
4. Make sure the hearing instrument is seated behind the ear.

The device will be easier to put on if you retract the tip of your ear with one hand while inserting the device with the other.

Putting on and removing of the device becomes easier with time. The more correct the position of the ear mold in the ear, the more comfortable the use of the device.

If the ear mold causes discomfort, burning sensation, bursting pain etc., contact your hearing care professional.

Attention! Never attempt to modify the ear mold, receiver or the housing of a hearing instrument.

How to put off a hearing instrument with custom ear molds

1. Hold the hook tube and pull it outward.
2. Consult your hearing care professional if you have difficulties removing the hearing instrument.

How to put on a hearing instrument with a thin tube

Hold the curve of the thin tube with fingers and gently insert the ear mold in the ear canal.


It is important to correctly position the ear mold in the ear.

If the ear mold is positioned correctly, the tube wire fits the ear flatly, without bulging.

Attention! Never attempt to modify the thin tube or the ear mold.

How to put off a hearing instrument with a thin tube

1. Hold the curve of the thin tube and pull it outward.
2. Consult your hearing care professional if you have difficulties removing the hearing instrument.

Sports lock


The sports lock is designed to fix the device on the ear during active movements of the user. The sports lock should be mounted by a specialist.

Volume control

The volume control is used to increase or decrease volume.

1. To increase volume, press the volume control up.
2. To decrease volume, press the volume control down.


You will hear a beep signal with any volume change. You will hear a long signal when the volume limit is reached (maximum or minimum).

If you have two hearing instruments with an active synchronization function, volume adjustment of the first device will be automatically duplicated in the second device.

Your hearing care professional may turn off the volume control if you do not need it or even remove it from the device housing and replace it with a special plug.

Program selector button

Depending on your experience with hearing instruments, individual hearing needs, and your way of life, your hearing care professional may install additional hearing programs. If additional programs have been added to your device, read the following instructions.


Program	Application
1	
2	
3	
4	

1. You can switch between programs by pushing the program selector button.
 2. You will hear a beep signal when switching from one program to another: the number of beeps indicates which program you have selected.
 3. When the hearing instrument is turned off and then back on, the hearing instrument always automatically returns to the first standard program.
- If you have two hearing instruments with an active synchronization function, program switch of the first device will be automatically duplicated in the second device.

Remote control device (purchased separately).

Remote control device (hereinafter RCD) will allow you to adjust settings of your hearing instrument easily and conveniently. Please read the following description and instructions for use of ReSound Unite Remote Control Device 2.

Description:

1. On/Off
2. Screen
3. Activation of wireless accessories
4. Return to the default settings
5. Mute of the hearing instrument microphones
6. Increase volume
7. Decrease volume
8. Select the left device
9. Select the right device
10. Program switch

11. Connection button

RCD On/Off

To turn the RCD on, press down the On/Off switch (1). You will see white indicator “I” and the RCD screen will light up.

To turn the RCD off, press up the On/Off switch.

Program switch

Button “P” (10) is used to switch programs that your hearing care professional has programmed into your hearing instrument. The hearing care professional can also assign a unique name to each program which will be displayed in the RCD screen. Press button “P” once to switch to the next program. Upon reaching the last program the next time you press the button the device will return to program 1. Thus, program switching is cyclic (e.g., 1, 2, 3, 1).

Each program switch is accompanied by an audible signal and the program number is displayed in the RCD screen. If you wish, your hearing care professional may disable acoustic indication.

For optimal hearing, it is recommended to set the same program on both devices while wearing.

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
Volume adjustment

Adjusting volume of two devices at the same time


The buttons labeled “+” (6) and “-” (7) are used to adjust volume. You can adjust volume of two devices at the same time and separately. By default, both hearing instruments (binaural hearing aids) are “active”; it is indicated by two arrows (< and >), located at the bottom of the RCD screen. Both devices can be simultaneously made louder by pressing “+” and quieter by pressing “-”. If this feature is not disabled by your hearing care professional, every change of volume is accompanied by an acoustic signal of hearing instruments and a change of the volume bar on the RCD screen. The volume level programmed by the hearing care professional by default is displayed as a green horizontal line on the volume bar.

Adjusting volume of the left hearing instrument only

To adjust volume in the left hearing instrument only, press button “<” (8). This activates the left side of the screen; the right side of the screen becomes inactive. An arrow is displayed at the bottom of the RCD screen on the left; the arrow on the right side of the screen will disappear and the volume control bar will become dim.

To increase or decrease volume of the left hearing instrument, use buttons “+” and button “-”.

Adjusting volume of the right hearing instrument only

To adjust volume in the right hearing instrument only, press button “>” (8). This activates the right side of the screen; the left side of the screen becomes inactive. An arrow is displayed at the bottom of the RCD screen on the right; the arrow on the left side of the screen will disappear and the volume control bar will become dim.

To increase or decrease volume of the right hearing instrument, use buttons “+” and button “-”.

To return to volume adjustment of two devices at the same time, press button “<” or “>” (of the inactive side). Both arrows will light up on the left and on the right again.

Important! If the binaural synchronization function is available in your hearing instruments and active, individual volume adjustment of one of the devices will automatically adjust volume of the another one.

Mute/Activation of hearing instrument microphones

You can mute or activate microphones of hearing instruments by pressing the mute button (5) on the RCD.

With binaural hearing aids, it is possible to mute/activate microphones only on both devices at the same time.

When microphones are muted, the microphone icon  changes to , and volume bars become dim.

To activate hearing instrument microphones again, press the mute button once more.

Activation of wireless accessories

If you use ReSound Unite wireless accessories, such as TV streamer and midget microphone, you can activate them with RCD by pressing button "Activation of wireless accessories" (3).

Return to the default settings

Default settings are the parameters of the hearing instrument set up by your hearing care professional and activated when you turn on the device.

To return to the default settings, click the appropriate button on the RCD (4). This will:

- activate Program 1 in both hearing instruments (binaural hearing aids);
- set programmed default volume;
- set the volume adjustment mode on both devices at the same time.

Button "Return to the default settings" will also stop transmission of the signal from the ReSound Unite wireless accessory.

Airplane mode

All wireless functions must be disabled when boarding the plane, because it is not allowed to use devices performing radio frequency transmission during the flight or some parts of the flight or when entering a zone where radio frequency transmission is prohibited.

Wireless connection in ReSound ENYA hearing instruments can be turned off by opening and closing the battery door of the device three times within 10 seconds (opening and closing, opening and closing, opening and closing).

To restore the wireless connection, open and close the battery door of the device once. Wireless communication will be restored in 10 seconds.

Important! If you open and close the battery door again during the 10 seconds when wireless connection is recovering, the airplane mode will be resumed. For this reason, do not open the battery door again until after 15 seconds to restore the wireless function.

Telephone use

During a phone call hold the phone to your ear as usual if your device has open ear mold or tulip ear mold. Try to find the best position when raising the phone to your ear if you have a powerful or custom ear mold. You might need to practice. The following advice may be helpful for using the phone.

1. Hold the telephone as you would normally.
2. Hold the telephone towards the top of the ear (closer to the microphones).
3. If whistling occurs, wait a few seconds before the hearing instrument adapts.
4. Any whistling may also be decreased by moving the telephone to the side or holding it slightly away from the ear.
5. Depending on your individual needs, your hearing care professional may install a program specifically for telephone use.


Listening to the radio/TV

While listening to the radio or TV, start by listening to the news programs, because anchors usually speak clearly, unlike speakers in other programs.

If it is difficult to understand what is being said on the radio or television, ask your hearing care professional about the available accessories that may help resolve this problem.

Mobile phone use

Your hearing instrument is designed in accordance with strict international standards of electromagnetic compatibility. However, not all mobile phones are compatible with hearing instruments. Various interferences may occur caused by the phone itself or the signal operator. If you have any difficulties while talking on the phone due to weak signal or insufficient speech intelligibility, ask your hearing care professional about the available accessories that may help resolve this problem.

PhoneNow

PhoneNow allows your hearing instrument to switch automatically to the phone conversation program when a telephone receiver is raised to the ear. When the telephone receiver is removed from the ear, the hearing instrument automatically returns to the previous program.

Placement of PhoneNow magnets

To properly attach the magnet, follow these instructions:

1. Clean the telephone receiver thoroughly (the part of the phone from where the sound comes). Use only a recommended cleaning agent.
2. Hold the telephone vertically, in a position similar to when making a telephone call.
3. Place the magnet just below the telephone receiver (just below the sound output). If necessary, move the magnet to another position to improve ease of use and comfort while speaking. You may also attach an additional magnet.

PhoneNow use

The telephone can be used in a normal manner. A short melody will indicate that the PhoneNow feature has automatically switched the hearing instrument to your phone program. You may need to move the telephone receiver slightly to find the best position for reliable PhoneNow activation and good signal quality.

If you have two hearing instruments with an active synchronization function, volume of the second device (not used with the phone) will be automatically reduced.

PhoneNow precautions

1. Keep magnets out of reach of pets, children and mentally challenged persons. If a magnet is swallowed, call a doctor immediately!
2. The magnet enforces the magnetic field of a telephone receiver. For this reason we recommend keeping the magnet at least 30 cm away from magnetically sensitive electronic devices, e.g. pacemakers, credit cards etc.
3. High distortion during a call may mean that the magnet is not in the optimal position relative to the telephone receiver. To improve quality of the connection, move the magnet.
4. Use only the magnets supplied by ReSound.

Telecoil (optional)

If your device features a telecoil function, then your hearing care professional can configure one additional program – Telephone (induction) coil. Telecoil catches the magnetic signal supplied by the phone and converts it to sound. This option enables better perception of the sound from the handset.

Induction loop

Many public spaces, such as theatres, churches, schools etc. are equipped with induction loop systems. The telecoil function enables reception of an excellent clean sound and impressive speech intelligibility.

If you do not hear the sound indoors using this program, the induction loop is not installed or is disabled. If the sound is not sufficiently clean and clear, go around the room to find the optimal position for perception.

FM receiver and Direct Audio Input (DAI) (optional)

FM system can significantly improve speech intelligibility in situations when the person you are talking to is at a large distance, or you are in a noisy or reverberating room.

Direct audio input (DAI) enables direct connection of your hearing instrument to such audio items as television, radio, audio player etc.; it helps to obtain excellent sound quality. The audio shoe is used to connect to an audio system's cable or a wireless FM system. Once you put the audio boot on, the hearing instrument will switch to the corresponding program automatically.

Connection/Disconnection of an audio shoe or an FM receiver

Connection of an audio shoe and an FM receiver

1. Align the audio shoe or the FM receiver with the groove located on the internal surface of the device's housing between the battery compartment and the serial number.
2. Having connected the audio shoe, insert it deeper towards the battery compartment.
3. Carefully snap the audio shoe.

Disconnection of an audio shoe and an FM receiver


4. Press and hold the button on the front surface of the audio shoe. For FM receiver: move the snap down with a nail.
5. Carefully remove the shoe from the hearing instrument.

DAI

FM

Important aspects of using FM systems

1. Do not use two transmitters in one FM channel
2. Do not use water or fluids for cleaning the FM receiver
3. Do not use the FM transmitter in places where the use of electronic devices is prohibited, e.g. in airplanes.
4. Keep in mind that the FM signal can be intercepted and heard by other receivers.
5. Before using the system in another country, contact an expert to make sure your radio channel is permitted in that country.
6. Only authorized service centers may repair your FM receiver and transmitter.

Battery door lock

Battery doors of your hearing instruments may be locked to prevent opening by a child. The lock system may vary depending on the model of your device.

Standard lock system

Your hearing care professional may use special labeling for your hearing instruments (to distinguish between the left-hand side device and the right-hand side device. This labeling may also be used to lock the battery compartment. After the battery door is locked, the hearing instrument may be turned off only after it is unlocked.

To lock the battery door:

1. Close the battery door.
2. Use a special lock tool (blue) to move the color marker (left – blue, right - red) to the right.

To unlock the battery door:

3. Slide the color marker to the left.


Hearing instrument care

Your hearing instrument is coated with a layer of protective hydrophobic material. Nevertheless, please follow the following recommendations to prolong the service life of your hearing instruments:

1. Keep your hearing instrument clean and dry. Wipe housing with a soft cloth or polish cloth after each use to remove moisture and dirt. Do not use water or solvents as these can damage the hearing instrument.
2. Never immerse hearing instruments in water or other liquids.
3. Avoid rough handling of hearing instruments or dropping them on hard surfaces or floor.
4. Do not leave hearing instruments exposed to or near the sources of heat or sunlight, such as in a hot parked car, because excessive heat may cause damage or deformation of the housing.
5. Do not wear hearing instruments while taking a shower, swimming in the pool, in heavy rain or in the sauna.
6. If your instrument gets wet, or if it has been exposed to high humidity, it should be left to dry out overnight. To do this, remove the battery and leave the battery door open. It is highly recommended to put the hearing instrument in a sealed container with an absorbent material for a night. Do not use the instrument until it is completely dry. Contact your hearing care professional regarding the selection of the absorbent material. You may also be offered special care sets for your hearing instrument.
7. Remove hearing instruments for a period of use (application) of cosmetic products, such as perfume, aftershave, hair spray, sunblock, etc., because they may penetrate the hearing instrument and cause damage.

Daily maintenance


It is important to keep your hearing instrument clean and dry. Every day, after use, the device should be cleaned with a soft dry cloth. Remove ear wax and other contaminants using a special brush. In order to avoid damage due to humidity, put your hearing instrument into a box with absorbent material. You may also use a special device care set.

Cleaning ear molds


1. To clean, remove the hook with the ear mold.
2. Clean the ear mold with warm soapy water, and rinse with tap water.
3. A wire loop can be used to remove the residual dirt from the tube.
4. After washing and blowing down (with or without a special bulb syringe), thoroughly dry out the ear mold.

With time, ear mold becomes hard and loses shape. Contact your hearing care professional regarding ear mold replacement.

Cleaning thin tubes and ear molds

1. To clean, remove the thin tube with the ear mold.
2. Clean the thin tube and the ear mold with warm soapy water, and rinse with tap water.
3. Use a special black wire to remove the clean the hook.
4. After washing the hook dry it thoroughly.

With time, hook becomes hard (after three months of use on the average). This negatively affects acoustic quality. Contact your hearing care professional regarding hook replacement.

Receiver care

Your hearing care professional will select a receiver with optimal power depending on the hearing loss level. The receiver is connected to the device with a wire, which transmits sound into the ear. It is important to ensure fitting of the receiver and the ear mold and their correct position within the ear. If the receiver causes irritation, redness or bursting pain, consult with your hearing care professional. Never attempt to modify shape of your ear mold or receiver. Regularly clean your receiver and ear mold. Please read the following recommendations.

Cleaning receivers and ear molds

Regularly clean the receiver tube and the ear mold. Use a soft dry cloth to wipe the surface of the receiver tube and the ear mold. Do not use water. Clean a UP (ultra powerful) receiver in the same manner. Please follow the wax filter replacement scheme (see below).

Cleaning custom ear molds

1. Remove the ear mold from the receiver tube.
2. Clean the ear mold with warm soapy water, and rinse with tap water.
3. After washing the ear mold dry it thoroughly.

Wax filter replacement

To replace the wax filter, follow these steps:


1. Remove the used wax filter. To do that, put the end of a rod designed to remove the filter in the used filter and carefully remove it.
2. To install a new filter, put the other end of the rod with a new filter in the sound output and press slightly. Remove the rod and ensure that the new filter remains in the receiver.

How to replace an ear mold

To avoid damage, we recommend that ear molds be replaced by a hearing care professional.

Standard ReSound ear mold


1. Put an ear mold on the thin tube (see the pic.).
2. Make sure that the ear mold is fit tightly.

ReSound tulip ear mold

To replace ear molds of this type, please keep in mind the following peculiarities. A tulip ear mold has two “petals”. It is important that the larger petal is above (over) the smaller one! To ensure that, follow these instructions:

1. Put an ear mold on the thin tube.
2. Turn the larger ear mold petal inside out and then expand it (see the pic.).

Use of ReSound hearing instruments with the mobile application

Purpose of the mobile application

Application GN ReSound for smartphones is designed for use with ReSound wireless hearing instruments. Mobile application GN ReSound sends and receives signals from ReSound wireless hearing instruments using mobile devices, for which this application has been developed.

How to use the mobile application

1. Mobile application update notifications should always be turned on. It is recommended to install all updates for proper operation of the application.
2. The application should be used only with GN ReSound devices, for which it is designed. GN ReSound is not responsible for any damages if the mobile application has been used with other devices.

Troubleshooting

Malfunction	Reason of malfunction	Remedy
No sound	Not turned on	Turn on
	Dead battery	Replace the battery
	Battery door not closed	Tightly close the battery door
	Clogged ear mold	Clean the ear mold
	Clogged hook	Clean the hook
	Standby	Press the top button of the volume control
Not loud enough	Incorrect ear mold placement	Insert the ear mold correctly
	Clogged ear mold	Clean the ear mold
	Clogged hook	Clean the hook
	Change in hearing sensitivity	Consult with your hearing care professional
	Excessive ear wax	Consult with your physician
	Minimum setting of the volume control	Increase volume or consult with your hearing care professional
	Minimum volume	Press the volume control up
Feedback (whistling)	Incorrect ear mold placement	Insert the ear mold correctly
	Incorrect ear mold placement on a hook	Place the ear mold correctly
	The hook lost elasticity and became hard	Replace the hook
	Hearing instrument hook damaged	Replace the hearing instrument hook
	Excessive ear wax	Consult with your physician
	Calibration of the hearing instrument required	Consult with your hearing care professional
	Settings of the hearing instrument not optimized	Consult with your hearing care professional
Distorted sound	Low battery	Replace the battery
	Improper ear mold	Consult with your hearing care professional
	Hearing instrument damaged	Consult with your hearing care professional
	Settings of the hearing instrument not optimized	Consult with your hearing care professional
Wireless connection	Possible cause – airplane	Open and close the battery door.


does not function	mode	Wireless connection shall restore in 10 seconds. If the issue persists, contact your hearing care professional
-------------------	------	--

Instructions for hearing care professionals

1. Additional hearing program installation instructions

Additional hearing programs can only be installed by a specialist when programming the hearing aid via Aventa 3 software. Only one program is available on the hearing aid by default. 1-3 additional programs may be installed depending on the aid's technological level.

To add a program, click the additional program bookmark located to the left of the active program bookmark.

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
Beeps and Volume Control
 Manage Accessories
 Select Sound
 Start
 Pre-Fit
 Fit
 Summary
 Right Ear
 Connect
 Left Ear
 Display Options
 Gain View
 Output View
 Gain Curves
 Target Curves
 Target Gain
 P1
 All-Around
 None
 Streamers
 TV / Mini Mic / Multi-
 Manage Programs
 Recalculate
 Autorelate
 Copy
 Paste
 Remote Control Program Name
 All
 Calibrate DFS
 Save

A list of available programs will be displayed upon clicking; choose the program you need.

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
Telecoil phone + Mic

Traffic

Programs are added in a sequential order: 1, 2, 3...


To remove additional programs, enter the "Manage programs" menu located below.

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In the pop-up window with the list of active programs choose the program you wish to remove and click the red cross next to it.

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Cancel

Apply

2. Instructions for switching off (at user's request), removing from the housing and replacing the volume control with a special plug


Switching the volume control off (at user's request)

A specialist may switch the volume control off via Aventa software at user's request. In that case pushing volume control buttons will not result in any changes, i.e. the volume will remain the same.

To switch the volume control off go to the "Beeps and Volume Control" bookmark located in the left menu of the "Fit" screen.

In the pop-up window choose the "Manual Volume Control" bookmark.

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
10 s

Copy

The volume control settings will appear below.

Put a marker (dot) next to the "Volume control" bookmark to switch the volume control off.

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The volume control is now switched off.

3. Removing the volume control from the housing

The volume may not only be switched off, but also completely removed from the housing if necessary. A special plug is put in the housing instead of the removed volume control. The volume control plug is supplied with the hearing aid and matches the housing in color and material.

Caution! All activities associated with removing the volume control must be conducted by a specialist.

To remove the volume control, use the Remover of the battery compartment.


There is a retaining Pin near the volume control on the hearing aid's housing.

Use Remover of the battery compartment to remove the Pin and, after that, the volume control.

Put a plug in instead of the removed volume control.

7. Precautions

1. Do not leave hearing instruments in the sun, near the open fire, or in a hot parked car.
2. Do not wear hearing instruments while taking a shower, in the pool, in the sauna, in heavy rain or snowfall.
3. If your hearing instrument becomes wet, remove the battery and place it with an open battery compartment in a sealed container with an absorbent. Your hearing care professional will help you in choosing the care products for your hearing instrument.
4. Remove hearing instruments for a period of use (application) of cosmetic products, such as perfume, aftershave, hair spray, sunblock, etc.
5. Low-power digital coded transmission is used during a wireless signal transfer that can affect the performance of nearby electronic equipment. To avoid this, do not use the instruments near devices that can be influenced by this type of transmission.
6. If an electromagnetic interference with other electronic devices occurs during the wireless communication, increase the distance between them.
7. Use only original ReSound accessories, for example, thin tubes, receivers and ear molds. Never attempt to modify the shape of the hearing instrument and an ear mold. This will lead to the withdrawal of the warranty and can cause damage to the device.
8. Connect ReSound devices only with the accessories that are supplied by ReSound and specially designed for this purpose.

8. Safety requirements

1. Consult with your hearing care professional if you find a foreign body in the ear canal, experience irritation or discomfort, or if you have excessive ear wax deposition associated with the use of hearing instruments.
2. Different types of radiation, such as NMR and MRI, CT scanners and devices for radiation therapy may damage hearing instruments. It is recommended not to wear hearing instruments during these diagnostic procedures. Other types of radiation, such as produced by alarm systems, fire sensors, radio equipment, mobile phones, etc. are not powerful enough to damage the hearing instrument, but may slightly affect sound quality.
3. Do not wear hearing instruments in mines, oil fields or other explosive areas, except for the cases when such areas are certified for the use of hearing instruments.
4. Do not allow others to use your hearing instrument. This can result in damage to the hearing instrument or hearing loss in individuals whom this device is not indicated.
5. For safety reasons, the use of the device by children or persons with mental disorders should be done under constant supervision.
6. Hearing instruments contain small parts that can be swallowed by a child. Please be careful, do not leave your child with a hearing instrument without supervision.
7. Hearing instruments should only be used as prescribed by your hearing care professional, otherwise it may result in hearing loss.
8. Always turn off wireless connection when boarding a plane and entering the places where radio communication use is prohibited.
9. Do not use a broken hearing instrument.

Safe use of batteries

Batteries of any size contain dangerous substances and must be disposed properly. The following information is for the safety of you and the environment.

1. Do not attempt to charge non-rechargeable batteries (zinc-air). They may explode or leak.
2. Do not attempt to dispose of the batteries by incineration. Please dispose of them according to the rules established in the Russian Federation, or return back to your hearing instrument center.
3. Do not put batteries in your mouth. If you accidentally swallow the battery, call your physician immediately because it may be harmful to your health.
4. Keep batteries out of reach of children, pets, and persons with mental disorders.
5. Remove the battery from the device if you do not use it for a long time.

9. Requirements to environmental protection when using the medical device

This medical device does not have a negative impact on people and the environment during use, transportation and storage.

Class A – epidemiologically safe waste similar to municipal solid waste in composition.

Danger class	Morphological composition features
(epidemiologically safe waste similar to MSW in composition)	Waste not contacting body fluids of patients and contagious patients. Stationery, packaging, furniture, equipment characterized by lost marketability, sweepings etc. Food scraps of central nutrition units and any divisions of an organization performing medical and/or pharmaceutical activities, apart from the departments treating contagious patients, including patients suffering from tuberculosis

10. Sterility

GN ReSound A/S hearing instruments are supplied non-sterile and should not be sterilized.

11. Transportation and storage

Transportation is carried out by all types of covered transport at the temperature from -20 to +60 °C and the relative humidity not more than 90% without condensation.

Store at the temperature from -20 to +60 °C and the relative humidity not more than 70% and the atmospheric pressure from 200 to 1,500 hPa.

12. Use

The device shall be used at the temperature from -10 to +40 °C and the relative humidity not more than 70% and the atmospheric pressure from 200 to 1,500 hPa.

13. Packaging

The device and accessories are packed in individual packages made of transparent polyethylene and placed in a shipping carton with cushioning layers of cardboard and polyethylene with manufacturer's labeling.

14. Labeling

The labeling of the medical device includes the following information:

- trademark of the manufacturer;
- addresses of the manufacturer and the authorized representative;
- name of the product;
- serial number.

Symbols used for labeling

	Item (serial) number
	Manufacturer, address of the manufacturing site of the medical device
	CE certified
	Keep dry
	See the instruction manual
	Type B applied part Degree of protection from electric shock (leak current)
	Class II equipment Protection from electric shock
	Temperature range
IP 58	Degree of protection from dust and water
	Special disposal Separate this object from regular waste and dispose of it in the safest possible manner, e.g., take it to a special disposal center to avoid damage to the environment
	This symbol indicates that the user should take note of the respective warnings listed in the user guide.

15. Service life

Average service life of the device is 5 years.

16. Disposal and destruction

Medical devices must be disposed of by local public law organizations.

17. Warranty

The warranty is valid for one year from the date of purchase. The warranty covers material defects and manufacturing flaws. The warranty does not cover accessories.

18. Applicable technical standards and guidelines

Documents/standards	Name
MDD 93/42/EEC	Council Directive 93/42/EEC of 14 June 1993 concerning medical devices
ANSI 63.19-2011	American National Standard Methods of Measurement of Compatibility between Wireless Communications Devices
IEC 60118-0:83+A1:94	Hearing aids – Part 0: Measurement of electroacoustic characteristics. Amendment 1 1994-01
IEC 60118-7:2005	Hearing aids – Part 7: Measurement of the performance characteristics of hearing aids for quality inspection for delivery purposes
IEC 60118-13:2004	Hearing aids – Part 13: Electromagnetic compatibility
ISO 10993-5:2009	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity
ISO 10993-10:2010	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization
ISO 14971:2007	Medical devices – Risk management
EN ISO 13485:2003+AC:2012	Medical devices – Quality management systems – Requirements for regulatory purposes

19. Maintenance and repair

Maintenance and repair is carried out by GN Hearing RUS, LLC (Limited Liability Company), registered at: Nizhniy Kislovskiy Pereulok, 7 building 1, premises 1, Moscow 125009, RF, 8 495 697 30 10\ 8 495 697 66 00, dganelin@gnresound.ru.

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