Blood glucose test-strips Medicare -Lis**User's manual****(Concerning the handling of the medical device on the territory of the Russian Federation)****Purpose**

Blood glucose test strips Medicare-Lis are intended for monitoring of capillary blood glucose level at home by people with diabetes or with a suspicion of diabetes and by health care professionals in medical facilities. They are used only with Blood glucose meter Medicare-Lis model BG-101. They are not intended for blood glucose testing of neonates.

Principal of measurements

Glucose contained in the blood sample when entering into a reaction with a special reagent in the test strip, produces a small electric current. The strength of current produced depends on the amount of glucose in the blood. The blood glucose meter Medicare-Lis measures the strength of the current and displays the results as a blood glucose level.

Characteristics

The blood sample volume is 0,3 μ L.

The time of measurement is 6 seconds.

The measurement range is 1,1 - 33,3 mmol/L (20-600 mg/dL) to a precision of 1 mg/dL.

Capillary blood from a finger or alternative testing site (palm, forearm) is used for the measurements.

Measurements accuracy: Blood glucose test-strips Medicare-Lis provide accurate measurement results if all the conditions described in this manual are observed.

Performance characteristics

Accuracy of measurements

> 95% in A zone of EGA vs YSI 2300 glucose analyzer:

Within ± 15 mg/dL (0.83 mmol/L) at glucose concentration < 75 mg/dL (or 4.2 mmol/L) and within 20% glucose concentration ≥ 75 mg/dL (or 4.2 mmol/L).

Precision of measurements

CVs(%) of intermediate precision and repeatability were less than 5%.

Standard deviation of measurement results ≤ 5 mg/dl at glucose concentration < 75 mg/dl

Coefficient of variation of measurement results ≤ 5 % at glucose concentration ≥ 75 mg/dl

Chemical composition of test strips

Each test strip contains:

1. glucose hydrogenase (FAD-dependent) -0,4%
2. mediators – 3,8%
3. buffer – 94,4%
4. non-reactive ingredients – 1,4%

Warning

Please, make sure that your blood glucose test strips Medicare-Lis are suitable for usage. Prior to the first use, make sure that the package is not damaged and closed.

Do not use test-strips if the blister or vial has a puncture or a tear in it.

Do not bend, cut and do not attempt in any way distort the test strip.

Do not cut the foil and do not use sharp instruments to take the test strip out of the blister.

If the drying agent has turned pink in a blister or in a vial, do not use the test strip.

Do not use expired test strips.

Avoid getting dirt, food or liquid on the test strip. Do not handle test strips with wet hands. All parts of the test strip can be touched but only with dry and clean hands. Remove the test strip from the package immediately before use.

CERTIFIED TRUE COPY

Victor Koh Koon Hor

Do not use for blood glucose testing of neonates.

Do not reuse. Test strips are intended for single use.

Store out of children's reach. A test strip may enter the respiratory system and cause choking. Do not swallow test strips. The test strip and its package contain substances that may be harmful if swallowed. If it happened, immediately see a doctor.

Use only capillary blood for blood glucose measurement, the use of other substances or plasma will lead to incorrect results.

Store the test-strips in accordance with the indicated temperature-humidity conditions, the violation of the indicated conditions may lead to obtaining of incorrect measurement results.

Store test strips only in their original package till the expiration date indicated on the package. To avoid damage or contamination, do not put test strips in any other container.

Hematocrit level below 30% or above 55% can lead to incorrect measurement results. Please, consult your doctor if you do not know your hematocrit level.

Body dehydration can lead to incorrect results. If you are suffering from severe dehydration, contact your doctor immediately.

Testing at altitudes up to 3,402 m (11161 feet) does not influence the result.

Do not change prescribed medication based on the test results without the consult of treating physician.

Additional information for healthcare professionals

1. Cholesterol content up to 500 mg/dL or triglycerides up to 2000 mg/dL do not affect test results significantly. However, glucose level at rates close to the indicated ones should be interpreted carefully.

2. Incorrect result may occur in hypotensive patients or patients in a state of shock. Inaccurate results are also possible in hyperosmolar state, with or without ketosis.

3. Interactions with other medical products: Reducing substances entering blood naturally (uric acid, bilirubin) or through therapeutic treatments (ascorbic acid, acetaminophen) don't affect the test results significantly. However, high concentration of these substances can affect the results. The substances and their acceptable concentrations are presented below:

Substance	The indicator below which the substance content in the blood is acceptable, and it does not affect the test result
Acetaminophen	$\leq 8,0$ mg/dL (0,53 mmol/L)
Ascorbic Acid	$\leq 5,0$ mg/dL (0,28 mmol/L)
Aspirin	≤ 60 mg/dL (3,33mmol/L)
Bilirubin	≤ 90 mg/dL (1,54 mmol/L)
Cholesterol	≤ 500 mg/dL (12,9mmol/L)
Creatinine	≤ 5 mg/dL (0,44 mmol/L)
Dopamine	$\leq 2,0$ mg/dL (0,11 mmol/L)
Galactose (Milk sugar)	≤ 900 mg/dL (50 mmol/L)
Gentisic Acid	$\leq 5,0$ mg/dL(0,32 mmol/L)
Hydroxyurea	$\leq 3,0$ mg/dL (0,39 mmol/L)
L-dopa	≤ 10 mg/dL (0,51 mmol/L)
Maltose (Malt sugar)	≤ 900 mg/dL (26,3 mmol/L)
Methyldopa	$\leq 3,0$ mg/dL (0,13 mmol/L)
Tolbutamide	≤ 400 mg/dL (14,8 mmol/L)
Triglycerides	≤ 2000 mg/dL (22,8 mmol/L)
Uric Acid	$\leq 8,0$ mg/dL (0,48 mmol/L)

Blood glucose measurement procedure

See your User's manual for the blood glucose meter Medicare-Lis model BG-101 for more detailed information concerning the test performance.

Test results

The result is shown in milligrams of glucose per deciliter of blood (mg/dL). The meter can display the range of 20-600 mg/dL.

Glucose level below 50 or above 250 mg/dL may indicate a serious disease. If your result is below 50 or above 250 mg/dL, please, contact a doctor immediately.

Conversion factor $1 \text{ mg/dL} \times 0,0555 \Rightarrow \text{mmol/L}$.

Incorrect results

If you have received the results which are not consistent with the way you feel:

1. Make sure that the blood sample is correctly applied to the test strip channel;
 2. Check the shelf life of the test strips;
 3. Check the performance of the meter and the test strips using control solution Medicare-Lis;
- Consult a doctor if the result is still high or low.

Control solution check

Use the control solution anytime you want to check performance of the meter and the test strip or practice the measurement performance procedure. Use only glucose control solution Medicare-Lis. The results of control measurements should be consistent with the results indicated on the test strip package.

Important: the control measurement result may vary with different control solution packages. Always refer to the information on the test strip package.

Possible forms of presentation

Blood glucose test strips Medicare-Lis are presented in following designs:

- in vials: 5, 10, 25, 50, 100, 200 pcs.
- in blisters: 5, 10, 15, 20, 50, 100, 200 pcs.

Shelf life: 24 months from the manufacture date.

IMPORTANT: the shelf life of test strips in vials from the first opening of the vial is 180 days.

Storage and transportation conditions

Store the packaging with test strips in a dry, cool place, away from sunlight at temperature 4-40 °C (38-104 °F) and relative humidity 10-85%. Do not refrigerate.

Transportation is carried out by all kinds of covered vehicles in accordance with the working regulations for this mode of transport at temperature 4-40 °C (38-104 °F) and relative humidity 10-85%.

Operation conditions

Use test-strips at temperature from 10 °C to 40 °C (39 °F – 104 °F) and relative humidity 10-85%.

Utilization conditions: The used products are Class B waste and disposed of in accordance with SanPiP 2.1.7.2790-10 "Sanitary - epidemiological requirements for the treatment of medical waste (approved by Decision of the Chief State Sanitary Physician of the Russian Federation of December 9, 2010 №163).

The description of symbols indicated on the test strip packaging is presented in the table below:

Symbol	Description
	For in vitro diagnostic

	Refer to the User's Manual
	Do not reuse
LOT	Lot number
	Caution, read accompanying documents.
	Temperature range / store at
	Expiration date
	Manufacturer
EC REP	EU representative
CE 0123	The product conforms to the requirements of Directive 98/79/EC for in vitro diagnostic medical devices.

The User's manual for test strips was compiled in April, 2016.

Registration certificate FSZ № dated

Manufacturer	Concerning the handling of the medical device on the territory of the Russian Federation contact the authorized representative
L-Tac MediCare Pte Ltd 25 Mandai Estate #07-11 Innovation Place Tower 1 Singapore 729930 Tel: +65 6334 2273 Fax: +65- 6863 8190 Manufacturer's address: OK Biotech Co., Ltd, № 91, Sec. 2, Gongdao, 5 Road, Hsinchu City, 30070, Taiwan ISO 13485: 2003 EC REP Medical Device Safety Service GmbH, Schiffgraben 41, D-30175 Hannover, Germany	105043 Formula-FR LLC, Moscow, Pervomayskaya street, 35/18 Tel.: 8 926 208 71 61

Amos Sim
Business manager
Management Representative

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