Control solution Medicare-Lis
User's manual

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(For the handling of the medical device on the territory of the Russian Federation)

Purpose: Control solution Medicare-Lis is intended for a quality check of a proper operation of a blood glucose meter Medicare-Lis model BG-101 and test-strips Medicare-Lis.

Control measurements should be performed in the following cases:

- when using the meter for the first time. It is necessary to perform the test with the control solution before using your own blood for the analysis. After you have made 3 tests and the result corresponds with the table value, you are ready to perform the analysis with you own blood;
- before using a new package of test-strips;
- if there is a suspicion that the meter and test-strips work improperly;
- if the result of the analysis is not consistent with how you feel or when you think that the results are incorrect;
- if test-strips were exposed to extreme conditions of the environment;
- if you dropped the meter;
- after battery replacement.

Package content:

User's manual, two 4 ml bottles with glucose control solutions (each bottle is intended for more than 100 control measurements).

Control solution 1 with a white cap (Level Normal) for a check in the normal glucose level concentrations interval.

Control solution 2 with a red cap (Level High) for a check in the high glucose level concentrations interval.

You need also: a blood glucose meter Medicare-Lis model BG-101 and test-strips Medicare-Lis, a tissue.

Preparation and control measurement performance: Please, follow the recommendations indicated in this manual and the user's manual for your blood glucose meter Medicare-Lis, model BG-101.

Control measurement procedure is carried out according to the following algorithm:

1. Insert test-strips Medicare-Lis to turn the meter on.
2. After the blood application indicator appears on the screen (●), press the central button and the indicator of CONTROL SOLUTION USAGE (⌚) will appear. The meter will not store these results in the memory if you preset the test as a control solution test. If you decide not to use the control solution mode, press the center button again and the indicator of CONTROL SOLUTION USAGE (⌚) will disappear.
3. Shake the bottle before usage. Squeeze one drop of the control solution out, wipe the dispenser with a tissue to remove remnants. Squeeze the second drop of the solution out on your finger.
4. Apply the drop of the solution from the finger on the absorbing channel of a test-strip Medicare-Lis till a beep sound appears.
5. The result will appear on a display within 6 second. Compare it with an optimal result which should fall within the range of the values of clause 7.3 of the present specification.

Glucose concentration in control solution	mg/dl	mmol/l
Normal level	65-115	3.6-6.4
High level	180-300	10.0-16.7

CERTIFIED TRUE COPY

Victor Koh

6. If the control measurement result does not conform to the optimal results, indicated in 7.3 of the present specification, check the shelf life of the control solution and repeat the measurement again. If the result is still incorrect, it means that a meter and/or test-strips work improperly. **DO NOT USE** a meter and test-strips to monitor your blood glucose level until the sampling results coincide with those from the table. Refer to the User's Manual for a blood glucose meter Medicare-Lis Model BG-101 for more detailed information. If you are unable to resolve the problem by yourself, please, contact the service provider.

IMPORTANT:

- control solution, test-strips and the meter should be at a room temperature of 20-30 °C during control measurements;
- do not add any liquid to the Control solution;
- not for consumption or injection;

Chemical composition of control solution: d-glucose (100 mg/dl (Level I) / 300 mg/dl (Level II)), polyvinyl acetate (aqueous emulsion) 10%, silicone dioxide 0,2%, sodium benzoate 0,2%, disodium EDTA 0,1%, food pigment red №6 0,05%, defoamer (polyethylene glycol 4000) 0,02%

System accuracy when measuring glucose concentration:

Standard deviation ≤ 5 mg/dl at glucose concentration < 75 mg/dl (4,2 mmol)

Coefficient of variation ≤ 5 % at glucose concentration ≥ 75 mg/dl (4,2 mmol)

Storage and transportation conditions:

Store and transport the control solution at temperature from +4 to +30 °C. Control solution should not be cooled, refrigerated, overheated, or exposed to direct sunlight.

Transportation of control solution is carried out by all kinds of covered vehicles in accordance with the working regulations for this mode of transport.

Shelf life: 18 months from the manufacture date.

Shelf life after the first opening of the bottle: 30 days from the first opening of the bottle.

It is recommended to indicate on the bottle when it has been opened for the first time and a monthly interval, during which the solution can be used.

The manufacturer guarantees the stability of the Control solution components till the expiration date if transportation, storage and usage conditions, indicated in the user's manual, are met.

The manufacturer guarantees the safety of the medical device, the absence of unacceptable risk of harm to human life, health and the environment when using the control solution according to its purpose under the conditions stipulated by the manufacturer.

Utilization: The used products are Class B waste and disposed of in accordance with SanPiP 2.1.7.2790-10 "Sanitary - epidemiological requirements for the treatment of medical waste (approved by Decision of the Chief State Sanitary Physician of the Russian Federation of December 9, 2010 №163).

The User's manual for control solution was compiled in April, 2016.

FOR IN VITRO DIAGNOSTICS

Symbol explanations


The following symbols are used on the package, on the device and in the user's manual.

ed in c eat the operly ling re ose n esolve	IVD	For in vitro diagnostic.		Temperature limitations / store at temperature from +4 to +30 °C
		Refer to the User's Manual.		Manufacturer
	LOT	Lot number.	REF	Reference number
uring		Caution, read accompanying documents.	CE 0123	The product conforms to the requirements of Directive 98/79/EC for in vitro diagnostic medical devices.
vel I odiu		Shelf life		

Registration certificate FSZ № dated

Manufacturer	Concerning the handling of the medical device on the territory of the Russian Federation contact the authorized representative
L-Tac MediCare Pte Ltd 25 Mandai Estate #07-11 Innovation Place Tower 1 Singapore 729930 Tel: +65 6334 2273 Fax: +65- 6863 8190 Manufacturer's address: OK Biotech Co., Ltd, № 91, Sec. 2, Gongdao, 5 Road, Hsinchu City, 30070, Taiwan	105043 Formula-FR LLC, Moscow, Pervomayskaya street, 35/18 Tel.: 8 926 208 71 61

Amos Sim
Business manager
Management Representative

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Label (bottle)

**Control Solution MediCare - Lis: Normal Level, High Level
(For the handling on the territory of the Russian Federation)**

Purpose: It is for a quality check of a proper operation of a blood glucose meter and test-strips Medicare-Lis. It is intended for usage as a part of blood glucose monitoring system Medicare-Lis model BG-101.

Precautionary measures: Do not take internally or inject! Use within 30 days after first opening the bottle.

Storage and transportation conditions: Store and transport at temperature from 4 to 30 °C (39 °F – 86 °F).

Utilization conditions: The used products are Class B waste and disposed of in accordance with SanPiP 2.1.7.2790-10 "Sanitary - epidemiological requirements for the treatment of medical waste (approved by Decision of the Chief State Sanitary Physician of the Russian Federation of December 9, 2010 №163).

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Date of opening: _____ (Discard within 30 days after first opening the bottle)

Volume: 4 ml

Registration certificate FSZ № _____ dated _____

CERTIFIED TRUE COPY

Julie Koh

Manufacturer	Concerning the handling of the medical device on the territory of the Russian Federation contact the authorized representative
 L-Tac MediCare Pte Ltd 25 Mandai Estate #07-11 Innovation Place Tower 1 Singapore 729930 Tel: +65 6334 2273 Fax: +65- 6863 8190 Manufacture address: OK Biotech Co., Ltd, № 91, Sec. 2, Gongdao, 5 Road, Hsinchu City, 30070, Taiwan	105043 Formula-FR LLC, Moscow, Pervomayskaya street , 35/18 Tel.: 8 926 208 71 61

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Label (case)
Control Solution MediCare - Lis: Normal Level, High Level
(For the handling on the territory of the Russian Federation)

Purpose: It is for a quality check of a proper operation of a blood glucose meter and test-strips Medicare-Lis. It is intended for usage as a part of blood glucose monitoring system Medicare-Lis model BG-101.

Package content:

User's manual.

Two 4 ml bottles with glucose control solutions (each bottle is intended for more than 100 control measurements).

Precautionary measures:

- The control solution should be held at a room temperature (20-30 °C/68-86 °F) prior to testing.
- Store in a closed bottle at a temperature from 4 °C (39 °F) to 30 °C (104 °F)! Do not refrigerate!
- Use within 30 days after first opening of the bottle!

Composition of control solution: d-glucose, polyvinyl acetate (aqueous emulsion) 10%, silicone dioxide 0,2%, sodium benzoate 0,2%, disodium EDTA 0,1%, food pigment red №6 0,05%, defoamer (polyethylene glycol 4000) 0,02%

Glucose concentration in control solution	mg/dl	mmol/l
Normal level	65-115	3.6-6.4
High level	180-300	10.0-16.7

Storage and transportation conditions: Store and transport at temperature from 4 to 30 °C (39 F – 86 F).

Utilization conditions: The used products are Class B waste and disposed of in accordance with SanPiP 2.1.7.2790-10 "Sanitary - epidemiological requirements for the treatment of medical waste (approved by Decision of the Chief State Sanitary Physician of the Russian Federation of December 9, 2010 №163).

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Manufacturer	Concerning the handling of the medical device on the territory of the Russian Federation contact the authorized representative
L-Tac MediCare Pte Ltd 25 Mandai Estate #07-11 Innovation Place Tower 1 Singapore 729930 Tel: +65 6334 2273 Fax: +65- 6863 8190 Manufacture address: OK Biotech Co., Ltd, № 91, Sec. 2, Gongdao, 5 Road, Hsinchu City, 30070, Taiwan	105043 Formula-FR LLC, Moscow, Pervomayskaya street , 35/18 Tel.: 8 926 208 71 61

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