**User's manual
on medical device
Blood Glucose Monitoring System MediCare-Lis**

**(Concerning the handling of the medical device on the territory of the Russian
Federation)**

Manufactured by L-Tac Medicare Pte., Ltd., Singapore

Victor Koh

[Signature]
Amos Sim
Business Director
25 January 2017

Thank you for choosing Blood Glucose Monitoring System MediCare-Lis. The system was produced in close collaboration with endocrinologists and, what is more important, with people suffering from diabetes.

This User's manual contains the information that you should know about the system. Please, read it carefully. The quick guide is also included for your convenience.

The system is designed to help you and your doctor to monitor your blood glucose levels. All the questions relating to the results of the analysis should be addressed to the attending physician.

This system is the latest development in the sphere of blood glucose monitoring. It is easy to use and it provides fast and accurate results using a minimal blood sample.

Blood Glucose Monitoring System MediCare-Lis:

- It does not require coding: it helps to save time and avoid subjective errors due to coding.
- It allows using an alternative site for the testing (palm).
- It is equipped with memory.

Page 1

IMPORTANT PRECAUTIONARY MEASURES

READ BEFORE USE

These precautions should be followed.

1. Strict control is necessary if the device is used near children, people with physical or mental disabilities and invalids.
2. Use the device only for the intended use described in the manual.
3. Do not use test strips supplied by other manufacturers.
4. Do not use the device if it works improperly or if it is damaged.
5. Before using any product to test blood glucose, read carefully all the instructions and perform the test. Perform all the control measurements as recommended by the endocrinologist.

KEEP THIS USER'S MANUAL WITH YOU

Page 2

Contents

Important information about health.....	5
Purpose of Blood Glucose Monitoring System MediCare-Lis.....	7
Section 1. Getting to know Blood Glucose Monitoring System MediCare-Lis	9
Blood glucose meter Medicare-Lis	9
Meter display.....	10
Lancing device and lancets.....	11
Preparing the lancing device for a test.....	12
Section 2. Performing blood glucose level measurement.....	13
Understanding test results.....	17
Analysis of blood obtained from alternative testing sites	18
Section 3. Meter setting, memory.....	19
Meter setting.....	19
Memory.....	22
Section 4. Control measurements.....	23
Ways to perform a control measurements.....	25

Page 3

Section 5. Maintenance, faults searching and troubleshooting	28
Battery replacement.....	28
Maintenance of the device, storage conditions and service life of the system.....	30
Display messages, faults searching and troubleshooting.....	31
Section 6. Technical information.....	35
Performance characteristics.....	35
Basic parameters and characteristics.....	36
Symbol explanations.....	37
Guarantee liabilities.....	38
Additional expandable materials	40
Guarantee card.....	41

Page 4

Important information about health

1. Apply only capillary blood for the blood glucose testing. The use of other substances, serum or plasma will lead to incorrect results.
2. Severe dehydration and critical water loss can lead to incorrect low results. If you think that you are suffering from severe dehydration, consult your doctor immediately.
3. If you have symptoms which are not consistent with the test results, but you have followed all the instructions described in the user's manual, consult your doctor.
4. Test results below 60 mg/dL (3,3 mmol/L)*1 indicate low blood glucose level (hypoglycemia). Test results higher than 240 mg/dL (13,3 mmol/L)*2 indicate high blood glucose level (hyperglycemia). If your results are below 60 mg/dL (3,3 mmol/L) or higher than 240 mg/dL (13,3 mmol/L), repeat the test and if your results are still below 60 mg/dL (3,3 mmol/L) or higher than 240 mg/dL (13,3 mmol/L), immediately consult your doctor.
5. The incorrect result can occur in severely hypotensive persons or patients in a state of shock. Incorrect low results can occur in persons in a hyperglycemic-hyperosmolar state with ketosis or without it. It's not recommended to test people having incurable diseases with this meter.
6. Abnormal red blood cell count (hematocrit level below 30 % or higher than 55%) can lead to incorrect testing results. Please, consult your doctor if you don't know your hematocrit level.

Page 5

7. Interactions with other medical products: Reducing substances entering blood naturally (uric acid, bilirubin) or through therapeutic treatments (ascorbic acid, acetaminophen) don't affect the test results significantly. However, high concentration of these substances can affect the results. The substances and their acceptable concentrations are presented below.

Substances	Uric Acid	Gentisic Acid	Acetaminophen	Hydroxyurea	Ascorbic Acid
The indicator below which the substance content in the blood is acceptable, and it does not affect the test result	≤8.0 mg/dL (0.48 mmol/L)	≤5.0 mg/dL (0.32 mmol/L)	≤8.0 mg/dL (0.53 mmol/L)	≤3.0 mg/dL (0.39 mmol/L)	≤5.0 mg/dL (0.28 mmol/L)

Substances	Dopamine	Dopa	Methyldopa	Tolbutamide	Triglycerides
The indicator below which the substance content in the blood is acceptable, and it does not affect the test result	≤2.0 mg/dL (0.11 mmol/L)	≤10 mg/dL (0.51 mmol/L)	≤3.0 mg/dL (0.13 mmol/L)	≤400 mg/dL (14,8 mmol/L л)	≤2000 mg/dL (22,8 mmol/L)

Substances	Aspirin	Bilirubin	Cholesterol	Creatinine	Galactose (milk sugar)	Maltose (malt sugar)
The indicator below which the substance content in the blood is acceptable, and it does not affect the test result	≤60 mg/dL (3,33 mmol/L)	≤90 mg/dL (1,54 mmol/L)	≤500 mg/dL (12,9 mmol/L)	≤5 mg/dL (0.44 mmol/L)	≤900 mg/dL (50 mmol/L)	≤900 mg/dL (26,3 mmol/L)

8. Please, look through the additional important information on the test-strips package insert/instruction.

Page 6

Purpose of Blood Glucose Monitoring System MediCare-Lis

Purpose

The blood glucose monitoring System MediCare-Lis is intended for monitoring of capillary blood glucose level outside the body (for in vitro diagnostic only) at home by people with diabetes or with a suspicion of diabetes.

Principals of test performance

Blood glucose level is measured by an electrical current that is produced by a blood sample mixing with the reagent (special chemicals) on the test strip. The electrical current changes depending on the amount of glucose in the blood sample. The meter measures the strength of the electrical current, calculates the blood glucose level and then displays your result in either mg/dL or mmol/L.

About the system

The system Medicare-I is uses the latest technology to provide you with easy and comfortable testing. It requires only a 0.3 μ L of blood sample and 6 seconds for measurement and getting the result.

Page 7

The composition of the system:

1. Blood Glucose Meter MediCare-Lis, model BG-101 – 1 pcs.
2. Alkaline battery - 1 pcs.
3. Lancing device MediCare - 1 pcs.
4. Sterile lancets for single use MediCare 30G – 10 pcs.
5. Carrying bag - 1 pcs.
6. User's manual for system - 1 pcs.
7. Quick guide – 1 pcs.

All the products have been tested and proven to work together for getting accurate measurement results.

CAUTION:


1. Use only test strips MediCare-Lis and control solution MediCare-Lis with the blood glucose meter MediCare-Lis, model BG-101. The use of test strips or control solutions of other brands may lead to incorrect results.
2. Users should not make medical decisions without first consulting with their doctor.
3. If you experience symptoms which are not consistent with the test results, see a doctor immediately.
4. The blood sample with a high content of reducing substances such as ascorbic acid and uric acid can give the result which is slightly higher than the actual glucose level.
5. The altitude above 3,402 meters may affect the accuracy of measurements.
6. The temperature beyond the range of 10°C to 40°C may affect the measurement results.
7. The magnetic fields, external electrical effects, electrostatic discharges, pressure or differential pressure, overload of thermal ignition sources can affect the measurement results.

Page 8

Section 1. Getting to know Blood Glucose Monitoring System MediCare-Lis

Blood glucose meter MediCare-Lis

Study the scheme and get to know all the components of your meter.

1. **LCD DISPLAY:** The large, easy to read display shows test results, messages, as well as test results stored in meter memory.
2. **CENTRAL BUTTON:** The main button is for turning the meter on or performing other functions.
3. **STRIP PORT:** Insert the test strip here. The meter will turn on automatically.
4. **SETTING BUTTON:** It is located on the battery compartment; press the button to enter the meter settings.
5. **BATTERY COMPARTMENT:** It contains ONE 3v Lithium battery, battery type is CR2032.
6. **METER MARKING:** Each meter has its unique number. It can be used when getting the warranty service. Please, do not alter or tear the label with the unique number off.

Page 9

METER DISPLAY

1. **DATE:** shows year, month and date.
2. **TIME:** the meter is programmed for a 12-hour format.
3. **AVERAGE VALUE:** appears when the meter shows 7/14/28 test average values in the memory mode.
4. **MEMORY SYMBOL:** appears in the memory mode.
5. **THERMOMETER SYMBOL:** appears when the ambient temperature is above or below the optimal range required for testing.
6. **CONTROL SOLUTION SYMBOL:** appears in Control Solution Mode. Your test result will not be stored in meter memory.
7. **MEASUREMENT UNITS:** appear with the test result either in mg/dL or mmol/L.
8. **BATTERY SYMBOL:** appears when battery is weak.
9. **BLOOD DROP SYMBOL:** flashes when it's necessary to apply a sample on the sensor.
10. **TEST STRIP SYMBOL:** appears when it's necessary to insert a test strip.

Page 10

Lancing device and lancets

The lancing device and lancets are used for getting samples of capillary blood from the puncture site.

Lancing device

1. Lancing device cap
2. Wheel of puncture depth selection
3. Lancet holder
4. Lancet's trigger button
5. Lancet ejector
6. Sliding barrel
7. Lancet
8. Lancet protective cap

Important Information about the lancets and lancing device

1. Ⓢ Lancets are for single use only.
2. Keep the lancing device and lancets clean during the shelf life period indicated on the package.
3. Use caution when removing and disposing the lancet from the device.
4. The meter and the lancing device should be used by one person only.
5. **DO NOT ALLOW** to use them to third parties, even if they are members of your family!
6. **NOT** for shared use!

Page 11

Preparing the lancing device for a test

1. Screw off the cap of the lancing device. Insert a lancet into the lancet holder and push down as far as it will go.

2. Twist off the protective cap.

3. Replace the cap on the lancing device and set the desired puncture depth

Important: the wheel of puncture depth selection provides 5 levels of skin penetration:

- 1-2 for soft or thin skin
- 3 – for average skin thickness
- 4-5 for thick or calloused skin

4. Cock the device - pull back the sliding barrel until a click and release it. If there is no a click, It may have been cocked when the lancet was inserted.

Page 12

Section 2. Performing blood glucose level measurement

Step 1. Wash your hands and the testing site

Wash your hands with water and soap. Rinse and dry completely. Warm your fingers for a rush of blood.

Step 2. Insert a Test Strip

Remove a new test strip from the package. Insert the test strip immediately into the strip port as indicated in the picture. The meter turns on automatically. When the blood symbol starts blinking, you can perform a test.

Step 3. Selection and puncture of the testing site

* For a Fingertip

Hold the prepared lancing device firmly against the side surface of your fingertip. Press the release button of a lancet.

* For alternative sites (palm or forearm)

Please, refer to the "Analysis of blood obtained from alternative testing sites" section. Consult your doctor before using the alternative site testing.

Page 13

Step 4. Obtain a Blood Sample

Massage carefully your finger or testing site to obtain the required blood volume. To perform the test, you will need only 0.3µL of blood. Do not smear the blood sample. To obtain the accurate result, wipe off the first blood drop and carefully squeeze another one.

Step 5. Apply blood on the test strip

Apply the blood to the tip of the test strip absorbent channel until the confirmation window is fully covered with blood. Then the meter begins to count down and shows up the test result in 6 seconds.

Step 6. Watch the result

The meter will display the blood glucose value in 6 seconds.

Page 14

Step 7. Remove the test strip to turn the meter off and dispose of it

Your test result is automatically stored in the meter memory. Turn the meter off by removing the test strip. Dispose of the used test strips in a proper way.

Step 8. Protecting the used lancet

Remove the lancing device cap and put a protective element on the sharp lancet tip.

Step 9. Dispose of the used lancet

Push the lancet ejector forward and remove the lancet. Close the lancet with a protective cap and place it in a sharp container for further disposal.

IMPORTANT:

1. The system needs only a small blood sample to perform a test. Choose different spots every time you perform a test. The test performance on the same spot may lead to soreness.
2. The first blood sample usually contains tissue and serum, which may affect the test result. It should be removed.
3. Blood should completely fill the confirmation window before the meter starts to count down. If you find out that the confirmation window is not filled with blood and the meter has started the countdown, never try to add more blood to the test strip. Discard the test strip and perform the test again.
4. If you do not apply a blood drop within three minutes, the meter will automatically turn off. You must remove the test strip, insert it again in the meter and restart the test performance.
5. If you have problems using test strips, contact your local distributor.

Understanding test results

Normal blood glucose level

Normal blood glucose level for fasting people, not suffering from diabetes, varies from 70 to 110 mg/dl (3.9 - 6.1 mmol/l) and is less than 140 mg/dl (7.8 mmol/L) when performing a test in 2 hours after a meal. Consult your doctor to determine your normal blood glucose level.

If the blood glucose level seems to be too high or low, or inconsistent with previous results, check the following:

1. Did you apply the blood drop on the test strip quickly enough after removing it out of the package?
 2. Was the volume of the blood sample sufficient?
 3. Was the test strip package damaged?
 4. Was the test strip used before the expiration date?
- Were the test strips stored away from extreme high or low temperatures or from areas of high humidity?

Perform a control test with the control solution Medicare-Lis (Medicare-Lis). If the control test provides results within acceptable values, repeat the analysis with your blood, using a new test strip. If the result still does not match your previous results, the normal blood glucose level or does not match your state of health, immediately consult a doctor.

Page 17

Analysis of blood obtained from alternative testing sites

Blood samples can be used not only from a fingertip, but also from other parts of the body (palms and forearms) for the analysis of blood glucose.

Attention! Please, consult your doctor before using the alternative site testing!

What are the advantages?

Fingertips are more sensitive because they are full of nerve endings (receptors). On other body parts, nerve endings are not so numerous and you will not feel such soreness as taking the sample from fingertips.

When to use alternative site testing?

- Immediately before meal or in a fasting state;
- in two hours or more after taking insulin.;
- in two hours or more after exercise;

*When it is **NOT ALLOWED** to use alternative site testing?*

- during a meal or immediately after a meal;
- after exercise;
- after a nervous breakdown;
- if the results of your regular measurements often vary;
- you have reason to believe that you have hypoglycemia or hyperglycemia.;
- you are pregnant

Important: It's necessary to rub the testing site for 20 minutes before performing the test to increase the accuracy of measurement results.

Page 18

Section 3. Meter setting, memory

Meter setting

It is necessary to insert the battery before use (see Section 5) and configure the device in the following sequence

Step 1. Setting


Open the battery compartment in the turned off mode. Press the SETTING BUTTON located at the right side of the battery compartment. The meter will turn on and will be ready for the setting.

Step 2. Set the year

The digit, denoting the year, appears and starts flashing. Press and release the central button to set the desired year, press the setting button to confirm and proceed to the month setting.

Step 3. Set the month

The digit, denoting the month, appears and starts flashing. Press and release the central button to set the desired month, press the setting button to confirm and proceed to the date setting.


Step 4. Set the date

The digit, denoting the date, appears and starts flashing. Press and release the central button to set the desired date, press the setting button to confirm and proceed to the hour setting.

Step 5. Set the hour

The digit, denoting the hour, appears and starts flashing. Press and release the central button to set the desired hour, press the setting button to confirm and proceed to the minutes setting.

Step 6. Set the minutes

The digit, denoting minutes, appears and starts flashing. Press and release the central button to set the desired minute, press the setting button to confirm and proceed to the following settings.

Step 7. Set the units of measurement

Default units (for example, mg/dL) will appear and will be flashing, press and hold the central button for 3 seconds to change the units to mmol/L. Press the setting button to confirm and proceed to the following settings.

mg dL

Step 8. Delete memory

When del. icon appears on the screen and blinks, you can choose the memory clearing. If you do not want clearing the memory, press the setting button again. If you want to delete all the previous results stored in the memory, press and hold the central button for 3 seconds. The icon "----" on the display means that the memory is cleared.

Step 9. Complete setting

Prior to the device's turning off, the icon «OFF» will appear and then the display will blink off. Installation is complete.

IMPORTANT:

1. The time, date and units of measurement can be changed **ONLY** in the setting mode. Therefore, it is not possible to change these parameters during the test performance.
2. The meter displays the average value for 7, 14 and 28 days, the access to this data is possible via meter memory.
3. The meter can display the results in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L). When you turn off the meter and replace the batteries, the units of measurement will not change.
4. The use of wrong units of measurement may lead to misreading of your analysis results, and, correspondingly, to incorrect treatment. Please, consult your doctor before you select the units of measurement.

Memory

The blood glucose meter Medicare-Lis automatically stores 120 results of previous measurements, letting you review them in order from the most recent to the oldest one. The meter also calculates and displays the average value for 7, 14 and 28 days. You can review the average test results or an individual result by entering the memory mode.

Memory mode

Step 1. Enter the memory review mode

Press the central button to turn the meter on. When the test strip blinking symbol appears, press the central button again to enter the memory mode.

Step 2. The average measurement result review.

After entering the memory mode, the average result for 7 days will appear. When pressing the central button again, the average result for 14 days and 28 days will appear, correspondingly.

Step 3. Individual result review.

After displaying the average value for 28 days, press the central button, the result and the time of the most recent measurement will be displayed. Continue pressing the central button for the consequent review of individual results. When the memory is full, the meter deletes the oldest result and adds a new one. When the latest value is displayed, the meter will display the average for 7 days again.

OFF

Step 4. Exit the memory mode.

If you want to exit the memory mode, press and hold the central button for three seconds to turn off the meter.

Note:

1. The control solution results are NOT stored in the meter memory if you set the meter correctly before the control test. (Refer to the "CONTROL MEASUREMENTS" section).
2. When using the meter for the first time, the icon "----" is displayed. When you recall test results or review the average result, this icon means that there are no results.
3. The average values are calculated from the latest results for 7, 14 and 28 days.
4. You can exit the memory mode at any time and turn off the meter by holding the central button for three (3) seconds.
5. If you leave the meter in the memory mode for more than one (1) minute, it will turn off automatically.

Section 4. Control measurements

Control solution MediCare-Lis

Control solution contains a known amount of glucose that reacts with test strips. Performing the test with the control solution, comparing the result with the table range on the control solution package, it's possible to check if the meter and the test strips work properly, as well as the accuracy of the presented results.

Control measurements should be performed in the following cases:

- when using the meter for the first time. It is necessary to perform the test with the control solution before using your own blood for the analysis. After you have made 3 tests and the result corresponds with the table value, you are ready to perform the analysis with you own blood;
- before using a new package of test-strips;
- if there is a suspicion that the meter and test-strips work improperly;
- if the result of the analysis is not consistent with how you feel or when you think that the results are incorrect.
- if test-strips were exposed to extreme conditions of the environment;
- if you dropped the meter;
- after battery replacement.

Important information!

1. Use only the control solution MediCare-Lis
2. Check the expiration date on the package. Do not use if expired.
3. The control solution, test strips and the meter should be at room temperature (20-30°C/68-86°F) when testing.
4. Shake the bottle before usage, remove the first drop of the control solution, and wipe off the dispenser to remove the remnants. It is necessary to get a good sample and an accurate result.
5. Store the control solution closed at temperature below 30°C (86°F). Do not refrigerate.

Ways to perform a control measurement

Step 1. Insert the test strip into the port

Insert the test strip into the port to turn on the meter.

Step 2. Proceed to Control solution check mode

After the blood application indicator appears on the screen (🩸), press the central button and the indicator of CONTROL SOLUTION USAGE (💧) will appear. The meter will not store these results in the memory if you preset the test as a control solution test. If you decide not to use the control solution mode, press the center button again and the indicator of CONTROL SOLUTION USAGE (💧) will disappear.

Step 3. Squeeze a drop of control solution

Shake the control solution bottle well. Remove the cap. Squeeze the first drop of the solution out and wipe it off from the dispenser with a clean paper or cloth. Squeeze the second drop out on your finger tip.

Page 25

Step 4. Apply the control solution

Apply the drop of the solution to the strip absorbent channel until the confirmation window is filled. The meter begins to count down.

Step 5. Check if the results are in the indicated range

After the meter counts down from 6 to 1, the result will appear on the screen. Compare it with the optimal result indicated on the test strip package.

Page 26

NOTE:

1. Make sure that the meter is in the Control Solution Mode so that the test result will not be stored in the meter memory.
2. To avoid contaminating the control solution with the content of the test strip, **DO NOT APPLY THE CONTROL SOLUTION DIRECTLY TO THE TEST STRIP!**
3. If the test result doesn't coincide with the information indicated on the test strip package, repeat the test again. If the result is still incorrect, it means that the meter and the test strip work improperly. **DO NOT** use the meter and the test-strips to monitor your blood glucose level until the sampling results coincide with those from the table. If you are unable to resolve the problem by yourself, please, contact your local distributor.

Page 27

Section 5. Maintenance, faults searching and troubleshooting

Battery replacement

The meter requires one lithium battery, type CR2032 (included with the purchase).

The battery service life depends on meter usage conditions. Always keep a spare battery at hand. The battery service life is 1000 tests or 12 months when testing 3 times a day.

If the battery indicator appears on the display, it means that the battery charge is low.

You can perform the tests with battery low charge but it will be necessary to replace it soon.

When the battery indicator flashes simultaneously with E-b, the meter will not give results and it is necessary to replace the battery immediately. Always have one spare battery with you to replace it anytime.

Ways of battery replacement

1. Make sure that the meter is turned off. Turn the device, press the buckle on the board covering the battery and lift it up. Remove the battery.

Page 28

2. Install the battery with the (+) side facing upward, push it until it is installed into the groove. In case of proper installation, you will hear a beep. Otherwise, remove the battery and install it again.

3. Close the battery cover.

NOTE:

1. The battery replacement does not affect the memory. However, the time and the date may need to be reset.
2. As with all small objects, the batteries should be stored out of children's reach. If the battery is swallowed, see a doctor immediately.
3. Batteries can spill chemicals if they not used for a long time. Remove the batteries of the meter if you are not going to use it for a long period of time (3 months or more.)
4. Discard the used or dead batteries according to the state regulations for electronic products recycling.

Page 29

Maintenance of the device, storage conditions and service life of the system

To avoid contamination of the device wash and dry your hands thoroughly before usage.

Cleaning

The device does not require any special maintenance. Any special methods of device cleaning are not provided, except cases when blood or Control Solution falls on it.

1. To clean the outer surface of the device wipe it with a cloth moistened with water or mild detergent, then wipe dry with a soft cloth. Do not flush with water.
2. Do not use organic solvents to clean the device, handle it carefully.

Storage and transportation conditions of the system

- Storage/transportation conditions of the system: 4°C~40°C (39°F~104°F) at relative humidity 10-85%.
- Storage/transportation conditions of the blood glucose meter: -20°C~50°C (-4°F~122°F) at relative humidity not higher than 90%.
- Always store and transport the device in its original packaging.
- Protect from direct sunlight and humidity.
- Protect from drops and strong shocks.
- The device is transported by all kinds of covered vehicles except unheated compartments of aircrafts in accordance with the rules approved in the established procedure for each mode of transport provided the storage conditions are observed.

System operational conditions: 10°C~40°C (50°F~104°F) at relative humidity 10-85%.

Page 30

Utilization conditions: The used products are Class B waste and disposed of in accordance with SanPiP 2.1.7.2790-10 "Sanitary - epidemiological requirements for the treatment of medical waste (approved by Decision of the Chief State Sanitary Physician of the Russian Federation of December 9, 2010 №163)

Service life of the device: 5 years

Display messages, faults searching and troubleshooting

The table of display messages and symbols is presented below. These messages help to identify the source of a problem but data do not always appear when a problem has occurred. Improper use may lead to incorrect results without displaying an error message. In the case of revealing a problem, refer to the column "actions".

Display	Description	Actions
	Display check	If some parts of the display are damaged, contact your local distributor.
	Moving strip	It is necessary to insert the test strip into the meter.
	Blinking blood drop	The meter is ready for a contact with a blood sample.
	Deleting memory	Deleting is complete.
	Test result is higher than 600 mg/dl (33.3 mmol/l)	If it is not consistent with the way you feel, perform a control test. Repeat the test with blood, if the result doesn't change, call a doctor immediately.


The result is lower than 20 mg/dl
(1,1 mmol/l)

If it is not consistent with the way
you feel, perform a control test.
Repeat the test with blood, if the
result doesn't change, call a doctor
immediately.

Low battery charge


Replace the battery

The battery is dead


Replace the battery immediately

Page 33

It is possible that:
1. Wet test-strip is used
2. The meter is broken

It is necessary:
1. Repeat the test with a new test-strip
2. Contact your local distributor

Temperature is out of the
operating range.

The meter is not working. Transfer the
meter to a room with temperature between
10°C - 40°C (50°F - 104°F). Wait for, at
least, 30 minutes. Do not heat or cool the
meter on purpose.

There are no
responses when a
test strip is
inserted into the
meter

It is possible that:
1. Battery is dead
2. Wrong strip is inserted
3. The meter is defective

It is necessary:
1. Replace the battery
2. Insert the test-strip correctly
3. Contact your local distributor

There are no
responses when a
blood sample is
applied to the test
strip

It is possible that:
1. Blood sample is not
sufficient
2. The meter is defective

It is necessary:
1. Repeat the test with a sufficient blood
sample.
2. Contact your local distributor

Page 34

Section 6. Technical information

Performance characteristics

Accuracy of measurements

> 95% in A zone of EGA vs YSI 2300 glucose analyzer:
Within 15 mg/dL (0.83 mmol/L) at glucose concentration < 75 mg/dL (or 4.2 mmol/L) and within 20% glucose concentration $\geq$ 75 mg/dL (or 4.2 mmol/L).

Precision of measurements

CVs(%) of intermediate precision and repeatability were less than 5%.
The device has been certified according to the following standards:

EN ISO 13485:2003, EN ISO 14971:2007-10-01, EN ISO 15197:2003, EN 980:2008, EN ISO 18113-4:2009, IEC/EN60601-1, IEC/EN60601-1-2, EN61010-1, EN61010-2-101:2002, EN 61326-1, Pr ISO 17511, EN 61326-2-6

Basic parameters and characteristics

Model name	BG-101
Analysis method	Electrochemical biosensor
Sample type	Capillary blood
Test result	Based on the measurement of glucose level
Alternative site testing	YES (palm, forearm)
Sample size	0.3 ml
Result processing time	6 seconds
Measuring range (sensitivity limit)	20-600 mg/dL (1.1-33.3 mmol/L)
Sensitivity in HCT	30~55%
Storage transportation conditions of the system	4°C~40°C (39°F~104°F) at relative humidity 10-85%.
Storage transportation conditions of the device	-20°C~50°C (-4°F~122°F) at relative humidity not higher than 90%.
Operational conditions of the system	10°C ~40°C (50°F ~104°F) at relative humidity10-85%.
Memory	120 results with date and time
Average value for	7, 14, 28 days
Power supply	One 3-volt Lithium Battery (battery type CR2032)
Battery service life	About 1000 measurements
Size of a model	80mm (length) x 48 mm (width) x 15mm (height)
Weight	39 g

Specifications are subject to change without prior notice to improve performance characteristics and qualities of the product.

Symbol explanations

The following symbols are used on the package, on the device and in the user's manual.

	For in vitro diagnostic.		Temperature limitations / store at +4 - +40°C
	Refer to the User's Manual.		Use by /Expiry date.
	Do not reuse.		Manufacturer.
	Lot number		Europe representative.
	Caution, read accompanying documents.		The product conforms to the requirements of Directive 98/79/EC for in vitro diagnostic medical devices.
	Do not dispose with household waste. Subject to special recycling.		Reference number

Registration certificate №

Declaration...

Last revision: 02.2017

Page 37

Guarantee liabilities

The manufacturer guarantees the compliance of the system with the stated parameters if operation, transportation and storage conditions, specified in this manual, are observed.

Guaranteed service life of the device: 24 months from the day of purchase (mark in the guarantee card).

To appeal to the guarantee service, it is necessary to provide the system in its original packaging or in packaging with the appropriate degree of protection, cashier's check or fully filled guarantee card with the stamp of commercial enterprise or technical service center, as well as the written notification of the manufacturing defect of quality or malfunction.

Repair or replacement under the guarantee does not give rise to any extension or renewal of the guarantee period.

Guarantee liabilities are valid, if the device is intact and was operated according to the user's manual and the descriptions given in the manual.

Guarantee liabilities does not cover consumables, device packaging, damage of display and fragile parts (namely, accidental fall of the device or drop of a heavy object on the device), defects caused by repair in the unofficial representative.

Page 38

The guarantee is not performed in the presence of traces of mechanical impact, dents, cracks, chips, etc on the device case (display included); in the presence of traces of opening the case, traces of attempts of repair in an unauthorized service center; in the presence of traces of ingress of moisture inside the casing, or exposure of the device to aggressive environment; in other cases of violation by a consumer of storage regulations, transportation and technical operation of the device, provided by the regulations, indicated in this user's manual.

Manufacturer and distributor/reseller has the right to deny in guarantee service if, during the guarantee period, the system was used not in accordance with the intended use or with violations of the operational requirements described in this user's manual, and if incorrect information was provided. In this case, the repair is made at the expense of the consumer.

WARNING: This device may not suit all users. It is recommended to consult your doctor for your safety.

Additional expandable materials

You will be able to purchase the following additional expandable materials in pharmacies or in medical equipment stores.

Test-strips

Blood Glucose Test-strips MediCare-Lis

Lancets

Sterile lancets for single use MediCare 30G

Control solution

Control solution MediCare-Lis

Guarantee card

Blood Glucose Monitoring System MediCare-Lis

Model number: BG-101

Device serial number: _____

User's name: _____

User's address: _____

User's phone number: _____

User's e-mail: _____

Failure/error description: _____

Place of purchase: _____

Date of purchase: _____

Signature: _____

Stamp of
an agent

Page 41

Manufacturer

L-Tac MediCare Pte Ltd

25 Mandai Estate #07-11

Innovation Place Tower 1

Singapore 729936

Tel: +65 6334 2273

Fax: +65- 6863 8190

Manufacture address:

OK Biotech Co., Ltd. № 91, Sec. 2,

Gongdao, 5 Road, Hsinchu City, 30070,

Taiwan

Concerning the handling of the medical device
on the territory of the Russian Federation
contact the authorized representative

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Amos Sim
Business Director
25 January 2017

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1 CERTIFIED TRUE COPY 1

Victor Koh Koon Hor

**Blood Glucose Monitoring System
MediCare-Lis**
(Concerning the handling of the medical device on the
territory of the Russian Federation)

Manufactured by L-Tac Medicare Pte., Ltd.,
Singapore

Amos Sim

Amos Sim
Business manager
Management Representative


Blood Glucose Monitoring System MediCare-Lis BG 101 Quick guide

Quick guide will help to set up the system and perform blood glucose measurements successively. You will also learn from the guide about additional functions of the system, about what the measurement results mean, about error messages, etc.

You need for measurements

Blood glucose meter

Test strips

Lancing device

How to perform measurements properly


1. Wash your hands, rinse them and dry completely before the beginning of the measurement.

2. Remove a new test strip from the package. Insert the test strip into the meter. The meter turns on automatically.

3. Use the lancing device to obtain the necessary blood amount which is 0,3 ml.

4. Massage carefully your finger or testing site to obtain the required blood volume. To perform the test, you will need only 0.3µL of blood. Do not smear the blood sample. To obtain the accurate result, wipe off the first blood drop and carefully squeeze another one.

5. Apply blood on the test strip.

6. The result will appear on the display in 6 seconds.

7. Remove the used test strip.

8. Remove the lancing device cap and put a tip of the used lancet on its protective cap.

9. Push the lancet ejector forward and remove the lancet. Dispose the lancet according to the safety regulations.

*Do not reuse the lancet.

The MediCare system is intended for capillary blood diagnostics outside the body (in vitro). Do not use plasma or serum.

2. Do not use for blood testing of neonates.
3. ⓧ Don't use the test strips more than once.
4. Do not use expired test strips.
5. Always store test strips in the original package.
6. Dehydration can lead to low results. If your body is severely dehydrated, contact your doctor immediately.
7. Inaccurate results are possible in patients with extremely low pressure or who are in a state of shock. Mistakenly low results are possible in persons in a hyperglycemic-hyperosmotic state with ketosis or without it. It is not allowed to perform tests to determine the glucose level in patients in critical condition.

Maintenance of the device, storage conditions and service life of the device

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1. To clean the outer surface of the device wipe it with a cloth moistened with water or mild detergent, then wipe dry with a soft cloth. Do not flush with water.

2. Do not use organic solvents to clean the device, handle it carefully.

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Service life of the device: 5 years

Guaranteed service life: 24 months from the day of purchase (mark in the guarantee card)

The time of day	Result in persons not suffering from diabetes
In the fasted state or before meals	Less than 6.1 mmol/l (110 mg/dl)
In 2 hours after meals	Less than 7.8 mmol/l (140 mg/dl)

Source: American Diabetes Association Clinical Practice Recommendations 2003

The time of day	My normal values
In the fasted state or before meals	
In 2 hours after meals	

Manufacturer

L-Tac MediCare Pte Ltd
25 Mandai Estate #07-11
Innovation Place Tower 1
Singapore 729930
Tel.: +65 63342273
Fax: +65- 68638190
Manufacturer's address:
OK Biotech Co., Ltd. 83,
Sec. 2, Gongdao, 5 Road,
Hsinchu City, 30070, 
Taiwan

Distributor in Russia, the organization authorized to receive claims

105043 Formula-FR LLC,
Moscow, Pervomayskaya
street, 35/18
Tel. **8 926 208 71 61**

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