APPROVED BY

Sophie REHAULT

Quality and Regulatory Affairs Director

BALT EXTRUSION

March 12th 2018**BALT EXTRUSION SAS**

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95160 MONTMORENCY


FRANCE

RCS Pontoise 315 633 081

INSTRUCTION FOR USELEO+ self-expandable intracranial stents
(BALT EXTRUSION, FRANCE)

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M. 

(seen exclusively to certify the above signature)

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2018

Name of the device:

LEO+ self-expandable intracranial stents

Design variants:

I. LEO+ self-expandable intracranial stent in the following versions:

1. LEO+ self-expandable intracranial stent in the following versions: LEO.3.5x12, LEO.3.5x18, LEO.3.5x25, LEO.3.5x30, LEO.3.5x35, LEO.3.5x50, LEO.4.5x15, LEO.4.5x20, LEO.4.5x25, LEO.4.5x30, LEO.4.5x40, LEO.4.5x50, LEO.4.5x75, LEO.5.5x25, LEO.5.5x30, LEO.5.5x35, LEO.5.5x50, LEO.5.5x60, LEO.5.5x75.
2. LEO+ Baby self-expandable intracranial stent, in the following versions: LEO.2.0x12, LEO.2.0x18, LEO.2.0x25, LEO.2.5x12, LEO.2.5x18, LEO.2.5x25, LEO.2.5x30, LEO.2.5x35, LEO.3.0x12, LEO.3.0x18, LEO.3.0x25, LEO.3.0x35.

Device delivery set:

LEO+ self-expandable intracranial stent

1. LEO.3.5x12, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x12) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F)) - 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

2. LEO.3.5x18, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x18) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F)) - 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

3. LEO.3.5x25, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x25) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F)) - 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

4. LEO.3.5x30, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x30) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F)) - 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

5. LEO.3.5x35, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x35) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F)) - 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

6. LEO.3.5x50, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.3.5x50) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+21 (2.4F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

7. LEO.4.5x15, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x15) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc (3 stickers).

8. LEO.4.5x20, comprised of:-

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x20) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

9. LEO.4.5x25, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x25) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

10. LEO.4.5x30, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x30) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

11. LEO.4.5x40, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x40) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

12. LEO.4.5x50, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x50) - 1 pc.
- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F))– 1 pc.
- Y-connector with hemostatic valve - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

13. LEO.4.5x75, comprised of:

- LEO+ self-expandable intracranial stent on the delivery system (LEO.4.5x75) - 1 pc.

- Reinforced microcatheter Vasco+ (Vasco+25 (3.0F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
14. LEO.5.5x25, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x25) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
15. LEO.5.5x30, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x30) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
16. LEO.5.5x35, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x35) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
17. LEO.5.5x50, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x50) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
18. LEO.5.5x60, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x60) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
19. LEO.5.5x75, comprised of:
- LEO+ self-expandable intracranial stent on the delivery system (LEO.5.5x75) - 1 pc.
 - Reinforced microcatheter Vasco+ (Vasco+28 (3.3F)) - 1 pc.
 - Y-connector with hemostatic valve - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).

LEO+ Baby self-expandable intracranial stent:

1. LEO.2.0x12, comprised of:
- LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.0x12) - 1 pc.
 - Rotating device (torquer) - 1 pc.

- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).
- 2. LEO.2.0x18, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.0x18) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 3. LEO.2.0x25, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.0x25) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 4. LEO.2.5x12, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.5x12) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 5. LEO.2.5x18, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.5x18) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 6. LEO.2.5x25, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.5x25) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 7. LEO.2.5x30, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.5x30) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 8. LEO.2.5x35, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.2.5x12) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 9. LEO.3.0x12, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.3.0x12) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 10. LEO.3.0x18, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.3.0x18) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.
 - Label with information about the implantable device - 1 pc. (3 stickers).
- 11. LEO.3.0x25, comprised of:
 - LEO+Baby self-expandable intracranial stent on the delivery system (LEO.3.0x25) - 1 pc.
 - Rotating device (torquer) - 1 pc.
 - User manual - 1 pc.

-Label with information about the implantable device - 1 pc. (3 stickers).

12. LEO.3.0x35 comprised of:

- LEO+Baby self-expandable intracranial stent on the delivery system (LEO.3.0x35) - 1 pc.
- Rotating device (torquer) - 1 pc.
- User manual - 1 pc.
- Label with information about the implantable device - 1 pc. (3 stickers).

Manufacturer: (BALT EXTRUSION, FRANCE)

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Development, manufacture, and sales of the device comply with ISO 13485: 2003/NF-EN ISO 13485 :2012 (certificate number: N° 9500 Rev. 4, issue date-03/03/2015, expiry date-02/03/2018, notified body: G-MED certification/LNE/G-MED Notified Body for Medical Devices)

Russian National Classification of Products by Economic Activities Code OK-034-2014:
32.50.22.190

Class depending on the degree of potential risk of medical use in accordance with the medical devices nomenclature classifier approved by the Order of the Ministry of Health of the Russian Federation No. 4H-3 dated June 6, 2012.

Device classification according to the Medical Device Nomenclature Classification stipulated by the Order of the Ministry of Health of the Russian Federation of June 6, 2012, No.4H –218150.

Type of body contact: Implanted stent - continuous contact with blood; stent delivery system - reinforced microcatheter – short-term contact with blood; torquer, Y-connector with a hemostatic valve, protective rim, terminal sterilization package (primary package of the device), individual package of torquer, individual package of reinforced microcatheter Vasco+, secondary package (cardboard box), transport package – short-term mediated contact with intact skin.

Conditions of use: medical and medical-preventive health care institutions (installation of the stent must only be exercised in healthcare institutions, where emergency medical aid can be provided). The product should only be used by experienced medical specialists in the field of interventional neuroradiology and the treatment of intracranial aneurysms. LEO+ and LEO+ Baby stents are intended for implantation to patients with cerebral aneurysms, regardless of gender and age (except for newborns, premature infants and infants) and who have no contraindications.

Purpose

The product is intended for implantation into the base or the parental artery of the intracranial aneurysm in order to ensure the delivery of embolizing agents to fill the aneurysm, provide blood coagulation within the aneurysm by slowing the blood flow within the aneurysm and / or providing support inside the aneurysm to the spiral for embolization of cerebral vessels.

Device special features

Description and design

LEO+ self-expandable intracranial stent (hereinafter referred to as LEO+ stent) and LEO+ Baby self-expandable intracranial stent (hereinafter referred to as LEO+ Baby stent) are self-expanding braided stents composed of 16 thin nitinol wires. Two longitudinal radiopaque platinum markers are placed along the entire length of LEO+ and LEO+ Baby, which improves its visualization and control. The LEO+ and LEO+ Baby stents have rounded ends with extended endings, which provide an improved fit to the wall and makes stent placement much easier.

When arranged, the stent transfers external radial force to the surface of the vessel lumen for self-expanding (the self-expanding stent has a mechanical spring-form property that enables it to expand to a predetermined diameter).

The LEO+ and LEO+ Baby stents allow repositioning if they are deployed to 90% of their deployment length.

The LEO+ Baby stent has a small profile, which allows access to vessels up to 1.5 mm in diameter. The LEO+ Baby stent is compatible with low profile microcatheter.


Figure 1 - Appearance of the stent

Stent installation system (delivery system)

LEO+ / LEO+ Baby stent installation system LEO includes:

Self-expanding nitinol stent (1) and delivery system. The delivery system consists of delivery conductor (2) and introducer (3). The LEO+/LEO+ Baby stent is pre-installed on a guide located inside the introducer (Fig. 2).


Figure 2 - Stent installation system

Installing LEO+ / LEO+ Baby requires the following components:

❖ Reinforced microcatheter Vasco+ (4) for installation (hereinafter referred to as delivering catheter Vasco+) (figure 3):

- Included in the package of the LEO+ stent delivery kit
- Not included with the LEO+ Baby stent (use only the microcatheter indicated on the stent label)


Figure 3 - Reinforced microcatheter for installation (delivering catheter Vasco+)

- ❖ Y-connector with hemostatic valve (5), which allows to install the introducer in the base of the microcatheter (Figure 4)


Figure 4 - Y-connector with hemostatic valve

- Included in the package of the LEO+ stent delivery kit
- Not included in the delivery kit of a LEO+ Baby stent
- ❖ Rotating device (torquer) (figure 5), which secures the guide firmly and assists to rotating and advancing the microcatheter.


Figure 5 - Rotating device (torquer)

Stents are supplied installed on the delivery guide according to the arrangement below (Figure 6):


Figure 6 – A stent installed on a delivery guide

The design of the Leo+ stent delivery system is presented in Figure 7.


Figure 7 - Design of Leo+ delivery system

The design of the Leo+ Baby stent delivery system is presented in Figure 8.


Figure 8 - Design of Leo+ Baby delivery system

Selecting the size of the stent

Successful implantation of the stent depends on the correct choice of the size of the stent. Incorrect stent size can lead to vascular damage or stent migration.

It is necessary to check the data on the device label regarding the maximum and minimum diameter of the artery, depending on the nominal size of the LEO+ / LEO Baby stent used. The stent should be chosen of such length that it would go beyond the neck of the aneurysm from both sides for a distance exceeding at least 1.5 times the diameter of the target artery (see Figure 9).


Figure 9. Determination of the size of LEO+/LEO+ Baby stent to be used

Technical specification

The device construction materials are shown in Table 1.

Table 1 – The device construction materials

Name of device component	Device materials	
	Leo+	Leo+ Baby
Stent (Implant)		
Braid	Nitinol (nickel-titanium alloy, 45% titanium CAS # 7440-32-6, 55% nickel CAS # 7440-02-0), produced by Fort Wayne Metals, USA	
X-ray contrast markers	Platinum 90% CAS # 7440-06-4 / iridium 10% CAS # 7439-88-5, Production Johnson Matthey, United Kingdom	Platinum 92% CAS # 7440-06-4 / Tungsten 8% CAS # 7439-88-5, Made By Johnson Matthey, United Kingdom
Delivery system		
Delivery guide	Nitinol (nickel-titanium alloy, 45% titanium CAS # 7440-32-6, 55% nickel CAS # 7440-02-0), produced by Fort Wayne Metals, USA	
Distal tip (platinum tip)	Platinum 92% CAS # 7440-06-4 / Tungsten 8% CAS # 7440-33-7, Production Johnson Matthey, United Kingdom	
Proximal part	Stainless steel (AISI SUS 304L) manufactured by Fort Wayne Metals, USA	
Introducer	Polyethylene (ALCUDIA 4810) manufactured by Repsol, Spain	
Fixing wings, radio-opaque marker of delivery guide	Platinum 90% CAS # 7440-06-4 / iridium 10% CAS # 7439-88-5, Production Fort Wayne Metals, USA	Platinum 90% CAS # 7440-06-4 / iridium 10% CAS # 7439-88-5, Production Fort Wayne Metals, USA
Stainless steel braid	Stainless steel (AISI SUS 304L) manufactured by Fort Wayne Metals, USA	
Protective ring	Polyethylene (Rigidex 5502 GA), manufactured by INEOS Olefins & Polymers, UK Clamps: Polyethylene (Rigidex 5502 GA), manufactured by INEOS Olefins & Polymers, UK, pigment Blue (B101)	
Terminal sterilization package for the stent (primary package)	Film 12/50 PET / PE STD (35709 G), Paper PS41 STD, manufactured by SUDPACK MEDICA, China	
Torquer	Polyamide Rilsan BESV0 A MED, by Arkema, France	
Individual package of torquer	Film 12/50 PET / PE STD (35709 G), Paper PS41 STD, manufactured by SUDPACK MEDICA, China	

Y-connector with hemostatic valve	Polyamide Rilsan BESV0 A MED, by Arkema, France Valve: silicon (NY- 600), manufactured by SF COMPOSITES, France Pigment: Red (red01)	Not included in the package
Individual packing of Y-connector with hemostatic valve	Film 12/50 PET / PE STD (35709 G), Paper PS41 STD, manufactured by SUDPACK MEDICA, China	-
Reinforced microcatheter Vasco+	Shaft: Stainless steel (AISI SUS 304L) manufactured by Fort Wayne Metals, USA Tubes polyamide Rilsan BESV0 A MED, by Arkema, France Hub: Polypropylene (Bormed RF825MO), Songhan Plastic Technology Co.,Ltd, China, CAS-No. 9010-79-1 green coloring substance (green C) CAS-No. 1328-53-6, by BASF – The Chemical Company, Germany X-ray markers: Platinum 90% CAS # 7440-06-4 / iridium 10% CAS # 7439-88-5, Production Fort Wayne Metals, USA	Not included in the package
Protective rim of microcatheter	Polyethylene (Rigidex 5502 GA), manufactured by INEOS Olefins & Polymers,	-
Individual package of reinforced microcatheter Vasco+	Film 12/50 PET / PE STD (35709 G), Paper PS41 STD, manufactured by SUDPACK MEDICA, China	-
Secondary package (cardboard box)	Cellulose board with clay coating (GU-1) with a density of 355 g / m ² , manufactured by CARTONNERIES DU VALOIS, France	
Transport package	Corrugated cardboard (AnT10) manufactured by CARTONNERIESDU VALOIS, France	

Type of body contact:

See Paragraph 1 of this document.

Device design:

Device technical specifications by device codes are shown in Table 2.

Table 2 – Technical characteristics of product versions (*the designation of the stent dimensions is shown in Figure 10)

Product	Product versions	Recommended artery diameter (mm)	Nominal diameter of stent (mm), $\pm 5\%$	Length of stent with nominal diameter (mm), $\pm 5\%$	Length of stent on delivery system (mm), $\pm 5\%$	Completely open stent		Delivering catheter	Stent cell size, (+0/-0,01)
						Diameter (mm), $\pm 5\%$	Length (mm), $\pm 5\%$		
Intracranial self-expanding LEO+ Baby stent	LEO.2,0x12	From 1.5 to 2.5	2.0	12	14	2.8	8	Vasco+10(1,9 F)	1.8 F (0.60 mm)
	LEO.2,0x18			18	21		12		
	LEO.2,0x25			25	30		16		
	LEO.2,5x12	From 2.0 to 3.1	2.5	12	17	3.2	8	Vasco+10(1,9 F)	1.9 F (0.63 mm)
	LEO.2,5x18			18	23		12		
	LEO.2,5x25			25	31		16		
	LEO.2,5x30			30	37		19		
	LEO.2,5x35			35	44		23		
	LEO.3,0x12	From 2.5 to 3.6	3.0	12	17	3.8	9	Vasco+21 (2,4F)	1.9 F (0.63 mm)
	LEO.3,0x18			18	24		13		
	LEO.3,0x25			25	32		16		
	LEO.3,0x35			35	47		24		
Intracranial self-expanding LEO+ stent	LEO.3,5x12	From 3.1 to 4.25	3.5	12	18	4.4	9	Vasco+21 (2,4F)	2.1 F (0.70 mm)
	LEO.3,5x18			18	24		12		
	LEO.3,5x25			25	34		17		
	LEO.3,5x30			30	40		19		
	LEO.3,5x35			35	48		24		
	LEO.3,5x50			50	70		36		
	LEO.4,5x15	From 4.25 to 5.30	4.5	15	24	5.4	12	Vasco+25 (3,0F)	2.4 F (0.80 mm)
	LEO.4,5x20			20	30		15		
	LEO.4,5x25			25	37		18		
	LEO.4,5x30			30	42		21		
	LEO.4,5x40			40	56		28		
	LEO.4,5x50			50	70		37		
	LEO.4,5x75			75	109		55		
	LEO.5,5x25	From 5.30 to 6.50	5.5	25	36	6.7	18	Vasco+28 (3,3F)	2.4 F (0.80 mm)
	LEO.5,5x30			30	43		21		

	LEO.5,5x35			35	51		26		
	LEO.5,5x50			50	72		35		
	LEO.5,5x60			60	86		42		
	LEO.5,5x75			75	104		51		



Figure 10 - Designation of stent dimensions

Table 3 – Stent technical specifications

Parameters	Specification
Description	Self-expanding bare metal stent
Shape	Cylindrical, with extended ends Extended ends contribute to better fitting of the implanted stent to the vessel wall. If the stent is installed into a vessel of a rather small diameter, such extension enables to avoid the "tulip" effect, which reduces the risk of thrombosis significantly and facilitates repeated passage through the stent.
Intended clinical use	Brain arteries
Raw material	Nitinol (nickel-titanium alloy, 45% titanium CAS # 7440-32-6, 55% nickel CAS # 7440-02-0), produced by Fort Wayne Metals, USA
Corrosion resistance	The value of corrosion potential difference before testing must be above 200 mV, that after testing, over 30 mV

The probability of contact, galvanic, pit, fretting corrosion, possible oxidation	The probability is low due to application of bioinert materials	
Hazard of issue perforation	Perforation, or rupture of issue may arise, if a stent is used not corresponding to the vessel diameter, or in case of non-compliance with precautions regarding use of the stent.	
Assessment of the stent displacement and migration	With a stent selected properly according to the vessel diameter, the probability of relative displacement of the stent and the tissue is negligible Using the stent is limited to brain arteries, which means that the stent is not affected by heartbeat.	
Design of stent	 16 braided nitinol wires (braiding at 60°(±5%)) 2 platinum nitinol wires all along the length of the stent	
Wire thickness	65 µm (±5%) (nitinol wire) 60 µm (±5%) (platinum wire)	
Stent surface/artery wall ratio (% of artery wall covered with outer surface of stent)	Stent designation	Stent surface/artery wall ratio, minimum, %
	LEO.2.0X12	14.40
	LEO.2.0X18	14.40
	LEO.2.0X25	14.45
	LEO.2.5X12	12.14
	LEO.2.5X18	12.34
	LEO.2.5X25	12.45
	LEO.2.5X30	12.13
	LEO.2.5X35	12.15

	LEO.3.0X12	11.52
	LEO.3.0X18	11.63
	LEO.3.0X25	11.67
	LEO.3.0X35	11.57
	LEO.3.5X12	12.24
	LEO.3.5X18	12.31
	LEO.3.5X25	12.31
	LEO.3.5X30	12.30
	LEO.3.5X35	12.34
	LEO.3.5X50	12.37
	LEO.4.5X15	12.42
	LEO.4.5X20	12.46
	LEO.4.5X25	12.47
	LEO.4.5X30	12.52
	LEO.4.5X40	12.57
	LEO.4.5X50	12.55
	LEO.4.5X75	12.58
	LEO.5.5X25	10.95
	LEO.5.5X30	11.00
	LEO.5.5X35	11.03
	LEO.5.5X50	11.06
	LEO.5.5X60	11.09
	LEO.5.5X75	11.10
Minimum bending/twisting radius	10 mm	
Resistance to failure of stent under plane-parallel load	stents are not compressed in clinical use	
Radial force produced by stent onto vessel wall in the open position at expansion and contraction.	Within 0.010 N/mm	

Extraction and replacement of implant	The stent is installed to be used indefinitely and is not supposed to be extracted from the patient's body and replaced.
Resistance to effects of biological fluids	The products are resistant to effects of biological fluids and secrets released by body tissues they contact with at operation: resistant to effects of aggressive biological fluids (blood) (the stent material is inert to effects of internal milieu)

Table 4 - Technical characteristics of stent delivery system

Parameters	Specification																						
Design	Delivery of stent along a guide.																						
Length of delivery system	2000 mm (+10/-100 mm)																						
Outer diameter of delivery system	0.38 mm (± 0.055 mm) for Leo+Baby 0.55 mm (± 0.055 mm) for Leo+																						
Length of the distal tip	9 mm (+1/-0 mm)																						
Diameter of the distal tip	0.22 mm (± 0.02 mm)																						
Length of introducer	600 mm (+5/-20 mm)																						
Length of stainless steel braid	4 mm (± 0.05 mm) for Leo+ For Leo+Baby: <table> <tr> <td>LEO.2.0X12</td><td>89 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.0X18</td><td>82 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.0X25</td><td>74 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.5X12</td><td>89 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.5X18</td><td>82 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.5X25</td><td>74 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.5X30</td><td>68 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.2.5X35</td><td>61 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.3.0X12</td><td>89 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.3.0X18</td><td>82 mm (± 0.1 mm)</td></tr> <tr> <td>LEO.3.0X25</td><td>74 mm (± 0.1 mm)</td></tr> </table>	LEO.2.0X12	89 mm (± 0.1 mm)	LEO.2.0X18	82 mm (± 0.1 mm)	LEO.2.0X25	74 mm (± 0.1 mm)	LEO.2.5X12	89 mm (± 0.1 mm)	LEO.2.5X18	82 mm (± 0.1 mm)	LEO.2.5X25	74 mm (± 0.1 mm)	LEO.2.5X30	68 mm (± 0.1 mm)	LEO.2.5X35	61 mm (± 0.1 mm)	LEO.3.0X12	89 mm (± 0.1 mm)	LEO.3.0X18	82 mm (± 0.1 mm)	LEO.3.0X25	74 mm (± 0.1 mm)
LEO.2.0X12	89 mm (± 0.1 mm)																						
LEO.2.0X18	82 mm (± 0.1 mm)																						
LEO.2.0X25	74 mm (± 0.1 mm)																						
LEO.2.5X12	89 mm (± 0.1 mm)																						
LEO.2.5X18	82 mm (± 0.1 mm)																						
LEO.2.5X25	74 mm (± 0.1 mm)																						
LEO.2.5X30	68 mm (± 0.1 mm)																						
LEO.2.5X35	61 mm (± 0.1 mm)																						
LEO.3.0X12	89 mm (± 0.1 mm)																						
LEO.3.0X18	82 mm (± 0.1 mm)																						
LEO.3.0X25	74 mm (± 0.1 mm)																						

	LEO.3.0X35 61 mm (± 0.1 mm)
Diameter of stainless steel braid	0.32 mm (± 0.1 mm) for Leo+ 0.21 mm (± 0.1 mm) for Leo+Baby
X-ray marker of delivery guide	$\varnothing 0.36$ mm x 0.42 mm x 0.5 mm (± 0.1 mm) for Leo+ $\varnothing 0.24$ mm x 0.45 mm x 0.5 mm (± 0.1 mm) for Leo+Baby
Fixing wings	$\varnothing 0.35$ mm x 0.63 mm x 0.6 mm (± 0.1 mm) for Leo+ $\varnothing 0.18$ mm x 0.28 mm x 0.65 mm (± 0.1 mm) for Leo+Baby
Effort required for the stent to open	No more than 5 N
Possibility of advancement	Advancement efforts under 3N
Pushability	enables the operator to arrange the stent freely, with no undesirable bending, or twisting, due to a smooth tip
Minimum radius of delivery system bending/twisting	10 mm
Rotatability	ensures free rotation due to a smooth tip
Strength of connection at rotation	At least 10 N•mm
Visualization	The stent is delivered to its destination under fluoroscopic control
Tensile strength of tubes	At least 10 N
Strength of connection of distal tip	At least 5 N
Strength of connection of guide marker	At least 10 N
Strength of connection of marker of fixing wings	At least 10 N
Strength of connection of distal part, stainless steel braiding, proximal part	At least 10 N


Technical characteristics of reinforced microcatheter Vasco+ (delivering catheter Vasco+) are presented in table 3.

Table 3 - Technical characteristics of microcatheter Vasco+ (delivering catheter Vasco+)

	Vasco+10(1,9F)	Vasco+21 (2,4F)	Vasco+25 (3,0F)	Vasco+28 (3,3F)
Overall length (mm)	1550±25	1550±25	1550±25	1550±25
Working length (mm)	1490±25	1490±25	1490±25	1490±25
Outer distal diameter, Fr (mm)	1.9 (0.63±0.05 mm)	2.4 (0.81±0.1mm)	3.0 (1.04±0.1 mm)	3.3 (1.08±0.1 mm)
Inner distal diameter " (mm)	0.017 (0.41±0.05 mm)	0.024 (0.60±0.05 mm)	0.029 (0.74±0.05 mm)	0.0323 (0.82±0.05 mm)
Outer proximal diameter Fr (mm)	2.2 (0.72±0.1mm)	2.7 (0.94±0.1mm)	3.3 (1.10±0.1mm)	3.4 (1.13±0.1mm)
Inner proximal diameter " (mm)	0.017 (0.41±0.05 mm)	0.024 (0.60±0.05 mm)	0.029 (0.74±0.05 mm)	0.0323 (0.82±0.05 mm)
Maximum diameter of compatible guide " (mm)	0.014 (0.36)	0.021 (0.53)	0.025 (0.64)	0.028 (0.71)
Number of X-ray contrast markers (pcs)	2	1	1	1
Width of X-ray contrast markers (mm)	0.900±0.005	0.900±0.005	0.900±0.005	0.900±0.005
Distance from distal tip to X-ray marker, mm	1.0 The second marker 30 mm away from the first one	1.0	1.0	1.0
Connecting joint	External diameter 4 mm, internal diameter 2 mm, 6% taper fitting, no screw thread	External diameter 4 mm, internal diameter 2 mm, 6% taper fitting, no screw thread	External diameter 4 mm, internal diameter 2 mm, 6% taper fitting, no screw thread	External diameter 4 mm, internal diameter 2 mm, 6% taper fitting, no screw thread
Strength of marker connection	At least 5 N	At least 5 N	At least 5 N	At least 5 N
Strength of connection of hub and rigid tube	At least 10 N	At least 10 N	At least 10 N	At least 10 N
Strength of connection of catheter tubes	At least 10 N	At least 10 N	At least 10 N	At least 10 N


* Unspecified deviations of dimensions ±5%

Vasco+10 design:


*All dimensions are indicated in mm. Size tolerances, which are not specified in Table 3, are $\pm 5\%$

Vasco+21 design:


*All dimensions are indicated in mm. Size tolerances, which are not specified in Table 3, are $\pm 5\%$

Vasco+25 (3.0F) design:

*All dimensions are indicated in mm. Size tolerances, which are not specified in Table 3, are $\pm 5\%$

Vasco+28 (3.3F) design:

*All dimensions are indicated in mm. Size tolerances, which are not specified in Table 3, are $\pm 5\%$

Y-connector with hemostatic valve:

It must pressurize well and prevent the blood loss.

Distal Luer tip is connected with the catheter or other devices.

Length in open position	68±1 mm
Length in closed position	63.3 ±1 mm
Internal diameter of valve	2.67 mm±5%
Connection joint	Luer lock with tapering of 6%
Valve air-tightness	Up to 300 kPa
Maximum weight	17 g

Torque device (torquer):

Parameter	Value
Length in fully open position	51.5 ±0.5 mm
Length in fully closed position	48 ±0.5 mm
External diameter of handle	0.7 ±0.1 mm
External diameter of cap	10 ±1 mm
Internal (working) diameter	1 ±0.01 mm
This product shall ensure gripping of the wire with diameter	From 0.36 to 0.89 mm
Maximum weight	4 g

Package:

	Overall dimensions	Specifications
Package for terminal sterilization (primary package): stent	290 x 255 mm ($\pm 5\%$)	Width of adhesive layer: 10 ± 1 mm, Tensile strength: ≥ 1.2 N Detachment/opening force: 2.5 ± 0.5 N, Weld joints shall be even – without burnt spots or folds; packing material surfaces shall be free from cracks, pressed folds, wrinkles, punctures, tears, material laminations, dirt, as well as sticking of internal surfaces of packing materials.
Package for terminal sterilization (primary package): Y-connector with hemostatic valve	184 x 100 mm ($\pm 5\%$)	Width of adhesive layer: 10 ± 1 mm, Tensile strength: ≥ 1.2 N Detachment/opening force: 2.5 ± 0.5 N, Weld joints shall be even – without burnt spots or folds; packing material surfaces shall be free from cracks, pressed folds, wrinkles, punctures, tears, material laminations, dirt, as well as sticking of internal surfaces of packing materials.
Package for terminal sterilization (primary package): torque device (torquer)	184 x 100 mm ($\pm 5\%$)	Width of adhesive layer: 10 ± 1 mm, Tensile strength: ≥ 1.2 N Detachment/opening force: 2.5 ± 0.5 N, Weld joints shall be even – without burnt spots or folds; packing material surfaces shall be free from cracks, pressed folds, wrinkles, punctures, tears, material laminations, dirt, as well as sticking of internal surfaces of packing materials.
Package for terminal sterilization (primary package): Vasco+ reinforced microcatheter	290 x 255 mm ($\pm 5\%$)	Width of adhesive layer: 10 ± 1 mm, Tensile strength: ≥ 1.2 N Detachment/opening force: 2.5 ± 0.5 N, Weld joints shall be even – without burnt spots or folds; packing material surfaces shall be free from cracks, pressed folds, wrinkles, punctures, tears, material laminations, dirt, as well as sticking of internal surfaces of packing materials.
Secondary package (cardboard box)	300x 260 mm ($\pm 5\%$)	-
Transportation container	526 x 380 x 306 mm ($\pm 5\%$)	-

Data on sterility

This product is supplied sterile. Sterilization parameters: gas sterilization method (ethylene oxide) in compliance with requirements of ISO 11135 (residual contamination level is 10^{-6}).

Repeated sterilization of this product is forbidden.

Repeated use is forbidden.

Transportation, storage and use

This product may be transported by any means of transportation, in covered vehicles and in compliance with cargo transportation regulations valid for the relevant means of transportation.

Transportation conditions: air temperature from $+5^{\circ}\text{C}$ to $+35^{\circ}\text{C}$, relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 – 1060 hPa.

Storage conditions: air temperature from $+5^{\circ}\text{C}$ to $+35^{\circ}\text{C}$, relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 – 1060 hPa.

The product shall be stored in the original package of the manufacturer in a clean, dry place, away from heat sources and direct sun light, at controlled temperature.

Use conditions: air temperature from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$, relative humidity from 60% to 100% (non-condensing) and atmospheric pressure of 500 – 1060 hPa.

Recommendations to medical specialist

Materials (Information for use for a medical specialist)

The product and its package do not contain any natural latex or polyvinylchloride.

Indications:

LEO+/LEO+ Baby self-expanding stent is used for treatment of intracranial aneurysms in combination with embolic coils.

Contraindications:

Use of LEO + stent is counter-indicative in the following cases:

- In patients with developing bacterial infection.
- In patients with contraindications for anticoagulant, anti-platelet therapy, or thrombolytics.
- In patients with known hypersensitivity to nickel-titanium.
- In patients with anatomic features, which prevent stent from advancing or expanding.
- In patients, who had not been administered anti-platelet drugs prior to the procedure.
- In patients resistant to the effect of anti-platelet and/or anticoagulant drugs.
- These products are counter-indicative for newborns, prematurely born children and infants.

Potential complications:


Complications and adverse manifestations may include without limitation:

- Death
- Renal insufficiency
- Sepsis / Infection
- Embolism
- Apoplexy
- Allergic reactions
- Perforation, rupture, dissection or any other arterial damage
- Diffuse intravascular coagulation
- Hemorrhage
- Arrhythmia
- Intima rupture or dissection
- Neurological disorders
- Aneurysm recanalization
- Migration / embolism
- Thrombosis
- Circulatory insufficiency
- Heart failure and ischemia
- AV fistula

- Necrosis of tissues
- Malposition of stent
- Vascular occlusion, or new stenosis, or new constriction at the place of implantation
- Hematoma
- Pain and tenderness
- Pyrogenic reaction

14. Safety precautions

- Do not use the product, if its package had been opened or damaged.
- Stents are apyrogenic.
- **This product is designed only** for single use. Do not reuse it. Any repeated use of this product will lead to a high risk of patient infection and product performance degradation.
- Do not re-sterilize it.
- **Keep this product in a dry place at a room** temperature and away from sunlight. Do not use the product, if expired.
- LEO+/LEO+ Baby stents shall be used only by medical specialists experienced in the field of interventional neuroradiology and treatment of intracranial aneurysms.
- Carefully determine the vessel diameter (± 0.2 mm) in the place chosen for stent implantation as well as the aneurysm neck length.
- **IMPORTANT: The stent ends shall be always positioned on straight-line portions of the artery.**

Stent positioning

- **Warning: From the perspective of currently existing literature, use of stents in bifurcation aneurysm is not recommended.**

Safety and effectiveness of stent-in-stent technique and implantation on the level of a collateral branch of a vessel have not been clinically assessed.

- Use of thrombocytic anti-platelet, anti-coagulating and vasodilating drugs is very important for successful implantation and follow-up monitoring.
- Do not use this product in patients resistant to the effect of anti-platelet drugs.
- Big (from 15 to 25 mm) and huge aneurysms (> 25 mm) are characterized by a high risk of hemorrhage. Implantation of LEO+ stent together with embolic coils into an aneurismal sack is very important, because it can help reduce this risk.
- Microcatheter and stent should be used under fluoroscopic control and with appropriate anti-coagulating agents.
- Do not attempt to retrack the stent into the introducer tube.
- Do not move the product, if it meets resistance, without finding out its causes first. If there is some resistance, the stent movement may damage the product and/or vessel.
- Follow the instructions for use of the products during the procedure.
- The catheter should be washed during the procedure.
- **Stent may be used only with the catheter, which is supplied in boxes** (LEO+ stents) and indicated on the label (LEO+/ LEO + Baby stents).
- Take care, while performing different actions with the product, including stent unpacking or use, to avoid any damages of the stent and its delivery system.

Information about MRI

LEO+ stent was found safe for use under certain parameters of magnetic resonance imaging (MRI) based on the terminology indicated by the American Society for Testing and Materials (ASTM).

The pre-clinical experience showed safe use of LEO+ stents under certain MRI conditions. Patients can be safely scanned provided that the following conditions are observed:

Static magnetic field

- static magnetic field of 3 T or less;
- maximum spatial gradient of magnetic field of 720 Gs/cm or less.

Heating associated with MRI

During the pre-clinical trial, LEO+ stent was causing a temperature rise during the MRI scanning within 15 minutes with the application of 1.5-Tesla (1.5 Tesla / 64 MHz, Magnetom, Siemens Medical Solutions, Malvern, PA) MRI systems. Numaris / 4, Version Syngo MMR 2002B DHHS) and 3-Tesla (3-Tesla / 128 MHz, Excite) software was used. Also, G3.0-052B, General Electric Healthcare, Milwaukee, WI software was used.

Maximum temperature change

MRI status

+ 2.0°C

1.5-T / 64-MHz

+ 2.1°C

3-T / 128-MHz

These temperature changes are not hazardous for human under the above-stated conditions.

Information about artifacts

The MR image quality can worsen, if the region of interest is within the same area or relatively close to the area, where the implant is positioned. That is why some optimization of MRI parameters may be needed to compensate for the presence of this implant.

Pulse sequence	T1-SE	T1-SE	GRE	GRE
Signal void size	486-mm ²	52-mm ²	831-mm ²	79-mm ²
Image plane orientation	parallel	perpendicular	parallel	perpendicular

A) Choice of stent size

Success of stent implantation depends on correct choice of the stent size. A stent of the wrong size may lead to damage of vessels or stent migration.

It is required to be guided by the data indicated on the product label regarding the permissible maximum and minimum artery diameter with respect to the nominal size of the applied LEO+/LEO+ Baby stent.

The stent length must be chosen such as to ensure the stent protrusion beyond the aneurysm neck from both sides by the length exceeding the target artery diameter, at least, 1.5 times (see Fig.1 below).


Fig. 1.

B) Catheter positioning

- Use only the microcatheter indicated on the label.
- Check up the state of the package to make sure that it has not been opened or damaged.
- Open the sterile package and retrieve the microcatheter using the sterility preservation method.
- Examine the microcatheter to make sure that it has not been damaged.


- Wash the microcatheter in sterile saline solution to ensure its patency.
- Insert the wire guide into the microcatheter (see the instructions to the wire guide) through a Y-connector with the valve, ensuring thus the infusion of the guiding catheter during the whole process.
- Advance gradually forward the assembled unit formed by the wire guide and the microcatheter, moving in turns first the wire guide and then the microcatheter.
- Direct the wire guide into the chosen vessel, applying carefully a torque force to it under the visualization control.

C) LEO+ stent preparation before implantation

- Check up the state of the package to make sure that it had not been opened or damaged. Open the package and, following the aseptic rules, retrieve the introducer tube and the Leo stent, fit on the delivery wire.
- Check up all the components to make sure that they are not damaged.
- Make sure that the end of the delivery wire is entirely within the introducer.
- **Attention: Do not let the stent expand (even partially) beyond the introducer.**
- **Attention: Do not give shape to the delivery wire tip.**
- Insert the introducer halfway into the Y-connector with the hemostatic valve assembled with the microcatheter.
- Close the Y-connector valve; do not connect the introducer to the microcatheter base, so that the liquid could flow back in the tube.
- A drop of liquid formed on the tube end opposite to the one located inside the catheter base signals that the introducer tube is filled with the liquid.
- Loosen the Y-connector valve.
- Push the introducer tube through the Y-connector, until it is connected to the microcatheter base.
- Close the valve so that the introducer tube remained in its place. **Attention: check up if there is a reliable contact and conformity between the introducer tube and the internal part of the microcatheter base.**
- Advance the delivery wire forward so that the stent was displaced from the introducer into the microcatheter. Advance the stent inside the microcatheter for, at least, 70 cm.
- Do not attempt to retrack the stent back into the introducer tube.
- Loosen the Y-connector valve and remove the introducer tube, retaining the delivery wire in its place.
- Tighten up the Y-connector valve.
- Then advance the stent in the microcatheter with the help of the delivery wire.
- **Attention: if any resistance is felt at any stage of the work with the stent, retrieve the system and replace it with a new one.**

D) Stent expansion

- Put the distal end of the stent (1) at the target place; the distal end shall protrude beyond the aneurysm neck for the length exceeding the target vessel diameter, at least, 1.5 times (see Fig.1).
- Carefully expand the stent, while advancing the delivery wire (2) **and** retracting the microcatheter (4).
- Note:** At this stage, the distal end of the delivery wire (2) will advance forward (as the stent will shorten, while getting expanded).
- **Attention: Make sure that this motion is done freely in the artery to avoid perforation.**
- At any moment of expanding (until the distal marker (7) of the microcatheter (4) coincides with the marker (6) of the delivery wire), the stent can be returned into the microcatheter by pulling the delivery wire (2).

Appearance of radio-opaque markers

Fig. 2

- If the stent position needs to be changed, pull carefully the microcatheter over the expanded stent, then shift the system and expand the stent at the new place.
- **Note:** While the microcatheter is being pulled over the stent to return it, the stent position may be stabilized, if needed, by pulling at the delivery wire.
- **Caution:** If any resistance is felt during the return of the stent into the microcatheter, stop the attempts to return the stent. Pull back slightly the microcatheter to expose the stent (trying not to exceed the limit, beyond which the return of the stent would be impossible) and then repeat again the attempt to return the stent.
- **Caution:** The full return of the stent is possible only once.
- The stent fully comes out of the microcatheter, when the marker (6) of the delivery wire (2) coincides with the marker (7) of the microcatheter (see Fig. 2 above).
- Keep pushing the delivery wire (2) until it comes onto the radio-opaque marker of the catheter (distal ring) (7) approximately for 2 mm.
- If at this moment the stent is not separated from the delivery wire (2), move the microcatheter (4) away by 2 cm in the proximal direction, trying not to move the delivery wire. Try to disconnect the stent by turning gently the delivery wire through 180° from right to left (Attention: do not turn it through more than 360°, otherwise the delivery wire may break!) and waiting until the distal part makes a turn during the revolving. Then again advance the microcatheter towards the stent, until it reaches the stent. Repeat these actions, if needed.
- **Attention:** To avoid the displacement of the stent after its expansion, do not move or retrieve the delivery wire, before you make sure first that the stent is completely disconnected.
- **Note:** Keep the expanded stent accessible by leaving the delivery wire inside the stent.
- Before retrieving the delivery wire, place the distal end of the microcatheter (4) inside the stent, to keep the distal end of the stent accessible.
- After making sure that the stent is correctly positioned and duly expanded, carefully remove the microcatheter and the delivery wire.

Technical repair and maintenance

The product shall not be subject to technical repair or maintenance. It is designed only for single use.

Labeling and packing

Each LEO+/LEO+ Baby stent on delivery system is enclosed in a protective plastic ring for additional protection against external effects during the transportation and storage, which has to be removed before the use and then sealed into the package for terminal sterilization (individual package). Afterwards, together with the instruction for use it shall be packed into a thick cardboard box. The box of LEO+ stents shall also contain the following items: a Y-connector with hemostatic valve, a torque device (torquer), Vasco+ reinforced microcatheter, each sealed in a separate package for terminal sterilization (individual package). The box of LEO+ Baby stents shall also contain a torque device (torque), which is sealed into a separate package for terminal sterilization (individual package).

For keeping stock and traceability of products, each package of the product is provided with adhesive labels – stickers – which shall be attached to the medical record and other records of the patient with the implanted stent. The sticker shall contain the following information:

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
c) production lot No. and/or sterilization lot No.: **LOT 00277548**

d) lot or model No. (based on reference No.): **REF LEO.3,5X25**

Sterile product. Its use is strictly forbidden, if the package is damaged.
The package conforms to ISO 11607.

The labeling of the secondary package (box) of LEO+ contains:

- Product name, - Lot No. (LOT),	
Reference No. (REF),	
Information about manufacturer	
- Text: "Made in France", - Expiry date,	
Manufacturing date,	
CE labeling,	
Symbol: "Single use only",	
Symbol: "Re-sterilization is forbidden",	
- Information about product supply package, - Symbol: "See instruction for use",	
Symbol: "Caution! See instruction for use",	

Symbol: "Keep dry",	
• Symbol: "Keep away from sunlight",	
• Symbol: "Do not use, if package is opened or damaged"	
• Symbol: "Sterilized with ethylene oxide",	
• Symbol: "MRI"	

• Symbol: "apyrogenic",	
• Non-toxic,	
• Compatible diameter of vessels,	For vessels with diameter from to
• Information on compatible microcatheter, including its basic dimensions,	Compatible only with
• Nominal diameter of stent, length of stent with nominal diameter,	
• Length of stent on delivery system,	
• Length and diameter of fully extended stent	

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Labeling on the individual package of LEO+ contains:

• Product name	
• Lot No. (LOT)	
• Reference No. (REF),	
• Information about manufacturer	
• Text: "Made in France",	
• Expiry date,	
• Manufacturing date,	
• CE labeling,	
• Symbol: "For single use only",	
• Symbol: "Re-sterilization is forbidden",	
• Information on stent components (stent and delivery system)	
• Symbol: "See instruction for use",	
• Symbol: "Caution! See instruction for use",	

• Symbol: “Keep dry”,	
• Symbol: “Keep away from sunlight”,	
• Symbol: “Do not use, if package is opened or damaged”	
• Symbol: “Sterilized with ethylene oxide”,	
• “MRI” symbol	

• Symbol: “apyrogenic”,	
• Non-toxic,	
• Compatible diameter of vessels,	For vessels with diameter from to
• Information on compatible microcatheter, including its basic dimensions,	Compatible only with
• Stent mesh size	
• Information on compatible microcatheter, including its basic dimensions,	
• Nominal diameter of stent, length of stent with nominal diameter	
• Length of stent on delivery system	
• Length and diameter of fully expanded stent	

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2016-10

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EAN13:

The labeling of the individual package of Vasco+ reinforced microcatheter contains:

• Product name,	
• Lot No. (LOT),	
• Reference No. (REF),	
• Information about manufacturer	
• Text: "Made in France",	
• Expiry date,	
• Manufacturing date,	
• CE labeling,	
• Symbol: "For single use only",	
• Symbol: "Re-sterilization is forbidden",	
• Information on stent components (stent and delivery system),	
• Symbol: "See instruction for use",	
• Symbol: "Caution! See instruction for use",	

• Symbol: "Keep dry",	
• Symbol: "Keep away from sunlight",	
• Symbol: "Do not use, if package is opened or damaged"	
• Symbol: "Sterilized with ethylene oxide",	
• "MRI" symbol	

• Symbol: "apyrogenic",	
• Non-toxic,	
• Contains 1 catheter	
• Overall dimensions	
• Text: "Not compatible with DMSO"	
• Max diameter of compatible wire	

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The labeling of the secondary package (box) of LEO+ Baby contains:

• Product name,	
• Lot No. (LOT),	
• Reference No. (REF),	
• Information about manufacturer	
• Text: "Made in France",	
• Expiry date,	
• Manufacturing date,	
• CE labeling,	
• Symbol: "For single use only",	
• Symbol: "Re-sterilization is forbidden",	
• Information on stent supply package,	
• Symbol: "See instruction for use",	
• Symbol: "Caution! See instruction for use",	

• Symbol: "Keep dry",	
• Symbol: "Keep away from sunlight",	
• Symbol: "Do not use, if package is opened or damaged"	
• Symbol: "Sterilized with ethylene oxide",	
• "MRI" symbol	

• Symbol: "apyrogenic",	
• Non-toxic,	
• Compatible diameter of vessels,	For vessels with diameter from to
• Stent mesh size,	
• Information on compatible microcatheter, its internal diameter (not included into the supply package),	Compatible with
• Nominal diameter of stent, length of stent with nominal diameter,	
• Length of stent on delivery system,	
• Length and diameter of fully expanded stent	

The labeling of the individual package for LEO+ Baby contains:


• Product name,	
• Lot No. (LOT),	
• Reference No. (REF),	
• Information about manufacturer	
• Text: "Made in France",	
• Expiry date,	
• Manufacturing date,	
• CE labeling,	
• Symbol: "For single use only",	
• Symbol: "Re-sterilization is forbidden",	
• Contains 1 product,	
• Information on stent components (stent and delivery system),	
• Symbol: "See instruction for use",	
• Symbol: "Caution! See instruction for use",	

• Symbol: "Keep dry",	
• Symbol: "Keep away from sunlight",	
• Symbol: "Do not use, if package is opened or damaged"	
• Symbol: "Sterilized with ethylene oxide",	
• "MRI" symbol	

• Symbol: "apyrogenic",	
• Non-toxic,	
• Compatible diameter of vessels,	For vessels with diameter from to
• Stent mesh size,	
• Information about a compatible microcatheter, its internal diameter (not included into the supply package),	Compatible with
• Nominal diameter of stent, length of stent with nominal diameter,	
• Length of stent on delivery system,	
• Length and diameter of fully expanded stent	

For transportation purposes the products shall be laid into a package made of corrugated cardboard. Each transportation container includes two products named "LEO+ self-expanding intracranial stent" or "LEO+ Baby self-expanding intracranial stent" enclosed in the secondary package.

Labeling of the transportation container

The labeling of the transportation container includes:

• Product name,	
• Lot No. (LOT),	
• Reference No. (REF),	
• Information about manufacturer	
• Expiry date,	

• Manufacturing date,	
• Symbol: "Keep dry",	
• Symbol: "Keep away from sunlight",	
• Symbol: "Do not use, if package is opened or damaged"	
• Symbol: "Storage temperature"	
• Number of products in a transportation container	
• Information about the authorized representative of the manufacturer	

Note: After the procedure of state marketing authorization on the RF territory, the primary and secondary packages as well as the transportation containers shall be provided with the stickers containing the following information: No. and date of Marketing Authorization.

Conformity to standards of the Russian Federation

GOST R 50444-92 "Medical instruments, apparatus and equipment. General specifications";
GOST ISO 10555-1-2011 "Sterile, single-use intravascular catheters. Part 1. General requirements";
GOST R ISO 25539-2-2012 "Cardiovascular implants. Endovascular implants. Part 2. Vascular stents";
GOST R ISO 14630-2011 "Non-active surgical implants. General requirements";
GOST R ISO 15223-1-2014 "Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied. Part 1. General requirements";
GOST ISO 11607-2011 "Packaging for terminally sterilized medical devices. General requirements";
GOST R ISO 11070-2010 "Sterile single-use intravascular catheter introducers. Technical requirements and test methods";
GOST R ISO 10993-2-2009 "Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements";
GOST ISO 10993-1-2011 "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing";
GOST ISO 10993-4-2011 "Medical devices. Biological evaluation of medical devices. Part 4. Selection of tests for interactions with blood";
GOST ISO 10993-5-2011 "Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity";
GOST ISO 10993-6-2011 "Medical devices. Biological evaluation of medical devices. Part 6. Tests for local effects after implantation";
GOST ISO 10993-7-2016 "Medical devices. Biological evaluation of medical devices. Part 7. Ethylene oxide sterilization results";
GOST ISO 10993-9-2015 "Medical devices. Biological evaluation of medical devices. Part 9. Framework identification for identification and quantification of potential degradation products";
GOST ISO 10993-10-2011 "Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity";
GOST ISO 10993-11-2011 "Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity";
GOST ISO 10993-12-2015 "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials";
GOST ISO 10993-13-2016 "Medical devices. Biological evaluation of medical devices. Part 13. Identification and quantification of degradation products from polymeric medical devices";
GOST ISO 10993-15-2011 "Medical devices. Biological evaluation of medical devices. Part 15. Identification and quantification of degradation products from metals and alloys";
GOST R 52770-2016 "Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests";
GOST 31214-2016 "Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity";
GOST 31209-2003 "Containers for blood and its components. Requirements for chemical and biological safety and methods of testing", par. 5.3
National Pharmacopoeia, Volume XIII, General Monograph 1.2.4.0005.15 "Pyrogenicity".
National Pharmacopoeia, Volume XIII, General Monograph 1.2.4.0003.15 "Sterility".

Data on disposal

Products shall be disposed of in compliance with SanPiN 2.1.7.2790-10 for Class B waste (epidemiologically hazardous waste). After hardware decontamination techniques with application of physical methods and changing of the waste appearance, thus ruling out the possibility of their repeated usage, Class B waste can be accumulated, temporarily stored, transported, eliminated and dumped together with Class A waste (epidemiologically safe waste brought more in line with solid municipal

waste in their composition). Packages with Class B decontaminated medical waste shall bear the labeling certifying that the waste was decontaminated.

Guarantees

Guaranteed shelf life is 5 years, sterilization validity time is 5 years.

BALT EXTRUSION, France, will not be held liable in the following cases:

- (I) If a product was misused, or used by untrained users;
- (II) If a product was used for indications other than it is designed for.
- (III) If a product was used together with the products not approved by BALT EXTRUSION;
- (IV) If a product was repeatedly used or sterilized.

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+33 (0)1 39 89 46 41

silvana.brown@balt-usa.com

Authorized representative of the Manufacturer on the RF territory:

"Perinatal Medical Center" LLC (PMC LLC)

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