
CERTIFICATE OF AUTHENTICITY AND SIGNATURE WITNESS

State of Massachusetts
County of Essex

On this 6st day of September 2017, before me, the undersigned Notary Public, personally appeared name of document designer, proved to me through satisfactory evidence of identification, which was their employee identification and/or driver's license, to be the person whose name is signed on the preceding or attached document in my presence.

(Official signature and seal of Notary)

Signature of Document Custodian (Nisarg Shah)

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, USA
Tel: +1 763 514 4000

Subscribed and sworn to before me this the
On this 6st day of September
2017

Heather P. Morose

My commission expires on June 22, 2023

Medtronic

DxTerity™

Diagnostic Catheter
Диагностичен катетър
Diagnosticky katetr
Diagnostisk kateter
Diagnostikkatheter
Διαγνωστικός καθετήρας
Cateter de diagnostic
Diagnostiline kateeter
Diagnostinen katetri
Catheter de diagnostic
Dijagnosticki kateter
Diagnosztikai katéter
Catelere diagnostico

Diagnostinis kateteris
Diagnostikas katetrs
Diagnostische katheter
Diagnostisk kateter
Cewnik diagnostyczny
Cateter de diagnostic
Cateter de diagnostic
Cateter de diagnostic
Диагностический катетер
Diagnosticky kateter
Dijagnostični kateter
Dijagnosticki kateter
Diagnostisk kateter
Dijagnostik Kateler
Діагностичний катетер

Instructions for Use • Инструкции за употреба •
Navod k použití • Bruksanvisning • Gebrauchsanweisung •
Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend •
Käyttöohjeet • Mode d'emploi • Upute za upotrebu •
Használati útmutató • Istruzioni per l'uso • 사용 지침 •
Naudojimo instrukcijos • Lietosanas pamācība •
Gebruiksaanwijzing • Bruksanvisning •
Instrukcja użytkowania • Instruções de utilização •
Instruções de utilização • Instrucţiuni de utilizare •
Инструкция по эксплуатации • Pokyny na používanie •
Navodila za uporabo • Uputstva za upotrebu •
Bruksanvisning • Kullanım Talimatları •
Інструкція з експлуатації • —

$$\frac{1}{\sqrt{\pi}} \int_{-\infty}^{\infty} f(x) e^{-x^2} dx = \frac{1}{\sqrt{\pi}}$$

USA.

```
master 20160505 alldocMod 20160504 alldocMod 20151009 alldocStyle_graph
20160504 alldocStyle 20160504 master_pageset 20150504 languageLookup 20160330
  20160113 characterJoin 20110211 3of9 20110211 3of9 svg 20110211
```

	Useable length: Исполняемая длина Vylna delka: Anvendelig længde: Nutzbare Länge: Deyilgiy uzunluq: Lõngitud ulil: Kasutatav pikkus: Käyttöpituus: Longueur utile: Iskonsuva dölina: Hasznos hossz: Lánghezza ulile: Naudngas lags: Lieldelgais ga: ms: Brukbare lengde: Nyttbar längde: Dugosc uzyteczna: Compimento ulil: Compimento ulil: Lungime utilizabila: Рабочая длина: Pouzitelna dĺžka: Uooraona dolona: Upotřebijiva dužina: Arbetslängd: Kulatilabir uzunluk: Корисна довжина
	Maximum guidewire diameter: Максимален диаметър на воден: Maximalni průměr vodičino drátu: Maksima guidewirediameter: Maximaler Durchmesser des Führungsdrahts: Μέγιστη διαμετρος: Duppuroq: Diametro máximo de la guía: Maksimaalne juhteliad: labimoot: Ohuainvajerni enimmisajapimite: Diametre maximum du fil-guide: Maksimaal: promjer žice vodičice: A vezetodróti legnagyobb átmérője: Diametro massimo del filo guida: Didžiausias kempiniosios vielos skersmuo: Maksimalni vaditlaistigas diametrs: Luax male diameter: voerdraad: Stierste diameter på ledetråder: Maksymalna średnica prowadnika: Diametre maximum do filo guida: Diametro maximum de filo guida: Diametre maximum de phidare: Максимальный диаметр проводника: Maksimaalnyi premer vodičeno drótu: Najveći premer vodične žice: Maksimaal: prenik vodič-zice: Ledarens maximala diameter: Maksimum kádvuz tel capi: Maksimaalnyi diameter: проводника: 最大導線直径
	Number of sideholes: Брой странични отвори: Pocer posirannich otvoru: Antal sidehulle: Anzahl der seitlichen Öffnungen: Apisúcs, trácsúcsúv otvora: Numero de orificios laterales: Kuligavade arv: Sivure-kyen maara: Nombre d'orifices lateraux: Broj bočnih otvora: Oldalnyílások száma: Numero di fori laterali: Sorknu anou skaicus: Sanu caurumi, skaits: Aantla: z yopeningen: Antal sidenuli: Lintza otvoro: bocznym: Numero de orificios laterais: Numero de orificios laterais: Numar de orificii laterale: Количество боковых отверстий: бо́чных отвори: Strvilo stranskih odprtij: Broj bočnih otvora: Anzahl seitlich: Yarı delikler: sayısı: Kúpsúcsúv, pótsúcsúv utvoro
	PSI / kPa: Maximum PSI/kPa: Максимум PSI/kPa: kPa: Maximum PSI/kPa: Maksimal PSI/kPa: Maximum psi/kPa: Mátyon: psi/kPa: Václ: maximum de presi on: PSI/kPa: Maksimaalne PSI/kPa: PSI/kPa: ampmassuun PSI/kPa: Maximum PSI/kPa: Maksimaal: PSI/kPa: Maximum PSI/kPa: Pressione massima: (PSI/kPa): Didžiausias slėgis: PSI/kPa: Maksimali PSI/kPa vertiba: Maximum psi/kPa: Maks: psi/kPa: Maksymalne ciśnienie w PSI/kPa: PSI/kPa maxima: PSI/kPa maxima: Valdeara: PSI/kPa maxima: Maksimaalnyi tlivnost na kú doim kPa: Maximum tlak psi/kPa: Najveća vrednost PSI/kPa: Maksimaal: pritisk u PSI/kPa: Högsta PSI/kPa: Maksimum PSI/kPa: Максимальный тиск (psi/kPa): 最大 PSI/kPa
	Do not bend: Не изгибать: Neohytteje: Maikke bajes: Nicht knicken: Mqy to kqumtut: No doblar: Mitte panutada: Ala la vuta: Ne püet: Ne savjati: Tilos meghajlítani: Noli peregati: Neleinkli: Nelodli: Niet buigen: Skat ikke bajes: Nie wyginac: Nao dobrar: Nao dobrar: A nu se îndoi: He cívdató: Neohybaje: Ne upogibati: Ne savijati: Far inte bajes: Bukmeyir: He arnitiati: 〇

DxTerity™

Diagnostic Catheter

1. Device description

DxTerity diagnostic catheters are available in various configurations and

2. Indications

The diagnostic catheter is indicated for cardiac and vascular procedures. It is designed to deliver radioopaque media, guidewires, and therapeutic agents to selected sites in the vascular system. The different configurations of the diagnostic catheter are designed to be used in arteries from access sites such as the radial, brachial, and femoral arteries.

3. Contraindications

None specific to this device.

4. Warnings and precautions

Caution: Read all instructions carefully. Failure to follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

4.1. General

- For single-patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a new site of contamination of the device, which could result in patient injury, illness, or death. Cleaning, disinfection, and resterilization may compromise the essential material and design characteristics of the device and lead to device failure.
- Do not use the diagnostic catheter if it has been damaged or any other anomaly is observed.
- Do not use after the use-by date.
- Use aseptic technique during the entire procedure.
- Warning:** Administer appropriate anticoagulation therapy to reduce potential thrombosis. If the patient is not appropriately anticoagulated, thrombus formation may occur. Testing without anticoagulation shows variable amounts of thrombus formation on the device surface.
- After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.
- Only physicians trained in intravascular angiographic procedures should use or supervise the use of the diagnostic catheter.
- Carefully inspect the diagnostic catheter before use to verify that the shape and condition of the catheter are suitable for the intended use.
- Store diagnostic catheters in a cool, dry, and dark area.
- Do not expose diagnostic catheters to light.
- The diagnostic catheter is supplied sterile and nonpyrogenic. Do not use the catheter if the catheter package has been opened or damaged before use.
- Do not heat or bend the diagnostic catheter before use.
- Select a diagnostic catheter with optimal tip shape and optimal size, taking into account the access site, the target site, and the patient's anatomy.

4.2. Interactions

- When using a drug or device with the diagnostic catheter, be aware of the properties and characteristics of the drug or device and exercise due caution to avoid damaging the diagnostic catheter.
- Avoid using alcohol, antiseptic solutions, and other solvents with the diagnostic catheter.

4.3. Preparation

- Apply appropriate anticoagulant therapy to the patient.
- To prevent or reduce clotting and air embolism, always fill the diagnostic catheter with either flushing solution or contrast medium before introducing it into the vascular system. Consider the use of systemic heparinization.
- Use a guidewire with the diagnostic catheter to ensure easy introduction of the catheter.
- The diagnostic catheter can accommodate the maximum guidewire size specified on the label.
- Use coated guidewires with the diagnostic catheter to enhance thromboresistance and lubricity.

4.4. Procedure

- Exercise caution when using 5 Fr and smaller diagnostic catheters and when preoperative imaging indicates that aortic disease is a reasonable expectation. Passing a 5 Fr or smaller diagnostic catheter through an aortic lesion may make detection of aortic disease difficult.
- To avoid damaging the diagnostic catheter after it has been advanced into the vessel, manipulate the guidewire carefully, particularly when negotiating a curve in the diagnostic catheter and when passing through the tip of the diagnostic catheter.
- Reinforce the diagnostic catheter shaft with a guidewire when torque is applied. Excessive torque without the use of a guidewire may result in kinking or twisting, which may make retrieval difficult.
- If resistance is met when advancing or withdrawing the diagnostic catheter, do not proceed until the cause has been determined using fluoroscopic examination.
- Exercise care when removing guidewires from multiple-curve diagnostic catheters.
- Before starting infusion, verify that the catheter has not been kinked or occluded. Failure to do so may cause the catheter to break, rupture, or separate, resulting in damage to the vessel.
- Do not exceed the pressure ratings specified on the label.
- Performance data for the diagnostic catheter is provided in Table 1. Do not exceed the maximum permissible injection pressure.

Instructions for Use English 3

Table 1. Diagnostic catheter

Flow rate (ml/sec)	Low rate	
	Maximum Flow Rate (ml/sec)	Maximum Permissible Pressure (psikPa)
1 (1.67 mm)	600 psi/412 kPa	200 psi/8274 kPa
	230 psi/8274 kPa	
	600 psi/412 kPa	200 psi/8274 kPa
	1200 psi/8274 kPa	

All measurements were taken using 100 cm length diagnostic catheter sample using a solution of 50% isovue-300/50% saline by volume at 27 °C

5. Adverse events

Adverse events that may occur or require intervention include, but are not limited to, the following:

- Thrombus formation
- Embolization
- Infection
- Arterial, venous, or heart wall damage or perforation
- Plaque dislodgement
- Hematoma at puncture site
- Cardiac arrhythmia
- Myocardial infarction
- Vascular occlusion
- Death
- Access site pain
- Allergic and neurological reactions
- Vasospasm
- Stroke and Transient ischemic Attack
- Aneurysm
- Arteriovenous fistula

6. Instructions for use

- Carefully remove the diagnostic catheter from the package
- Flush the diagnostic catheter with heparinized saline solution
- The diagnostic catheter may be introduced into the selected vessel by the Seldinger technique, Seldinger technique, or introducer methods. When using an introducer, ensure that it is a French size greater than or equal to the selected diagnostic catheter
- Introduce the diagnostic catheter into the vascular system over a guidewire and carefully advance it under fluoroscopic control
- To remove the diagnostic catheter, reverse the guidewire to straighten the tip of the diagnostic catheter. Slowly withdraw the diagnostic catheter and guidewire under fluoroscopic control. Medtronic recommends using aspiration during withdrawal

7. Warranty disclaimer

ALTHOUGH THE DIAGNOSTIC CATHETER, HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC, INC. AND AFFILIATES (COLLECTIVELY, "MEDTRONIC"), HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CAUSE FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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
[Handwritten signature]

CERTIFICATE OF AUTHENTICITY AND COPY

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On this 16th day of March 2018, I certify that the attached document is a true, exact, complete, and unaltered copy made by me of Appendix to IFU,
presented to me by Varsha Kamat

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Subscribed and sworn to before me this the
16th day of March 2018

Heather P Morose, Notary Public
Commonwealth of Massachusetts
My commission expires on June 22, 2023

Appendix to the Instruction for use (version #M052478T001 Rev 1D)

DxTerity Diagnostic Angiographic Catheters

DxTerity diagnostic angiographic catheters in the set of versions:

1. DxTerity diagnostic angiographic catheters Angio-Kit kit, consisting of:

- DxTerity diagnostic catheter, length of catheter 100 cm – 2.PC
- DxTerity diagnostic catheter, length of catheter 110 cm – 1 PC.

models:

AKSTAN501, AKSTAN502, AKSTAN601, AKSTAN602:

2. DxTerityTRA diagnostic angiographic catheters, consisting of:

- DxTerity TRA diagnostic catheter, the length of the catheter 100 cm – 1 PC.
- DxTerity diagnostic catheter, length of catheter 110 cm – 1 PC.

models:

AKUNIV501, AKUNIV502, AKUNIV503, AKUNIV504, AKUNIV601, AKUNIV602, AKUNIV603, AKUNIV604

1.1. With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity, location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;

Developer:

Medtronic Vascular,
3576 Unocal Place,
Santa Rosa,
CA 95403,
USA.

Manufacturer:

Medtronic, Inc.
710 Medtronic Parkway, Minneapolis, MN 55432, USA
www.medtronic.com
Tel: 763-514-4000

Addresses of Manufacturing Sites for medical device:

Availmed S.A. de C.V.,
C. Industrial Lt. 001,
Mz. 105 No 20905 Int. A,
Col. Cd. Industrial,
Tijuana,
Baja California,
22444 Mexico.

Available upon request

DxTerity Diagnostic Angiographic Catheters

1.2. With respect to authorized representative of the manufacturer of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity

Authorized Representative of Manufacturer:

Medtronic Ireland, Parkmore Business Park West, Galway, Ireland.

Tel: +353 91 708000

Fax: +353 91 757524

Holder of Registration Certificate in Russia

"Medtronic" LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10

Tel: (495) 580-73-77, Fax: (495) 580-73-78

www.medtronic.ru

2. Technical characteristics of the medical device

DxTerity Angio-Kit	AKSTAN501	DxTerity™ Angio-Kit	100cm 100cm 110cm	5	• Judkins Left 4.0 • Judkins Right 4.0 • Straight Pigtail	• 5F 100cm DxTerity Judkins Left 4.0 Catheter • 5F 100cm DxTerity Judkins Right 4.0 Catheter • 5F 110cm DxTerity Straight Pigtail Catheter
	AKSTAN502	DxTerity™ Angio-Kit	100cm 100cm 110cm	5	• Judkins Left 4.0 • Judkins Right 4.0 • 145° Pigtail	• 5F 100cm DxTerity Judkins Left 4.0 Catheter • 5F 100cm DxTerity Judkins Right 4.0 Catheter • 5F 110cm DxTerity 145° Pigtail Catheter
	AKSTAN601	DxTerity™ Angio-Kit	100cm 100cm 110cm	6	• Judkins Left 4.0 • Judkins Right 4.0 • Straight Pigtail	• 6F 100cm DxTerity Judkins Left 4.0 Catheter • 6F 100cm DxTerity Judkins Right 4.0 Catheter • 6F 110cm DxTerity Straight Pigtail Catheter
	AKSTAN602	DxTerity™ Angio-Kit	100cm 100cm 110cm	6	• Judkins Left 4.0 • Judkins Right 4.0 • 145° Pigtail	• 6F 100cm DxTerity Judkins Left 4.0 Catheter • 6F 100cm DxTerity Judkins Right 4.0 Catheter

Available upon request

DxTerity Diagnostic Angiographic Catheters

DxTerity TRA Angio Kit						• 6F 110cm DxTerity 145° Pigtail Catheter
	AKUNIV501	DxTerity™ TRA Angio-Kit	100cm 110cm	5	• TRapease 4.0 • Straight Pigtail	• 5F 100cm DxTerity TRA TRapease 4.0 Catheter • 5F 110cm DxTerity Straight Pigtail Catheter
	AKUNIV502	DxTerity™ TRA Angio-Kit	100cm 110cm	5	• TRapease 4.0 • 145° Pigtail	• 5F 100cm DxTerity TRA TRapease 4.0 Catheter • 5F 110cm DxTerity 145° Pigtail Catheter
	AKUNIV503	DxTerity™ TRA Angio-Kit	100cm 110cm	5	• TRAnsformer 3.5 • Straight Pigtail	• 5F 100cm DxTerity TRA TRAnsformer 3.5 Catheter • 5F 110cm DxTerity Straight Pigtail Catheter
	AKUNIV504	DxTerity™ TRA Angio-Kit	100cm 110cm	5	• TRAnsformer 3.5 • 145° Pigtail	• 5F 100cm DxTerity TRA TRAnsformer 3.5 Catheter • 5F 110cm DxTerity 145° Pigtail Catheter
	AKUNIV601	DxTerity™ TRA Angio-Kit	100cm 110cm	6	• TRapease 4.0 • Straight Pigtail	• 6F 100cm DxTerity TRA TRapease 4.0 Catheter • 6F 110cm DxTerity Straight Pigtail Catheter
	AKUNIV602	DxTerity™ TRA Angio-Kit	100cm 110cm	6	• TRapease 4.0 • 145° Pigtail	• 6F 100cm DxTerity TRA TRapease 4.0 Catheter • 6F 110cm DxTerity™ 145° Pigtail Catheter
	AKUNIV603	DxTerity™ TRA Angio-Kit	100cm 110cm	6	• TRAnsformer 3.5 • Straight Pigtail	• 6F 100cm DxTerity™ TRA TRAnsformer 3.5 Catheter

Available upon request

DxTerity Diagnostic Angiographic Catheters

						• 6F 110cm DxTerity Straight Pigtail Catheter
	AKUNIV604	DxTerity™ TRA Angio-Kit	100cm 110cm	6	• TRAnsformer 3.5 • 145° Pigtail	• 6F 100cm DxTerity TRA TRAnsformer 3.5 Catheter • 6F 110cm DxTerity 145° Pigtail Catheter

3. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)

Operation condition:

Exercise caution when using 5 Fr and smaller diagnostic catheters and when preoperative imaging indicates that ostial disease is a reasonable expectation. Passing a 5 Fr or smaller diagnostic catheter through an ostial lesion may make detection of ostial disease difficult.

To avoid damaging the diagnostic catheter after it has been advanced into the vessel, manipulate the guidewire carefully, particularly when negotiating a curve in the diagnostic catheter and when passing through the tip of the diagnostic catheter.

Reinforce the diagnostic catheter shaft with a guidewire when torque is applied. Excessive torque without the use of a guidewire may result in kinking or twisting, which may make retrieval difficult.

If resistance is met when advancing or withdrawing the diagnostic catheter, do not proceed until the cause has been determined using fluoroscopic examination.

Exercise care when removing guidewires from multiple-curve diagnostic catheters.

Before starting infusion, verify that the catheter has not been kinked or occluded. Failure to do so may cause the catheter to break, rupture, or separate, resulting in damage to the vessel.

Do not exceed the pressure ratings specified on the label.

Temperature, °C From +32 to + 42

Humidity, % Up To 100

Relative pressure, kPa From 700 to 1060

Store diagnostic catheters in a cool, dry, and dark area.

Do not expose diagnostic catheters to UV light.

Do not heat or bend the diagnostic catheter before use.

Temperature from +20 °C to +25 °C

Relative Humidity of 20% -80%

Pressure from 100kPa to 103 kPa

4. List of national standards applied by the manufacturer

DxTerity diagnostic angiographic catheters in the set of versions is compliant with active national standards mentioned in the Declaration of Conformity.

5. Warnings and/or precautions taken by the consumer when disposing of the medical device, accessories and consumables used together with it (if any), including information about the infectious, microbial, environmental or physical hazard of the medical device

When using a drug or device with the diagnostic catheter, be aware of the properties and characteristics of the drug or device and exercise due caution to avoid damaging the diagnostic catheter.

Available upon request

Appendix to the Instruction for use

DxTerity Diagnostic Angiographic Catheters

Avoid using alcohol, antiseptic solutions, and other solvents with the diagnostic catheter.

6. Procedure and conditions for the disposal or destruction of a medical device.

After use, this product may be a potential biohazard (hazard class B). Handle and dispose of all such device is accordance with accepted medical practice and applicable local, regional and national laws and regulations.

8. Shelf life

The Use by Date is indicated on the product label for each CFN.

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© 2016 Medtronic