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
I, Celine Curran, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Finisklin Business Park, Sligo, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: CLL

Date: 29 May 2019

Celine Curran,
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Finisklin,
Sligo,
Ireland

**Abbott Ireland Diagnostics Division
Finisklin Business Park
Sligo
County Sligo
Ireland**

**SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH**

May Hart 29/5/2019

Authorised Official

Created November 2016.

NAME

Alinity i CMV IgG Avidity Controls (also referred to as CMV Avidity Ctrls)

INTENDED USE

The Alinity i CMV IgG Avidity Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer.

For additional information, refer to the Alinity i CMV IgG Avidity reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CONTROL L and **CONTROL H** contain recalcified human plasma.

Preservatives: ProClin 950 and sodium azide.

The controls are at the following ranges:

Control	Quantity	Color	CMV IgG Avidity
			RANGE (%Avi)
CONTROL L	1 x 8.0 mL	Natural	< 50.0
CONTROL H	1 x 8.0 mL	Natural	≥ 60.0

NOTE: The controls have been designed exclusively for use with this assay and assay format and therefore are not suitable for use with other assays. The controls cannot be run as patient specimens.

If multiple control replicates are required, note that only one avidity determination can be ordered per sample cup.

The controls are manufactured and released within specific ranges and not to a specified target value. Therefore lot-to-lot differences are expected. Levey-Jennings may be utilized for data tracking, however application of Westgard Rules is not recommended due to the stated lot-to-lot variability.

STANDARDIZATION

The Alinity i CMV IgG Avidity assay is standardized against internal reference material.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

Safety Precautions

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in the Low and High Control is reactive for CMV IgG, and nonreactive for HBsAg, HIV-1 Ag or HIV-1 RNA, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to:

CONTROL L and **CONTROL H**

WARNING	Contains methylisothiazolones and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.

Prevention

P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response

P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

Disposal

P501	Dispose of contents / container in accordance with local regulations.
------	---

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 8 drops of the low control and 8 drops of the high control into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Control Number
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Control (Low, High)
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Ireland
	Range
	List Number

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Finlekin Business Park
Sligo
Ireland
+353-71-9171712

0843

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Created November 2016.

©2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i CMV IgG Avidity Controls for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the maximum shelf life of

Alinity i CMV IgG Controls Avidity – 12 months from date of manufacturing

The dating of each Alinity i CMV IgG Avidity Reagent Kit and Controls master lot is determined by the kit component with the shortest expiration dating.

Device mustn't be used if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

LIST OF APPLICABLE STANDARDS

EN ISO 9001

ISO 13485

EN ISO 14971

EN 13641

EN 13612

EN ISO 15223-1

EN 980

EN ISO 18113-1

EN ISO 18113-2

EN ISO 18113-3

EN ISO 23640

EN ISO 20776-1

EN 13975

EN ISO 17511

POTENTIAL USER

The Alinity i CMV IgG assay is intended for use in a variety of hospital clinical laboratory environments including, but not limited to; hospitals and commercial laboratories, whose personnel have been trained to use the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

TECHNICAL CHARACTERISTICS

Density at 20 °C:

- Low Control, High Control: 1,028 g/cm³

pH-value at 20 °C:

- Low Control, High Control: 7,3

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

MATERIALS OF ANIMAL/HUMAN ORIGIN

Materials of animal/human origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS

Not applicable.

CONTRAINDICATIONS

Not applicable.

QUALITY CONTROL PROCEDURES

The recommended control requirement for the Alinity i CMV IgG Avidity assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

- If multiple control replicates are required, note that only one avidity determination can be ordered per sample cup.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect.

Follow the established quality control procedures for your laboratory. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

- Review quality control results and acceptance criteria following a change of reagent lot.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

MATERIALS REQUIRED BUT NOT PROVIDED

- Alinity i CMV IgG Avidity assay file package
- 07P4310 Alinity i CMV IgG Avidity Reagent Kit
- 07P42 Alinity i CMV IgG Reagent Kit
- 07P4201 Alinity i CMV IgG Calibrators
- 07P4210 Alinity i CMV IgG Controls or other control material
- 07P44 Alinity i CMV IgM Reagent Kit (optional)
- 07P4401 Alinity i CMV IgM Calibrator (optional)
- 07P4410 Alinity i CMV IgM Controls or other control material (optional)
- 06P1160 Alinity Trigger Solution
- 06P1265 Alinity Pre-Trigger Solution
- 06P1368 Alinity i-series Concentrated Wash Buffer

ANALYTICAL CHARACTERISTICS

DESCRIPTION	PRODUCT REQUIREMENT
Within Laboratory (Total) Precision	• $\leq 14\%$ CV for samples representing the grayzone avidity range (50.0% Avi to 59.9% Avi)

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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