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
I, Celine Curran, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Finisklin Business Park, Sligo, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Celine Curran

Date: 29 May 2019

Celine Curran,
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Finisklin,
Sligo,
Ireland

Abbott Ireland Diagnostics Division
Finisklin Business Park
Sligo
County Sligo
Ireland

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Mary Harte 29/5/2019 **Authorised Official**

Revised February 2017.

REF 07P4322

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity i CMV IgG Avidity Reagent Kit (also referred to as CMV Avidity)

INTENDED USE

The Alinity i CMV IgG Avidity assay is a chemiluminescent microparticle immunoassay (CMIA) used for the determination of the avidity of IgG antibodies to Cytomegalovirus in human serum and plasma on the Alinity i analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Infections with Cytomegalovirus (CMV), a member of the herpesvirus family, are common in man and are usually mild and asymptomatic.¹ However, in pregnant women,² newborns,³ and immunocompromised individuals⁴ CMV infections may pose a significant medical risk. CMV infection remains difficult to diagnose on symptoms alone since a high percentage of infections remains asymptomatic. *In utero* infection may result in sequelae of varying degrees including mental retardation, chorioretinitis, hearing loss and neurologic problems.⁵ Since the risk of *in utero* virus transmission and CMV related damage of the fetus is markedly increased during primary infection, reliable recognition of primary CMV infection is of high importance for pregnant women.⁵ An individual may undergo primary infection with CMV, reinfection with exogenous virus or reactivation of latent virus. Although presence of anti-CMV IgG reduces the likelihood of CMV related complications, it does not assure complete protection from disease.

CMV acquisition in infants can occur transplacentally following maternal infection, during birth by contact with the virus excreted from the cervix⁶ or following birth through the ingestion of infected maternal breast milk.⁷ Both seronegative individuals and infants may acquire CMV through infected blood products^{3, 8} or contact with an infected individual.^{9, 10}

The functional binding affinity or avidity of IgG antibodies increases progressively over time after immunization, also known as maturation of the humoral immune response. High percentage of low avidity IgG antibodies may indicate a primary infection whereas high percentage of high avidity IgG antibodies may indicate a recurrent infection.

The Alinity i CMV IgG Avidity assay is a qualitative method for the determination of the avidity of IgG antibodies to Cytomegalovirus in human serum and plasma. It is used as an aid in the differentiation between primary and non-primary infection.

If primary infection needs to be excluded, CMV IgG reactive samples should be tested for CMV IgM and CMV IgG Avidity. A positive CMV IgM result in combination with a low avidity result is a strong indicator for a primary CMV infection within the last 4 months.

CMV IgG	CMV IgM	CMV IgG Avidity	Indication for...
nonreactive	nonreactive	N/A	no infection
reactive	nonreactive	high avidity	past infection; low risk for <i>in utero</i> transmission
reactive	reactive	low avidity	primary infection; high risk for <i>in utero</i> transmission
reactive	reactive	high avidity	non-primary infection; low risk for <i>in utero</i> transmission

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay consists of 2 single tests that are both two-step immunoassays for the determination of Cytomegalovirus in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology. The avidity of anti-CMV IgG in the sample is calculated using relative light units (RLUs) of both tests.¹¹ An aliquot of the sample and Pre-Treatment Reagent 2 (blocking agent) are combined. A second aliquot of sample is combined with Pre-Treatment Reagent 1 (buffer). CMV virus lysate (strain AD169) coated paramagnetic microparticles are combined with each aliquot of pretreated sample and incubated. The mixture is washed. Murine anti-human IgG acridinium-labeled conjugate is added to create a reaction mixture and incubated. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). The % Avidity is calculated from the RLUs obtained from the sample pretreated with the blocking agent and the RLUs obtained from the unblocked sample.

To determine the avidity of a sample, a reactive Alinity i CMV IgG result of the respective sample is required to enable the Alinity i analyzer to select the correct dilution protocol for the Alinity i CMV IgG Avidity assay. This is ensured by ordering automated assay panels consisting of at least Alinity i CMV IgG and Alinity i CMV IgG Avidity. For details refer to the INSTRUMENT PROCEDURE section.

It is strongly recommended to test only CMV IgM reactive samples with the Alinity i CMV IgG Avidity assay. This is imperative for samples with a low Alinity i CMV IgG concentration value between 6 to 15 AU/mL.

Refer also to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS, Specimen Types section.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity i CMV IgG Avidity Reagent Kit 07P43

NOTE: This product is composed of 4 components, which are packaged as a 2 cartridge reagent set. Both cartridges are required to perform the assay.

Volumes (mL) listed in the table below indicate the volume per cartridge set.

REF	07P4322
Tests per cartridge set	100 (50 avidity determinations)
Number of cartridge sets per kit	2
Tests per kit	200 (100 avidity determinations)
MICROPARTICLES	6.6 mL
CONJUGATE	6.1 mL
PRE-TREATMENT 1	5.9 mL
PRE-TREATMENT 2	4.2 mL
MICROPARTICLES CMV virus lysate (strain AD169) coated microparticles in TRIS buffered saline. Minimum concentration: 0.08% solids. Preservatives: ProClin 300 and antimicrobial agents.	
CONJUGATE Murine anti-human IgG acridinium-labeled conjugate in MES buffer. Minimum concentration: 44 ng/mL. Preservatives: sodium azide and antimicrobial agents.	
PRE-TREATMENT 1 Pre-Treatment reagent 1 containing TRIS buffer. Preservatives: ProClin 300 and antimicrobial agents.	
PRE-TREATMENT 2 Pre-Treatment reagent 2 containing CMV virus lysate (strain AD169) in TRIS buffer. Preservatives: ProClin 300 and antimicrobial agents.	

Warnings and Precautions

- IVD**
- For In Vitro Diagnostic Use

Safety Precautions

CAUTION: This product contains human-sourced and/or potentially infectious components. Refer to the **REAGENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹²⁻¹⁵

The following warnings and precautions apply to: MICROPARTICLES / PRE-TREATMENT 1 / PRE-TREATMENT 2	
WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.
The following warnings and precautions apply to: CONJUGATE	
WARNING	Contains polyethylene glycol octylphenyl ether and sodium azide.
H319	Causes serious eye irritation.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
 - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.
- After mixing, place reagent cartridges in an upright position for 8 hours before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.
- Do not have more than one lot of reagents on board.**

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 8 hours before use.
Onboard	System Temperature	30 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a control value is out of the specified range. Associated test results are invalid, and samples must be retested.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity i CMV IgG Avidity assay files must be installed on the Alinity i analyzer prior to performing the assay.

Assay panels, calculated assay files, and retest rules will be installed automatically when installing CMV Avidity (assay file 3000) together with the following required assay files:

- CMV IgG R
- CMV IgM
- CMV Avi1
- CMV Avi2

Do not modify the assay parameters for assays starting with "zz" or "CMV Avi".

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay on the ARCHITECT i System.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum Serum separator
Plasma	Lithium heparin Lithium heparin (Plasma separator) Potassium EDTA Sodium citrate Sodium heparin

- Performance has not been established for the use of cadaveric specimens or the use of bodily fluids other than human serum and plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- To determine the avidity of a sample, a reactive Alinity i CMV IgG result of the respective sample is required to enable the Alinity i analyzer to select the correct dilution protocol for the Alinity i CMV IgG Avidity assay. This is ensured by ordering automated assay panels consisting of at least Alinity i CMV IgG and Alinity i CMV IgG Avidity.
- **It is strongly recommended to test only CMV IgM reactive samples with the Alinity i CMV IgG Avidity assay. This is imperative for samples with a low Alinity i CMV IgG concentration value between 6 to 15 AU/mL.**

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.
- All samples (controls and patient specimens) should be tested within 5 hours of being placed on board the Alinity i analyzer.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter
- they require repeat testing.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge at a minimum of 100 000 g-minutes.
- Examples of acceptable time and force ranges that meet this criterion are listed in the table below.

Centrifugation time using alternate RCF values can be calculated using the following formula:

$$\text{Minimum Centrifugation time (minutes)} = \frac{100\,000 \text{ g-minutes}}{\text{RCF}}$$

Recentrifugation Time (Minutes)	RCF (x g)	g-Minutes
10	10 000	100 000
20	5000	100 000
40	2500	100 000

$$\text{RCF} = 1.12 \times r_{\text{max}} (\text{rpm}/1000)^2$$

RCF -	The relative centrifugal force generated during centrifugation.
rpm -	The revolutions per minute of the rotor on which the specimens are being spun (usually the digital readout on the centrifuge will indicate the rpm).
Centrifugation Time -	The time should be measured from the time the rotor reaches the required RCF or rpm to the time it begins decelerating.
r_{max} -	Radius of the rotor in millimeters. NOTE: If custom tube adapters (i.e., adapters not defined by the centrifuge manufacturer) are used, then the radius (r_{max}) should be manually measured in millimeters and the RCF calculated.
g-minutes -	The unit of measure for the product of RCF (x g) and centrifugation time (minutes).

- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen storage conditions were verified on the ARCHITECT i System.

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	Room temperature (15 to 30°C)	3 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.
	2 to 8°C	14 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.

Remove serum or plasma from the clot, red blood cells, or separator gel if stored longer than the maximum 2-8°C or room temperature storage time and store frozen at -10°C or colder.

Avoid more than 6 freeze/thaw cycles.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P43 Alinity i CMV IgG Avidity Reagent Kit

Materials Required but not Provided

- Alinity i CMV IgG Avidity assay file package
- 07P4310 Alinity i CMV IgG Avidity Controls
- 07P42 Alinity i CMV IgG Reagent Kit
- 07P4201 Alinity i CMV IgG Calibrators
- 07P4210 Alinity i CMV IgG Controls or other control material
- 07P44 Alinity i CMV IgM Reagent Kit (optional)
- 07P4401 Alinity i CMV IgM Calibrator (optional)
- 07P4410 Alinity i CMV IgM Controls or other control material (optional)
- 06P1160 Alinity Trigger Solution
- 06P1265 Alinity Pre-Trigger Solution
- 06P1368 Alinity i-series Concentrated Wash Buffer


For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- It is strongly recommended to use the automated assay panels for Alinity i CMV IgG Avidity testing since this will ensure the correct selection of CMV assay testing steps and dilution protocols required for avidity determinations. Refer to figure below. Please contact Abbott for information on ordering the CMV assays using an LIS (Host).

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- Alinity i CMV IgG, CMV IgM, and one determination of CMV IgG Avidity can be sampled from the same sample cup. If multiple avidity determinations are required for the same sample, place the respective number of differently barcoded sample cups with aliquots of the same sample on the instrument.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 1
 - Priority:
 - Sample volume for each test: 300 µL for the complete Alinity i CMV assay panel testing using one sample cup
 - ≤ 3 hours on the reagent and sample manager:
 - Sample volume for each test: 300 µL for the complete Alinity i CMV assay panel testing using one sample cup
 - > 3 hours on the reagent and sample manager:
 - Additional sample volume is required. For information on onboard sample storage, refer to the Alinity ci-series Operations Manual, Section 4.

- If the result of the initial test meets the retest criteria of the assay panel, the Alinity i analyzer generates a rerun order. Tests scheduled for rerun are displayed on the Orders screen (All Orders tab and Rerun tab) and Sample Status screen. The R (rerun) code is assigned to the test. The system can be configured to automatically reposition the samples for retest. If the system is not configured for automatic repositioning, samples to be retested may need to be reloaded.
- A rerun order for the "CMV IgG R" assay will be placed for samples that have an Alinity i CMV IgG result of > 250 AU/mL. These samples will be tested automatically in a 1:10 dilution performed by the instrument. The D (dilution) and R (rerun) code are assigned to the test.
- Refer to the Alinity i CMV IgG Avidity control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples must not be diluted manually prior to running the Alinity i CMV IgG Avidity assay.

Calibration

Calibration is not required for the Alinity i CMV IgG Avidity assay.

Quality Control Procedures

The recommended control requirement for the Alinity i CMV IgG Avidity assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

- If multiple control replicates are required, note that only one avidity determination can be ordered per sample cup.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent lot.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

The following results are reported depending on the panel chosen:

- Alinity i CMV IgG result in AU/mL.
- Alinity i CMV IgM result in Index (or S/CO).
- Alinity i CMV IgG Avidity ("CMV Avi") result in %Avi.

NOTE: Results with the unit of measurement "AdjRLU" or results from assay files starting with "zz" may be reported. These results are intermediate and not relevant for sample evaluation.

Calculation

The Alinity i analyzer calculates the avidity for the Alinity i CMV IgG Avidity assay using the following calculation and stores the result. Avidity [%] = 100 x (1- (CMV Avi1/CMV Avi2)).

Interpretation of Results

Result	Interpretation
< 50.0 %Avi	low avidity
50.0 - 59.9 %Avi	grayzone
≥ 60.0 %Avi	high avidity

No clinical interpretation can be drawn from a grayzone result. It is recommended to take a second sample within an appropriate period of time (e.g. 2 weeks) and repeat testing.

For details on configuring the Alinity i analyzer to use grayzone and high avidity interpretations, refer to the Alinity ci-series Operations Manual, Section 2.

The grayzone and high avidity interpretations are editable parameters and should be utilized per end user requirements.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, results should be used in conjunction with other data; e.g., results of other tests (CMV IgG, CMV IgM), clinical impressions, etc.
- If the CMV IgG Avidity results are inconsistent with clinical evidence, additional testing is recommended to confirm the result.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.¹⁶ Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA).^{17, 18} Specimens containing HAMA may produce anomalous values when tested with assay kits that employ mouse monoclonal antibodies.¹⁶

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity i analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.¹⁹ Testing was conducted using 3 lots of the Alinity i CMV IgG Avidity Reagent Kit and 3 lots of the Alinity i CMV IgG Avidity Controls and 1 instrument. Two controls and 3 human plasma panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Sample	n	Mean (%Avi)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD (Range ^b)	%CV (Range ^b)
Low Avidity Control	355	22.2	2.45	11.1	2.47 (1.97-3.19)	11.2 (7.6-14.7)
High Avidity Control	355	87.8	1.13	1.3	1.16 (0.60-1.81)	1.3 (0.7-2.1)
Panel 1	356	12.2	2.94	24.2	2.97 (2.48-3.77)	24.4 (19.3-31.9)
Panel 2	355	54.9	1.96	3.6	2.25 (1.44-3.27)	4.1 (2.6-6.0)
Panel 3	355	89.0	1.33	1.5	1.37 (0.44-2.28)	1.5 (0.5-2.6)

^a Includes within-run, between-run, and between-day variability.

^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Specificity

This study was performed on the ARCHITECT i System.

A total of 492 specimens from pregnant females, diagnostic/hospitalized patients and randomly selected volunteer blood donors were evaluated. 33 specimens grayzone on either assay were not included in the calculation of specificity. Data for specificity are summarized in the following table.

Sample Type	ARCHITECT CMV IgG Avidity		Commercially-Available Diagnostic Kit	
	Observed Specificity	Lower 95% Confidence Limit	Observed Specificity	Lower 95% Confidence Limit
Pregnant Females	100.0% (399/399)	99.1%	99.5% (397/399)	98.2%
Diagnostic/ Hospitalized Patients	96.8% (61/63)	89.0%	98.4% (62/63)	91.5%
Blood Donors	96.7% (29/30)	82.8%	100.0% (30/30)	88.4%
Total	99.4% (489/492)	98.2%	99.4% (489/492)	98.2%

Sensitivity

This study was performed on the ARCHITECT i System.

The study included specimens from patients with documented primary CMV infection. Three commercially-available seroconversion panels with a total number of 49 bleeds, of which 37 bleeds represent the acute CMV infection phase, were tested. The sensitivity was assessed to be 97.3% (36/37) on ARCHITECT CMV IgG Avidity (lower 95% confidence limit: 85.8%) compared to 75.7% (28/37) sensitivity on a commercially-available diagnostic kit (lower 95% confidence limit: 58.8%). The following table shows data from the seroconversion panels tested.

Panel Member	Days after first draw	ARCHITECT CMV IgG [AU/mL]	ARCHITECT CMV IgM [Index]	ARCHITECT CMV IgG Avidity [%Avi]	Commercially- Available Diagnostic Kit [% Avidity]
Low Avidity/ Nonreactive		< 6.0	< 0.85	< 50.0	< 35
Grayzone		N/A	0.85-0.99	50.0-59.9	35-45
High Avidity/ Reactive		≥ 6.0	≥ 1.00	≥ 60.0	> 45
RP-003-01	1	8.1	1.60	1.0	N/A
RP-003-02	4	21.8	4.71	2.4	N/A
RP-003-03	8	51.1	5.70	10.5	N/A
RP-003-04	51	86.2	1.74	25.4	24.1
RP-003-05	55	94.1	1.44	24.7	31.2
RP-003-06	59	91.2	1.29	30.4	27.9
RP-003-07	65	84.0	1.20	31.8	29.6
RP-003-08	67	85.1	1.18	32.7	30.3
RP-003-09	72	80.8	1.08	37.9	50.2
RP-003-10	74	78.4	1.05	39.9	28.4
RP-003-11	79	80.3	1.07	40.5	31.0
RP-003-12	84	85.1	1.07	39.4	43.3
RP-003-13	88	94.7	1.15	40.3	26.7
RP-003-14	95	87.1	1.05	45.1	45.7
RP-003-15	99	89.4	1.02	44.9	44.1
RP-019-01	1	1.2	0.26	N/A	N/A
RP-019-02	5	1.1	0.24	N/A	N/A
RP-019-03	8	1.1	0.25	N/A	N/A
RP-019-04	12	1.1	0.27	N/A	N/A
RP-019-05	15	1.1	0.24	N/A	N/A
RP-019-06	21	1.4	0.35	N/A	N/A
RP-019-07	26	1.8	0.41	N/A	N/A
RP-019-08	29	3.3	0.99	N/A	N/A
RP-019-09	33	10.9	4.22	1.2	N/A
RP-019-10	36	18.6	6.30	8.3	N/A
RP-019-11	43	76.5	7.05	12.9	14.0
RP-019-12	50	88.1	6.06	9.2	12.9
RP-019-13	57	98.0	5.30	18.8	13.0
RP-019-14	68	107.4	4.11	30.1	13.7
RP-019-15	75	109.9	3.68	21.9	18.1
RP-019-16	82	108.0	2.98	22.8	18.5
RP-019-17	86	107.8	3.19	28.3	17.8
RP-019-18	89	102.7	2.87	29.5	18.2
RP-019-19	96	106.4	2.66	28.3	24.3
RP-019-20	104	100.1	2.50	29.9	22.1
RP-019-21	109	93.8	2.24	34.1	24.5
RP-019-22	113	101.4	2.43	31.9	21.6
RP-019-23	116	108.9	2.66	32.4	20.5
RP-019-24	121	101.8	2.12	39.1	21.3
RP-019-25	124	104.2	2.18	37.6	18.2
PTC901-01	0	3.2	0.15	N/A	N/A
PTC901-02	2	2.0	0.12	N/A	N/A
PTC901-03	7	1.9	0.13	N/A	N/A
PTC901-04	10	1.7	0.13	N/A	N/A
PTC901-05	17	3.6	3.09	N/A	16.1
PTC901-06	24	42.8	5.29	15.8	16.8
PTC901-07	32	77.0	4.87	25.3	17.7
PTC901-08	51	101.1	4.50	22.7	17.3
PTC901-09	67	126.1	4.07	29.7	19.3

N/A in the table indicates that anti-CMV IgG was either not detectable (on a commercially-available diagnostic kit) or the result was nonreactive (on the ARCHITECT CMV IgG assay). Therefore no avidity result could be generated.

Percent Agreement

A study was performed to determine the percent agreement between the Alinity i CMV IgG Avidity and a commercially-available CMV IgG Avidity assay by testing 80 specimens from a population consisting of 35 low avidity and 44 high avidity specimens. One grayzone specimen was excluded from the evaluation. The specimens were tested in singlicate on 3 Alinity i analyzers using 3 lots of the Alinity i CMV IgG Avidity Reagent Kit and 1 lot of Alinity i CMV IgG Avidity Controls. The specimens were also tested in singlicate using a commercially-available CMV IgG Avidity assay.

For each reagent lot, the positive percent agreement was 100% (35/35) and the negative percent agreement was 100% (44/44). The lower 95% confidence interval (of the assumed comparator assay specificity of 100%) was 91.96%. The lower 95% confidence interval (of the assumed comparator assay sensitivity of 100%) was 90.00%.

Alinity i CMV IgG Avidity	Commercially-Available CMV IgG Avidity Assay	
	Low Avidity	High Avidity
Low Avidity	35	0
High Avidity	0	44

Positive % agreement = 100.00% (35/35) 95% CI (90.00 - 100.00%)
Negative % agreement = 100.00% (44/44) 95% CI (91.96 - 100.00%)

Interference

Potentially Interfering Endogenous Substances

This study was performed on the ARCHITECT i System.

No interference was observed between experimental controls and low avidity or high avidity specimens tested with elevated levels of the following potentially interfering substances:

Potentially Interfering Substance	Interferent Level
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 3000 mg/dL
Protein	≤ 12 g/dL
Hemoglobin	≤ 500 mg/dL

Potentially Interfering Other Conditions

This study was performed on the ARCHITECT i System.

Potential interference was further assessed on specimens from patients with diseases unrelated to CMV infection. High avidity results can be expected since all specimens were anti-CMV IgG reactive and anti-CMV IgM nonreactive thus indicating the absence of a primary infection. Data are summarized in the table below.

Interfering Substance	Total per Category (N)	High Avidity [N]	Grayzone Avidity [N]	Low Avidity [N]
Anti-nuclear antibody	9	9	0	0
Epstein-Barr virus	8	8	0	0
Influenza vaccinees	2	2	0	0
Human anti-mouse antibody	5	4	0	1
Herpes simplex virus type 1	9	8	1	0
Herpes simplex virus type 2	9	8	1	0
Hyperpolyclonal IgG	9	9	0	0
Hyperpolyclonal IgM	8	8	0	0
Measles	6	6	0	0
Rheumatoid factor	6	5	1	0
Systemic lupus erythematosus	6	5	1	0
Varicella zoster virus	7	7	0	0

BIBLIOGRAPHY

- Hecker M, Qui D, Marquardt K, Bein G, Hackstein H. Continuous cytomegalovirus seroconversion in a large group of healthy blood donors. *Vox Sang.* 2004 Jan;86(1):41-44.
- Munro SC, Hall B, Whybin LR, Leader L, Robertson P, Maine GT, Rawlinson WD. Diagnosis of and screening for cytomegalovirus infection in pregnant women. *J Clin Microbiol.* 2005 Sep;43(9):4713-4718.
- Cannon MJ, Davis KF. Washing our hands of the congenital cytomegalovirus disease epidemic BMC Public Health. 2005 Jun 20;5:70.
- Pass R, Griffiths C, August A. Antibody Response to Cytomegalovirus after Renal Transplantation: Comparison of Patients with Primary and Recurrent Infections. *J Infect Dis* 1983;147:40-46.
- Lazzarotto T, Gabrielli L, Lanari M, Guerra B, Bellucci T, Sassi M, Landini MP. Congenital cytomegalovirus infection: recent advances in the diagnosis of maternal infection. *Hum Immunol.* 2004 May;65(5):410-415.
- Reynolds D, Stagno S, Hosty T, Tiller M, Alford C. Maternal Cytomegalovirus Excretion and Perinatal Infection. *N Engl J Med* 1973;289:1-5.
- Stagno S, Reynolds D, Pass R, et al. Breast milk and the risk of cytomegalovirus infection. *N Engl J Med* 1980;302:1073-1076.
- Bowden R, Sayers M, Flournoy N, Newton R, Banaji M, Thomas E, et al. Cytomegalovirus Immune Globulin and Seronegative Blood Products to Prevent Primary Cytomegalovirus Infection After Bone Marrow Transplantation. *N Engl J Med* 1986;314(16):1006-1010.
- Yow M, White N, Taber L, et al. Acquisition of cytomegalovirus infection from birth to 10 years: a longitudinal serologic study. *J Pediatr* 1987;10:37-42.
- Adler S. Molecular epidemiology of cytomegalovirus: viral transmission among children attending a day care center, their parents, and caretakers. *J Pediatr* 1988;122:366-371.
- Maine, G.T., Hsu S.T. et al. Methods of the determination of antibody IgG avidity. Patent PCT/US2006/042818.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

■ Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONJUGATE	Conjugate
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
INVERSIONS PERFORMED	Inversions Performed
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TREATMENT 1	Pre-Treatment 1
PRE-TREATMENT 2	Pre-Treatment 2
PRODUCT OF IRELAND	Product of Ireland
REF	List Number
SN	Serial number

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Finisklin Business Park
Sligo
Ireland
+353-71-9171712

0843

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised February 2017.

©2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i CMV IgG Avidity Reagent Kit for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the maximum shelf life of

Alinity i CMV IgG Avidity Reagent Kit – 12 months from date of manufacturing

The dating of each Alinity i CMV IgG Reagent Kit, Calibrators, and Controls master lot is determined by the kit component with the shortest expiration dating.

Device mustn't be used if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

LIST OF APPLICABLE STANDARDS

EN ISO 9001

ISO 13485

EN ISO 14971

EN 13641

EN 13612

EN ISO 15223-1

EN 980

EN ISO 18113-1

EN ISO 18113-2

EN ISO 18113-3

EN ISO 23640

EN ISO 20776-1

EN 13975

EN ISO 17511

POTENTIAL USER

The Alinity i CMV IgG assay is intended for use in a variety of hospital clinical laboratory environments including, but not limited to; hospitals and commercial laboratories, whose personnel have been trained to use the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

TECHNICAL CHARACTERISTICS

Density at 20 °C:

- Microparticles: 1,047 g/cm³
- Conjugate: 1,031 g/cm³
- Pretreatment 1,2: 1,047 g/cm³

pH-value at 20 °C:

- Microparticles: 8,4
- Conjugate: 6,6
- Pretreatment 1,2: 8,4

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

MATERIALS OF ANIMAL/HUMAN ORIGIN

Materials of animal/human origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

INDICATIONS

Not applicable.

CONTRAINDICATIONS

Not applicable.

ANALYTICAL CHARACTERISTICS

ASSAY NAME	ASSAY LIST NUMBER	DESCRIPTION	PRODUCT REQUIREMENT
CMV IgG Avidity Serum/Plasma	07P43	Within Laboratory (Total) Precision	• $\leq 14\%$ CV for samples representing the grayzone avidity range (50.0% Avi to 59.9% Avi)
		Sensitivity	• The sensitivity shall be equal to or better than that of a commercially available IVD CMV IgG Avidity assay when tested on the same population with documented primary CMV infection, defined as time between the last seronegative and seropositive bleeds is less than or equal to 4 months
		Specificity	• The specificity shall be equal to or better than that of a commercially available IVD CMV IgG Avidity assay when tested on the same population consisting of CMV immune blood donors, diagnostic/hospitalized patients and pregnant women that are all negative for CMV IgM antibody and positive for CMV IgG
		Interference	The assay shall have no qualitative performance difference for the tested potentially interfering substances within the ranges of less than 50.0% Avi and equal to or greater than 60.0% Avi.

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- Hecker M, Qui D, Marquardt K, Bein G, Hackstein H. Continuous cytomegalovirus seroconversion in a large group of healthy blood donors. *Vox Sang.* 2004 Jan;86(1):41-44.
- Munro SC, Hall B, Whybin LR, Leader L, Robertson P, Maine GT, Rawlinson WD. Diagnosis of and screening for cytomegalovirus infection in pregnant women. *J Clin Microbiol.* 2005 Sep;43(9):4713-4718.
- Cannon MJ, Davis KF. Washing our hands of the congenital cytomegalovirus disease epidemic BMC Public Health. 2005 Jun 20;5:70.
- Pass R, Griffiths C, August A. Antibody Response to Cytomegalovirus after Renal Transplantation: Comparison of Patients with Primary and Recurrent Infections. *J Infect Dis* 1983;147:40-46.
- Lazarotto T, Gabrielli L, Lanari M, Guerra B, Bellucci T, Sassi M, Landini MP. Congenital cytomegalovirus infection: recent advances in the diagnosis of maternal infection. *Hum Immunol.* 2004 May;65(5):410-415.
- Reynolds D, Stagno S, Hosty T, Tiller M, Alford C. Maternal Cytomegalovirus Excretion and Perinatal Infection. *N Engl J Med* 1973;289:1-5.
- Stagno S, Reynolds D, Pass R, et al. Breast milk and the risk of cytomegalovirus infection. *N Engl J Med* 1980;302:1073-1076.
- Bowden R, Sayers M, Flournoy N, Newton R, Banaji M, Thomas E, et al. Cytomegalovirus Immune Globulin and Seronegative Blood Products to Prevent Primary Cytomegalovirus Infection After Bone Marrow Transplantation. *N Engl J Med* 1986;314(16):1006-1010.
- Yow M, White N, Taber L, et al. Acquisition of cytomegalovirus infection from birth to 10 years: a longitudinal serologic study. *J Pediatr* 1987;10:37-42.
- Adler S. Molecular epidemiology of cytomegalovirus: viral transmission among children attending a day care center, their parents, and caretakers. *J Pediatr* 1988;122:366-371.
- Maine. G.T., Hsu S.T. et al. Methods of the determination of antibody IgG avidity. Patent PCT/US2006/042818.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

16. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
17. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
18. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
19. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.

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