Unet J. al

APPROVED


August 3, 2018

L.S.

INSTRUCTIONS
FOR USE OF A MEDICAL DEVICE

Dental aspiration device for removal of liquids
and particles, reusable, with accessories.

II Edition

refle sale utano
forampt

EURONDA S.p.A.

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Tel: 0445.329811 - Fax 0445.865246

Partita IVA 00595740242

1. Name of the medical device:

Dental aspiration device for removal of liquids and particles, reusable, with accessories, in the following versions:

- cannula MONOART® EM19
- cannula MONOART® EM21
- cannula MONOART® EM21 Aspirator Tip
- cannula MONOART® EM40 Aspirator Tip

Accessories:

- adapter for cannula MONOART® EM40 Aspirator Tip

hereinafter referred to as the **devices**, the cannula, cannula EM19, cannula EM21, cannula EM21 Aspirator Tip and cannula EM40 Aspirator Tip.

2. Manufacturer, authorized representative of the manufacturer

2.1 Name and address of the manufacturer:

EURONDA S.P.A., Italy

Via dell' Artigianato - 36030 Montebelluna Precalcino (VI), Italy

tel. +39 0445 329811, fax: +39 0445 865246

web-site: info@euronda.com

2.2 Authorized representative of the manufacturer

PROTECO, JSC

Registered office: 57 11 Salova Str., Saint Petersburg 192102, Russia

tel./fax: +7-812-635-8890

web-site: www.protecodent.ru

e-mail: info@protecodent.ru

3. Purpose of the medical device:

The cannula is designed to remove liquids (saliva) and waste (particles) from the patient's mouth cavity to ensure maximum effectiveness of any dental procedures.

The device may be reused if duly sterilized (autoclaved).

The dental aspiration device for removal of liquids and particles, reusable, with accessories is recommended to the following customers:

healthcare professionals performing various dental manipulations or procedures

4. Indications, contraindications and precautions associated with the intended use of the medical device

The dental aspiration device for removal of liquids and particles, reusable, with accessories is used for removing a limited amount of liquid (saliva) and suspended solid particles from the patient's mouth during any types of dental manipulations or procedures.

Contraindications: The use of the devices is contraindicated in the event of individual intolerance to the materials of the device.

Use the devices only as intended (see the application scope).

Read the instructions for use before use. Do not use defective devices or if the package is damaged!

Sterilize the device before reuse. Disinfect the device and perform the pre-sterilization cleaning after use.

Adverse effects associated with the intended use of the medical device: irritation or allergic reaction caused by individual intolerance to the materials of the device.

5. Application scope:

Application scope - dental office at various healthcare facilities, private dental practices and dental clinics, including private dental clinics.

formation on the procedure of use

1 Application

Application: to ensure maximum effectiveness of any dental procedures. The use of cannulas ensures suction (aspiration) of **liquids** and waste (particles) from any parts of the patient's mouth cavity, including the most difficult to access parts.

5.2 Rules of application

5.2.1 Mode of application: open the package and remove the cannula or the adapter. Install the device into the aspiration and/or separation system of any dental treatment unit, diameter 10 mm. It is approved for use in the Russian Federation, by connecting the device with the outlet port of the saliva aspirator's hitch on the dental treatment unit. While working, direct the cannula at the processed surface to fully remove liquids or particles.

5.2.2. Cannulas and adapters are shipped non-sterile and thus shall be fully sterilized according to the instructions for use before the first and each subsequent use.

After use, place the cannula or the adapter in a special container filled with a disinfecting solution. Maintain the concentration and the duration specified by the manufacturer. Use only disinfectants and detergents compatible with polypropylene (PP) or **ethylene-propylene copolymers**. Observe instructions for the concentration, the temperature and the duration specified by the manufacturer of disinfectants and detergents.

5.2.3 Pre-sterilization cleaning (PSC). Remove the largest impurities from the cannula or the adapter and then rinse it with flowing water.

5.2.4. Manual cleaning:

- Clean the cannula with a nylon brush and a detergent-disinfectant compatible with polypropylene (PP) or ethylene-propylene copolymer until all impurities have been removed. Do not use metal brushes to clean cannulas or adapters. Use a laboratory brush to clean the interior.
- Thoroughly **rinse the** cannula or the adapter with flowing water to fully remove residual detergent-disinfectant.
- Dry with compressed air.

5.2.5 Ultrasonic cleaning:

- Put the **cannula or** the adapter in a mesh basket and place it in an ultrasonic bath filled with a detergent or a disinfectant compatible with polypropylene (PP) or ethylene-propylene copolymer.
- Set the cleaning cycle according to the instructions provided by the manufacturer. The temperature **should** not exceed +45°C.
- Thoroughly rinse the device with flowing water to fully remove residual detergent-disinfectant.
- **Dry with compressed air.**

5.2.6 Sterilization. Choose a sterilization Kraft paper bag fit with the cannula or the adapter and leave enough space to prevent stress when sealing the bag. The devices shall be autoclaved with water-saturated steam in fractional vacuum observing the following conditions:

- use a type B cycle.
- sterilization temperature 134°C.
- sterilization duration 5 minutes.
- drying 10 minutes.

5.3 Special instructions and precautions for use

Read the instructions for use thoroughly. Use the cannula only as **intended** (see the application scope). Store new or unused cannulas in a dry, clean and secure place. Visually **examine the** cannula before each use to detect signs of wear and tear. If the cannula is damaged, do not use it and dispose of it immediately.

For professional dental use.

EURONDA S.p.A.

Via del Annunziato, 7 - 36100 MONTECCHIO P.
Tel 0445 329811 ra - Fax 0445 865246
Partita IVA 00505740242

6 Description of the medical device and primary technical characteristics

6.1 Description

The cannula is a bent tube without sharp edges to prevent the risk of irritating or injuring the oral mucosa which helps to smoothly remove saliva and particles.

6.2 Description of the accessories intended for use in combination with the medical device

The cannula MONOART® EM 40 Aspirator Tip features an accessory adapter:
Adapter for cannula MONOART® EM 40 Aspirator Tip
Adapters are designed to connect aspirator tips and dental saliva aspirators with aspiration system tubing (diameter 16 mm) of dental treatment units.
Adapter for cannula MONOART® EM40 Aspirator Tip, inner diameter 11.2-10.6 mm, outer diameter 16.1-15.0 mm

6.3 Primary technical characteristics

Various sizing options, smooth cannula interior and aerodynamic form aid in the patient's comfort.

Dimensions of cannulas:

cannula MONOART® EM19 - diameter 16 mm, length 114 mm

cannula MONOART® EM21 - diameter 16 mm, length 129 mm

cannula MONOART® EM21 Aspirator Tip - diameter 11 mm, length 150 mm

cannula MONOART® EM40 Aspirator Tip - diameter 11 mm, length 200 mm

6.3.1 Cannulas MONOART® EM19, MONOART® EM21 and MONOART® EM21 Aspirator Tip are manufactured in the following colors: blue, gray, lilac, orange, yellow, green, lime, pink, burgundy.

6.3.2 Cannula MONOART® EM40 Aspirator Tip and cannula adapter are manufactured in gray color.

6.3.3 Ergonomic non-slip handles facilitate operation, while linear form helps to maintain hygiene and perform cleaning.

6.4 Recommendations:

Cannula MONOART® EM19 provides great comfort to adult patients owing to round shape. Cannula MONOART® EM21 is recommended for procedures in difficult-to-access cavities. Cannula MONOART® EM40 Aspirator Tip is a universal shaped tip for surgical aspirators. The device is shorter and bent at the end for better control during operations with patients and quieter aspiration. Furthermore, smaller size helps to reduce the risk of swallowing reflexes. 6.5 The cannula has an aerodynamic form, does not have sharp edges, defects, nicks, cracks or burrs.

7 Installation and commissioning procedure

Prevent the device from mechanical damage and contamination during use.
Use personal protection equipment (PPE) according to the activities performed during use.

8 Peculiarities of interaction with the environment (requirements to premises and personnel qualification)

The devices are safe for dental personnel of various healthcare facilities.

9 Conditions of transportation and storage

9.1 Transportation

Cannulas may be transported by any enclosed means of transportation according to the applicable shipping rules.

A carton package or a bag ensures integrity of the device during transportation.

- Ambient temperature: -50°C – +40°C

- Relative humidity: 45-80%

EURODENT S.P.A.
Via S. Rocco 10, 20139 Milano (MI)
Tel. +39 02 5749 1111
Fax +39 02 5749 1112
www.eurodent.it

EURODENT S.P.A.

Storage

The package must be protected from dust, moisture and the risk of recontamination.
Store new or unused cannulas in a dry, clean and secure place.

- Ambient temperature: +15 °C – +25 °C
- Relative humidity: 45-60%

10 Shelf life

The manufacturer ensures compliance of the medical device to the requirements affecting quality, effectiveness, safety and functional purpose provided that the customer does not violate conditions of use, transportation and storage.

The warranty storage life is 5 years provided that conditions of storage and transportation are met.
Shelf life is unlimited.

11 Packaging

11.1 Cannulas are packaged in consumer packaging: carton packages or ziplock (gripper) bags made of polyethylene film. 10 pcs. per package


Adapters are packaged in consumer packaging: ziplock (gripper) bags made of polyethylene film.

11.2 Packaged cannulas and adapters are placed in corrugated boxes ensuring integrity of devices during transportation.

12 Labeling

12.1 Consumer (group) packaging labeling:

- manufacturer's trademark;
- name and address of the manufacturer;
- name and address of the authorized representative of the manufacturer;
- name of the device;
- article number, color, size;
- amount;
- date of production;
- lot number (LOT);
- shelf life;
- storage conditions;
- method and mode of sterilization, number of cycles;
- registration certificate number and date;
- CE mark;
- statement: "For professional dental use"

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12.2 Explanation of symbols used in group (consumer) packaging:

Symbol	Explanation	Symbol	Explanation
	Date of production		Manufacturer
	Lot code		Temperature range (if required)
	Article number, catalog number		Keep dry (if required)
	Non-sterile		Consult the instructions for use
	CE mark		GOST R compliance
	Keep away from sunlight (if required)		Do not dispose in a household waste bin

13 Scope of supply

13.1 The scope of supply of the device must conform to table No. 3

		Table No. 3
No.	Name	Amount (pcs.)
1	Cannula* MONOART*	10
2	Instructions for use	1

* versions and color of the device are selected from the available range

13.2 Adapter for cannula MONOART® EM40 Aspirator Tip is not included in the scope, it is supplied as per request.

14 Applied standards

The devices comply with the essential requirements of EU Council Directive 93/42/EEC dated June 14, 1993 (Medical Devices Directive) and Directive 2007/47 CE.

The cannulas comply with requirements of ISO 9001:2015, ISO 13485:2016, ISO 14971:2012, DIN EN 13060:2015, DIN EN 285:2016, EN 1041:2013, ISO 15223-1:2016, ISO 10993-1:2009, ISO 10993-5:2009, ISO 10993-10:2010, ISO 10993-12:2012 and ISO 17664:2004.

15 Procedure for disposal and destruction of the device

15.1 Disposal shall be performed in compliance with norms and requirements in force in the Russian Federation.

In terms of biological, toxicological and radiation hazards, the medical devices "Dental aspiration device for removal of liquids and particles, reusable, with accessories" are classified as biologically hazardous medical waste: materials and tools, items contaminated with blood and/or other body fluids.

Consumer packaging and defective goods are considered biologically safe waste similar in composition to the municipal solid waste (MSW).

15.2 After use, place the cannula or the adapter in a special container filled with a disinfecting solution. Maintain the concentration and the duration specified by the manufacturer. Use only disinfectants and detergents compatible with polypropylene (PP). Observe instructions for the

EURONDA S.p.A.

Via dell'Industria 7 - 36100 MONTECCHIO PAV.
Tel. 0445 929811 - Fax 0445 865246

... IVA 00595740442

EURONDA S.p.A.

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concentration, the temperature and the duration specified by the manufacturer of disinfectants and detergents.

Disinfected cannulas or adapters shall be collected and transferred to special licensed medical waste disposal organizations.

Do not dispose of the used devices with regular household waste!

16 Interaction with other medical devices (or general-purpose devices)

This medical device is designed to be used in combination with other devices or equipment. The entire combination, including the connection system, shall be safe and not impair characteristics of the device.

17 Radiation protection

The device does not emit potentially hazardous levels of radiation.

18 Rules of procedure for complaints




If adverse effects not specified in the instructions for use of the device, adverse reactions to the use of the device, peculiarities of interaction of medical devices, facts and circumstances threatening life and health of patients and healthcare professionals during application and use of this device have been detected, address a message containing this information to the competent healthcare authority organization and the authorized representative of the manufacturer:

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