ReSound Aventa 3

ReSound Aventa 3 Installation Guide

ReSound

rediscover hearing

Table of Contents

PC Requirements for Aventa 3	3
Installation FAQ	4
NOAH vs. Operating systems	5
Installing Aventa 3 software	6
Installing .NET Framework	12
Airlink installation	16
Speedlink installation	17
Hi-PRO USB installation	17
Noahlink Installation	17
AutoPlay function	18
Requirement when running 64-bit	19
Windows 8, 7 and Vista on 64 & 32-bit	20

PC Requirements for Aventa 3

Before installing Aventa 3 it is necessary to make sure that the machine meets the requirements for running this fitting software application. These requirements are summarized in the tables below.

ReSound Aventa 3 Minimum PC requirements	ReSound Aventa 3 Recommended PC requirements
1 GHz or higher	2.93 GHz or higher
2 MB RAM	4 GB RAM
3 GB Free Hard Disk Space	6 GB Free Hard Disk Space
USB Port	USB Port
DVD-ROM Drive	DVD-ROM Drive
Sound Card: 16 bit, stereo, Microsoft DirectSound and DirectSound 3D compatible	Sound Card: 16 bit, stereo, Microsoft DirectSound and DirectSound 3D compatible
1400x1050 Screen Resolution (Full) 1920x1080 Screen Resolution (Wide)	1400x1050 Screen Resolution (Full) 1920x1080 Screen Resolution (Wide)
Printer	Printer
XP Professional (SP3 or higher)	Windows 7 32-bit (SP1 or higher)
* NOAH 3.1.2 Patch E or higher	* NOAH 4.4

* NOAH 2 is not supported

Slower performance should be expected on systems that meet minimum PC requirements.
Performance can be optimized by adherence to recommended PC requirements.

Supported Operating systems		
Windows XP	32-bit	SP3 or higher
Windows Vista	32 & 64-bit*	SP1 or higher
Windows 7	32 & 64-bit*	All
Windows 8 (Except Windows 8 RT)	32 & 64-bit**	All

* NOAH 3.7 is needed to run 64-bit

** NOAH 4 is needed to run on Windows

Old ReSound fitting software Not working on 64-bit
Resource II
Aventa 1.x
Danafit
Auditfit

Installation FAQ

Before installing Aventa 3

We recommend that you browse through these FAQ's before installing. You should be able to find the answer to most of your questions here.

What happens if my PC doesn't meet the recommended requirements for Aventa 3?

If you for example have the recommended Aventa 2.9 PC requirements but you are running Aventa 3 then the software will have a slower launch, it will take longer to get to initial fit etc. ReSound Aventa 3 will simply not perform as intended.

Will the installation of Aventa 3 affect Aventa 2.9 in any way if the computer meets the recommended requirements?

No. Aventa 3 and Aventa 2.9 can coexist on the same computer.

Can Aventa 3 run on my Mac computer?

No. Aventa is not compatible with Mac computers and will not be able to run on these.

Is Aventa 3 compatible with 64-bit operating systems?

Yes. Aventa 3 can be installed on 64-bit operating systems. Additionally, a 64-bit Speedlink driver will automatically be loaded during such Aventa 3 installations.

Use of NOAH 3.6.1 Patch G is recommended. Does this mean that I can no longer use NOAH 3.1.2 Patch E (or other earlier version of Noah)?

No. But, Noah 3.6.1 Patch G and subsequent version of Noah utilize a more robust database which greatly reduces the risk of data corruption.

Does Aventa 3 need to be installed before plugging in the Airlink?

Yes. The Airlink driver is installed during the Aventa 3 installation. Once Aventa 3 is installed, Airlink can be plugged in and recognized by the computer.

Is .NET is required to run Aventa 3?

Yes. Aventa will check the version of .NET on your computer during the Aventa 3 installation. If your computer already has the proper version of .NET, the Aventa 3 installation will proceed. If it is determined that .NET must be installed or updated, a .NET installation screen will appear and you will be guided through the .NET install.

NOAH vs. Operating systems

	Windows XP Pro (32-bit)	Windows Vista (32-bit)	Windows Vista (64-bit)	Windows 7 (32/64-bit)	Windows 8 (32/64-bit)
Noah 2	NO	NO	NO	NO	NO
Noah 3.1.2 E	OK	NO	NO	NO	NO
Noah 3.5.2 C	OK	NO	NO	NO	NO
Noah 3.6.1 G	OK	OK	NO	NO	NO
Noah 3.7	OK	OK	OK	OK	NO
Noah 4.4	NO	OK	OK	OK	OK



Installing Aventa 3 Software

If the PC requires a .NET Framework update you will be notified about this during the Aventa 3 installation. Should you wish to follow the .NET installation process go to page 12.

If installing Aventa 3 on Windows Vista or Windows 7 or 8 please see page 20 before inserting the DVD.

1. Close all applications that may be running on the PC, including NOAH.
2. Any previously installed version of Aventa 3 does not need to be manually un-installed prior to installation of a newer Aventa 3 version. Previous installation is automatically removed.
3. Insert **Aventa 3 DVD** and follow the on-screen instructions as prompted.

Aventa will automatically start up the installation process.


Occurs only when AutoPlay is enabled


Installing Aventa 3 Software

The Aventa 3 InstallShield Wizard will now prepare for installation.

4. Aventa 3 is now ready to be installed on your PC. Click **Next** to start the installation.



Installing Aventa 3 Software

5. System Requirements screen inform you if the system requirements are met. Click **Next**.

If the recommended requirements are not met you will see the yellow warning icon displayed. This does not mean that you cannot install Aventa 3 – this is for your information, that only the minimum system requirements are met.


6. Accept the License Agreement and click **Next**.


Installing Aventa 3 Software

7. Select your country in the dropdown. If a Standalone version is required, check the Standalone box. If Web Updates functionality is preferred, make sure the Allow Web Updates box is checked. Click **Next**.

8. This screen displays the sound files enable for the specific market version. If sound file language or availability need to be changed then do so and click **Next**.

Note: Versions for markets without language specific speech files available will have the non-speech sound files set as default.


Installing Aventa 3 Software

9. Ready to install Aventa 3. Click the **Install** button.


The installation of Aventa 3 now continues. This may take several minutes.

Installing Aventa 3 Software

11. A prompt will occur, indicating that Aventa 3 was successfully installed. Click **Finish** to exit the wizard.

Note: If .NET Framework 4 was installed during this installation and requires a reboot, the installation will prompt you to restart the machine. Click Yes on the Restart prompt to complete the installation and restart your machine.


Installing .NET Framework

For Windows Vista or XP users

It is required to have both **.NET Framework 3.5 SP1** and the **.NET Framework 4** installed. If you begin installation of Aventa 3, without both **.NET Framework 3.5 SP1** and the **.NET Framework 4** installed, here is the procedure to be followed for proper .NET Framework installation.


1. Insert the Aventa 3 DVD

If you have AutoRun enabled, the **.NET Framework 3.5 SP1** will automatically start up the installation process. The following screen will appear.


2. If AutoRun is not enabled, go to the CD/DVD drive and click the **Setup.exe** file. The **.NET Framework 3.5 SP1** installations will now begin.

The process of extracting files may be running for a few minutes.


The setup process will take a couple of minutes. The screen to the right will be displayed during this timeframe.


Installing .NET Framework

4. Accept the License Agreement and click **Install**.


Download & Install progress screen will be displayed.

Installing .NET Framework

5. A message will appear, indicating that Microsoft .NET Framework 3.5 SP1 has been installed successfully. Click **Exit**.


6. Go to the CD/DVD drive again and run **Setup.exe**. If .NET Framework 4 or higher is not installed on the PC, the following message is shown. Click **OK** to install the .NET Framework 4.


No license agreement is required for installation of .NET Framework 4.

The setup process will take a couple of minutes. The screen below will be displayed during this timeframe.

7. When completing the .NET Framework 4 installation, the Aventa 3 installer is automatically initiated and prompts for language selection.

Installing .NET Framework

For Windows 7 users





If you install the **Aventa 3 DVD** without the **.NET Framework 4** installed, a message will appear.

8. Click **OK** and **.NET Framework 4** will automatically start the installation process.

The process of extracting files may be running for a few minutes. Do not click **Cancel** during **.NET Framework 4** installation.

The setup process will take a couple of minutes. The screen to the right will be displayed during this timeframe.

9. When completing the **.NET Framework 4** installation, the Aventa 3 installer is automatically initiated and prompts for language selection. For installation of Aventa 3, go back to page 6.

Airlink installation

ReSound Airlink is a USB plug and play programming interface enabling you to fit wireless ReSound hearing instruments.

1. Before you insert Airlink, Aventa 3 must be installed, as the Airlink driver is installed during fitting software installation.
2. a. Airlink: For the best possible wireless performance it is strongly recommended to place the Airlink upright in a USB cradle. Please note that the USB hub must be USB 2.0 (or higher) compliant.



The illustrations display the optimal setup for wireless fitting with Airlink.

- b. Airlink 2: Place Airlink 2 on the table and insert USB cable to the USB port.

3. After inserting Airlink you can launch Aventa 3 and start to fit wireless ReSound products.

Do not attempt to use Airlink with any other fitting software unless explicitly instructed. Airlink can co-exist on the same PC together with programming interfaces (Speed-link, Hi-PRO, NOAHlink). It is possible, however, that some wireless programming interfaces from other hearing instrument manufacturers might interfere with the Airlink while it is in use.

Speedlink installation


Speedlink is ReSound's high-speed hearing instrument programming interface.

To use Speedlink with your Aventa 3 fitting software make sure Aventa 3 is installed.


1. Plug the Speedlink cable into any USB port on your PC.
2. Launch Aventa 3
3. After installation: Make sure Speedlink programming interface is chosen in the Connect dialog or in Edit preferences (Machine preferences)

For further information on Speedlink compatibility and technical specification please see the Speedlink user guide that came with you Speedlink.

In the Connect dialog:

Edit preferences (Machine preferences):

Hi-PRO USB installation

Please follow installation guide delivered with your Hi-PRO USB. Instructions can also be downloaded from Otometrics website.

<http://www.otometrics.com/>

Instructions can be found under:

Products > Solution for: > Fitting > Hi-PRO USB

Noahlink Installation

Please follow installation guide delivered with your NOAHlink. Instructions can also be downloaded from Himsa website.

<http://www.himsa.com>

Instructions can be found under:

Support > NOAHlink Knowledgebase > Installing > Installing

AutoPlay function

If you have AutoPlay enabled on your PC, the Aventa 3 software will automatically start up when you insert the disk.


The AutoPlay function will take you directly to the installation process.

If you wish to view or print content from the disk, please following procedure below:

How to see content on the Aventa 3 DVD

Go to your CD/DVD drive and right click.

Click on open.

Requirements when running 64-bit

Before installing Aventa 3 it is necessary to make sure that the machine and programming interface meets the requirements for running this fitting software application on 64-bit versions of Windows. This information is summarized below.

	Windows Vista (64-bit) Aventa 2.9	Windows 7 & 8 (64-bit) Aventa 2.9	Windows Vista (64-bit) + Aventa 3.x	Windows 7 (64-bit) Aventa 3.x	Windows 8 (64-bit) Aventa 3.x
Airlink	Not supported	Not supported	Works fine on 64-bit No new driver is required Remember to insert Airlink before you open Aventa 3.x	Works fine on 64-bit No new driver is required Remember to insert Airlink before you open Aventa 3.x	Works fine on 64-bit No new driver is required Remember to insert Airlink before you open Aventa 3.x
Speedlink	Driver for 64-bit required Driver Zip file GN ReSound Speedlink Instruction Speedlink 64-bit install doc. If Aventa 3.x is already installed you don't need separate 64-bit driver	Driver for 64-bit required Driver Zip file GN ReSound Speedlink Instruction Speedlink 64-bit install doc. If Aventa 3.x is already installed you don't need separate 64-bit driver	64-bit Driver is automatically installed when installing Aventa 3.x	64-bit Driver is automatically installed when installing Aventa 3.x	64-bit Driver is automatically installed when installing Aventa 3.x
NOAHlink	Works fine on 64-bit Minimum version 01.54.01 is needed	Works fine on 64-bit Minimum version 01.54.01 is needed	Works fine on 64-bit Minimum version 01.54.01 is needed	Works fine on 64-bit Minimum version 01.54.01 is needed	Works fine on 64-bit Minimum version 01.54.01 is needed
Hi-Pro (Serial)	Driver for 64-bit required Driver Zip file Hi_pro_ installation_8-49-8080_01	Driver for 64-bit required Driver Zip file Hi_pro_ installation_8-49-8080_01	Driver for 64-bit required Driver Zip file Hi_pro_ installation_8-49-8080_01	Driver for 64-bit required Driver Zip file Hi_pro_ installation_8-49-8080_01	Driver for 64-bit required Driver Zip file Hi_pro_ installation_8-49-8080_01
Hi-Pro (USB)	"Only" if using Aurical See instruction: 64-bit - Aurical with built-in Hi-PRO	"Only" if using Aurical See instruction: 64-bit - Aurical with built-in Hi-PRO	"Only" if using Aurical See instruction: 64-bit - Aurical with built-in Hi-PRO	"Only" if using Aurical See instruction: 64-bit - Aurical with built-in Hi-PRO	"Only" if using Aurical See instruction: 64-bit - Aurical with built-in Hi-PRO
Noah 3.7	Works fine on 64-bit No patch is needed	Works fine on 64-bit No patch is needed	Works fine on 64-bit No patch is needed	Works fine on 64-bit No patch is needed	Not Supported Minimum requirement is Noah 4.0

Administrator Privileges on Windows 8, 7 and Vista on 64 & 32-bit

Before installing Aventa 3 on Windows 8, 7 or Vista with 64 or 32-bit, it's recommended that you install Aventa with Administrator privileges. This information is summarized below.

The recommendation is to install Aventa with Administrator privileges.
If you're not sure of your user account configuration, or you are using Network configuration, please ask your local IT support for guidance.

On new 64 & 32-bit machines you cannot start the Setup in the regular way.

Right click on Setup and choose **"Run as Administrator"**
Now follow normal procedure, when prompted on the screen.

If you observe a problem with Administrator rights, it might be because **"UAC"**
(User Account control) is enabled within the operating system.
Please change/disable the function.

How to change UAC in System configuration on Windows 8 and 7


Go to start, Click Run, Type: **"msconfig"**, go to Tools, Scroll down to Change UAC Settings, Click Launch and drag the slider down to 'Never notify' and click 'OK'.
Remember to restart the machine afterwards.

How to disable UAC in System configuration on Windows Vista

Go to start, Click Run, Type: **"msconfig"**, go to Tools, Scroll down to Disable UAC and click OK.
Remember to restart the machine afterwards.

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The translation of this text was performed by me, the translator Mamedov Timur Djavanshirovich. The requirements regarding the text of the translation (maximum correctness, literate statement) have been explained to me.

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I, Martinova Natalia Andreevna, deputy notary of the city of Moscow for Akimov Gleb Borisovich, certify the authenticity of signature, made by the translator Mamedov Timur Djavanshirovich in my presence. His identity is established.

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It is registered in the register under №
Fee collected: 100 rubles

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(signature)

Seal: Notary of Moscow Akimov
ITN 770400047857

Stitched up, numbered and sealed () pages

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In this guide you will find:

HOW TO START


Airlink	3
Other programming interfaces	4

GOOD TO KNOW

Start screen	5
Product selection	5
Pre-Fit screen	5
Fit screen	6
- Gain Adjustment	6
- Advanced Features	6
- Beeps and Volume Control	6
Summary screen	7
Feature Demo Animations	7








PRECAUTIONS AND WARNINGS

General Precautions.	8
General Warnings.	9

Every time you fit a client, it's an opportunity to build your reputation. When you choose ReSound, we understand the trust you've placed in us. We strive to provide you with flexible, simple-to-fit software so you can make a real difference in the lives of people with hearing loss.

Airlink

ReSound Aventa 3 supports two generations of Airlink; the Airlink dongle and the latest Airlink 2.


Both generations offer the great client experience of wireless fitting and the comfort it brings.

Being the latest generation, Airlink 2 is the recommended option of the two. It's designed to ensure optimal usage and better usability.

- Place Airlink 2 on the table with a clear line-of-sight to the hearing instruments, which should be within a range of 10 feet (3 metres).
- Avoid placing the Airlink 2 in a USB hub with other USB devices (e.g. Bluetooth dongle), as this can decrease the efficiency of the Airlink 2.
- When fitting hearing instruments inside a sound booth place the Airlink 2 inside or close to the booth.
- It is recommended not to use USB cables between the Airlink 2 and the PC exceeding a length of 10 feet (3 metres).

Note: If using the Airlink dongle it's recommended to use a USB hub that places the Airlink in an upright/vertical position as this provides for a better signal broadcast. Refer to the pictures to the left for more guidance.

Other programming interfaces


Aventa 3 works with:

- Airlink 2
- Airlink dongle
- Speedlink
- Hi-PRO
- NOAHlink

You can select programming interface on top of the Start Screen.

Note: Animations in Connection Assistance will show you how to connect specific models of hearing aids with the selected programming interface.


Start screen

After launching the Aventa software you will see the Start screen.


- 1 Select the patient's experience with hearing aids for an optimized gain prescription.
- 2 Let Aventa recommend the program defaults or select them based on your patient interview.
- 3 Connect to the hearing aids chosen for this fitting.
- 4 Use AutoFit as the easiest most direct path to a fitting.

Product selection

- 1 Product family selection - here you can view all available products per family.
- 2 Use filters to help you easily find the most relevant products. For example, "Audiogram Suitable" is a filter that shows only products with fitting ranges that cover your patient's individual audiogram.
- 3 Fitting range and patient audiogram.
- 4 Technical Information on the selected hearing aid.


Pre-Fit screen

- 1 Counsel on the audiogram using the hearing loss simulator for a 3rd party present during the fitting session.
- 2 Perform an in-situ audiogram by presenting the pure tones from the hearing aids.
- 3 Reconfigure the hearing aid software to match the actual hardware.
- 4 Display audiogram overlays such as the speech banana and familiar sounds.
- 5 Have a closer look at how the hearing aids looks and how it sits on the ear.


Fit screen

Gain Adjustment


- 1 Adjust gain changes in increments of 1 dB, 2 dB or 3 dB.
- 2 Choose and fine-tune environmental programs and Wireless program.
- 3 Muting and linking the hearing aids.
- 4 Pair the hearing instruments with wireless accessories.

Advanced Features


- 1 Feature settings are accessed by clicking "Advanced Features" in the "tools" section. These feature settings apply to the currently selected program.
- 2 Adjust the features of the hearing aids as needed.

Beeps and Volume Control


- 1 Click here to adjust beeps and volume settings.
- 2 In the "Standard Beeps" tab, all the acoustic indicators can be demonstrated, changed, added and removed as needed. Low or high frequency beeps can also be selected.
- 3 "Manual Controls" tab contains options such as customizing push button function, enabling and disabling the program button and adjusting the volume control range.

Summary screen


- 1 Confirm number of programs.
- 2 Customize the Clinician report and the Patient handout, which provides the patient with useful information.
- 3 Type in notes and they will be printed on patient's handout.
- 4 Save options let you select whether to save data in both the hearing aids and database, only in one place, or not at all.

Feature Demo Animations

Demo animations illustrate the features available in the hearing aids. You can find the animations in Help menu.

General Precautions ⚠

- This Aventa fitting software will apply recommended amplification settings based on available audiometrical information. Settings will be specific to each fitting.
- Manual modification of the fitting parameters will impact and change the amplification level prescribed. This change will be audible to the patient when the instruments are connected.
- Caution should be exercised with patients sensitive to sounds, for example suffering from tinnitus or hyperacusis.
- Hearing device performance may drift and decline over time, please ensure devices are performing within specifications prior to performing in situ threshold testing.
- Feature accuracy is reliant on ear canal being completely occluded.

General Warnings

- When connected, the Aventa Fitting Software controls the acoustic amplification levels in a hearing instrument. In some acoustic environments, over amplification can cause discomfort and damage to the patients hearing.
- The Fitting System provides initial default amplification settings based on hearing threshold levels. These default settings could be higher than stable levels and cause feedback when amplification is first applied.
- The Fitting System initiates the feedback calibration procedure. Feedback calibration uses broadband noise to measure the amount of sound leaking from the hearing instrument. The output level is designed to be on the border of "uncomfortable" level based on hearing threshold level at a specific frequency. The sound will be ramped up in volume and cease when the calibration data is received. It is possible for the level to exceed a patients comfort level, but it needs to be determined if it can obtain hazardous risk levels.
- The Fitting System initiates feedback calibration to measure the receiver to microphone transfer function. This is used primarily by the hearing instrument to manage feedback suppression but it is also used to display the limits of stable gain.
- The maximum stable gain estimates are based on feedback calibration data, the presence of active feedback suppression, and a headroom estimate. The headroom value is meant to be conservative however the presence of directionality has shown to give inaccurate estimates of max. stable gain. The risk is that, under some situations, the Fitting System could show that the hearing instrument is stable when it is actually close to unstable and in risk of feedback.
- If feedback calibration has not been performed during fitting the Max Stable Gain is not known and the device could cause feedback without warning.
- The Fitting System uses the Audiogram+ algorithm to interpret audio gram data for determining optimal gain settings. The parameter used is "first time user". This may give a less than optimal initial fitting but should not pose a safety risk. The algorithm has been independently validated with the specification and verified to be consistent with Fitting Software.
- The Fitting System uses feedback calibration measurements to compute Max Stable Gain values. These values are an estimate of the amount of gain that can be safely applied to a hearing instrument before it begins to cause feedback. There is a warning when this gain is reached, and the over gain values are highlighted with bold, red text. The safety margins give a "close" approximation of the actual feedback border. However, it's an estimation and feedback can occur before the warning is given. Sustained feedback on high power devices can damage residual hearing.
- Warning to hearing care practitioners: Special care should be exercised in selecting and fitting hearing instrument(s) whose maximum sound pressure level exceeds 132 dB SPL with an IEC 60711: 1981 occluded ear simulator, because there may be a risk of impairing the remaining hearing of the hearing instrument user.
- Warning to hearing care practitioners: Special care should be exercised in selecting and fitting hearing instrument(s) utilizing Tinnitus Sound Generator. The maximum output of the tinnitus sound generator feature falls into the range that can cause hearing loss according to OSHA regulations. For further details please consult the user guide of the relevant hearing instrument that includes the Tinnitus Sound Generator feature. In accordance with NIOSH recommendations the user should not use the sound generator for more than eight (8) hours a day when this is set to a level of 85db SPL or above. When the sound generator is set to levels of 90db SPL or above the user should not use the sound generator for more than two (2) hours per day. In no case should the sound generator be worn at uncomfortable levels.
- Children and physically or mentally challenged users will require guardian supervision while wearing the device.

WARNING points out a situation that could lead to serious injuries, CAUTION indicates a situation that could lead to minor and moderate injuries.

0297 Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to ReSound A/S.

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Country: Denmark

Dette offentlige dokument / This public document

2. er underskrevet af / has been signed by
Pernille Hellesøe

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6. den 7. juni 2016

the 07 June 2016

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
by the Ministry of Foreign Affairs of Denmark

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
ReSound Unite

USER GUIDE
ReSound Unite™ TV Streamer 2

ReSound
rediscover hearing

Welcome

Thank you for selecting ReSound.

ReSound Unite TV Streamer 2 will provide you with high quality streamed audio from your TV, sound system, personal computer or similar device directly to your hearing instruments. Please read this manual carefully to avoid any problems in the operation and fully benefit from the TV Streamer 2. Consult your hearing care professional if you have any questions.


GN ReSound A/S
Lautrupbjerg 7
2750 Ballerup Denmark

Contents



Description.....	3
Getting Started	3
Pairing to your hearing instruments.....	6
Start and stop streaming.....	7
Controls.....	8
Adjustable delay mode.....	9
Important Information.....	9
Troubleshooting Guide	10
Technical specifications.....	12
Warranty and Repair	13
Acceptance and sale certificate.....	13

Description

1. Micro-USB power socket
2. Pairing button
3. Mono/stereo switch
4. Analog mini-jack audio input
5. Optical (Toslink) audio input
6. Coaxial audio input

7. Increase volume
8. Decrease volume
9. Power light (green)
10. Activity light (yellow)
11. Charger wall plug
12. Charger mini-USB plug

Getting Started

Installing the TV Streamer 2

Depending on your TV, sound system, computer or other audio setup, the installation of TV Streamer 2 may be done in several ways described in the following section. The following steps will guide you through the pairing and connecting process, which can be easily accomplished by almost anyone. Start by connecting cables between the TV streamer 2 and the audio source and finish by pairing the TV Streamer 2 to your hearing instruments.


 TV Streamer 2 supports both analog and digital audio input. Please read the following sections for information about supported audio formats and how to connect the TV Streamer 2 to audio sources.

Connecting the TV Streamer 2 to an analog audio source

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer 2.
2. Insert the male stereo mini-jack on the audio cable into the female stereo mini-jack input on the back of the TV Streamer 2 (4).
3. Plug the other end audio cord (with red and white plugs) into the "AUDIO OUT" terminals on the audio source (while matching connector colors).
4. The "AUDIO OUT" is usually marked "L – R" for left and right, and can usually be found on the back of the audio source.

Note!

With this installation, most audio sources will send sound to both hearing instruments and additional speakers connected to the audio source simultaneously. On some audio sources you can adjust the balance between the AUDIO OUT and the speakers. Consult your audio source manual for instructions.


Alternative analog audio cable installations:

If it is not possible to connect the TV Streamer 2 to the audio output of the audio source using the red and white plugs, the following installation alternatives can be used:

Mini-jack adaptor:

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer 2.
2. Insert the male stereo mini-jack on the audio cable into the female stereo mini-jack input on the back of the TV Streamer (4).
3. Connect the stereo mini-jack adaptor to the red and white plugs at the other end of the audio cable.
4. Plug the stereo mini-jack adaptor into the audio sources ("HEADPHONE") output terminal, or alternatively into your computers' audio terminal.


Note: The HEADPHONE output terminal is usually located on the side or front of the audio source. When using the HEADPHONE terminal, usually only the person wearing the hearing instruments will hear the sound, and audio source speakers are muted. If muted speakers are not desired, please refer to the first installation, or the one on the following page. Your audio source may support speakers and headphone being on simultaneously. Check your audio source manual for instructions.

SCART adapter:

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer 2.
2. Insert the male stereo mini-jack on the audio cable into the female stereo mini-jack input on the back of the TV Streamer 2 (4).
3. Connect the SCART adapter to the red and white plugs at the other end of the audio cable.
4. Plug the SCART adapter into the audio source's SCART terminal.

Note: The SCART output terminal is usually located on the back or side of the audio source. Please consult the manual of your audio source if additional configuration of the SCART output is needed.

Connecting the TV Streamer 2 to a digital audio source:


Using the optical (Toslink) audio input

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer 2.

2. Insert one of the plugs on the optical audio cable into the optical audio input socket (5) on TV streamer 2.


Note: Most optical cables are delivered with a small protective cap to protect the tip of the cable. Please remove this protective cap before inserting the cable.

3. Plug the other end of the optical audio cable into the "DIGITAL AUDIO OUT/OPTICAL" on your audio source.

Using the coaxial audio input

1. Connect the power supply to the wall socket and insert the cable into the micro-USB socket on the TV Streamer 2.
2. Connect one of the plugs on the coaxial audio cable to the coaxial plug (6 on page 8) on TV Streamer 2.
3. Plug the other end of the coaxial audio cable into the "DIGITAL AUDIO OUT/COAXIAL" on your audio source.

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i Note! Note! The following information is intended primarily for your hearing care professional, who in most cases will complete the necessary tasks before the TV Streamer 2 can be used with your hearing instruments. However, the process is designed to be simple and most people can do it themselves.

It is possible to pair up to three streaming devices (e.g. two TV Streamers and one Clip-On Microphone) with a given pair of hearing instruments. Therefore channel 1, 2 or 3 must be selected when starting the pairing process.

Pairing the TV Streamer 2 in channel 1

1. Make sure that the TV Streamer 2 power supply is plugged into a wall socket and that the green light indicator on the front of the TV Streamer 2 (9) is ON. Your hearing instruments must be OFF - open the battery door.

2. Press the pairing button on the TV Streamer 2 (2) once. The yellow light indicator on the front of the TV Streamer 2 will now blink once every 2 seconds to indicate that the TV Streamer 2 is ready to be paired in channel 1. Your TV Streamer 2 will remain in pairing mode for 20 seconds.
3. While pairing mode is activated, close the battery doors of both hearing instruments and make sure they are ON. A successful pairing will be indicated by an audible melody being played in both hearing instruments. You are now ready to use the TV Streamer 2.

Pairing the TV Streamer 2 in channel 2

1. Press the pairing button twice. The yellow light indicator on the front of the TV Streamer 2 will now blink twice every 2 seconds to indicate that it is ready to be paired in channel 2.
2. While in pairing mode close the battery doors on the hearing instruments.

Pairing the TV Streamer 2 in channel 3

1. Press the pairing button three times. The yellow light indicator on the front of the TV Streamer 2 will now blink three times every 2 seconds to indicate that it is ready to be paired in channel 3.
2. While in pairing mode close the battery doors on the hearing instruments.

Pairing blinking pattern

Pairing button pushes	Light indicator pattern	Channel
1		1
2		2
3		3


Start and stop streaming

Turn on the TV, sound system, computer or other audio device connected to TV Streamer 2. Make sure your hearing instruments are turned on as well. Streaming may be started in one of two ways.

1. Press and hold the push button on either hearing instrument for about three seconds. This must be done separately on both hearing instruments, in order to hear streamed audio in both ears, unless your hearing instruments have been set up with device-to-device synchronization. Once TV Streamer 2 has been connected, an audible melody will play in your hearing instruments and will be followed by the streamed sound.
2. If you have a Remote Control (optional), simply press the streaming button on this Remote Control. Once the hearing instruments connect to TV Streamer 2, an audible melody will play in both ears to confirm the presence of streamed audio and a streaming icon will appear in the display. The streamed sound can now be heard in your hearing instruments.

Start streaming on a 2nd or 3rd TV Streamer 2 unit

If more than one (up to three) TV Streamer 2 device is available, you may press the streaming button on the Remote Control a 2nd or 3rd time, depending on which TV Streamer 2 device you want to listen to.

Up to 3 TV Streamer 2 units are initially labeled by letter designator A, B or C, and the streaming device in use is shown on the Remote Control display.

If you do not have a Remote Control, press and hold the push button for about 3 seconds on the hearing instrument a 2nd or 3rd time to access the 2nd or 3rd streaming unit.

Stop streaming

Shortly press the push button on either hearing instrument. This will cancel streaming on the side pressed and return to the last program used.

Note: If your hearing instruments have device-to-device program synchronization, press the Program button only on one of the instruments. The second hearing instrument will switch to the same program automatically.

With the Remote Control (optional) you just have to press the "P" button located in the center of the key pad. This will return the hearing instruments to the last program used.

You may also press the home button to return to your default settings.

If you are using the telephone with Phone Clip+ or Auto Phone, streaming from the TV Streamer 2 will be temporarily disrupted. When you have ended your phone call, streaming will resume.

i *Note! The stable operating range of TV Streamer 2 is about 7 meters. You may get out of range from your TV Streamer 2. The connection with the audio source may be then lost. If you get back in range within 5 minutes, your hearing instruments will automatically reconnect and continue streaming.*

Controls

2" Mute or

The TV Streamer 2 is designed to provide a comfortable listening level from the moment it is activated. However, for some devices or listening environments it may be desirable to adjust the listening volume which can be accomplished in a variety of ways. Keep in mind that volume adjustments on the TV Streamer 2 apply to the streamed signal only, not the volume of normal environmental sounds. Make sure that your hearing instruments are comfortably adjusted before you change the settings with the volume control.

Use the "+" and "-" keys (7 and 8) on TV Streamer 2 to adjust volume to a comfortable level. In most cases this adjustment should only be made one time, as TV Streamer 2 will store the latest set volume level. Use ReSound Unite Remote Control 2 (optional) to adjust both the streaming volume level and the hearing instrument microphone volume level in the hearing instrument.* Ask your hearing care professional for more information about this accessory.

The volume can, in some cases, be further adjusted at the signal source, for example by turning the TV volume down. This will adjust the volume of the streamed sound transmitted

wirelessly to the hearing instruments. Adjusting the TV volume will affect the listening volume in the room as well.

Mono/stereo control

The two-position switch (3) may be used to set Stereo (two channels) or Mono (one channel) audio streaming. For binaural fittings, the switch should preferably be set to the Stereo position. For monaural fitting it should be set to the Mono position.

**Adjustment of both hearing instrument and streaming volume is only available with ReSound Unite Remote Control 2.*

Adjustable delay mode

In rare cases there may be a misalignment between the audio streamed to the hearing instruments and the sound coming from the TV speakers (echo) or even a misalignment between the streamed audio and the TV picture (lip-sync). This may be evident in more complex audio setups where the TV Streamer 2 is not connected directly to the TV (e.g. through a set-top box or other auxiliary unit). To overcome such misalignments the streaming delay from the TV Streamer 2 can be adjusted as follows:

1. Make sure that the TV Streamer 2 is turned on and that your hearing instruments are receiving the streamed audio.
2. Press and hold the pairing button (2) for 3 seconds. Now the TV Streamer 2 will enter into the adjustable delay mode, where the green light indicator on the front of the TV Streamer 2 will blink once every 2 seconds. The yellow light indicator will be on whenever an audio signal is present.
3. Use the "+" and "-" key (7 and 8) on the top of the TV Streamer 2 to increase or decrease the delay of the streamed signal in 10 ms steps. The streamed signal can be delayed in the range of 0-250 ms.
4. When the optimal delay has been reached, press and hold the pairing button for 3 seconds. Now the TV Streamer 2 will resume normal operation, indicated by the green light indicator being constantly on.

The offset set in adjustable delay mode will be stored in the TV Streamer 2.

If you do not exit the adjustable delay mode by pressing and holding the pairing button down for 3 seconds, the TV Streamer 2 will automatically resume normal operation after 30 minutes.

To reset the delay to 0 ms (factory setting), press and hold the pairing button while simultaneously powering on the TV Streamer 2.

Important Information

Maintenance and Care

- Clean the TV Streamer 2 using a damp cloth.
- Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.
- When the TV Streamer 2 is not being used for longer periods of time, unplug the power supply from the wall socket.

General warnings

- Keep the TV Streamer 2 out of reach of children under 3 years.
- ReSound Unite™ devices use low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected. In that case, move the TV Streamer 2 away from the affected electronic device.
- Use the TV Streamer 2 only with the power supply unit that has been supplied together with the product. The use of any other types may damage or destroy the product and could be dangerous. Use of other power supplies may invalidate the warranty.
- When you disconnect the power cord, grasp and pull the plug, not the cord.
- Never use a power supply that is damaged.
- Do not attempt to disassemble or make modifications to the power supply as it may expose you to dangerous electric shock
- Do not use the power supply outdoors or in damp areas.

General precautions

- Your hearing instrument and the TV Streamer 2 were given a unique communication network code during the fitting. This ensures that TV Streamer 2 will not affect hearing instruments worn by others.
- High-powered electronic equipment, larger electronic installations and metallic structures may significantly reduce the operating range.
- Use only original ReSound accessories.
- Do not make any changes or modifications to the TV Streamer 2.
- Opening the TV Streamer 2 might damage it. If problems occur which cannot be resolved using recommendations given in the Troubleshooting Guide, consult your hearing care professional.
- The TV Streamer 2 may only be repaired by an authorized service center.
- Do not use the TV Streamer 2 in areas where RF transmission is prohibited, e.g. airplanes, etc.
- Only connect the TV Streamer 2 to connections for which it is explicitly intended.

Troubleshooting Guide

SYMPTOM	CAUSE
<i>The streaming button on the Remote Control is pressed but with no effect at all</i>	The TV Streamer 2 and hearing instruments have not been paired. The Remote Control and hearing instruments are not within wireless range.
<i>The push button on the hearing instrument has been pressed for more than 3 seconds but with no effect at all</i>	The TV Streamer 2 and hearing instruments have not been paired.
<i>The pairing button on the TV Streamer 2 is pressed but no pairing melody is played in the hearing instrument</i>	The TV Streamer 2 and hearing instrument are not within wireless range. The TV Streamer 2 and hearing instrument have not been simultaneously in pairing mode.

POSSIBLE REMEDY

Carry out pairing process.

Ensure that Remote Control and hearing instruments are within wireless range and activate streaming again.

Carry out pairing process.

Ensure that the TV Streamer 2 and hearing instruments are within wireless range, and repeat pairing process.

Repeat pairing process and ensure that the battery doors on both hearing instruments are closed within 20 seconds after the pairing button has been pressed on the TV Streamer 2.

SYMPTOM	CAUSE
<i>The streaming button on the Remote Control is pressed but no streamed audio signal is found when searched for (indicated by an 'X' icon on the display)</i>	The TV Streamer 2 and hearing instrument are not within wireless range. The TV Streamer 2 is not cabled correctly to the TV. The TV Streamer 2 is powered off or the sound from the TV output has been muted.
<i>There is no sound in the hearing instruments although still being in the streaming program</i>	The TV Streamer 2 and hearing instrument are not within wireless range. The TV might have been powered off or the sound from the TV has been muted.
<i>The sound from the TV Streamer 2 is distorted</i>	The audio input level from the TV is too high.

POSSIBLE REMEDY

Ensure that Remote Control and hearing instruments are within wireless range and activate streaming again.

Connect the cable between the TV Streamer 2 and the TV correctly and activate streaming again.

Power on the TV Streamer 2 or unmute the TV and activate streaming again.

Ensure that the TV Streamer 2 and hearing instruments are within wireless range.

Power on or unmute the TV.

Adjust the volume using the volume button on the TV Streamer 2 until the sound is no longer distorted.

SYMPTOM	CAUSE
<i>The sound level from the TV Streamer 2 is very low</i>	The audio input level from the TV is too low.
<i>The sound from the TV Streamer 2 is distorted or drop-outs occur frequently</i>	The TV Streamer 2 and hearing instrument are on the edge of the wireless range. The TV Streamer 2 and hearing instruments are not within sufficient "line-of-sight".
<i>The sound from the TV Streamer 2 is not synchronized with the TV picture</i>	Your TV is not able to synchronize the sound from the selected audio outputs and the picture.

POSSIBLE REMEDY

Adjust the volume using the "+" and "-" keys on the TV Streamer 2 until the sound is

sufficiently loud.

Move a little closer to the TV Streamer 2.

Ensure that the TV Streamer 2 is placed in an adequate position, and that you are within normal reach without any significant obstacles obstructing the streaming.

If possible, try using another audio output from your TV. Alternatively, follow the instructions in the “Adjustable delay mode” section of this user guide. If this does not help, contact your TV dealer.

The TV Streamer 2 audio streaming has almost no latency and does not contribute itself to any lip sync error.

SYMPTOM	CAUSE
<i>The sound from the TV Streamer 2 is not synchronized with the sound from the TV loudspeakers</i>	Your TV is not able to synchronize the sound from the selected audio outputs with the sound from the TV loudspeakers.
<i>The sound from the hearing instruments is either too low or too high</i>	The audio volume level of the hearing instruments is not suitable for listening.
<i>A hearing instrument unintentionally drops out of the streaming program</i>	The TV Streamer 2 and hearing instruments have been out of wireless range for more than 5 minutes. The battery in the hearing instrument is so depleted that it no longer supports audio streaming.

POSSIBLE REMEDY

If possible, try using another audio output from your TV. Alternatively, follow the instructions in the “Adjustable delay mode” section of this user guide. If this does not help, contact your TV dealer.

The TV Streamer 2 audio streaming has almost no latency and does not contribute itself to any echo effects.

Adjust the volume using the “+” and “-” keys on the TV Streamer 2 until the sound is suitable..

Alternatively, use the “+” and “-” keys on the Remote Control (optional) for this operation.

Ensure that TV Streamer 2 and hearing instruments are within wireless range and activate streaming again.

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Technical specifications

Operating range:	up to 7 meters
Power supply:	external wall adapter
Wireless:	2.4 GHz
Operating temperature:	0 to 55°C
Storage temperature:	-20 to 60°C

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Date of sale _____ Seal

Manufacturer:
GN ReSound A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark

12-month warranty leaflet No. 2 for ReSound Unite TV Streamer 2

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12-month warranty No. 2 for ReSound Unite TV Streamer 2

Serial No. _____
Serial number _____
Date of audit _____ Stamp
Date of sale _____
Representative of repair contractor _____
_____ Seal

12-month warranty leaflet No. 2 for ReSound Unite TV Streamer 1

Serial number _____
Date of audit _____
Received _____
Representative of repair contractor _____
_____ Seal

12-month warranty No. 2 for ReSound Unite TV Streamer 1

Serial number _____
Date of audit _____ Stamp
Date of sale _____
Representative of repair contractor _____
_____ Seal

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I, Martynova Natalia Andreevna, deputy notary of the city of Moscow for Akimov Gleb Borisovich, certify the authenticity of signature, made by the translator Mamedov Timur Djavanshirovich in my presence. His identity is established.

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It is registered in the register under
Fee collected: 100 rubles

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Stitched up, numbered and sealed (1) pages

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
Welcome

Thank you for selecting ReSound.

ReSound Unite Mini Microphone is a high-tech and robust product that streams speech and audio signals directly to your hearing aids.

The Mini Microphone can also be connected to a computer, cell phone or other audio device and can transform your hearing aids into wireless headphones.

Please read this manual carefully to avoid any problems in the operation and fully benefit from the Mini Microphone.

Consult your hearing care professional if you have any questions.


GN ReSound A/S
Lautrupbjerg 7
2750 Ballerup - Denmark

Contents


Description.....	3
Getting Started	4
Charging the Battery	4
Turning your Mini Microphone ON and OFF	4
Pairing to your hearing instruments.....	4
Start streaming	6
Stop streaming	6
Wearing Your ReSound Unite™ Mini Microphone.....	7
Light indicator.....	7
Using the Line-in	7
Volume control	8
Important Information.....	8
Maintenance and Care.....	8
General warnings	8
General precautions	9
Built-in battery care	9
Charger care.....	10
Troubleshooting Guide	11
Technical specifications.....	13
Warranty and Repair.....	13
Acceptance and sale certificate	14

Description

1. Mini-USB connector for charging
2. Microphone inlet
3. Line-in input
4. ON/OFF switch

5. Volume control
6. Light indicator
7. Clip
8. Pairing button


Getting Started

Charging the Battery

Before using your ReSound Unite™ Mini Microphone, you need to fully charge the rechargeable battery. To charge it, plug the charging cable as shown in the picture to the right. Then plug the other end of the charging cable into a power outlet.

During the charging process, the LED indicator will show a solid red light. When fully charged and still connected to the charger, the LED will turn solid green interrupted by a red blink every 2 seconds.

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- When charging your ReSound Unite™ Mini Microphone for the first time, leave it to charge for at least 3 hours, even if the battery indicator indicates that the battery is fully charged.
- For safety reasons, recharge the battery only with the charger supplied by ReSound together with the Mini Microphone.

Do not charge the Unite Mini Microphone while it is on your body.

Turning your Mini Microphone ON and OFF

To turn ON your Mini Microphone, slide the ON/ OFF switch (4) down to the ON position. A white "I" icon will be visible on the switch and the LED indicator will now blink green every 2.5 seconds.

To turn it OFF, slide the ON/OFF switch up to the OFF position. The LED will now stop blinking.

If battery level is low the LED will blink yellow every 2 seconds.

Pairing to your hearing instruments

Note! The following information is intended primarily for your hearing care professional, who in most cases will complete the tasks that are necessary before the ReSound Unite™


Mini Microphone can be used with your hearing instruments. However, the process is designed to be simple and most people can do it themselves.

It is possible to pair up to three Unite™ streaming devices, e.g. one Mini Microphone and two TV streamers, with a given pair of hearing instruments. Therefore, channel 1, 2 or 3 must be selected when starting the pairing process.

Pairing the Mini Microphone in channel 1

1. Make sure ReSound Unite™ Mini Microphone is turned ON. Your hearing instruments must be OFF - open the battery door.
2. Remove the silver cap from the Mini Microphone. With the tip of a pen or similar object, press the small pairing button one time on the front of your ReSound Unite™ Mini Microphone (as shown in the picture). The LED will now blink yellow once every 2 seconds to indicate that is ready to be paired in channel 1.

Your ReSound Unite™ Mini Microphone will remain in pairing mode for 20 seconds.

3. While pairing mode is activated (20 sec), close the battery doors of both hearing instruments and make sure they are ON. A successful pairing will be indicated by an audible melody played in both hearing instruments. You are now ready to use the Mini Microphone!

Pairing the Mini Microphone in channel 2 or 3

To pair the Mini Microphone in channel 2, press the pairing button twice. The LED will now blink yellow twice every 2 seconds to indicate that it is ready to be paired in channel 2. Close the battery doors on the hearing instruments.

To pair the Mini Microphone in channel 3, press the pairing button three times. The LED will now blink yellow three times every 2 seconds to indicate that it is ready to be paired in channel 3. Close the battery doors on the hearing instruments.

Pairing blinking pattern:

Pairing button pushes	LED blinking pattern	Channel
1		1
2		2
3		3

Start streaming


Turn on the Mini Microphone and make sure your hearing instruments are turned on as well. Streaming may be started in one of two ways.

1. Press and hold the push button on both hearing instruments for about three seconds in order to listen to the streamed audio signal. If your hearing instruments have been set up with device-to-device synchronization, press the button only on one of the instruments. The second instrument will switch to the same program automatically. Once ReSound Unite™ Mini Microphone has been connected an audible melody will play in your hearing instruments and will be followed by the streamed sound.
2. If you have a ReSound Unite™ Remote Control, simply press the Unite™ streaming button. Doing so will briefly show a search icon in the remote control display. Once the hearing instruments connect to ReSound Unite™ Mini Microphone, a short melody will be played in both ears to confirm the presence of streamed audio and a streaming icon will appear in the display. The streamed sound can now be heard in your hearing instruments.

Start streaming on a 2nd or 3rd ReSound Unite™ Mini Microphone unit

If your hearing instrument is paired to more than one streaming device, you can start streaming in two ways.

1. Press and hold the hearing instrument push button for about 3 seconds a 2nd or 3rd time to access the 2nd or 3rd ReSound Unite™ streaming device.
2. If you use a Unite™ Remote Control, press the streaming button on the Unite™ Remote Control a 2nd or 3rd time to access the 2nd or 3rd ReSound Unite™ streaming device.

Searching for Mini
Microphone

Connected to Mini
Microphone

Streaming button

The up to 3 ReSound Unite™ streaming devices are initially labeled by letter designator A, B or C, and the device in use is shown on the Remote Control display.

Stop streaming

- By briefly pressing the Program button on the ReSound Unite™ Remote Control, you will return to the previous program of your hearing instrument.
- With the ReSound Unite™ Remote Control (optional) you just have to press the "P" button, located in the center of the key pad. This will return the hearing instruments to the last program used. You may also press the home button to return to your default settings.
- If you are using the telephone with ReSound Unite™ Phone Clip or PhoneNow™, streaming of ReSound Unite™ Mini Microphone will be temporarily disrupted. When you have ended your phone call, streaming will resume.

Note! You may get out of range (>7 m) from your ReSound Unite™ Mini Microphone. The connection will be then interrupted. If you get back in range within 5 minutes, your hearing instruments will automatically reconnect and continue streaming.

Wearing Your ReSound Unite™ Mini Microphone

Clip your ReSound Unite™ Mini Microphone on the speaker's jacket or other clothing. Place the device within a range of 10-20 cm from the speaker's mouth.

Light indicator


The Light indicator (6) provides information on the status of your device. The LED patterns and their meanings are given in the table below.

Blinks	Meaning
	Normal operation
	Battery low
	Charging
	Fully charged - still connected to charger
	Pairing in channel 1
	Pairing in channel 2
	Pairing in channel 3

Using the Line-in

The Line-in (3) can be used to stream audio signals from most auxiliary audio devices through the ReSound Unite™ Mini Microphone's line input directly to your hearing instrument.

1. Insert the included jack cable into the Line-in plug on the ReSound Unite™ Mini Microphone.
2. Plug the other end of the jack cable into the AUDIO OUT on the audio device you want to stream from.

i Note: When using the ReSound Unite™ Mini Microphone with a TV, the HEADPHONE output terminal is usually located on the side or front of the TV. When using the HEADPHONE output terminal on the TV, usually only the person wearing the hearing instruments will hear the sound, and the TV speakers are muted. Your TV may support speakers and HEADPHONE being on simultaneously. Check your TV manual for instructions.

Volume control

ReSound Unite™ Mini Microphone has been designed to provide a comfortable listening level from the moment it is activated. However, for some devices or listening environments it may be desirable to adjust the listening volume. This can be accomplished in a variety of ways. Keep in mind that volume adjustments apply to the streamed signal only, not the volume of normal environmental sounds.

- Use the “+” and “-” keys (5) to adjust volume to a comfortable level. Make sure that your hearing instruments are comfortably adjusted before you change the settings with the volume control.

Note: When using the Mini Microphone as a microphone the last volume setting will be stored. When using the Line-in, the volume setting will revert to the default level when the Mini Microphone is turned OFF and ON again to avoid an uncomfortable loud initial sound level when used with different audio sources.

- Use the ReSound Unite™ Remote Control to adjust the streaming volume in the hearing instrument. Ask your hearing care professional for more information about this product.
- The volume can in some cases be further adjusted at the signal source, for example by turning the TV itself, stereo system, computer, etc. down. In many cases this will adjust the volume of both streamed audio (transmitted wirelessly to your hearing instruments) and non-streamed audio (sounds reaching the hearing instrument microphones in the normal fashion). Of course, this will also affect listening volume for others in the room.

Important Information

Maintenance and Care

- Clean the ReSound Unite™ Mini Microphone using a damp cloth. Never use household cleaning products (washing powder, soap, etc.) or alcohol to clean the device.
- When ReSound Unite™ Mini Microphone is not in use, turn it OFF and store it safely.
- Protect ReSound Unite™ Mini Microphone from excessive moisture and heat. Protect it from excessive shock and vibration.
- Protect the microphone and charger inlets from debris and dirt. If necessary, use the cleaning brush provided with your hearing instruments to clean these areas.

General warnings

- Young children can wrap cords and lanyards / key chains around their neck causing the risk of strangulation. Therefore, always keep cords and lanyards / key chains out of reach of children.
- ReSound Unite™ devices use low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected.
- In that case, move the ReSound Unite™ Mini Microphone away from the affected electronic device.
- Avoid exposing your product to rain, moisture or other liquids to protect against damage to the product or injury to you.
- Keep the product and any cords and cables away from operating machinery.

- If the product overheats, if the product has been dropped or damaged, if the product has a damaged cord or plug, or if the product has been dropped in a liquid, discontinue use and contact your hearing care professional.
- Dispose of the product according to applicable regulations of the Russian Federation.
- Do not use the product in environments where there is a danger of ignition of flammable gasses.

General precautions

- Do not use this device in locations where it is forbidden to use electronic devices (e.g. airplane).
- X-ray radiation (e.g. CT scans, MRI scans) may adversely affect the correct functioning of this device. We recommend that you switch the device OFF before undergoing X-ray procedures and keep it outside the room.
- High-powered electronic equipment, larger electronic installations and metallic structures may impair and significantly reduce the operating range.
- Your hearing instruments and this device were given a unique communication code during the fitting. This ensures that the device will not affect hearing instruments worn by others.
- Do not make any changes or modifications to this device.
- Opening the ReSound Unite™ Mini Microphone might damage it. If problems occur which cannot be resolved using recommendations given in Troubleshooting Guide, consult your hearing care professional.
- ReSound Unite™ Mini Microphone may only be repaired by an authorized service center.
- Do only connect ReSound Unite™ Mini Microphone to connections which it is explicitly intended for.

Built-in battery care

- Your wireless device is powered by a rechargeable battery.
- The full performance of a new battery is achieved only after two or three complete charge and discharge cycles.
- The battery can be charged and discharged hundreds of times, but will eventually wear out.
- Unplug the charger from the electrical plug and the product when not in use, to save energy and minimize the risk of fire.
- If left unused, a fully charged battery will lose its charge over time.
- Leaving the product in hot places, such as in a closed car in the summer conditions, will cause the battery to heat, reduce its capacity and lifetime.
- A product with a cold battery may not work temporarily, even when the battery is fully charged. Battery performance is particularly limited in temperatures well below freezing.
- Caution! The battery used in this product may present a risk of fire or chemical burn if mistreated.
- Do not attempt to open the product or replace the battery. It is built-in and not changeable. Use of other batteries may present a risk of fire or explosion and injury.
- Recharge the ReSound Unite™ Mini Microphone only with a power supply unit that has been supplied together with the product.

- Dispose of the product with the built-in battery according to applicable regulations of the Russian Federation.

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- Do not attempt to charge your product with anything other than the ReSound charger provided. The use of any other charger may damage or destroy the device and the battery. Use of other chargers may invalidate the warranty.
- Charge the product according to the instructions supplied with the product.
- When you disconnect the power cord, grasp and pull the plug, not the cord.
- Never use a charger that is damaged.
- Do not attempt to disassemble the charger as it may be damaged and you may be exposed to electric shock.
- Avoid charging your product in extremely high or low temperatures.
- Do not charge the ReSound Unite™ Mini Microphone while it is on your body.

Troubleshooting Guide

<i>SYMPTOM</i>	<i>CAUSE</i>
<i>The streaming button on the Remote Control is pressed but with no effect at all</i>	The Mini Microphone and hearing instruments have not been paired. The Remote Control and hearing instruments are not within wireless range.
<i>The push button on the hearing instrument has been pressed for more than 3 seconds but with no effect at all</i>	The Mini Microphone and hearing instruments have not been paired.
<i>The pairing button on the Mini Microphone is pressed but no pairing melody is played in the hearing instrument</i>	The Mini Microphone and hearing instrument are not within wireless range. The Mini Microphone and hearing instrument have not been simultaneously in pairing mode.

POSSIBLE REMEDY

Carry out pairing process.
Assure that Remote Control and hearing instruments are within wireless range and activate streaming again.
Carry out pairing process.
Assure that Mini Microphone and hearing instruments are within wireless range and repeat pairing process.
repeat pairing process and assure that the battery doors on both hearing instruments are closed within 20 seconds after the pairing button has been pressed on the Mini Microphone.

Troubleshooting Guide

<i>SYMPTOM</i>	<i>CAUSE</i>
<i>The streaming button on the Remote Control is pressed but no streamed audio signal is found when searched for (indicated by an 'X' icon on the display)</i>	The Mini Microphone and hearing instrument are not within wireless range. The Unite Mini Microphone is not powered on. The battery level of the Mini Microphone is too low to support streaming. An audio device has been connected to the Line-in plug but the device is powered off or the output has been muted.

POSSIBLE REMEDY

Assure that the Mini Microphone and hearing instruments are within wireless range and activate streaming again.
Power on the unite Mini Microphone and activate streaming again.
Charge the battery.
Power on or unmute the external audio device.

Troubleshooting Guide

<i>SYMPTOM</i>	<i>CAUSE</i>
-----------------------	---------------------

There is no sound in the hearing instruments although still being in the streaming program	The Mini Microphone and hearing instrument are not within wireless range. The Mini Microphone is no longer powered on. The battery level of the Mini Microphone is too low to support streaming. An audio device has been connected to the Line-in plug but the device is powered off or the output has been muted.
--	---

POSSIBLE REMEDY

Assure that the Mini Microphone and hearing instruments are within wireless range again.
 Power on the unite Mini Microphone and activate streaming again.
 Charge the battery.
 Power on or unmute the external audio device.

Troubleshooting Guide

SYMPTOM

CAUSE

A hearing instrument unintentionally drops out of the streaming program	The Mini Microphone and hearing instrument have not been within wireless range for more than 5 minutes. The battery in the hearing instrument is so depleted that it no longer supports audio streaming.
The sound from the Mini Microphone is too high or is distorted	An audio device has been connected to the Line-in plug and the audio level from this device is too high.

POSSIBLE REMEDY

Take the necessary corrective action as described for the 'no sound' symptom above and activate streaming again.
 Replace the battery in the hearing instrument with a new one, and activate streaming again.
 Adjust the volume using the "+" and "-" keys on the Unite Mini Microphone, or using the volume control on the external audio device itself, until the best sound level and quality is obtained.

Troubleshooting Guide

SYMPTOM

CAUSE

The sound level from the Mini Microphone is very low	The Mini Microphone is not worn correctly by the speaker. The sound received by the Mini Microphone is very low. An audio device has been connected to the Line-in plug and the audio level from this device is too low.
--	---

POSSIBLE REMEDY

Clip the unite Mini Microphone on the speaker's jacket or other clothing or hang it around

the speaker's neck using the lanyard, so that the device is within 10-20 cm from the speaker's mouth

Adjust the volume using the "+" and "-" keys on the Mini Microphone.

Adjust the volume using the "+" and "-" keys on the Unite Mini Microphone, or using the volume control on the external audio device itself, until the best sound level and quality is obtained.

Troubleshooting Guide

SYMPTOM

The sound from the Mini Microphone is still either too low or too high

The sound from the Unite Mini Microphone is distorted or drop-outs occur frequently

CAUSE

The audio volume level of the hearing instruments is not suitable for listening.

The Mini Microphone and the hearing instruments are on the edge of the wireless range.

The Unite Mini Microphone and hearing instruments are not within sufficient "line-of-sight".

POSSIBLE REMEDY

Adjust the volume either on the hearing instruments themselves (if provided with a volume control), or alternatively (option), using the "+" and "-" keys on the Remote Control.

Assure that the Mini Microphone and hearing instruments are within wireless range.

Assure that no significant obstacles will obstruct the hassle-free streaming from the Mini Microphone to the hearing instruments.

Technical specifications

Operating range:	up to 7 meters or more
Power supply:	built-in battery
Wireless:	2.4 GHz
Operating temperature:	0 to 45°C
Storage temperature:	-20 to 45°C

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The product is covered by a free warranty for a period of one year from the day of purchase (the repair service will mark the type of performed *repair on the warranty card*):

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DK-2750 Ballerup

Denmark

12-month warranty leaflet No. 2 for ReSound Unite Mini Microphone

Serial number _____

Date of audit _____

Received _____

Representative of repair contractor
_____ Seal

12-month warranty card No. 2 for ReSound Unite Mini Microphone

Serial number _____

Date of audit _____ Stamp

Date of sale _____

Representative of repair contractor
_____ Seal

12-month warranty leaflet No. 1 for ReSound Unite Mini Microphone

Serial number _____

Date of audit _____

Received _____

Representative of repair contractor
_____ Seal

12-month warranty card No. 1 for ReSound Unite Mini Microphone

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Translation of this text was performed by me, the translator Mamedov Timur Djavanshirovich. The requirements of the text of the translation (maximum correctness, literate statement) have been explained to me.

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Stitched up, numbered and sealed (24) pages

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
USER GUIDE
ReSound Unite™ Phone Clip+

ReSound
rediscover hearing

Welcome

Thank you for selecting ReSound. Our products feature innovative technology, attractive design and reliable quality. ReSound Unite™ products are specifically designed to allow you to stay in touch with the outside world.

Phone Clip+ Bluetooth clip will allow you to connect your wireless ReSound hearing instruments to Bluetooth-enabled mobile phones. This enables clear and hassle-free phone connectivity, whether you are at home, in the office or on the move.

Phone Clip+ also doubles as a basic remote control making it possible for you to adjust to volume of your hearing instruments and toggle through programs, when you are not talking on the phone.

You can also chose to mute your hearing instrument microphone to shut out environmental sound while you are talking on the phone (Verso wireless family only).

Please read this manual carefully to avoid any damage or problems in the operation of the Phone Clip+.


GN ReSound A/S
Lautrupbjerg 7
2750 Ballerup - Denmark

Contents


Description	3
Charging the Battery	3
Turning your Phone Clip+ ON and OFF	4
Pairing to hearing instruments	4
What is Bluetooth?	5
Pairing to your mobile phone	6
Connecting	6
Phone operation	6
Remote Control operation	7
Volume Control	7
Mute or unmute the microphone(s) on you hearing instruments	8
Toggle through hearing instrument programs	8
Operating Time	8
Light indicators	8
Important Information	8
TROUBLESHOOTING GUIDE	11
Warranty and Repair	13
Acceptance and sale certificate	13

Description

1. ON/OFF switch
2. Call pick-up/hang-up button
3. Hearing instrument program toggle
4. Hearing instrument microphone mute button
5. Micro-USB for charging


6. Volume control
7. Microphone inlet
8. Mounting clip
9. Hearing instrument pairing button
10. Bluetooth pairing button

Getting Started

Charging the Battery

Before using your Phone Clip+, you need to fully charge the rechargeable battery. To charge it, plug the charging cable into the Phone Clip+ as shown in the picture to the right.

Then plug the other end of the charging cable into a power outlet. During the charging process, the LED indicator will show a solid red light.

When fully charged and still connected to the charger, the LED will turn solid green interrupted by a red blink every 2 seconds.

Charging an empty battery typically takes about 3 hours. It is safe to leave the charger connected overnight - the battery cannot be overcharged.

- When charging your Phone Clip+ for the first time, leave it to charge for at least 3 hours, even if the battery indicator indicates that the battery is fully charged.
- For safety reasons, recharge the battery only with chargers supplied by ReSound.

Do not charge the ReSound Unite Phone Clip+ while it is on your body.

Turning your Phone Clip+ ON and OFF

To turn ON your Phone Clip+ slide the ON/ OFF switch (1) down to the ON position. A white "I" icon will be visible on the base of the switch and the LED indicator will now blink green every 2.5 seconds.

To turn it OFF, slide the ON/OFF switch up to the OFF position.

The LED will now stop blinking. If battery level is low the LED will blink yellow every 2 seconds indicating that the unit should be recharged.

Pairing to hearing instruments

Note! The following information is intended primarily for your Hearing Care Professional, who in most cases will complete the tasks that are necessary before the Phone Clip+ can be used with your hearing instruments. However, the process is designed to be simple and should be easily accomplished by almost anyone.

1. Make sure Phone Clip+ is turned on. Your hearing instruments must be OFF. This is done by opening the battery door on the hearing instrument(s).
2. Remove the silver cap from the Phone Clip+. With the tip of a pen or similar object, press the small white pairing button one time on the back of your Phone Clip+ (see picture to the right). The LED will now blink yellow every 2 seconds to indicate that it is ready to be paired to your hearing instruments. Your Phone Clip+ will remain in pairing mode for 20 seconds.
3. While pairing mode is active, close the battery doors on both hearing instruments and make sure they are ON. A successful pairing will be indicated by an audible melody played in both hearing instruments. Also the LED on the Phone Clip will remain lit (yellow) for 2 seconds.

What is Bluetooth?

Bluetooth provides a way for devices like mobile phones, computers, etc. to communicate wirelessly. First, check whether the phone you wish to use offers Bluetooth connectivity. Look for the following symbol on the device or its user guide.

Second, the range of possible applications on a Bluetooth-enabled device is defined by “profiles”. The device you wish to use with your Phone Clip+ must support the appropriate Bluetooth profiles.

To receive or place phone calls your mobile phone must support either the “Handsfree” or the “Headset” profile. Many, but not all mobile phones support these Bluetooth profiles. Look for this information in your mobile phone user guide.

Bluetooth Operating Range

In free air, your Phone Clip+ supports Bluetooth operation up to a range of approximately 10 meters. Line-of-sight between the devices is not required. However, the following factors may affect the maximum range:

- interference in the environment
- the mobile phone which you are connecting to may support a more limited range than the 10 meters supported by the Phone Clip+

If your mobile phone moves out of range, the connection will be terminated. If the device is then moved again within the range, the Phone Clip+ will attempt to re-establish connection. Whether it achieves to do so will depend on delay since the link was broken and on the mobile phone properties.

Bluetooth Pairing

If your mobile phone is Bluetooth-enabled and supports the appropriate profiles, it can communicate with the Phone Clip+ only after completing two simple procedures: pairing and connecting.

The pairing procedure allows you to control which Bluetooth devices are permitted to communicate with each other. It needs to be completed only once for each device you wish to use with your Phone Clip+.

You can pair up to 8 Bluetooth devices with the Phone Clip+, but only 2 Bluetooth devices can be connected/active at a time. The pairing and connecting procedure is described in the following section.

Pairing to your mobile phone

The pairing procedure is controlled by your mobile phone. Different mobile phones have different menu structures. You should therefore refer to your mobile phone's manual if you cannot complete the pairing procedure with the generic steps described here.

1. Turn you Phone Clip+ ON and remove the silver cap. Place your mobile phone next to the Phone Clip+, and make sure your phone is ON.
2. Find the connectivity settings on your phone. Look for the "Bluetooth" feature in the menu.
3. On your mobile phone, ensure that Bluetooth is turned ON.
4. With the tip of a pen or similar object, press the small blue Bluetooth pairing button one time on the back of your Phone Clip+ (see picture below). Now the Phone Clip+ will be in Bluetooth pairing mode for two minutes.
5. On the mobile phone choose to search for Bluetooth devices.
6. Your mobile phone should present a list of Bluetooth devices discovered. Select "Hearing Aid Phone" from this list.
7. Your mobile phone may then ask you to enter a passcode. If so, enter "0000" (four zeros). Your Phone Clip+ should now be paired to your mobile phone. Some mobile phones may ask which Bluetooth service you wish to enable. Select "Headset".

Connecting

Once your mobile phone and Phone Clip+ are paired, the mobile phone may need to be connected before it will actually send the appropriate audio signals. This function should also appear in your mobile phone's "Bluetooth" menu, usually under a list of "Paired Devices". In that list, highlight "ReSound Unite™" and select "Connect". When your Phone Clip+ and mobile phone are connected, you may see a headset symbol on the main screen of your mobile phone, similar to this: 

The connection between the Phone Clip+ and the mobile phone will be maintained for as long as both devices are turned ON and remain within range of each other. If either device is switched OFF or gets out of range, you may have to perform the connection procedure again.

On some phones you can configure "Automatic Reconnection". In this case, you will not need to connect your devices again. Refer to the user guide of your mobile phone for instructions on how to do this if supported.

Phone clip operation

In the below table you will find an overview of how you operate the Phone Clip+ when using it in connection with phone calls.

Answering and ending a call	Push the call pick-up/hang-up button once
Rejecting a call	Double click the call pick-up/hang-up button

Transferring a call initiated on the mobile phone to Phone Clip+	Push the call pick-up/hang-up button once
Put an active call on hold	Press and hold the call pick-up/hang-up button for 2 seconds
Last number redial	Double click the call pick-up/hang-up button while no call is active nor pending

Note:

Except for the answer and end call operations one or more of the below described operations may be phone dependent, i.e. they may either differ slightly or be absent altogether.

If another call is already active it will be terminated and the new incoming call - indicated by knocking tones - will become active. If another call is currently on hold it will become the active call instead.

If another call is already active it will remain active, only the new incoming call - indicated by knocking tones - will be rejected

If another call is currently on hold or an incoming call is pending - this will now become the active call instead.

If two mobile phones are connected this operation applies to the last paired mobile phone.

Voice dialing	Press and hold the call pick-up/hang-up button for 2 seconds.
Increase or decrease the volume of the speaker at the other end of the phone call	Push the + sign on the side of the Phone Clip+ to increase the volume. Push the - sign on the side of the Phone Clip+ to decrease the volume.
Mute or unmute your hearing instrument microphones during a phone call	Push the microphone mute/unmute button  on the front of the Phone Clip+

If two mobile phones are connected this operation applies to the last paired mobile phone. This functionality requires that the Voice dialing functionality is available and activated on your Bluetooth phone.

This functionality is only available for the wireless Verso hearing instrument family.

Remote Control operation

When not on the phone, the Phone Clip+ can be used as a remote control for your hearing instruments*.

Volume Control

Push the + sign on the volume control (6) to increase the volume of your hearing instruments. Push the - sign on the volume control (6) to decrease the volume of your hearing instruments. Depending on the phone, you can also utilize the volume control that is built in your mobile phone.

Please note that if you have two hearing instruments changing the volume will apply to both hearing instruments in unison.

*For the family of wireless Alera hearing instruments only: The Alera hearing instrument(s) must be initially paired once with a ReSound Unite Remote Control (either a customer or dispenser remote control) to enable the remote control functionality in the Phone Clip+.

Mute or unmute the microphone(s) on you hearing instruments)

(Only available for the wireless Resound Verso hearing instrument family)

To mute/unmute your hearing instrument microphones push the microphone mute/unmute button (4) once on the front of the Phone Clip+.

Toggle through hearing instrument programs

Push the program button (3) to toggle through available programs in your hearing instrument(s).

Operating Time

Battery consumption and the operating time of your Phone Clip+ depend greatly on the usage of the device. Your Phone Clip+ battery can be recharged hundreds of times. If the operating time reduces substantially, please contact your hearing care professional.

Wearing Your Phone Clip+

Clip your Phone Clip+ on your sweater, jacket or other clothing so that the microphone inlet faces upwards. For optimal voice pickup keep the device within a range of 10- 30 cm from your mouth. Your Phone Clip+ can be also attached to a car's sun visor.

Light indicators

Light indicator at the top of your Phone Clip+ serve as a multipurpose user interface, providing information on the status of your device.

Blinks	Meaning
	Normal operation
	Battery low
	Charging
	Fully charged, still connected to charger
	Ready for pairing with hearing instrument(s)
	Ready for pairing with Bluetooth device
	Bluetooth call active

Important Information

Maintenance and Care

- Clean the ReSound Unite Phone Clip+ using a soft, dry cloth. Never use household cleaning products (washing powder, soap, gel, etc...) or alcohol to clean the device.
- Protect ReSound Unite Phone Clip+ from excessive moisture and heat. Protect it from excessive shock and vibration.
- Protect the microphone and charger inlets from debris, dirt and moisture during cleaning and operation.
- When ReSound Unite Phone Clip+ is not in use, turn it OFF and store it safely

General precautions

- To prevent damage or injury, keep this device out of reach of children or people with mental disorders.
- ReSound Unite Phone Clip+ uses low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected. In that case, move the ReSound Unite Phone Clip+ away from the affected electronic device.
- When using ReSound Unite Phone Clip+ and the device is affected by electromagnetic interference, move away from the source.
- Avoid exposing your product to rain, moisture or other liquids to protect against damage to the product or injury to you.
- Keep the product and any cords and cables away from operating machinery.
- If the product overheats, if the product has been dropped or damaged, if the product has a damaged cord or plug, or if the product has been dropped in a liquid, discontinue use and contact your hearing care professional.
- Dispose of the product according to applicable regulations of the Russian Federation.
- Do not use in areas where there are explosive hazards or danger of ignition of flammable gasses.

General warnings

- Do not use this device in locations where it is forbidden to use electronic devices and infrared transmission, e.g. airplanes.
- High-powered electronic equipment, larger electronic installations and metallic structures may impair and significantly reduce the operating range.
- If the hearing instruments do not respond to the device because of an unusual field disturbance, move away from the disturbing field.
- Your hearing instruments and this device were given a unique communication code during the fitting. This ensures that the device will not affect hearing instruments worn by others.
- Do not open or attempt to repair or make any modifications to this device. This may result in damage and your warranty will be terminated
- If any problems occur, use the troubleshooting table at the end of this manual. If problems occur which cannot be resolved, consult your hearing care professional. ReSound Unite Phone Clip+ may only be repaired by an authorized service center.

Built-in battery care

- The ReSound Unite Phone Clip+ is powered by a rechargeable battery.
- The full performance of a new battery is achieved only after two or three complete charge and discharge cycles.
- The battery can be charged and discharged hundreds of times, but will eventually wear out.
- Unplug the charger from the electrical plug and the product when not in use. to save energy and minimize the risk of fire.
- If left unused, a fully charged battery will lose its charge over time.
- Leaving the product in hot places, such as in a closed car in the summer conditions, will cause the battery to heat, reduce its capacity and lifetime.
- A product with a cold battery may not work temporarily, even when the battery is fully charged. Battery performance is particularly limited in temperatures well below freezing.
- Do not attempt to open the product and replace the battery. It is built-in and not changeable. Use of other batteries may present a risk of fire or explosion and injury.

- Recharge the ReSound Unite Phone Clip+ only with a power supply unit that has been supplied together with the product.
- The battery may explode if damaged.
- Do not dispose the ReSound Unite Phone Clip+ as household waste as it may explode or result in a fire.
- Dispose of the ReSound Unite Phone Clip + as a device with a built-in battery according to the applicable regulations of the Russian Federation

Charger care

Make sure to read this carefully before use.

- Do not attempt to charge your product with anything other than the ReSound charger provided together with your ReSound Unite Phone Clip+. The use of any other charger may damage or destroy the device and the battery. Use of other chargers may invalidate the warranty.
- Charge the product according to the instructions supplied with the product.
- When you disconnect the power cord, grasp and pull the plug, not the cord.
- Never use a charger that is damaged.
- Do not attempt to disassemble the charger as it may be damaged and you may be exposed to electric shock.
- Avoid charging your product in extremely high or low temperatures.
- Do not charge the ReSound Unite Phone Clip+ while it is on your body.

TROUBLESHOOTING GUIDE

SYMPTOM	CAUSE
<i>The Phone Clip+ does not respond to any (long or short) presses on the call pick-up/hang-up button</i>	The battery in the Phone Clip+ is depleted
<i>The pairing button on the Phone Clip+ is pressed but no pairing melody is played in the hearing instrument after the battery doors have been closed</i>	The Phone Clip+ and hearing instrument are not within wireless range The Phone Clip+ and hearing instrument have not been simultaneously in pairing mode
<i>The Phone Clip+ cannot be found by the mobile phone when searching for it in the Bluetooth menu</i>	The Phone Clip+ and hearing instrument are not within wireless range The Phone Clip+ and mobile phone are not within wireless range

POSSIBLE REMEDY

Charge the battery.

Assure that Phone Clip+ and hearing instruments are within wireless range, and repeat pairing process

Repeat pairing process, and assure that the battery doors on both hearing instruments are closed within 20 seconds after the pairing button has been pressed on the Phone Clip+

Assure that Phone Clip+ and mobile phone are within wireless range, and repeat Bluetooth pairing process

The Phone Clip+ is only in Bluetooth pairing mode for 120 seconds after having been activated. Repeat the Bluetooth pairing process.

TROUBLESHOOTING GUIDE

SYMPTOM	CAUSE
<i>The Bluetooth pairing failed</i>	The Phone Clip+ left Bluetooth pairing mode before mobile phone completed the pairing session. The wrong 4-digit PIN number has been entered.
<i>The mobile phone is ringing but the light indicator on the Phone Clip+ does not flash, and no ringing tones will sound in the hearing instruments</i>	The Phone Clip+ is powered off. The Phone Clip+ and mobile phone have not been paired. The Bluetooth connection between the Phone Clip+ and the mobile phone is lost.

POSSIBLE REMEDY

Repeat the Bluetooth pairing process, and assure that the process is completed within 120 seconds after having been activated in the Phone Clip+.

Repeat the Bluetooth pairing process and if a PIN number is requested, enter 0000 (four zeros).

Power on the Phone Clip+.

Perform the Bluetooth pairing process and make a Bluetooth connection between the Phone Clip+ and the mobile phone.

Assure that Bluetooth is enabled in the mobile phone and make a Bluetooth connection between the Phone Clip+ and the mobile phone again, e.g. by pressing the multifunction button briefly on the Phone Clip+.

As for any details of the Bluetooth operation of the mobile phone, please refer to its user guide.

TROUBLESHOOTING GUIDE

SYMPTOM	CAUSE
<i>The mobile phone is ringing and the light indicator on the Phone Clip+ does flash, but no ringing tones will sound in one or both hearing instruments</i>	The Phone Clip+ and hearing instrument are not within wireless range. The battery in the hearing instruments is so depleted that it no longer supports audio streaming.
<i>The remote user does not hear my voice clearly</i>	The Phone Clip+ microphone does not pick up your voice correctly.
<i>The sound from the hearing instruments is either to low or too high</i>	The audio input level is not suitable for listening.
<i>There is no sound in the hearing instruments although being in the phone program</i>	The Phone Clip+ and hearing instrument are not within wireless range.

POSSIBLE REMEDY

Assure that Phone Clip+ and hearing instruments are within wireless range.

Replace the battery in the hearing instrument with a new one.

Assure that the Phone Clip+ microphone is kept close to your mouth (at a distance of 10-30 cm), and that the Phone Clip+ does not rub against any clothing during the call.

Adjust the volume on the mobile phone until the sound is suitable, or adjust the volume using the "+" and "-" keys on the Phone Clip+ until the sound is suitable. Alternatively (option), use the "+" and "-" keys on the Remote Control for doing this operation.

Assure that Phone Clip+ and hearing instruments are within wireless range.

Technical specifications

Operating range:	up to 3 meters
Power supply:	built-in battery
Wireless:	2.4 GHz
Operating temperature:	0 to 45°C
Storage temperature:	-20 to 45°C

Temperature test, transport and storage information

The product is subject to various tests in temperature and damp heating cycling between -25°C and +70°C according to internal standards of the manufacturer.

Warranty and Repair

The warranty maintenance of ReSound Unite™ devices is provided by

The product is covered by a free warranty for a period of one year from the day of purchase (the repair service will mark the type of performed repair on the warranty card):

- provided that the date of purchase is specified in this user guide and warranty card and that they bear the seal of the trading company;
- provided that the device is clean. If the date of purchase is not specified, the warranty shall come into force from the date of audit.

The following ReSound Unite™ devices are not covered by warranty:

- devices with any signs of mechanical damage;
- devices with any signs of chemical exposure;
- devices subjected to unauthorized disassembly and misuse in violation of operation conditions specified in this user guide.

In the above cases, the repair shall be performed at the buyer's expense.

Acceptance and sale certificate

ReSound Unite™ products are registered in the Russian Federation.

Registration Certificate No. RZN 2015/2727 dated June 4, 2015.

Declaration of Conformity ROSS DK.IM28.D02086 dated June 25, 2015.

ReSound Unite™ Phone Clip+

Serial number _____

Date of audit _____ Stamp

Date of sale _____ Seal

Manufacturer:

GN ReSound A/S

Lautrupbjerg 7

DK-2750 Ballerup

Denmark

12-month warranty leaflet No. 1 for the ReSound Unite Phone Clip+

Serial number _____

Date of audit _____

Received _____

Representative of repair contractor

_____ Seal

12-month warranty No. 1 for the ReSound Unite Phone Clip+

Serial number _____
Date of audit _____ Stamp
Date of sale _____
Representative of repair contractor _____ Seal

Please consult your local hearing care professional concerning disposal of your accessory.

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Stitched up, numbered and sealed 29 pages

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

APPENDIX TO INSTRUCTION MANUAL

Digital Hearing Instrument LINX2 with Accessories

APPROVE

GN ReSound A/S, Denmark

(position)
MARIA I. JENSEN
(name)

(signature)

" " _____ 2016
"day" month (numerals)

Stamp

GN ReSound A/S
Lautrupbjerg 7
2750 Ballerup - Denmark

The Danish Chamber of Commerce do hereby
confirm that the company is a member of our organization
and know to us as worthy of confidence.

09 AUG 2016

Danish Chamber of Commerce
Secretary: Pernille Hellesøe

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Iram Ilyas

20. Earwax filter (8 pcs).
21. Tube cleaning wire
22. Meter
23. Box for thin tubes
24. Set of thin tubes in a box
25. Thin tubes of various sizes for the hearing instruments, at most 100 pcs.
26. Tulip mold, at most 50 pcs.
27. Thick mold, small, at most 50 pcs.
28. Thick mold, medium, at most 50 pcs.
29. Thick mold, large, at most 50 pcs.
30. Open mold, small, at most 50 pcs.
31. Open mold, medium, at most 50 pcs.
32. Open mold, large, at most 50 pcs.
33. Hearing instrument packaging:
 - Carton box
 - Plastic case
 - Protective cover
 - Polishing cloth
 - Cleaning brush
34. Remote control device RC-2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
35. Streamer TV2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
36. Midget microphone RS UNITE, MINI MIC (EU, US, AUS, UK, CAN) c/w packaging
37. Hands-free phone kit PHONE CLIP PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
38. Instructions for the behind-the-ear hearing aids
39. Set of batteries, size 13

Name of medical device: Digital Hearing Instrument LINX2 with Accessories

Description of medical device:

Digital hearing instrument LINX2 with accessories

Digital hearing instrument LINX2, models: HI,LS977-DW, HI,LS988-DW, HI,LS777-DW, HI,LS788-DW, HI,LS577-DW, HI,LS588-DW

Accessories:

1. Color elements for housing
2. Volume control blank plug
3. Tool for battery compartment lock
4. Child earhook
5. Adolescent earhook
6. Red color marker
7. Blue color marker
8. Hearing instrument connecting adapter
9. Extension cable for connection with hearing instrument
10. 500-mm connecting cable for devices, left
11. 500-mm connecting cable for devices, right
12. 2500-mm connecting cable for devices, left
13. 2500-mm connecting cable for devices, right
14. Magnet "Autophone", at most 10 pcs.
15. Remover of battery compartment, at most 10 pcs.
16. Integrated FM adaptor, NA
17. Integrated FM adaptor, NB
18. Audio boot
19. Thin tube adaptor
20. Earwax filter (8 pcs).
21. Tube cleaning wire
22. Meter
23. Box for thin tubes
24. Set of thin tubes in a box
25. Thin tubes of various sizes for the hearing instruments, at most 100 pcs.
26. Tulip mold, at most 50 pcs.
27. Thick mold, small, at most 50 pcs.
28. Thick mold, medium, at most 50 pcs.
29. Thick mold, large, at most 50 pcs.
30. Open mold, small, at most 50 pcs.
31. Open mold, medium, at most 50 pcs.
32. Open mold, large, at most 50 pcs.
33. Hearing instrument packaging:
 - Carton box
 - Plastic case
 - Protective cover
 - Polishing cloth
 - Cleaning brush
34. Remote control device RC-2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
35. Streamer TV2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
36. Midget microphone RS UNITE, MINI MIC (EU, US, AUS, UK, CAN) c/w packaging
37. Hands-free phone kit PHONE CLIP PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging
38. Instructions for the behind-the-ear hearing aids
39. Set of batteries, size 13

Manufacturing

Primary manufacturer

GN ReSound A/S, Lautrupbjerg 7, DK-2750 Ballerup, Denmark

Manufacturing sites

GN ReSound A/S, Lautrupbjerg 7, DK-2750 Ballerup, Denmark

GN ReSound, Component Mfg. Danmark, Industrivej 6, 4720 Praesto, Denmark

GN ReSound China Ltd., 15, Chuang Xin Road, CTP Area, Xiamen 361006, P.R. of China

GN ReSound A/S digital hearing instrument LINX2 with accessories must comply with the in-house technical standard used by GN ReSound A/S, Denmark on the development and manufacture of digital hearing instruments LINX2 with accessories (this standard is confidential proprietary information because it contains know-how technologies).

Development rights on this medical device are owned by GN ReSound A/S, Denmark.

The final inspection of the device is performed by the quality control department in accordance with the requirements approved by GN ReSound A/S.

Contact with human body

Items	Materials contacting with human body
Color elements for housing	Terluran GP-22 ABS plastic Colorants: Pyrazolone Orange (3520-72-7) – 90%; Barium sulfate BaSO ₄ (7727-43-7) – 10%
Streamer TV2, PLUS, UNITE, RS	Cycoloy C1200HF-701 ABS plastic, black (top part); Cycoloy C1200HF ABS plastic, gray (bottom part); PEBAX 6333 SP 01 thermoplastic elastomer (buttons)
Midget microphone RS UNITE, MINI MIC	Cycoloy C1200HF-701 ABS plastic, black (top part); Cycoloy C1200HF ABS plastic, gray (bottom part); PEBAX 6333 SP 01 thermoplastic elastomer (buttons)
Hands-free phone kit PHONE CLIP PLUS, UNITE, RS	Cycoloy C1200HF-701 ABS plastic, black (top part); Cycoloy C1200HF ABS plastic, gray (bottom part); PEBAX 6333 SP 01 thermoplastic elastomer (buttons)
Earwax filter	Eltex Med PH23T630 low-density polyethylene (LDPE)

Minimal package contents:

- digital hearing instrument LINX2,
- instruction for use
- plastic case.

Appendix to the Technical Description

HI,LS977-DW, HI,LS777-DW, HI,LS577-DW

Hearing compensation: grade II-III hearing loss

The **Smart Start** feature delays turning on the hearing instrument for 5 or 10 seconds after closing the battery door.

The **PhoneNow** feature automatically enables the phone conversation program when a phone with an attached PhoneNow magnet is brought to the ear.

The **Comfort Phone** feature is used with binaural hearing aids. When a phone with an attached PhoneNow magnet is brought to one ear, the instrument in this ear automatically activates the phone conversation program, and the opposite instrument automatically reduces the gain by 10 dB to prevent the surrounding sounds from interfering with the conversation.

The **Softswitching** feature automatically switches the instrument from the omnidirectional to the directional mode and vice versa depending on the presence of speech and noise in the environment.

The **Environmental Optimizer** feature classifies acoustic situations into seven categories (silence, low speech, loud speech, speech in a moderate noise, speech in a loud noise, low noise and loud noise) and sets the amplification parameters depending on the situation. When a user moves from one situation to another the instrument automatically adjusts the amplification parameters.

Noise Tracker II is a multilevel noise suppression system that detects the noise and subtracts it from the overall signal.

Expansion is a method of suppression of low level noises, like noises from refrigerator, computer fans etc.

Windguard is a method of suppression of wind noise.

Sound Shaper is a method of compression of high frequency sounds into the medium frequency range. It is used when the patient does not respond to high frequencies.

DFS Ultra II with music mode is a technology that suppresses the feedback (whistling) and provides a special music mode that does not suppress high tones of a music composition.

The **Auto DFS** feature is enabled if the feedback cannot be calibrated. During the use of the hearing instrument, this feature automatically calculates critical frequencies and gains that cause the feedback.

Acceptance Manager is a technology that automatically gradually increases the gain of the hearing instrument over the given period of time. It is used to remotely bring the client from inexperienced to the experienced category.

Tinnitus Sound Generator is a method of tinnitus therapy (ringing in the ears). This module emits various sounds that distract the user from the tinnitus sensation. These can be various relaxing sounds from the nature or a constant or modulated white noise.

Low-frequency gain is a technology of additional amplification of bass frequencies.

The **Onboard Analyzer II** feature records statistical data on the use of the hearing instrument by the patient, including the average wear time during the day, the intensity of volume control and program switching, the duration of various acoustic situations etc.

In Situ Audiometry makes it possible to perform an audiometry through the hearing instrument. The tones are presented via the hearing instrument and the patient's response is recorded. This feature can be used if no audiometer is available.

Wireless 1
and wires.

* All type
omnidirec
mode. The
the fitting

Adaptive
depending

Multiada
omnidirec
noise and

Automat
automatic

Binaural
the binau

Natural
omnidire
mode. Th

Binaura
on the ac
omnidire
mode wh
Spatial S
bilateral
sound ty

Model d
Battery :
Function
Fully fle
Synchro
Synchro
SmartSt
PhoneN
Comfort
2.4 GHz
Direct il
ReSoun
ReSoun
required
ReSoun
Audio
WARP
Environ

Wireless fitting with Airlink 2 can be used to fit a hearing instrument without additional boots and wires. To do so, the Airlink 2 programmer unit is connected to a computer of the specialist.

* All types of directionality represent a directional mix: low frequencies are processed in the omnidirectional mode, while medium and high frequencies are processed in the directional mode. The threshold frequency below which the omnidirectional mode is used is calculated by the fitting program individually for each patient.

Adaptive Directionality is an automatic change of the directional diagram in the back depending on the intensity and location of the noise source.

Multiadaptive Directionality is an automatic change of the directional diagram from the omnidirectional pattern to the directional pattern depending on the intensity and location of the noise and speech sources.

Automatic Adaptive Directionality is a variant of the multiadaptive directionality that automatically changes the beamwidth depending on the intensity of the signal in the front.

Binaural Directionality, Natural Directionality, and Spatial Sense features are available only in the binaural configuration (use of two instruments).

Natural Directionality is a configuration when one instrument (in one ear) is set in the omnidirectional mode and the other instrument (in the opposite ear) is set in the directional mode. This allows the patient to perceive all sounds and focus on interesting information.

Binaural Directionality is an automatic change of directionality in both instruments depending on the acoustic situation. The following options are available: both ears are in the omnidirectional mode; both ears are in the directional mode; the left ear is in the directional mode while the right ear is in the omnidirectional mode and vice versa.

Spatial Sense is the use of the omnidirectional mode in both ears with the Auricle Imitation and bilateral compression. This technology allows the patient to get a spatial sense, to distinguish sound types and to locate sound sources.

Model overview	HI,LS977-DW	HI,LS777-DW	HI,LS577-DW
Model description			
Battery size	13		
Functional technologies			
Fully flexible programs	4	4	4
Synchronized program push button	•	•	•
Synchronized volume control	•	•	•
SmartStart	•	•	•
PhoneNow	•	•	•
Comfort Phone	•	•	•
2.4 GHz binaural link	•	⊙	○
Direct iPhone connectivity	•	•	•
ReSound Unite wireless accessories	•	•	•
ReSound Control application (Phone Clip+ required)	•	•	•
ReSound Smart application	•	•	•
Audiological features			
WARP signal processing channels	17	14	12
Environment Classifier	•	•	•

Binaural Directionality II	•		
Spatial Sense	•		
Binaural Directionality		•	
Mixed Directionality	•	•	•
- Adjustable directional mix	•		
Natural Directionality II		•	•
Synchronized SoftSwitching	•	•	
SoftSwitching			•
Automatic Adaptive Directionality	•		
Multiadaptive Directionality		•	
Adaptive Directionality			•
Binaural Environmental Optimizer II	•		
Environmental Optimizer		•	
NoiseTracker II	•	⊙	○
Expansion	•	⊙	○
Windguard	•	⊙	○
Sound Shaper	•	•	•
DFS Ultra II	•	•	•
- Music Mode	•	•	•
Auto DFS	•	•	•
Synchronized Acceptance Manager	•	•	•
Tinnitus Sound Generator	•	•	•
Fitting features			
Aventa 3.8 software or higher	•	•	•
Onboard Analyzer II	•	•	•
In Situ Audiometry	•	•	•
AIRLINK2* wireless fitting	•	•	•


• - Main ⊙ – Advanced ○ – Basic

* Programmer (programming interface is AIRLINK2) – Registration Certificate No. P3H 2015/2727 dated June 4, 2015.

Appendix to the Specification

Parameter	
Compression ratio	1:0 3:0
AGC attack time (ms)	10
AGC recovery time (ms)	50
Electrical input sensitivity by 70 dB at 1600 GHz (mV)	2.6

HI,LS988-DW, HI,LS788-DW, HI,LS588-DW

Hearing compensation: grade III-IV hearing loss

Model overview	HI,LS988-DW	HI,LS788-DW	HI,LS588-DW
Model description			
Battery size	13		
Functional technologies			
Fully flexible programs	4	4	4
Synchronized program push button	•	•	•
Synchronized volume control	•	•	•
SmartStart	•	•	•
PhoneNow	•	•	•
Comfort Phone	•	•	•
2.4 GHz binaural link	•	⊙	○
Direct iPhone connectivity	•	•	•
ReSound Unite wireless accessories	•	•	•
ReSound Control application (Phone Clip+ required)	•	•	•
ReSound Smart application	•	•	•
Audiological features			
WARP signal processing channels	17	14	12
Environment Classifier	•	•	•
Binaural Directionality II	•		
Spatial Sense	•		
Binaural Directionality		•	
Mixed Directionality	•	•	•
- Adjustable directional mix	•		
Natural Directionality II		•	•
Synchronized SoftSwitching	•	•	
SoftSwitching			•
Automatic Adaptive Directionality	•		
Multiadaptive Directionality		•	
Adaptive Directionality			•
Binaural Environmental Optimizer II	•		
Environmental Optimizer		•	
NoiseTracker II	•	⊙	○
Expansion	•	⊙	○
Windguard	•	⊙	○
Sound Shaper	•	•	•
DFS Ultra II	•	•	•

- Music Mode	•	•	•
Auto DFS	•	•	•
Synchronized Acceptance Manager	•	•	•
Tinnitus Sound Generator	•	•	•
Gain strategy (WDRC/Semi-linear/Linear)	•	•	⊙
Low-frequency gain	•	⊙	○
Fitting features			
Aventa 3.8 software or higher	•	•	•
Onboard Analyzer II	•	•	•
In Situ Audiometry	•	•	•
AIRLINK2* wireless fitting	•	•	•

• Main ⊙ Advanced ○ Basic

* Programmer (programming interface is AIRLINK2) – Registration Certificate No. P3H 2015/2727 dated June 4, 2015.

Appendix to the Specification

Parameter	
Compression ratio	0:1 3:0
AGC attack and recovery time (ms)	10
AGC recovery time (ms)	50
Electrical input sensitivity by 70 dB at 1600 GHz (mV)	1.8

Specification of the magnet that automatically switches the hearing instrument to the phone conversation mode.

Magnet power that automatically enables the phone conversation mode – 1370–1436 mT







Weight of the magnet – not more than 2.8 g





Dimensions – diameter 14 mm, thickness 1.1 mm

Directionality Specification







HI,LS977-DW, HI,LS777-DW, HI,LS577-DW





1. Adaptive Directionality



2. Fixed Directionality









HI,LS988-DW, HI,LS788-DW, HI,LS588-DW





1. Adaptive Directionality










2. Fixed Directionality









ACCESSORIES:

Name	Photo	Weight (g)**	Size (mm)
1. Color elements for housing Intended use: These elements change the color of the instrument housing.		Weight of the set – 14.7 Weight of the battery door cover – 0.3 Weight of the housing side cover – 0.2 Weight of the housing top cover – 0.2	Size of the 110×53 Size of the door cover 18×9 Size of the side cover 33×17 Size of the top cover 33×12
2. Volume control blank plug Intended use: This plug fits in the slot that appears after removing the volume control from the housing.		0.1	8.0×2.0
3. Tool for battery compartment lock Intended use: This tool opens the battery door.		1.9	5.6×0.9
4. Child earhook 5. Adolescent earhook Intended use: This part conducts the sound from the hearing instrument to the earmold and into the ear.		0.20 0.22	17.0×16.0 19.0×16.0
6. Red color marker 7. Blue color marker Intended use: Red and blue plastic inserts are placed onto the hearing instruments to mark the side (blue for the left ear and red for the right ear).		0.09	32.0×5.0
8. Hearing instrument connecting adapter Intended use: This adaptor connects the hearing instrument to a cable and to a computer.		0.02	23.0×4.0

14. Magnet "Autophone", at most 10 pcs. Intended use: This magnet automatically switches the hearing instrument to the phone conversation mode.		2.5	12.0x5.0
15. Remover of battery compartment, at most 10 pcs. Intended use: This tool is used to remove the battery compartment from the hearing instrument housing.		5.0	50.0x10.0
16. Integrated FM adaptor, NA 17. Integrated FM adaptor, NB Intended use: This adaptor is used to connect the hearing instrument to FM systems.		1.10	19.0x12.0
18. Audio boot Intended use: The audio boot is used to connect the hearing instrument to audio systems.		0.03	22.0x9.0x10.0
19. Thin tube adaptor Intended use: This adaptor is used to connect a thin tube to the hearing instrument.		0.15	9.0x7.0x10.0
20. Earwax filter (8 pcs) Intended use: This filter protects the hearing instrument from earwax ingress.		1 piece – 0.3 The set – 5.1	41x4.0x10.0 50.0x35.0
21. Tube cleaning wire Intended use: This wire is used to clean a tube from earwax.		0.01	100.0x0.5
22. Meter Intended use: This round-shaped ruler is used to measure the size of the tube connecting the instrument to an earmold (a necessary		1.7	60.0x40.0

	accessory for prescribing an open-fit hearing aid).			
2.0×5.0	23. Thin tube box Intended use: This box is used to store thin tubes.		495.0	360.0×225.0×45.0
0.0×07.0	24. Set of thin tubes in a box Items: 0 (right and left), 1 (right and left), 2 (right and left) and 3 (right and left). Intended use: This is a complete set of thin tubes of all available sizes.		Weight of the set – 520.0 Weights of thin tubes see below.	Size of the set box – 360.0×225.0×45.0 Sizes of thin tubes see below.
0×12.0×	25. Thin tubes of various sizes for the hearing instruments, at most 100 pcs. Items: a thin tube with the sports clip: 0 (right and left), 1 (right and left), 2 (right and left) and 3 (right and left). Intended use: These are thin tubes for the open-fit hearing aid technology. They conduct the sound from the hearing instrument into the ear.		0 – 0.11 1 – 0.13 2 – 0.15 3 – 0.17	Outer diameter – 13.0 Length of the tube: 0 – 59.0 1 – 62.0 2 – 68.0 3 – 74.0 Length of the sports clip – 40.0
×9.0×19	26. Tulip mold, at most 50 pcs. Intended use: This dome is connected to the hearing instrument and used by patients with impaired hearing.		0.15	10.0×12.0×4.0
×7.0×7	27. Thick mold, small, at most 50 pcs. 28. Thick mold, medium, at most 50 pcs. 29. Thick mold, large, at most 50 pcs. Intended use: This dome is connected to the hearing instrument and used by patients with impaired hearing.		0.07 0.09 0.1	7.0×7.0×7.0 10.0×9.0×9.0 12.0×10.0×12.0

30. Open mold, small, at most 50 pcs. 31. Open mold, medium, at most 50 pcs. 32. Open mold, large, at most 50 pcs. Intended use: This dome is connected to the hearing instrument and used by patients with impaired hearing.		0.01 0.01 0.01	5.0×5.0× 6.0×8.0× 8.0×10.0×
33. Hearing instrument packaging: - Carton box - Plastic case - Protective cover - Polishing cloth - Cleaning brush Intended use: These materials are used to store the hearing instrument.	    	85.0 48.0 0,50 0.13 0.8	160.0×120.0 0 90.0×43.0× 90.0×70.0× 145.0×145.0 62.0×9.0
34. Remote control device RC-2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging Items: USB cable, AC/DC adaptor, transport bag, transport bag belt. Intended use: The remote control for the hearing aid with accessories uses the 2.4 GHz wireless protocol. Adaptor: Input voltage: 100–240		Weight of the retail package 210.0 Weight of the remote control 32.0	Dimension retail package 160.0×120.0 0 Dimension remote control 95×41×

V, current: 350 mA Output voltage: ± 5 V, 1000 mA Discharge time is up to 80 hours depending on the intensity of use. Capacity: 240 mAh			
35. Streamer TV2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging Items: USB cable, TV cable, optical audio cable, mini jack plug, and AC/DC adaptor. Intended use: This device is connected to a TV audio output to transfer the TV audio signal directly to the hearing instrument via the 2.4 GHz wireless network. Adaptor: Input voltage: 100–240 V, current: 350 mA Output voltage: ± 5 V, 1000 mA		Weight of the retail package 340.0 Weight of the streamer 128.0	Dimensions of the retail package 175.0×110.0×60.0 Dimensions of the streamer 82.0×92.0×45.0
36. Midget microphone RS UNITE, MINI MIC (EU, US, AUS, UK, CAN) c/w packaging Items: USB cable, charger, interface cable (3.5 mm jack plug), transport bag, transport bag belt. Intended use: This is an outer wireless microphone for the hearing instrument with a wireless range up to 7 m. Adaptor: Input voltage: 100–240 V, current: 80 mA Output voltage: ± 5 V, 550 mA Discharge time: 8.5 hours in the speech mode. Capacity: 170 mAh		Weight of the retail package 240.0 Weight of the microphone 21.5	Dimensions of the retail package 175.0×110.0×60.0 Dimensions of the microphone 62.0×33.0×18.0
37. Hands-free phone kit PHONE CLIP PLUS, UNITE, RS (EU, US, AUS, UK, CAN) c/w packaging Items: USB cable, AC/DC adaptor, transport bag, transport bag belt. Intended use: This device transfers the audio signal from a telephone to the hearing instrument via the wireless network. Adaptor: Input voltage: 100–240 V, current: 350 mA Output voltage: ± 5 V, 1000 mA		Weight of the retail package 195.0 Weight of the phone clip 26.5	Dimensions of the retail package 160×120×45 Dimensions of the phone clip 68.0×33.0×22.0

Discharge time: 6 hours in the speech mode and 80 hours in the stand-by mode. Capacity: 240 mAh			
38. Instructions for the behind-the-ear hearing aids Intended use: This document contains a description of the hearing instrument, recommendations for use and a warranty card.			
Set of size 13 batteries Intended use: A battery powers the hearing instrument. Battery voltage: 1.4 V Capacity: 0,385 Ah		Weight of the set – 12.5 Weight of a battery – 1.6	Dimensions set: 8.2×4.8× Diameter battery –

**** Note:**

The sizes are indicated with a tolerance of $\pm 2\%$. The weights are indicated as a minimum weight.

User guides for:

- the streamer TV2, PLUS, UNITE, RS (EU, US, AUS, UK, CAN),
 - the midget microphone RS UNITE, MINI MIC (EU, US, AUS, UK, CAN), and
 - the hands-free phone kit PHONE CLIP PLUS, UNITE, RS (EU, US, AUS, UK, CAN)
- can be found in the Internet on the following web pages:

<http://www.resound.com/en-US/hearing-aids/accessories/tv-streamer2>

<http://www.resound.com/en-US/hearing-aids/accessories/mini-microphone>

<http://www.resound.com/en-US/hearing-aids/accessories/phoneclip-plus>

Environment protection requirements for medical devices

This medical device does not cause any negative effect on human or the environment during its use, transportation, and storage.

This device is a class A waste, i.e. epidemiologically safe waste similar to domestic solid waste by its composition.

Hazard class	Morphological composition characteristics
Epidemiologically safe waste similar to domestic solid waste by its composition.	Waste products that had no contact with patient's biological fluids or with infectious disease patients. Stationery, packaging, furniture, and equipment that lost consumer appeal. Sweepings from outdoor areas and other such waste. Food waste from central nutrition units and from any unit of a medical and/or pharmaceutical institution, except infectious disease institutions (including phthisiatrical institutions).

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