Instruction for Use of the Medical Device
Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE

“APPROVE”

O&O mdc S.r.l.
Italy

Business Director

(position)

(name)

(signature)

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«day» month (numerals)

O&O mdc S.r.l.

S. Leg.: Via Orazio, 30 - 00193 Roma
S. Op.va: Via Grotte Portella, 28/34 - 00044 Fregene (RM)
C.F.P. IVA: 09539461005

1 Medical device name (designated by the manufacturer) with the indication of code OKPD (Russian Classification of Products by Economic Activities)-2:

Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE.

OKPD 2: 32.50.13.120. "Ophthalmic Devices and Accessories".

Versions:

1. Easy 2.4
2. Easy 2.2
3. Easy 2.0
4. Easy 1.8
5. One 2.4
6. One 2.0
7. TwinHFA 2.2
8. TwinHFA 2.0
9. Swing 2.0
10. Swing 1.8

2 Medical Device Classification

Depending on the potential risk degree related to the medical use, the medical device "Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE" refers to Class 2a.

Medical device type in accordance with the Medical Devices Nomenclature Classification: 120460 "Disposable Manual Intraocular Lens Injector".

3 Main Technical Characteristics

The Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE represents a set, which consists of an injector and a cartridge.

The injector of the surgical ophthalmic device for intraocular lens (IOL) implantation INSIDE is made of polycarbonate or polymethylmethacrylate for medical use. This is a specialized system that utilizes a low pressure to transfer the lens from the product cartridge by hydraulic effect and safely introduce the IOL through small surgical incisions. The surgical ophthalmic device for intraocular lens implantation INSIDE provides for a smooth, precise movement of the IOL, prevents haptics scratching and stretching, and provides for a controlled lens opening during injection into the eye cavity.

The cartridge is made of polypropylene for medical use and represents a specialized tip for the surgical ophthalmic device for intraocular lens implantation INSIDE, which allows accurate and gentle IOL placement into the eye capsule bag, prevents IOL haptics from scratching and stretching.

The surgical ophthalmic device for intraocular lens implantation INSIDE is used to fold an intraocular lens and introduce it into the eye cavity.

Scope of Delivery

The scope of delivery of the medical device includes:

1. The ophthalmic surgical device for intraocular lens implantation INSIDE – 8 pcs. or 10 pcs.
2. Sticker with information on the product – 8 pcs. or 10 pcs.
3. Instruction for use – 1 pc.

The surgical ophthalmic device for intraocular lens implantation INSIDE consists of the following parts:


- injector case;
- push rod with a special soft plunger;
- cartridge;
- stopper.

Table 1. Main Dimensions and Parameters of Injector

Parameter Name	Easy 2.4	Easy 2.2	Easy 2.0	Easy 1.8	One 2.4	One 2.0	Twin HFA 2.2	Twin HFA 2.0	Swing 2.0	Swing 1.8
<i>Main dimensions of injector</i>										
- overall length with an extended rod and inserted cartridge (mm)	147.0	147.0	147.0	147.0	146.0	146.0	166.0	166.0	147.0	147.0
- overall length with a retracted rod without the cartridge (mm)	92.0	92.0	92.0	92.0	92.0	92.0	81.0	81.0	92.0	92.0
- support length (mm)	33.8	33.8	33.8	33.8	33.8	33.8	36.7	36.7	33.8	33.8
- support width (mm)	13.9	13.9	13.9	13.9	13.9	13.9	17.0	17.0	13.9	13.9
- support thickness (mm)	2.0	2.0	2.0	2.0	2.0	2.0	5.9	5.9	2.0	2.0
- minimal applied force required for a full indentation of the piston rod (H)	8.9	8.9	8.9	8.9	8.9	8.9	2.2	2.2	8.9	8.9

Table 2. Main Parameters of Cartridge

Parameter Name	Easy 2.4	Easy 2.2	Easy 2.0	Easy 1.8	One 2.4	One 2.0	Twin HFA 2.2	Twin HFA 2.0	Swing 2.0	Swing 1.8
<i>Main Parameters of Cartridge</i>										
- overall length (mm)	32.0	32.0	32.0	32.0	32.0	32.0	49.0	49.0	39.0	39.0
- internal tip diameter (mm)	1.6	1.5	1.4	1.4	1.6	1.4	1.4	1.6	1.5	1.4
- external tip diameter (mm)	2.0	1.9	1.8	1.7	2.1	1.9	1.8	2.0	1.9	1.8
- fold width (mm)	8.0	8.0	8.0	8.0	8.0	8.0	13.5	13.5	8.5	8.5
- fold length (mm)	16.0	16.0	16.0	16.0	16.0	16.0	16.5	16.5	14.6	14.6

Handwritten signature in blue ink.

4 Indications and Contraindications, Side Effects

The surgical ophthalmic device for intraocular lens implantation INSIDE is used in course of surgical operations to replace the removed human lens affected by the cataract.

The use of the surgical ophthalmic device for intraocular lens implantation INSIDE is recommended in the following cases:

- during or after cataract extraction in the elderly (especially in case of single-eye cataracts);
- patients who may face problems with contact lenses and/or whose vision is poorly adapted to the glasses (this group may include patients with occupational restrictions, the elderly, with diseases such as yellow spots of the retina - to preserve peripheral vision).

Contraindications:

Contraindications to implantation into the eye capsule bag have not been found.

Potential pre-, intra- and postsurgical complications and adverse events:

During any surgical operation there is a risk of complications. The most common complications associated with cataract surgery or IOL implantation (some of which requiring a repeated surgery) include, but are not limited to:

- corneal endothelial damage and/or edema;
- flat anterior chamber after the lens removal;
- detachment of the posterior limiting membrane;
- dehiscence;
- thermal burn;
- astigmatism, edematic / bullous keratopathy;
- uveitis;
- bleeding at one or several eye sites;
- radial ruptures of the anterior capsule wall;
- capsule rupture;
- capsular phimosis and capsule blockade syndrome;
- delayed capsule rupture with IOL displacement to the posterior eye chamber;
- turbidity of the posterior capsule wall;
- Zinn's zonules damage with a further IOL displacement (dislocation, decentering);
- iris prolapse, iris damage, seclusio pupillae, iris capture, epithelium intussusception, pupillary block;
- IOL damage during implantation;
- post-surgical turbidity/IOL calcification;
- incorrect IOL allocation during the surgery;
- IOL displacement;
- retinal detachment;
- vitreous body departure;

- increased intraocular pressure (closed-/open-angle glaucoma);
- cystoid macular edema;
- round pathological membrane;
- hypopyon;
- entophthalmia;
- iris inflammation (iritis);
- blindness.

The following complications may lead to a repeated surgical intervention:

- poor visual acuity, including due to an incorrect calculation of the IOL optical power;
- IOL displacement (decentering, inclination, axial displacement);
- pupillary block, iris capture;
- wound bleeding;
- retinal detachment.

Provision of Information on Adverse Side Reactions

Claims submitted by customers, including quality claims, as well as information about adverse reactions and other observations potentially related to the medical device, shall be reported to the manufacturer immediately. The required information shall contain the information about the adverse reaction, therapy method used, version and serial number of the used product.

5 Procedural Requirements

This medical device shall be used only in the clinical settings.

6 Special Recommendations

A successful IOL implantation requires a high surgeon's expertise level. The surgeon shall have an extensive experience in assisting or supervising such surgeries and successfully complete one or more respective trainings before performing the implantation on his/her own.

For the procedure, the operating surgeon shall have knowledge of all materials provided by the manufacturer for a proper preparation and intraocular lens implantation.

7 Product Labeling

	Vaporized hydrogen peroxide (VHP) sterilization		Lot No.
	Shelf life		Attention

Patrick Meunier

	Disposable		Re-sterilization prohibited
	Not to be used if the package is damaged		Transportation temperature
	Storage temperature		Store in a dry place
	Not to be exposed to sunlight		See instructions for use
	Manufacturer		Maximum and minimum incision size

8 Medical Device Operating Procedure

Easy Version

Not to be used if the package is damaged.

Not to be used after shelf time expiration.

Intended for professional use.

Single use (disposable) medical device.

Reuse may lead to severe infections in the patient.

Easy medical device with a cartridge allows for an accurate IOL placement during implantation into the eye, thus minimizing any complications which may occur in course of implantation.

- Remove the injector and its cartridge from the opened contour cell package and place them in a sterile area (Fig. 2);
- Open and hold the cartridge folds in such a way that it is possible to easily position the lens;
- Fill the inner groove, cartridge end and tip with a viscous elastic solution (Fig.3);
- Always lubricate the cartridge with a viscous elastic substance;
- Place the IOL so that the lens is located in the center (Fig.4). For lenses with "C loop" haptic shape (loop in the form of C) make sure that the leading haptic element is oriented to the left;
- Fold the closing haptic element in the direction of the optical system and carefully inspect the lens location in the cartridge (Fig. 5);
- Make sure that the back haptic element is not the closing one and is not located on the outer side of the cartridge;
- Make sure that neither the haptic element nor the optical system is jammed between the folds (Fig. 6);
- In case of jamming, inject the viscous elastic substance to the back of the cartridge ;
- Make sure that the elastomeric plug is smoothly and properly secured at the pin tip;
- For products of the version Easy, install the cartridge into the injector case by inserting it into the hole located on the frontal part of the injector (Fig. 7);
- Slowly push the stock until the leading haptic element reaches the tip of the cartridge;
- Place the cartridge tip into the incision;
- To deliver the lens to the anterior eye chamber, press the stock (Fig. 8);
- Surgeons independently determine the applied force as may be necessary and depending on the surgical intervention type.

Botchko N. A.

Fig.1

Fig.2


Fig.3


Fig.4


Fig.5


Fig.6

Fig.7



Fig.8

Diagram 1. Operating Procedure for Easy Version

For the New

[Handwritten signature]

One Version

Not to be used if the package is damaged.

Not to be used after shelf time expiration.


Intended for professional use.

Single use (disposable) medical device.

Reuse may lead to severe infections in the patient.

One medical device with a cartridge allows for an accurate IOL placement during implantation into the eye, thus minimizing any complications which may occur in course of implantation.

- Remove the injector and its cartridge from the opened contour cell package and place them in a sterile area (Fig. 2);
- Open and hold the cartridge folds in such a way that it is possible to easily position the lens;
- Fill the inner groove, cartridge end and tip with a viscous elastic solution (Fig.3);
- Always lubricate the cartridge with a viscous elastic substance;
- Place the IOL so that the lens is located in the center. For lenses with "C loop" haptic shape (loop in the form of C) make sure that the leading haptic element is oriented to the left (Fig.4);
- Fold the closing haptic element in the direction of the optical system and carefully inspect the lens location in the cartridge (Fig. 5);
- Make sure that the back haptic element is not the closing one and is not located on the outer side of the cartridge;
- Make sure that neither the haptic element nor the optical system is jammed between the folds (Fig. 6);
- In case of jamming, inject the viscous elastic substance to the back of the cartridge;
- Make sure that the elastomeric plug is smoothly and properly secured at the pin tip;
- Slowly push the stock until the leading haptic element reaches the tip of the cartridge;
- Place the cartridge tip into the incision;
- To deliver the lens to the anterior eye chamber, press the stock (Fig. 7);
- Surgeons independently determine the applied force as may be necessary and depending on the surgical intervention type.

Polina Mennikov

Fig.1


Fig.2

Fig.3

Fig.4


Fig.5


Fig.6


Fig.7

Diagram 2. Operating Procedure for One Version

Handwritten signature

Swing Version

Not to be used if the package is damaged.

Not to be used after shelf time expiration.

Intended for professional use.

Single use (disposable) medical device.


Reuse may lead to severe infections in the patient.

Swing medical device with a pre-installed cartridge allows for an accurate IOL placement during implantation into the eye, thus minimizing any complications which may occur in course of implantation.

- Remove the injector and its cartridge from the opened contour cell package and place them in a sterile area (Fig. 2);
- Hold the injector and fill the inner groove, faucet and cartridge tip with a viscous elastic solution (Fig.3);
- Always lubricate the cartridge with a viscous elastic substance;
- Place the IOL so that the lens is located in the center. For lenses with "C loop" haptic shape (loop in the form of C) make sure that the leading haptic element is oriented to the left (Fig.4);
- Close the folds and inject more viscous elastic solution through the respective hole in the cartridge (Fig. 5,6);
- For products of the version Swing, when used with hydrophilic IOL, make sure, that the elastomeric plug is smoothly and properly secured at the pin tip;
- When the Inside delivery systems are used, slowly push the stock until the leading haptic element reaches the tip of the cartridge;
- Place the cartridge tip into the wound;

When Swing products are used. to deliver the lens to the anterior eye chamber, press the stock (Fig. 7);

- Surgeons independently determine the applied force as may be necessary and depending on the surgical intervention type.

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
Fig.6

Fig.7

Diagram 3. Operating Procedure for Swing Version

Patent Memo

Twin Version

Not to be used if the package is damaged.

Not to be used after shelf time expiration.


Intended for professional use.

Single use (disposable) medical device.

Reuse may lead to severe infections in the patient.

Twin medical device with a pre-installed cartridge allows for an accurate IOL placement during implantation into the eye, thus minimizing any complications which may occur in course of implantation.

- Remove the injector and its cartridge from the opened contour cell package and place them in a sterile area (Fig. 2);
- Hold the injector and fill the inner groove, faucet and cartridge tip with a viscous elastic solution (Fig.3);
- Always lubricate the cartridge with a viscous elastic substance;
- Place the IOL so that the lens is located in the center. For lenses with "C loop" haptic shape (loop in the form of C) make sure that the leading haptic element is oriented to the left (Fig.4);
- Close the folds and inject more viscous elastic solution through the respective hole in the cartridge (Fig. 5,6);
- For products of version Twin, use a disconnecting device for capturing and/or pushing depending on the used delivery method;
- For products of version Twin, when used with hydrophilic IOL, make sure, that the elastomeric plug is smoothly and properly secured at the stock tip;
- When any Twin delivery systems are used, slowly push the stock until the leading haptic element reaches the tip of the cartridge;
- Place the cartridge tip into the incision;
- When a Twin product is used to deliver the lens to the anterior eye chamber, grip the stock (Fig. 7);
- Surgeons independently determine the applied force as may be necessary and depending on the surgical intervention type.

A handwritten signature in blue ink is written over a circular official stamp. The stamp contains a star and some text, including "NOTARIAT". Below the signature, the page number "- 13 -" is printed.



Diagram 4. Operating Procedure for Twin Version

9 Product Sterilization

This medical device is supplied sterile.

Sterilization method – vaporized hydrogen peroxide (VHP).

The sterilization method has been validated in accordance with the requirements of ISO 14937:2009 (GOST R ISO 14937-2012) and ISO 11135 (GOST ISO 11135-2012).

10 Warranty Obligations

The company O&O MDC S.R.L., Italy, assumes the warranty obligation toward the end user in respect of the medical device provided there are no defects and manufacturer's specifications are complied with, within 12 months from the date of sale, subject to the terms and conditions set out in the instructions for use.

11 Transportation

Transportation may be performed out by all transportation types in the closed vehicles in accordance with the rules of transportation, applicable to the respective transportation type at the temperature from 5°C to 45°C.

12 Shelf Life

The Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE is supplied sterile and remains sterile provided the package is intact and not opened, for 36 months.

This Medical Device shall not be used after the expiration date specified on the package.


Storage: from 10 °C and to 30 °C, to be stored in a dry, dark place.

13 Disposal and Destruction of the Medical Device

This medical device: "Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE" and its parts are disposed of as medical waste referred to Class B "Epidemiologically Hazardous Waste" in accordance with SanPiN (Russian Sanitary Rules and Norms) 2.1.7.2790, as in the course of a surgery they contact and get contaminated with biological fluids and secretions of body tissues (blood and/or other biological fluids).

The packaging materials, a cardboard box and blister package of this product shall be disposed in accordance with SanPiN 2.1.7.2790, as Class A Waste (Epidemiologically Safe Waste, similar by composition to solid household waste).

The Surgical Ophthalmic Device for Intraocular Lens Implantation INSIDE is to be disposed if the consumer package integrity violation or shelf life expiration in the sterile package.

Rachel Mennin

The product and its package shall be disposed in accordance with the instruction (or another document) regulating medical waste treatment at specific medical institution. This document shall be developed in accordance with the requirements of SanPiN 2.1.7.2790.10 "Sanitary and Epidemiological Requirements to Medical Waste Treatment".

14 Responsibility

The manufacturer, O&O MDC S.R.L. shall not be held liable for any improper selection of the medical device by the attending physician, improper handling and use, as well as for the used equipment and other medical malpractice.

15 Information about the Medical Device Manufacturer

O&O MDC S.R.L.

Via Orazio, 30, 00193, Roma, Italia

+39 06 940 16099 ono-assist@oo-mdc.com

Medical Device manufacturing site:

Via di Grotte Portella 28-34, 00044 Frascati (RM) – Italy

tel: +39 06 9401 6099

e-mail: ono-assist@oo-mdc.com

Official Representative in the Russian Federation: MK-Business Ltd., 620130, Russia, Ekaterinburg, Aviationnaya str., bldg. 48, apt. 80, 8 (343)286 23 92, info@mk-business.ru

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