APPROVED BY

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9. Sceau :

10. Signature : **MICHEL LERNOUT**
PREMIER AVOCAT GENERAL

This instruction for use has been prepared to provide information for purposes of registration of a medical device in the territory of the Russian Federation. The information is provided in the amount required in accordance with the Rules for Registration of Medical Devices (Russian Government Decree No. 1416 dated 27.12.2012). The technical document contents comply with the requirements of the Order of the Ministry of Health of the Russian Federation No. 11n dated 19.01.2017 "On approval of requirements to contents of technical and operational documentation of the medical product manufacturer (producer)".

1. Name of the device:

"Sterile single-use intravascular occlusal balloon catheters with double lumen, in versions"

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quality.assistant@balt.fr

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Development, production and sale of the device complies with the requirements of EN ISO 13485:2012 + AC : 2012 (Certificate No. 513975 MP2012, issue date February 6, 2018, valid until March 31, 2019, notified body for medical devices: DQS Medizinprodukte GmbH); the device bears CE-marking (EC certificate number: 513975MR2; date of issue – 01.12.2017, date of expiry – 31.07.2022, notified body: DQS Medizinprodukte GmbH).

Device code according to the All - Russian Classifier of Products OK-005-93: 94 3630;

OKPD2 - 32.50.13.110 - Syringes, needles, catheters, cannulas and similar instruments;

Class depending on the degree of potential risk of medical use in accordance with the medical devices nomenclature classifier approved by the Order of the Ministry of Health of the Russian Federation No. 4n-3 dated June 6, 2012;

Class depending on the degree of potential risk of medical use in accordance with the medical devices nomenclature classifier approved by the Order of the Ministry of Health of the Russian Federation No. 4n dated June 6, 2012 - 3

Device classification according to the Medical Device Nomenclature Classification stipulated by the Order of the Ministry of Health of the Russian Federation of June 6, 2012, No.4n – 131560;

This device is sterile. For single use.

Conditions of use: medical and medical-preventive health care institutions, only as prescribed by a medical specialist.

2. Purpose

The device is intended for use in peripheral and neurovascular system where temporal occlusion is required. The balloon catheters provide temporary vascular occlusion which is useful in selective stopping or controlling blood flow as well as in vasospasm treatment. The balloon catheters also provide balloon embolization of intracranial aneurysms.

3. Device special features

Intravascular double lumen balloon catheters ECLIPSE2L and COPENIC2L are balloon catheters with an internal lumen for the guide wire and a separate lumen for inflation and deflation of the balloon. Inflation and deflation of the balloon catheters is allowed with any position of the guide wire. The intravascular double lumen balloon catheters are of oval shape. They are produced with two different diameters as well as with various balloon lengths (see label). The Eclipse2L and Copernic2L catheters are equipped with 3 radio opaque markers, one of which enables visualization of the distal end of the catheter while the two other markers indicate borders of the cylindrical part of the balloon. Intravascular double-lumen balloon catheters ECLIPSE2L and COPENIC2L are intended for use in peripheral and neurovascular system where temporal occlusion is required. The balloon catheters provide temporary vascular occlusion which is useful in selective stopping or controlling blood flow as well as in vasospasm treatment. The balloon catheters also provide balloon embolization of intracranial aneurysms. They are also intended for use in the peripheral vasculature for administration of diagnostic agents (contrast solutions) and therapeutic agents (embolization material, compatible with the inner diameter of the balloon catheter). The catheters labeled ECL2L6x7SN and COP2L6x10SN have an additional proximal marker located 3 cm away from the distal marker. It enables precise deployment of detaching embolizing coils (see user manual for use of the detaching coils). Intravascular double-lumen balloon catheters ECLIPSE2L and COPENIC2L are also equipped with a distal hatch to deflate the balloon before use. Technical characteristics of the devices are given in clause 4. Each device is packed in an individual package, each individual package is packed in a cardboard box, which also contains user manual. For more information on packaging and labeling, see clause 5. Appearance of the device is shown in Figures 1-4.


Figure 1. Appearance of the device without package

Figure 2: Appearance of the device ECLIPSE2L in the final package for sterilization

Figure 3: Appearance of the device COPERNIC2L in the final package for sterilization

Figure 4: Device structure

Note. All components of the product are a single unit and are not used separately from each other.

4. Technical specification

Table 1 – Materials of the device

COMPONENT	RAW MATERIAL, USED IN THE CATHETERS ECLIPSE2L / COPERNIC2L	SUPPLIER
Connector (Hub)	Polyamide, mark Grilamid TR90	EMS-Chimie Ag, Switzerland
	Pigment green (green C), CAS # 1328-53-6	PLASTICOULEUR, France
Rigid tube	Polyamide, mark Grilamid L20G	EMS-Chimie Ag, Switzerland
	Pigment green (green C), CAS # 1328-53-6	PLASTICOULEUR, France
Intermediate tube	Polyetherblockamide PEBAX 7033 SA	ARKEMA, France
	Pigment blue CAS # 57455-37-5	PLASTICOULEUR, France
	Pigment White (WHITE TiO2) CAS # 13463-67-7	PLASTICOULEUR, France
Flexible tube	Polyetherblockamide PEBAX4033SP01	ARKEMA, France
	Polyetherblockamide PEBAX2533SD02	ARKEMA, France
Super flexible tube	Thermoplastic elastomer mark THERMOPLASTIC TF6ST	KRAIBURG, Germany
	Pigment blue CAS # 57455-37-5	PLASTICOULEUR, France
	Pigment White (WHITE TiO2) CAS # 13463-67-7	PLASTICOULEUR, France
Extra flexible tube	Thermoplastic elastomer mark THERMOPLASTIC TF6ST	KRAIBURG, Germany
	Thermoplastic polyurethane mark ESTANE 58562	LUBRIZOL, Germany
Internal tube	Polytetrafluorethylene PTFE	Lake Region Medical, USA
Hydrophilic coating	Polyvinyl pyrrolidone (PVP)	Accros Organics, Belgium
Balloon	CHRONOPRENE 25A	AdvanSource Biomaterials, USA
Braid ‡	Nitinol (nickel-titanium alloy, 45% titanium CAS # 7440-32-6, 55% nickel CAS # 7440-02-0)	FORT WAYNE Metals, Ireland
Adhesive ‡	Brand LOCTITE® 4011 Composition: Epoxy resin CAS # 25068-38-6 - 35-50% Butanediol diglycidyl ether CAS # 2425-79-8 - 25-35% Isophorone diamine CAS # 2855-13-2 - 7-10% Trimethylhexamethylenediamine CAS # 25620-58-0 - 5-7% Araldite® 2020	HENKEL, Germany
	Ethyl-2-cyanoacrylate CAS # 7085-85-0 - 99.9% - 100% Hydroquinone CAS # 123-31-9 - 0,01 - <0,1%	HUNTSMAN, USA
Radiopaque marker ‡	Gold 99% brand z199.99, impurities 1% (not more) Lead 0.003, Iron 0.004, Antimony 0.001, Bismuth 0.002, Copper 0.007, Silver 0.008	Johnson Matthey, United Kingdom
Mandren *	Stainless steel (AISI SUS 304L)	Fort Wayne Metals, USA
Protective rim *	Polyethylene (Rigidex 5502 GA)	INEOS Olefins & Polymers, UK
Rim holders *	Polyethylene (Rigidex 5502 GA) Pigment blue CAS # 57455-37-5	INEOS Olefins & Polymers, UK

* Components of the product that do not come into contact with circulating blood, but contact directly or indirectly with the human body

(‡): The components are not intended for direct or indirect body contact.

Note: The pigment ratio is generally from 1 g of dye per 5 kg of polymer.

Table 3: The device technical characteristics

	ECL2L6 x7SN	ECL2L6 x9SN	ECL2L6 x12SN	ECL2L6 x15SN	ECL2L6 x20SN	ECL2L6 x7	ECL2L6 x9	ECL2L6 x12	ECL2L6 x15	ECL2L6 x20	COP2L6 x10SN	COP2L6 x15SN	COP2L6 x20SN	COP2L6 x30SN	COP2L6 x10	COP2L6 x15	COP2L6 x20	COP2L6 x30		
Total length, mm	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)	1600 (+33/-20)		
Effective length, mm	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)	1500 (+33/-20)		
Diameter and length of the distal part, mm	0,68 x3	0,68x3	0,68x3	0,68 x3	0,68x3	0,68x3	0,68x3	0,68x3	0,68x3	0,68x3	0,85x3	0,85x3	0,85x3	0,85x3	0,85x15	0,85x15	0,85x15	0,85x15		
Overall dimensions of the balloon, (Max. Ø x length) mm	6x7	6x9	6x12	6x15	6x20	6x7	6x9	6x12	6x15	6x20	6x10	6x15	6x20	6x30	6x10	6x15	6x20	6x30		
The maximum volume of balloon (the balloon is hypercompatible with no tolerance required), ml	0.20	0.25	0.30	0.40	0.45	0.20	0.25	0.30	0.40	0.45	0.25	0.40	0.45	0.65	0.25	0.40	0.45	0.65		
Maximum filling time of the balloon, sec	not more than 1 minute The filling time of the balloon depends on the filling rate used by the surgeon with a syringe																			
Maximum time of balloon discharge, sec	6	8	10	14		16	6	8	10	14		16	8	14	16	22	8	14	16	22
Maximum allowable pressure in the balloon:	a balloon with a low pressure of 50 kPa (0.5 atm)																			
Balloon burst pressure	Not more than 200 kPa (2 atm)									Not more than 350 kPa (3.5 atm)										
Tensile strength, N	2±0.1																			
Inner diameter (Main lumen), mm	0.41	0.41	0.41	0.41	0.41	0.41	0.41	0.41	0.41	0.41	0.60	0.60	0.60	0.60	0.60	0.60	0.60	0.60	0.60	
Double lumen thickness, mm	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	0.9	
Double lumen width, mm	1.15	1.15	1.15	1.15	1.15	1.15	1.15	1.15	1.15	1.15	1.26	1.26	1.26	1.26	1.26	1.26	1.26	1.26	1.26	
Weight, g	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.3±0.038	3.6±0.004	3.6±0.004	3.6±0.004	3.6±0.004	3.6±0.004	3.6±0.004	3.6±0.004	3.6±0.004	
Connector type	Luer lock according to ISO 594																			
Hydrophilic coating	Polyvinylpyrrolidone (PVP) is freely soluble in alcohol and in water, transforming into a colloidal solution. This is a slightly yellow or white amorphous odorless powder.																			

Unspecified permissible values: ± 1%

Appearance of the device: The device surface has no foreign inclusions, bubbles or streaks. The surface is smooth, has no technological and surface defects and provides minimal trauma to vessels when being used. The catheter is biologically safe and compatible with biological body fluids. Cannula and other parts of the catheter do not allow fluid leakage.

The device is not implantable and is removed from the human body after being used.

5. The device range:

Sterile single-use intravascular occlusal balloon catheters with double lumen, in versions:

- I. Sterile single-use intravascular occlusal balloon catheters with double lumen ECLIPSE2L: ECL2L6x7SN, ECL2L6x9SN, ECL2L6x12SN, ECL2L6x15SN, ECL2L6x20SN, ECL2L6x7, ECL2L6x9, ECL2L6x12, ECL2L6x15, ECL2L6x20.
 - II. Sterile single-use intravascular occlusal balloon catheters with double lumen COPERNIC2L: COP2L6x10SN, COP2L6x15SN, COP2L6x20SN, COP2L6x30SN, COP2L6x10, COP2L6x15, COP2L6x20, COP2L6x30.
- user manual (to be supplied with each device) - 1 pc.

6. Marking and packaging

Intravascular double-lumen sterile balloon catheter for single use is supplied packaged in an individual Tyvek pack which is the final package for sterilization and in an individual cardboard package. User manual is attached to each device.

6.1. Technical characteristics of the final package for sterilization:

One side of the packs is made of Tyvek® 2 FS UNCo 7502918M coating and the other side is made of PET/PE 65 in accordance with standard NF EN ISO11607-1 (2009). The final package for sterilization is glued together.

Table 4: Technical characteristics of the pack

Basic properties	Results
Pack overall dimensions (LxWxT), mm	290±5x255±5x0.3
Width of adhesive layer on the long side, mm	8
Width of adhesive layer on the wide side, mm	8
Thickness of adhesive layer, mm	0.3
Tyvek paper - Unit weight	59.5 g/m² (56.5-62.5)
Tyvek paper - Thickness	155 µm (75- 235)
Tyvek paper – Porosity (The Gurley Method)	22 s (9-35)
Tyvek paper – Porosity (The Bendtsen Method)	529 ml/min (155-885)
Tyvek paper - Mullenburst test	925 kPa (760-1,110)
Tyvek paper - MD breaking strength test	156 N/2.54 cm (135 to 177)
Tyvek paper - XD breaking strength test	157 N/2.54 cm (127 - 187)
Tyvek paper - Elmendorf MD breaking test	2.8 N (1.7-4.0)
Tyvek paper - Elmendorf XD breaking test	3.7 N (2.8-4.6)
Tyvek paper - Stratification	2.70 N/2.54 cm (2.05-3.25)
Tyvek paper - Hydrostatic flow pressure	145 cm H ₂ O (125- 165)
PET film - Thickness	65 µm ± 5%
PET film - Total weight	63g/m²
PET film - Breaking test	30- 55 N/mm²
PET film - Elongation at break	50- 100 %
PET film - Elongation at break	80- 160 %
PET film - O ₂ permeability	90m³/m²/bar
PET film - CO ₂ permeability	345 m³/m²/bar
PET film - N ₂ permeability	50 m³/m²/bar
PET film - Water vapor permeability	< 3.0 g/m²
PET film - Sealing parameters	115-150°C
Sealing	> 1.5 N / 15 mm
Manufacturer of final package for sterilization	SUDPACK MEDICA, France

Unspecified permissible values: ±1%

Table 5: Technical characteristics of the individual cardboard pack

Basic properties	Results
Box overall dimensions (LxWxT), mm	260±4x300±3x14±1
Cardboard density	365 g/m²±1%
Material	Cardboard
Manufacturer	IMPRIMERIE LLP/BUIGNET Address: 29, avenue des Chataigniers – 95150 Taverny

6.2. Appearance of the final package marking:

Figure 5: ECL2L6x7SN package and box marking

Figure 6: ECL2L6x9SN package and box marking

Figure 7: ECL2L6x12SN package and box marking

Figure 8: ECL2L6x15SN package and box marking

Figure 9: ECL2L6x20SN package and box marking

Figure 10: ECL2L6x7 package and box marking

Figure 11: ECL2L6x9 package and box marking

Figure 12: ECL2L6x12 package and box marking

Figure 13: ECL2L6x15 package and box marking

Figure 14: ECL2L6x20 package and box marking

Figure 15: COP2L6x10SN package and box marking

Figure 16: COP2L6x15SN package and box marking


Figure 17: COP2L6x20SN package and box marking


Figure 18: COP2L6x30SN package and box marking

Figure 19: COP2L6x10 package and box marking


Figure 20: COP2L6x15 package and box marking


Figure 21: COP2L6x20 package and box marking


Figure 22: COP2L6x30 package and box marking

Figure 23: Box layout


Figure 24: Draft of the marking for medical product "Sterile single-use intravascular occlusal balloon catheters with double lumen, in versions" designed for medical product distribution in the territory of the Russian Federation. The marking requirements are applied for all products.

Marking on the final sterilization package is identical to marking on the individual cardboard package.

The marking contains:

- Device name and model;
- Information on the manufacturer;
- Schematic representation of the device;
- Quantity;

- Diameter of the compatible catheter;
- Diameter of the compatible guidewire;
- Lot number;
- Serial number;
- Date of manufacture;
- Expiration date;

Technical characteristics of the product:

- Ø - Diameter of the balloon
- L - Length of the balloon
- Volume - Maximum balloon volume
- A - Overall length, mm
- B - Effective length, mm
- C - Length of distal part, mm
- D - Diameter of distal part, mm
- E - Internal diameter (main lumen), mm
- F - Thickness of double lumen, mm
- G - Width of double lumen, mm
- "Do not use if package is damaged" symbol;
- "Do not reuse" symbol;
- "Do not resterilize" symbol,
- "Do not expose to sunlight" symbol;
- "Refer to user manual" symbol;
- "Caution! Refer to accompanying documentation" symbol.
- Information on the method of sterilization;

Symbol legends list

	Reference number		Refer to user manual
	Serial number		Keep dry!
	Use before		Prohibition on resterilization
	Date of manufacture		Does not contain diethylhexyl phthalate, may contain particles of polyhexahydrotriazine
	CE label		European article number
	EO-sterilized		Manufacturer
	Attention! Refer to user manual		"Caution! Fragile!" symbol
	Prohibition on reuse		Non-pyrogenic
	Do not expose to sunlight		Non-toxic
	Do not use if package is damaged		

Cardboard packages with the device are packed in a corrugated cardboard transport box.

One transport box contains – no more than 50 cardboard packages.

Transport box marking contains:

- Device name and model;
- Information on the manufacturer;
- Quantity;
- Lot number;
- Serial number;
- Date of manufacture;
- Expiration date;
- "Do not expose to sunlight" symbol;
- "Keep dry!" symbol;

7. Compliance with the Russian Federation standards & Control methods in the Russian Federation standards

GOST ISO 10555-1-2011 Sterile, single-use intravascular catheters. Part 1. General technical requirements;

GOST ISO 10555-4-2012 Sterile, single-use intravascular catheters. Part 4. Catheters for balloon expansion;

GOST R 50444-92 "Medical devices, apparatuses and equipment. General technical specifications";
 GOST 31814-2012. "Conformity assessment. General sampling rules for products testing upon confirmation of conformity".
 GOST 31214-2016 "Medical devices. Requirements for samples and documentation provided for toxicological, sanitary and chemical tests, sterility and pyrogenicity tests".
 GOST ISO 10993-1-2011. "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and research".
 GOST ISO 10993-4-2011. "Medical devices. Biological evaluation of medical devices. Part 4. Selection of tests for interaction with blood".
 GOST ISO 10993-5-2011 Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
 GOST ISO 10993-7-2016. Medical devices. Biological evaluation of medical devices. Part 7. Ethylene oxide sterilization residuals.
 GOST ISO 10993-10-2011. Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity.
 GOST ISO 10993-11-2011. "Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity".
 GOST ISO 10993-12-2015. "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and control materials".
 GOST ISO 10993-17-2011. Medical devices. Biological evaluation of medical devices. Part 17. Establishment of allowable limits for leachable substance
 GOST ISO 10993-18-2011. Medical devices. Biological evaluation of medical devices. Part 18. Chemical characterization of materials.
 GOST R ISO 15223-1-2014 "Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements";
 GOST ISO 11607-2011 "Packaging for terminally sterilized medical devices";
 GOST ISO 11135-2012 Medical devices. Validation and routine control of ethylene oxide sterilization.
 GOST EN 556-1-2011. Sterilization of medical devices. Requirements for medical devices to be designated "sterile". Part 1. Requirements for terminally sterilized medical devices.
 GOST ISO 11138-2-2012 Sterilization of health care products. Biological indicators. Part 2. Biological indicators for ethylene oxide sterilization.
 GOST ISO 11737-1-2012 Sterilization of medical devices. Microbiological methods. Part 1. Estimation of population of microorganisms on products.
 GOST R 52770-2016. Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests.
 GOST 25375-82 "Methods, means and conditions of sterilization and disinfection of medical use products. Terms and definitions".
 GOST ISO 14971-2011 Medical devices. Application of risk management to medical devices.
 GOST R IEC 62366-2013 Medical devices. Application of usability engineering to medical devices.
 GOST ISO 14698-1-2005. Cleanrooms and associated controlled environments. Biocontamination control. Part 1. General principles and methods

9. *Description of accessories and medical products used in conjunction with the registered product:*

Medical product "Sterile single-use intravascular occlusal balloon catheters with double lumen, in versions" is suitable for use with any guides with a maximum diameter of 0.35 mm, embolization devices, and also used in conjunction with a liquid embolic agent (Squid) produced by BALT EXTRUSION SAS, as well as for the installation of self-expanding stents:

ECLIPSE 2L is suitable for the installation of LEO+Baby stent manufactured by BALT EXTRUSION SAS;

COPERNIC 2L is suitable for installation of LEO+3.5 or SILK+2.0 to 4.5 stents manufactured by BALT EXTRUSION SAS.

10. **Transportation and storage**

Transportation conditions: air temperature from -10°C to +35°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 - 1,060 hPa away from direct sunlight.

Storage conditions: store the device in the manufacturer's package in a dry, clean place, away from heat sources and direct sunlight, at controlled temperature. Shelf life is indicated on the label of each device. Air temperature from +15°C to +27°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 - 1,060 hPa.

Operating conditions: air temperature from +20°C to +25°C at relative humidity from 35% to 65% (non-condensing) and atmospheric pressure of 500 - 1,060 hPa away from direct sunlight. The device is temporary operated in the human body at the temperature from +36°C to +42°C. The device is operated in the human body for less than 24 hours.

11. **Information on the device sterility**

The device is supplied sterilized with Ethylene Oxide sterilization method.

12. **Disposal information**

The device shall be disposed in accordance with SanPiN 2.1.7.2790-10 for class B waste (epidemiologically hazardous waste). After the apparatus decontamination methods using physical methods and changing the waste appearance, excluding the possibility of their reuse, class B waste may be accumulated, temporarily stored, transported, destroyed and buried together with class A waste (epidemiologically safe waste approximated to solid household waste). A package of disinfected class B medical wastes shall have a mark indicating the completed waste disinfection.

13. **Warranties**

Shelf life - 5 years. Sterilization validity - 5 years.

Address for complaints submission:

Perinatal Medicine Center Limited Liability Company

Registered/Business address: 109004, Moscow, Tovarishchesky per., 20, bldg. 4

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+33 (0)1 39 89 46 41

+33 (0)1 37 17 03 46

14. General indications, contraindications and side effects

Indications:

Balloon catheters ECLIPSE2L and COPENIC2L are intended for use in the central nervous system vessels and peripheral vessels where their temporal occlusion is required. The balloon catheters provide temporary vascular occlusion which is useful in selective stopping of blood flow as well as in vasospasm treatment. The balloon catheters are also applied to provide balloon embolization of intracranial aneurysms with balloon assistance.

They are also recommended for:

use in peripheral vessels and in the central nervous system vessels and peripheral vessels to inject diagnostic (radiopaque materials) and medical products (embolization materials, compatible with the inner diameters of ECLIPSE 2L and COPENIC 2L balloon catheters) and stent placement.

Contraindications:

Application of ECLIPSE 2L and COPENIC 2L catheters is contraindicated if the doctor considers the procedure with application of ECLIPSE 2L and COPENIC 2L catheters can aggravate the patient's condition.

- Not intended for embolectomy or angioplasty.
- Not intended for application in coronary vessels.
- Not intended for application in children or neonates.

Potential complications:

Possible complications, connected with use of balloon catheters or with intravascular procedures include, without limitation, the following disorders:

- complications in the access site area;
- allergic reaction;
- aneurysm rupture;
- neurological disorders;
- transitory ischaemic attack;
- vascular occlusion;
- puncture of the vessel or aneurysm;
- vascular spasm;
- hematoma at the access site;
- embolism, ischemia;
- hemorrhage;
- pseudoaneurysm;
- epileptic seizure;
- cerebral accident;
- infection;
- vessel dissection;
- thrombus formation;
- death.

15. Application method.

15.1. Preparation before use

Choose a balloon catheter according to the vessel's size.

Before extraction of the balloon catheter from the dispensing tube, damp the guidewire lumen with physiological saline.

Extract the balloon catheter, pulling it out of the dispensing tube. In case of resistance, repeat the damping procedure until the balloon catheter is wet enough to be easily extracted out of the dispensing tube.

Perform thorough visual examination of the balloon catheter to ensure it is not damaged.

Do not allow the balloon catheter to get dry before insertion it in the sheath introducer.

Do not repack the moisturized balloon catheter.

15.2. Balloon preparation

Wash the guidewire lumen of the balloon catheter with physiological saline. Insert the damped guidewire carefully in the respective lumen of the balloon catheter.

CAUTION: while placing the distal end of the balloon catheter on the guidewire, ensure the distal end of the balloon catheter is not damaged.

prepare the mixture of radiopaque material (at least 50%) with physiological saline.

CAUTION: viscosity and concentration of the radiopaque material will influence the time of inflation and deflation of the balloon.

Fill the balloon inflation port at its base with the mixture of radiopaque material and physiological saline by means of the introducing needle set on a syringe of 1 ml volume.

Remove the introducing needle and connect the tap to the 1 ml volume syringe, filled with the mixture of radiopaque material and physiological saline.

Wash the tap with the mixture of radiopaque material and physiological saline and connect it to the inflation port at the balloon catheter's base.

Moisturize the balloon well enough by immersing it in the physiological saline for at least 30 s. Take the balloon out of the physiological saline and, holding it upright with its distal end up, slowly inject the mixture of radiopaque material and physiological saline. Initially, the balloon is inflated with air from the catheter and then it is gradually filled with the liquid while the air passes outflows from distal micro holes for removal of air. When the balloon is filled with the liquid completely, the procedure of air removal is over; after it, the balloon can be immersed in the physiological saline again.

CAUTION: the liquid should not reach distal holes for air removal until complete air removal therefore the balloon should be held upright with its distal end up.

CAUTION: during the air removal procedure, the liquid should be injected slowly in order to avoid bursting of the balloon.

After air removal from the balloon catheter inflating lumen, deflate the balloon completely holding it in the state of immersion. Inflate the balloon with the mixture of radiopaque material and physiological saline again to ensure absence of air in it. If the balloon still contains air, hold it upright with its distal end down to enable air bubbles rise to the proximal end. Deflate the balloon: the air will be driven to the syringe by the liquid inside the balloon.

CAUTION: during preparation, the balloon catheter should remain moisturized therefore it should be immersed in physiological saline from time to time when necessary.

Deflate the balloon and then inflate it again up to its nominal volume to ensure the lack of any defects or damages in the balloon catheter before use, do not apply the catheter in case of detection of any malfunctions. Deflate the balloon again, holding its distal end in physiological saline, and allow the pressure in the catheter to stabilize. The balloon catheter is ready for use upon complete filling of the catheter and the balloon.

CAUTION: do not connect any high pressure devices to the port to inflate the balloon as it can result in its bursting.

CAUTION: do not inflate the balloon with air or any other gas.

CAUTION: improper preparation can lead to penetration of air into the system, which may prevent due X-ray visualization.

15.3. Usage

Connect the rotating hemostatic valve to the balloon catheter lumen for the guidewire. Install the line of continuous washing with physiological saline and connect it to the side sleeve of the rotating hemostatic valve.

Choose an appropriate sheath catheter or diagnostic catheter. Fix the rotating hemostatic valve at the proximal base of the sheath catheter or diagnostic catheter. To prevent blood reflux in the lumen of the catheter, connect the line of continuous washing with physiological saline to the side sleeve of the rotating hemostatic valve.

Open the rotating hemostatic valve at the proximal base of the sheath catheter or diagnostic catheter. Insert the balloon catheter via the rotating hemostatic valve.

Deploy the balloon catheter and guidewire to the desired place in the vascular system under X-ray monitoring. Fix the valve of the rotating hemostatic valve around the balloon catheter carefully to prevent any leakage from the rotating hemostatic valve. The hemostatic valve should ensure advance of the balloon catheter after fixation.

CAUTION: do not tighten the rotating hemostatic valve too hard around the balloon catheter. Too tight fixation of the valve can make the balloon inflation/deflation process slower.

CAUTION: do not advance the balloon catheter or guidewire in case of resistance. In case of resistance, it is necessary to establish the reason for such resistance by means of X-ray monitoring.

By means of the 1 ml syringe, filled with the appropriate solution of radiopaque material, inflate the balloon slowly up to the required diameter.

CAUTION: do not exceed the recommended inflation volume, indicated on the label, in order to avoid bursting of the balloon.

CAUTION: to ensure safety of the patient, inflation and deflation of the balloon should always be performed under X-ray visualization.

Close the stopcock after inflation if necessary.

Deflation of the balloon should be made under X-ray monitoring to ensure complete deflation before extraction of the balloon catheter. After the procedure is complete, remove the balloon catheter and guidewire slowly.

15.4. PRECAUTIONS FOR USE

Do not use if the package is opened or damaged. If the package is intact, the devices are sterile.

A single-use medical device. Do not reuse. Any reuse of the device results in a high risk of infection of the patient and deterioration of the device performance.

Do not resterilize these devices.

Keep in a dry, dark place at room temperature.

Observe the shelf life.

These devices shall only be used by medical specialists in the sphere of interventional neuroradiology and/or interventional radiology.

Check the vessel's size by means of X-ray examination. Ensure the balloon catheter complies with the vessel's size. Do not exceed the recommended inflation volume, indicated on the label, in order to avoid bursting of the balloon. No test for compatibility of use of the balloon catheter with liquid embolizing agents was performed with the exception of embolism systems Onyx (® Ev 3) and Squid (® Emboflu A.G.).

Viscosity and concentration of the radiopaque material will influence the time of inflation and deflation of the balloon. During preparation, deflation of the balloon is only allowed with the balloon distal tip immersed in physiological saline or radiopaque material to avoid penetration of air inside the balloon.

Do not connect any high pressure devices to the port to inflate the balloon as it can result in its bursting.

Do not inflate the balloon with air or any other gas while it is inside the body.

Improper preparation can lead to penetration of air in the system. It can prevent performance of due X-ray visualization. While placing the distal end of the balloon catheter on the guidewire, ensure the distal end of the balloon catheter is not damaged.

Do not tighten the rotating hemostatic valve too hard around the balloon catheter. Too tight fixation of the valve can make the balloon inflation/deflation process slower.

Do not advance the balloon catheter or guidewire in case of resistance. In case of resistance, it is necessary to establish the reason for such resistance by means of X-ray monitoring.

To ensure safety of the patient, inflation and deflation of the balloon should always be performed under x-ray visualization.

After preparation of the balloon and before its use inflate it again up to its nominal volume to inspect the balloon catheter for any defects or damages before use. Do not use the catheter in case of detection of any malfunctions.

Check compatibility of the balloon catheter in case of application of other accessory units, generally used in intravascular interventional procedures. The doctor should specialize in percutaneous intravascular procedures and be aware of potential complications, connected with the procedure.

The balloon catheter has a hydrophilic surface and should be moisturized before use. Do not allow the balloon catheter to dry after moisturizing.

The balloon catheter should be manipulated carefully to avoid accidental damage of the device. Do not expose the balloon catheter face to organic solvents, which can damage the balloon catheter and/or its coating.

Before use, it is essential to ensure the inner diameter of the guidewire or accessory application is compatible with the inner diameter of the balloon catheter, specified on the label.

The balloon catheter should be manipulated with special care when inserted in winding vessels to avoid any damage. In case of detection of unusual resistance, do not advance or remove the device until the reason for such resistance is established.

Presence of calcification, uneven surface or earlier inserted devices can result in damage of the balloon catheter as well as hinder its insertion or removal.

The balloon catheter does not contain latex of PVC materials.

15.6. Legal provisions and conditions of use - no warranty / disclaimer

The use of this device is equated with express and unconditional consent to these legal terms and conditions. If you do not agree with these legal terms and conditions, you may not use this device.

This user manual is not a contractual document and may be changed by BALT without notice. Thus, the user shall not be only informed by BALT about the updates of this manual, but also use the device according to the current manual. This user manual is subject to French laws, any disputes or controversial issues related to this user manual and holding BALT liable will be submitted to the competent judicial authorities of Pontoise (France) for consideration. The user irrevocably agrees with this provision on the exclusive jurisdiction. BALT may not be held liable in the following cases:

- Incorrect use of the device and/or use thereof by untrained users,
- If the device is used over the limits of the indications for use for which it is intended according to the applicable manual,
- In the event of changes in the device without the prior written consent of BALT (including the change of labels and package),
- If the device is used in combination with accessory components, products or materials not approved by BALT,
- In case of repeated use or re-processing of the device,
- If the device seems to be damaged or if the device defect is found when being used (in this case it is necessary to save the device and package and contact BALT or your local distributor, the devices should be returned in their original package upon receipt the consent BALT in advance).

The device is delivered on an "as is" basis. BALT waives all warranties, express or implied, with respect to any device, its compliance with a particular purpose, its quality, its merchantability, etc. Purposes, descriptions or technical specifications specified in BALT's printed documents, including in this user manual, are intended solely for the general description of the device at the time of its release and do not constitute any express warranty. In no event BALT may be held liable based on a claim for contractual liability, non-contractual liability or any other claim for any direct or indirect, incidental or collateral damages incurred by the user or by a third party due to (a) the value of the goods purchased by way of replacement (b) special, indirect, punitive, incidental or collateral damages, or reduction in profit or turnover.

The liability of BALT for defective products will comply with the provisions of the applicable legislation of the European Union with respect to defective products.

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