INDIVIDUAL ACKNOWLEDGMENT

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On this the 25 day of June, 2018, before me,
Danielle Jean Rasmussen the undersigned Notary Public,
Name of Notary Public
personally appeared Upuli Meemaduma
Name(s) of Signer(s)

☒ personally known to me – OR –

☐ proved to me on the basis of satisfactory evidence

to be the person(s) whose name(s) is/are subscribed to the within instrument, and acknowledged to me that he/she/they executed the same for the purposes therein stated.

WITNESS my hand and official seal.

Place Notary Seal and/or Stamp Above

Danielle Jean Rasmussen
Signature of Notary Public
Plymouth, MN

Other Required Information (Printed Name of Notary, Residence, etc.)

OPTIONAL

Although the information in this section is not required by law, it may prove valuable to persons relying on the document and could prevent fraudulent removal and reattachment of this form to another document.

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Title or Type of Document: IFU, Nitrex

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above: _____

Right Thumbprint of Signer

Top of thumb here

INSTRUCTIONS FOR USE

FOR ONE-TIME USE ONLY - DO NOT RESTERILIZE

READ INSTRUCTIONS PRIOR TO USE

CONTENTS OF UNOPENED, UNDAMAGED PACKAGE

ARE STERILE

PRODUCT DESCRIPTION

The NITREX™ Nitinol Guidewire is constructed of a Nitinol core wire with a gold-plated tungsten coil. The wire has a proprietary silicone coating for ease in guidewire placement. A torque device is packaged with 0.014 in. (0.36 mm) through 0.018 in. (0.46 mm) diameter NITREX Nitinol guidewires.

INDICATIONS FOR USE

The 0.014 in. (0.36 mm) and 0.018 in. (0.46 mm) diameter NITREX Nitinol Guidewires are intended for use in the peripheral and coronary vasculature. The 0.025 in. (0.64 mm) and 0.035 in. (0.89 mm) diameter NITREX Nitinol Guidewires are indicated for use in the peripheral vasculature.

CONTRAINDICATIONS

Nitinol Guidewires are not intended for use in the neurovasculature or coronary vasculature.

Do not use substances that are incompatible with any of the product's components (see Product Components).

POTENTIAL COMPLICATIONS

Potential complications associated with guidewires include, but are not limited to:

- Vessel Perforation
- Irritation to vessel causing vessel spasm
- Infection
- Embolization
- Cardiac Tamponade
- Septis

WARNINGS/CAUTIONS

WARNING: This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-wrapping could increase the risk of patient infection and risk of compromised device performance.

WARNING: Handle with care to reduce the possibility of coil separation, uncoiling and breakage.

WARNING: Do not use a wire that has been damaged in any way. DO NOT USE ANY DEVICE IF THE PACKAGE IS OPEN OR DAMAGED.

WARNING: Do not reshape the guidewire by any means. Attempting to reshape the wire may cause damage resulting in embolization of wire fragments into the vessel.

WARNING: Do not advance or withdraw guidewires against resistance until the cause of the resistance has been determined. Excess force against resistance may result in damage to guidewire, catheter or vessel perforation.

WARNING: Do not withdraw NITREX Guidewires through a metal trocar or a metal needle.

WARNING: Maintain fluoroscopic view of the guidewire to ensure the tip rotates freely when torque is applied.

PRODUCT COMPONENTS

NITREX Nitinol Guidewires contain the following components that may come into contact with the skin or blood:

- Guidewire
- Nitinol
- Gold-Plated Tungsten Coil
- Proprietary Silicone Coating
- Epoxy Adhesive
- PTFE Coating

INSTRUCTIONS FOR USE

(Microtome is not included with 0.025 in. (0.64 mm) or 0.035 in. (0.89 mm) diameter NITREX Nitinol Guidewires.)

NOTE: Prior to use, flush the hoop with sterile saline.

1. Remove wire retention clip.

2. Carefully remove the guidewire from the dispenser and inspect for any defects.

3. Introduce the guidewire, flexible end first, into an appropriate access device. (Reference product specific instructions for use for access device.)

4. Advance the wire slowly under continuous fluoroscopic observation. It may be advantageous to rotate the wire during advancement or withdrawal.

5. Directional control of the guidewire tip can be aided by using a microtome pin vise, if desired.

STORAGE CONDITIONS

- Do not store in direct sunlight.
- Keep dry.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

WARRANTY DISCLAIMER

Although the guidewire hereafter referred to as "product" has been manufactured under carefully controlled conditions, ev3 Inc. has no control over the conditions under which this product is used. ev3 Inc., therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ev3 Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind ev3 Inc. to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

ev3 Inc.,
4600 Nathan Lane North
Plymouth, MN 55442, USA
763 398 7000
1 800 333 3333
Uguti Mameduma
Associate RA Specialist

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INDIVIDUAL ACKNOWLEDGMENT

State/Commonwealth of Minnesota }
County of Hennepin } ss.

On this the 25 day of June, 2018


Danielle Jean Rasmussen, the undersigned Notary Public,
Name of Notary Public

Upulil Meemaduma personally appeared
Name(s) of Signer(s)

☒ personally known to me - OR -

☐ proved to me on the basis of satisfactory
evidence

to be the person(s) whose name(s) I am
acknowledged to me that he/she/they expected
the same for the purposes therein stated.
WITNESS my hand and official seal

Place "Notary Seal" and/or Stamp Above

Signature of Notary Public

Plymouth, MN

Other Required Information (Printed Name of Municipality)

OPTIONAL

Although the information in this section is not required by law, it may prove valuable to persons
relying on the document and could prevent fraudulent removal and reattachment of this page
to another document.

Description of Any Attached Document

Title or Type of Document: IFU, Enter

Document Date: _____
Number of Pages: _____

Signer(s) Other Than Named Above: _____

Enteer™

Guidewire

CE
0120

English

DEVICE DESCRIPTION

The Enteer™ Guidewire is available in a nominal diameter of 0.014" (0.36 mm). The Enteer Guidewire has a coiled distal end where the coil is radiopaque to facilitate advancement into and withdrawal from the vasculature using fluoroscopy. The proximal shaft is coated with polytetrafluoroethylene (PTFE) while the distal shaft in the region of the coil is coated with hydrophilic coating. The narrowed distal tip is intended to aid in guidewire control. Refer to the product label for product specifications (e.g. wire length, length of tip radiopacity, tip diameter).

INDICATIONS

The Enteer Guidewire is intended to facilitate placement of balloon dilatation catheters or other intravascular devices during percutaneous transluminal angioplasty (PTA). The Enteer Guidewire is not to be used in cerebral blood vessels.

CONTRAINDICATIONS

None known

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the Enteer Guidewire.
- The Enteer Guidewire is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing or re-sterilizing may increase the risk of patient infection and risk of compromised device performance.
- Do not use the Enteer Guidewire if package is open and sterile barrier is broken.
- Do not push the guidewire if guidewire tip prolapse is observed.
- If guidewire tip prolapse is observed, do not allow tip to remain in a prolapsed position. Prolonged guidewire tip prolapse may cause guidewire damage and/or guidewire tip separation. Determine the cause of the guidewire tip prolapse under fluoroscopy and take necessary remedial action.

PRECAUTIONS

- The Enteer Guidewire should be stored in a cool and dry place.
- Use prior to the "Use By" date.
- Before insertion of the guidewire, administer appropriate anticoagulant and vasodilator therapy.
- The Enteer Guidewire should be handled with care. Prior to use and during the procedure, inspect the packaging and guidewire for bends, kinks, or other damage. Discontinue use if the guidewire becomes damaged.
- Exercise care during a procedure to reduce the possibility of accidental breakage, kinking, bending or coil separation.
- Do not attempt to straighten a guidewire that has been kinked.
- To reduce the potential for guidewire breakage, do not advance a kinked guidewire into a vessel or endovascular catheter.
- Do not rotate the guidewire if significant resistance is felt.
- The Enteer Guidewire should only be manipulated under fluoroscopic observation.

POTENTIAL ADVERSE EFFECTS

Potential adverse effects include, but are not limited to, the following:

- Acute myocardial infarction
- Hemorrhage or hematoma
- Vessel trauma requiring surgical repair or intervention
- Artery spasm
- Embolism
- Stroke
- Neurological deficit
- Infection

SUPPLIES

The following supplies are not provided and need to be available and prepped prior to use of the Enteer Guidewire:

- Guiding catheter
- Angiographic imaging supplies (i.e. radiopaque contrast, manifold, tubing, etc.)
- Sterile heparinized saline
- Other supplies as needed to complete the procedure

PRODUCT PREPARATION

1. Inspect the device package prior to opening. Do not use if package is opened or damaged.
2. Remove the Enteer Guidewire from the sterile packaging.
3. After removing the guidewire from the holder tube, inspect to make sure it is not damaged.
4. Ensure all devices are exposed to sterile heparinized saline.
5. If desired and according to standard practice, gently shape the flexible tip of the guidewire.
6. After pulling the guidewire out of the body, wipe it with gauze soaked in heparinized saline.
7. Upon completion of the use of the guidewire, discard according to hospital standards and national legal requirements.

DIRECTIONS FOR USE

Back-load the guidewire into an endovascular catheter or "wire exchange" technique

1. Prep the Enteer Guidewire per preparation instructions.
2. Insert a guidewire insertion tool per manufacturer's instructions as to communicate with the guidewire port of the endovascular catheter.
3. Carefully insert distal guidewire tip through the insertion tool and into the endovascular catheter.
4. Remove the insertion tool by withdrawing it over the proximal end of the guidewire.
5. Advance the Enteer Guidewire through the endovascular catheter under fluoroscopy to the intended target vessel. Make sure the guidewire moves smoothly and torque the wire if necessary.
6. A guidewire torque device may be used per manufacturer's instructions.

Front-load endovascular catheter onto the guidewire or "wire first" technique

1. Access the vasculature with a guiding catheter and hemostatic Y-adapter per established procedure.
2. Assure the Enteer Guidewire has been prepped per instructions.
3. Advance the guidewire into the vasculature to the target vessel using fluoroscopic guidance. Make sure the wire moves smoothly and torque the wire if necessary.
4. A guidewire torque device may be used per manufacturer's instructions.
5. Once the target location has been accessed, load the distal tip of the endovascular catheter onto the proximal portion of the guidewire.
6. Loosen the knob on the hemostatic Y-adapter.
7. Track the endovascular catheter over the Enteer Guidewire and into the proximal guiding catheter lumen.
8. Tighten the hemostatic Y-adapter to create a seal around the endovascular catheter. Do not over-tighten.
9. Advance the endovascular catheter over the guidewire using fluoroscopic guidance per established procedure.

ev3 Inc.,
4600 Nathan Lane North
Plymouth, MN 55442, USA
Tel: +1 763,398 7000
www.ev3.net

Upul Meemaduria
Associate RA Specialist 25/06/2018

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, ev3 Inc. has no control over the conditions under which this product is used. ev3 Inc. therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ev3 Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind ev3 Inc. to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

DEFINITION OF SYMBOLS

	Manufacturer		Sterilized using ethylene oxide
	Catalogue Number		Batch Code
	Authorized representative in the European community		For prescription use only
	Keep Dry		Keep away from sunlight
	Use By		Do not reuse
	Do not use if package is damaged		Consult instructions for use
	Telephone		Facsimile
	Package Quantity		

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©2015 Covidien.
NITREX is a trademark of a Covidien
Enteer is a trademark of a Covidien company.