MILESTONE
H E L P I N G
P A T I E N T S

To the Federal Service on Surveillance in
Healthcare of the Russian Federation

Diego CORTESI - C.E.O.
Milestone s.r.l.

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

Manual for use
Vacuum Bags for closing of biological specimens
on TissueSAFE and SealSAFE

Thank you for having selected our system and welcome to the ever growing world club of users for Milestone laboratory instrumentation.

We are sure that you will be completely satisfied with this new tool entering your laboratory.

We invite you to read carefully this operator manual and to keep it in reach for convenient and fast consulting. For any possible clarification or any request for assistance please contact either our Representative in your country:

Authorized representative in Russia:

BioLine LLC

197101, Saint-Petersburg,

Pinsky alleyway, 3, lit.A

Phone: +7(812)320 49 49

Hotline: 8 800 333 00 49

main@bioline.ru

Or contact to the manufacturer

Milestone s.r.l.

Via Fatebenefratelli, 1/5

24010 Sorisole (BG) Italy

Tel. +39.035.412 8264

Fax +39.035.575498

web site www.milestonemedsrl.com

e-mail medical@milestonemedsrl.com

Please read the user manual before using the device.

MM – VACUUM BAGS – Manual for use

2018-01

SUMMARY

SUMMARY	4
1. INTRODUCTION	5
1.1. Symbols used in this manual	5
1.2. Designated uses	6
1.3. Technical specifications	7
1.3.1. Stand up vacuum bags	8
1.3.2. Non-Stand up Vacuum bags	9
1.4. Transportation and storage conditions	10
1.5. Compatible reagents for the fixation	10
1.6. Warning information and precautions	10
2. USE THE VACUUM BAGS WITH SEALSAFE	11
2.1. Run a program	11
2.2. Vacuum Only	11
2.3. Vacuum + Fixative	20
2.3.1. Reagent for the fixation management during a "Vacuum + Fixative" program	30
2.3.2. Sealing a vacuum bag more times for the same specimen (single use)	30
2.4. Specimens transportation guidelines	34
3. USE THE VACUUM BAGS WITH TISSUESAFE	35
3.1. Control panel	35
3.2. Select a program	35
3.3. Run a program	36
3.4. Human Organs: Vacuum Selection Guideline	40
3.5. Specimens transportation guidelines	41
4. REPAIR AND SERVICE INFORMATION	42
5. WARRANTY	42
6. SHELF LIFE	42
7. DISPOSAL	42

1.

1.1.

1. INTRODUCTION

1.1. Symbols used in this manual

Please read carefully this instruction.

An instruction accompanied by this symbol provides a cautionary statement: **WARNING**. Failure to follow the instruction may endanger the user or cause damage to the instrument.

Biohazard.

CE logo: This instrument complies with European community directives.

IVD medical device according to 98/79/EC directive.

Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Indicates the manufacturer's batch code so that the batch or lot can be identified.

Catalog number.

Manufacturing date.

Manufacturer.

The number of tests, which are performed using medical devices.

Storage limits of temperature.

Storage humidity limits.

Keep away from sunlight.

1.2. Designated uses

The Vacuum bags (hereinafter "vacuum bags", "plastic bags", "bags") should be used for vacuum packaging of biological specimens in histopathology laboratories for purpose of the following transportation and the storage as follows:

- Preservation and transportation of biological specimens and their supports and/or containers in fresh conditions or in fixative;
- Long term specimen preservation (if required with a defined quantity of fixative into the vacuum bag).

Individual catalog number for one bag/REF	Catalog number of the group package / REF	Name/ type of Vacuum bag	Quantity in package
68047	68047SS*	Vacuum bags SMALL	200 pcs.in package
68055	68055SS*	Vacuum bags MEDIUM	200 pcs.in package
68054	68054SS*	Vacuum bags LARGE	200 pcs.in package
68046	68046SS*	Vacuum bags EXTRA LARGE	200 pcs.in package
68545	68545SS*	Vacuum bags SU* SMALL	200 pcs.in package
68409	68409SS*	Vacuum bags SU* LARGE	200 pcs.in package

* SS - additional marking symbol to the catalog number, indicating the group package, which contains 200 pieces of Bags.

On the individual labeling of the Vacuum bag, the catalog number is indicated without the additional marking symbol SS.

*SU – English abbreviation Stand Up, it means to stand upright.

In-vitro medical device professional use only (European directive 98/79/EC).

Contraindications: not applicable.

Type of samples intended for the preservation and the transportation in vacuum bags: any bio-specimens, human organs and tissues.

The Vacuum bags must be operated exclusively with the reagents listed in section "Compatible reagents for the fixation" 1.5.

The Vacuum bags must be operated exclusively with TissueSAFE, TissueSAFE PLUS and SealSAFE systems produced by Milestone s.r.l.

The Vacuum bags must be operated exclusively as reported in this manual.

Material

Features

- Specimen
- With do
- Applied
- "Single

Each bag
the sealin
the chapter

Milestone
paragraph

Contents

The each p

1.3. Technical specifications

The particular design and the special multi-layer of Vacuum bags, in combination with a dedicated sealing bar (supplied by an electrical intensity specific for thickness and material), allow a perfect double sealing, that resists even in case of accidental fall or external compression. Therefore the use of Vacuum bags is to be considered exclusive to the SealSAFE, TissueSAFE and TissueSAFE PLUS instrumentations. The use of such bags for a different type of instrumentation or the use of vacuum bags not provided by MILESTONE for the above-mentioned instruments preclude the safety of use, resulting in biological risk for the operator and risk of impairment of quality for the specimen, making void the manufacturer's warranty and the CE-IVD.


Material

- BOPA/PE 145 µm (Bioriented Polyamide / Polyethylene)
- Vacuum bags do not contain neither chlorine nor any of its derivatives. Therefore they can be disposed by incineration without dioxine emission.
- High mechanical strength

Features of Bags

- Specimen-carrier bags with double welding for enhanced vacuum
- With dotted lines that suggest where to cut
- Applied sealable pocket for documents
- "Single use only"

Each bag is marked with the symbol "Single use only" (see the circle in figure below); if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient, as explained in the chapter 2.3.2


Milestone S.r.l. provides two types of stand up vacuum bags and four types not stand up. See the following paragraphs for details.

Contents of delivery

The each package is contained 200 pcs of Bags and Manual for Use.

1.3.1. Stand up vacuum bags

Milestone supplies two models differ only in size.

Pic.1 Stand up Vacuum bags

Characteristics:

INDIVIDUAL CATALOG NUMBER	DESCRIPTION
68409	VACUUM BAG SU* LARGE Thickness of the side of Vacuum Bag 145 μ m External dimensions: 330mmx440mm Internal dimensions: 304mm x 430 (365) mm
68545	VACUUM BAG SU* SMALL Thickness of the side of Vacuum Bag 145 μ m External dimensions: 220mmx310mm Internal dimensions: 194mmx300 (280)mm

*SU – English abbreviation Stand Up, it means to stand upright.

Tolerances on the compositions: absolute tolerance for each measure ± 7 mm
(Thicknesses: $\pm 10\%$)

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the specimen from the same patient, as shown in the chapter 2.3.2)
- Due to the special design with the stand up bottom, the Vacuum bag is Stand Up and can stand alone on a table or other horizontal surface
- High mechanical strength
- Applied sealable pocket for documents that can be opened and closed more times
- Continuous lines that suggest where sealing
- Dotted lines that suggest where cutting

- Double barrier layer of material BOPA/PE (polyamide and polyethylene)
- Indication of lot number and expiration.

1.3.2. Non-Stand up Vacuum bags


Four models differ only in size are available:

REF 68047

REF 68055

REF 68054

REF 68046

Pic.2 Not stand up Vacuum bags

Characteristics:

INDIVIDUAL CATALOG NUMBER	DESCRIPTION
68047	VACUUM BAG SMALL Thickness of the side of Vacuum Bag 145 µm Dimensions: External 250mm x 300mm Internal 224mmx170 mm
68055	VACUUM BAG MEDIUM Thickness of the side of Vacuum Bag 145 µm Dimensions: External 330mmx385mm Internal 304mmx195 mm
68054	VACUUM BAG LARGE Thickness of the side of Vacuum Bag 145 µm Dimensions: External 330mmx555mm Internal 304mmx 365 mm
68046	VACUUM BAG EXTRA LARGE Thickness of the side of Vacuum Bag 145 µm Dimensions: External 330mmx555mm Internal 304mmx535mm

Tolerances on the compositions: absolute tolerance for each measure ± 7 mm
(Thicknesses: $\pm 10\%$)

2.**2.1.**

After th

It's poss

• Va

• Va

When s
insertedWhen s
on the clIn additi
possible**2.2. Va**

A program

1. at th
Anat
2. in A
be w

There are tv

1. If two c
without
2. If only
in the p

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient see the chapter 2.3.2)
- High mechanical strength
- Applied sealable pocket for documents that can be opened and closed more times
- Continuous lines that suggest where sealing
- Double barrier layer of material BOPA/PE (polyamide and polyethylene)
- Indication of lot number and expiration.

1.4. Transportation and storage conditions

Temperature: –20°C up to +30°C

Humidity: up to 65%, non–condensing

Avoid exposure to sunlight.

1.5. Compatible reagents for the fixation

1. Formalin, buffered
2. Alcohol-based fixatives

DO NOT USE REAGENTS DIFFERENT FROM THOSE LISTED ABOVE.
In case of doubt please contact: application@milestonemedsrl.com.

Reagents other than those listed above may damage the bag and some components of the instrument.

1.6. Warning information and precautions


In the diagnostic laboratory the work is associated with the direct contact of personnel with various types of biological materials, which can lead to a risk of infection by **infectious agents**. For your safety during the work with vacuum bags, you should to wear gloves and work clothes.

After using Vacuum bags should be treated as biohazard material.

2. USE THE VACUUM BAGS WITH SEALS SAFE

2.1. Run a program

After the user log in, the following screen appears.


It's possible to run two different types of programs by selecting the corresponding icon in the screen:

- **Vacuum Only**
- **Vacuum + Fixative**

When selecting **"Vacuum Only"**, a pre-set vacuum is applied to the vacuum bag where the specimen has been inserted (the vacuum level depends on the type of specimen).

When selecting **"Vacuum + Fixative"**, a defined quantity of fixative is loaded into the vacuum bag depending on the chosen ratio and then a pre-set vacuum is applied.

In addition, you can press the icon **"Favorites"** where the favorite programs can be saved. From this page, it's possible to directly run a program by selecting the corresponding icon.

2.2. Vacuum Only

A program "Vacuum Only" can be used in 2 different cases:

1. at the end of a surgical operation, to put under vacuum a fresh specimen before being sent to the Anatomic Pathology,
2. in Anatomic Pathology, to put under vacuum a specimen to be preserved in stock (this specimen must be well fixed).

When running this program after an organ excision, remove the excess of blood (or any other organic fluid) with a medical cloth. This avoids any putrefactive effect caused by the blood/organic fluid stagnation.


There are two ways of operation to correctly insert a specimen into a bag:

1. If two operators are available: one can keep the bag opened and the second one can insert the specimen without making dirty the external and internal upper part of a bag.
2. If only one operator is available: before inserting a specimen into a bag, reverse the upper part (as shown in the picture below) in a way to do not make dirty the external and internal upper part of a bag.

Follow the instructions below to run a **"Vacuum Only"** program.

- Press the icon **"Vacuum only"** from the main menu and a keyboard named **"Bag ID"** appears, as shown below.


- The barcode data (Bag ID) is automatically entered into the system after being scanned, as shown below. The code of the bag is positioned on its back side, on the left, as shown in the left figure below. Below in the right figure the label containing the barcode of the bag is shown in detail: on the top there is the bag type (for example VACUUM BAG SU LARGE) and below the barcode.

When the barcode is correctly read, the reader emits a confirmation sound.

- When the bag ID is read, automatically another keyboard appears, named **"Case ID"**. Scan the ID of patient as shown below, or enter it manually by the keyboard and press **"Enter"**.


The maximum number of characters allowed for CASE IDs is 50.

as shown

If you enter the barcode of the bag and/or of the patient by the barcode reader (not manually) you don't have to press any button on the screen but the code is automatically read and included in the details of the program.

The following screen appears, showing the Bag ID previously scanned the Case ID (if present) and the default vacuum protocols P1-P4.

own below
elow in the
e bag type

Additional information shown in the screen is the weight measured in the cavity; it should be zero because the bag has not been inserted yet.

When the Bag ID is scanned, automatically a TARE of the weight measured into the cavity is performed.

In case it is necessary to perform the TARE again, it is possible to press the button TARE before inserting the bag into the cavity and then choosing the program of interest.

- Insert the bag into the cavity. The vacuum bag with the specimen must be placed with the open side directed toward the operator.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).

Always position the topmost lines on the bag over the sealing bar.


he ID of th

The base of the bag should rest on the base of the metallic scale; if not, it's required to use the movable metallic tool to keep the base of the bag in position and so as to prevent it from slipping downwards. See the image below. When you don't use it leave it flat on the metallic slope.

- When the bag has been correctly positioned (with the metallic sealing bar between the two white lines) lower the two sidebars that help to keep the bag correctly positioned.

- Above it's shown the procedure to lower the sidebar on the left side of the cavity, repeat the same on the right one.


In the picture below you can see an example of a bag correctly positioned.

It's very important to lower the sidebars before running a process to keep the bag in position and seal it correctly.

- Now you can read the weight of the specimen before starting the vacuum program. For example below is shown a weight of 180g.


- Choose the proper vacuum program from P1 to P4.

The right program selection must be done according to the type of specimen.

Under the icons P1-P4 there is the button "i" (as indicated by the arrow below); if pressed, the list of all the types of specimens to be used for each vacuum program is shown.

Press the icon "i" and follow the indications to choose the vacuum program according to the type of specimen that you have **ONLY** if it is fresh. If the specimen is already fixed and you want to put it under vacuum for archive, it is strongly suggested to use the program P2 for ALL types of specimens.

When the specimens are **FRESH**, you can also refer to the "Vacuum Only" Selection Guideline of the Operator Manual to choose the proper vacuum program. It contains the same information that you get by pressing the info icons.

If you want to run a "Vacuum Only" program of a bag containing a specimen immersed in fixative, do not use the P1 program.

If you want to run a "Vacuum Only" program to put under vacuum a specimen for archive, do not use the P1 program. It is strongly suggested to use the program P2 for ALL types of specimens, as explained above.

The programs P1-P4 differ only in vacuum level; P1 applies the highest level of vacuum and P4 the lowest.

- When you select the vacuum program (for example P1) the following screen appears.

In the lower part of the screen appears the message **"Close cover and keep pressed for few seconds"** that indicates the operator to close the lid and hold it for a few seconds to facilitate the creation of vacuum in the chamber.

- Now the operator has to manually close the cover and to keep it pressed for few seconds, as shown below. When the cover is closed, the vacuum pump will start creating vacuum inside the chamber.

It's very important to keep the cover pressed for few seconds when the vacuum pump starts (right figure above).


The vacuum progression in the chamber is shown on the control panel display under **"Sealing time left"**.


While the vacuum process is running, if you press **"Abort"** (see the arrow in figure above), the program will immediately stopped and the cover will be able to be opened.

While the vacuum process is running, if you press **"Stop & Seal"** (see the arrow in figure above), the program will be immediately stopped, the plastic bag will be sealed and the cover will be able to be opened.

- When the desired vacuum level is achieved, the vacuum bag is automatically sealed.
- On completion of the sealing procedure, an audible notice sounds and a message advises the user that is possible to open the cover.

- Automatically a label showing the process data is printed (if the printer is activated)


If a CASE ID has a number of characters >14, it will be automatically printed vertically, otherwise it will be printed horizontally.

It is possible to choose the details to be printed.


The data printed are (if enabled):

- ✓ PROCESSING (Name of the program)
- ✓ Case ID
- ✓ Bag ID
- ✓ Category ("VACUUM ONLY")
- ✓ Specimen weight (g)
- ✓ User
- ✓ Date
- ✓ Time
- ✓ Outcome

- Open the cover and extract the sealed bag.
- Tear off the label using the jagged edge.

- Remove the protective film on the back of the label, as shown below. There is a cutting line in the centre of the film that should facilitate its removal, so it is suggested to remove first one part of the film, then the second one.

- Attach the label to the bag.


Do not attach the label on the top edge of the bag because it may need to be opened and sealed again. In total it's possible to open the bag and seal it up to 4 times (maximum).

It is required to leave bags and labels away from direct sunlight.

When working with Ethanol based fixatives, it is required to not attach the label to the bag because it may discolor if in contact with that kind of fixative. Please in this case don't remove the protective film and insert the label into the plastic pocket of the bag.

If you need to print the label more times, you can press the icon "Print the last process" on the main menu

- Now the operator can annotate the patient data, by writing on the label of the bag (if desired).
- When running this program at the end of a surgical operation, the operator requesting the histological examination will fill the "Requisition Form". This document needs to be placed in the external pocket on the left of the label where the first annotation is written. At this point, the operator can seal the requisition form in order to avoid any accidental separation from the specimen.

When running this program at the end of a surgical operation, the plastic bag containing a fresh specimen needs to be placed into an appropriate refrigerator under controlled temperature conditions (at 4°C).

At the end of a program, it is required to verify that the sealing is uniform all along the side of the bag.

In case the sealing is not correct:

- Open the bag by using some scissors.
- Clean the upper internal part of the bag with a dry clean cloth. Thoroughly remove any residual liquid inside of the bag must be clean and dry to achieve of complete sealing.
- Repeat the vacuum program.


2.3. Vacuum + Fixative

A program **"Vacuum + Fixative"** can be used in 2 different cases:

1. at the end of a surgical operation, to put under vacuum a fresh specimen with a desired quantity of fixative before being sent to Anatomic Pathology (when they don't want to receive fresh specimens),
2. in the Anatomic Pathology, to put under vacuum a specimen with a desired quantity of fixative to run, for example, the overnight fixation into the bag.

Follow the instructions below to run a **"Vacuum+Fixative"** program.

- Press the icon **"Vacuum+Fixative"** from the main menu and a keyboard named **"Bag ID"** appears, as shown below.



- The barcode data (Bag ID) is automatically entered into the system after being scanned, as shown below. The code of the bag is positioned on its back side, on the left, as shown in the left figure below. Below in the right figure the label containing the barcode of the bag is shown in detail: on the top there is the bag type (for example VACUUM BAG SU LARGE) and below the barcode.

When the barcode is correctly read, the reader emits a confirmation sound.


- If the function **"enable Case ID"** is active, when the bag ID is read, automatically another keyboard appears named **"Case ID"**. Scan the ID of the patient as shown below, or enter it manually by the keyboard and press **"Enter"**.

The maximum number of characters allowed for CASE IDs is 50.

If you enter the barcode of the bag and/or of the patient by the barcode reader (not manually) you don't have to press any button on the screen but the code is automatically read and memorized on the details of the program.

- The following screen appears, showing the Bag ID previously scanned, the Case ID (if present) and the default protocols. The default ratios "weight of specimen: weight of fixative" are: "1:1", "1:2", "1:2.5" and "1:3", as shown below.


If you want to apply a different ratio, you can customize a default program under "Settings" function.

Milestone s.r.l. does not suggest any ratio. An internal protocol in the laboratory/Institute should be validated.

When the Bag ID is scanned, automatically a TARE of the weight measured into the cavity is performed.

In case it is necessary to perform the TARE again, it is possible to press the button TARE before inserting the bag into the cavity and then choosing the program of interest.

The internal part of the bag should be clean and dry to achieve a good sealing. It is also very important to correctly position the bag over the sealing bar, as explained in this section

- Then insert the bag into the cavity. The vacuum bag with the specimen must be placed with the open side directed toward the operator.


When running this program after an organ excision, in case of bleeding organ, remove the excess of blood (or any other organic fluid) with a medical cloth. This to avoid any putrefactive effect caused by the blood/organic fluid stagnation.

Before running a "Vacuum + Fixative" program, it's required to grossly dissect large

organs and to open natural cavities for proper fixative penetration.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).

Always position the topmost lines on the bag over the sealing bar.


The base of the bag should rest on the base of the metallic scale; if not, it's required to use the movable metallic tool to keep the base of the bag in position and so as to prevent it from slipping downwards. See the image below. When you don't use it leave it flat on the metallic slope.

Metallic tool is UP to support the base of the bag.

- When the bag has been correctly positioned (with the metallic sealing bar between the two white lines) lower the two sidebars that help to keep the bag correctly positioned.

- Above it's shown the procedure to lower the sidebar on the left side of the cavity, repeat the same on the right one.


In the picture below you can see an example of a bag correctly positioned.

It's very important to lower the sidebars before running a process to keep the bag in position and seal it correctly.

- Now read the weight of the specimen before to start the program.
- Select the "Vacuum+Fixative" program depending on the desired ratio "weight of specimen: weight of fixative", for example "1:2". The following screen appears.

- In the lower part of the screen appears the message "Arm injection nozzle".

Before arming the injection nozzle, check that the open side of the bag is quite opened. To correctly load the fixative into the bag, **it is important to check that the two sides of the bag are not stuck together and that the bottom side of the bag adheres well to the metallic plate.** See below.

- Now the operator can insert the injection nozzle into the bag.

- When inserted, check again that the open side of the bag is positioned on the sealing bar.

- Now the operator has to manually close the cover and keep it pressed for few seconds.

At the end of a program, it is required to verify that the sealing is uniform all along the side of the bag.

In case the sealing of the bag is not correct:

1. Open the bag by using some scissors.
2. Clean the upper internal part of the bag with a dry clean cloth. Thoroughly remove any residual liquid. The inside of the bag must be clean and dry to achieve of complete sealing.
3. Run immediately a program **"Vacuum Only"** to seal again the bag (leaving the reagent for the fixation inside the bag). Choose the default program P2. See the chapter 2.2 to run a **"Vacuum Only"** program.

In the example reported in this chapter, the ratio "1:2" is selected and it means that if the weight of specimen is 460g, the unit will load into the bag 920g of fixative.


By default each type of bag has a minimum and maximum limit that controls the filling of reagent.

1. **MINIMUM LIMIT:** a minimum quantity of fixative will be loaded into the bag, even if the quantity of fixative should be LOWER (if calculated on the base of the weight of the specimen and the selected ratio); this value avoids cases when the specimen is too little and the fixative could not be enough to cover it.
2. **MAXIMUM LIMIT:** a maximum quantity of fixative will be loaded into the bag, even if the quantity of fixative should be HIGHER (calculated on the base of the weight of the specimen and the selected ratio); this value avoids cases when the fixative to be loaded is too high for the volume of the bag and it could leak out from the bag.

The MINIMUM and the MAXIMUM limit values depend on the type of bag; before running a program when the user scans the BAG ID of the bag, the software will identify which is the type of bag in use and consequently its limits.


If the weight of the reagent to be loaded is less than the MINIMUM quantity established for the bag in use, a message appears advising the user that the minimum quantity will be loaded.

Press "OK" to continue (in the figure below the value "40g" is only an example because the minimum quantity depends on the size of the bag).

The fixative weight declared does not include to quantity of liquid (more or less 60g) loaded into the bag at the end of a procedure, when the fluid injector is cleaned before going out of the bag.

If the sum between the weight of the fixative to be loaded (depending on the ratio chosen) and the weight of the specimen is higher than the MAXIMUM limit defined for the bag, a message warns the user that it is not possible to reach the selected ratio. The new ratio is shown. Press "Yes" to continue or "No" to abort the program. Below an example.


In the table below there are the MINIMUM and MAXIMUM limits for all types of bags:

DESCRIPTION	MINIMUM LIMIT (Fixative weight)	MAXIMUM LIMIT (Fixative weight + Specimen weight)
VACUUM BAG SU LARGE	150g	3000g
VACUUM BAG SU SMALL	40g	300g
VACUUM BAG EXTRA LARGE	150g	3000g
VACUUM BAG LARGE	100g	2000g
VACUUM BAG MEDIUM	40g	150g
VACUUM BAG SMALL	40g	40g

If the weight of the reagent for the fixation to be loaded is higher than the MAXIMUM quantity established for the bag in use and the maximum ratio that is possible to reach is less than 1:0.2, a message appears advising the user that the specimen is too heavy and that it's not possible to continue.

In case the user doesn't want to run a process by using a ratio lower than the selected or when the specimen is too heavy to continue, it's required to manage the specimen in a different way, as:


- using a bigger bag (if available),
- splitting the specimen into 2 bags,
- inserting the specimen in an external plastic container to start fixation.

When the cover is closed, the vacuum pump will start creating vacuum inside the chamber.

The vacuum progression in the chamber is shown on the control panel display under **"Sealing time left"**.


It's very important to keep the cover pressed for few seconds when the vacuum pump starts (see below).

What happens during a "Vacuum + Fixative" program, when the cavity is closed:

- a fixed level of vacuum is applied, that is necessary to load the reagent for the fixation into the bag,
- after a few seconds the vacuum pump stops, some air is pumped into the bag to open it a bit and the reagent for the fixation is loaded into the bag,
- during the loading step, the "Sealing time left" bar disappears and both the Fixative weight and message "Loading reagent, do not open the cover" are shown (as shown below),

the loading stops when the desired RATIO "weight of specimen: weight of fixative" is achieved (in example below the specimen weight is 460g, the ratio is 1:2 and so the loading stops when the fixative weight is 920g),

- for a few seconds the fluid injector sprinkles the remaining reagent into the bag,
- then it comes back automatically to the initial position (out of the bag),
- now the vacuum pump restarts to reach a fixed level of vacuum.


It is possible

The data p

- PRO
- Case
- Bag I
- Categ
- Spec
- Fixati
- Spec
- the M
- Reag
- User
- Date
- Time
- Outco

- Open th

- Tear off

If you press **"Abort"** (see the arrow in figure above), the program will be immediately stopped and the cover will be able to be opened.

If you press **"Stop & Seal"** (see the arrow in figure above), the program will be immediately stopped, the plastic bag will be sealed and the cover will be able to be opened.

- finally the bag is automatically sealed.
- On completion of the sealing procedure, an audible notice sounds and a message advises the user that it's possible to open the cover.
- Automatically a label showing the process data is printed (if the printer is activated).


If a CASE ID has a number of characters >14, it will be automatically printed vertically, otherwise it will be printed horizontally.

It is possible to choose the details to be printed.

The data printed are (if enabled):

- PROCESSING
- Case ID
- Bag ID
- Category ("VACUUM+FIXATIVE")
- Specimen weight (g)
- Fixative weight (g)
- Spec/Fix ratio (Specimen weight / Fixative weight ratio. If the selected ratio cannot be reached because of the MIN or the MAX limits, the "Set ratio" and the "Effective ratio" will be shown)
- Reagent for the fixation
- User
- Date
- Time
- Outcome

- Open the cover and extract the sealed bag.
- Tear off the label using the jagged edge.

- Remove the protective film on the back of the label, as shown below. There is a cutting line in the center of the film that should facilitate its removal, so it is suggested to remove first one part of the film, then the second one.

- Attach the label to the bag.

Do not attach the label on the top edge of the bag because it may need to be opened and sealed again. In total it's possible to open the bag and seal it up to 4 times (maximum).


In case the
Only” progr
See the cha


It is required to leave bags and labels away from direct sunlight.


When working with Ethanol based fixatives, it is required to not attach the label to the bag because it may discolor if in contact with that kind of fixative. Please in this case don't remove the protective film and insert the label into the plastic pocket of the bag.

If you need to print the label more times, you can press the icon "Print the last process" on the main menu.

- Open the cover and extract the sealed bag.

At the end of a program, it is required to verify that the sealing is uniform all along the side of the bag.

In case the sealing of the bag has not been performed correctly, open the bag and run immediately a "**Vacuum Only**" program to seal it again (leave the specimen and the formalin inside). Choose the default program P2. See the chapter 2.2 to get instructions on how to run a "Vacuum Only" program.


the bag
se don't

2.3.1. Reagent for the fixation management during a “Vacuum + Fixative” program


The reagent for the fixation is always loaded from the oldest tank (in the reagent page when a tank is set as full, the loading date is saved).

During a process, if a tank becomes empty before the completion of the loading step, the unit will try to load reagent for 30 seconds. On completion of this time, when the fixative weight doesn't increase for 30 seconds, if the loading is not completed (because the tank is empty):

- **when another FULL tank exists**, the instrument starts automatically to load reagent for the fixation from the next oldest tank available and the old tank is set up as empty;
- **when all tanks are empty**, an audio alarm and message appears advising the user that he can enter in the reagent page, set up minimum one tank as FULL and after that the loading step will be automatically completed by the unit. See the message below.

Press “Yes” to enter directly to the Reagent screen. Set up as FULL at least one tank. Then press the arrow (circle below) to exit from the reagent screen and the loading procedure will restart immediately.

2.3.2. Sealing a vacuum bag more times for the same specimen (single use)

If you have a TissueSAFE, TissueSAFE PLUS or SealSAFE next to the operating room and a SealSAFE Anatomic Pathology, the SAME bag for the SAME specimen can be used from the operation room to the storage.

The software of instruments allows for sealing each bag up to 4 times.

Below an example on how to use a bag for all life of a specimen is reported:

- Next to the Surgery room, the same bag can be sealed the first time by a program “Vacuum only” if SealSAFE is present or by a vacuum program if a TissueSAFE is present; in this way, the specimen arrives to the pathology lab in fresh conditions.


Wa
con

- Then th
- The par
(by a pr
sure.

set as

to load
onds, if

on from


enter in
natically

Then in Anatomic Pathology, if another SealSAFE is present, the same bag can be opened and sealed again by a program "Vacuum+Fixative"; in this way an overnight fixation of the specimen can be executed (as you can see in the picture the bag filled by fixative is stand up).

the arrow

Finally the bag can be emptied of fixative before extraction of the specimen.
See below how to empty the fixative before opening the bag. First cut only a lateral side of the bag to empty the fixative and then open completely the bag.

SealSAFE in
room to them only" if a
e specimen

Waste reagent for the fixation should be disposed in an approved incinerator in compliance with all federal, provincial and local government regulations.


- Then the specimen can be extracted from the bag and a biopsy is taken to be processed.
- The part of the specimen not to be processed can be sealed for the third time to archive, without reagent (by a program "Vacuum Only") or with new fixative (by a program "Vacuum+Fixative") if the fixation is not sure.

In the image below the archiving without reagent for the fixation.

If the fixation is not sure, it's required to execute the archiving with fixative.

The top edge of the bags contains helpful indications for areas to seal and cut. This process can occur up to three times for sealing and twice for cutting. See below.


As shown in figure above, the lines where to cut are dashed.

In case of incorrect sealing the bag can be sealed up to 4 times (maximum).

VACUUM BAG SU LARGE:

When the bag has been sealed and cut 2 times, before to insert the bag in the cavity to seal it for the third time it's required to use the tool of the metal plate to keep the bag in position. See below. On the metal plate there are numbered positions etched on the metal plate: you have to set the tool at the position number 6.

It's important to use the bar of the metallic plate to keep the bag correctly positioned over the sealing bar because it's shorter after two cuts.

For the first and second sealing, do not use the bar: leave it flat on the metallic plate, as shown below.

VACUUM
When using
position.


You have
the second

Other type
evaluate if

VACUUM BAG SU SMALL:

When using the vacuum bag small size it's always required to use the tool of the metal plate to keep the bag in position. See below.


You have to set the tool at the position number 5 before running the first sealing, at the position number 4 at the second sealing and at the position number 3 at the third sealing.

Other types of bags (the not stand up), it's required to check their size before running a program and to evaluate if using or not using the metallic tool is needed.

the sealing

ow.

2.4. Specimens transportation guidelines

If the time required by the courier for the specimens' transportation to the Anatomic Pathology is:

- **less than 30 minutes:** put the specimens in the Thermal bag with the cool packs at +4°C,
- **less than 1 hour:** put the specimens into the Refrigerator, close to the SealSAFE, for at least 2 Hr (pre-cooling procedure) and later place them in the Thermal bag with the cool packs at +4°C,
- **more than 3 hours:** use of the Mobile Refrigerator connected to an external power supply (12Vdc/115/230V~).

For details on Refrigerators models available from Milestone s.r.l., or contact directly the distributor or the local support.

With the system at controlled temperature designed by Milestone s.r.l. is ensured a perfect preservation of the morphology, antigenicity, DNA/RNA and proteins are preserved up to 48/72 hours in an optimal manner (depending on the initial conditions and the type of sample).

In none of the three cases described above the total time of the transport must exceed 72 hours (including the time spent in the Surgery Room/Refrigerator of the Laboratory of Histology). For proper preservation of the morphology, the samples must be packed and stored at 4°C ± 2°C.

3.

3.1.

W
Th
inv

Cor

1.
2.
3.
4.
5.
6.

3.2.

Whe
displ

Tiss

To m
+ (P
- (M

3. USE THE VACUUM BAGS WITH TISSUESAFE


Specimens under vacuum have to be placed and maintained at a temperature of +4°C.

3.1. Control panel

When the TissueSAFE unit is switched on, the control panel will appear as shown below.

The operator can access the functionality of the unit through the panel, which controls all the parameters involved in the bio-specimens preservation process.

Control panel description:

1. Forced sealing
2. Interrupt a vacuum process
3. Display
4. Programs selection
5. Minus
6. Plus

3.2. Select a program

When the unit is switched on, the operator can run a default protocol (from P01 to P04) appearing on the display (3) of the control panel.

TissueSAFE is able to memorize up to 30 different protocols.

To move from P01 to another protocol, press the button:


- + (Plus) to move forward,
- (Minus) to go back to previous program.


The right program selection must be done according to the Vacuum Selection Guide (see 3.4).

3.3. Run a program

After the user log in, the following screen appears:

In the above picture the display is showing P01. It means that the unit is ready to start the vacuum program n.1. After the excision of the organ by the surgeon, it should be placed in the appropriate (by dimension) Milestone plastic bag. The operator can use the designed shelf (if available) on the right of the vacuum cavity (left figure below).

The vacuum bag with the specimen has to be placed with the open side directed toward the operator, exactly as shown in the picture above.


Make sure of using plastic bags of the appropriate size according to the specimen dimensions. Use of inappropriate plastic bag size (too small) may cause spill out of the specimen and a non correct vacuum phase and/or sealing.

After the organ excision, in case of bleeding organ, remove the excess of blood (or any other organic fluid) with a medical cloth. This to avoid any putrefactive effect caused by the blood/organic fluid stagnation.


The plastic bag must be placed with the metallic sealing bar between the two white lines (see picture below).


If the median axis of the sample does not fit with the height of the sealing bar, remove one or both of the spacers from the cavity (see pictures below).


The operator has to manually close the cover, as shown below. When the cover is closed, the vacuum pump will start creating the vacuum inside the chamber.


The vacuum progression inside the chamber is shown on the control panel display.

After the vacuum step, the vacuum bag is sealed.

On completion of the sealing procedure, an audible alarm will sound, the following screen will appear and the cover will automatically open.


Now the operator can annotate the patient data, by writing on the label of the plastic bag.

At the end of the surgical operation, the operator requesting the histological examination will fill the "Requisition Form". This document needs to be placed in the external pocket on the right of the label where the first annotation is written.

At this point, the operator can seal the requisition form in order to avoid any accidental separation from the specimen.

The plastic bag containing the specimen now needs to be placed into an appropriate refrigerator under controlled temperature conditions (at 4°C).

3.4. Human Organs: Vacuum Selection Guideline

<i>P1</i>	<i>P2</i>	<i>P3</i>	<i>P4</i>
Kidney	Lung	Brain	Cyst > 10cm in diameter
Breast	Placenta	Eye	
Uterus	Prostate	Cyst <10cm in diameter	
Liver	Gall bladder		
Muscle	Bowel		
Thyroid	Testis		
Thymus	Appendix		
Skin			
Tongue			
Spleen			
Colon			
Bone			
Pancreas			
Bladder			

This guide is intended for specimens directly packaged with vacuum sealed plastic bags.

When sealing very soft specimens (For ex: ovarian cyst and lymph nodes), we suggest the use of plastic containers with perforated lids. In this case the specimen dimension must not exceed 2/3 of the container inner volume.

When using plastic containers with perforated lids, choose the protocol depending on the type of organ following the table above.

3.5.

SE

SE

SE

If the t
Depar

1. les

2. les


(pre

3. mo

115

3.5. Specimens transportation guidelines

If the time required by the courier for the specimens' transportation to the Anatomical Pathology Department is:

- less than 30 minutes:** put the specimens in the Thermal bag with the cool packs at +4°C.
- less than 1 hour:** put the specimens into the refrigerator, close to the TissueSAFE, for at least 2 Hr (pre-cooling procedure) and later place them in the Thermal bag with the cool packs at +4°C,
- more than 3 hours:** use of the Mobile refrigerator connected to an external power supply (12/24Vdc, 115/230V~).

With the system at controlled temperature designed by Milestone s.r.l. is ensured a perfect preservation of the morphology, antigenicity, DNA/RNA and proteins are preserved up to 48/72 hours in an optimal manner.

In none of the three cases described above the total time of the transport must exceed 72 hours (including the time spent in the Operating Room/Refrigerator of the Laboratory of Histology). For proper preservation of the morphology, the samples must be packed and stored at $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

4. REPAIR AND SERVICE INFORMATION

The medical product is not a subject of technical service (maintenance) and repair.
To receive additional service or support, please contact with Customer's Support or your local distributor.

Customer's Support

Tel.: +39 035 572735

E-mail: customersupport@milestonemedsrl.com

Authorized representative in Russia:

BioLine LLC

Pinsky alleyway, 3, lit.A
197101, Saint-Petersburg,
Russia

E-mail: service@bioline.ru

Phone: +7(812)320 49 49

Fax: +7(812) 320 49 40

Hotline: 8 800 333 00 49

In case of any undesired reactions or events during the use of the product, please immediately contact with the Customer's Support or your local distributor.

5. WARRANTY

Milestone s.r.l. guarantees that the product has been subjected to a comprehensive quality control procedure based on the our in-house testing standards, and that the product is faultless and complies with all technical specifications and/or agreed characteristics warranted.

The manufacturer guarantees high quality of product operation during all term until expiry date under condition of storage and methodic recommendations observance.

6. SHELF LIFE

The shelf life of Vacuum Bags is 10 years from the manufacturing date when stored as recommended. The manufacturing date is indicated on the group package.

7. DISPOSAL

Vacuum bags do not contain neither chlorine nor any of its derivatives. Therefore they can be disposed by incineration without dioxine emission.

The Vacuum bags supplied by Milestone s.r.l. should be treated as biohazard material after their use.

Do not use Vacuum Bags beyond their expiry date indicated on the package. Waste material should be disposed in an approved incinerator in compliance with all federal, provincial and local government regulations. For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by "СанПиН 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class "Б" (epidemiologically hazardous waste).

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