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1. Product purpose.

Tests to determine hormones in urine "FRAUTEST, versions: Test to determine pregnancy FRAUTEST express; Test to determine pregnancy FRAUTEST double control; Test to determine pregnancy FRAUTEST expert; Test to determine pregnancy FRAUTEST comfort; Test to determine pregnancy FRAUTEST exclusive; Test to determine ovulation FRAUTEST ovulation; Set of tests to determine ovulation and pregnancy FRAUTEST planning" (hereinafter - Tests) are intended for the qualitative determination of the quantitative content of human chorionic gonadotropin (HCG) and luteinizing hormone (LH) in urine samples.

2. Scope.

Self-test.

3. General indications for use.

Tests are used for early diagnosis of pregnancy and the calculation of expected date of ovulation.

4. Conditions of use.

Non-invasive disposable medical device for self-diagnosis at home and clinical settings.

5. The main technical characteristics.

Tests for self-diagnosis are non-sterile single-use product. Tests represent – strip, cartridge and stream test strip which consists of glass fiber strips with monoclonal antibody against the luteinizing hormone (LH) or chorionic gonadotropin (HCG), solid anti-mouse immunoglobine antibody of nitrocellulose membrane, colloidal gold devices- monoclonal antibody against LH or HCG.

Tests should not have external defects.

Technical specifications shall comply with Addendum 1 of this document.

The general view of the test is presented in Figures 1, 2, 3.

Figure 1.


Figure 2.


Figure 3.

6. Used materials.

Used materials shall comply with Annex 1 of this document.

7. Contraindications.

None.

8. Side effects

None.

9. The principle of operation.

Pregnancy tests FRAUTEST are independent immunoassay devices which are intended for the qualitative detection of human chorionic gonadotropin (HCG) in urine for early detection of pregnancy. Human chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. During normal pregnancy HCG may be detected in urine in 7 days after conception and doubles every 2 days. During the passage of the last menstrual period the level of HCG in the urine is approximately 100 mIU /ml with peak levels of 100 000-200 000 mIU /mL at the end of the first trimester. The presence of HCG soon after conception and its subsequent increase in concentration during early gestational growth makes it an ideal marker for the early detection of pregnancy.

When this method is used a unique combination of colored and solid phase monoclonal antibody composition for selective determination of HCG in the sample with extremely high sensitivity is applied. Less than in 5 minutes HCG may be detected at such a low level as 15mIU/ml.

As soon as the test sample passes through the absorbent device the labeled antibody-dye conjugate binds to HCG to form an antibody-antigen complex. This complex binds to antibodies to HCG positive reaction zone («T» zone) and produces a pinkish purple bands at the concentration of HCG more than 15mIU/ml. In the absence of HCG in the zone of positive response no line is available. Unbound conjugate binds the reagents in the control zone (zone "C"), in this case, pinkish-purple streak that demonstrates that the reagents are functioning properly shall appear.

Ovulation tests FRAUTEST are independent immunoassay for the qualitative detection of luteinizing hormone (LH) in urine to determine ovulation.

Luteinizing hormone (LH) is a glycoprotein hormone that is released by basophil cells of the anterior pituitary, the main function is to stimulate ovulation, the release of luteinizing hormone. All healthy women ovulate regularly, every month, so there is the problem of achieving control of pregnancy by ovulation prediction, especially for infertile women who hope that the artificial insemination in the near future will be paid more attention.

This method uses a unique combination of monoclonal dye complex compounds and polyclonal solid phase antibodies to selective identification of the LH in the control sample with an extremely high level of sensitivity. Less than within 10 minutes the LH level, which is at a low level of 25 mIU / mL may be determined.

As soon as the test sample passes through the absorbent device the labeled antibody-dye of conjugate binds to LH to form an antibody-antigen complex. In the event of presence of LH in the sample, on the reading portion of the test results appears purplish-pink strip by an immunological method of coloring with gold special composition with antibodies LH. On the other hand, the strip of light color will always appear in the control area. The control strip is a reference to the color intensity of 25 mIU / mL LH. If the intensity of the test strip is equal to or greater than the control / reference strip the test result is positive, indicating that the surge of LH is probably still growing.

10. Preparation for the application.

Before the application prepare the container for urine collection.

11. Method of application and evaluation of the results.

In accordance with Annex 2.

12. Precautions

Pregnancy tests FRAUTEST:

- Do not use after expiry date.
- Do not use the test to retest. Discard the test after a single application.
- Do not use the test if the individual packaging is damaged or opened.
- Do not touch the membrane area.
- After opening the foil pouch the test shall be used immediately. Prolonged exposure to ambient humidity will damage the product.

- Carefully read the instructions for use before carrying out the test.
- The final result of the diagnosis shall not be based on the results of a single test, but should be determined only by a physician after evaluating all clinical and laboratory findings.

Ovulation Test FRAUTEST:

- Do not take medications of HCG or LH which may affect the test result, moreover, during the pregnancy, menopause or using contraceptives the result will not be normal.
- Do not use after expiry date.
- Do not use the test to retest. Discard the test after a single application.
- Do not use the test if the individual packaging is damaged or opened.
- Do not touch the area of the membrane.
- After opening the foil pouch the test shall be used immediately. Prolonged exposure to ambient humidity will damage the product.
- Carefully read the instructions for use before carrying out the test.
- The final result of the diagnosis shall not be based on the results of a single test but should be determined only by a physician after evaluating all clinical and laboratory findings.

13. Marking.

The marking is applied directly on the individual and consumer packing. The marking is applied by any means (printing, stamping, die), allowing its brightness, clarity and readability.

Marking of individual packing:

- reference designation "HCG", "Ovulation"
- number of batch «LOT»:
- valid until "⌚:" or «EXP:».

Marking of consumer packing:

- name of the producer (manufacturer), the name of the country of origin, address:
- name of the importing company/authorized organization, address:
- product name, trademark:
- date of the product expiration:
- valid until «⌚:»:
- batch number «LOT: ...»:
- disposal instructions: "Do not throw down the drain":
- method of storage:
- product description:
- product barcode:
- information about the certification (if any):
- registration number.

Transport marking:

- name of the producer (manufacturer), the name of the country of origin, address:
- name of the importing company/authorized organization, address:
- product name, trademark:
- quantity:
- shelf life:
- batch number «LOT: ...»:
- method of storage:
- product barcode.

14. Storage conditions and shelf life.

Shelf life: 3 years from the date of manufacture. Do not use after expiry date. Store in sealed containers at temperatures between + 4°C to + 30°C in a dry place, protected from direct sunlight and out of reach of children.

15. Terms and conditions of transportation.

It is transported by any mode of transport.

It is transported at a temperature from -50°C to + 50°C in a dry place, protected from direct sunlight.

16. Maintenance and repair.

The product is disposable and does not require maintenance and carrying out routine repairs.

17. Warranty.

The manufacturer is not responsible for any consequences resulting from the improper use. Warranty is 3 years.

18. The method of disposal.

It is forbidden to lower the test in the toilet. The used test shall be disposed as a household waste.

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20. Authorised representative in Russia.

"BOLEAR" Ltd, jur. address: Russia, 115230, Moscow, Kashirskoye sh. 12, bldg. 1, act. address: 115201, Moscow, Kashirskiy passage, 13, tel.: +7 (495) 787-92-68, www.bolear.ru.

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I hereby certify the above signature of Wolfgang Hypko, residing at Am Jahnplatz 5 68642 Bürstadt, born on 09.03.1943, personally known to me.

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- In the Companies' Register of the Local court of Darmstadt the Axiom Gesellschaft für Diagnostica und Biochemica mbH with its headquarters located in Bürstadt is registered under registration number HRB 61822*
- In the Companies' Register, the above mentioned Wolfgang Hypko is registered as director.*
- According to the register; he has the power to represent the company individually and he is not bound by the restrictions of article 181 BGB (German Civil Code), which limits a representative's power to enter into contracts with himself or with himself representing also a third party at the same time.*

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
“06” June 2016

ANNEX 2 TO MAINTENANCE DOCUMENTATION

Tests to determine hormones in urine FRAUTEST, versions: Test to determine pregnancy FRAUTEST express; Test to determine pregnancy FRAUTEST double control; Test to determine pregnancy FRAUTEST expert; Test to determine pregnancy FRAUTEST comfort; Test to determine pregnancy FRAUTEST exclusive; Test to determine ovulation FRAUTEST ovulation; Set of tests to determine ovulation and pregnancy FRAUTEST planning.


Produced by company AXIOM Gesellschaft für Diagnostica und Biochemica mbH, Am Jahnplatz 5, 68642 Bürstadt, Germany.

Test to determine pregnancy FRAUTEST express

Directions for use.

1. Open the individual package and remove the test device, while holding its colored tip;
2. Immerse the test device vertically into the container with urine so that the tip of the arrow points toward urine, the level of urine should not exceed the MAX line.
3. Remove the test device from urine in 15 seconds.
4. Place the test device on a flat surface, assess the result after 3-5 minutes, but not later than within minutes.
5. Discard the test device after use.

Assessment of results.

- a - absorbent cotton
- b - test line
- c - control line

Negative result

Positive result

Invalid

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test was carried out incorrectly, if the purple-pink line does not appear in the region of the test and control line. This happens very rarely. The reason for this may be failure to comply with the instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:


- The test was carried out before the term of positive result obtainment (in particular, this result may be observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG at normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

Test to determine pregnancy FRAUTEST double control

Directions for use.

1. Open the individual package and remove the test device, while holding its colored tip;
2. Immerse the test device vertically into the container with urine so that the tip of the arrow points toward urine, the level of urine should not exceed the MAX line.
3. Remove the test device from urine in 15 seconds.
4. Place the test device on a flat surface, assess the result after 3-5 minutes, but not later than within 10 minutes.
5. Discard the test device after use.

Assessment of results.

- a - absorbent cotton
b - test line
c - control line

Negative result

Positive result

Invalid

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test was carried out incorrectly, if the purple-pink line does not appear in the region of the test and control line. This happens very rarely. The reason for this may be failure to comply with the instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Pregnil, Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:


- The test was carried out before the term of positive result obtainment (in particular, this result may be observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG within normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

Test to determine pregnancy FRAUTEST expert

Directions for use.

1. Open the individual package and remove the test device.
2. Draw urine out of the container for urine collection using a dropper.
3. Transfer urine to the specimen well using a disposable dropper.
4. Place the test device on a flat surface with its window facing upward, assess the result after 3-5 minutes, but not later than within 10 minutes.
5. Discard the test device after use.

Assessment of results.

a – specimen well
b – test line
c – control line

Negative result

Positive result

Invalid

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test has been carried out incorrectly, if the purple-pink line does not appear in the region of test and control line. This happens very rarely. The reason for this may be failure to comply with instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:


- The test was carried out before the term of positive result obtainment (in particular, this result is observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG at normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

Test to determine pregnancy FRAUTEST comfort

Directions for use.

1. Open the individual package and remove the test device.
2. Remove the cap.
3. The test device should be held with its absorbent tip facing downward. Place the absorbent tip under the urine stream, until it gets completely wet
4. Place the test device on a flat surface with its window facing upward, assess the result after 3-5 minutes, but not later than within 10 minutes.
5. Put on the cap.
6. Discard the test device after use.

Assessment of results.

- a - cap
b - absorbent tip
c - test line
d - control line

Negative result

Positive result

Invalid

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test has been carried out incorrectly, if the purple-pink line does not appear in the region of the test and control line. This happens very rarely. The reason for this may be failure to comply with the instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Pregnil, Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:


- The test was carried out before the term of positive result obtainment (in particular, this result may be observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG within normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

Test to determine pregnancy FRAUTEST exclusive

Directions for use.

1. Open the individual package and remove the test device.
2. Remove the cap.
3. The test device should be held with its absorbent tip facing downward. Place the absorbent tip in the urine stream, until it gets completely wet.
4. Place the test device on a flat surface with its window facing upward, assess the result after 3-5 minutes, but not later than within 10 minutes.
5. Put on the cap.
6. Discard the test device after use.

Assessment of results.

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test has been carried out incorrectly, if the purple-pink line does not appear in the region test and control line. This happens very rarely. The reason for this may be failure to comply with instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:

- The test was carried out before the term of positive result obtainment (in particular, this result is observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG in normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

Test to determine ovulation FRAUTEST ovulation

When the testing should be performed.


In order to determine the day of testing initiation, women of reproductive age should deduct 17 days from the number of days in their menstrual period. Most women have a regular period that continues for 28 days. In case of irregular period, choose the shortest period for counting. For example, if your period continues for 28 days, testing should be initiated at Day 11 after the start of the period. Use this table to determine the day of testing initiation. In order to determine the day of increase in LH production, you should carry out urine analysis approximately at the same time every day. Generally, in most women with a regular menstrual period, 5 days of testing is enough. Two most favorable days occur upon sharp rise in production of the hormone needed for ovulation. The highest probability of pregnancy occurs with sexual relations in favorable days, namely, at the days of ovulation.

Menstrual cycle	Day for test initiation
21-22 days	Day 6
23-24 days	Day 7
25 days	Day 8
26 days	Day 9
27 days	Day 10
28 days	Day 11
29 days	Day 12
30 days	Day 13
31 days	Day 14
32 days	Day 15
33 days	Day 16
34 days	Day 17
35 days	Day 18
36 days	Day 19
37 days	Day 20
38 days	Day 21

Directions for use.

1. Open the individual package and remove the test device, while holding its colored tip;
2. Immerse the test device vertically into the container with urine attached so that the tip of the arrow points toward urine, the level of urine should not exceed the MAX line.
3. Remove the test device from urine in 10 seconds and place on a dry flat surface.
4. Wait for colored lines to appear. Depending on LH content in urine, the result can be observed at least in 40 seconds. However, result confirmation requires full testing time of 10 minutes. Do not interpret the result after 30 minutes.
5. Discard the test device after use.

Assessment of results.

- a - absorbent cotton
- b - test line
- c - control line

Negative result

Positive result

Invalid

Negative result: only one purple-pink line appears in the result control region, and the line in the test region is absent, or the test line is lighter than the control line. If no increase in luteinizing hormone content is detected, you need to continue daily testing (once a day) until detection of increase in hormone content.

Positive result: two purple-pink lines appear, and the color of the test line is the same or more intense than the color of the control line. This result indicates the level of LH surge and most probably your ovulation will occur within 24-48 hours. The greatest chance to get pregnant occurs with sexual relations during these 2 favorable days. Once you have determined that the LH content is increased, there is no need to continue testing.

Invalid: there are no lines in the control region, the result is that the test has been carried out incorrectly. This happens very rarely. Repeat the test procedure with a new test device, while strictly following requirements of this instruction.

The test is not intended for prevention of pregnancy!

Set of tests to determine ovulation and pregnancy FRAUTEST planning

Test to determine ovulation

When the testing should be performed.

In order to determine the day of testing initiation, women of reproductive age should deduct 17 days from the number of days in their menstrual period. Most women have a regular period that continues for 28 days. In case of irregular period, choose the shortest period for counting.

For example, if your period continues for 28 days, testing should be initiated at Day 11 after the start of the period. Use this table to determine the day of testing initiation.

In order to determine the day of increase in LH production, you should carry out urine analysis approximately at the same time every day. Generally, in most women with a regular menstrual period, 5 days of testing is enough.


Two most favorable days occur upon sharp rise in production of the hormone needed for ovulation. The highest probability of pregnancy occurs with sexual relations in favorable days, namely, at the days of ovulation.

Menstrual cycle	Day for test initiation
21-22 days	Day 6
23-24 days	Day 7
25 days	Day 8
26 days	Day 9
27 days	Day 10
28 days	Day 11
29 days	Day 12
30 days	Day 13
31 days	Day 14
32 days	Day 15
33 days	Day 16
34 days	Day 17
35 days	Day 18
36 days	Day 19
37 days	Day 20
38 days	Day 21

Directions for use.

1. Open the individual package and remove the test device, while holding its colored tip;
2. Immerse the test device vertically into the container with urine attached so that the tip of the arrow points toward urine, the level of urine should not exceed the MAX line.
3. Remove the test device from urine in 10 seconds and place on a dry flat surface.
4. Wait for colored lines to appear. Depending on LH content in urine, the result can be observed at least in 40 seconds. However, result confirmation requires full testing time of 10 minutes. Do not interpret the result after 30 minutes.
5. Discard the test device after use.

Assessment of results.


- a - absorbent cotton
- b - test line
- c - control line

Negative result

Positive result

Invalid

Negative result: only one purple-pink line appears in the result control region, and the line in the test region is absent, or the test line is lighter than the control line. If no increase in luteinizing hormone content is detected, you need to continue daily testing (once a day) until detection of increase in hormone content.

Positive result: two purple-pink lines appear, and the color of the test line is the same or more intense than the color of the control line. This result indicates the level of LH surge and most probably your ovulation will occur within 24-48 hours. The greatest chance to get pregnant occurs with sexual relations during these 2 favorable days. Once you have determined that the LH content is increased, there is no need to continue testing.

Invalid: there are no lines in the control region, the result is that the test has been carried out incorrectly. This happens very rarely. Repeat the test procedure with a new test device, while strictly following requirements of this instruction.


The test is not intended for prevention of pregnancy!

Test to determine pregnancy

Directions for use.

1. Open the individual package and remove the test device, while holding its colored tip;
2. Immerse the test device vertically into the container with urine so that the tip of the arrow points toward urine, the level of urine should not exceed the MAX line.
3. Remove the test device from urine in 15 seconds.
4. Place the test device on a flat surface, assess the result after 3-5 minutes, but not later than within 10 minutes.
5. Discard the test device after use.

Assessment of results.

- a - absorbent cotton
- b - test line
- c - control line

Negative result

Positive result

Invalid

Negative result: You are not pregnant, if only one pink line appears in the control region.

Positive result: You are pregnant, if two pink lines appear: the control and test line.

Invalid: The test has been carried out incorrectly, if the purple-pink line does not appear in the region of the test and control line. This happens very rarely. The reason for this may be failure to comply with instructions, for instance:

- Insufficient amount of urine on the absorbent cotton;
- During the test, the test device was held horizontally rather than vertically.

During the test, it is possible to obtain false-positive and false-negative result.

False-positive (there is no pregnancy, but the test shows a positive result). Possible reasons:

- The test was performed during administration of drugs for infertility, containing HCG (Pregnil, Profasi), or within 10 days after the end of such treatment course;
- Presence of tumors (chorionepithelioma of the uterus or hydatiform mole);
- Incomplete removal of fetal tissue after previous abortion.

False-negative (there is pregnancy, but the test shows a negative result). Possible reasons:

- The test was carried out before the term of positive result obtainment (in particular, this result may be observed in women with an irregular menstrual cycle);
- Pathology of the kidneys and cardio-vascular system, which prevents urinary excretion of HCG within normal concentrations;
- Increased fluid intake or administration of diuretics before the test.

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I hereby certify too that, according to the Companies' Register which I have inspected on July 13th 2016:

- In the Companies' Register of the Local court of Darmstadt the Axiom Gesellschaft für Diagnostica und Biochemica mbH with its headquarters located in Bürstadt is registered under registration number HRB 61822*
- In the Companies' Register, the above mentioned Wolfgang Hypko is registered as director.*
- According to the register, he has the power to represent the company individually and he is not bound by the restrictions of article 181 BGB (German Civil Code), which limits a representative's power to enter into contracts with himself or with himself representing also a third party at the same time.*

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1. Used Materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110—160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha HCG antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

Table 4

Product	Characteristics
Strip	Sizes: 0.1 x 0.4 x 7.5 cm

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1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110— 160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha HCG antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

**Table 4.**

Product	Characteristics
Strip	Sizes: 0.1 x 0.4 x 7.5 cm

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1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110— 160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
5	Housing	---	ABS	750 4D002	Kumho Petrochemical Co., Ltd.	Korea
6	Pipette	---	Polyethylene	Lupolen 1800U	LyondellBasel	Germany
7	Container for urine	---	PET	TexPet	TK Chemical Corp.	Korea

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha hCG antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4


Table 4.

Product	Characteristics
Cassette	Sizes: 2.5x7.5 cm

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1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

S. No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110—160 µm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 µm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
5	Housing	---	ABS	750 4D002	Kumho Petrochemical Co., Ltd.	Korea
6	Cap	---	ABS	750 4D002	Kumho Petrochemical Co., Ltd.	Korea

Table 2. Membrane

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha hCG antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

Table 4.

Product	Characteristics
Cassette	Sizes: 1.3 x 14.4 cm

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1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110— 160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
5	Housing	---	ABS	750 4D002	Kumho Petrochemical Co., Ltd.	Korea
6	Cap	---	ABS	750 4D002	Kumho Petrochemical Co., Ltd.	Korea

1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110—160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
5	Container for urine	---	PET	TexPet	TK Chemical Corp.	Korea

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: 4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha LH antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 1 mg/mL Preservative: 0.1 % NaN_3 Storage: 4°C Elisa test: passed

Table 3. Colloidal gold

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 1 mg/mL Preservative: 0.1 % NaN ₃ Storage: 4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

Product	Characteristics
Strip	Sizes: 0.1 x 0.4 x 7.5 cm

Set of tests to determine ovulation and pregnancy FRAUTEST planning:
Test to determine pregnancy

1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110— 160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
7	Container for urine	---	PET	TexPet	TK Chemical Corp.	Korea

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C
Test line	Anti-alpha hCG antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: +4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 2 mg/mL Preservative: 0.1 % NaN ₃ Storage: +4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

**Table 4**

Product	Characteristics
Strip	Sizes: 0.1 x 0.4 x 7.5 cm

Set of tests to determine ovulation and pregnancy FRAUTEST planning:

Test to determine ovulation

1. Used materials

Materials employed must comply with Table 1, 2 and 3

Table 1. Test strip

No	Component	Specification	Material	Trade name	Manufacturer	Country of origin
1	Membrane	Membrane with high throughput capacity. Thickness: 110—160 μm Absorption time: 1.5 – 3.0 mm/sec Intensity: ≥ 1 lb/in	See Table 2	-----	AXIOM GmbH	Germany
2	Absorbent paper (Glass fiber)	Absorbent cotton Thickness: 15 mil Pore size: 39 μm	Glass fiber	AX	AXIOM GmbH	Germany
3	Colloidal gold	Particle size: 45 nm	See Table 3	-----	AXIOM GmbH	Germany
4	Plastic plate	Size: 8.0 x 32 cm Thickness: 0.5 mm Color: white Weight: 19 g	Polypropylene	Moplen HP 648T	LyondellBasel	Germany
7	Container for urine	---	PET	TexPet	TK Chemical Corp.	Korea

Table 2. Membrane

Membrane lines	Material name	Specification
Control line	Goat anti-mouse immunoglobulin G	Classification: immunoglobulin G Activity: double diffusion: 1.8 Concentration: 2 mg/mL Preservative: 0.1 % NaN_3 Storage: 4°C Buffer: PH 7.4 PBS(phosphate buffered saline)
Test line	Anti-alpha LH antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 1 mg/mL Preservative: 0.1 % NaN_3 Storage: 4°C Elisa test: passed

Table 3. Colloidal gold

	Material name	Specification
Antibodies in colloidal gold	Anti-beta human chorionic gonadotropin antibodies	Classification: immunoglobulin G ₁ Purity: electrophoresis shows a purified band Concentration: 1 mg/mL Preservative: 0.1 % NaN ₃ Storage: 4°C Elisa test: passed

2. Quality indicators

Quality indicators must comply with Table 4

Product	Characteristics
Strip	Sizes: 0.1 x 0.4 x 7.5 cm

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